

June 2010

This distribution contains change pages for patch MD*1.0*21 of the Clinical Procedures 1.0 Technical Manual.

The change pages for CP Patch 11, should be inserted before the change pages for CP Patch 21:

File Name:
MD_1_P21_TM.PDF

Patch:
MD*1.0*21

Patch MD*1.0*21 pages:

Replace Pages:
Title page
Revision History
Table of Contents
4-1 to 4-4
6-17 to 6-26

With Pages:
Title page
Revision History
Table of Contents
4-1 to 4-4
6-17 to 6-28



CLINICAL PROCEDURES

TECHNICAL MANUAL AND PACKAGE SECURITY GUIDE

Version 1.0

April 2004

Revised June 2010

Department of Veterans Affairs
Office of Information & Technology
Office of Enterprise Development

Revision History

Description	Date	
Originally released.	April 2004	
¹ Patch MD*1.0*1 released.	July 2004	
Patch MD*1.0*2 released.	August 2006	
² Patch MD*1.0*5 released August 2006. Updated File List, Package Default Definition, Parameter Definitions, and menu options.	Documented February 2008	Shirley Ackerman, Alfred Bustamante
³ Patch MD*1.0*14 released. Updated Routine Descriptions, File List, Parameter Definitions, Protocols, menu options, and Cross References. Deleted bad references to Sample Reports in Ch. 15.	March 2008	Shirley Ackerman, Alfred Bustamante
⁴ Patch MD*1.0*6 released. Added description of Hemodialysis module and 508 Compliance to Introduction; updated Routine Descriptions, File List, Package Default Definition, Remote Procedure Calls, Parameter Definitions, menu options, Cross References, Callable Routines, External Relations, Internal Relations, and Glossary. Removed individual vendor contact information from Ch.15.	May 2008	Shirley Ackerman, Alfred Bustamante
⁵ Patch MD*1.0*11 released. Updated Routine Description, File and Field Description, Parameter Definition, and Menu Options By Name.	June 2009	Shirley Ackerman, Alfred Bustamante
⁶ Patch MD*1.0*21 released. Updated Routine Description, Parameter Definition, and Menu Options By Name.	June 2010	Shirley Ackerman, Rachel Wilder

¹ Patch MD*1.0*1 and MD*1.0*2 July 2004 Patch 2 release added.

² Patch MD*1.0*5 August 2006 Patch 5 release added.

³ Patch MD*1.0*14 March 2008 Patch 14 release added.

⁴ Patch MD*1.0*6 May 2008 Patch 6 release added.

⁵ Patch MD*1.0*11 June 2009 Patch 11 release added.

⁶ Patch MD*1.0*21 June 2010 Patch 21 release added.

Table of Contents

1. Introduction	1-1
Benefits	1-1
508 Compliance	1-2
2. Implementation and Maintenance	2-1
3. Clinical Instrument Interface Specifications	3-1
4. Routine Descriptions	4-1
5. File List and Related Information	5-1
File and Field Descriptions	5-1
Package Default Definition	5-18
6. Exported Options	6-1
Delphi Components	6-1
Remote Procedure Calls (RPC)	6-3
Parameter Definitions	6-14
Protocols	6-21
HL7 Application Parameters	6-22
HL Logical Links	6-23
Menu Options by Name	6-24
7. Cross-References	7-1
8. Archiving and Purging	8-1
9. Callable Routines	9-1
10. External Relations	10-1
11. Internal Relations	11-1
12. Package-wide Variables	12-1
13. SAC Exemptions	13-1
14. Software Product Security	14-1
Security Management	14-1
Security Features	14-1
15. Vendor Interfaces	15-1
List of Vendor Interfaces	15-1
Device Setup Instructions	15-1
Clinivision.....	15-2
Endoworks	15-6
Muse.....	15-8
Sensormedics V-MAX.....	15-9
B. Braun	15-12
Fresenius Medical Care	15-13
Gambro	15-14

16. Glossary 16-1

4. Routine Descriptions

```
MDAPI ; HOIFO/DP/NCA - CP API Calls ; [05-05-2003 10:28]
; ;1.0;CLINICAL PROCEDURES;**6**;Apr 01, 2004;Build 103
MDAPI1 ; HOIFO/NCA - Electrocardiogram Data Extraction ;12/4/02 12:32
; ;1.0;CLINICAL PROCEDURES;**1**;Apr 01, 2004;Build 4
MDAR7M ; HOIFO/NCA - Get Text Impression ;2/27/09 12:38
; ;1.0;CLINICAL PROCEDURES;**21**;Apr 01, 2004;Build 24
MDARP3 ; HOIFO/NCA - Get Procedures for Medicine ;1/13/04 14:35
; ;1.0;CLINICAL PROCEDURES;**10,13**;Apr 01, 2004;Build 22
MDARSET ; HOIFO/NCA - High Volume Check-In Setup ;6/30/09 10:00
; ;1.0;CLINICAL PROCEDURES;**21**;Apr 01, 2004;Build 24
MDCLN ; HOIFO/NCA - Cleanup Disabled Studies ;4/19/01 11:52
; ;1.0;CLINICAL PROCEDURES;**21**;Apr 01, 2004;Build 24
MDCVT ; HOIFO/DP/NCA - Medicine Package Conversion ;10/20/04 12:49
; ;1.0;CLINICAL PROCEDURES;**5**;Apr 01, 2004;Build 4
MDCVT1 ; HOIFO/NCA - Medicine Package Conversion (Cont.) ;1/6/05 15:12
; ;1.0;CLINICAL PROCEDURES;**5**;Apr 01, 2004;Build 4
MDCVTU ; HOIFO/NCA - Medicine Conversion Verification Utility ; [08-28-2003
11:34]
; ;1.0;CLINICAL PROCEDURES;**5**;Apr 01, 2004;Build 4
MDESPRT ; HOIFO/NCA - ELECTRONIC SIGNATURE PRINT ;12/21/04 09:24
; ;1.0;CLINICAL PROCEDURES;**5**;Apr 01, 2004;Build 4
MDHL7A ; HOIFO/WAA - Routine to Decode HL7 for CP ;05/21/09 15:57
; ;1.0;CLINICAL PROCEDURES;**6,11,21**;Apr 01, 2004;Build 24
MDHL7B ; HOIFO/WAA -Bi-directional interface routine ;7/23/01 11:41
; ;1.0;CLINICAL PROCEDURES;;Apr 01, 2004;Build 4
MDHL7BH ; HOIFO/WAA -Bi-directional interface (HL7) routine ;10/26/09 09:21
; ;1.0;CLINICAL PROCEDURES;**11,21**;Apr 01, 2004;Build 24
MDHL7D ; HOIFO/WAA -B-Braun, Fresenius Dialysis ; 06/08/00
; ;1.0;CLINICAL PROCEDURES;**6**;Apr 01, 2004;Build 103
MDHL7E ; HOIFO/WAA -Olympus/CMore/Pentax Endoscopy ; 06/08/00
; ;1.0;CLINICAL PROCEDURES;;Apr 01, 2004;Build 4
MDHL7K1 ; HOIFO/WAA-KenitDx Interface ; 06/08/00
; ;1.0;CLINICAL PROCEDURES;**21**;Apr 01, 2004;Build 24
MDHL7K2 ; HOIFO/WAA -HP EnConcert Echo ; 06/08/00
; ;1.0;CLINICAL PROCEDURES;;Apr 01, 2004;Build 4
MDHL7M1 ; HOIFO/WAA - Muse EKG ; [02-06-2002 16:13]
; ;1.0;CLINICAL PROCEDURES;;Apr 01, 2004;Build 4
MDHL7MCA ; HOIFO/REL-Routine to Decode HL7 for MEDICINE ; [05-07-2001 10:38]
; ;1.0;CLINICAL PROCEDURES;**6**;Apr 01, 2004;Build 103
MDHL7MCX ; HIRMF/WAA - Generate HL7 Error Message for MEDICINE ; [05-07-2001
10:38]
; ;1.0;CLINICAL PROCEDURES;;Apr 01, 2004;Build 4
MDHL7P1 ; HOIFO/WAA-Sensormedics,Jaeger Pulmonary ; 06/08/00
; ;1.0;CLINICAL PROCEDURES;;Apr 01, 2004;Build 4
MDHL7R1 ; HOIFO/WAA -Clinivision Resporatory ; 06/13/02
; ;1.0;CLINICAL PROCEDURES;;Apr 01, 2004;Build 4
MDHL7U ; HOIFO/WAA -Routine utilities for CP ;7/23/01 11:41
; ;1.0;CLINICAL PROCEDURES;;Apr 01, 2004;Build 4
MDHL7U1 ; HOIFO/WAA -Routine utilities for CP PROCESSING OBX ; 7/26/00
; ;1.0;CLINICAL PROCEDURES;**11**;Apr 01, 2004;Build 68
MDHL7U2 ; HOIFO/WAA -Utilities for CP PROCESSING OBX text ; 7/26/00
; ;1.0;CLINICAL PROCEDURES;;Apr 01, 2004;Build 4
MDHL7U3 ; HOIFO/WAA -Utilities for CP to process HL7 messages ;02/17/10 15
:59
; ;1.0;CLINICAL PROCEDURES;**6,21**;Apr 01, 2004;Build 24
MDHL7X ; HOIFO/WAA -Generate HL7 Error Message ; 06/08/00
; ;1.0;CLINICAL PROCEDURES;**6**;Apr 01, 2004;Build 103
MDHL7XXX ; HOIFO/DP - Loopback device for CP ;4/10/09 09:20
; ;1.0;CLINICAL PROCEDURES;**21**;Apr 01, 2004;Build 24
MDKRPC1 ; HOIFO/FT-RPC to return patient data ;2/19/08 13:13
; ;1.0;CLINICAL PROCEDURES;**6**;Apr 01, 2004;Build 103
MDKRPC2 ; HOIFO/DP - RPC Calls (Cont.) ;11/27/07 09:42
; ;1.0;CLINICAL PROCEDURES;**6**;Apr 01, 2004;Build 103
MDKUTL ; HOIFO/DP - Renal Utilities ;11/29/07 14:45
; ;1.0;CLINICAL PROCEDURES;**14**;Apr 01, 2004;Build 22
MDKUTLR ; HOIFO/DP - Renal Utilities RPC;11/29/07 14:45
; ;1.0;CLINICAL PROCEDURES;**6**;Apr 01, 2004;Build 103
MDNCHK ; HOIFO/NCA - CP Multiple Result Check ;4/26/05 15:17
```

¹ Patch MD*1.0*21 June 2010 Update routine list with new routines and patch history changes.

Routine Descriptions

```

MDOUTOR ;1.0;CLINICAL PROCEDURES;**11,21**;Apr 01, 2004;Build 68
; HOIFO/NCA - Post Conversion Routine ; [04-14-2003 10:51]
;1.0;CLINICAL PROCEDURES;**5**;Apr 01, 2004;Build 4
MDPCE ; HIRMF/NCA - Routine For Data Extract ;6/9/08 13:29
;1.0;CLINICAL PROCEDURES;**5,21**;Apr 01, 2004;Build 24
MDPCE1 ; HOIFO/NCA - Updated Routine For Data Extract ; [05-28-2002 12:55]
;1.0;CLINICAL PROCEDURES;**6**;Apr 01, 2004;Build 103
MDPCE2 ; HOIFO/NCA - Routine For Data Extract For Hemo Dialysis;9/10/04 11
:23 ;1/20/10 10:00
;1.0;CLINICAL PROCEDURES;**6,21**;Apr 01, 2004;Build 24
MDPFTP1 ;HOIFO/NCA - PFT REPORT-DEMO INFO ;3/15/04 11:55
;1.0;CLINICAL PROCEDURES;**2**;Apr 01, 2004;Build 4
MDPFTP2 ; HOIFO/NCA - PFT REPORT-VOLUMES ;3/15/04 10:00
;1.0;CLINICAL PROCEDURES;**2**;Apr 01, 2004;Build 4
MDPFTP2A ; HOIFO/NCA - PFT REPORT-FLOWS ;3/17/04 08:22
;1.0;CLINICAL PROCEDURES;**2**;Apr 01, 2004;Build 4
MDPFTP3 ; HOIFO/NCA - PFT REPORT-SPECIAL STUDIES (PT 2) ;3/17/04 12:48
;1.0;CLINICAL PROCEDURES;**2**;Apr 01, 2004;Build 4
MDPOST ; HOIFO/DP - Post Init ;2/18/04 11:39
;1.0;CLINICAL PROCEDURES;Apr 01, 2004;Build 4
MDPOST04 ; HOIFO/DP - Post Init ; 2/18/04 11:39
;1.0;CLINICAL PROCEDURES;**4**;Apr 01, 2004;Build 6
MDPOST06 ; HOIFO/DP - Post Init ;2/7/07 16:15
;1.0;CLINICAL PROCEDURES;**6**;Apr 01, 2004;Build 103
MDPOST1 ; HOIFO/NCA/DP - Build CP DEFINITION file (#702.01) - Optional Post
Init ; [12-04-2002 13:06]
;1.0;CLINICAL PROCEDURES;Apr 01, 2004;Build 4
MDPOST21 ; HOIFO/NCA - Post Init ;2/7/07 16:15
;1.0;CLINICAL PROCEDURES;**21**;Apr 01, 2004;Build 24
MDPOST6A ;HOIFO/NCA-Convert Existing Notes to New File ;11/28/07 14:31
;1.0;CLINICAL PROCEDURES;**6**;Apr 01, 2004;Build 103
MDPS1 ; HOIFO/NCA - CP/Medicine Report Generator ;5/18/04 09:48
;1.0;CLINICAL PROCEDURES;**2,10,13,21**;Apr 01, 2004;Build 24
MDPS2 ; HOIFO/NCA - CP/Medicine Report Generator (Cont.) ;5/18/04 09:41
;1.0;CLINICAL PROCEDURES;**2**;Apr 01, 2004;Build 4
MDPS3 ; HOIFO/NCA - Remote Data View Data Retriever for CP ;8/26/05 14:37
;1.0;CLINICAL PROCEDURES;**2,5,13**;Apr 01, 2004;Build 22
MDPS4 ; HOIFO/NCA - Retrieve List of Consult Procedures ;1/26/06 12:45
;1.0;CLINICAL PROCEDURES;**13**;Apr 01, 2004;Build 22
MDPS5 ; HOIFO/NCA - Retrieve List of Consult Procedures for RDV ;3/4/05 1
3:29
;1.0;CLINICAL PROCEDURES;**13**;Apr 01, 2004;Build 22
MDPSU ; HOIFO/NCA - CP/Medicine Report Generator Utility;5/18/04 09:48
;1.0;CLINICAL PROCEDURES;**21**;Apr 01, 2004;Build 24
MDPSUL ; HOIFO/NCA - HS Component Utility;5/18/04 09:48 ;10/5/09 09:33
;1.0;CLINICAL PROCEDURES;**21**;Apr 01, 2004;Build 24
MDPURGE ;HOIFO/NCA - Study Clean-Up process ;6/18/08 10:15
;1.0;CLINICAL PROCEDURES;**11**;Apr 01, 2004;Build 68
MDRPCNT ; HOIFO/NCA - Document Handler Object (TMDNOTE) ;5/23/05 15:50
;1.0;CLINICAL PROCEDURES;**6**;Apr 01, 2004;Build 103
MDRPCNT1 ; HOIFO/NCA - Object RPCs (TMDNOTE) Continued 2;10/29/04 12:20 ;2/2
5/09 16:08
;1.0;CLINICAL PROCEDURES;**6,21**;Apr 01, 2004;Build 24
MDRPCOD ; HOIFO/DP - Object RPCs (TMDProcedureDef) ; [01-09-2003 15:20]
;1.0;CLINICAL PROCEDURES;Apr 01, 2004;Build 4
MDRPCOG ; HOIFO/DP - CP Gateway ; [01-09-2003 15:20]
;1.0;CLINICAL PROCEDURES;**6**;Apr 01, 2004;Build 103
MDRPCOL ; HOIFO/DP - Object RPCs (Logfile) ; [02-11-2002 13:41]
;1.0;CLINICAL PROCEDURES;Apr 01, 2004;Build 4
MDRPCOO ; HOIFO/DP - Object RPCs (TMDOutput) ; [03-24-2003 15:44]
;1.0;CLINICAL PROCEDURES;Apr 01, 2004;Build 4
MDRPCOP ; HOIFO/DP - Object RPCs (TMDPatient) ;3/12/08 09:16
;1.0;CLINICAL PROCEDURES;**4,6,11**;Apr 01, 2004;Build 68
MDRPCOP1 ; HOIFO/DP - Object RPCs (TMDPatient) - Cont. ; 01-09-2003 15:21
;1.0;CLINICAL PROCEDURES;**6**;Apr 01, 2004;Build 103
MDRPCOR ; HOIFO/DP - Object RPCs (TMDRecordId) ; [01-10-2003 09:14]
;1.0;CLINICAL PROCEDURES;**17**;Apr 01, 2004;Build 4
MDRPCOT ; HOIFO/DP/NCA - Object RPCs (TMDTransaction) ;10/26/09 10:23
;1.0;CLINICAL PROCEDURES;**5,6,11,21**;Apr 01, 2004;Build 24
MDRPCOT1 ; HOIFO/NCA/DP - Object RPCs (TMDTransaction) - Continued ;3/13/09
11:18
;1.0;CLINICAL PROCEDURES;**5,11,21**;Apr 01, 2004;Build 24
MDRPCOT2 ; HOIFO/NCA - Object RPCs (TMDTransaction) Continued 2;10/29/04 12:
20 ;3/12/08 09:18
;1.0;CLINICAL PROCEDURES;**6,21**;Apr 01, 2004;Build 24

```

```

MDRPCOTH ; HOIFO/NCA - Process High Volume Procedure Results ;2/27/09 10:08
;;1.0;CLINICAL PROCEDURES;**21**;Apr 01, 2004;Build 24
MDRPCOU ; HOIFO/DP - Object RPCs (TMDUser) ; [01-09-2003 15:21]
;;1.0;CLINICAL PROCEDURES;;Apr 01, 2004;Build 4
MDRPCOV ; HOIFO/DP - Object RPCs (TMDParameter) ; [04-15-2003 12:42]
;;1.0;CLINICAL PROCEDURES;;Apr 01, 2004;Build 4
MDRPCOW ; HOIFO/DP/NCA - Billing Widget ;10/3/05 12:17
;;1.0;CLINICAL PROCEDURES;**6**;Apr 01, 2004;Build 103
MDRPCU ; HOIFO/DP - Object RPC Utilities ; [05-23-2003 10:16]
;;1.0;CLINICAL PROCEDURES;**4**;Apr 01, 2004;Build 6
MDRPCW ; HOIFO/NCA - Calls to AICS;04/01/2003 ;01/21/10 11:51
;;1.0;CLINICAL PROCEDURES;**6,21**;Apr 01, 2004;Build 24
MDRPCW1 ; HOIFO/NCA - MD TMDENCOUNTER Object; [05-28-2002 12:55] ;2/16/10 1
6:17
;;1.0;CLINICAL PROCEDURES;**6,21**;Apr 01, 2004;Build 24
MDRPCWU ; HOIFO/NCA - CPT Code Query; [05-28-2002 12:55] ;2/16/10 16:17
;;1.0;CLINICAL PROCEDURES;**21**;Apr 01, 2004;Build 24
MDSTATU ; HOIFO/NCA - Print List of Document Titles Needed ;10/21/04 13:44
;;1.0;CLINICAL PROCEDURES;**5**;Apr 01, 2004;Build 4
MDSTUDL ; HOIFO/NCA - Clinical Procedures Studies List ;10/26/05 11:46
;;1.0;CLINICAL PROCEDURES;**6**;Apr 01, 2004;Build 103
MDSTUDW ; HOIFO/NCA - Print a List of Procedures With Incomplete Workload ;3
/2/09 10:00
;;1.0;CLINICAL PROCEDURES;**21**;Apr 01, 2004;Build 24
MDUXML ; HOIFO/WAA -Utilities for XML text ; 7/26/00
;;1.0;CLINICAL PROCEDURES;**6**;Apr 01, 2004;Build 103
MDUXMLM ; HOIFO/WAA -Utilities for XML text ; 7/26/00
;;1.0;CLINICAL PROCEDURES;**6**;Apr 01, 2004;Build 103
MDUXMLOX ; HOIFO/WAA -OBX converter XML text ; 7/26/00
;;1.0;CLINICAL PROCEDURES;**6**;Apr 01, 2004;Build 103
MDUXMLU1 ; HOIFO/WAA -Utilities for XML text ; 7/26/00
;;1.0;CLINICAL PROCEDURES;**6**;Apr 01, 2004;Build 103
MDWCAN ;HOIFO/NCA - Process No-Shows and Cancels ;7/29/08 09:50
;;1.0;CLINICAL PROCEDURES;**11,21**;Apr 01, 2004;Build 24
MDWCHK ; HOIFO/NCA - Create CP Studies for Existing Procedures ;12/13/07 1
5:52
;;1.0;CLINICAL PROCEDURES;**14**;Apr 01,2004;Build 22
MDWOR ; HOIFO/NCA - Main Routine to Decode HL7 ;9/8/08 15:20
;;1.0;CLINICAL PROCEDURES;**14,11,21**;Apr 01,2004;Build 24
MDWORC ; HOIFO/NCA - Main Routine to Decode HL7 from Consult ;1/8/08 15:00
;;1.0;CLINICAL PROCEDURES;**14**;Apr 01,2004;Build 22
MDWORSR ; HOIFO/NCA - Daily Schedule Studies;7/2/04 12:39 ;10/15/08 13:39
;;1.0;CLINICAL PROCEDURES;**14,11,21**;Apr 01,2004;Build 24
MDWSETUP ; HOIFO/NCA - Auto Study Check-In Setup ;3/18/08 14:14
;;1.0;CLINICAL PROCEDURES;**14,11**;Apr 01, 2004;Build 68
MDXMLFM ; HOIFO/DP - Fileman -> XML Utilities ; [01-10-2003 09:14]
;;1.0;CLINICAL PROCEDURES;**6**;Apr 01, 2004;Build 103
MDXMLFM1 ; HOIFO/DP/NCA - Data -> XML Utilities ; [01-10-2003 09:14]
;;1.0;CLINICAL PROCEDURES;**6**;Apr 01, 2004;Build 103

```

Routine Descriptions

|

DESCRIPTION:

This parameter stores a list of valid file types and the associated extensions of these files.

PRECEDENCE: 1 ENTITY FILE: SYSTEM

NAME: **MD GATEWAY**

MULTIPLE VALUED: Yes

VALUE TERM: Parameter Value

VALUE DOMAIN: 1:255

INSTANCE DOMAIN: 1:255

PRECEDENCE: 1

DISPLAY TEXT: CP Gateway Parameters

INSTANCE TERM: Parameter Name

VALUE DATA TYPE: free text

INSTANCE DATA TYPE: free text

ENTITY FILE: SYSTEM

¹NAME: **MD GET HIGH VOLUME**

MULTIPLE VALUED: Yes

VALUE TERM: Get String

INSTANCE DATA TYPE: pointer

INSTANCE HELP: Enter a high volume procedure.

INSTANCE SCREEN CODE: I

+\$(^MDS(702.01,+Y,0),"^",6)'=2&(+\$(^MDS(702.01,+Y,0),"^",11)'=2)&(\$(^MDS(702.01,+Y,0),"^",9)>0)

DESCRIPTION:

This parameter will contain a free text string that contains two pieces of data delimited by a semicolon ';'. The two pieces of data are: 1) 1/0 (Yes/No) to indicate whether or not the text of the result should be added to the note, 2) 1/0 (Yes/No) to enter the text of the result as the significant finding of the Consult. (If you enter a 0, the note will be auto closed with the text inside.)

Example string: 1;0

PRECEDENCE: 1

ENTITY FILE: SYSTEM

NAME: **MD HFS SCRATCH**

DISPLAY TEXT: Vista Scratch HFS Directory

MULTIPLE VALUED: No

VALUE DATA TYPE: free text

VALUE HELP: Enter in an OS level directory

DESCRIPTION:

Contains the directory specification for the Kernel OPEN^%ZISH call. This directory should be accessible for read/write operations by all CP users.

PRECEDENCE: 1

ENTITY FILE: SYSTEM

NAME: **MD IMAGING XFER**

MULTIPLE VALUED: No

VALUE DATA TYPE: free text

DESCRIPTION:

This parameter contains the name of a network server, share, and path (UNC) to a location where Clinical Procedures can put files for pick-up by the Imaging background processor for archiving.

PRECEDENCE: 1

ENTITY FILE: SYSTEM

¹ Patch MD*1.0*21 June 2010 Parameter Definition added.

Exported Options

¹ NAME: **MDK APPLICATION INSTALL** DISPLAY TEXT: MDK Application
Install
MULTIPLE VALUED: Yes
INSTANCE TERM: Installation Distribution Info
VALUE TERM: Distribution Info Value PROHIBIT EDITING: No
VALUE DATA TYPE: free text VALUE DOMAIN: 1:250
INSTANCE DATA TYPE: free text INSTANCE DOMAIN: 1:250
DESCRIPTION:
This parameter is used to store the Hemodialysis application
distribution information. The information includes the following:
1) Date/Time when application first launched.
2) User Name
3) System Option Loaded (Y/N)
4) Workstation of where the application was launched.
PRECEDENCE: 1 ENTITY FILE: SYSTEM

NAME: **MDK GUI VERSION**
DISPLAY TEXT: Hemodialysis Version Compatibility
MULTIPLE VALUED: Yes INSTANCE TERM: Application:Version
VALUE TERM: Compatible with current server version
PROHIBIT EDITING: No VALUE DATA TYPE: yes/no
INSTANCE DATA TYPE: free text INSTANCE DOMAIN: 1:40
DESCRIPTION:
This parameter is used to store the application:versions that are compatible
with the current server version of Hemodialysis. Instance format
of APPLICATION:VERSION (example: HEMODIALYSIS.EXE:0.0.0.0).
PRECEDENCE: 1 ENTITY FILE: SYSTEM

² NAME: **MD MEDICINE CONVERTED** DISPLAY TEXT: Medicine Package
Converted
MULTIPLE VALUED: No VALUE TERM: Yes/No
PROHIBIT EDITING: No VALUE DATA TYPE: yes/no
DESCRIPTION:
Used to determine if the Medicine Package has been converted.
PRECEDENCE: 1 ENTITY FILE: SYSTEM

³ NAME: **MD NOT ADMN CLOSE MUSE NOTE** DISPLAY TEXT: NOT ADMN Close Muse
Note
MULTIPLE VALUED: No VALUE TERM: Yes/No
PROHIBIT EDITING: No VALUE DATA TYPE: yes/no
DESCRIPTION:
This parameter is used to indicate the note should not be
administratively closed with the proxy user CLINICAL, DEVICE PROXY
SERVICE but the interpreter of the procedure for the MUSE device.
The default is "No".
PRECEDENCE: 1 ENTITY FILE: SYSTEM

NAME: **MD OFFLINE MESSAGE** DISPLAY TEXT: Offline message
MULTIPLE VALUED: No VALUE TERM: Offline Message
VALUE DATA TYPE: word processing
DESCRIPTION:

¹ Patch MD*1.0*6 May 2008 Parameter Definition added.

² Patch MD*1.0*5 August 2006 Parameter Definition added.

³ Patch MD*1.0*21 June 2010 Parameter Definition added.

This parameter contains a message to display to the users when the Clinical Procedures application is offline.

PRECEDENCE: 1 ENTITY FILE: SYSTEM

¹NAME: **MD OLYMPUS 7** DISPLAY TEXT: MD OLYMPUS 7
 MULTIPLE VALUED: No VALUE TERM: Yes/No
 PROHIBIT EDITING: No VALUE DATA TYPE: yes/no
 VALUE HELP: Enter Yes/No whether you have Olympus version 7.3.7.
 DESCRIPTION:
 This parameter definition indicates whether the Olympus device is version 7.3.7. The value is Yes/No. The default value is "No".
 PRECEDENCE: 1 ENTITY FILE: SYSTEM

NAME: **MD ONLINE**
 DISPLAY TEXT: Clinical Procedure Online/Offline
 MULTIPLE VALUED: No
 VALUE TERM: Is Clinical Procedures Online
 PROHIBIT EDITING: No VALUE DATA TYPE: yes/no
 VALUE HELP: Enter 'Yes' to allow access to CP
 DESCRIPTION:
 This parameter controls access to the Clinical Procedures package.
 PRECEDENCE: 1 ENTITY FILE: SYSTEM

²NAME: **MD USE APPOINTMENT** DISPLAY TEXT: Use Appointment
 Location
 MULTIPLE VALUED: No VALUE TERM: Use Appointment location
 VALUE DATA TYPE: yes/no
 DESCRIPTION:
 Set this value to Yes to allow CPUser to use the location of the appointment selected during CP study check-in for the workload. Otherwise, the hospital location of the CP Definition will be used.
 Enter RETURN to continue or '^' to exit:
 If no value is entered, the default value is No.
 PRECEDENCE: 1 ENTITY FILE: SYSTEM

³NAME: **MD USE APPT WITH PROCEDURE**
 DISPLAY TEXT: Use Appointment With Procedure
 MULTIPLE VALUED: No
 VALUE TERM: Use appointment with procedure
 PROHIBIT EDITING: No VALUE DATA TYPE: yes/no
 DESCRIPTION:
 Enter "Y" or "N" for Yes/No on whether your site selects the appointment scheduled for outpatients during the procedure request in CPRS.
 PRECEDENCE: 1 ENTITY FILE: SYSTEM

⁴NAME: **MD USE NOTE** DISPLAY TEXT: Use Note
 VALUE TERM: Yes/No VALUE DATA TYPE: yes/no

¹ Patch MD*1.0*11 June 2009 Parameter Definition added.

² Patch MD*1.0*11 June 2009 Parameter Definition added.

³ Patch MD*1.0*14 March 2008 Parameter Definition added.

⁴ Patch MD*1.0*21 June 2010 Parameter Definition Added.

Exported Options

DESCRIPTION:

This parameter indicates that Clinical Procedures will use the note for the text of the result instead of the Significant Finding field in Consult.

PRECEDENCE: 1

ENTITY FILE: SYSTEM

NAME: MD USER DEFAULTS

MULTIPLE VALUED: Yes

VALUE TERM: Parameter value

VALUE DATA TYPE: free text

INSTANCE DATA TYPE: free text

DESCRIPTION:

This parameter is used to store a users default parameter settings. Each setting is defined on the client.

PRECEDENCE: 1

DISPLAY TEXT: CP User Defaults

INSTANCE TERM: Parameter setting

PROHIBIT EDITING: No

VALUE DOMAIN: 1:250

INSTANCE DOMAIN: 1:250

ENTITY FILE: USER

NAME: MD VERSION CHK

MULTIPLE VALUED: Yes

VALUE TERM: Compatible with current server version

PROHIBIT EDITING: No

INSTANCE DATA TYPE: free text

DESCRIPTION:

This parameter is used to store the application:versions that are compatible with the current server version of Clinical Procedures. Instance format of APPLICATION:VERSION (example: CPMANAGER.EXE:0.0.0.0).

PRECEDENCE: 1

DISPLAY TEXT: Version Compatibility

INSTANCE TERM: Application:Version

VALUE DATA TYPE: yes/no

INSTANCE DOMAIN: 1:30

ENTITY FILE: SYSTEM

NAME: MD WEBLINK

DISPLAY TEXT: Clinical Procedures Home Page

MULTIPLE VALUED: No

VALUE DATA TYPE: free text

DESCRIPTION:

This parameter contains the web address for the Clinical Procedures home page. This can be modified to a local address in the event that the pages are downloaded to be displayed from a local server location.

PRECEDENCE: 1

VALUE TERM: Web Address

VALUE DOMAIN: 1:250

ENTITY FILE: SYSTEM

Protocols

NAME: **MCAR Device Client** ITEM TEXT: Instrument Device Client
 TYPE: subscriber CREATOR: ACKERMAN,NIEN-CHIN
 PACKAGE: MEDICINE
 DESCRIPTION: Subscriber protocol for sending data to Vista from clinical instruments.
 TIMESTAMP: 59276,54156 RECEIVING APPLICATION: MCAR-INST
 TRANSACTION MESSAGE TYPE: ORU EVENT TYPE: R01
 PROCESSING ID: P LOGICAL LINK: MCAR INST
 VERSION ID: 2.3 RESPONSE MESSAGE TYPE: ACK
 PROCESSING ROUTINE: D ^MDHL7A SENDING FACILITY REQUIRED?: NO
 RECEIVING FACILITY REQUIRED?: NO

NAME: **MCAR Device Server** ITEM TEXT: Instrument HL7 Event Driver
 TYPE: event driver CREATOR: ACKERMAN,NIEN-CHIN
 PACKAGE: MEDICINE
 DESCRIPTION: This protocol is used by the HL7 package to send results to Vista from various clinical instrumentation.
 TIMESTAMP: 59276,54156 SENDING APPLICATION: INST-MCAR
 TRANSACTION MESSAGE TYPE: ORU EVENT TYPE: R01
 PROCESSING ID: P VERSION ID: 2.3
 SENDING FACILITY REQUIRED?: NO RECEIVING FACILITY REQUIRED?: NO
 SUBSCRIBERS: MCAR Device Client

NAME: **MCAR ORM CLIENT** TYPE: subscriber
 CREATOR: ACKERMAN,NIEN-CHIN RECEIVING APPLICATION: INST-MCAR
 EVENT TYPE: O02 RESPONSE MESSAGE TYPE: ORR
 SENDING FACILITY REQUIRED?: NO RECEIVING FACILITY REQUIRED?: NO
 SECURITY REQUIRED?: NO ROUTING LOGIC: Q

NAME: **MCAR ORM SERVER**
 ITEM TEXT: Clinical Procedures ORM Protocol Server
 TYPE: event driver CREATOR: ACKERMAN,NIEN-CHIN
 TIMESTAMP: 59276,54156 SENDING APPLICATION: MCAR-INST
 TRANSACTION MESSAGE TYPE: ORM EVENT TYPE: O01
 VERSION ID: 2.3
 SUBSCRIBERS: MCAR ORM CLIENT

¹NAME: **MD RECEIVE GMRC**
 ITEM TEXT: Clinical Procedures receives messages from Consult
 TYPE: action CREATOR: ACKERMAN,NIEN-CHIN
 PACKAGE: CLINICAL PROCEDURES
 DESCRIPTION: This protocol receives messages from Consult. (IA 3140)
 ENTRY ACTION: D EN^MDWORC(.XQORMSG) TIMESTAMP: 60934,38793

NAME: **MD RECEIVE OR**
 ITEM TEXT: Clinical Procedures receives order msgs from CPRS
 TYPE: action CREATOR: ACKERMAN,NIEN-CHIN
 PACKAGE: CLINICAL PROCEDURES
 DESCRIPTION: This protocol receives order messages from CPRS. (IA 3135)
 ENTRY ACTION: D EN^MDWOR(.XQORMSG) TIMESTAMP: 60934,38793

¹ Patch MD*1.0*14 March 2008 Protocols added to support the auto study check-in.

HL7 Application Parameters

NAME: **INST-MCAR**

COUNTRY CODE: USA

HL7 FIELD SEPARATOR: |

ACTIVE/INACTIVE: ACTIVE

HL7 ENCODING CHARACTERS: ^~\&

NAME: **MCAR-INST**

FACILITY NAME: VISTA

COUNTRY CODE: USA

HL7 FIELD SEPARATOR: |

ACTIVE/INACTIVE: ACTIVE

MAIL GROUP: POSTMASTER

HL7 ENCODING CHARACTERS: ^~\&

HL Logical Links

NODE: **MCAR INST**

DEVICE TYPE: Single-threaded Server
 AUTOSTART: Enabled
 TASK NUMBER: 526320
 QUEUE SIZE: 100
 RE-TRANSMISSION ATTEMPTS: 3
 ACK TIMEOUT: 60
 TCP/IP PORT: 9026
 PERSISTENT: NO
 IN QUEUE BACK POINTER: 331
 OUT QUEUE BACK POINTER: 220

LLP TYPE: TCP

STATE: Reading
 TIME STARTED: MAR 04, 2004@06:46:17
 SHUTDOWN LLP ?: NO

READ TIMEOUT: 60
 EXCEED RE-TRANSMIT ACTION: ignore
 TCP/IP SERVICE TYPE: SINGLE LISTENER
 STARTUP NODE: DEV:ISC4A1
 IN QUEUE FRONT POINTER: 331
 OUT QUEUE FRONT POINTER: 210

NODE: **MCAR OUT**

DEVICE TYPE: Non-Persistent Client
 AUTOSTART: Enabled
 TASK NUMBER: 529066
 QUEUE SIZE: 100
 READ TIMEOUT: 60
 EXCEED RE-TRANSMIT ACTION: ignore
 TCP/IP PORT: 9028
 PERSISTENT: NO
 IN QUEUE BACK POINTER: 202
 OUT QUEUE BACK POINTER: 206

LLP TYPE: TCP

STATE: Openfail
 TIME STARTED: MAR 04, 2004@06:45:47
 SHUTDOWN LLP ?: NO

RE-TRANSMISSION ATTEMPTS: 3
 ACK TIMEOUT: 60
 TCP/IP ADDRESS: 10.3.17.157
 TCP/IP SERVICE TYPE: CLIENT (SENDER)
 STARTUP NODE: DEV:ISC4A1
 IN QUEUE FRONT POINTER: 202
 OUT QUEUE FRONT POINTER: 202

Menu Options by Name

NAME: **MD GUI USER** MENU TEXT: MD GUI USER
 TYPE: Broker (Client/Server) CREATOR: ACKERMAN,NIEN-CHIN
 TIMESTAMP OF PRIMARY MENU: 59331,44145
 RPC: MD TMDOUTPUT
 RPC: MD TMDPARAMETER
 RPC: MD TMDPATIENT
 RPC: MD TMDPROCEDURE
 RPC: MD TMDRECORDID
 RPC: MD TMDTRANSACTION
 RPC: MD TMDUSER
 RPC: MD UTILITIES
 UPPERCASE MENU TEXT: MD GUI USER

NAME: **MD GUI MANAGER** MENU TEXT: MD GUI MANAGER
 TYPE: Broker (Client/Server) CREATOR: ACKERMAN,NIEN-CHIN
 TIMESTAMP OF PRIMARY MENU: 59385,45622
 RPC: MD TMDOUTPUT
 RPC: MD TMDPARAMETER
 RPC: MD TMDPATIENT
 RPC: MD TMDPROCEDURE
 RPC: MD TMDRECORDID
 RPC: MD TMDTRANSACTION
 RPC: MD TMDUSER
 RPC: MD UTILITIES
 RPC: MD GATEWAY
 UPPERCASE MENU TEXT: MD GUI MANAGER

¹NAME: **MD AUTO CHECK-IN SETUP** MENU TEXT: Auto Study Check-In Setup
 TYPE: run routine CREATOR: ACKERMAN,NIEN-CHIN
 PACKAGE: CLINICAL PROCEDURES
 DESCRIPTION: This option is used to populate the XPAR parameters MD USE
 APPT WITH PROCEDURE, MD CHECK-IN PROCEDURE LIST, MD CLINIC QUICK LIST, and MD
 CLINICS WITH MULT PROC. The four XPAR parameters are used for the auto study
 check-in. Users can use the option to indicate whether their site use and
 schedule appointments. They can populate a list of procedures and associated
 clinics that need a CP study checked-in.
 ROUTINE: EN1^MDWSETUP
 UPPERCASE MENU TEXT: AUTO STUDY CHECK-IN SETUP

NAME: **MD SCHEDULED STUDIES** MENU TEXT: Scheduled Studies
 TYPE: run routine CREATOR: ACKERMAN,NIEN-CHIN
 PACKAGE: CLINICAL PROCEDURES
 DESCRIPTION: This option is tasked to run daily. It will process the HL7
 messages that need to be sent to the device on a daily basis for CP studies.
 ROUTINE: EN1^MDWORSR SCHEDULING RECOMMENDED: YES
 UPPERCASE MENU TEXT: SCHEDULED STUDIES

NAME: **MD STUDY CHECK-IN** MENU TEXT: Study Check-in
 TYPE: run routine CREATOR: ACKERMAN,NIEN-CHIN
 PACKAGE: CLINICAL PROCEDURES
 DESCRIPTION: This option is tasked to run daily. It checks-in CP studies

¹ Patch MD*1.0*14 March 2008 Options added to support the auto study check-in.

for procedures that require multiple encounters such as Hemodialysis, Respiratory Therapy, and Sleep Studies.

ROUTINE: CLINICPT^MDWORSR SCHEDULING RECOMMENDED: YES
UPPERCASE MENU TEXT: STUDY CHECK-IN

¹ NAME: **MD HIGH VOLUME PROCEDURE SETUP** MENU TEXT: High Volume Procedure Setup
TYPE: run routine CREATOR: ACKERMAN,NIEN-CHIN
PACKAGE: CLINICAL PROCEDURES
DESCRIPTION: This option will populate the XPAR Parameters MD GET HIGH VOLUME and MD NOT ADMN CLOSE MUSE NOTE. It will let the user populate a list of Clinical Procedures procedures set it up for high volume procedure process.
ROUTINE: EN1^MDARSET
UPPERCASE MENU TEXT: HIGH VOLUME PROCEDURE SETUP

NAME: **MD PROC W/INCOMPLETE WORKLOAD**
MENU TEXT: Print list of Procedure with incomplete workload
TYPE: run routine CREATOR: ACKERMAN,NIEN-CHIN
PACKAGE: CLINICAL PROCEDURES
DESCRIPTION: This option prints a list of procedures that has incomplete workload for the visit.
ROUTINE: E1^MDSTUDW
UPPERCASE MENU TEXT: PRINT LIST OF PROCEDURE WITH I

NAME: **MD PROCESS RESULTS** MENU TEXT: MD Process Results
TYPE: run routine CREATOR: ACKERMAN,NIEN-CHIN
PACKAGE: CLINICAL PROCEDURES E ACTION PRESENT: YES
DESCRIPTION: This task is ran daily for every hour to process results and update Consults package.
ENTRY ACTION: N ZTSAVE S ZTSAVE("DUZ")=DUZ,ZTSAVE("DUZ(")=""
ROUTINE: PROCESS^MDHL7XXX UPPERCASE MENU TEXT: MD PROCESS RESULTS

² NAME: **MD HEMODIALYSIS USER** MENU TEXT: HEMODIALYSIS USER
TYPE: Broker (Client/Server) CREATOR: ACKERMAN,NIEN-CHIN
TIMESTAMP OF PRIMARY MENU: 60387,39853
RPC: MDK GET VISTA DATA
RPC: MDK GET/SET RENAL DATA
RPC: MDK UTILITY
RPC: VAFCTFU CONVERT DFN TO ICN
RPC: VAFCTFU CONVERT ICN TO DFN
RPC: MD TMDWIDGET
RPC: MD TMDNOTE
RPC: MD TMDCIDC
RPC: MD TMDLEX
RPC: MD TMDENCOUNTER
RPC: GMV MANAGER
RPC: MD GATEWAY
RPC: MD TMDSUBMITU
RPC: ORWPT PTINQ
RPC: GMV PTSELECT
RPC: DG SENSITIVE RECORD ACCESS
RPC: DG SENSITIVE RECORD BULLETIN
RPC: MD TMDRECORDID

¹ Patch MD*1.0*21 June 2010 Options added to support high volume procedures enhancement.

² Patch MD*1.0*6 May 2008 Hemodialysis User menu option added.

Exported Options

UPPERCASE MENU TEXT: HEMODIALYSIS USER

NAME: MD STUDIES LIST

MENU TEXT: Clinical Procedures Studies List

TYPE: run routine

CREATOR: ACKERMAN,NIEN-CHIN

PACKAGE: CLINICAL PROCEDURES

DESCRIPTION: This option will generate a list of Clinical Procedures studies.

ROUTINE: EN2^MDSTUDL

UPPERCASE MENU TEXT: CLINICAL PROCEDURES STUDIES LI

¹NAME: MDCVT MANAGER

MENU TEXT: Medicine to CP Conversion Manager

TYPE: menu

CREATOR: ACKERMAN,NIEN-CHIN

PACKAGE: CLINICAL PROCEDURES

DESCRIPTION: This is the Medicine to CP Manager menu option. This menu option consists of options to assist the site in converting the Medicine reports to Clinical Procedures text reports.

ITEM: MDCVT SETUP SYNONYM: 1

DISPLAY ORDER: 1

ITEM: MDCVT RUN SYNONYM: 3

DISPLAY ORDER: 3

ITEM: MDCVT SUMMARY SYNONYM: 4

DISPLAY ORDER: 4

ITEM: MDCVT DISK SPACE SYNONYM: 5

DISPLAY ORDER: 5

ITEM: MDCVT LIST OF TIU TITLES SYNONYM: 6

DISPLAY ORDER: 6

ITEM: MDCVT TOTALS SYNONYM: 7

DISPLAY ORDER: 7

ITEM: MDCVT ERROR LOG SYNONYM: 8

DISPLAY ORDER: 8

ITEM: MDCVT CONVERSION LOCKOUT SYNONYM: 9

DISPLAY ORDER: 9

ITEM: MDCVT BUILD CONVERSION LIST SYNONYM: 2

DISPLAY ORDER: 2

TIMESTAMP: 60459,53192

TIMESTAMP OF PRIMARY MENU: 59904,24363

UPPERCASE MENU TEXT: MEDICINE TO CP CONVERSION MANA

NAME: MDCVT SETUP

MENU TEXT: Conversion Setup

TYPE: run routine

CREATOR: ACKERMAN,NIEN-CHIN

PACKAGE: CLINICAL PROCEDURES

X ACTION PRESENT: YES

DESCRIPTION: This option will bring up a setup screen for the site to setup the Medicine Report Conversion parameter setup. This parameter setup allows the site to control which Medicine reports will be converted and which CP Definition and TIU title to link to.

EXIT ACTION: K DDSFILE,DR,DA

ROUTINE: SETUP^MDCVT

UPPERCASE MENU TEXT: CONVERSION SETUP

NAME: MDCVT RUN

MENU TEXT: Run the Conversion Process

TYPE: run routine

CREATOR: ACKERMAN,NIEN-CHIN

PACKAGE: CLINICAL PROCEDURES

DESCRIPTION: This option will start the Medicine Report conversion to

¹ Patch MD*1.0*5 August 2006 Patch 5 menu options added.

Clinical Procedures. This option will only convert reports for procedures that have the "CONVERT Y/N" field set to "Yes" under the MEDICINE FILE PARAMETERS in the CP CONVERSION file (#703.9).

ROUTINE: EN^MDCVT

UPPERCASE MENU TEXT: RUN THE CONVERSION PROCESS

NAME: **MDCVT SUMMARY**

MENU TEXT: Summary of Conversion

Process

TYPE: print

CREATOR: ACKERMAN,NIEN-CHIN

PACKAGE: CLINICAL PROCEDURES

DESCRIPTION: This option will generate a Medicine Report Conversion report. This report consists of a listing of all Medicine records that were processed in the conversion in variable pointer format and the status of the conversion whether the record was converted, skipped, or errored. If the record was converted, the total number of lines and bytes that the record was converted to in a TIU document will be displayed. If the record errored, the reason why

it errored will be displayed. If the record was skipped, the reason why it was skipped will be displayed.

DIC {DIP}: MDD(703.9,

L.: 0

FLDS: [MD CONVERSION SUMMARY]

BY: [MD CONVERSION SUMMARY]

UPPERCASE MENU TEXT: SUMMARY OF CONVERSION PROCESS

NAME: **MDCVT DISK SPACE**

MENU TEXT: Disk Space Requirements

TYPE: run routine

CREATOR: ACKERMAN,NIEN-CHIN

PACKAGE: CLINICAL PROCEDURES

DESCRIPTION: This option will generate a summary of the Medicine report conversion. This summary consists of a list of the files converted to Clinical Procedures, the count of records converted, the total lines and Bytes

the records were converted in each file.

ROUTINE: SUMMARY^MDCVT

UPPERCASE MENU TEXT: DISK SPACE REQUIREMENTS

NAME: **MDCVT LIST OF TIU TITLES**

MENU TEXT: List of TIU Titles Needed

TYPE: run routine

CREATOR: ACKERMAN,NIEN-CHIN

PACKAGE: CLINICAL PROCEDURES

DESCRIPTION: This option will allow the user to generate a list of Medicine procedures and the TIU titles needed to be created for the procedures that will be used for the Medicine report conversion. The PRINT NAME of the procedures in the PROCEDURE/SUBSPECIALTY file (#697.2) will be used in the display. This list will list the procedures and titles for a Medicine Package

Procedure, if the "Convert Y/N" parameter is set to "Yes" and the "Use TIU Note Title" parameter is blank in the Conversion Setup option.

ROUTINE: DISP^MDSTATU

UPPERCASE MENU TEXT: LIST OF TIU TITLES NEEDED

NAME: **MDCVT TOTALS**

MENU TEXT: Conversion Totals By Status

TYPE: run routine

CREATOR: ACKERMAN,NIEN-CHIN

PACKAGE: CLINICAL PROCEDURES

DESCRIPTION: This option will verify that the Medicine reports conversion is

complete and are in appropriate statuses.

ROUTINE: TOTALS^MDCVT

Exported Options

UPPERCASE MENU TEXT: CONVERSION TOTALS BY STATUS

NAME: **MDCVT ERROR LOG** MENU TEXT: Error Log
TYPE: print CREATOR: ACKERMAN,NIEN-CHIN
PACKAGE: CLINICAL PROCEDURES
DESCRIPTION: This option generates a log of all the errors that occurred with each Medicine report during the conversion. The listing consists of the CONVERSION ID and ERROR MESSAGE. The CONVERSION ID consists of the record # concatenated with a ";" and the global location (e.g., "345;MCAR(699,)").
DIC {DIP}: MDD(703.9, L.: 0
FLDS: [MD CONVERSION ERRORS] BY: [MD CONVERSION ERRORS]
UPPERCASE MENU TEXT: ERROR LOG

NAME: **MDCVT CONVERSION LOCKOUT** MENU TEXT: Conversion Lockout
TYPE: run routine CREATOR: ACKERMAN,NIEN-CHIN
PACKAGE: CLINICAL PROCEDURES
DESCRIPTION: This option will let the user place a specialty/procedure or ALL specialty/procedures Enter/Edit and Report options 'OUT OF SERVICE' in the
Medicine package. It will also set Kernel site parameter MD MEDICINE CONVERTED to "YES" when all specialties/procedures enter/edit and report options are disabled or when the user indicated that all Medicine reports has been converted.
ROUTINE: LOCKOUT^MDCVT LOCKOUT UPPERCASE MENU TEXT: CONVERSION

NAME: **MDCVT BUILD CONVERSION LIST** MENU TEXT: Build Conversion List
TYPE: action CREATOR: ACKERMAN,NIEN-CHIN
PACKAGE: CLINICAL PROCEDURES E ACTION PRESENT: YES
X ACTION PRESENT: YES
DESCRIPTION: The user will need to run this option before using the [MDCVT RUN], Run the Conversion Process, option. This option will let the user build
the conversion list of the Medicine file records for the CP CONVERSION file (#703.9). It will populate the CONVERSION LOG sub-file (#703.92) with all entries in the "AC" cross reference in the MEDICAL PATIENT file (#690) and set
the STATUS field as "Ready to Convert" for each entry. This option can be queued. Once the conversion list is built, this option can also be used to add new additional entries in the Medicine file into the conversion list. This option will not overwrite the existing entries in the CONVERSION LOG but add to the list.
EXIT ACTION: K MDS ENTRY ACTION: S MDS=\$\$BLD^MDCVT1()
UPPERCASE MENU TEXT: BUILD CONVERSION LIST

¹NAME: **MD PROCESS NOSHOW/CANCEL**
MENU TEXT: Process No Show/Cancel Studies
TYPE: run routine CREATOR: ACKERMAN,NIEN-CHIN
PACKAGE: CLINICAL PROCEDURES
DESCRIPTION: This option is tasked to run daily. It will check for any appointment that is No Show or Cancelled for CP studies in the "Pending Instrument Data" status.
ROUTINE: EN1^MDWCAN
UPPERCASE MENU TEXT: PROCESS NO SHOW/CANCEL STUDIES

¹ Patch MD*1.0*11 June 2009 Add new exported option.