

ClO Terminology Dictionary

and

Clinical Data Model

Date/Time Generated: JUL 25,2014@10:08

M Environment: DEV,DEV

Abbreviation..... AB BOW S
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Abdominal Bowel Sounds
Display Name..... Abdomen Bowel Sounds
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 6 entries

- 1. NORMAL (n) {553C9A97-A3D5-C52F-0E4E-C7D168C99A2B}
2. HYPOACTIVE (HYPOACTIVE) {2E0C516D-3858-4A1F-A2F3-BF0AB9E3A7FC}
3. HYPERACTIVE (HYPERACTIV) {1A548FFE-8645-4C0D-899D-5D91430CC1B1}
4. FLATULENCE (FLATULENCE) {F6ECE7D8-D9C8-485C-8FAC-03DA9D7E0B66}
5. ABSENT (a) {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0}
6. BRUIT (BRUIT) {2627E78A-B8B3-44EE-BC6F-4DAB16394F8C}

LOCATION: 5 entries

- 1. ABDOMEN LUQ (ABD LUQ) {3B31E1EC-CE74-4C76-BB61-10C678431FCA}
2. ABDOMEN RUQ (ABD RUQ) {F29C9FA1-717A-49B7-8000-95070E9AACA0}
3. ABDOMEN RLQ (ABD RLQ) {2C48394F-6877-4738-856C-E933644D8067}
4. ABDOMEN LLQ (ABD LLQ) {39831F59-D7B0-4156-AC57-4C622A93CA58}
5. GENERAL UMBILICUS (GEN UMBILI) {34CEEF50-14E8-44B2-9FEA-3BAB36D08DFA}

UNITS:



=====

TERM: ABDOMINAL GIRTH  
ID: {F70E6642-2719-22BE-BE87-DEF0A884F177}

=====

Page: 3

Abbreviation..... AG  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Girth of patient abdomen  
Display Name..... Abdominal Girth  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid..... 4536404

UNITS:

-----

CENTIMETER            0 > 76 \* 152 < 300      Dec Prec: 0

Abbreviation..... Wan Pre Me
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description.....
Display Name..... Activities of Daily Living: Wandering
Prevention Measures
Help Text..... For patients at risk for wandering
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Void.....

DATA VALUE: 6 entries

- 10. ASSIGNED ONE TO ONE SITTER (One to One)
{07D48368-101F-4974-9EC5-7CB182A39281}
20. NO WANDERING PREVENTION MEASURES APPLIED (No Wan App)
{6967F640-223E-46CE-82F9-0E4D87745D94}
30. NO WANDERING PREVENTION MEASURES NEEDED (No Wan Nee)
{8D39BA2B-24C6-4355-9AC2-4AAAD2E66D9B}
40. PATIENT LOCATOR DEVICE IN PLACE (Loc Device)
{CAF7B45E-0143-472D-BABF-9A92857419E9}
50. SECURE TRACK IN PLACE (Secure Trk)
{33F72416-BAD4-44AD-825E-93BCBFFB9670}
60. WANDER GUARD IN PLACE (Wander Grd)
{B49D735F-A6ED-4E8F-9BD0-F8514245B262}

UNITS:
-----

Abbreviation..... ADL ADL
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Activities of daily living assessment
Display Name..... Activities of Daily Living
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 25 entries

- 1. COMPLETE BATH (COMP BATH) {A7CB1F56-F88E-403B-AABC-105BE129FBB0}
2. PARTIAL BATH (PART. BATH) {60B9C7AD-9E6E-4412-9BB3-57AFF956C1FC}
3. SELF AM CARE BY PATIENT (SELF AM CR)
{C48A46C1-0E68-479C-B803-08C80617ADFF}
4. TURN AND POSITION (TURN POSIT)
{44A4F470-1760-4073-AB3A-D59CA9EFE7F3}
5. UP TO CHAIR (UP TO CHAI) {42E598E5-82EC-428B-B5FA-1129E3095C78}
6. MEAL (MEAL) {85D38DCD-6A4A-4908-ABE2-4E92ABEF1321}
7. AMBULATED (AMBULATED) {017D99FE-91D5-49F1-AC15-50620879C466}
8. SNACK (SNACK) {60B77974-4102-4511-A105-E24850B924D8}
9. SELF CARE OF BEDSIDE ENVIRONMENT (S CARE BED)
{A1BEC7F4-D404-4824-8BB7-108E0A6C43A2}
10. SHOWER (SHOWER) {C259929A-44DB-48FD-A2EC-BC6ECF9D0476}
11. SELF SHOWER (SELF SHOVE) {C35CCAE8-27D9-46F2-8971-2C4CF4924CCB}
12. SHAMPOO (SHAMPOO) {95CCB23C-E614-40CD-84F0-EE4EFBC95271}
13. DRY SHAMPOO (DRY SHAMP) {8910A186-2590-4352-A884-7D692A9D2B3C}
14. NAIL CARE (NAIL CARE) {B783A171-1F0A-4135-8535-85776B81AB68}
15. TOENAIL CARE (TOEN CARE) {6CF3ABCE-FF74-43E2-81E5-D23339017CA9}
16. FOOT CARE (FOOT CARE) {332E8338-6DAE-4BC9-B195-63DAD57937D0}
17. SELF FOOT CARE (SELF FOOT)
{F2DCAF0E-389D-4987-8C73-66E35266CC78}
18. TURNING SELF INDEPENDENTLY (TURN SELF)
{DD804CEB-F7CA-4D42-BF6A-447E38B55E98}
19. SHAVE (SHAVE) {A5CB33DF-C33B-4DD3-B00D-78D78F7966C7}
20. SELF SHAVE (SELF SHAVE) {90F1760C-4C12-4976-ADCB-DBF456A34023}
21. FACE WASH (FACE WASH) {E1162A8B-D26D-4EB3-988F-50583B1849EA}
22. SELF FACE WASH (S FACE WAS)
{77457D46-5676-46E6-B38F-FE3AF70B5CBC}
23. PERINEAL CARE (PERIN CARE)
{D88736E4-E0EB-4FF7-9199-288A4D051A20}
24. SELF TRACH CARE (S TRACH CA)
{70BE79EC-107F-4CD3-B964-7919B3DDE222}
25. CLOTHED SELF (CLOTH S) {33802A30-FAA7-426B-ADAC-5297818DA7DE}

UNITS:

Abbreviation..... ADL AMB AS  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Ambulation Assistance  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 8 entries

- 
1. INDEPENDENT (INDEPENDEN) {11BC0CF7-1C29-49C7-B7C8-2BCF8FF82250}
  2. WITHOUT ASSISTANCE (WO ASSIS)  
{41CE081D-CB0E-49AF-936A-4FD78E61681B}
  3. ONE PERSON ASSIST (1 PERS AS)  
{A3D14AA8-D658-412F-BBE8-D8AD195ADA3D}
  4. TWO PERSON ASSIST (2 PERS AST)  
{0547070B-C9BC-4406-BDC0-BF9C6C23D41A}
  5. WALKER (WALKER) {91247001-5E89-46A0-A841-B4FAE9226E63}
  6. CANE (CANE) {562B2E0B-0298-44C6-B33B-BC57DE6CBC3C}
  7. CRUTCHES (CRUTCHES) {DDA81269-B25C-4F7D-A059-BB1602B2B179}
  8. GAIT BELT (GAIT BELT) {AD9AB19B-92B6-40BE-B906-FBFC46C39E0D}

UNITS:  
-----

=====

TERM: ADL - AMBULATION DISTANCE  
ID: {49097AB5-ACC4-4FC5-91FB-B62DFF5EBA51}

=====

Page: 7

Abbreviation..... ADL AM D  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Ambulation Distance  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

FOOT                    0 > 0 \* 9999 < 9999      Dec Prec: 0



=====

TERM: ADL - BED MOBILITY  
ID: {ADB040D9-C136-4AF6-809D-CDD868DA3DCA}

=====

Page: 8

Abbreviation..... ADL BED MO  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Bed Mobility  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 9 entries

- 
1. INDEPENDENT (INDEPENDEN) {11BC0CF7-1C29-49C7-B7C8-2BCF8FF82250}
  2. SUPERVISION/SET-UP HELP (SUPERVISED)  
{596D028B-B60D-4244-A1B0-882177C8D25F}
  3. LIMITED ASSISTANCE (LIM ASSIST)  
{9A8C1310-C840-431E-BA21-065B7EA7B768}
  4. EXTENSIVE ASSISTANCE (EXT ASSIST)  
{C402B659-E880-4070-8345-13E5508B4D29}
  5. TOTAL DEPENDENCE (TOT DEPEND)  
{38702CDD-E440-4EFA-9583-36EEE719F7C3}
  6. ONE PERSON ASSIST (1 PERS AS)  
{A3D14AA8-D658-412F-BBE8-D8AD195ADA3D}
  7. TWILL TAPE (Twill Tape) {C52BE250-732D-03BB-B8F9-72BE4D546807}
  8. TWO PERSON ASSIST (2 PERS AST)  
{0547070B-C9BC-4406-BDC0-BF9C6C23D41A}
  9. NOT DONE (NOT DONE) {C3FE84AC-EEA1-4153-AE6C-69755C2E6BB1}

UNITS:  
-----

Abbreviation..... ADL DRESS  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Dressing  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 8 entries

- 
1. INDEPENDENT (INDEPENDEN) {11BC0CF7-1C29-49C7-B7C8-2BCF8FF82250}
  2. SUPERVISION/SET-UP HELP (SUPERVISED)  
{596D028B-B60D-4244-A1B0-882177C8D25F}
  3. LIMITED ASSISTANCE (LIM ASSIST)  
{9A8C1310-C840-431E-BA21-065B7EA7B768}
  4. EXTENSIVE ASSISTANCE (EXT ASSIST)  
{C402B659-E880-4070-8345-13E5508B4D29}
  5. TOTAL DEPENDENCE (TOT DEPEND)  
{38702CDD-E440-4EFA-9583-36EEE719F7C3}
  6. ONE PERSON ASSIST (1 PERS AS)  
{A3D14AA8-D658-412F-BBE8-D8AD195ADA3D}
  7. TWO PERSON ASSIST (2 PERS AST)  
{0547070B-C9BC-4406-BDC0-BF9C6C23D41A}
  8. NOT DONE (NOT DONE) {C3FE84AC-EEA1-4153-AE6C-69755C2E6BB1}

UNITS:  
-----

Abbreviation..... ADL EATING  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Eating  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 9 entries

-----

1. INDEPENDENT (INDEPENDEN) {11BC0CF7-1C29-49C7-B7C8-2BCF8FF82250}
2. SUPERVISION/SET-UP HELP (SUPERVISED)  
{596D028B-B60D-4244-A1B0-882177C8D25F}
3. SET UP HELP ONLY (SETUP ONLY)  
{DACC77AB-5AE7-4F1B-A3E1-C347D3DD7219}
4. LIMITED ASSISTANCE (LIM ASSIST)  
{9A8C1310-C840-431E-BA21-065B7EA7B768}
5. EXTENSIVE ASSISTANCE (EXT ASSIST)  
{C402B659-E880-4070-8345-13E5508B4D29}
6. TOTAL DEPENDENCE (TOT DEPEND)  
{38702CDD-E440-4EFA-9583-36EEE719F7C3}
7. ONE PERSON ASSIST (1 PERS AS)  
{A3D14AA8-D658-412F-BBE8-D8AD195ADA3D}
8. TWO PERSON ASSIST (2 PERS AST)  
{0547070B-C9BC-4406-BDC0-BF9C6C23D41A}
9. NOT DONE (NOT DONE) {C3FE84AC-EEA1-4153-AE6C-69755C2E6BB1}

UNITS:

-----

=====

TERM: ADL - EMOTIONAL SUPPORT  
ID: {8E6573BB-DED8-47D3-AA52-C56404C5D1D9}

=====

Page: 11

Abbreviation..... EMO SUPPT  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Emotional Support  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. PATIENT (PATIENT) {7881F205-BCD1-4948-AE66-AD572F5A02F2}
  2. FAMILY (FAMILY) {1233C4BB-5AC5-4883-AF2A-F5C93CD8814F}
  3. FRIEND (FRIEND) {03902E48-AF86-432A-8372-37796A55F625}
  4. SIGNIFICANT OTHER (SO) {082C42E8-9E8A-4A96-B3AD-8B0E0A62AF00}
  5. CAREGIVER (CAREGIVER) {F9B8BABE-2F05-437D-85C8-585FC1FFC07F}

UNITS:  
-----

Abbreviation..... ADL OTHER
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Other Activity
Help Text..... Other ADL Activities
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 14 entries

- 1. BEDREST (BEDREST) {8CD28057-132E-4CEE-9EA1-F03C8AF50CE8}
2. DANGLED (DANGLED) {40665750-112F-401D-AEFA-8E19F45EC0C9}
3. CHAIR (CHAIR) {5B8A5DBD-7462-4653-8D1F-CA8B1226DE38}
4. COMMUNE (COMMUNE) {75C896CA-F105-4EF0-B641-3F2B1F984F1B}
5. PASSIVE RANGE OF MOTION (PASS ROM) {DD43B9EF-E18F-4203-8564-89EBBAD963A4}
6. ACTIVE RANGE OF MOTION (ACTIVE ROM) {F55E68F3-6788-47AE-B962-119413B4A089}
7. GERI-CHAIR (g) {167F9786-8F24-054F-020A-658A6EE5B0B4}
8. GURNEY CHAIR (GURN CHAIR) {B6EA0944-8414-42FD-88C0-50155DE82896}
9. BATHROOM PRIVILEGES (BATHR PRIV) {29C1B6DD-3308-4DED-9970-1322A97C633B}
10. TRANSFER WITH ASSIST (TRANS ASST) {62A5BFF5-9C7C-4412-83A8-F872095DE005}
11. FULL SPINAL PRECAUTIONS (FULL S P) {42ADB74D-12CC-4111-B181-626C0F01BC7C}
12. THORACIC SPINAL PRECAUTIONS (THOR S P) {D55AEF23-7F5E-45DD-A178-F217E6AF3873}
13. CERVICAL SPINAL PRECAUTIONS (CERV S P) {A15F7127-2367-464A-AB7D-1E6A7C09AA21}
14. LUMBAR SPINAL PRECAUTIONS (LUMB S P) {5AFBBBAB-F919-436A-9C5A-71DCF5D79D18}

UNITS:
-----

Abbreviation..... ADL REPOSI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Reposition  
Display Name..... Reposition  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 9 entries

- 
1. LEFT (LEFT) {81063CEE-2E83-4C95-AAEE-76B7DCB828D4}
  2. RIGHT (RIGHT) {AB74C579-DF14-4D7B-A9EF-5F119089F976}
  3. SUPINE (SUPINE) {A18F9534-17AF-47D7-8D14-2AA457A889BD}
  4. PRONE (PRONE) {41180E09-D9D6-4903-BAE9-28CD5AAC4991}
  5. SELF CARE (SELF CARE) {0E068188-75B8-46C5-8947-A6CC5C7849B5}
  6. ASSISTED (ASSISTED) {65069804-39BF-41ED-9EA3-BF2C41389F10}
  7. CARDIAC REPOSITION (CARD REPOS)  
{E133FE5A-20E1-418F-BD90-D2DFE14908C0}
  8. FOWLER'S (FOWLER'S) {8E691813-A9AC-48BE-9401-528A7D762A32}
  9. HIGH FOWLER'S (HIGH FOWLE)  
{EA23C94D-EE91-4552-9F74-E051FAF3E194}

UNITS:  
-----

=====

TERM: ADL - VISITORS

Page: 14

ID: {C0C8630F-6535-4C89-9556-D97CBCB41952}

=====

Abbreviation..... ADL VISITO  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Visitors  
Display Name..... Visitors  
Help Text..... Has the patient had any visitors today.  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

-----

1. NO (n)	{FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y)	{483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

-----

Abbreviation..... ADL ORAL C
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Oral Care
Display Name..... Oral Care
Help Text..... Oral Care
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 12 entries

- 1. DONE (DONE) {60412C8D-132F-46EA-8DDF-FCC282C4CC8B}
2. ASSISTED (ASSISTED) {65069804-39BF-41ED-9EA3-BF2C41389F10}
3. SELF CARE (SELF CARE) {0E068188-75B8-46C5-8947-A6CC5C7849B5}
4. OFFERED (OFFERED) {715091C2-9668-4678-900E-6D55D3335A30}
5. REFUSED (REFUSED) {D5B3E135-4D68-4528-BBC9-6E8A7A997104}
6. NOT DONE (NOT DONE) {C3FE84AC-EEA1-4153-AE6C-69755C2E6BB1}
7. SUCTIONED ONLY (SUCTION ONLY) {7A8B23FD-69C5-4931-A6D4-E92F8FD70817}
8. LIP SALVE (LIP SALVE) {612F78AB-7386-4176-9896-02DE5A8E836D}
9. TEETH BRUSHED (TEETH BRUSH) {20841794-6E6D-4262-B1FE-9E3F20F8CC49}
10. MOUTHWASH USED (MOUTHWASH) {9F7EF355-12A4-4EBA-92ED-0F38273E411B}
11. ETT REPOSITIONED (ETT REPOSITION) {95B35EDB-76B6-433E-B526-01C543A7CCEA}
12. ETT RETAPED/SECURED (ETT RETAPE) {6D8A2EDA-CAC9-44D2-BDA9-3B7F20908F1E}

UNITS:



=====

TERM: ADI - SPIRITUAL CARE  
ID: {42608477-6B3D-4036-BE05-F0DA3ED218B9}

=====

Page: 16

Abbreviation..... SPIRIT CA  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Spiritual care  
Help Text..... Has Patient Had Spiritual Care  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----



=====

TERM: AIRWAY ARTIFICIAL CUFF STATUS  
ID: {18958F14-76F1-B4F3-DB01-F5FB4C8A70BC}

=====

Page: 18

Abbreviation..... ART AIR C  
Active..... Yes  
Boolean Value False..... Deflated  
Boolean Value True..... Inflated  
Data Type..... Boolean  
Description..... AIRWAY CUFF INFLATED OR DEFLATED  
Display Name..... Artificial Airway Cuff Status  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

=====

TERM: AIRWAY ARTIFICIAL POSITION AT LIP/NOSE  
ID: {A9032CCC-5479-87C2-27F8-46D8B4F68272}

=====

Page: 19

Abbreviation..... ART AIR PO  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Artificial Airway Position at Lip/Nose  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
CENTIMETER            11.9 > 12.0 \* 28.5 < 28.6 Dec Prec: 1

Abbreviation..... ART AIR F  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Artifical Airway Size -French  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

FRENCH                    3.9 > 4.0 \* 14.0 < 14.1    Dec Prec: 1

Abbreviation..... ART AIR mm  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Artificial Airway Size - mm  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 13 entries

- 
- |               |  |
|---------------|--|
| 1. 4 (4)      | {D2280C7E-DE50-5A8D-C3FD-5862B69C0E1A} |
| 2. 4.5 (4.5)  | {EC6C0C72-8C49-4535-92B2-7E543007EE32} |
| 3. 5 (5)      | {188EDCE7-0726-4234-D8D2-92C8052CD221} |
| 4. 5.5 (5.5)  | {4878F059-EC9D-44FC-9A61-734C291530E7} |
| 5. 6 (6)      | {640CF2F9-3BBF-AB13-185B-CD02A609FB5A} |
| 6. 6.5 (6.5)  | {CAFCBDB1-8533-4C10-8DD1-2BD07A4CF6D3} |
| 7. 7 (7)      | {7A010324-D2C2-5D4B-848A-7B8878666D20} |
| 8. 7.5 (7.5)  | {7517EE9C-7C8B-4AF0-B60A-88F4E24BFA5E} |
| 9. 8 (8)      | {518E7CFC-102E-351E-C7AF-FB575E3C8F4D} |
| 10. 8.5 (8.5) | {96723567-151D-4B5C-9C65-B3EB6D584272} |
| 11. 9 (9)     | {58FF2FD4-569F-F454-9F5B-FE6F9E34D8CC} |
| 12. 9.5 (9.5) | {252555A1-ADB2-4948-8764-C952F2587DC5} |
| 13. 10 (10)   | {0C827160-18DB-816F-7FAA-A3C175C6D33B} |

UNITS:

-----

Abbreviation..... ART AIRW T  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... TYPE OF ARTIFICIAL AIRWAY  
Display Name..... Artificial Airway Type  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 8 entries

- 
1. NASAL ENDOTRACHEAL (Nasal Endo)  
{C631352E-CAB5-0409-6133-E5D11ECC3293}
  2. ORAL ENDOTRACHEAL (ORAL ENDO)  
{466CAD16-D7F3-9D73-6500-959C4DA50D9A}
  3. TRACHEOSTOMY (TRACHEOSTO) {2B0B7599-847E-852F-D16F-091CAAA4B205}
  4. LARYNGECTOMY (LARYNGECTO) {225CB76A-BAF5-827E-729B-C3D60AD5B50A}
  5. NASAL TRUMPET (NASAL TRUM)  
{2FFB8510-793C-D6A4-838E-20FF7CD26533}
  6. BITE BLOCK (BITE BLOCK) {E7CC51F3-08A8-5FEE-38B2-16B74E1F542E}
  7. ORAL AIRWAY (ORAL AIRWA) {B0CF4270-A68D-1322-0621-1DFF8D677B40}
  8. CASS ENDOTRACHEAL TUBE (CASS)  
{02BF805C-DF35-4DAB-A99F-2ED673174E49}

UNITS:  
-----

Abbreviation..... ALARM PARA  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Alarm Parameters Verified  
Display Name..... Alarm Parameters Verified  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----



Abbreviation..... ALARMS  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Alarms  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 7 entries

- 
1. CENTRAL ECG (cecg) {26D36175-E4E3-C364-3018-1B0DB4DC4793}
  2. BED EXIT ALARM (BED EXIT A)  
{AC53759E-FFAE-5E0F-2EBA-E56A2712357C}
  3. WHEELCHAIR EXIT ALARM (WHEEL EXIT)  
{7D06564E-ECDE-4A11-BEFE-92AB7E448AA3}
  4. IV PUMP (IV PUMP) {EFB92C16-CFA3-18BB-9D54-590089DA3962}
  5. VENTILATOR (v) {A5A647ED-A56D-5012-AD5B-788B3DCB3498}
  6. UROMETER (u) {13074A3A-19D7-300C-7BFF-8DC18DDE2951}
  7. IABP (IABP) {763FE975-F5E5-80A6-AB16-E903A829E665}

UNITS:  
-----











=====

TERM: ALDRETE - TOTAL

Page: 30

ID: {A947CF6C-094E-49B9-A856-E4103162D875}

=====

Abbreviation..... ALD TOTAL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Modified Aldrete Scoring System Total  
Display Name..... Aldrete - Total Score  
Help Text..... Total value of Aldrete scores  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MODIFIED ALDRETE SCORE     -1 > 0 \* 10 < 11     Dec Prec: 0

Abbreviation..... Anes Plan
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description.....
Display Name..... Anesthesia Planned
Help Text.....
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 9 entries

- 1. LOCAL (Local) {64A0FD20-4584-446C-8A2A-57F2FEF3DD10}
2. REGIONAL BLOCK (Region Blk) {3F829DFC-D3F6-4D08-AACD-A9A29ED82F19}
3. BIER BLOCK (Bier Blk) {241C28E7-6040-4ABC-B76C-4255452C66FD}
4. GENERAL (General) {6ED377F7-55E4-4743-B781-124267567200}
5. IV SEDATION (IV Sed) {789B61EB-5E83-4024-850D-D92370C4EC9A}
6. SPINAL (Spinal) {AC02181C-0A2C-4F7A-B983-B32EF1BCD150}
7. EPIDURAL (Epidural) {63317A63-16FF-4D55-A034-420C65B17AEB}
8. MODERATE SEDATION (Mod Sed) {C874EDC7-24BF-4F39-8FDA-E8D1193A8428}
9. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}

UNITS:



Abbreviation..... APPEARAN  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Appearance  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
1. NEAT (NEAT) {8B01F7B8-B626-4143-A092-1E732EEB4A05}
  2. APPROPRIATE (a) {1CFB285F-34DF-73A2-B89C-F6EE1F40E7E9}
  3. INAPPROPRIATE (Inappropri)  
{7E3B54C2-8A49-28EC-055E-FEA87599421C}
  4. DISHEVELED (DISHEVEL) {0455549E-5CCB-470D-9D7F-F3959E6D55FE}

UNITS:

-----

=====

TERM: ARM BAND - ALLERGIES  
ID: {D72DC694-6192-8D57-FD8B-02E49EF45DE5}

=====

Page: 33

Abbreviation..... AB ALLERGY  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Is the arm band for allergies present  
Display Name..... Arm Band - Allergies  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

Abbreviation..... AB BL TR  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Arm Band - Blood Transfusion IV  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

Abbreviation..... AB DNR  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Is patient armband for DNR present  
Display Name..... Arm Band - DNR  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:

-----

=====

TERM: ARM BAND - ELOPMENT  
ID: {950320E8-0C68-29C0-E39A-355CEB0119C6}

=====

Page: 36

Abbreviation..... AB ELOPE  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Arm Band - Elopment  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEA649CD} |

UNITS:  
-----

=====

TERM: ARM BAND - HIGH FALL RISK  
ID: {48BCE5FD-3703-4F3B-619B-8E5FD2E25D2E}

=====

Page: 37

Abbreviation..... AB FALL RI  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Arm Band - High Fall Risk  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEA649CD} |

UNITS:  
-----

=====

TERM: ARM BAND - PERSONAL MEDICAL ALERT  
ID: {4EF7E9C0-4B52-0EA5-B5ED-D29E655089FE}

=====

Page: 38

Abbreviation..... AB PM ALER  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Arm Band - Personal Medical Alert  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

=====

TERM: ARM BAND - WANDERER  
ID: {F0D2FBB6-0843-239A-356E-D903542E9117}

=====

Page: 39

Abbreviation..... AB WANDER  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Arm Band - Wanderer  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----



Abbreviation..... AB PRESENT  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Arm Bands Present  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 9 entries

- 
1. N/A (n) {B3DBE3FD-C5CB-3894-607A-BE13319594D0}
  2. IDENTIFICATION (i) {39274F69-556C-7F05-6D17-CE8A03D1867D}
  3. DO NOT RESUSCITATE (dnr) {0BD8A134-C563-7E00-A6B1-F7C543CEF6FF}
  4. ALLERGY (a) {EBF87F6F-1BDB-6147-C444-2FF787DEAC7E}
  5. HIGH FALL RISK (hfr) {276C15F3-5BED-0FE0-5649-337A23AB1386}
  6. ELOPEMENT (e) {EB030907-FB90-CA9D-CF50-E26C5D9FE671}
  7. PERSONAL MEDICAL ALERT BRACELET (pmab)  
{E63D5FDE-4509-5751-A11C-C023E296858A}
  8. BLOOD TRANSFUSION IV (BI TRAN IV)  
{27A0BE4C-433C-F920-BAC2-066F304BDF59}
  9. WANDERER (w) {C5BD737B-F9A6-5800-B45C-2585507C3EDC}

UNITS:

-----

=====

TERM: ARTERIAL LINE ALARMS ON  
ID: {9267950F-51B3-4C38-B280-C62D7B076F6B}

=====

Page: 41

Abbreviation..... ALINE ALAR  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Arterial Lines Alarms On  
Display Name..... Arterial Line Alarms On  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----



Abbreviation..... ART LIN DC
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Arterial Line Dressing Change
Display Name..... Arterial Line Dressing Change
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

- 1. DRY (d) {D8F4F0FD-1B19-2515-3815-EC88D175B305}
2. INTACT (i) {935C9CAF-6C88-6610-C6E9-4C1B822923AF}
3. CHANGED (CHANGED) {74B14023-D74D-41F6-AA48-B0E5137623CD}
4. REINFORCED (REINFORCED) {EC524C50-0CDE-4AA1-8059-E4323850575E}

UNITS:

Abbreviation..... A LINE DRE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Arterial Line Dressing  
Help Text..... Type of dressing used on arterial line  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. TRANSPARENT (TRANSPAREN) {0E73D24D-0EFC-08F9-196A-D30DD1BB3228}
  2. PRESSURE DRESSING (PRESS DRES)  
{6D9DB9EF-85EB-8321-03FC-879AA4C5615B}
  3. GAUZE AND TAPE (GAU & TAPE)  
{5BC40405-B5C4-47CF-81C3-520732E100AC}

UNITS:  
-----

=====

TERM: ARTERIAL LINE SITE CARE  
ID: {2CCDDAA0-9563-8B09-11C1-51E1E2823A48}

=====

Page: 45

Abbreviation..... A LINE CAR  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Has Arterial Line Site Care Been Done?  
Display Name..... Arterial Line Site Care  
Help Text..... Date arterial line inserted  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

Abbreviation..... ART LI CON
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Arterial Line Site Condition
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 8 entries

- 2. RED (RED) {515D7CFB-3527-30E7-EC4B-B607C272DFA0}
3. SWOLLEN (SWOLLEN) {DCE86B10-C393-8DA9-6518-F2943F0E1A4E}
4. DISCOMFORT (DISCOMFORT) {77D058C9-1EB4-96C0-C9E3-666C840B884C}
50. PATENT (Patent) {FOC38BCB-1D11-436C-8B8C-4297C27DDE04}
60. CLEAN (Clean) {1F187B99-9963-4AD5-98F6-D8D24D2DDD30}
70. BLANCHED (Blanched) {EBB2EA8D-04BC-49E7-812A-62B46F0F4678}
80. OOZING (Oozing) {C23549D7-3238-4523-BA31-AB71609F01C2}
90. OCCLUDED (Occluded) {25509BE3-F946-4AD7-80E8-5CE8A2E61AD6}

LOCATION: 8 entries

- 1. L RADIAL (L RAD) {302F78AD-4067-8B6B-BD15-993EC8561AE6}
2. R RADIAL (R RADIAL) {BA10E6E5-DDF4-345E-68AF-21D32E9D972D}
3. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
4. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
5. L AXILLARY (L AXILLARY) {21AF62CA-BEED-2A96-2CD6-F2350A01995E}
6. R AXILLARY (R AXILLARY) {3044F28B-FED2-4069-7E3B-2A691C5771F9}
7. L BRACHIAL (LBRA) {75F6CE33-1AB3-A1DA-70B8-04825290B577}
8. R BRACHIAL (RBRA) {A88E69A8-F6AE-80B0-2CD4-687B121F9DDE}

UNITS:

=====

TERM: ARTERIAL LINE SIZE  
ID: {9D96EBB7-19E0-4832-B61F-D2A42CE15737}

=====

Page: 47

Abbreviation..... ART LIN SI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Arterial Line in size Gauge  
Display Name..... Arterial Line Size  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

LOCATION: 8 entries

- 
1. L RADIAL (L RAD)            {302F78AD-4067-8B6B-BD15-993EC8561AE6}
  2. R RADIAL (R RADIAL)        {BA10E6E5-DDF4-345E-68AF-21D32E9D972D}
  3. L FEMORAL (LFEM)           {50AF6A08-CE37-6859-F4E6-762660A982F7}
  4. R FEMORAL (R FEMORAL)     {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
  5. L AXILLARY (L AXILLARY)    {21AF62CA-BEED-2A96-2CD6-F2350A01995E}
  6. R AXILLARY (R AXILLARY)    {3044F28B-FED2-4069-7E3B-2A691C5771F9}
  7. L BRACHIAL (LBRA)          {75F6CE33-1AB3-A1DA-70B8-04825290B577}
  8. R BRACHIAL (RBRA)          {A88E69A8-F6AE-80B0-2CD4-687B121F9DDE}

UNITS:

-----

GAUGE                    14 > 18 \* 26 < 30        Dec Prec: 0



=====

TERM: ARTERIAL LINE WAVEFORM  
ID: {1B11B2EF-C905-40F6-9F6F-195FEC0CD8E6}

=====

Page: 48

Abbreviation..... ALINE WAVE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Arterial Line Waveform  
Display Name..... Arterial Line Waveform  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. NORMAL (n) {553C9A97-A3D5-C52F-0E4E-C7D168C99A2B}
  2. DAMPENED (DAMPENED) {77C777B3-7E56-41E3-8E70-DE144E6A3424}
  3. HYPERDYNAMIC (HYPERDYNAM) {818A8A99-04D9-4989-A9BA-ECFACED74605}
  4. POSITIONAL (POSITIONAL) {9A74D16B-5E95-4EA8-88AF-E53C80311A5D}
  5. NOT MONITORED (NOT MONITO)  
{0E6A8A07-48E0-428F-8D86-22572CABCF48}

UNITS:  
-----

=====

TERM: ARTERIAL LINE ZEROED  
ID: {A500DBBF-6E69-4CB3-926F-6B9988BB8B68}

=====

Page: 49

Abbreviation..... ALINE ZERO  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Arterial Line Zeroed  
Display Name..... Arterial Line Zeroed  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

Abbreviation..... ASA Phys  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... ASA Physical Status  
Help Text.....  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 6 entries

- 
- |                          |  |
|--------------------------|--|
| 1. ASA PS 1 (ASA PS 1)   | {BAF1730E-C358-17F1-8DB8-4B6097AFF066} |
| 2. ASA PS 2 (ASA PS 2)   | {070C0F56-3D9C-5430-EF9E-BF3B96DFB595} |
| 3. ASA PS 3 (ASA PS 3)   | {B041CE25-82D0-EDAD-C673-A64BCD6647F4} |
| 4. ASA PS 4 (ASA PS 4)   | {040841C3-D933-E687-C7E0-29F439C08778} |
| 5. ASA PS 5 (ASA PS 5)   | {1F7B73A6-3F8E-E230-1C23-C7DF92032E58} |
| 6. EMERGENCY (Emergency) | {3C7B5E35-63CA-32EE-3D64-AE4A052F55E1} |

UNITS:

-----

Abbreviation..... BABINSKI  
Active..... Yes  
Boolean Value False..... Negative  
Boolean Value True..... Positive  
Data Type..... Boolean  
Description..... Babinski Reflex Test  
Display Name..... Babinski  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
  - 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEE649CD}

UNITS:  
-----

Abbreviation..... Barriers  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... Barriers to Learning  
Help Text.....  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. PHYSICAL IMPAIRMENT (Phys Impai)  
    {4D744081-134E-420D-9712-AF114977E490}
  2. COGNITIVE IMPAIRMENT (Cog Impair)  
    {6D051F44-2124-4914-A415-B98B2AC516E9}
  3. OTHER (OTHER)            {19094E5B-AEC3-44B0-9BD5-44E79E7EB105}

UNITS:  
-----

=====

TERM: BED POSITION

Page: 53

ID: {D65573FE-6CAB-3A81-4572-ACC9488D8BAF}

=====

Abbreviation..... BEDPOS  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Bed Position  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. HIGH (h) {D88CB5EC-2AD4-2F80-2D43-40ECB0B7D390}
  2. LOW (l) {B6CE52A4-C1DD-FDE3-F87E-8E0C98FCC4AC}
  3. SEMIFOWLERS HIGH (sh) {1BBE7175-2BD2-BFB5-1E1B-E02443433072}
  4. SEMIFOWLERS LOW (sl) {A265AF02-53FC-B182-8AE3-A06E95718D7E}
  5. CARDIAC CHAIR (cc) {D4A2CFC1-D3A4-696B-4062-AFB2DAFFD76A}

UNITS:

-----

Abbreviation..... BEHAVIOR  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Is Patients Behavior Appropriate?  
Display Name..... Behavior  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
1. APPROPRIATE (a) {1CFB285F-34DF-73A2-B89C-F6EE1F40E7E9}
  2. INAPPROPRIATE (Inappropri) {7E3B54C2-8A49-28EC-055E-FEA87599421C}
  3. RESTING (Resting) {D750373F-12E2-4D95-AF0D-F7731905F723}
  4. APPEARS ASLEEP (App Asleep) {8D1436B8-ADB0-4092-A0CE-383D74131AA8}

UNITS:

-----

Abbreviation..... CUR LOC  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Behavior Current Location  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 10 entries

- 
1. PATIENT'S ROOM (pr) {0F13E11D-9F78-14A3-0568-AF26393F39F1}
  2. BATHROOM (b) {D20D91A7-E9F5-D000-5E0B-772BA9F5F72E}
  3. HALL (h) {BFCEA10F-2663-199E-88D8-E263122D7D26}
  4. DAY ROOM (dr) {269125EF-FA82-B17A-6B56-0CEAB4788899}
  5. QUIET ROOM (qr) {CCD96D41-C273-F6F3-DA18-8818B1CC9980}
  6. GROUP ROOM (gr) {D41C3D9F-8F10-1A0F-D571-20DE0EA36328}
  7. SMOKING PORCH (sp) {5AD3F05B-2A5B-E2B4-FD47-66712CFA44E6}
  8. SMOKING AREA (sa) {8AFDBCB0-9FC8-523F-5D07-4D8CCBA7036F}
  9. MEDICATION ROOM (mr) {9571AA9F-F538-55EC-37F6-99FE5677B970}
  10. NURSING STATION (ns) {B6E48F92-B9BA-488F-3A2E-D74568A6B855}

UNITS:

-----



Abbreviation..... BEH OBS  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Behavior  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 15 entries

- 
- 1. COOPERATIVE (c) {C176406E-EBDE-4ADE-53DE-CF191D8AFC85}
  - 2. CALM (c) {E206C3F5-6977-302C-51FB-AF2155B5FFD9}
  - 3. ANXIOUS (a) {381BDB41-7E8A-6FDB-841B-2ECC1889B9EF}
  - 4. HOSTILE (h) {29EF2C04-7DB0-4491-2B5E-D92A6A558721}
  - 5. COMBATIVE (c) {DF65D42D-6CAA-657E-AEE1-C070C1FDF177}
  - 6. INEFFECTIVE COPING (ic) {704B1A55-B639-60F4-E3A5-16F942ABC259}
  - 7. WITHDRAWN (w) {E74B0F04-0E3F-632F-9409-C4953D95B1BA}
  - 8. UNCOOPERATIVE (u) {2D737EFB-F4CA-7AEA-9EB2-6FF133742815}
  - 9. SLEEPS POORLY (sp) {AB9424B4-D0B1-148B-C800-FF535F3FAC84}
  - 10. AGITATED (a) {BAD2CCCE-9083-6F61-1196-C0C9A8DDBFFC}
  - 11. SEDATED (s) {29D9ADD6-2C1F-F6D5-F4D0-B2BB86416997}
  - 12. RESTLESS (r) {8FCCAD26-DBA0-7067-DA72-30D342BFA8E8}
  - 13. ALERT (a) {BC7E2630-179A-5C80-E6C7-598D42F0BFD6}
  - 14. DROWSY (d) {6B350814-A88A-0F05-C382-9F27830FBC09}
  - 15. LETHARGIC (l) {4ECA2E44-093A-9354-848E-CDC72B165072}

UNITS:

-----

Abbreviation..... PRE LVL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Behavior Precaution Level  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
1. CLOSE OBSERVATION (co) {2679BCAC-8753-FE73-56F8-3AB632161AF4}
  2. 1:1 OBSERVATION (1o) {B1A7AF4B-9BD2-5F48-79B3-F35861AA5C45}
  3. PATIENT STAYS WITHIN LINE OF SITE OF NURSING STATION (pswlosons)  
{CE3E3531-7D9B-A2B4-00B3-00EBCDD09D3D}
  4. PATIENT SIGNS SAFETY CONTRACT (PT CONTRAC)  
{C65734EC-2CAF-4547-8FD8-CBCD9C706003}

UNITS:  
-----

Abbreviation..... PRE REA  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Behavior Precaution Level  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. SUICIDAL (s) {B104E9C7-A78E-82DC-49A4-19B94A639AFE}
  2. ESCAPE (e) {965F3201-9D63-66AB-FA25-99479AE7145A}
  3. ELOPEMENT (e) {EB030907-FB90-CA9D-CF50-E26C5D9FE671}
  4. ASSAULTIVE (a) {99ECD8C7-191A-3B9C-2A7A-803A921065FA}
  5. MEDICAL DEVICE SAFETY (mds)  
{1985CABE-4F53-9FCF-F444-F9AC53C04DBD}

UNITS:  
-----





=====

TERM: BLADDER ASSESSMENT  
ID: {BEDA282C-C39A-4134-DA6D-7BFCC567A058}

=====

Page: 61

Abbreviation..... BLAD ASSMT  
Active..... Yes  
Boolean Value False..... Nondistended  
Boolean Value True..... Distended  
Data Type..... Boolean  
Description..... Assessment of bladder  
Display Name..... Bladder Assessment  
Help Text..... Is bladder nondistended or distended?  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----



Abbreviation..... BP
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Complex
Description..... Complex observation of Systolic and Diastolic pressures.
Display Name..... Blood Pressure
Help Text..... Enter in the blood pressure values as ###/###.
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid..... 4500634

LOCATION: 7 entries

- 1. BRACHIAL (Bra) {67426127-09E1-E566-F31E-15714A5B2B7C}
2. FEMORAL (Fem) {E3442C9E-8677-7296-D314-7365319C3716}
3. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
4. L LEG (LL) {B088C809-E528-B785-21C2-14E45E1FD8AF}
5. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
6. R LEG (RL) {C5E40A7B-6FC4-267D-8CE8-3571D21D81FD}
7. RADIAL (RADIAL) {D52C4103-2AB3-A0C7-16E3-CF95A1B3B789}

METHOD: 11 entries

- 1. AFTER EXERCISE (AE) {D129EF4B-3A99-9966-54F1-2C01A36814EE}
2. ARTERIAL LINE (AL) {D4B3F4ED-8B33-431E-4164-BD165E823776}
3. AT REST (AR) {E8C289A6-0C3D-A7B2-BB85-7F1A23F22B50}
4. CUFF-AUTOMATED (CA) {088AAD2B-2E7B-F835-376F-3D4340790976}
5. CUFF-MANUAL (CM) {5F56C163-2B04-112F-55DE-E31374E0FB58}
6. DOPPLER (Dop) {52683C19-F9DC-876B-CA3F-53E7C864CF04}
7. INVASIVE (I) {5DB8746E-5446-F47A-883C-A44B91275A6E}
8. MONITOR (Mo) {63FC24EA-50B8-B302-D948-ABE3517CF6FC}
9. NON-INVASIVE (N-I) {F50ABFED-5FD7-FC95-20FE-1E3038C9F299}
10. PALPATED (Pal) {295D8E72-12A1-08AF-6BE8-0F9CD238F5EF}
11. WITH ACTIVITY (Act) {544A262B-EA34-F82A-46A9-F80740515D65}

POSITION: 6 entries

- 1. LYING (Ly) {DE28D2BF-6F8D-64C3-0092-4CEB5E0902B9}
2. SEMIFOWLERS (Sfl) {C2F6A945-3829-AABA-3D67-54D1829D3EB9}
3. SITTING (Sit) {DD4DCDF1-0A22-2CC5-404F-0E472C73961C}
4. STANDING (St) {0776AE42-C370-B7B3-9B4D-9C914BA8F231}
5. TRENDELENBURG (Tr) {42F3353B-7041-EEE1-A87E-941D3BE1CBCD}
6. SUPINE (SUPINE) {DEA8A744-B445-4A9E-875D-485E060F2ACC}

UNITS:



=====

TERM: BLOOD PRESSURE CUFF SIZE  
ID: {46DCF900-6F2C-46D0-ADB8-2EC17E7659FF}

=====

Page: 64

Abbreviation..... BP CUFF SI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Blood Pressure Cuff Size  
Help Text..... Blood Pressure Cuff Size  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. ADULT (ADULT) {37B557BE-B932-4757-8E08-5EA315E8AD1F}
  2. ADULT LARGE (AD LARGE) {38E64D72-2806-47FA-B82F-51844574F850}
  3. ADULT SMALL (AD SMALL) {0161A4E5-713A-406B-B3FD-9ED13BEB5366}
  4. THIGH (THIGH) {65AE788B-52F5-44D3-B9C1-49DAE46FF5AB}
  5. PEDIATRIC (PEDIATRIC) {06865E22-3ED9-4EDD-9749-6F8641BF29DE}

UNITS:  
-----

=====

TERM: BODY SURFACE AREA  
ID: {5C9A522B-C5ED-4D59-7BF8-3D88ADEBAB07}

=====

Page: 65

Abbreviation..... M2  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Body Surface Area  
Display Name..... BSA  
Help Text..... Body Surface Area  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

SQUARE METERS            0.50 > 1.00 \* 1.75 < 4.00 Dec Prec: 2

=====

TERM: BOWEL MOVEMENT EVENT(S)  
ID: {2EB1DAFB-1E56-40A8-B371-D46EF7AF2435}

=====

Page: 66

Abbreviation..... BM EVENT  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Number of bowel movements observed  
Display Name..... Bowel Movement Event  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
EVENT(S)                    -1 > 0 \* 8 < 30                    Dec Prec: 0



Abbreviation..... BR SC FRIC  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Braden Scale Friction and Shear  
Display Name..... Braden Scale Friction and Shear  
Help Text..... 1. Problem: Requires moderate to  
maximum assistance in moving. 2  
Potential Problem: Moves feebly or  
requires minimum assistance.3 No  
Apparent Problem: Moves in bed and  
in chair independently and has  
sufficient muscle strength to lift up  
completely  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. BRADEN FRICTION AND SHEAR 1 PROBLEM (BS F/S 1)  
{AF8D2952-E0DA-411E-91CC-0114F4EB1C33}
  2. BRADEN SCALE FRICTION AND SHEAR 2 POTENTIAL PROBLEM (BS F/S 2)  
{B56338A9-8956-4D27-8BA0-DA50E07359E0}
  3. BRADEN SCALE FRICTION AND SHEAR 3 NO APPARENT PROBLEM (BS F/S 3)  
{9850E2BF-6600-4E66-8FCD-1517D903BC0F}

UNITS:  
-----

=====

TERM: BRADEN SCALE MOBILITY  
ID: {85F9FB55-26A3-6203-C1AE-6EFED5952B3D}

=====

Page: 69

Abbreviation..... BR SC MOBI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Braden Scale Mobility - ability to  
change and control body position  
Display Name..... Braden Scale Mobility  
Help Text..... 1=Completely Immobile, 2=Very  
Limited, 3=Slightly Limited, 4=No  
Impairment  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
1. BRADEN SCALE MOBILITY 1 COMPLETELY IMMOBILE (BS MOB 1)  
{7607574E-D4C2-4E64-8DC5-ABE7D00917A6}
  2. BRADEN SCALE MOBILITY 2 VERY LIMITED (BS MOB 2)  
{868EFA33-877E-4B91-B4FD-9655F66A7160}
  3. BRADEN SCALE MOBILITY 3 SLIGHTLY LIMITED (BS MOB 3)  
{C54B0A8F-46F1-4DFA-9C24-AB6C2A29CA82}
  4. BRADEN SCALE MOBILITY 4 NO LIMITATION (BS MOB 4)  
{84BFB8D0-D6E4-4384-88EE-28B7C4EDB275}

UNITS:  
-----

Abbreviation..... BR SC MOIS  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Braden Scale Moisture - degree to  
which skin is exposed to moisture  
Display Name..... Braden Scale - Moisture  
Help Text..... 1=Constantly Moist, 2=Very Moist,  
3=Occasionally Moist, 4=Rarely Moist  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
1. BRADEN SCALE MOISTURE 1 CONSTANTLY MOIST (BS MOI 1)  
{3E2569BA-ECF0-4622-8F4A-26378A2CB466}
  2. BRADEN SCALE MOISTURE 2 VERY MOIST (BS MOI 2)  
{2C6D15A3-6166-4D7D-B200-B58EC9C2CABE}
  3. BRADEN SCALE MOISTURE 3 OCCASIONALLY MOIST (BS MOI 3)  
{D7F2714B-B0C7-4000-8588-8582C1B32ECE}
  4. BRADEN SCALE MOISTURE 4 RARELY MOIST (BS MOI 4)  
{745EA3B9-18C9-4CF3-9DB8-9806A02656E8}

UNITS:  
-----







=====

TERM: BRADEN SCALE TOTAL  
ID: {42F65882-FAE1-45E1-84E9-A2A71052BD0C}

=====

Page: 73

Abbreviation..... BRADEN TOT  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Braden Scale Total  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

BRADEN SCALE            5 > 6 \* 23 < 24            Dec Prec: 0

Abbreviation..... BREATH SOU  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... SOUND OBSERVATIONS FROM PATIENT'S LUNGS  
Display Name..... Breath Sounds  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 7 entries

- 
1. RHONCHI (RHONCHI) {850A9D2A-2F7F-A812-3D56-EE0948B0A0CC}
  2. CRACKLES (c) {6B64C108-4424-00D9-525F-D89D549476D6}
  3. CLEAR (c) {8981C365-2B46-704F-C9F4-975DF1C23C2D}
  4. ABSENT (a) {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0}
  5. DIMINISHED (d) {49DE419E-DC3C-F7CD-771D-742031E0DEBC}
  6. INSPIRATORY WHEEZES (INSP WHEEZ)  
{3FA7F54C-DE17-4DFF-B5F6-C2CD05EC0B24}
  7. EXPIRATORY WHEEZES (EXP WHEEZE)  
{CFC103A8-8B02-4942-BBDD-8D24C930883E}

LOCATION: 5 entries

- 
1. LEFT LOWER LOBE (LLL) {14155C13-874D-6E26-01D8-2C03B349F913}
  2. LEFT UPPER LOBE (LUL) {21BC3BA2-B22E-6A23-A377-5C2D5C1575C7}
  3. R LOWER LOBE (RLL) {418D350A-3602-1918-22C8-53D093AF49B0}
  4. R MIDDLE LOBE (RML) {57879159-350C-2055-34D3-C6EFA39CA98F}
  5. R UPPER LOBE (RUL) {098C11B4-B556-FBBE-7E21-AEC9D7E4408A}

UNITS:

-----

=====

TERM: CALL LIGHT WITHIN REACH  
ID: {ABAC4D24-C449-117E-BE4C-BFBC1F599A57}

=====

Page: 75

Abbreviation..... CALL LT  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Call Light Within Reach  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEA649CD} |

UNITS:  
-----

=====

TERM: CAPILLARY REFILL  
ID: {951BBDAA-883A-E4F7-AD05-B086A71A2FEC}

=====

Page: 76

Abbreviation..... CAP REFILL  
Active..... Yes  
Boolean Value False..... Greater than 3 seconds  
Boolean Value True..... Less than 3 seconds  
Data Type..... Boolean  
Description..... RATE OF CAPILLARY REFILL  
Display Name..... Capillary Refill  
Help Text..... RATE OF CAPILLARY REFILL  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

=====

TERM: CARDIAC ABNORMAL BPM  
ID: {1B168030-C39C-8120-915A-36073D044B17}

=====

Page: 77

Abbreviation..... ABN  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Cardiac Abnormal BPM  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

BEATS PER MINUTE      0 > 1 \* 200 < 201      Dec Prec: 0

=====

TERM: CARDIAC ALARM PARAMETERS VERIFIED  
ID: {564B52F5-824E-3DBF-5EDB-493DDFFC30A3}

=====

Page: 78

Abbreviation..... CARD ALARM  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Cardiac Alarm Parameters Verified  
Help Text..... Cardiac Alarm Parameters Verified  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

Abbreviation..... BRUIT  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Bruit  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |                |  |
|----------------|--|
| 1. PRESENT (p) | {384DC20F-1625-CE1E-1C85-A3B3924BB2A0} |
| 2. ABSENT (a)  | {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0} |

LOCATION: 7 entries

- 
- |                                |  |
|--------------------------------|--|
| 1. L FEMORAL (LFEM)            | {50AF6A08-CE37-6859-F4E6-762660A982F7} |
| 2. R FEMORAL (R FEMORAL)       | {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC} |
| 3. L CAROTID (L CAROTID)       | {FCACD1EB-C12B-2515-4A68-4E9E28A72DBC} |
| 4. R CAROTID (R CAROTID)       | {8280ACD2-72D6-202A-89AF-86BCA7736427} |
| 5. R TEMPORAL (R TEMPOR)       | {1F50461D-5ECF-25B2-A29A-E02803CB7814} |
| 6. L TEMPORAL (L TEMPO)        | {375E65CF-A59E-A1E9-EB6D-D53CD818939E} |
| 7. ABDOMINAL AORTA (ABD AORTA) | {C4BD42DD-1D6B-4B2C-B6D5-19E1185E06E3} |

UNITS:

-----



Abbreviation..... CARD E  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Cardiac Ectopy  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 13 entries

- 
- 1. PJC (p) {1D2D1FEA-2006-3C6C-FAD4-72384F2033C9}
  - 2. PAC (p) {848AB18A-598D-8E4E-8B1B-ED7478CE340B}
  - 3. PVC (p) {33FF92F1-84CC-7CFA-41CF-B271A26C59B8}
  - 4. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
  - 6. MULTIFOCAL PVC (Mult PVC) {34AAE3E0-C0AA-4262-90A2-FCAD14E8023A}
  - 7. UNIFOCAL PVC (Uni PVC) {5FE7D7FD-048B-4948-A4DC-C9F453AB31C5}
  - 8. PAIRS (Pairs) {668B0D85-E507-4E12-B07A-A90CB98730DC}
  - 9. RUNS OF VT (Runs of VT) {3BA2FF31-1A56-4B4B-AC13-E8608BC3FCC2}
  - 10. BIGEMINY (Bigeminy) {729FCAD7-EDA1-4CFD-A6CB-E291DDCF0255}
  - 11. TRIGEMINY (Trigeminy) {13D65F9F-B51C-4CEC-8002-037C8323626C}
  - 12. QUADRAGEMINY (Quadrage mi) {899BBEBB-D103-4864-BBE5-C9DF60F1C5BB}
  - 13. JUNCTIONAL ESCAPE BEAT (Junctional)  
{54F8E569-801E-4512-A9BD-C04CDA5CEE7C}
  - 14. REMNANT P-WAVE (Remnant P)  
{3BBAA695-337E-41C4-BB83-8452917C07DF}

UNITS:

-----

=====

TERM: CARDIAC HEART SOUNDS  
ID: {476105AC-C83C-ECEC-6C72-7C7C8EE87853}

=====

Page: 81

Abbreviation..... CARD HS  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Cardiac Heart Sounds  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 7 entries

- 
1. S1 (s) {1AB9DF0A-8238-E8FB-8520-4749E13DDE1D}
  2. S2 (s) {57C043D3-B227-3E7A-C1E6-5129A2800718}
  3. S3 (s) {9374B57F-DE1F-B9C4-26B5-2DD6B194DC26}
  4. S4 (s) {9C6CA302-85B6-8FFC-EAB4-8FED62C535A5}
  5. RUB (r) {0040835D-1CF0-7D55-EEBB-A94A1E0D1206}
  6. CLICK (c) {F1F3A1A8-35E1-61B4-C783-3314C8247380}
  7. MURMUR (MURMUR) {591058E0-4B71-E432-B7D7-EBD42A8C4225}

UNITS:  
-----

=====

TERM: CARDIAC INDEX

Page: 82

ID: {BFC29A9A-C666-68B5-C5CF-36B809CF66C2}

=====

Abbreviation..... CI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Cardiac Index  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----  
1. MONITOR (Mo) {63FC24EA-50B8-B302-D948-ABE3517CF6FC}

UNITS:

-----  
LITERS PER MINUTE PER METER2 (BSA)  
-0.1 > 0.0 \* 8.0 < 8.1 Dec Prec: 1

=====

TERM: CARDIAC OUTPUT  
ID: {0B7C4991-12E5-10DE-F25D-CB81BFA1ADAD}

=====

Page: 83

Abbreviation..... CO  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Cardiac Output  
Display Name..... Cardiac Output  
Help Text..... Cardiac Output  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----

1. MONITOR (Mo) {63FC24EA-50B8-B302-D948-ABE3517CF6FC}

UNITS:

-----

LITERS PER MINUTE 0.00 > 4.00 \* 8.00 < 16.00  
Dec Prec: 2

Abbreviation..... CARDR  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... HEART RHYTHM  
Display Name..... Cardiac Rhythm  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 32 entries

- 
1. SINUS RHYTHM (sr) {A34F80EA-3D6A-076A-DB2A-E4E5C5AACA11}
  2. SINUS BRADYCARDIA (sb) {D6EF4940-63AE-22AF-1A80-D8E6F0F07839}
  3. AGONAL (a) {1621B411-33E3-3A6A-78A6-40B7BDC3E8CB}
  4. ASYSTOLE (a) {6E87ECBC-3160-2CEF-D2FF-EEBD1DAA96F2}
  5. ATRIAL FIB/FLUTTER (af) {9C1BC2F6-3BF0-BEBB-BE87-50E9F7984544}
  6. ATRIAL FIBRILLATION (a fib)  
{7B5EBAF2-55B6-2905-3D89-B28CE6A509CB}
  7. ATRIAL FLUTTER (af) {9CB8D413-6762-D78E-F40E-25DB66CD23B3}
  8. L BUNDLE BRANCH BLOCK (lbbb)  
{C3F2E19E-495B-7BF6-99EF-B8236D6B9C3D}
  9. R BUNDLE BRANCH BLOCK (rbbb)  
{FDE0D174-5220-F069-D284-20399D9489B0}
  10. IDIOVENTRICULAR (i) {8BD08966-F536-8E74-0BDA-A95D0D7F8B6C}
  11. JUNCTIONAL (j) {4A339D84-7051-1258-877D-85AFC1CA20FE}
  12. PACED: AV SEQUENTIAL (pas)  
{3C485674-224C-47DA-510D-92F1A0CA54A1}
  13. PACED: ATRIAL (pa) {88837DD4-098D-B77F-A754-89D5C4F5870F}
  14. PACED: VENTRICULAR (pv) {A63F6E6E-A7C4-0928-81FB-FE2F70F713C6}
  15. PAROXYSMAL ATRIAL TACHYCARDIA (pat)  
{27BC1FAF-A403-6EA8-044C-F1271F12A62A}
  16. SINUS ARREST (sa) {3CD8F7AC-3FE0-F342-089C-787232D2A5C2}
  17. SINUS ARRHYTHMIA (sa) {07A24AC9-B486-A4B8-BADB-04D2729542CE}
  18. SUPRAVENTRICULAR TACHYCARDIA (st)  
{F712412F-D39A-58F7-E4F1-6AC8B77445D7}
  19. VENTRICULAR FIBRILLATION (vf)  
{400B70C2-6E15-458F-70C6-FA2597DE2556}
  20. VENTRICULAR TACHYCARDIA (vt)  
{138258AD-3A2C-2A6C-9DF8-615CCC303DC2}
  21. FIRST DEGREE AV BLOCK (fdab)  
{4FF046C9-B638-F86C-5D51-006B3D5225E3}
  22. SECOND DEGREE AV BLOCK: TYPE 1 (sdabt1)  
{508BF663-96FB-EEBE-4F06-D2F1A02BD099}
  23. SECOND DEGREE AV BLOCK: TYPE 2 (sdabt2)  
{3367754A-08D1-31A9-E850-1B36358E358A}
  24. THIRD DEGREE AV BLOCK (tdab)  
{1AC7E31B-830E-2309-3AF7-8FB03003DD86}
  25. SINUS PAUSE (sp) {7F8F89A9-6B16-B159-55B3-DEFDF5B4B4F6}
  26. SINUS BLOCK (sb) {94D6EFED-3D08-E724-BC40-9C398919A142}
  27. SINUS TACHYCARDIA (SINUS T)  
{0AC41AF3-7656-A35D-E397-21D35405812A}
  28. PAROXYSMAL SUPRAVENTRICULAR TACHYCARDIA (PST)  
{9C910A04-95D2-4EB7-A4CB-B4A74FE7CCEE}
  29. MUTIFOCAL ATRIAL TACHYCARDIA (mat)  
{562C6BBB-C210-4467-93B3-8E50640E7CCE}

- 30. JUNCTIONAL TACHYCARDIA (jt)  
{D6EF01C1-A8F9-40FD-AE24-785041526041}
- 31. TORSADES DE POINTES (Tor de Poi)  
{E936B787-7912-4CAE-A465-AD1A1FEA68CA}
- 32. PULSELESS ELECTRICAL ACTIVITY (Pulseless)  
{8E21E500-78D4-462D-9483-DD81D3887EE5}

UNITS:

-----

Abbreviation..... CAR RH LD  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Cardiac Rhythm Lead  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 12 entries

- 
- |              |  |
|--------------|--|
| 1. I (I)     | {BEEDFF9A-9401-4D2A-8345-126D8CFFC8D8} |
| 2. II (II)   | {F165C13B-397E-4182-8D72-BD1D74CFBA00} |
| 3. III (III) | {D8FE8FA9-9BB5-416E-AE2C-60E627F76E5A} |
| 4. AVR (AVR) | {04AF2551-1185-4BA9-AEED-977D562D188A} |
| 5. AVL (AVL) | {9124A193-5F50-448D-9E60-C5869EC4029C} |
| 6. AVF (AVF) | {C8E97206-58AA-47C0-BD2E-7A30E56DB46A} |
| 7. V1 (V1)   | {4CAA3A92-F3DB-4847-8523-380EFEE6CFD5} |
| 8. V2 (V2)   | {05A3232E-A30D-4771-A1B5-E7705B717646} |
| 9. V3 (V3)   | {E150126E-4AA7-4F96-961E-FD1169C77158} |
| 10. V4 (V4)  | {255260F9-5CC7-41B1-831C-F9319FA27506} |
| 11. V5 (V5)  | {D5C10744-0EB5-46FA-AE6C-CCF571318F8A} |
| 12. V6 (V6)  | {B86F72D7-7884-4470-92A8-EF9277C5B584} |

UNITS:  
-----

Abbreviation..... THRILL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Thrill  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- 1. PRESENT (p)                    {384DC20F-1625-CE1E-1C85-A3B3924BB2A0}
  - 2. ABSENT (a)                    {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0}

LOCATION: 8 entries

- 
- 1. R FEMORAL (R FEMORAL)       {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
  - 2. L FEMORAL (LFEM)             {50AF6A08-CE37-6859-F4E6-762660A982F7}
  - 3. R CAROTID (R CAROTID)       {8280ACD2-72D6-202A-89AF-86BCA7736427}
  - 4. L CAROTID (L CAROTID)       {FCACD1EB-C12B-2515-4A68-4E9E28A72DBC}
  - 5. R TEMPORAL (R TEMPOR)       {1F50461D-5ECF-25B2-A29A-E02803CB7814}
  - 6. L TEMPORAL (L TEMPO)        {375E65CF-A59E-A1E9-EB6D-D53CD818939E}
  - 7. ABDOMINAL AORTA (ABD AORTA)  
                                  {C4BD42DD-1D6B-4B2C-B6D5-19E1185E06E3}
  - 8. HEART (HEART)                {F7D48A36-B512-42DD-B92F-07ACF2B95E08}

UNITS:

-----



=====

TERM: CARDIAC THRILL TIMING  
ID: {D52A1728-1692-4742-886A-B119BB7EE598}

=====

Page: 87

Abbreviation..... THRILL TIM  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Thrill Timing  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
1. CONTINUOUS (CONTINUOUS) {676CD2DA-F882-4BE5-99BB-510EA7669749}
  2. SYSTOLIC (SYSTOLIC) {859B5700-0990-42FF-8129-8DB24EFDA8D1}
  3. DIASTOLIC (DIASTOLIC) {3E9D1DCF-6B3B-482C-AE5F-442814F4B5C0}
  4. PRESYSTOLIC (PRESYSTOLI) {BB85B8EE-36F8-4AC2-9AB0-0BABCFA9FCBB}

UNITS:

-----

Abbreviation..... CAUTI: Alt  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... CAUTI: alternatives attempted prior  
to cath  
Help Text.....  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

-----

- 10. BLADDER SCAN (Blad Scan)  
{6C66BC97-1042-4E53-A49F-69AB7ABEE551}
- 20. EXTERNAL CONDOM CATH (Ex Con Cat)  
{747E493A-C535-44EA-AD6D-27CE3E8E8CAC}
- 30. BEDSIDE COMMODORE (Bedside Cm)  
{E2B0C66F-335A-418E-AEB8-C7B666AE2F8D}
- 40. TOILETING BY SCHEDULE (Toilet Sch)  
{D47E6118-E5AD-46A0-8DCE-B28CDAFD6154}

UNITS:

-----

=====

TERM: CAUTI: CATHETER IN PLACE  
ID: {A1A9B834-08E3-461E-8E52-03C49D1B63B0}

=====

Page: 89

Abbreviation..... CAUTI: Cat  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... CAUTI: Catheter in place  
Help Text.....  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 10 entries

- 
- 10. FROM OPERATING ROOM (From Op Rm) {2B26BD84-3D4D-4931-99C7-8204E9D49827}
  - 20. FROM EMERGENCY DEPARTMENT (Fr Emergen) {9F4AA078-FFA2-4F04-A171-BA0486FCCA0E}
  - 30. FROM INTERVENTIONAL RADIOLOGY (Fr Int Rad) {617DA086-B121-4E2F-8B02-41133EEDD50D}
  - 40. FROM ANOTHER VA FACILITY (Fr Oth VA) {66AF9A88-6404-4E05-9640-4FB0C6513BAF}
  - 50. FROM A NON-VA OUTSIDE FACILITY (Fr Non-VA) {8090F49E-5056-4F7C-932C-538581996F8F}
  - 60. FROM ANOTHER UNIT (Fr An Unit) {D94B5252-7E8A-48EB-A481-0B18C68CAEF7}
  - 70. FROM MICU (From MICU) {971D83AE-B75C-420B-B23D-844397CEE968}
  - 80. FROM SICU (From SICU) {892002CC-E07A-49A5-B909-1F88665B3D65}
  - 90. FROM CLC (From CLC) {8B7B82C4-2307-47B7-87D7-88F70CF3B9FB}
  - 100. FROM HOME (From Home) {2BB78383-F1D2-4FF8-AAE2-C16DC3BD4E8F}

UNITS:

-----

Abbreviation..... CAUTI: In  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... CAUTI: Indications  
Help Text.....  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 6 entries

- 
- 10. COMFORT CARE (Comfort Ca) {51669E41-DF78-4E49-BA5F-DE6F1B5F1148}
  - 20. NEUROGENIC BLADDER (Neuro Blad)  
{D386E29C-C4B7-4A9D-8092-CA978F9CC4D6}
  - 30. STAGE 3 OR 4 PRESSURE ULCER (Press Ulce)  
{883029EF-BB3A-4146-95DD-3C6F1C0D097D}
  - 40. SELECTED SURGERY (Sel Surg)  
{6DBC6356-0F00-4001-AF94-20A479995A6C}
  - 50. UNABLE TO VOID (Un to Void)  
{1F2A844D-4E14-4759-B8F7-82B58E157720}
  - 60. OTHER (OTHER) {19094E5B-AEC3-44B0-9BD5-44E79E7EB105}

UNITS:  
-----

Abbreviation..... CAUTI: Req  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... CAUTI: Insertion Requirements  
Help Text.....  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 9 entries

- 
- 10. HAND HYGIENE BEFORE AND AFTER (Hand Hyg)  
{6005023D-3C02-4C93-99E8-2A8167574E93}
  - 20. STERILE GLOVES-DRAPES USED (Ster Glove)  
{F26BBB1B-1483-40CD-AC47-828F71B9187C}
  - 30. ASEPTIC TECHNIQUE FOLLOWED (Aseptic)  
{2DE69DDB-F455-4570-A821-C81F39FD2574}
  - 40. SINGLE USE LUBRICANT PACKET (Single Use)  
{7F748EDB-77D4-42D2-BE75-54F73598D4B0}
  - 50. CATH SECURED WITH DRAIN BAG BELOW PT'S BLADDER LEVEL (Cath Secur)  
{A2D2ABF1-88D4-46A2-B935-33211292FF25}
  - 60. DRAINAGE FLOW IS UNOBSTRUCTED (Unobstruct)  
{DB66B2DD-50FE-468B-8892-1A8FF157BCB0}
  - 70. CATH SIZE DOCUMENTED (Cath Size)  
{FB8BE4D4-41C7-4F8C-AD17-74F8FE170BC3}
  - 80. CATHETERIZATION TIME DOCUMENTED (Cath Time)  
{4EF58FCA-A5E1-454C-B435-5BFB487618F4}
  - 90. PROVIDED PATIENT EDUCATION (Pat Ed)  
{25D0998B-3600-40C9-8608-9A0D91CCBEE0}

UNITS:  
-----

=====

TERM: CAUTION-VITALS MAY NOT DISPLAY IN CORRECT SEQUENCE. REVIEW TIMES CAREFULLY

Page: 92

ID: {130B6426-2DE1-4FBE-ADDD-AA8C66F24A37}

=====

Abbreviation..... Caution  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... String  
Description.....  
Display Name..... Caution-Vitals may not display in  
correct sequence. Review times  
carefully  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: CENTRAL LINE ALARMS ON  
ID: {152A037A-A6CC-4F12-8B85-25CEE78772C4}

=====

Page: 93

Abbreviation..... C-LINE ALA  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Central Line Alarms On  
Display Name..... Central Line Alarms On  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

Abbreviation..... C LINE CHA  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Central Line Changed  
Help Text..... Central Line changed to  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
1. NEW SITE (NEW SITE) {10B1FF61-E58E-4608-9CCB-9232BC964418}
  2. OVER WIRE (OVER WIRE) {157CD456-1020-43AE-9AAA-1686F1574B6F}
  3. INSERTED THROUGH EXISTING INTRODUCER (INS EXS IN)  
{20BC183A-197F-4464-A27C-47F63C5B842A}
  4. INSERTED THROUGH NEW INTRODUCER (INS NEW IN)  
{7FCA403D-41E5-4BAA-B155-B58F35EA49B3}

UNITS:  
-----



=====

TERM: CENTRAL LINE DRESSING CHANGE  
ID: {DDEA3089-9246-40F3-9A67-56A05B393309}

=====

Page: 95

Abbreviation..... C LI DR CH  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Central Line Dressing Change  
Display Name..... Central Line Dressing Change  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
1. DRY (d) {D8F4F0FD-1B19-2515-3815-EC88D175B305}
  2. INTACT (i) {935C9CAF-6C88-6610-C6E9-4C1B822923AF}
  3. CHANGED (CHANGED) {74B14023-D74D-41F6-AA48-B0E5137623CD}
  4. REINFORCED (REINFORCED) {EC524C50-0CDE-4AA1-8059-E4323850575E}

UNITS:  
-----

=====

TERM: CENTRAL LINE DRESSING TYPE  
ID: {AF276338-5F2C-1CEF-F009-AC4BAED0DDE3}

=====

Page: 96

Abbreviation..... C LINE DRE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Central Line Dressing Type  
Help Text..... Type of dressing on central line  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. ANTI-MICROBIAL DRESSING (AB DRESSIN)  
{9734C5E9-D099-47EB-93B1-D82201F64331}
  2. PRESSURE DRESSING (PRESS DRES)  
{6D9DB9EF-85EB-8321-03FC-879AA4C5615B}
  3. TRANSPARENT (TRANSPAREN) {0E73D24D-0EFC-08F9-196A-D30DD1BB3228}

UNITS:  
-----



=====

TERM: CENTRAL LINE SITE CARE  
ID: {4620B5B3-22F7-45DA-A52E-6D3EAE1C4202}

=====

Page: 98

Abbreviation..... C-LINE CAR  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Central Line Site Care  
Display Name..... Central Line Site Care  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

Abbreviation..... C LINE CON
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Central Line Site Condition
Help Text..... Condition of central line
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 7 entries

- 2. RED (RED) {515D7CFB-3527-30E7-EC4B-B607C272DFA0}
3. SWOLLEN (SWOLLEN) {DCE86B10-C393-8DA9-6518-F2943F0E1A4E}
4. DISCOMFORT (DISCOMFORT) {77D058C9-1EB4-96C0-C9E3-666C840B884C}
50. PATENT (Patent) {F0C38BCB-1D11-436C-8B8C-4297C27DDE04}
60. CLEAN (Clean) {1F187B99-9963-4AD5-98F6-D8D24D2DDD30}
70. OCCLUDED (Occluded) {25509BE3-F946-4AD7-80E8-5CE8A2E61AD6}
80. OOZING (Oozing) {C23549D7-3238-4523-BA31-AB71609F01C2}

UNITS:
-----

=====

TERM: CENTRAL LINE SIZE FRENCH  
ID: {234D330A-C51B-4D56-8631-9418BB081E96}

=====

Page: 100

Abbreviation..... CLINE FREN  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Central Line Size (French)  
Display Name..... Central Line Size (French)  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

FRENCH                    0.0 > 6.0 \* 10.0 < 50.0    Dec Prec: 1

=====

TERM: CENTRAL LINE SIZE GAUGE  
ID: {870FAA87-69C9-4E00-97A8-CD5CC79B41D6}

=====

Page: 101

Abbreviation..... C-LINE GAU  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Central Line Size (Gauge)  
Display Name..... Central Line Size (Gauge)  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
GAUGE                    0 > 12 \* 21 < 40        Dec Prec: 0

Abbreviation..... C LINE TYP
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Central Line Type
Help Text..... Type of central line
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 6 entries

- 1. SINGLE (SINGLE) {8D305E0F-9035-E396-E768-40D7E6455BE4}
2. DOUBLE (DOUBLE) {A486852F-150B-24E3-8EB3-8CAABB3D0520}
3. TRIPLE (TRIPLE) {077BE8CF-F2AB-56E1-1F25-51DF80F8054D}
4. INTRODUCER (INTRODUCER) {B46F7CA4-C94E-68FD-2892-52866763B5A0}
5. PICC (PICC) {3CFE33D0-5A35-88CC-152C-D7501D1858D9}
6. DIALYSIS (DIALYSIS) {D218E23C-5A34-09B0-DE43-DC792DCD0F1C}

LOCATION: 16 entries

- 1. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
2. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
3. L INTERNAL JUGULAR (LIJ) {AFBB4FE5-8A11-2047-7900-DADF2FC4430D}
4. R INTERNAL JUGULAR (R INT JUGU) {7ABF6916-D379-198D-6454-F0518A74844A}
5. L EXTERNAL JUGULAR (L EXT JUGU) {6F7B331B-4F83-0B79-1755-F6F673D647D6}
6. R EXTERNAL JUGULAR (R EXT JUGU) {4D78F704-9E3B-907E-A792-4DC4F4CFF170}
7. L SUBCLAVIAN (LSCL) {71B93594-15F9-AA4A-AEF2-E222EC28C00D}
8. R SUBCLAVIAN (R SUBCLAVI) {53DEAE83-93AE-123B-14E5-4DE38CD4F752}
9. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
10. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
11. L ANTECUBITAL (L ANTECUBI) {C21838CC-BAB5-1D7F-E2A4-A7C663FB28C3}
12. R ANTECUBITAL (R ANTECUBI) {7FB00A65-B4B2-06A6-BD89-98AC112D8D4E}
13. L UPPER ARM (L UPP ARM) {5F915696-1B08-32E9-6B95-49246117D660}
14. R UPPER ARM (R UPP ARM) {3EF20E8E-F68C-E4F2-8F8C-FAD3B44DDAA4}
15. L WRIST (L WRIST) {7BCAFACE-FF3B-46DA-99F9-6437A90B7126}
16. R WRIST (R WRIST) {FDD9FE1D-F8A2-09AD-82D7-AA478693DD01}

UNITS:



=====

TERM: CENTRAL LINE WAVEFORM  
ID: {4F847833-9A36-4EAE-BC73-B7AB069D39E4}

=====

Page: 103

Abbreviation..... C-LINE WAV  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Central Line Waveform  
Display Name..... Central Line Waveform  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. NORMAL (n) {553C9A97-A3D5-C52F-0E4E-C7D168C99A2B}
  2. DAMPENED (DAMPENED) {77C777B3-7E56-41E3-8E70-DE144E6A3424}
  3. HYPERDYNAMIC (HYPERDYNAM) {818A8A99-04D9-4989-A9BA-ECFACED74605}
  4. POSITIONAL (POSITIONAL) {9A74D16B-5E95-4EA8-88AF-E53C80311A5D}
  5. NOT MONITORED (NOT MONITO)  
{0E6A8A07-48E0-428F-8D86-22572CABCF48}

UNITS:  
-----

=====

TERM: CENTRAL LINE ZEROED  
ID: {7868021A-72D4-451F-A2E3-E3106DEC52B4}

=====

Page: 104

Abbreviation..... C-LINE ZER  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Central Line Zeroed  
Display Name..... Central Line Zeroed  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

=====

TERM: CENTRAL VENOUS PRESSURE  
ID: {D30F98A7-4C5D-12E8-AB4D-9C85A4332EC3}

=====

Page: 105

Abbreviation..... CVP  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Central Venous Pressure  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid..... 4688719

LOCATION: 8 entries

- 
1. L SUBCLAVIAN (LSCL) {71B93594-15F9-AA4A-AEF2-E222EC28C00D}
  2. R SUBCLAVIAN (R SUBCLAVI) {53DEAE83-93AE-123B-14E5-4DE38CD4F752}
  3. L INTERNAL JUGULAR (LIJ) {AFBB4FE5-8A11-2047-7900-DADF2FC4430D}
  4. R INTERNAL JUGULAR (R INT JUGU)  
{7ABF6916-D379-198D-6454-F0518A74844A}
  5. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
  6. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
  7. L BRACHIAL (LBRA) {75F6CE33-1AB3-A1DA-70B8-04825290B577}
  8. R BRACHIAL (RBRA) {A88E69A8-F6AE-80B0-2CD4-687B121F9DDE}

METHOD: 2 entries

- 
1. CARDIAC OUTPUT MONITOR (COM)  
{CCBE5744-7CB9-9DFA-938F-FDDBC73B45B3}
  2. MONITOR (Mo) {63FC24EA-50B8-B302-D948-ABE3517CF6FC}

UNITS:

-----

MILLIMETERS OF MERCURY 0 > 2 \* 6 < 19 Dec Prec: 0

=====

TERM: CEREBRAL PERFUSION PRESSURE - DO NOT CHECK TOTAL OR SUBTOTALS Page: 106  
ID: {166BD1C4-D983-464A-9A6E-B3A7AC7885EB}

=====

Abbreviation..... CPP  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Cerebral Perfusion Pressure  
Display Name..... CPP - Do Not Check Total or Subtotals  
Help Text..... Cerebral Perfusion Pressure  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
MILLIMETERS OF MERCURY      0 > 50 \* 150 < 300      Dec Prec: 0

=====

TERM: CHEST EXPANSION  
ID: {BAF272C6-1F0A-40E1-918B-3BABC5F22407}

=====

Page: 107

Abbreviation..... CHEST EXPA  
Active..... Yes  
Boolean Value False..... Asymmetrical  
Boolean Value True..... Symmetrical  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Chest Expansion  
Help Text..... Symetrical or Asymetrical Chest  
Expansion  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

ID: {57AE16AF-8EC2-5800-131F-0AF254C3FA76}

=====

Abbreviation..... CHEST PN  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... EVALUATION OF PATIENTS CHEST DISCOMFORT  
Display Name..... Chest Pain  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 11 entries

-----

1. 0 (0)	{2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
2. 1 (1)	{8F5D10A8-FE17-47BC-DACA-6972D12CC00A}
3. 2 (2)	{38E22B79-A2B6-3401-B195-C955F24E1278}
4. 3 (3)	{E714454D-6A67-91F6-2FB7-D8F24F184328}
5. 4 (4)	{D2280C7E-DE50-5A8D-C3FD-5862B69C0E1A}
6. 5 (5)	{188EDCE7-0726-4234-D8D2-92C8052CD221}
7. 6 (6)	{640CF2F9-3BBF-AB13-185B-CD02A609FB5A}
8. 7 (7)	{7A010324-D2C2-5D4B-848A-7B8878666D20}
9. 8 (8)	{518E7CFC-102E-351E-C7AF-FB575E3C8F4D}
10. 9 (9)	{58FF2FD4-569F-F454-9F5B-FE6F9E34D8CC}
11. 10 (10)	{0C827160-18DB-816F-7FAA-A3C175C6D33B}

LOCATION: 3 entries

-----

1. L CHEST (L CHEST)	{E332E497-FE0C-47CB-A0A7-217FBF52DC97}
2. R CHEST (R CHEST)	{7E30EC0E-815A-4C2E-A9FC-3C2B2D6E3CD4}
3. MID-STERNAL (MID-STERN)	{3F5168F1-9082-4FB2-A364-3CDE33B1214B}

UNITS:

-----

Abbreviation..... DRAIN CHAR  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... DESCRIPTION OF PATIENT'S DRAINAGE  
Display Name..... Chest Tube Characteristics  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
1. SEROUS (SEROUS) {9BDA347E-29E9-80C7-91FA-016192D4257E}
  2. SEROSANGUINOUS (SEROSANGUI)  
{19A8D332-899D-177A-2247-ABCDAAEE3BFA8}
  3. PURULENT (PURULENT) {86C01494-743A-C6A2-8883-01CA5200C79B}
  4. BLOODY (BLOODY) {D7904E84-8331-D6DD-F2C0-7F3039C5746A}

UNITS:  
-----

Abbreviation..... CH TU SU T  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Chest Tube Drainage Suction Type  
Display Name..... Chest Tube Drainage Suction Type  
Help Text..... Type of Chest Tube Suction  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

-----

1. WATERSEAL (WATERSEAL) {6D91B99D-F199-09ED-78BE-B54F2CAC4E06}
2. GRAVITY (GRAVITY) {C0233F70-376D-41CE-9F14-714C7B35244F}
3. CONTINUOUS SUCTION (CONT SUCTI)  
{8C10DAC2-9AE5-40B9-99C7-25D47F4F13BA}

LOCATION: 12 entries

-----

1. L UPPER ANTERIOR CHEST (L U A CHES)  
{FE63B58A-DE39-4C03-AE61-E34A5487C767}
2. R UPPER ANTERIOR CHEST (R U A CHES)  
{67890049-A4A8-4241-BDF3-573729BD6EB5}
3. L MIDDLE ANTERIOR CHEST (L M A CHES)  
{36A83F77-117E-474B-AE34-0B2C02E131AE}
4. R MIDDLE ANTERIOR CHEST (R M A CHES)  
{2F751A9C-66CA-4C22-ACC9-6E81FA63FA2B}
5. L LOWER ANTERIOR CHEST (L L A CHES)  
{6073B3C6-AA4D-4BBA-9FF2-32ED40CF92AD}
6. R LOWER ANTERIOR CHEST (R L A CHES)  
{96BCF731-B661-4455-B5E7-C4AE2509505F}
7. L UPPER POSTERIOR CHEST (L U P CHES)  
{6D5674E9-EF04-4EF1-84A0-3F5ED2FB532A}
8. R UPPER POSTERIOR CHEST (R U P CHES)  
{D8F03FA5-1FD7-46EE-8FA4-17B3B7E62141}
9. L MIDDLE POSTERIOR CHEST (L M P CHES)  
{5A656939-F82B-4DEC-9F1D-B5FA644E3DA9}
10. R MIDDLE POSTERIOR CHEST (R M P CHES)  
{2BC487D6-6CE8-499E-94A0-A79EAA8E199B}
11. L LOWER POSTERIOR CHEST (L L P CHES)  
{2C19F10B-7546-4A1F-A5DC-3867359ED971}
12. R LOWER POSTERIOR CHEST (R L P CHES)  
{5163180E-FFD2-4BFB-BF8B-4C6B18EFE5D4}

UNITS:

-----





=====

TERM: CHEST TUBE DRESSING CHANGE  
ID: {8DE63ADE-D3AC-49EE-BB24-ED44AB9971F5}

=====

Page: 112

Abbreviation..... CH T DR CH  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Chest Tube Dressing Change  
Display Name..... Chest Tube Dressing Change  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

Abbreviation..... CH TU ASSE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Chest Tube Site Assessment  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. NORMAL (n) {553C9A97-A3D5-C52F-0E4E-C7D168C99A2B}
  2. AIR LEAK (AIR LEAK) {52E40E25-6D28-53EC-BC05-D1FEAA4616EF}
  3. FLUCTUATES (Fluctuates) {C2697650-EF53-5F2E-EEE3-7B2E298F2FA0}
  4. CREPITUS (CREPITUS) {56A97506-05F0-4F37-8FD9-3AA82191408A}
  5. SQ EMPHYSEMA (SQ EMPHYS.) {B8D03538-0168-4CDC-FF41-BA1185EBCC87}

UNITS:  
-----

=====

TERM: CHEST TUBE SUCTION (CM H2O)  
ID: {BD9724EC-B5BB-4F95-8C24-B0A1451C33A6}

=====

Page: 114

Abbreviation..... CHST TU SU  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Chest Tube Suction (cm H2O)  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

CENTIMETERS H2O            0 > 10 \* 20 < 30            Dec Prec: 0

=====

TERM: CHIEF CONCERN/COMPLAINT  
ID: {2A6683ED-9258-4CE9-8793-2FD0B2A56E6F}

=====

Page: 115

Abbreviation..... Chief Com  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... String  
Description.....  
Display Name..... Chief Concern/Complaint  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

Abbreviation..... CG
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... circumference
Display Name..... Circumference/Girth
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid..... 4688720

LOCATION: 11 entries

- 1. ABDOMINAL (Abd) {AE1DA31F-2C28-98C2-501A-F3D5ADB2501F}
2. ANKLE (Ank) {F2A2CF53-BFA5-9279-DB22-EA2627988E5B}
3. CALF (Cif) {7923E496-90D6-B6D5-726A-A5B8DF3436CF}
4. HEAD (Hd) {C3209F44-90C6-FF1B-3169-265B294F276A}
5. HIP (Hip) {D2C8BD11-7FB2-7217-A663-EBFC8C9428E7}
6. KNEE (Kn) {B577308A-4C37-47BD-AD84-38E6342B6DF7}
7. LOWER ARM (LrA) {D7627720-D6B0-AC00-8FA3-F7B569ED06FC}
8. NECK (Ne) {D4812C3B-AAC2-1B61-679B-3CD8B5CB947C}
9. THIGH (THIGH) {2DA93C27-12B0-C76C-85D5-B7D6375F3AD3}
10. UPPER ARM (UpA) {B43A2B31-1F35-3F4D-5D23-882541ABE477}
11. WRIST (Wri) {C4F1453C-A20A-E48C-0BD8-FAF5E5568ADA}

UNITS:

INCH 29 > 30 \* 60 < 61 Dec Prec: 0
CENTIMETER 75 > 76 \* 156 < 157 Dec Prec: 0

Abbreviation..... CIWA AG  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Agitation  
Help Text..... CIWA Agitation Observation  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 8 entries

- 
- 1. 0 - NORMAL ACTIVITY (CIWA AG 0)  
{9FF32084-A432-467C-B0A2-C8BC9F3D6D13}
  - 2. 1- SOMEWHAT MORE THAN NORMAL ACTIVITY (CIWA AG 1)  
{497F7A3D-72C8-4986-A1E6-74EA0E5704F9}
  - 3. 2 (2) {38E22B79-A2B6-3401-B195-C955F24E1278}
  - 4. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}
  - 5. 4 - MODERATELY FIDGETY AND RESTLESS (CIWA AG 4)  
{B1515E1F-4305-40D5-84B4-5E37AEF0AD3D}
  - 6. 5 (5) {188EDCE7-0726-4234-D8D2-92C8052CD221}
  - 7. 6 (6) {640CF2F9-3BBF-AB13-185B-CD02A609FB5A}
  - 8. 7 - PACING BACK AND FORTH MOST OF THE INTERVIEW (CIWA AG 7)  
{9F03A167-CC40-4D65-9FD8-CFB3FEEAF833}

UNITS:  
-----

=====

TERM: CIWA: ANXIETY  
ID: {9552124C-8650-406F-86FA-AB31D04AC8EF}

=====

Page: 118

Abbreviation..... CIWA AX  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Anxiety  
Help Text..... Ask "Do You Feel Nervous?"  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 8 entries

- 
- 1. 0 - NO ANXIETY AT EASE (CIWA AX 0)  
{DA61A13E-6443-4E8A-9C76-762FEACDB4E7}
  - 2. 1 - MILDLY ANXIOUS (CIWA AX 1)  
{401113EE-5CE7-4AE9-80C1-1C3252E706D0}
  - 3. 2 (2) {38E22B79-A2B6-3401-B195-C955F24E1278}
  - 4. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}
  - 5. 4 - MODERATELY ANXIOUS (CIWA AX 4)  
{9A01FEFD-844C-4FAB-BEC3-96333495FA62}
  - 6. 5 (5) {188EDCE7-0726-4234-D8D2-92C8052CD221}
  - 7. 6 (6) {640CF2F9-3BBF-AB13-185B-CD02A609FB5A}
  - 8. 7 - ACUTE PANIC (CIWA AX 7)  
{3F5E6786-FD67-405F-A535-EC1609910F14}

UNITS:  
-----





Abbreviation..... CIWA HA  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Headache, Fullness In Head  
Help Text..... Ask "Does your head feel different?  
Does it feel like there is a band  
around your head? Do not rate for  
dizziness or lightheadedness."  
Otherwise, rate severity.  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Void.....

DATA VALUE: 8 entries

- 
1. 0 - ABSENT (0-a) {75F14311-F113-963C-4F05-B8AE17F044A9}
  2. 1 - VERY MILD (CIWA HA 1) {F789C1EA-CC7A-4368-87BA-08AFFE492EE6}
  3. 2 - MILD (CIWA HA 2) {6B752397-0A4F-4176-90E7-E220DB2BF1C4}
  4. 3 - MODERATE (CIWA HA 3) {F26EB828-6864-4384-9EBB-35A6CDCACC59}
  5. 4 - MODERATELY SEVERE (CIWA HA 4)  
{83B3DE0F-3BA4-4B0A-A222-F372BD403094}
  6. 5 - SEVERE (CIWA HA 5) {7BD7B91F-DF6F-47B9-BAA6-3EF249340578}
  7. 6 - VERY SEVERE (CIWA HA 6)  
{241C1EF4-2E2A-4D65-9991-607BB6B29E01}
  8. 7 - EXTREMELY SEVERE (CIWA HA 7)  
{F8DCC910-D36D-480E-867E-B45221170C49}

UNITS:

-----



Abbreviation..... CIWA OR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Orientation And Clouding Of Sensorium
Help Text..... Ask "What day is this? Where are you?
Who am I?"
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

- 1. 0 - ORIENTED AND CAN DO SERIAL ADDITIONS (CIWA OR 0)
{D8015F54-6E9E-457A-955C-38565BDDDBB6}
2. 1 - CANNOT DO SERIAL ADDITIONS (CIWA OR 1)
{F4446E01-E326-420E-8EA8-04451AA9D963}
3. 2 - DISORIENTED FOR DATE < 2 CALENDAR DAYS (CIWA OR 2)
{01B894B0-3E70-47FA-820A-83CA48B93EF9}
4. 3 - DISORIENTED FOR DATE > 2 CALENDAR DAYS (CIWA OR 3)
{68C351F9-C8DE-49C5-B3BD-D69E97F3E516}
5. 4 - DISORIENTED FOR PLACE OR PERSON (CIWA OR 4)
{03BBFFA2-233A-41A1-8212-F988889C5E72}

UNITS:

Abbreviation..... CIWA PS  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Paroxysmal Sweats  
Help Text..... Paroxysmal Sweats  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 8 entries

- 
- 1. 0 - NO SWEAT VISIBLE (CIWA PS 0)  
{7C83E616-4414-4766-AAF5-5F5442FFBF09}
  - 2. 1 - BARELY PERCEPTIBLE SWEATING / PALMS MOIST (CIWA PS 1)  
{93B5B833-14B5-4ADD-97E8-7BF14A8BA01A}
  - 3. 2 (2) {38E22B79-A2B6-3401-B195-C955F24E1278}
  - 4. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}
  - 5. 4 - BEADS OF SWEAT OBVIOUS ON FOREHEAD (CIWA PS 4)  
{8388B1A0-46EA-4003-A32B-A2C708A6F056}
  - 6. 5 (5) {188EDCE7-0726-4234-D8D2-92C8052CD221}
  - 7. 6 (6) {640CF2F9-3BBF-AB13-185B-CD02A609FB5A}
  - 8. 7 - DRENCHING SWEAT (CIWA PS 7)  
{15D2E2CC-BAED-403B-A141-B848FFE4A9D0}

UNITS:  
-----



=====

TERM: CIWA: TOTAL

Page: 125

ID: {BD47158C-B56B-4B92-AD2E-33B7149EF6C0}

=====

Abbreviation..... CIWA TOTAL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Total CIWA Score  
Help Text..... Add All CIWA Scores Together  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

CIWA SCALE                    -1 > 0 \* 9 < 68                    Dec Prec: 0

=====

TERM: CIWA: TREMOR  
ID: {E87CAD13-134C-4556-8D55-49450A0D544C}

=====

Page: 126

Abbreviation..... CIWA TR  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Tremor  
Help Text..... Arms extended and fingers spread apart.  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 8 entries

- 
1. 0 - NO TREMOR (CIWA TR 0) {2D1A5B39-4C5D-4345-84B4-DDCBC8E25B4D}
  2. 1 - NOT VISIBLE BUT CAN BE FELT FINGERTIP (CIWA TR 1)  
{31833D55-9018-47EE-8341-A0863744CEF8}
  3. 2 (2) {38E22B79-A2B6-3401-B195-C955F24E1278}
  4. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}
  5. 4 - MODERATE WITH PATIENTS ARMS EXTENDED (CIWA TR 4)  
{79CE1486-151C-4492-AFD2-7845DD69DA9E}
  6. 5 (5) {188EDCE7-0726-4234-D8D2-92C8052CD221}
  7. 6 (6) {640CF2F9-3BBF-AB13-185B-CD02A609FB5A}
  8. 7 - SEVERE EVEN WITH ARMS NOT EXTENDED (CIWA TR 7)  
{FE0BE72E-913A-48CD-BABF-C5B742AFED93}

UNITS:  
-----





Abbreviation..... CO BEH
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Close Observation Behavior
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 7 entries

- 1. WALKING OR PACING (wop) {CE7AABCB-84B7-A310-3A65-0AF17FB7B0E7}
2. LYING OR SITTING (los) {52A940B1-55AF-AB4E-9E59-F93515A1296C}
3. QUIET OR RESTING (qor) {76E4908F-014C-2088-FF7F-33EBF223DE40}
4. SLEEPING (s) {87338420-F351-F98C-7A80-4990B87ACE2B}
5. WITH PHYSICIAN (wp) {1673C53F-DA8C-2CEC-679B-B01D914FE533}
6. WITH THERAPIST (wt) {795A1300-1D60-9D67-692D-A2AB2B190596}
7. GERI-CHAIR (g) {167F9786-8F24-054F-020A-658A6EE5B0B4}

UNITS:
-----

Abbreviation..... CO CUR LOC  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Close Observation Current Location  
Display Name..... Close Observation Current Location  
Help Text..... Close Observation Current Location  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 10 entries

- 
1. PATIENT'S ROOM (pr) {0F13E11D-9F78-14A3-0568-AF26393F39F1}
  2. BATHROOM (b) {D20D91A7-E9F5-D000-5E0B-772BA9F5F72E}
  3. HALL (h) {BFCEA10F-2663-199E-88D8-E263122D7D26}
  4. DAY ROOM (dr) {269125EF-FA82-B17A-6B56-0CEAB4788899}
  5. QUIET ROOM (qr) {CCD96D41-C273-F6F3-DA18-8818B1CC9980}
  6. GROUP ROOM (gr) {D41C3D9F-8F10-1A0F-D571-20DE0EA36328}
  7. SMOKING AREA (sa) {8AFDBCB0-9FC8-523F-5D07-4D8CCBA7036F}
  8. SMOKING PORCH (sp) {5AD3F05B-2A5B-E2B4-FD47-66712CFA44E6}
  9. MEDICATION ROOM (mr) {9571AA9F-F538-55EC-37F6-99FE5677B970}
  10. NURSING STATION (ns) {B6E48F92-B9BA-488F-3A2E-D74568A6B855}

UNITS:

-----

=====

TERM: CLOSE OBSERVATION REASON FOR PRECAUTION  
ID: {1B131F5A-CD32-38CA-19A5-1038A07C9867}

=====

Page: 130

Abbreviation..... CO RFP  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Close Observation Reason for Precaution  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. SUICIDAL (s) {B104E9C7-A78E-82DC-49A4-19B94A639AFE}
  2. ESCAPE (e) {965F3201-9D63-66AB-FA25-99479AE7145A}
  3. ASSAULTIVE (a) {99ECD8C7-191A-3B9C-2A7A-803A921065FA}
  4. ELOPEMENT (e) {EB030907-FB90-CA9D-CF50-E26C5D9FE671}
  5. MEDICAL DEVICE SAFETY (mds)  
{1985CABE-4F53-9FCF-F444-F9AC53C04DBD}

UNITS:  
-----

=====

TERM: CODE STATUS

Page: 131

ID: {95248A8C-E05A-94F3-CA15-6F113997A621}

=====

Abbreviation..... CS  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Code Status  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
1. FULL CODE (fc) {1DED65E6-2F2E-8C5F-D009-BE6467A9A2AA}
  2. DO NOT RESUSCITATE (dnr) {0BD8A134-C563-7E00-A6B1-F7C543CEF6FF}

UNITS:  
-----

=====

TERM: COLLOID TYPE

Page: 132

ID: {C22CB099-BC3F-42AE-9A11-67B99093D959}

=====

Abbreviation..... COLLOID TY  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Type of Colloid  
Display Name..... Colloid Type  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. 5 PERCENT ALBUMIN (5 PERC ALB)  
    {38CE47CB-736B-475E-8E89-C4C13F2DF184}
  2. 25 PERCENT ALBUMIN (25 PERC AL)  
    {2FCE19CA-B5AB-4F15-8FEA-CF156D1E4CBC}
  3. 6 PERCENT HETASTARCH (6 PERC HET)  
    {D94072C9-87E8-46CA-929F-4904F3BD964C}

UNITS:  
-----

=====

TERM: COMA SCORE EYE  
ID: {EAAF4E71-3BE6-5058-07FF-1A3D977BC569}

=====

Page: 133

Abbreviation..... CSE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... 0: CLOSED, 1: NONE, 2: TO PAIN, 3: TO  
SPEECH, 4: OPEN SPONTANEOUSLY  
Display Name..... Coma Score  
Help Text..... 0: CLOSED, 1: NONE, 2: TO PAIN, 3: TO  
SPEECH, 4: OPEN SPONTANEOUSLY  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Void.....

DATA VALUE: 5 entries

- 
- 1. CLOSED (CLOSED) {9BA6EE03-F562-4DE2-99A7-337F5A5E34F2}
  - 2. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
  - 3. TO PAIN (TO PAIN) {AA1F818D-DCFC-4FBC-97D3-BE86D309119E}
  - 4. TO SPEECH (TO SPEECH) {2DDD94FA-D4D0-48B0-B4F4-F25C42E4AF96}
  - 5. OPEN SPONTANEOUSLY (OPEN SPONT)  
{013F1473-BDC9-4E7C-9709-C577357709F6}

LOCATION: 2 entries

- 
- 1. L EYE (L EYE) {D714A0CA-EAD6-EC35-DA54-FF027C23FF45}
  - 2. R EYE (R EYE) {258F0C1D-67BD-ECF8-3152-EB289DEF4681}

UNITS:

-----

Abbreviation..... CSM  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Coma Score Motor  
Help Text..... 1: None, 2: Extends to pain, 3:  
Flexion to pain, 4: Withdraws to  
pain, 5: Localizes to pain, 6: Obeys  
Commands  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 6 entries

- 
1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
  2. EXTENDS TO PAIN (EXT TO PAI)  
{A62B2C4D-BDB7-49FA-9669-DF6732063DA6}
  3. FLEXION TO PAIN (FLEX TO PA)  
{B9F5046F-4F8A-4A20-9D47-90432DC8FEC2}
  4. WITHDRAWS TO PAIN (WITHD TO P)  
{013527C8-FF6B-49CD-BF45-9338F7413E81}
  5. LOCALIZES TO PAIN (LOC TO PAI)  
{4072DFF9-2D90-423F-B6F0-0F2B924025D9}
  6. OBEYS COMMANDS (OBEYS COMM)  
{EE88EE08-4C71-4C70-B93E-65F9091544C1}

UNITS:  
-----



Abbreviation..... CSV
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... 0: Tubed/Dysphagic, 1: None, 2:
Incomprehensible, 3: Inappropriate,
4: Confused, 5: Oriented
Display Name..... Coma Score Verbal
Help Text..... 0: Tubed/Dysphagic, 1: None, 2:
Incomprehensible, 3: Inappropriate,
4: Confused, 5: Oriented
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 6 entries

- 1. TUBED/DYSPHAGIC (TUBE/DYSPH) {82AAFAE3-50B0-4F2A-91CF-7900849A14EB}
2. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
3. INCOMPREHENSIBLE (INCOMPREHE) {916CD781-563D-4F54-8752-800D37BFAA3B}
4. INAPPROPRIATE (Inappropri) {7E3B54C2-8A49-28EC-055E-FEA87599421C}
5. CONFUSED (CONFUSED) {31DAAA60-BE89-49A7-A10D-EB0DA82F12C0}
6. ORIENTED (ORIENTED) {B8DC9245-6D81-40E4-8B47-B6413E49F152}

UNITS:



Abbreviation..... COUGH  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Patients cough character  
Display Name..... Cough  
Help Text..... Patients cough character  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. PRODUCTIVE (p) {22BE0499-F9FE-3A27-2605-61EBEC0ED6AB}
  2. NON-PRODUCTIVE (np) {9AB01716-E395-18B9-6C5D-9BE6AEE41D5B}
  3. ABSENT (a) {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0}

UNITS:

-----

=====

TERM: CPAP Page: 138  
ID: {11AD0BDB-8E33-C04E-29B2-CCB22FDF9EEC}

=====

Abbreviation..... CPAP  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Continuous Positive Airway Pressure  
Display Name..... CPAP  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid..... 4688690

UNITS:

-----

CENTIMETERS H2O            -1 > 0 \* 15 < 100      Dec Prec: 0



=====

TERM: DENTURES

Page: 140

ID: {96B79B0B-3D26-42E9-A082-A601E5832633}

=====

Abbreviation..... Dentures  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... Dentures  
Help Text.....  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 6 entries

- 
1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
  2. FULL UPPER (Full Up) {1808E638-62BB-4978-8706-34BD74B427DE}
  3. FULL LOWER (Full Low) {D38022EF-F4F0-44A7-82FE-B277A2A0A165}
  4. PARTIAL UPPER (Part Up) {D7B30583-C22C-4B7E-BE7F-36FA9A577DAB}
  5. PARTIAL LOWER (Part Low) {2556CE7E-F211-4992-BCAF-09677C20BEFE}
  6. DENTAL IMPLANTS (Implants)  
{D3AF96DE-6920-4C93-8105-99EE574FC092}

UNITS:  
-----

=====

TERM: DIAGNOSIS

Page: 141

ID: {025498CA-716F-440A-9638-CF34CDE29413}

=====

Abbreviation..... Diagnosis  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... String  
Description.....  
Display Name..... Diagnosis  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: DIAGNOSIS TYPE

Page: 142

ID: {3289986D-87C6-4796-AA19-76010682BBAE}

=====

Abbreviation..... DIAG TYPE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Type of Diagnosis  
Display Name..... Diagnosis Type  
Help Text..... Type of Diagnosis  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. POST-OPERATIVE (POST-OP) {9FA19656-849F-4A65-92FC-755A4C9E8056}
  2. ADMISSION (ADMISSION) {483C2830-813A-46C3-A012-CCD70044B2A2}
  3. POST PROCEDURE DIAGNOSIS (PO PRO DIA)  
{933C7D2C-88F8-4C18-ACBB-A16A295D1B3A}
  4. DISCHARGE (DISCHARGE) {ADF67CC2-E05B-4B8D-A233-4AE1C35B2D20}
  5. TRANSFER (TRANSFER) {DFDB4783-25BE-485B-A9BA-621F1004AAF2}

UNITS:

-----



Abbreviation..... D CATH CON
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Dialysis Catheter Site Condition
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 7 entries

- 2. RED (RED) {515D7CFB-3527-30E7-EC4B-B607C272DFA0}
3. SWOLLEN (SWOLLEN) {DCE86B10-C393-8DA9-6518-F2943F0E1A4E}
4. DISCOMFORT (DISCOMFORT) {77D058C9-1EB4-96C0-C9E3-666C840B884C}
170. PATENT (Patent) {F0C38BCB-1D11-436C-8B8C-4297C27DDE04}
180. CLEAN (Clean) {1F187B99-9963-4AD5-98F6-D8D24D2DDD30}
190. OCCLUDED (Occluded) {25509BE3-F946-4AD7-80E8-5CE8A2E61AD6}
200. OOZING (Oozing) {C23549D7-3238-4523-BA31-AB71609F01C2}

LOCATION: 16 entries

- 1. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
2. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
3. L INTERNAL JUGULAR (LIJ) {AFBB4FE5-8A11-2047-7900-DADF2FC4430D}
4. R INTERNAL JUGULAR (R INT JUGU) {7ABF6916-D379-198D-6454-F0518A74844A}
5. L EXTERNAL JUGULAR (L EXT JUGU) {6F7B331B-4F83-0B79-1755-F6F673D647D6}
6. R EXTERNAL JUGULAR (R EXT JUGU) {4D78F704-9E3B-907E-A792-4DC4F4CFF170}
7. L SUBCLAVIAN (LSCL) {71B93594-15F9-AA4A-AEF2-E222EC28C00D}
8. R SUBCLAVIAN (R SUBCLAVI) {53DEAE83-93AE-123B-14E5-4DE38CD4F752}
9. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
10. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
11. L ANTECUBITAL (L ANTECUBI) {C21838CC-BAB5-1D7F-E2A4-A7C663FB28C3}
12. R ANTECUBITAL (R ANTECUBI) {7FB00A65-B4B2-06A6-BD89-98AC112D8D4E}
13. L UPPER ARM (L UPP ARM) {5F915696-1B08-32E9-6B95-49246117D660}
14. R UPPER ARM (R UPP ARM) {3EF20E8E-F68C-E4F2-8F8C-FAD3B44DDAA4}
15. L WRIST (L WRIST) {7BCAFACE-FF3B-46DA-99F9-6437A90B7126}
16. R WRIST (R WRIST) {FDD9FE1D-F8A2-09AD-82D7-AA478693DD01}

UNITS:

=====

TERM: DIALYSIS DRESSING TYPE  
ID: {A07D4532-93DF-178D-EB55-070FCF1C570F}

=====

Page: 144

Abbreviation..... DIAL DRESS  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Dialysis Dressing Type  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
1. TRANSPARENT (TRANSPAREN) {0E73D24D-0EFC-08F9-196A-D30DD1BB3228}
  2. PRESSURE DRESSING (PRESS DRES)  
{6D9DB9EF-85EB-8321-03FC-879AA4C5615B}

UNITS:  
-----

=====

TERM: DIALYSIS RATIO

Page: 145

ID: {F0CDCBBA-E65B-FE01-2E15-640FADF9DF20}

=====

Abbreviation..... DIAL RATIO  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Dialysis Ratio  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

-----

1. 1:1 (1:1)	{5C07BBD0-FDB9-0CD7-F2C2-616913E0281B}
2. 1:2 (1:2)	{287CDCE0-36F3-65F6-E6A5-9EF7CF8976AE}
3. 1:3 (1:3)	{062B4F93-8F48-7373-F679-11F57CE98BB3}
4. 1:4 (1:4)	{3F1B75FF-C168-D413-0BDF-6581E7A4EEAD}
5. 1:8 (1:8)	{94790C4C-3AF2-6F4C-58B8-87021CA97A53}

UNITS:

-----







Abbreviation..... DIASTOL PR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Diastolic Pressure
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

LOCATION: 4 entries

- 1. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
2. L LEG (LL) {B088C809-E528-B785-21C2-14E45E1FD8AF}
3. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
4. R LEG (RL) {C5E40A7B-6FC4-267D-8CE8-3571D21D81FD}

METHOD: 5 entries

- 1. CUFF (Cu) {E59A2543-1860-1E1C-BEF1-6E7F6B6FBC5E}
2. DOPPLER (Dop) {52683C19-F9DC-876B-CA3F-53E7C864CF04}
3. ARTERIAL LINE (AL) {D4B3F4ED-8B33-431E-4164-BD165E823776}
4. NON-INVASIVE (N-I) {F50ABFED-5FD7-FC95-20FE-1E3038C9F299}
5. INVASIVE (I) {5DB8746E-5446-F47A-883C-A44B91275A6E}

POSITION: 3 entries

- 1. LYING (Ly) {DE28D2BF-6F8D-64C3-0092-4CEB5E0902B9}
2. SITTING (Sit) {DD4DCDF1-0A22-2CC5-404F-0E472C73961C}
3. STANDING (St) {0776AE42-C370-B7B3-9B4D-9C914BA8F231}

UNITS:

MILLIMETERS OF MERCURY 0 > 40 \* 90 < 300 Dec Prec: 0

Abbreviation..... DIET  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Type of diet ordered  
Display Name..... DIET  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 46 entries

-----

1. 1800 CALORIE CONTROLLED (1800 CAL C)  
{C564D7FE-EC2C-401E-A77D-79CC5FC2B2F8}
2. 1800 DIABETIC (1800 DIAB) {FC71C7D7-F24E-4316-A455-8674804BE092}
3. 2000 CALORIE CONTROLLED (2000 CAL C)  
{DB2FF9DC-298E-4881-94D4-47A7AF0E2BD7}
4. 2000 DIABETIC (2000 DIAB) {DF3056C9-3AEA-4F4E-B758-E35590080AFB}
5. 2200 CALORIE CONTROLLED (2200 CAL C)  
{77F7C138-B2CB-4D76-8C0B-D68F317D4DE0}
6. 2200 DIABETIC (2200 DIAB) {1FB3EFB6-B92F-44FB-831E-4641C91C7A8D}
7. 2400 CALORIE CONTROLLED (2400 CAL C)  
{F3EA28B5-FC96-4F34-B2BE-A96A196B38A2}
8. 2400 DIABETIC (2400 DIAB) {2018FB65-3D4E-47E0-90E9-205A45EDF130}
9. BLAND (BLAND) {38930AFE-ECF7-4548-857C-854E1B5BE148}
10. CARDIAC DIET (CARD DIET) {218A883D-32A0-438D-B453-7CBCE02693A3}
11. CLEAR LIQUID DIET (CL LI DIET)  
{D540348D-3E1D-4E42-9EDC-43C56CF56206}
12. DUMPING SYNDROME (DUMP SYND)  
{0173C971-63CB-45A5-9E93-D966EB3B957D}
13. DYSPHAGIA PHASE I (DYSP PH I)  
{8011795B-77A5-4FF4-8EB9-EB6E1AA1FD4D}
14. DYSPHAGIA STAGE II (DYSP ST II)  
{294BF58A-47FD-43A1-8A1E-9E348987B004}
15. DYSPHAGIA STAGE III (DYS ST III)  
{DCF91830-C5CA-47B2-A7C6-4DD971F953EA}
16. DYSPHAGIA STAGE IV (DYSP ST IV)  
{E738569A-5104-4C93-8EFC-47D80765F739}
17. DYSPHAGIA STAGE V (DYSPH ST V)  
{ABE7D85D-1821-4592-B4E8-84670EE484E4}
18. FAT RESTRICTED 50 GM DIET (F RES 50gm)  
{8E97A639-B7A6-4B3B-9BEF-0A46F3B46E43}
19. FECAL FAT TEST DIET (FEC F TEST)  
{BA1E8E4A-7698-4EFB-B1DD-18C444EE82D5}
20. FULL LIQUID DIET (FULL LIQ D)  
{CAC22234-7B69-44C3-8A0E-9A84A684E429}
21. FULL LIQUID DIABETIC (FULL L DIA)  
{B89FAE3C-07EB-4AF6-B123-1153ADADC79F}
22. GASTRIC BYPASS CLEAR LIQUID (GAS BYP LI)  
{6A87E820-7BEA-4636-B102-F9CB83BDD4D9}
23. GASTRIC BYPASS PUREED (GAS BY PUR)  
{9F9EEB7F-7211-45A2-BCDE-B8398596CBF6}
24. GASTRIC BYPASS REGULAR (GAS BY REG)  
{EE831906-EA51-4EE4-AAF2-583DF366A942}
25. GLUTEN RESTRICTED (GLUT RESTR)  
{0805A164-D387-4E67-9D42-BAF1F2543337}



26. HIGH FIBER (HIGH FIBER) {3CAF7EE0-6A2D-41A3-BEDC-C57ECBBA72A2}
27. HIGH PROTEIN (HIGH PROTE) {1595A29D-CA91-45AE-893B-856EB909E4DD}
28. HIGH CALORIE (HIGH CALOR) {FC453BF0-1E7E-4733-9605-AC3425F16205}
29. LACTOSE RESTRICTED (LACT RESTR)  
{78A9066A-48E2-4285-88E3-A180CE771B5A}
30. LIBERAL DIABETIC (LIBRL DIAB)  
{9B66B47B-E77E-4F69-BDFA-5323B28A283D}
31. LISTED DIET (LIST DIET) {24FBBFD9-B426-4DD9-9843-3428A1C91676}
32. LOW CHOLESTEROL (LOW CHOLE)  
{88F0170A-2560-48C6-9604-982453BCFC3B}
33. LOW FIBER (LOW FIBER) {2834D9A0-B6F7-4C21-96ED-1721D32B3B2B}
34. LOW PROTEIN 60 GM (L PROT 60g)  
{BF5B6D8B-01A4-42D5-8551-6E78A17C43F7}
35. LOW TYRAMINE (LOW TYRAMI) {DB067A83-9D0C-4FC2-BAC8-6EC5641CAE42}
36. MECHANICAL (MECHANICAL) {01EA793C-B524-4685-A3EC-60A317CBACC5}
37. NEUTROPENIC (NEUTROPENI) {5CB9620C-34D9-48F7-AAF0-C595D06FC8AE}
38. NPO DIET (NPO DIET) {BEC8F771-B082-40A2-9EE8-553222A25FBF}
39. POST CATH DIET (POST CA DT)  
{B18FEC50-18E1-4EED-8B37-E7DBAA7D4D3B}
40. POTASSIUM MODIFIED 60 MEQ (K MOD 60M)  
{6EC2F10D-CC44-4872-8473-2445A835ABB0}
41. PUREED (PUREED) {99F8AE5F-E390-4583-8A92-CAC13E21BCF4}
42. REGULAR (REGULAR) {3D3F3D4C-ECD0-2056-A0D1-43CC2F2BA2C8}
43. RENAL DIET (REN DIET) {089DCCC2-B12B-4E3B-8D12-FCA5E974A55A}
44. SODIUM CONTROLLED 2 GM (SOD CONT 2)  
{8A1311AF-6C7D-407A-962F-0FAE6C341F9A}
45. SODIUM CONTROLLED 4 GM (SOD CON 4G)  
{26538A2F-5502-4D94-A2A7-D9E0BF825C47}
46. VEGETARIAN DIET (VEG DIET)  
{6B8143BE-843F-4944-9ABB-7175AFE2DF17}

UNITS:

-----

Abbreviation..... DIS ACTI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Discharge Activities
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 12 entries

- 1. ALERT AND ORIENTED (AL AN OR) {C031F776-36DF-4570-AFAF-A30B8BFA635B}
2. VITAL SIGNS STABLE (VI SI ST) {E6ECD067-54A0-4214-A5F8-43D08C9FEAD6}
3. STEADY GATE OR BASELINE (ST GA BA) {0D47B2CE-88A7-4FE5-964B-7A8107A5DB91}
4. TRANSPORT FACILITATOR PRESENT (TR FA PR) {645BF269-6B23-4377-95ED-71291F075B25}
5. TOLERATING FLUIDS/NOURISHMENT (TO OR FL) {E107F272-CD27-46B8-8C50-F9FF776DDCD9}
6. IV DISCONTINUED (IV DISC) {FDF951E7-5B72-4166-969B-073D5EC74800}
7. VOIDED (VOIDED) {7F23E601-A1E5-4016-8A12-06B1099041CF}
8. DISCHARGE MEDS PLANNED (DI ME PL) {04E9DF8B-8807-47A5-AE88-180F431F3BDE}
9. MED RECONCILIATION DOCUMENTED (ME RE DO) {9ABF3E2E-4D63-43EF-A80D-9DCF5C820226}
10. DISCHARGE ORDER ON CHART (DI OR CH) {F7267E73-07BC-44FB-AB0F-FAA07F6E2C1B}
11. GIVEN DISCHARGE INSTRUCTION AND VERBALIZES UNDERSTANDING (DI IN UN) {C9A61418-88B2-42C4-B5F4-FC3CF76AB39D}
12. DRESSINGS CHECKED (DRESS CH) {21014F53-29CE-4914-8760-6976D8DBA742}

UNITS:

Abbreviation..... Dis Vis
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description.....
Display Name..... Distance Vision
Help Text.....
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 21 entries

- 10. 20/15 (20/15) {0DB4E45B-A96C-4D11-B59E-A89C9F531146}
20. 20/20 (20/20) {FE086353-09B2-4314-84E9-FA909253AAED}
30. 20/25 (20/25) {4A5CBDA4-00D2-479E-A6C8-CCD8B40AECFB}
40. 20/30 (20/30) {667B5FCD-A85D-4054-BFA9-3DA9DDD21306}
50. 20/40 (20/40) {048C6D5D-B7BF-4F2B-BFCC-BB4328A333CF}
60. 20/50 (20/50) {C679C160-D93A-431C-8F28-800D64A019D5}
70. 20/60 (20/60) {AEAA52C2-1450-45B4-A512-73E339718B02}
80. 20/70 (20/70) {B0279832-605C-4441-B568-C131636EB634}
90. 20/80 (20/80) {60F51F22-C739-4BD1-A77E-D224517F861F}
100. 20/100 (20/100) {ECB85701-8B21-407A-9E0B-C015834E984B}
110. 20/200 (20/200) {C3CAF0D6-22A6-4FD5-B8BB-55771F416D4E}
120. 20/300 (20/300) {E6976753-667B-418B-A6E3-0A7CD010CAC0}
130. 20/400 (20/400) {AE8E4C3E-D031-4047-B44E-06F4EB29C73B}
140. 20/800 (20/800) {1303A571-B57C-4757-AC5B-286E6B09343A}
150. CF: COUNTING FINGERS (CF) {3AFD9EDB-98A2-498A-AD95-F5285130CC03}
160. CSM: CENTRAL STEADY MAINTAINED (CSM) {43085BCB-4264-4E4A-BC17-A9C252A89BA8}
170. LP: LIGHT PERCEPTION (LP) {3B5C4FA1-6770-4347-B5FB-2B0554D71339}
180. LPP: LIGHT PERCEPTION WITH PROJECTION (LPP) {29ED35AA-49C5-4835-AF43-17173333A2BB}
190. NLP: NO LIGHT PERCEPTION (NLP) {F710E3C5-9B5A-4826-876D-B5675F9AA67E}
200. HM: HAND MOTIONS (HM) {5180AD3D-4377-435B-9F47-FE59F3E3501B}
210. PROSTHESIS (PROSTHESIS) {CCD0EF23-FE9C-4B20-942E-9CF1CD7CF6E9}

LOCATION: 2 entries

- 10. RIGHT EYE (Right Eye) {56F81421-F1C7-4FD4-9757-E0EDEF26626E}
20. LEFT EYE (Left Eye) {D1B20983-6076-47E6-9578-2EFB3CB6A47B}

UNITS:

=====

TERM: DO NOT PROCESS

Page: 153

ID: {4615254C-EC67-46D0-BF70-9A54BEB4B32D}

=====

Abbreviation..... DNP  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... String  
Description..... Do not process  
Display Name..... Do Not Process  
Help Text..... Do Not Process  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

Abbreviation..... DORFL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Dorsi Flexion  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |                |  |
|----------------|--|
| 1. ABSENT (a)  | {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0} |
| 2. PRESENT (p) | {384DC20F-1625-CE1E-1C85-A3B3924BB2A0} |

UNITS:  
-----

=====

TERM: DRAIN SUCTION CM H20  
ID: {CB2C5928-71E5-429F-9A0E-3FFBF21F62FF}

=====

Page: 155

Abbreviation..... Drain Suct  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description.....  
Display Name..... Drain Suction cm H20  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

CENTIMETERS H2O            0 > 10 \* 20 < 30            Dec Prec: 0

=====

TERM: DRAIN SUCTION TYPE  
ID: {550360E9-AFFD-43DE-9142-F3AF7B95185D}

=====

Page: 156

Abbreviation..... Drain Type  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... Drain Suction Type  
Help Text..... When "Number" is selected, enter  
numeric value as a "Comment"  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 8 entries

- 
1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
  2. INTERMITTENT (INTERMITTE) {DD09ED7C-B977-412F-A678-5130F77D239E}
  3. CONTINUOUS (CONTINUOUS) {676CD2DA-F882-4BE5-99BB-510EA7669749}
  4. BULB (Bulb) {2390453E-3353-40EB-B2D7-EF770B67F037}
  5. LOW (l) {B6CE52A4-C1DD-FDE3-F87E-8E0C98FCC4AC}
  6. MEDIUM (MEDIUM) {62C94066-A7B6-CAA8-064B-2AF7C0DA6323}
  7. HIGH (h) {D88CB5EC-2AD4-2F80-2D43-40ECB0B7D390}
  8. NUMBER: SPECIFY (Specify) {397D773F-78F9-4E69-9A84-B518109A95AB}

UNITS:  
-----

Abbreviation..... DRAIN CHAR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Drainage Characteristics
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 15 entries

- 1. THIN (THIN) {CB8D8B5C-E824-D030-1241-D8EB814671E9}
2. THICK (THICK) {1554A49E-1468-46B2-F851-531DDB1D1A98}
3. YELLOW (YELLOW) {90DC31AD-6211-9F25-4237-79F7122A2763}
4. AMBER (AMBER) {90A13048-4C35-E01D-3482-64D782001D4A}
5. WHITE (WHITE) {AB04EF54-5246-EF3A-E1D9-02E9E24369DD}
6. BROWN (BROWN) {A1C32FBC-7B4D-707C-CA58-92A2E201368E}
7. GREEN (GREEN) {3D1F1553-56FC-68CE-79A9-E4F6357B6AAF}
8. RED (RED) {515D7CFB-3527-30E7-EC4B-B607C272DFA0}
9. BLACK (BLACK) {7A55AFB9-6861-4918-0746-3D67D47FE5C3}
10. CLAY (CLAY) {A4E4D16C-B45C-0934-ED22-8BF640080108}
11. ORANGE (ORANGE) {77B1B667-212D-88BA-D8E6-5C22F6ACDC9C}
12. TAN (TAN) {AB994ED7-8916-FAE0-0E97-CD4E0EF80ED1}
13. PURULENT (PURULENT) {86C01494-743A-C6A2-8883-01CA5200C79B}
14. PINK (PINK) {1A7FAF84-D69C-4EEB-7BC5-D21A598A1F4D}
15. CLEAR (c) {8981C365-2B46-704F-C9F4-975DF1C23C2D}

UNITS:



Abbreviation..... DVT PROPHY  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... DVT Prophylaxis devices  
Display Name..... DVT Prophylaxis  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
1. ANTI-EMBOLUS STOCKINGS KNEE HIGH (ANT EMB KN)  
{4FBFF7A8-D719-98B9-4A46-D76D7BB5213E}
  2. ANTI-EMBOLUS STOCKINGS THIGH HIGH (ANT EMB ST)  
{9B32D3F3-47F1-886F-A59A-D8AB38A34754}
  3. SEQUENTIAL COMPRESSION DEVICE (SCD)  
{29D657DB-6679-8A55-0E33-235B5DD349F2}
  4. ELASTIC WRAP (ELAST WRAP) {9EB52FA5-32D8-823F-2E3C-83375DA62342}

UNITS:  
-----

ID: {64C54E39-9541-4FD2-8F3F-A32053F6889F}

=====

Abbreviation..... EDEMA AMOU  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Edema Amount  
Display Name..... Edema Amount  
Help Text..... 0=Absent, +1=Slight (0-2 mm),  
+2=Moderate (3-5 mm), +3=Marked (6-8  
mm), +4=Very Marked (more than 8 mm)  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Void.....

DATA VALUE: 5 entries

- 
1. 0 - ABSENT (0-a) {75F14311-F113-963C-4F05-B8AE17F044A9}
  2. +1 = SLIGHT (+1 SLIGHT) {7C0245BA-B261-4B68-A84C-D737CA634230}
  3. +2 = MODERATE (+2 = MODER)  
{CCF23A60-068C-443F-AFD0-2D56EEE0446C}
  4. +3 = MARKED (+3 = MARKE) {41ACC586-0C68-4780-8B2D-F36B21FCEDB2}
  5. +4 = VERY MARKED (+4 = VY MA)  
{55EF9D42-B4F8-4468-BA1B-98BCD4075E58}

UNITS:

-----

Abbreviation..... EDEMA
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Observation Of Patients Edema
Display Name..... Edema Severity
Help Text..... 0=Absent, +1=Slight (0-2 mm),
+2=Moderate (3-5 mm), +3=Marked (6-8
mm), +4=Very Marked (more than 8 mm)
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

- 1. 0 - ABSENT (0-a) {75F14311-F113-963C-4F05-B8AE17F044A9}
2. +1 = SLIGHT (+1 SLIGHT) {7C0245BA-B261-4B68-A84C-D737CA634230}
3. +2 = MODERATE (+2 = MODER) {CCF23A60-068C-443F-AFD0-2D56EEE0446C}
4. +3 = MARKED (+3 = MARKE) {41ACC586-0C68-4780-8B2D-F36B21FCEDB2}
5. +4 = VERY MARKED (+4 = VY MA) {55EF9D42-B4F8-4468-BA1B-98BCD4075E58}

LOCATION: 28 entries

- 1. GENERALIZED (GENERIZED) {17073C6C-8691-CADC-2238-CBAAB0A5C648}
2. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
3. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
4. L HAND (L HAND) {639DF25A-19E3-B504-D323-4B3F6ED39A24}
5. R HAND (R HAND) {A54D3A80-8537-D079-9CCE-7DC56346D132}
6. BOTH HANDS (BOTH HANDS) {63AA2685-CB1D-4680-B2F4-2199A22D1EC6}
7. TORSO (TORSO) {EA1AD477-8BFF-E575-23C1-B9CFAF2B9986}
8. ABDOMEN (ABDOMEN) {109FE041-2521-6EF2-7F95-1FD1587796AB}
9. LEFT NECK (L NECK) {40522B5F-F23E-EADB-AE64-8EDE3B9E1155}
10. R NECK (R NECK) {7E4D84D3-535B-6164-C6C1-837C6F997A0D}
11. BILATERAL NECK (BILAT NECK) {744BF4A2-C55E-4D55-8711-58700DD1C446}
12. L FACE (L FACE) {06F00A70-F475-CD40-665D-AF0569A1E75F}
13. R FACE (R FACE) {D1B15677-EF6F-0125-3601-68C1FE6386CE}
14. FACE (FACE) {C638EA14-7903-6A94-496A-C6560D51BEE3}
15. L FOOT (L FOOT) {2EC3FBEC-4EA2-9E7A-17D4-FDC2F9D28046}
16. R FOOT (R FOOT) {7D9F5FF6-232C-2488-4504-A721FF06928F}
17. FEET (FEET) {944CAA0E-5165-A6AF-BE3D-85455440BF40}
18. L KNEE (L KNEE) {3A9E30B4-09A9-967A-58EF-69A01EC4F859}
19. R KNEE (R KNEE) {BCD18D53-07B6-2C41-64FC-3A07A7E979E7}
20. KNEES (KNEES) {247D6621-3A85-F1F7-59BF-1FB0361B3C21}
21. SACRAL (SACRAL) {0CF8A8A7-AB35-B803-4544-0C91645D9D1C}
22. PERIORBITAL (PERIORBITA) {DBAC362B-6278-EDB4-7E28-3180321A38B2}
23. L EYE (L EYE) {D714A0CA-EAD6-EC35-DA54-FF027C23FF45}
24. R EYE (R EYE) {258F0C1D-67BD-ECF8-3152-EB289DEF4681}
25. L LEG (LL) {B088C809-E528-B785-21C2-14E45E1FD8AF}
26. R LEG (RL) {C5E40A7B-6FC4-267D-8CE8-3571D21D81FD}
27. INCISIONAL SITE (INCISON SI) {EB36A66D-751D-4031-9F68-0BB03E267321}
28. BOTH LEGS (BOTH LEGS) {5A15A231-4695-4CF3-AD0F-2E7950208967}

UNITS:

-----

Abbreviation..... EDU BARR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Barriers
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 15 entries

- 1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
2. HEARING (HEARING) {D1A8AC66-1A5F-4F7D-AA89-529AF7A4AA11}
3. LANGUAGE (LANGUAGE) {0666538E-4F44-4D34-A5EA-69CF2F1F86FA}
4. LIMITED ATTENION SPAN (LIM ATT S) {51553020-0CB9-483E-AC1C-15CA6D9380A3}
5. MEMORY (MEMORY) {E30C1930-B98B-4382-B892-34E9C6008D1F}
6. PAIN (PAIN) {F2F23EE7-0A06-46B8-A496-B7A4651A551F}
7. SEDATED (s) {29D9ADD6-2C1F-F6D5-F4D0-B2BB86416997}
8. VISUAL IMPAIRMENT (VIS IMPAIR) {7B2107CE-002D-43BE-9102-FCD5B9FAD617}
9. ANXIOUS (a) {381BDB41-7E8A-6FDB-841B-2ECC1889B9EF}
10. CONFUSED (CONFUSED) {31DAAA60-BE89-49A7-A10D-EB0DA82F12C0}
11. COMMUNCATION BARRIER (COMM BARRI) {1A24946D-D966-4FAB-B6EA-940FA9EF60D3}
12. IMPAIRED LISTENING SKILLS (IMP LIS SK) {9CE0E7E7-79E5-4910-918C-2022C9830171}
13. EMOTIONAL (EMOTIONAL) {626C22E6-EF2A-489F-AA79-D27541A0BC45}
14. UNABLE TO ASSESS (UNA TO AS) {4D11A5E0-1CE8-4390-9EDD-ECEACCAA991B}
15. UNABLE TO READ (ILLITERATE) {6D0C2C48-BDD2-489B-94C8-49E6DAC93F2E}

UNITS:





Abbreviation..... EDU TOPIC
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Topic
Help Text..... Patient Education
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 50 entries

- 1. ACTIVITIES OF DAILY LIVING (ADL) {12F2CA66-3DCE-443D-883C-119CE4838B4A}
2. ADVANCED CARE PLANNING (ADV CAR PL) {F0ED0F8A-1AD6-4F3E-A4C8-4FB915474FE7}
3. ANALGESIA (ANALGESIA) {5452ED11-25E1-437B-8048-3BB3E03BE3B9}
4. CLINICAL ASSESSMENT (CLIN ASSME) {803BEA4B-E22F-48C2-B20E-D5B3C8C7084C}
5. COMPLIANCE (COMPLIANC) {EF43367E-9898-44C3-8AC8-F5BD28051F1D}
6. DAILY PLAN (DAILY PLAN) {B29357A5-4294-4EA0-A0D3-59182F3C94FB}
7. DISCHARGE (DISCHARGE) {ADF67CC2-E05B-4B8D-A233-4AE1C35B2D20}
8. DISEASE PROCESS (DISEA PROC) {E80807AE-6208-421D-BE4C-445C470D32FA}
9. EQUIPMENT (EQUIPMENT) {04C0F805-C7B3-4A9A-91D8-0F32F3265874}
10. HOME HEALTH (HOME HEALT) {F94FB70A-AFFD-4A5E-A829-0AF933B044BB}
11. INFECTION (INFECTION) {DA00DB58-640D-4294-9468-EEF070DBAE8E}
12. INFECTION PREVENTION AND SCREENING (INF PRE SC) {BA58329F-59EC-49F0-BACE-75AA5EC57F52}
13. MEDICATION (MEDICATION) {22E6DEBB-9662-466B-990C-8DA500C0A183}
14. MENTAL HEALTH (MEN HEAL) {ED55FD2C-E2EC-4458-9AE4-FD7F0E6BD7F6}
15. MENTAL HEALTH HOLD (MEN HEA HO) {31560285-7BD5-4E53-BB77-93B244367F42}
16. NEW ADMISSION EDUCATION (N ADM EDU) {1890BC97-FEEB-4BED-A054-3C754C01817B}
17. NUTRITION (NUTRITION) {429A28EF-8A1F-4247-85E8-68D4EFC189C}
18. ORIENTATION (ORIENTA) {10FA5B0E-F4A7-44E4-A482-7C80F145C4D2}
19. PAIN (PAIN) {F2F23EE7-0A06-46B8-A496-B7A4651A551F}
20. POST-OPERATIVE (POST-OP) {9FA19656-849F-4A65-92FC-755A4C9E8056}
21. POST-PROCEDURE (POST-PROC) {0F943BD7-BD13-412E-B4EA-9B5526587340}
22. PRE-OPERATIVE (PRE-OP) {80C418D5-54B1-4CC3-9B19-5653B8D73C45}
23. PRE-PROCEDURE (PRE-PROC) {D2A7B025-9AB9-4C69-B163-D5EBB78B82BE}
24. SAFETY (SAFETY) {AD905A8F-33BD-4EFE-929E-A2B575797402}
25. SPIRITUAL (SPIRITUAL) {9670E2B9-6938-4487-A857-25546775B7A2}
26. TREATMENTS (TREATMENTS) {7A6879B5-F781-498E-9D71-3FA0243C3F21}
27. PRE-SEDATION (PRE SEDA) {684C764E-2BE9-4C99-A2B2-C92D75C73751}
28. POST-SEDATION (POST SED) {B829A249-CB14-4919-86F8-4A8C8D6B7753}
29. ENVIRONMENT (ENVIRO) {C61AA262-F6AF-46AC-A215-0D56FD01B234}
30. ENTERAL TUBE/STOMA CARE (Tube/Stoma) {B91B94EF-CFE6-4628-BEF2-D06E5870D574}
31. ENTERAL TUBE MEDICATION ADMINISTRATION (Med Admin) {68D124A7-69CF-414B-8A47-BC666A373E24}
32. ENTERAL FEEDING FORMULA/WATER ADMINISTRATION (Feed Form) {D2EA8840-1123-4660-AD25-587458BF1E89}
33. ENTERAL FEEDING PUMP (Feed Pump)



- {31FF069C-0060-45C3-898F-20F1A9D1552B}
34. PCA PUMP (PCA Pump) {C6F139AF-C42B-4577-83E1-A12675C31A21}
35. EPIDURAL PUMP (Epid Pump) {26C83111-5FA8-4398-8D3E-787B31A2A9E9}
360. BLOOD GLUCOSE MONITORING (BI Glu Mon)  
{525C8F3C-F8F5-4158-9CDE-3A7B306BADC5}
370. COMMUNITY RESOURCES (Comm Res)  
{0088604F-6A99-4D99-B451-5262734E0D17}
380. DIAGNOSIS/HEALTH STATUS (Diag/Healt)  
{42F128CE-63B7-4AA9-A121-B08C507F9941}
390. FINANCIAL CONCERNS (Fin Concer)  
{A7D1159D-DCA1-47AF-8A07-3EB05493CEFD}
400. GENERAL HEALTH (Gen Healt)  
{D1B589E7-2F59-4A76-8E7B-945F345FA1A0}
410. OSTOMY CARE (Ostomy Car) {BC834BC9-5861-4B94-BA34-FDF5C9FDF979}
420. PATIENT RIGHTS AND RESPONSIBILITIES (Pat Rights)  
{DA69AAEE-DBDC-4CAA-849A-7242F0854385}
430. PERSONAL HYGIENE AND/OR GROOMING (Per Hyg)  
{0DDC2293-90E8-4E0B-A6B1-9C5C946AC831}
440. PREVENTIVE HEALTH REHAB TECHNIQUES (Rehab Tech)  
{2BCABB3F-6688-4BA8-A9E2-D6CA5E38BB7B}
450. SELF-CARE SKILLS (Self Skill)  
{39145E84-85F2-486C-9FF8-27787EC8026A}
460. SUPPORT GROUPS (Supp Group)  
{7D67F542-4210-42FE-B083-7DA9D29055BD}
470. CAREGIVER STRESS (Care Stres)  
{A6A201E7-259C-485E-A35E-9F2F3489B24F}
480. FAMILY NEEDS (Fam Needs) {CEB07929-0A6E-47E9-99CD-BD9333304C9C}
490. OTHER (Other) {4CBA16B3-6BBA-4B8D-AA12-D721AECFE286}
500. NONE IDENTIFIED AT THIS VISIT (None)  
{1B1A62CC-FB5C-4051-95E1-448CA3D0B5F8}

UNITS:

-----

Abbreviation..... EKG  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... EKG Lead  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

LOCATION: 5 entries

- 
- 1. CHEST (CHEST) {8C3C30FE-6E06-C35E-8226-CBF0C94014CC}
  - 2. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
  - 3. L LEG (LL) {B088C809-E528-B785-21C2-14E45E1FD8AF}
  - 4. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
  - 5. R LEG (RL) {C5E40A7B-6FC4-267D-8CE8-3571D21D81FD}

METHOD: 1 entries

- 
- 1. SKIN (ADHESIVE INDICATOR) (SKIN) {C7CE7D84-6FC7-406A-7E56-B080637509DF}

QUALITY: 1 entries

- 
- 1. ACTUAL (A) {5842A4D9-56A1-7A5A-13C9-F481699C3A09}

UNITS:

-----

MILLIMETERS -21.00 > -20.00 \* 20.00 < 21.00  
Dec Prec: 2

=====

TERM: ENTERAL FEEDING TUBE TYPE  
ID: {C0633B19-939C-4213-B7C7-C04DC2ADEF0C}

=====

Page: 166

Abbreviation..... EF TUBE TY  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Type of Enteral Feeding Tube used  
Display Name..... E.F. Tube Type  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 8 entries

- 
1. NASOGASTRIC LARGE BORE TUBE (NG LG BORE)  
{C3FEE16A-9998-455B-BEB8-CC55AF890202}
  2. NASAL INTESTINAL SMALL BORE (N INT SM B)  
{03DCB9EB-9B11-4C3B-869B-B1842BBCC844}
  3. PEG TUBE (PEG TUBE) {970C3755-2B3F-462F-AA87-EA33F846843F}
  4. JEJUNOSTOMY (JEJUNOSTOM) {1ED3EC41-3ACF-42C1-B64A-280300EDBCFC}
  5. GASTRIC TUBE (G TUBE) {65F774BA-F977-46D3-AE41-7F7C74CFF47E}
  6. ORAL GASTRIC TUBE (OG TUBE)  
{C6995D6A-939C-4C31-A3C2-4504CB5597CD}
  7. DUODENOSTOMY (DUODENOSTO) {FA44135D-6806-42E7-B371-8C3912B3582E}
  8. TRANSPYLORIC (TRANSPYOLO) {78FA5E0A-EA21-40F3-8E8D-3B3FCC1ED060}

UNITS:  
-----

Abbreviation..... En Deliver  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... Enteral Nutrition Delivery Method  
Help Text.....  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 7 entries

- 
1. PUMP - CONTINUOUS (Pump- cont)  
{EEE9625C-7EED-4A31-8C13-5E3F0996EDC8}
  2. PUMP - INTERMITTENT (Pump - int)  
{1D0AB4E0-5DF7-446C-BF2C-8336E7C28FCA}
  3. BAG - CONTINUOUS (Bag - cont)  
{F072EF8A-B38B-4D75-82BF-D5EF851784B0}
  4. BAG - INTERMITTENT (Bag - inte)  
{3761CDD3-AA2A-4F2B-9954-865E0A3D8C9A}
  5. SYRINGE (Syringe ) {8ED1EF3E-79AC-4507-A35B-473F31510FB7}
  6. OPEN SYSTEM (Open syste) {2047EDE8-5D72-4541-A23A-81F097A6CA2A}
  7. CLOSED SYSTEM (Closed sys)  
{D076BAB6-C4D1-42F1-9036-80BF4B72E695}

UNITS:  
-----



=====

TERM: ENTERAL NUTRITION RATE  
ID: {75A4F240-8D8E-4F7D-84D3-3DBCE1C20F00}

=====

Page: 169

Abbreviation..... ENT NUT RA  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Rate of Enteral Nutrition  
Display Name..... Enteral Nutrition Rate  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MILLILITERS PER HOUR      -1 > 0 \* 500 < 1000      Dec Prec: 0

=====

TERM: ENTERAL NUTRITION SYSTEM MAINTENANCE  
ID: {9B0AB028-886A-49A5-83BF-A9E4BC787EBF}

=====

Page: 170

Abbreviation..... En Sys Mnt  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... Enteral Nutrition System Maintenance  
Help Text.....  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
1. BAG CHANGE (Bag Chg) {BC58D518-A742-410B-BE33-732C04D303BD}
  2. SYRINGE CHANGE (Syr Chg) {1F052EF1-F62C-4442-8B2B-D15280226C69}
  3. TUBING CHANGE (Tube Chg) {6418C2A2-D244-4C5B-A398-904B99AAACEA4}
  4. STOPCOCK CHANGE (Stop Chg)  
{C4A6E58E-DA4B-478F-823B-D2CFC9860072}

UNITS:  
-----

=====

TERM: ENTERAL NUTRITION TUBE PLACEMENT VERIFIED  
ID: {45513946-76B4-4941-87A4-2737720A16E4}

=====

Page: 171

Abbreviation..... Tube Place  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description.....  
Display Name..... Enteral Nutrition Tube Placement  
Verified  
Help Text..... Can specify method of verification in  
comments  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Void.....

UNITS:  
-----



Abbreviation..... ENT NUT TY
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Type of Enteral Nutrition Used
Display Name..... Enteral Nutrition Type
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 43 entries

- 1. MAGNACAL (MAGNACAL) {C810A618-FB1A-45BB-A540-EADF1414E343}
3. OSMOLITE HN (OSMOLIT HN) {E4186CA6-302F-4513-B125-30B094151423}
4. ULTRACAL (ULTRACAL) {179AE9A5-154D-4CD2-B70B-66E84F5D0932}
6. VITAL HN ENTERAL NUTRITION (VITAL HN) {E45849D4-05EA-4DC1-BD9C-3FF66DD4D292}
7. D5W ENTERAL NUTRITION (D5W) {D6664641-BF43-4157-9A30-8AD8CFE96518}
8. D10W ENTERAL NUTRITION (D10W ENT N) {5049BCEF-CAE5-490D-9EB5-4A204A4D85E1}
9. D50W ENTERAL NUTRITION (D50W ENT N) {528A4339-0261-4EC5-821C-7CE5592666DC}
17. BOOST (Boost) {2D15A3E9-41D2-4099-ADDC-7A1FAD770FF5}
18. BOOST PLUS (Boost Plus) {432C8151-10CE-4804-939D-6B6B609BE54A}
19. ENSURE PLUS (Ensure Plu) {519D1249-B6C0-49ED-BA86-1B27C0F964C9}
20. FIBER STAT (Fiber Stat) {201C8A8B-4417-4BD1-AAE1-FE8CBABF2511}
21. GLUCERNA 1.0 (Glu 1.0) {1C465710-80F2-4C13-B8D5-B1D084151104}
22. GLUCERNA 1.2 (Glu 1.2) {F58F1014-D1C4-44EB-9B01-FB88F2B3D751}
23. ISOSOURCE VHN (Iso VHN) {D68525BF-F922-464E-97DE-1AB767974483}
24. JEVITY 1.0 (Jevity 1.0) {637F4523-0417-4EA4-9991-379C8362EC00}
25. JEVITY 1.2 (Jevity 1.2) {0B3ADECE-FC19-4D90-912B-6EDC3D35C845}
26. JEVITY 1.5 (Jevity 1.5) {ED69D4F6-7EFA-4954-85F0-CCED2BA80A39}
27. JUVEN (Juven) {2112F76A-D581-4C55-94CA-4B4CC68813AD}
28. NEPRO (Nepro) {00FFC7A1-5CDA-45AF-BD50-D77666A23851}
29. NEPRO CARB STEADY (Nepro Car) {6B2B132A-ABCF-4EAC-8B5C-4AFCE885915D}
30. NUTRIHEP (Nutrihep) {ABD98675-EB60-40C7-BCE6-73FD787F6D17}
31. OSMOLITE 1.0 (Osmo 1.0) {D5EE7DA5-707B-4000-A133-0C571824A376}
32. OSMOLITE 1.2 (Osmo 1.2) {CEF57C44-82E5-43F5-8646-F8FBCD4609A0}
33. OSMOLITE 1.5 (Osmo 1.5) {8AF4CD1A-631A-4E95-A631-ACFF775284BA}
34. OXEPA (Oxepa) {1CE54D02-5FBF-43A0-93F4-E5D5DB86FEF6}
35. PERATIVE (Perative) {9F387502-6FE8-4AFA-AF6E-BE4D0CF60947}
36. PIVOT (Pivot) {2E68ED54-57F5-490D-A982-668121F0CF5A}
37. PIVOT 1.5 (Pivot 1.5) {9126F0C6-9716-4C33-8B7C-178D46899D44}
38. PRO STAT SUGAR FREE 64 (ProSta 64) {3622E0A3-8BDF-4ADA-AC4D-20762C10C798}
39. PROMOTE (Promote) {6377C8D3-F6DD-4BD2-877F-DB52525DD65D}
40. PROMOTE WITH FIBER (Pro W/Fib) {33E04B6F-2976-402C-B3B7-B75705430340}
41. PROPASS PROTEIN POWDER (Pro Prot) {2307C41A-92D2-4DD4-BC1F-2E78F7477D98}
42. PRO-STAT (Pro-Stat) {E39E36C7-0A47-459D-BA75-CF6186247B40}
43. PROTENEX WOUND CARE (Protenex) {3B5ED9D2-3DBC-42FF-ACC4-3594EFE33627}
44. PULMOCARE (Pulmocare) {8D683106-B5EE-46E1-B1C6-0C7843EDD15C}

- 45. RESOURCE BREEZE (Re Breeze) {D3A38687-AD62-4E1D-A945-1201168E1CF6}
- 46. RESOURCE DM (Re Dm) {C9D4A461-81F9-4392-98DA-129B00E5807F}
- 47. SUPLENA (Suplena) {F606979A-083E-40C1-9EBB-EA56053EF5FD}
- 48. TWO CAL HN (Two Cal HN) {3CAF568F-4844-4B3F-97A2-86EDE50BB2DE}
- 49. VITAL 1.0 (Vital 1.0) {D3CA653C-419D-4186-8704-E171650062DF}
- 50. VITAL 1.5 (Vital 1.5) {3665F539-7BBD-4AD9-A4D0-0A33506036D7}
- 51. VITAL AF (Vital AF) {D7FCB7B2-EA11-4A75-AD23-1BEDEE03DC58}
- 52. VIVONEX RTF (Viv RTF) {461E7430-CC94-47D1-A90D-9694314C7729}

PRODUCT: 15 entries

- 
- 1. FORMULA 2.0 CAL (2.0 C FORM) {E5FB9590-25AA-435A-B9B4-4812EA15C217}
  - 2. FORMULA D5W ENTERAL NUTRITION (D5W ENT NU) {9EA15A1B-DA1C-49EC-B463-82DD5F125F82}
  - 3. FORMULA D10W ENTERAL NUTRITION (D10W ENT N) {F8FB55C4-6229-427F-987D-3D469E307D0D}
  - 4. FORMULA D50W ENTERAL NUTRITION (D50W ENT N) {19ED9C83-42C1-42AF-AB4D-1A34B621EFB9}
  - 5. FORMULA - PROBALANCE (PROBALANCE) {4EA729FE-08EA-490B-822D-696061E5B55F}
  - 6. GLUCERNA DIABETIC (GLUCERNA D) {C2F2CCFE-6D77-49B4-BF5E-F1A285326440}
  - 7. HIGH PROTEIN FORMULA (HIGH PROTE) {3618C874-3EAE-41CF-BE2B-3DFB960E034F}
  - 8. MAGNACAL (MAGNACAL) {86688633-DFCA-4FBE-99F5-8E0B3FAFE695}
  - 9. MODULAR CUSTOMIZED FORMULA (CUST FORMU) {249D095E-9676-4CC5-8A39-097F96A4C032}
  - 10. NEPRO RENAL FORMULA (NEPR RENAL) {59C9B396-A705-445B-B566-F030DA4B2B79}
  - 11. NUTREN 1.5 (NUTREN 1.5) {B93EF897-B9A1-45C6-AC9A-B386554A3186}
  - 13. OSMOLITE HN (OSMOLIT HN) {3622D5C6-1D7B-4377-94B7-5E33C642C4C5}
  - 14. ULTRACAL (ULTRACAL) {B18E6151-9DC6-4380-899A-E84715BFC61B}
  - 15. VITAL ENTERAL NUTRITION TYPE (VIT ENT NU) {D77B607C-85C4-428E-AAE6-B319F33F4E9F}
  - 16. VITAL HN ENTERAL NUTRITION (VIT NUT HN) {3973F52F-D53F-4D80-928A-BD6154D5C434}

UNITS:

---

Abbreviation..... Asp Color  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... Enteral Tube Aspirate Color  
Help Text.....  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 14 entries

- 
- |                         |  |
|-------------------------|--|
| 1. YELLOW (YELLOW)      | {90DC31AD-6211-9F25-4237-79F7122A2763} |
| 2. AMBER (AMBER)        | {90A13048-4C35-E01D-3482-64D782001D4A} |
| 3. WHITE (WHITE)        | {AB04EF54-5246-EF3A-E1D9-02E9E24369DD} |
| 4. BROWN (BROWN)        | {A1C32FBC-7B4D-707C-CA58-92A2E201368E} |
| 5. GREEN (GREEN)        | {3D1F1553-56FC-68CE-79A9-E4F6357B6AAF} |
| 6. RED (RED)            | {515D7CFB-3527-30E7-EC4B-B607C272DFA0} |
| 7. BLACK (BLACK)        | {7A55AFB9-6861-4918-0746-3D67D47FE5C3} |
| 8. CLAY (CLAY)          | {A4E4D16C-B45C-0934-ED22-8BF640080108} |
| 9. ORANGE (ORANGE)      | {77B1B667-212D-88BA-D8E6-5C22F6ACDC9C} |
| 10. TAN (TAN)           | {AB994ED7-8916-FAE0-0E97-CD4E0EF80ED1} |
| 11. PURULENT (PURULENT) | {86C01494-743A-C6A2-8883-01CA5200C79B} |
| 12. PINK (PINK)         | {1A7FAF84-D69C-4EEB-7BC5-D21A598A1F4D} |
| 13. CLEAR (c)           | {8981C365-2B46-704F-C9F4-975DF1C23C2D} |
| 140. GOLD (Gold)        | {1D1A2586-C13C-41FA-A3AD-39294BEBA83E} |

UNITS:

-----

=====

TERM: ENTERAL TUBE ASPIRATE PH  
ID: {C6DD77B7-2501-4B20-AAF8-FF40554257CF}

=====

Page: 174

Abbreviation..... Asp PH  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description.....  
Display Name..... Enteral Tube Aspirate pH  
Help Text..... pH of enteral tube output  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

pH                    -1.0 > 0.0 \* 14.0 < 15.0   Dec Prec: 1

=====

TERM: ENTERAL TUBE RESIDUAL VOLUME  
ID: {22EBBA65-D612-43B8-A97F-E750C2ADC9D0}

=====

Page: 175

Abbreviation..... Res Vol  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description.....  
Display Name..... Enteral Tube Residual Volume  
Help Text..... volume of enteral tube output  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MILLILITER            -1 > 0 \* 300 < 301      Dec Prec: 0





=====

TERM: EPIDURAL: CONTINUOUS RATE (ML/HR)  
ID: {45F9DF95-EF21-4BE4-8F43-B680A3CDDC47}

=====

Page: 178

Abbreviation..... C Rate Hr  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description.....  
Display Name..... Epidural: Continuous Rate (ml/hr)  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MILLILITERS PER HOUR      0.01 > 0.01 \* 999.99 < 999.99  
Dec Prec: 2



Abbreviation..... Derm Left
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description.....
Display Name..... Epidural: Dermatome Sensory Level -
Left
Help Text.....
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 22 entries

- 1. C1 (C1) {F942693D-F272-4A70-81FF-E750BAF082B9}
2. C2 (C2) {EE1C5D26-12DE-42ED-A46C-63D66BC145C5}
3. C3 (C3) {B53D0E0E-7BF0-4D26-B4DD-D072035E307E}
4. C4 (C4) {A7B5DC47-8E3E-4B23-826C-7F7D5FE0C6BF}
5. C5 (C5) {30EA5569-8E29-4B5E-ABA3-598B7C39A8D4}
6. T1 (T1) {78DB9C11-B0FE-4A08-A3AE-8BEA2CA4D91F}
7. T2 (T2) {7536B154-F3CA-4B21-8CEC-301F204B19A9}
8. T3 (T3) {F8BDC550-D9B9-46CD-9BD1-52E32AA513FD}
9. T4 (T4) {E95970BB-3FDA-4027-AAF0-EE5670642321}
10. T5 (T5) {01B568D4-F9F4-44A7-94A4-EFF830691940}
11. T6 (T6) {0C2C1A00-7CCB-4608-AECB-30BFB52A5342}
12. T7 (T7) {385B254A-F256-4626-9C12-7D8984F4B605}
13. T8 (T8) {68337EB5-902D-4DCB-AFDD-8AC59849F58B}
14. T9 (T9) {06F10D8D-2E9D-454D-AD6B-B750300523DD}
15. T10 (T10) {EB6E695A-95AF-4E9D-AEFA-B1484BABFC4E}
16. T11 (T11) {3C40B21A-965B-48C4-A65D-9A01CF2327F1}
17. T12 (T12) {36FA36A9-FBB2-4F16-A3BA-7E3260DA768D}
18. L1 (L1) {E6FA1F4F-B32B-4523-966C-9AFAD9C44316}
19. L2 (L2) {0E5E7F8E-BA82-499B-A57A-8D32C2113006}
20. L3 (L3) {0CBAF5CD-AABA-4FF0-A25D-01D20CDD53C3}
21. L4 (L4) {94632F85-0152-47C2-9D82-B49C36ABF495}
22. L5 (L5) {87E93C0D-B607-4DD2-83EC-71131C072019}

UNITS:

Abbreviation..... Derm Right
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description.....
Display Name..... Epidural: Dermatome Sensory Level -
Right
Help Text.....
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 22 entries

- 1. C1 (C1) {F942693D-F272-4A70-81FF-E750BAF082B9}
2. C2 (C2) {EE1C5D26-12DE-42ED-A46C-63D66BC145C5}
3. C3 (C3) {B53D0E0E-7BF0-4D26-B4DD-D072035E307E}
4. C4 (C4) {A7B5DC47-8E3E-4B23-826C-7F7D5FE0C6BF}
5. C5 (C5) {30EA5569-8E29-4B5E-ABA3-598B7C39A8D4}
6. T1 (T1) {78DB9C11-B0FE-4A08-A3AE-8BEA2CA4D91F}
7. T2 (T2) {7536B154-F3CA-4B21-8CEC-301F204B19A9}
8. T3 (T3) {F8BDC550-D9B9-46CD-9BD1-52E32AA513FD}
9. T4 (T4) {E95970BB-3FDA-4027-AAF0-EE5670642321}
10. T5 (T5) {01B568D4-F9F4-44A7-94A4-EFF830691940}
11. T6 (T6) {0C2C1A00-7CCB-4608-AECB-30BFB52A5342}
12. T7 (T7) {385B254A-F256-4626-9C12-7D8984F4B605}
13. T8 (T8) {68337EB5-902D-4DCB-AFDD-8AC59849F58B}
14. T9 (T9) {06F10D8D-2E9D-454D-AD6B-B750300523DD}
15. T10 (T10) {EB6E695A-95AF-4E9D-AEFA-B1484BABFC4E}
16. T11 (T11) {3C40B21A-965B-48C4-A65D-9A01CF2327F1}
17. T12 (T12) {36FA36A9-FBB2-4F16-A3BA-7E3260DA768D}
18. L1 (L1) {E6FA1F4F-B32B-4523-966C-9AFAD9C44316}
19. L2 (L2) {0E5E7F8E-BA82-499B-A57A-8D32C2113006}
20. L3 (L3) {0CBAF5CD-AABA-4FF0-A25D-01D20CDD53C3}
21. L4 (L4) {94632F85-0152-47C2-9D82-B49C36ABF495}
22. L5 (L5) {87E93C0D-B607-4DD2-83EC-71131C072019}

UNITS:

Abbreviation..... Epi Meds  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... Epidural: Medication  
Help Text.....  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
- 10. {8BC7F294-7126-4843-91A0-14F805EC4578} (Morphine)  
      {8BC7F294-7126-4843-91A0-14F805EC4578}
  - 20. HYDROMORPHONE (Hydromorph)  
      {8E6E9DB1-C969-4274-9761-C9830A9FCF67}
  - 30. ROPIVICAINE (Ropivicain) {911BE8BC-6737-4DBF-96AF-A501D3D54493}
  - 40. BUPIVICAINE (Bupivicain) {03B908B3-365A-47F6-8D71-8F8FE1D2C16B}
  - 50. FENTANYL (Fentanyl) {C0740659-05FB-4513-BD0F-4F4111B32215}

UNITS:  
-----



Abbreviation..... Side Eff
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description.....
Display Name..... Epidural: Side Effects
Help Text.....
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 19 entries

- 1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
2. NAUSEA (NAUSEA) {6BF33DC3-0300-4643-BB14-27C35692B806}
3. VOMITING (VOMITING) {921F463C-1639-4D13-8E4E-6D3A1B6927F9}
4. ITCHING (Itching) {33304859-7DF4-4D89-86E8-9AA3395CC80A}
5. URINARY RETENTION (Urin Reten) {BDD3A568-00E4-4148-8A5D-00319728022A}
6. CONSTIPATION (CONSTIPAT.) {4A3B0D47-88BF-4FEB-AA0D-7AD7FED2C5C3}
7. HYPOTENSION (Hypotensio) {1AE4E3BB-3B78-42D2-BDAC-83563A0546E1}
8. HEADACHE (Headache) {9338E67C-201C-4066-A734-8BD4B47462BD}
9. BLURRED VISION (Blurred Vi) {1CD10C9B-28E5-49BA-AA6A-0FA3686A9E85}
10. SLURRED SPEECH (Slurred Sp) {6C3D295A-71C9-4DA8-A172-995F431D94C9}
11. CONFUSION (Confusion) {C577D53B-2924-41C6-AAF0-426F5259B36F}
12. POOR BAL (Poor Balan) {959AFDB4-6984-452D-9235-B154409DC0ED}
13. TINGLING MOUTH/LIPS (Ting Mou/L) {E603F979-310D-438B-80A2-1F52E4DE71C6}
14. UNEXPECTED NUMBNESS (Unexp Numb) {6418B135-56AF-4A76-A207-55D7AE0C5364}
15. EXTENDED PERIOD OF IMMOBILITY (Ex Pd Immo) {3DD00E20-6B5E-47EF-B23C-ABA2C1ABAF0A}
16. SYMPTOMS OF MEDICATION ALLERGY (Sym Med Al) {8E73F22D-7583-4F24-B946-091CF8508CED}
17. REPIRATORY DEPRESSION (Resp Depr) {7FE6D54C-1D7E-4028-B84C-E8EC416F9104}
18. CNS DEPRESSION (CNS Depres) {EE349460-8B60-427B-AE98-72007DFD2F6A}
19. OTHER (OTHER) {19094E5B-AEC3-44B0-9BD5-44E79E7EB105}

UNITS:



Abbreviation..... ESI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... EMERGENCY SEVERITY INDEX  
Display Name..... ESI  
Help Text..... EMERGENCY SEVERITY INDEX  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
- |          |  |
|----------|--|
| 1. 1 (1) | {8F5D10A8-FE17-47BC-DACA-6972D12CC00A} |
| 2. 2 (2) | {38E22B79-A2B6-3401-B195-C955F24E1278} |
| 3. 3 (3) | {E714454D-6A67-91F6-2FB7-D8F24F184328} |
| 4. 4 (4) | {D2280C7E-DE50-5A8D-C3FD-5862B69C0E1A} |
| 5. 5 (5) | {188EDCE7-0726-4234-D8D2-92C8052CD221} |

UNITS:  
-----

=====

TERM: ER TRIAGE CONDITION  
ID: {E4848AF5-F7CA-4B23-B436-B902B0F5555A}

=====

Page: 186

Abbreviation..... TRIAGE CON  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Triage Condition  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
- |                        |  |
|------------------------|--|
| 1. STABLE (STABLE)     | {A474E7BC-47EF-4472-A19A-172C441099EC} |
| 2. UNSTABLE (UNSTABLE) | {5841C521-8070-4034-A69D-27DAB2990EB3} |
| 3. OTHER (OTHER)       | {19094E5B-AEC3-44B0-9BD5-44E79E7EB105} |

UNITS:

-----



=====

TERM: ER TRIAGE DISPOSITION  
ID: {4D5AB8A3-0883-4D48-AAF7-73933D411038}

=====

Page: 187

Abbreviation..... TRIAGE DIS  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Triage Disposition  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. EMERGENCY DEPARTMENT (ER DEPT)  
    {D2E511F7-8B80-4928-A64B-3F1EE17395F7}
  2. WAITING ROOM (WAIT ROOM) {13476486-F895-4FFC-AB42-B274827B47DA}
  3. OTHER (OTHER) {19094E5B-AEC3-44B0-9BD5-44E79E7EB105}

UNITS:  
-----

Abbreviation..... ET LOC
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... LOCATION OF ENDOTRACHEAL TUBE
Display Name..... Endotracheal Tube Location
Help Text..... Measured in centimeters from lip or
nose
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 31 entries

- 1. 13 (13) {030242DB-3F82-7517-3DEE-4D50E3E8AE23}
2. 13.5 (13.5) {34D24648-B2D9-9C6C-6D71-F80E977CE185}
3. 14 (14) {EDF1C2B1-123B-011F-9B02-7D71088971C1}
4. 14.5 (14.5) {373C8F8C-5D36-4AD6-6A76-513064669CB5}
5. 15 (15) {FE8C049F-B64C-4991-2F50-3C561BA4825C}
6. 15.5 (15.5) {1BBF5122-2FFA-6B27-445B-3291C3F83A36}
7. 16 (16) {9E4B84FF-430A-39FB-E6AA-E1031814796D}
8. 16.5 (16.5) {7D4EA1D4-5233-760F-FD2B-A70FB0969815}
9. 17 (17) {2E89B9C8-5E80-F3C9-FE66-4396CB63201A}
10. 17.5 (17.5) {FDE6FFD9-00D0-9C55-1A7F-4AA70CE313BA}
11. 18 (18) {F2815231-0FCE-4988-95E9-D9AF46C57086}
12. 18.5 (18.5) {75ABCB85-4D6E-0BBC-157F-CA69B7737402}
13. 19 (19) {478FFB89-4941-1474-04E1-C838B2467439}
14. 19.5 (19.5) {58F0F98D-299A-8C4E-F50F-AE21B6E696D6}
15. 20.5 (20.5) {DB1F50D2-AF99-7BB3-88F5-B58C0BDE7987}
16. 21 (21) {C6188159-9EF1-DA5C-C015-2142A1205C4C}
17. 21.5 (21.5) {A259E466-773F-519C-7B5B-1570D1B03B70}
18. 22 (22) {DA839156-C7A9-76D7-58CD-7EC2B016EF40}
19. 22.5 (22.5) {D402D24F-1EF2-1AA0-DD1C-FA8320AA929A}
20. 23 (23) {8FB6A8B1-B692-BDC7-43B9-0332F33D826A}
21. 23.5 (23.5) {DDB76295-7E4A-FEAC-324D-37878F62771A}
22. 24 (24) {2BE934CF-40B3-2263-9FA7-6AC9DE0C7C8A}
23. 24.5 (24.5) {D815830B-5D94-0975-D007-E4FDD966155B}
24. 25 (25) {42DD0655-142D-6901-D3FE-66D295BBFCD1}
25. 25.5 (25.5) {9D704E9E-FB1D-FB83-2142-EDBF3805E8F7}
26. 26 (26) {181395F6-29FD-5A41-5AD7-BFC9CF564F3B}
27. 26.5 (26.5) {1D464967-6CFA-E29E-548D-EF68AFCDEC61}
28. 27 (27) {64D25A86-6516-0910-B4B5-644B85378876}
29. 27.5 (27.5) {8C7BE365-AFD0-3FFD-B8EC-6EC914763CFC}
30. 28 (28) {66B5A1F5-5396-10C4-643E-3E8CAB6D6CC8}
31. 28.5 (28.5) {C99C9298-D4D2-BC02-36D3-F612B6080A87}

UNITS:

=====

TERM: FALL PRECAUTIONS  
ID: {7912B52F-20BC-9128-883B-08AF0EC48045}

=====

Page: 189

Abbreviation..... FALL PRE  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Fall Precautions  
Help Text..... Fall Precautions Taken  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

Abbreviation..... Fall Prev  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... Fall Prevention  
Help Text.....  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 17 entries

-----

- 10. ORIENTED TO SURROUNDINGS/ROUTINE (Oriented)  
{A8F3F756-E5B5-4B8D-AF21-E45746A64E46}
- 20. INCREASE ROUNDS (Inc Rds) {873180AF-9BB9-44F4-B9E4-80F8BDDF9769}
- 30. NON-SKID SLIPPERS (Non-Skid )  
{D7641B31-2CB9-4284-B505-74B131BF4789}
- 40. PLACED NEAR TO NURSE STATION (Near N/Sta)  
{7790D901-36F0-4DFD-9BF5-8C3A913F30E1}
- 50. NIGHT LIGHT PROVIDED (Night Lite)  
{A6F3E41A-92A1-4AAE-9EFC-84D21E8A32BD}
- 60. PERSONAL ARTICLES AT REACH (At Reach)  
{166DE5BA-67F1-4BDB-9378-062862519BE4}
- 70. ORIENTED TO CALL FOR ASSISTANCE (Call/Asst)  
{1F8390CF-2F88-486A-8583-AF7C7A64F6B5}
- 80. COMFORT MEASURES (Com Measu)  
{A5F09A1F-921E-4469-9652-0A136BCD9C07}
- 90. ESTABLISHED ROUTINE TOILETING SCHEDULE (Toil Sched)  
{A66F5259-5C01-4597-B881-9267DEBC017C}
- 100. AVOIDED LIGHT AND NOISE OVERSTIMULATION (Avoid Lght)  
{135BE5DC-B93B-4A72-97AC-B423A6E81538}
- 110. ELIMINATED ENVIRONMENTAL HAZARDS (Eli Env Hz)  
{78E77129-F53F-491D-8C86-36AE948FCCC6}
- 120. ALARMS ON (Alarms On) {A07F739F-3E78-405D-9012-E43CD84B2CBC}
- 130. HELMET (Helmet) {E65974F8-BCE7-4361-868D-C0B4D780EA21}
- 140. HIP-PROTECTORS (Hip-Prot)  
{EA08514E-84C1-4CF4-B5B1-E5E14EDA82C1}
- 150. FALL RISK SIGN (Fall Sign)  
{E9670720-957E-487E-B0F0-7F7A8AEF69DE}
- 160. FALL RISK ID (Fall Risk) {EF1BE81B-656B-474C-B19E-CCDAD8BB4100}
- 170. OTHER (OTHER) {19094E5B-AEC3-44B0-9BD5-44E79E7EB105}

UNITS:

-----

=====

TERM: FAMILY ORIENTED TO UNIT  
ID: {E1CB7972-95FE-FE7B-9E6E-ECF1FF4699EB}

=====

Page: 191

Abbreviation..... FOTU  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Family Oriented to Unit  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEA649CD} |

UNITS:  
-----

=====

TERM: FETAL HEART TONES

Page: 192

ID: {A2E22A44-E924-ADDE-2B8E-0251666B4DE6}

=====

Abbreviation..... FT  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... From Vitals  
Display Name..... Fetal Heart Tones  
Help Text..... From Vitals  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid..... 4688721

UNITS:

-----

BEATS PER MINUTE      50 > 120 \* 160 < 225      Dec Prec: 0

=====

TERM: FINGERSTICK GLUCOSE  
ID: {805F184D-EF11-446F-A272-B08249A477E9}

=====

Page: 193

Abbreviation..... FINGER GLU  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Fingerstick Glucose Test Result  
Display Name..... Fingerstick Glucose  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

mg/dl                    29 > 65 \* 110 < 601      Dec Prec: 0

=====

TERM: FLATUS

Page: 194

ID: {9999776E-D262-4F64-BDEC-16A7E5885185}

=====

Abbreviation..... FLATUS  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Flatus  
Display Name..... Flatus  
Help Text..... Flatus / Flatulence  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
1. UNKNOWN (u) {1F872397-8635-18A3-E252-AE9B20B88C7D}
  2. PATIENT REPORTS NOT PRESENT (PT RP N PR)  
{AD0EF44D-7231-4E8C-A3AE-0E43DD0C44F2}
  3. PATIENT REPORTS PRESENT (PT REP PRE)  
{1726BC7E-5A82-41AC-B77A-BF62263E1EF6}
  4. PRESENT (p) {384DC20F-1625-CE1E-1C85-A3B3924BB2A0}

UNITS:  
-----



=====

TERM: FLOWSHEET ADDENDUM  
ID: {72F2398C-E7DB-418F-89A4-79E07895295F}

=====

Page: 195

Abbreviation..... FLOW ADDEN  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... String  
Description..... Flowsheet Addendum  
Display Name..... Flowsheet Addendum  
Help Text..... 250 Characters ONLY: Place to enter  
status updates, notes, significant  
events, etc. That are relevant to the  
patient and entire flowsheet.  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

Abbreviation..... FOOTS  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Foot Strength  
Help Text..... Evaluation of foot strength  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
- |               |  |
|---------------|--|
| 1. ABSENT (a) | {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0} |
| 2. WEAK (w)   | {BE2D9305-E2DA-F492-66CF-164EE58B8179} |
| 3. STRONG (s) | {A12154E1-544E-F372-3839-0D44E5BC476B} |
| 4. REFLEX (r) | {F21585F9-5206-A6B3-2BF6-F416E5842366} |
| 6. N/A (n)    | {B3DBE3FD-C5CB-3894-607A-BE13319594D0} |

LOCATION: 3 entries

- 
- |                                |  |
|--------------------------------|--|
| 1. L FOOT (L FOOT)             | {2EC3FBEC-4EA2-9E7A-17D4-FDC2F9D28046} |
| 2. R FOOT (R FOOT)             | {7D9F5FF6-232C-2488-4504-A721FF06928F} |
| 3. BILATERAL FOOT (BILAT FOOT) | {C758847D-468F-479A-B12D-F4E46F0E5E3A} |

UNITS:

-----

=====

TERM: FUNDAL HEIGHT

Page: 197

ID: {EEAB8762-624F-7BA3-4001-114FD229BA69}

=====

Abbreviation..... FH  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... From Vitals  
Display Name..... Fundal Height  
Help Text..... Starting at 24 weeks  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid..... 4688722

UNITS:

-----

CENTIMETER                    15 > 21 \* 43 < 60            Dec Prec: 0

ID: {21F9FD70-0E30-2772-2E3C-D36C9D3ABA9D}

=====

Abbreviation..... GAG  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Gag Reflex  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |                |  |
|----------------|--|
| 1. ABSENT (a)  | {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0} |
| 2. PRESENT (p) | {384DC20F-1625-CE1E-1C85-A3B3924BB2A0} |

UNITS:  
-----

Abbreviation..... Gait  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... Gait  
Help Text..... Assessment of Gait  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
- |                           |  |
|---------------------------|--|
| 10. STEADY (Steady)       | {42CF2FDA-E774-4784-99D3-24545E30C46F} |
| 20. UNSTEADY (Unsteady)   | {7AC26DBA-B7F1-401C-9D4C-5C0F06C87152} |
| 30. SHUFFLING (Shuffling) | {7EA3A898-1D90-472C-BE90-2110330A1ADA} |
| 40. SCOOTER (Scooter)     | {663BC475-C20E-4B8B-B8B5-B16CC79A93C6} |

UNITS:

-----

=====

TERM: GASTRIC RESIDUAL AMOUNT DISCARDED  
ID: {E74268D1-6DB1-40BB-810B-4A33D65F73A4}

=====

Page: 200

Abbreviation..... GAS DISCAR  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Gastric Residual Amount Discarded  
Display Name..... Gastric Residual Amount Discarded  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
MILLILITER            -1 > 0 \* 750 < 10000    Dec Prec: 0

=====

TERM: GASTRIC RESIDUAL CHECKED  
ID: {7CFB9F8C-FDDB-4CB2-B977-E6570D84BAFE}

=====

Page: 201

Abbreviation..... GAS RES CH  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Was the Gastric Residual Checked?  
Display Name..... Gastric Residual Checked  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

=====

TERM: GASTRIC RESIDUAL NOT DISCARDED  
ID: {D93EDA7E-706C-4914-A774-5BE46D4A3085}

=====

Page: 202

Abbreviation..... GAS NOT DI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Gastric Residual Amount Not Discarded  
Display Name..... Gastric Residual Amount Not Discarded  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
MILLILITER            -1 > 0 \* 500 < 10000    Dec Prec: 0



=====

TERM: GASTRIC pH Page: 203  
ID: {5ECE05DB-47C5-43C9-A3CD-0CD31E285BA8}

=====

Abbreviation..... GASTRIC pH  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Gastric pH  
Display Name..... Gastric pH  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

pH                    -1 > 0 \* 7 < 15            Dec Prec: 0

=====

TERM: GCS - EYE OPENING  
ID: {11D9DC0B-FE5C-9374-0904-96896A69E41B}

=====

Page: 204

Abbreviation..... GCS EYE OP  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Eye Opening  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
- 1. 4 - SPONTANEOUSLY (4-s) {4612FC84-FCA8-8C85-50DC-A9A42C195F51}
  - 2. 3 - TO VERBAL COMMAND (3-tvc)  
{EBF4B748-2861-C81F-C0D0-B9FE99FC7A8F}
  - 3. 2 - TO PAIN (2-tp) {49CE3C40-5FC2-002E-1C89-3E5DFA4CCF0F}
  - 4. 1 - NO RESPONSE (1-nr) {26D3168D-40EA-C918-D321-54BED4837311}
  - 5. 0 - ABSENT (0-a) {75F14311-F113-963C-4F05-B8AE17F044A9}

UNITS:  
-----

Abbreviation..... GCS MOTOR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Motor Response
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 6 entries

- 1. 6 - OBEYS COMMANDS (6-oc) {0D4956F4-7E79-8D17-EDF3-274B96BBC89A}
2. 5 - LOCALIZES PAIN (5-lp) {CF13650F-899B-6AE1-5913-DC1F4776E5E4}
3. 4 - WITHDRAWS FROM PAIN (4-wfp) {D808ACF1-30CB-29AB-4957-E91D63C4DF18}
4. 3 - ABNORMAL FLEXION FROM PAIN (3-affp) {9AED8954-5F5D-D2C2-3CAF-4F8EF92009C2}
5. 2 - ABNORMAL EXTENSION FROM PAIN (2-aefp) {9B4951D8-B271-DFC3-0A00-1DC9984039E3}
6. 1 - NO RESPONSE (1-nr) {26D3168D-40EA-C918-D321-54BED4837311}

UNITS:

=====

TERM: GCS - TOTAL

Page: 206

ID: {09118FCE-11BC-2679-68FC-67913B3EAAD0}

=====

Abbreviation..... GCS TOTAL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Total GCS  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

GLASGOW COMA SCALE	2 > 3 * 15 < 16	Dec Prec: 0
--------------------	-----------------	-------------

=====

TERM: GCS - VERBAL RESPONSE  
ID: {B4BA8135-48CB-2930-B5F2-5B93688C07FE}

=====

Page: 207

Abbreviation..... GCS VERBAL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Verbal Response  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. 5 - ORIENTED AND CONVERSES (5-oac)  
    {53AD29B1-BA51-4BDE-4F4C-4268619C93A4}
  2. 4 - DISORIENTED AND CONVERSES (4-dac)  
    {D43F646D-A130-5D06-6DC6-DDB971027231}
  3. 3 - INAPPROPRIATE WORDS (3-iw)  
    {FE000191-243E-4D32-1E1E-546EE2CF88A5}
  4. 2 - INCOMPREHENSIBLE SOUNDS (2-is)  
    {F2261C3D-79A3-055B-702D-446B7A6A3113}
  5. 1 - NO RESPONSE (1-nr)    {26D3168D-40EA-C918-D321-54BED4837311}

UNITS:  
-----

=====

TERM: GENITAL ASSESSMENT  
ID: {139D5B8A-4DF7-4764-A1C9-835C73A13E41}

=====

Page: 208

Abbreviation..... GEN. ASSMT  
Active..... Yes  
Boolean Value False..... Normal  
Boolean Value True..... Abnormal  
Data Type..... Boolean  
Description..... General condition of genitals  
Display Name..... Genital Assessment  
Help Text..... Are genitals normal or abnormal?  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
  - 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEE649CD}

UNITS:  
-----

=====

TERM: GI SYMPTOMS

Page: 209

ID: {2566D442-22EA-4523-8D77-007B858C1EA0}

=====

Abbreviation..... GI SYMPTOM  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... GI Symptoms  
Display Name..... GI Symptoms  
Help Text..... GI Symptoms  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 7 entries

- 
1. NAUSEA (NAUSEA) {6BF33DC3-0300-4643-BB14-27C35692B806}
  2. VOMITING (VOMITING) {921F463C-1639-4D13-8E4E-6D3A1B6927F9}
  3. PROJECTILE VOMITING (PROJ VOMIT)  
{0D2CB23B-B266-40E1-A6FF-20B733655B69}
  4. CRAMPING (CRAMPING) {3EC3E352-5469-4A85-80E4-3AFEB61224C4}
  5. HEARTBURN (HEARTBURN) {5ECEAF5E-6B18-4177-A18D-2DEE87796CC2}
  6. DIARRHEA (DIARRHEA) {353472EF-3499-46B1-AD7C-554FE01D8CE2}
  7. CONSTIPATION (CONSTIPAT.) {4A3B0D47-88BF-4FEB-AA0D-7AD7FED2C5C3}

UNITS:  
-----

Abbreviation..... GI DRAIN C  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... GI Tube Drainage Characteristics  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |                  |  |
|------------------|--|
| 1. THICK (THICK) | {1554A49E-1468-46B2-F851-531DDB1D1A98} |
| 2. THIN (THIN)   | {CB8D8B5C-E824-D030-1241-D8EB814671E9} |

UNITS:  
-----



=====

TERM: GI TUBE DRAINAGE COLOR  
ID: {4D0CB409-25EA-40CB-ABE5-1BB2E0FDE1C3}

=====

Page: 211

Abbreviation..... GI DRA COL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... GI Tube Drainage Color  
Display Name..... GI Tube Drainage Color  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 13 entries

- 
- |                         |  |
|-------------------------|--|
| 1. YELLOW (YELLOW)      | {90DC31AD-6211-9F25-4237-79F7122A2763} |
| 2. AMBER (AMBER)        | {90A13048-4C35-E01D-3482-64D782001D4A} |
| 3. WHITE (WHITE)        | {AB04EF54-5246-EF3A-E1D9-02E9E24369DD} |
| 4. BROWN (BROWN)        | {A1C32FBC-7B4D-707C-CA58-92A2E201368E} |
| 5. GREEN (GREEN)        | {3D1F1553-56FC-68CE-79A9-E4F6357B6AAF} |
| 6. RED (RED)            | {515D7CFB-3527-30E7-EC4B-B607C272DFA0} |
| 7. BLACK (BLACK)        | {7A55AFB9-6861-4918-0746-3D67D47FE5C3} |
| 8. CLAY (CLAY)          | {A4E4D16C-B45C-0934-ED22-8BF640080108} |
| 9. ORANGE (ORANGE)      | {77B1B667-212D-88BA-D8E6-5C22F6ACDC9C} |
| 10. TAN (TAN)           | {AB994ED7-8916-FAE0-0E97-CD4E0EF80ED1} |
| 11. PURULENT (PURULENT) | {86C01494-743A-C6A2-8883-01CA5200C79B} |
| 12. PINK (PINK)         | {1A7FAF84-D69C-4EEB-7BC5-D21A598A1F4D} |
| 13. CLEAR (c)           | {8981C365-2B46-704F-C9F4-975DF1C23C2D} |

UNITS:

-----

=====

TERM: GI TUBE PLACEMENT CONFIRMATION  
ID: {BFA579B5-CEA7-43EF-84C0-C4893EFD0F3A}

=====

Page: 212

Abbreviation..... GI T PLACE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... GI Tube Placement Confirmation  
Display Name..... GI Tube Placement Confirmation  
Help Text..... GI Tube Placement Confirmation  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
1. X-RAY (X-RAY) {C82CB7F8-2456-4A08-9FD2-4E498EB1DB20}
  2. ASPIRATION (ASPIRATION) {68B91085-8E3C-4B94-9F97-6305F5DFD900}
  3. AUSCULTATION (AUSCULTATI) {B9D756DF-74F5-4604-B781-D7BC512FC453}
  4. pH (pH) {104FE39C-B0A2-40B5-A523-B11E5F2BA272}

UNITS:  
-----

=====

TERM: GI TUBE SIZE

Page: 213

ID: {9F0E347C-C292-41A6-B213-340B2ECD716E}

=====

Abbreviation..... GI TU SIZE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... GI Tube Size  
Display Name..... GI Tube Size  
Help Text..... GI Tube Size  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

FRENCH                      7 > 8 \* 18 < 25              Dec Prec: 0

=====

TERM: GI TUBE SUCTION TYPE  
ID: {884A773C-187A-B910-AEC0-59815DEE579E}

=====

Page: 214

Abbreviation..... TUBE SUC T  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... TYPE OF SUCTION APPLIED TO GI TUBE  
Display Name..... GI Tube Suction Type  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 6 entries

- 
1. HIGH (h) {D88CB5EC-2AD4-2F80-2D43-40ECB0B7D390}
  2. LOW (l) {B6CE52A4-C1DD-FDE3-F87E-8E0C98FCC4AC}
  3. CONTINUOUS (CONTINUOUS) {676CD2DA-F882-4BE5-99BB-510EA7669749}
  4. INTERMITTENT (INTERMITTE) {DD09ED7C-B977-412F-A678-5130F77D239E}
  5. CLAMPED (CLAMPED) {B86865B8-078A-4996-B5EC-46A4A99D5A7D}
  6. GRAVITY (GRAVITY) {C0233F70-376D-41CE-9F14-714C7B35244F}

UNITS:  
-----

Abbreviation..... GI TUBE TY  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... TYPE OF GI TUBE  
Display Name..... GI Tube Type  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 10 entries

-----

1. NASOGASTRIC TUBE (NG TUBE)  
{6475E343-A310-48A5-B543-D29573E7684D}
2. GASTRIC TUBE (G TUBE) {65F774BA-F977-46D3-AE41-7F7C74CFF47E}
3. JEJUNOSTOMY (JEJUNOSTOM) {1ED3EC41-3ACF-42C1-B64A-280300EDBCFC}
4. PEG TUBE (PEG TUBE) {970C3755-2B3F-462F-AA87-EA33F846843F}
5. NASAL SMALL BOWEL FEEDING TUBE (NAS SM BOW)  
{A2968C8A-DDBD-4C3B-AFE6-3E2FED41C413}
6. ORAL GASTRIC TUBE (OG TUBE)  
{C6995D6A-939C-4C31-A3C2-4504CB5597CD}
7. SMALL BORE FEEDING TUBE (SM BO TUBE)  
{7E347B8A-14E7-4093-9A5B-56E6AC0D78DF}
8. GASTROSTOMY TUBE (GAST TUBE)  
{EF017A39-745E-4C29-821A-5B25F103DF99}
9. GASTROSTOMY FEEDING BUTTON (GAST F BUT)  
{1627B53B-374B-42B4-BBFE-463EE37B6110}
10. EWALD TUBE (EWALD TUBE) {B83AB4C1-A098-4FAB-A645-61FE1FEADBBE}

UNITS:

-----

=====

TERM: GU BLADDER IRRIGATION FREQUENCY  
ID: {D16099EA-DF97-4481-8B8F-7A25E88D15B1}

=====

Page: 216

Abbreviation..... GU BL IR F  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Bladder Irrigation Frequency  
Display Name..... Bladder Irrigation Frequency  
Help Text..... Bladder Irrigation Frequency  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. CONTINUOUS (CONTINUOUS) {676CD2DA-F882-4BE5-99BB-510EA7669749}
  2. INTERMITTENT (INTERMITTE) {DD09ED7C-B977-412F-A678-5130F77D239E}
  3. PRN (PRN) {36493D81-954B-4C56-B0A2-9988EA785EBD}

UNITS:  
-----

=====

TERM: GU CATHETER CARE PERFORMED  
ID: {F99A31A1-6D6D-4532-81F4-35EA957233D7}

=====

Page: 217

Abbreviation..... GU CA CARE  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Catheter Care Performed  
Display Name..... Catheter Care Performed  
Help Text..... Catheter Care Performed  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

=====

TERM: GU CATHETER SIZE  
ID: {222FF9DB-7D48-4004-B9B9-49B41D8E510B}

=====

Page: 218

Abbreviation..... GU CA SIZE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Catheter Size  
Display Name..... Catheter Size  
Help Text..... Catheter Size  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

FRENCH	11 > 12 * 30 < 31	Dec Prec: 0
--------	-------------------	-------------



Abbreviation..... GU CATH TY  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Type of catheter used  
Display Name..... Catheter Type  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 13 entries

- 
1. INDWELLING CATHETER (INDW CATH)  
{EC16724C-D1EA-5650-58D4-38A89D4456F9}
  2. EXTERNAL CATHETER (EXT. CATHE)  
{C8AEE936-DD91-FE91-4DD6-6F210FB59726}
  3. SUPRAPUBIC (SUPRAPUBIC) {69E0176F-52FC-A3C4-7A4E-AD57DBD6E572}
  4. L NEPHROSTOMY (L NEPHROST)  
{D8FC4C18-081E-A7B9-E394-25F194BE9757}
  5. R NEPHROSTOMY (R NEPHROST)  
{33661137-688D-5611-4368-B1DFD8DB1469}
  6. UROSTOMY (UROSTOMY) {1F9E528D-6D78-54EE-94E6-649C6F2056A3}
  7. STRAIGHT CATHETER (STR. CATH)  
{5F4A8A7F-E720-9458-2360-1DC14E7D05F3}
  8. ILEAL CONDUIT (ILEAL COND)  
{032CB6FE-77D2-702E-A0AE-46176F5CAA57}
  9. L STENT (LEFT STENT) {FD917F0B-CC96-0CB0-BED3-1D7C51A77B69}
  10. R STENT (RT. STENT) {6D0EED3B-1F08-11D2-165D-97A9F58327D5}
  11. FOLEY (FOLEY) {9FE7D3A7-67CF-CF54-F721-2F8035CE2BF1}
  12. 3 WAY FOLEY (3WAY FOLEY) {32408760-F517-4530-AEC0-54FE0A6A062B}
  13. MUSHROOM CATHETER (MUSH CATHE)  
{4FFC451E-14B5-46E4-99C9-99A17BF90AE4}

LOCATION: 16 entries

- 
1. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
  2. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
  3. L INTERNAL JUGULAR (LIJ) {AFBB4FE5-8A11-2047-7900-DADF2FC4430D}
  4. R INTERNAL JUGULAR (R INT JUGU)  
{7ABF6916-D379-198D-6454-F0518A74844A}
  5. L EXTERNAL JUGULAR (L EXT JUGU)  
{6F7B331B-4F83-0B79-1755-F6F673D647D6}
  6. R EXTERNAL JUGULAR (R EXT JUGU)  
{4D78F704-9E3B-907E-A792-4DC4F4CFF170}
  7. L SUBCLAVIAN (LSCL) {71B93594-15F9-AA4A-AEF2-E222EC28C00D}
  8. R SUBCLAVIAN (R SUBCLAVI) {53DEAE83-93AE-123B-14E5-4DE38CD4F752}
  9. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
  10. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
  11. L ANTECUBITAL (L ANTECUBI)  
{C21838CC-BAB5-1D7F-E2A4-A7C663FB28C3}
  12. R ANTECUBITAL (R ANTECUBI)  
{7FB00A65-B4B2-06A6-BD89-98AC112D8D4E}
  13. L UPPER ARM (L UPP ARM) {5F915696-1B08-32E9-6B95-49246117D660}
  14. R UPPER ARM (R UPP ARM) {3EF20E8E-F68C-E4F2-8F8C-FAD3B44DDAA4}
  15. L WRIST (L WRIST) {7BCAFACE-FF3B-46DA-99F9-6437A90B7126}
  16. R WRIST (R WRIST) {FDD9FE1D-F8A2-09AD-82D7-AA478693DD01}

UNITS:

-----

=====

TERM: GU VOIDING Page: 220  
ID: {BE9D4DD5-CF6E-C973-C1F1-D70D9F554E38}

=====

Abbreviation..... GU VOIDING  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Method of urine elimination  
Display Name..... Voiding  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
- 1. CONTINENT (CONTINENT) {3B658506-F42E-BAF5-04C4-5AD49699E0D2}
  - 2. INCONTINENT (INCONTINEN) {BA2E9C32-8D03-8FC9-5E02-EAE2026FDC63}
  - 3. CATHETER (CATHETER) {B0D1AAD7-012D-D3E8-3216-1F71A1EC685E}

UNITS:  
-----

Abbreviation..... HAIR CONDI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Hair Condition  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
1. NORMAL (n) {553C9A97-A3D5-C52F-0E4E-C7D168C99A2B}
  2. ALOPECIA PARTIAL (ALOPE PART)  
{98EDBD8D-205A-4630-A09D-6B4222F5B2A2}
  3. ALOPECIA TOTAL (ALOPE TOTA)  
{ECDCE3C4-81F6-4FA1-BADA-97A20BD9A826}
  4. HIRSUTISM (HIRSUTISM) {A909A4FC-080A-4E9A-867A-0B30D4633602}

UNITS:

-----

Abbreviation..... HAND GRIP  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Hand Grip  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
- |               |  |
|---------------|--|
| 1. ABSENT (a) | {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0} |
| 2. WEAK (w)   | {BE2D9305-E2DA-F492-66CF-164EE58B8179} |
| 3. STRONG (s) | {A12154E1-544E-F372-3839-0D44E5BC476B} |
| 4. REFLEX (r) | {F21585F9-5206-A6B3-2BF6-F416E5842366} |
| 7. N/A (n)    | {B3DBE3FD-C5CB-3894-607A-BE13319594D0} |

LOCATION: 3 entries

- 
- |                                |  |
|--------------------------------|--|
| 1. L HAND (L HAND)             | {639DF25A-19E3-B504-D323-4B3F6ED39A24} |
| 2. R HAND (R HAND)             | {A54D3A80-8537-D079-9CCE-7DC56346D132} |
| 3. BILATERAL HAND (BILAT HAND) | {C782FE56-6117-493F-847F-83173407F0B5} |

UNITS:

-----

Abbreviation..... Handoff Rp  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... String  
Description.....  
Display Name..... Handoff Report  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: HEAD CIRCUMFERENCE  
ID: {33827E3C-5DBB-083C-D8BE-4DFD7D42071F}

=====

Page: 224

Abbreviation..... HC  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... From Vitals  
Display Name..... Head Circumference  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid..... 4536405

UNITS:

-----

CENTIMETER            0 > 1 \* 100 < 101      Dec Prec: 0

=====

TERM: HEAD OF BED ELEVATION  
ID: {97D6AA11-1F1A-1853-23E5-25F2F41E07CF}

=====

Page: 225

Abbreviation..... HOB E  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Head of Bed Elevation  
Help Text..... Value in degrees  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

DEGREES (ANGLE)            -1 > 0 \* 90 < 99            Dec Prec: 0



=====

TERM: HEAT/MOISTURE EXCHANGER PRESENT  
ID: {1AAC7785-C695-44AB-9E0E-96257D66F432}

=====

Page: 226

Abbreviation..... HME PRESEN  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Heat & Moisture Exchanger Present  
Display Name..... Heat & Moisture Exchanger Present  
Help Text..... Heat & Moisture Exchanger Present  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
  - 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEE649CD}

UNITS:  
-----

=====

TERM: HEIGHT

Page: 227

ID: {B440216B-0FB3-1950-7859-7C1BE398FE4A}

=====

Abbreviation..... HT  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Height of Patient  
Display Name..... Height  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid..... 4688724

QUALITY: 4 entries

-----

1. ACTUAL (A) {5842A4D9-56A1-7A5A-13C9-F481699C3A09}
2. ESTIMATED (E) {9EAF26B-7D87-DA66-F991-022D095FB5A0}
3. ESTIMATED BY ARM SPAN (EAS)  
{486CAE5A-BE56-06BA-C65E-1F8131DB1C2D}
4. STATED (Sta) {6362C194-B7F2-D38C-5B72-3BE5DE09438E}

UNITS:

-----

INCH 0.00 > 1.00 \* 96.00 < 100.00  
Dec Prec: 2

CENTIMETER 0.0 > 1.0 \* 2999.0 < 3000.0  
Dec Prec: 1

=====  
TERM: HOMANS SIGN

Page: 228

ID: {996F1DCC-A6EE-4D52-B889-6E0E6185C74C}  
=====

Abbreviation..... HOMAN SIGN  
Active..... Yes  
Boolean Value False..... Negative  
Boolean Value True..... Positive  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Homan's Sign  
Help Text..... Positive is calf pain with  
                        Dorsiflexion and Plantar Flexion  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries  
-----

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
- 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:  
-----

Abbreviation..... I:E RATIO  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Complex  
Description..... Inspiratory Time to Expiratory Time  
Ratio  
Display Name..... I:E Ratio  
Help Text..... Inspiratory Time to Expiratory Time  
Ratio  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Void.....

UNITS:  
-----

=====

TERM: IABP ASSISTED SYSTOLIC PRESSURE  
ID: {21DA5247-9369-4E7E-8089-29087553F2AD}

=====

Page: 230

Abbreviation..... IABP A S P  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... IABP Assisted Systolic Pressure  
Display Name..... IABP Assisted Systolic Pressure  
Help Text..... IABP Assisted Systolic Pressure  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
MILLIMETERS OF MERCURY    0 > 80 \* 150 < 250    Dec Prec: 0

Abbreviation..... IABP BAEDP  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... IABP Balloon Assisted Aortic End  
Diastolic Pressure  
Display Name..... IABP BAEDP  
Help Text..... IABP Balloon Assisted Aortic End  
Diastolic Pressure  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Void.....

UNITS:

-----  
MILLIMETERS OF MERCURY 0 > 40 \* 80 < 200 Dec Prec: 0

=====

TERM: IABP BALLOON SIZE  
ID: {FBA3FB62-E2E2-4246-8D39-D58E5811B1E2}

=====

Page: 232

Abbreviation..... IABP B SIZ  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Size of the IABP  
Display Name..... Balloon Size  
Help Text..... Size of the IABP  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |                  |  |
|------------------|--|
| 1. 34 ml (34 ml) | {25952BB5-50B1-4CC0-9483-8AF891767045} |
| 2. 40 ml (40 ml) | {37A2998F-AB1C-4493-9C7A-FFCB05ED2302} |

UNITS:  
-----

=====

TERM: IABP CATHETER FLUSH  
ID: {D034E8DA-1C70-44D2-B556-F21860827AF9}

=====

Page: 233

Abbreviation..... IABP CAT F  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Was the Intra-aortic Balloon Pump  
Flushed  
Display Name..... IABP Catheter Flush  
Help Text..... Was the Intra-aortic Balloon Pump  
Flushed  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Void.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEA649CD} |

UNITS:

-----



=====

TERM: IABP CIRCULATION ASSESSED  
ID: {59A9ADB9-6630-41B3-AE5E-495AB232E287}

=====

Page: 234

Abbreviation..... IABP CIR A  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... IABP Circulation Assessed  
Display Name..... IABP Circulation Assessed  
Help Text..... IABP Circulation Assessed  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----



=====

TERM: IABP PEAK DIASTOLIC AUGMENTED PRESSURE  
ID: {24E18FB7-94D9-4F2D-9D6D-3CAF3C1F9951}

=====

Page: 236

Abbreviation..... IABP PDAP  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... BP Peak Diastolic Augmented Pressure  
Display Name..... IABP Peak Diastolic Augmented Pressure  
Help Text..... BP Peak Diastolic Augmented Pressure  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
MILLIMETERS OF MERCURY    0 > 80 \* 120 < 200    Dec Prec: 0

Abbreviation..... IABP RATIO  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... IABP Ratio  
Display Name..... IABP Ratio  
Help Text..... IABP Ratio  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
- |              |  |
|--------------|--|
| 1. 1:1 (1:1) | {5C07BBD0-FDB9-0CD7-F2C2-616913E0281B} |
| 2. 1:2 (1:2) | {287CDCE0-36F3-65F6-E6A5-9EF7CF8976AE} |
| 3. 1:3 (1:3) | {062B4F93-8F48-7373-F679-11F57CE98BB3} |
| 4. 1:4 (1:4) | {3F1B75FF-C168-D413-0BDF-6581E7A4EEAD} |
| 5. 1:8 (1:8) | {94790C4C-3AF2-6F4C-58B8-87021CA97A53} |

UNITS:

-----

=====

TERM: IABP SITE CONDITION  
ID: {AA81E246-45C3-258F-4FC8-49CDD592B2D4}

=====

Page: 238

Abbreviation..... IABP CONDI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... IABP Site Condition  
Help Text..... Condition of IABP Site  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 7 entries

-----

2. RED (RED)	{515D7CFB-3527-30E7-EC4B-B607C272DFA0}
3. SWOLLEN (SWOLLEN)	{DCE86B10-C393-8DA9-6518-F2943F0E1A4E}
4. DISCOMFORT (DISCOMFORT)	{77D058C9-1EB4-96C0-C9E3-666C840B884C}
50. PATENT (Patent)	{F0C38BCB-1D11-436C-8B8C-4297C27DDE04}
60. CLEAN (Clean)	{1F187B99-9963-4AD5-98F6-D8D24D2DDD30}
70. OCCLUDED (Occluded)	{25509BE3-F946-4AD7-80E8-5CE8A2E61AD6}
80. OOZING (Oozing)	{C23549D7-3238-4523-BA31-AB71609F01C2}

UNITS:  
-----

Abbreviation..... IABP SIZE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Intra-aortic Balloon Pump Size  
Display Name..... IABP Size  
Help Text..... Intra-aortic Balloon Pump Size  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
- 1. 7.5 (7.5) {7517EE9C-7C8B-4AF0-B60A-88F4E24BFA5E}
  - 2. 8 (8) {518E7CFC-102E-351E-C7AF-FB575E3C8F4D}
  - 3. 8.5 (8.5) {96723567-151D-4B5C-9C65-B3EB6D584272}
  - 4. 9 (9) {58FF2FD4-569F-F454-9F5B-FE6F9E34D8CC}
  - 5. 9.5 (9.5) {252555A1-ADB2-4948-8764-C952F2587DC5}

LOCATION: 5 entries

- 
- 1. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
  - 2. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
  - 3. L AXILLARY (L AXILLARY) {21AF62CA-BEED-2A96-2CD6-F2350A01995E}
  - 4. R AXILLARY (R AXILLARY) {3044F28B-FED2-4069-7E3B-2A691C5771F9}
  - 5. MEDIASTINAL (MEDIASTINA) {D99997D2-61EF-CAC2-A854-D906E3A97447}

UNITS:

-----

Abbreviation..... IABP TIMIN  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... IABP Timing  
Display Name..... IABP Timing  
Help Text..... IABP Timing  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
- |              |  |
|--------------|--|
| 1. 1:1 (1:1) | {5C07BBD0-FDB9-0CD7-F2C2-616913E0281B} |
| 2. 1:2 (1:2) | {287CDCE0-36F3-65F6-E6A5-9EF7CF8976AE} |
| 3. 1:3 (1:3) | {062B4F93-8F48-7373-F679-11F57CE98BB3} |
| 4. 1:4 (1:4) | {3F1B75FF-C168-D413-0BDF-6581E7A4EEAD} |
| 5. 1:8 (1:8) | {94790C4C-3AF2-6F4C-58B8-87021CA97A53} |

UNITS:

-----

Abbreviation..... IABP TI CH  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Was The IABP Timing Checked?  
Display Name..... IABP Timing Check  
Help Text..... Was The IABP Timing Checked?  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
  - 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEE649CD}

UNITS:  
-----





=====

TERM: IABP UNASSISTED SYSTOLIC PRESSURE  
ID: {BA1A966E-51BE-42D7-84F6-2C119AD043F8}

=====

Page: 243

Abbreviation..... IABP U S P  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... IABP Unassisted Systolic Pressure  
Display Name..... IABP Unassisted Systolic Pressure  
Help Text..... IABP Unassisted Systolic Pressure  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
MILLIMETERS OF MERCURY    0 > 80 \* 150 < 250    Dec Prec: 0

=====

TERM: INACTIVE - ARTERIAL LINE SITE  
ID: {F7924F9E-16C8-0C2C-9461-FBC45435191A}

=====

Page: 244

Abbreviation..... A LINE SIT  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Arterial Line Site  
Help Text..... Site of arterial line  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - ARTERIAL LINE TYPE  
ID: {708908E0-9BCD-4DCA-A33D-CB5C6A729526}

=====

Page: 245

Abbreviation..... ALINE TYPE  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Arterial Line Type  
Display Name..... Arterial Line Type  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - AUDIOMETRY  
ID: {FFD29134-8BB2-248E-0412-93C2C08B076F}

=====

Page: 246

Abbreviation..... AUD  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Audiometry  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid..... 4688718

UNITS:

-----

=====

TERM: INACTIVE - BISPECTRAL INDEX  
ID: {4819B466-D0D0-4851-A434-E1691ECA1E65}

=====

Page: 247

Abbreviation..... BIS  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Bispectral Index  
Display Name..... Bispectral Index  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:  
-----

=====

TERM: INACTIVE - BODY MASS INDEX  
ID: {54398789-EC32-4F1D-8E02-913FC96FD198}

=====

Page: 248

Abbreviation..... BMI  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... BMI  
Help Text..... Body Mass Index  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - CARDIAC AV FISTULA BRUIT  
ID: {A0E99499-8021-47FF-AEA2-EB33B9648B0E}

=====

Page: 249

Abbreviation..... AV FIST BR  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... AV Fistula Bruit  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----



=====

TERM: INACTIVE - CATHETER TYPE  
ID: {01941D20-3FD7-A4BF-6108-712C7EBC33A1}

=====

Page: 250

Abbreviation..... CATHETER T  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Catheter Type  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - CENTRAL LINE SITE  
ID: {D6FC875B-CC24-2981-B85E-07FF9CD798AB}

=====

Page: 251

Abbreviation..... C LINE SIT  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Central Line Site  
Help Text..... Site of Central Line  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - CHEST TUBE DRESSING INTACT  
ID: {525B2489-75C2-4AB3-9441-A457BCD7B257}

=====

Page: 252

Abbreviation..... CH TU DR I  
Active..... No  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Chest Tube Dressing Intact  
Display Name..... Chest Tube Dressing Intact  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - CHEST TUBE SITE  
ID: {C5224B54-65A5-2010-F9E3-1C4570500A11}

=====

Page: 253

Abbreviation..... CHEST TU S  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... SITE OF CHEST TUBE  
Display Name..... Chest Tube Site  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:  
-----

=====

TERM: INACTIVE - COUGH AND DEEP BREATH EVENTS  
ID: {5F7D2087-018E-4DA4-8B36-02743F65A769}

=====

Page: 254

Abbreviation..... C&DB EVENT  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Cough and Deep Breath Events  
Display Name..... Cough and Deep Breath Events  
Help Text..... Cough and Deep Breath Events  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:  
-----

=====

TERM: INACTIVE - DATE

Page: 255

ID: {5BB40FB6-173F-9656-8D5B-7A2F5A5F3B82}

=====

Abbreviation..... DATE  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Complex  
Description..... Todays Date  
Display Name..... Date  
Help Text..... Date: Year, Date: Month, Date: Day  
                  are all available  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - DATE DAY  
ID: {CD91474C-1EB3-304B-4E3A-B85F6976C1E7}

=====

Page: 256

Abbreviation..... DATE DAY  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... The day of the month  
Display Name..... Date: Day  
Help Text..... The day of the month  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - DATE MONTH  
ID: {D4E1590F-7B12-FD7C-B545-9CF83BA68BAC}

=====

Page: 257

Abbreviation..... DATE MONTH  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Month of the year  
Display Name..... Date: Month  
Help Text..... Month of the year  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----



=====

TERM: INACTIVE - DATE YEAR  
ID: {57029BBB-5904-DE1D-F58D-AF13A890C6ED}

=====

Page: 258

Abbreviation..... DATE YEAR  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... The Year  
Display Name..... Date: Year  
Help Text..... The Year  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - DIALYSIS CATHETER SITE  
ID: {FC099337-E02E-84E6-7338-E423ECD3DDCB}

=====

Page: 259

Abbreviation..... DIAL CTH S  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Dialysis Catheter Site  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - EDEMA LOCATION  
ID: {8C4F762B-FD4A-CB97-326F-74C4A2ED5548}

=====

Page: 260

Abbreviation..... EDEMA LOCA  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Location of identified edema  
Display Name..... Edema Location  
Help Text..... Location of identified edema  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - ENTERAL NUTRITION FORMULA  
ID: {1C7EF762-4D93-438B-8CBD-2CB98E6C8686}

=====

Page: 261

Abbreviation..... ENT NUT FO  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Enteral Nutrition Formula Used  
Display Name..... Enteral Nutrition Formula  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - GI EMESIS APPEARANCE  
ID: {1B314FCB-6B36-46E0-8306-60C361377BAB}

=====

Page: 262

Abbreviation..... GI EME APP  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Emesis Appearance  
Display Name..... Emesis Appearance  
Help Text..... Emesis Appearance  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - HEARING  
ID: {813CCC94-3D64-5093-BC6C-053EFD9948F9}

=====

Page: 263

Abbreviation..... HEARING  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Hearing  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid..... 4688723

UNITS:

-----



=====

TERM: INACTIVE - IABP SITE  
ID: {EA276ABE-52D0-6D4D-1B0A-7E6443BB264F}

=====

Page: 265

Abbreviation..... IABP SITE  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... IABP Site  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----



=====

TERM: INACTIVE - INCISION LOCATION  
ID: {8A9D3611-6984-4FA1-9C16-4DC215F9B4D6}

=====

Page: 266

Abbreviation..... IN LOCATIO  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Incision Location  
Display Name..... Incision Location  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - INCISION SITE  
ID: {9DBE601B-4D61-D280-1F13-0FDE3533B0A4}

=====

Page: 267

Abbreviation..... WOUND INCS  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... SITE OF WOUND INCISION  
Display Name..... Incision Site  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----



=====

TERM: INACTIVE - IV FLUID

Page: 269

ID: {60B33951-8ABD-41B1-A9FA-5EE1F85CAB97}

=====

Abbreviation..... IV FLUID  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... IV FLUID  
Display Name..... IV FLUID  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - IV INSERTION DATE  
ID: {D01454A0-E430-8B54-C7A6-11B04EBAA417}

=====

Page: 270

Abbreviation..... IV DATE  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... IV Insertion Date  
Help Text..... Date IV Inserted  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - IV SITE

Page: 271

ID: {7F252F53-B53D-D8BD-097B-78197EE99FBA}

=====

Abbreviation..... IV SITE  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Site of IV  
Display Name..... IV Site  
Help Text..... Site of IV  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - OR OUTPUT - CHEST TUBE  
ID: {F6ECA901-354B-4C5B-B03F-C17F7F6448B6}

=====

Page: 272

Abbreviation..... OR OUT CHE  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... OR Output - Chest Tube  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - OR OUTPUT - DRAIN  
ID: {1A99C3C7-9FE1-46D3-9A32-57A2C735DFC6}

=====

Page: 273

Abbreviation..... OR OUT DRA  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... OR Output - Drain  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----



=====

TERM: INACTIVE - OR OUTPUT - PUMP LOSS  
ID: {6A424E02-67FF-49F6-9B33-853BE8545B9B}

=====

Page: 274

Abbreviation..... OR OUT PUM  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... OR Output - Pump Loss  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - OR OUTPUT - URINE  
ID: {3ECAAD17-1BA5-4D0A-BEDC-EC512D32EFF9}

=====

Page: 275

Abbreviation..... OR OUT URI  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... OR Ouput - Urine  
Help Text..... contributes to output total value  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - OSTOMY TYPE  
ID: {975E1FD8-521C-4818-9975-98A30B12C1EF}

=====

Page: 276

Abbreviation..... OSTOMY TYP  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Ostomy Type  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - OXYGEN DELIVERED  
ID: {CFD876D1-AC28-F63F-3B6C-680F5B66A1BD}

=====

Page: 277

Abbreviation..... D02  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Oxygen Delivered  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----



=====

TERM: INACTIVE - OXYGEN EXTRACTION RATIO  
ID: {D6C6849F-6134-417F-AD8B-7853F5AE3451}

=====

Page: 279

Abbreviation..... O2ER  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Complex  
Description..... Ratio of Oxygen Consumed (VO2) to  
Oxygen delivered (DO2)  
Display Name..... O2ER  
Help Text..... Ratio of Oxygen Consumed (VO2) to  
Oxygen delivered (DO2)  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Void.....

UNITS:  
-----

=====

TERM: INACTIVE - PA CATHETER INSERTION DATE  
ID: {64ADD596-950E-EDB5-D7DE-38F5F8601D9E}

=====

Page: 280

Abbreviation..... PA INS DAT  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... PA Catheter Insertion Date  
Help Text..... Date PA Catheter Inserted  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - PA CATHETER SITE  
ID: {1E2B676E-24D4-28D6-9EC4-2200F3DA31B6}

=====

Page: 281

Abbreviation..... PA SITE  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... PA Catheter Site  
Help Text..... Site of PA Catheter  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----



=====

TERM: INACTIVE - PACEMAKER DATE INSERTED  
ID: {B395D1F6-2788-CA17-618C-62719C1CAF4A}

=====

Page: 282

Abbreviation..... PDI  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Pacemaker Date Inserted  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - PACEMAKER SITES  
ID: {668C0AE7-061F-C573-A119-8ED8D006F018}

=====

Page: 283

Abbreviation..... PMSI  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Pacemaker Site  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - PAIN LOCATION  
ID: {20836431-E2C3-41CA-BE14-A8FB1F6A8637}

=====

Page: 284

Abbreviation..... PAIN LOC  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Location of patients Pain  
Display Name..... Pain Location  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

Abbreviation..... POD#  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... NUMBER OF DAYS SINCE OPERATION  
Display Name..... Post Op Day #  
Help Text..... Number of days since operation.  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:  
-----

=====

TERM: INACTIVE - PRESSURE ULCER SITE  
ID: {580747C6-F7B1-4CCA-AE24-E82C436746B6}

=====

Page: 286

Abbreviation..... PR UL SITE  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Pressure Ulcer Site  
Display Name..... Pressure Ulcer Site  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----



=====

TERM: INACTIVE - PULSE OXIMETRY  
ID: {9B830A66-0278-20E1-DB38-089A952DF381}

=====

Page: 288

Abbreviation..... PO2  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... O2 saturation  
Display Name..... Pulse Oximetry  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid..... 4500637

UNITS:

-----

=====

TERM: INACTIVE - SKIN DISCHARGE COLOR  
ID: {C4633087-CBC5-4924-ADF6-846D13624D3E}

=====

Page: 289

Abbreviation..... SKI DIS CO  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Skin Discharge Color  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----



=====

TERM: INACTIVE - SKIN DISCHARGE ODOR  
ID: {4E9309BD-6F1E-4D2C-B500-CE724A92A0F6}

=====

Page: 290

Abbreviation..... SKI DIS OD  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Skin Discharge Odor  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - SKIN DRAINAGE CHARACTERISTICS  
ID: {9F12C134-BA23-40C9-8CC9-0AC08D8891CA}

=====

Page: 291

Abbreviation..... SKI DRA CH  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Skin Drainage Characteristics  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - VENTILATOR IN-LINE NEBULIZER  
ID: {12804E1D-9B65-4027-BCD4-1A28C9C57730}

=====

Page: 292

Abbreviation..... VENT I-L N  
Active..... No  
Boolean Value False..... Off  
Boolean Value True..... On  
Data Type..... Boolean  
Description..... Ventilator In-Line Nebulizer  
Display Name..... Ventilator In-Line Nebulizer  
Help Text..... Ventilator In-Line Nebulizer  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - WOUND SITE  
ID: {7BFBDAAF-779E-4C7A-82F3-DA9E0E7DD597}

=====

Page: 293

Abbreviation..... WOU SITE  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Wound Site  
Display Name..... Wound Site  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INACTIVE - WOUND TOPICAL TREATMENT TYPE  
ID: {29E59078-A908-496E-ABC7-B2C48A4E6400}

=====

Page: 294

Abbreviation..... WOU TOP TR  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Topical Treatment  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: INCENTIVE SPIROMETRY  
ID: {C1E6301C-CB9C-4909-9500-9AA01C42583D}

=====

Page: 295

Abbreviation..... INC SPIROM  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Incentive Spirometry  
Display Name..... Incentive Spirometry  
Help Text..... Incentive Spirometry  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
- |            |  |
|------------|--|
| 1. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |
| 2. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 3. N/A (n) | {B3DBE3FD-C5CB-3894-607A-BE13319594D0} |

UNITS:  
-----

Abbreviation..... INC DIS CH  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Incision Discharge Characteristics  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 8 entries

- 
1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
  2. BLOODY (BLOODY) {D7904E84-8331-D6DD-F2C0-7F3039C5746A}
  3. SEROUS (SEROUS) {9BDA347E-29E9-80C7-91FA-016192D4257E}
  4. SEROSANGUINOUS (SEROSANGUI)  
{19A8D332-899D-177A-2247-ABCDAAEE3BFA8}
  5. PURULENT (PURULENT) {86C01494-743A-C6A2-8883-01CA5200C79B}
  6. CLEAR (c) {8981C365-2B46-704F-C9F4-975DF1C23C2D}
  7. CLOUDY (CLOUDY) {3CCFCCCC-00FC-6B00-1D58-A77C02FDB401}
  8. VISCOUS (VISCOUS) {0E8019D5-F2DB-46B8-B750-41568DD0179C}

UNITS:

-----

=====

TERM: INCISION DRAINAGE AMOUNT  
ID: {77AAAB12-979C-2624-5CFE-7FD5F2310AD5}

=====

Page: 297

Abbreviation..... I DR AMOUN  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Amount of incision drainage  
Display Name..... Incision Drainage Amount  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. INCISION DRAINAGE AMOUNT NONE (I DR A NON)  
{5D0D9027-4B20-43A4-B437-EA7B77CB2105}
  2. INCISION DRAINAGE AMOUNT SCANT (I DR A SCA)  
{A476E395-46AA-44E7-8A00-CE22C7066BB1}
  3. INCISION DRAINAGE AMOUNT SMALL (I DR A SMA)  
{4479479D-8ADA-4FBF-A5A2-426D49FC94D3}
  4. INCISION DRAINAGE AMOUNT MODERATE (I DR A MOD)  
{9B0487A9-DCB4-42D3-92E6-805C412B5D34}
  5. INCISION DRAINAGE AMOUNT HEAVY (I DR A HEA)  
{D77FCEA0-D247-4DD5-873E-2F95C7FF48A6}

QUALITY: 1 entries

- 
1. ESTIMATED (E) {9EAF26B-7D87-DA66-F991-022D095FB5A0}

UNITS:  
-----





Abbreviation..... INC DRE TY  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... OBSERVATION OF INCISION WOUND DRESSING  
Display Name..... Incision Dressing Type  
Help Text..... Incision Dressing Type  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 15 entries

- 
1. INCISION DRESSING TYPE NONE (INC DR NON)  
{81F07C8D-A463-48A0-A64A-A143CECB5DD6}
  2. INCISION DRESSING TYPE PRESSURE (INC DR PRE)  
{B5823373-E5D8-40A5-954C-6A193D7E7511}
  3. INCISION DRESSING TYPE ANTIMICROBIAL (INC DR ANT)  
{35F5F5C1-B76A-44D2-8576-9D7D3ECA3197}
  4. INCISION DRESSING TYPE GAUZE (INC DR GAU)  
{E4AC9307-0512-4F9D-A254-FE19E654461F}
  5. INCISION DRESSING TYPE TRANSPARENT (INC DR TRA)  
{E198EDAD-AA67-49F9-A065-6BB11FE39ED6}
  6. INCISION DRESSING TYPE ANTIMICROBIAL DISK (INC DR ADI)  
{55A33062-7D6B-4CDF-9131-D4ADC38F6CE3}
  7. INCISION DRESSING TYPE GAUZE WRAP (INC DR GWR)  
{E1DBEC29-1430-4DA0-99DC-49B9389D2440}
  8. INCISION DRESSING TYPE FOAM DRESSING (INC DR FOA)  
{5C2C2D24-6022-4BBF-9065-03FB5971CB3F}
  9. INCISION DRESSING TYPE GAUZE PACKING (INC DR GPA)  
{3561DAEB-8A30-4E16-B920-F7DFE022EBB1}
  10. INCISION DRESSING TYPE HYDROCOLLOIDAL (INC DR HYD)  
{424CCDCA-5196-4BBD-AE8B-0F68DDCC3FCF}
  11. INCISION DRESSING TYPE CALCIUM ALGINATE (INC DR CAL)  
{B20558BA-FB9E-499E-82DF-49616ACE9375}
  12. INCISION DRESSING TYPE HYPERTONIC HYDROGEL (INC DR HYP)  
{25820A63-AD76-4198-BEB1-EF95DDEC5B35}
  13. INCISION DRESSING TYPE HYDROGEL (INC DR HYD)  
{8CBAC6DF-F51E-456B-8F9A-A68E4452812C}
  14. INCISION DRESSING TYPE BALSAM OF PERU (INC DR BAL)  
{A1FB45C6-6315-4282-8D89-DCB0C018E357}
  15. INCISION DRESSING TYPE BISMUTH GAUZE (INC DR BIS)  
{CA88C5DC-3890-49A2-800A-0528CE194BA8}

UNITS:  
-----

=====

TERM: INCISION EDGE  
ID: {3FEA28A9-1C53-47D7-95F9-BFFFB25576C6}

=====

Page: 300

Abbreviation..... INC ED  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Incision Edge  
Display Name..... Incision Edge  
Help Text..... Incision Edge  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. INCISION EDGE INDISTINCT/ATTACHED (INC ED ATT)  
    {AF7E7150-525E-49A4-AC23-042B788E8805}
  2. INCISION EDGE NON-ATTACHED (INC ED NON)  
    {F457F951-F1DC-47AC-9336-AC093E149862}
  3. INCISION EDGE TUNNELING (INC ED TUN)  
    {492CD783-3C32-40F0-AE35-5A89E90409A7}
  4. INCISION EDGE UNDERMINING (INC ED UND)  
    {0F3E968A-B867-4123-8C05-ED6B12D364CD}
  5. INCISION EDGE ROLLED (INC ED ROL)  
    {FFA33CAE-5D8B-4501-8120-A536BE53296F}

UNITS:  
-----

=====

TERM: INCISION ERYTHEMA  
ID: {D87F3C97-622C-43D3-9C7F-55C2C848941B}

=====

Page: 301

Abbreviation..... INC ERYTHE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Incision Erythema  
Display Name..... Incision Erythema  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
  2. LIMITED TO INCISION (LIM TO INC)  
{81A61254-55B6-47E1-8BE2-F4A47F624BC6}
  3. EXTENDS INTO PERIWOUND <1CM (EX IN PER<)  
{0EBCEAE0-E17A-484F-BFE2-51F17AA761CB}
  4. EXTENDS INTO PERIWOUND>1CM (EX IN PER>)  
{1F31FE9B-DA19-434C-8C58-67CDC5385231}

UNITS:  
-----

=====

TERM: INCISION ODOR  
ID: {2584F1F9-CD15-450C-BDDB-D0B2B1924303}

=====

Page: 302

Abbreviation..... INC OD  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Incision Odor  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
1. INCISION ODOR NONE (INC OD NON)  
    {CED90890-6FBC-4A4A-B0B8-6F75F0DD567A}
  2. INCISION ODOR FOUL (INC OD FOU)  
    {ECDC0988-BE82-41E3-8104-A66A5F8A1093}

UNITS:  
-----

Abbreviation..... INCIS PRES
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... Incision Present
Display Name..... Incision Present
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

LOCATION: 85 entries

- 1. L THORACOTOMY (L THORACOT) {313C11A1-559A-4F68-8718-4012416EC189}
2. R THORACOTOMY (R THORACOT) {A18E4BB3-0718-48D6-BD49-720A70C0EF60}
3. L TRAPEZIUS (L TRAPEZIU) {B74CF943-C10F-4A1B-A19A-4975E674A97A}
4. R TRAPEZIUS (R TRAPEZIU) {69753B67-AFA1-4C93-BA78-E266F9894E46}
5. DONOR SITE (DONOR SITE) {F0A65974-7650-4A4B-AC88-0945E37A0AD2}
6. L GLUTEUS (L GLUTEUS) {88C4E547-04CD-4D52-9416-CAEB6D898D21}
7. R GLUTEUS (R GLUTEUS) {B4D32045-68D1-4457-9106-7937887FA6F6}
8. L LOWER ANTERIOR CHEST (L L A CHES) {6073B3C6-AA4D-4BBA-9FF2-32ED40CF92AD}
9. R LOWER ANTERIOR CHEST (R L A CHES) {96BCF731-B661-4455-B5E7-C4AE2509505F}
10. L MIDDLE ANTERIOR CHEST (L M A CHES) {36A83F77-117E-474B-AE34-0B2C02E131AE}
11. R MIDDLE ANTERIOR CHEST (R M A CHES) {2F751A9C-66CA-4C22-ACC9-6E81FA63FA2B}
12. L UPPER ANTERIOR CHEST (L U A CHES) {FE63B58A-DE39-4C03-AE61-E34A5487C767}
13. R UPPER ANTERIOR CHEST (R U A CHES) {67890049-A4A8-4241-BDF3-573729BD6EB5}
14. L LOWER POSTERIOR CHEST (L L P CHES) {2C19F10B-7546-4A1F-A5DC-3867359ED971}
15. R LOWER POSTERIOR CHEST (R L P CHES) {5163180E-FFD2-4BFB-BF8B-4C6B18EFE5D4}
16. L MIDDLE POSTERIOR CHEST (L M P CHES) {5A656939-F82B-4DEC-9F1D-B5FA644E3DA9}
17. R MIDDLE POSTERIOR CHEST (R M P CHES) {2BC487D6-6CE8-499E-94A0-A79EAA8E199B}
18. L UPPER POSTERIOR CHEST (L U P CHES) {6D5674E9-EF04-4EF1-84A0-3F5ED2FB532A}
19. R UPPER POSTERIOR CHEST (R U P CHES) {D8F03FA5-1FD7-46EE-8FA4-17B3B7E62141}
20. MIDLINE ABDOMINAL (MID ABDOMI) {37ED2947-9BDA-4207-94F1-6B8C0FE12ECF}
21. UMBILICAL (UMBILICAL) {8874B300-0751-49AF-8211-CF4C009C68B5}
22. TRANSVERSE ABDOMINAL (TRANS ABDO) {1C10FF98-7B20-45D0-B420-714981832B3D}
23. L ANKLE (L ANKLE) {D5A692EF-9D7C-473C-AF96-7EF1C01FA2AF}

24. R ANKLE (R ANKLE) {4E6DF352-F399-4CE9-908C-48132C3D856B}  
 25. L LOWER ARM (L LOW ARM) {63AA95E8-256D-4390-8CB3-F00A1B3F91CD}  
 26. R LOWER ARM (R LOW ARM) {8DE6E277-EFFF-40DD-9E83-D7295AACB6AA}  
 27. L UPPER ARM (L UPP ARM) {5F915696-1B08-32E9-6B95-49246117D660}  
 28. R UPPER ARM (R UPP ARM) {3EF20E8E-F68C-E4F2-8F8C-FAD3B44DDAA4}  
 29. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}  
 30. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}  
 31. L INGUINAL (L INGUINAL) {627D172E-8F33-48C9-BE89-D036E07D46BF}  
 32. R INGUINAL (R INGUINAL) {395233DF-3B89-49E0-946B-EB5C02E51178}  
 33. L FLANK (L FLANK) {6E453E8D-B456-4BFD-A015-F9C013830F30}  
 34. R FLANK (R FLANK) {F5700ED2-646F-4B75-A55B-111C27FB939C}  
 35. L FOOT (L FOOT) {2EC3FBEC-4EA2-9E7A-17D4-FDC2F9D28046}  
 36. R FOOT (R FOOT) {7D9F5FF6-232C-2488-4504-A721FF06928F}  
 37. L HIP (L HIP) {786CD881-27AD-4136-B9CC-8423CB38A431}  
 38. R HIP (R HIP) {A48CBCAE-5404-4934-9321-ECEF660A5350}  
 39. L KNEE (L KNEE) {3A9E30B4-09A9-967A-58EF-69A01EC4F859}  
 40. R KNEE (R KNEE) {BCD18D53-07B6-2C41-64FC-3A07A7E979E7}  
 41. L THIGH (L THIGH) {16C2312F-1E0E-45E0-B675-160F40C23F57}  
 42. R THIGH (R THIGH) {B7A2D4A0-801E-4578-967A-F1CFB86FB96B}  
 43. L CALF (L CALF) {3DAA7817-936B-444F-8A2D-D3AB09144C10}  
 44. R CALF (R CALF) {0298A452-4901-44BF-A7A8-F42DD3BA8B63}  
 45. L LOWER LEG (L LOW LEG) {7B8E92FA-E065-4D06-9FB6-4234656CFD61}  
 46. R LOWER LEG (R LOW LEG) {C8EA4D71-D97E-46BD-AE9D-C10954C525CD}  
 47. MIDLINE NECK (MID NECK) {FF666713-278D-46D6-8D25-E5D57E78EA08}  
 48. L NECK (L NECK) {C4EA6279-05E2-4195-A638-3794DCBA7B78}  
 49. R NECK (R NECK) {7E4D84D3-535B-6164-C6C1-837C6F997A0D}  
 50. PERINEAL (PERINEAL) {1801441C-D5EF-4775-9E58-E4985496F7E0}  
 51. SPINE CERVICAL (SPINE CERV)  
     {9D28E1BF-D014-4DA4-A832-1F3C25C148E9}  
 52. SPINE THORACIC (SPINE THOR)  
     {2CACBB4F-3EF3-4FBF-AC3E-41D72F0178CB}  
 53. SPINE LUMBAR (SPINE LUMB) {E7E411AB-5C29-4906-B0A9-E8658AFFFF91}  
 54. L SUBCLAVIAN (LSCL) {71B93594-15F9-AA4A-AEF2-E222EC28C00D}  
 55. R SUBCLAVIAN (R SUBCLAVI) {53DEAE83-93AE-123B-14E5-4DE38CD4F752}  
 56. L THORACIC ANTERIOR (L THOR ANT)  
     {C85B5FB9-153B-4A8B-9BE6-5524A572ACAC}  
 57. R THORACIC ANTERIOR (R THOR ANT)  
     {BD3DF757-2C65-4D6A-9C80-5DC95BB41CFB}  
 58. L THORACIC POSTERIOR (L THOR POS)  
     {68E3381F-7605-4786-85A8-523D215A2B97}  
 59. R THORACIC POSTERIOR (R THOR POS)  
     {E91085D2-E665-4EB6-A086-8CF573DECD8F}  
 60. L THORACIC LATERAL (L THOR LAT)  
     {48D26CEB-2E33-4F57-AB3A-F5A46108FFBF}  
 61. R THORACIC LATERAL (R THOR LAT)  
     {55C48294-D3F3-4013-B552-BA3344D9513B}  
 62. L THORACOTOMY (L THORACOT)  
     {313C11A1-559A-4F68-8718-4012416EC189}  
 63. R THORACOTOMY (R THORACOT)  
     {A18E4BB3-0718-48D6-BD49-720A70C0EF60}  
 64. THYROID (THYROID) {6AC68004-FE58-46D8-9A1E-4EEF72AAD71B}  
 65. L PARIETAL (L PARIETAL) {2687ED24-F35E-456C-99C9-F6A8EC4364E0}  
 66. R PARIETAL (R PARIETAL) {28C927CD-1342-483B-B366-0B7EA8ACABAD}  
 67. L ORBIT (L ORBIT) {1D798E35-5095-4BA0-8F9F-FA96668A5239}  
 68. R ORBIT (R ORBIT) {DFFD974A-1684-49AE-8230-8E12411EA62E}  
 69. L MANDIBULAR (L MANDIBUL) {CB9DB373-5F0A-446B-909F-F76392CB3A04}  
 70. R MANDIBULAR (R MANDIBUL) {303A1DE2-CB9A-498F-9CC5-DFA8FBFF10C1}  
 71. L ORAL (L ORAL) {057B787A-8750-4563-80DB-25110CCC1146}  
 72. R ORAL (R ORAL) {9CEEF934-AA3D-4BDE-BC3B-9E275C51CC9C}  
 73. L FACE (L FACE) {06F00A70-F475-CD40-665D-AF0569A1E75F}  
 74. R FACE (R FACE) {D1B15677-EF6F-0125-3601-68C1FE6386CE}  
 75. L EAR (L EAR) {DA9C2DD1-D6FF-4899-A62B-1096FAA36C59}

- 76. R EAR (R EAR) {AE31327D-A03C-488C-8A8E-084E7E6E7FDE}
- 77. L THUMB (L THUMB) {F82F4731-2CC5-4705-91BE-05E40C25B419}
- 78. R THUMB (R THUMB) {B414BD68-4155-45B3-985D-B2F8252A20FE}
- 79. L FINGER (L FINGER) {BB825645-1630-498C-A31B-6728D422EEFC}
- 80. R FINGER (R FINGER) {E12CA0EF-3267-4B93-8643-230E2226C3D6}
- 81. L TOE (L TOE) {C666DE86-E7AE-4AB3-8AE1-91496219E63D}
- 82. R TOE (R TOE) {B65727AB-3205-4DCE-833F-6DA8554636F9}
- 83. L OCCIPITAL (L OCCIPTAL) {A4234F41-CD6F-4CE1-9210-316DB847BDE2}
- 84. R OCCIPITAL (R OCCIPITA) {5DB951DB-CB03-482E-B0DA-CC5EE8FB1A67}
- 85. OCCIPITAL (OCCIPITAL) {5A3E384B-1C87-4A19-A449-DB91AB443FFC}

UNITS:

-----



=====

TERM: INCISION STAPLES INTACT  
ID: {FA75C02B-3698-4D85-A293-0ACC8C1D3252}

=====

Page: 304

Abbreviation..... IN STA INT  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Incision Stapes Intact  
Display Name..... Incision Stapes Intact  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

=====

TERM: INCISION SUTURES INTACT  
ID: {A5C84A03-633B-4A9B-891B-292449B8C409}

=====

Page: 305

Abbreviation..... IN SUT INT  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Incision Sutures Intact  
Display Name..... Incision Sutures Intact  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

=====

TERM: INCISION TENDER TO PALPATION  
ID: {C6ACAC1F-0AEF-4F1B-A487-35A92207AE22}

=====

Page: 306

Abbreviation..... IN TENDER  
Active..... Yes  
Boolean Value False..... Patient Denies  
Boolean Value True..... Appropriate  
Data Type..... Boolean  
Description..... Incision Tender To Palpation  
Display Name..... Incision Tender To Palpation  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

=====

TERM: INFUSION SET CHANGED  
ID: {BA138A48-8B25-40CF-8345-2AD9E12CC840}

=====

Page: 307

Abbreviation..... INF SET CH  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Was the infusion set changed?  
Display Name..... Infusion Set Changed  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEA649CD} |

UNITS:  
-----

=====

TERM: INSULIN PUMP SQ INFUSION ON?  
ID: {0B8DD5F9-E685-4D44-81D1-04D17DE1E5DA}

=====

Page: 308

Abbreviation..... In Pump on  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description.....  
Display Name..... Insulin Pump SQ Infusion On?  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----



- {65CBDCBD-99DB-44FE-9230-D62886208EBF}
19. PLATELETS PHERESIS IRRADIATED (PLATEL P I) {70C7155A-5AB6-4733-B43F-606F8BBB3B4A}
  20. PPN (PPN) {BF402283-3E6F-4DF1-9CC6-4DE35D755C7C}
  21. PRBC (PRBC) {BCE37EFD-78D6-4A13-A5E4-EA278BF18E02}
  22. PRBC - AUTOLOGOUS (PRBC-A) {1580990C-0B35-4DDD-84AA-682B341C9B76}
  23. PRBC - LEUKOREduced (PRBC-L) {13DB3170-6C26-4178-B429-C02AFC541AC7}
  24. PRBC - LEUKOREduced - DIRECTED (PRBC-L-D) {9C54B7F7-FBD9-4FA8-B2A1-893A25A9CD2C}
  25. PRBC - LEUKOREduced - IRRADIATED (PRBC-L-I) {B4431519-0102-49ED-94D6-1738E375C5A4}
  26. PRBC - LEUKOREduced - WASHED (PRBC-L-W) {BB58F35A-1086-4061-B485-FA35DFB79A0D}
  27. SANDOSOURCE PEPTIDE (SAND PEPTI) {B543461C-7316-4750-9918-B30E96AB3C95}
  28. TPN (TPN) {55341B62-E356-4C38-A52D-4F196751E9EB}
  29. WATER (WATER) {ACCC998E-0F64-510D-B843-3D30B2AD67CF}
  30. WBC - WASHED (WBC WASHED) {F1460F94-203C-473D-AC7C-92E8ADFFC0F5}
  31. CRYSTALLOID (CRYSTALLOI) {F8B0C9B5-F386-4C88-896A-092C4A95A031}
  32. PLATELETS (PLATELETS) {3F4FA2EF-DF05-4A6D-839F-F97D3C2AFCF7}
  33. WHOLE BLOOD (WHOL BLOOD) {144EE01A-16B5-4F9B-9A1C-D420D0A9A2D0}

UNITS:

-----

MILLILITER	-1 > 0 * 9999 < 10000	Dec Prec: 0
OUNCE	-1 > 0 * 64 < 350	Dec Prec: 0
CUBIC CENTIMETER	-1 > 0 * 9999 < 10000	Dec Prec: 0

=====

TERM: INTAKE - CBI IRRIGATION TO SUBTRACT (NEGATIVE NUMBER)  
ID: {182D805C-CF2D-4829-8925-6916D027020D}

=====

Page: 310

Abbreviation..... CBI Irriga  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description.....  
Display Name..... Intake - CBI Irrigation to Subtract  
(negative number)  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
MILLILITER                    -99999 > -99999 \* 0 < 0    Dec Prec: 0





=====

TERM: INTAKE - EF/SUPPLEMENT CALORIES  
ID: {58506F9F-9D10-47F6-A195-5CE7D749FFBC}

=====

Page: 312

Abbreviation..... INT EF CAL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Enteral Feeding and Supplement Calories  
Display Name..... Intake - EF/Supplement Calories  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

CALORIE                    -1 > 0 \* 500 < 2001      Dec Prec: 0

=====

TERM: INTAKE - ENTERAL FEEDING TUBE  
ID: {769017E8-B98E-45DC-9706-FE6475B7FF8D}

=====

Page: 313

Abbreviation..... INT EF TU  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Intake - EF Tube  
Display Name..... Intake - EF Tube  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MILLILITER	-1 > 0 * 9999 < 10000	Dec Prec: 0
CUBIC CENTIMETER	-1 > 0 * 9999 < 10000	Dec Prec: 0





Abbreviation..... INT IRRIGA  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Irrigation Intakes  
Display Name..... Intake - Irrigation  
Help Text..... Fluid intake from irrigations  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 2 entries

- 
- 1. BLADDER (BLAD) {E3B8D8AC-2846-B84B-75D9-7F50B55B4ADD}
  - 2. OTHER (OTHER) {AED275BD-8A42-46BE-BCDF-472EF4C31D65}

PRODUCT: 5 entries

- 
- 1. SALINE (SALINE) {19D84001-70DB-47AA-9149-A37B7E0BAEF1}
  - 2. BACTERIOSTATIC WATER (BACTER H2O)  
{EF13BCA6-3AF2-47E9-87A0-82642B529368}
  - 3. DAKINS SOLUTION (DAKINS SOL)  
{51651D24-B909-427B-9FF0-3C2859A63D73}
  - 4. AMPHOTERICIN B (AMPHOTER B)  
{E4ED1488-6A27-4622-A8EE-AB6B53DB86D0}
  - 5. ANTIBIOTIC (AB) {BC1045A3-FF58-41C0-B48F-C2CFD6E522A7}

UNITS:

-----

MILLILITER -1 > 0 \* 9999 < 10000 Dec Prec: 0  
CUBIC CENTIMETER -1 > 0 \* 9999 < 10000 Dec Prec: 0

Abbreviation..... INT BLOOD
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Blood Products given to patient
Display Name..... Intake - IV Blood Products
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

PRODUCT: 18 entries

- 1. CRYOPRECIPITATE (CRYOPRECIP) {5329E737-E059-480D-83F8-1146B091C14E}
2. FFP (FFP) {9DD0ECF5-634F-49CF-A529-A2657240A80C}
3. FFP - THAWED (FFP-THAWED) {AAF452A3-0DCA-4FF4-BF96-08E60D2D2875}
4. GRANULOCYTES PHERESIS (GRAN PHERE) {2B90878D-9C29-4731-9525-D7DD7FCB513B}
5. PLATELETS PHERESIS (PLATELET P) {65CBDCBD-99DB-44FE-9230-D62886208EBF}
6. PLATELETS PHERESIS IRRADIATED (PLATEL P I) {70C7155A-5AB6-4733-B43F-606F8BBB3B4A}
7. PRBC - AUTOLOGOUS (PRBC-A) {1580990C-0B35-4DDD-84AA-682B341C9B76}
8. PRBC - LEUKOREduced (PRBC-L) {13DB3170-6C26-4178-B429-C02AFC541AC7}
9. PRBC - LEUKOREduced - DIRECTED (PRBC-L-D) {9C54B7F7-FBD9-4FA8-B2A1-893A25A9CD2C}
10. PRBC - LEUKOREduced - IRRADIATED (PRBC-L-I) {B4431519-0102-49ED-94D6-1738E375C5A4}
11. PRBC - LEUKOREduced - WASHED (PRBC-L-W) {BB58F35A-1086-4061-B485-FA35DFB79A0D}
12. WBC - WASHED (WBC WASHED) {F1460F94-203C-473D-AC7C-92E8ADFFC0F5}
13. PRBC (PRBC) {BCE37EFD-78D6-4A13-A5E4-EA278BF18E02}
14. PLATELETS (PLATELETS) {3F4FA2EF-DF05-4A6D-839F-F97D3C2AFCF7}
15. PLASMAPHERESIS (PLASMAPHER) {E1214367-8BB0-4FAA-A916-EC707569A9C4}
16. 25 PERCENT ALBUMIN (25 PERC AL) {CBC7E32F-E457-4A7F-8F27-030EBC6ECBD5}
17. 5 PERCENT ALBUMIN (5 PERC ALB) {A2D0F9A0-A9EC-4540-9663-B8B19A189198}
18. Rh IMMUNE GLOBULIN (Rh IM GLOB) {1FE7BD0C-7252-427A-8216-36FCA4637089}

UNITS:

MILLILITER -1 > 0 \* 9999 < 10000 Dec Prec: 0
CUBIC CENTIMETER -1 > 0 \* 9999 < 10000 Dec Prec: 0

=====

TERM: INTAKE - IV CALORIES  
ID: {868852E2-3D29-4007-BB36-62F6D90B1685}

=====

Page: 318

Abbreviation..... INT IV CAL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Calories taken in by IV  
Display Name..... Itake - IV Calories  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

CALORIE                    -1 > 0 \* 9999 < 20000    Dec Prec: 0



=====

TERM: INTAKE - IV CARDIAC OUTPUT INJECTATE  
ID: {FE023C43-32D0-4937-A100-0F14CDCACEF2}

=====

Page: 319

Abbreviation..... INT CO INJ  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Cardiac Output Injectate  
Display Name..... Cardiac Output Injectate  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MILLILITER	-1 > 0 * 9999 < 10000	Dec Prec: 0
CUBIC CENTIMETER	-1 > 0 * 9999 < 10000	Dec Prec: 0

=====

TERM: INTAKE - IV COLLOID  
ID: {7CD55D76-7D00-414E-BF24-315A5A92A998}

=====

Page: 320

Abbreviation..... INT IV COL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Intake of IV Colloids  
Display Name..... Intake IV Colloids  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

PRODUCT: 4 entries

- 
1. ALBUMIN (ALBUMIN) {EF4CC3AE-32B2-4EBC-BBDE-7DE7CCEAF081}
  2. 5 PERCENT ALBUMIN (5 PERC ALB)  
   {A2D0F9A0-A9EC-4540-9663-B8B19A189198}
  3. 25 PERCENT ALBUMIN (25 PERC AL)  
   {CBC7E32F-E457-4A7F-8F27-030EBC6ECBD5}
  4. 6 PERCENT HETASTARCH (6 PERC HET)  
   {8A82F108-D2FE-4C21-93BB-344502235AE8}

UNITS:

-----

MILLILITER	-1 > 0 * 9999 < 20000	Dec Prec: 0
CUBIC CENTIMETER	-1 > 0 * 9999 < 10000	Dec Prec: 0

=====

TERM: INTAKE - IV CRYSTALLOIDS  
ID: {DD597777-861C-44FD-82E6-1FFF925D0095}

=====

Page: 321

Abbreviation..... INT IV CRY  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Intake of IV Crystalloids  
Display Name..... Intake - IV Crystalloids  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MILLILITER	-1 > 0 * 9999 < 10000	Dec Prec: 0
CUBIC CENTIMETER	-1 > 0 * 9999 < 10000	Dec Prec: 0





=====

TERM: INTAKE - IV NUTRITION  
ID: {DDD9CA13-29FA-4BDB-9558-B0DFEF76FA15}

=====

Page: 324

Abbreviation..... IN IV NUTR  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... IV Nutrition  
Display Name..... IV Nutrition  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

PRODUCT: 3 entries

- 
1. CRYSTALLOID (CRYSTALLOI) {F8B0C9B5-F386-4C88-896A-092C4A95A031}
  2. COLLOID (COLLOID) {5B912A05-9ECA-43C1-A117-BE0843313CDD}
  3. LIPIDS (LIPIDS) {80E456A6-CA3D-4C2D-8DCD-06C63AE1F320}

UNITS:

-----

MILLILITER -1 > 0 \* 9999 < 10000 Dec Prec: 0  
CUBIC CENTIMETER -1 > 0 \* 9999 < 10000 Dec Prec: 0

=====

TERM: INTAKE - LINE FLUSH  
ID: {0DAF17DD-4039-406C-9F53-5AB46E172648}

=====

Page: 325

Abbreviation..... IN LINE FL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Intake from line flush  
Display Name..... Intake - Line Flush  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

PRODUCT: 2 entries

- 
1. SALINE (SALINE) {19D84001-70DB-47AA-9149-A37B7E0BAEF1}
  2. LINE FLUSH HEPARIN (LN FL HEPA)  
{E8B68E5A-6B9A-401F-9D1A-E3C907C7B81C}

UNITS:

-----

MILLILITER            -1 > 0 \* 9999 < 10000    Dec Prec: 0  
CUBIC CENTIMETER    -1 > 0 \* 9999 < 10000    Dec Prec: 0





=====

TERM: INTAKE - PERCENT MEAL CONSUMED  
ID: {A3C27D46-477C-476E-8D56-A652B8BC3933}

=====

Page: 327

Abbreviation..... INT % MEAL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Percentage of meal consumed by patient  
Display Name..... Percent Meal Consumed  
Help Text..... Percentage of meal consumed by patient  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
PERCENTAGE                    -1 > 0 \* 100 < 101      Dec Prec: 0

=====

TERM: INTAKE - PO CALORIES

Page: 328

ID: {DDCB1A84-6EA1-4ECC-AF57-EEA6754E8540}

=====

Abbreviation..... INT PO CAL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Pharyngeal Oral caloric intake  
Display Name..... Intake - PO Calories  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

CALORIE                    -1 > 0 \* 9999 < 10000    Dec Prec: 0

Abbreviation..... INT PO FLU
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Oral intake of fluids
Display Name..... Intake PO Fluids
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

PRODUCT: 10 entries

- 1. NUTRITIONAL SUPPLEMENT (NUTR SUPPL) {17F8F1DA-4FC6-48E5-9E6E-1F67CD06D200}
2. WATER (WATER) {ACCC998E-0F64-510D-B843-3D30B2AD67CF}
3. JUICE (JUICE) {91FFF519-1677-BB1D-D959-65B1BFA10285}
4. FULL LIQUID (FL LIQUID) {C8405270-CCA3-4BFF-B9D0-398590FBB1D9}
6. OTHER PO FLUID (PO OTHER) {1EA44424-4DE2-42B8-A50F-60DC56204502}
70. MILK (Milk) {6C4DF9C1-D7BA-429A-BDAF-00FBF35F6B09}
80. COFFEE (Coffee) {BEFE282A-F07F-4B3A-A1FC-5FBC0DE75AE1}
90. SODA: DIET (Soda: Diet) {0EFA0EC7-EB52-4F84-B863-5CD8BC5509DA}
100. SODA: NON-DIET (Soda: Non) {B15D88DC-07B6-4980-98B2-A59B5662797A}
110. CLEAR LIQUIDS (Clear Liq) {BBE3A652-ACDC-4797-B140-D4DB4F474B13}

UNITS:

MILLILITER 0 > 1 \* 10000 < 10001 Dec Prec: 0
OUNCE 0 > 1 \* 338 < 339 Dec Prec: 0
CUBIC CENTIMETER 0 > 1 \* 10000 < 10001 Dec Prec: 0

Abbreviation..... INT TPN CA  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Total parenteral nutrition calories  
Display Name..... Intake - TPN Calories  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

CALORIE                    -1 > 0 \* 2000 < 2001    Dec Prec: 0

=====

TERM: INTAKE - TRANSFUSION  
ID: {DCD2263B-0390-4754-950A-6A410C0C59EB}

=====

Page: 331

Abbreviation..... INT TRANSF  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Transfusion Type  
Display Name..... Transfusion  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

CUBIC CENTIMETER        -1 > 0 \* 9999 < 10000    Dec Prec: 0

Abbreviation..... ICP
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Intracranial Pressure - Do Not Check
Total or Subtotals
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

- 1. VENTRICULAR DRAIN (VENT DRAIN)
{5FF035D1-9C97-415D-BBEB-943506395A0C}

UNITS:

MILLIMETERS OF MERCURY 0 > 8 \* 20 < 30 Dec Prec: 0



Abbreviation..... DR UN HR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Only to be used with products with
dosage measured in "Units"
Display Name..... IV DOSE RATE - UNITS/HOUR - Do Not
Check Total or Subtotals
Help Text..... Only to be used with products with
dosage measured in "Units"
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

-----
UNITS PER HOUR 0 > 1 \* 1200 < 2001 Dec Prec: 0



=====

TERM: IV DOSE RATE - UNITS/kg/HOUR - DO NOT CHECK TOTAL OR SUBTOTALS Page: 335  
ID: {D1749D90-3EEE-41BE-BBF5-B3B00C18E8D2}

=====

Abbreviation..... IVDR UKH  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Only to be used with products with  
dosage measured in "Units"  
Display Name..... IV DOSE RATE - UNITS/kg/HOUR - Do Not  
Check Total or Subtotals  
Help Text..... Only to be used with products with  
dosage measured in "Units"  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Void.....

UNITS:

-----  
UNITS/kg/HOUR            0 > 10 \* 51 < 101            Dec Prec: 0



Abbreviation..... IV DR mc
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... IV DOSE RATE - mcg/TIME - Do Not
Check Total or Subtotals
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

-----
MICROGRAMS PER HOUR -1.00 > 0.00 \* 9999.00 < 10000.00
Dec Prec: 2
MICROGRAMS PER MINUTE -1.00 > 0.00 \* 9999.00 < 10000.00
Dec Prec: 2







Abbreviation..... IV FLOW R
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... IV Flow Rate in ml / min or hour
Display Name..... IV Flow Rate - Do Not Check Total or
Subtotals
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

-----
MILLILITERS PER HOUR -1.00 > 0.00 \* 9999.00 < 10000.00
Dec Prec: 2
MILLILITERS PER MINUTE -1.00 > 0.00 \* 9999.00 < 10000.00
Dec Prec: 2

=====

TERM: IV FLOW RATE - ml/kg/TIME - DO NOT CHECK TOTAL OR SUBTOTALS    Page: 342  
ID: {5BDF4F11-9713-4B20-8DF1-49BB33DEAB67}

=====

Abbreviation..... IV FL kg  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... IV FLOW RATE - ml/kg/Time - Do Not  
                    Check Total or Subtotals  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Void.....

UNITS:

-----

MILLILITERS PER KILOGRAM PER HOUR  
                    0.00 > 1.00 \* 40.00 < 100.00  
                                    Dec Prec: 2

MILLILITERS PER KILOGRAM PER MINUTE  
                    0.00 > 0.33 \* 1.33 < 1.70 Dec Prec: 2



Abbreviation..... IV LI CARE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Care of IV Lines  
Display Name..... IV Line Care  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 20 entries

-----

1. PATENT CLEAN (PAT CLEAN) {68688246-AE58-2EFB-BF9D-5E00D8CFDA8E}
2. DRESSING CHANGED (DRESS CHAN)  
{E8CF4B7A-D1AB-49A2-8E7C-3375C3BC436F}
3. STARTED (STARTED) {D6E105C4-5EAD-424E-8E4C-195EB7783A3D}
4. TUBING CHANGED/LABELED (TUBE CHANG)  
{19BA337C-6E28-411E-9702-48B262C593FA}
5. HEPLOCK (HEPLOCK) {C827FEDE-59D6-4C45-BDDA-4884FA0882F3}
6. BAG CHANGED (BAG CHANGE) {818990DE-4F95-4F8E-92D8-B4D3B975CB36}
7. LEVELED (LEVELED) {D72523B2-E0CD-40E1-AC14-2FA43636A456}
8. PRESSURE CALIBRATION (PRESS CALI)  
{38C12E69-2CD2-4F61-98D9-06AB857DA51A}
9. ARTERIAL LINE CARE (ART LINE C)  
{787D232D-DFB8-44E2-8A0C-9D7D6B80764E}
10. CENTRAL VENOUS CATHETER (CENT VEN C)  
{EE91E841-C09B-4432-A121-8EBAB92E9CD2}
11. PICC (PICC) {3CFE33D0-5A35-88CC-152C-D7501D1858D9}
12. FLUSHED (nt) {C7C54C1A-6E75-4BCF-4C57-0F8622A5BD16}
13. CHANGED (CHANGED) {74B14023-D74D-41F6-AA48-B0E5137623CD}
14. DISCONTINUED (DISCONTINU) {15E289FB-09DF-4731-9CF3-48A4E0031C31}
15. SALINE LOCK (SALIN LOCK) {A94D7639-2A45-4319-A6C7-5A892AA6B1F4}
16. ALCOHOL PORT PROTECTOR (Al Port Pr)  
{785F47CF-A279-46F9-9073-B6A6A70E8888}
17. NEUTRAL PRESSURE CONNECTOR (N Pres Con)  
{263F1A2C-4384-4B6A-B930-B9AF85B7F965}
18. POSITIVE PRESSURE CONNECTOR (P Pres Con)  
{1A662A2F-4666-4CFD-9E5C-7C5115BB7C2D}
19. NEEDLELESS HUB DISINFECTED AND ACCESSED PER POLICY (Needleless)  
{5D9D8D90-9B39-4DF3-8095-755235687F45}
20. MEDICAL NECESSITY CONFIRMED (med nec)  
{F380D431-196B-46C9-8A1C-C57E7ECBFC8A}

UNITS:

-----

Abbreviation..... IV PORT #
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... IV Port descriptor or number
Display Name..... IV Port / Number
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 15 entries

- 1. PROXIMAL PORT (PROX PORT) {3244EF4F-F18A-4250-96F5-ABCEA7BFA578}
2. MEDIAL PORT (MEDI PORT) {D1734D73-2E1A-4538-8884-05D8A78CF556}
3. DISTAL PORT (DIST PORT) {E493DF34-AD0E-4951-9715-5AD4EFA14036}
4. 1 (1) {8F5D10A8-FE17-47BC-DACA-6972D12CC00A}
5. 2 (2) {38E22B79-A2B6-3401-B195-C955F24E1278}
6. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}
7. 4 (4) {D2280C7E-DE50-5A8D-C3FD-5862B69C0E1A}
8. 5 (5) {188EDCE7-0726-4234-D8D2-92C8052CD221}
9. 6 (6) {640CF2F9-3BBF-AB13-185B-CD02A609FB5A}
10. 7 (7) {7A010324-D2C2-5D4B-848A-7B8878666D20}
11. 8 (8) {518E7CFC-102E-351E-C7AF-FB575E3C8F4D}
12. 9 (9) {58FF2FD4-569F-F454-9F5B-FE6F9E34D8CC}
13. 10 (10) {0C827160-18DB-816F-7FAA-A3C175C6D33B}
14. 11 (11) {BD51A7B3-FACF-890C-472E-E3B28519CC16}
15. 12 (12) {ACC5CA06-3ABC-476D-EFB0-D82FEEE129C7}

LOCATION: 20 entries

- 1. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
2. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
3. L ANTECUBITAL (L ANTECUBI) {C21838CC-BAB5-1D7F-E2A4-A7C663FB28C3}
4. R ANTECUBITAL (R ANTECUBI) {7FB00A65-B4B2-06A6-BD89-98AC112D8D4E}
5. L UPPER ARM (L UPP ARM) {5F915696-1B08-32E9-6B95-49246117D660}
6. R UPPER ARM (R UPP ARM) {3EF20E8E-F68C-E4F2-8F8C-FAD3B44DDAA4}
7. L WRIST (L WRIST) {7BCAFACE-FF3B-46DA-99F9-6437A90B7126}
8. R WRIST (R WRIST) {FDD9FE1D-F8A2-09AD-82D7-AA478693DD01}
9. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
10. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
11. L LEG (LL) {B088C809-E528-B785-21C2-14E45E1FD8AF}
12. R LEG (RL) {C5E40A7B-6FC4-267D-8CE8-3571D21D81FD}
13. L FOOT (L FOOT) {2EC3FBEC-4EA2-9E7A-17D4-FDC2F9D28046}
14. R FOOT (R FOOT) {7D9F5FF6-232C-2488-4504-A721FF06928F}
15. L INTERNAL JUGULAR (LIJ) {AFBB4FE5-8A11-2047-7900-DADF2FC4430D}
16. R INTERNAL JUGULAR (R INT JUGU) {7ABF6916-D379-198D-6454-F0518A74844A}
17. L EXTERNAL JUGULAR (L EXT JUGU) {6F7B331B-4F83-0B79-1755-F6F673D647D6}
18. R EXTERNAL JUGULAR (R EXT JUGU) {4D78F704-9E3B-907E-A792-4DC4F4CFF170}
19. L SUBCLAVIAN (LSCL) {71B93594-15F9-AA4A-AEF2-E222EC28C00D}
20. R SUBCLAVIAN (R SUBCLAVI) {53DEAE83-93AE-123B-14E5-4DE38CD4F752}

UNITS:

-----

Abbreviation..... IV ST COND
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... IV Site Condition
Help Text..... Condition of IV Site
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 7 entries

- 2. RED (RED) {515D7CFB-3527-30E7-EC4B-B607C272DFA0}
3. SWOLLEN (SWOLLEN) {DCE86B10-C393-8DA9-6518-F2943F0E1A4E}
4. DISCOMFORT (DISCOMFORT) {77D058C9-1EB4-96C0-C9E3-666C840B884C}
50. PATENT (Patent) {FOC38BCB-1D11-436C-8B8C-4297C27DDE04}
60. CLEAN (Clean) {1F187B99-9963-4AD5-98F6-D8D24D2DDD30}
70. OCCLUDED (Occluded) {25509BE3-F946-4AD7-80E8-5CE8A2E61AD6}
80. OOZING (Oozing) {C23549D7-3238-4523-BA31-AB71609F01C2}

LOCATION: 20 entries

- 1. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
2. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
3. L ANTECUBITAL (L ANTECUBI) {C21838CC-BAB5-1D7F-E2A4-A7C663FB28C3}
4. R ANTECUBITAL (R ANTECUBI) {7FB00A65-B4B2-06A6-BD89-98AC112D8D4E}
5. L UPPER ARM (L UPP ARM) {5F915696-1B08-32E9-6B95-49246117D660}
6. R UPPER ARM (R UPP ARM) {3EF20E8E-F68C-E4F2-8F8C-FAD3B44DDAA4}
7. L WRIST (L WRIST) {7BCAFACE-FF3B-46DA-99F9-6437A90B7126}
8. R WRIST (R WRIST) {FDD9FE1D-F8A2-09AD-82D7-AA478693DD01}
9. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
10. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
11. L LEG (LL) {B088C809-E528-B785-21C2-14E45E1FD8AF}
12. R LEG (RL) {C5E40A7B-6FC4-267D-8CE8-3571D21D81FD}
13. L FOOT (L FOOT) {2EC3FBEC-4EA2-9E7A-17D4-FDC2F9D28046}
14. R FOOT (R FOOT) {7D9F5FF6-232C-2488-4504-A721FF06928F}
15. L INTERNAL JUGULAR (LIJ) {AFBB4FE5-8A11-2047-7900-DADF2FC4430D}
16. R INTERNAL JUGULAR (R INT JUGU) {7ABF6916-D379-198D-6454-F0518A74844A}
17. L EXTERNAL JUGULAR (L EXT JUGU) {6F7B331B-4F83-0B79-1755-F6F673D647D6}
18. R EXTERNAL JUGULAR (R EXT JUGU) {4D78F704-9E3B-907E-A792-4DC4F4CFF170}
19. L SUBCLAVIAN (LSCL) {71B93594-15F9-AA4A-AEF2-E222EC28C00D}
20. R SUBCLAVIAN (R SUBCLAVI) {53DEAE83-93AE-123B-14E5-4DE38CD4F752}

UNITS:

=====

TERM: LABS PER ORDERS

Page: 346

ID: {6D66F3BC-9A35-4CE2-8677-00938428CD19}

=====

Abbreviation..... Labs p/Ord  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description.....  
Display Name..... Labs per Orders  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----



Abbreviation..... LVEF - CAT
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Categorical Left Ventricular Ejection
Fraction
Display Name..... LVEF - Categorical
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 8 entries

- 1. HYPERDYNAMIC (HYPERDYNAM) {818A8A99-04D9-4989-A9BA-ECFACED74605}
2. NORMAL (n) {553C9A97-A3D5-C52F-0E4E-C7D168C99A2B}
3. BORDERLINE REDUCED (BORD REDUC)
{63703AB1-6249-403B-BDE8-FECA9FA36CF5}
4. MILD-MODERATELY REDUCED (MI-MOD-RED)
{41DA2518-0E49-4E92-A733-3ACF3D2E8199}
5. MODERATELY REDUCED (MOD REDUCE)
{91F73300-7EE0-428C-A838-5B3BE062D8BC}
6. MODERATE-SEVERELY REDUCED (MO-SEV RED)
{6D031C70-2F95-4111-A9B8-BC22CEC1B17D}
7. SEVERELY REDUCED (SEV REDUCE)
{A895BACE-85AC-45B2-A190-4A6816DEDC1D}
8. LESS THAN 30 (<30) {5AA98A0D-9B05-425D-97CA-B8B12425AACA}

METHOD: 7 entries

- 1. ECHOCARDIOGRAPHY (ECG) {CA44697E-3BBA-4767-A30D-6548A72F9639}
2. NUCLEAR PERFUSION IMAGING (NPI)
{D0F7DB21-0D66-4161-8214-DACCAF41ADAC}
3. MULTIGATED NUCLEAR VENTRICULOGRAPHY (MUGA)
{4D73168B-4DF4-4EA2-83B9-ACE853D5DD42}
4. CONTRAST LEFT VENTRICULOGRAPHY (CON L VENT)
{35BEA63F-9696-4C03-AE1D-D76E4DDBB017}
5. MAGNETIC RESONANCE IMAGING (MRI)
{50EBFFF8-DF0D-41F2-B11B-33A26AA40DBA}
6. COMPUTED TOMOGRAPHY (COM TOMOGR)
{B4676AE0-BFAF-4649-AB08-3710E39F9B8D}
7. OTHER (OTHER) {AED275BD-8A42-46BE-BCDF-472EF4C31D65}

UNITS:

=====

TERM: LEFT VENTRICULAR STROKE WORK  
ID: {E9D668F0-920C-7E29-A6A1-618FC8ECD19B}

=====

Page: 349

Abbreviation..... LVSW  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Left Ventricular Stroke Work  
Display Name..... Left Ventricular Stroke Work  
Help Text..... Left Ventricular Stroke Work  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----

1. CARDIAC OUTPUT MONITOR (COM)  
{CCBE5744-7CB9-9DFA-938F-FDDBC73B45B3}

UNITS:

-----

STROKE WORK PER BEAT 29.9 > 30.0 \* 150.0 < 150.1  
Dec Prec: 1



ID: {49E9E897-E19A-8FE9-416C-F98CA66EFBC9}

=====

Abbreviation..... LEGMOV  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Leg Movements  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |                |  |
|----------------|--|
| 1. ABSENT (a)  | {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0} |
| 2. PRESENT (p) | {384DC20F-1625-CE1E-1C85-A3B3924BB2A0} |

UNITS:  
-----

Abbreviation..... LOC
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Level of Consciousness
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

- 1. ALERT (a) {BC7E2630-179A-5C80-E6C7-598D42F0BFD6}
2. DROWSY (d) {6B350814-A88A-0F05-C382-9F27830FBC09}
3. RESTLESS (r) {8FCCAD26-DBA0-7067-DA72-30D342BFA8E8}
4. AGITATED (a) {BAD2CCCE-9083-6F61-1196-C0C9A8DDBFFC}
5. LETHARGIC (l) {4ECA2E44-093A-9354-848E-CDC72B165072}

UNITS:

Abbreviation..... LLES  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Lower Limb Extremities Strength  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 6 entries

- 
- |                      |  |
|----------------------|--|
| 1. ABSENT (a)        | {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0} |
| 2. WEAK (w)          | {BE2D9305-E2DA-F492-66CF-164EE58B8179} |
| 3. STRONG (s)        | {A12154E1-544E-F372-3839-0D44E5BC476B} |
| 4. REFLEX (r)        | {F21585F9-5206-A6B3-2BF6-F416E5842366} |
| 5. EQUAL (Equal)     | {AC079584-158D-EBF2-E13A-C2E919375E4D} |
| 6. UNEQUAL (Unequal) | {73319200-164F-B39F-649C-B7FB1AFBDBB8} |

LOCATION: 3 entries

- 
- |                              |  |
|------------------------------|--|
| 1. L LOWER LIMB (LLLI)       | {4838CB05-A1DD-F4FB-587C-2D0C17A2484E} |
| 2. R LOWER LIMB (R LOW LIMB) | {A8E26A05-311A-F196-40B5-B5202D9F8058} |
| 3. BILATERAL LEG (BILAT LEG) | {0713761D-E0FC-DA03-02A1-F20EFEB99F45} |

UNITS:

-----

=====

TERM: LVSWI Page: 353  
ID: {7945FB00-2FA8-6D4C-BC85-0110EA75AC5C}

=====

Abbreviation..... LVSWI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Left Ventricular Stroke Work Index  
Display Name..... LVSWI  
Help Text..... Left Ventricular Stroke Work Index  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----

1. CARDIAC OUTPUT MONITOR (COM)  
{CCBE5744-7CB9-9DFA-938F-FDDBC73B45B3}

UNITS:

-----

STROKE WORK PER BEAT PER M2 19 > 20 \* 150 < 151 Dec Prec: 0  
GM-MM2/BEAT 19 > 20 \* 150 < 151 Dec Prec: 0

=====

TERM: MALLAMPATI SCORE  
ID: {33452BD5-3E25-4B92-8849-F31EE04879AD}

=====

Page: 354

Abbreviation..... Mallampati  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... Mallampati Score  
Help Text.....  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
- |                          |  |
|--------------------------|--|
| 1. CLASS I (Class I)     | {4A1B7DD6-3846-4455-A06D-0E3C0FB3752A} |
| 2. CLASS II (Class II)   | {2F5F3798-3EFD-4DEA-9FAB-883F9E326EA7} |
| 3. CLASS III (Class III) | {C0FFDE12-241C-4087-9404-DB128F6B793A} |
| 4. CLASS IV (Class IV)   | {176030D0-2576-4D2D-A5AB-E5DA38DADEB5} |

UNITS:

-----











Abbreviation..... Dress Self  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... MDS: Dressing-Self Performance  
Help Text..... How resident puts on, fastens and  
takes off all items of clothing,  
including donning/removing a  
prosthesis or TED hose. Dressing  
includes putting on and changing  
pajamas and housedresses.  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Void.....

DATA VALUE: 6 entries

- 
- 10. 0 - INDEPENDENT - NO HELP OR STAFF OVERSIGHT (0 Ind)  
{4A125BF9-A2C1-4503-9182-BF3E2A708547}
  - 20. 1 - SUPERVISION - OVERSIGHT (1 Superv)  
{7400DA9E-E89A-44F4-8516-CB6D5931D97D}
  - 30. 2 - LIMITED ASSISTANCE (2 Limited)  
{35934E53-08EB-46E0-8136-96FD6174ACF1}
  - 40. 3 - EXTENSIVE ASSISTANCE (3 Extensiv)  
{4B1FCD6E-7211-4151-ACD8-DF2F614A4337}
  - 50. 4 - TOTAL DEPENDENCE - FULL STAFF PERFORMANCE (4 Full Stf)  
{E8B7D02B-E80A-4851-A0E2-EF27CF75A5BB}
  - 60. 8 - ACTIVITY DID NOT OCCUR (8 d/n occ)  
{8B586CC9-2360-4060-B660-021C25E28F0F}

UNITS:  
-----

Abbreviation..... Dress Supp  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... MDS: Dressing-Support  
Help Text..... How resident puts on, fastens and  
takes off all items of clothing,  
including donning/removing a  
prosthesis or TED hose. Dressing  
includes putting on and changing  
pajamas and housedresses.  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Void.....

DATA VALUE: 5 entries

- 
- 10. 0 - NO SETUP OR PHYSICAL HELP FROM STAFF (0 No Setup)  
{730049CC-4CD5-449E-A0E8-389AB61A02BD}
  - 20. 1 - SETUP HELP ONLY (1 Setup )  
{5494773C-99FE-4A96-B12E-D734188517CE}
  - 30. 2 - ONE PERSON PHYSICAL ASSIST (2 one asst)  
{C675F08C-1BEE-4CDC-AE01-5CE429AB2105}
  - 40. 3 - TWO + PERSONS PHYSICAL ASSIST (3 + asst)  
{AE401096-DF23-4A6E-8842-E1F58829902A}
  - 50. 8 - ACTIVITY DID NOT OCCUR (8 d/n occ)  
{B016EA6E-B60A-4075-9B27-17B7C86D8E9F}

UNITS:  
-----

Abbreviation..... Eat Self
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description.....
Display Name..... MDS: Eating-Self Performance
Help Text..... How resident eats and drinks,
regardless of skill. Do not include
eating/drinking during medication
pass. Includes intakes of nourishment
by other means (eg. Tube feeding,
total parenteral nutrition, IV fluids
administered for nutrition or
hydration)
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 6 entries

- 10. 0 - INDEPENDENT - NO HELP OR STAFF OVERSIGHT (0 Ind)
{4A125BF9-A2C1-4503-9182-BF3E2A708547}
20. 1 - SUPERVISION - OVERSIGHT (1 Superv)
{7400DA9E-E89A-44F4-8516-CB6D5931D97D}
30. 2 - LIMITED ASSISTANCE (2 Limited)
{35934E53-08EB-46E0-8136-96FD6174ACF1}
40. 3 - EXTENSIVE ASSISTANCE (3 Extensiv)
{4B1FCD6E-7211-4151-ACD8-DF2F614A4337}
50. 4 - TOTAL DEPENDENCE - FULL STAFF PERFORMANCE (4 Full Stf)
{E8B7D02B-E80A-4851-A0E2-EF27CF75A5BB}
60. 8 - ACTIVITY DID NOT OCCUR (8 d/n occ)
{8B586CC9-2360-4060-B660-021C25E28F0F}

UNITS:

Abbreviation..... Eat Supp
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description.....
Display Name..... MDS: Eating-Support
Help Text..... How resident eats and drinks,
regardless of skill. Do not include
eating/drinking during medication
pass. Includes intakes of nourishment
by other means (eg. Tube feeding,
total parenteral nutrition, IV fluids
administered for nutrition or
hydration)
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

- 10. 0 - NO SETUP OR PHYSICAL HELP FROM STAFF (0 No Setup)
20. 1 - SETUP HELP ONLY (1 Setup )
30. 2 - ONE PERSON PHYSICAL ASSIST (2 one asst)
40. 3 - TWO + PERSONS PHYSICAL ASSIST (3 + asst)
50. 8 - ACTIVITY DID NOT OCCUR (8 d/n occ)

UNITS:
=====

Abbreviation..... Locom Self
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description.....
Display Name..... MDS: MDS: Locomotion Off Unit-Self
Performance
Help Text..... How resident moves to and returns
from off-unit locations (eg. areas
set aside for dining, activities,
treatments). If facility only has 1
floor, how resident moves to/from
distant areas on floor. If in
wheelchair, self-sufficiency once in
chair
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 6 entries

- 10. 0 - INDEPENDENT - NO HELP OR STAFF OVERSIGHT (0 Ind)
{4A125BF9-A2C1-4503-9182-BF3E2A708547}
20. 1 - SUPERVISION - OVERSIGHT (1 Superv)
{7400DA9E-E89A-44F4-8516-CB6D5931D97D}
30. 2 - LIMITED ASSISTANCE (2 Limited)
{35934E53-08EB-46E0-8136-96FD6174ACF1}
40. 3 - EXTENSIVE ASSISTANCE (3 Extensiv)
{4B1FCD6E-7211-4151-ACD8-DF2F614A4337}
50. 4 - TOTAL DEPENDENCE - FULL STAFF PERFORMANCE (4 Full Stf)
{E8B7D02B-E80A-4851-A0E2-EF27CF75A5BB}
60. 8 - ACTIVITY DID NOT OCCUR (8 d/n occ)
{8B586CC9-2360-4060-B660-021C25E28F0F}

UNITS:

Abbreviation..... Locom Supp  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... MDS: Locomotion Off Unit-Support  
Help Text..... How resident moves to and returns  
from off-unit locations (eg. areas  
set aside for dining, activities,  
treatments). If facility only has 1  
floor, how resident moves to/from  
distant areas on floor. If in  
wheelchair, self-sufficiency once in  
chair  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
- 10. 0 - NO SETUP OR PHYSICAL HELP FROM STAFF (0 No Setup)  
{730049CC-4CD5-449E-A0E8-389AB61A02BD}
  - 20. 1 - SETUP HELP ONLY (1 Setup )  
{5494773C-99FE-4A96-B12E-D734188517CE}
  - 30. 2 - ONE PERSON PHYSICAL ASSIST (2 one asst)  
{C675F08C-1BEE-4CDC-AE01-5CE429AB2105}
  - 40. 3 - TWO + PERSONS PHYSICAL ASSIST (3 + asst)  
{AE401096-DF23-4A6E-8842-E1F58829902A}
  - 50. 8 - ACTIVITY DID NOT OCCUR (8 d/n occ)  
{B016EA6E-B60A-4075-9B27-17B7C86D8E9F}

UNITS:  
-----

Abbreviation..... Locom Unt  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... MDS: Locomotion on Unit-Self  
Performance  
Help Text..... How resident moves between locations  
in his/her room and adjacent corridor  
on same floor. If in wheelchair,  
self-sufficiency once in chair  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 6 entries

- 
- 10. 0 - INDEPENDENT - NO HELP OR STAFF OVERSIGHT (0 Ind)  
{4A125BF9-A2C1-4503-9182-BF3E2A708547}
  - 20. 1 - SUPERVISION - OVERSIGHT (1 Superv)  
{7400DA9E-E89A-44F4-8516-CB6D5931D97D}
  - 30. 2 - LIMITED ASSISTANCE (2 Limited)  
{35934E53-08EB-46E0-8136-96FD6174ACF1}
  - 40. 3 - EXTENSIVE ASSISTANCE (3 Extensiv)  
{4B1FCD6E-7211-4151-ACD8-DF2F614A4337}
  - 50. 4 - TOTAL DEPENDENCE - FULL STAFF PERFORMANCE (4 Full Stf)  
{E8B7D02B-E80A-4851-A0E2-EF27CF75A5BB}
  - 60. 8 - ACTIVITY DID NOT OCCUR (8 d/n occ)  
{8B586CC9-2360-4060-B660-021C25E28F0F}

UNITS:  
-----





Abbreviation..... Hygi Self
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description.....
Display Name..... MDS: Personal Hygiene-Self Performance
Help Text..... How resident maintains personal
hygiene, including combing hair,
brushing teeth, shaving, applying
makeup, washing/drying face and hands
(excludes baths and showers)
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 6 entries

- 10. 0 - INDEPENDENT - NO HELP OR STAFF OVERSIGHT (0 Ind)
20. 1 - SUPERVISION - OVERSIGHT (1 Superv)
30. 2 - LIMITED ASSISTANCE (2 Limited)
40. 3 - EXTENSIVE ASSISTANCE (3 Extensiv)
50. 4 - TOTAL DEPENDENCE - FULL STAFF PERFORMANCE (4 Full Stf)
60. 8 - ACTIVITY DID NOT OCCUR (8 d/n occ)

UNITS:

Abbreviation..... Hygi Supp
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description.....
Display Name..... MDS: Personal Hygiene-Support
Help Text..... How resident maintains personal
hygiene, including combing hair,
brushing teeth, shaving, applying
makeup, washing/drying face and hands
(excludes baths and showers)
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

- 10. 0 - NO SETUP OR PHYSICAL HELP FROM STAFF (0 No Setup)
20. 1 - SETUP HELP ONLY (1 Setup )
30. 2 - ONE PERSON PHYSICAL ASSIST (2 one asst)
40. 3 - TWO + PERSONS PHYSICAL ASSIST (3 + asst)
50. 8 - ACTIVITY DID NOT OCCUR (8 d/n occ)

UNITS:
=====

Abbreviation..... Toil Self
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description.....
Display Name..... MDS: Toilet Use-Self Performance
Help Text..... How resident uses the toilet room,
commode, bedpan or urinal; transfers
on/off toilet; cleanses after
elimination; changes pad; manages
ostomy or catheter; adjusts clothes.
Do not include emptying bedpan,
urinal, commode, catheter or ostomy
bag.
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 6 entries

- 10. 0 - INDEPENDENT - NO HELP OR STAFF OVERSIGHT (0 Ind)
{4A125BF9-A2C1-4503-9182-BF3E2A708547}
20. 1 - SUPERVISION - OVERSIGHT (1 Superv)
{7400DA9E-E89A-44F4-8516-CB6D5931D97D}
30. 2 - LIMITED ASSISTANCE (2 Limited)
{35934E53-08EB-46E0-8136-96FD6174ACF1}
40. 3 - EXTENSIVE ASSISTANCE (3 Extensiv)
{4B1FCD6E-7211-4151-ACD8-DF2F614A4337}
50. 4 - TOTAL DEPENDENCE - FULL STAFF PERFORMANCE (4 Full Stf)
{E8B7D02B-E80A-4851-A0E2-EF27CF75A5BB}
60. 8 - ACTIVITY DID NOT OCCUR (8 d/n occ)
{8B586CC9-2360-4060-B660-021C25E28F0F}

UNITS:

Abbreviation..... Toil Supp
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description.....
Display Name..... MDS: Toilet Use-Support
Help Text..... How resident uses the toilet room,
commode, bedpan or urinal; transfers
on/off toilet; cleanses after
elimination; changes pad; manages
ostomy or catheter; adjusts clothes.
Do not include emptying bedpan,
urinal, commode, catheter or ostomy
bag.
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

- 10. 0 - NO SETUP OR PHYSICAL HELP FROM STAFF (0 No Setup)
20. 1 - SETUP HELP ONLY (1 Setup )
30. 2 - ONE PERSON PHYSICAL ASSIST (2 one asst)
40. 3 - TWO + PERSONS PHYSICAL ASSIST (3 + asst)
50. 8 - ACTIVITY DID NOT OCCUR (8 d/n occ)

UNITS:





Abbreviation..... Wlk Cor Sf  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... MDS: Walk in Corridor-Self Performance  
Help Text..... How resident walks in corridor on unit  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 6 entries

- 
- 10. 0 - INDEPENDENT - NO HELP OR STAFF OVERSIGHT (0 Ind)  
{4A125BF9-A2C1-4503-9182-BF3E2A708547}
  - 20. 1 - SUPERVISION - OVERSIGHT (1 Superv)  
{7400DA9E-E89A-44F4-8516-CB6D5931D97D}
  - 30. 2 - LIMITED ASSISTANCE (2 Limited)  
{35934E53-08EB-46E0-8136-96FD6174ACF1}
  - 40. 3 - EXTENSIVE ASSISTANCE (3 Extensiv)  
{4B1FCD6E-7211-4151-ACD8-DF2F614A4337}
  - 50. 4 - TOTAL DEPENDENCE - FULL STAFF PERFORMANCE (4 Full Stf)  
{E8B7D02B-E80A-4851-A0E2-EF27CF75A5BB}
  - 60. 8 - ACTIVITY DID NOT OCCUR (8 d/n occ)  
{8B586CC9-2360-4060-B660-021C25E28F0F}

UNITS:  
-----



Abbreviation..... Wlk Cor Sp  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... MDS: Walk in Corridor-Support  
Help Text..... How resident walks in corridor on unit  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
- 10. 0 - NO SETUP OR PHYSICAL HELP FROM STAFF (0 No Setup)  
{730049CC-4CD5-449E-A0E8-389AB61A02BD}
  - 20. 1 - SETUP HELP ONLY (1 Setup )  
{5494773C-99FE-4A96-B12E-D734188517CE}
  - 30. 2 - ONE PERSON PHYSICAL ASSIST (2 one asst)  
{C675F08C-1BEE-4CDC-AE01-5CE429AB2105}
  - 40. 3 - TWO + PERSONS PHYSICAL ASSIST (3 + asst)  
{AE401096-DF23-4A6E-8842-E1F58829902A}
  - 50. 8 - ACTIVITY DID NOT OCCUR (8 d/n occ)  
{B016EA6E-B60A-4075-9B27-17B7C86D8E9F}

UNITS:  
-----

Abbreviation..... Wlk Rm Sel
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description.....
Display Name..... MDS: Walk in room-Self Performance
Help Text..... How resident walks between locations
in his/her room
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 6 entries

- 10. 0 - INDEPENDENT - NO HELP OR STAFF OVERSIGHT (0 Ind)
{4A125BF9-A2C1-4503-9182-BF3E2A708547}
20. 1 - SUPERVISION - OVERSIGHT (1 Superv)
{7400DA9E-E89A-44F4-8516-CB6D5931D97D}
30. 2 - LIMITED ASSISTANCE (2 Limited)
{35934E53-08EB-46E0-8136-96FD6174ACF1}
40. 3 - EXTENSIVE ASSISTANCE (3 Extensiv)
{4B1FCD6E-7211-4151-ACD8-DF2F614A4337}
50. 4 - TOTAL DEPENDENCE - FULL STAFF PERFORMANCE (4 Full Stf)
{E8B7D02B-E80A-4851-A0E2-EF27CF75A5BB}
60. 8 - ACTIVITY DID NOT OCCUR (8 d/n occ)
{8B586CC9-2360-4060-B660-021C25E28F0F}

UNITS:

Abbreviation..... Wlk Rm Sp  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... MDS: Walk in Room-Support  
Help Text..... How resident walks between locations  
                                in his/her room  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE:                                5 entries

- 
- 10. 0 - NO SETUP OR PHYSICAL HELP FROM STAFF (0 No Setup)  
                                {730049CC-4CD5-449E-A0E8-389AB61A02BD}
  - 20. 1 - SETUP HELP ONLY (1 Setup )  
                                {5494773C-99FE-4A96-B12E-D734188517CE}
  - 30. 2 - ONE PERSON PHYSICAL ASSIST (2 one asst)  
                                {C675F08C-1BEE-4CDC-AE01-5CE429AB2105}
  - 40. 3 - TWO + PERSONS PHYSICAL ASSIST (3 + asst)  
                                {AE401096-DF23-4A6E-8842-E1F58829902A}
  - 50. 8 - ACTIVITY DID NOT OCCUR (8 d/n occ)  
                                {B016EA6E-B60A-4075-9B27-17B7C86D8E9F}

UNITS:  
-----

=====

TERM: MEAN ARTERIAL PRESSURE  
ID: {9D9D9FD8-C352-6434-E796-C583407EBD45}

=====

Page: 377

Abbreviation..... MAP  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Mean Arterial Pressure  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 5 entries

- 
1. CUFF (Cu) {E59A2543-1860-1E1C-BEF1-6E7F6B6FBC5E}
  2. ARTERIAL LINE (AL) {D4B3F4ED-8B33-431E-4164-BD165E823776}
  3. INVASIVE (I) {5DB8746E-5446-F47A-883C-A44B91275A6E}
  4. CARDIAC OUTPUT MONITOR (COM)  
{CCBE5744-7CB9-9DFA-938F-FDDBC73B45B3}
  5. NON-INVASIVE (N-I) {F50ABFED-5FD7-FC95-20FE-1E3038C9F299}

UNITS:

-----

MILLIMETERS OF MERCURY    25 > 50 \* 100 < 300    Dec Prec: 0

=====

TERM: MEAN BLOOD PRESSURE  
ID: {9038EF74-932A-42EE-992F-A20BD8784BB7}

=====

Page: 378

Abbreviation..... BP MEAN  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Mean Blood Pressure  
Display Name..... Mean Blood Pressure  
Help Text..... Mean Blood Pressure  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 5 entries

- 
1. ARTERIAL LINE (AL) {D4B3F4ED-8B33-431E-4164-BD165E823776}
  2. NON-INVASIVE (N-I) {F50ABFED-5FD7-FC95-20FE-1E3038C9F299}
  3. CUFF (Cu) {E59A2543-1860-1E1C-BEF1-6E7F6B6FBC5E}
  4. IABP (IABP) {4352A41F-AAF3-4544-9590-BAFD7977DF49}
  5. INVASIVE (I) {5DB8746E-5446-F47A-883C-A44B91275A6E}

UNITS:

-----

MILLIMETERS OF MERCURY    20 > 50 \* 100 < 200    Dec Prec: 0

=====

TERM: MEAN BLOOD PRESSURE ARTERIAL LINE  
ID: {4493FB14-8C5E-9995-C08F-D0BC9A3B5D0B}

=====

Page: 379

Abbreviation..... MBPAL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Mean Blood Pressure Arterial Line  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
MILLIMETERS OF MERCURY      24 > 25 \* 200 < 201      Dec Prec: 0

=====

TERM: MEAN BLOOD PRESSURE NON-INVASIVE  
ID: {A89EC390-60BB-FE79-6F40-F9E2BAFE5BB5}

=====

Page: 380

Abbreviation..... MBPNI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Mean Blood Pressure Non-Invasive  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
MILLIMETERS OF MERCURY    24 > 25 \* 200 < 201    Dec Prec: 0

=====

TERM: MEAN BP CUFF

Page: 381

ID: {1F7DAE88-DE53-4807-F51F-47AB8A672EA1}

=====

Abbreviation..... MBPC  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Mean BP Cuff  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MILLIMETERS OF MERCURY      24 > 25 \* 200 < 201      Dec Prec: 0



=====

TERM: MED LAST DOSE: ANTIBIOTIC  
ID: {F521A3A6-27C9-48F7-8203-6E5634783198}

=====

Page: 382

Abbreviation..... MLD AnBio  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... String  
Description.....  
Display Name..... Med Last Dose: Antibiotic  
Help Text..... Date and time of Last Dose  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: MED LAST DOSE: ANTICOAGULANT  
ID: {4404F3B5-B163-4F86-B2B1-E1E9E0687937}

=====

Page: 383

Abbreviation..... MLD AnCoag  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... String  
Description.....  
Display Name..... Med Last Dose: Anticoagulant  
Help Text..... Date and Time of Last Dose  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:  
-----

=====

TERM: MED LAST DOSE: ASPIRIN  
ID: {40F5DD1D-4534-463F-B4C8-5FF0D589C878}

=====

Page: 384

Abbreviation..... MLD Aspi  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... String  
Description.....  
Display Name..... Med Last Dose: Aspirin  
Help Text..... Date and Time of Last Dose  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: MED LAST DOSE: ESA ERYTHROPOETIN  
ID: {E63799D0-44A6-40A0-9968-FD0974BEDD17}

=====

Page: 385

Abbreviation..... MLD: ESA  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... String  
Description.....  
Display Name..... Med Last Dose: ESA Erythropoetin  
Help Text..... Enter last time/dose of medications  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:  
-----

=====

TERM: MED LAST DOSE: HYPOGLYCEMIC AGENT  
ID: {26E3B7EF-95A6-4ED3-A483-59DF78E4EB94}

=====

Page: 386

Abbreviation..... MLD: Hypog  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... String  
Description.....  
Display Name..... Med Last Dose: Hypoglycemic Agent  
Help Text..... Enter last time/dose of medications  
                                  for Diabetes  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Void.....

UNITS:  
-----





=====

TERM: MEWS: HEART RATE  
ID: {6D3266EC-E57D-4299-9AAC-80D0842C0D3F}

=====

Page: 389

Abbreviation..... MEWS: HR  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... MEWS: Heart Rate  
Help Text..... (Modified Early Warning System) Heart  
rate per minute  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 6 entries

-----

10. 2 - HR 0-39 (HR 0-39) {1243C1C9-C72F-4765-8D19-D689E1FC44F5}  
20. 1 - HR 40-50 (HR 40-50 ) {EC915B1F-BED1-497D-AD10-4FD0C2ADEDCE}  
30. 0 - HR 51-100 (HR 51-100) {6A262124-7FE8-4520-9762-87D0A3DB8AE9}  
40. 1 - HR 101-110 (HR 101-110)  
{FCB6C37E-1BB4-418F-8203-0916B4056DEB}  
50. 2 - HR 111-129 (HR 111-129)  
{6E1B9533-ECBE-4E2E-90D8-50C9C72AE0D3}  
60. 3 - HR 130 OR MORE (HR 130 +)  
{CBDB9BE1-832D-4F4A-9CF3-8768F054DDD2}

UNITS:  
-----



=====

TERM: MEWS: HOURLY URINE  
ID: {11CAC7DF-6CF0-404C-ADB5-747B1AF6F48E}

=====

Page: 390

Abbreviation..... MEWS: Urin  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... MEWS: Hourly Urine  
Help Text..... (Modified Early Warning System)  
Hourly urine for 2 hours  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Void.....

DATA VALUE: 3 entries

- 
- 10. 3 - LESS THAN 10MLS/HR (< 10mls/hr)  
{E013F8AC-7DBE-4F13-BF48-0317A46332CC}
  - 20. 2 - LESS THAN 30MLS/HR (< 30mls/hr)  
{A85A05F7-3443-4E62-ACF4-3334E3EF5675}
  - 30. 1 - LESS THAN 45MLS/HR (< 45mls/hr)  
{3AED380E-00F9-4C29-B303-52E08F8911B3}

UNITS:  
-----

=====

TERM: MEWS: RESPIRATORY  
ID: {BBB4C1F6-B418-4960-A3C7-87BAF9DB7595}

=====

Page: 391

Abbreviation..... MEWS: Resp  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... MEWS: Respiratory  
Help Text..... (Modified Early Warning System)  
Respiratory rate per minute  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

-----

10. 2 - RR 0-8 (RR 0-8) {73DC19E6-865C-4161-9CBE-3A8D44BC4D84}  
20. 0 - RR 9-14 (RR 9-14) {F32D69EF-7AC9-4A59-A63B-6E941BF5F06B}  
30. 1 - RR 15-20 (RR 15-20) {1CF2487B-30C6-46AB-A051-B4499F825DC2}  
40. 2 - RR 21-29 (RR 21-29) {72DF48C2-B8F3-4ED4-91D1-87157707C3EC}  
50. 3 - RR 30 OR MORE (RR 30 or +)  
{6BFA9D84-AAD7-4DC0-BA75-530CE5C12E77}

UNITS:  
-----

Abbreviation..... MEWS-Score  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... (Modified Early Warning System) MEWS  
Score  
Display Name..... MEWS: Score  
Help Text..... (Modified Early Warning System) MEWS  
Score  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Void.....

UNITS:

-----  
MEWS: SCORE UNIT            0 > 0 \* 15 < 15            Dec Prec: 0





Abbreviation..... Appli Remo
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description.....
Display Name..... Moderate sedation: Pre-assessment
appliances removed
Help Text.....
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 6 entries

- 1. DENTURES (Dentures) {ECABE4D0-97BE-418D-AC7D-29E40D166581}
2. EYEGASSES (Eyeglasses) {E4E5B291-CE6D-409F-A0A0-33EAA47C0C27}
3. CONTACTS (Contacts) {6BB5B9A8-6AB0-4EF3-B8B6-C4243AD1B806}
4. PROSTHESIS (PROSTHESIS) {CCD0EF23-FE9C-4B20-942E-9CF1CD7CF6E9}
5. HEARING AIDS (Hear Aids) {504FF149-8889-444E-925B-C952314FDF98}
6. OTHER (SEE COMMENTS) (Other) {7908D30B-8BD4-4A02-964D-8E6DBCA24A81}

UNITS:

Abbreviation..... Health HX
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description.....
Display Name..... Moderate Sedation: Pre-assessment
CPRS Health HX Reviewed
Help Text.....
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 15 entries

- 1. AICD (AICD) {EDB0E730-A1DC-4639-888B-30E956CCA958}
2. ARTIFICIAL VALVES (Art Valves) {088BB4A5-A93E-47E2-B21F-BC2957D52693}
3. BMI>35KG/M2 (BMI>35KG/M) {FD3D1D68-DD58-44B5-BEB2-0D4324F8353E}
4. DIABETES (Diabetes) {E295E970-759D-44CD-AD43-6047382A790D}
5. GRAFTS (Grafts) {E2AEE65A-7387-4A58-80E6-0212A1F3D1C7}
6. IMPLANTS (Implants) {14998852-14FD-4C6A-BB8A-5073D6094B11}
7. PACER (Pacer) {F8C0D639-D7B2-4342-A9FC-8A0F54BEB7DE}
8. KNOWN/SUSPECTED OBSTRUCTIVE SLEEP APNEA (Known/Susp) {BD47DD8C-7057-403B-B393-EBCBF9907B68}
10. METAL IMPLANTS (Metal Impl) {E7424C97-0A81-4C4E-916E-0AE812FFDB96}
11. PREGNANT, NO (Preg, No) {D08F5D48-965E-4537-83BC-14722E28B0AE}
12. PREGNANT, YES (Preg, Yes) {CF4323D7-9C8C-4FA5-9AB5-1AE4FFD1B535}
13. SEIZURES (Seizures) {DE974E49-2AB6-4781-9930-468F2B781BA9}
14. SHUNTS (Shunts) {8D9CB03B-6D63-4205-86C3-869421BEBD04}
15. SNORES (Snores) {69FE4ECF-E1FE-4235-AD44-5AF9E634E000}
16. STROKE (Stroke) {F0EE98A8-040F-4923-8A21-059DA7B17995}

UNITS:
-----





=====
TERM: MODERATE SEDATION: PREASSESS: PREVIOUS MOD SED. COMPLICATIONS Page: 398
ID: {6242F2E9-8402-44B6-B7C3-74B3106314AD}
=====

Abbreviation..... Prev Compl
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description.....
Display Name..... Moderate Sedation: PreAssess:
Previous Mod Sed. Complications
Help Text.....
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 6 entries

- 1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
2. HYPOTENSION (Hypotensio) {1AE4E3BB-3B78-42D2-BDAC-83563A0546E1}
3. HYPOXIA (Hypoxia) {7599629D-1206-4A1A-9D02-EEE384E5FE51}
4. APNEA (Apnea) {BF6A89F6-146B-4BE9-AD36-FF2E2C79022E}
5. ARRHYTHMIA (Arrhythmi) {F12E4852-27E7-4C2C-989F-E973EBDF2C62}
6. OTHER (SEE COMMENTS) (Other)
{7908D30B-8BD4-4A02-964D-8E6DBCA24A81}

UNITS:
-----



Abbreviation..... Mod Sed Pr  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... Moderate Sedation: Preparation  
Help Text.....  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 6 entries

- 
1. CODE CART IS ASSESSIBLE (Code Cart)  
{DB670E39-CF90-4068-8E38-44DCA5F365B5}
  2. BAG/MASK & O2 SET UP (Bag/Mask)  
{248ECF6D-D89B-4318-BFB7-D28CF5C9C5B0}
  3. FUNCTIONAL SECTION SET UP (Func Sect)  
{E156C048-11D5-4F62-8B1A-0CEA931E03EF}
  4. CARDIO-RESP MONITOR (Cardio-Re)  
{4402EC2C-2096-4DF7-BD4A-AE66C67035E0}
  5. PULSE OXIMETER (Pulse Ox) {773A956B-0B08-4D72-AAD4-630DF8268694}
  6. END-TIDAL CO2 MONITOR (End-Tidal)  
{609E3607-E441-465D-B1FC-C6986C58C11C}

UNITS:  
-----



=====

TERM: MONITOR PLACED  
ID: {FA699021-BD76-4DCA-B680-6E0F523B67FF}

=====

Page: 402

Abbreviation..... Monitor  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... Monitor Placed  
Help Text.....  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
1. CARDIAC (Cardiac) {2A4308E5-C843-4844-B908-32DA0A0B50E3}
  2. RESPIRATION (Respiratio) {5D6198F0-FE26-4F94-9118-D49A7EC44017}
  3. BLOOD PRESSURE (Blood Pres)  
{D6E2D19A-0B74-476A-90BF-B7C98CC7361A}
  4. PULSE OXIMETER (Pulse Ox) {D681E1E3-E914-421C-A19E-37C72F5E869F}

UNITS:  
-----

=====

TERM: MONITORING CAMERA  
ID: {3C959476-0C22-4C00-88D8-99C944D2D4A5}

=====

Page: 403

Abbreviation..... MON CAMERA  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Monitoring Camera  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. PATIENT INFORMED (PAT INFORM)  
    {074A1D19-3920-49F1-8958-2EB6F9CF05BA}
  2. FAMILY INFORMED (FAM INFORM)  
    {77135D32-9CDD-4D3F-9B62-6E8189B77844}
  3. UNABLE TO INFORM (UNABLE INF)  
    {77428C87-B3F9-4BBD-B781-F7D89BCF9492}

UNITS:  
-----

=====

TERM: MORSE: AMBULATORY AID DURING WALKING  
ID: {472ED5D6-2DFB-4EAF-9CCB-78D77C2D685B}

=====

Page: 404

Abbreviation..... MORSE AID  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Morse: Ambulatory Aid During Walking  
Display Name..... Morse: Ambulatory Aid During Walking  
Help Text..... 0: None/bedrest/Nurse Assis./  
Wheelchair, 15: Crutches/Cane/Walker,  
30: Furniture  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Void.....

DATA VALUE: 3 entries

-----

1. 0 (0)	{2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
2. 15 (15)	{FE8C049F-B64C-4991-2F50-3C561BA4825C}
3. 30 (30)	{1CF40F33-4B6B-FE03-98FB-55DFA8FB4145}

UNITS:

-----

=====

TERM: MORSE: GAIT

Page: 405

ID: {7CCE584F-DDAB-090B-9E73-41D5A8696F98}

=====

Abbreviation..... MORSE GAIT  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Morse: Gait  
Help Text..... 0: Normal/Bedrest/Wheelchair, 10:  
Weak, 20: Impaired  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

-----

1. 0 (0)	{2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
2. 10 (10)	{0C827160-18DB-816F-7FAA-A3C175C6D33B}
3. 20 (20)	{532FBDF8-7F7B-9551-300C-F2CFEF656163}

UNITS:

-----





=====

TERM: MORSE: IV THERAPY SALINE LOCK  
ID: {5FB26C06-AF8D-ECBE-FC10-442510EC52F1}

=====

Page: 407

Abbreviation..... MORSE IV  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Morse: IV Therapy Saline Lock  
Help Text..... 0 if NO, 25 if YES  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- 1. 0 (0) {2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
  - 2. 20 (20) {532FBDF8-7F7B-9551-300C-F2CFEF656163}

UNITS:  
-----

=====

TERM: MORSE: MENTAL STATUS OF PT OWN ABILITY  
ID: {E15F905B-5D1C-47CA-A23A-5EC47302D9A3}

=====

Page: 408

Abbreviation..... MORSE PT S  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Morse: Mental Status of Patients Own  
Ability  
Display Name..... Morse: Mental Status of Patients Own  
Ability  
Help Text..... 0: Oriented to own ability, 15  
Overestimates, forgets limitations  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

-----

1. 0 (0)	{2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
2. 15 (15)	{FE8C049F-B64C-4991-2F50-3C561BA4825C}

UNITS:  
-----

=====

TERM: MORSE: SCORE TOTAL  
ID: {80F2ACB5-5E42-48B8-8EB8-CADAD149E7BB}

=====

Page: 409

Abbreviation..... MORSE TOTA  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Morse Score Total  
Display Name..... Morse: Score Total  
Help Text..... Morse Score Total  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MORSE SCORE                    -1 > 0 \* 135 < 136      Dec Prec: 0

=====

TERM: MORSE: SECONDARY DIAGNOSIS LISTED ON CHART  
ID: {305B5629-D6B2-440A-B46F-E0535AE843F9}

=====

Page: 410

Abbreviation..... MORSE 2NDA  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Morse: Secondary Diagnosis Listed On C  
Help Text..... 0 = No, 25 = Yes  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- 1. 0 (0) {2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
  - 2. 15 (15) {FE8C049F-B64C-4991-2F50-3C561BA4825C}

UNITS:  
-----

Abbreviation..... MOTOR ACTI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Motor Activity  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 6 entries

- 
1. NORMAL POWER (np) {0E45C061-0744-8A21-7256-003702A5005B}
  2. NO MOVEMENT (nm) {2618E301-379B-BAB9-9EDE-98FD4B60D311}
  3. FLICKER OF CONTRACTION (foc)  
{218B322B-84AA-F6F6-5207-D86F9A4311F0}
  4. MOVEMENT WITH GRAVITY ELIMINATED (mwge)  
{B86E1131-1976-FA8E-74FE-B1BC2F0FDE3F}
  5. MOVE AGAINST GRAVITY BUT NO RESISTANCE (magbnr)  
{F995A165-8536-C493-2269-286E8EFF5E43}
  6. MOVE AGAINST GRAVITY AND RESISTANCE, EASILY OVERCOME BY THE EXAMINER (magareobte)  
{2ECCC197-4A81-025F-3EF3-6214004887A9}

LOCATION: 4 entries

- 
1. L LOWER EXTREMITY (LLEX) {FC6F5056-779A-BF76-86D0-A237936C1ECC}
  2. L UPPER EXTREMITY (LUEX) {48806E79-B1DA-129B-EB52-453292A3DDC6}
  3. R LOWER EXTREMITY (R LOW EXTR)  
{7FC3CECC-1E6E-A686-FF04-2BEF8A9B9DCA}
  4. R UPPER EXTREMITY (R UPP EXTR)  
{C65BA269-8F05-6FB9-62CB-466A286CF2C0}

UNITS:

-----

=====

TERM: MOTOR CHECK

Page: 412

ID: {521171CF-BD9D-2F28-4042-EA885C40E556}

=====

Abbreviation..... MOCH  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Motor Check  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- 1. INTENTIONAL (i) {D5BBBFFC-F739-1E69-FED4-288D8E803D04}
  - 2. REFLEX (r) {F21585F9-5206-A6B3-2BF6-F416E5842366}

UNITS:

-----

=====

TERM: MUCOUS MEMBRANE CONDITION  
ID: {4D9B0E42-4A15-4699-A80E-31588D8C13FD}

=====

Page: 413

Abbreviation..... MUC MEMBRA  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Mucous Membrane Condition  
Display Name..... Mucous Membrane Condition  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
- |                 |  |
|-----------------|--|
| 1. PINK (PINK)  | {1A7FAF84-D69C-4EEB-7BC5-D21A598A1F4D} |
| 2. CYANOTIC (c) | {06BD77DB-1849-D731-0F27-F8CE48755EA9} |
| 3. PALE (p)     | {0EB39D28-7ABE-A4AB-D334-3F795E65618D} |

UNITS:  
-----



Abbreviation..... NAIL COLOR  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Nail Color  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. NORMAL FOR ETHNIC GROUP (ETH. NORM.)  
{3B2CD951-58CE-FF53-AD5D-6C65C07658D3}
  2. CYANOTIC (c) {06BD77DB-1849-D731-0F27-F8CE48755EA9}
  3. DISCOLORATION (DISCOLORAT)  
{BBE3E910-CC49-430C-BD90-FB0A97F1B183}
  4. PALE (p) {0EB39D28-7ABE-A4AB-D334-3F795E65618D}
  5. HEMORRHAGE (HEMORRHAGE) {E00AC828-83F3-4F69-B97B-D6AD1D87A440}

UNITS:  
-----

=====

TERM: NAIL SHAPE

Page: 415

ID: {6351A78B-67BC-46A4-B7E4-72D63C8C2917}

=====

Abbreviation..... NAIL SHAPE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Nail Shape  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. NORMAL (n) {553C9A97-A3D5-C52F-0E4E-C7D168C99A2B}
  2. CLUBBING (CLUBBING) {AE77D344-6F0F-452B-B1AA-EC8D9F092CFC}
  3. SPOONING (SPOONING) {CA384C9F-D87D-4E2D-B4C0-89BD85692B16}

UNITS:

-----

=====

TERM: NAIL TEXTURE

Page: 416

ID: {AD3F008F-6AA8-483F-AAA6-5F797109F4A6}

=====

Abbreviation..... NAIL TEXTU  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Nail Texture  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
- |                    |  |
|--------------------|--|
| 1. NORMAL (n)      | {553C9A97-A3D5-C52F-0E4E-C7D168C99A2B} |
| 2. THICK (THICK)   | {1554A49E-1468-46B2-F851-531DDB1D1A98} |
| 3. PITTED (PITTED) | {BEA42FD8-91C7-4BF6-936E-726D78064F11} |

UNITS:

-----

=====

TERM: NECK VEIN CONDITION  
ID: {D76459ED-8572-409E-85A4-F64228AE9A36}

=====

Page: 417

Abbreviation..... NECK VEIN  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Neck Vein Condition  
Display Name..... Neck Vein Condition  
Help Text..... Condition of Neck Veins  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. NORMAL (n) {553C9A97-A3D5-C52F-0E4E-C7D168C99A2B}
  2. FLAT (Flat) {EA4F1980-EDF3-82BC-B677-D747F48F1163}
  3. DISTENDED (DISTENDED) {B5C529CA-2BE0-4698-A153-718F4D58E643}

UNITS:

-----

=====

TERM: NEED FOR FLUIDS/TOILETING ASSESSED  
ID: {816D5DB3-7AEA-E73F-08D8-68C58D4FB4FE}

=====

Page: 418

Abbreviation..... NFFTA  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Need for Fluids/Toileting Assessed  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEA649CD} |

UNITS:  
-----

=====

TERM: NOSE SIZE

Page: 419

ID: {0A7CD784-279C-9F56-67CC-1C0624B26F4A}

=====

Abbreviation..... ns  
Active..... No  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Describe the size of the honker  
Display Name..... Nose Size  
Help Text..... Honker etiquette  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: NUMBER OF DAYS  
ID: {C2D0AF1E-C88B-48ED-A8BE-E3C258B4317D}

=====

Page: 420

Abbreviation..... # of Days  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description.....  
Display Name..... Number of Days  
Help Text..... Number of days since event or device  
                          days  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
DAY                   -1 > -1 \* 9999 < 9999   Dec Prec: 0

ID: {656AA278-6C1E-C1B8-6876-A85780A1608A}

=====

Abbreviation..... ORIENT  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Orientation  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 8 entries

-----

1. ORIENTED TO TIME (ott) {C79EBE25-C51B-478C-CB39-25943792F191}
2. ORIENTED TO PLACE (otp) {8E5933E2-2DF0-C8A4-7A92-20E792F42EA6}
3. ORIENTED TO PERSON (otp) {9912CFB9-AC55-6445-3EC1-0F1DF98C0561}
4. ORIENTED TO SITUATION (ots)  
{955A2F3B-E444-99DF-81DD-5E6BF9D22E49}
5. DISORIENTED TO TIME (dtt) {2AD1A28A-4121-4B82-CFB8-E5A8D3095704}
6. DISORIENTED TO PLACE (dtp)  
{00C98809-C0D0-C7B4-0615-8F7D31017FE8}
7. DISORIENTED TO PERSON (DISOR TO P)  
{E5D83C03-7EF5-770B-930F-7B3B3345048C}
8. DISORIENTED TO SITUATION (dts)  
{7143B3D7-5569-4378-D3AA-10C4675F0BA5}

UNITS:

-----





=====

TERM: OSTOMY CARE

Page: 423

ID: {BB127CE7-04B9-48BA-85CB-AD895D3DA32D}

=====

Abbreviation..... Ostomy Car  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... Ostomy Care  
Help Text.....  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
- 10. APPLIANCE CHANGED (app change)  
{CB4BC887-F738-4EE3-9A5A-39C9CDF7488C}
  - 20. BRIDGE PRESENT (bridge pre)  
{5C474C17-33C9-492D-833E-64F5EF049BF4}
  - 30. CONNECTED TO DRAINAGE BAG (drain bag)  
{1901C4A7-671C-43BD-BEDD-1D48BCB43465}
  - 40. PATIENT INSTRUCTED (pat instru)  
{D89A50A6-352D-4EB8-ABD7-5D732A060D9A}
  - 50. FAMILY INSTRUCTED (fam instru)  
{0FD5988C-2C4D-4786-99B3-3287A355E813}

UNITS:  
-----

=====

TERM: OSTOMY TYPE

Page: 424

ID: {CA5CC473-EE95-4E81-8E24-DC23B47A7F42}

=====

Abbreviation..... OSTOMY TYP  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Type of Ostomy  
Display Name..... OSTOMY TYPE  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. COLOSTOMY (COLOSTOMY) {53F939AB-6AF4-4DD5-8D64-17AC1D1736B7}
  2. JEJUNOSTOMY (JEJUNOSTOM) {1ED3EC41-3ACF-42C1-B64A-280300EDBCFC}
  3. SIGMOIDOSTOMY (SIGMOIDOST)  
{0C84B317-59B3-433B-A9AD-6588C1361E0C}
  4. ILEOSTOMY (ILEOSTOMY) {D5DF49C8-C6BA-44C1-9840-7765907730A3}
  70. ILEAL-CONDUIT (Ileal-Con) {C91C4E4A-2D0D-438A-8EF5-8717D703A777}

UNITS:

-----

=====

TERM: OUTPUT - BLOOD  
ID: {F90805CD-A0E1-4943-B871-D93B445EA639}

=====

Page: 425

Abbreviation..... OUT BLOOD  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Output of blood lost  
Display Name..... Output - Blood  
Help Text..... Output of blood lost  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MILLILITER	-1 > 0 * 9999 < 10000	Dec Prec: 0
CUBIC CENTIMETER	-1 > 0 * 9999 < 10000	Dec Prec: 0

=====

TERM: OUTPUT - CBI DRAINAGE  
ID: {08B2CC2A-AE5A-47C3-9FFE-CD0619F8B3BE}

=====

Page: 426

Abbreviation..... CBI Drain  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description.....  
Display Name..... Output - CBI Drainage  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MILLILITER            0 > 0 \* 99999 < 99999    Dec Prec: 0

=====

TERM: OUTPUT - CHEST TUBE  
ID: {478FE098-139E-B9EB-9F87-C03A847B6789}

=====

Page: 427

Abbreviation..... OUT CHE TU  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... VOLUME OF CHEST TUBE DRAINAGE  
Display Name..... Chest Tube Drainage Volume (ml)  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MILLILITER	-1.0 > 0.0 * 9999.0 < 10000.0
	Dec Prec: 1
CUBIC CENTIMETER	-1.0 > 0.0 * 9999.0 < 10000.0
	Dec Prec: 1

Abbreviation..... OUT CHTU E  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... AMOUNT OF CHEST TUBE DRAINAGE  
Display Name..... Chest Tube Drainage Amount Estimated  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
  2. SCANT (SCANT) {8EF6C62D-0D94-37BB-60AD-86E306E8120D}
  3. SMALL (s) {5991F892-B28C-3998-E733-D8FFA26D7E21}
  4. MODERATE (MODERATE) {C3C25064-C9FC-EA1F-3E74-32A25F0FC5E3}
  5. LARGE (l) {E60D5CF9-920D-44E2-8A83-136DAB1FE251}

UNITS:  
-----

ID: {C7954A6C-4214-4CB2-B80E-8F4128D05847}

=====

Abbreviation..... OUT DIALYS  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Dialysis Output  
Display Name..... Dialysis Output  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 4 entries

- 
1. HEMODIALYSIS (HEMODIALYS) {DE7CD394-9E69-46B0-A207-5A18F7C5A1A7}
  2. PERITONEAL DIALYSIS (PERI DIALY)  
{6B4AA38B-AFCC-42CE-9F32-59B32E8A6BBB}
  3. CVVHD (CVVHD) {4CC9AB0D-0B6C-43F9-A93E-A8228EC3A77F}
  4. CRRT (CRRT) {AB37FC1D-DD7A-4103-A777-7B8D0F876889}

UNITS:

-----

MILLILITER -1 > 0 \* 9999 < 10000 Dec Prec: 0  
CUBIC CENTIMETER -1 > 0 \* 9999 < 10000 Dec Prec: 0



Abbreviation..... OUTPUT
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Output - Do Not Check Total or
Subtotals
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

PRODUCT: 10 entries

- 1. CARDIAC OUTPUT INJECTATE (CARD OUT I)
{2128B705-DDF4-4776-87C3-AA82F854E313}
2. EF RESIDUAL (EF RESID) {35E269DF-F246-4DF6-9D68-7FD61F8FE088}
3. NG RESIDUAL (NG RESID) {B03C9557-84B5-4B05-B784-C87FA90D27A7}
4. HEME (HEME) {A0ADA229-E779-4562-84B5-67121B88C443}
5. HEMODIALYSIS (HEMODIALYS) {8409B357-C41E-4F4A-BB52-F09AD9704C30}
6. MEDIASTINAL (MEDIASTINA) {A0550B31-B937-4FC7-8F7E-071F75215311}
7. BLADDER IRRIGATION (BLAD IRRIG)
{BCBB106B-ECE4-4232-87C7-D90EDB47F80D}
8. WOUND VACUUM (WOUND VAC) {9F4983A6-B622-4E4D-A7F7-ADA7A4B7CEBC}
9. GASTRIC DRAINAGE (GAST DRAIN)
{C01D18F9-191A-4D34-B45D-D8D70ACA67BE}
10. WOUND DRAINAGE (WOUN DRAIN)
{73038CBA-1F31-47E2-A78E-50ED1067CBEA}

UNITS:

MILLILITER -1 > 0 \* 9999 < 10000 Dec Prec: 0
OUNCE -1 > 0 \* 338 < 339 Dec Prec: 0

=====

TERM: OUTPUT - DRAIN

Page: 431

ID: {B8FD2F72-FCDD-484A-98E1-14B3637B797A}

=====

Abbreviation..... OUT DRAIN  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Drain Output  
Display Name..... Drain Output  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 5 entries

- 
1. JP DRAIN (JP DRAIN) {94A32247-6190-4047-B581-A212FE4C9453}
  2. HEMOVAC (HEMOVAC) {79DA0646-B811-4954-8BE1-C133398CD574}
  3. PENROSE (PENROSE) {CBC3C6DE-E5DE-41D5-83AB-30A4F3CFA9E1}
  4. CHEST TUBE (CH TUBE M) {F4F467D3-BB57-4E6F-9EA0-C0FBFB66F5B1}
  5. OTHER (OTHER) {AED275BD-8A42-46BE-BCDF-472EF4C31D65}

UNITS:

-----

MILLILITER            -1 > 0 \* 9999 < 10000    Dec Prec: 0  
CUBIC CENTIMETER    -1 > 0 \* 9999 < 10000    Dec Prec: 0

=====

TERM: OUTPUT - EMESIS

Page: 432

ID: {56BD11AB-FE61-4785-B3CA-C528A5F4EBD6}

=====

Abbreviation..... OUT EMESIS  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Emesis Output  
Display Name..... Emesis Output  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

PRODUCT: 1 entries

-----

1. EMESIS (EMESIS)            {A7B1605E-1DFE-4FC6-A69B-0E18C1A291FF}

UNITS:

-----

MILLILITER            -1 > 0 \* 9999 < 10000    Dec Prec: 0  
CUBIC CENTIMETER    -1 > 0 \* 9999 < 10000    Dec Prec: 0

=====

TERM: OUTPUT - GI TUBE

Page: 433

ID: {CB54872E-CF17-4BA5-BF3B-80B0E6C745F8}

=====

Abbreviation..... OUT GI TU  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... GI Tube Output  
Display Name..... GI Tube Output  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MILLILITER	-1 > 0 * 9999 < 10000	Dec Prec: 0
CUBIC CENTIMETER	-1 > 0 * 9999 < 10000	Dec Prec: 0

=====

TERM: OUTPUT - STOOL  
ID: {8B7CD3AF-5441-40A7-BBAA-11B822B529FD}

=====

Page: 434

Abbreviation..... OUT STOOL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Measured amount of stool  
Display Name..... Stool Output  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MILLILITER	-1 > 0 * 1000 < 2000	Dec Prec: 0
CUBIC CENTIMETER	-1 > 0 * 1000 < 2000	Dec Prec: 0

Abbreviation..... OUT ST EST
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Stool Amount Estimated
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

- 1. SMEAR (SMEAR) {4B99BBD3-443E-4954-8ADE-5FB3EE50540A}
2. SMALL (s) {5991F892-B28C-3998-E733-D8FFA26D7E21}
3. MODERATE (MODERATE) {C3C25064-C9FC-EA1F-3E74-32A25F0FC5E3}
4. LARGE (l) {E60D5CF9-920D-44E2-8A83-136DAB1FE251}
5. COPIUS (COPIUS) {C74A6979-0E35-AA12-F992-69ACA3B0FA23}

QUALITY: 1 entries

- 1. ESTIMATED (E) {9EAF26B-7D87-DA66-F991-022D095FB5A0}

UNITS:

=====

TERM: OUTPUT - URINE

Page: 436

ID: {F4998CFC-19CA-402A-AB0F-C32FA194DAA5}

=====

Abbreviation..... OUT URINE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Urine Output  
Display Name..... Urine Output  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MILLILITER	-1 > 0 * 9999 < 10000	Dec Prec: 0
CUBIC CENTIMETER	-1 > 0 * 9999 < 10000	Dec Prec: 0

=====

TERM: OUTPUT - URINE ESTIMATED  
ID: {C5116BFA-332E-4C05-8FF7-4C828D8999BF}

=====

Page: 437

Abbreviation..... OUT UR EST  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Output - Urine Estimated  
Display Name..... Output - Urine Estimated  
Help Text..... Output - Urine Estimated  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
1. SLIGHT (SLIGHT) {805549C1-5DC5-4A2D-BC1F-046E44928843}
  2. SMALL (s) {5991F892-B28C-3998-E733-D8FFA26D7E21}
  3. MODERATE (MODERATE) {C3C25064-C9FC-EA1F-3E74-32A25F0FC5E3}
  4. LARGE (l) {E60D5CF9-920D-44E2-8A83-136DAB1FE251}

UNITS:  
-----



=====

TERM: OXYGEN CONSUMPTION  
ID: {5EC773D2-D17D-EE6B-5079-BED57318658D}

=====

Page: 438

Abbreviation..... V02  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Oxygen Consumption  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MILLILITERS PER MINUTE PER METER SQUARE  
0 > 1 \* 300 < 301      Dec Prec: 0



=====

TERM: PA CATHETER ALARMS ON  
ID: {8129EB35-E85A-4C3B-86AF-34FB2ABE168C}

=====

Page: 440

Abbreviation..... PA ALARM O  
Active..... Yes  
Boolean Value False..... Off  
Boolean Value True..... On  
Data Type..... Boolean  
Description..... PA Catheter Alarms On  
Display Name..... PA Catheter Alarms On  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

=====

TERM: PA CATHETER DRESSING TYPE  
ID: {16F7EF6A-E0B6-9931-1F74-039E9D18AC57}

=====

Page: 441

Abbreviation..... PA DRESSIN  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... PA Catheter Dressing Type  
Help Text..... Type of dressing on PA Catheter site  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
1. TRANSPARENT (TRANSPAREN) {0E73D24D-0EFC-08F9-196A-D30DD1BB3228}
  2. PRESSURE DRESSING (PRESS DRES)  
{6D9DB9EF-85EB-8321-03FC-879AA4C5615B}

UNITS:  
-----

=====

TERM: PA CATHETER DYNAMIC RESPONSE TEST  
ID: {7ED7B728-5683-4660-B301-79E72C3FD1C9}

=====

Page: 442

Abbreviation..... PA CATH DR  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... PA Catheter Dynamic Response Test  
Display Name..... PA Catheter Dynamic Response Test  
Help Text..... PA Catheter Dynamic Response Test  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. NORMAL (n) {553C9A97-A3D5-C52F-0E4E-C7D168C99A2B}
  2. DAMPENED (DAMPENED) {77C777B3-7E56-41E3-8E70-DE144E6A3424}
  3. HYPERDYNAMIC (HYPERDYNAM) {818A8A99-04D9-4989-A9BA-ECFACED74605}

UNITS:  
-----

=====

TERM: PA CATHETER INSERTION MARK  
ID: {7E62F0A7-9233-4BD3-870B-189BC86A179E}

=====

Page: 443

Abbreviation..... PA INS MAR  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... PA Catheter Insertion Mark  
Display Name..... PA Catheter Insertion Mark (cm)  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
CENTIMETER            -1.0 > 0.0 \* 10.0 < 20.0    Dec Prec: 1

=====

TERM: PA CATHETER INSERTION TECHNIQUE  
ID: {30F7B681-FCDB-4349-9082-92BDF3618CA3}

=====

Page: 444

Abbreviation..... PACATH INS  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... PA Catheter Insertion Technique  
Display Name..... PA Catheter Insertion Technique  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
1. NEW SITE (NEW SITE) {10B1FF61-E58E-4608-9CCB-9232BC964418}
  2. OVER WIRE (OVER WIRE) {157CD456-1020-43AE-9AAA-1686F1574B6F}
  3. INSERTED THROUGH EXISTING INTRODUCER (INS EXS IN)  
{20BC183A-197F-4464-A27C-47F63C5B842A}
  4. INSERTED THROUGH NEW INTRODUCER (INS NEW IN)  
{7FCA403D-41E5-4BAA-B155-B58F35EA49B3}

UNITS:  
-----

=====

TERM: PA CATHETER REPOSITIONED/DISCONTINUED REASON  
ID: {959199C6-C6DE-4F32-9303-F06CD0C80A28}

=====

Page: 445

Abbreviation..... PA CATH RE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... PA Catheter Repositioned/Discontinued  
Reason  
Display Name..... PA Catheter Repositioned/Discontinued  
Reason  
Help Text..... PA Catheter Repositioned/Discontinued  
Reason  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 6 entries

- 
1. VENTRICULAR WHIP PRESENT (VENT WHIP)  
{15CB9645-690D-446B-A380-F14CEEC9C6AA}
  2. DAMPENED (DAMPENED) {77C777B3-7E56-41E3-8E70-DE144E6A3424}
  3. WEDGED (WEDGED) {03D7F6FD-399D-4A6B-9900-59655D76EE95}
  4. UNABLE TO WEDGE (UNABLE WED)  
{5D35F546-F1EA-49C6-8B1D-22D24526AF8B}
  5. WAVEFORM REFLECTS ATRIAL PLACEMENT (WAVE ATRIA)  
{A1249745-41C2-4A4F-A86D-7E5E17477B06}
  6. BY PROVIDER ORDER (BY ORDER)  
{C038BC16-44CF-4933-BECC-B0F1A1B12603}

UNITS:  
-----



=====

TERM: PA CATHETER SITE CARE  
ID: {F5F81475-6135-4B0C-A144-658AB8C075C7}

=====

Page: 446

Abbreviation..... PA SITE CA  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... PA Catheter Site Care  
Display Name..... PA Catheter Site Care  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

=====

TERM: PA CATHETER SITE CONDITION  
ID: {00A7358A-DBE8-CB1D-FFC0-046F3C5E266F}

=====

Page: 447

Abbreviation..... PA CONDITI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... PA Catheter Site Condition  
Help Text..... Condition of PA catheter  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 7 entries

-----

2. RED (RED)	{515D7CFB-3527-30E7-EC4B-B607C272DFA0}
3. SWOLLEN (SWOLLEN)	{DCE86B10-C393-8DA9-6518-F2943F0E1A4E}
4. DISCOMFORT (DISCOMFORT)	{77D058C9-1EB4-96C0-C9E3-666C840B884C}
50. PATENT (Patent)	{F0C38BCB-1D11-436C-8B8C-4297C27DDE04}
60. CLEAN (Clean)	{1F187B99-9963-4AD5-98F6-D8D24D2DDD30}
70. OCCLUDED (Occluded)	{25509BE3-F946-4AD7-80E8-5CE8A2E61AD6}
80. OOZING (Oozing)	{C23549D7-3238-4523-BA31-AB71609F01C2}

UNITS:  
-----

=====

TERM: PA CATHETER SIZE

Page: 448

ID: {706045DE-B1D4-46DB-B6CE-A171BAFB0CB9}

=====

Abbreviation..... PA CATH SI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... PA Catheter Size  
Display Name..... PA Catheter Size  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
FRENCH                    -1.0 > 5.0 \* 10.0 < 20.0   Dec Prec: 1

=====

TERM: PA CATHETER TYPE  
ID: {A1F96B2C-AA70-0D43-8F1E-8DD2BDDFA8F3}

=====

Page: 449

Abbreviation..... PA TYPE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... PA Catheter Type  
Help Text..... Type of PA catheter  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
- 1. PLAIN (PLAIN) {483EE48F-F72B-A917-DCD0-E8B69A3FA36F}
  - 2. OXIMETRIC (OXIMETRIC) {DA1F5AAA-A349-C17B-74CD-1794492E0A13}
  - 3. PACER (PACER) {79499460-ACAB-33B3-52CF-D2D2E5BCD0EE}
  - 4. VIP PORT (VIP PORT) {3D0426B3-F8F3-BAA4-71C7-8A09F911EF08}

LOCATION: 8 entries

- 
- 1. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
  - 2. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
  - 3. L INTERNAL JUGULAR (LIJ) {AFBB4FE5-8A11-2047-7900-DADF2FC4430D}
  - 4. R INTERNAL JUGULAR (R INT JUGU)  
{7ABF6916-D379-198D-6454-F0518A74844A}
  - 5. L EXTERNAL JUGULAR (L EXT JUGU)  
{6F7B331B-4F83-0B79-1755-F6F673D647D6}
  - 6. R EXTERNAL JUGULAR (R EXT JUGU)  
{4D78F704-9E3B-907E-A792-4DC4F4CFF170}
  - 7. L SUBCLAVIAN (LSCL) {71B93594-15F9-AA4A-AEF2-E222EC28C00D}
  - 8. R SUBCLAVIAN (R SUBCLAVI) {53DEAE83-93AE-123B-14E5-4DE38CD4F752}

UNITS:  
-----

=====

TERM: PA CATHETER WAVEFORM  
ID: {9EF25DC1-2B10-4248-807E-DDED1B76ABC}

=====

Page: 450

Abbreviation..... PA WAVEFOR  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... PA Catheter Waveform  
Display Name..... PA Catheter Waveform  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. NORMAL (n) {553C9A97-A3D5-C52F-0E4E-C7D168C99A2B}
  2. DAMPENED (DAMPENED) {77C777B3-7E56-41E3-8E70-DE144E6A3424}
  3. HYPERDYNAMIC (HYPERDYNAM) {818A8A99-04D9-4989-A9BA-ECFACED74605}
  4. POSITIONAL (POSITIONAL) {9A74D16B-5E95-4EA8-88AF-E53C80311A5D}
  5. NOT MONITORED (NOT MONITO)  
{0E6A8A07-48E0-428F-8D86-22572CABCF48}

UNITS:  
-----

=====

TERM: PA CATHETER ZEROED  
ID: {E8B0247F-BD12-4902-87F7-012AA4C1D73D}

=====

Page: 451

Abbreviation..... PA ZEROED  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... PA Catheter Zeroed  
Display Name..... PA Catheter Zeroed  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

=====

TERM: PACEMAKER ATRIAL mA  
ID: {2551E3D7-CD46-909E-1D79-082B92AEDF7A}

=====

Page: 452

Abbreviation..... PAmA  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Pacemaker Atrial mA  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MILLIAMPS            -0.01 > 0.00 \* 24.00 < 48.00  
                              Dec Prec: 2

=====

TERM: PACEMAKER AV INTERVAL  
ID: {118070EB-2F23-9C3D-68D1-B57D423449E1}

=====

Page: 453

Abbreviation..... PAVI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Pacemaker AV Interval  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

SECOND                    0.01 > 0.12 \* 0.20 < 0.50 Dec Prec: 2



=====

TERM: PACEMAKER CAPTURE  
ID: {4D93B7C5-5302-085D-DFC5-EBEED03A6B31}

=====

Page: 454

Abbreviation..... PMC  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Pacemaker Capture  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEA649CD} |

UNITS:  
-----

Abbreviation..... PMM
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Pacemaker Mode
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 13 entries

- 1. AOO (AOO) {86F25E48-AADB-5B30-4B89-E8BE63DAE7D9}
2. AAO (AAO) {E87E0670-2C0F-6E6A-DEC9-F52E5C530A2A}
3. AAI (AAI) {2DD18ACA-B98E-24D3-E760-6D082A76074B}
4. AAIR (AAIR) {3B4B6F37-B631-BE32-F1C1-65B99D0C5018}
5. AVO (AVO) {02EA752A-D156-96E6-3D28-3D389FFC7B58}
6. AVI (AVI) {8BB49456-FB80-4916-6EBE-06837BF12B0C}
7. VVI (VVI) {2C87F8F3-2B99-80E7-04EB-E6FE88168627}
8. VDD (VDD) {E733E525-EF36-2F93-AECC-BAD1C3E2C633}
9. DVI (DVI) {4D6E3081-7BB5-45DB-E8E6-0D3D7869E6BD}
10. DOO (DOO) {749BCE3D-FF09-B796-FF2A-ECD2D116C96F}
11. DDI (DDI) {0C1E359F-B1F4-7994-D337-37622CB14BC8}
12. DDD (DDD) {7AA7152E-48F0-8550-AD0D-9F298E90B765}
13. DDDR (DDDR) {A3CC0D9D-D329-1A52-2DD1-911A1E04B871}

UNITS:

ID: {A7D28B32-845F-0E16-FDAF-F654FA6EE8A4}

=====

Abbreviation..... PMO  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Pacemaker On  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
  - 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

LOCATION: 8 entries

- 
- 1. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
  - 2. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
  - 3. L BRACHIAL (LBRA) {75F6CE33-1AB3-A1DA-70B8-04825290B577}
  - 4. R BRACHIAL (RBRA) {A88E69A8-F6AE-80B0-2CD4-687B121F9DDE}
  - 5. L INTERNAL JUGULAR (LIJ) {AFBB4FE5-8A11-2047-7900-DADF2FC4430D}
  - 6. R INTERNAL JUGULAR (R INT JUGU)  
{7ABF6916-D379-198D-6454-F0518A74844A}
  - 7. L SUBCLAVIAN (LSCL) {71B93594-15F9-AA4A-AEF2-E222EC28C00D}
  - 8. R SUBCLAVIAN (R SUBCLAVI) {53DEAE83-93AE-123B-14E5-4DE38CD4F752}

UNITS:

-----



=====

TERM: PACEMAKER SENSE  
ID: {7F716D0A-ACA4-7E33-62D7-988BA22F0D7F}

=====

Page: 458

Abbreviation..... PMS  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Pacemaker Sense  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

=====

TERM: PACEMAKER SENSITIVITY ATRIAL  
ID: {6CBB234D-694E-CB88-3844-6B9E9E262320}

=====

Page: 459

Abbreviation..... PSA  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Pacemaker Sensitivity Atrial  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
MILLIVOLTS            0.01 > 0.05 \* 2.00 < 6.00 Dec Prec: 2

=====

TERM: PACEMAKER SENSITIVITY VENTRICULAR  
ID: {305A8823-01D3-F3FE-D66E-3CACE736CD30}

=====

Page: 460

Abbreviation..... PSV  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Pacemaker Sensitivity Ventricular  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MILLIVOLTS            0.01 > 0.05 \* 2.00 < 6.00 Dec Prec: 2

=====

TERM: PACEMAKER TYPE  
ID: {3B5B64D7-9091-8511-64CE-848456CAD2CD}

=====

Page: 461

Abbreviation..... PACEM TYPE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Pacemaker Type  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 6 entries

- 
1. TRANSCUTANEOUS (t) {F4315F92-8804-1CB4-4D03-CCFDC95C6E71}
  2. PERMANENT (p) {A31FF319-F8C3-7AAC-D7FC-765C89732D09}
  3. TEMPORARY (TRANSVENOUS) (tt)  
{F6A74188-3D79-D749-540F-4EA32201B0E6}
  4. EPICARDIAL LEADS (el) {2B551575-3D50-96BD-7BAA-50CC9796361E}
  5. ATRIAL (a) {266D2FDA-5D6A-203D-6FF0-8282D49BFD67}
  6. VENTRICULAR (v) {307A4572-3CF0-7875-7313-7FD4AB23F8E7}

LOCATION: 8 entries

- 
1. L BRACHIAL (LBRA) {75F6CE33-1AB3-A1DA-70B8-04825290B577}
  2. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
  3. L INTERNAL JUGULAR (LIJ) {AFBB4FE5-8A11-2047-7900-DADF2FC4430D}
  4. L SUBCLAVIAN (LSCL) {71B93594-15F9-AA4A-AEF2-E222EC28C00D}
  5. R BRACHIAL (RBRA) {A88E69A8-F6AE-80B0-2CD4-687B121F9DDE}
  6. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
  7. R INTERNAL JUGULAR (R INT JUGU)  
{7ABF6916-D379-198D-6454-F0518A74844A}
  8. R SUBCLAVIAN (R SUBCLAVI) {53DEAE83-93AE-123B-14E5-4DE38CD4F752}

UNITS:

-----



=====

TERM: PACEMAKER UNIT ID  
ID: {CDA2992F-349C-9904-45E3-20D01DC87917}

=====

Page: 462

Abbreviation..... PUID  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... String  
Description..... No Description  
Display Name..... Pacemaker Unit ID  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: PACEMAKER VENTRICULAR mA  
ID: {D64D8FED-93B8-FABD-5361-96CD28F99472}

=====

Page: 463

Abbreviation..... PVmA  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Pacemaker Ventricular mA  
Help Text..... Pacemaker Ventricular mA  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MILLIAMPS            -0.01 > 0.00 \* 24.00 < 48.00  
                              Dec Prec: 2

Abbreviation..... PMWI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Pacemaker Wires
Help Text..... Pacemaker Wires
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 7 entries

- 1. 1A +1V (1+) {31711946-9F1B-0AD9-BDBE-846EDBC86AE2}
2. 1A + 2V (1+2) {FE47CEC3-C3A5-0F5C-3EE8-50F30FD4E0C4}
3. 2A (2) {C75BDC9F-C756-9051-CAA3-A27DDD94B35E}
4. 2V (2) {D24B5FA4-1DA6-36DE-49DB-627A9B82E968}
5. 2A + 1V (2+1) {1849AA0D-1572-55E0-A780-ACE602E2C72D}
6. 2A + 2V (2+2) {17203251-68E8-A00A-CF54-E2CC55B8EEDD}
7. 1A + 3V (1+3) {CEDF220D-E593-B8F9-DB08-D32800AF3F5E}

UNITS:
-----

Abbreviation..... P DIS C
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... PACU - Discharge Criteria
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 22 entries

- 1. TOLERATING FLUIDS/NOURISHMENT (TO OR FL)
2. ABLE TO AMBULATE (AMBULAT)
3. ALERT AND ORIENTED (AL AN OR)
4. PT DENIES NAUSEA / NO VOMITING (DEN NAUS)
5. ALDRETE SCORE 8 OR GREATER (ALD 8+)
6. AIRWAY REFLEXES INTACT / SECRETIONS (AIR REFL)
7. PAIN CONTROLLED PER PT (PAIN CON)
8. PT PAIN AT ACCEPTABLE LEVEL (PT PA OK)
9. VOIDED (VOIDED)
10. ABSENT TO MINIMAL BLEEDING (ABS BLEE)
11. CMS INTACT DISTAL TO SURGICAL SITE (CMS INTA)
12. SKIN WARM AND DRY COLOR APPROPRIATE (S W A DR)
13. MAINTAINING EFFECTIVE VENTILATION (EFF VENT)
14. ALDRETE SCORE AT BASELINE (ALD BASE)
15. STABLE BP AND PULSE (Stable BP)
16. SaO2% AT ACCEPTABLE LEVEL BY ORDERS (SaO2 by Or)
17. SaO2% AT ACCEPTABLE LEVEL WITH PRESCRIBED O2 LEVEL (SaO2 w/Pre)
18. 20 MINUTES SINCE LAST IV NARCOTIC/SEDATIVE (20 Min Nar)
19. RESTING RESPIRATORY RATE BASELINE OR ABOVE (Rest Base)
20. MENTAL STATUS - AWAKE ORIENTED (M/S Awake)
21. MENTAL STATUS - BASELINE FOR KNOWN DEMENTIA (M/S Demen)
22. 60 MINUTES SINCE ADMINISTRATION OF REVERSAL AGENT (60 Min Rev)

UNITS:

-----

=====

TERM: PACU - DISCHARGE CRITERIA MET  
ID: {3B2F659B-ABA4-41E4-9861-B4B0FE3C873D}

=====

Page: 466

Abbreviation..... P DIS ME  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Discharge Criteria Met  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEA649CD} |
| 2. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |

UNITS:  
-----

Abbreviation..... PAC DISP  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... PACU - Disposition  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 12 entries

- 
1. ICU (ICU) {E7FA69FE-1088-4853-AFBD-2D815E3C2EB0}
  2. MICU (MICU) {28244B62-644C-4DA0-A350-EAFE0D2FA7C0}
  3. SICU (SICU) {EAEB2DE6-B057-4143-B44C-99DAB9B7B4AC}
  4. STEP DOWN (STEP DOW) {F6AC5CA7-0953-4E21-B1A4-6BFE1983D515}
  5. WARD (WARD) {A16CB5CB-3401-4F6C-8467-7F78FC929872}
  6. AMBULATORY CARE UNIT (ACU)  
{7FB7FBC5-42A6-4BA4-9CA6-494D77C4C549}
  7. PROCEDURE ROOM (PROC ROO) {CC9A2AB0-7A73-4DE7-9B7F-DC2C160AE462}
  8. IMAGING (IMAGING) {791A322E-F45E-40EB-B95D-E5B62E1D0A75}
  9. OPERATING ROOM (OR) {37DACE14-2A28-45F4-9DEA-35A1A433F4E6}
  10. EMERGENCY DEPARTMENT (ER DEPT)  
{D2E511F7-8B80-4928-A64B-3F1EE17395F7}
  11. WAITING ROOM (WAIT ROOM) {13476486-F895-4FFC-AB42-B274827B47DA}
  12. HOME (HOME) {C8E5E320-361B-421C-A6A9-13DC16125440}

UNITS:

-----

=====

TERM: PACU - PT RECEIVED REVERSAL AGENT  
ID: {08A85E74-4365-41AA-8C26-4D712C9364DE}

=====

Page: 468

Abbreviation..... REV AGEN  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Patient Received Reversal Agent  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEA649CD} |
| 2. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |

UNITS:  
-----



=====

TERM: PACU - TIME OF REVERSAL AGENT  
ID: {43AD957E-55D2-44BA-B591-C4B31EB64E55}

=====

Page: 469

Abbreviation..... Time of Rv  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... String  
Description.....  
Display Name..... PACU - Time of Reversal Agent  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

Abbreviation..... PAIN CHARA
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Description of pain patient is
experiencing
Display Name..... Pain Character
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 28 entries

- 1. ACHING (ACHING) {736BC1A4-2DA7-4A56-ACC0-28086EA103E0}
2. AGONIZING (AGONIZING) {F43DF9E7-555B-49FF-987C-35F083D547AA}
3. ANGINAL (ANGINAL) {ECF02415-4124-46D3-AA01-F21AC40FEFF7}
4. ANNOYING (ANNOYING) {53437530-57C9-4EA8-AC0D-36B5F2BD051D}
5. BURNING (BURNING) {D75E2D5D-9CF9-4E70-A2EB-0632EBDB82C9}
6. CONSTANT (CONSTANT) {0FB8E43C-2106-4C37-9430-38675953A85B}
7. CRAMPING (CRAMPING) {3EC3E352-5469-4A85-80E4-3AFEB61224C4}
8. CRUSHING (CRUSHING) {92D4BAF1-63E3-4FC0-B096-B5CBC4736245}
9. DEEP (DEEP) {78AD569C-16D4-CFC3-6B6F-41FB9FFAB9D7}
10. DULL (DULL) {F0F80976-F6A0-4363-9C75-AAAE9937DB49}
11. ELECTRIC (ELECTRIC) {0C221385-61AF-4660-A129-35E610EA71AB}
12. GNAWING (GNAWING) {8AC8EDB0-D7B7-4368-BE95-C35CBADD43E0}
13. HEAVY (HEAVY) {1F10CD98-03FA-4705-BF7D-941E881C4E3D}
14. INTERMITTENT (INTERMITTE) {DD09ED7C-B977-412F-A678-5130F77D239E}
15. KNIFE-LIKE (KNIFE-LIKE) {FBA0EE3F-134B-431B-879A-9B17D65CD464}
16. NUMB (NUMB) {652634C0-5334-4D56-BE57-5F58F952FF0D}
17. POSITIONAL (POSITIONAL) {9A74D16B-5E95-4EA8-88AF-E53C80311A5D}
18. PRESSURE (PRESSURE) {387AFB7F-3973-4D77-B2E6-0724F344D8F6}
19. RADIATING (RADIATING) {B6807939-F809-45C1-912A-ADEF2E0DB0DC}
20. SHARP (SHARP) {20608864-E5CF-4FE3-8B6D-AE2FEE086503}
21. SHOOTING (SHOOTING) {73C48DB1-6FE6-46B7-B1F6-E93ACBE44C27}
22. SQUEEZING (SQUEEZING) {A148587A-1B99-4D0D-B11E-87ADC3D2A811}
23. STABBING (STABBING) {684640D5-B3AB-47DA-98D8-5434ABF22B7E}
24. STINGING (STINGING) {DF235170-F501-4E8D-89DD-5A3C76F4B13B}
25. TENDER (TENDER) {30486CC7-3097-4F95-9238-FD6C2126CD10}
26. THROBBING (THROBBING) {3ACE2B10-1418-47F8-A4F3-0C41A699B5D9}
27. TINGLING (t) {AE8E3101-AC3D-DEC3-0293-58F901F82EB8}
28. UNBEARABLE (UNBEARABLE) {96611993-45F3-4D1D-BF7F-5E5C46E5569E}

UNITS:

=====

TERM: PAIN DURATION

Page: 471

ID: {DF714694-E547-4D39-9355-ADF3B419F9CD}

=====

Abbreviation..... PAIN DURAT  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Pain Duration  
Display Name..... Pain Duration  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
1. CONSTANT (CONSTANT) {0FB8E43C-2106-4C37-9430-38675953A85B}
  2. INTERMITTENT (INTERMITTE) {DD09ED7C-B977-412F-A678-5130F77D239E}

UNITS:  
-----

ID: {709DFACE-B211-4160-B291-359777ED8D0F}

=====

Abbreviation..... PAIN GOAL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Pain Goal  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 11 entries

-----

1. 0 (0)	{2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
2. 1 (1)	{8F5D10A8-FE17-47BC-DACA-6972D12CC00A}
3. 2 (2)	{38E22B79-A2B6-3401-B195-C955F24E1278}
4. 3 (3)	{E714454D-6A67-91F6-2FB7-D8F24F184328}
5. 4 (4)	{D2280C7E-DE50-5A8D-C3FD-5862B69C0E1A}
6. 5 (5)	{188EDCE7-0726-4234-D8D2-92C8052CD221}
7. 6 (6)	{640CF2F9-3BBF-AB13-185B-CD02A609FB5A}
8. 7 (7)	{7A010324-D2C2-5D4B-848A-7B8878666D20}
9. 8 (8)	{518E7CFC-102E-351E-C7AF-FB575E3C8F4D}
10. 9 (9)	{58FF2FD4-569F-F454-9F5B-FE6F9E34D8CC}
11. 10 (10)	{0C827160-18DB-816F-7FAA-A3C175C6D33B}

UNITS:

-----

Abbreviation..... PAIN LEVEL
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Standard pain scale
Display Name..... Pain Level
Help Text..... Enter 0-10, 99 or Unable to respond
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid..... 4500635

DATA VALUE: 12 entries

- 1. 0 (0) {2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
2. 1 (1) {8F5D10A8-FE17-47BC-DACA-6972D12CC00A}
3. 2 (2) {38E22B79-A2B6-3401-B195-C955F24E1278}
4. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}
5. 4 (4) {D2280C7E-DE50-5A8D-C3FD-5862B69C0E1A}
6. 5 (5) {188EDCE7-0726-4234-D8D2-92C8052CD221}
7. 6 (6) {640CF2F9-3BBF-AB13-185B-CD02A609FB5A}
8. 7 (7) {7A010324-D2C2-5D4B-848A-7B8878666D20}
9. 8 (8) {518E7CFC-102E-351E-C7AF-FB575E3C8F4D}
10. 9 (9) {58FF2FD4-569F-F454-9F5B-FE6F9E34D8CC}
11. 10 (10) {0C827160-18DB-816F-7FAA-A3C175C6D33B}
12. 99 (9) {B013C05B-1AF0-1547-4F03-378883F887C5}

LOCATION: 58 entries

- 1. HEAD (Hd) {C3209F44-90C6-FF1B-3169-265B294F276A}
2. FACE (FACE) {C638EA14-7903-6A94-496A-C6560D51BEE3}
3. NECK (Ne) {D4812C3B-AAC2-1B61-679B-3CD8B5CB947C}
4. STOMACH (STOMACH) {613C20B9-CC9C-4F5D-B5B1-69B66A663540}
5. ABDOMEN (ABDOMEN) {109FE041-2521-6EF2-7F95-1FD1587796AB}
6. BACK (BACK) {22C26239-C692-4E93-B522-E5E8392BA3B1}
7. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
8. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
9. L HAND (L HAND) {639DF25A-19E3-B504-D323-4B3F6ED39A24}
10. R HAND (R HAND) {A54D3A80-8537-D079-9CCE-7DC56346D132}
11. L HIP (L HIP) {786CD881-27AD-4136-B9CC-8423CB38A431}
12. R HIP (R HIP) {A48CBCAE-5404-4934-9321-ECEF660A5350}
13. L GROIN (L GROIN) {5488E4EE-C1A2-4B02-8763-A53C412117A7}
14. R GROIN (R GROIN) {BAE1C889-8CCF-4CEB-9EBD-0243BC6F096A}
15. L THIGH (L THIGH) {16C2312F-1E0E-45E0-B675-160F40C23F57}
16. R THIGH (R THIGH) {B7A2D4A0-801E-4578-967A-F1CFB86FB96B}
17. L KNEE (L KNEE) {3A9E30B4-09A9-967A-58EF-69A01EC4F859}
18. R KNEE (R KNEE) {BCD18D53-07B6-2C41-64FC-3A07A7E979E7}
19. L FOOT (L FOOT) {2EC3FBEC-4EA2-9E7A-17D4-FDC2F9D28046}
20. R FOOT (R FOOT) {7D9F5FF6-232C-2488-4504-A721FF06928F}
21. L LEG (LL) {B088C809-E528-B785-21C2-14E45E1FD8AF}
22. R LEG (RL) {C5E40A7B-6FC4-267D-8CE8-3571D21D81FD}
23. COCCYX (COCCYX) {49D65E9D-3D5B-40E3-B7CF-1003F80B852C}
24. GASTROINTESTINAL (GI) {C27C66DB-3A6A-487E-9981-FA7FD323C0C3}
25. L EAR (L EAR) {DA9C2DD1-D6FF-4899-A62B-1096FAA36C59}
26. R EAR (R EAR) {AE31327D-A03C-488C-8A8E-084E7E6E7FDE}
27. L SHOULDER (L SHOULDER) {91CBD2FA-2743-42EB-946B-6933980D9337}
28. R SHOULDER (R SHOULDER) {D27AE130-442E-4DD6-B391-8D7E8AD494A2}
29. L ELBOW (L ELBOW) {BF4317E7-2024-47EF-9E5F-3D21D0547885}

- 30. R ELBOW (R ELBOW) {5BAFF87B-6993-4E76-862F-9F386D9B71C5}
- 31. L ANKLE (L ANKLE) {D5A692EF-9D7C-473C-AF96-7EF1C01FA2AF}
- 32. R ANKLE (R ANKLE) {4E6DF352-F399-4CE9-908C-48132C3D856B}
- 33. PERINEAL (PERINEAL) {1801441C-D5EF-4775-9E58-E4985496F7E0}
- 34. BLADDER (BLADDER) {18E083F2-C4BB-407E-B0A3-20C4942FC97F}
- 35. L JAW (L JAW) {65EBD3CF-76D8-4F45-9B1B-3575D988BE8E}
- 36. R JAW (R JAW) {7846AEE2-B6F4-4B5A-909F-56E386FD82A8}
- 37. TOOTH (TOOTH) {8407F780-7C7F-49BF-B1DD-4FDD578CEAB3}
- 38. TEETH (TEETH) {38B36D8E-A9C3-4EE6-8727-EECDA4C6034E}
- 39. L NOSE (L NOSE) {C7707847-5135-4749-B196-002E3E59BB54}
- 40. R NOSE (R NOSE) {D80BCB19-E71F-43EA-8990-6C121EE7F7A1}
- 41. NOSE (NOSE) {85F7003E-6D4E-4A79-A534-8A57FBC2D106}
- 42. PENIS (PENIS) {2D82CE55-83B6-4160-A754-07C71FD388E8}
- 43. ABDOMEN RUQ (ABD RUQ) {F29C9FA1-717A-49B7-8000-95070E9AACA0}
- 44. ABDOMEN LUQ (ABD LUQ) {3B31E1EC-CE74-4C76-BB61-10C678431FCA}
- 45. ABDOMEN RLQ (ABD RLQ) {2C48394F-6877-4738-856C-E933644D8067}
- 46. ABDOMEN LLQ (ABD LLQ) {39831F59-D7B0-4156-AC57-4C622A93CA58}
- 47. INCISIONAL SITE (INCISON SI)  
{EB36A66D-751D-4031-9F68-0BB03E267321}
- 48. CHEST (CHEST) {8C3C30FE-6E06-C35E-8226-CBF0C94014CC}
- 49. L CHEST (L CHEST) {E332E497-FE0C-47CB-A0A7-217FBF52DC97}
- 50. R CHEST (R CHEST) {7E30EC0E-815A-4C2E-A9FC-3C2B2D6E3CD4}
- 51. L SCAPULA (L SCAPULA) {3381EC50-771B-41C8-85C8-6EB0D208D54E}
- 52. R SCAPULA (R SCAPULA) {EA657261-4D72-40C2-B70D-2A2E27556A01}
- 53. L UPPER ARM (L UPP ARM) {5F915696-1B08-32E9-6B95-49246117D660}
- 54. R UPPER ARM (R UPP ARM) {3EF20E8E-F68C-E4F2-8F8C-FAD3B44DDAA4}
- 55. SUBSTERNAL (SUBSTERN) {3DE7E0F1-BCA4-43EC-82CB-A87C6506EF46}
- 56. L LOWER ARM (L LOW ARM) {63AA95E8-256D-4390-8CB3-F00A1B3F91CD}
- 57. R LOWER ARM (R LOW ARM) {8DE6E277-EFFF-40DD-9E83-D7295AACB6AA}
- 58. SCROTUM (SCROTUM) {BDC109F5-0D30-4914-A62C-CEFEB448B547}

UNITS:

-----

Abbreviation..... PAI NONVER
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Pain Non-Verbal Cues
Display Name..... Pain Non-Verbal Cues
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 10 entries

- 1. CRYING (CRYING) {D87AB50C-0B18-40DD-884C-6AAEE4DD97FC}
2. GRIMACING (GRIMACING) {69EFDB2D-A0C4-48D4-A7E0-A862F34545C8}
3. MOANING (MOANING) {A71E3515-AA20-4B0E-9CB8-1612F76BA406}
4. RESTLESS (r) {8FCCAD26-DBA0-7067-DA72-30D342BFA8E8}
5. SPLINTING (SPLINTING) {D7FC5AA9-E63A-4A88-BE03-A50D0A8751F9}
6. GUARDING (GUARDING) {9EFD10A7-B1B7-43DA-B289-646AF5C0A25C}
7. PAINFUL EXPRESSION (PAINFUL EX) {608554F9-29F9-4B27-BE63-486AF81810BD}
8. IRRITABLE (Irritable) {10651C11-2736-4056-BA84-E544358BA877}
9. SLEEPING (s) {87338420-F351-F98C-7A80-4990B87ACE2B}
10. NO SIGNS AND SYMPTOMS (NO SIG&SYM) {91285942-6C82-4CB1-AAFB-8D3DA15A0B8A}

UNITS:

=====

TERM: PAIN RELIEVING FACTOR INTERVENTION  
ID: {325DB732-E996-4DF5-9F5F-EF3B55D80935}

=====

Page: 475

Abbreviation..... PN FACT IN  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Pain Relieving Factor Intervention  
Display Name..... Pain Relieving Factor Intervention  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 9 entries

- 
1. COLD (COLD) {B7B422BA-22EC-117D-0A0C-8419E9447D46}
  2. HEAT (HEAT) {E0FCA383-3AB8-4E03-9ABC-00DE3FE01431}
  3. QUIET ROOM (qr) {CCD96D41-C273-F6F3-DA18-8818B1CC9980}
  4. REPOSITION (REPOSITION) {DF3039CA-5977-4454-81D4-DB0AF966F85D}
  5. MEDICATION (MEDICATION) {22E6DEBB-9662-466B-990C-8DA500C0A183}
  6. REFUSE MEDICATION (REFUSE MED)  
{17D6956B-4CF6-492C-B113-B5A93D9E215E}
  7. REQUEST FOR MORE MEDICATION (REQ MOR ME)  
{0639081C-BEEC-4CAD-B420-5F2C4DBA38D0}
  8. PRE-MEDICATED FOR ACTIVITY (PRE-MED AC)  
{D6DA772A-FA22-4170-B928-855A2F68711F}
  9. PHYSICIAN NOTIFIED (PHYS NOTIF)  
{4BD68949-74F6-42A6-AEF5-E9A33DC31FFC}

UNITS:  
-----



=====

TERM: PAINAD BODY LANGUAGE  
ID: {6DF1DD0D-C735-4D34-ABF6-A274621EC177}

=====

Page: 476

Abbreviation..... PAINAD BL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... PAINAD Body Language  
Help Text..... Pain Assessment in Advanced Dementia  
(PAINAD)  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. PAINAD BL 0 (PAINAD BL0) {C77929DB-A59A-4A24-B57B-BF111770C2AB}
  2. PAINAD BL 1 (PAINAD BL1) {CC1E1516-086D-4B79-8264-0C3BD67D8552}
  3. PAINAD BL 2 (PAINAD BL2) {34C3E304-DF66-40A7-906F-22AEA3F76FAD}

UNITS:

-----

=====

TERM: PAINAD BREATHING

Page: 477

ID: {EBF3C2FD-5906-4FC8-8E75-CF998FA797C3}

=====

Abbreviation..... PAINAD BR  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... PAINAD Breathing  
Help Text..... Pain Assessment in Advanced Dementia  
(PAINAD)  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

-----

1. PAINAD BR 0 (PAINAD B0) {58E71B2E-F5C3-419F-8B71-A9FD701E0202}
2. PAINAD BR 1 (PAINAD BR1) {727D1B04-A702-4806-857B-78BCBACB0FE1}
3. PAINAD BR 2 (PAINAD BR2) {37978012-ADAC-416C-8920-48CDAAA5D59D}

UNITS:

-----

=====

TERM: PAINAD CONSOLABILITY  
ID: {869E7E2E-04C5-46A5-8811-C7F94E066A2C}

=====

Page: 478

Abbreviation..... PAINAD CO  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... PAINAD Consolability  
Help Text..... Pain Assessment in Advanced Dementia  
(PAINAD)  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. PAINAD CO 0 (PAINAD CO0) {06A89B55-5822-42F7-A689-DBD91B62979F}
  2. PAINAD CO 1 (PAINAD CO1) {8E5DDB47-1955-40AD-A008-A9CE4646A707}
  3. PAINAD CO 2 (PAINAD CO2) {1DFFB735-6272-4B60-8F42-CCCAA701C7A7}

UNITS:

-----

=====

TERM: PAINAD FACIAL EXPRESSION  
ID: {CB17C4F3-2E4B-4355-9E63-DE165CE47A25}

=====

Page: 479

Abbreviation..... PAINAD FE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Pain Assessment in Advanced Dementia  
(PAINAD)  
Display Name..... Painad Facial Expression  
Help Text..... Pain Assessment in Advanced Dementia  
(PAINAD)  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. PAINAD FE 0 (PAINAD FE0) {B3DE1BC6-E7A0-48C5-8C8E-FBE033AD1046}
  2. PAINAD FE 1 (PAINAD FE1) {66C01EE9-8F36-40B5-AED2-DBC2E94F860D}
  3. PAINAD FE 2 (PAINAD FE2) {35688D06-94F2-4B36-8940-529803D4BC43}

UNITS:  
-----

=====

TERM: PAINAD NEGATIVE VOCALIZATION  
ID: {815392EF-49E0-48D9-BCC0-39FD3333B255}

=====

Page: 480

Abbreviation..... PAINAD NV  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Pain Assessment in Advanced Dementia  
(PAINAD)  
Display Name..... PAINAD Negative Vocalization  
Help Text..... Pain Assessment in Advanced Dementia  
(PAINAD)  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. PAINAD NV 0 (PAINAD NV0) {63F7F990-2DEA-4B35-AC52-D23B8E13B1FC}
  2. PAINAD NV 1 (PAINAD NV1) {D517E6C9-9D98-4D0F-A8F4-261BED4812D1}
  3. PAINAD NV 2 (PAINAD NV2) {3C8D6B1E-53E6-4D2C-BD8A-9C548466E1D4}

UNITS:  
-----

=====

TERM: PAINAD TOTAL

Page: 481

ID: {4E35BBC8-DD3C-4EEA-A0BD-DADBD1252F28}

=====

Abbreviation..... PAINAD TO  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... PAINAD Total Score  
Help Text..... Pain Assessment in Advanced Dementia  
(PAINAD)  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
PAINAD SCALE            0 > 1 \* 10 < 11            Dec Prec: 0

=====

TERM: PATIENT ABLE TO RESPOND  
ID: {08A622B7-3939-427E-A532-7F12F8790EF7}

=====

Page: 482

Abbreviation..... PAT RESPON  
Active..... Yes  
Boolean Value False..... Unable to respond  
Boolean Value True..... Able to repond  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Patient Able To Respond  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

=====

TERM: PATIENT EDUCATION PROVIDED PRIOR TO SEDATION  
ID: {59683A04-FA78-40FA-8DDF-A728BA34C72D}

=====

Page: 483

Abbreviation..... Pat Ed Pro  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description.....  
Display Name..... Patient Education Provided Prior to  
Sedation  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----



=====

TERM: PCA - CLINICAL BOLUS DOSE  
ID: {3705805D-6ABB-412B-8DF2-648E634237F9}

=====

Page: 484

Abbreviation..... PCA BOL DO  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... PCA - Clinical Bolus Dose  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MICROGRAMS	-1.00 > 0.00 * 1000000.00 < 1001000.00 Dec Prec: 2
MILLIGRAMS	-1.00 > 0.00 * 1000.00 < 1001.00 Dec Prec: 2

=====

TERM: PCA - CONTINUOUS MCG/HR (BASAL)  
ID: {3E22F87D-E696-4401-90D0-2BE623E75DE9}

=====

Page: 485

Abbreviation..... PCA - MCG  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description.....  
Display Name..... PCA - Continuous mcg/hr (Basal)  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
MICROGRAMS PER HOUR      0.0 > 6.2 \* 100.0 < 300.0 Dec Prec: 1

=====

TERM: PCA - CONTINUOUS MG/HR (BASAL)  
ID: {39A3FFFD-8AB1-4691-A5CE-7FD268F6F5CB}

=====

Page: 486

Abbreviation..... PCA - MG  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description.....  
Display Name..... PCA - Continuous mg/hr (Basal)  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
MILLIGRAMS PER HOUR      0.0 > 0.1 \* 8.0 < 20.0    Dec Prec: 1

=====

TERM: PCA - DEMAND DOSE - DO NOT CHECK TOTAL OR SUBTOTALS  
ID: {48E49A30-E83F-4A5B-A1EA-C58D92D08288}

=====

Page: 487

Abbreviation..... PCA DEM DO  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... PCA Demand Dose - Do Not Check Total  
or Subtotals  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Void.....

UNITS:

-----

MICROGRAMS	-1.00 > 0.00 * 1000000.00 < 1001000.00 Dec Prec: 2
MILLIGRAMS	-1.00 > 0.00 * 1000.00 < 1001.00 Dec Prec: 2

=====

TERM: PCA - DOSAGE DELIVERED  
ID: {BBD517C0-EDFE-452C-9B2C-5B21400E173A}

=====

Page: 488

Abbreviation..... PCA DAG DE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... PCA - Dosage Delivered  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MICROGRAMS	-0.01 > 0.00 * 8000000.00 < 8001000.00 Dec Prec: 2
MILLIGRAMS	-1.00 > 0.00 * 8000.00 < 8001.00 Dec Prec: 2



Abbreviation..... Lev of Con  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... PCA - Level of Consciousness  
Help Text.....  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 9 entries

- 
1. ALERT (a) {BC7E2630-179A-5C80-E6C7-598D42F0BFD6}
  2. AWAKE (Awake) {D71412D4-5782-46BA-AE6C-24EDC586AE55}
  3. APPEARS ASLEEP (Asleep) {245A1E78-B5E0-4440-963C-A6ABB416D9B9}
  4. EASY TO AROUSE (Easy to Ar)  
{19BB5A04-F979-449C-AA93-D5A21BD368F6}
  5. DROWSY (d) {6B350814-A88A-0F05-C382-9F27830FBC09}
  6. SEDATED (s) {29D9ADD6-2C1F-F6D5-F4D0-B2BB86416997}
  7. STUPOROUS (Stuporous) {4C25FF58-3346-4B06-ACE5-6ADB942877A8}
  8. DIFFICULT TO AROUSE (Diff to Ar)  
{5EF640D1-CE7F-49AB-80DC-6BCFF47BC86C}
  9. UNABLE TO AROUSE (Un to Ar)  
{05B32DD9-E634-4B60-A58B-2E915EA2834D}

UNITS:  
-----







Abbreviation..... PCA LINE C  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... PCA Line Care  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 10 entries

- 
1. BAG CHANGED (BAG CHANGE) {818990DE-4F95-4F8E-92D8-B4D3B975CB36}
  2. CENTRAL VENOUS CATHETER (CENT VEN C)  
{EE91E841-C09B-4432-A121-8EBAB92E9CD2}
  3. DISCONTINUED (DISCONTINU) {15E289FB-09DF-4731-9CF3-48A4E0031C31}
  4. DRESSING CHANGED (DRESS CHAN)  
{E8CF4B7A-D1AB-49A2-8E7C-3375C3BC436F}
  5. HEPLOCK (HEPLOCK) {C827FEDE-59D6-4C45-BDDA-4884FA0882F3}
  6. PATENT CLEAN (PAT CLEAN) {68688246-AE58-2EFB-BF9D-5E00D8CFDA8E}
  7. SALINE LOCK (SALIN LOCK) {A94D7639-2A45-4319-A6C7-5A892AA6B1F4}
  8. STARTED (STARTED) {D6E105C4-5EAD-424E-8E4C-195EB7783A3D}
  9. SYRINGE CHANGED (SYR CHANGE)  
{BE84F060-7E03-49DD-A6BC-787375C8E1F1}
  10. TUBING CHANGED/LABELED (TUBE CHANG)  
{19BA337C-6E28-411E-9702-48B262C593FA}

UNITS:  
-----



=====

TERM: PCA - MODE

Page: 495

ID: {7CAE5B29-5855-4E16-85DF-4886AF9A75E6}

=====

Abbreviation..... PCA MODE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... PCA Mode  
Help Text..... Patient Controlled Analgesia Mode  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. BASAL RATE ONLY (BASAL ONL)  
    {4BFFF2BC-837C-4D62-B811-0E547247D4D0}
  2. PCA AND BASAL RATE (PCA/BASAL)  
    {D711BB51-0662-4DD0-8654-C737F6A6686A}
  3. PCA ONLY (PCA ONLY)      {30E0441A-1098-4FE9-B277-BF2BBD6B1667}

UNITS:

-----

=====

TERM: PCA - PUMP CLEARED

Page: 496

ID: {9CD3B8C5-40C0-4EA4-8325-0B4399A91F53}

=====

Abbreviation..... PCA PUMP C  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... PCA - Pump Cleared  
Help Text..... PCA Pump Usage History Cleared  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

-----

1. YES (y)	{483DE806-2EDB-34B0-FE11-C29EEEA649CD}
2. NO (n)	{FC275445-1225-9912-18A3-DF36245893D3}

UNITS:  
-----



Abbreviation..... PCA: Side  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... PCA - Side Effects  
Help Text.....  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 10 entries

- 
1. NAUSEA (NAUSEA) {6BF33DC3-0300-4643-BB14-27C35692B806}
  2. VOMITING (VOMITING) {921F463C-1639-4D13-8E4E-6D3A1B6927F9}
  3. URINARY RETENTION (Urinary Re)  
{959EBBDC-7F67-42EC-8CF0-295F709BB433}
  4. WEAKNESS (Weakness) {E4C6D569-7961-4CB4-8E90-A02A31F98094}
  5. HYPOTENSION (Hypotensio) {1AE4E3BB-3B78-42D2-BDAC-83563A0546E1}
  6. ITCHING (Itching) {3A1FB47B-F274-458D-BDF6-E6882D6CB600}
  7. TINNITUS (Tinnitus) {846F5DD2-F8AC-44BA-95C6-BF63A07E7E32}
  8. CONFUSION (Confusion) {4C1B8247-75C4-4C60-970B-3D7F07BD1164}
  9. SOMNOLENT (Somnolent) {008A3F20-DD3D-4775-AA18-CCA7B817A77D}
  10. RESPIRATORY DEPRESSION (Resp Depr)  
{A4801B08-0A98-4FDB-B075-C8D84DEAD5A0}

UNITS:  
-----

=====

TERM: PCA - TOTAL DOSE 24 HR  
ID: {C18B5B5C-B1D5-440C-964B-FD7651AF2F32}

=====

Page: 499

Abbreviation..... Tot 24 hr  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description.....  
Display Name..... PCA - Total Dose 24 hr  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MICROGRAMS	-1.00 > -1.00 * 9999.00 < 9999.00 Dec Prec: 2
MILLIGRAMS	-1.00 > -1.00 * 9999.00 < 9999.00 Dec Prec: 2



=====

TERM: PCA - TOTAL MCG PER SHIFT  
ID: {5EEA95D7-D6D8-4FA4-AD1F-C40303A145A5}

=====

Page: 500

Abbreviation..... PCA mcg sh  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Micrograms per Shift  
Display Name..... PCA - Total mcg Per Shift  
Help Text..... Micrograms per Shift  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MICROGRAMS PER SHIFT      0.0 > 0.0 \* 800.0 < 5000.0  
Dec Prec: 1

=====

TERM: PCA - TOTAL MG PER SHIFT  
ID: {58CCC88B-404B-4D47-A131-472950622B2B}

=====

Page: 501

Abbreviation..... PCA mg shf  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Milligrams per Shift  
Display Name..... PCA - Total mg Per Shift  
Help Text..... Milligrams per Shift  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
MILLIGRAMS PER SHIFT      0.0 > 0.0 \* 64.0 < 300.0    Dec Prec: 1





Abbreviation..... PK INSP PR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Peak Inspiratory Pressure - Do Not
Check Total or Subtotals
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

CENTIMETERS H2O 0 > 5 \* 100 < 200 Dec Prec: 0

=====

TERM: PLAN OF CARE - REVISION INDICATED  
ID: {DE8BBB70-21B4-4955-BEC4-C1E71BBA3E31}

=====

Page: 505

Abbreviation..... POC REV I  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Plan Of Care Revision Indicated  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

=====

TERM: PLANTAR FLEXION  
ID: {4BBDBB76-053B-CA20-908E-017F5956DF5E}

=====

Page: 506

Abbreviation..... PLANFL  
Active..... Yes  
Boolean Value False..... Absent  
Boolean Value True..... Present  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Plantar Flexion  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

=====

TERM: PR INTERVAL

Page: 507

ID: {3657545A-7138-94C0-03A5-3136C51B0369}

=====

Abbreviation..... PR INTERVA  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... The PR interval is measured from the  
beginning of the P wave to the  
beginning of the QRS complex  
Display Name..... PRI  
Help Text..... The PR interval is measured from the  
beginning of the P wave to the  
beginning of the QRS complex  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
SECOND                    0.00 > 0.12 \* 0.20 < 1.00 Dec Prec: 2



=====

TERM: PREPROCEDURE: 12 LEAD EKG DONE?  
ID: {98DB7D48-ED01-4F73-BD71-82C62F0EBBD2}

=====

Page: 508

Abbreviation..... Pre-Prc 12  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description.....  
Display Name..... Preprocedure: 12 lead EKG done?  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

Abbreviation..... Pre-Allerg  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... Preprocedure: Allergies Reassessed  
and Documented  
Help Text.....  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

-----

1. PATIENT STATED (Pat Stated)  
{843654E1-4BC9-4277-9EFC-0B388A70A907}
2. FAMILY STATED (Fam Stated)  
{07A4B749-65F2-48E8-AD93-D142B04055B4}
3. CHART VERIFIED CORRECT (Chart Ver)  
{319E702F-C7A9-409D-881F-3F4D1D9AF393}
4. ALLERGY WRISTBAND (Wristband)  
{BB7F6746-102D-4CBA-8E53-EA465176812B}

UNITS:

-----

=====

TERM: PREPROCEDURE: ANTIBIOTIC GIVEN?  
ID: {07D521A2-01EA-404F-96C9-40EBE772E514}

=====

Page: 510

Abbreviation..... Pre-Prc Gv  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description.....  
Display Name..... Preprocedure: Antibiotic given?  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: PREPROCEDURE: ANTIBIOTIC ORDERED?  
ID: {5FBF230D-E096-4980-B96D-B7BB264AE733}

=====

Page: 511

Abbreviation..... Pre-Prc An  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description.....  
Display Name..... Preprocedure: Antibiotic ordered?  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----



=====

TERM: PREPROCEDURE: CHEST XRAY DONE?  
ID: {37CC251D-A494-4E4C-AECA-D90C4AEAEF4F}

=====

Page: 513

Abbreviation..... Pre-Prc Xr  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description.....  
Display Name..... Preprocedure: Chest Xray done?  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: PREPROCEDURE: DENTURES REMOVED?  
ID: {BA505699-D88A-41FB-AAC7-DCBE13D35D0C}

=====

Page: 514

Abbreviation..... Pre-Dentur  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description.....  
Display Name..... Preprocedure: Dentures removed?  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:  
-----

Abbreviation..... Clothing  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... Preprocedure: Disposition of Clothing  
Help Text.....  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 7 entries

- 
1. CLOTHING ROOM (Cloth Room)  
{8F116C60-10B1-4033-93F2-9121629CD608}
  2. LOCKER (Locker) {36AE4F3F-E2BC-4DD3-8DD0-95BD0E38B1B9}
  3. SPOUSE (Spouse) {2C601365-23F9-4939-9B63-F73720A17199}
  4. FAMILY (FAMILY) {1233C4BB-5AC5-4883-AF2A-F5C93CD8814F}
  5. FRIEND (FRIEND) {03902E48-AF86-432A-8372-37796A55F625}
  6. UNDER STRETCHER/BED (Stretcher)  
{0B945D0B-FE87-46A5-A2A9-A3FC8CC578AC}
  7. N/A (n) {B3DBE3FD-C5CB-3894-607A-BE13319594D0}

UNITS:  
-----



Abbreviation..... Pre-Val  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... Preprocedure: Disposition of Valuables  
Help Text.....  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 8 entries

- 
1. CLOTHING ROOM (Cloth Room)  
{8F116C60-10B1-4033-93F2-9121629CD608}
  2. AGENT CASHIER (Cashier) {3849A887-F353-4215-88FD-DDD0188BE1BB}
  3. LOCKER (Locker) {36AE4F3F-E2BC-4DD3-8DD0-95BD0E38B1B9}
  4. SPOUSE (Spouse) {2C601365-23F9-4939-9B63-F73720A17199}
  5. FAMILY (FAMILY) {1233C4BB-5AC5-4883-AF2A-F5C93CD8814F}
  6. FRIEND (FRIEND) {03902E48-AF86-432A-8372-37796A55F625}
  7. LOCK BOX (Lock Box) {974A2F44-C918-4E85-ABAF-03A21498F98B}
  8. N/A (n) {B3DBE3FD-C5CB-3894-607A-BE13319594D0}

UNITS:  
-----

Abbreviation..... Pre-driver
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description.....
Display Name..... Preprocedure: Driver/Escort to Home
Confirmed
Help Text.....
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

- 1. SPOUSE (Spouse) {2C601365-23F9-4939-9B63-F73720A17199}
2. FRIEND (FRIEND) {03902E48-AF86-432A-8372-37796A55F625}
3. RELATIVE (Relative) {F4EEB0A5-D25A-45DC-A5AA-BBB7ABA08750}
4. VOLUNTEER (VOLUNTEE) {F5406412-DE48-4F84-8A31-0B1D79DB576D}

UNITS:

=====

TERM: PREPROCEDURE: FINGERSTICK GLUCOSE DONE?  
ID: {98B3541C-372E-4876-ABD2-74E04AF23DF0}

=====

Page: 518

Abbreviation..... Pre-Finger  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description.....  
Display Name..... Preprocedure: Fingerstick Glucose done?  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:  
-----

=====

TERM: PREPROCEDURE: GLASSES / CONTACTS REMOVED?  
ID: {FBE31829-FEFC-4C11-8D62-EDCF94EA28EF}

=====

Page: 519

Abbreviation..... Pre-Glass  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description.....  
Display Name..... Preprocedure: Glasses / Contacts  
removed?  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:  
-----

=====

TERM: PREPROCEDURE: H AND P CURRENT?  
ID: {CF191D5F-A3FF-4169-8A47-1C799BA22C48}

=====

Page: 520

Abbreviation..... Pre-Prc HP  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description.....  
Display Name..... Preprocedure: H and P current?  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: PREPROCEDURE: HEARING AIDS IN PLACE?  
ID: {173DAC5E-4332-4D0B-A232-98DA65305234}

=====

Page: 521

Abbreviation..... Pre-Hear  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description.....  
Display Name..... Preprocedure: Hearing Aids in place?  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:  
-----



=====

TERM: PREPROCEDURE: HOSPITAL GOWN ON?  
ID: {62B89BA3-009B-4A2A-A6E4-E82B527FBC40}

=====

Page: 523

Abbreviation..... Pre-Gown  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description.....  
Display Name..... Preprocedure: Hospital Gown on?  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----



=====

TERM: PREPROCEDURE: INFORMED CONSENT CURRENT  
ID: {068445ED-622A-4549-842D-0D88922EF8CE}

=====

Page: 524

Abbreviation..... Pre-Prc In  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description.....  
Display Name..... Preprocedure: Informed Consent current?  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: PREPROCEDURE: JEWELRY REMOVED?  
ID: {8220E017-8816-452A-B2E1-07209A167781}

=====

Page: 525

Abbreviation..... Pre-Jewel  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description.....  
Display Name..... Preprocedure: Jewelry Removed?  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: PREPROCEDURE: LABS DONE?  
ID: {88D2EB1F-C5E5-4190-A697-17C1E078A12B}

=====

Page: 526

Abbreviation..... Pre-Prc La  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description.....  
Display Name..... Preprocedure: Labs Done?  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: PREPROCEDURE: NPO SINCE (DATE/TIME)  
ID: {64F57F8D-DDB7-4CB7-9C71-F9C6C334F717}

=====

Page: 527

Abbreviation..... Pre-NPO  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... String  
Description.....  
Display Name..... Preprocedure: NPO Since (Date/Time)  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: PREPROCEDURE: PLANNED PROCEDURE  
ID: {696F0ABC-ADA5-4BDF-8EE6-CF9FCA84581C}

=====

Page: 528

Abbreviation..... Pre-Prc PI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... String  
Description.....  
Display Name..... Preprocedure: Planned Procedure  
Help Text..... Enter name of procedure to be performed  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:  
-----



=====

TERM: PREPROCEDURE: PREMEDICATION ORDERED?  
ID: {13AD97AD-C093-4341-BC4E-834A5B016727}

=====

Page: 530

Abbreviation..... Pre-Meds  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description.....  
Display Name..... Preprocedure: Premedication ordered?  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----







=====

TERM: PREPROCEDURE: TYPE AND CROSS MATCH DONE?  
ID: {6CC89C30-4F13-4ECA-9D80-2A80BEB13120}

=====

Page: 533

Abbreviation..... Pre-Prc Cr  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description.....  
Display Name..... Preprocedure: Type and Cross Match  
done?  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:  
-----

Abbreviation..... Pre-Prc TS  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description.....  
Display Name..... Preprocedure: Type and Screen done?  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:  
-----

=====

TERM: PREPROCEDURE: VITALS SIGNS DONE?  
ID: {749E6863-8E6B-43F2-A521-84959FF0260B}

=====

Page: 535

Abbreviation..... Pre-Vitals  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description.....  
Display Name..... Preprocedure: Vitals Signs done?  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:  
-----

Abbreviation..... Pre-Voided  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description.....  
Display Name..... Preprocedure: Voided?  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: PRESSURE ULCER COLOR  
ID: {CC98B260-63A7-43DF-956C-AFE2EE8974D4}

=====

Page: 537

Abbreviation..... PR UL COLO  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Pressure Ulcer Color  
Display Name..... Pressure Ulcer Color  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 6 entries

- 
- 2. BLACK/NECROTIC (BLACK/NECR)  
{53519F9A-B844-41C7-A013-F9D2EA8E4BC0}
  - 3. YELLOW/FIBRINOUS (YELLOW/FIB)  
{47CF52FE-0374-4690-ADA1-643E0954E035}
  - 40. RED (RED) {515D7CFB-3527-30E7-EC4B-B607C272DFA0}
  - 50. PINK / GRANULATING (Pink/Gran)  
{05245B27-3292-48B2-9F37-6AB9BB92395B}
  - 60. SLOUGH (SLOUGH) {548EB536-1AC4-4234-8C4D-DC54E4CAA6F6}
  - 70. OTHER (Other) {4CBA16B3-6BBA-4B8D-AA12-D721AECFE286}

UNITS:  
-----

=====

TERM: PRESSURE ULCER DEPTH  
ID: {A829C7B5-B0E5-4B7B-882A-BA63337C65F8}

=====

Page: 538

Abbreviation..... PR UL DEPT  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Pressure Ulcer Depth  
Display Name..... Pressure Ulcer Depth  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

CENTIMETER            -1.00 > 0.00 \* 20.00 < 40.00  
                              Dec Prec: 2

Abbreviation..... PR UL DRAI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Pressure Ulcer Drainage Character  
Display Name..... Pressure Ulcer Drainage  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 8 entries

- 
1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
  2. BLOODY (BLOODY) {D7904E84-8331-D6DD-F2C0-7F3039C5746A}
  3. SEROUS (SEROUS) {9BDA347E-29E9-80C7-91FA-016192D4257E}
  4. SEROSANGUINOUS (SEROSANGUI)  
{19A8D332-899D-177A-2247-ABCDAAEE3BFA8}
  5. PURULENT (PURULENT) {86C01494-743A-C6A2-8883-01CA5200C79B}
  6. CLEAR (c) {8981C365-2B46-704F-C9F4-975DF1C23C2D}
  7. CLOUDY (CLOUDY) {3CCFCCCC-00FC-6B00-1D58-A77C02FDB401}
  8. VISCOUS (VISCOUS) {0E8019D5-F2DB-46B8-B750-41568DD0179C}

UNITS:

-----







=====

TERM: PRESSURE ULCER ODOR  
ID: {48FE49CA-D1FE-4812-ACB2-D4540B0EC2B9}

=====

Page: 542

Abbreviation..... PR UL ODOR  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Pressure Ulcer Odor  
Display Name..... Pressure Ulcer Odor  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
  2. FOUL SMELLING (FOUL SMELL)  
{761346A9-329D-4F42-8D80-0E9D8D40653C}

UNITS:  
-----

=====

TERM: PRESSURE ULCER PHOTOGRAPHED  
ID: {453C53C4-516D-4FE6-9032-496517880D3F}

=====

Page: 543

Abbreviation..... PR UL PHOT  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Was the pressure ulcer photographed?  
Display Name..... Pressure Ulcer Photographed  
Help Text..... Was the pressure ulcer photographed?  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

Abbreviation..... PR ULC BOO
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Pressure Ulcer Prevention Boots
Display Name..... Pressure Ulcer Prevention Devices
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 9 entries

- 1. ELBOW PAD (ELBOW PAD) {26942F81-1E33-4962-B772-A96FFBCED814}
2. FOAM AIR HEEL ELEVATOR (FO AI HEEL) {4EC26AA3-4577-4808-8728-F360C5C5A0D4}
3. HEEL LIFT (HEEL LIFT) {974F90E6-382D-421A-A4CE-94A74CBFC32}
5. INFLATABLE HEEL LIFT (INF HEEL L) {793B34A4-2F21-42DA-881D-1E6E0B4EAD6B}
6. RIGID HEEL SPLINT (RIG HEEL S) {5A231387-9FAC-4B2D-94D6-5DF4312E9AF8}
7. MULTIPODOUS (MULTIPODOU) {ED098F9D-B204-4310-813A-98D91FF9DD25}
8. WAFFLE (WAFFLE) {29A0FA73-CAA7-4660-9861-79A09F6DD98C}
90. SPECIALTY MATTRESS/BED (Specialty) {6B84BE9C-9AAC-430B-8F0D-0001E2F2AAC7}
100. OTHER (ADD COMMENTS) (Other) {264AB0B9-64B6-4EC4-BBAC-8A879CCDB23B}

UNITS:
-----

Abbreviation..... PRE ULC RF
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Pressure Ulcer Risk Factors
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 9 entries

- 1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
2. BARIATRIC PATIENT (BARIAC PAT) {6F614334-685E-4533-A475-BDF5440855EA}
3. DEVICE RELATED PRESSURE (DEV REL PR) {3AED6107-A173-4B53-8EB5-2FB2B3336FC5}
4. END OF LIFE CARE (END LIFE) {047D8977-0567-4B90-AD6A-3AAF8DB35238}
5. HYPOALBUMINEMIA (HYPOALBUM) {A3397F38-2491-48E7-8366-A6FA2D91E866}
6. VASOPRESSOR (VASOPRESS) {591FCB09-1287-48C6-8007-4CE961296BA9}
7. REFUSE TO MOVE (REF TO MOV) {E14588CE-2498-4FE8-AB0B-E954DC1696AC}
8. UNSTABLE (UNSTABLE) {5841C521-8070-4034-A69D-27DAB2990EB3}
9. VERY LOW BMI (LOW BMI) {78B3A37B-938B-42E5-A403-75E206722BC3}

UNITS:
-----

=====

TERM: PRESSURE ULCER SINUS/UNDERMINING  
ID: {CAE9260D-3B9E-433E-8FCD-F5D98A0EB02D}

=====

Page: 546

Abbreviation..... PR UL SINU  
Active..... Yes  
Boolean Value False..... Not Present  
Boolean Value True..... Present  
Data Type..... Boolean  
Description..... Pressure Ulcer Sinus/Undermining  
Display Name..... Pressure Ulcer Sinus/Undermining  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

Abbreviation..... PR UL STAG  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Pressure Ulcer Stage  
Display Name..... Pressure Ulcer Stage  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 12 entries

- 
1. NONBLANCHABLE ERYTHEMA (N BLAN ERY)  
{30273ACE-03F3-4818-995B-26586AAB5499}
  2. REDDENING OF THE SKIN (RED O SKIN)  
{716AEC57-0835-4D07-A78C-D8453599A3D2}
  3. NO BREAK IN EPIDERMIS (NO BREAK E)  
{615D96CD-6096-446B-BDDF-3F7242936AF9}
  4. PART THICKNESS SKIN LOSS INVOLVE EPIDERMIS/DERMIS (PART EPIDE)  
{EAD3A260-CC46-4B47-AA49-91CFE13EEE4B}
  5. FULL THICK SKIN LOSS INVOLVING SUBCUTANEOUS TISSUE (FULL SKINL)  
{14C31B14-1D35-4CCB-81A3-C6E30801950D}
  6. FULL THICK SKIN LOSS INVOLVING MUSCLE, BONE, SUPPO (F SK L MBS)  
{E027FBEF-1C83-41B7-B848-26D9D62F4266}
  7. SUSPECTED DEEP TISSUE INJURY (SUS DEEP T)  
{E82B4082-646F-4B5A-B68E-BB0872A4F6D2}
  8. STAGE I (STAGE I) {B376AAAF-7143-4422-BCA2-CACF1201817C}
  9. STAGE II (STAGE II) {B520BEA3-841E-490A-A3CD-63D1027B561F}
  10. STAGE III (STAGE III) {A2F0866C-ECEE-4FCD-A34B-498232662D7E}
  11. STAGE IV (STAGE IV) {284334AA-E488-4D45-93EF-9B93B2086A63}
  12. UNSTAGEABLE (UNSTAGEABL) {5C507F64-2EC0-4FD8-A786-94D76D90C0D1}

LOCATION: 51 entries

- 
1. COCCYX (COCCYX) {49D65E9D-3D5B-40E3-B7CF-1003F80B852C}
  2. HEAD OCCIPITAL (HEAD OCCIP)  
{2066CA47-5A16-447F-9919-24ED435B627F}
  3. L ANKLE (L ANKLE) {D5A692EF-9D7C-473C-AF96-7EF1C01FA2AF}
  4. L ANTECUBITAL (L ANTECUBI)  
{C21838CC-BAB5-1D7F-E2A4-A7C663FB28C3}
  5. L EAR (L EAR) {DA9C2DD1-D6FF-4899-A62B-1096FAA36C59}
  6. L ELBOW (L ELBOW) {BF4317E7-2024-47EF-9E5F-3D21D0547885}
  7. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
  8. L FOOT (L FOOT) {2EC3FBEC-4EA2-9E7A-17D4-FDC2F9D28046}
  9. L GLUTEUS (L GLUTEUS) {88C4E547-04CD-4D52-9416-CAEB6D898D21}
  10. L HEAD (L HEAD) {3737DC2D-DC9E-4FC8-839D-0355987B62D1}
  11. L HEEL (L HEEL) {C9154881-17C7-4921-A0EE-BFB962C59CEF}
  12. L INNER KNEE (L INN KNEE) {79853581-CCF0-49E9-9D91-D16DE02A836D}
  13. L ISCHIAL (L ISCHIAL) {DCDAD741-9599-4479-9D55-3279F800DBE6}
  14. L LOWER ARM (L LOW ARM) {63AA95E8-256D-4390-8CB3-F00A1B3F91CD}
  15. L OUTER KNEE (L OUT KNEE) {B143C935-2687-49C3-9F65-8D69C1775ABC}
  16. L RIBS (L RIBS) {BDFEE74F-B3CC-4AD9-AE49-78B0BF19A219}
  17. L SCAPULA (L SCAPULA) {3381EC50-771B-41C8-85C8-6EB0D208D54E}
  18. L SHOULDER (L SHOULDER) {91CBD2FA-2743-42EB-946B-6933980D9337}
  19. L TROCHANTER (L TROCHANT) {C91647A4-255D-40E4-A24D-01274E5D799F}
  20. L UPPER ARM (L UPP ARM) {5F915696-1B08-32E9-6B95-49246117D660}



- 21. PERINEAL (PERINEAL) {1801441C-D5EF-4775-9E58-E4985496F7E0}
- 22. R ANKLE (R ANKLE) {4E6DF352-F399-4CE9-908C-48132C3D856B}
- 23. R ANTECUBITAL (R ANTECUBI)  
{7FB00A65-B4B2-06A6-BD89-98AC112D8D4E}
- 24. R EAR (R EAR) {AE31327D-A03C-488C-8A8E-084E7E6E7FDE}
- 25. R ELBOW (R ELBOW) {5BAFF87B-6993-4E76-862F-9F386D9B71C5}
- 26. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
- 27. R FOOT (R FOOT) {7D9F5FF6-232C-2488-4504-A721FF06928F}
- 28. R GLUTEUS (R GLUTEUS) {B4D32045-68D1-4457-9106-7937887FA6F6}
- 29. R HEAD (R HEAD) {42A95AAE-3DA9-4C96-A4E7-A15C9840217A}
- 30. R HEEL (R HEEL) {0020DD20-2EA4-4362-ADF0-07BFA220C60E}
- 31. R INNER KNEE (R INN KNEE) {2AE24BF1-EA1B-490D-9DB7-EEFC9314C560}
- 32. R ISCHIAL (R ISCHIAL) {97E50F6C-CF75-4E59-AFDE-80FC5CE6D9E3}
- 33. R LOWER ARM (R LOW ARM) {8DE6E277-EFFF-40DD-9E83-D7295AACB6AA}
- 34. R OUTER KNEE (R OUT KNEE) {330F1DF1-8C99-4984-9965-7404ABC9562D}
- 35. R RIBS (R RIBS) {260FD564-5E58-419B-9F73-62B6037D1572}
- 36. R SCAPULA (R SCAPULA) {EA657261-4D72-40C2-B70D-2A2E27556A01}
- 37. R SHOULDER (R SHOULDER) {D27AE130-442E-4DD6-B391-8D7E8AD494A2}
- 38. R TROCHANTER (R TROCHANT) {1A3E21B2-C391-40D2-8059-12290E667D9E}
- 39. R UPPER ARM (R UPP ARM) {3EF20E8E-F68C-E4F2-8F8C-FAD3B44DDAA4}
- 40. SACRAL (SACRAL) {0CF8A8A7-AB35-B803-4544-0C91645D9D1C}
- 41. FOREHEAD (FOREHEAD) {5F892722-956E-474C-BA08-4074142D49F6}
- 42. L LATERAL HEEL (L LAT HEEL)  
{1A10FB20-7BAD-4D2B-9CBD-1A1A01DDCBF1}
- 43. R LATERAL HEEL (R LAT HEEL)  
{A1896CF9-0246-4323-8DED-E4B5FB12CACE}
- 44. L MEDIAL HEEL (L MED HEEL)  
{47AD7BA6-4EBB-4A08-B03C-FE73EB240725}
- 45. R MEDIAL HEEL (R MED HEEL)  
{C0526527-0AA9-4ECF-B768-99B2A1F5F3BF}
- 46. L LATERAL MALLEOLUS (L LAT MALL)  
{9A523B7C-64DD-4E47-9543-165788108D4C}
- 47. R LATERAL MALLEOLUS (R LAT MALL)  
{950FC45B-93E8-424F-A5D4-DE646166D7B0}
- 48. L MEDIAL MALLEOLUS (L MED MALL)  
{330E7682-407F-43CD-84FC-1CDA41F08C68}
- 49. R MEDIAL MALLEOLUS (R MED MALL)  
{C70FB4FA-6516-4567-8E26-46BA439F1DF4}
- 50. SCROTUM (SCROTUM) {BDC109F5-0D30-4914-A62C-CEFEB448B547}
- 51. TRACHEOSTOMY (TRACHEOSTO) {4D12329A-56C8-4175-986F-2F0B68EB7E6D}

UNITS:

-----

Abbreviation..... PR UL TREA
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Methods used for treatment of
Pressure Ulcer
Display Name..... Pressure Ulcer Treatment
Help Text..... Methods used for treatment of
Pressure Ulcer
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 24 entries

- 1. ANTIBIOTIC OINTMENT (AB OINTMEN) {352D9690-77BB-44DC-919E-5357876DEAC4}
2. BACTERICIDAL (BACTERIOCI) {680E4A3C-A40A-461E-818D-003AE8A1EFC8}
3. BACTERIOSTATIC (BACTERIOST) {AB74D109-8ED7-4A6C-B028-7A041781DF67}
4. CADEXOMER IODINE (CADEXO IOD) {C83E41AA-6A1F-43CC-AB2E-7F6555FEC9D2}
5. CALCIUM ALGINATE (CA ALGINAT) {2B2ACD25-A77B-4D39-8598-B26FD9BAFFC2}
6. COLLAGEN (COLLAGEN) {8D734495-2FAF-41C0-9146-55B22583A830}
7. DILUTE DAKINS SOLUTION (DL DAKINS) {66DEF4D0-54C0-4ADF-8FFD-F4E022C30F43}
8. ENZYMATIC OINTMENT (ENZY OINTM) {E781DA75-8777-4884-A093-FB623774DB8F}
9. GAUZE (GAUZE) {CF10C658-4334-4F3B-A38E-E98894458333}
10. GAUZE PACKING (GAUZE PACK) {6C093717-82CB-4966-929E-C909310649C1}
11. GAUZE WRAP (GAUZE WRAP) {040A9648-E6E1-4AC5-93C7-53A49537B15C}
12. HYDROCOLLOIDAL DRESSING (HCOLL DRES) {D85741EF-0DB3-48CF-BC9A-88C834A72F36}
13. HYDROFIBER (H-FIBER) {42740A7F-A831-4FBD-9E59-88F315176674}
14. HYDROFIBER WITH SILVER (HYFIB W SI) {04BD9512-DE5D-421A-9818-8A33D1F0693D}
15. HYDROGEL DRESSING (H-GEL DRES) {F6192D10-D6EA-42C2-AB15-270659D9925C}
16. OPEN TO AIR (OPE TO AIR) {EFB715D0-5666-4281-9034-A0326A14B2BB}
17. NEGATIVE PRESSURE THERAPY DEVICE (NEG PRE TH) {8AF0A069-66AD-4114-99D8-238AD987F0FF}
18. SALINE MOIST DRESSING (SAL MOI DR) {19F98514-8ABF-461C-9969-C31177FE7DE7}
19. TRANSPARENT DRESSING (TRANSP DRE) {777C7722-AEA8-4A0B-9ADD-1157D37AE539}
20. PETROLATUM (PETROLATUM) {DA297C75-4776-4C8D-B0D7-A0FF0E4D1B16}
21. PORCINE SUBMUCOSA (PORC SUBMU) {6F3917DA-74BD-4790-9D4D-B284E27B7623}
22. SILVADENE (SILVADENE) {0889134F-C419-4366-B12E-EAB775719410}
23. XENADERM (XENADERM) {F181D70E-A231-4D2B-A1E6-DDF4DBED31CD}
24. XEROFORM (XEROFORM) {7D7C3CDC-253E-4262-B8D4-3B70DE9AEBE4}

UNITS:

=====

TERM: PRESSURE ULCER TUNNEL  
ID: {584775CB-EB35-47BB-A7D4-62966A23D0DF}

=====

Page: 549

Abbreviation..... PRE ULC TU  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Pressure Ulcer Tunnelling  
Help Text..... Is the Pressure Ulcer Tunnelling?  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

=====

TERM: PRESSURE ULCER TUNNEL LENGTH  
ID: {D7FC23D9-8C2B-42F0-A545-5102AFC22335}

=====

Page: 550

Abbreviation..... PR UL TU L  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Pressure Ulcer Tunnel Length  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
CENTIMETER            -1.0 > 0.0 \* 15.9 < 16.0    Dec Prec: 1

=====

TERM: PRESSURE ULCER UNDERMINING  
ID: {CF9D2160-4E70-452C-B0D1-D14A4E5662CE}

=====

Page: 551

Abbreviation..... PRE ULC UM  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Pressure Ulcer Undermining  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----



Abbreviation..... Proc Loc  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... Procedure Location  
Help Text..... Document physical environment in  
which the procedure is occurring.  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 12 entries

- 
1. PULMONARY PROCEDURE ROOM (Pul Proc Rm)  
{D082F661-EBA4-4EAA-987C-67E2BDFBE94D}
  2. RADIOLOGY PROCEDURE ROOM (Radio Rm)  
{744142F5-D301-4509-A96D-06AB8C1D2C81}
  3. SPECIAL PROCEDURE ROOM (Special Rm)  
{60AA3B10-B3EE-42EF-9CEC-0CAE66DFE164}
  4. EMERGENCY DEPT (Emer Dept)  
{E8490DF5-278C-4102-9675-C9B40D596A10}
  5. GI PROCEDURE ROOM (GI Proc Rm)  
{924532F4-B52F-4EA9-8C7D-59168B340118}
  6. GU PROCEDURE ROOM (GU Proc Rm)  
{D4B6678F-C1DD-4A23-9045-574251F6CD74}
  7. DENTAL (Dental) {0A60F06A-ED63-4B06-8FF9-EEFB12567D65}
  8. RADIOLOGY (Radiology) {AF29AA13-EAA2-4B72-AA9E-994B5CCC4705}
  9. INTERVENTIONAL RADIOLOGY (Interv Rad)  
{322B66E6-814F-43AF-B669-E503519DB72F}
  10. ICU (ICU) {E7FA69FE-1088-4853-AFBD-2D815E3C2EB0}
  11. BEDSIDE (Bedside) {481BBBB3-A1F0-4F49-96DA-06BF3884E9BD}
  12. OTHER (Other) {4CBA16B3-6BBA-4B8D-AA12-D721AECFE286}

UNITS:

-----

Abbreviation..... Proc Pers
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description.....
Display Name..... Procedure Personnel
Help Text..... Identify personnel present
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 17 entries

- 10. PHYSICIAN (PHYSICIA) {DB0C5AE0-459A-4E04-A202-D1CF54B84A6A}
20. NP (NP) {E7EA4017-4374-41A1-9D3D-54051CEDF63C}
30. PA (PA) {C057C47D-780A-4230-B577-2C528AD09DEB}
40. DENTIST (Dentist) {DE4A47EB-35D3-4992-B8DA-3E5D51E0EA3C}
50. CNS (CNS) {487A4DE0-A822-437D-87A2-75FA86D7BBDE}
60. ANESTHESIOLOGIST (Anesthesio) {44BAD970-DD8F-4015-9C41-AC7868CEBE79}
70. CRNA (CRNA) {41D71B84-D87B-4219-8E1F-AF0FF5CF5EDF}
80. LVN/LPN (LVN/LPN) {DC89CD99-80C5-4455-A280-3A19E74DCE7D}
90. HEALTH TECH (Health Tec) {7C5F34C7-DBE6-4FCB-B46D-F802FBBE4CEF}
100. RESP TECH (Resp Tech) {8A7BA130-3EDC-4D07-A852-846B37399FC9}
110. RAD TECH (Rad Tech) {EEA2A91A-DD5A-4062-9DCD-6A909220ABD9}
120. STUDENT (STUDENT) {613AFF67-9121-4F6C-A44F-1ACFEF10E612}
130. RN (RN) {A103E506-EAA5-4777-8126-268CD723AD87}
140. NSG ASSIST (Nsg Assist) {14FA9638-3786-4FEA-BE8F-2345C1F455C3}
150. VENDOR REP (Vendor Rep) {4D90C5A3-7395-4BDF-9B5B-BC85F1BE28F7}
160. TECH (Tech) {13AAF5D5-FEBC-4E09-B972-5D86BA02F307}
170. OTHER (SEE COMMENTS) (Other) {7908D30B-8BD4-4A02-964D-8E6DBCA24A81}

UNITS:



=====

TERM: PROCEDURE TIME  
ID: {709E6959-A58C-454C-A2A4-2AF69483F21B}

=====

Page: 555

Abbreviation..... Proc Time  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... Procedure Time  
Help Text..... When entering data, make sure the  
save date/time reflects the start or  
stop time.  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Void.....

DATA VALUE: 2 entries

- 
1. START TIME: (Start Time) {4D0190FA-FFE4-4D5F-8A15-0517E8AB9B14}
  2. STOP TIME: (Stop Time:) {784CA4D8-988F-4A59-8AB5-48AFC9CA210B}

UNITS:  
-----

=====

TERM: PROCEDURE, PLANNED  
ID: {87C9E441-B4EA-4F72-9CDF-B9DCC072F86A}

=====

Page: 556

Abbreviation..... Proc, Plan  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... String  
Description.....  
Display Name..... Procedure, Planned  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: PULMONARY ARTERY DIASTOLIC  
ID: {A7126912-49DD-A952-8731-D06A323D2C3E}

=====

Page: 557

Abbreviation..... PAD  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Pulmonary Artery Diastolic Pressure  
Display Name..... Pulmonary Artery Diastolic Pressure  
Help Text..... Pulmonary Artery Diastolic Pressure  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MILLIMETERS OF MERCURY      0 > 8 \* 15 < 30      Dec Prec: 0

=====

TERM: PULMONARY ARTERY MEAN PRESSURE  
ID: {BE486939-CA23-58EF-ADE2-A34D19C4F41F}

=====

Page: 558

Abbreviation..... PAM  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Pulmonary Artery Mean Pressure  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----  
1. CARDIAC OUTPUT MONITOR (COM)  
{CCBE5744-7CB9-9DFA-938F-FDDBC73B45B3}

UNITS:

-----  
MILLIMETERS OF MERCURY 0 > 9 \* 19 < 50 Dec Prec: 0

=====

TERM: PULMONARY ARTERY SYSTOLIC  
ID: {4DFDF10B-3501-87E1-0EB8-AD3772E60012}

=====

Page: 559

Abbreviation..... PAS  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Pulmonary Artery Systolic  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MILLIMETERS OF MERCURY      0 > 15 \* 25 < 50      Dec Prec: 0

=====

TERM: PULMONARY CAPILLARY WEDGE PRESSURE  
ID: {9D3EB2BE-9492-CDE3-03EC-07FE580E22E0}

=====

Page: 560

Abbreviation..... PCWP WEDGE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Pulmonary Capillary Wedge Pressure  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----  
1. CARDIAC OUTPUT MONITOR (COM)  
{CCBE5744-7CB9-9DFA-938F-FDDBC73B45B3}

UNITS:

-----  
MILLIMETERS OF MERCURY 0 > 8 \* 12 < 30 Dec Prec: 0

=====

TERM: PULMONARY CARE

Page: 561

ID: {A7906E15-CBF1-A881-F3C7-4DF41A6A85B6}

=====

Abbreviation..... PUL CARE  
Active..... Yes  
Boolean Value False..... False  
Boolean Value True..... True  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Pulmonary Care  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

=====

TERM: PULMONARY SECRETION AMOUNT  
ID: {FEED2CE3-F697-D1AB-AB1B-72A7968F0EB0}

=====

Page: 562

Abbreviation..... PUL SEC AM  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... AMOUNT OF SPUTUM FROM PATIENT  
Display Name..... Pulmonary Secretion Amount  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. SCANT (SCANT) {8EF6C62D-0D94-37BB-60AD-86E306E8120D}
  2. SMALL (s) {5991F892-B28C-3998-E733-D8FFA26D7E21}
  3. MODERATE (MODERATE) {C3C25064-C9FC-EA1F-3E74-32A25F0FC5E3}
  4. LARGE (l) {E60D5CF9-920D-44E2-8A83-136DAB1FE251}
  5. COPIUS (COPIUS) {C74A6979-0E35-AA12-F992-69ACA3B0FA23}

UNITS:  
-----



Abbreviation..... SPUTUM COL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Color of sputum  
Display Name..... Pulmonary Secretion Color  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 8 entries

- 
- |                    |  |
|--------------------|--|
| 1. CLEAR (c)       | {8981C365-2B46-704F-C9F4-975DF1C23C2D} |
| 2. WHITE (WHITE)   | {AB04EF54-5246-EF3A-E1D9-02E9E24369DD} |
| 3. YELLOW (YELLOW) | {90DC31AD-6211-9F25-4237-79F7122A2763} |
| 4. TAN (TAN)       | {AB994ED7-8916-FAE0-0E97-CD4E0EF80ED1} |
| 5. BROWN (BROWN)   | {A1C32FBC-7B4D-707C-CA58-92A2E201368E} |
| 6. GREEN (GREEN)   | {3D1F1553-56FC-68CE-79A9-E4F6357B6AAF} |
| 7. RED (RED)       | {515D7CFB-3527-30E7-EC4B-B607C272DFA0} |
| 8. PINK (PINK)     | {1A7FAF84-D69C-4EEB-7BC5-D21A598A1F4D} |

UNITS:

-----

=====

TERM: PULMONARY SECRETION CONSISTENCY  
ID: {633A949F-FED1-7A86-5DEA-C62DAC802B7B}

=====

Page: 564

Abbreviation..... P SECR CON  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Sputum Consistency  
Display Name..... Pulmonary Secretion Consistency  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. THICK (THICK) {1554A49E-1468-46B2-F851-531DDB1D1A98}
  2. TENACIOUS (TENACIOUS) {EE90B04C-F74B-43F8-F3F2-8C2D4359D0B3}
  3. THIN (THIN) {CB8D8B5C-E824-D030-1241-D8EB814671E9}

UNITS:  
-----





=====

TERM: PULMONARY SUCTIONING SITE  
ID: {D011DC21-02E6-463B-DF02-6E4ADD94D703}

=====

Page: 567

Abbreviation..... PUL SUC SI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... SITE OF PATIENTS PULMONARY SUCTIONING  
Display Name..... Pulmonary Suctioning Site  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. ORAL (ORAL) {BBAF0DD0-1F52-49EE-9A47-FA9323D6C529}
  2. NASAL (NASAL) {0F73F0CB-E25E-4611-9625-62E6F1596A9B}
  3. ARTIFICIAL AIRWAY (ART AIRWAY)  
{8DEE2358-8892-4BFB-AFAE-6C134D1CCAD4}

UNITS:  
-----

=====

TERM: PULMONARY VASCULAR RESISTANCE  
ID: {56E4D7F3-EC49-B76A-B223-3C9E62E8B0ED}

=====

Page: 568

Abbreviation..... PVR  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Pulmonary Vascular Resistance  
Display Name..... Pulmonary Vascular Resistance  
Help Text..... Pulmonary Vascular Resistance  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
DYNE \* SECONDS / CM(-5) \* M2 50 > 150 \* 250 < 500      Dec Prec: 0

=====

TERM: PULMONARY VASCULAR RESISTANCE INDEX  
ID: {F98F291B-14F4-AA32-F218-E146364D9754}

=====

Page: 569

Abbreviation..... PVRI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Pulmonary Vascular Resistance Index  
Display Name..... Pulmonary Vascular Resistance Index  
Help Text..... Pulmonary Vascular Resistance Index  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----  
1. CARDIAC OUTPUT MONITOR (COM)  
    {CCBE5744-7CB9-9DFA-938F-FDDBC73B45B3}

UNITS:

-----  
DYNE \* SECONDS / CM(-5) / M2 20 > 80 \* 240 < 350      Dec Prec: 0

Abbreviation..... P
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Beats per minute
Display Name..... Pulse - Do Not Check Total or Subtotals
Help Text..... Heart beats per minute,
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid..... 4500636

LOCATION: 12 entries

- 1. APICAL (Ap) {C5CF5C51-63CC-33FB-9BAE-FB1DD9900911}
2. AV FISTULA (AVF) {F977EEF6-EA6A-D56D-D1F7-29B91510A02F}
3. BILATERAL PERIPHERALS (BILAT PERI) {7F134529-36C5-E9D4-ABD5-64B404F9E706}
4. BRACHIAL (Bra) {67426127-09E1-E566-F31E-15714A5B2B7C}
5. CAROTID (Car) {9536F203-1797-9336-6F69-604249E1B3A0}
6. DORSALIS PEDIS (DP) {9657D743-1632-2044-B4C5-9EF291A2F6D0}
7. FEMORAL (Fem) {E3442C9E-8677-7296-D314-7365319C3716}
8. PERIPHERAL (Per) {36B31698-E84A-5792-8EAE-F374DFD1B55A}
9. POPLITEAL (Pop) {77030DF2-CE06-D9E2-0684-EE70209CBE73}
10. POSTERIOR TIBIAL (PT) {010F2FFE-DEBC-DF9D-9CF2-AD192B16FFD5}
11. RADIAL (RADIAL) {D52C4103-2AB3-A0C7-16E3-CF95A1B3B789}
12. ULNAR (Uln) {81FC2EA0-472D-06D7-F5A2-E71B0B7DE091}

METHOD: 13 entries

- 1. AFTER EXERCISE (AE) {D129EF4B-3A99-9966-54F1-2C01A36814EE}
2. AT REST (AR) {E8C289A6-0C3D-A7B2-BB85-7F1A23F22B50}
3. AUSCULTATE (Aus) {0C629C1F-92D3-C641-B02C-9EA43B401DC9}
4. DOPPLER (Dop) {52683C19-F9DC-876B-CA3F-53E7C864CF04}
5. INVASIVE (I) {5DB8746E-5446-F47A-883C-A44B91275A6E}
6. LEAD (L) {B628FDCE-05E0-DACF-E66E-0F8380AB465A}
7. MONITOR (Mo) {63FC24EA-50B8-B302-D948-ABE3517CF6FC}
8. NON-INVASIVE (N-I) {F50ABFED-5FD7-FC95-20FE-1E3038C9F299}
9. PALPATED (Pal) {295D8E72-12A1-08AF-6BE8-0F9CD238F5EF}
10. WITH ACTIVITY (Act) {544A262B-EA34-F82A-46A9-F80740515D65}
11. ELECTROCARDIOGRAM (ECG) {CF5F1AB9-3B74-7FBA-E38B-C9EF54FE4DC5}
12. ELECTROCARDIOGRAM TRACE 1 (ECG TR 1) {E0D449BA-9F06-4A27-8FE9-24C1BA0746DB}
13. ELECTROCARDIOGRAM TRACE 2 (ECG TR 2) {73EACBC1-C1FB-4B2C-96AF-C5C1CA07DA72}

POSITION: 4 entries

- 1. LYING (Ly) {DE28D2BF-6F8D-64C3-0092-4CEB5E0902B9}
2. SEMIFOWLERS (Sfl) {C2F6A945-3829-AABA-3D67-54D1829D3EB9}
3. SITTING (Sit) {DD4DCDF1-0A22-2CC5-404F-0E472C73961C}
4. STANDING (St) {0776AE42-C370-B7B3-9B4D-9C914BA8F231}

QUALITY: 2 entries

- 1. ACTUAL (A) {5842A4D9-56A1-7A5A-13C9-F481699C3A09}
2. ESTIMATED (E) {9EAF26B-7D87-DA66-F991-022D095FB5A0}



UNITS:

-----  
BEATS PER MINUTE

0 > 60 \* 100 < 300

Dec Prec: 0

Abbreviation..... PULSE STR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Quality of pulse
Display Name..... Pulse Strength
Help Text..... Select one of the measurements for
pulse quality
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 6 entries

- 1. 3 - STRONG (3-s) {C55EAB72-7816-5440-84E9-6E7E190ED197}
2. 2 - NORMAL (2-n) {25FAADFC-D280-78AC-4F01-CFCF34ECD3FF}
3. 1 - WEAK (1 - WEAK) {A41570DC-2F35-DDF6-977B-7C262150DFA2}
4. 0 - ABSENT (0-a) {75F14311-F113-963C-4F05-B8AE17F044A9}
5. D - DOPPLER (d-d) {2B77F099-0520-474D-A4BC-22BF824C8364}
6. AMPUTEE (AMPUTEE) {5804286F-2B23-91BD-FCB5-6AD18C918203}

LOCATION: 8 entries

- 1. R RADIAL (R RADIAL) {BA10E6E5-DDF4-345E-68AF-21D32E9D972D}
2. L RADIAL (L RAD) {302F78AD-4067-8B6B-BD15-993EC8561AE6}
3. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
4. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
5. R POST TIBIAL (R PO TIB) {88CDA53A-5153-4950-701A-3F4F2007CC6B}
6. L POST TIBIAL (L PO TIB) {BED130EE-31CC-C7E9-6597-78685BCE91AC}
7. R PEDAL (R PEDAL) {9532B614-F4E1-316C-1322-905EB58E58C4}
8. L PEDAL (LEFT PEDAL) {D84E9953-8C17-3130-8F21-A30D2A99FDE4}

UNITS:

ID: {9FCBF7CF-39B8-9710-6F41-9B7DDDC3BDD2}

=====

Abbreviation..... PREA  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Pupil Reaction  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 9 entries

-----

- 1. FIXED (f) {697D8C5A-976F-F423-30EC-A04102CF85F0}
- 2. SLUGGISH (s) {3215A82E-8DCC-C7C0-8AE1-2D4C9D4A4EFB}
- 3. EYE CLOSED (ec) {140B7D2F-F9C2-0B4C-3CDC-7355EDDB5030}
- 4. DRESSING (d) {FFEF0714-3B5B-E2B4-D760-3C03D45E59EA}
- 5. REACTIVE (Reactive) {31B44ADB-0000-A5B0-5F56-59D4A08C05D2}
- 6. NONREACTIVE (nonreact) {AB6B7D92-B94C-F3AA-E32F-6035D9CF2EF9}
- 7. BRISK (brisk) {980AE687-0A78-458B-BA7B-5D646C38A5F5}
- 8. PINPOINT (pinpoint) {4E53DDAC-0BF4-66D5-D614-9763BBF1FC8E}
- 9. DILATED (Dilated) {BF98ADC4-C69E-B17A-9BDE-D426D7F4C395}

LOCATION: 2 entries

-----

- 1. L EYE (L EYE) {D714A0CA-EAD6-EC35-DA54-FF027C23FF45}
- 2. R EYE (R EYE) {258F0C1D-67BD-ECF8-3152-EB289DEF4681}

UNITS:

-----

Abbreviation..... PS  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Pupil Size  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

LOCATION: 2 entries

- 
- 1. L EYE (L EYE) {D714A0CA-EAD6-EC35-DA54-FF027C23FF45}
  - 2. R EYE (R EYE) {258F0C1D-67BD-ECF8-3152-EB289DEF4681}

UNITS:

-----

MILLIMETERS 0 > 1 \* 9 < 11 Dec Prec: 0

Abbreviation..... QRS  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... The QRS complex is a structure on the  
ECG that corresponds to the  
depolarization of the ventricles.  
Display Name..... QRS Complex  
Help Text..... The QRS complex is a structure on the  
ECG that corresponds to the  
depolarization of the ventricles.  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
SECOND                    0.02 > 0.04 \* 0.12 < 0.24 Dec Prec: 2

Abbreviation..... QT  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... The QT interval is measured from the  
beginning of the QRS complex to the  
end of the T wave.  
Display Name..... QT Interval  
Help Text..... The QT interval is measured from the  
beginning of the QRS complex to the  
end of the T wave.  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
SECOND                    0.10 > 0.30 \* 0.46 < 1.00 Dec Prec: 2

=====

TERM: RECTAL TUBE PRESENT  
ID: {D05DB527-8EDD-4324-B7B5-989AC004486E}

=====

Page: 576

Abbreviation..... REC TU PRE  
Active..... Yes  
Boolean Value False..... Not Present  
Boolean Value True..... Present  
Data Type..... Boolean  
Description..... Rectal Tube Present  
Display Name..... Rectal Tube Present  
Help Text..... Rectal Tube Present  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

Abbreviation..... RESP PATT
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Respiratory Pattern
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 11 entries

- 1. REGULAR (REGULAR) {3D3F3D4C-ECD0-2056-A0D1-43CC2F2BA2C8}
2. IRREGULAR (IRREGULAR) {6CCE27A2-F1EB-3EC8-B140-CDA241D159F2}
3. DEEP (DEEP) {78AD569C-16D4-CFC3-6B6F-41FB9FFAB9D7}
4. SHALLOW (SHALLOW) {CD205617-BF98-4177-971A-056AFECC823A}
5. DYSPNEIC (DYSPENIC) {4C6B42D9-694F-5EB2-6396-7ED85D5B154E}
6. ORTHOPNEIC (ORTHOPNIC) {FF373BD7-FD84-DF43-7A44-1B8E62917116}
7. ABDOMINAL (ABDOMINAL) {29E279F4-5D2E-4644-A2FA-CF85BEDFD480}
8. CHEYNE-STOKES (CHEYNE-STO) {58DD2B80-B522-499C-2C60-C48A654BAD84}
9. KUSSMAUL (KUSSMAUL) {0C5ED066-F445-ABFE-F7B1-D158EC7C2976}
10. AGONAL (a) {1621B411-33E3-3A6A-78A6-40B7BDC3E8CB}
11. VENTILATOR (v) {A5A647ED-A56D-5012-AD5B-788B3DCB3498}

UNITS:



Abbreviation..... RR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... RESPIRATIONS PER MINUTE
Display Name..... Respiratory Rate - Do Not Check Total
or Subtotals
Help Text..... Breaths per minute
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid..... 4688725

METHOD: 8 entries

- 1. AFTER EXERCISE (AE) {D129EF4B-3A99-9966-54F1-2C01A36814EE}
2. ASSISTED VENTILATOR (AV) {A540F333-5A09-5F21-7B51-2B5DF05E32C0}
3. AT REST (AR) {E8C289A6-0C3D-A7B2-BB85-7F1A23F22B50}
4. CONTROLLED VENTILATOR (CV) {0CFC2183-5245-E7AB-96ED-93AF4A3FEC64}
5. MONITOR (Mo) {63FC24EA-50B8-B302-D948-ABE3517CF6FC}
6. SPONTANEOUS (S) {5E9463BC-C871-54FD-F800-DCFDEDB8A93B}
7. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}
8. WITH ACTIVITY (Act) {544A262B-EA34-F82A-46A9-F80740515D65}

POSITION: 6 entries

- 1. LYING (Ly) {DE28D2BF-6F8D-64C3-0092-4CEB5E0902B9}
2. SITTING (Sit) {DD4DCDF1-0A22-2CC5-404F-0E472C73961C}
3. STANDING (St) {0776AE42-C370-B7B3-9B4D-9C914BA8F231}
4. SEMIFOWLERS (Sfl) {C2F6A945-3829-AABA-3D67-54D1829D3EB9}
5. SUPINE (SUPINE) {DEA8A744-B445-4A9E-875D-485E060F2ACC}
6. TRENDELENBURG (Tr) {42F3353B-7041-EEE1-A87E-941D3BE1CBCD}

UNITS:

BREATHS PER MINUTE -1 > 12 \* 20 < 201 Dec Prec: 0

Abbreviation..... RES CRITER  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Restraint Criteria  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 7 entries

- 
1. RES CRIT: LINES (R.C. LINE)  
{27F545F2-5338-4BC8-B5CE-1D321D28F8D6}
  2. RES CRIT: EQUIPMENT (R.C. EQUI)  
{38C2C56C-6D57-4750-8E68-F1193B16190A}
  3. RES CRIT: AMBULATE (R.C. AMBU)  
{54B72608-2C13-46CD-8034-00179669762C}
  4. RES CRIT: UNABLE TO FOLLOW DIRECTIONS (R.C. DIRE)  
{41358EAE-9912-48F8-8931-D3C635406D21}
  5. RES CRIT: COMBATIVE (R.C. COMB)  
{38DBF61A-D364-4215-B4AA-36B660F96462}
  6. RES CRIT: AGGRESSIVE (R.C. AGGR)  
{7A47B361-85EA-4C88-A1AF-CD8FF7D35D97}
  7. RES CRIT: HARM TO SELF (R.C. HARM)  
{4DC33341-8F53-4D30-BC28-1E32DB54F30D}

UNITS:  
-----

=====

TERM: RESTRAINT EDUCATION  
ID: {08B0A919-200B-4874-993B-E9924E9A439F}

=====

Page: 580

Abbreviation..... RES EDUCAT  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Restraint Education  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
1. FAMILY INFORMED (FAM INFORM)  
    {77135D32-9CDD-4D3F-9B62-6E8189B77844}
  2. PATIENT INFORMED (PAT INFORM)  
    {074A1D19-3920-49F1-8958-2EB6F9CF05BA}

UNITS:  
-----

=====

TERM: RESTRAINT EDUCATION TYPE  
ID: {05693D33-977E-4089-988F-4FE029D64BE4}

=====

Page: 581

Abbreviation..... RES EDU  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Restraint Education Type  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. RESTRAINTS EXPLAINED (RES EXPL)  
    {1C063DD2-C4A3-4027-A245-6DD51DF5A3E9}
  2. RESTRAINTS ED TYPE (RES ED TY)  
    {5560A332-C365-4DAF-BA33-610C4AE301F5}
  3. RESTRAINTS ED REMOVAL (RES ED REM)  
    {B19316C8-687B-48FB-9EBC-9857662B9F85}
  4. RESTRAINTS BROCHURE (RES BRO P)  
    {1E54C2A1-8B8D-4F0D-B247-91E19722CDED}
  5. RESTRAINTS TIME LIMIT (RES TIME L)  
    {513B483E-AD7B-4B1A-BD01-91133C605C3C}

UNITS:  
-----

Abbreviation..... R.E. ALT  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Alternative Measures Attempted  
Help Text..... Restraint Evaluation: Alternative  
Measures Attempted  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 11 entries

- 
1. RES EVAL: DECREASE STIMULATION (R.E. DECR)  
{0841626E-994E-4346-B701-1662A00209EB}
  2. RES EVAL: REALITY ORIENTATION (R.E. REALI)  
{221C363C-2261-40FC-8400-8B1D2E693F69}
  3. RES EVAL: RELOCATE (R.E. RELO)  
{D8165678-ECEF-4F4D-A1C4-C616349E1873}
  4. RES EVAL: INVOLVEMENT WITH FAMILY (R.E. INVO)  
{0EF31047-70A5-4072-953A-B62E1ACDF84B}
  5. RES EVAL: AMBULATION/EXERCISE (R.E. EXER)  
{46F4D1EB-3439-4BE4-A78C-453BDBF6C230}
  6. RES EVAL: EDUCATION (R.E. EDUC)  
{76341E6E-7F0A-429C-B768-9CA83C452FD7}
  7. RES EVAL: DIVERSIONAL ACTIVITIES (R.E. DIVER)  
{569F6D49-C564-4910-835A-B05179086E1F}
  8. RES EVAL: MEDICATION REVIEW (R.E. REVI)  
{914A0F82-5863-4C84-B3AB-D5A22BD340D8}
  9. RES EVAL: UNABLE TO ATTEMPT (R.E. UNAB)  
{979B0477-5789-4319-9562-6752048C4005}
  10. BED EXIT ALARM (BED EXIT A)  
{AC53759E-FFAE-5E0F-2EBA-E56A2712357C}
  11. WHEELCHAIR EXIT ALARM (WHEEL EXIT)  
{7D06564E-ECDE-4A11-BEFE-92AB7E448AA3}

UNITS:  
-----









Abbreviation..... RES EV LRM  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Least Restrictive Measures Evaluated  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 10 entries

- 
1. RESTRAINTS LIMIT SETTING (RES LIM SE)  
{F97BA57E-ED33-D6FE-5E35-E6EFB0625773}
  2. REDIRECT PATIENT ACTIVITY (REDIR ACT)  
{BFF0ABB1-953A-4275-9EF2-08BFF660AA8D}
  3. REMOVE FROM STIM. SETTING (REM SET)  
{2010BA89-B399-4A30-9ACF-C65AB57B6351}
  4. OFFER PRN MEDS (OFFER PRN)  
{B7AF4903-913C-4848-A255-0B4767B81AE2}
  5. VERBAL DE-ESCALATION (VERB DE-ES)  
{8AB46F8C-B4D5-4800-BCF7-52CA3F32E99B}
  6. ONE TO ONE SUPERVISION (1 SUPERV)  
{CFC9C653-BFE9-40B2-BA7A-944116C7FCB7}
  7. FREQUENT VERBAL REMINDERS (FR VER RE)  
{F30D791D-BB80-4AB3-B37F-22B6B547EA65}
  8. PREVENT ACCESS TO DEVICE (PREV ACCE)  
{E7455083-DD2C-4562-A6E2-9626A4CDC4DD}
  9. REORIENTED TO ENVIRONMENT (REORIENT)  
{4C47B764-4BA7-4D71-9B0B-D42F120E9A4A}
  10. REPOSITION (REPOSITION) {DF3039CA-5977-4454-81D4-DB0AF966F85D}

UNITS:  
-----





=====

TERM: RESTRAINT EVAL - REPOSITION/RANGE OF MOTION  
ID: {54C82D71-B32F-48A2-87A9-683D6023D905}

=====

Page: 589

Abbreviation..... R.E. REPOS  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Reposition / Range of Motion Evaluated  
Help Text..... Restraint Evaluation: Reposition /  
Range of Motion Evaluated  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEA649CD} |

UNITS:  
-----





=====

TERM: RESTRAINT INITIATED BY  
ID: {6939E11B-F414-435F-93B2-D69334A50FA1}

=====

Page: 592

Abbreviation..... RES INIT  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Seclusion/Restraint Initiated By  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
1. REGISTERED NURSE (RN) {524F344F-FDCD-423F-AD34-291D5DB4271B}
  2. PHYSICIAN ORDER (PHY ORDER)  
    {AB89E640-83B7-4E6C-8C8D-2FA8C2564CCB}

UNITS:  
-----

Abbreviation..... RES LRME
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Least Restrictive Measures Evaluated
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 10 entries

- 1. RES LRM: LIMIT SETTING (RLRM SETT) {96417DDF-8562-4EFF-B4F6-0933FD91CD36}
2. RES LRM: REDIRECT (RLRM REDI) {07DE67AF-E226-4794-ADF5-17E9C56C9FDF}
3. RES LRM: REMOVE FROM SETTING (RLRM REMO) {857DBF6E-7EFF-4F07-9AC1-7FDB37CC7618}
4. RES LRM: PRN MEDS (RLRM: MEDS) {CDF02FA3-FC72-4CCC-8197-41710A69A2CB}
5. RES LRM: VERBAL DE-ESCALATION (RLRM VERB) {CE8618C8-CAD7-4C74-8566-ACB5E66B082C}
6. RES LRM: SUPERVISION (RLRM SUPR) {18FA8CB0-0155-4D1C-9B14-AE51A0C06280}
7. RES LRM: VERBAL REMINDERS (RLRM REMI) {987B5A32-DE2E-47E5-9B06-1AEE86096AB2}
8. RES LRM: DEVICE ACCESS (RLRM ACCE) {1662F5E1-30E7-4E8A-BCC5-1738868435B8}
9. RES LRM: REORIENT (RLRM REOR) {DFCC8245-180D-4F3F-AEFC-5D55AFDDC46D}
10. REPOSITION (REPOSITION) {DF3039CA-5977-4454-81D4-DB0AF966F85D}

UNITS:



Abbreviation..... RES MENT A
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Mental Assessment Patient Evaluation
Display Name..... Mental Assessment
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

- 1. ORIENTED (ORIENTED) {B8DC9245-6D81-40E4-8B47-B6413E49F152}
2. DISORIENTED (DISORIENTE) {A7D1F2EA-CD7C-4100-8552-44834684C6DE}
3. RESTLESS (r) {8FCCAD26-DBA0-7067-DA72-30D342BFA8E8}
4. SLEEPING (s) {87338420-F351-F98C-7A80-4990B87ACE2B}
5. SEDATED (s) {29D9ADD6-2C1F-F6D5-F4D0-B2BB86416997}

UNITS:

=====

TERM: RESTRAINT ORDER TYPE  
ID: {843976A3-7C4A-4E09-A273-4589923DDDC2}

=====

Page: 595

Abbreviation..... RES ORD TY  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Order Type  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
1. SECLUSION (SECLUSION)      {65B57F6C-C849-4BF7-8705-175361BC900F}
  2. RESTRAINT (RESTRAINT)      {0429322D-A987-4913-B64D-A6414645B214}

UNITS:  
-----



=====

TERM: RESTRAINT SAFETY CARE  
ID: {1450E52B-CBF3-BA02-F4E0-4AE793AE9614}

=====

Page: 597

Abbreviation..... RESSC  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Restraint Safety Care  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
1. CIRCULATION CHECK (ccq) {BA32C2A8-F15C-84EF-4859-6277EFA6463A}
  2. DEVICE RELEASE (drq) {9F97947F-CFFE-3285-F4C3-A8DF294DA46C}

UNITS:  
-----

Abbreviation..... RES TYPE A
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Restraint Type
Help Text..... Location and type of restraint used
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 19 entries

- 1. RESTRAINT SOFT (RESTR SOFT) {C02CAC4B-5E9C-4060-BECD-0C0820EDA39B}
2. RESTRAINT SOFT 2 POINT (RES S2POIN) {66B71CFE-2603-40C9-9625-24448012BAB5}
3. RESTRAINT SOFT 4 POINT (RES S4POIN) {6820E997-60A0-4802-BE26-CA29B4F48F32}
4. RESTRAINT SOFT ARMS + WAIST (RES S A+W) {4A322BB6-0F63-4BB7-9290-C0655BD7D772}
5. RESTRAINT SOFT ARMS + LOCKED WAIST (RES S A+LW) {C480D67C-958C-42F9-9203-98E2F2BDBEC8}
6. RESTRAINT SINGLE HAND MITT (RES S HA M) {A5F1D037-FB5E-402B-9F25-22DE48404D57}
7. RESTRAINT BILATERAL HAND MITT (RES BI H M) {5E293110-1235-4307-80C2-0FFFA007E1DA}
8. RESTRAINT VEST (RES VEST) {643378B2-86D7-4DB7-93ED-9DCE16D85E27}
9. RESTRAINT BILATERAL ANKLE (RES BI ANK) {8A0B0DE5-E652-40E7-9B67-2F4D27685A68}
10. RESTRAINT WRIST (RES WRIST) {03038E52-A66F-4EC0-8DC1-7B0BB596FFC0}
11. RESTRAINT GERI CHAIR (RES GERI C) {6B8F87BB-C2AB-4D92-AE04-24F19B26D2F1}
12. RESTRAINT HARD PLASTIC/LEATHER 2 POINT (RES P/L 2P) {17CB8739-612A-4319-946A-2062C5245C90}
13. RESTRAINT HARD PLASTIC/LEATHER 4 POINT (RES P/L 4P) {94DFD0BA-E941-4B09-8FD7-87E9AD7C7E6F}
14. RESTRAINT ONE ARM (RES ONE AR) {839914DA-6D63-43EB-961E-70DA56737A92}
15. RESTRAINT ONE LEG (RES ONE LE) {2790F0FC-B49A-4C99-9968-BCADD7CC63AC}
16. RESTRAINT WAIST LOCKED (RES WAI LO) {1EE7AB42-CE17-4D49-849B-CA550966E58C}
17. RESTRAINT WAIST (RES WAIST) {253E2757-5A57-4410-9C0E-F4A907A1EC16}
18. RESTRAINT HARD PLASTIC/LEATHER (RES H PL/L) {9BB528F4-7016-4030-B73F-4037ADF2621D}
19. RESTRAINT STRAP (RES STRAP) {642FDACE-B320-42BC-9514-F5FCD6583B19}

LOCATION: 11 entries

- 1. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
2. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
3. L WRIST (L WRIST) {7BCAFACE-FF3B-46DA-99F9-6437A90B7126}
4. BILATERAL WRIST (BILAT WRIS)

- {89B4D30F-9E46-F8D6-14F6-1CC0843DDC21}
- 5. BILATERAL LEG (BILAT LEG) {0713761D-E0FC-DA03-02A1-F20EFEB99F45}
- 6. BILATERAL ANKLES (BILA ANKLE)  
{E56A03A5-5CE0-47BF-A14D-F6DF66F4E769}
- 7. CHEST (CHEST) {8C3C30FE-6E06-C35E-8226-CBF0C94014CC}
- 8. L LEG (LL) {B088C809-E528-B785-21C2-14E45E1FD8AF}
- 9. R LEG (RL) {C5E40A7B-6FC4-267D-8CE8-3571D21D81FD}
- 10. R WRIST (R WRIST) {FDD9FE1D-F8A2-09AD-82D7-AA478693DD01}
- 11. WAIST (WAIST) {784DD311-8B48-4FE3-96B6-471CCE453CD5}

UNITS:

-----

Abbreviation..... RES TY ARM  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Arm Restraints  
Display Name..... Arm Restraints  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. RESTRAINT SINGLE HAND MITT (RES S HA M)  
{A5F1D037-FB5E-402B-9F25-22DE48404D57}
  2. RESTRAINT HARD PLASTIC/LEATHER (RES H PL/L)  
{9BB528F4-7016-4030-B73F-4037ADF2621D}
  3. RESTRAINT SOFT (RESTR SOFT)  
{C02CAC4B-5E9C-4060-BECD-0C0820EDA39B}

LOCATION: 6 entries

- 
1. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
  2. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
  3. BOTH ARMS (BOTH ARMS) {A493AB32-9DE9-40AC-8C46-7EFF395B63B5}
  4. L WRIST (L WRIST) {7BCAFACE-FF3B-46DA-99F9-6437A90B7126}
  5. R WRIST (R WRIST) {FDD9FE1D-F8A2-09AD-82D7-AA478693DD01}
  6. BILATERAL WRIST (BILAT WRIS)  
{89B4D30F-9E46-F8D6-14F6-1CC0843DDC21}

UNITS:

-----

Abbreviation..... RES TY LEG  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Leg Restraints  
Display Name..... Leg Restraints  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. RESTRAINT HARD PLASTIC/LEATHER (RES H PL/L)  
{9BB528F4-7016-4030-B73F-4037ADF2621D}
  2. RESTRAINT SOFT (RESTR SOFT)  
{C02CAC4B-5E9C-4060-BECD-0C0820EDA39B}
  3. RESTRAINT STRAP (RES STRAP)  
{642FDACE-B320-42BC-9514-F5FCD6583B19}

LOCATION: 6 entries

- 
1. L LEG (LL) {B088C809-E528-B785-21C2-14E45E1FD8AF}
  2. R LEG (RL) {C5E40A7B-6FC4-267D-8CE8-3571D21D81FD}
  3. BOTH LEGS (BOTH LEGS) {5A15A231-4695-4CF3-AD0F-2E7950208967}
  4. L ANKLE (L ANKLE) {D5A692EF-9D7C-473C-AF96-7EF1C01FA2AF}
  5. R ANKLE (R ANKLE) {4E6DF352-F399-4CE9-908C-48132C3D856B}
  6. BILATERAL ANKLES (BILA ANKLE)  
{E56A03A5-5CE0-47BF-A14D-F6DF66F4E769}

UNITS:

-----



=====

TERM: RESTRAINT TYPE OTHER  
ID: {40448F75-66C2-4FA1-AA02-6B827C9A06F8}

=====

Page: 601

Abbreviation..... RES TYP OT  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Restraint Type Other  
Display Name..... Restraint Type Other  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. GERI-CHAIR (g) {167F9786-8F24-054F-020A-658A6EE5B0B4}
  2. RESTRAINT VEST (RES VEST) {643378B2-86D7-4DB7-93ED-9DCE16D85E27}
  3. RESTRAINT WAIST (RES WAIST)  
{253E2757-5A57-4410-9C0E-F4A907A1EC16}

LOCATION: 3 entries

- 
1. TORSO (TORSO) {EA1AD477-8BFF-E575-23C1-B9CFAF2B9986}
  2. CHEST (CHEST) {8C3C30FE-6E06-C35E-8226-CBF0C94014CC}
  3. WAIST (WAIST) {784DD311-8B48-4FE3-96B6-471CCE453CD5}

UNITS:  
-----





=====

TERM: RIGHT VENTRICULAR STROKE WORK  
ID: {3A72DA75-5112-7980-4F4F-E1E34FFCD1A6}

=====

Page: 604

Abbreviation..... RVSW  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Right Ventricular Stroke Work  
Display Name..... Right Ventricular Stroke Work  
Help Text..... Right Ventricular Stroke Work  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----

1. CARDIAC OUTPUT MONITOR (COM)  
{CCBE5744-7CB9-9DFA-938F-FDDBC73B45B3}

UNITS:

-----

STROKE WORK PER BEAT 0.0 > 5.0 \* 30.0 < 50.0 Dec Prec: 1

=====

TERM: RIGHT VENTRICULAR STROKE WORK INDEX  
ID: {EE450A40-E87C-1F5C-FDEF-BAAD030B8149}

=====

Page: 605

Abbreviation..... RVSWI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Right Ventricular Stroke Work Index  
Display Name..... Right Ventricular Stroke Work Index  
Help Text..... Right Ventricular Stroke Work Index  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----

1. CARDIAC OUTPUT MONITOR (COM)  
{CCBE5744-7CB9-9DFA-938F-FDDBC73B45B3}

UNITS:

-----

STROKE WORK PER BEAT PER M2 1.9 > 2.0 \* 18.0 < 18.1 Dec Prec: 1  
GM-MM2/BEAT 3.9 > 4.0 \* 18.0 < 18.1 Dec Prec: 1

=====

TERM: SAFETY DEVICE WAIST BELT APPLIED  
ID: {27E75399-989F-4598-A972-646FBFF51A20}

=====

Page: 606

Abbreviation..... WAIST BE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Safety Device Waist Belt Applied  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEA649CD} |
| 2. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |

UNITS:  
-----



Abbreviation..... 2nd Verf  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... Second Verifier  
Help Text.....  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 8 entries

- 
- |               |  |
|---------------|--|
| 10. RN (RN)   | {A103E506-EAA5-4777-8126-268CD723AD87} |
| 20. LVN (LVN) | {8F8300FD-5A33-4B01-9F0B-901E2A3D96AB} |
| 30. NA (NA)   | {1993F8D9-D536-4ECA-AFAD-366F79BA0CBD} |
| 40. MD (MD)   | {9AE2F01F-A8EF-4423-BD1C-55E9A1CD50A8} |
| 50. DO (DO)   | {F70738A6-7BFA-4D76-881C-6A4B4F4E41E4} |
| 60. RT (RT)   | {6E374EFA-B7CE-44C2-ADCE-DF9C132E951E} |
| 70. NP (NP)   | {E7EA4017-4374-41A1-9D3D-54051CEDF63C} |
| 80. PA (PA)   | {C057C47D-780A-4230-B577-2C528AD09DEB} |

UNITS:

-----



=====

TERM: SEDATION SCALE - RAMSEY  
ID: {CF2E7E29-8C96-461B-9074-CA94455F62D6}

=====

Page: 609

Abbreviation..... SED RAMSEY  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Ramsey Sedation Scale  
Display Name..... RAMSEY  
Help Text..... Ramsey Sedation Scale  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 6 entries

- 
1. RAMSEY ANXIOUS (RAM ANXIOUS)  
{88205A0C-7065-4616-8FCF-75059462DD15}
  2. RAMSEY COOPERATIVE (RAM COOPER)  
{56846ADA-0F0E-44AA-9779-A567B9D24C4F}
  3. RAMSEY SEDATED (RAM SEDATED)  
{CC12405F-60DB-485F-94E5-78C1EB32AB74}
  4. RAMSEY ASLEEP BRISK RESPONSE (RAM ASL BR)  
{00F4A117-F9E7-4B60-A521-F2E8725AECD3}
  5. RAMSEY ASLEEP SLUGGISH (RAM ASL SL)  
{AA61604F-00F9-4103-B27F-E817FBC9577B}
  6. RAMSEY ASLEEP NO RESPONSE (RAM ASL NO)  
{3FF5BA0A-2D3C-4935-A432-23DAA4761AB7}

UNITS:  
-----

Abbreviation..... RASS
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Richmond Agitation Sedation Scale
Display Name..... RASS
Help Text..... Richmond Agitation Sedation Scale
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 10 entries

- 1. RASS COMBATIVE (RASS COMBA) {A97E064A-7709-42F9-9A88-4AAFF1020883}
2. RASS VERY AGITATED (RASS V AGI) {7528CA1D-A9A4-4F4E-8943-E69F3B42E189}
3. RASS AGITATED (RASS AGITA) {D81A5D5D-E7FC-445A-A814-495731ABA0F8}
4. RASS RESTLESS (RASS RESTL) {2CCFC1ED-D807-45F0-AE69-980E214AD5C3}
5. RASS ALERT AND CALM (RASS AL CA) {F1834BF6-F8F1-4E13-BC25-B13E756E0FF0}
6. RASS DROWSY (RASS DROWS) {5E6BDAC6-4F81-4110-A0AF-FED746D9FF0D}
7. RASS LIGHT SEDATION (RASS LT SE) {FDEEC084-3922-4065-A863-4893E42F8F93}
8. RASS MODERATE SEDATION (RASS MD SE) {D12B0A3E-B80E-494F-B89F-A79E8B81B080}
9. RASS DEEP SEDATION (RASS DP SE) {83CB4E64-3581-452F-987A-D31CB3AB38D2}
10. RASS UNAROUSABLE (RASS UNARO) {1578E147-7D2C-4469-828E-FCA97109B27C}

UNITS:

=====

TERM: SEDATION SCALE - SAS  
ID: {8281E806-CEC9-4236-B984-495B8A3973E5}

=====

Page: 611

Abbreviation..... SED SC SAS  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Riker Sedation-Agitation Scale  
Display Name..... SAS  
Help Text..... Riker Sedation-Agitation Scale  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 7 entries

- 
1. SAS DANGEROUS AGITATION (SAS DA AGI)  
    {0701FB25-1DF0-4771-8C81-AC19A985A3B6}
  2. SAS VERY AGITATED (SAS VY AGI)  
    {B501606A-EEA7-470E-81BE-F51D422272F1}
  3. SAS AGITATED (SAS AGITAT) {107143D4-9527-4A26-8CB7-EEBB656B68B3}
  4. SAS CALM AND COOPERATIVE (SAS C COOP)  
    {F0633778-280C-43C3-9E6E-9C35C18179F6}
  5. SAS SEDATED (SAS SEDATE) {AFEC3F29-52D3-4A49-91F8-05AA3891FEC8}
  6. SAS VERY SEDATED (SAS VY SED)  
    {4B13529B-61C5-4B68-857A-3A5351DC4BE2}
  7. SAS UNAROUSABLE (SAS UNAROU)  
    {3DDFF9FE-A88A-471F-BF36-E5B569DEBDAA}

UNITS:  
-----

Abbreviation..... CPRS note  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... String  
Description.....  
Display Name..... See CPRS Progress Note for Additional  
Data and Assessments  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:  
-----

=====

TERM: SEIZURE CHARACTER  
ID: {15B23878-6E4C-4FCE-9F49-1C728A4A65D2}

=====

Page: 613

Abbreviation..... SEIZU CHAR  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Seizure Character  
Display Name..... Seizure Character  
Help Text..... Character of Seizure  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. TONIC/CLONIC (TON/CLON) {81932BAA-6C79-47A3-85B3-592824B6C09C}
  2. PARTIAL IMPAIRED (PART IMPAI)  
{2F1B8C52-1C94-4EA7-AB82-117694FD4D52}
  3. PARTIAL SIMPLE (PART SIMP)  
{F4A5E7ED-92AB-4021-A63C-E1D948EFE6C1}

UNITS:  
-----

=====

TERM: SEIZURE DURATION  
ID: {DAABDF54-F64C-6306-87EA-B8F5DE034DA0}

=====

Page: 614

Abbreviation..... SEID  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Seizure Duration  
Help Text..... Value in minutes  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MINUTE                    0 > 1 \* 15 < 30                    Dec Prec: 0

=====

TERM: SEIZURE FOCI

Page: 615

ID: {719641C6-BA10-0BA3-65F1-17DC4B1B92E0}

=====

Abbreviation..... SEIF  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Seizure FOCI  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

-----

1. STARTED ON LEFT (sol) {AA6767EE-57F6-F565-76E9-EB219644492F}
2. STARTED ON RIGHT (sor) {0C690E36-E710-244D-44FC-5DB5BCC0FF7E}
3. UNKNOWN (u) {1F872397-8635-18A3-E252-AE9B20B88C7D}

UNITS:

-----

=====

TERM: SEIZURE PRECAUTIONS  
ID: {9BD54461-23E5-F317-E5AC-5FAF87EECA0C}

=====

Page: 616

Abbreviation..... SEI PREC  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Seizure Precautions  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:

---



=====

TERM: SEIZURES

Page: 617

ID: {BB006255-D69A-E29C-1286-D71E6FF15830}

=====

Abbreviation..... SEI  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Seizures  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

=====

TERM: SENSATION FOOT  
ID: {CCF961DC-8590-2498-3836-0F010461F69A}

=====

Page: 618

Abbreviation..... SENSF  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Sensation Foot  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
- |                 |  |
|-----------------|--|
| 1. INTACT (i)   | {935C9CAF-6C88-6610-C6E9-4C1B822923AF} |
| 2. NUMBNESS (n) | {037F2F79-36D9-A0B3-C23E-8EE525932AF3} |
| 3. TINGLING (t) | {AE8E3101-AC3D-DEC3-0293-58F901F82EB8} |
| 4. ABSENT (a)   | {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0} |

LOCATION: 2 entries

- 
- |                    |  |
|--------------------|--|
| 1. L FOOT (L FOOT) | {2EC3FBEC-4EA2-9E7A-17D4-FDC2F9D28046} |
| 2. R FOOT (R FOOT) | {7D9F5FF6-232C-2488-4504-A721FF06928F} |

UNITS:  
-----

=====

TERM: SENSATION HAND  
ID: {7B880F9C-83F8-3395-9983-F10A1EF2F97A}

=====

Page: 619

Abbreviation..... SENSH  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Sensation Hand  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
- 1. INTACT (i) {935C9CAF-6C88-6610-C6E9-4C1B822923AF}
  - 2. NUMBNESS (n) {037F2F79-36D9-A0B3-C23E-8EE525932AF3}
  - 3. TINGLING (t) {AE8E3101-AC3D-DEC3-0293-58F901F82EB8}
  - 4. ABSENT (a) {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0}

LOCATION: 2 entries

- 
- 1. L HAND (L HAND) {639DF25A-19E3-B504-D323-4B3F6ED39A24}
  - 2. R HAND (R HAND) {A54D3A80-8537-D079-9CCE-7DC56346D132}

UNITS:

-----

Abbreviation..... SLEX  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Sensation Lower Extremity  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
- 1. INTACT (i) {935C9CAF-6C88-6610-C6E9-4C1B822923AF}
  - 2. NUMBNESS (n) {037F2F79-36D9-A0B3-C23E-8EE525932AF3}
  - 3. TINGLING (t) {AE8E3101-AC3D-DEC3-0293-58F901F82EB8}
  - 4. ABSENT (a) {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0}

LOCATION: 2 entries

- 
- 1. L LOWER EXTREMITY (LLEX) {FC6F5056-779A-BF76-86D0-A237936C1ECC}
  - 2. R LOWER EXTREMITY (R LOW EXTR)  
{7FC3CECC-1E6E-A686-FF04-2BEF8A9B9DCA}

UNITS:

-----

=====

TERM: SENSATION UPPER EXTREMITY  
ID: {25313956-8B3E-FA21-0AA7-8E9DD4ABFF8D}

=====

Page: 621

Abbreviation..... SUEX  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Sensation Upper Extremity  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
- |                 |  |
|-----------------|--|
| 1. INTACT (i)   | {935C9CAF-6C88-6610-C6E9-4C1B822923AF} |
| 2. NUMBNESS (n) | {037F2F79-36D9-A0B3-C23E-8EE525932AF3} |
| 3. TINGLING (t) | {AE8E3101-AC3D-DEC3-0293-58F901F82EB8} |
| 4. ABSENT (a)   | {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0} |

LOCATION: 2 entries

- 
- |                                   |  |
|-----------------------------------|--|
| 1. L UPPER EXTREMITY (LUEX)       | {48806E79-B1DA-129B-EB52-453292A3DDC6} |
| 2. R UPPER EXTREMITY (R UPP EXTR) | {C65BA269-8F05-6FB9-62CB-466A286CF2C0} |

UNITS:

-----

Abbreviation..... SEQ DEV AP
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Sequential Compression Device
Application
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 8 entries

- 1. INITIATED (INITIATED) {562E5FFB-4E53-4614-ACAC-E7D9B6E1B781}
2. CONTINUOUS (CONTINUOUS) {676CD2DA-F882-4BE5-99BB-510EA7669749}
3. INTERRUPTED (INTERRUPTED) {AB53B9AB-DA3F-42D3-BEB2-E02A452ED777}
4. DISCONTINUED (DISCONTINUED) {15E289FB-09DF-4731-9CF3-48A4E0031C31}
5. REPOSITION (REPOSITION) {DF3039CA-5977-4454-81D4-DB0AF966F85D}
6. OFF WHILE AMBULATORY (OFF AMBULATORY) {FF4610BF-8C92-408F-882A-1D6B2C2F146F}
7. ON WHILE SITTING IN CHAIR (ON IN CHAIR) {4139F4C3-BB68-4477-BF63-9BE7813AADA7}
8. OFF DURING BATHING (OFF BATHING) {5D0EFDD9-F170-4F34-9ABE-CE2D376E65F4}

UNITS:

=====

TERM: SEWS - AGITATION

Page: 623

ID: {012BB6B1-0CF8-4F51-832C-D85EF58C0831}

=====

Abbreviation..... SEWS AGITA  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Severity Of Ethanol Withdrawl Scale  
Display Name..... Agitation  
Help Text..... Agitation: fidgety, restless, pacing?  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

-----

1. 0 (0)	{2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
2. 3 (3)	{E714454D-6A67-91F6-2FB7-D8F24F184328}

UNITS:  
-----

=====

TERM: SEWS - ANXIETY

Page: 624

ID: {1F7184F3-45F0-4D93-808C-FB797589814F}

=====

Abbreviation..... SEWS ANXIE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Severity of Ethanol Withdrawl Scale  
Display Name..... Anxiety  
Help Text..... Do you feel that something bad is  
about to happen to you right now?  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

-----

1. 0 (0)	{2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
2. 3 (3)	{E714454D-6A67-91F6-2FB7-D8F24F184328}

UNITS:

-----









=====

TERM: SEWS - NAME ONLY  
ID: {7F010DD2-E6C5-45EF-800D-37C238ED69A1}

=====

Page: 628

Abbreviation..... SEWS N ONL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Severity Of Ethanol Withdrawl Scale  
Display Name..... Disorientation: Knows Name Only  
Help Text..... Disorientation: Knows Name Only  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- 1. 0 (0) {2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
  - 2. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}

UNITS:  
-----

=====

TERM: SEWS - NAUSEA

Page: 629

ID: {168E97DD-48FA-45BF-A3C0-8FCB5B316B6D}

=====

Abbreviation..... SEWS NAUSE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Severity Of Ethanol Withdrawl Scale  
Display Name..... Nausea and Dry Heaves or Vomitting?  
Help Text..... Nausea and Dry Heaves or Vomitting?  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

-----

1. 0 (0)	{2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
2. 3 (3)	{E714454D-6A67-91F6-2FB7-D8F24F184328}

UNITS:  
-----

ID: {4B7D9035-79F8-4736-9A20-87671D76E12B}

=====

Abbreviation..... SEWS SWEAT  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Severity Of Ethanol Withdrawl Scale  
Display Name..... Sweating  
Help Text..... Sweating: Includes moist palms,  
sweating now?  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

-----

- 1. 0 (0) {2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
- 2. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}

UNITS:

-----



=====

TERM: SEWS - TREMOR

Page: 632

ID: {8BC41692-A9A8-4900-BB2D-5E60E5625798}

=====

Abbreviation..... SEWS TREMO  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Severity Of Ethanol Withdrawl Scale  
Display Name..... Tremor  
Help Text..... Tremor: with arms extended, eyes  
closed?  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

-----

1. 0 (0)	{2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
2. 3 (3)	{E714454D-6A67-91F6-2FB7-D8F24F184328}

UNITS:

-----



ID: {657A9E0B-12CC-4D33-B4A8-99827CAE5CFF}

=====

Abbreviation..... SEWS VITAL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Severity Of Ethanol Withdrawl Scale  
Display Name..... Vitals: Pulse > 110, Diastolic BP >  
90, Temp > 99.6  
Help Text..... Vitals: Pulse > 110, Diastolic BP >  
90, Temp > 99.6  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

-----

1. 0 (0)	{2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
2. 3 (3)	{E714454D-6A67-91F6-2FB7-D8F24F184328}

UNITS:

-----

Abbreviation..... SD RAIL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Side Rails  
Help Text..... Side Rail Positions  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 14 entries

-----

1. ALL UP (ALL UP) {48199083-ACC2-408C-9242-7A6313CD7A7F}
2. ALL DOWN (ALL DOWN) {76C14E79-95E2-49B7-B847-3B114653E98C}
3. HEAD UP (hu) {BB1C79C0-80F1-7CDE-CC10-7CA0C704AE35}
4. R HEAD UP (R HEAD UP) {D03DBEC6-F4DD-441A-AAB2-653CC34FCDBF}
5. L HEAD UP (L HEAD UP) {BB39792D-B572-4A10-BA81-8685809E7437}
6. HEAD DOWN (hd) {80BF3425-D252-F917-052E-DC54D39E464E}
7. R HEAD DOWN (R HEAD DOW) {FDFE0638-444C-4023-B586-985E84EEA1FA}
8. L HEAD DOWN (L HEAD DOW) {3C919CB8-1C0A-498D-95A6-C46D25744CEC}
9. FOOT UP (fu) {49CCC59E-002E-8387-E8D6-D57AAC089EF3}
10. R FOOT UP (R FOOT UP) {5928B710-6EB1-4183-93CB-6C72B37A97C3}
11. L FOOT UP (L FOOT UP) {A15CD8C5-5CD1-4B90-B79C-5D3D331C7466}
12. FOOT DOWN (fd) {13F73BE9-2FA8-FE88-443B-9BA3FE6EDFD7}
13. R FOOT DOWN (R FOOT DOW) {DBCC0BA2-3097-42C7-AE9E-DF71FA779162}
14. L FOOT DOWN (L FOOT DOW) {5B18895B-DAD5-4874-87DB-DACFFF88B170}

UNITS:

-----

=====

TERM: SIGNIFICANT EVENT  
ID: {0DBEF78C-8CFC-0173-2A39-02C24EB10A5F}

=====

Page: 635

Abbreviation..... se  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... String  
Description..... No Description  
Display Name..... Significant Event  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: SINUS DRAINAGE CHARACTER  
ID: {05B77F59-77DE-424F-ACD1-413E70D52429}

=====

Page: 636

Abbreviation..... SIN DRA CH  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Sinus Drainage Character  
Help Text..... Describe The Sinus Drainage  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
1. SEROUS (SEROUS) {9BDA347E-29E9-80C7-91FA-016192D4257E}
  2. SEROSANGUINOUS (SEROSANGUI)  
{19A8D332-899D-177A-2247-ABCDAAEE3BFA8}
  3. CLEAR (c) {8981C365-2B46-704F-C9F4-975DF1C23C2D}
  4. CLOUDY (CLOUDY) {3CCFCCCC-00FC-6B00-1D58-A77C02FDB401}

UNITS:  
-----

=====

TERM: SINUS DRAINAGE PRESENT  
ID: {0F9CDF0D-A7A8-408D-BBFD-A6F74E5B8D6F}

=====

Page: 637

Abbreviation..... SIN DRAI P  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Sinus Drainage Present  
Help Text..... Is Sinus Drainage Present?  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEA649CD} |

UNITS:  
-----

=====

TERM: SINUS NECROTIC TISSUE APPEARANCE  
ID: {0578E21E-CFBE-440A-A532-8135AD414085}

=====

Page: 638

Abbreviation..... SNEC TI AP  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Sinus Necrotic Tissue Appearance  
Help Text..... Sinus Necrotic Tissue Appearance  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
1. SOFT (SOFT) {E68B4D4B-B86B-4004-A861-986645267F6E}
  2. YELLOW (YELLOW) {90DC31AD-6211-9F25-4237-79F7122A2763}
  3. HARD (HARD) {7934B543-7028-4CBC-9308-99F5A952A1DE}
  4. DRY (d) {D8F4F0FD-1B19-2515-3815-EC88D175B305}

UNITS:  
-----

=====

TERM: SINUS NECROTIC TISSUE PRESENT  
ID: {299921C9-5A13-4D9E-BC0D-8F777ED490D4}

=====

Page: 639

Abbreviation..... SIN NEC TI  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Sinus Necrotic Tissue Present  
Help Text..... Is Sinus Necrotic Tissue Present?  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

=====

TERM: SINUS TRACT SURROUNDING TISSUE CHARACTER  
ID: {59D298AC-1BE0-4C13-B75B-1CE331454151}

=====

Page: 640

Abbreviation..... SIN TI CHA  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Sinus Tract Surrounding Tissue  
Character  
Help Text..... Sinus Tract Surrounding Tissue  
Character  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Void.....

DATA VALUE: 4 entries

- 
1. ERYTHEMA (ERYTHEMA) {0083E1B3-AE2B-4CCF-8794-86E612B9DCEA}
  2. DISCOLORATION (DISCOLORAT)  
{BBE3E910-CC49-430C-BD90-FB0A97F1B183}
  3. INDURATION (INDURAITON) {A1669576-45D9-435C-87B9-488E9CF6BBA4}
  4. MACERATION (MACERATION) {B7C63E05-E193-40AA-A572-F1ACA0162B8B}

UNITS:  
-----



=====

TERM: SINUS TRACT ULCER  
ID: {D89C6750-E1F0-4817-9EDD-205EE310A09A}

=====

Page: 641

Abbreviation..... SIN TRA UL  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Sinus Tract Ulcer  
Help Text..... Is a Sinus Tract Ulcer Present?  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----



=====

TERM: SINUS TRACT ULCER GRANULIZATION  
ID: {64BE9437-2EAB-40F2-93BC-2D0C9DBDC1EB}

=====

Page: 643

Abbreviation..... SIN ULC GR  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Sinus Tract Ulcer Granulization  
Help Text..... Is Sinus Tract Ulcer Granulization  
Present?  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

=====

TERM: SKIN ALTERATION

Page: 644

ID: {6B4325F4-F776-48C1-902D-2CFFC73A22BB}

=====

Abbreviation..... SKIN ALTER  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Skin Alteration  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
  2. WOUND (WOUND) {1EC40149-051E-4750-8651-2B2D0DD2DEAC}

UNITS:  
-----

Abbreviation..... SKIN COLOR  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Observation of Patients Skin Color  
Display Name..... Skin Color  
Help Text..... Observation of Patients Skin Color  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 13 entries

-----

1. PALE (p) {0EB39D28-7ABE-A4AB-D334-3F795E65618D}
2. MOTTLED (m) {4B5FFB10-2CFF-2064-1553-4CBB7133B0C5}
3. ICTERIC (ICTERIC) {274ED5BB-5A34-20EE-A9B0-81AB7803780E}
4. JAUNDICE (j) {B17F4603-7C68-660D-96A7-119140A998D3}
5. FLUSHED (nt) {C7C54C1A-6E75-4BCF-4C57-0F8622A5BD16}
6. DUSKY (d) {2ABF845E-4E2F-883F-3FCC-F72B8507FF85}
7. CYANOTIC (c) {06BD77DB-1849-D731-0F27-F8CE48755EA9}
8. ASHEN (ASHEN) {7431E4EC-DCF1-4966-85CB-29BF65311855}
9. NORMAL FOR ETHNIC GROUP (ETH. NORM.)  
{3B2CD951-58CE-FF53-AD5D-6C65C07658D3}
10. PINK (PINK) {1A7FAF84-D69C-4EEB-7BC5-D21A598A1F4D}
11. BRUISED (BRUISED) {FA5B38A7-7744-186F-055A-8458D696AF33}
12. HEMATOMA (HEMATOMA) {6A24B7EE-62F3-C952-C1A3-086CDBB1FE68}
13. ECCHYMOSES (ECCHYMOSES) {DF076906-CC17-6AD1-76E7-B07678A126DE}

UNITS:

-----

=====

TERM: SKIN CONDITION

Page: 646

ID: {831BEB82-92CE-07A0-F36B-09C94030D088}

=====

Abbreviation..... SKIN COND.  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Skin Condition  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
1. INTACT (i) {935C9CAF-6C88-6610-C6E9-4C1B822923AF}
  2. BREAKDOWN (BREAKDOWN) {9CFD21C0-A805-47C4-E559-8033C754DE80}

UNITS:  
-----

=====

TERM: SKIN DRAINAGE AMOUNT  
ID: {B973CF56-27AD-447B-9B12-A4326997E027}

=====

Page: 647

Abbreviation..... SKI DRA AM  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Skin Device Drainage Amount  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
  2. SCANT (SCANT) {8EF6C62D-0D94-37BB-60AD-86E306E8120D}
  3. SMALL (s) {5991F892-B28C-3998-E733-D8FFA26D7E21}
  4. MODERATE (MODERATE) {C3C25064-C9FC-EA1F-3E74-32A25F0FC5E3}
  5. HEAVY (HEAVY) {1F10CD98-03FA-4705-BF7D-941E881C4E3D}

UNITS:  
-----

=====

TERM: SKIN ERYTHEMA EXTENDS AROUND INCISION  
ID: {A88D0DBE-7CC6-4983-B9AC-AF9CCF9653AB}

=====

Page: 648

Abbreviation..... ER EXI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Skin Erythema Extension Around Incision  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
1. EXTENDS AROUND INCISION NONE (ER EXI NON)  
{91F1A7C5-885C-4B13-9F7F-189C93E505CC}
  2. EXTENDS AROUND INCISION EXTENDS INTO PERIWOUND LESS THAN 1 CM (ER EXI <  
1) {095709D0-6190-49A3-AB08-BF4575DC8396}
  3. EXTENDS AROUND INCISION MORE THAN 1 CM (ER EXI >1)  
{11F13A1F-7AC2-4629-A960-CD92655EE101}
  4. LIMITED TO INCISION (ER EXI LIM)  
{C9ED5E2C-8ED1-445A-8A3E-C1207CE3843D}

UNITS:  
-----



=====

TERM: SKIN ERYTHEMA EXTENSION AROUND WOUND  
ID: {6BABC7F1-2C62-41A7-83E2-F5F4B5ACE087}

=====

Page: 649

Abbreviation..... ERY EX  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Skin Erythema Extension Around Wound  
Display Name..... Erythema Extension Around Wound  
Help Text..... Skin Erythema Extension Around Wound  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
- 1. EXTENSION AROUND WOUND NONE (ERY EX NON)  
{815C7BDF-C9A8-4AFE-B2BE-2C2FF24AA495}
  - ) 2. EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND MORE THAN 1 CM (ERY EX >1)  
{175B2505-C596-4036-8FC3-DFEA8EFB76F9}
  - ) 3. EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND LESS THAN 1 CM (ERY EX <1)  
{48C7FA69-D31A-4AAA-84F3-79AC863B71FD}
  - 4. EXTENSION AROUND WOUND LIMITED TO WOUND (ERY EX LIM)  
{2DC83B85-B3B6-4679-9F48-7DD26EFC4D4}

UNITS:  
-----

Abbreviation..... INCI CLOSU  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Skin Incision Closure  
Help Text..... Type of material used to close incision  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 8 entries

- 
1. SKIN INCISION WOUND CLOSURE SUTURES (INC CL SUT)  
{DED8DEDA-82C9-4D3A-A81E-641583314919}
  2. INCISION WOUND CLOSURE GLUE (IN CL GLU)  
{9EEB436E-2991-4E9A-9F8A-BFEC196784CB}
  3. SKIN INCISION WOUND CLOSURE STAPLES (INC CL STA)  
{4F48D42D-C206-423C-970B-02C03B2EE0B6}
  4. SKIN INCISION WOUND CLOSURE STERI-STRIPS (INC CL STE)  
{DF52B6F3-54EE-46B3-AFDF-4A701B065E04}
  5. INCISION WOUND CLOSURE NOT CLOSED (INC CL)  
{D63A05BF-1205-45F6-A78D-15D198BCB1C6}
  6. SKIN INCISION WOUND CLOSURE OPEN WOUND (INC CL OPE)  
{AD018B6F-13E5-41C0-BDBC-306CA0B1B475}
  7. INCISION WOUND CLOSURE APPROXIMATED (INC CL APP)  
{2DD5713E-D81C-484B-A0DF-6A5730B396DF}
  8. INCISION WOUND CLOSURE DEHISCENCE (INC CL DEH)  
{DCA1D137-0D9A-42D2-9270-7B3A82DCB6EB}

UNITS:  
-----

=====

TERM: SKIN INDURATION EXTENSION AROUND INCISION  
ID: {36580CA5-4FB0-4562-9652-0B3BD7D04307}

=====

Page: 651

Abbreviation..... INC IN  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Induration Extension Around Incision  
Display Name..... Induration Extension Around Incision  
Help Text..... Induration Extension Around Incision  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
1. SKIN INDURATION AROUND INCISION NONE (INC IN NON)  
{A9B95680-F9A0-44C0-929C-70C00B9AA5B7}
  2. SKIN INDURATION AROUND INCISION EXTENDS INTO PERIWOUND LESS THAN 1 cm (  
INC IN <1) {2A5097BF-955C-4F79-B483-7EE14AF05786}
  3. SKIN INDURATION AROUND INCISION EXTENDS INTO PERIWOUND MORE THAN 1 cm (  
INC IN >1) {D8EE448D-00EB-412D-95AC-2E2E92C53EEB}
  4. SKIN INDURATION AROUND INCISION LIMITED TO INCISION (INC IN LIM)  
{88BCE7F9-F21F-427A-A50B-C1BA831639AB}

UNITS:  
-----

Abbreviation..... WOU IN  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Skin Induration Extension Around Wound  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
1. SKIN INDURATION EXTENSION AROUND WOUND NONE (WOU IN NON)  
{2C57CD27-4FDD-4DE8-88A0-D438CBE535D8}
  2. SKIN INDURATION EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND < 1 cm (W  
OU IN <1)  
{3594CCCD-FA51-4971-A61E-FFA3B4AC7EE0}
  3. SKIN INDURATION EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND > 1 cm (W  
OU IN >1)  
{7DF7D08D-3E5E-4D18-94E4-A0DFDA5DB72F}
  4. SKIN INDURATION EXTENSION AROUND WOUND LIMITED TO WOUND (WOU IN LIM)  
{1828A9EA-FF9D-4A49-9543-15B85FBF211E}

UNITS:  
-----

=====

TERM: SKIN MOISTURE

Page: 653

ID: {BE934424-A7CF-EEF0-8098-26B3AFC4B2FE}

=====

Abbreviation..... SKIN MOIST  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Skin Moisture  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
- |                    |  |
|--------------------|--|
| 1. DRY (d)         | {D8F4F0FD-1B19-2515-3815-EC88D175B305} |
| 2. MOIST (m)       | {AA73A9D6-6645-5DA7-C602-A537A2DFBFA0} |
| 3. DIAPHORETIC (d) | {F5F1E644-A1BD-0137-54B5-1C02678C1D5A} |
| 4. NORMAL (n)      | {553C9A97-A3D5-C52F-0E4E-C7D168C99A2B} |

UNITS:

-----

=====

TERM: SKIN PATCHES PRESENT  
ID: {D48AF8EF-A439-4D1B-9987-A04D9E030B45}

=====

Page: 654

Abbreviation..... SKIN PAT P  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Are Skin Patches Present  
Display Name..... Skin Patches Present  
Help Text..... Are Skin Patches Present? (EKG,  
Medication...)  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEA649CD} |

UNITS:  
-----

Abbreviation..... SK RISK PR  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Are Skin Risk Factors Present?  
Display Name..... Skin Risk Factors Present  
Help Text..... Does The Patient Have Spinal Cord  
                                  Injury, Paralysis, Or Neurological  
                                  Disease?  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Void.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

=====

TERM: SKIN SWELLING EXTENSION AROUND INCISION  
ID: {AB9716EF-235C-4434-A96D-5BC289C9458A}

=====

Page: 656

Abbreviation..... INC SW  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Skin Swelling Extension Around Incision  
Display Name..... Skin Swelling Extension Around Incision  
Help Text..... Skin Swelling Extension Around Incision  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. SKIN SWELLING EXTENSION AROUND INCISION NONE (INC SW NON)  
{54C9FC46-899C-460B-9524-D6F196B9FF31}
  2. SKIN SWELLING EXTENSION AROUND INCISION EXTENDS INTO PERIWOUND < 1 cm (  
INC SW <1) {ABA0527C-8FAD-4E34-9F6F-2AA97EA06E20}
  3. SKIN SWELLING EXTENSION AROUND INCISION EXTENDS INTO PERIWOUND > 1 cm (  
INC SW >1) {106A2AD3-F482-4DEE-967E-6F06C6026944}
  4. SKIN SWELLING EXTENSION AROUND INCISION LIMITED TO WOUND (INC SW LIM)  
{B8503D54-6925-430E-8AB6-692412347D5D}
  5. SKIN SWELLING EXTENSION AROUND INCISION FLUCTUANT (INC SW FLU)  
{5D38DC7A-8E40-4D27-A039-2D6A557B4FFD}

UNITS:  
-----



=====

TERM: SKIN SWELLING EXTENSION AROUND WOUND  
ID: {FE9F8966-B828-473C-9285-E189A448E304}

=====

Page: 657

Abbreviation..... WOU SW  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Skin Swelling Extension Around Wound  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. SKIN SWELLING EXTENSION AROUND WOUND NONE (WOU SW NON)  
{550E3B92-2C69-403B-A997-2DBE61374ABF}
  2. SKIN SWELLING EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND < 1 cm (WOU SW <1)  
{15F314A0-6B77-4EB5-92C0-FFEC4D82C75C}
  3. SKIN SWELLING EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND > 1 cm (WOU SW >1)  
{EA6AC933-0737-4C28-9F31-D36F3452DA96}
  4. SKIN SWELLING EXTENSION AROUND WOUND LIMITED TO WOUND (WOU SW LIM)  
{413CD1EC-309C-45D0-A1E0-E829EA615E2A}
  5. SKIN SWELLING EXTENSION AROUND WOUND FLUCTANT (WOU SW FLU)  
{DB5077E7-D339-4C87-B8CA-D16269DB2935}

UNITS:  
-----

Abbreviation..... SKIN TEMP  
 Active..... Yes  
 Boolean Value False.....  
 Boolean Value True.....  
 Data Type..... PickList  
 Description..... Skin Temperature  
 Display Name..... Skin Temperature  
 Help Text..... Skin Temperature  
 Multi Select Picklist..... No  
 Term Type..... OBSERVATION  
 Vuid.....

DATA VALUE: 9 entries

- 1. COOL (COOL) {48A1FEA8-FD75-DAF8-174A-9C6AB2044439}
- 2. WARM (WARM) {E03BBEB4-5643-968F-190F-309A1293FCF9}
- 3. COLD (COLD) {B7B422BA-22EC-117D-0A0C-8419E9447D46}
- 4. HOT (HOT) {401EA6CD-12DE-AB0E-2B52-8B2E62F7C34E}
- 5. CONSISTENT WITH BODY TEMPERATURE (CON W TEMP)  
{941FC5EF-6AF8-4DA1-A43C-95575AC74B57}
- 6. WARMER THAN GENERAL BODY TEMP (WARM B TEM)  
{AE43F2B4-B585-43F4-9E97-D61E225A803B}
- 7. WARMER THAN OPPOSITE BODY PART (WARM OPP P)  
{5DBF8614-5067-49C5-BD02-C1B21C9DFECC}
- 8. COOLER THAN OPPOSITE BODY PART (COOL OPP P)  
{AFA2635C-8E80-4C71-A1B4-4AB117096FD2}
- 9. COOLER THAN GENERAL BODY TEMPERATURE (CoolGenT)  
{03D7D4CF-9B3C-49CA-9650-213E9BABE064}

LOCATION: 49 entries

- 1. HEAD (Hd) {C3209F44-90C6-FF1B-3169-265B294F276A}
- 2. FACE (FACE) {C638EA14-7903-6A94-496A-C6560D51BEE3}
- 3. NECK (Ne) {D4812C3B-AAC2-1B61-679B-3CD8B5CB947C}
- 4. ABDOMEN (ABDOMEN) {109FE041-2521-6EF2-7F95-1FD1587796AB}
- 5. BACK (BACK) {22C26239-C692-4E93-B522-E5E8392BA3B1}
- 6. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
- 7. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
- 8. L HAND (L HAND) {639DF25A-19E3-B504-D323-4B3F6ED39A24}
- 9. R HAND (R HAND) {A54D3A80-8537-D079-9CCE-7DC56346D132}
- 10. L HIP (L HIP) {786CD881-27AD-4136-B9CC-8423CB38A431}
- 11. R HIP (R HIP) {A48CBCAE-5404-4934-9321-ECEF660A5350}
- 12. L GROIN (L GROIN) {5488E4EE-C1A2-4B02-8763-A53C412117A7}
- 13. R GROIN (R GROIN) {BAE1C889-8CCF-4CEB-9EBD-0243BC6F096A}
- 14. L THIGH (L THIGH) {16C2312F-1E0E-45E0-B675-160F40C23F57}
- 15. R THIGH (R THIGH) {B7A2D4A0-801E-4578-967A-F1CFB86FB96B}
- 16. L KNEE (L KNEE) {3A9E30B4-09A9-967A-58EF-69A01EC4F859}
- 17. R KNEE (R KNEE) {BCD18D53-07B6-2C41-64FC-3A07A7E979E7}
- 18. L FOOT (L FOOT) {2EC3FBEC-4EA2-9E7A-17D4-FDC2F9D28046}
- 19. R FOOT (R FOOT) {7D9F5FF6-232C-2488-4504-A721FF06928F}
- 20. L LEG (LL) {B088C809-E528-B785-21C2-14E45E1FD8AF}
- 21. R LEG (RL) {C5E40A7B-6FC4-267D-8CE8-3571D21D81FD}
- 22. L EAR (L EAR) {DA9C2DD1-D6FF-4899-A62B-1096FAA36C59}
- 23. L SHOULDER (L SHOULDER) {91CBD2FA-2743-42EB-946B-6933980D9337}
- 24. R SHOULDER (R SHOULDER) {D27AE130-442E-4DD6-B391-8D7E8AD494A2}
- 25. L ELBOW (L ELBOW) {BF4317E7-2024-47EF-9E5F-3D21D0547885}
- 26. R ELBOW (R ELBOW) {5BAFF87B-6993-4E76-862F-9F386D9B71C5}
- 27. L ANKLE (L ANKLE) {D5A692EF-9D7C-473C-AF96-7EF1C01FA2AF}

- 28. R ANKLE (R ANKLE) {4E6DF352-F399-4CE9-908C-48132C3D856B}
- 29. COCCYX (COCCYX) {49D65E9D-3D5B-40E3-B7CF-1003F80B852C}
- 30. L JAW (L JAW) {65EBD3CF-76D8-4F45-9B1B-3575D988BE8E}
- 31. R JAW (R JAW) {7846AEE2-B6F4-4B5A-909F-56E386FD82A8}
- 32. L NOSE (L NOSE) {C7707847-5135-4749-B196-002E3E59BB54}
- 33. R NOSE (R NOSE) {D80BCB19-E71F-43EA-8990-6C121EE7F7A1}
- 34. NOSE (NOSE) {85F7003E-6D4E-4A79-A534-8A57FBC2D106}
- 35. PENIS (PENIS) {2D82CE55-83B6-4160-A754-07C71FD388E8}
- 36. ABDOMEN LUQ (ABD LUQ) {3B31E1EC-CE74-4C76-BB61-10C678431FCA}
- 37. ABDOMEN RUQ (ABD RUQ) {F29C9FA1-717A-49B7-8000-95070E9AACAO}
- 38. ABDOMEN LLQ (ABD LLQ) {39831F59-D7B0-4156-AC57-4C622A93CA58}
- 39. ABDOMEN RLQ (ABD RLQ) {2C48394F-6877-4738-856C-E933644D8067}
- 40. L CHEST (L CHEST) {E332E497-FE0C-47CB-A0A7-217FBF52DC97}
- 41. R CHEST (R CHEST) {7E30EC0E-815A-4C2E-A9FC-3C2B2D6E3CD4}
- 42. CHEST (CHEST) {8C3C30FE-6E06-C35E-8226-CBF0C94014CC}
- 43. L SCAPULA (L SCAPULA) {3381EC50-771B-41C8-85C8-6EB0D208D54E}
- 44. R SCAPULA (R SCAPULA) {EA657261-4D72-40C2-B70D-2A2E27556A01}
- 45. L UPPER ARM (L UPP ARM) {5F915696-1B08-32E9-6B95-49246117D660}
- 46. R UPPER ARM (R UPP ARM) {3EF20E8E-F68C-E4F2-8F8C-FAD3B44DDAA4}
- 47. L LOWER ARM (L LOW ARM) {63AA95E8-256D-4390-8CB3-F00A1B3F91CD}
- 48. R LOWER ARM (R LOW ARM) {8DE6E277-EFFF-40DD-9E83-D7295AACB6AA}
- 49. SCROTUM (SCROTUM) {BDC109F5-0D30-4914-A62C-CEFEB448B547}

UNITS:

-----

Abbreviation..... INC TP  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Skin Temperature at Incision Site  
Display Name..... Skin Temperature at Incision Site  
Help Text..... Skin Temperature at Incision Site  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

-----

1. SKIN TEMPERATURE AT INCISION SITE CONSISTENT WITH BODY TEMPERATURE (INC TP CON) {D349175A-4B20-4243-BE40-E01E50C5583F}
2. SKIN TEMPERATURE AT INCISION SITE WARMER THAN GENERAL BODY TEMPERATURE (WARM BD OP) {5F888E9D-8F65-49ED-BC5D-F6F287A9D478}
3. SKIN TEMPERATURE AT INCISION SITE WARMER THAN OPPOSITE BODY PART (WARME R OP) {F518787A-325A-45F3-82B4-73C856802073}
4. SKIN TEMPERATURE AT INCISION SITE COOLER THAN GENERAL BODY TEMPERATURE (INC TP CBT) {D357C975-857D-4293-B3EF-EFA7A4586357}
5. SKIN TEMPERATURE AT INCISION SITE COOLER THAN OPPOSITE BODY PART (INC T P COP) {85307AD8-F730-43F0-BCB3-882594D9DDF8}

UNITS:

-----

Abbreviation..... WOU TP  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Skin Temperature at Wound Site  
Display Name..... Skin Temperature at Wound Site  
Help Text..... Skin Temperature at Wound Site  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

-----

1. SKIN TEMPERATURE AT WOUND SITE CONSISTENT WITH BODY TEMPERATURE (WOU TP CBT) {672D56C1-D1A3-4DDF-8815-F05DEE197061}
2. SKIN TEMPERATURE AT WOUND SITE WARMER THAN GENERAL BODY TEMPERATURE (WOU TP WGB) {03BB8F9A-3DAC-4994-9F09-AA1DEEF48090}
3. SKIN TEMPERATURE AT WOUND SITE WARMER THAN OPPOSITE BODY PART (WOU TP WOB) {6290BD23-24BF-4A3A-8555-DCE52D2274B3}
4. SKIN TEMPERATURE AT WOUND SITE COOLER THAN OPPOSITE BODY PART (WOU TP COB) {A904C666-2810-4211-B53A-2961DCE6894E}
5. SKIN TEMPERATURE AT WOUND SITE COOLER THAN GENERAL BODY TEMPERATURE (WOU TP CGB) {7203993E-3334-44EC-9FD9-14224253283D}

UNITS:

-----

=====

TERM: SKIN TURGOR

Page: 661

ID: {4D90E9F2-5E61-67AA-1ACA-874AED4E6C79}

=====

Abbreviation..... SKIN TURG  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Skin Turgor  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. GOOD (GOOD) {39ECA304-49F9-44B7-9735-B33CEC9D8119}
  2. DECREASED (DECREASED) {D910FE54-A6AC-4BE3-AE35-C9E6679A93BD}
  3. POOR (POOR) {A3D75A12-A791-446D-AF63-9435054DFC4C}

UNITS:  
-----

Abbreviation..... WOU CLOSUR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Wound Closure
Help Text..... Type of material used to close wound
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 8 entries

- 1. SKIN WOUND CLOSURE APPROXIMATED (WOU CL APP)
{83D48E0F-8643-436B-9ECC-80795857C322}
2. SKIN WOUND CLOSURE DEHISCENCE (WOU CL DEH)
{425A460F-19C0-4A4A-8E4E-712C66968C38}
3. SKIN WOUND CLOSURE GLUE (WOU CL GLU)
{29FE4AF7-39EE-400A-BB3D-E07786F11BFD}
4. SKIN WOUND CLOSURE NOT CLOSED (WOU CL NOT)
{0E69F1FA-35F6-4F34-8547-CC3602670A7B}
5. SKIN WOUND CLOSURE OPEN WOUND (WOU CL OPE)
{C6F2AD2C-932F-434A-8A7A-513029092FFD}
6. SKIN WOUND CLOSURE STAPLES (WOU CL STA)
{FE04BFE3-DF81-48E7-9FA1-1CF4ECAD71}
7. SKIN WOUND CLOSURE STERI-STRIPS (WOU CL STE)
{AD451791-6B49-4A51-ABC3-796D08BEF910}
8. SKIN WOUND CLOSURE SUTURES (WOU CL SUT)
{F2840AE4-216B-477C-9FE6-DC91E8CE5364}

UNITS:

=====

TERM: SLEEP AMOUNT

Page: 663

ID: {E68F85F4-35A1-4FAC-BA86-90BF1FD5D63F}

=====

Abbreviation..... SLEEP AMOU  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Sleep Amount  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

HOUR                    -0.1 > 0.0 \* 24.0 < 24.1    Dec Prec: 1



ID: {7727D9C8-C3F2-4676-A2B8-23F4DAC72B90}

=====

Abbreviation..... SLEEP QUAL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Sleep Qualifiers  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 6 entries

-----

1. NAPS ONLY (NAPS ONLY) {CF749F5C-2E4B-419A-B020-C81B30A7E482}
2. HAS NOT SLEPT (NO SLEEP) {3BF9DB54-5D1E-422D-90AF-918F81C9EA85}
3. ASLEEP AT ROUNDS (SLEEP ROUN)  
{24E1D9F0-D8F1-4DC8-911F-6170F343C087}
4. DENIES HAVING SLEPT (DENIES SL)  
{1AC4B0CF-EC41-426B-89AB-68D5A9BAFADA}
5. SLEPT WELL (SLEPT WELL) {18573E82-568F-4D2E-B9E1-6E4589EBEE2D}
6. INTERMITTENT SLEEP (INT SLEEP)  
{5005C433-5892-4F4D-A64D-4C94AD1853E2}

UNITS:

-----

=====

TERM: SPECIALTY BED/MATTRESS  
ID: {274B7CFD-D091-47F2-9F35-602F14D75340}

=====

Page: 665

Abbreviation..... SPEC BED/M  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Type of specialty Bed or Mattress  
Display Name..... Specialty Bed/Mattress  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 9 entries

- 
1. BARIATRIC FRAME (BARI FRAME)  
{CF95B426-F280-451B-86F2-CEF7A1A2B967}
  2. BARIATRIC LOW AIRLOSS MATTRESS REPLACEMENT (BAR L AIR)  
{220C6E92-1D0A-4513-B3A9-8065A6EC7DB1}
  3. LOW AIRLOSS/ALTERNATING PRESSURE COMBO (L AIRL ALT)  
{ADDD1260-C95B-4D7B-AB89-3749D76A3700}
  4. LOW AIRLOSS OVERLAY (L AIR OVRL)  
{282C439C-A697-43E6-85C1-786352EB56BC}
  5. PERCUSSION BED (PERC BED) {BE0B086B-A2F3-446A-ADE0-9589E92689FE}
  6. ROTATIONAL THERAPY (ROT THERAP)  
{BB795594-720C-40D7-9B7C-6792BC849DA0}
  7. STANDARD PRESSURE REDUCTION (STAN P RED)  
{556537CD-A84D-4493-81F6-51CBD844B5B7}
  8. STRYKER FRAME (STRYKER FR)  
{71A1DB72-C860-42D8-BA4F-B71AFA2164B1}
  9. FLOTATION BED (FLOT BED) {FA6329AB-85A6-4301-BC14-D0D089249356}

UNITS:  
-----

=====

TERM: SPECIFIC GRAVITY  
ID: {CFB61154-54BD-4DE2-B83B-BC497AEC0F33}

=====

Page: 666

Abbreviation..... SPEC GRAVI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Specific Gravity  
Display Name..... Specific Gravity  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

PRODUCT: 1 entries

-----

1. URINE (u) {1F61B077-1D34-0CBA-E55B-65A8B6A7CAE7}

UNITS:

-----

SI UNIT 0.500 > 1.000 \* 1.030 < 2.000  
Dec Prec: 3

=====

TERM: ST SEGMENT AVF

Page: 667

ID: {F5ECC376-8003-62AF-64E7-08D52B244095}

=====

Abbreviation..... STAVF  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... ST Segment Average Vector Foot  
Display Name..... ST Segment AVF  
Help Text..... The part of an electrocardiogram  
immediately following the QRS complex  
and merging into the T wave  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Void.....

METHOD: 1 entries

-----

1. ELECTROCARDIOGRAM (ECG) {CF5F1AB9-3B74-7FBA-E38B-C9EF54FE4DC5}

UNITS:

-----

MILLIVOLTS -0.50 > -0.05 \* 0.10 < 1.00  
Dec Prec: 2  
MILLIMETERS -5.00 > -0.50 \* 1.00 < 10.00  
Dec Prec: 2

Abbreviation..... STAVL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... ST Segment Average Vector Left  
Display Name..... ST Segment AVL  
Help Text..... The part of an electrocardiogram  
immediately following the QRS complex  
and merging into the T wave  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Void.....

METHOD: 1 entries

-----

1. ELECTROCARDIOGRAM (ECG) {CF5F1AB9-3B74-7FBA-E38B-C9EF54FE4DC5}

UNITS:

-----

MILLIVOLTS	-0.50 > -0.05 * 0.10 < 1.00 Dec Prec: 2
MILLIMETERS	-5.00 > -0.50 * 1.00 < 10.00 Dec Prec: 2



Abbreviation..... ST1  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... ST Segment 1  
Display Name..... ST Segment 1  
Help Text..... The part of an electrocardiogram  
immediately following the QRS complex  
and merging into the T wave  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----  
1. ELECTROCARDIOGRAM (ECG) {CF5F1AB9-3B74-7FBA-E38B-C9EF54FE4DC5}

UNITS:

-----  
MILLIVOLTS -0.50 > -0.05 \* 0.10 < 1.00  
Dec Prec: 2  
MILLIMETERS -5.00 > -0.50 \* 1.00 < 10.00  
Dec Prec: 2

=====

TERM: ST SEGMENT II  
ID: {6F42CE65-5097-918C-93B6-7D3000CA402E}

=====

Page: 671

Abbreviation..... ST2  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... ST Segment 2  
Display Name..... ST Segment 2  
Help Text..... The part of an electrocardiogram  
immediately following the QRS complex  
and merging into the T wave  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Void.....

METHOD: 1 entries

-----

1. ELECTROCARDIOGRAM (ECG) {CF5F1AB9-3B74-7FBA-E38B-C9EF54FE4DC5}

UNITS:

-----

MILLIVOLTS -0.50 > -0.05 \* 0.10 < 1.00  
Dec Prec: 2  
MILLIMETERS -5.00 > -0.50 \* 1.00 < 10.00  
Dec Prec: 2



=====

TERM: ST SEGMENT III  
ID: {C2D177EB-A435-9EAA-E518-2039E8592E69}

=====

Page: 672

Abbreviation..... ST 3  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... ST Segment 3  
Display Name..... ST Segment 3  
Help Text..... The part of an electrocardiogram  
                  immediately following the QRS complex  
                  and merging into the T wave  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Void.....

METHOD: 1 entries

-----

1. ELECTROCARDIOGRAM (ECG) {CF5F1AB9-3B74-7FBA-E38B-C9EF54FE4DC5}

UNITS:

-----

MILLIVOLTS	-0.50 > -0.05 * 0.10 < 1.00
	Dec Prec: 2
MILLIMETERS	-5.00 > -0.50 * 1.00 < 10.00
	Dec Prec: 2



Abbreviation..... STV2  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... ST Segment from Lead V2  
Display Name..... ST Segment V2  
Help Text..... The part of an electrocardiogram  
immediately following the QRS complex  
and merging into the T wave from Lead  
V2  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----

1. ELECTROCARDIOGRAM (ECG) {CF5F1AB9-3B74-7FBA-E38B-C9EF54FE4DC5}

UNITS:

-----

MILLIVOLTS -0.50 > -0.05 \* 0.10 < 1.00  
Dec Prec: 2  
MILLIMETERS -5.00 > -0.50 \* 1.00 < 10.00  
Dec Prec: 2



=====  
TERM: ST SEGMENT V4

Page: 676

ID: {A48D67C4-49D9-4EFD-DF99-7088D8E2DFBA}

=====  
Abbreviation..... STV4  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... ST Segment from Lead 4  
Display Name..... ST Segment V4  
Help Text..... The part of an electrocardiogram  
immediately following the QRS complex  
and merging into the T wave from Lead  
V4  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----  
1. ELECTROCARDIOGRAM (ECG) {CF5F1AB9-3B74-7FBA-E38B-C9EF54FE4DC5}

UNITS:

-----  
MILLIVOLTS -0.50 > -0.05 \* 0.10 < 1.00  
Dec Prec: 2  
MILLIMETERS -5.00 > -0.50 \* 1.00 < 10.00  
Dec Prec: 2

=====

TERM: ST SEGMENT V5  
ID: {0B6F3CF6-9DDB-17EC-EFC0-37610387E274}

=====

Page: 677

Abbreviation..... STV5  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... ST Segment from Lead 5  
Display Name..... ST Segment V5  
Help Text..... The part of an electrocardiogram  
immediately following the QRS complex  
and merging into the T wave from Lead  
V5  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----  
1. ELECTROCARDIOGRAM (ECG) {CF5F1AB9-3B74-7FBA-E38B-C9EF54FE4DC5}

UNITS:

-----  
MILLIVOLTS -0.50 > -0.05 \* 0.10 < 1.00  
Dec Prec: 2  
MILLIMETERS -5.00 > -0.50 \* 1.00 < 10.00  
Dec Prec: 2

Abbreviation..... STV6
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... ST Segment from Lead 6
Display Name..... ST Segment V6
Help Text..... The part of an electrocardiogram
immediately following the QRS complex
and merging into the T wave from Lead
V6
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

- 1. ELECTROCARDIOGRAM (ECG) {CF5F1AB9-3B74-7FBA-E38B-C9EF54FE4DC5}

UNITS:

MILLIVOLTS -0.50 > -0.05 \* 0.10 < 1.00
Dec Prec: 2
MILLIMETERS -5.00 > -0.50 \* 1.00 < 10.00
Dec Prec: 2

=====

TERM: ST/T WAVE MEASUREMENT  
ID: {F825C20C-82DE-4DD6-BC78-4E25030D0DB1}

=====

Page: 679

Abbreviation..... ST/T measu  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description.....  
Display Name..... ST/T Wave Measurement  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

MILLIMETERS            -9.0 > -9.0 \* 9.0 < 9.0    Dec Prec: 1



Abbreviation..... STOOL CHAR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Stool Characteristics
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 10 entries

- 1. NORMAL STOOL (NORM STOOL) {486B0693-CAC5-48EE-B2E6-1F73C651418A}
2. LOOSE (LOOSE) {1346A74E-3AA1-4D81-876B-91201F9717A5}
3. SOFT (SOFT) {E68B4D4B-B86B-4004-A861-986645267F6E}
4. DIARRHEA (DIARRHEA) {353472EF-3499-46B1-AD7C-554FE01D8CE2}
5. CONSTIPATION (CONSTIPAT.) {4A3B0D47-88BF-4FEB-AA0D-7AD7FED2C5C3}
6. TARRY (TARRY) {09364271-6307-471F-A334-A8F7E051EF46}
7. MUCOID (MUCOID) {073E5514-DACE-4F43-ADE5-91A1A87D6601}
8. BLOODY (BLOODY) {D7904E84-8331-D6DD-F2C0-7F3039C5746A}
9. FRANK BLOOD (FRANK BLD) {CFAE6E6B-87B2-4CBA-A280-C5233B7995ED}
10. CLOTS (CLOTS) {B52B88E0-6D23-4657-97CD-81E3EF64991B}

UNITS:

=====

TERM: STOOL COLOR

Page: 681

ID: {85012C5D-09F3-4078-9A4D-82564BF57235}

=====

Abbreviation..... STOOL COLO  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Stool Color  
Display Name..... Stool Color  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 9 entries

- 
1. BROWN (BROWN) {A1C32FBC-7B4D-707C-CA58-92A2E201368E}
  2. YELLOW/GREEN (YELLOW/GRN) {1D0E7E81-3331-4D07-9B94-3D4BF72A0F0D}
  3. BLACK (BLACK) {7A55AFB9-6861-4918-0746-3D67D47FE5C3}
  4. RED (RED) {515D7CFB-3527-30E7-EC4B-B607C272DFA0}
  5. MAROON (MAROON) {1F518B23-C79F-4F06-9DF3-ADC5C96D3029}
  6. CLEAR (c) {8981C365-2B46-704F-C9F4-975DF1C23C2D}
  7. GREEN (GREEN) {3D1F1553-56FC-68CE-79A9-E4F6357B6AAF}
  8. RASH/EXFOLIATION (RASH/EXFOL)  
{AB83ABA9-28A7-6484-67E7-4E66B5526942}
  9. YELLOW (YELLOW) {90DC31AD-6211-9F25-4237-79F7122A2763}

UNITS:

-----

Abbreviation..... NUM STOOLS
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Number of stools
Display Name..... Number of stools
Help Text..... Number of stools
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 11 entries

- 1. 1 (1) {8F5D10A8-FE17-47BC-DACA-6972D12CC00A}
2. 2 (2) {38E22B79-A2B6-3401-B195-C955F24E1278}
3. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}
4. 4 (4) {D2280C7E-DE50-5A8D-C3FD-5862B69C0E1A}
5. 5 (5) {188EDCE7-0726-4234-D8D2-92C8052CD221}
6. 6 (6) {640CF2F9-3BBF-AB13-185B-CD02A609FB5A}
7. 7 (7) {7A010324-D2C2-5D4B-848A-7B8878666D20}
8. 8 (8) {518E7CFC-102E-351E-C7AF-FB575E3C8F4D}
9. 9 (9) {58FF2FD4-569F-F454-9F5B-FE6F9E34D8CC}
10. 10 (10) {0C827160-18DB-816F-7FAA-A3C175C6D33B}
11. N/A (n) {B3DBE3FD-C5CB-3894-607A-BE13319594D0}

UNITS:



=====

TERM: STROKE VOLUME

Page: 684

ID: {97938D03-0837-552B-5BC2-3F8D1A38E386}

=====

Abbreviation..... SV  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Stroke Volume  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----

1. CARDIAC OUTPUT MONITOR (COM)  
{CCBE5744-7CB9-9DFA-938F-FDDBC73B45B3}

UNITS:

-----

MILLILITERS PER BEAT 30 > 60 \* 130 < 200 Dec Prec: 0

=====

TERM: STROKE VOLUME INDEX  
ID: {0FFED66D-43F3-E4CA-18B7-DBFCAA2F5EC5}

=====

Page: 685

Abbreviation..... SVI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Stroke Volume Index  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----

1. CARDIAC OUTPUT MONITOR (COM)  
{CCBE5744-7CB9-9DFA-938F-FDDBC73B45B3}

UNITS:

-----

MILLILITERS PER BEAT PER METER SQUARE (BSA)  
20 > 41 \* 51 < 100 Dec Prec: 0



Abbreviation..... O2FLOWRATE
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Supplemental Oxygen Flow Rate
Display Name..... Supplemental Oxygen Flow Rate
Help Text..... Supplemental Oxygen Flow Rate
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 12 entries

- 1. AMBU BAG (AB) {E9D82DB9-D8A3-10E4-CE51-34BFF48F4603}
2. BIPAP (BIP) {3CFD1D11-EF8D-58BA-6381-B24956344D2E}
3. BIPAP-CPAP (PAP) {2BCBE562-8F44-662F-469B-13AC0E863A3D}
4. CPAP (CP) {A041FB71-08EF-F3FE-82FA-051FF74C824F}
5. MASK (M) {4FB6B24A-D56E-4EE3-6CBC-E54C9734424A}
6. NASAL CANNULA (NC) {EFFD38AE-6AF5-77FD-72E1-E0F7E3DEBA29}
7. NEBULIZER (NEBULIZER) {3CF184A3-2688-4988-B956-D5B0F303E8E8}
8. NON RE-BREATHING (NRB) {1B9AA055-FD7D-7DCC-4EDB-BDE3C821AA8E}
9. PARTIAL RE-BREATHING (PRB) {6F1DF170-5128-2E3F-1831-12FB0BB5E95C}
10. PENDANT CANNULA (PC) {A1451228-4ABB-DBCB-31AD-578BCD93012F}
11. RESERVOIR CANNULA (RC) {B332C44D-29BE-765C-932C-63EBFC9631FC}
12. VENTURI MASK (VM) {E0401675-9C4E-2AFA-2B93-F352B781E963}

UNITS:

LITERS PER MINUTE -1 > 0 \* 20 < 21 Dec Prec: 0



=====

TERM: SYSTEMIC VASCULAR RESISTANCE  
ID: {BD2C3F43-186B-A312-2514-186EFE223B65}

=====

Page: 688

Abbreviation..... SVR  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Systemic Vascular Resistance  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
DYNE \* SECONDS / CM(-5) \* M2 400 > 800 \* 1200 < 3000 Dec Prec: 0

=====

TERM: SYSTEMIC VASCULAR RESISTANCE INDEX  
ID: {62C5A9F1-C0D2-DD0C-5A01-459C833375EE}

=====

Page: 689

Abbreviation..... SVRI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Systemic Vascular Resistance Index  
Display Name..... Systemic Vascular Resistance Index  
Help Text..... Systemic Vascular Resistance Index  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----  
1. CARDIAC OUTPUT MONITOR (COM)  
    {CCBE5744-7CB9-9DFA-938F-FDDBC73B45B3}

UNITS:

-----  
DYNE \* SECONDS / CM(-5) / M2 1500 > 1760 \* 2600 < 4000 Dec Prec: 0

Abbreviation..... SYSTOL PRE
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Systolic Pressure
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

LOCATION: 4 entries

- 1. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
2. L LEG (LL) {B088C809-E528-B785-21C2-14E45E1FD8AF}
3. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
4. R LEG (RL) {C5E40A7B-6FC4-267D-8CE8-3571D21D81FD}

METHOD: 5 entries

- 1. CUFF (Cu) {E59A2543-1860-1E1C-BEF1-6E7F6B6FBC5E}
2. DOPPLER (Dop) {52683C19-F9DC-876B-CA3F-53E7C864CF04}
3. ARTERIAL LINE (AL) {D4B3F4ED-8B33-431E-4164-BD165E823776}
4. INVASIVE (I) {5DB8746E-5446-F47A-883C-A44B91275A6E}
5. NON-INVASIVE (N-I) {F50ABFED-5FD7-FC95-20FE-1E3038C9F299}

POSITION: 3 entries

- 1. LYING (Ly) {DE28D2BF-6F8D-64C3-0092-4CEB5E0902B9}
2. SITTING (Sit) {DD4DCDF1-0A22-2CC5-404F-0E472C73961C}
3. STANDING (St) {0776AE42-C370-B7B3-9B4D-9C914BA8F231}

UNITS:

MILLIMETERS OF MERCURY 0 > 80 \* 150 < 300 Dec Prec: 0

=====

TERM: SaO2% Page: 691  
ID: {D2D98D5C-552B-9913-31FB-E4FE5BB32912}

=====

Abbreviation..... SaO2  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Arterial Oxygen Saturation  
Display Name..... SaO2%  
Help Text..... Arterial Oxygen Saturation  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

PERCENTAGE                    50 > 92 \* 100 < 101      Dec Prec: 0

Abbreviation..... SpO2  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Spot Oxygen Saturation  
Display Name..... SpO2%  
Help Text..... Oxygen Saturation  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid..... 4500637

METHOD: 2 entries

- 
- 13. MONITOR (Mo) {63FC24EA-50B8-B302-D948-ABE3517CF6FC}
  - 19. PULSE OX MONITOR (POM) {A7E3E90C-4C8B-87D7-C2F8-49D9F3846CE8}

UNITS:

-----

PERCENTAGE 50 > 92 \* 100 < 101 Dec Prec: 0

Abbreviation..... SvO2%  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Venous Oxygen Saturation  
Display Name..... SvO2%  
Help Text..... Venous Oxygen Saturation  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----  
PERCENTAGE                    30 > 60 \* 75 < 100      Dec Prec: 0



=====

TERM: TIME OUT: PATIENT IDENTIFICATION METHOD  
ID: {339EAE62-5A1A-4F43-8D83-F9152EE5AB94}

=====

Page: 695

Abbreviation..... ID Method  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description.....  
Display Name..... Time Out: Patient Identification Method  
Help Text..... Method of Patient Identification  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
- 1. VERBALIZED BY PATIENT (Verbal)  
    {6AF901B8-A19D-4C89-869F-2ECBF2D105A7}
  - 2. ID BAND (ID Band)        {5526C5AB-64C8-4AB1-B77F-DAF45D26B439}
  - 3. ID CARD (ID Card )      {CFD540F0-3D73-498F-B590-CCEF11499078}

UNITS:  
-----



=====

TERM: TIME OUT: VERBAL CONFIRMATION BY TEAM OF PROCEDURE TO BE PERFORMED

Page: 696

ID: {AD9203C2-342A-4397-B71D-1AE9F36D4582}

=====

Abbreviation..... T/O Confir  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description.....  
Display Name..... Time Out: Verbal confirmation by team  
of procedure to be performed  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

=====

TERM: TIME OUT: VERBAL CONFIRMATION BY TEAM OF PROCEDURE AND LATERALITY

Page: 697

ID: {C90F1CBF-5CC7-49DB-8643-D220B20D954C}

=====

Abbreviation..... T/O Lat  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description.....  
Display Name..... Time Out: Verbal confirmation by team  
of procedure and laterality  
Help Text.....  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Void.....

UNITS:

-----

=====

TERM: TONOMETRY

Page: 698

ID: {C06989EF-4B0F-4941-B1A7-FA9D81A480FF}

=====

Abbreviation..... TON  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Tonometry  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid..... 4688726

UNITS:

-----  
MILLIMETERS OF MERCURY      0 > 0 \* 80 < 80      Dec Prec: 0





=====

TERM: TRACHEA MIDLINE DEVIATED  
ID: {E952DF23-0FC7-5809-2135-C8AA60FE159B}

=====

Page: 701

Abbreviation..... TRACH MIDL  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... No Description  
Display Name..... Trachea Midline Deviated  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:

---

=====

TERM: TRACTION TYPE

Page: 702

ID: {5A9CB3E0-7141-45E9-BE82-AA026A7889AB}

=====

Abbreviation..... TRACT TYPE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Traction Type  
Display Name..... Traction Type  
Help Text..... Type of traction used  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. BUCKS TRACTION (BUCK TRACT)  
{5E63F920-02D7-4729-8F30-A53A0A6E1EA0}
  2. SLING TRACTION (SLING TRAC)  
{304E924C-C95A-4DBE-B30B-191BC4D26EF7}
  3. PIN TRACTION (PIN TRACT) {C75104A5-3B6B-4CFE-8E2E-05877522C811}

UNITS:

-----





=====

TERM: TRAIN OF FOUR RESPONSE  
ID: {86E4352A-A72E-4287-8C17-ABCF97DCAF94}

=====

Page: 704

Abbreviation..... TOF RESP  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Train of Four response  
Display Name..... Train of Four Response  
Help Text..... Record the number of responses  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

-----

1. 1:4 (1:4)	{3F1B75FF-C168-D413-0BDF-6581E7A4EEAD}
2. 2:4 (2:4)	{7D423C9A-42B3-4FDD-999B-889193492F5C}
3. 3:4 (3:4)	{E4EBD0E2-FE8C-40FA-9DCA-475DCE43B1A2}
4. 4:4 (4:4)	{E9FCD5A0-995A-401D-8E85-6F0C753609B6}

UNITS:  
-----



=====

TERM: TRANSPORT FACILITATOR  
ID: {9F4001FD-C520-4360-87CC-1420803FF72A}

=====

Page: 706

Abbreviation..... TRANS FA  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Transport Facilitator  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. VOLUNTEER (VOLUNTEE) {F5406412-DE48-4F84-8A31-0B1D79DB576D}
  2. STAFF (STAFF) {BE8D678F-CE19-41FC-8C12-7A5DE7F0C7E2}
  3. PHYSICIAN (PHYSICIA) {DB0C5AE0-459A-4E04-A202-D1CF54B84A6A}
  4. STUDENT (STUDENT) {613AFF67-9121-4F6C-A44F-1ACFEF10E612}
  5. FAMILY/FRIEND (FAM/FRIE) {FA8FE272-E659-41EF-A394-F907A5F55ED8}

UNITS:  
-----

Abbreviation..... TRANS MO  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Transportation Mode/Device  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 10 entries

- 
1. ABLE TO AMBULATE (AMBULAT)  
{4E1D40CA-F5C4-442E-8BBA-C6E4C4781F57}
  2. CRUTCH (CRUTCH) {2F040C2C-7A95-4FFA-8DAC-8B2EA9A90259}
  3. CRUTCHES (CRUTCHES) {DDA81269-B25C-4F7D-A059-BB1602B2B179}
  4. CANE (CANE) {562B2E0B-0298-44C6-B33B-BC57DE6CBC3C}
  5. CANES (CANES) {DF37C541-1C7D-4527-A4D2-31867B5ED438}
  6. MOTORIZED WHEELCHAIR (MO WH CH)  
{4CA9CF41-F52C-41A8-B8A7-100A56097FDE}
  7. WHEELCHAIR (WH CHAIR) {EC1B82CE-4037-466F-90AC-52CD5C3B26BE}
  8. STRETCHER (STRETCHER) {93CB521A-3A67-41FB-A59C-3DD54D810B11}
  9. WALKER (WALKER) {91247001-5E89-46A0-A841-B4FAE9226E63}
  10. CHAIR (CHAIR) {5B8A5DBD-7462-4653-8D1F-CA8B1226DE38}

UNITS:  
-----



=====

TERM: TUBE POSITION VERIFIED  
ID: {1B065DF3-2D96-4F90-8CA2-7F377A0C2D28}

=====

Page: 709

Abbreviation..... TUBE POS V  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Was the tube position verified  
Display Name..... Tube Position Verified  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

=====

TERM: TURN AND SKIN CARE

Page: 710

ID: {42B46B96-C43B-4F36-B7F9-A3413F02EEA7}

=====

Abbreviation..... TURN SKIN  
Active..... Yes  
Boolean Value False..... Not Done  
Boolean Value True..... Done  
Data Type..... Boolean  
Description..... Turn and Skin Care  
Display Name..... Turn and Skin Care  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEA649CD} |

UNITS:  
-----





Abbreviation..... UP. LIM. S  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Upper Limb Extremities Strength  
Help Text..... No Help Text  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 6 entries

- 
- |                      |  |
|----------------------|--|
| 1. ABSENT (a)        | {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0} |
| 2. WEAK (w)          | {BE2D9305-E2DA-F492-66CF-164EE58B8179} |
| 3. STRONG (s)        | {A12154E1-544E-F372-3839-0D44E5BC476B} |
| 4. REFLEX (r)        | {F21585F9-5206-A6B3-2BF6-F416E5842366} |
| 5. EQUAL (Equal)     | {AC079584-158D-EBF2-E13A-C2E919375E4D} |
| 6. UNEQUAL (Unequal) | {73319200-164F-B39F-649C-B7FB1AFBDBB8} |

LOCATION: 3 entries

- 
1. L UPPER LIMB (L Upp Limb) {E6453716-62D3-DFCB-F84D-F6F2FDFC20B7}
  2. R UPPER LIMB (R UPP LIMB) {5E028BFF-ED9B-014F-74F0-C4137DB074FB}
  3. BILATERAL WRIST (BILAT WRIS)  
{89B4D30F-9E46-F8D6-14F6-1CC0843DDC21}

UNITS:

-----

Abbreviation..... URINE CHAR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Urine Characteristics
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 8 entries

- 1. ANURIC (ANURIC) {0CCA4B44-7A59-4042-BAB9-EB2395A63C51}
2. CLEAR (c) {8981C365-2B46-704F-C9F4-975DF1C23C2D}
3. CLOUDY (CLOUDY) {3CCFCCCC-00FC-6B00-1D58-A77C02FDB401}
4. HEMATURIA (HEMATURIA) {ABD21F4E-BC48-D43A-38CF-75DB7704C4B8}
5. PURULENT (PURULENT) {86C01494-743A-C6A2-8883-01CA5200C79B}
6. SEDIMENT (SEDIMENT) {4D4E4236-55CE-7ECD-FC55-6AD0CFD4D9A5}
7. CLOTS (CLOTS) {B52B88E0-6D23-4657-97CD-81E3EF64991B}
8. THREADS (Threads) {029BB303-B428-41D0-978B-6A8176108D7D}

UNITS:

=====

TERM: URINE COLOR  
ID: {80D67904-3FE2-98BA-6810-F15D88A637B3}

=====

Page: 714

Abbreviation..... URINE COLO  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Color of urine  
Display Name..... Urine Color  
Help Text..... Color of urine  
Multi Select Picklist..... Yes  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 14 entries

- 
- |                              |  |
|------------------------------|--|
| 1. YELLOW (YELLOW)           | {90DC31AD-6211-9F25-4237-79F7122A2763} |
| 2. WHITE (WHITE)             | {AB04EF54-5246-EF3A-E1D9-02E9E24369DD} |
| 3. BROWN (BROWN)             | {A1C32FBC-7B4D-707C-CA58-92A2E201368E} |
| 4. GREEN (GREEN)             | {3D1F1553-56FC-68CE-79A9-E4F6357B6AAF} |
| 5. RED (RED)                 | {515D7CFB-3527-30E7-EC4B-B607C272DFA0} |
| 6. TAN (TAN)                 | {AB994ED7-8916-FAE0-0E97-CD4E0EF80ED1} |
| 7. PINK (PINK)               | {1A7FAF84-D69C-4EEB-7BC5-D21A598A1F4D} |
| 8. CLEAR (c)                 | {8981C365-2B46-704F-C9F4-975DF1C23C2D} |
| 9. AMBER (AMBER)             | {90A13048-4C35-E01D-3482-64D782001D4A} |
| 10. BLACK (BLACK)            | {7A55AFB9-6861-4918-0746-3D67D47FE5C3} |
| 11. CLAY (CLAY)              | {A4E4D16C-B45C-0934-ED22-8BF640080108} |
| 12. ORANGE (ORANGE)          | {77B1B667-212D-88BA-D8E6-5C22F6ACDC9C} |
| 13. DYE PRESENT (DYE PRESEN) | {E865AD3E-6DAD-C3F9-6F93-7AD6F71234A4} |
| 14. TEA (TEA)                | {E2C0C287-DFDC-4F7B-8A67-7A339B8E8669} |

UNITS:  
-----

Abbreviation..... URINE VOID
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Urine Voids
Help Text..... Number of urine voids since last entry
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 11 entries

- 1. 1 (1) {8F5D10A8-FE17-47BC-DACA-6972D12CC00A}
2. 2 (2) {38E22B79-A2B6-3401-B195-C955F24E1278}
3. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}
4. 4 (4) {D2280C7E-DE50-5A8D-C3FD-5862B69C0E1A}
5. 5 (5) {188EDCE7-0726-4234-D8D2-92C8052CD221}
6. 6 (6) {640CF2F9-3BBF-AB13-185B-CD02A609FB5A}
7. 7 (7) {7A010324-D2C2-5D4B-848A-7B8878666D20}
8. 8 (8) {518E7CFC-102E-351E-C7AF-FB575E3C8F4D}
9. 9 (9) {58FF2FD4-569F-F454-9F5B-FE6F9E34D8CC}
10. 10 (10) {0C827160-18DB-816F-7FAA-A3C175C6D33B}
11. N/A (n) {B3DBE3FD-C5CB-3894-607A-BE13319594D0}

UNITS:

Abbreviation..... VOIDING
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Urine Voiding Method
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 12 entries

- 1. CONTINENT (CONTINENT) {3B658506-F42E-BAF5-04C4-5AD49699E0D2}
2. DYSURIA (DYSURIA) {75AC5ECD-8975-406A-30C7-E2EB23F10F84}
3. INCONTINENT (INCONTINEN) {BA2E9C32-8D03-8FC9-5E02-EAE2026FDC63}
4. DIVERSION (DIVERSION) {087AAE9F-7FE1-5717-6423-07574C637875}
5. INDWELLING CATHETER (INDW CATH) {EC16724C-D1EA-5650-58D4-38A89D4456F9}
6. EXTERNAL CATHETER (EXT. CATHE) {C8AEE936-DD91-FE91-4DD6-6F210FB59726}
7. STRAIGHT CATHETER (STR. CATH) {5F4A8A7F-E720-9458-2360-1DC14E7D05F3}
8. ILEAL CONDUIT (ILEAL COND) {032CB6FE-77D2-702E-A0AE-46176F5CAA57}
9. SUPRAPUBIC (SUPRAPUBIC) {69E0176F-52FC-A3C4-7A4E-AD57DBD6E572}
10. L NEPHROSTOMY (L NEPHROST) {D8FC4C18-081E-A7B9-E394-25F194BE9757}
11. R NEPHROSTOMY (R NEPHROST) {33661137-688D-5611-4368-B1DFD8DB1469}
12. REPORTED (REPORTED) {1D13BAC5-D8F9-455C-B815-8E256943492E}

QUALITY: 2 entries

- 1. ACTUAL (A) {5842A4D9-56A1-7A5A-13C9-F481699C3A09}
2. STATED (Sta) {6362C194-B7F2-D38C-5B72-3BE5DE09438E}

UNITS:



=====

TERM: VENTILATOR CIRCUIT CHANGE  
ID: {C0494BB3-5896-41BD-89B4-86C684078604}

=====

Page: 718

Abbreviation..... VENT CIR C  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Ventilator Circuit Change  
Display Name..... Ventilator Circuit Change  
Help Text..... Ventilator Circuit Change  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
  - 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEEA649CD}

METHOD: 1 entries

- 
- 1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:  
-----

=====

TERM: VENTILATOR EXHALED VOLUME  
ID: {404A3850-382D-4E86-B7F3-29DDAF1616BC}

=====

Page: 719

Abbreviation..... VENT EXHAL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Ventilator Exhaled Volume  
Display Name..... Ventilator Exhaled Volume  
Help Text..... Ventilator Exhaled Volume  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

-----

LITER 0.0 > 1.0 \* 99.0 < 100.0 Dec Prec: 1



=====

TERM: VENTILATOR EXPIRATORY TIME  
ID: {004C6736-FFB2-4002-B92C-E2F17D03AF4E}

=====

Page: 720

Abbreviation..... VENT EXP T  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Ventilator Expiratory Time  
Display Name..... Ventilator Expiratory Time  
Help Text..... Ventilator Expiratory Time  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

-----

SECOND 0.00 > 0.01 \* 9.99 < 10.00  
Dec Prec: 2

=====

TERM: VENTILATOR FLOW RATE  
ID: {548E1F72-B501-49DD-B22B-1366CE06C2BE}

=====

Page: 721

Abbreviation..... VENT FL RA  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Ventilator Flow Rate  
Display Name..... Ventilator Flow Rate  
Help Text..... Ventilator Flow Rate  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

-----

LITERS PER MINUTE -1 > 0 \* 150 < 151 Dec Prec: 0

=====

TERM: VENTILATOR INSPIRATORY TIME  
ID: {35948B65-6D8F-4BD4-8F4E-CE785D8D829D}

=====

Page: 722

Abbreviation..... VENT INS T  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Ventilator Inspiratory Time  
Display Name..... Ventilator Inspiratory Time  
Help Text..... Ventilator Inspiratory Time  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

-----

SECOND 0.00 > 0.01 \* 9.99 < 10.00  
Dec Prec: 2

=====

TERM: VENTILATOR MEAN AIR PRESSURE  
ID: {BD28B791-7E79-4116-8DF1-6296B60FCC57}

=====

Page: 723

Abbreviation..... VENT MAP  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Ventilator Mean Air Pressure  
Display Name..... Ventilator Mean Air Pressure  
Help Text..... Ventilator Mean Air Pressure  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

-----

CENTIMETERS H2O	0 > 1 * 99 < 100	Dec Prec: 0
MILLIBAR	0 > 20 * 80 < 100	Dec Prec: 0

=====

TERM: VENTILATOR MINUTE VOLUME  
ID: {1C718F1D-EE59-4743-888D-7767864C4EB0}

=====

Page: 724

Abbreviation..... VENT MIN V  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Ventilator Minute Volume  
Display Name..... Ventilator Minute Volume  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

-----

LITERS PER MINUTE -1 > -1 \* 30 < 30 Dec Prec: 0

=====

TERM: VENTILATOR MODE  
ID: {21CF6705-7597-53D0-F5EA-CBB704128A4B}

=====

Page: 725

Abbreviation..... VENT MODE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Ventilator Mode  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 7 entries

- 
1. CONTROL MODE VENTILATION (CMV)  
{2E79351C-BA31-FC04-B636-F39728B23EA6}
  2. ASSIST CONTROL (AC) {862ED94C-9DC3-DFAF-F824-0F7DE348C2CB}
  3. SYNCHRONOUS INTERMITTENT MANDATORY VENTILATION (SIMV)  
{A606BB2B-5C00-FEA1-127F-1D3455AE5996}
  4. CONTINUOUS POSITIVE AIR PRESSURE (CPAP)  
{387B42F0-F587-AAFC-ECA0-66067F18BEC6}
  5. PRESSURE CONTROL (PRESS CONT)  
{C3CDF756-484F-1483-9883-9B63DEF3B012}
  6. BIPHASIC POSITIVE AIRWAY PRESSURE (BiPAP)  
{4802AB1E-A0AD-BB9F-E8BA-33863BCD3466}
  7. SPONTANEOUS (SPONTANEOU) {C137C4A7-08A6-4064-8F95-4A4476E2E3A8}

METHOD: 1 entries

- 
1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:  
-----







Abbreviation..... VT PT RATE
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Rate that the patient actually achieved
Display Name..... Ventilator Patient Rate - Do Not
Check Total or Subtotals
Help Text..... Value that is determined by taking
the number of breaths total per
minute minus the breaths set to be
delivered by the machine
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

BREATHS PER MINUTE 0 > 1 \* 119 < 120 Dec Prec: 0



=====

TERM: VENTILATOR PCV SET PRESSURE  
ID: {BC521E68-B7AA-4B6E-AAC1-0594EE27D94E}

=====

Page: 730

Abbreviation..... VENT PCV S  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Pressure Control Ventilator Pressure  
Display Name..... Ventilator PCV Set Pressure  
Help Text..... Pressure Control Ventilator Pressure  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

-----

CENTIMETERS H2O 0 > 1 \* 99 < 100 Dec Prec: 0



=====

TERM: VENTILATOR PEEP  
ID: {4847847A-7EB7-9A8F-D565-1DF26438A616}

=====

Page: 732

Abbreviation..... VENT PEEP  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... PEEP: Positive End Expiratory Pressure  
Display Name..... Ventilator PEEP  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

-----

CENTIMETERS H2O -1 > 0 \* 20 < 21 Dec Prec: 0

=====

TERM: VENTILATOR PRESSURE SUPPORT  
ID: {B6ABD066-7B39-126F-C4C2-A4243C2F93E9}

=====

Page: 733

Abbreviation..... V PRES SUP  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Ventilator Pressure Support  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

-----

CENTIMETERS H2O -1 > 0 \* 20 < 40 Dec Prec: 0

=====

TERM: VENTILATOR PSV SET PRESSURE  
ID: {B2494FE0-0B78-4BD9-97B6-DE604B711801}

=====

Page: 734

Abbreviation..... VENT PSV S  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Pressure Support Ventilator Pressure  
Display Name..... Ventilator PSV Set Pressure  
Help Text..... Pressure Support Ventilator Pressure  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

-----

CENTIMETERS H2O 0 > 1 \* 99 < 100 Dec Prec: 0

=====

TERM: VENTILATOR RATE  
ID: {001D591B-6CE7-F51F-098C-79AC6B3E0E09}

=====

Page: 735

Abbreviation..... VENT RATE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Ventilator Rate Set  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----

1. VENTILATOR (Ven)      {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

-----

BREATHS PER MINUTE      0 > 8 \* 40 < 120      Dec Prec: 0



=====

TERM: VENTILATOR SENSITIVITY  
ID: {52F9AF8D-CF31-430B-8A11-A8115BFF4B2E}

=====

Page: 736

Abbreviation..... VENT SENSI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Ventilator Sensitivity  
Display Name..... Ventilator Sensitivity  
Help Text..... Ventilator Sensitivity  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

-----

CENTIMETERS H2O -10.0 > -9.9 \* -0.1 < 0.0 Dec Prec: 1



=====

TERM: VENTILATOR TIDAL VOLUME SET  
ID: {90D4A5CA-E21B-06F6-7AAC-FC6D4B11E7AC}

=====

Page: 738

Abbreviation..... TIDL VOL S  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Tidal Volume Set  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

METHOD: 1 entries

-----

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

-----

MILLILITER	0 > 800 * 1200 < 5000	Dec Prec: 0
LITER	0.000 > 0.800 * 1.200 < 5.000	Dec Prec: 3

=====

TERM: VENTILATOR TYPE  
ID: {58151FCB-646B-4CC5-A89F-6D34DFB7E622}

=====

Page: 739

Abbreviation..... VENT TYPE  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Type of Ventilator  
Display Name..... Ventilator Type  
Help Text..... Type of Ventilator  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 13 entries

- 
1. VENTILATOR PB 840 (PB 840)  
{25A7B0DA-E437-4201-B8F3-3663B1B91C1A}
  2. VENTILATOR ACHIEVA P5O2 (ACHIEVA P5)  
{1690EEA5-3C3B-4D30-AC5C-214C94626699}
  3. VENTILATOR ADULT STAR (ADULT STAR)  
{879B13DA-4874-42D4-99EC-8E4ADB31B14E}
  4. VENTILATOR BEAR 1000 (BEAR 1000)  
{03C30F99-A27C-427C-B0ED-872F0DE646BD}
  5. VENTILATOR DRAEGER (DRAEGER)  
{DD18C42F-EF13-460B-8676-9ED1E429B54B}
  6. VENTILATOR PB 7200ae (PB7200ae)  
{86AE485F-ECD1-4B8E-820F-71CF159F318B}
  7. VENTILATOR SERVO (SERVO) {D6F06715-1ADD-4941-B026-B9A5E463943E}
  8. BIPHASIC POSITIVE AIRWAY PRESSURE (BiPAP)  
{4802AB1E-A0AD-BB9F-E8BA-33863BCD3466}
  9. CONTINUOUS POSITIVE AIR PRESSURE (CPAP)  
{387B42F0-F587-AAFC-ECA0-66067F18BEC6}
  10. VENTILATOR OTHER (VENT OTHER)  
{138B7B3D-9F8B-447A-9785-12AB396AE6DE}
  11. BIPAP - RESPIRONICS (BIPAP RESP)  
{AA94A723-4C68-489F-88AC-2C6F276EAB76}
  12. BIPAP - VISION (BIPAP VISI)  
{AA39C765-BA1F-4289-B459-200A57C07FFF}
  13. BIPAP/CPAP (BIPAP/CPAP) {8D767C67-3712-410F-A96B-69C731E03C5D}

METHOD: 1 entries

- 
1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:  
-----

=====

TERM: VENTILATOR WAVE FORM  
ID: {17F52B3E-B4C9-483A-B02F-0F4DD809E673}

=====

Page: 740

Abbreviation..... VENT WA FO  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Ventilator Wave Form Setting  
Display Name..... Ventilator Wave Form  
Help Text..... Ventilator Wave Form Setting  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 3 entries

- 
1. SQUARE (SQUARE) {63F70678-2CAF-48E8-964D-2805D8B9A43B}
  2. ACCELERATING (ACCELERATI) {514FE300-87EE-4C2C-AC70-5E57DBECAFF7}
  3. DECELERATING (DECELERATI) {B259D4D4-82B8-454A-B53D-F12F4952BB61}

METHOD: 1 entries

- 
1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:  
-----

=====

TERM: VISION CORRECTED  
ID: {ED022AC1-EBE4-E708-684D-63D00628A94C}

=====

Page: 741

Abbreviation..... VC  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Vision Corrected  
Help Text..... Enter Acuity value XX/20  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid..... 4688727

LOCATION: 2 entries

- 
- 1. L EYE (L EYE)            {D714A0CA-EAD6-EC35-DA54-FF027C23FF45}
  - 2. R EYE (R EYE)            {258F0C1D-67BD-ECF8-3152-EB289DEF4681}

UNITS:

-----

FOOT                    -1 > 0 \* 200 < 201      Dec Prec: 0

=====

TERM: VISION UNCORRECTED  
ID: {BEA5E565-D728-F5B3-0A3A-0528C42A45C4}

=====

Page: 742

Abbreviation..... VU  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Vision Uncorrected  
Help Text..... Enter Acuity value XX/20  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid..... 4688728

UNITS:

-----

FOOT                    -1 > 0 \* 200 < 201      Dec Prec: 0

=====

TERM: VT/VC - DO NOT CHECK TOTAL OR SUBTOTALS  
ID: {FF8D8ABA-E89E-429A-AE65-7C25E22EF199}

=====

Page: 743

Abbreviation..... VT/VC  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Tidal Volume / Vital Capacity  
Display Name..... Vt/VC - Do Not Check Total or Subtotals  
Help Text..... Tidal Volume / Vital Capacity  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

LITER	0.00 > 1.00 * 99.00 < 100.00
	Dec Prec: 2



=====

TERM: WARMING COOLING DEVICE  
ID: {34F4BAF5-8388-4431-8AE1-F04670A25FE4}

=====

Page: 744

Abbreviation..... WM CL DEV  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Warming, Cooling Device  
Help Text..... Warming or Cooling Device  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 11 entries

-----

1. AIR FLOW DEVICE (AIR FLOW D)  
{60AF6215-CE2A-4EEE-9596-909E3E50BE71}
2. K-PAD (K-PAD) {10A0C893-F279-48ED-ACAA-1DCA4AD9DB24}
3. WARMING BLANKET (WARM BLANK)  
{23D20AE5-DE25-4456-8F77-BF0411FB19A6}
4. HUMIDIFIER (HUMIDIFIER) {A4CE67D1-F885-423B-B914-392D8CAC697E}
5. HEAT LAMP (HEAT LAMP) {4E9DE229-F810-4AB0-9FD7-4596759CEA11}
6. HEAT CHEM PACK (HT CHEM PA)  
{D759F891-2535-40C2-A68C-678170E56F7A}
7. COOL CHEM PACK (CO CHEM PA)  
{6DE5E644-28C2-4190-A2EF-1A766E014A65}
8. ICE PACK (ICE PACK) {E46ED309-D60D-45DF-BFA3-F74DD13FC4AE}
9. WARMING BOTTLE (WARM BOTTL)  
{AD1BEFE5-4CE6-4076-8095-9BE7387CFE25}
10. WARM HUMIDITY WITH O2 (WM HUM O2)  
{E62E07C3-4CA2-426B-8EBA-48C406837D9C}
11. COOLING BLANKET (COOL BLANK)  
{13A42B12-4959-4674-AF97-2D1242CFD934}

UNITS:

-----

=====

TERM: WARMING COOLING DEVICE APPLICATION  
ID: {17908A29-AF84-4217-A035-1CE6CA0F7316}

=====

Page: 745

Abbreviation..... WM CL APPL  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Warming/Cooling Device Application  
Help Text..... Current State of Warming/Cooling Device  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 4 entries

- 
1. WARMING ON (WARMING ON) {6FFDE928-DE1E-4B25-BA3C-EB29E0CDEB35}
  2. WARMING OFF (WARMING OF) {98FAC977-8E2E-406A-9217-A77E1DDCDF50}
  3. COOLING ON (COOLING ON) {911CA17D-614E-445A-969E-4F5220E469EB}
  4. COOLING OFF (COOLING OF) {14F573ED-0D10-4E0F-A1F3-A11F29D6DB11}

UNITS:  
-----

Abbreviation..... WT  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Weight of patient  
Display Name..... Weight  
Help Text..... Patients weight in lb or kg.  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid..... 4500639

METHOD: 8 entries

-----

1. BED (B) {20188D9F-3815-1F7F-34EC-BA77F75B310E}
2. CHAIR (Ch) {BBF5C05C-220D-D052-C7E8-236E791A4FC9}
3. LIFT SCALE (LS) {8D741461-25C5-3955-BE2C-CBD27317CBA1}
4. STANDING WEIGHT (SW) {BE16FA09-16DA-510B-D235-B7F0F93AB105}
5. WHEELCHAIR SCALE (WcS) {8C496753-9F42-2A36-E5CA-750188306ED4}
6. WITH CAST OR BRACE (WCB) {47709C05-DF21-FB7D-DE7F-0699D1B7E457}
7. WITH PROSTHESIS (WP) {73C93D2D-934B-CA55-3557-E756A89710F5}
8. WITHOUT PROSTHESIS (WOP) {9AAC852-E70D-1048-FCE7-112D46C12619}

POSITION: 3 entries

-----

1. STANDING (St) {0776AE42-C370-B7B3-9B4D-9C914BA8F231}
2. SITTING (Sit) {DD4DCDF1-0A22-2CC5-404F-0E472C73961C}
3. LYING (Ly) {DE28D2BF-6F8D-64C3-0092-4CEB5E0902B9}

QUALITY: 4 entries

-----

1. ACTUAL (A) {5842A4D9-56A1-7A5A-13C9-F481699C3A09}
2. DRY (D) {B319EFF4-BA4A-85B1-2617-E00E46A54ABE}
3. ESTIMATED (E) {9EAF26B-7D87-DA66-F991-022D095FB5A0}
4. STATED (Sta) {6362C194-B7F2-D38C-5B72-3BE5DE09438E}

UNITS:

-----

POUND 0.0 > 1.0 \* 1499.0 < 1500.0  
Dec Prec: 1

KILOGRAM 0.0 > 0.5 \* 680.3 < 680.4 Dec Prec: 1

=====

TERM: WEIGHT AT ADMISSION  
ID: {F9C859E2-ACD7-44F5-8E0D-3166C3D0A61F}

=====

Page: 747

Abbreviation..... WT AT ADMI  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Weight at Admission  
Help Text..... Weight of patient at admission time  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

POUND	0.0 > 1.0 * 1499.0 < 1500.0 Dec Prec: 1
KILOGRAM	0.0 > 0.5 * 680.3 < 680.4 Dec Prec: 1

=====

TERM: WEIGHT DRY

Page: 748

ID: {8594A12F-97BD-41A5-802A-76E5317448C6}

=====

Abbreviation..... WEIGHT DRY  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... No Description  
Display Name..... Dry Weight  
Help Text..... Patient Dry Weight  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

QUALITY: 3 entries

-----

1. ACTUAL (A) {5842A4D9-56A1-7A5A-13C9-F481699C3A09}
2. ESTIMATED (E) {9EAF26B-7D87-DA66-F991-022D095FB5A0}
3. STATED (Sta) {6362C194-B7F2-D38C-5B72-3BE5DE09438E}

UNITS:

-----

POUND 0.0 > 1.0 \* 1499.0 < 1500.0  
Dec Prec: 1

KILOGRAM 0.0 > 0.5 \* 680.3 < 680.4 Dec Prec: 1

=====

TERM: WOUND BED COLOR  
ID: {2B890957-8886-4BC8-84AC-D46DDAE733A0}

=====

Page: 749

Abbreviation..... WOU BED CO  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Wound Bed Color  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 8 entries

- 
- |                    |  |
|--------------------|--|
| 1. BLACK (BLACK)   | {7A55AFB9-6861-4918-0746-3D67D47FE5C3} |
| 2. GRAY (GRAY)     | {D38ABEE0-52D7-4316-A713-40ACF4250599} |
| 3. GREEN (GREEN)   | {3D1F1553-56FC-68CE-79A9-E4F6357B6AAF} |
| 4. PINK (PINK)     | {1A7FAF84-D69C-4EEB-7BC5-D21A598A1F4D} |
| 5. RED (RED)       | {515D7CFB-3527-30E7-EC4B-B607C272DFA0} |
| 6. TAN (TAN)       | {AB994ED7-8916-FAE0-0E97-CD4E0EF80ED1} |
| 7. WHITE (WHITE)   | {AB04EF54-5246-EF3A-E1D9-02E9E24369DD} |
| 8. YELLOW (YELLOW) | {90DC31AD-6211-9F25-4237-79F7122A2763} |

UNITS:  
-----



=====

TERM: WOUND DEPTH

Page: 751

ID: {B8FDF2BE-E6ED-49E3-8F42-4CEF94E29C47}

=====

Abbreviation..... WOU DEPTH  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Wound Depth  
Display Name..... Wound Depth  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

CENTIMETER	-1.00 > 0.00 * 200.00 < 400.00
	Dec Prec: 2



Abbreviation..... WOU DR CHA  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Wound Drainage  
Display Name..... Wound Discharge Characteristics  
Help Text..... Wound Discharge Characteristics Pick  
One  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Void.....

DATA VALUE: 7 entries

-----

1. WOUND DISCHARGE CHARACTERISTICS NONE (WOU DI NON)  
{929F2FA5-923F-4D88-8EDF-E0EBEE9CC7A8}
2. WOUND DISCHARGE CHARACTERISTICS BLOODY (WOU DI BLO)  
{628DBFBC-9B77-48EA-99FE-0DF000376506}
3. WOUND DISCHARGE CHARACTERISTICS CLOUDY (WOU DI CLO)  
{159C9993-A2E1-4133-9D5F-315EF9985EC4}
4. WOUND DISCHARGE CHARACTERISTICS PURULENT (WOU DI PUR)  
{C1E10410-B6D8-4429-BD94-0B179992EAAB}
5. WOUND DISCHARGE CHARACTERISTICS SEROSANGUINOUS (WOU DI SSA)  
{B431D8EF-5B53-4DCF-A390-431EB5ED3704}
6. WOUND DISCHARGE CHARACTERISTICS SEROUS (WOU DI SER)  
{2124005E-5C36-48AC-AA6B-03E06510D70A}
7. WOUND DISCHARGE CHARACTERISTICS VISCOUS (WOU DI VIS)  
{382FA6D9-1CB2-4ECC-A458-8598ABA708F8}

UNITS:

-----

Abbreviation..... WOU DIS CO
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Wound Discharge Color
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 15 entries

- 1. AMBER (AMBER) {90A13048-4C35-E01D-3482-64D782001D4A}
2. BLACK (BLACK) {7A55AFB9-6861-4918-0746-3D67D47FE5C3}
3. BLUE (BLUE) {8FA78575-E1DD-4E6E-BC43-9E5C92E2B7E5}
4. BROWN (BROWN) {A1C32FBC-7B4D-707C-CA58-92A2E201368E}
5. COLORLESS (COLORLESS) {D093D8C5-E284-47D4-8FA6-DCD5B030B10E}
6. CLAY (CLAY) {A4E4D16C-B45C-0934-ED22-8BF640080108}
7. GREEN (GREEN) {3D1F1553-56FC-68CE-79A9-E4F6357B6AAF}
8. MAROON (MAROON) {1F518B23-C79F-4F06-9DF3-ADC5C96D3029}
9. ORANGE (ORANGE) {77B1B667-212D-88BA-D8E6-5C22F6ACDC9C}
10. PINK (PINK) {1A7FAF84-D69C-4EEB-7BC5-D21A598A1F4D}
11. RED (RED) {515D7CFB-3527-30E7-EC4B-B607C272DFA0}
12. TAN (TAN) {AB994ED7-8916-FAE0-0E97-CD4E0EF80ED1}
13. VIOLET (VIOLET) {1249AE97-CF1F-4604-A17A-5FC85F6D0185}
14. WHITE (WHITE) {AB04EF54-5246-EF3A-E1D9-02E9E24369DD}
15. YELLOW (YELLOW) {90DC31AD-6211-9F25-4237-79F7122A2763}

UNITS:
=====

=====

TERM: WOUND DISCHARGE ODOR  
ID: {67D752F7-CC83-4D76-96EB-ECE6DBA40210}

=====

Page: 754

Abbreviation..... WOU DIS OD  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Wound Discharge Odor  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
  2. FOUL SMELLING (FOUL SMELL)  
{761346A9-329D-4F42-8D80-0E9D8D40653C}

UNITS:  
-----

ID: {8EC56666-A337-AA7F-5109-A10155D4FD58}

=====

Abbreviation..... WOU DR TY  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Type of drain used  
Display Name..... Wound Drain Type  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 10 entries

-----

1. BULB SUCTION (Bulb Sucti) {B0839A82-BB79-BB9F-5E4E-0ED67F716A50}
2. HEMOVAC (HEMOVAC) {82DA2DA7-E20B-E571-659C-457DF42472C2}
3. JP DRAIN (JP Drain) {20F3B3D5-D3A4-96FD-5446-253FD031DF2B}
4. LUMBAR DRAIN (LUM DRAIN) {FF36A046-1799-0734-286B-2272C34E3056}
5. PENROSE DRAIN (Penrose D) {DE508E95-D044-ACCA-EE69-33B9A8454E6C}
6. SUMP (Sump) {DAD881B9-F1EE-FDAD-AEB0-68C754FDC24E}
7. TWILL TAPE (Twill Tape) {C52BE250-732D-03BB-B8F9-72BE4D546807}
8. VENTRICULOSTOMY (Ventriculo)  
{ECF20317-4708-7661-D4C2-9BDB4CAAC4FA}
9. WOUND VAC DEVICE (WOUND VAC)  
{A03F4978-969D-DCB8-DCE2-26EB33F097B3}
10. VACUUM ASSISTED CLOSURE (VAC ASS CL)  
{9C3DD561-FECF-47FA-BBA5-4C08FF9A4CA1}

UNITS:

-----

=====

TERM: WOUND DRAINAGE AMOUNT  
ID: {3BABF5D4-BD30-4814-A397-EEEEAE0EE1D13}

=====

Page: 756

Abbreviation..... WOU DRA AM  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Amount of drainage from wound  
Display Name..... Wound Drainage Amount  
Help Text..... Amount of drainage from wound  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. WOUND DRAINAGE AMOUNT NONE (WOU DR A N)  
{92910398-311F-D4D7-8BD5-AB66C3DC5E72}
  2. WOUND DRAINAGE AMOUNT SCANT (W DR SCANT)  
{9C3CD0BD-F960-4098-B14C-17B20703B242}
  3. WOUND DRAINAGE AMOUNT SMALL (W DR A SMA)  
{3D8EDF1A-EC4B-41EF-A194-73FF60074A57}
  4. WOUND DRAINAGE AMOUNT MODERATE (W DR A MOD)  
{E5DBFD9C-A2AC-4009-98CB-7246F6AA702D}
  5. WOUND DRAINAGE AMOUNT HEAVY (W DR A HEA)  
{BF959172-6672-4FEB-8FB7-0568405F3C83}

UNITS:  
-----

=====

TERM: WOUND DRESSING CHANGE  
ID: {01A45F49-9EBB-44E7-8B9A-6920886F51BB}

=====

Page: 757

Abbreviation..... WOU DR CHA  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Has Wound Dressing Been Changed?  
Display Name..... Wound Dressing Changed  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

Abbreviation..... WOU DR TY
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Wound Dressing Type
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 15 entries

- 1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
2. ANTI-MICROBIAL DRESSING (AB DRESSIN) {9734C5E9-D099-47EB-93B1-D82201F64331}
3. ANTIMICROBIAL DISK (WOU DR ADI) {1F1745EA-D763-4AFD-9CED-82E974E57692}
4. BALSAM OF PERU (WOU DR PER) {DD412BB1-B1C1-4A48-A38F-832629448E06}
5. BISMUTH GAUZE (WOU DR BIS) {010DED59-C808-4D92-A7F8-87AEB86B4439}
6. CALCIUM ALGINATE (CA ALGINAT) {2B2ACD25-A77B-4D39-8598-B26FD9BAFFC2}
7. FOAM DRESSING (WOU DR FOA) {D18A5EDD-0EF6-470E-8370-4EDB6D322477}
8. GAUZE (GAUZE) {CF10C658-4334-4F3B-A38E-E98894458333}
9. GAUZE PACKING (GAUZE PACK) {6C093717-82CB-4966-929E-C909310649C1}
10. GAUZE WRAP (GAUZE WRAP) {040A9648-E6E1-4AC5-93C7-53A49537B15C}
11. HYDROCOLLOIDAL DRESSING (HCOLL DRES) {D85741EF-0DB3-48CF-BC9A-88C834A72F36}
12. HYDROGEL DRESSING (H-GEL DRES) {F6192D10-D6EA-42C2-AB15-270659D9925C}
13. HYPERTONIC HYDROGEL (WOU DR HYP) {533D6FC3-5183-4417-AB89-9FAD9BAD0F44}
14. PRESSURE (PRESSURE) {387AFB7F-3973-4D77-B2E6-0724F344D8F6}
15. TRANSPARENT (TRANSPAREN) {0E73D24D-0EFC-08F9-196A-D30DD1BB3228}

UNITS:

=====

TERM: WOUND EDGE

Page: 759

ID: {BC1B996C-35E6-4A05-BDF6-00DF84EA90EF}

=====

Abbreviation..... WOU ED  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Wound Edge  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. WOUND EDGE INDISTINCT/ATTACHED (WOU ED IND)  
{C899235C-4A07-4420-92BC-8F818B8E0A7F}
  2. WOUND EDGE NON-ATTACHED (WOU ED NON)  
{1D5C1317-B4FF-421C-96A0-C20C1089E672}
  3. WOUND EDGE TUNNELING (WOU ED TUN)  
{97425206-5261-46C1-905F-E03E2F071381}
  4. WOUND EDGE UNDERMINING (WOU ED UND)  
{AAA3BB3A-3A3B-4F39-9AC7-3BE1A82BDC02}
  5. WOUND EDGE ROLLED (WOU ED ROL)  
{24D9CFEC-1ED6-4069-9715-1F9FCB45BBD3}

UNITS:  
-----



Abbreviation..... WOU HEA ST  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Wound Healing Status  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 8 entries

- 
1. EPITHELIALIZATION (EPITHELIAL) {5FC43A6D-921D-4736-A234-9993DB6499B7}
  2. ESCHAR (ESCHAR) {555E03C5-FBFE-4EB9-816F-25CE34476B3C}
  3. FIBRINOUS (FIBRINOUS) {D286EAE2-6902-4B7E-845D-BBB316941B2C}
  4. GRANULATING (GRANULATI) {F8930276-3954-4BFE-82B0-6C44A8E9CAE0}
  5. HYPERGRANULIZATION (HYPERGRANU) {C14B9FB0-6F33-4ACB-88B3-4B799594F845}
  6. NECROTIC (NECROTIC) {86885A34-E4EC-4BB5-BD94-59DD4B0F1A0E}
  7. SLOUGH (SLOUGH) {548EB536-1AC4-4234-8C4D-DC54E4CAA6F6}
  8. SMOOTH (SMOOTH) {5D71A0D5-52ED-4B8D-B359-AC7D500660D8}

UNITS:  
-----

Abbreviation..... WOU IDENTI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Wound Identified
Display Name..... Wound Identified
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

- 1. ADMISSION (ADMISSION) {483C2830-813A-46C3-A012-CCD70044B2A2}
2. DISCHARGE (DISCHARGE) {ADF67CC2-E05B-4B8D-A233-4AE1C35B2D20}
3. TRANSFER (TRANSFER) {DFDB4783-25BE-485B-A9BA-621F1004AAF2}
4. POST-OPERATIVE (POST-OP) {9FA19656-849F-4A65-92FC-755A4C9E8056}
5. POST-FALL (POST-FALL) {9BC82B93-CA4A-44E3-A3BA-7AF72BAFD44B}

LOCATION: 70 entries

- 1. R WRIST (R WRIST) {FDD9FE1D-F8A2-09AD-82D7-AA478693DD01}
2. L WRIST (L WRIST) {7BCAFACE-FF3B-46DA-99F9-6437A90B7126}
3. R TESTICLE (R TESTICLE) {0947D277-50AD-47F8-A76D-0EC34EB3B168}
4. L TESTICLE (L TESTICLE) {1E376539-E89F-4472-8274-CFD11275E11F}
5. R SHOULDER (R SHOULDER) {D27AE130-442E-4DD6-B391-8D7E8AD494A2}
6. L SHOULDER (L SHOULDER) {91CBD2FA-2743-42EB-946B-6933980D9337}
7. RECTAL (RECTAL) {0E7069F4-D3E6-67A1-E821-EF1879C3C267}
8. NECK (Ne) {D4812C3B-AAC2-1B61-679B-3CD8B5CB947C}
9. ORAL (O) {7A7ACA1B-ACBD-B08D-D40C-54BA8E724B8C}
10. R UPPER LEG (R UPPER LE) {113BD5ED-5A0A-47CA-B100-91B244BFD04A}
11. L UPPER LEG (L UPPER LE) {1DB35895-F74C-4B88-A3EE-3E556CB6559A}
12. R LOWER LEG (R LOW LEG) {C8EA4D71-D97E-46BD-AE9D-C10954C525CD}
13. L LOWER LEG (L LOW LEG) {7B8E92FA-E065-4D06-9FB6-4234656CFD61}
14. R INNER KNEE (R INN KNEE) {2AE24BF1-EA1B-490D-9DB7-EEFC9314C560}
15. R OUTER KNEE (R OUT KNEE) {330F1DF1-8C99-4984-9965-7404ABC9562D}
16. L INNER KNEE (L INN KNEE) {79853581-CCF0-49E9-9D91-D16DE02A836D}
17. L OUTER KNEE (L OUT KNEE) {B143C935-2687-49C3-9F65-8D69C1775ABC}
18. R HIP (R HIP) {A48CBCAE-5404-4934-9321-ECEF660A5350}
19. L HIP (L HIP) {786CD881-27AD-4136-B9CC-8423CB38A431}
20. L HEAD (L HEAD) {3737DC2D-DC9E-4FC8-839D-0355987B62D1}
21. R HEAD (R HEAD) {42A95AAE-3DA9-4C96-A4E7-A15C9840217A}
22. R HAND (R HAND) {A54D3A80-8537-D079-9CCE-7DC56346D132}
23. L HAND (L HAND) {639DF25A-19E3-B504-D323-4B3F6ED39A24}
24. GENITALS (GENITALS) {83437106-3336-4925-92C8-F959561687FC}
25. GENERALIZED (GENERIZED) {17073C6C-8691-CADC-2238-CBAAB0A5C648}
26. R FOOT (R FOOT) {7D9F5FF6-232C-2488-4504-A721FF06928F}
27. L FOOT (L FOOT) {2EC3FBEC-4EA2-9E7A-17D4-FDC2F9D28046}
28. R FACE (R FACE) {D1B15677-EF6F-0125-3601-68C1FE6386CE}
29. L FACE (L FACE) {06F00A70-F475-CD40-665D-AF0569A1E75F}
30. R EYE (R EYE) {258F0C1D-67BD-ECF8-3152-EB289DEF4681}
31. L EYE (L EYE) {D714A0CA-EAD6-EC35-DA54-FF027C23FF45}
32. R ELBOW (R ELBOW) {5BAFF87B-6993-4E76-862F-9F386D9B71C5}
33. L ELBOW (L ELBOW) {BF4317E7-2024-47EF-9E5F-3D21D0547885}
34. R EAR (R EAR) {AE31327D-A03C-488C-8A8E-084E7E6E7FDE}
35. L EAR (L EAR) {DA9C2DD1-D6FF-4899-A62B-1096FAA36C59}
36. R CHEST (R CHEST) {7E30EC0E-815A-4C2E-A9FC-3C2B2D6E3CD4}

- 37. L CHEST (L CHEST) {E332E497-FE0C-47CB-A0A7-217FBF52DC97}
- 38. R GLUTEUS (R GLUTEUS) {B4D32045-68D1-4457-9106-7937887FA6F6}
- 39. L GLUTEUS (L GLUTEUS) {88C4E547-04CD-4D52-9416-CAEB6D898D21}
- 40. UPPER BACK (UPPER BACK) {C1DA4294-E0DC-438B-810D-7F91E170244A}
- 41. MID BACK (MID BACK) {F8C515E0-F437-4596-B030-DA1811038CDF}
- 42. LOWER BACK (LOWER BACK) {7E761680-61BC-450F-9905-C45CE760EC45}
- 43. R UPPER ARM (R UPP ARM) {3EF20E8E-F68C-E4F2-8F8C-FAD3B44DDAA4}
- 44. L UPPER ARM (L UPP ARM) {5F915696-1B08-32E9-6B95-49246117D660}
- 45. R LOWER ARM (R LOW ARM) {8DE6E277-EFFF-40DD-9E83-D7295AACB6AA}
- 46. L LOWER ARM (L LOW ARM) {63AA95E8-256D-4390-8CB3-F00A1B3F91CD}
- 47. R ANKLE (R ANKLE) {4E6DF352-F399-4CE9-908C-48132C3D856B}
- 48. L ANKLE (L ANKLE) {D5A692EF-9D7C-473C-AF96-7EF1C01FA2AF}
- 49. UPPER ABDOMEN (UP ABDOM) {234C677B-1CE6-498B-9375-BB5DC3958842}
- 50. R ABDOMEN (R ABDOMEN) {5C45EAA8-96F5-4EC7-8AB5-A9543CD6C414}
- 51. LOWER ABDOMEN (LO ABDOM) {077C9CF5-DCA7-4D94-A7A3-F3C153CEA47D}
- 52. L ABDOMEN (L ABDOMEN) {46B37481-4105-492D-96F0-E2D3895379C1}
- 53. COCCYX (COCCYX) {49D65E9D-3D5B-40E3-B7CF-1003F80B852C}
- 54. HEAD OCCIPITAL (HEAD OCCIP)  
{2066CA47-5A16-447F-9919-24ED435B627F}
- 55. L ANTECUBITAL (L ANTECUBI)  
{C21838CC-BAB5-1D7F-E2A4-A7C663FB28C3}
- 56. R ANTECUBITAL (R ANTECUBI)  
{7FB00A65-B4B2-06A6-BD89-98AC112D8D4E}
- 57. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
- 58. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
- 59. L HEEL (L HEEL) {C9154881-17C7-4921-A0EE-BFB962C59CEF}
- 60. R HEEL (R HEEL) {0020DD20-2EA4-4362-ADF0-07BFA220C60E}
- 61. L ISCHIAL (L ISCHIAL) {DCDAD741-9599-4479-9D55-3279F800DBE6}
- 62. R ISCHIAL (R ISCHIAL) {97E50F6C-CF75-4E59-AFDE-80FC5CE6D9E3}
- 63. L RIBS (L RIBS) {BDFEE74F-B3CC-4AD9-AE49-78B0BF19A219}
- 64. R RIBS (R RIBS) {260FD564-5E58-419B-9F73-62B6037D1572}
- 65. L SCAPULA (L SCAPULA) {3381EC50-771B-41C8-85C8-6EB0D208D54E}
- 66. R SCAPULA (R SCAPULA) {EA657261-4D72-40C2-B70D-2A2E27556A01}
- 67. L TROCHANTER (L TROCHANT) {C91647A4-255D-40E4-A24D-01274E5D799F}
- 68. R TROCHANTER (R TROCHANT) {1A3E21B2-C391-40D2-8059-12290E667D9E}
- 69. PERINEAL (PERINEAL) {1801441C-D5EF-4775-9E58-E4985496F7E0}
- 70. SACRAL (SACRAL) {0CF8A8A7-AB35-B803-4544-0C91645D9D1C}

UNITS:

-----

Abbreviation..... WOU LEN  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Wound Length  
Display Name..... Wound Length  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

CENTIMETER	-1.00 > 0.00 * 200.00 < 400.00
	Dec Prec: 2

=====

TERM: WOUND ODOR

Page: 763

ID: {F71814E1-2D36-4074-B8F4-253AC6BDA1AA}

=====

Abbreviation..... WOU ODOR  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Wound Odor  
Display Name..... Wound Odor  
Help Text..... Wound Odor  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

-----

1. WOUND ODOR NONE (WOU OD NON)  
    {D2E388F7-B531-4D60-A2BE-830F11B9AA75}
2. WOUND ODOR FOUL (WOU OD FOU)  
    {99A74112-EA71-4896-A4DD-911379E5B453}

UNITS:

-----

=====

TERM: WOUND PHOTOGRAPHED  
ID: {398E171E-93A6-4DB3-82E7-5C0D118A80EE}

=====

Page: 764

Abbreviation..... WOU PHOTOG  
Active..... Yes  
Boolean Value False..... No  
Boolean Value True..... Yes  
Data Type..... Boolean  
Description..... Has Wound Been Photographed  
Display Name..... Wound Photographed  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
- |            |  |
|------------|--|
| 1. NO (n)  | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEE649CD} |

UNITS:  
-----

=====

TERM: WOUND TREATMENT  
ID: {4883E272-E8D4-48D4-A7C5-B00F5315B6DC}

=====

Page: 765

Abbreviation..... WOU TREATM  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... Wound Treatment  
Display Name..... Wound Treatment  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 5 entries

- 
1. HYDROCOLLOIDAL DRESSING (HCOLL DRES)  
{D85741EF-0DB3-48CF-BC9A-88C834A72F36}
  2. TRANSPARENT DRESSING (TRANSP DRE)  
{777C7722-AEA8-4A0B-9ADD-1157D37AE539}
  3. SALINE MOIST DRESSING (SAL MOI DR)  
{19F98514-8ABF-461C-9969-C31177FE7DE7}
  4. OPEN TO AIR (OPE TO AIR) {EFB715D0-5666-4281-9034-A0326A14B2BB}
  5. XEROFORM (XEROFORM) {7D7C3CDC-253E-4262-B8D4-3B70DE9AEBE4}

UNITS:  
-----

Abbreviation..... WOU TYPE
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Type of wound
Display Name..... Wound Type
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 29 entries

- 1. ABRASION (ABRASION) {7E962B43-24B1-46D2-A79A-54811117271A}
2. BLISTERING (BLISTERING) {3A4212D5-F2A8-4636-BB7F-F36FA8CED192}
3. BRUISED (BRUISED) {FA5B38A7-7744-186F-055A-8458D696AF33}
4. BURN (BURN) {97B37D7A-28C6-4ABF-8D09-1FC49D5F160F}
5. CELLULITIS (CELLULITIS) {709F8795-0D4A-4E98-AB82-DD3E2AA7EA7B}
6. CRUSHING (CRUSHING) {92D4BAF1-63E3-4FC0-B096-B5CBC4736245}
7. DIABETIC/VASCULAR ULCER (DIA/VAS UL) {E209E53E-EFF1-48F5-9581-75236C89D9B9}
8. HEMATOMA (HEMATOMA) {6A24B7EE-62F3-C952-C1A3-086CDBB1FE68}
9. PREVIOUS PERCUTANEOUS SITE (PRE PER SI) {F2ECCF5C-D91A-40F5-9C09-452C2C81DA66}
10. RASH/EXFOLIATION (RASH/EXFOL) {AB83ABA9-28A7-6484-67E7-4E66B5526942}
11. SKIN TEAR (SKIN TEAR) {2E230C6E-119F-4F9F-8797-2FFA9D76F337}
12. SKIN TRAUMA (SKIN TRAUM) {91C8EE13-2C4A-468C-94EB-C53E0D153EBA}
13. SURGICAL DONOR (SURG DONOR) {F3EE10EA-5776-479B-A58C-F59D6526F90B}
14. SURGICAL GRAFT (SURG GRAFT) {328BA3EA-409B-4C5E-92E7-C6ECE8F983E1}
15. SURGICAL (SURGICAL) {2C6CFEF1-6122-45BA-A9F8-F33B5B3C0E90}
16. VASCULAR (VASCULAR) {60523BF8-6372-4760-881C-A62F85D2A0F6}
17. VASCULAR/NEUROPATHIC (VASC/NEURO) {6B91CEED-4863-46EA-A9AE-3B0D522F1BDC}
18. BITE (Bite) {FC33E0AE-3130-4881-B9EF-9B81FBEC6AC0}
19. BULLET (BULLET) {1570D8E2-D936-494D-AF6E-13838EA34EFD}
20. CONTUSION (Contusion) {5285197A-CC0A-4D12-A1F8-966765C687F1}
21. INCONTENENCE DERMATITIS (INCON DERM) {9F6C5BA2-79C7-4B09-B19A-6DDC444936E5}
22. LACERATION (LACERATION) {F4548BB2-39AF-4E0D-B52F-472EA080D6FF}
23. PENETRATING (PENETRATIN) {00A2D4C8-EDDE-4593-A1F8-A87A836B3379}
24. PRESSURE ULCER (PRE ULCER) {7279639E-498D-4EA3-A3FA-4A6D611E6F48}
25. PUNCTURE (PUNCTURE) {28CF34AF-BD37-4BA4-A59B-C4AFCE2EB8DD}
26. TAPE BURN (TAPE BURN) {DDD75AFD-DA8D-4E73-873F-86D4DF4AF47A}
27. VASCULAR LESION ARTERIAL (VASLES ART) {EA565D3F-D984-4553-B597-DFE844554CE9}
28. VASCULAR LESION MIXED (VASLES MIX) {8684AB2D-9162-4DC2-A5FA-A14558678689}
29. VASCULAR LESION VENOUS (VASLES VEN) {3F8CEC69-5E13-4F26-AB86-34FE2805E3C6}

UNITS:



=====

TERM: WOUND VACUUM ASSISTED CLOSURE  
ID: {04748983-B61B-4121-AF49-FA1FDA1D9646}

=====

Page: 767

Abbreviation..... WOU VAC AS  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... PickList  
Description..... No Description  
Display Name..... Vacuum Assisted Closure Type  
Help Text..... No Help Text  
Multi Select Picklist..... No  
Term Type..... OBSERVATION  
Vuid.....

DATA VALUE: 2 entries

- 
1. INTERMITTENT (INTERMITTE) {DD09ED7C-B977-412F-A678-5130F77D239E}
  2. CONSTANT (CONSTANT) {0FB8E43C-2106-4C37-9430-38675953A85B}

UNITS:  
-----

=====

TERM: WOUND WIDTH

Page: 768

ID: {7EED882C-3CAE-4A54-88FA-D0B380D81658}

=====

Abbreviation..... WOU WIDTH  
Active..... Yes  
Boolean Value False.....  
Boolean Value True.....  
Data Type..... Numeric  
Description..... Wound Width  
Display Name..... Wound Width  
Help Text..... No Help Text  
Multi Select Picklist.....  
Term Type..... OBSERVATION  
Vuid.....

UNITS:

-----

CENTIMETER	-1.00 > 0.00 * 200.00 < 400.00
	Dec Prec: 2



DATA VALUE: 0

-----  
ABBREVIATION.....0  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....0  
HELP\_TEXT.....No Help Text  
ID.....{2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}  
TERM.....0

DATA VALUE: 1

-----  
ABBREVIATION.....1  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....1  
HELP\_TEXT.....No Help Text  
ID.....{8F5D10A8-FE17-47BC-DACA-6972D12CC00A}  
TERM.....1

DATA VALUE: 2

-----  
ABBREVIATION.....2  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....2  
HELP\_TEXT.....No Help Text  
ID.....{38E22B79-A2B6-3401-B195-C955F24E1278}  
TERM.....2

DATA VALUE: 3

-----  
ABBREVIATION.....3  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....3  
HELP\_TEXT.....No Help Text  
ID.....{E714454D-6A67-91F6-2FB7-D8F24F184328}  
TERM.....3

DATA VALUE: 4

-----  
ABBREVIATION.....4  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....4  
HELP\_TEXT.....No Help Text  
ID.....{D2280C7E-DE50-5A8D-C3FD-5862B69C0E1A}  
TERM.....4

DATA VALUE: 4.5

-----  
ABBREVIATION.....4.5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....4.5  
HELP\_TEXT.....No Help Text  
ID.....{EC6C0C72-8C49-4535-92B2-7E543007EE32}  
TERM.....4.5

DATA VALUE: 5

-----  
ABBREVIATION.....5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....5  
HELP\_TEXT.....No Help Text  
ID.....{188EDCE7-0726-4234-D8D2-92C8052CD221}  
TERM.....5

DATA VALUE: 5.5

-----  
ABBREVIATION.....5.5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....5.5  
HELP\_TEXT.....No Help Text  
ID.....{4878F059-EC9D-44FC-9A61-734C291530E7}  
TERM.....5.5

DATA VALUE: 6

-----  
ABBREVIATION.....6  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....6  
HELP\_TEXT.....No Help Text  
ID.....{640CF2F9-3BBF-AB13-185B-CD02A609FB5A}  
TERM.....6

DATA VALUE: 6.5

-----  
ABBREVIATION.....6.5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....6.5  
HELP\_TEXT.....No Help Text  
ID.....{CAFCBDB1-8533-4C10-8DD1-2BD07A4CF6D3}  
TERM.....6.5

DATA VALUE: 7

-----  
ABBREVIATION.....7  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....7  
HELP\_TEXT.....No Help Text  
ID.....{7A010324-D2C2-5D4B-848A-7B8878666D20}  
TERM.....7

DATA VALUE: 7.5

-----  
ABBREVIATION.....7.5  
ACTIVE.....Yes  
DESCRIPTION.....7.5  
DISPLAY\_NAME.....7.5  
HELP\_TEXT.....No Help Text  
ID.....{7517EE9C-7C8B-4AF0-B60A-88F4E24BFA5E}  
TERM.....7.5

DATA VALUE: 8

-----  
ABBREVIATION.....8  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....8  
HELP\_TEXT.....No Help Text  
ID.....{518E7CFC-102E-351E-C7AF-FB575E3C8F4D}  
TERM.....8

DATA VALUE: 8.5

-----  
ABBREVIATION.....8.5  
ACTIVE.....Yes  
DESCRIPTION.....8.5  
DISPLAY\_NAME.....8.5  
HELP\_TEXT.....8.5  
ID.....{96723567-151D-4B5C-9C65-B3EB6D584272}  
TERM.....8.5

DATA VALUE: 9

-----  
ABBREVIATION.....9  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....9  
HELP\_TEXT.....No Help Text  
ID.....{58FF2FD4-569F-F454-9F5B-FE6F9E34D8CC}  
TERM.....9

DATA VALUE: 9.5

-----  
ABBREVIATION.....9.5  
ACTIVE.....Yes  
DESCRIPTION.....9.5  
DISPLAY\_NAME.....9.5  
HELP\_TEXT.....9.5  
ID.....{252555A1-ADB2-4948-8764-C952F2587DC5}  
TERM.....9.5

DATA VALUE: 10

-----  
ABBREVIATION.....10  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....10  
HELP\_TEXT.....No Help-Text  
ID.....{0C827160-18DB-816F-7FAA-A3C175C6D33B}  
TERM.....10

DATA VALUE: 11

-----  
ABBREVIATION.....11  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....11  
HELP\_TEXT.....11  
ID.....{BD51A7B3-FACF-890C-472E-E3B28519CC16}  
TERM.....11

DATA VALUE: 12

-----  
ABBREVIATION.....12  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....12  
HELP\_TEXT.....No Help Text  
ID.....{ACC5CA06-3ABC-476D-EFB0-D82FEEE129C7}  
TERM.....12

DATA VALUE: 12.5

-----  
ABBREVIATION.....12.5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....12.5  
HELP\_TEXT.....No Help Text  
ID.....{97A947AE-2BCC-F75D-1794-6AF9BAC1497B}  
TERM.....12.5

DATA VALUE: 13

-----  
ABBREVIATION.....13  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....13  
HELP\_TEXT.....No Help Text  
ID.....{030242DB-3F82-7517-3DEE-4D50E3E8AE23}  
TERM.....13

DATA VALUE: 13.5

-----  
ABBREVIATION.....13.5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....13.5  
HELP\_TEXT.....No Help Text  
ID.....{34D24648-B2D9-9C6C-6D71-F80E977CE185}  
TERM.....13.5

DATA VALUE: 14

-----  
ABBREVIATION.....14  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....14  
HELP\_TEXT.....No Help Text  
ID.....{EDF1C2B1-123B-011F-9B02-7D71088971C1}  
TERM.....14

DATA VALUE: 14.5

-----  
ABBREVIATION.....14.5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....14.5  
HELP\_TEXT.....No Help Text  
ID.....{373C8F8C-5D36-4AD6-6A76-513064669CB5}  
TERM.....14.5

DATA VALUE: 15

-----  
ABBREVIATION.....15  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....15  
HELP\_TEXT.....No Help Text  
ID.....{FE8C049F-B64C-4991-2F50-3C561BA4825C}  
TERM.....15



DATA VALUE: 15.5

-----  
ABBREVIATION.....15.5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....15.5  
HELP\_TEXT.....No Help Text  
ID.....{1BBF5122-2FFA-6B27-445B-3291C3F83A36}  
TERM.....15.5

DATA VALUE: 16

-----  
ABBREVIATION.....16  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....16  
HELP\_TEXT.....No Help Text  
ID.....{9E4B84FF-430A-39FB-E6AA-E1031814796D}  
TERM.....16

DATA VALUE: 16.5

-----  
ABBREVIATION.....16.5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....16.5  
HELP\_TEXT.....No Help Text  
ID.....{7D4EA1D4-5233-760F-FD2B-A70FB0969815}  
TERM.....16.5

DATA VALUE: 17

-----  
ABBREVIATION.....17  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....17  
HELP\_TEXT.....No Help Text  
ID.....{2E89B9C8-5E80-F3C9-FE66-4396CB63201A}  
TERM.....17

DATA VALUE: 17.5

-----  
ABBREVIATION.....17.5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....17.5  
HELP\_TEXT.....No Help Text  
ID.....{FDE6FFD9-00D0-9C55-1A7F-4AA70CE313BA}  
TERM.....17.5

DATA VALUE: 18

-----  
ABBREVIATION.....18  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....18  
HELP\_TEXT.....No Help Text  
ID.....{F2815231-0FCE-4988-95E9-D9AF46C57086}  
TERM.....18

DATA VALUE: 18.5

-----  
ABBREVIATION.....18.5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....18.5  
HELP\_TEXT.....No Help Text  
ID.....{75ABCB85-4D6E-0BBC-157F-CA69B7737402}  
TERM.....18.5

DATA VALUE: 19

-----  
ABBREVIATION.....19  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....19  
HELP\_TEXT.....No Help Text  
ID.....{478FFB89-4941-1474-04E1-C838B2467439}  
TERM.....19

DATA VALUE: 19.5

-----  
ABBREVIATION.....19.5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....19.5  
HELP\_TEXT.....No Help Text  
ID.....{58F0F98D-299A-8C4E-F50F-AE21B6E696D6}  
TERM.....19.5

DATA VALUE: 20

-----  
ABBREVIATION.....20  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....20  
HELP\_TEXT.....No Help Text  
ID.....{532FBDF8-7F7B-9551-300C-F2CFEF656163}  
TERM.....20

DATA VALUE: 20.5

-----  
ABBREVIATION.....20.5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....20.5  
HELP\_TEXT.....No Help Text  
ID.....{DB1F50D2-AF99-7BB3-88F5-B58C0BDE7987}  
TERM.....20.5

DATA VALUE: 21

-----  
ABBREVIATION.....21  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....21  
HELP\_TEXT.....No Help Text  
ID.....{C6188159-9EF1-DA5C-C015-2142A1205C4C}  
TERM.....21

DATA VALUE: 21.5

-----  
ABBREVIATION.....21.5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....21.5  
HELP\_TEXT.....No Help Text  
ID.....{A259E466-773F-519C-7B5B-1570D1B03B70}  
TERM.....21.5

DATA VALUE: 22

-----  
ABBREVIATION.....22  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....22  
HELP\_TEXT.....No Help Text  
ID.....{DA839156-C7A9-76D7-58CD-7EC2B016EF40}  
TERM.....22

DATA VALUE: 22.5

-----  
ABBREVIATION.....22.5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....22.5  
HELP\_TEXT.....No Help Text  
ID.....{D402D24F-1EF2-1AA0-DD1C-FA8320AA929A}  
TERM.....22.5

DATA VALUE: 23

-----  
ABBREVIATION.....23  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....23  
HELP\_TEXT.....No Help Text  
ID.....{8FB6A8B1-B692-BDC7-43B9-0332F33D826A}  
TERM.....23

DATA VALUE: 23.5

-----  
ABBREVIATION.....23.5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....23.5  
HELP\_TEXT.....No Help Text  
ID.....{DDB76295-7E4A-FEAC-324D-37878F62771A}  
TERM.....23.5

DATA VALUE: 24

-----  
ABBREVIATION.....24  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....24  
HELP\_TEXT.....No Help Text  
ID.....{2BE934CF-40B3-2263-9FA7-6AC9DE0C7C8A}  
TERM.....24

DATA VALUE: 24.5

-----  
ABBREVIATION.....24.5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....24.5  
HELP\_TEXT.....No Help Text  
ID.....{D815830B-5D94-0975-D007-E4FDD966155B}  
TERM.....24.5

DATA VALUE: 25

-----  
ABBREVIATION.....25  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....25  
HELP\_TEXT.....No Help Text  
ID.....{42DD0655-142D-6901-D3FE-66D295BBFCD1}  
TERM.....25

DATA VALUE: 25.5

-----  
ABBREVIATION.....25.5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....25.5  
HELP\_TEXT.....No Help Text  
ID.....{9D704E9E-FB1D-FB83-2142-EDBF3805E8F7}  
TERM.....25.5

DATA VALUE: 26

-----  
ABBREVIATION.....26  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....26  
HELP\_TEXT.....No Help Text  
ID.....{181395F6-29FD-5A41-5AD7-BFC9CF564F3B}  
TERM.....26

DATA VALUE: 26.5

-----  
ABBREVIATION.....26.5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....26.6  
HELP\_TEXT.....No Help Text  
ID.....{1D464967-6CFA-E29E-548D-EF68AFCDEC61}  
TERM.....26.5

DATA VALUE: 27

-----  
ABBREVIATION.....27  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....27  
HELP\_TEXT.....No Help Text  
ID.....{64D25A86-6516-0910-B4B5-644B85378876}  
TERM.....27

DATA VALUE: 27.5

-----  
ABBREVIATION.....27.5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....27.5  
HELP\_TEXT.....No Help Text  
ID.....{8C7BE365-AFD0-3FFD-B8EC-6EC914763CFC}  
TERM.....27.5

DATA VALUE: 28

-----  
ABBREVIATION.....28  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....28  
HELP\_TEXT.....No Help Text  
ID.....{66B5A1F5-5396-10C4-643E-3E8CAB6D6CC8}  
TERM.....28

DATA VALUE: 28.5

-----  
ABBREVIATION.....28.5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....28.5  
HELP\_TEXT.....No Help Text  
ID.....{C99C9298-D4D2-BC02-36D3-F612B6080A87}  
TERM.....28.5

DATA VALUE: 29

-----  
ABBREVIATION.....29  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....29  
HELP\_TEXT.....29  
ID.....{93F62F6F-5EFE-C451-5908-47C895B7177C}  
TERM.....29

DATA VALUE: 30

-----  
ABBREVIATION.....30  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....30  
HELP\_TEXT.....30  
ID.....{1CF40F33-4B6B-FE03-98FB-55DFA8FB4145}  
TERM.....30

DATA VALUE: 31

-----  
ABBREVIATION.....31  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....31  
HELP\_TEXT.....31  
ID.....{5EEDDCF5-7D59-1526-2522-FA59099C8C21}  
TERM.....31

DATA VALUE: 99

-----  
ABBREVIATION.....9  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....99  
HELP\_TEXT.....No Help Text  
ID.....{B013C05B-1AF0-1547-4F03-378883F887C5}  
TERM.....99

DATA VALUE: +1 = SLIGHT

-----  
ABBREVIATION.....+1 SLIGHT  
ACTIVE.....Yes  
DESCRIPTION.....+1 = Slight  
DISPLAY\_NAME.....+1 = Slight  
HELP\_TEXT.....No Help Text  
ID.....{7C0245BA-B261-4B68-A84C-D737CA634230}  
TERM.....+1 = SLIGHT

DATA VALUE: +2 = MODERATE

-----  
ABBREVIATION.....+2 = MODER  
ACTIVE.....Yes  
DESCRIPTION.....+2 = Moderate  
DISPLAY\_NAME.....+2 = Moderate  
HELP\_TEXT.....No Help Text  
ID.....{CCF23A60-068C-443F-AFD0-2D56EEE0446C}  
TERM.....+2 = MODERATE

DATA VALUE: +3 = MARKED

-----  
ABBREVIATION.....+3 = MARKE  
ACTIVE.....Yes  
DESCRIPTION.....+3 = Marked  
DISPLAY\_NAME.....+3 = Marked  
HELP\_TEXT.....No Help Text  
ID.....{41ACC586-0C68-4780-8B2D-F36B21FCEDB2}  
TERM.....+3 = MARKED

DATA VALUE: +4 = VERY MARKED

-----  
ABBREVIATION.....+4 = VY MA  
ACTIVE.....Yes  
DESCRIPTION.....+4 = Very Marked  
DISPLAY\_NAME.....+4 = Very Marked  
HELP\_TEXT.....No Help Text  
ID.....{55EF9D42-B4F8-4468-BA1B-98BCD4075E58}  
TERM.....+4 = VERY MARKED

DATA VALUE: 0 - ABSENT

-----  
ABBREVIATION.....0-a  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....0 - Absent  
HELP\_TEXT.....No Help Text  
ID.....{75F14311-F113-963C-4F05-B8AE17F044A9}  
TERM.....0 - ABSENT

DATA VALUE: 0 - ALERT

-----  
ABBREVIATION.....Alert  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....0 - Alert  
HELP\_TEXT.....  
ID.....{A144FB6A-D6FA-4E0F-9F3C-3F1756464D05}  
TERM.....0 - ALERT

DATA VALUE: 0 - HR 51-100

-----  
ABBREVIATION.....HR 51-100  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....0 - HR 51-100  
HELP\_TEXT.....  
ID.....{6A262124-7FE8-4520-9762-87D0A3DB8AE9}  
TERM.....0 - HR 51-100

DATA VALUE: 0 - INDEPENDENT - NO HELP OR STAFF OVERSIGHT

-----  
ABBREVIATION.....0 Ind  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....0 - Independent - No Help or Staff  
Oversight  
HELP\_TEXT.....  
ID.....{4A125BF9-A2C1-4503-9182-BF3E2A708547}  
TERM.....0 - INDEPENDENT - NO HELP OR STAFF  
OVERSIGHT

DATA VALUE: 0 - INDEPENDENT - NO HELP PROVIDED

-----  
ABBREVIATION.....0 ind  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....0 - Independent - No Help Provided  
HELP\_TEXT.....  
ID.....{07AB4004-1793-410A-9872-00BB0B322DFE}  
TERM.....0 - INDEPENDENT - NO HELP PROVIDED



DATA VALUE: 0 - NO ANXIETY AT EASE

-----  
ABBREVIATION.....CIWA AX 0  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....0 - No Anxiety, At Ease  
HELP\_TEXT.....No Help Text  
ID.....{DA61A13E-6443-4E8A-9C76-762FEACDB4E7}  
TERM.....0 - NO ANXIETY AT EASE

DATA VALUE: 0 - NO NAUSEA/VOMITING

-----  
ABBREVIATION.....CIWA NV 0  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....0 - No Nausea / Vomiting  
HELP\_TEXT.....No Help Text  
ID.....{0569555E-E7B3-49E3-9A18-68D23CFB2B8D}  
TERM.....0 - NO NAUSEA/VOMITING

DATA VALUE: 0 - NO SETUP OR PHYSICAL HELP FROM STAFF

-----  
ABBREVIATION.....0 No Setup  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....0 - No Setup or Physical Help from  
Staff  
HELP\_TEXT.....  
ID.....{730049CC-4CD5-449E-A0E8-389AB61A02BD}  
TERM.....0 - NO SETUP OR PHYSICAL HELP FROM  
STAFF

DATA VALUE: 0 - NO SWEAT VISIBLE

-----  
ABBREVIATION.....CIWA PS 0  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....0 - No Sweat Visible  
HELP\_TEXT.....No Help Text  
ID.....{7C83E616-4414-4766-AAF5-5F5442FFBF09}  
TERM.....0 - NO SWEAT VISIBLE

DATA VALUE: 0 - NO TREMOR

-----  
ABBREVIATION.....CIWA TR 0  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....0 - No Tremor  
HELP\_TEXT.....No Help Text  
ID.....{2D1A5B39-4C5D-4345-84B4-DDCBC8E25B4D}  
TERM.....0 - NO TREMOR

DATA VALUE: 0 - NONE

-----  
ABBREVIATION.....CIWA TC 0  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....0 - None  
HELP\_TEXT.....No Help Text  
ID.....{5F9F33C7-EFCE-4AC1-BCA3-1967F68EB597}  
TERM.....0 - NONE

DATA VALUE: 0 - NORMAL ACTIVITY

-----  
ABBREVIATION.....CIWA AG 0  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....0 - Normal Activity  
HELP\_TEXT.....No Help Text  
ID.....{9FF32084-A432-467C-B0A2-C8BC9F3D6D13}  
TERM.....0 - NORMAL ACTIVITY

DATA VALUE: 0 - ORIENTED AND CAN DO SERIAL ADDITIONS

-----  
ABBREVIATION.....CIWA OR 0  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....0 - Oriented And Can Do Serial  
Additions  
HELP\_TEXT.....No Help Text  
ID.....{D8015F54-6E9E-457A-955C-38565BDDDBB6}  
TERM.....0 - ORIENTED AND CAN DO SERIAL  
ADDITIONS

DATA VALUE: 0 - RR 9-14

-----  
ABBREVIATION.....RR 9-14  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....0 - RR 9-14  
HELP\_TEXT.....  
ID.....{F32D69EF-7AC9-4A59-A63B-6E941BF5F06B}  
TERM.....0 - RR 9-14

DATA VALUE: 0 - SBP 101-199

-----  
ABBREVIATION.....SBP101-199  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....0 - SBP 101-199  
HELP\_TEXT.....  
ID.....{ED9B73BA-8B9E-49C0-83A2-5B2397028564}  
TERM.....0 - SBP 101-199

DATA VALUE: 0 - TEMP 36.1-38 C

-----  
ABBREVIATION.....36.1-38 C  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....0 - Temp 36.1-38 C  
HELP\_TEXT.....  
ID.....{87CCCC92-0761-47FC-BD44-C37590832EA7}  
TERM.....0 - TEMP 36.1-38 C

DATA VALUE: 0 = 0 EXTREMITIES

-----  
ABBREVIATION.....ALD ACT 0  
ACTIVE.....Yes  
DESCRIPTION.....0 = 0 Extremities  
DISPLAY\_NAME.....0 = Moves 0 Extremities  
HELP\_TEXT.....0 = 0 Extremities  
ID.....{6E248A9A-B509-4BD0-A01D-9ED5CE5B8FC9}  
TERM.....0 = 0 EXTREMITIES

DATA VALUE: 0 = <90%

-----  
ABBREVIATION.....ALD O2 0  
ACTIVE.....Yes  
DESCRIPTION.....0 = <90% with O2  
DISPLAY\_NAME.....0 = <90% with O2  
HELP\_TEXT.....0 = <90% with O2  
ID.....{30042730-F600-44FD-A3D0-53172B70082B}  
TERM.....0 = <90%

DATA VALUE: 0 = ABSENT PULSE

-----  
ABBREVIATION.....0p  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....0 = Absent Pulse  
HELP\_TEXT.....No Help Text  
ID.....{4253CB74-51AE-BBFA-214F-937762F7638C}  
TERM.....0 = ABSENT PULSE

DATA VALUE: 0 = APNEIC

-----  
ABBREVIATION.....ALD RESP 0  
ACTIVE.....Yes  
DESCRIPTION.....0 = Apneic  
DISPLAY\_NAME.....0 = Apneic  
HELP\_TEXT.....0 = Apneic  
ID.....{649C42B0-6E0A-4E86-AAC5-70A8C20EC511}  
TERM.....0 = APNEIC

DATA VALUE: 0 = BP +/- 50% PREOP

-----  
ABBREVIATION.....ALD CIR 0  
ACTIVE.....Yes  
DESCRIPTION.....0 = BP +/- 50% Preop  
DISPLAY\_NAME.....0 = BP +/- 50% Preop  
HELP\_TEXT.....0 = BP +/- 50% Preop  
ID.....{675F2CD1-DD7B-473B-B8EC-C2BF4AC3B2}  
TERM.....0 = BP +/- 50% PREOP

DATA VALUE: 0 = NOT RESPONDING

-----  
ABBREVIATION.....ALD CON 0  
ACTIVE.....Yes  
DESCRIPTION.....0 = Not Responding  
DISPLAY\_NAME.....0 = Not Responding  
HELP\_TEXT.....0 = Not Responding  
ID.....{E8765215-B0B0-430C-8C89-AB969A6AD12E}  
TERM.....0 = NOT RESPONDING

DATA VALUE: 1 - HR 101-110

-----  
ABBREVIATION.....HR 101-110  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....1 - HR 101-110  
HELP\_TEXT.....  
ID.....{FCB6C37E-1BB4-418F-8203-0916B4056DEB}  
TERM.....1 - HR 101-110

DATA VALUE: 1 - BARELY PERCEPTIBLE SWEATING / PALMS MOIST

-----  
ABBREVIATION.....CIWA PS 1  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....1 - Barely Perceptible Sweating /  
Palms Moist  
HELP\_TEXT.....No Help Text  
ID.....{93B5B833-14B5-4ADD-97E8-7BF14A8BA01A}  
TERM.....1 - BARELY PERCEPTIBLE SWEATING /  
PALMS MOIST

DATA VALUE: 1 - CANNOT DO SERIAL ADDITIONS

-----  
ABBREVIATION.....CIWA OR 1  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....1 - Cannot Do Serial Additions Or Is  
Uncertain About Date  
HELP\_TEXT.....No Help Text  
ID.....{F4446E01-E326-420E-8EA8-04451AA9D963}  
TERM.....1 - CANNOT DO SERIAL ADDITIONS

DATA VALUE: 1 - HR 40-50

-----  
ABBREVIATION.....HR 40-50  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....1 - HR 40-50  
HELP\_TEXT.....  
ID.....{EC915B1F-BED1-497D-AD10-4FD0C2ADEDCE}  
TERM.....1 - HR 40-50

DATA VALUE: 1 - LESS THAN 45MLS/HR

-----  
ABBREVIATION.....< 45mls/hr  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....1 - Less than 45mls/hr  
HELP\_TEXT.....  
ID.....{3AED380E-00F9-4C29-B303-52E08F8911B3}  
TERM.....1 - LESS THAN 45MLS/HR

DATA VALUE: 1 - MILD NAUSEA / NO VOMITING

-----  
ABBREVIATION.....CIWA NV 1  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....1 - Mild Nausea / No Vomiting  
HELP\_TEXT.....No Help Text  
ID.....{582CBFB6-60A0-4F6D-936A-9D2D6F871B4A}  
TERM.....1 - MILD NAUSEA / NO VOMITING

DATA VALUE: 1 - MILDLY ANXIOUS

-----  
ABBREVIATION.....CIWA AX 1  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....1 - Mildly Anxious  
HELP\_TEXT.....No Help Text  
ID.....{401113EE-5CE7-4AE9-80C1-1C3252E706D0}  
TERM.....1 - MILDLY ANXIOUS

DATA VALUE: 1 - NEW AGITATION CONFUSION

-----  
ABBREVIATION.....New Agitat  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....1 - New Agitation Confusion  
HELP\_TEXT.....  
ID.....{F458199B-A547-419B-A691-F15E7CD61F8C}  
TERM.....1 - NEW AGITATION CONFUSION

DATA VALUE: 1 - NO RESPONSE

-----  
ABBREVIATION.....1-nr  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....1 - No Response  
HELP\_TEXT.....No Help Text  
ID.....{26D3168D-40EA-C918-D321-54BED4837311}  
TERM.....1 - NO RESPONSE

DATA VALUE: 1 - NOT VISIBLE BUT CAN BE FELT FINGERTIP

-----  
ABBREVIATION.....CIWA TR 1  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....1 - Not Visible But Can Be Felt  
                  Fingertip To Fingertip  
HELP\_TEXT.....No Help Text  
ID.....{31833D55-9018-47EE-8341-A0863744CEF8}  
TERM.....1 - NOT VISIBLE BUT CAN BE FELT  
                  FINGERTIP

DATA VALUE: 1 - RESPONDS TO VOICE

-----  
ABBREVIATION.....Voice  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....1 - Responds to Voice  
HELP\_TEXT.....  
ID.....{69C7DA2C-893A-40AB-969C-E003B5D2F05D}  
TERM.....1 - RESPONDS TO VOICE

DATA VALUE: 1 - RR 15-20

-----  
ABBREVIATION.....RR 15-20  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....1 - RR 15-20  
HELP\_TEXT.....  
ID.....{1CF2487B-30C6-46AB-A051-B4499F825DC2}  
TERM.....1 - RR 15-20

DATA VALUE: 1 - SBP 81-100

-----  
ABBREVIATION.....SBP 81-100  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....1 - SBP 81-100  
HELP\_TEXT.....  
ID.....{93D627C1-2FF3-4FCC-B0AF-83A54183AE3D}  
TERM.....1 - SBP 81-100

DATA VALUE: 1 - SETUP HELP ONLY

-----  
ABBREVIATION.....1 Setup  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....1 - Setup Help Only  
HELP\_TEXT.....  
ID.....{5494773C-99FE-4A96-B12E-D734188517CE}  
TERM.....1 - SETUP HELP ONLY

DATA VALUE: 1 - SUPERVISION - OVERSIGHT

-----  
ABBREVIATION.....1 Superv  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....1 - Supervision - Oversight  
HELP\_TEXT.....  
ID.....{7400DA9E-E89A-44F4-8516-CB6D5931D97D}  
TERM.....1 - SUPERVISION - OVERSIGHT

DATA VALUE: 1 - SUPERVISION - OVERSIGHT HELP ONLY

-----  
ABBREVIATION.....1 Supervi  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....1 - Supervision - Oversight Help Only  
HELP\_TEXT.....  
ID.....{F7DFB6E8-F550-4072-B693-93BF96A23834}  
TERM.....1 - SUPERVISION - OVERSIGHT HELP ONLY

DATA VALUE: 1 - TEMP 35.1-36 C

-----  
ABBREVIATION.....35.1-36 C  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....1 - Temp 35.1-36 C  
HELP\_TEXT.....  
ID.....{AF9BC2E6-9C52-4030-8422-28DCBDD43DD9}  
TERM.....1 - TEMP 35.1-36 C

DATA VALUE: 1 - TEMP 38.1-38.5 C

-----  
ABBREVIATION.....38.1-38.5C  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....1 - Temp 38.1-38.5 C  
HELP\_TEXT.....  
ID.....{21326CB3-A848-4984-BB34-5594AD7E5D32}  
TERM.....1 - TEMP 38.1-38.5 C

DATA VALUE: 1 - VERY MILD

-----  
ABBREVIATION.....CIWA HA 1  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....1 - Very Mild  
HELP\_TEXT.....No Help Text  
ID.....{F789C1EA-CC7A-4368-87BA-08AFFE492EE6}  
TERM.....1 - VERY MILD

DATA VALUE: 1 - VERY MILD HARSHNESS OR ABILITY TO FRIGHTEN

-----  
ABBREVIATION.....CIWA AU 1  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....1 - Very Mild Harshness Or Ability To  
Frighten  
HELP\_TEXT.....No Help Text  
ID.....{E7D03A9C-B975-48B2-8CA6-2FC6616C6AF6}  
TERM.....1 - VERY MILD HARSHNESS OR ABILITY TO  
FRIGHTEN

DATA VALUE: 1 - VERY MILD ITCHING PINS/NEEDLES BURNING

-----  
ABBREVIATION.....CIWA TC 1  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....1 - Very Mild Itching, Pins/Needles,  
Burning, Numbness  
HELP\_TEXT.....No Help Text  
ID.....{00642E0D-E1CF-4EE0-91F4-CBAF209AA54C}  
TERM.....1 - VERY MILD ITCHING PINS/NEEDLES  
BURNING

DATA VALUE: 1 - VERY MILD SENSITIVITY

-----  
ABBREVIATION.....CIWA VS 1  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....1 - Very Mild Sensitivity  
HELP\_TEXT.....No Help Text  
ID.....{8BF96CF1-3EA7-468C-838A-284430EA850D}  
TERM.....1 - VERY MILD SENSITIVITY

DATA VALUE: 1 - WEAK

-----  
ABBREVIATION.....1 - WEAK  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....1 - Weak  
HELP\_TEXT.....No Help Text  
ID.....{A41570DC-2F35-DDF6-977B-7C262150DFA2}  
TERM.....1 - WEAK



DATA VALUE: 1 = 2 EXTREMITIES

-----  
ABBREVIATION.....ALD ACT 1  
ACTIVE.....Yes  
DESCRIPTION.....1 = 2 Extremities  
DISPLAY\_NAME.....1 = Moves 2 Extremities  
HELP\_TEXT.....1 = 2 Extremities  
ID.....{6DCB4762-72CE-48C2-9AFE-B19969ED1DE0}  
TERM.....1 = 2 EXTREMITIES

DATA VALUE: 1 = AROUSABLE

-----  
ABBREVIATION.....ALD CON 1  
ACTIVE.....Yes  
DESCRIPTION.....1 = Arousable with minimal stimulation  
DISPLAY\_NAME.....1 = Arousable with minimal stimulation  
HELP\_TEXT.....1 = Arousable with minimal stimulation  
ID.....{B2A07804-9829-42D0-9C49-A61FD51EEAA0}  
TERM.....1 = AROUSABLE

DATA VALUE: 1 = BP +/- 20-49% PREOP

-----  
ABBREVIATION.....ALD CIR 1  
ACTIVE.....Yes  
DESCRIPTION.....1 = BP +/- 20-49% Preop  
DISPLAY\_NAME.....1 = BP +/- 20-49% Preop  
HELP\_TEXT.....1 = BP +/- 20-49% Preop  
ID.....{D4A6100D-1228-4860-9BD6-B43F30D4FF94}  
TERM.....1 = BP +/- 20-49% PREOP

DATA VALUE: 1 = DYSPNEA

-----  
ABBREVIATION.....ALD RESP 1  
ACTIVE.....Yes  
DESCRIPTION.....1 = Dyspnea, shallow or limited  
DISPLAY\_NAME.....1 = Dyspnea, shallow or limited  
HELP\_TEXT.....1 = Dyspnea, shallow or limited  
ID.....{3DAE60A2-57DA-4121-A507-BC43497C7B67}  
TERM.....1 = DYSPNEA

DATA VALUE: 1 = NEEDS O2

-----  
ABBREVIATION.....ALD O2 1  
ACTIVE.....Yes  
DESCRIPTION.....1 = Needs O2 to maintain >90%  
DISPLAY\_NAME.....1 = Needs O2 to maintain >90%  
HELP\_TEXT.....1 = Needs O2 to maintain >90%  
ID.....{5D97F6BB-7B02-4644-8CE5-014C89826731}  
TERM.....1 = NEEDS O2

DATA VALUE: 1 HOUR

-----  
ABBREVIATION.....1 HOUR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....1 Hour  
HELP\_TEXT.....No Help Text  
ID.....{319BA174-1624-483C-B666-4998127B0800}  
TERM.....1 HOUR

DATA VALUE: 1+=PALPABLE, THREADY, EASILY OBLITERATED

-----  
ABBREVIATION.....1teo  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....1+=Palpable, Thready, Easily  
Obliterated  
HELP\_TEXT.....No Help Text  
ID.....{3AE92556-1337-BFE2-C3EE-9676327935C3}  
TERM.....1+=PALPABLE, THREADY, EASILY  
OBLITERATED

DATA VALUE: 1- MILD ANXIOUS

-----  
ABBREVIATION.....CIWA AX 1  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....1 - Mild Anxious  
HELP\_TEXT.....No Help Text  
ID.....{CAB6FB3A-4088-4D2C-AE1E-623352511D73}  
TERM.....1- MILD ANXIOUS

DATA VALUE: 1- SOMEWHAT MORE THAN NORMAL ACTIVITY

-----  
ABBREVIATION.....CIWA AG 1  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....1 - Somewhat More Than Normal Activity  
HELP\_TEXT.....No Help Text  
ID.....{497F7A3D-72C8-4986-A1E6-74EA0E5704F9}  
TERM.....1- SOMEWHAT MORE THAN NORMAL ACTIVITY

DATA VALUE: 1800 CALORIE CONTROLLED

-----  
ABBREVIATION.....1800 CAL C  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....1800 Calorie Controlled Diet  
HELP\_TEXT.....No Help Text  
ID.....{C564D7FE-EC2C-401E-A77D-79CC5FC2B2F8}  
TERM.....1800 CALORIE CONTROLLED

DATA VALUE: 1800 DIABETIC

-----  
ABBREVIATION.....1800 DIAB  
ACTIVE.....Yes  
DESCRIPTION.....1800 Calorie Diabetic Diet  
DISPLAY\_NAME.....1800 Diabetic  
HELP\_TEXT.....No Help Text  
ID.....{FC71C7D7-F24E-4316-A455-8674804BE092}  
TERM.....1800 DIABETIC

DATA VALUE: 1:1

-----  
ABBREVIATION.....1:1  
ACTIVE.....Yes  
DESCRIPTION.....1 to 1 ratio  
DISPLAY\_NAME.....1:1  
HELP\_TEXT.....No Help Text  
ID.....{5C07BBD0-FDB9-0CD7-F2C2-616913E0281B}  
TERM.....1:1

DATA VALUE: 1:1 OBSERVATION

-----  
ABBREVIATION.....1o  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....1:1 Observation  
HELP\_TEXT.....No Help Text  
ID.....{B1A7AF4B-9BD2-5F48-79B3-F35861AA5C45}  
TERM.....1:1 OBSERVATION

DATA VALUE: 1:2

-----  
ABBREVIATION.....1:2  
ACTIVE.....Yes  
DESCRIPTION.....1 to 2 ratio  
DISPLAY\_NAME.....1:2  
HELP\_TEXT.....No Help Text  
ID.....{287CDCE0-36F3-65F6-E6A5-9EF7CF8976AE}  
TERM.....1:2

DATA VALUE: 1:3

-----  
ABBREVIATION.....1:3  
ACTIVE.....Yes  
DESCRIPTION.....1 to 3 ratio  
DISPLAY\_NAME.....1:3  
HELP\_TEXT.....No Help Text  
ID.....{062B4F93-8F48-7373-F679-11F57CE98BB3}  
TERM.....1:3

DATA VALUE: 1:4

-----  
ABBREVIATION.....1:4  
ACTIVE.....Yes  
DESCRIPTION.....1 to 4 ratio  
DISPLAY\_NAME.....1:4  
HELP\_TEXT.....No Help Text  
ID.....{3F1B75FF-C168-D413-0BDF-6581E7A4EEAD}  
TERM.....1:4

DATA VALUE: 1:8

-----  
ABBREVIATION.....1:8  
ACTIVE.....Yes  
DESCRIPTION.....1 to 8 ratio  
DISPLAY\_NAME.....1:8  
HELP\_TEXT.....No Help Text  
ID.....{94790C4C-3AF2-6F4C-58B8-87021CA97A53}  
TERM.....1:8

DATA VALUE: 1A + 2V

-----  
ABBREVIATION.....1+2  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....1A + 2V  
HELP\_TEXT.....No Help Text  
ID.....{FE47CEC3-C3A5-0F5C-3EE8-50F30FD4E0C4}  
TERM.....1A + 2V

DATA VALUE: 1A + 3V

-----  
ABBREVIATION.....1+3  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....1A + 3V  
HELP\_TEXT.....No Help Text  
ID.....{CEDF220D-E593-B8F9-DB08-D32800AF3F5E}  
TERM.....1A + 3V

DATA VALUE: 1A +1V

-----  
ABBREVIATION.....1+  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....1A +1V  
HELP\_TEXT.....No Help Text  
ID.....{31711946-9F1B-0AD9-BDBE-846EDBC86AE2}  
TERM.....1A +1V

DATA VALUE: 2 - ABNORMAL EXTENSION FROM PAIN

-----  
ABBREVIATION.....2-aefp  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....2 - Abnormal Extension From Pain  
HELP\_TEXT.....No Help Text  
ID.....{9B4951D8-B271-DFC3-0A00-1DC9984039E3}  
TERM.....2 - ABNORMAL EXTENSION FROM PAIN

DATA VALUE: 2 - BEADS OF SWEAT OBVIOUS ON FOREHEAD

-----  
ABBREVIATION.....CIWA PS 2  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....2 - Beads Of Sweat Obvious On Forehead  
HELP\_TEXT.....No Help Text  
ID.....{38EA4194-F5F4-42CB-8DE2-715167888F86}  
TERM.....2 - BEADS OF SWEAT OBVIOUS ON FOREHEAD

DATA VALUE: 2 - DISORIENTED FOR DATE < 2 CALENDAR DAYS

-----  
ABBREVIATION.....CIWA OR 2  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....2 - Disoriented For Date By No More  
Than 2 Calendar Days  
HELP\_TEXT.....No Help Text  
ID.....{01B894B0-3E70-47FA-820A-83CA48B93EF9}  
TERM.....2 - DISORIENTED FOR DATE < 2 CALENDAR  
DAYS

DATA VALUE: 2 - HR 0-39

-----  
ABBREVIATION.....HR 0-39  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....2 - HR 0-39  
HELP\_TEXT.....  
ID.....{1243C1C9-C72F-4765-8D19-D689E1FC44F5}  
TERM.....2 - HR 0-39

DATA VALUE: 2 - HR 111-129

-----  
ABBREVIATION.....HR 111-129  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....2 - HR 111-129  
HELP\_TEXT.....  
ID.....{6E1B9533-ECBE-4E2E-90D8-50C9C72AE0D3}  
TERM.....2 - HR 111-129

DATA VALUE: 2 - INCOMPREHENSIBLE SOUNDS

-----  
ABBREVIATION.....2-is  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....2 - Incomprehensible Sounds  
HELP\_TEXT.....No Help Text  
ID.....{F2261C3D-79A3-055B-702D-446B7A6A3113}  
TERM.....2 - INCOMPREHENSIBLE SOUNDS

DATA VALUE: 2 - INTERMITTENT NAUSEA WITH DRY HEAVES

-----  
ABBREVIATION.....CIWA NV 2  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....2 - Intermittent Nausea With Dry Heaves  
HELP\_TEXT.....No Help Text  
ID.....{E92221EE-1298-4D9D-9EFD-35F615118C96}  
TERM.....2 - INTERMITTENT NAUSEA WITH DRY HEAVES

DATA VALUE: 2 - LESS THAN 30MLS/HR

-----  
ABBREVIATION.....< 30mls/hr  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....2 - Less than 30mls/hr  
HELP\_TEXT.....  
ID.....{A85A05F7-3443-4E62-ACF4-3334E3EF5675}  
TERM.....2 - LESS THAN 30MLS/HR

DATA VALUE: 2 - LIMITED ASSISTANCE

-----  
ABBREVIATION.....2 Limited  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....2 - Limited Assistance  
HELP\_TEXT.....  
ID.....{35934E53-08EB-46E0-8136-96FD6174ACF1}  
TERM.....2 - LIMITED ASSISTANCE

DATA VALUE: 2 - MILD

-----  
ABBREVIATION.....CIWA HA 2  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....2 - Mild  
HELP\_TEXT.....No Help Text  
ID.....{6B752397-0A4F-4176-90E7-E220DB2BF1C4}  
TERM.....2 - MILD

DATA VALUE: 2 - MILD HARSHNESS OR ABILITY TO FRIGHTEN

-----  
ABBREVIATION.....CIWA AU 2  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....2 - Mild Harshness Or Ability To  
                  Frighten  
HELP\_TEXT.....No Help-Text  
ID.....{71204969-E614-46A4-ADA5-39FAF81A79B5}  
TERM.....2 - MILD HARSHNESS OR ABILITY TO  
                  FRIGHTEN

DATA VALUE: 2 - MILD ITCHING PINS/NEEDLES BURNING NUMBNESS

-----  
ABBREVIATION.....CIWA TC 2  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....2 - Mild Itching, Pins/Needles,  
                  Burning, Numbness  
HELP\_TEXT.....No Help Text  
ID.....{1A591F17-CFB3-4B24-8AD6-ECF09318E2B6}  
TERM.....2 - MILD ITCHING PINS/NEEDLES BURNING  
                  NUMBNESS

DATA VALUE: 2 - MILD SENSITIVITY

-----  
ABBREVIATION.....CIWA VS 2  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....2 - Mild Sensitivity  
HELP\_TEXT.....No Help Text  
ID.....{ED600249-E5E8-4020-955E-A5CDA18B5DE9}  
TERM.....2 - MILD SENSITIVITY

DATA VALUE: 2 - MODERATE WITH PATIENTS ARMS EXTENDED

-----  
ABBREVIATION.....CIWA TR 2  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....2 - Moderate With Patient's Arms  
                  Extended  
HELP\_TEXT.....No Help Text  
ID.....{B4E7B7FB-B88E-435A-904E-C613AC937323}  
TERM.....2 - MODERATE WITH PATIENTS ARMS  
                  EXTENDED

DATA VALUE: 2 - NORMAL

-----  
ABBREVIATION.....2-n  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....2 - Normal  
HELP\_TEXT.....No Help Text  
ID.....{25FAADFC-D280-78AC-4F01-CFCF34ECD3FF}  
TERM.....2 - NORMAL

DATA VALUE: 2 - ONE PERSON PHYSICAL ASSIST

-----  
ABBREVIATION.....2 one asst  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....2 - One Person Physical Assist  
HELP\_TEXT.....  
ID.....{C675F08C-1BEE-4CDC-AE01-5CE429AB2105}  
TERM.....2 - ONE PERSON PHYSICAL ASSIST

DATA VALUE: 2 - PHYSICAL HELP LIMITED TO TRANSFER ONLY

-----  
ABBREVIATION.....2 Phys Tra  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....2 - Physical Help Limited to Transfer  
Only  
HELP\_TEXT.....  
ID.....{0D087E2B-ACB7-47DD-A1F8-30C0F89CE2DE}  
TERM.....2 - PHYSICAL HELP LIMITED TO TRANSFER  
ONLY

DATA VALUE: 2 - RESPONDS TO PAIN

-----  
ABBREVIATION.....Pain  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....2 - Responds to Pain  
HELP\_TEXT.....  
ID.....{AE6B5DEC-6510-4FEE-8461-D2E930534ACE}  
TERM.....2 - RESPONDS TO PAIN

DATA VALUE: 2 - RR 0-8

-----  
ABBREVIATION.....RR 0-8  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....2 - RR 0-8  
HELP\_TEXT.....  
ID.....{73DC19E6-865C-4161-9CBE-3A8D44BC4D84}  
TERM.....2 - RR 0-8

DATA VALUE: 2 - RR 21-29

-----  
ABBREVIATION.....RR 21-29  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....2 - RR 21-29  
HELP\_TEXT.....  
ID.....{72DF48C2-B8F3-4ED4-91D1-87157707C3EC}  
TERM.....2 - RR 21-29



DATA VALUE: 2 - SBP 200 OR MORE

-----  
ABBREVIATION.....SBP 200+  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....2 - SBP 200 or more  
HELP\_TEXT.....  
ID.....{C3B80866-749C-402A-B7CA-A1B797D36885}  
TERM.....2 - SBP 200 OR MORE

DATA VALUE: 2 - SBP 71-80

-----  
ABBREVIATION.....SBP 71-80  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....2 - SBP 71-80  
HELP\_TEXT.....  
ID.....{063CE20D-435C-4176-9A8B-3FC3DB16FF11}  
TERM.....2 - SBP 71-80

DATA VALUE: 2 - TEMP 35.0 C OR LESS

-----  
ABBREVIATION.....35.0 C or<  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....2 - Temp 35.0 C or less  
HELP\_TEXT.....  
ID.....{D7987CD3-1866-4AE9-A2EA-C96F2CFE2850}  
TERM.....2 - TEMP 35.0 C OR LESS

DATA VALUE: 2 - TEMP MORE THAN 38.6 C

-----  
ABBREVIATION.....> 38.6 C  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....2 - Temp more than 38.6 C  
HELP\_TEXT.....  
ID.....{B45FAD7F-9438-488D-BDE1-8A64A84E51F6}  
TERM.....2 - TEMP MORE THAN 38.6 C

DATA VALUE: 2 - TO PAIN

-----  
ABBREVIATION.....2-tp  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....2 - To Pain  
HELP\_TEXT.....No Help Text  
ID.....{49CE3C40-5FC2-002E-1C89-3E5DFA4CCF0F}  
TERM.....2 - TO PAIN

DATA VALUE: 2 = 4 EXTREMITIES

-----  
ABBREVIATION.....ALD ACT 2  
ACTIVE.....Yes  
DESCRIPTION.....2 = 4 Extremities  
DISPLAY\_NAME.....2 = Moves 4 Extremities  
HELP\_TEXT.....2 = 4 Extremities  
ID.....{3434DF4B-5CB6-417F-B5BE-CEF2D63FE8E0}  
TERM.....2 = 4 EXTREMITIES

DATA VALUE: 2 = BP +/- 20% PREOP

-----  
ABBREVIATION.....ALD CIR 2  
ACTIVE.....Yes  
DESCRIPTION.....2 = BP +/- 20% Preop  
DISPLAY\_NAME.....2 = BP +/- 20% Preop  
HELP\_TEXT.....2 = BP +/- 20% Preop  
ID.....{04B67413-7A88-409B-8005-471A7C84B879}  
TERM.....2 = BP +/- 20% PREOP

DATA VALUE: 2 = BREATHES FREELY

-----  
ABBREVIATION.....ALD RESP 2  
ACTIVE.....Yes  
DESCRIPTION.....2 = Breathes freely, coughs freely  
DISPLAY\_NAME.....2 = Breathes freely, coughs freely  
HELP\_TEXT.....2 = Breathes freely, coughs freely  
ID.....{278CD44C-538E-4015-9D98-D9312A7039BA}  
TERM.....2 = BREATHES FREELY

DATA VALUE: 2 = FULLY AWAKE

-----  
ABBREVIATION.....ALD CON 2  
ACTIVE.....Yes  
DESCRIPTION.....2 = Fully awake  
DISPLAY\_NAME.....2 = Fully awake  
HELP\_TEXT.....2 = Fully awake  
ID.....{903D0F56-ADEC-4527-BA66-C8A568E3EAC0}  
TERM.....2 = FULLY AWAKE

DATA VALUE: 2 = ON ROOM AIR

-----  
ABBREVIATION.....ALD O2 2  
ACTIVE.....Yes  
DESCRIPTION.....2 = On room air  
DISPLAY\_NAME.....2 = >92% on room air  
HELP\_TEXT.....2 = On room air  
ID.....{D58C5B46-CCD3-4E6D-86E8-7A1CFA92BC1C}  
TERM.....2 = ON ROOM AIR

DATA VALUE: 2+=NORMAL, NOT EASILY OBLITERATED

-----  
ABBREVIATION.....2neo  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....2+=Normal, Not Easily Obliterated  
HELP\_TEXT.....No Help Text  
ID.....{50484DA4-EC15-47E5-3923-268C4B5D6428}  
TERM.....2+=NORMAL, NOT EASILY OBLITERATED

DATA VALUE: 2- MILD HARSHNESS OR ABILITY TO FRIGHTEN

-----  
ABBREVIATION.....CIWA AU 2  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....2 - Mild Harshness Or Ability To  
                  Frighten  
HELP\_TEXT.....No Help Text  
ID.....{42DAFB5B-E9AF-4558-95D4-2FFF7BE779F1}  
TERM.....2- MILD HARSHNESS OR ABILITY TO  
                  FRIGHTEN

DATA VALUE: 2.0 CAL FORMULA

-----  
ABBREVIATION.....2.0 C FORM  
ACTIVE.....Yes  
DESCRIPTION.....2.0 Cal Enteral Nutrition Formula  
DISPLAY\_NAME.....2.0 Cal Formula  
HELP\_TEXT.....No Help Text  
ID.....{4A595230-45D1-40F0-AEF3-11E75DF418CB}  
TERM.....2.0 CAL FORMULA

DATA VALUE: 20 MINUTES SINCE LAST IV NARCOTIC/SEDATIVE

-----  
ABBREVIATION.....20 Min Nar  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....20 Minutes Since Last IV  
                  Narcotic/Sedative  
HELP\_TEXT.....  
ID.....{FBDB50B7-19B5-472C-80E5-F4758E141BB1}  
TERM.....20 MINUTES SINCE LAST IV  
                  NARCOTIC/SEDATIVE

DATA VALUE: 20/100

-----  
ABBREVIATION.....20/100  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....20/100  
HELP\_TEXT.....  
ID.....{ECB85701-8B21-407A-9E0B-C015834E984B}  
TERM.....20/100

DATA VALUE: 20/15

-----  
ABBREVIATION.....20/15  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....20/15  
HELP\_TEXT.....  
ID.....{0DB4E45B-A96C-4D11-B59E-A89C9F531146}  
TERM.....20/15

DATA VALUE: 20/20

-----  
ABBREVIATION.....20/20  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....20/20  
HELP\_TEXT.....  
ID.....{FE086353-09B2-4314-84E9-FA909253AAED}  
TERM.....20/20

DATA VALUE: 20/200

-----  
ABBREVIATION.....20/200  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....20/200  
HELP\_TEXT.....  
ID.....{C3CAF0D6-22A6-4FD5-B8BB-55771F416D4E}  
TERM.....20/200

DATA VALUE: 20/25

-----  
ABBREVIATION.....20/25  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....20/25  
HELP\_TEXT.....  
ID.....{4A5CBDA4-00D2-479E-A6C8-CCD8B40AECFB}  
TERM.....20/25

DATA VALUE: 20/30

-----  
ABBREVIATION.....20/30  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....20/30  
HELP\_TEXT.....  
ID.....{667B5FCD-A85D-4054-BFA9-3DA9DDD21306}  
TERM.....20/30

DATA VALUE: 20/300

-----  
ABBREVIATION.....20/300  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....20/300  
HELP\_TEXT.....  
ID.....{E6976753-667B-418B-A6E3-0A7CD010CAC0}  
TERM.....20/300

DATA VALUE: 20/40

-----  
ABBREVIATION.....20/40  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....20/40  
HELP\_TEXT.....  
ID.....{048C6D5D-B7BF-4F2B-BFCC-BB4328A333CF}  
TERM.....20/40

DATA VALUE: 20/400

-----  
ABBREVIATION.....20/400  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....20/400  
HELP\_TEXT.....  
ID.....{AE8E4C3E-D031-4047-B44E-06F4EB29C73B}  
TERM.....20/400

DATA VALUE: 20/50

-----  
ABBREVIATION.....20/50  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....20/50  
HELP\_TEXT.....  
ID.....{C679C160-D93A-431C-8F28-800D64A019D5}  
TERM.....20/50

DATA VALUE: 20/60

-----  
ABBREVIATION.....20/60  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....20/60  
HELP\_TEXT.....  
ID.....{AEAA52C2-1450-45B4-A512-73E339718B02}  
TERM.....20/60

DATA VALUE: 20/70

-----  
ABBREVIATION.....20/70  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....20/70  
HELP\_TEXT.....  
ID.....{B0279832-605C-4441-B568-C131636EB634}  
TERM.....20/70

DATA VALUE: 20/80

-----  
ABBREVIATION.....20/80  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....20/80  
HELP\_TEXT.....  
ID.....{60F51F22-C739-4BD1-A77E-D224517F861F}  
TERM.....20/80

DATA VALUE: 20/800

-----  
ABBREVIATION.....20/800  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....20/800  
HELP\_TEXT.....  
ID.....{1303A571-B57C-4757-AC5B-286E6B09343A}  
TERM.....20/800

DATA VALUE: 2000 CALORIE CONTROLLED

-----  
ABBREVIATION.....2000 CAL C  
ACTIVE.....Yes  
DESCRIPTION.....2000 Calorie Controlled Diet  
DISPLAY\_NAME.....2000 Calorie Controlled Diet  
HELP\_TEXT.....No Help Text  
ID.....{DB2FF9DC-298E-4881-94D4-47A7AF0E2BD7}  
TERM.....2000 CALORIE CONTROLLED

DATA VALUE: 2000 DIABETIC

-----  
ABBREVIATION.....2000 DIAB  
ACTIVE.....Yes  
DESCRIPTION.....2000 Diabetic Diet  
DISPLAY\_NAME.....2000 Diabetic Diet  
HELP\_TEXT.....No Help Text  
ID.....{DF3056C9-3AEA-4F4E-B758-E35590080AFB}  
TERM.....2000 DIABETIC

DATA VALUE: 2200 CALORIE CONTROLLED

-----  
ABBREVIATION.....2200 CAL C  
ACTIVE.....Yes  
DESCRIPTION.....2200 Calorie Controlled Diet  
DISPLAY\_NAME.....2200 Calorie Controlled Diet  
HELP\_TEXT.....No Help Text  
ID.....{77F7C138-B2CB-4D76-8C0B-D68F317D4DE0}  
TERM.....2200 CALORIE CONTROLLED

DATA VALUE: 2200 DIABETIC

-----  
ABBREVIATION.....2200 DIAB  
ACTIVE.....Yes  
DESCRIPTION.....2200 Diabetic Diet  
DISPLAY\_NAME.....2200 Diabetic Diet  
HELP\_TEXT.....No Help Text  
ID.....{1FB3EFB6-B92F-44FB-831E-4641C91C7A8D}  
TERM.....2200 DIABETIC

DATA VALUE: 2400 CALORIE CONTROLLED

-----  
ABBREVIATION.....2400 CAL C  
ACTIVE.....Yes  
DESCRIPTION.....2400 Calorie Controlled Diet  
DISPLAY\_NAME.....2400 Calorie Controlled Diet  
HELP\_TEXT.....No Help Text  
ID.....{F3EA28B5-FC96-4F34-B2BE-A96A196B38A2}  
TERM.....2400 CALORIE CONTROLLED

DATA VALUE: 2400 DIABETIC

-----  
ABBREVIATION.....2400 DIAB  
ACTIVE.....Yes  
DESCRIPTION.....2400 Diabetic Diet  
DISPLAY\_NAME.....2400 Diabetic Diet  
HELP\_TEXT.....No Help Text  
ID.....{2018FB65-3D4E-47E0-90E9-205A45EDF130}  
TERM.....2400 DIABETIC

DATA VALUE: 25 PERCENT ALBUMIN

-----  
ABBREVIATION.....25 PERC AL  
ACTIVE.....Yes  
DESCRIPTION.....25 Percent Albumin  
DISPLAY\_NAME.....25 Percent Albumin  
HELP\_TEXT.....No Help Text  
ID.....{2FCE19CA-B5AB-4F15-8FEA-CF156D1E4CBC}  
TERM.....25 PERCENT ALBUMIN

DATA VALUE: 2:4

-----  
ABBREVIATION.....2:4  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....2:4  
HELP\_TEXT.....No Help Text  
ID.....{7D423C9A-42B3-4FDD-999B-889193492F5C}  
TERM.....2:4

DATA VALUE: 2A

-----  
ABBREVIATION.....2  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....2A  
HELP\_TEXT.....No Help Text  
ID.....{C75BDC9F-C756-9051-CAA3-A27DDD94B35E}  
TERM.....2A

DATA VALUE: 2A + 1V

-----  
ABBREVIATION.....2+1  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....2A + 1V  
HELP\_TEXT.....No Help Text  
ID.....{1849AA0D-1572-55E0-A780-ACE602E2C72D}  
TERM.....2A + 1V

DATA VALUE: 2A + 2V

-----  
ABBREVIATION.....2+2  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....2A + 2V  
HELP\_TEXT.....No Help Text  
ID.....{17203251-68E8-A00A-CF54-E2CC55B8EEDD}  
TERM.....2A + 2V

DATA VALUE: 2V

-----  
ABBREVIATION.....2  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....2V  
HELP\_TEXT.....No Help Text  
ID.....{D24B5FA4-1DA6-36DE-49DB-627A9B82E968}  
TERM.....2V



DATA VALUE: 3 - HR 130 OR MORE

-----  
ABBREVIATION.....HR 130 +  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....3 - HR 130 or More  
HELP\_TEXT.....  
ID.....{CBDB9BE1-832D-4F4A-9CF3-8768F054DDD2}  
TERM.....3 - HR 130 OR MORE

DATA VALUE: 3 - ABNORMAL FLEXION FROM PAIN

-----  
ABBREVIATION.....3-affp  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....3 - Abnormal Flexion From Pain  
HELP\_TEXT.....No Help Text  
ID.....{9AED8954-5F5D-D2C2-3CAF-4F8EF92009C2}  
TERM.....3 - ABNORMAL FLEXION FROM PAIN

DATA VALUE: 3 - BEADS OF SWEAT OBVIOUS ON FOREHEAD

-----  
ABBREVIATION.....CIWA PS 3  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....3 - Beads Of Sweat Obvious On Forehead  
HELP\_TEXT.....No Help Text  
ID.....{663591B6-3FB4-4060-8A49-18747F857B16}  
TERM.....3 - BEADS OF SWEAT OBVIOUS ON FOREHEAD

DATA VALUE: 3 - DISORIENTED FOR DATE > 2 CALENDAR DAYS

-----  
ABBREVIATION.....CIWA OR 3  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....3 - Disoriented For Date By More Than  
2 Calendar Days  
HELP\_TEXT.....No Help Text  
ID.....{68C351F9-C8DE-49C5-B3BD-D69E97F3E516}  
TERM.....3 - DISORIENTED FOR DATE > 2 CALENDAR  
DAYS

DATA VALUE: 3 - EXTENSIVE ASSISTANCE

-----  
ABBREVIATION.....3 Extensiv  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....3 - Extensive Assistance  
HELP\_TEXT.....  
ID.....{4B1FCD6E-7211-4151-ACD8-DF2F614A4337}  
TERM.....3 - EXTENSIVE ASSISTANCE



DATA VALUE: 3 - MODERATE ITCHING PINS/NEEDLES BURNING

-----  
ABBREVIATION.....CIWA TC 3  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....3 - Moderate Itching, Pins/Needles,  
                  Burning, Numbness  
HELP\_TEXT.....No Help Text  
ID.....{B7E2BFDD-7D85-4BCB-BE1D-0A220C1CEDAE}  
TERM.....3 - MODERATE ITCHING PINS/NEEDLES  
                  BURNING

DATA VALUE: 3 - MODERATE SENSITIVITY

-----  
ABBREVIATION.....CIWA VS 3  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....3 - Moderate Sensitivity  
HELP\_TEXT.....No Help Text  
ID.....{D673A630-1F1E-4734-AADE-17586C48E8B9}  
TERM.....3 - MODERATE SENSITIVITY

DATA VALUE: 3 - MODERATE WITH PATIENTS ARMS EXTENDED

-----  
ABBREVIATION.....CIWA TR 3  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....3 - Moderate With Patient's Arms  
                  Extended  
HELP\_TEXT.....No Help Text  
ID.....{36FF0AB3-F5ED-4CBA-A6F3-110DFCB65244}  
TERM.....3 - MODERATE WITH PATIENTS ARMS  
                  EXTENDED

DATA VALUE: 3 - PHYSICAL HELP IN PART OF BATHING ACTIVITY

-----  
ABBREVIATION.....3 Phys Bat  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....3 - Physical Help in Part of Bathing  
                  Activity  
HELP\_TEXT.....  
ID.....{1F24B08B-A5C1-4430-AC56-B40C347E1DDE}  
TERM.....3 - PHYSICAL HELP IN PART OF BATHING  
                  ACTIVITY

DATA VALUE: 3 - RR 30 OR MORE

-----  
ABBREVIATION.....RR 30 or +  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....3 - RR 30 or More  
HELP\_TEXT.....  
ID.....{6BFA9D84-AAD7-4DC0-BA75-530CE5C12E77}  
TERM.....3 - RR 30 OR MORE

DATA VALUE: 3 - SBP 70 OR LESS

-----  
ABBREVIATION.....SBP 70 or<  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....3 - SBP 70 or less  
HELP\_TEXT.....  
ID.....{84DE09CE-FCB0-4E07-A72F-C1064CECF094}  
TERM.....3 - SBP 70 OR LESS

DATA VALUE: 3 - STRONG

-----  
ABBREVIATION.....3-s  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....3 - Strong  
HELP\_TEXT.....No Help Text  
ID.....{C55EAB72-7816-5440-84E9-6E7E190ED197}  
TERM.....3 - STRONG

DATA VALUE: 3 - TO VERBAL COMMAND

-----  
ABBREVIATION.....3-tvc  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....3 - To Verbal Command  
HELP\_TEXT.....No Help Text  
ID.....{EBF4B748-2861-C81F-C0D0-B9FE99FC7A8F}  
TERM.....3 - TO VERBAL COMMAND

DATA VALUE: 3 - TWO + PERSONS PHYSICAL ASSIST

-----  
ABBREVIATION.....3 + asst  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....3 - Two + Persons Physical Assist  
HELP\_TEXT.....  
ID.....{AE401096-DF23-4A6E-8842-E1F58829902A}  
TERM.....3 - TWO + PERSONS PHYSICAL ASSIST

DATA VALUE: 3 - UNRESPONSIVE

-----  
ABBREVIATION.....Unresponsi  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....3 - Unresponsive  
HELP\_TEXT.....  
ID.....{F40F8A1E-ED15-424A-9481-3527B0E69DF1}  
TERM.....3 - UNRESPONSIVE

DATA VALUE: 3 WAY FOLEY

-----  
ABBREVIATION.....3WAY FOLEY  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....3 Way Foley  
HELP\_TEXT.....No Help Text  
ID.....{32408760-F517-4530-AEC0-54FE0A6A062B}  
TERM.....3 WAY FOLEY

DATA VALUE: 3+=BOUNDING, EASILY PALPABLE, CANNOT BE OBLITERATED

-----  
ABBREVIATION.....3epcbo  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....3+=Bounding, Easily Palpable, Cannot  
Be Obliterated  
HELP\_TEXT.....No Help Text  
ID.....{7A7A07ED-26B6-A3BF-C564-94BD4F6D4E1C}  
TERM.....3+=BOUNDING, EASILY PALPABLE, CANNOT  
BE OBLITERATED

DATA VALUE: 34 ml

-----  
ABBREVIATION.....34 ml  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....34 ml  
HELP\_TEXT.....No Help Text  
ID.....{25952BB5-50B1-4CC0-9483-8AF891767045}  
TERM.....34 ml

DATA VALUE: 3:4

-----  
ABBREVIATION.....3:4  
ACTIVE.....Yes  
DESCRIPTION.....3:4  
DISPLAY\_NAME.....3:4  
HELP\_TEXT.....No Help Text  
ID.....{E4EBD0E2-FE8C-40FA-9DCA-475DCE43B1A2}  
TERM.....3:4

DATA VALUE: 4 - BEADS OF SWEAT OBVIOUS ON FOREHEAD

-----  
ABBREVIATION.....CIWA PS 4  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....4 - Beads Of Sweat Obvious On Forehead  
HELP\_TEXT.....No Help Text  
ID.....{8388B1A0-46EA-4003-A32B-A2C708A6F056}  
TERM.....4 - BEADS OF SWEAT OBVIOUS ON FOREHEAD

DATA VALUE: 4 - DISORIENTED AND CONVERSES

-----  
ABBREVIATION.....4-dac  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....4 - Disoriented And Converses  
HELP\_TEXT.....No Help Text  
ID.....{D43F646D-A130-5D06-6DC6-DDB971027231}  
TERM.....4 - DISORIENTED AND CONVERSES

DATA VALUE: 4 - DISORIENTED FOR PLACE OR PERSON

-----  
ABBREVIATION.....CIWA OR 4  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....4 - Disoriented For Place/Or Person  
HELP\_TEXT.....No Help Text  
ID.....{03BBFFA2-233A-41A1-8212-F988889C5E72}  
TERM.....4 - DISORIENTED FOR PLACE OR PERSON

DATA VALUE: 4 - INTERMITTENT NAUSEA WITH DRY HEAVES

-----  
ABBREVIATION.....CIWA NV 4  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....4 - Intermittent Nausea With Dry Heaves  
HELP\_TEXT.....No Help Text  
ID.....{7FBE10F0-472F-4F21-BE55-0892259141B5}  
TERM.....4 - INTERMITTENT NAUSEA WITH DRY HEAVES

DATA VALUE: 4 - MODERATE WITH PATIENTS ARMS EXTENDED

-----  
ABBREVIATION.....CIWA TR 4  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....4 - Moderate With Patient's Arms  
Extended  
HELP\_TEXT.....No Help Text  
ID.....{79CE1486-151C-4492-AFD2-7845DD69DA9E}  
TERM.....4 - MODERATE WITH PATIENTS ARMS  
EXTENDED

DATA VALUE: 4 - MODERATELY ANXIOUS

-----  
ABBREVIATION.....CIWA AX 4  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....4 - Moderately Anxious, Or Guarded,  
Anxiety Is Inferred  
HELP\_TEXT.....No Help Text  
ID.....{9A01FEFD-844C-4FAB-BEC3-96333495FA62}  
TERM.....4 - MODERATELY ANXIOUS

DATA VALUE: 4 - MODERATELY FIDGETY AND RESTLESS

-----  
ABBREVIATION.....CIWA AG 4  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....4 - Moderately Fidgety And Restless  
HELP\_TEXT.....No Help Text  
ID.....{B1515E1F-4305-40D5-84B4-5E37AEF0AD3D}  
TERM.....4 - MODERATELY FIDGETY AND RESTLESS

DATA VALUE: 4 - MODERATELY SEVERE

-----  
ABBREVIATION.....CIWA HA 4  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....4 - Moderately Severe  
HELP\_TEXT.....No Help Text  
ID.....{83B3DE0F-3BA4-4B0A-A222-F372BD403094}  
TERM.....4 - MODERATELY SEVERE

DATA VALUE: 4 - MODERATELY SEVERE HALLUCINATIONS

-----  
ABBREVIATION.....CIWA 4  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....4 - Moderately Severe Hallucinations  
HELP\_TEXT.....No Help Text  
ID.....{13FCDA97-A967-4993-AED5-13EBA61F3102}  
TERM.....4 - MODERATELY SEVERE HALLUCINATIONS

DATA VALUE: 4 - SPONTANEOUSLY

-----  
ABBREVIATION.....4-s  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....4 - Spontaneously  
HELP\_TEXT.....No Help Text  
ID.....{4612FC84-FCA8-8C85-50DC-A9A42C195F51}  
TERM.....4 - SPONTANEOUSLY

DATA VALUE: 4 - TOTAL DEPENDENCE

-----  
ABBREVIATION.....4 Tot Dep  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....4 - Total Dependence  
HELP\_TEXT.....  
ID.....{062C8F50-294A-48CF-8E81-0C6D7787A4B1}  
TERM.....4 - TOTAL DEPENDENCE

DATA VALUE: 4 - TOTAL DEPENDENCE - FULL STAFF PERFORMANCE

-----  
ABBREVIATION.....4 Full Stf  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....4 - Total Dependence - Full Staff  
Performance  
HELP\_TEXT.....  
ID.....{E8B7D02B-E80A-4851-A0E2-EF27CF75A5BB}  
TERM.....4 - TOTAL DEPENDENCE - FULL STAFF  
PERFORMANCE

DATA VALUE: 4 - WITHDRAWS FROM PAIN

-----  
ABBREVIATION.....4-wfp  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....4 - Withdraws From Pain  
HELP\_TEXT.....No Help Text  
ID.....{D808ACF1-30CB-29AB-4957-E91D63C4DF18}  
TERM.....4 - WITHDRAWS FROM PAIN

DATA VALUE: 4 HOUR

-----  
ABBREVIATION.....4 HOUR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....4 Hour  
HELP\_TEXT.....No Help Text  
ID.....{91B2A49D-0308-47C9-B220-63A0E318A201}  
TERM.....4 HOUR

DATA VALUE: 40 ml

-----  
ABBREVIATION.....40 ml  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....40 ml  
HELP\_TEXT.....No Help Text  
ID.....{37A2998F-AB1C-4493-9C7A-FFCB05ED2302}  
TERM.....40 ml

DATA VALUE: 4:4

-----  
ABBREVIATION.....4:4  
ACTIVE.....Yes  
DESCRIPTION.....4:4  
DISPLAY\_NAME.....4:4  
HELP\_TEXT.....No Help Text  
ID.....{E9FCD5A0-995A-401D-8E85-6F0C753609B6}  
TERM.....4:4



DATA VALUE: 5 - CONSTANT NAUSEA W FREQUENT HEAVES AND VOMITING

-----  
ABBREVIATION.....CIWA NV 5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....5 - Constant Nausea W Frequent Heaves  
and Vomiting  
HELP\_TEXT.....No Help Text  
ID.....{3A90064F-0B58-48AB-ABEE-188F215F4E33}  
TERM.....5 - CONSTANT NAUSEA W FREQUENT HEAVES  
AND VOMITING

DATA VALUE: 5 - DRENCHING SWEAT

-----  
ABBREVIATION.....CIWA PS 5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....5 - Drenching Sweat  
HELP\_TEXT.....No Help Text  
ID.....{4CAA2CFF-6DA5-4238-A25B-4D3A960808A9}  
TERM.....5 - DRENCHING SWEAT

DATA VALUE: 5 - LOCALIZES PAIN

-----  
ABBREVIATION.....5-lp  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....5 - Localizes Pain  
HELP\_TEXT.....No Help Text  
ID.....{CF13650F-899B-6AE1-5913-DC1F4776E5E4}  
TERM.....5 - LOCALIZES PAIN

DATA VALUE: 5 - ORIENTED AND CONVERSES

-----  
ABBREVIATION.....5-oac  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....5 - Oriented And Converses  
HELP\_TEXT.....No Help Text  
ID.....{53AD29B1-BA51-4BDE-4F4C-4268619C93A4}  
TERM.....5 - ORIENTED AND CONVERSES

DATA VALUE: 5 - SEVERE

-----  
ABBREVIATION.....CIWA HA 5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....5 - Severe  
HELP\_TEXT.....No Help Text  
ID.....{7BD7B91F-DF6F-47B9-BAA6-3EF249340578}  
TERM.....5 - SEVERE

DATA VALUE: 5 - SEVERE EVEN WITH ARMS NOT EXTENDED

-----  
ABBREVIATION.....CIWA TR 5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....5 - Severe Even With Arms Not Extended  
HELP\_TEXT.....No Help Text  
ID.....{611A344A-C540-4861-81D7-681DAF1AA473}  
TERM.....5 - SEVERE EVEN WITH ARMS NOT EXTENDED

DATA VALUE: 5 - SEVERE HALLUCINATIONS

-----  
ABBREVIATION.....CIWA 5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....5 - Severe Hallucinations  
HELP\_TEXT.....No Help Text  
ID.....{14DB90F7-0182-44AC-8B5B-CE86A0E5A9E6}  
TERM.....5 - SEVERE HALLUCINATIONS

DATA VALUE: 5 PERCENT ALBUMIN

-----  
ABBREVIATION.....5 PERC ALB  
ACTIVE.....Yes  
DESCRIPTION.....5 Percent Albumin  
DISPLAY\_NAME.....5 Percent Albumin  
HELP\_TEXT.....No Help Text  
ID.....{38CE47CB-736B-475E-8E89-C4C13F2DF184}  
TERM.....5 PERCENT ALBUMIN

DATA VALUE: 6 - CONSTANT NAUSEA W FREQUENT HEAVES AND VOMITING

-----  
ABBREVIATION.....CIWA NV 6  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....6 - Constant Nausea W Frequent Heaves  
and Vomiting  
HELP\_TEXT.....No Help Text  
ID.....{446FE39A-EB53-4D6B-B5B8-1370A5D55EEB}  
TERM.....6 - CONSTANT NAUSEA W FREQUENT HEAVES  
AND VOMITING

DATA VALUE: 6 - DRENCHING SWEAT

-----  
ABBREVIATION.....CIWA PS 6  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....6 - Drenching Sweat  
HELP\_TEXT.....No Help Text  
ID.....{7294B311-9032-4FE1-9779-363D19FD1FD3}  
TERM.....6 - DRENCHING SWEAT

DATA VALUE: 6 - EXTREMELY SEVERE HALLUCINATIONS

-----  
ABBREVIATION.....CIWA 6  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....6 - Extremely Severe Hallucinations  
HELP\_TEXT.....No Help Text  
ID.....{9C63ED29-2B47-4E59-A026-1FDB0666AFCB}  
TERM.....6 - EXTREMELY SEVERE HALLUCINATIONS

DATA VALUE: 6 - OBEYS COMMANDS

-----  
ABBREVIATION.....6-oc  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....6 - Obeys Commands  
HELP\_TEXT.....No Help Text  
ID.....{0D4956F4-7E79-8D17-EDF3-274B96BBC89A}  
TERM.....6 - OBEYS COMMANDS

DATA VALUE: 6 - SEVERE EVEN WITH ARMS NOT EXTENDED

-----  
ABBREVIATION.....CIWA TR 6  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....6 - Severe Even With Arms Not Extended  
HELP\_TEXT.....No Help Text  
ID.....{9BCF893C-76D0-4A3E-80D0-281A917F1A08}  
TERM.....6 - SEVERE EVEN WITH ARMS NOT EXTENDED

DATA VALUE: 6 - VERY SEVERE

-----  
ABBREVIATION.....CIWA HA 6  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....6 - Very Severe  
HELP\_TEXT.....No Help Text  
ID.....{241C1EF4-2E2A-4D65-9991-607BB6B29E01}  
TERM.....6 - VERY SEVERE

DATA VALUE: 6 HOUR

-----  
ABBREVIATION.....6 HOUR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....6 Hour  
HELP\_TEXT.....No Help Text  
ID.....{0D1F000A-01F0-422C-9CB8-2A04CD479444}  
TERM.....6 HOUR

DATA VALUE: 6 PERCENT HETASTARCH

-----  
ABBREVIATION.....6 PERC HET  
ACTIVE.....Yes  
DESCRIPTION.....6 Percent Hetastarch  
DISPLAY\_NAME.....6 Percent Hetastarch  
HELP\_TEXT.....No Help Text  
ID.....{D94072C9-87E8-46CA-929F-4904F3BD964C}  
TERM.....6 PERCENT HETASTARCH

DATA VALUE: 60 MINUTES SINCE ADMINISTRATION OF REVERSAL AGENT

-----  
ABBREVIATION.....60 Min Rev  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....60 Minutes Since Administration of  
Reversal Agent  
HELP\_TEXT.....  
ID.....{87627CD7-2DA9-4BDA-BB74-68C4F03AA3D9}  
TERM.....60 MINUTES SINCE ADMINISTRATION OF  
REVERSAL AGENT

DATA VALUE: 7 - ACUTE PANIC

-----  
ABBREVIATION.....CIWA AX 7  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....7 - Equivalent To Acute Panic States  
HELP\_TEXT.....No Help Text  
ID.....{3F5E6786-FD67-405F-A535-EC1609910F14}  
TERM.....7 - ACUTE PANIC

DATA VALUE: 7 - CONSTANT NAUSEA W FREQUENT HEAVES AND VOMITING

-----  
ABBREVIATION.....CIWA NV 7  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....7 - Constant Nausea W Frequent Heaves  
and Vomiting  
HELP\_TEXT.....No Help Text  
ID.....{7F1076C4-C655-4BF2-8BB1-036F5169B619}  
TERM.....7 - CONSTANT NAUSEA W FREQUENT HEAVES  
AND VOMITING

DATA VALUE: 7 - CONTINUOUS HALLUCINATIONS

-----  
ABBREVIATION.....CIWA 7  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....7 - Continuous Hallucinations  
HELP\_TEXT.....No Help Text  
ID.....{19ECC57C-26A6-4CAE-97C4-2A860F337BE2}  
TERM.....7 - CONTINUOUS HALLUCINATIONS

DATA VALUE: 7 - DRENCHING SWEAT

-----  
ABBREVIATION.....CIWA PS 7  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....7 - Drenching Sweat  
HELP\_TEXT.....No Help Text  
ID.....{15D2E2CC-BAED-403B-A141-B848FFE4A9D0}  
TERM.....7 - DRENCHING SWEAT

DATA VALUE: 7 - EXTREMELY SEVERE

-----  
ABBREVIATION.....CIWA HA 7  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....7 - Extremely Severe  
HELP\_TEXT.....No Help Text  
ID.....{F8DCC910-D36D-480E-867E-B45221170C49}  
TERM.....7 - EXTREMELY SEVERE

DATA VALUE: 7 - PACING BACK AND FORTH MOST OF THE INTERVIEW

-----  
ABBREVIATION.....CIWA AG 7  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....7 - Pacing Back And Forth The  
Interview/Thrashing  
HELP\_TEXT.....No Help Text  
ID.....{9F03A167-CC40-4D65-9FD8-CFB3FEEAF833}  
TERM.....7 - PACING BACK AND FORTH MOST OF THE  
INTERVIEW

DATA VALUE: 7 - SEVERE EVEN WITH ARMS NOT EXTENDED

-----  
ABBREVIATION.....CIWA TR 7  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....7 - Severe Even With Arms Not Extended  
HELP\_TEXT.....No Help Text  
ID.....{FE0BE72E-913A-48CD-BABF-C5B742AFED93}  
TERM.....7 - SEVERE EVEN WITH ARMS NOT EXTENDED

DATA VALUE: 8 - ACTIVITY DID NOT OCCUR

-----  
ABBREVIATION.....8 d/n occ  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....8 - Activity Did Not Occur  
HELP\_TEXT.....  
ID.....{8B586CC9-2360-4060-B660-021C25E28F0F}  
TERM.....8 - ACTIVITY DID NOT OCCUR

DATA VALUE: 8 HOUR

-----  
ABBREVIATION.....8 HOUR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....8 Hour  
HELP\_TEXT.....No Help Text  
ID.....{05E46D8F-C165-4552-8063-B735DECD4AF3}  
TERM.....8 HOUR

DATA VALUE: AAI

-----  
ABBREVIATION.....AAI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....AAI  
HELP\_TEXT.....No Help Text  
ID.....{2DD18ACA-B98E-24D3-E760-6D082A76074B}  
TERM.....AAI

DATA VALUE: AAIR

-----  
ABBREVIATION.....AAIR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....AAIR  
HELP\_TEXT.....No Help Text  
ID.....{3B4B6F37-B631-BE32-F1C1-65B99D0C5018}  
TERM.....AAIR

DATA VALUE: AAO

-----  
ABBREVIATION.....AAO  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....AAO  
HELP\_TEXT.....No Help Text  
ID.....{E87E0670-2C0F-6E6A-DEC9-F52E5C530A2A}  
TERM.....AAO

DATA VALUE: ABDOMINAL

-----  
ABBREVIATION.....ABDOMINAL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Abdominal  
HELP\_TEXT.....No Help Text  
ID.....{29E279F4-5D2E-4644-A2FA-CF85BEDFD480}  
TERM.....ABDOMINAL

DATA VALUE: ABLE TO AMBULATE

-----  
ABBREVIATION.....AMBULAT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Able to Ambulate  
HELP\_TEXT.....No Help Text  
ID.....{4E1D40CA-F5C4-442E-8BBA-C6E4C4781F57}  
TERM.....ABLE TO AMBULATE

DATA VALUE: ABRASION

-----  
ABBREVIATION.....ABRASION  
ACTIVE.....Yes  
DESCRIPTION.....Abrasion  
DISPLAY\_NAME.....Abrasion  
HELP\_TEXT.....No Help Text  
ID.....{7E962B43-24B1-46D2-A79A-5481117271A}  
TERM.....ABRASION

DATA VALUE: ABSENT

-----  
ABBREVIATION.....a  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Absent  
HELP\_TEXT.....No Help Text  
ID.....{6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0}  
TERM.....ABSENT

DATA VALUE: ABSENT TO MINIMAL BLEEDING

-----  
ABBREVIATION.....ABS BLEE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Absent to Minimal Bleeding  
HELP\_TEXT.....No Help Text  
ID.....{123460FB-1A8D-45AB-B0A2-489EAB482F67}  
TERM.....ABSENT TO MINIMAL BLEEDING

DATA VALUE: ACCELERATING

-----  
ABBREVIATION.....ACCELERATI  
ACTIVE.....Yes  
DESCRIPTION.....Accelerating  
DISPLAY\_NAME.....Accelerating  
HELP\_TEXT.....Accelerating  
ID.....{514FE300-87EE-4C2C-AC70-5E57DBECAFF7}  
TERM.....ACCELERATING

DATA VALUE: ACHING

-----  
ABBREVIATION.....ACHING  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Aching  
HELP\_TEXT.....No Help Text  
ID.....{736BC1A4-2DA7-4A56-ACC0-28086EA103E0}  
TERM.....ACHING

DATA VALUE: ACTIVE RANGE OF MOTION

-----  
ABBREVIATION.....ACTIVE ROM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Active Range Of Motion  
HELP\_TEXT.....No Help Text  
ID.....{F55E68F3-6788-47AE-B962-119413B4A089}  
TERM.....ACTIVE RANGE OF MOTION

DATA VALUE: ACTIVITIES OF DAILY LIVING

-----  
ABBREVIATION.....ADL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Activities of Daily Living  
HELP\_TEXT.....No Help Text  
ID.....{12F2CA66-3DCE-443D-883C-119CE4838B4A}  
TERM.....ACTIVITIES OF DAILY LIVING

DATA VALUE: ADMISSION

-----  
ABBREVIATION.....ADMISSION  
ACTIVE.....Yes  
DESCRIPTION.....Admission  
DISPLAY\_NAME.....Admission  
HELP\_TEXT.....No Help Text  
ID.....{483C2830-813A-46C3-A012-CCD70044B2A2}  
TERM.....ADMISSION

DATA VALUE: ADMISSION DIAGNOSIS

-----  
ABBREVIATION.....ADM DIAG  
ACTIVE.....Yes  
DESCRIPTION.....Admission Diagnosis  
DISPLAY\_NAME.....Admission Diagnosis  
HELP\_TEXT.....Admission Diagnosis  
ID.....{2EEC3F2F-9802-46C5-AB47-008AFA093A7F}  
TERM.....ADMISSION DIAGNOSIS



DATA VALUE: ADULT

-----  
ABBREVIATION.....ADULT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Adult  
HELP\_TEXT.....No Help Text  
ID.....{37B557BE-B932-4757-8E08-5EA315E8AD1F}  
TERM.....ADULT

DATA VALUE: ADULT LARGE

-----  
ABBREVIATION.....AD LARGE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Adult Large  
HELP\_TEXT.....No Help Text  
ID.....{38E64D72-2806-47FA-B82F-51844574F850}  
TERM.....ADULT LARGE

DATA VALUE: ADULT SMALL

-----  
ABBREVIATION.....AD SMALL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Adult Small  
HELP\_TEXT.....No Help Text  
ID.....{0161A4E5-713A-406B-B3FD-9ED13BEB5366}  
TERM.....ADULT SMALL

DATA VALUE: ADVANCED CARE PLANNING

-----  
ABBREVIATION.....ADV CAR PL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Advanced Care Planning  
HELP\_TEXT.....No Help Text  
ID.....{F0ED0F8A-1AD6-4F3E-A4C8-4FB915474FE7}  
TERM.....ADVANCED CARE PLANNING

DATA VALUE: AGENT CASHIER

-----  
ABBREVIATION.....Cashier  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Agent Cashier  
HELP\_TEXT.....  
ID.....{3849A887-F353-4215-88FD-DDD0188BE1BB}  
TERM.....AGENT CASHIER

DATA VALUE: AGITATED

-----  
ABBREVIATION.....a  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Agitated  
HELP\_TEXT.....No Help Text  
ID.....{BAD2CCCE-9083-6F61-1196-C0C9A8DDBFFC}  
TERM.....AGITATED

DATA VALUE: AGONAL

-----  
ABBREVIATION.....a  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Agonal  
HELP\_TEXT.....No Help Text  
ID.....{1621B411-33E3-3A6A-78A6-40B7BDC3E8CB}  
TERM.....AGONAL

DATA VALUE: AGONIZING

-----  
ABBREVIATION.....AGONIZING  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Agonizing  
HELP\_TEXT.....No Help Text  
ID.....{F43DF9E7-555B-49FF-987C-35F083D547AA}  
TERM.....AGONIZING

DATA VALUE: AICD

-----  
ABBREVIATION.....AICD  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....AICD  
HELP\_TEXT.....  
ID.....{EDB0E730-A1DC-4639-888B-30E956CCA958}  
TERM.....AICD

DATA VALUE: AIR FLOW DEVICE

-----  
ABBREVIATION.....AIR FLOW D  
ACTIVE.....Yes  
DESCRIPTION.....Air Flow Device  
DISPLAY\_NAME.....Air Flow Device  
HELP\_TEXT.....Air Flow Device  
ID.....{60AF6215-CE2A-4EEE-9596-909E3E50BE71}  
TERM.....AIR FLOW DEVICE



DATA VALUE: ALCOHOL PORT PROTECTOR

-----  
ABBREVIATION.....AI Port Pr  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Alcohol Port Protector  
HELP\_TEXT.....  
ID.....{785F47CF-A279-46F9-9073-B6A6A70E8888}  
TERM.....ALCOHOL PORT PROTECTOR

DATA VALUE: ALDRETE SCORE 8 OR GREATER

-----  
ABBREVIATION.....ALD 8+  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Aldrete Score 8 or Greater  
HELP\_TEXT.....No Help Text  
ID.....{1FACE6D7-BC92-4579-8665-C5E1D3475461}  
TERM.....ALDRETE SCORE 8 OR GREATER

DATA VALUE: ALDRETE SCORE AT BASELINE

-----  
ABBREVIATION.....ALD BASE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Aldrete Score at Baseline  
HELP\_TEXT.....No Help Text  
ID.....{CC86D48F-9473-4AD8-AA5B-3BA0D1CCA8FD}  
TERM.....ALDRETE SCORE AT BASELINE

DATA VALUE: ALERT

-----  
ABBREVIATION.....a  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Alert  
HELP\_TEXT.....No Help Text  
ID.....{BC7E2630-179A-5C80-E6C7-598D42F0BFD6}  
TERM.....ALERT

DATA VALUE: ALERT AND ORIENTED

-----  
ABBREVIATION.....AL AN OR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Alert and Oriented  
HELP\_TEXT.....No Help Text  
ID.....{C031F776-36DF-4570-AFAF-A30B8BFA635B}  
TERM.....ALERT AND ORIENTED

DATA VALUE: ALL DOWN

-----  
ABBREVIATION.....ALL DOWN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....All Down  
HELP\_TEXT.....No Help Text  
ID.....{76C14E79-95E2-49B7-B847-3B114653E98C}  
TERM.....ALL DOWN

DATA VALUE: ALL UP

-----  
ABBREVIATION.....ALL UP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....All Up  
HELP\_TEXT.....No Help Text  
ID.....{48199083-ACC2-408C-9242-7A6313CD7A7F}  
TERM.....ALL UP

DATA VALUE: ALLERGY

-----  
ABBREVIATION.....a  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Allergy  
HELP\_TEXT.....No Help Text  
ID.....{EBF87F6F-1BDB-6147-C444-2FF787DEAC7E}  
TERM.....ALLERGY

DATA VALUE: ALLERGY WRISTBAND

-----  
ABBREVIATION.....Wristband  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Allergy Wristband  
HELP\_TEXT.....  
ID.....{BB7F6746-102D-4CBA-8E53-EA465176812B}  
TERM.....ALLERGY WRISTBAND

DATA VALUE: ALOPECIA PARTIAL

-----  
ABBREVIATION.....ALOPE PART  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Alopecia Partial  
HELP\_TEXT.....No Help Text  
ID.....{98EDBD8D-205A-4630-A09D-6B4222F5B2A2}  
TERM.....ALOPECIA PARTIAL

DATA VALUE: ALOPECIA TOTAL

-----  
ABBREVIATION.....ALOPE TOTA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Alopecia Total  
HELP\_TEXT.....No Help Text  
ID.....{ECDCE3C4-81F6-4FA1-BADA-97A20BD9A826}  
TERM.....ALOPECIA TOTAL

DATA VALUE: AMBER

-----  
ABBREVIATION.....AMBER  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Amber  
HELP\_TEXT.....No Help Text  
ID.....{90A13048-4C35-E01D-3482-64D782001D4A}  
TERM.....AMBER

DATA VALUE: AMBULATE AND ASSIST

-----  
ABBREVIATION.....AMB ASSIS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Ambulate and Assist  
HELP\_TEXT.....No Help Text  
ID.....{EB53ADD9-B5B0-4174-95C7-C4B33E80DBDC}  
TERM.....AMBULATE AND ASSIST

DATA VALUE: AMBULATED

-----  
ABBREVIATION.....AMBULATED  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Ambulated  
HELP\_TEXT.....No Help Text  
ID.....{017D99FE-91D5-49F1-AC15-50620879C466}  
TERM.....AMBULATED

DATA VALUE: AMBULATORY CARE UNIT

-----  
ABBREVIATION.....ACU  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Ambulatory Care Unit  
HELP\_TEXT.....No Help Text  
ID.....{7FB7FBC5-42A6-4BA4-9CA6-494D77C4C549}  
TERM.....AMBULATORY CARE UNIT

DATA VALUE: AMPUTEE

-----  
ABBREVIATION.....AMPUTEE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Amputee  
HELP\_TEXT.....No Help-Text  
ID.....{5804286F-2B23-91BD-FCB5-6AD18C918203}  
TERM.....AMPUTEE

DATA VALUE: ANALGESIA

-----  
ABBREVIATION.....ANALGESIA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Analgesia  
HELP\_TEXT.....No Help Text  
ID.....{5452ED11-25E1-437B-8048-3BB3E03BE3B9}  
TERM.....ANALGESIA

DATA VALUE: ANESTHESIOLOGIST

-----  
ABBREVIATION.....Anesthesio  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Anesthesiologist  
HELP\_TEXT.....  
ID.....{44BAD970-DD8F-4015-9C41-AC7868CEBE79}  
TERM.....ANESTHESIOLOGIST

DATA VALUE: ANGINAL

-----  
ABBREVIATION.....ANGINAL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....ANGINAL  
HELP\_TEXT.....No Help Text  
ID.....{ECF02415-4124-46D3-AA01-F21AC40FEFF7}  
TERM.....ANGINAL

DATA VALUE: ANNOYING

-----  
ABBREVIATION.....ANNOYING  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Annoying  
HELP\_TEXT.....No Help Text  
ID.....{53437530-57C9-4EA8-AC0D-36B5F2BD051D}  
TERM.....ANNOYING

DATA VALUE: ANTI-EMBOLUS STOCKINGS KNEE HIGH

-----  
ABBREVIATION.....ANT EMB KN  
ACTIVE.....Yes  
DESCRIPTION.....ANTI-EMBOLUS STOCKINGS KNEE HIGH  
DISPLAY\_NAME.....Anti-Embolus Stockings Knee High  
HELP\_TEXT.....ANTI-EMBOLUS STOCKINGS KNEE HIGH  
ID.....{4FBFF7A8-D719-98B9-4A46-D76D7BB5213E}  
TERM.....ANTI-EMBOLUS STOCKINGS KNEE HIGH

DATA VALUE: ANTI-EMBOLUS STOCKINGS THIGH HIGH

-----  
ABBREVIATION.....ANT EMB ST  
ACTIVE.....Yes  
DESCRIPTION.....ANTI-EMBOLUS STOCKINGS THIGH HIGH  
DISPLAY\_NAME.....Anti-Embolus Stockings Thigh High  
HELP\_TEXT.....ANTI-EMBOLUS STOCKINGS THIGH HIGH  
ID.....{9B32D3F3-47F1-886F-A59A-D8AB38A34754}  
TERM.....ANTI-EMBOLUS STOCKINGS THIGH HIGH

DATA VALUE: ANTI-MICROBIAL DRESSING

-----  
ABBREVIATION.....AB DRESSIN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Anti-Microbial Dressing  
HELP\_TEXT.....No Help Text  
ID.....{9734C5E9-D099-47EB-93B1-D82201F64331}  
TERM.....ANTI-MICROBIAL DRESSING

DATA VALUE: ANTIBIOTIC OINTMENT

-----  
ABBREVIATION.....AB OINTMEN  
ACTIVE.....Yes  
DESCRIPTION.....Antibiotic Ointment  
DISPLAY\_NAME.....Antibiotic Ointment  
HELP\_TEXT.....No Help Text  
ID.....{352D9690-77BB-44DC-919E-5357876DEAC4}  
TERM.....ANTIBIOTIC OINTMENT

DATA VALUE: ANTIMICROBIAL DISK

-----  
ABBREVIATION.....WOU DR ADI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Antimicrobial Disk  
HELP\_TEXT.....No Help-Text  
ID.....{1F1745EA-D763-4AFD-9CED-82E974E57692}  
TERM.....ANTIMICROBIAL DISK



DATA VALUE: ANURIC

-----  
ABBREVIATION.....ANURIC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Anuric  
HELP\_TEXT.....No Help Text  
ID.....{0CCA4B44-7A59-4042-BAB9-EB2395A63C51}  
TERM.....ANURIC

DATA VALUE: ANXIOUS

-----  
ABBREVIATION.....a  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Anxious  
HELP\_TEXT.....No Help Text  
ID.....{381BDB41-7E8A-6FDB-841B-2ECC1889B9EF}  
TERM.....ANXIOUS

DATA VALUE: AOO

-----  
ABBREVIATION.....AOO  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....AOO  
HELP\_TEXT.....No Help Text  
ID.....{86F25E48-AADB-5B30-4B89-E8BE63DAE7D9}  
TERM.....AOO

DATA VALUE: APNEA

-----  
ABBREVIATION.....Apnea  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Apnea  
HELP\_TEXT.....  
ID.....{BF6A89F6-146B-4BE9-AD36-FF2E2C79022E}  
TERM.....APNEA

DATA VALUE: APPEARS ASLEEP

-----  
ABBREVIATION.....App Asleep  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Appears Asleep  
HELP\_TEXT.....  
ID.....{8D1436B8-ADB0-4092-A0CE-383D74131AA8}  
TERM.....APPEARS ASLEEP

DATA VALUE: APPLIANCE CHANGED

-----  
ABBREVIATION.....app change  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Appliance Changed  
HELP\_TEXT.....  
ID.....{CB4BC887-F738-4EE3-9A5A-39C9CDF7488C}  
TERM.....APPLIANCE CHANGED

DATA VALUE: APPROPRIATE

-----  
ABBREVIATION.....a  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Appropriate  
HELP\_TEXT.....No Help Text  
ID.....{1CFB285F-34DF-73A2-B89C-F6EE1F40E7E9}  
TERM.....APPROPRIATE

DATA VALUE: ARRHYTHMIA

-----  
ABBREVIATION.....Arrhythmi  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Arrhythmia  
HELP\_TEXT.....  
ID.....{F12E4852-27E7-4C2C-989F-E973EBDF2C62}  
TERM.....ARRHYTHMIA

DATA VALUE: ARTERIAL LINE CARE

-----  
ABBREVIATION.....ART LINE C  
ACTIVE.....Yes  
DESCRIPTION.....Arterial Line Care  
DISPLAY\_NAME.....Arterial Line Care  
HELP\_TEXT.....No Help Text  
ID.....{787D232D-DFB8-44E2-8A0C-9D7D6B80764E}  
TERM.....ARTERIAL LINE CARE

DATA VALUE: ARTIFICIAL AIRWAY

-----  
ABBREVIATION.....ART AIRWAY  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Artificial Airway  
HELP\_TEXT.....No Help Text  
ID.....{8DEE2358-8892-4BFB-AFAE-6C134D1CCAD4}  
TERM.....ARTIFICIAL AIRWAY

DATA VALUE: ARTIFICIAL VALVES

-----  
ABBREVIATION.....Art Valves  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Artificial Valves  
HELP\_TEXT.....  
ID.....{088BB4A5-A93E-47E2-B21F-BC2957D52693}  
TERM.....ARTIFICIAL VALVES

DATA VALUE: ASA PS 1

-----  
ABBREVIATION.....ASA PS 1  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....ASA PS 1  
HELP\_TEXT.....  
ID.....{BAF1730E-C358-17F1-8DB8-4B6097AFF066}  
TERM.....ASA PS 1

DATA VALUE: ASA PS 2

-----  
ABBREVIATION.....ASA PS 2  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....ASA PS 2  
HELP\_TEXT.....  
ID.....{070C0F56-3D9C-5430-EF9E-BF3B96DFB595}  
TERM.....ASA PS 2

DATA VALUE: ASA PS 3

-----  
ABBREVIATION.....ASA PS 3  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....ASA PS 3  
HELP\_TEXT.....  
ID.....{B041CE25-82D0-EDAD-C673-A64BCD6647F4}  
TERM.....ASA PS 3

DATA VALUE: ASA PS 4

-----  
ABBREVIATION.....ASA PS 4  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....ASA PS 4  
HELP\_TEXT.....  
ID.....{040841C3-D933-E687-C7E0-29F439C08778}  
TERM.....ASA PS 4

DATA VALUE: ASA PS 5

-----  
ABBREVIATION.....ASA PS 5  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....ASA PS 5  
HELP\_TEXT.....  
ID.....{1F7B73A6-3F8E-E230-1C23-C7DF92032E58}  
TERM.....ASA PS 5

DATA VALUE: ASEPTIC TECHNIQUE FOLLOWED

-----  
ABBREVIATION.....Aseptic  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Aseptic Technique Followed  
HELP\_TEXT.....  
ID.....{2DE69DDB-F455-4570-A821-C81F39FD2574}  
TERM.....ASEPTIC TECHNIQUE FOLLOWED

DATA VALUE: ASHEN

-----  
ABBREVIATION.....ASHEN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Ashen  
HELP\_TEXT.....No Help-Text  
ID.....{7431E4EC-DCF1-4966-85CB-29BF65311855}  
TERM.....ASHEN

DATA VALUE: ASLEEP AT ROUNDS

-----  
ABBREVIATION.....SLEEP ROUN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Appeared Asleep At All Rounds  
HELP\_TEXT.....No Help Text  
ID.....{24E1D9F0-D8F1-4DC8-911F-6170F343C087}  
TERM.....ASLEEP AT ROUNDS

DATA VALUE: ASPIRATION

-----  
ABBREVIATION.....ASPIRATION  
ACTIVE.....Yes  
DESCRIPTION.....Aspiration  
DISPLAY\_NAME.....Aspiration  
HELP\_TEXT.....No Help Text  
ID.....{68B91085-8E3C-4B94-9F97-6305F5DFD900}  
TERM.....ASPIRATION

DATA VALUE: ASSAULTIVE

-----  
ABBREVIATION.....a  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Assaultive  
HELP\_TEXT.....No Help Text  
ID.....{99ECD8C7-191A-3B9C-2A7A-803A921065FA}  
TERM.....ASSAULTIVE

DATA VALUE: ASSIGNED ONE TO ONE SITTER

-----  
ABBREVIATION.....One to One  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Assigned One to One Sitter  
HELP\_TEXT.....  
ID.....{07D48368-101F-4974-9EC5-7CB182A39281}  
TERM.....ASSIGNED ONE TO ONE SITTER

DATA VALUE: ASSIST CONTROL

-----  
ABBREVIATION.....AC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Assist Control  
HELP\_TEXT.....No Help Text  
ID.....{862ED94C-9DC3-DFAF-F824-0F7DE348C2CB}  
TERM.....ASSIST CONTROL

DATA VALUE: ASSISTED

-----  
ABBREVIATION.....ASSISTED  
ACTIVE.....Yes  
DESCRIPTION.....Assisted  
DISPLAY\_NAME.....Assisted  
HELP\_TEXT.....No Help Text  
ID.....{65069804-39BF-41ED-9EA3-BF2C41389F10}  
TERM.....ASSISTED

DATA VALUE: ASYSTOLE

-----  
ABBREVIATION.....a  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Asystole  
HELP\_TEXT.....No Help Text  
ID.....{6E87ECBC-3160-2CEF-D2FF-EEBD1DAA96F2}  
TERM.....ASYSTOLE

DATA VALUE: ATRIAL

-----  
ABBREVIATION.....a  
ACTIVE.....Yes  
DESCRIPTION.....ATRIAL  
DISPLAY\_NAME.....Atrial  
HELP\_TEXT.....No Help Text  
ID.....{266D2FDA-5D6A-203D-6FF0-8282D49BFD67}  
TERM.....ATRIAL

DATA VALUE: ATRIAL FIB/FLUTTER

-----  
ABBREVIATION.....af  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Atrial Fib/Flutter  
HELP\_TEXT.....No Help Text  
ID.....{9C1BC2F6-3BF0-BEBB-BE87-50E9F7984544}  
TERM.....ATRIAL FIB/FLUTTER

DATA VALUE: ATRIAL FIBRILLATION

-----  
ABBREVIATION.....a fib  
ACTIVE.....Yes  
DESCRIPTION.....Atrial Fibrillation  
DISPLAY\_NAME.....Atrial Fibrillation  
HELP\_TEXT.....Atrial Fibrillation  
ID.....{7B5EBAF2-55B6-2905-3D89-B28CE6A509CB}  
TERM.....ATRIAL FIBRILLATION

DATA VALUE: ATRIAL FLUTTER

-----  
ABBREVIATION.....af  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Atrial Flutter  
HELP\_TEXT.....No Help Text  
ID.....{9CB8D413-6762-D78E-F40E-25DB66CD23B3}  
TERM.....ATRIAL FLUTTER

DATA VALUE: AUSCULTATION

-----  
ABBREVIATION.....AUSCULTATI  
ACTIVE.....Yes  
DESCRIPTION.....Auscultation  
DISPLAY\_NAME.....Auscultation  
HELP\_TEXT.....Auscultation  
ID.....{B9D756DF-74F5-4604-B781-D7BC512FC453}  
TERM.....AUSCULTATION

DATA VALUE: AVF

-----  
ABBREVIATION.....AVF  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....AVF  
HELP\_TEXT.....No Help Text  
ID.....{C8E97206-58AA-47C0-BD2E-7A30E56DB46A}  
TERM.....AVF

DATA VALUE: AVI

-----  
ABBREVIATION.....AVI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....AVI  
HELP\_TEXT.....No Help Text  
ID.....{8BB49456-FB80-4916-6EBE-06837BF12B0C}  
TERM.....AVI

DATA VALUE: AVL

-----  
ABBREVIATION.....AVL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....AVL  
HELP\_TEXT.....No Help Text  
ID.....{9124A193-5F50-448D-9E60-C5869EC4029C}  
TERM.....AVL

DATA VALUE: AVO

-----  
ABBREVIATION.....AVO  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....AVO  
HELP\_TEXT.....No Help Text  
ID.....{02EA752A-D156-96E6-3D28-3D389FFC7B58}  
TERM.....AVO

DATA VALUE: AVOIDED LIGHT AND NOISE OVERSTIMULATION

-----  
ABBREVIATION.....Avoid Lght  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Avoided Light and Noise Overstimulation  
HELP\_TEXT.....  
ID.....{135BE5DC-B93B-4A72-97AC-B423A6E81538}  
TERM.....AVOIDED LIGHT AND NOISE OVERSTIMULATION

DATA VALUE: AVR

-----  
ABBREVIATION.....AVR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....AVR  
HELP\_TEXT.....No Help Text  
ID.....{04AF2551-1185-4BA9-AEED-977D562D188A}  
TERM.....AVR

DATA VALUE: AWAKE

-----  
ABBREVIATION.....Awake  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Awake  
HELP\_TEXT.....  
ID.....{D71412D4-5782-46BA-AE6C-24EDC586AE55}  
TERM.....AWAKE

DATA VALUE: BACTERICIDAL

-----  
ABBREVIATION.....BACTERIOCI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bactericidal  
HELP\_TEXT.....No Help-Text  
ID.....{680E4A3C-A40A-461E-818D-003AE8A1EFC8}  
TERM.....BACTERICIDAL

DATA VALUE: BACTERIOSTATIC

-----  
ABBREVIATION.....BACTERIOST  
ACTIVE.....Yes  
DESCRIPTION.....Bacteriostatic  
DISPLAY\_NAME.....Bacteriostatic  
HELP\_TEXT.....No Help Text  
ID.....{AB74D109-8ED7-4A6C-B028-7A041781DF67}  
TERM.....BACTERIOSTATIC

DATA VALUE: BAG - CONTINUOUS

-----  
ABBREVIATION.....Bag - cont  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Bag - Continuous  
HELP\_TEXT.....  
ID.....{F072EF8A-B38B-4D75-82BF-D5EF851784B0}  
TERM.....BAG - CONTINUOUS



DATA VALUE: BAG - INTERMITTENT

-----  
ABBREVIATION.....Bag - inte  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Bag - Intermittent  
HELP\_TEXT.....  
ID.....{3761CDD3-AA2A-4F2B-9954-865E0A3D8C9A}  
TERM.....BAG - INTERMITTENT

DATA VALUE: BAG CHANGE

-----  
ABBREVIATION.....Bag Chg  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Bag Change  
HELP\_TEXT.....  
ID.....{BC58D518-A742-410B-BE33-732C04D303BD}  
TERM.....BAG CHANGE

DATA VALUE: BAG CHANGED

-----  
ABBREVIATION.....BAG CHANGE  
ACTIVE.....Yes  
DESCRIPTION.....Changed Bag  
DISPLAY\_NAME.....Bag Changed  
HELP\_TEXT.....No Help Text  
ID.....{818990DE-4F95-4F8E-92D8-B4D3B975CB36}  
TERM.....BAG CHANGED

DATA VALUE: BAG/MASK & O2 SET UP

-----  
ABBREVIATION.....Bag/Mask  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Bag/Mask & O2 Set Up  
HELP\_TEXT.....  
ID.....{248ECF6D-D89B-4318-BFB7-D28CF5C9C5B0}  
TERM.....BAG/MASK & O2 SET UP

DATA VALUE: BALSAM OF PERU

-----  
ABBREVIATION.....WOU DR PER  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Balsam Of Peru  
HELP\_TEXT.....No Help-Text  
ID.....{DD412BB1-B1C1-4A48-A38F-832629448E06}  
TERM.....BALSAM OF PERU

DATA VALUE: BARIATRIC FRAME

-----  
ABBREVIATION.....BARI FRAME  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bariatric Frame  
HELP\_TEXT.....No Help Text  
ID.....{CF95B426-F280-451B-86F2-CEF7A1A2B967}  
TERM.....BARIATRIC FRAME

DATA VALUE: BARIATRIC LOW AIRLOSS MATTRESS REPLACEMENT

-----  
ABBREVIATION.....BAR L AIR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bariatric Low Airloss Mattress  
Replacement  
HELP\_TEXT.....No Help Text  
ID.....{220C6E92-1D0A-4513-B3A9-8065A6EC7DB1}  
TERM.....BARIATRIC LOW AIRLOSS MATTRESS  
REPLACEMENT

DATA VALUE: BARIATRIC PATIENT

-----  
ABBREVIATION.....BARIAC PAT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bariatric Patient  
HELP\_TEXT.....No Help Text  
ID.....{6F614334-685E-4533-A475-BDF5440855EA}  
TERM.....BARIATRIC PATIENT

DATA VALUE: BASAL RATE ONLY

-----  
ABBREVIATION.....BASAL ONL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Basal Rate Only  
HELP\_TEXT.....No Help Text  
ID.....{4BFFF2BC-837C-4D62-B811-0E547247D4D0}  
TERM.....BASAL RATE ONLY

DATA VALUE: BATHROOM

-----  
ABBREVIATION.....b  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bathroom  
HELP\_TEXT.....No Help Text  
ID.....{D20D91A7-E9F5-D000-5E0B-772BA9F5F72E}  
TERM.....BATHROOM

DATA VALUE: BATHROOM PRIVILEGES

-----  
ABBREVIATION.....BATHR PRIV  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bathroom Privileges  
HELP\_TEXT.....No Help Text  
ID.....{29C1B6DD-3308-4DED-9970-1322A97C633B}  
TERM.....BATHROOM PRIVILEGES

DATA VALUE: BED

-----  
ABBREVIATION.....BED  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bed  
HELP\_TEXT.....No Help Text  
ID.....{7E447C9C-41A0-43C9-AF63-38D6EA68279C}  
TERM.....BED

DATA VALUE: BED EXIT ALARM

-----  
ABBREVIATION.....BED EXIT A  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bed Exit Alarm  
HELP\_TEXT.....No Help Text  
ID.....{AC53759E-FFAE-5E0F-2EBA-E56A2712357C}  
TERM.....BED EXIT ALARM

DATA VALUE: BEDREST

-----  
ABBREVIATION.....BEDREST  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bedrest  
HELP\_TEXT.....No Help Text  
ID.....{8CD28057-132E-4CEE-9EA1-F03C8AF50CE8}  
TERM.....BEDREST

DATA VALUE: BEDSIDE

-----  
ABBREVIATION.....Bedside  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Bedside  
HELP\_TEXT.....  
ID.....{481BBBB3-A1F0-4F49-96DA-06BF3884E9BD}  
TERM.....BEDSIDE

DATA VALUE: BEDSIDE COMMUNE

-----  
ABBREVIATION.....Bedside Cm  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Bedside Commode  
HELP\_TEXT.....  
ID.....{E2B0C66F-335A-418E-AEB8-C7B666AE2F8D}  
TERM.....BEDSIDE COMMUNE

DATA VALUE: BIER BLOCK

-----  
ABBREVIATION.....Bier Blk  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Bier Block  
HELP\_TEXT.....  
ID.....{241C28E7-6040-4ABC-B76C-4255452C66FD}  
TERM.....BIER BLOCK

DATA VALUE: BIGEMINY

-----  
ABBREVIATION.....Bigeminy  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Bigeminy  
HELP\_TEXT.....  
ID.....{729FCAD7-EDA1-4CFD-A6CB-E291DDCF0255}  
TERM.....BIGEMINY

DATA VALUE: BILATERAL

-----  
ABBREVIATION.....Bilateral  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bilateral  
HELP\_TEXT.....No Help Text  
ID.....{FCD5361B-C674-473C-8F29-7B08650BA338}  
TERM.....BILATERAL

DATA VALUE: BIPAP - RESPIRONICS

-----  
ABBREVIATION.....BIPAP RESP  
ACTIVE.....Yes  
DESCRIPTION.....BiPap machine by Respironics  
DISPLAY\_NAME.....Bipap - Respironics  
HELP\_TEXT.....No Help Text  
ID.....{AA94A723-4C68-489F-88AC-2C6F276EAB76}  
TERM.....BIPAP - RESPIRONICS

DATA VALUE: BIPAP - VISION

-----  
ABBREVIATION.....BIPAP VISI  
ACTIVE.....Yes  
DESCRIPTION.....BiPap ventilator made by Vision  
DISPLAY\_NAME.....Bipap - Vision  
HELP\_TEXT.....No Help Text  
ID.....{AA39C765-BA1F-4289-B459-200A57C07FFF}  
TERM.....BIPAP - VISION

DATA VALUE: BIPAP/CPAP

-----  
ABBREVIATION.....BIPAP/CPAP  
ACTIVE.....Yes  
DESCRIPTION.....BiPap / CPAP machine  
DISPLAY\_NAME.....BiPap / CPAP machine  
HELP\_TEXT.....No Help Text  
ID.....{8D767C67-3712-410F-A96B-69C731E03C5D}  
TERM.....BIPAP/CPAP

DATA VALUE: BIPHASIC POSITIVE AIRWAY PRESSURE

-----  
ABBREVIATION.....BiPAP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....BiPap  
HELP\_TEXT.....No Help Text  
ID.....{4802AB1E-A0AD-BB9F-E8BA-33863BCD3466}  
TERM.....BIPHASIC POSITIVE AIRWAY PRESSURE

DATA VALUE: BISMUTH GAUZE

-----  
ABBREVIATION.....WOU DR BIS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bismuth Gauze  
HELP\_TEXT.....No Help-Text  
ID.....{010DED59-C808-4D92-A7F8-87AEB86B4439}  
TERM.....BISMUTH GAUZE

DATA VALUE: BITE

-----  
ABBREVIATION.....Bite  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bite  
HELP\_TEXT.....No Help Text  
ID.....{FC33E0AE-3130-4881-B9EF-9B81FBEC6AC0}  
TERM.....BITE

DATA VALUE: BITE BLOCK

-----  
ABBREVIATION.....BITE BLOCK  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bite Block  
HELP\_TEXT.....No Help Text  
ID.....{E7CC51F3-08A8-5FEE-38B2-16B74E1F542E}  
TERM.....BITE BLOCK

DATA VALUE: BLACK

-----  
ABBREVIATION.....BLACK  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Black  
HELP\_TEXT.....No Help Text  
ID.....{7A55AFB9-6861-4918-0746-3D67D47FE5C3}  
TERM.....BLACK

DATA VALUE: BLACK/NECROTIC

-----  
ABBREVIATION.....BLACK/NECR  
ACTIVE.....Yes  
DESCRIPTION.....Black / Necrotic  
DISPLAY\_NAME.....Black / Necrotic  
HELP\_TEXT.....No Help Text  
ID.....{53519F9A-B844-41C7-A013-F9D2EA8E4BC0}  
TERM.....BLACK/NECROTIC

DATA VALUE: BLADDER SCAN

-----  
ABBREVIATION.....Blad Scan  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Bladder Scan  
HELP\_TEXT.....  
ID.....{6C66BC97-1042-4E53-A49F-69AB7ABEE551}  
TERM.....BLADDER SCAN

DATA VALUE: BLANCHED

-----  
ABBREVIATION.....Blanched  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Blanched  
HELP\_TEXT.....  
ID.....{EBB2EA8D-04BC-49E7-812A-62B46F0F4678}  
TERM.....BLANCHED

DATA VALUE: BLAND

-----  
ABBREVIATION.....BLAND  
ACTIVE.....Yes  
DESCRIPTION.....BLAND DIET  
DISPLAY\_NAME.....BLAND  
HELP\_TEXT.....No Help Text  
ID.....{38930AFE-ECF7-4548-857C-854E1B5BE148}  
TERM.....BLAND

DATA VALUE: BLISTERING

-----  
ABBREVIATION.....BLISTERING  
ACTIVE.....Yes  
DESCRIPTION.....Blistering  
DISPLAY\_NAME.....Blistering  
HELP\_TEXT.....No Help Text  
ID.....{3A4212D5-F2A8-4636-BB7F-F36FA8CED192}  
TERM.....BLISTERING

DATA VALUE: BLOOD GLUCOSE MONITORING

-----  
ABBREVIATION.....BI Glu Mon  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Blood Glucose Monitoring  
HELP\_TEXT.....  
ID.....{525C8F3C-F8F5-4158-9CDE-3A7B306BAD5}  
TERM.....BLOOD GLUCOSE MONITORING

DATA VALUE: BLOOD PRESSURE

-----  
ABBREVIATION.....Blood Pres  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Blood Pressure  
HELP\_TEXT.....  
ID.....{D6E2D19A-0B74-476A-90BF-B7C98CC7361A}  
TERM.....BLOOD PRESSURE

DATA VALUE: BLOOD STREAKED

-----  
ABBREVIATION.....BLOOD STRE  
ACTIVE.....Yes  
DESCRIPTION.....Blood Streaked  
DISPLAY\_NAME.....Blood Streaked  
HELP\_TEXT.....Blood Streaked  
ID.....{AEC5607D-806F-49AE-8436-DC47F49F39BC}  
TERM.....BLOOD STREAKED

DATA VALUE: BLOOD TRANSFUSION IV

-----  
ABBREVIATION.....BI TRAN IV  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Blood Transfusion IV  
HELP\_TEXT.....No Help Text  
ID.....{27A0BE4C-433C-F920-BAC2-066F304BDF59}  
TERM.....BLOOD TRANSFUSION IV

DATA VALUE: BLOODY

-----  
ABBREVIATION.....BLOODY  
ACTIVE.....Yes  
DESCRIPTION.....Containing blood  
DISPLAY\_NAME.....Bloody  
HELP\_TEXT.....No Help Text  
ID.....{D7904E84-8331-D6DD-F2C0-7F3039C5746A}  
TERM.....BLOODY

DATA VALUE: BLUE

-----  
ABBREVIATION.....BLUE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Blue  
HELP\_TEXT.....No Help Text  
ID.....{8FA78575-E1DD-4E6E-BC43-9E5C92E2B7E5}  
TERM.....BLUE

DATA VALUE: BLUNTED

-----  
ABBREVIATION.....Blunted  
ACTIVE.....Yes  
DESCRIPTION.....Blunted  
DISPLAY\_NAME.....Blunted  
HELP\_TEXT.....No Help Text  
ID.....{F119F2DA-9182-2480-BA75-02D39F4FF69B}  
TERM.....BLUNTED

DATA VALUE: BLURRED VISION

-----  
ABBREVIATION.....Blurred Vi  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Blurred Vision  
HELP\_TEXT.....  
ID.....{1CD10C9B-28E5-49BA-AA6A-0FA3686A9E85}  
TERM.....BLURRED VISION



DATA VALUE: BMI>35KG/M2

-----  
ABBREVIATION.....BMI>35KG/M  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....BMI > 35kg/m2  
HELP\_TEXT.....  
ID.....{FD3D1D68-DD58-44B5-BEB2-0D4324F8353E}  
TERM.....BMI>35KG/M2

DATA VALUE: BONE

-----  
ABBREVIATION.....BONE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bone  
HELP\_TEXT.....No Help Text  
ID.....{3729DEB9-0967-4B32-86B7-409CF27C8625}  
TERM.....BONE

DATA VALUE: BOOST

-----  
ABBREVIATION.....Boost  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Boost  
HELP\_TEXT.....  
ID.....{2D15A3E9-41D2-4099-ADDC-7A1FAD770FF5}  
TERM.....BOOST

DATA VALUE: BOOST PLUS

-----  
ABBREVIATION.....Boost Plus  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Boost Plus  
HELP\_TEXT.....  
ID.....{432C8151-10CE-4804-939D-6B6B609BE54A}  
TERM.....BOOST PLUS

DATA VALUE: BORDERLINE REDUCED

-----  
ABBREVIATION.....BORD REDUC  
ACTIVE.....Yes  
DESCRIPTION.....Borderline Reduced  
DISPLAY\_NAME.....Borderline Reduced  
HELP\_TEXT.....No Help Text  
ID.....{63703AB1-6249-403B-BDE8-FECA9FA36CF5}  
TERM.....BORDERLINE REDUCED

DATA VALUE: BRADEN FRICTION AND SHEAR 1 PROBLEM

-----  
ABBREVIATION.....BS F/S 1  
ACTIVE.....Yes  
DESCRIPTION.....Braden Scale Friction and Shear 1  
                  Problem  
DISPLAY\_NAME.....1 - Problem  
HELP\_TEXT.....Requires moderate to maximum  
                  assistance in moving. Complete  
                  lifting without sliding against  
                  sheets is impossible. Frequently  
                  slides down in bed or chair,  
                  requiring frequent repositioning with  
                  maximum assistance.  
ID.....{AF8D2952-E0DA-411E-91CC-0114F4EB1C33}  
TERM.....BRADEN FRICTION AND SHEAR 1 PROBLEM

DATA VALUE: BRADEN SCALE ACTIVITY 1 BEDFAST

-----  
ABBREVIATION.....BS ACT 1  
ACTIVE.....Yes  
DESCRIPTION.....Braden Scale Activity 1 Bedfast  
DISPLAY\_NAME.....1 - Bedfast  
HELP\_TEXT.....Confined to bed  
ID.....{B9382966-5DE9-4D1C-BEC7-4DEBC5D6DCF5}  
TERM.....BRADEN SCALE ACTIVITY 1 BEDFAST

DATA VALUE: BRADEN SCALE ACTIVITY 2 CHAIRFAST

-----  
ABBREVIATION.....BS ACT 2  
ACTIVE.....Yes  
DESCRIPTION.....Braden Scale Activity 2 Chairfast  
DISPLAY\_NAME.....2 - Chairfast  
HELP\_TEXT.....Ability to walk severely limited or  
                  non-existent. Cannot bear own weight  
                  and/or must be assisted into chair or  
                  wheelchair.  
ID.....{CDBB8C3B-B53A-4A0D-91ED-8A53CFE1485F}  
TERM.....BRADEN SCALE ACTIVITY 2 CHAIRFAST

DATA VALUE: BRADEN SCALE ACTIVITY 3 WALKS OCCASIONALLY

-----  
ABBREVIATION.....BS ACT 3  
ACTIVE.....Yes  
DESCRIPTION.....Braden Scale Activity 3 Walks  
                  Occasionally  
DISPLAY\_NAME.....3 - Walks Occasionally  
HELP\_TEXT.....Walks occasionally during day, but  
                  for very short distances, with or  
                  without assistance. Spends majority  
                  of each shift in bed or chair.  
ID.....{E47B47BC-6097-4602-9184-BF5A4F6D686A}  
TERM.....BRADEN SCALE ACTIVITY 3 WALKS  
                  OCCASIONALLY

DATA VALUE: BRADEN SCALE ACTIVITY 4 WALKS FREQUENTLY

-----  
ABBREVIATION.....BS ACT 4  
ACTIVE.....Yes  
DESCRIPTION.....Braden Scale Activity 4 Walks

Frequently  
DISPLAY\_NAME.....4 - Walks Frequently  
HELP\_TEXT.....Walks outside room at least twice a  
day and inside room at least once  
every two hours during waking hours.  
ID.....{C5A8C2EA-E042-4804-9812-77603C47FC6A}  
TERM.....BRADEN SCALE ACTIVITY 4 WALKS  
FREQUENTLY

DATA VALUE: BRADEN SCALE FRICTION AND SHEAR 2 POTENTIAL PROBLEM

-----  
ABBREVIATION.....BS F/S 2  
ACTIVE.....Yes  
DESCRIPTION.....Braden Scale Friction and Shear 2  
Potential Problem  
DISPLAY\_NAME.....2 - Potential Problem  
HELP\_TEXT.....Moves feebly or requires minimum  
assistance. During a move skin  
probably slides to some extent  
against sheets, chair, restraints or  
other devices. Maintains relatively  
good position in chair or bed most of  
the time but occasionally slides down  
ID.....{B56338A9-8956-4D27-8BA0-DA50E07359E0}  
TERM.....BRADEN SCALE FRICTION AND SHEAR 2  
POTENTIAL PROBLEM

DATA VALUE: BRADEN SCALE FRICTION AND SHEAR 3 NO APPARENT PROBLEM

-----  
ABBREVIATION.....BS F/S 3  
ACTIVE.....Yes  
DESCRIPTION.....Braden Scale Friction and Shear 3 No  
Apparent Problem  
DISPLAY\_NAME.....3 - No Apparent Problem  
HELP\_TEXT.....No Apparent Problem  
ID.....{9850E2BF-6600-4E66-8FCD-1517D903BC0F}  
TERM.....BRADEN SCALE FRICTION AND SHEAR 3 NO  
APPARENT PROBLEM

DATA VALUE: BRADEN SCALE MOBILITY 1 COMPLETELY IMMOBILE

-----  
ABBREVIATION.....BS MOB 1  
ACTIVE.....Yes  
DESCRIPTION.....Braden Scale Mobility 1 Completely  
Immobile  
DISPLAY\_NAME.....1 - Completely Immobile  
HELP\_TEXT.....Does not make even slight changes in  
body or extremity position without  
assistance  
ID.....{7607574E-D4C2-4E64-8DC5-ABE7D00917A6}  
TERM.....BRADEN SCALE MOBILITY 1 COMPLETELY  
IMMOBILE

DATA VALUE: BRADEN SCALE MOBILITY 2 VERY LIMITED

-----  
ABBREVIATION.....BS MOB 2  
ACTIVE.....Yes  
DESCRIPTION.....Braden Scale Mobility 2 Very Limited  
DISPLAY\_NAME.....2 - Very Limited  
HELP\_TEXT.....Makes occasional slight changes in  
body or extremity position but unable  
to make frequent or significant  
changes independently  
ID.....{868EFA33-877E-4B91-B4FD-9655F66A7160}  
TERM.....BRADEN SCALE MOBILITY 2 VERY LIMITED

DATA VALUE: BRADEN SCALE MOBILITY 3 SLIGHTLY LIMITED

-----  
ABBREVIATION.....BS MOB 3

ACTIVE.....Yes  
DESCRIPTION.....Braden Scale Mobility 3 Slightly  
Limited  
DISPLAY\_NAME.....3 - Slightly Limited  
HELP\_TEXT.....Makes frequent though slight changes  
in body or extremity position  
independently  
ID.....{C54B0A8F-46F1-4DFA-9C24-AB6C2A29CA82}  
TERM.....BRADEN SCALE MOBILITY 3 SLIGHTLY  
LIMITED

DATA VALUE: BRADEN SCALE MOBILITY 4 NO LIMITATION

-----  
ABBREVIATION.....BS MOB 4  
ACTIVE.....Yes  
DESCRIPTION.....Braden Scale Mobility 4 No Limitation  
DISPLAY\_NAME.....4 - No Limitation  
HELP\_TEXT.....Makes major and frequent changes in  
                  position without assistance  
ID.....{84BFB8D0-D6E4-4384-88EE-28B7C4EDB275}  
TERM.....BRADEN SCALE MOBILITY 4 NO LIMITATION

DATA VALUE: BRADEN SCALE MOISTURE 1 CONSTANTLY MOIST

-----  
ABBREVIATION.....BS MOI 1  
ACTIVE.....Yes  
DESCRIPTION.....Braden Scale Moisture 1 Constantly  
                  Moist  
DISPLAY\_NAME.....1 - Constantly Moist  
HELP\_TEXT.....Skin is kept moist almost constantly  
                  by perspiration, urine, etc. Dampness  
                  is detected every time patient is  
                  moved or turned  
ID.....{3E2569BA-ECF0-4622-8F4A-26378A2CB466}  
TERM.....BRADEN SCALE MOISTURE 1 CONSTANTLY  
                  MOIST

DATA VALUE: BRADEN SCALE MOISTURE 2 VERY MOIST

-----  
ABBREVIATION.....BS MOI 2  
ACTIVE.....Yes  
DESCRIPTION.....Braden Scale Moisture 2 Very Moist  
DISPLAY\_NAME.....2 - Very Moist  
HELP\_TEXT.....Skin is often, but not always moist.  
                  Linen must be changed at least once a  
                  shift  
ID.....{2C6D15A3-6166-4D7D-B200-B58EC9C2CABE}  
TERM.....BRADEN SCALE MOISTURE 2 VERY MOIST

DATA VALUE: BRADEN SCALE MOISTURE 3 OCCASIONALLY MOIST

-----  
ABBREVIATION.....BS MOI 3  
ACTIVE.....Yes  
DESCRIPTION.....Braden Scale Moisture 3 Occasionally  
                  Moist  
DISPLAY\_NAME.....3 - Occasionally Moist  
HELP\_TEXT.....Skin is occasionally moist, requiring  
                  an extra linen change approximately  
                  once a day  
ID.....{D7F2714B-B0C7-4000-8588-8582C1B32ECE}  
TERM.....BRADEN SCALE MOISTURE 3 OCCASIONALLY  
                  MOIST

DATA VALUE: BRADEN SCALE MOISTURE 4 RARELY MOIST

-----  
ABBREVIATION.....BS MOI 4  
ACTIVE.....Yes  
DESCRIPTION.....Braden Scale Moisture 4 Rarely Moist  
DISPLAY\_NAME.....4 - Rarely Moist  
HELP\_TEXT.....Skin is usually dry, linen only  
                  requires changing at routine

intervals

ID.....{745EA3B9-18C9-4CF3-9DB8-9806A02656E8}

TERM.....BRADEN SCALE MOISTURE 4 RARELY MOIST

DATA VALUE: BRADEN SCALE NUTRITION 1 VERY POOR

-----  
ABBREVIATION.....BS NUT 1  
ACTIVE.....Yes  
DESCRIPTION.....Braden Scale Nutrition 1 Very Poor  
DISPLAY\_NAME.....1 - Very Poor  
HELP\_TEXT.....Never eats a complete meal. Rarely  
                  eats more than 1/2 of any food  
                  offered. Eats 2 servings or less of  
                  protein per day. Takes fluids poorly,  
                  does not take a liquid supplement. OR  
                  is NPO maintained on clear liquids or  
                  IV's for more than 5 days  
ID.....{28A86074-A6B2-455E-ABA4-19195C2E3679}  
TERM.....BRADEN SCALE NUTRITION 1 VERY POOR

DATA VALUE: BRADEN SCALE NUTRITION 2 PROBABLY INADEQUATE

-----  
ABBREVIATION.....BS NUT 2  
ACTIVE.....Yes  
DESCRIPTION.....Braden Scale Nutrition 2 Probably  
                  Inadequate  
DISPLAY\_NAME.....2 - Probably Inadequate  
HELP\_TEXT.....Rarely eats a complete meal and  
                  generally eats only about 2 of any  
                  food offered. Protein intake includes  
                  only 3 servings of meat or dairy  
                  products per day. Occasionally will  
                  take a dietary supplement. OR  
                  receives less than optimum amount of  
                  liquid.  
ID.....{7DAB22F6-22E4-4A8B-A6A5-CD1277EB6B87}  
TERM.....BRADEN SCALE NUTRITION 2 PROBABLY  
                  INADEQUATE

DATA VALUE: BRADEN SCALE NUTRITION 3 ADEQUATE

-----  
ABBREVIATION.....BS NUT 3  
ACTIVE.....Yes  
DESCRIPTION.....Braden Scale Nutrition 3 Adequate  
DISPLAY\_NAME.....3 - Adequate  
HELP\_TEXT.....Eats over half of most meals. Eats a  
                  total of 4 servings of protein (meat,  
                  dairy products) per day. Occasionally  
                  will refuse a meal, but will usually  
                  take a supplement when offered. OR is  
                  on a tube feeding or TPN regimen.  
ID.....{3A8A9C59-8DC5-49A9-BAC9-882CB75D8FD5}  
TERM.....BRADEN SCALE NUTRITION 3 ADEQUATE

DATA VALUE: BRADEN SCALE NUTRITION 4 EXCELLENT

-----  
ABBREVIATION.....BS NUT 4  
ACTIVE.....Yes  
DESCRIPTION.....Braden Scale Nutrition 4 Excellent  
DISPLAY\_NAME.....4 - Excellent  
HELP\_TEXT.....Eats most of every meal. Never  
                  refuses a meal. Usually eats a total  
                  of 4 or more servings of meat and  
                  dairy products. Occasionally eats



between meals. Does NOT require  
supplements  
ID.....{F09E4D35-1AA8-49A5-9B38-FBFFAFCFC181}  
TERM.....BRADEN SCALE NUTRITION 4 EXCELLENT

DATA VALUE: BRADEN SCALE SENSORY PERCEPTION 1 COMPLETELY LIMITED

-----  
ABBREVIATION.....BS SENC 1  
ACTIVE.....Yes  
DESCRIPTION.....Braden Scale Sensory Perception 1  
                  Completely Limited  
DISPLAY\_NAME.....1 - Completely Limited  
HELP\_TEXT.....Unresponsive (does not moan, flinch  
                  or grasp) to painful stimuli, due to  
                  diminished level of consciousness or  
                  sedation. OR limited ability to feel  
                  pain over most of body  
ID.....{60F646B3-790F-4C5A-B2A0-845915E120E6}  
TERM.....BRADEN SCALE SENSORY PERCEPTION 1  
                  COMPLETELY LIMITED

DATA VALUE: BRADEN SCALE SENSORY PERCEPTION 2 VERY LIMITED

-----  
ABBREVIATION.....BS SENC 2  
ACTIVE.....Yes  
DESCRIPTION.....Braden Scale Sensory Perception 2  
Very Limited  
DISPLAY\_NAME.....2 - Very Limited  
HELP\_TEXT.....Responds only to painful stimuli.  
Cannot Communicate discomfort except  
by moaning or restlessness. OR has a  
sensory impairment which limits the  
ability to feel pain or discomfort  
over 1/2 of the body  
ID.....{98362E63-7245-4E7F-8C30-896A9B923199}  
TERM.....BRADEN SCALE SENSORY PERCEPTION 2  
VERY LIMITED

DATA VALUE: BRADEN SCALE SENSORY PERCEPTION 3 SLIGHTLY LIMITED

-----  
ABBREVIATION.....BS SENC 3  
ACTIVE.....Yes  
DESCRIPTION.....Braden Scale Sensory Perception 3  
Slightly Limited  
DISPLAY\_NAME.....3 - Slightly Limited  
HELP\_TEXT.....Responds to verbal commands, but  
cannot always communicate discomfort  
or the need to be turned. OR has some  
sensory impairment which limits  
ability to feel pain or discomfort in  
one or two extremities.  
ID.....{1D01D47D-CDC2-40EE-837B-5417D6EE47D6}  
TERM.....BRADEN SCALE SENSORY PERCEPTION 3  
SLIGHTLY LIMITED

DATA VALUE: BRADEN SCALE SENSORY PERCEPTION 4 NO IMPAIRMENT

-----  
ABBREVIATION.....BS SENC 4  
ACTIVE.....Yes  
DESCRIPTION.....Braden Scale Sensory Perception 4 No  
Impairment  
DISPLAY\_NAME.....4 - No Impairment  
HELP\_TEXT.....Responds to verbal commands. Has no  
sensory deficit which would limit  
ability to feel or voice pain or  
discomfort.  
ID.....{10304656-E5E8-4F32-8358-A9F80511EFDA}  
TERM.....BRADEN SCALE SENSORY PERCEPTION 4 NO  
IMPAIRMENT

DATA VALUE: BREAKDOWN

-----  
ABBREVIATION.....BREAKDOWN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Breakdown  
HELP\_TEXT.....Use comment to describe breakdown  
ID.....{9CFD21C0-A805-47C4-E559-8033C754DE80}  
TERM.....BREAKDOWN

DATA VALUE: BRIDGE PRESENT

-----  
ABBREVIATION.....bridge pre  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Bridge Present  
HELP\_TEXT.....  
ID.....{5C474C17-33C9-492D-833E-64F5EF049BF4}  
TERM.....BRIDGE PRESENT

DATA VALUE: BRIGHT

-----  
ABBREVIATION.....Bright  
ACTIVE.....Yes  
DESCRIPTION.....Bright  
DISPLAY\_NAME.....Bright  
HELP\_TEXT.....No Help Text  
ID.....{D16542CE-0403-2C68-A361-125CB3C75ABA}  
TERM.....BRIGHT

DATA VALUE: BRISK

-----  
ABBREVIATION.....brisk  
ACTIVE.....Yes  
DESCRIPTION.....Brisk  
DISPLAY\_NAME.....Brisk  
HELP\_TEXT.....No Help Text  
ID.....{980AE687-0A78-458B-BA7B-5D646C38A5F5}  
TERM.....BRISK

DATA VALUE: BROWN

-----  
ABBREVIATION.....BROWN  
ACTIVE.....Yes  
DESCRIPTION.....The color brown  
DISPLAY\_NAME.....Brown  
HELP\_TEXT.....The color brown  
ID.....{A1C32FBC-7B4D-707C-CA58-92A2E201368E}  
TERM.....BROWN

DATA VALUE: BRUISED

-----  
ABBREVIATION.....BRUISED  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bruised  
HELP\_TEXT.....No Help Text  
ID.....{FA5B38A7-7744-186F-055A-8458D696AF33}  
TERM.....BRUISED

DATA VALUE: BRUIT

-----  
ABBREVIATION.....BRUIT  
ACTIVE.....Yes  
DESCRIPTION.....Bruit  
DISPLAY\_NAME.....Bruit  
HELP\_TEXT.....No Help Text  
ID.....{2627E78A-B8B3-44EE-BC6F-4DAB16394F8C}  
TERM.....BRUIT

DATA VALUE: BUCKS TRACTION

-----  
ABBREVIATION.....BUCK TRACT  
ACTIVE.....Yes  
DESCRIPTION.....Buck's Traction  
DISPLAY\_NAME.....Buck's Traction  
HELP\_TEXT.....Buck's Traction  
ID.....{5E63F920-02D7-4729-8F30-A53A0A6E1EA0}  
TERM.....BUCKS TRACTION

DATA VALUE: BULB

-----  
ABBREVIATION.....Bulb  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Bulb  
HELP\_TEXT.....  
ID.....{2390453E-3353-40EB-B2D7-EF770B67F037}  
TERM.....BULB

DATA VALUE: BULB SUCTION

-----  
ABBREVIATION.....Bulb Sucti  
ACTIVE.....Yes  
DESCRIPTION.....Bulb Suction Drain  
DISPLAY\_NAME.....Bulb Suction Drain  
HELP\_TEXT.....No Help Text  
ID.....{B0839A82-BB79-BB9F-5E4E-0ED67F716A50}  
TERM.....BULB SUCTION

DATA VALUE: BULLET

-----  
ABBREVIATION.....BULLET  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bullet  
HELP\_TEXT.....No Help Text  
ID.....{1570D8E2-D936-494D-AF6E-13838EA34EFD}  
TERM.....BULLET

DATA VALUE: BUPIVICAINE

-----  
ABBREVIATION.....Bupivicain  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Bupivicaïne  
HELP\_TEXT.....  
ID.....{03B908B3-365A-47F6-8D71-8F8FE1D2C16B}  
TERM.....BUPIVICAINE

DATA VALUE: BURN

-----  
ABBREVIATION.....BURN  
ACTIVE.....Yes  
DESCRIPTION.....Burn  
DISPLAY\_NAME.....Burn  
HELP\_TEXT.....No Help Text  
ID.....{97B37D7A-28C6-4ABF-8D09-1FC49D5F160F}  
TERM.....BURN

DATA VALUE: BURNING

-----  
ABBREVIATION.....BURNING  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Burning  
HELP\_TEXT.....No Help Text  
ID.....{D75E2D5D-9CF9-4E70-A2EB-0632EBDB82C9}  
TERM.....BURNING

DATA VALUE: BY PROVIDER ORDER

-----  
ABBREVIATION.....BY ORDER  
ACTIVE.....Yes  
DESCRIPTION.....By Provider Order  
DISPLAY\_NAME.....By Provider Order  
HELP\_TEXT.....No Help Text  
ID.....{C038BC16-44CF-4933-BECC-B0F1A1B12603}  
TERM.....BY PROVIDER ORDER

DATA VALUE: C1

-----  
ABBREVIATION.....C1  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....C1  
HELP\_TEXT.....  
ID.....{F942693D-F272-4A70-81FF-E750BAF082B9}  
TERM.....C1

DATA VALUE: C2

-----  
ABBREVIATION.....C2  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....C2  
HELP\_TEXT.....  
ID.....{EE1C5D26-12DE-42ED-A46C-63D66BC145C5}  
TERM.....C2

DATA VALUE: C3

-----  
ABBREVIATION.....C3  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....C3  
HELP\_TEXT.....  
ID.....{B53D0E0E-7BF0-4D26-B4DD-D072035E307E}  
TERM.....C3

DATA VALUE: C4

-----  
ABBREVIATION.....C4  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....C4  
HELP\_TEXT.....  
ID.....{A7B5DC47-8E3E-4B23-826C-7F7D5FE0C6BF}  
TERM.....C4

DATA VALUE: C5

-----  
ABBREVIATION.....C5  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....C5  
HELP\_TEXT.....  
ID.....{30EA5569-8E29-4B5E-ABA3-598B7C39A8D4}  
TERM.....C5

DATA VALUE: CADEXOMER IODINE

-----  
ABBREVIATION.....CADEXO IOD  
ACTIVE.....Yes  
DESCRIPTION.....Cadexomer Iodine  
DISPLAY\_NAME.....Cadexomer Iodine  
HELP\_TEXT.....No Help Text  
ID.....{C83E41AA-6A1F-43CC-AB2E-7F6555FECD92}  
TERM.....CADEXOMER IODINE

DATA VALUE: CALCIUM ALGINATE

-----  
ABBREVIATION.....CA ALGINAT  
ACTIVE.....Yes  
DESCRIPTION.....Calcium Alginate  
DISPLAY\_NAME.....Calcium Alginate  
HELP\_TEXT.....No Help Text  
ID.....{2B2ACD25-A77B-4D39-8598-B26FD9BAFFC2}  
TERM.....CALCIUM ALGINATE

DATA VALUE: CALM

-----  
ABBREVIATION.....c  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Calm  
HELP\_TEXT.....No Help Text  
ID.....{E206C3F5-6977-302C-51FB-AF2155B5FFD9}  
TERM.....CALM

DATA VALUE: CANE

-----  
ABBREVIATION.....CANE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Cane  
HELP\_TEXT.....No Help Text  
ID.....{562B2E0B-0298-44C6-B33B-BC57DE6CBC3C}  
TERM.....CANE

DATA VALUE: CANES

-----  
ABBREVIATION.....CANES  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Canes  
HELP\_TEXT.....No Help Text  
ID.....{DF37C541-1C7D-4527-A4D2-31867B5ED438}  
TERM.....CANES

DATA VALUE: CARDIAC

-----  
ABBREVIATION.....Cardiac  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Cardiac  
HELP\_TEXT.....  
ID.....{2A4308E5-C843-4844-B908-32DA0A0B50E3}  
TERM.....CARDIAC

DATA VALUE: CARDIAC CHAIR

-----  
ABBREVIATION.....cc  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Cardiac Chair  
HELP\_TEXT.....No Help Text  
ID.....{D4A2CFC1-D3A4-696B-4062-AFB2DAFFD76A}  
TERM.....CARDIAC CHAIR



DATA VALUE: CARDIAC DIET

-----  
ABBREVIATION.....CARD DIET  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....CARDIAC DIET  
HELP\_TEXT.....No Help Text  
ID.....{218A883D-32A0-438D-B453-7CBCE02693A3}  
TERM.....CARDIAC DIET

DATA VALUE: CARDIAC OUTPUT INJECTATE

-----  
ABBREVIATION.....CARD OUT I  
ACTIVE.....Yes  
DESCRIPTION.....Cardiac Output Injectate  
DISPLAY\_NAME.....Cardiac Output Injectate  
HELP\_TEXT.....No Help Text  
ID.....{19605CAB-8D54-4704-819B-20ABDAB6C604}  
TERM.....CARDIAC OUTPUT INJECTATE

DATA VALUE: CARDIAC REPOSITION

-----  
ABBREVIATION.....CARD REPOS  
ACTIVE.....Yes  
DESCRIPTION.....Cardiac Reposition  
DISPLAY\_NAME.....Cardiac Reposition  
HELP\_TEXT.....No Help Text  
ID.....{E133FE5A-20E1-418F-BD90-D2DFE14908C0}  
TERM.....CARDIAC REPOSITION

DATA VALUE: CARDIO-RESP MONITOR

-----  
ABBREVIATION.....Cardio-Re  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Cardio-Resp Monitor  
HELP\_TEXT.....  
ID.....{4402EC2C-2096-4DF7-BD4A-AE66C67035E0}  
TERM.....CARDIO-RESP MONITOR

DATA VALUE: CARDIOLOGY PROCEDURE ROOM

-----  
ABBREVIATION.....Cardio Rm  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Cardiology Procedure Room  
HELP\_TEXT.....  
ID.....{744142F5-D301-4509-A96D-06AB8C1D2C81}  
TERM.....CARDIOLOGY PROCEDURE ROOM

DATA VALUE: CAREGIVER

-----  
ABBREVIATION.....CAREGIVER  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Caregiver  
HELP\_TEXT.....No Help Text  
ID.....{F9B8BABE-2F05-437D-85C8-585FC1FFC07F}  
TERM.....CAREGIVER

DATA VALUE: CAREGIVER STRESS

-----  
ABBREVIATION.....Care Stres  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Caregiver Stress  
HELP\_TEXT.....  
ID.....{A6A201E7-259C-485E-A35E-9F2F3489B24F}  
TERM.....CAREGIVER STRESS

DATA VALUE: CASS ENDOTRACHEAL TUBE

-----  
ABBREVIATION.....CASS  
ACTIVE.....Yes  
DESCRIPTION.....CASS Endotracheal Tube  
DISPLAY\_NAME.....CASS Endotracheal Tube  
HELP\_TEXT.....Mucous Slurping Endotracheal Tube  
ID.....{02BF805C-DF35-4DAB-A99F-2ED673174E49}  
TERM.....CASS ENDOTRACHEAL TUBE

DATA VALUE: CATH SECURED WITH DRAIN BAG BELOW PT'S BLADDER LEVEL

-----  
ABBREVIATION.....Cath Secur  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Cath Secured with Drain Bag Below  
Pt's Bladder Level  
HELP\_TEXT.....  
ID.....{A2D2ABF1-88D4-46A2-B935-33211292FF25}  
TERM.....CATH SECURED WITH DRAIN BAG BELOW  
PT'S BLADDER LEVEL

DATA VALUE: CATH SIZE DOCUMENTED

-----  
ABBREVIATION.....Cath Size  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Cath Size Documented  
HELP\_TEXT.....  
ID.....{FB8BE4D4-41C7-4F8C-AD17-74F8FE170BC3}  
TERM.....CATH SIZE DOCUMENTED

DATA VALUE: CATHETER

-----  
ABBREVIATION.....CATHETER  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Catheter  
HELP\_TEXT.....No Help Text  
ID.....{B0D1AAD7-012D-D3E8-3216-1F71A1EC685E}  
TERM.....CATHETER

DATA VALUE: CATHETERIZATION TIME DOCUMENTED

-----  
ABBREVIATION.....Cath Time  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME..... Catheterization Time Documented  
HELP\_TEXT.....  
ID.....{4EF58FCA-A5E1-454C-B435-5BFB487618F4}  
TERM.....CATHETERIZATION TIME DOCUMENTED

DATA VALUE: CELLULITIS

-----  
ABBREVIATION.....CELLULITIS  
ACTIVE.....Yes  
DESCRIPTION.....Cellulitis  
DISPLAY\_NAME.....Cellulitis  
HELP\_TEXT.....No Help Text  
ID.....{709F8795-0D4A-4E98-AB82-DD3E2AA7EA7B}  
TERM.....CELLULITIS

DATA VALUE: CENTRAL ECG

-----  
ABBREVIATION.....cecg  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Central ECG  
HELP\_TEXT.....No Help Text  
ID.....{26D36175-E4E3-C364-3018-1B0DB4DC4793}  
TERM.....CENTRAL ECG

DATA VALUE: CENTRAL VENOUS CATHETER

-----  
ABBREVIATION.....CENT VEN C  
ACTIVE.....Yes  
DESCRIPTION.....Central Venous Catheter  
DISPLAY\_NAME.....Central Venous Catheter  
HELP\_TEXT.....No Help Text  
ID.....{EE91E841-C09B-4432-A121-8EBAB92E9CD2}  
TERM.....CENTRAL VENOUS CATHETER

DATA VALUE: CERVICAL SPINAL PRECAUTIONS

-----  
ABBREVIATION.....CERV S P  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Cervical Spinal Precautions  
HELP\_TEXT.....No Help Text  
ID.....{A15F7127-2367-464A-AB7D-1E6A7C09AA21}  
TERM.....CERVICAL SPINAL PRECAUTIONS

DATA VALUE: CF: COUNTING FINGERS

-----  
ABBREVIATION.....CF  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....CF: Counting Fingers  
HELP\_TEXT.....  
ID.....{3AFD9EDB-98A2-498A-AD95-F5285130CC03}  
TERM.....CF: COUNTING FINGERS

DATA VALUE: CHAIR

-----  
ABBREVIATION.....CHAIR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Chair  
HELP\_TEXT.....No Help Text  
ID.....{5B8A5DBD-7462-4653-8D1F-CA8B1226DE38}  
TERM.....CHAIR

DATA VALUE: CHANGED

-----  
ABBREVIATION.....CHANGED  
ACTIVE.....Yes  
DESCRIPTION.....Changed  
DISPLAY\_NAME.....Changed  
HELP\_TEXT.....No Help Text  
ID.....{74B14023-D74D-41F6-AA48-B0E5137623CD}  
TERM.....CHANGED

DATA VALUE: CHART VERIFIED CORRECT

-----  
ABBREVIATION.....Chart Ver  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Chart Verified Correct  
HELP\_TEXT.....  
ID.....{319E702F-C7A9-409D-881F-3F4D1D9AF393}  
TERM.....CHART VERIFIED CORRECT

DATA VALUE: CHEYNE-STOKES

-----  
ABBREVIATION.....CHEYNE-STO  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Cheyne-Stokes  
HELP\_TEXT.....No Help Text  
ID.....{58DD2B80-B522-499C-2C60-C48A654BAD84}  
TERM.....CHEYNE-STOKES

DATA VALUE: CIRCULATION CHECK

-----  
ABBREVIATION.....ccq  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Circulation Check  
HELP\_TEXT.....No Help Text  
ID.....{BA32C2A8-F15C-84EF-4859-6277EFA6463A}  
TERM.....CIRCULATION CHECK

DATA VALUE: CLAMPED

-----  
ABBREVIATION.....CLAMPED  
ACTIVE.....Yes  
DESCRIPTION.....Clamped  
DISPLAY\_NAME.....Clamped  
HELP\_TEXT.....No Help Text  
ID.....{B86865B8-078A-4996-B5EC-46A4A99D5A7D}  
TERM.....CLAMPED

DATA VALUE: CLASS I

-----  
ABBREVIATION.....Class I  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Class I  
HELP\_TEXT.....  
ID.....{4A1B7DD6-3846-4455-A06D-0E3C0FB3752A}  
TERM.....CLASS I

DATA VALUE: CLASS II

-----  
ABBREVIATION.....Class II  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Class II  
HELP\_TEXT.....  
ID.....{2F5F3798-3EFD-4DEA-9FAB-883F9E326EA7}  
TERM.....CLASS II

DATA VALUE: CLASS III

-----  
ABBREVIATION.....Class III  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Class III  
HELP\_TEXT.....  
ID.....{C0FFDE12-241C-4087-9404-DB128F6B793A}  
TERM.....CLASS III

DATA VALUE: CLASS IV

-----  
ABBREVIATION.....Class IV  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Class IV  
HELP\_TEXT.....  
ID.....{176030D0-2576-4D2D-A5AB-E5DA38DADEB5}  
TERM.....CLASS IV

DATA VALUE: CLAY

-----  
ABBREVIATION.....CLAY  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Clay  
HELP\_TEXT.....No Help Text  
ID.....{A4E4D16C-B45C-0934-ED22-8BF640080108}  
TERM.....CLAY

DATA VALUE: CLEAN

-----  
ABBREVIATION.....Clean  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Clean  
HELP\_TEXT.....  
ID.....{1F187B99-9963-4AD5-98F6-D8D24D2DDD30}  
TERM.....CLEAN

DATA VALUE: CLEAR

-----  
ABBREVIATION.....c  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Clear  
HELP\_TEXT.....No Help Text  
ID.....{8981C365-2B46-704F-C9F4-975DF1C23C2D}  
TERM.....CLEAR

DATA VALUE: CLEAR LIQUID DIET

-----  
ABBREVIATION.....CL LI DIET  
ACTIVE.....Yes  
DESCRIPTION.....Clear Liquid Diet  
DISPLAY\_NAME.....Clear Liquid Diet  
HELP\_TEXT.....No Help Text  
ID.....{D540348D-3E1D-4E42-9EDC-43C56CF56206}  
TERM.....CLEAR LIQUID DIET

DATA VALUE: CLICK

-----  
ABBREVIATION.....c  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Click  
HELP\_TEXT.....No Help Text  
ID.....{F1F3A1A8-35E1-61B4-C783-3314C8247380}  
TERM.....CLICK

DATA VALUE: CLINICAL ASSESSMENT

-----  
ABBREVIATION.....CLIN ASSME  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Clinical Assessment  
HELP\_TEXT.....No Help Text  
ID.....{803BEA4B-E22F-48C2-B20E-D5B3C8C7084C}  
TERM.....CLINICAL ASSESSMENT

DATA VALUE: CLOSE OBSERVATION

-----  
ABBREVIATION.....co  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Close Observation  
HELP\_TEXT.....No Help Text  
ID.....{2679BCAC-8753-FE73-56F8-3AB632161AF4}  
TERM.....CLOSE OBSERVATION

DATA VALUE: CLOSED

-----  
ABBREVIATION.....CLOSED  
ACTIVE.....Yes  
DESCRIPTION.....CLOSED  
DISPLAY\_NAME.....CLOSED  
HELP\_TEXT.....No Help Text  
ID.....{9BA6EE03-F562-4DE2-99A7-337F5A5E34F2}  
TERM.....CLOSED

DATA VALUE: CLOSED SYSTEM

-----  
ABBREVIATION.....Closed sys  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Closed System  
HELP\_TEXT.....  
ID.....{D076BAB6-C4D1-42F1-9036-80BF4B72E695}  
TERM.....CLOSED SYSTEM

DATA VALUE: CLOTHED SELF

-----  
ABBREVIATION.....CLOTH S  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Clothed Self  
HELP\_TEXT.....No Help Text  
ID.....{33802A30-FAA7-426B-ADAC-5297818DA7DE}  
TERM.....CLOTHED SELF

DATA VALUE: CLOTHING ROOM

-----  
ABBREVIATION.....Cloth Room  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Clothing Room  
HELP\_TEXT.....  
ID.....{8F116C60-10B1-4033-93F2-9121629CD608}  
TERM.....CLOTHING ROOM

DATA VALUE: CLOTS

-----  
ABBREVIATION.....CLOTS  
ACTIVE.....Yes  
DESCRIPTION.....Clots Present  
DISPLAY\_NAME.....Clots  
HELP\_TEXT.....No Help Text  
ID.....{B52B88E0-6D23-4657-97CD-81E3EF64991B}  
TERM.....CLOTS

DATA VALUE: CLOUDY

-----  
ABBREVIATION.....CLOUDY  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Cloudy  
HELP\_TEXT.....No Help Text  
ID.....{3CCFCCCC-00FC-6B00-1D58-A77C02FDB401}  
TERM.....CLOUDY



DATA VALUE: CLUBBING

-----  
ABBREVIATION.....CLUBBING  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Clubbing  
HELP\_TEXT.....No Help Text  
ID.....{AE77D344-6F0F-452B-B1AA-EC8D9F092CFC}  
TERM.....CLUBBING

DATA VALUE: CMS INTACT DISTAL TO SURGICAL SITE

-----  
ABBREVIATION.....CMS INTA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....CMS Intact Distal to Surgical Site  
HELP\_TEXT.....No Help Text  
ID.....{CE2D2CDE-71C5-4993-8E46-C4799301A851}  
TERM.....CMS INTACT DISTAL TO SURGICAL SITE

DATA VALUE: CNS

-----  
ABBREVIATION.....CNS  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....CNS  
HELP\_TEXT.....  
ID.....{487A4DE0-A822-437D-87A2-75FA86D7BBDE}  
TERM.....CNS

DATA VALUE: CNS DEPRESSION

-----  
ABBREVIATION.....CNS Depres  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....CNS Depression  
HELP\_TEXT.....  
ID.....{EE349460-8B60-427B-AE98-72007DFD2F6A}  
TERM.....CNS DEPRESSION

DATA VALUE: CODE CART IS ASSESSIBLE

-----  
ABBREVIATION.....Code Cart  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Code Cart is Accessible  
HELP\_TEXT.....  
ID.....{DB670E39-CF90-4068-8E38-44DCA5F365B5}  
TERM.....CODE CART IS ASSESSIBLE

DATA VALUE: COFFEE GROUND

-----  
ABBREVIATION.....COFFEE GRO  
ACTIVE.....Yes  
DESCRIPTION.....Coffee Ground  
DISPLAY\_NAME.....Coffee Ground  
HELP\_TEXT.....Coffee Ground  
ID.....{A5B3119A-C739-4B8D-87C2-7D858B94F380}  
TERM.....COFFEE GROUND

DATA VALUE: COGNITIVE IMPAIRMENT

-----  
ABBREVIATION.....Cog Impair  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Cognitive Impairment  
HELP\_TEXT.....  
ID.....{6D051F44-2124-4914-A415-B98B2AC516E9}  
TERM.....COGNITIVE IMPAIRMENT

DATA VALUE: COLD

-----  
ABBREVIATION.....COLD  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Cold  
HELP\_TEXT.....No Help Text  
ID.....{B7B422BA-22EC-117D-0A0C-8419E9447D46}  
TERM.....COLD

DATA VALUE: COLLAGEN

-----  
ABBREVIATION.....COLLAGEN  
ACTIVE.....Yes  
DESCRIPTION.....Collagen  
DISPLAY\_NAME.....Collagen  
HELP\_TEXT.....No Help Text  
ID.....{8D734495-2FAF-41C0-9146-55B22583A830}  
TERM.....COLLAGEN

DATA VALUE: COLORLESS

-----  
ABBREVIATION.....COLORLESS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Colorless  
HELP\_TEXT.....No Help Text  
ID.....{D093D8C5-E284-47D4-8FA6-DCD5B030B10E}  
TERM.....COLORLESS

DATA VALUE: COLOSTOMY

-----  
ABBREVIATION.....COLOSTOMY  
ACTIVE.....Yes  
DESCRIPTION.....Colostomy  
DISPLAY\_NAME.....Colostomy  
HELP\_TEXT.....No Help Text  
ID.....{53F939AB-6AF4-4DD5-8D64-17AC1D1736B7}  
TERM.....COLOSTOMY

DATA VALUE: COMBATIVE

-----  
ABBREVIATION.....c  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Combative  
HELP\_TEXT.....No Help Text  
ID.....{DF65D42D-6CAA-657E-AEE1-C070C1FDF177}  
TERM.....COMBATIVE

DATA VALUE: COMFORT CARE

-----  
ABBREVIATION.....Comfort Ca  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Comfort Care  
HELP\_TEXT.....  
ID.....{51669E41-DF78-4E49-BA5F-DE6F1B5F1148}  
TERM.....COMFORT CARE

DATA VALUE: COMFORT MEASURES

-----  
ABBREVIATION.....Com Measu  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Comfort Measures  
HELP\_TEXT.....  
ID.....{A5F09A1F-921E-4469-9652-0A136BCD9C07}  
TERM.....COMFORT MEASURES

DATA VALUE: COMMUNE

-----  
ABBREVIATION.....COMMUNE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Commune  
HELP\_TEXT.....No Help-Text  
ID.....{75C896CA-F105-4EF0-B641-3F2B1F984F1B}  
TERM.....COMMUNE

DATA VALUE: COMMUNCATION BARRIER

-----  
ABBREVIATION.....COMM BARRI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Communication Barrier  
HELP\_TEXT.....No Help Text  
ID.....{1A24946D-D966-4FAB-B6EA-940FA9EF60D3}  
TERM.....COMMUNCATION BARRIER

DATA VALUE: COMMUNITY RESOURCES

-----  
ABBREVIATION.....Comm Res  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Community Resources  
HELP\_TEXT.....  
ID.....{0088604F-6A99-4D99-B451-5262734E0D17}  
TERM.....COMMUNITY RESOURCES

DATA VALUE: COMPLETE BATH

-----  
ABBREVIATION.....COMP BATH  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Complete Bath  
HELP\_TEXT.....No Help Text  
ID.....{A7CB1F56-F88E-403B-AABC-105BE129FBB0}  
TERM.....COMPLETE BATH

DATA VALUE: COMPLIANCE

-----  
ABBREVIATION.....COMPLIANC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Compliance  
HELP\_TEXT.....No Help Text  
ID.....{EF43367E-9898-44C3-8AC8-F5BD28051F1D}  
TERM.....COMPLIANCE

DATA VALUE: CONFUSED

-----  
ABBREVIATION.....CONFUSED  
ACTIVE.....Yes  
DESCRIPTION.....Confused  
DISPLAY\_NAME.....Confused  
HELP\_TEXT.....No Help Text  
ID.....{31DAAA60-BE89-49A7-A10D-EB0DA82F12C0}  
TERM.....CONFUSED

DATA VALUE: CONFUSION

-----  
ABBREVIATION.....Confusion  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Confusion  
HELP\_TEXT.....  
ID.....{4C1B8247-75C4-4C60-970B-3D7F07BD1164}  
TERM.....CONFUSION

DATA VALUE: CONNECTED TO DRAINAGE BAG

-----  
ABBREVIATION.....drain bag  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Connected to Drainage Bag  
HELP\_TEXT.....  
ID.....{1901C4A7-671C-43BD-BEDD-1D48BCB43465}  
TERM.....CONNECTED TO DRAINAGE BAG

DATA VALUE: CONSISTENT WITH BODY TEMPERATURE

-----  
ABBREVIATION.....CON W TEMP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Consistent With Body Temperature  
HELP\_TEXT.....No Help Text  
ID.....{941FC5EF-6AF8-4DA1-A43C-95575AC74B57}  
TERM.....CONSISTENT WITH BODY TEMPERATURE

DATA VALUE: CONSTANT

-----  
ABBREVIATION.....CONSTANT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Constant  
HELP\_TEXT.....No Help Text  
ID.....{0FB8E43C-2106-4C37-9430-38675953A85B}  
TERM.....CONSTANT

DATA VALUE: CONSTIPATION

-----  
ABBREVIATION.....CONSTIPAT.  
ACTIVE.....Yes  
DESCRIPTION.....Constipation  
DISPLAY\_NAME.....Constipation  
HELP\_TEXT.....No Help Text  
ID.....{4A3B0D47-88BF-4FEB-AA0D-7AD7FED2C5C3}  
TERM.....CONSTIPATION

DATA VALUE: CONSTRICTED

-----  
ABBREVIATION.....Constricte  
ACTIVE.....Yes  
DESCRIPTION.....Constricted  
DISPLAY\_NAME.....Constricted  
HELP\_TEXT.....No Help Text  
ID.....{A7F66C5B-89EF-DEEE-A6C0-0DCDFEB3AEB4}  
TERM.....CONSTRICTED

DATA VALUE: CONTACT

-----  
ABBREVIATION.....CONTACT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Contact  
HELP\_TEXT.....No Help Text  
ID.....{958F2B5E-F980-47B6-A7F7-EE7D5520B61F}  
TERM.....CONTACT

DATA VALUE: CONTACTS

-----  
ABBREVIATION.....Contacts  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Contacts  
HELP\_TEXT.....  
ID.....{6BB5B9A8-6AB0-4EF3-B8B6-C4243AD1B806}  
TERM.....CONTACTS

DATA VALUE: CONTEMPLATIVE

-----  
ABBREVIATION.....CONTEMPLAT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Contemplative  
HELP\_TEXT.....No Help Text  
ID.....{5194DB04-7A68-491E-B117-BD089C912DFB}  
TERM.....CONTEMPLATIVE

DATA VALUE: CONTINENT

-----  
ABBREVIATION.....CONTINENT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Continent  
HELP\_TEXT.....No Help Text  
ID.....{3B658506-F42E-BAF5-04C4-5AD49699E0D2}  
TERM.....CONTINENT

DATA VALUE: CONTINUOUS

-----  
ABBREVIATION.....CONTINUOUS  
ACTIVE.....Yes  
DESCRIPTION.....Continuous  
DISPLAY\_NAME.....Continuous  
HELP\_TEXT.....No Help Text  
ID.....{676CD2DA-F882-4BE5-99BB-510EA7669749}  
TERM.....CONTINUOUS

DATA VALUE: CONTINUOUS POSITIVE AIR PRESSURE

-----  
ABBREVIATION.....CPAP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....CPAP  
HELP\_TEXT.....No Help Text  
ID.....{387B42F0-F587-AAFC-ECA0-66067F18BEC6}  
TERM.....CONTINUOUS POSITIVE AIR PRESSURE

DATA VALUE: CONTINUOUS SUCTION

-----  
ABBREVIATION.....CONT SUCTI  
ACTIVE.....Yes  
DESCRIPTION.....Continuous Suction  
DISPLAY\_NAME.....Continuous Suction  
HELP\_TEXT.....No Help Text  
ID.....{8C10DAC2-9AE5-40B9-99C7-25D47F4F13BA}  
TERM.....CONTINUOUS SUCTION

DATA VALUE: CONTROL MODE VENTILATION

-----  
ABBREVIATION.....CMV  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Control Mode Ventilation  
HELP\_TEXT.....No Help Text  
ID.....{2E79351C-BA31-FC04-B636-F39728B23EA6}  
TERM.....CONTROL MODE VENTILATION

DATA VALUE: CONTUSION

-----  
ABBREVIATION.....Contusion  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Contusion  
HELP\_TEXT.....No Help Text  
ID.....{5285197A-CC0A-4D12-A1F8-966765C687F1}  
TERM.....CONTUSION

DATA VALUE: COOL

-----  
ABBREVIATION.....COOL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Cool  
HELP\_TEXT.....No Help Text  
ID.....{48A1FEA8-FD75-DAF8-174A-9C6AB2044439}  
TERM.....COOL

DATA VALUE: COOL CHEM PACK

-----  
ABBREVIATION.....CO CHEM PA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Cool Chem Pack  
HELP\_TEXT.....No Help Text  
ID.....{6DE5E644-28C2-4190-A2EF-1A766E014A65}  
TERM.....COOL CHEM PACK

DATA VALUE: COOLER THAN GENERAL BODY TEMPERATURE

-----  
ABBREVIATION.....CoolGenT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Cooler Than General Body Temperature  
HELP\_TEXT.....No Help Text  
ID.....{03D7D4CF-9B3C-49CA-9650-213E9BABE064}  
TERM.....COOLER THAN GENERAL BODY TEMPERATURE

DATA VALUE: COOLER THAN OPPOSITE BODY PART

-----  
ABBREVIATION.....COOL OPP P  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Cooler Than Opposite Body Part  
HELP\_TEXT.....No Help Text  
ID.....{AFA2635C-8E80-4C71-A1B4-4AB117096FD2}  
TERM.....COOLER THAN OPPOSITE BODY PART

DATA VALUE: COOLING BLANKET

-----  
ABBREVIATION.....COOL BLANK  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Cooling Blanket  
HELP\_TEXT.....No Help Text  
ID.....{13A42B12-4959-4674-AF97-2D1242CFD934}  
TERM.....COOLING BLANKET



DATA VALUE: COOLING OFF

-----  
ABBREVIATION.....COOLING OF  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Cooling Off  
HELP\_TEXT.....No Help Text  
ID.....{14F573ED-0D10-4E0F-A1F3-A11F29D6DB11}  
TERM.....COOLING OFF

DATA VALUE: COOLING ON

-----  
ABBREVIATION.....COOLING ON  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Cooling On  
HELP\_TEXT.....No Help Text  
ID.....{911CA17D-614E-445A-969E-4F5220E469EB}  
TERM.....COOLING ON

DATA VALUE: COOPERATIVE

-----  
ABBREVIATION.....c  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Cooperative  
HELP\_TEXT.....No Help Text  
ID.....{C176406E-EBDE-4ADE-53DE-CF191D8AFC85}  
TERM.....COOPERATIVE

DATA VALUE: COPIUS

-----  
ABBREVIATION.....COPIUS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Copius  
HELP\_TEXT.....No Help Text  
ID.....{C74A6979-0E35-AA12-F992-69ACA3B0FA23}  
TERM.....COPIUS

DATA VALUE: CRACKLES

-----  
ABBREVIATION.....c  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Crackles  
HELP\_TEXT.....No Help Text  
ID.....{6B64C108-4424-00D9-525F-D89D549476D6}  
TERM.....CRACKLES

DATA VALUE: CRAMPING

-----  
ABBREVIATION.....CRAMPING  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Cramping  
HELP\_TEXT.....No Help Text  
ID.....{3EC3E352-5469-4A85-80E4-3AFEB61224C4}  
TERM.....CRAMPING

DATA VALUE: CREPITUS

-----  
ABBREVIATION.....CREPITUS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Crepitus  
HELP\_TEXT.....No Help Text  
ID.....{56A97506-05F0-4F37-8FD9-3AA82191408A}  
TERM.....CREPITUS

DATA VALUE: CRNA

-----  
ABBREVIATION.....CRNA  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....CRNA  
HELP\_TEXT.....  
ID.....{41D71B84-D87B-4219-8E1F-AF0FF5CF5EDF}  
TERM.....CRNA

DATA VALUE: CRUSHING

-----  
ABBREVIATION.....CRUSHING  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Crushing  
HELP\_TEXT.....No Help Text  
ID.....{92D4BAF1-63E3-4FC0-B096-B5CBC4736245}  
TERM.....CRUSHING

DATA VALUE: CRUTCH

-----  
ABBREVIATION.....CRUTCH  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Crutch  
HELP\_TEXT.....No Help Text  
ID.....{2F040C2C-7A95-4FFA-8DAC-8B2EA9A90259}  
TERM.....CRUTCH

DATA VALUE: CRUTCHES

-----  
ABBREVIATION.....CRUTCHES  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Crutches  
HELP\_TEXT.....No Help Text  
ID.....{DDA81269-B25C-4F7D-A059-BB1602B2B179}  
TERM.....CRUTCHES

DATA VALUE: CRYING

-----  
ABBREVIATION.....CRYING  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Crying  
HELP\_TEXT.....No Help Text  
ID.....{D87AB50C-0B18-40DD-884C-6AAEE4DD97FC}  
TERM.....CRYING

DATA VALUE: CRYSTALLOID

-----  
ABBREVIATION.....CRYSTALLOI  
ACTIVE.....Yes  
DESCRIPTION.....Crystalloid  
DISPLAY\_NAME.....Crystalloid  
HELP\_TEXT.....No Help Text  
ID.....{4EC07E19-88B2-4C13-99B6-8233DF612997}  
TERM.....CRYSTALLOID

DATA VALUE: CSM: CENTRAL STEADY MAINTAINED

-----  
ABBREVIATION.....CSM  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....CSM: Central Steady Maintained  
HELP\_TEXT.....  
ID.....{43085BCB-4264-4E4A-BC17-A9C252A89BA8}  
TERM.....CSM: CENTRAL STEADY MAINTAINED

DATA VALUE: CULTURAL

-----  
ABBREVIATION.....CULTURAL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Cultural  
HELP\_TEXT.....No Help Text  
ID.....{530D166D-D7C1-4204-94DC-83F4E140F11C}  
TERM.....CULTURAL

DATA VALUE: CYANOTIC

-----  
ABBREVIATION.....c  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Cyanotic  
HELP\_TEXT.....No Help Text  
ID.....{06BD77DB-1849-D731-0F27-F8CE48755EA9}  
TERM.....CYANOTIC

DATA VALUE: D - DOPPLER

-----  
ABBREVIATION.....d-d  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....D - Doppler  
HELP\_TEXT.....No Help Text  
ID.....{2B77F099-0520-474D-A4BC-22BF824C8364}  
TERM.....D - DOPPLER

DATA VALUE: D10W ENTERAL NUTRITION

-----  
ABBREVIATION.....D10W ENT N  
ACTIVE.....Yes  
DESCRIPTION.....D10W Enteral Nutrition Type  
DISPLAY\_NAME.....D10W Enteral Nutrition  
HELP\_TEXT.....No Help Text  
ID.....{5049BCEF-CAE5-490D-9EB5-4A204A4D85E1}  
TERM.....D10W ENTERAL NUTRITION

DATA VALUE: D50W ENTERAL NUTRITION

-----  
ABBREVIATION.....D50W ENT N  
ACTIVE.....Yes  
DESCRIPTION.....D50W Enteral Nutrition Type  
DISPLAY\_NAME.....D50W Enteral Nutrition  
HELP\_TEXT.....No Help Text  
ID.....{528A4339-0261-4EC5-821C-7CE5592666DC}  
TERM.....D50W ENTERAL NUTRITION

DATA VALUE: D5W ENTERAL NUTRITION

-----  
ABBREVIATION.....D5W  
ACTIVE.....Yes  
DESCRIPTION.....D5W Enteral Nutrition Type  
DISPLAY\_NAME.....D5W Enteral Nutrition  
HELP\_TEXT.....No Help Text  
ID.....{D6664641-BF43-4157-9A30-8AD8CFE96518}  
TERM.....D5W ENTERAL NUTRITION

DATA VALUE: D=DOPPLER

-----  
ABBREVIATION.....d  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....D=Doppler  
HELP\_TEXT.....No Help Text  
ID.....{BCF37855-1E6E-9544-B316-D80A651697E4}  
TERM.....D=DOPPLER

DATA VALUE: DAILY PLAN

-----  
ABBREVIATION.....DAILY PLAN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Daily Plan  
HELP\_TEXT.....No Help Text  
ID.....{B29357A5-4294-4EA0-A0D3-59182F3C94FB}  
TERM.....DAILY PLAN

DATA VALUE: DAMPENED

-----  
ABBREVIATION.....DAMPENED  
ACTIVE.....Yes  
DESCRIPTION.....Dampened  
DISPLAY\_NAME.....Dampened  
HELP\_TEXT.....No Help Text  
ID.....{77C777B3-7E56-41E3-8E70-DE144E6A3424}  
TERM.....DAMPENED

DATA VALUE: DANGLED

-----  
ABBREVIATION.....DANGLED  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Dangled  
HELP\_TEXT.....No Help Text  
ID.....{40665750-112F-401D-AEFA-8E19F45EC0C9}  
TERM.....DANGLED

DATA VALUE: DATE OF BIRTH

-----  
ABBREVIATION.....DOB  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Date of Birth  
HELP\_TEXT.....  
ID.....{F5EA7F14-946E-4B4A-A7AC-613D0CBC3D6F}  
TERM.....DATE OF BIRTH

DATA VALUE: DAY ROOM

-----  
ABBREVIATION.....dr  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Day Room  
HELP\_TEXT.....No Help Text  
ID.....{269125EF-FA82-B17A-6B56-0CEAB4788899}  
TERM.....DAY ROOM

DATA VALUE: DDD

-----  
ABBREVIATION.....DDD  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....DDD  
HELP\_TEXT.....No Help Text  
ID.....{7AA7152E-48F0-8550-AD0D-9F298E90B765}  
TERM.....DDD

DATA VALUE: DDDR

-----  
ABBREVIATION.....DDDR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....DDDR  
HELP\_TEXT.....No Help Text  
ID.....{A3CC0D9D-D329-1A52-2DD1-911A1E04B871}  
TERM.....DDDR

DATA VALUE: DDI

-----  
ABBREVIATION.....DDI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....DDI  
HELP\_TEXT.....No Help Text  
ID.....{0C1E359F-B1F4-7994-D337-37622CB14BC8}  
TERM.....DDI

DATA VALUE: DECELERATING

-----  
ABBREVIATION.....DECELERATI  
ACTIVE.....Yes  
DESCRIPTION.....Decelerating  
DISPLAY\_NAME.....Decelerating  
HELP\_TEXT.....Decelerating  
ID.....{B259D4D4-82B8-454A-B53D-F12F4952BB61}  
TERM.....DECELERATING

DATA VALUE: DECREASED

-----  
ABBREVIATION.....DECREASED  
ACTIVE.....Yes  
DESCRIPTION.....Decreased  
DISPLAY\_NAME.....Decreased  
HELP\_TEXT.....No Help Text  
ID.....{D910FE54-A6AC-4BE3-AE35-C9E6679A93BD}  
TERM.....DECREASED

DATA VALUE: DEEP

-----  
ABBREVIATION.....DEEP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Deep  
HELP\_TEXT.....No Help Text  
ID.....{78AD569C-16D4-CFC3-6B6F-41FB9FFAB9D7}  
TERM.....DEEP

DATA VALUE: DENIES HAVING SLEPT

-----  
ABBREVIATION.....DENIES SL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Denies Having Slept  
HELP\_TEXT.....No Help Text  
ID.....{1AC4B0CF-EC41-426B-89AB-68D5A9BAFADA}  
TERM.....DENIES HAVING SLEPT

DATA VALUE: DENIES NEED

-----  
ABBREVIATION.....DENIES NE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Denies Need  
HELP\_TEXT.....No Help Text  
ID.....{28C8C733-5A28-4DDC-8D2E-F9D0E6CD9B28}  
TERM.....DENIES NEED

DATA VALUE: DENTAL

-----  
ABBREVIATION.....Dental  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Dental  
HELP\_TEXT.....  
ID.....{0A60F06A-ED63-4B06-8FF9-EEFB12567D65}  
TERM.....DENTAL

DATA VALUE: DENTAL IMPLANTS

-----  
ABBREVIATION.....Implants  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Dental Implants  
HELP\_TEXT.....  
ID.....{D3AF96DE-6920-4C93-8105-99EE574FC092}  
TERM.....DENTAL IMPLANTS

DATA VALUE: DENTIST

-----  
ABBREVIATION.....Dentist  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Dentist  
HELP\_TEXT.....  
ID.....{DE4A47EB-35D3-4992-B8DA-3E5D51E0EA3C}  
TERM.....DENTIST

DATA VALUE: DENTURES

-----  
ABBREVIATION.....Dentures  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Dentures  
HELP\_TEXT.....  
ID.....{ECABE4D0-97BE-418D-AC7D-29E40D166581}  
TERM.....DENTURES

DATA VALUE: DESCRIP

-----  
ABBREVIATION.....N  
ACTIVE.....Yes  
DESCRIPTION.....1.5 Cal Enteral Nutrition Formula  
DISPLAY\_NAME.....N  
HELP\_TEXT.....No Help Text  
ID.....{47FD2D75-2C45-476F-82AF-7B40108E9C07}  
TERM.....DESCRIP

DATA VALUE: DEVICE RELATED PRESSURE

-----  
ABBREVIATION.....DEV REL PR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Device Related Pressure  
HELP\_TEXT.....No Help Text  
ID.....{3AED6107-A173-4B53-8EB5-2FB2B3336FC5}  
TERM.....DEVICE RELATED PRESSURE



DATA VALUE: DEVICE RELEASE

-----  
ABBREVIATION.....drq  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Device Release  
HELP\_TEXT.....No Help Text  
ID.....{9F97947F-CFFE-3285-F4C3-A8DF294DA46C}  
TERM.....DEVICE RELEASE

DATA VALUE: DIABETES

-----  
ABBREVIATION.....Diabetes  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Diabetes  
HELP\_TEXT.....  
ID.....{E295E970-759D-44CD-AD43-6047382A790D}  
TERM.....DIABETES

DATA VALUE: DIABETIC/VASCULAR ULCER

-----  
ABBREVIATION.....DIA/VAS UL  
ACTIVE.....Yes  
DESCRIPTION.....Type of wound  
DISPLAY\_NAME.....Diabetic / Vascular Ulcer  
HELP\_TEXT.....No Help Text  
ID.....{E209E53E-EFF1-48F5-9581-75236C89D9B9}  
TERM.....DIABETIC/VASCULAR ULCER

DATA VALUE: DIAGNOSIS/HEALTH STATUS

-----  
ABBREVIATION.....Diag/Healt  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Diagnosis/Health Status  
HELP\_TEXT.....  
ID.....{42F128CE-63B7-4AA9-A121-B08C507F9941}  
TERM.....DIAGNOSIS/HEALTH STATUS

DATA VALUE: DIALYSIS

-----  
ABBREVIATION.....DIALYSIS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Dialysis Line  
HELP\_TEXT.....No Help Text  
ID.....{D218E23C-5A34-09B0-DE43-DC792DCD0F1C}  
TERM.....DIALYSIS

DATA VALUE: DIAPHORETIC

-----  
ABBREVIATION.....d  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Diaphoretic  
HELP\_TEXT.....No Help Text  
ID.....{F5F1E644-A1BD-0137-54B5-1C02678C1D5A}  
TERM.....DIAPHORETIC

DATA VALUE: DIARRHEA

-----  
ABBREVIATION.....DIARRHEA  
ACTIVE.....Yes  
DESCRIPTION.....Diarrhea  
DISPLAY\_NAME.....Diarrhea  
HELP\_TEXT.....No Help Text  
ID.....{353472EF-3499-46B1-AD7C-554FE01D8CE2}  
TERM.....DIARRHEA

DATA VALUE: DIASTOLIC

-----  
ABBREVIATION.....DIASTOLIC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Diastolic  
HELP\_TEXT.....No Help Text  
ID.....{3E9D1DCF-6B3B-482C-AE5F-442814F4B5C0}  
TERM.....DIASTOLIC

DATA VALUE: DIFFICULT TO AROUSE

-----  
ABBREVIATION.....Diff to Ar  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Difficult to Arouse  
HELP\_TEXT.....  
ID.....{5EF640D1-CE7F-49AB-80DC-6BCFF47BC86C}  
TERM.....DIFFICULT TO AROUSE

DATA VALUE: DILATED

-----  
ABBREVIATION.....Dilated  
ACTIVE.....Yes  
DESCRIPTION.....Dilated  
DISPLAY\_NAME.....Dilated  
HELP\_TEXT.....No Help Text  
ID.....{BF98ADC4-C69E-B17A-9BDE-D426D7F4C395}  
TERM.....DILATED

DATA VALUE: DILUTE DAKINS SOLUTION

-----  
ABBREVIATION.....DL DAKINS  
ACTIVE.....Yes  
DESCRIPTION.....Dilute Dakins Solution  
DISPLAY\_NAME.....Dilute Dakins Solution  
HELP\_TEXT.....No Help Text  
ID.....{66DEF4D0-54C0-4ADF-8FFD-F4E022C30F43}  
TERM.....DILUTE DAKINS SOLUTION

DATA VALUE: DILUTE SODIUM HYPOCHLORITE SOLUTION

-----  
ABBREVIATION.....D NA HYPO  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Dilute Sodium Hypochlorite Solution  
HELP\_TEXT.....No Help Text  
ID.....{14D2B1FB-1A5D-42AA-ACCF-A5033741DB61}  
TERM.....DILUTE SODIUM HYPOCHLORITE SOLUTION

DATA VALUE: DIMINISHED

-----  
ABBREVIATION.....d  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Diminished  
HELP\_TEXT.....No Help Text  
ID.....{49DE419E-DC3C-F7CD-771D-742031E0DEBC}  
TERM.....DIMINISHED

DATA VALUE: DISCHARGE

-----  
ABBREVIATION.....DISCHARGE  
ACTIVE.....Yes  
DESCRIPTION.....Discharge  
DISPLAY\_NAME.....Discharge  
HELP\_TEXT.....No Help Text  
ID.....{ADF67CC2-E05B-4B8D-A233-4AE1C35B2D20}  
TERM.....DISCHARGE

DATA VALUE: DISCHARGE MEDS PLANNED

-----  
ABBREVIATION.....DI ME PL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Discharge Meds Planned  
HELP\_TEXT.....No Help Text  
ID.....{04E9DF8B-8807-47A5-AE88-180F431F3BDE}  
TERM.....DISCHARGE MEDS PLANNED

DATA VALUE: DISCHARGE ORDER ON CHART

-----  
ABBREVIATION.....DI OR CH  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Discharge Order On Chart  
HELP\_TEXT.....No Help Text  
ID.....{F7267E73-07BC-44FB-AB0F-FAA07F6E2C1B}  
TERM.....DISCHARGE ORDER ON CHART

DATA VALUE: DISCOLORATION

-----  
ABBREVIATION.....DISCOLORAT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Discoloration  
HELP\_TEXT.....No Help Text  
ID.....{BBE3E910-CC49-430C-BD90-FB0A97F1B183}  
TERM.....DISCOLORATION

DATA VALUE: DISCOMFORT

-----  
ABBREVIATION.....DISCOMFORT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Discomfort  
HELP\_TEXT.....No Help Text  
ID.....{77D058C9-1EB4-96C0-C9E3-666C840B884C}  
TERM.....DISCOMFORT

DATA VALUE: DISCONTINUED

-----  
ABBREVIATION.....DISCONTINU  
ACTIVE.....Yes  
DESCRIPTION.....Discontinued  
DISPLAY\_NAME.....Discontinued  
HELP\_TEXT.....No Help Text  
ID.....{15E289FB-09DF-4731-9CF3-48A4E0031C31}  
TERM.....DISCONTINUED

DATA VALUE: DISEASE PROCESS

-----  
ABBREVIATION.....DISEA PROC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Disease Process  
HELP\_TEXT.....No Help Text  
ID.....{E80807AE-6208-421D-BE4C-445C470D32FA}  
TERM.....DISEASE PROCESS

DATA VALUE: DISHEVELED

-----  
ABBREVIATION.....DISHEVEL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Disheveled  
HELP\_TEXT.....No Help Text  
ID.....{0455549E-5CCB-470D-9D7F-F3959E6D55FE}  
TERM.....DISHEVELED

DATA VALUE: DISORIENTED

-----  
ABBREVIATION.....DISORIENTE  
ACTIVE.....Yes  
DESCRIPTION.....Disoriented  
DISPLAY\_NAME.....Disoriented  
HELP\_TEXT.....No Help Text  
ID.....{A7D1F2EA-CD7C-4100-8552-44834684C6DE}  
TERM.....DISORIENTED

DATA VALUE: DISORIENTED TO PERSON

-----  
ABBREVIATION.....DISOR TO P  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Disoriented To Person  
HELP\_TEXT.....No Help Text  
ID.....{E5D83C03-7EF5-770B-930F-7B3B3345048C}  
TERM.....DISORIENTED TO PERSON

DATA VALUE: DISORIENTED TO PLACE

-----  
ABBREVIATION.....dtp  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Disoriented To Place  
HELP\_TEXT.....No Help Text  
ID.....{00C98809-C0D0-C7B4-0615-8F7D31017FE8}  
TERM.....DISORIENTED TO PLACE

DATA VALUE: DISORIENTED TO SITUATION

-----  
ABBREVIATION.....dts  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Disoriented To Situation  
HELP\_TEXT.....No Help Text  
ID.....{7143B3D7-5569-4378-D3AA-10C4675F0BA5}  
TERM.....DISORIENTED TO SITUATION

DATA VALUE: DISORIENTED TO TIME

-----  
ABBREVIATION.....dt  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Disoriented To Time  
HELP\_TEXT.....No Help Text  
ID.....{2AD1A28A-4121-4B82-CFB8-E5A8D3095704}  
TERM.....DISORIENTED TO TIME

DATA VALUE: DISTAL PORT

-----  
ABBREVIATION.....DIST PORT  
ACTIVE.....Yes  
DESCRIPTION.....Distal Port  
DISPLAY\_NAME.....Distal Port  
HELP\_TEXT.....No Help Text  
ID.....{E493DF34-AD0E-4951-9715-5AD4EFA14036}  
TERM.....DISTAL PORT

DATA VALUE: DISTENDED

-----  
ABBREVIATION.....DISTENDED  
ACTIVE.....Yes  
DESCRIPTION.....Distended  
DISPLAY\_NAME.....Distended  
HELP\_TEXT.....No Help Text  
ID.....{B5C529CA-2BE0-4698-A153-718F4D58E643}  
TERM.....DISTENDED

DATA VALUE: DIVERSION

-----  
ABBREVIATION.....DIVERSION  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Diversion  
HELP\_TEXT.....No Help Text  
ID.....{087AAE9F-7FE1-5717-6423-07574C637875}  
TERM.....DIVERSION

DATA VALUE: DO

-----  
ABBREVIATION.....DO  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....DO  
HELP\_TEXT.....  
ID.....{F70738A6-7BFA-4D76-881C-6A4B4F4E41E4}  
TERM.....DO

DATA VALUE: DO NOT RESUSCITATE

-----  
ABBREVIATION.....dnr  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Do Not Resuscitate  
HELP\_TEXT.....No Help Text  
ID.....{0BD8A134-C563-7E00-A6B1-F7C543CEF6FF}  
TERM.....DO NOT RESUSCITATE

DATA VALUE: DONE

-----  
ABBREVIATION.....DONE  
ACTIVE.....Yes  
DESCRIPTION.....Done  
DISPLAY\_NAME.....Done  
HELP\_TEXT.....No Help Text  
ID.....{60412C8D-132F-46EA-8DDF-FCC282C4CC8B}  
TERM.....DONE

DATA VALUE: DOO

-----  
ABBREVIATION.....DOO  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....DOO  
HELP\_TEXT.....No Help Text  
ID.....{749BCE3D-FF09-B796-FF2A-ECD2D116C96F}  
TERM.....DOO

DATA VALUE: DOUBLE

-----  
ABBREVIATION.....DOUBLE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Double  
HELP\_TEXT.....No Help Text  
ID.....{A486852F-150B-24E3-8EB3-8CAABB3D0520}  
TERM.....DOUBLE

DATA VALUE: DOWN

-----  
ABBREVIATION.....d  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Down  
HELP\_TEXT.....No Help Text  
ID.....{36ECC86E-2D73-D4A9-5E13-9972E2ECF2E2}  
TERM.....DOWN

DATA VALUE: DRAINAGE FLOW IS UNOBSTRUCTED

-----  
ABBREVIATION.....Unobstruct  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Drainage Flow is Unobstructed  
HELP\_TEXT.....  
ID.....{DB66B2DD-50FE-468B-8892-1A8FF157BCB0}  
TERM.....DRAINAGE FLOW IS UNOBSTRUCTED

DATA VALUE: DRESSING

-----  
ABBREVIATION.....d  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Dressing  
HELP\_TEXT.....No Help Text  
ID.....{FFEF0714-3B5B-E2B4-D760-3C03D45E59EA}  
TERM.....DRESSING

DATA VALUE: DRESSING CHANGED

-----  
ABBREVIATION.....DRESS CHAN  
ACTIVE.....Yes  
DESCRIPTION.....Changed Dressing  
DISPLAY\_NAME.....Dressing Changed  
HELP\_TEXT.....No Help Text  
ID.....{E8CF4B7A-D1AB-49A2-8E7C-3375C3BC436F}  
TERM.....DRESSING CHANGED

DATA VALUE: DRESSING WITH PRESCRIPTION RX

-----  
ABBREVIATION.....DR W RX  
ACTIVE.....Yes  
DESCRIPTION.....Dressing With Prescription Rx  
DISPLAY\_NAME.....Dressing With Prescription Rx  
HELP\_TEXT.....No Help Text  
ID.....{AD5F9FDC-F280-4F1C-99E9-2BE0F9FFB15C}  
TERM.....DRESSING WITH PRESCRIPTION RX

DATA VALUE: DRESSINGS CHECKED

-----  
ABBREVIATION.....DRESS CH  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Dressings Checked  
HELP\_TEXT.....No Help Text  
ID.....{21014F53-29CE-4914-8760-6976D8DBA742}  
TERM.....DRESSINGS CHECKED



DATA VALUE: DROPLET

-----  
ABBREVIATION.....DROPLET  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Droplet  
HELP\_TEXT.....No Help Text  
ID.....{CE1DD314-8A30-4785-BD56-1540CDA85146}  
TERM.....DROPLET

DATA VALUE: DROWSY

-----  
ABBREVIATION.....d  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Drowsy  
HELP\_TEXT.....No Help Text  
ID.....{6B350814-A88A-0F05-C382-9F27830FBC09}  
TERM.....DROWSY

DATA VALUE: DRY

-----  
ABBREVIATION.....d  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Dry  
HELP\_TEXT.....No Help Text  
ID.....{D8F4F0FD-1B19-2515-3815-EC88D175B305}  
TERM.....DRY

DATA VALUE: DRY AND INTACT

-----  
ABBREVIATION.....DRY INTACT  
ACTIVE.....Yes  
DESCRIPTION.....Dry and Intact  
DISPLAY\_NAME.....Dry and Intact  
HELP\_TEXT.....No Help Text  
ID.....{77B62A18-F098-4764-86FA-66373E08E8C3}  
TERM.....DRY AND INTACT

DATA VALUE: DRY SHAMPOO

-----  
ABBREVIATION.....DRY SHAMP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Dry Shampoo  
HELP\_TEXT.....No Help Text  
ID.....{8910A186-2590-4352-A884-7D692A9D2B3C}  
TERM.....DRY SHAMPOO

DATA VALUE: DRY STERILE DRESSING

-----  
ABBREVIATION.....DRY ST DRE  
ACTIVE.....Yes  
DESCRIPTION.....Dry Sterile Dressing  
DISPLAY\_NAME.....Dry Sterile Dressing  
HELP\_TEXT.....No Help Text  
ID.....{3069DFCE-1B59-467C-97CE-9DBAF9AF42CC}  
TERM.....DRY STERILE DRESSING

DATA VALUE: DULL

-----  
ABBREVIATION.....DULL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Dull  
HELP\_TEXT.....No Help Text  
ID.....{F0F80976-F6A0-4363-9C75-AAAE9937DB49}  
TERM.....DULL

DATA VALUE: DUMPING SYNDROME

-----  
ABBREVIATION.....DUMP SYND  
ACTIVE.....Yes  
DESCRIPTION.....Dumping Syndrome Diet  
DISPLAY\_NAME.....Dumping Syndrome Diet  
HELP\_TEXT.....No Help Text  
ID.....{0173C971-63CB-45A5-9E93-D966EB3B957D}  
TERM.....DUMPING SYNDROME

DATA VALUE: DUODENOSTOMY

-----  
ABBREVIATION.....DUODENOSTO  
ACTIVE.....Yes  
DESCRIPTION.....Duodenostomy  
DISPLAY\_NAME.....Duodenostomy  
HELP\_TEXT.....No Help Text  
ID.....{FA44135D-6806-42E7-B371-8C3912B3582E}  
TERM.....DUODENOSTOMY

DATA VALUE: DUSKY

-----  
ABBREVIATION.....d  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Dusky  
HELP\_TEXT.....No Help Text  
ID.....{2ABF845E-4E2F-883F-3FCC-F72B8507FF85}  
TERM.....DUSKY

DATA VALUE: DVI

-----  
ABBREVIATION.....DVI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....DVI  
HELP\_TEXT.....No Help Text  
ID.....{4D6E3081-7BB5-45DB-E8E6-0D3D7869E6BD}  
TERM.....DVI

DATA VALUE: DYE PRESENT

-----  
ABBREVIATION.....DYE PRESEN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Dye Present  
HELP\_TEXT.....No Help Text  
ID.....{E865AD3E-6DAD-C3F9-6F93-7AD6F71234A4}  
TERM.....DYE PRESENT

DATA VALUE: DYSPHAGIA PHASE I

-----  
ABBREVIATION.....DYSP PH I  
ACTIVE.....Yes  
DESCRIPTION.....Dysphagia Stage I Diet  
DISPLAY\_NAME.....Dysphagia Stage I Diet  
HELP\_TEXT.....No Help Text  
ID.....{8011795B-77A5-4FF4-8EB9-EB6E1AA1FD4D}  
TERM.....DYSPHAGIA PHASE I

DATA VALUE: DYSPHAGIA STAGE II

-----  
ABBREVIATION.....DYSP ST II  
ACTIVE.....Yes  
DESCRIPTION.....Dysphagia Stage II Diet  
DISPLAY\_NAME.....Dysphagia Stage II Diet  
HELP\_TEXT.....No Help Text  
ID.....{294BF58A-47FD-43A1-8A1E-9E348987B004}  
TERM.....DYSPHAGIA STAGE II

DATA VALUE: DYSPHAGIA STAGE III

-----  
ABBREVIATION.....DYS ST III  
ACTIVE.....Yes  
DESCRIPTION.....Dysphagia Stage III Diet  
DISPLAY\_NAME.....Dysphagia Stage III Diet  
HELP\_TEXT.....No Help Text  
ID.....{DCF91830-C5CA-47B2-A7C6-4DD971F953EA}  
TERM.....DYSPHAGIA STAGE III

DATA VALUE: DYSPHAGIA STAGE IV

-----  
ABBREVIATION.....DYSP ST IV  
ACTIVE.....Yes  
DESCRIPTION.....Dysphagia Stage IV Diet  
DISPLAY\_NAME.....Dysphagia Stage IV Diet  
HELP\_TEXT.....No Help Text  
ID.....{E738569A-5104-4C93-8EFC-47D80765F739}  
TERM.....DYSPHAGIA STAGE IV

DATA VALUE: DYSPHAGIA STAGE V

-----  
ABBREVIATION.....DYSPH ST V  
ACTIVE.....Yes  
DESCRIPTION.....Dysphagia Stage V Diet  
DISPLAY\_NAME.....Dysphagia Stage V Diet  
HELP\_TEXT.....No Help Text  
ID.....{ABE7D85D-1821-4592-B4E8-84670EE484E4}  
TERM.....DYSPHAGIA STAGE V

DATA VALUE: DYSPNEIC

-----  
ABBREVIATION.....DYSPENIC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Dyspenic  
HELP\_TEXT.....No Help Text  
ID.....{4C6B42D9-694F-5EB2-6396-7ED85D5B154E}  
TERM.....DYSPNEIC

DATA VALUE: DYSURIA

-----  
ABBREVIATION.....DYSURIA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Dysuria  
HELP\_TEXT.....No Help Text  
ID.....{75AC5ECD-8975-406A-30C7-E2EB23F10F84}  
TERM.....DYSURIA

DATA VALUE: EAGER TO LEARN

-----  
ABBREVIATION.....EAGER  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Eager To Learn  
HELP\_TEXT.....No Help Text  
ID.....{A20AF51B-C89C-4308-AE17-C66144FE68BE}  
TERM.....EAGER TO LEARN

DATA VALUE: EASY TO AROUSE

-----  
ABBREVIATION.....Easy to Ar  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Easy to Arouse  
HELP\_TEXT.....  
ID.....{19BB5A04-F979-449C-AA93-D5A21BD368F6}  
TERM.....EASY TO AROUSE

DATA VALUE: ECCHYMOSIS

-----  
ABBREVIATION.....ECCHYMOSIS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Ecchymosis  
HELP\_TEXT.....No Help Text  
ID.....{DF076906-CC17-6AD1-76E7-B07678A126DE}  
TERM.....ECCHYMOSIS

DATA VALUE: ELASTIC WRAP

-----  
ABBREVIATION.....ELAST WRAP  
ACTIVE.....Yes  
DESCRIPTION.....ELASTIC WRAP  
DISPLAY\_NAME.....Elastic Wrap  
HELP\_TEXT.....ELASTIC WRAP  
ID.....{9EB52FA5-32D8-823F-2E3C-83375DA62342}  
TERM.....ELASTIC WRAP

DATA VALUE: ELBOW PAD

-----  
ABBREVIATION.....ELBOW PAD  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Elbow Pad  
HELP\_TEXT.....No Help Text  
ID.....{26942F81-1E33-4962-B772-A96FFBCED814}  
TERM.....ELBOW PAD

DATA VALUE: ELECTRIC

-----  
ABBREVIATION.....ELECTRIC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....ELECTRIC  
HELP\_TEXT.....No Help Text  
ID.....{0C221385-61AF-4660-A129-35E610EA71AB}  
TERM.....ELECTRIC

DATA VALUE: ELEMENTAL FORMULA

-----  
ABBREVIATION.....ELE FORM  
ACTIVE.....Yes  
DESCRIPTION.....Elemental Enteral Nutrition Formula  
DISPLAY\_NAME.....Elemental Enteral Nutrition Formula  
HELP\_TEXT.....No Help Text  
ID.....{DCE1E81C-1AC7-48F4-97AA-F1FC4D0601D8}  
TERM.....ELEMENTAL FORMULA

DATA VALUE: ELIMINATED ENVIRONMENTAL HAZARDS

-----  
ABBREVIATION.....Eli Env Hz  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Eliminated Environmental Hazards  
HELP\_TEXT.....  
ID.....{78E77129-F53F-491D-8C86-36AE948FCCC6}  
TERM.....ELIMINATED ENVIRONMENTAL HAZARDS

DATA VALUE: ELOPEMENT

-----  
ABBREVIATION.....e  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Elopement  
HELP\_TEXT.....No Help Text  
ID.....{EB030907-FB90-CA9D-CF50-E26C5D9FE671}  
TERM.....ELOPEMENT

DATA VALUE: EMERGENCY

-----  
ABBREVIATION.....Emergency  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Emergency  
HELP\_TEXT.....  
ID.....{3C7B5E35-63CA-32EE-3D64-AE4A052F55E1}  
TERM.....EMERGENCY

DATA VALUE: EMERGENCY DEPARTMENT

-----  
ABBREVIATION.....ER DEPT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Emergency Department  
HELP\_TEXT.....No Help Text  
ID.....{D2E511F7-8B80-4928-A64B-3F1EE17395F7}  
TERM.....EMERGENCY DEPARTMENT

DATA VALUE: EMERGENCY DEPT

-----  
ABBREVIATION.....Emer Dept  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Emergency Dept  
HELP\_TEXT.....  
ID.....{E8490DF5-278C-4102-9675-C9B40D596A10}  
TERM.....EMERGENCY DEPT

DATA VALUE: EMOTIONAL

-----  
ABBREVIATION.....EMOTIONAL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Emotional  
HELP\_TEXT.....No Help Text  
ID.....{626C22E6-EF2A-489F-AA79-D27541A0BC45}  
TERM.....EMOTIONAL

DATA VALUE: END OF LIFE CARE

-----  
ABBREVIATION.....END LIFE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....End Of Life Care  
HELP\_TEXT.....No Help Text  
ID.....{047D8977-0567-4B90-AD6A-3AAF8DB35238}  
TERM.....END OF LIFE CARE

DATA VALUE: END-TIDAL CO2 MONITOR

-----  
ABBREVIATION.....End-Tidal  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....End-Tidal CO2 Monitor  
HELP\_TEXT.....  
ID.....{609E3607-E441-465D-B1FC-C6986C58C11C}  
TERM.....END-TIDAL CO2 MONITOR

DATA VALUE: ENSURE PLUS

-----  
ABBREVIATION.....Ensure Plu  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Ensure Plus  
HELP\_TEXT.....  
ID.....{519D1249-B6C0-49ED-BA86-1B27C0F964C9}  
TERM.....ENSURE PLUS

DATA VALUE: ENTERAL FEEDING FORMULA/WATER ADMINISTRATION

-----  
ABBREVIATION.....Feed Form  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Enteral Feeding Formula/Water  
Administration  
HELP\_TEXT.....  
ID.....{D2EA8840-1123-4660-AD25-587458BF1E89}  
TERM.....ENTERAL FEEDING FORMULA/WATER  
ADMINISTRATION

DATA VALUE: ENTERAL FEEDING PUMP

-----  
ABBREVIATION.....Feed Pump  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Enteral Feeding Pump  
HELP\_TEXT.....  
ID.....{31FF069C-0060-45C3-898F-20F1A9D1552B}  
TERM.....ENTERAL FEEDING PUMP

DATA VALUE: ENTERAL TUBE MEDICATION ADMINISTRATION

-----  
ABBREVIATION.....Med Admin  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Enteral Tube Medication Administration  
HELP\_TEXT.....  
ID.....{68D124A7-69CF-414B-8A47-BC666A373E24}  
TERM.....ENTERAL TUBE MEDICATION ADMINISTRATION

DATA VALUE: ENTERAL TUBE/STOMA CARE

-----  
ABBREVIATION.....Tube/Stoma  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Enteral Tube/Stoma Care  
HELP\_TEXT.....  
ID.....{B91B94EF-CFE6-4628-BEF2-D06E5870D574}  
TERM.....ENTERAL TUBE/STOMA CARE

DATA VALUE: ENVIRONMENT

-----  
ABBREVIATION.....ENVIRO  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Environment  
HELP\_TEXT.....No Help Text  
ID.....{C61AA262-F6AF-46AC-A215-0D56FD01B234}  
TERM.....ENVIRONMENT



DATA VALUE: ENZYMATIC OINTMENT

-----  
ABBREVIATION.....ENZY OINTM  
ACTIVE.....Yes  
DESCRIPTION.....Enzymatic Ointment  
DISPLAY\_NAME.....Enzymatic Ointment  
HELP\_TEXT.....No Help Text  
ID.....{E781DA75-8777-4884-A093-FB623774DB8F}  
TERM.....ENZYMATIC OINTMENT

DATA VALUE: EPICARDIAL LEADS

-----  
ABBREVIATION.....el  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Epicardial Leads  
HELP\_TEXT.....No Help Text  
ID.....{2B551575-3D50-96BD-7BAA-50CC9796361E}  
TERM.....EPICARDIAL LEADS

DATA VALUE: EPIDURAL

-----  
ABBREVIATION.....Epidural  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Epidural  
HELP\_TEXT.....  
ID.....{63317A63-16FF-4D55-A034-420C65B17AEB}  
TERM.....EPIDURAL

DATA VALUE: EPIDURAL PUMP

-----  
ABBREVIATION.....Epid Pump  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Epidural Pump  
HELP\_TEXT.....  
ID.....{26C83111-5FA8-4398-8D3E-787B31A2A9E9}  
TERM.....EPIDURAL PUMP

DATA VALUE: EPITHELIALIZATION

-----  
ABBREVIATION.....EPITHELIAL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Epithelialization  
HELP\_TEXT.....No Help Text  
ID.....{5FC43A6D-921D-4736-A234-9993DB6499B7}  
TERM.....EPITHELIALIZATION

DATA VALUE: EQUAL

-----  
ABBREVIATION.....Equal  
ACTIVE.....Yes  
DESCRIPTION.....Equal  
DISPLAY\_NAME.....Equal  
HELP\_TEXT.....No Help Text  
ID.....{AC079584-158D-EBF2-E13A-C2E919375E4D}  
TERM.....EQUAL

DATA VALUE: EQUIPMENT

-----  
ABBREVIATION.....EQUIPMENT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Equipment  
HELP\_TEXT.....No Help Text  
ID.....{04C0F805-C7B3-4A9A-91D8-0F32F3265874}  
TERM.....EQUIPMENT

DATA VALUE: ERYTHEMA

-----  
ABBREVIATION.....ERYTHEMA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Erythema  
HELP\_TEXT.....No Help Text  
ID.....{0083E1B3-AE2B-4CCF-8794-86E612B9DCEA}  
TERM.....ERYTHEMA

DATA VALUE: ESCAPE

-----  
ABBREVIATION.....e  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Escape  
HELP\_TEXT.....No Help Text  
ID.....{965F3201-9D63-66AB-FA25-99479AE7145A}  
TERM.....ESCAPE

DATA VALUE: ESCHAR

-----  
ABBREVIATION.....ESCHAR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Eschar  
HELP\_TEXT.....No Help Text  
ID.....{555E03C5-FBFE-4EB9-816F-25CE34476B3C}  
TERM.....ESCHAR

DATA VALUE: ESTABLISHED ROUTINE TOILETING SCHEDULE

-----  
ABBREVIATION.....Toil Sched  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Established Routine Toileting Schedule  
HELP\_TEXT.....  
ID.....{A66F5259-5C01-4597-B881-9267DEBC017C}  
TERM.....ESTABLISHED ROUTINE TOILETING SCHEDULE

DATA VALUE: ETT REPOSITIONED

-----  
ABBREVIATION.....ETT REPOSI  
ACTIVE.....Yes  
DESCRIPTION.....ETT Repositioned  
DISPLAY\_NAME.....ETT Repositioned  
HELP\_TEXT.....No Help Text  
ID.....{95B35EDB-76B6-433E-B526-01C543A7CCEA}  
TERM.....ETT REPOSITIONED

DATA VALUE: ETT RETAPED/SECURED

-----  
ABBREVIATION.....ETT RETAPE  
ACTIVE.....Yes  
DESCRIPTION.....ETT Retaped / Secured  
DISPLAY\_NAME.....ETT Retaped / Secured  
HELP\_TEXT.....No Help Text  
ID.....{6D8A2EDA-CAC9-44D2-BDA9-3B7F20908F1E}  
TERM.....ETT RETAPED/SECURED

DATA VALUE: EWALD TUBE

-----  
ABBREVIATION.....EWALD TUBE  
ACTIVE.....Yes  
DESCRIPTION.....Ewald Tube  
DISPLAY\_NAME.....Ewald Tube  
HELP\_TEXT.....No Help Text  
ID.....{B83AB4C1-A098-4FAB-A645-61FE1FEADBBE}  
TERM.....EWALD TUBE

DATA VALUE: EXPIRATORY WHEEZES

-----  
ABBREVIATION.....EXP WHEEZE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Expiratory Wheezes  
HELP\_TEXT.....No Help Text  
ID.....{CFC103A8-8B02-4942-BBDD-8D24C930883E}  
TERM.....EXPIRATORY WHEEZES

DATA VALUE: EXTENDED PERIOD OF IMMOBILITY

-----  
ABBREVIATION.....Ex Pd Immo  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Extended Period of Immobility  
HELP\_TEXT.....  
ID.....{3DD00E20-6B5E-47EF-B23C-ABA2C1ABAF0A}  
TERM.....EXTENDED PERIOD OF IMMOBILITY

DATA VALUE: EXTENDS AROUND INCISION EXTENDS INTO PERIWOUND LESS THAN 1 CM

-----  
ABBREVIATION.....ER EXI <1  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Less Than 1 cm  
HELP\_TEXT.....No Help-Text  
ID.....{095709D0-6190-49A3-AB08-BF4575DC8396}  
TERM.....EXTENDS AROUND INCISION EXTENDS INTO  
PERIWOUND LESS THAN 1 CM

DATA VALUE: EXTENDS AROUND INCISION MORE THAN 1 CM

-----  
ABBREVIATION.....ER EXI >1  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Extends Into Periwound More Than 1 cm  
HELP\_TEXT.....No Help-Text  
ID.....{11F13A1F-7AC2-4629-A960-CD92655EE101}  
TERM.....EXTENDS AROUND INCISION MORE THAN 1 CM

DATA VALUE: EXTENDS AROUND INCISION NONE

-----  
ABBREVIATION.....ER EXI NON  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....None  
HELP\_TEXT.....No Help-Text  
ID.....{91F1A7C5-885C-4B13-9F7F-189C93E505CC}  
TERM.....EXTENDS AROUND INCISION NONE

DATA VALUE: EXTENDS INTO PERIWOUND <1CM

-----  
ABBREVIATION.....EX IN PER<  
ACTIVE.....Yes  
DESCRIPTION.....Extends Into Periwound <1cm  
DISPLAY\_NAME.....Extends Into Periwound <1cm  
HELP\_TEXT.....No Help Text  
ID.....{0EBCEAE0-E17A-484F-BFE2-51F17AA761CB}  
TERM.....EXTENDS INTO PERIWOUND <1CM

DATA VALUE: EXTENDS INTO PERIWOUND>1CM

-----  
ABBREVIATION.....EX IN PER>  
ACTIVE.....Yes  
DESCRIPTION.....Extends Into Periwound > 1cm  
DISPLAY\_NAME.....Extends Into Periwound > 1cm  
HELP\_TEXT.....No Help Text  
ID.....{1F31FEBB-DA19-434C-8C58-67CDC5385231}  
TERM.....EXTENDS INTO PERIWOUND>1CM

DATA VALUE: EXTENDS TO PAIN

-----  
ABBREVIATION.....EXT TO PAI  
ACTIVE.....Yes  
DESCRIPTION.....Extends To Pain  
DISPLAY\_NAME.....Extends To Pain  
HELP\_TEXT.....No Help Text  
ID.....{A62B2C4D-BDB7-49FA-9669-DF6732063DA6}  
TERM.....EXTENDS TO PAIN

DATA VALUE: EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND LESS THAN 1 CM

-----  
ABBREVIATION.....ERY EX <1  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Extends Into Periwound Less Than 1 cm  
HELP\_TEXT.....No Help-Text  
ID.....{48C7FA69-D31A-4AAA-84F3-79AC863B71FD}  
TERM.....EXTENSION AROUND WOUND EXTENDS INTO  
PERIWOUND LESS THAN 1 CM

DATA VALUE: EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND MORE THAN 1 CM

-----  
ABBREVIATION.....ERY EX >1  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Extends Into Periwound Greater Than 1  
cm  
HELP\_TEXT.....No Help-Text  
ID.....{175B2505-C596-4036-8FC3-DFEA8EFB76F9}  
TERM.....EXTENSION AROUND WOUND EXTENDS INTO  
PERIWOUND MORE THAN 1 CM

DATA VALUE: EXTENSION AROUND WOUND LIMITED TO WOUND

-----  
ABBREVIATION.....ERY EX LIM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Limited To Wound  
HELP\_TEXT.....No Help-Text  
ID.....{2DC83B85-B3B6-4679-9F48-7DD26EFFC4D4}  
TERM.....EXTENSION AROUND WOUND LIMITED TO WOUND

DATA VALUE: EXTENSION AROUND WOUND NONE

-----  
ABBREVIATION.....ERY EX NON  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....None  
HELP\_TEXT.....No Help-Text  
ID.....{815C7BDF-C9A8-4AFE-B2BE-2C2FF24AA495}  
TERM.....EXTENSION AROUND WOUND NONE

DATA VALUE: EXTENSIVE ASSISTANCE

-----  
ABBREVIATION.....EXT ASSIST  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Extensive Assistance  
HELP\_TEXT.....No Help Text  
ID.....{C402B659-E880-4070-8345-13E5508B4D29}  
TERM.....EXTENSIVE ASSISTANCE

DATA VALUE: EXTERNAL CATHETER

-----  
ABBREVIATION.....EXT. CATHE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....External Catheter  
HELP\_TEXT.....No Help Text  
ID.....{C8AEE936-DD91-FE91-4DD6-6F210FB59726}  
TERM.....EXTERNAL CATHETER

DATA VALUE: EXTERNAL CONDOM CATH

-----  
ABBREVIATION.....Ex Con Cat  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....External Condom Cath  
HELP\_TEXT.....  
ID.....{747E493A-C535-44EA-AD6D-27CE3E8E8CAC}  
TERM.....EXTERNAL CONDOM CATH

DATA VALUE: EXTREMELY DRY

-----  
ABBREVIATION.....EXTR DRY  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Extremely Dry  
HELP\_TEXT.....No Help Text  
ID.....{51677B19-80B9-4787-B6E2-E0BCF6601AA1}  
TERM.....EXTREMELY DRY

DATA VALUE: EYE CLOSED

-----  
ABBREVIATION.....ec  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Eye Closed  
HELP\_TEXT.....No Help Text  
ID.....{140B7D2F-F9C2-0B4C-3CDC-7355EDDB5030}  
TERM.....EYE CLOSED

DATA VALUE: EYEGLASSES

-----  
ABBREVIATION.....Eyeglasses  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Eyeglasses  
HELP\_TEXT.....  
ID.....{E4E5B291-CE6D-409F-A0A0-33EAA47C0C27}  
TERM.....EYEGLASSES

DATA VALUE: FACE WASH

-----  
ABBREVIATION.....FACE WASH  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Face Wash  
HELP\_TEXT.....No Help Text  
ID.....{E1162A8B-D26D-4EB3-988F-50583B1849EA}  
TERM.....FACE WASH

DATA VALUE: FALL RISK ID

-----  
ABBREVIATION.....Fall Risk  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Fall Risk ID  
HELP\_TEXT.....  
ID.....{EF1BE81B-656B-474C-B19E-CCDAD8BB4100}  
TERM.....FALL RISK ID

DATA VALUE: FALL RISK SIGN

-----  
ABBREVIATION.....Fall Sign  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Fall Risk Sign  
HELP\_TEXT.....  
ID.....{E9670720-957E-487E-B0F0-7F7A8AEF69DE}  
TERM.....FALL RISK SIGN

DATA VALUE: FAMILY

-----  
ABBREVIATION.....FAMILY  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Family  
HELP\_TEXT.....No Help Text  
ID.....{1233C4BB-5AC5-4883-AF2A-F5C93CD8814F}  
TERM.....FAMILY

DATA VALUE: FAMILY INFORMED

-----  
ABBREVIATION.....FAM INFORM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Family Informed  
HELP\_TEXT.....No Help Text  
ID.....{77135D32-9CDD-4D3F-9B62-6E8189B77844}  
TERM.....FAMILY INFORMED

DATA VALUE: FAMILY INSTRUCTED

-----  
ABBREVIATION.....fam instru  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Family Instructed  
HELP\_TEXT.....  
ID.....{0FD5988C-2C4D-4786-99B3-3287A355E813}  
TERM.....FAMILY INSTRUCTED

DATA VALUE: FAMILY NEEDS

-----  
ABBREVIATION.....Fam Needs  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Family Needs  
HELP\_TEXT.....  
ID.....{CEB07929-0A6E-47E9-99CD-BD9333304C9C}  
TERM.....FAMILY NEEDS

DATA VALUE: FAMILY STATED

-----  
ABBREVIATION.....Fam Stated  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Family Stated  
HELP\_TEXT.....  
ID.....{07A4B749-65F2-48E8-AD93-D142B04055B4}  
TERM.....FAMILY STATED



DATA VALUE: FAMILY/FRIEND

-----  
ABBREVIATION.....FAM/FRIE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Family/Friend  
HELP\_TEXT.....No Help Text  
ID.....{FA8FE272-E659-41EF-A394-F907A5F55ED8}  
TERM.....FAMILY/FRIEND

DATA VALUE: FASCIA

-----  
ABBREVIATION.....FASCIA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Fascia  
HELP\_TEXT.....No Help Text  
ID.....{EEE8535C-B887-4A34-A93F-C1C4AFC1296F}  
TERM.....FASCIA

DATA VALUE: FAT RESTRICTED 50 GM DIET

-----  
ABBREVIATION.....F RES 50gm  
ACTIVE.....Yes  
DESCRIPTION.....Fat Restricted 50 gm Diet  
DISPLAY\_NAME.....Fat Restricted 50 gm Diet  
HELP\_TEXT.....No Help Text  
ID.....{8E97A639-B7A6-4B3B-9BEF-0A46F3B46E43}  
TERM.....FAT RESTRICTED 50 GM DIET

DATA VALUE: FECAL FAT TEST DIET

-----  
ABBREVIATION.....FEC F TEST  
ACTIVE.....Yes  
DESCRIPTION.....Fecal Fat Test Diet  
DISPLAY\_NAME.....Fecal Fat Test Diet  
HELP\_TEXT.....No Help Text  
ID.....{BA1E8E4A-7698-4EFB-B1DD-18C444EE82D5}  
TERM.....FECAL FAT TEST DIET

DATA VALUE: FENTANYL

-----  
ABBREVIATION.....Fentanyl  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Fentanyl  
HELP\_TEXT.....  
ID.....{C0740659-05FB-4513-BD0F-4F4111B32215}  
TERM.....FENTANYL

DATA VALUE: FHWN1

-----  
ABBREVIATION.....FHWN1  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....FHWN1  
HELP\_TEXT.....No Help Text  
ID.....{93B28875-51BF-4294-B2F9-58BE9365DC81}  
TERM.....FHWN1

DATA VALUE: FIBER FORMULA

-----  
ABBREVIATION.....FIBER FORM  
ACTIVE.....Yes  
DESCRIPTION.....Fiber Enteral Nutrition Formula  
DISPLAY\_NAME.....Fiber Enteral Nutrition Formula  
HELP\_TEXT.....No Help Text  
ID.....{1AF14DA4-8D59-420B-9895-4256E7A818BE}  
TERM.....FIBER FORMULA

DATA VALUE: FIBER STAT

-----  
ABBREVIATION.....Fiber Stat  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Fiber Stat  
HELP\_TEXT.....  
ID.....{201C8A8B-4417-4BD1-AAE1-FE8CBABF2511}  
TERM.....FIBER STAT

DATA VALUE: FIBRINOUS

-----  
ABBREVIATION.....FIBRINOUS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Fibrinous  
HELP\_TEXT.....No Help Text  
ID.....{D286EAE2-6902-4B7E-845D-BBB316941B2C}  
TERM.....FIBRINOUS

DATA VALUE: FINANCIAL CONCERNS

-----  
ABBREVIATION.....Fin Concer  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Financial Concerns  
HELP\_TEXT.....  
ID.....{A7D1159D-DCA1-47AF-8A07-3EB05493CEFD}  
TERM.....FINANCIAL CONCERNS

DATA VALUE: FIRM

-----  
ABBREVIATION.....FIRM  
ACTIVE.....Yes  
DESCRIPTION.....Firm  
DISPLAY\_NAME.....Firm  
HELP\_TEXT.....No Help Text  
ID.....{11E23F20-39A0-4A6D-B198-04BB36B0C909}  
TERM.....FIRM

DATA VALUE: FIRST DEGREE AV BLOCK

-----  
ABBREVIATION.....fdab  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....First Degree Av Block  
HELP\_TEXT.....No Help Text  
ID.....{4FF046C9-B638-F86C-5D51-006B3D5225E3}  
TERM.....FIRST DEGREE AV BLOCK

DATA VALUE: FIXED

-----  
ABBREVIATION.....f  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Fixed  
HELP\_TEXT.....No Help Text  
ID.....{697D8C5A-976F-F423-30EC-A04102CF85F0}  
TERM.....FIXED

DATA VALUE: FLAT

-----  
ABBREVIATION.....Flat  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Flat  
HELP\_TEXT.....No Help Text  
ID.....{EA4F1980-EDF3-82BC-B677-D747F48F1163}  
TERM.....FLAT

DATA VALUE: FLATULENCE

-----  
ABBREVIATION.....FLATULENCE  
ACTIVE.....Yes  
DESCRIPTION.....Flatulence  
DISPLAY\_NAME.....Flatulence  
HELP\_TEXT.....No Help Text  
ID.....{F6ECE7D8-D9C8-485C-8FAC-03DA9D7E0B66}  
TERM.....FLATULENCE

DATA VALUE: FLEXION TO PAIN

-----  
ABBREVIATION.....FLEX TO PA  
ACTIVE.....Yes  
DESCRIPTION.....Flexion To Pain  
DISPLAY\_NAME.....Flexion To Pain  
HELP\_TEXT.....No Help Text  
ID.....{B9F5046F-4F8A-4A20-9D47-90432DC8FEC2}  
TERM.....FLEXION TO PAIN

DATA VALUE: FLICKER OF CONTRACTION

-----  
ABBREVIATION.....foc  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Flicker Of Contraction  
HELP\_TEXT.....No Help Text  
ID.....{218B322B-84AA-F6F6-5207-D86F9A4311F0}  
TERM.....FLICKER OF CONTRACTION

DATA VALUE: FLOTATION BED

-----  
ABBREVIATION.....FLOT BED  
ACTIVE.....Yes  
DESCRIPTION.....Flotation Bed  
DISPLAY\_NAME.....Flotation Bed  
HELP\_TEXT.....No Help Text  
ID.....{FA6329AB-85A6-4301-BC14-D0D089249356}  
TERM.....FLOTATION BED

DATA VALUE: FLUCTUANCE

-----  
ABBREVIATION.....FLUCTUANCE  
ACTIVE.....Yes  
DESCRIPTION.....Fluctuance  
DISPLAY\_NAME.....Fluctuance  
HELP\_TEXT.....No Help Text  
ID.....{A6AD2660-621D-4F39-8485-A013D7491566}  
TERM.....FLUCTUANCE

DATA VALUE: FLUCTUATES

-----  
ABBREVIATION.....Fluctuates  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Fluctuates  
HELP\_TEXT.....No Help Text  
ID.....{C2697650-EF53-5F2E-EEEE3-7B2E298F2FA0}  
TERM.....FLUCTUATES

DATA VALUE: FLUSHED

-----  
ABBREVIATION.....nt  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Flushed  
HELP\_TEXT.....No Help Text  
ID.....{C7C54C1A-6E75-4BCF-4C57-0F8622A5BD16}  
TERM.....FLUSHED

DATA VALUE: FOAM AIR HEEL ELEVATOR

-----  
ABBREVIATION.....FO AI HEEL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Foam Air Heel Elevator  
HELP\_TEXT.....No Help Text  
ID.....{4EC26AA3-4577-4808-8728-F360C5C5A0D4}  
TERM.....FOAM AIR HEEL ELEVATOR

DATA VALUE: FOAM DRESSING

-----  
ABBREVIATION.....WOU DR FOA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Foam Dressing  
HELP\_TEXT.....No Help-Text  
ID.....{D18A5EDD-0EF6-470E-8370-4EDB6D322477}  
TERM.....FOAM DRESSING

DATA VALUE: FOLEY

-----  
ABBREVIATION.....FOLEY  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Foley  
HELP\_TEXT.....No Help Text  
ID.....{9FE7D3A7-67CF-CF54-F721-2F8035CE2BF1}  
TERM.....FOLEY

DATA VALUE: FOOT CARE

-----  
ABBREVIATION.....FOOT CARE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Foot Care  
HELP\_TEXT.....No Help Text  
ID.....{332E8338-6DAE-4BC9-B195-63DAD57937D0}  
TERM.....FOOT CARE

DATA VALUE: FOOT DOWN

-----  
ABBREVIATION.....fd  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Foot Down  
HELP\_TEXT.....No Help Text  
ID.....{13F73BE9-2FA8-FE88-443B-9BA3FE6EDFD7}  
TERM.....FOOT DOWN

DATA VALUE: FOOT UP

-----  
ABBREVIATION.....fu  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Foot Up  
HELP\_TEXT.....No Help Text  
ID.....{49CCC59E-002E-8387-E8D6-D57AAC089EF3}  
TERM.....FOOT UP

DATA VALUE: FOUL SMELLING

-----  
ABBREVIATION.....FOUL SMELL  
ACTIVE.....Yes  
DESCRIPTION.....Foul Smelling  
DISPLAY\_NAME.....Foul Smelling  
HELP\_TEXT.....No Help Text  
ID.....{761346A9-329D-4F42-8D80-0E9D8D40653C}  
TERM.....FOUL SMELLING

DATA VALUE: FOWLER'S

-----  
ABBREVIATION.....FOWLER'S  
ACTIVE.....Yes  
DESCRIPTION.....Fowler's  
DISPLAY\_NAME.....Fowler's  
HELP\_TEXT.....No Help Text  
ID.....{8E691813-A9AC-48BE-9401-528A7D762A32}  
TERM.....FOWLER'S

DATA VALUE: FRANK BLOOD

-----  
ABBREVIATION.....FRANK BLD  
ACTIVE.....Yes  
DESCRIPTION.....Frank Blood  
DISPLAY\_NAME.....Frank Blood  
HELP\_TEXT.....No Help Text  
ID.....{CFAE6E6B-87B2-4CBA-A280-C5233B7995ED}  
TERM.....FRANK BLOOD

DATA VALUE: FREQUENT VERBAL REMINDERS

-----  
ABBREVIATION.....FR VER RE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Frequent Verbal Reminders  
HELP\_TEXT.....No Help Text  
ID.....{F30D791D-BB80-4AB3-B37F-22B6B547EA65}  
TERM.....FREQUENT VERBAL REMINDERS

DATA VALUE: FRIEND

-----  
ABBREVIATION.....FRIEND  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Friend  
HELP\_TEXT.....No Help Text  
ID.....{03902E48-AF86-432A-8372-37796A55F625}  
TERM.....FRIEND

DATA VALUE: FROM A NON-VA OUTSIDE FACILITY

-----  
ABBREVIATION.....Fr Non-VA  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....From a Non-VA Outside Facility  
HELP\_TEXT.....  
ID.....{8090F49E-5056-4F7C-932C-538581996F8F}  
TERM.....FROM A NON-VA OUTSIDE FACILITY

DATA VALUE: FROM ANOTHER UNIT

-----  
ABBREVIATION.....Fr An Unit  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....From Another Unit  
HELP\_TEXT.....  
ID.....{D94B5252-7E8A-48EB-A481-0B18C68CAEF7}  
TERM.....FROM ANOTHER UNIT

DATA VALUE: FROM ANOTHER VA FACILITY

-----  
ABBREVIATION.....Fr Oth VA  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....From Another VA Facility  
HELP\_TEXT.....  
ID.....{66AF9A88-6404-4E05-9640-4FB0C6513BAF}  
TERM.....FROM ANOTHER VA FACILITY

DATA VALUE: FROM CLC

-----  
ABBREVIATION.....From CLC  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....From CLC  
HELP\_TEXT.....  
ID.....{8B7B82C4-2307-47B7-87D7-88F70CF3B9FB}  
TERM.....FROM CLC

DATA VALUE: FROM EMERGENCY DEPARTMENT

-----  
ABBREVIATION.....Fr Emergen  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....From Emergency Department  
HELP\_TEXT.....  
ID.....{9F4AA078-FFA2-4F04-A171-BA0486FCCA0E}  
TERM.....FROM EMERGENCY DEPARTMENT

DATA VALUE: FROM HOME

-----  
ABBREVIATION.....From Home  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....From Home  
HELP\_TEXT.....  
ID.....{2BB78383-F1D2-4FF8-AAE2-C16DC3BD4E8F}  
TERM.....FROM HOME

DATA VALUE: FROM INTERVENTIONAL RADIOLOGY

-----  
ABBREVIATION.....Fr Int Rad  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....From Interventional Radiology  
HELP\_TEXT.....  
ID.....{617DA086-B121-4E2F-8B02-41133EEDD50D}  
TERM.....FROM INTERVENTIONAL RADIOLOGY

DATA VALUE: FROM MICU

-----  
ABBREVIATION.....From MICU  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....From MICU  
HELP\_TEXT.....  
ID.....{971D83AE-B75C-420B-B23D-844397CEE968}  
TERM.....FROM MICU



DATA VALUE: FROM OPERATING ROOM

-----  
ABBREVIATION.....From Op Rm  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....From Operating Room  
HELP\_TEXT.....  
ID.....{2B26BD84-3D4D-4931-99C7-8204E9D49827}  
TERM.....FROM OPERATING ROOM

DATA VALUE: FROM SICU

-----  
ABBREVIATION.....From SICU  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....From SICU  
HELP\_TEXT.....  
ID.....{892002CC-E07A-49A5-B909-1F88665B3D65}  
TERM.....FROM SICU

DATA VALUE: FULL CODE

-----  
ABBREVIATION.....fc  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Full Code  
HELP\_TEXT.....No Help Text  
ID.....{1DED65E6-2F2E-8C5F-D009-BE6467A9A2AA}  
TERM.....FULL CODE

DATA VALUE: FULL LIQUID DIABETIC

-----  
ABBREVIATION.....FULL L DIA  
ACTIVE.....Yes  
DESCRIPTION.....Full Liquid Diabetic  
DISPLAY\_NAME.....Full Liquid Diabetic  
HELP\_TEXT.....No Help Text  
ID.....{B89FAE3C-07EB-4AF6-B123-1153ADADC79F}  
TERM.....FULL LIQUID DIABETIC

DATA VALUE: FULL LIQUID DIET

-----  
ABBREVIATION.....FULL LIQ D  
ACTIVE.....Yes  
DESCRIPTION.....Full Liquid Diet  
DISPLAY\_NAME.....Full Liquid Diet  
HELP\_TEXT.....No Help Text  
ID.....{CAC22234-7B69-44C3-8A0E-9A84A684E429}  
TERM.....FULL LIQUID DIET

DATA VALUE: FULL LOWER

-----  
ABBREVIATION.....Full Low  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Full Lower  
HELP\_TEXT.....  
ID.....{D38022EF-F4F0-44A7-82FE-B277A2A0A165}  
TERM.....FULL LOWER

DATA VALUE: FULL NAME

-----  
ABBREVIATION.....Full Name  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Full Name  
HELP\_TEXT.....  
ID.....{A915EA98-D18D-4887-8C50-4B6C0DF44D1D}  
TERM.....FULL NAME

DATA VALUE: FULL SPINAL PRECAUTIONS

-----  
ABBREVIATION.....FULL S P  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Full Spinal Precautions  
HELP\_TEXT.....No Help Text  
ID.....{42ADB74D-12CC-4111-B181-626C0F01BC7C}  
TERM.....FULL SPINAL PRECAUTIONS

DATA VALUE: FULL SSN

-----  
ABBREVIATION.....Full SSN  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Full SSN  
HELP\_TEXT.....  
ID.....{3E2887AF-2DC5-40FA-9052-7B04347CC802}  
TERM.....FULL SSN

DATA VALUE: FULL THICK SKIN LOSS INVOLVING MUSCLE, BONE, SUPPO

-----  
ABBREVIATION.....F SK L MBS  
ACTIVE.....Yes  
DESCRIPTION.....Full Thickness Skin Loss Involving  
Muscle, Bone, or Supporting Tissue  
DISPLAY\_NAME.....Full Thickness Skin Loss Involving  
Muscle, Bone, or Supporting Tissue  
HELP\_TEXT.....No Help Text  
ID.....{E027FBEF-1C83-41B7-B848-26D9D62F4266}  
TERM.....FULL THICK SKIN LOSS INVOLVING  
MUSCLE, BONE, SUPPO

DATA VALUE: FULL THICK SKIN LOSS INVOLVING SUBCUTANEOUS TISSUE

-----  
ABBREVIATION.....FULL SKINL  
ACTIVE.....Yes  
DESCRIPTION.....Full Thickness Skin Loss Involving  
                  Subcutaneous Tissue  
DISPLAY\_NAME.....Full Thickness Skin Loss Involving  
                  Subcutaneous Tissue  
HELP\_TEXT.....No Help Text  
ID.....{14C31B14-1D35-4CCB-81A3-C6E30801950D}  
TERM.....FULL THICK SKIN LOSS INVOLVING  
                  SUBCUTANEOUS TISSUE

DATA VALUE: FULL UPPER

-----  
ABBREVIATION.....Full Up  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Full Upper  
HELP\_TEXT.....  
ID.....{1808E638-62BB-4978-8706-34BD74B427DE}  
TERM.....FULL UPPER

DATA VALUE: FUNCTIONAL SECTION SET UP

-----  
ABBREVIATION.....Func Sect  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Functional Suction Set Up  
HELP\_TEXT.....  
ID.....{E156C048-11D5-4F62-8B1A-0CEA931E03EF}  
TERM.....FUNCTIONAL SECTION SET UP

DATA VALUE: GAIT BELT

-----  
ABBREVIATION.....GAIT BELT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Gait Belt  
HELP\_TEXT.....No Help Text  
ID.....{AD9AB19B-92B6-40BE-B906-FBFC46C39E0D}  
TERM.....GAIT BELT

DATA VALUE: GASTRIC BYPASS CLEAR LIQUID

-----  
ABBREVIATION.....GAS BYP LI  
ACTIVE.....Yes  
DESCRIPTION.....Gastric Bypass Clear Liquid  
DISPLAY\_NAME.....Gastric Bypass Clear Liquid  
HELP\_TEXT.....No Help Text  
ID.....{6A87E820-7BEA-4636-B102-F9CB83BDD4D9}  
TERM.....GASTRIC BYPASS CLEAR LIQUID

DATA VALUE: GASTRIC BYPASS PUREED

-----  
ABBREVIATION.....GAS BY PUR  
ACTIVE.....Yes  
DESCRIPTION.....Gastric Bypass Pureed  
DISPLAY\_NAME.....Gastric Bypass Pureed  
HELP\_TEXT.....No Help Text  
ID.....{9F9EEB7F-7211-45A2-BCDE-B8398596CBF6}  
TERM.....GASTRIC BYPASS PUREED

DATA VALUE: GASTRIC BYPASS REGULAR

-----  
ABBREVIATION.....GAS BY REG  
ACTIVE.....Yes  
DESCRIPTION.....Gastric Bypass Regular  
DISPLAY\_NAME.....Gastric Bypass Regular  
HELP\_TEXT.....No Help Text  
ID.....{EE831906-EA51-4EE4-AAF2-583DF366A942}  
TERM.....GASTRIC BYPASS REGULAR

DATA VALUE: GASTRIC TUBE

-----  
ABBREVIATION.....G TUBE  
ACTIVE.....Yes  
DESCRIPTION.....Gastric Feeding Tube  
DISPLAY\_NAME.....G Tube  
HELP\_TEXT.....No Help Text  
ID.....{65F774BA-F977-46D3-AE41-7F7C74CFF47E}  
TERM.....GASTRIC TUBE

DATA VALUE: GASTROSTOMY FEEDING BUTTON

-----  
ABBREVIATION.....GAST F BUT  
ACTIVE.....Yes  
DESCRIPTION.....Gastrostomy Feeding Button  
DISPLAY\_NAME.....Gastrostomy Feeding Button  
HELP\_TEXT.....No Help Text  
ID.....{1627B53B-374B-42B4-BBFE-463EE37B6110}  
TERM.....GASTROSTOMY FEEDING BUTTON

DATA VALUE: GASTROSTOMY TUBE

-----  
ABBREVIATION.....GAST TUBE  
ACTIVE.....Yes  
DESCRIPTION.....Gastrostomy Tube  
DISPLAY\_NAME.....Gastrostomy Tube  
HELP\_TEXT.....No Help Text  
ID.....{EF017A39-745E-4C29-821A-5B25F103DF99}  
TERM.....GASTROSTOMY TUBE

DATA VALUE: GAUZE

-----  
ABBREVIATION.....GAUZE  
ACTIVE.....Yes  
DESCRIPTION.....Gauze  
DISPLAY\_NAME.....Gauze  
HELP\_TEXT.....No Help Text  
ID.....{CF10C658-4334-4F3B-A38E-E98894458333}  
TERM.....GAUZE

DATA VALUE: GAUZE AND TAPE

-----  
ABBREVIATION.....GAU & TAPE  
ACTIVE.....Yes  
DESCRIPTION.....Gauze And Tape Dressing  
DISPLAY\_NAME.....Gauze and Tape  
HELP\_TEXT.....No Help Text  
ID.....{5BC40405-B5C4-47CF-81C3-520732E100AC}  
TERM.....GAUZE AND TAPE

DATA VALUE: GAUZE PACKING

-----  
ABBREVIATION.....GAUZE PACK  
ACTIVE.....Yes  
DESCRIPTION.....Gauze Packing  
DISPLAY\_NAME.....Gauze Packing  
HELP\_TEXT.....No Help Text  
ID.....{6C093717-82CB-4966-929E-C909310649C1}  
TERM.....GAUZE PACKING

DATA VALUE: GAUZE WRAP

-----  
ABBREVIATION.....GAUZE WRAP  
ACTIVE.....Yes  
DESCRIPTION.....Gauze Wrap  
DISPLAY\_NAME.....Gauze Wrap  
HELP\_TEXT.....No Help Text  
ID.....{040A9648-E6E1-4AC5-93C7-53A49537B15C}  
TERM.....GAUZE WRAP

DATA VALUE: GENERAL

-----  
ABBREVIATION.....General  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....General  
HELP\_TEXT.....  
ID.....{6ED377F7-55E4-4743-B781-124267567200}  
TERM.....GENERAL

DATA VALUE: GENERAL HEALTH

-----  
ABBREVIATION.....Gen Healt  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....General Health  
HELP\_TEXT.....  
ID.....{D1B589E7-2F59-4A76-8E7B-945F345FA1A0}  
TERM.....GENERAL HEALTH

DATA VALUE: GERI-CHAIR

-----  
ABBREVIATION.....g  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Geri-Chair  
HELP\_TEXT.....No Help Text  
ID.....{167F9786-8F24-054F-020A-658A6EE5B0B4}  
TERM.....GERI-CHAIR

DATA VALUE: GI PROCEDURE ROOM

-----  
ABBREVIATION.....GI Proc Rm  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....GI Procedure Room  
HELP\_TEXT.....  
ID.....{924532F4-B52F-4EA9-8C7D-59168B340118}  
TERM.....GI PROCEDURE ROOM

DATA VALUE: GIVEN DISCHARGE INSTRUCTION AND VERBALIZES UNDERSTANDING

-----  
ABBREVIATION.....DI IN UN  
ACTIVE.....No  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Given Discharge Instructions and  
Verbalizes Understanding  
HELP\_TEXT.....No Help Text  
ID.....{C9A61418-88B2-42C4-B5F4-FC3CF76AB39D}  
TERM.....GIVEN DISCHARGE INSTRUCTION AND  
VERBALIZES UNDERSTANDING

DATA VALUE: GLUCERNA 1.0

-----  
ABBREVIATION.....Glu 1.0  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Glucerna 1.0  
HELP\_TEXT.....  
ID.....{1C465710-80F2-4C13-B8D5-B1D084151104}  
TERM.....GLUCERNA 1.0

DATA VALUE: GLUCERNA 1.2

-----  
ABBREVIATION.....Glu 1.2  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Glucerna 1.2  
HELP\_TEXT.....  
ID.....{F58F1014-D1C4-44EB-9B01-FB88F2B3D751}  
TERM.....GLUCERNA 1.2

DATA VALUE: GLUTEN RESTRICTED

-----  
ABBREVIATION.....GLUT RESTR  
ACTIVE.....Yes  
DESCRIPTION.....Gluten Restricted  
DISPLAY\_NAME.....Gluten Restricted  
HELP\_TEXT.....No Help Text  
ID.....{0805A164-D387-4E67-9D42-BAF1F2543337}  
TERM.....GLUTEN RESTRICTED

DATA VALUE: GNAWING

-----  
ABBREVIATION.....GNAWING  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Gnawing  
HELP\_TEXT.....No Help Text  
ID.....{8AC8EDB0-D7B7-4368-BE95-C35CBADD43E0}  
TERM.....GNAWING

DATA VALUE: GOLD

-----  
ABBREVIATION.....Gold  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Gold  
HELP\_TEXT.....  
ID.....{1D1A2586-C13C-41FA-A3AD-39294BEBA83E}  
TERM.....GOLD

DATA VALUE: GOOD

-----  
ABBREVIATION.....GOOD  
ACTIVE.....Yes  
DESCRIPTION.....Good  
DISPLAY\_NAME.....Good  
HELP\_TEXT.....No Help Text  
ID.....{39ECA304-49F9-44B7-9735-B33CEC9D8119}  
TERM.....GOOD

DATA VALUE: GRAFTS

-----  
ABBREVIATION.....Grafts  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Grafts  
HELP\_TEXT.....  
ID.....{E2AEE65A-7387-4A58-80E6-0212A1F3D1C7}  
TERM.....GRAFTS

DATA VALUE: GRANULATING

-----  
ABBREVIATION.....GRANULATI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Granulating  
HELP\_TEXT.....No Help Text  
ID.....{F8930276-3954-4BFE-82B0-6C44A8E9CAE0}  
TERM.....GRANULATING

DATA VALUE: GRAVITY

-----  
ABBREVIATION.....GRAVITY  
ACTIVE.....Yes  
DESCRIPTION.....Gravity  
DISPLAY\_NAME.....Gravity  
HELP\_TEXT.....No Help Text  
ID.....{C0233F70-376D-41CE-9F14-714C7B35244F}  
TERM.....GRAVITY

DATA VALUE: GRAY

-----  
ABBREVIATION.....GRAY  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Gray  
HELP\_TEXT.....No Help Text  
ID.....{D38ABEE0-52D7-4316-A713-40ACF4250599}  
TERM.....GRAY

DATA VALUE: GREEN

-----  
ABBREVIATION.....GREEN  
ACTIVE.....Yes  
DESCRIPTION.....The color green  
DISPLAY\_NAME.....Green  
HELP\_TEXT.....The color green  
ID.....{3D1F1553-56FC-68CE-79A9-E4F6357B6AAF}  
TERM.....GREEN



DATA VALUE: GRIMACING

-----  
ABBREVIATION.....GRIMACING  
ACTIVE.....Yes  
DESCRIPTION.....Grimacing  
DISPLAY\_NAME.....Grimacing  
HELP\_TEXT.....No Help Text  
ID.....{69EFDB2D-A0C4-48D4-A7E0-A862F34545C8}  
TERM.....GRIMACING

DATA VALUE: GROUP ROOM

-----  
ABBREVIATION.....gr  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Group Room  
HELP\_TEXT.....No Help Text  
ID.....{D41C3D9F-8F10-1A0F-D571-20DE0EA36328}  
TERM.....GROUP ROOM

DATA VALUE: GU PROCEDURE ROOM

-----  
ABBREVIATION.....GU Proc Rm  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....GU Procedure Room  
HELP\_TEXT.....  
ID.....{D4B6678F-C1DD-4A23-9045-574251F6CD74}  
TERM.....GU PROCEDURE ROOM

DATA VALUE: GUARDING

-----  
ABBREVIATION.....GUARDING  
ACTIVE.....Yes  
DESCRIPTION.....Guarding  
DISPLAY\_NAME.....Guarding  
HELP\_TEXT.....No Help Text  
ID.....{9EFD10A7-B1B7-43DA-B289-646AF5C0A25C}  
TERM.....GUARDING

DATA VALUE: GURNEY CHAIR

-----  
ABBREVIATION.....GURN CHAIR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Gurney Chair  
HELP\_TEXT.....No Help Text  
ID.....{B6EA0944-8414-42FD-88C0-50155DE82896}  
TERM.....GURNEY CHAIR

DATA VALUE: HALL

-----  
ABBREVIATION.....h  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Hall  
HELP\_TEXT.....No Help Text  
ID.....{BFCEA10F-2663-199E-88D8-E263122D7D26}  
TERM.....HALL

DATA VALUE: HAND HYGIENE BEFORE AND AFTER

-----  
ABBREVIATION.....Hand Hyg  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Hand Hygiene Before and After  
HELP\_TEXT.....  
ID.....{6005023D-3C02-4C93-99E8-2A8167574E93}  
TERM.....HAND HYGIENE BEFORE AND AFTER

DATA VALUE: HARD

-----  
ABBREVIATION.....HARD  
ACTIVE.....Yes  
DESCRIPTION.....Hard  
DISPLAY\_NAME.....Hard  
HELP\_TEXT.....No Help Text  
ID.....{7934B543-7028-4CBC-9308-99F5A952A1DE}  
TERM.....HARD

DATA VALUE: HAS NOT SLEPT

-----  
ABBREVIATION.....NO SLEEP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Has Not Slept  
HELP\_TEXT.....No Help Text  
ID.....{3BF9DB54-5D1E-422D-90AF-918F81C9EA85}  
TERM.....HAS NOT SLEPT

DATA VALUE: HEAD DOWN

-----  
ABBREVIATION.....hd  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Head Down  
HELP\_TEXT.....No Help Text  
ID.....{80BF3425-D252-F917-052E-DC54D39E464E}  
TERM.....HEAD DOWN

DATA VALUE: HEAD UP

-----  
ABBREVIATION.....hu  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Head Up  
HELP\_TEXT.....No Help Text  
ID.....{BB1C79C0-80F1-7CDE-CC10-7CA0C704AE35}  
TERM.....HEAD UP

DATA VALUE: HEADACHE

-----  
ABBREVIATION.....Headache  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Headache  
HELP\_TEXT.....  
ID.....{9338E67C-201C-4066-A734-8BD4B47462BD}  
TERM.....HEADACHE

DATA VALUE: HEALTH TECH

-----  
ABBREVIATION.....Health Tec  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Health Tech  
HELP\_TEXT.....  
ID.....{7C5F34C7-DBE6-4FCB-B46D-F802FBBE4CEF}  
TERM.....HEALTH TECH

DATA VALUE: HEARING

-----  
ABBREVIATION.....HEARING  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Hearing  
HELP\_TEXT.....No Help Text  
ID.....{D1A8AC66-1A5F-4F7D-AA89-529AF7A4AA11}  
TERM.....HEARING

DATA VALUE: HEARING AIDS

-----  
ABBREVIATION.....Hear Aids  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Hearing Aids  
HELP\_TEXT.....  
ID.....{504FF149-8889-444E-925B-C952314FDF98}  
TERM.....HEARING AIDS

DATA VALUE: HEARTBURN

-----  
ABBREVIATION.....HEARTBURN  
ACTIVE.....Yes  
DESCRIPTION.....Heartburn  
DISPLAY\_NAME.....Heartburn  
HELP\_TEXT.....Heartburn  
ID.....{5ECEAF5E-6B18-4177-A18D-2DEE87796CC2}  
TERM.....HEARTBURN

DATA VALUE: HEAT

-----  
ABBREVIATION.....HEAT  
ACTIVE.....Yes  
DESCRIPTION.....Heat  
DISPLAY\_NAME.....Heat  
HELP\_TEXT.....No Help Text  
ID.....{E0FCA383-3AB8-4E03-9ABC-00DE3FE01431}  
TERM.....HEAT

DATA VALUE: HEAT CHEM PACK

-----  
ABBREVIATION.....HT CHEM PA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Heat Chemi Pack  
HELP\_TEXT.....No Help Text  
ID.....{D759F891-2535-40C2-A68C-678170E56F7A}  
TERM.....HEAT CHEM PACK

DATA VALUE: HEAT LAMP

-----  
ABBREVIATION.....HEAT LAMP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Heat Lamp  
HELP\_TEXT.....Heat Lamp  
ID.....{4E9DE229-F810-4AB0-9FD7-4596759CEA11}  
TERM.....HEAT LAMP

DATA VALUE: HEAVY

-----  
ABBREVIATION.....HEAVY  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Heavy  
HELP\_TEXT.....No Help Text  
ID.....{1F10CD98-03FA-4705-BF7D-941E881C4E3D}  
TERM.....HEAVY

DATA VALUE: HEEL LIFT

-----  
ABBREVIATION.....HEEL LIFT  
ACTIVE.....Yes  
DESCRIPTION.....Heel Lift  
DISPLAY\_NAME.....Heel Lift  
HELP\_TEXT.....No Help Text  
ID.....{974F90E6-382D-421A-A4CE-94A74CBEFC32}  
TERM.....HEEL LIFT

DATA VALUE: HEEL PILLOW

-----  
ABBREVIATION.....HEEL PILLO  
ACTIVE.....Yes  
DESCRIPTION.....Heel Pillow  
DISPLAY\_NAME.....Heel Pillow  
HELP\_TEXT.....No Help Text  
ID.....{F0F1D997-C780-40BF-AA88-78262E81BFAA}  
TERM.....HEEL PILLOW

DATA VALUE: HELMET

-----  
ABBREVIATION.....Helmet  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Helmet  
HELP\_TEXT.....  
ID.....{E65974F8-BCE7-4361-868D-C0B4D780EA21}  
TERM.....HELMET

DATA VALUE: HEMATOMA

-----  
ABBREVIATION.....HEMATOMA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Hematoma  
HELP\_TEXT.....No Help Text  
ID.....{6A24B7EE-62F3-C952-C1A3-086CDBB1FE68}  
TERM.....HEMATOMA

DATA VALUE: HEMATURIA

-----  
ABBREVIATION.....HEMATURIA  
ACTIVE.....Yes  
DESCRIPTION.....Blood present in urine  
DISPLAY\_NAME.....Hematuria  
HELP\_TEXT.....Blood present in urine  
ID.....{ABD21F4E-BC48-D43A-38CF-75DB7704C4B8}  
TERM.....HEMATURIA

DATA VALUE: HEMORRHAGE

-----  
ABBREVIATION.....HEMORRHAGE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Hemorrhage  
HELP\_TEXT.....No Help Text  
ID.....{E00AC828-83F3-4F69-B97B-D6AD1D87A440}  
TERM.....HEMORRHAGE

DATA VALUE: HEMOVAC

-----  
ABBREVIATION.....HEMOVAC  
ACTIVE.....Yes  
DESCRIPTION.....Hemovac  
DISPLAY\_NAME.....Hemovac  
HELP\_TEXT.....No Help Text  
ID.....{82DA2DA7-E20B-E571-659C-457DF42472C2}  
TERM.....HEMOVAC

DATA VALUE: HEPATIC FORMULA

-----  
ABBREVIATION.....HEPAT FORM  
ACTIVE.....Yes  
DESCRIPTION.....Hepatic Enteral Nutrition Formula  
DISPLAY\_NAME.....Hepatic Enteral Nutrition Formula  
HELP\_TEXT.....No Help Text  
ID.....{1BC0A09A-CCD7-457F-BC08-9B6FC5144BF3}  
TERM.....HEPATIC FORMULA

DATA VALUE: HEPLOCK

-----  
ABBREVIATION.....HEPLOCK  
ACTIVE.....Yes  
DESCRIPTION.....Heplock  
DISPLAY\_NAME.....Heplock  
HELP\_TEXT.....No Help Text  
ID.....{C827FEDE-59D6-4C45-BDDA-4884FA0882F3}  
TERM.....HEPLOCK

DATA VALUE: HIGH

-----  
ABBREVIATION.....h  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....High  
HELP\_TEXT.....No Help Text  
ID.....{D88CB5EC-2AD4-2F80-2D43-40ECB0B7D390}  
TERM.....HIGH

DATA VALUE: HIGH CALORIE

-----  
ABBREVIATION.....HIGH CALOR  
ACTIVE.....Yes  
DESCRIPTION.....High Calorie  
DISPLAY\_NAME.....High Calorie  
HELP\_TEXT.....No Help Text  
ID.....{FC453BF0-1E7E-4733-9605-AC3425F16205}  
TERM.....HIGH CALORIE

DATA VALUE: HIGH FALL RISK

-----  
ABBREVIATION.....hfr  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....High Fall Risk  
HELP\_TEXT.....No Help Text  
ID.....{276C15F3-5BED-0FE0-5649-337A23AB1386}  
TERM.....HIGH FALL RISK

DATA VALUE: HIGH FIBER

-----  
ABBREVIATION.....HIGH FIBER  
ACTIVE.....Yes  
DESCRIPTION.....High Fiber  
DISPLAY\_NAME.....High Fiber  
HELP\_TEXT.....No Help Text  
ID.....{3CAF7EE0-6A2D-41A3-BEDC-C57ECBBA72A2}  
TERM.....HIGH FIBER

DATA VALUE: HIGH FOWLER'S

-----  
ABBREVIATION.....HIGH FOWLE  
ACTIVE.....Yes  
DESCRIPTION.....High Fowler's  
DISPLAY\_NAME.....High Fowler's  
HELP\_TEXT.....No Help Text  
ID.....{EA23C94D-EE91-4552-9F74-E051FAF3E194}  
TERM.....HIGH FOWLER'S

DATA VALUE: HIGH PROTEIN

-----  
ABBREVIATION.....HIGH PROTE  
ACTIVE.....Yes  
DESCRIPTION.....High Protein  
DISPLAY\_NAME.....High Protein  
HELP\_TEXT.....No Help Text  
ID.....{1595A29D-CA91-45AE-893B-856EB909E4DD}  
TERM.....HIGH PROTEIN

DATA VALUE: HIP-PROTECTORS

-----  
ABBREVIATION.....Hip-Prot  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Hip-Protectors  
HELP\_TEXT.....  
ID.....{EA08514E-84C1-4CF4-B5B1-E5E14EDA82C1}  
TERM.....HIP-PROTECTORS

DATA VALUE: HIRSUTISM

-----  
ABBREVIATION.....HIRSUTISM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Hirsutism  
HELP\_TEXT.....No Help Text  
ID.....{A909A4FC-080A-4E9A-867A-0B30D4633602}  
TERM.....HIRSUTISM

DATA VALUE: HM: HAND MOTIONS

-----  
ABBREVIATION.....HM  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....HM: Hand Motions  
HELP\_TEXT.....  
ID.....{5180AD3D-4377-435B-9F47-FE59F3E3501B}  
TERM.....HM: HAND MOTIONS

DATA VALUE: HOME

-----  
ABBREVIATION.....HOME  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Home  
HELP\_TEXT.....No Help Text  
ID.....{C8E5E320-361B-421C-A6A9-13DC16125440}  
TERM.....HOME

DATA VALUE: HOME HEALTH

-----  
ABBREVIATION.....HOME HEALT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Home Health  
HELP\_TEXT.....No Help Text  
ID.....{F94FB70A-AFFD-4A5E-A829-0AF933B044BB}  
TERM.....HOME HEALTH



DATA VALUE: HOSTILE

-----  
ABBREVIATION.....h  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Hostile  
HELP\_TEXT.....No Help Text  
ID.....{29EF2C04-7DB0-4491-2B5E-D92A6A558721}  
TERM.....HOSTILE

DATA VALUE: HOT

-----  
ABBREVIATION.....HOT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Hot  
HELP\_TEXT.....No Help Text  
ID.....{401EA6CD-12DE-AB0E-2B52-8B2E62F7C34E}  
TERM.....HOT

DATA VALUE: HUMIDIFIER

-----  
ABBREVIATION.....HUMIDIFIER  
ACTIVE.....Yes  
DESCRIPTION.....Humidifier  
DISPLAY\_NAME.....Humidifier  
HELP\_TEXT.....Humidifier  
ID.....{A4CE67D1-F885-423B-B914-392D8CAC697E}  
TERM.....HUMIDIFIER

DATA VALUE: HYDROCOLLOIDAL DRESSING

-----  
ABBREVIATION.....HCOLL DRES  
ACTIVE.....Yes  
DESCRIPTION.....Hydrocolloidal Dressing  
DISPLAY\_NAME.....Hydrocolloidal Dressing  
HELP\_TEXT.....No Help Text  
ID.....{D85741EF-0DB3-48CF-BC9A-88C834A72F36}  
TERM.....HYDROCOLLOIDAL DRESSING

DATA VALUE: HYDROFIBER

-----  
ABBREVIATION.....H-FIBER  
ACTIVE.....Yes  
DESCRIPTION.....Hydrofiber  
DISPLAY\_NAME.....Hydrofiber  
HELP\_TEXT.....No Help Text  
ID.....{42740A7F-A831-4FBD-9E59-88F315176674}  
TERM.....HYDROFIBER

DATA VALUE: HYDROFIBER WITH SILVER

-----  
ABBREVIATION.....HYFIB W SI  
ACTIVE.....Yes  
DESCRIPTION.....Hydrofiber With Silver  
DISPLAY\_NAME.....Hydrofiber With Silver  
HELP\_TEXT.....No Help Text  
ID.....{04BD9512-DE5D-421A-9818-8A33D1F0693D}  
TERM.....HYDROFIBER WITH SILVER

DATA VALUE: HYDROGEL DRESSING

-----  
ABBREVIATION.....H-GEL DRES  
ACTIVE.....Yes  
DESCRIPTION.....Hydrogel Dressing  
DISPLAY\_NAME.....Hydrogel Dressing  
HELP\_TEXT.....No Help Text  
ID.....{F6192D10-D6EA-42C2-AB15-270659D9925C}  
TERM.....HYDROGEL DRESSING

DATA VALUE: HYDROMORPHONE

-----  
ABBREVIATION.....Hydromorph  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Hydromorphone  
HELP\_TEXT.....  
ID.....{8E6E9DB1-C969-4274-9761-C9830A9FCF67}  
TERM.....HYDROMORPHONE

DATA VALUE: HYPERACTIVE

-----  
ABBREVIATION.....HYPERACTIV  
ACTIVE.....Yes  
DESCRIPTION.....Hyperactive  
DISPLAY\_NAME.....Hyperactive  
HELP\_TEXT.....No Help Text  
ID.....{1A548FFE-8645-4C0D-899D-5D91430CC1B1}  
TERM.....HYPERACTIVE

DATA VALUE: HYPERDYNAMIC

-----  
ABBREVIATION.....HYPERDYNAM  
ACTIVE.....Yes  
DESCRIPTION.....Hyperdynamic  
DISPLAY\_NAME.....Hyperdynamic  
HELP\_TEXT.....No Help Text  
ID.....{818A8A99-04D9-4989-A9BA-ECFACED74605}  
TERM.....HYPERDYNAMIC

DATA VALUE: HYPERGRANULIZATION

-----  
ABBREVIATION.....HYPERGRANU  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Hypergranulization  
HELP\_TEXT.....No Help Text  
ID.....{C14B9FB0-6F33-4ACB-88B3-4B799594F845}  
TERM.....HYPERGRANULIZATION

DATA VALUE: HYPERTONIC HYDROGEL

-----  
ABBREVIATION.....WOU DR HYP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Hypertonic Hydrogel  
HELP\_TEXT.....No Help-Text  
ID.....{533D6FC3-5183-4417-AB89-9FAD9BAD0F44}  
TERM.....HYPERTONIC HYDROGEL

DATA VALUE: HYPOACTIVE

-----  
ABBREVIATION.....HYPOACTIVE  
ACTIVE.....Yes  
DESCRIPTION.....Hypoactive  
DISPLAY\_NAME.....Hypoactive  
HELP\_TEXT.....No Help Text  
ID.....{2E0C516D-3858-4A1F-A2F3-BF0AB9E3A7FC}  
TERM.....HYPOACTIVE

DATA VALUE: HYPOALBUMINEMIA

-----  
ABBREVIATION.....HYPOALBUM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Hypoalbuminemia  
HELP\_TEXT.....No Help Text  
ID.....{A3397F38-2491-48E7-8366-A6FA2D91E866}  
TERM.....HYPOALBUMINEMIA

DATA VALUE: HYPOTENSION

-----  
ABBREVIATION.....Hypotensio  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Hypotension  
HELP\_TEXT.....  
ID.....{1AE4E3BB-3B78-42D2-BDAC-83563A0546E1}  
TERM.....HYPOTENSION

DATA VALUE: HYPOXIA

-----  
ABBREVIATION.....Hypoxia  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Hypoxia  
HELP\_TEXT.....  
ID.....{7599629D-1206-4A1A-9D02-EEE384E5FE51}  
TERM.....HYPOXIA

DATA VALUE: I

-----  
ABBREVIATION.....I  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....I  
HELP\_TEXT.....No Help Text  
ID.....{BEEDFF9A-9401-4D2A-8345-126D8CFFC8D8}  
TERM.....I

DATA VALUE: IABP

-----  
ABBREVIATION.....IABP  
ACTIVE.....Yes  
DESCRIPTION.....Intra-Aortic Balloon Pump  
DISPLAY\_NAME.....IABP  
HELP\_TEXT.....No Help Text  
ID.....{763FE975-F5E5-80A6-AB16-E903A829E665}  
TERM.....IABP

DATA VALUE: ICE PACK

-----  
ABBREVIATION.....ICE PACK  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Ice Pack  
HELP\_TEXT.....No Help Text  
ID.....{E46ED309-D60D-45DF-BFA3-F74DD13FC4AE}  
TERM.....ICE PACK

DATA VALUE: ICTERIC

-----  
ABBREVIATION.....ICTERIC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Icteric  
HELP\_TEXT.....Relating to or affected with jaundice  
ID.....{274ED5BB-5A34-20EE-A9B0-81AB7803780E}  
TERM.....ICTERIC

DATA VALUE: ICU

-----  
ABBREVIATION.....ICU  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....ICU  
HELP\_TEXT.....No Help Text  
ID.....{E7FA69FE-1088-4853-AFBD-2D815E3C2EB0}  
TERM.....ICU

DATA VALUE: ID BAND

-----  
ABBREVIATION.....ID Band  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....ID Band  
HELP\_TEXT.....  
ID.....{5526C5AB-64C8-4AB1-B77F-DAF45D26B439}  
TERM.....ID BAND

DATA VALUE: ID CARD

-----  
ABBREVIATION.....ID Card  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....ID Card  
HELP\_TEXT.....  
ID.....{CFD540F0-3D73-498F-B590-CCEF11499078}  
TERM.....ID CARD

DATA VALUE: IDENTIFICATION

-----  
ABBREVIATION.....i  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Identification  
HELP\_TEXT.....No Help Text  
ID.....{39274F69-556C-7F05-6D17-CE8A03D1867D}  
TERM.....IDENTIFICATION

DATA VALUE: IDIOVENTRICULAR

-----  
ABBREVIATION.....i  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Idioventricular  
HELP\_TEXT.....No Help Text  
ID.....{8BD08966-F536-8E74-0BDA-A95D0D7F8B6C}  
TERM.....IDIOVENTRICULAR

DATA VALUE: II

-----  
ABBREVIATION.....II  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....II  
HELP\_TEXT.....No Help Text  
ID.....{F165C13B-397E-4182-8D72-BD1D74CFBA00}  
TERM.....II

DATA VALUE: III

-----  
ABBREVIATION.....III  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....III  
HELP\_TEXT.....No Help Text  
ID.....{D8FE8FA9-9BB5-416E-AE2C-60E627F76E5A}  
TERM.....III

DATA VALUE: ILEAL CONDUIT

-----  
ABBREVIATION.....ILEAL COND  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Ileal Conduit  
HELP\_TEXT.....No Help Text  
ID.....{032CB6FE-77D2-702E-A0AE-46176F5CAA57}  
TERM.....ILEAL CONDUIT

DATA VALUE: ILEAL-CONDUIT

-----  
ABBREVIATION.....Ileal-Con  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Ileal-Conduit  
HELP\_TEXT.....  
ID.....{C91C4E4A-2D0D-438A-8EF5-8717D703A777}  
TERM.....ILEAL-CONDUIT

DATA VALUE: ILEOSTOMY

-----  
ABBREVIATION.....ILEOSTOMY  
ACTIVE.....Yes  
DESCRIPTION.....Ileostomy  
DISPLAY\_NAME.....Ileostomy  
HELP\_TEXT.....No Help Text  
ID.....{D5DF49C8-C6BA-44C1-9840-7765907730A3}  
TERM.....ILEOSTOMY

DATA VALUE: IMAGING

-----  
ABBREVIATION.....IMAGING  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Imaging  
HELP\_TEXT.....No Help Text  
ID.....{791A322E-F45E-40EB-B95D-E5B62E1D0A75}  
TERM.....IMAGING

DATA VALUE: IMPAIRED LISTENING SKILLS

-----  
ABBREVIATION.....IMP LIS SK  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Impaired Listening Skills  
HELP\_TEXT.....No Help Text  
ID.....{9CE0E7E7-79E5-4910-918C-2022C9830171}  
TERM.....IMPAIRED LISTENING SKILLS

DATA VALUE: IMPLANTS

-----  
ABBREVIATION.....Implants  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Implants  
HELP\_TEXT.....  
ID.....{14998852-14FD-4C6A-BB8A-5073D6094B11}  
TERM.....IMPLANTS

DATA VALUE: INAPPROPRIATE

-----  
ABBREVIATION.....Inappropri  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Inappropriate  
HELP\_TEXT.....No Help Text  
ID.....{7E3B54C2-8A49-28EC-055E-FEA87599421C}  
TERM.....INAPPROPRIATE

DATA VALUE: INCISION DISCHARGE CHARACTERISTICS BLOODY

-----  
ABBREVIATION.....INC DI BLO  
ACTIVE.....Yes  
DESCRIPTION.....Incision Discharge Characteristics  
                  Bloody  
DISPLAY\_NAME.....Bloody  
HELP\_TEXT.....Incision Discharge Characteristics  
                  Bloody  
ID.....{D185E55B-EC21-49B6-A040-5422382C3B90}  
TERM.....INCISION DISCHARGE CHARACTERISTICS  
                  BLOODY

DATA VALUE: INCISION DISCHARGE CHARACTERISTICS CLOUDY

-----  
ABBREVIATION.....INC DI CLO  
ACTIVE.....Yes  
DESCRIPTION.....Incision Discharge Characteristics  
                  Cloudy  
DISPLAY\_NAME.....Cloudy  
HELP\_TEXT.....Incision Discharge Characteristics  
                  Cloudy  
ID.....{3EBBB25E-1B90-415B-BCB5-3C3DAC33C118}  
TERM.....INCISION DISCHARGE CHARACTERISTICS  
                  CLOUDY

DATA VALUE: INCISION DISCHARGE CHARACTERISTICS NONE

-----  
ABBREVIATION.....INC DI NON  
ACTIVE.....Yes  
DESCRIPTION.....Incision Discharge Characteristics None  
DISPLAY\_NAME.....None  
HELP\_TEXT.....Incision Discharge Characteristics None  
ID.....{4F118F41-8C09-400E-A303-0955DFDF54E2}  
TERM.....INCISION DISCHARGE CHARACTERISTICS NONE

DATA VALUE: INCISION DISCHARGE CHARACTERISTICS PURULENT

-----  
ABBREVIATION.....INC DI PUR  
ACTIVE.....Yes  
DESCRIPTION.....Incision Discharge Characteristics  
                  Purulent  
DISPLAY\_NAME.....Purulent  
HELP\_TEXT.....Incision Discharge Characteristics  
                  Purulent  
ID.....{B2D71BDF-41AB-4F10-ABE3-A9D9F338F6F1}  
TERM.....INCISION DISCHARGE CHARACTERISTICS  
                  PURULENT

DATA VALUE: INCISION DISCHARGE CHARACTERISTICS SEROSANGUINOUS

-----  
ABBREVIATION.....INC DI SSA  
ACTIVE.....Yes  
DESCRIPTION.....Incision Discharge Characteristics  
                  Serosanguinous  
DISPLAY\_NAME.....Serosanguinous  
HELP\_TEXT.....Incision Discharge Characteristics  
                  Serosanguinous  
ID.....{102F19A9-CC28-4156-88DE-1B4C8A682AC2}  
TERM.....INCISION DISCHARGE CHARACTERISTICS  
                  SEROSANGUINOUS

DATA VALUE: INCISION DISCHARGE CHARACTERISTICS SEROUS

-----  
ABBREVIATION.....INC DI SER  
ACTIVE.....Yes  
DESCRIPTION.....Incision Discharge Characteristics  
                  Serous  
DISPLAY\_NAME.....Serous  
HELP\_TEXT.....Incision Discharge Characteristics  
                  Serous  
ID.....{CC4F3BA1-AD0D-42BA-9D76-1FE66F864FC8}  
TERM.....INCISION DISCHARGE CHARACTERISTICS



SEROUS

DATA VALUE: INCISION DISCHARGE CHARACTERISTICS VISCOUS

-----  
ABBREVIATION.....INC DI VIS  
ACTIVE.....Yes  
DESCRIPTION.....Incision Discharge Characteristics  
                  Viscous  
DISPLAY\_NAME.....Viscous  
HELP\_TEXT.....Incision Discharge Characteristics  
                  Viscous  
ID.....{A38C2D8C-25B1-41E5-A69F-7A5C439C6EB9}  
TERM.....INCISION DISCHARGE CHARACTERISTICS  
                  VISCOUS

DATA VALUE: INCISION DRAINAGE AMOUNT HEAVY

-----  
ABBREVIATION.....I DR A HEA  
ACTIVE.....Yes  
DESCRIPTION.....Incision Drainage Amount Heavy  
DISPLAY\_NAME.....Heavy  
HELP\_TEXT.....Incision Drainage Amount Heavy  
ID.....{D77FCEA0-D247-4DD5-873E-2F95C7FF48A6}  
TERM.....INCISION DRAINAGE AMOUNT HEAVY

DATA VALUE: INCISION DRAINAGE AMOUNT MODERATE

-----  
ABBREVIATION.....I DR A MOD  
ACTIVE.....Yes  
DESCRIPTION.....Incision Drainage Amount Moderate  
DISPLAY\_NAME.....Moderate  
HELP\_TEXT.....Incision Drainage Amount Moderate  
ID.....{9B0487A9-DCB4-42D3-92E6-805C412B5D34}  
TERM.....INCISION DRAINAGE AMOUNT MODERATE

DATA VALUE: INCISION DRAINAGE AMOUNT NONE

-----  
ABBREVIATION.....I DR A NON  
ACTIVE.....Yes  
DESCRIPTION.....Incision Drainage Amount None  
DISPLAY\_NAME.....None  
HELP\_TEXT.....Incision Drainage Amount None  
ID.....{5D0D9027-4B20-43A4-B437-EA7B77CB2105}  
TERM.....INCISION DRAINAGE AMOUNT NONE

DATA VALUE: INCISION DRAINAGE AMOUNT SCANT

-----  
ABBREVIATION.....I DR A SCA  
ACTIVE.....Yes  
DESCRIPTION.....Incision Drainage Amount Scant  
DISPLAY\_NAME.....Scant  
HELP\_TEXT.....Incision Drainage Amount Scant  
ID.....{A476E395-46AA-44E7-8A00-CE22C7066BB1}  
TERM.....INCISION DRAINAGE AMOUNT SCANT

DATA VALUE: INCISION DRAINAGE AMOUNT SMALL

-----  
ABBREVIATION.....I DR A SMA  
ACTIVE.....Yes  
DESCRIPTION.....Incision Drainage Amount Small  
DISPLAY\_NAME.....Small  
HELP\_TEXT.....Incision Drainage Amount Small  
ID.....{4479479D-8ADA-4FBF-A5A2-426D49FC94D3}  
TERM.....INCISION DRAINAGE AMOUNT SMALL

DATA VALUE: INCISION DRESSING TYPE ANTIMICROBIAL

-----  
ABBREVIATION.....INC DR ANT  
ACTIVE.....Yes  
DESCRIPTION.....Incision Dressing Type Antimicrobial  
DISPLAY\_NAME.....Antimicrobial  
HELP\_TEXT.....Incision Dressing Type Antimicrobial  
ID.....{35F5F5C1-B76A-44D2-8576-9D7D3ECA3197}  
TERM.....INCISION DRESSING TYPE ANTIMICROBIAL

DATA VALUE: INCISION DRESSING TYPE ANTIMICROBIAL DISK

-----  
ABBREVIATION.....INC DR ADI  
ACTIVE.....Yes  
DESCRIPTION.....Incision Dressing Type Antimicrobial  
Disk  
DISPLAY\_NAME.....Antimicrobial Disk  
HELP\_TEXT.....Incision Dressing Type Antimicrobial  
Disk  
ID.....{55A33062-7D6B-4CDF-9131-D4ADC38F6CE3}  
TERM.....INCISION DRESSING TYPE ANTIMICROBIAL  
DISK

DATA VALUE: INCISION DRESSING TYPE BALSAM OF PERU

-----  
ABBREVIATION.....INC DR BAL  
ACTIVE.....Yes  
DESCRIPTION.....Incision Dressing Type Balsam of Peru  
DISPLAY\_NAME.....Balsam Of Peru  
HELP\_TEXT.....Incision Dressing Type Balsam of Peru  
ID.....{A1FB45C6-6315-4282-8D89-DCB0C018E357}  
TERM.....INCISION DRESSING TYPE BALSAM OF PERU

DATA VALUE: INCISION DRESSING TYPE BISMUTH GAUZE

-----  
ABBREVIATION.....INC DR BIS  
ACTIVE.....Yes  
DESCRIPTION.....Incision Dressing Type Bismuth Gauze  
DISPLAY\_NAME.....Bismuth Gauze  
HELP\_TEXT.....Incision Dressing Type Bismuth Gauze  
ID.....{CA88C5DC-3890-49A2-800A-0528CE194BA8}  
TERM.....INCISION DRESSING TYPE BISMUTH GAUZE

DATA VALUE: INCISION DRESSING TYPE CALCIUM ALGINATE

-----  
ABBREVIATION.....INC DR CAL  
ACTIVE.....Yes  
DESCRIPTION.....Incision Dressing Type Calcium Alginate  
DISPLAY\_NAME.....Calcium Alginate  
HELP\_TEXT.....Incision Dressing Type Calcium Alginate  
ID.....{B20558BA-FB9E-499E-82DF-49616ACE9375}  
TERM.....INCISION DRESSING TYPE CALCIUM ALGINATE

DATA VALUE: INCISION DRESSING TYPE FOAM DRESSING

-----  
ABBREVIATION.....INC DR FOA  
ACTIVE.....Yes  
DESCRIPTION.....Incision Dressing Type Foam dressing  
DISPLAY\_NAME.....Foam Dressing  
HELP\_TEXT.....Incision Dressing Type Foam dressing  
ID.....{5C2C2D24-6022-4BBF-9065-03FB5971CB3F}  
TERM.....INCISION DRESSING TYPE FOAM DRESSING

DATA VALUE: INCISION DRESSING TYPE GAUZE

-----  
ABBREVIATION.....INC DR GAU  
ACTIVE.....Yes  
DESCRIPTION.....Incision Dressing Type Gauze  
DISPLAY\_NAME.....Gauze  
HELP\_TEXT.....Incision Dressing Type Gauze  
ID.....{E4AC9307-0512-4F9D-A254-FE19E654461F}  
TERM.....INCISION DRESSING TYPE GAUZE

DATA VALUE: INCISION DRESSING TYPE GAUZE PACKING

-----  
ABBREVIATION.....INC DR GPA  
ACTIVE.....Yes  
DESCRIPTION.....Incision Dressing Type Gauze Packing  
DISPLAY\_NAME.....Gauze Packing  
HELP\_TEXT.....Incision Dressing Type Gauze Packing  
ID.....{3561DAEB-8A30-4E16-B920-F7DFE022EBB1}  
TERM.....INCISION DRESSING TYPE GAUZE PACKING

DATA VALUE: INCISION DRESSING TYPE GAUZE WRAP

-----  
ABBREVIATION.....INC DR GWR  
ACTIVE.....Yes  
DESCRIPTION.....Incision Dressing Type Gauze Wrap  
DISPLAY\_NAME.....Gauze Wrap  
HELP\_TEXT.....Incision Dressing Type Gauze Wrap  
ID.....{E1DBEC29-1430-4DA0-99DC-49B9389D2440}  
TERM.....INCISION DRESSING TYPE GAUZE WRAP

DATA VALUE: INCISION DRESSING TYPE HYDROCOLLOIDAL

-----  
ABBREVIATION.....INC DR HYD  
ACTIVE.....Yes  
DESCRIPTION.....Incision Dressing Type Hydrocolloidal  
DISPLAY\_NAME.....Hydrocolloidal  
HELP\_TEXT.....Incision Dressing Type Hydrocolloidal  
ID.....{424CCDCA-5196-4BBD-AE8B-0F68DDCC3FCF}  
TERM.....INCISION DRESSING TYPE HYDROCOLLOIDAL

DATA VALUE: INCISION DRESSING TYPE HYDROGEL

-----  
ABBREVIATION.....INC DR HYD  
ACTIVE.....Yes  
DESCRIPTION.....Incision Dressing Type Hydrogel  
DISPLAY\_NAME.....Hydrogel  
HELP\_TEXT.....Incision Dressing Type Hydrogel  
ID.....{8CBAC6DF-F51E-456B-8F9A-A68E4452812C}  
TERM.....INCISION DRESSING TYPE HYDROGEL

DATA VALUE: INCISION DRESSING TYPE HYPERTONIC HYDROGEL

-----  
ABBREVIATION.....INC DR HYP  
ACTIVE.....Yes  
DESCRIPTION.....Incision Dressing Type Hypertonic  
Hydrogel  
DISPLAY\_NAME.....Hypertonic Hydrogel  
HELP\_TEXT.....Incision Dressing Type Hypertonic  
Hydrogel  
ID.....{25820A63-AD76-4198-BEB1-EF95DDEC5B35}  
TERM.....INCISION DRESSING TYPE HYPERTONIC  
HYDROGEL

DATA VALUE: INCISION DRESSING TYPE NONE

-----  
ABBREVIATION.....INC DR NON  
ACTIVE.....Yes  
DESCRIPTION.....Incision Dressing Type None  
DISPLAY\_NAME.....None  
HELP\_TEXT.....Incision Dressing Type None  
ID.....{81F07C8D-A463-48A0-A64A-A143CECB5DD6}  
TERM.....INCISION DRESSING TYPE NONE

DATA VALUE: INCISION DRESSING TYPE PRESSURE

-----  
ABBREVIATION.....INC DR PRE  
ACTIVE.....Yes  
DESCRIPTION.....Incision Dressing Type Pressure  
DISPLAY\_NAME.....Pressure  
HELP\_TEXT.....Incision Dressing Type Pressure  
ID.....{B5823373-E5D8-40A5-954C-6A193D7E7511}  
TERM.....INCISION DRESSING TYPE PRESSURE

DATA VALUE: INCISION DRESSING TYPE TRANSPARENT

-----  
ABBREVIATION.....INC DR TRA  
ACTIVE.....Yes  
DESCRIPTION.....Incision Dressing Type Gauze  
DISPLAY\_NAME.....Transparent  
HELP\_TEXT.....Incision Dressing Type Gauze  
ID.....{E198EDAD-AA67-49F9-A065-6BB11FE39ED6}  
TERM.....INCISION DRESSING TYPE TRANSPARENT

DATA VALUE: INCISION EDGE INDISTINCT/ATTACHED

-----  
ABBREVIATION.....INC ED ATT  
ACTIVE.....Yes  
DESCRIPTION.....Incision Edge Indistinct, Attached  
DISPLAY\_NAME.....Indistinct, Attached  
HELP\_TEXT.....Incision Edge Indistinct, Attached  
ID.....{AF7E7150-525E-49A4-AC23-042B788E8805}  
TERM.....INCISION EDGE INDISTINCT/ATTACHED

DATA VALUE: INCISION EDGE NON-ATTACHED

-----  
ABBREVIATION.....INC ED NON  
ACTIVE.....Yes  
DESCRIPTION.....Incision Edge Non Attached  
DISPLAY\_NAME.....Non-Attached  
HELP\_TEXT.....Incision Edge Non Attached  
ID.....{F457F951-F1DC-47AC-9336-AC093E149862}  
TERM.....INCISION EDGE NON-ATTACHED

DATA VALUE: INCISION EDGE ROLLED

-----  
ABBREVIATION.....INC ED ROL  
ACTIVE.....Yes  
DESCRIPTION.....Incision Edge Rolled  
DISPLAY\_NAME.....Rolled  
HELP\_TEXT.....Incision Edge Rolled  
ID.....{FFA33CAE-5D8B-4501-8120-A536BE53296F}  
TERM.....INCISION EDGE ROLLED

DATA VALUE: INCISION EDGE TUNNELING

-----  
ABBREVIATION.....INC ED TUN  
ACTIVE.....Yes  
DESCRIPTION.....Incision Edge Tunneling  
DISPLAY\_NAME.....Tunneling  
HELP\_TEXT.....Incision Edge Tunneling  
ID.....{492CD783-3C32-40F0-AE35-5A89E90409A7}  
TERM.....INCISION EDGE TUNNELING

DATA VALUE: INCISION EDGE UNDERMINING

-----  
ABBREVIATION.....INC ED UND  
ACTIVE.....Yes  
DESCRIPTION.....Incision Edge Undermining  
DISPLAY\_NAME.....Undermining  
HELP\_TEXT.....Incision Edge Undermining  
ID.....{0F3E968A-B867-4123-8C05-ED6B12D364CD}  
TERM.....INCISION EDGE UNDERMINING

DATA VALUE: INCISION ODOR FOUL

-----  
ABBREVIATION.....INC OD FOU  
ACTIVE.....Yes  
DESCRIPTION.....Incision Odor Foul  
DISPLAY\_NAME.....Foul  
HELP\_TEXT.....Incision Odor Foul  
ID.....{ECDC0988-BE82-41E3-8104-A66A5F8A1093}  
TERM.....INCISION ODOR FOUL

DATA VALUE: INCISION ODOR NONE

-----  
ABBREVIATION.....INC OD NON  
ACTIVE.....Yes  
DESCRIPTION.....Incision Odor None  
DISPLAY\_NAME.....None  
HELP\_TEXT.....Incision Odor None  
ID.....{CED90890-6FBC-4A4A-B0B8-6F75F0DD567A}  
TERM.....INCISION ODOR NONE

DATA VALUE: INCISION WOUND CLOSURE APPROXIMATED

-----  
ABBREVIATION.....INC CL APP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Approximated  
HELP\_TEXT.....No Help-Text  
ID.....{2DD5713E-D81C-484B-A0DF-6A5730B396DF}  
TERM.....INCISION WOUND CLOSURE APPROXIMATED

DATA VALUE: INCISION WOUND CLOSURE DEHISCENCE

-----  
ABBREVIATION.....INC CL DEH  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Dehiscence  
HELP\_TEXT.....No Help-Text  
ID.....{DCA1D137-0D9A-42D2-9270-7B3A82DCB6EB}  
TERM.....INCISION WOUND CLOSURE DEHISCENCE

DATA VALUE: INCISION WOUND CLOSURE GLUE

-----  
ABBREVIATION.....IN CL GLU  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Glue  
HELP\_TEXT.....No Help-Text  
ID.....{9EEB436E-2991-4E9A-9F8A-BFEC196784CB}  
TERM.....INCISION WOUND CLOSURE GLUE

DATA VALUE: INCISION WOUND CLOSURE NOT CLOSED

-----  
ABBREVIATION.....INC CL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Not Closed  
HELP\_TEXT.....No Help-Text  
ID.....{D63A05BF-1205-45F6-A78D-15D198BCB1C6}  
TERM.....INCISION WOUND CLOSURE NOT CLOSED

DATA VALUE: INCISIONAL SITE

-----  
ABBREVIATION.....INCIS SITE  
ACTIVE.....Yes  
DESCRIPTION.....INCISIONAL SITE  
DISPLAY\_NAME.....Incisional Site  
HELP\_TEXT.....INCISIONAL SITE  
ID.....{19E69DE2-097B-0FA5-449D-B2E9161A3E05}  
TERM.....INCISIONAL SITE

DATA VALUE: INCOMPREHENSIBLE

-----  
ABBREVIATION.....INCOMPREHE  
ACTIVE.....Yes  
DESCRIPTION.....Incomprehensible  
DISPLAY\_NAME.....Incomprehensible  
HELP\_TEXT.....No Help Text  
ID.....{916CD781-563D-4F54-8752-800D37BFAA3B}  
TERM.....INCOMPREHENSIBLE

DATA VALUE: INCONTENENCE DERMATITIS

-----  
ABBREVIATION.....INCON DERM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Incontenance Dermatitis  
HELP\_TEXT.....No Help Text  
ID.....{9F6C5BA2-79C7-4B09-B19A-6DDC444936E5}  
TERM.....INCONTENENCE DERMATITIS



DATA VALUE: INCONTINENT

-----  
ABBREVIATION.....INCONTINEN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Incontinent  
HELP\_TEXT.....No Help Text  
ID.....{BA2E9C32-8D03-8FC9-5E02-EAE2026FDC63}  
TERM.....INCONTINENT

DATA VALUE: INCREASE ROUNDS

-----  
ABBREVIATION.....Inc Rds  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Increase Rounds  
HELP\_TEXT.....  
ID.....{873180AF-9BB9-44F4-B9E4-80F8BBD9769}  
TERM.....INCREASE ROUNDS

DATA VALUE: INDEPENDENT

-----  
ABBREVIATION.....INDEPENDEN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Independent  
HELP\_TEXT.....No Help-Text  
ID.....{11BC0CF7-1C29-49C7-B7C8-2BCF8FF82250}  
TERM.....INDEPENDENT

DATA VALUE: INDICATES UNDERSTANDING

-----  
ABBREVIATION.....UNDERSTAND  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Indicates Understanding  
HELP\_TEXT.....No Help Text  
ID.....{C4AEFE5F-C388-4E00-9663-6B308934284B}  
TERM.....INDICATES UNDERSTANDING

DATA VALUE: INDIFFERENT

-----  
ABBREVIATION.....INDIFF  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Indifferent  
HELP\_TEXT.....No Help Text  
ID.....{7C3B96CF-C3D0-4117-A952-94BD9FBE05A5}  
TERM.....INDIFFERENT

DATA VALUE: INDURATION

-----  
ABBREVIATION.....INDURAITON  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Induration (Hardness)  
HELP\_TEXT.....No Help Text  
ID.....{A1669576-45D9-435C-87B9-488E9CF6BBA4}  
TERM.....INDURATION

DATA VALUE: INDWELLING CATHETER

-----  
ABBREVIATION.....INDW CATH  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Indwelling Catheter  
HELP\_TEXT.....No Help Text  
ID.....{EC16724C-D1EA-5650-58D4-38A89D4456F9}  
TERM.....INDWELLING CATHETER

DATA VALUE: INEFFECTIVE COPING

-----  
ABBREVIATION.....ic  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Ineffective Coping  
HELP\_TEXT.....No Help Text  
ID.....{704B1A55-B639-60F4-E3A5-16F942ABC259}  
TERM.....INEFFECTIVE COPING

DATA VALUE: INFECTION

-----  
ABBREVIATION.....INFECTION  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Infection  
HELP\_TEXT.....No Help Text  
ID.....{DA00DB58-640D-4294-9468-EEF070DBAE8E}  
TERM.....INFECTION

DATA VALUE: INFECTION PREVENTION AND SCREENING

-----  
ABBREVIATION.....INF PRE SC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Infection Prevention and Screening  
HELP\_TEXT.....No Help Text  
ID.....{BA58329F-59EC-49F0-BACE-75AA5EC57F52}  
TERM.....INFECTION PREVENTION AND SCREENING

DATA VALUE: INFLATABLE HEEL LIFT

-----  
ABBREVIATION.....INF HEEL L  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Inflatable Heel Lift  
HELP\_TEXT.....No Help Text  
ID.....{793B34A4-2F21-42DA-881D-1E6E0B4EAD6B}  
TERM.....INFLATABLE HEEL LIFT

DATA VALUE: INITIATED

-----  
ABBREVIATION.....INITIATED  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Initiated  
HELP\_TEXT.....No Help Text  
ID.....{562E5FFB-4E53-4614-ACAC-E7D9B6E1B781}  
TERM.....INITIATED

DATA VALUE: INSERTED THROUGH EXISTING INTRODUCER

-----  
ABBREVIATION.....INS EXS IN  
ACTIVE.....Yes  
DESCRIPTION.....Inserted Through Existing Introducer  
DISPLAY\_NAME.....Inserted Through Existing Introducer  
HELP\_TEXT.....No Help Text  
ID.....{20BC183A-197F-4464-A27C-47F63C5B842A}  
TERM.....INSERTED THROUGH EXISTING INTRODUCER

DATA VALUE: INSERTED THROUGH NEW INTRODUCER

-----  
ABBREVIATION.....INS NEW IN  
ACTIVE.....Yes  
DESCRIPTION.....Inserted Through New Introducer  
DISPLAY\_NAME.....Inserted Through New Introducer  
HELP\_TEXT.....No Help Text  
ID.....{7FCA403D-41E5-4BAA-B155-B58F35EA49B3}  
TERM.....INSERTED THROUGH NEW INTRODUCER

DATA VALUE: INSPIRATORY WHEEZES

-----  
ABBREVIATION.....INSP WHEEZ  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Inspiratory Wheezes  
HELP\_TEXT.....No Help Text  
ID.....{3FA7F54C-DE17-4DFF-B5F6-C2CD05EC0B24}  
TERM.....INSPIRATORY WHEEZES

DATA VALUE: INTACT

-----  
ABBREVIATION.....i  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Intact  
HELP\_TEXT.....No apparent breakdown  
ID.....{935C9CAF-6C88-6610-C6E9-4C1B822923AF}  
TERM.....INTACT

DATA VALUE: INTENTIONAL

-----  
ABBREVIATION.....i  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Intentional  
HELP\_TEXT.....No Help Text  
ID.....{D5BBBFFC-F739-1E69-FED4-288D8E803D04}  
TERM.....INTENTIONAL

DATA VALUE: INTERMITTENT

-----  
ABBREVIATION.....INTERMITTE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Intermittent  
HELP\_TEXT.....No Help Text  
ID.....{DD09ED7C-B977-412F-A678-5130F77D239E}  
TERM.....INTERMITTENT

DATA VALUE: INTERMITTENT SLEEP

-----  
ABBREVIATION.....INT SLEEP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....States Has Slept Intermittently  
HELP\_TEXT.....No Help Text  
ID.....{5005C433-5892-4F4D-A64D-4C94AD1853E2}  
TERM.....INTERMITTENT SLEEP

DATA VALUE: INTERRUPTED

-----  
ABBREVIATION.....INTERRUPTED  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Interrupted  
HELP\_TEXT.....No Help Text  
ID.....{AB53B9AB-DA3F-42D3-BEB2-E02A452ED777}  
TERM.....INTERRUPTED

DATA VALUE: INTERVENTIONAL RADIOLOGY

-----  
ABBREVIATION.....Interv Rad  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Interventional Radiology  
HELP\_TEXT.....  
ID.....{322B66E6-814F-43AF-B669-E503519DB72F}  
TERM.....INTERVENTIONAL RADIOLOGY

DATA VALUE: INTRODUCER

-----  
ABBREVIATION.....INTRODUCER  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Introducer  
HELP\_TEXT.....No Help Text  
ID.....{B46F7CA4-C94E-68FD-2892-52866763B5A0}  
TERM.....INTRODUCER

DATA VALUE: IRREGULAR

-----  
ABBREVIATION.....IRREGULAR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Irregular  
HELP\_TEXT.....No Help Text  
ID.....{6CCE27A2-F1EB-3EC8-B140-CDA241D159F2}  
TERM.....IRREGULAR

DATA VALUE: IRRITABLE

-----  
ABBREVIATION.....Irritable  
ACTIVE.....Yes  
DESCRIPTION.....Irritable  
DISPLAY\_NAME.....Irritable  
HELP\_TEXT.....No Help Text  
ID.....{10651C11-2736-4056-BA84-E544358BA877}  
TERM.....IRRITABLE

DATA VALUE: ISOSOURCE VHN

-----  
ABBREVIATION.....Iso VHN  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Isosource VHN  
HELP\_TEXT.....  
ID.....{D68525BF-F922-464E-97DE-1AB767974483}  
TERM.....ISOSOURCE VHN

DATA VALUE: ITCHING

-----  
ABBREVIATION.....Itching  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Itching  
HELP\_TEXT.....  
ID.....{3A1FB47B-F274-458D-BDF6-E6882D6CB600}  
TERM.....ITCHING

DATA VALUE: IV DISCONTINUED

-----  
ABBREVIATION.....IV DISC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....IV Discontinued  
HELP\_TEXT.....No Help Text  
ID.....{FDF951E7-5B72-4166-969B-073D5EC74800}  
TERM.....IV DISCONTINUED

DATA VALUE: IV PUMP

-----  
ABBREVIATION.....IV PUMP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....IV Pump  
HELP\_TEXT.....No Help Text  
ID.....{EFB92C16-CFA3-18BB-9D54-590089DA3962}  
TERM.....IV PUMP

DATA VALUE: IV SEDATION

-----  
ABBREVIATION.....IV Sed  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....IV Sedation  
HELP\_TEXT.....  
ID.....{789B61EB-5E83-4024-850D-D92370C4EC9A}  
TERM.....IV SEDATION

DATA VALUE: IVPB

-----  
ABBREVIATION.....IVPB  
ACTIVE.....Yes  
DESCRIPTION.....IV Piggy Back  
DISPLAY\_NAME.....IV Piggy Back  
HELP\_TEXT.....No Help Text  
ID.....{50EE76CF-5A9C-4EEB-9201-37CD90509355}  
TERM.....IVPB

DATA VALUE: JAUNDICE

-----  
ABBREVIATION.....j  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Jaundice  
HELP\_TEXT.....No Help Text  
ID.....{B17F4603-7C68-660D-96A7-119140A998D3}  
TERM.....JAUNDICE

DATA VALUE: JEJUNOSTOMY

-----  
ABBREVIATION.....JEJUNOSTOM  
ACTIVE.....Yes  
DESCRIPTION.....Jejunostomy EF Tube  
DISPLAY\_NAME.....Jejunostomy  
HELP\_TEXT.....No Help Text  
ID.....{1ED3EC41-3ACF-42C1-B64A-280300EDBCFC}  
TERM.....JEJUNOSTOMY

DATA VALUE: JEVITY 1.0

-----  
ABBREVIATION.....Jevity 1.0  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Jevity 1.0  
HELP\_TEXT.....  
ID.....{637F4523-0417-4EA4-9991-379C8362EC00}  
TERM.....JEVITY 1.0

DATA VALUE: JEVITY 1.2

-----  
ABBREVIATION.....Jevity 1.2  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Jevity 1.2  
HELP\_TEXT.....  
ID.....{0B3ADECE-FC19-4D90-912B-6EDC3D35C845}  
TERM.....JEVITY 1.2

DATA VALUE: JEVITY 1.5

-----  
ABBREVIATION.....Jevity 1.5  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Jevity 1.5  
HELP\_TEXT.....  
ID.....{ED69D4F6-7EFA-4954-85F0-CCED2BA80A39}  
TERM.....JEVITY 1.5

DATA VALUE: JOINT CAPSULE

-----  
ABBREVIATION.....JOINT CAPS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Joint Capsule  
HELP\_TEXT.....No Help Text  
ID.....{F10E4DA8-DF84-4625-A317-724E7C9E028C}  
TERM.....JOINT CAPSULE

DATA VALUE: JP DRAIN

-----  
ABBREVIATION.....JP Drain  
ACTIVE.....Yes  
DESCRIPTION.....Jackson Pratt Drain  
DISPLAY\_NAME.....JP Drain  
HELP\_TEXT.....No Help Text  
ID.....{20F3B3D5-D3A4-96FD-5446-253FD031DF2B}  
TERM.....JP DRAIN

DATA VALUE: JUNCTIONAL

-----  
ABBREVIATION.....j  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Junctional  
HELP\_TEXT.....No Help Text  
ID.....{4A339D84-7051-1258-877D-85AFC1CA20FE}  
TERM.....JUNCTIONAL

DATA VALUE: JUNCTIONAL ESCAPE BEAT

-----  
ABBREVIATION.....Junctional  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Junctional Escape Beat  
HELP\_TEXT.....  
ID.....{54F8E569-801E-4512-A9BD-C04CDA5CEE7C}  
TERM.....JUNCTIONAL ESCAPE BEAT

DATA VALUE: JUNCTIONAL TACHYCARDIA

-----  
ABBREVIATION.....jt  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Junctional Tachycardia  
HELP\_TEXT.....  
ID.....{D6EF01C1-A8F9-40FD-AE24-785041526041}  
TERM.....JUNCTIONAL TACHYCARDIA



DATA VALUE: JUVEN

-----  
ABBREVIATION.....Juven  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Juven  
HELP\_TEXT.....  
ID.....{2112F76A-D581-4C55-94CA-4B4CC68813AD}  
TERM.....JUVEN

DATA VALUE: K-PAD

-----  
ABBREVIATION.....K-PAD  
ACTIVE.....Yes  
DESCRIPTION.....K-pad  
DISPLAY\_NAME.....K-pad  
HELP\_TEXT.....K-pad  
ID.....{10A0C893-F279-48ED-ACAA-1DCA4AD9DB24}  
TERM.....K-PAD

DATA VALUE: KNIFE-LIKE

-----  
ABBREVIATION.....KNIFE-LIKE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Knife-Like  
HELP\_TEXT.....No Help Text  
ID.....{FBA0EE3F-134B-431B-879A-9B17D65CD464}  
TERM.....KNIFE-LIKE

DATA VALUE: KNOWN/SUSPECTED OBSTRUCTIVE SLEEP APNEA

-----  
ABBREVIATION.....Known/Susp  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Known/Suspected Obstructive Sleep Apnea  
HELP\_TEXT.....  
ID.....{BD47DD8C-7057-403B-B393-EBCBF9907B68}  
TERM.....KNOWN/SUSPECTED OBSTRUCTIVE SLEEP APNEA

DATA VALUE: KUSSMAUL

-----  
ABBREVIATION.....KUSSMAUL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Kussmaul  
HELP\_TEXT.....No Help Text  
ID.....{0C5ED066-F445-ABFE-F7B1-D158EC7C2976}  
TERM.....KUSSMAUL

DATA VALUE: L BUNDLE BRANCH BLOCK

-----  
ABBREVIATION.....lbbb  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Bundle Branch Block  
HELP\_TEXT.....No Help Text  
ID.....{C3F2E19E-495B-7BF6-99EF-B8236D6B9C3D}  
TERM.....L BUNDLE BRANCH BLOCK

DATA VALUE: L FOOT DOWN

-----  
ABBREVIATION.....L FOOT DOW  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Foot Down  
HELP\_TEXT.....No Help Text  
ID.....{5B18895B-DAD5-4874-87DB-DACFFF88B170}  
TERM.....L FOOT DOWN

DATA VALUE: L FOOT UP

-----  
ABBREVIATION.....L FOOT UP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Foot Up  
HELP\_TEXT.....No Help Text  
ID.....{A15CD8C5-5CD1-4B90-B79C-5D3D331C7466}  
TERM.....L FOOT UP

DATA VALUE: L HEAD DOWN

-----  
ABBREVIATION.....L HEAD DOW  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Head Down  
HELP\_TEXT.....No Help Text  
ID.....{3C919CB8-1C0A-498D-95A6-C46D25744CEC}  
TERM.....L HEAD DOWN

DATA VALUE: L HEAD UP

-----  
ABBREVIATION.....L HEAD UP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Head Up  
HELP\_TEXT.....No Help Text  
ID.....{BB39792D-B572-4A10-BA81-8685809E7437}  
TERM.....L HEAD UP

DATA VALUE: L NECK

-----  
ABBREVIATION.....LEFT NECK  
ACTIVE.....Yes  
DESCRIPTION.....Left Neck  
DISPLAY\_NAME.....L Neck  
HELP\_TEXT.....Left Neck  
ID.....{6A53228A-5BD1-6AAF-029F-701C7EFE7C48}  
TERM.....L NECK

DATA VALUE: L NEPHROSTOMY

-----  
ABBREVIATION.....L NEPHROST  
ACTIVE.....Yes  
DESCRIPTION.....Left Nephrostomy  
DISPLAY\_NAME.....L Nephrostomy  
HELP\_TEXT.....Left Nephrostomy  
ID.....{D8FC4C18-081E-A7B9-E394-25F194BE9757}  
TERM.....L NEPHROSTOMY

DATA VALUE: L STENT

-----  
ABBREVIATION.....LEFT STENT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Stent  
HELP\_TEXT.....No Help Text  
ID.....{FD917F0B-CC96-0CB0-BED3-1D7C51A77B69}  
TERM.....L STENT

DATA VALUE: L1

-----  
ABBREVIATION.....L1  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....L1  
HELP\_TEXT.....  
ID.....{E6FA1F4F-B32B-4523-966C-9AFAD9C44316}  
TERM.....L1

DATA VALUE: L2

-----  
ABBREVIATION.....L2  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....L2  
HELP\_TEXT.....  
ID.....{0E5E7F8E-BA82-499B-A57A-8D32C2113006}  
TERM.....L2

DATA VALUE: L3

-----  
ABBREVIATION.....L3  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....L3  
HELP\_TEXT.....  
ID.....{0CBAF5CD-AABA-4FF0-A25D-01D20CDD53C3}  
TERM.....L3

DATA VALUE: L4

-----  
ABBREVIATION.....L4  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....L4  
HELP\_TEXT.....  
ID.....{94632F85-0152-47C2-9D82-B49C36ABF495}  
TERM.....L4

DATA VALUE: L5

-----  
ABBREVIATION.....L5  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....L5  
HELP\_TEXT.....  
ID.....{87E93C0D-B607-4DD2-83EC-71131C072019}  
TERM.....L5

DATA VALUE: LACERATION

-----  
ABBREVIATION.....LACERATION  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Laceration  
HELP\_TEXT.....No Help Text  
ID.....{F4548BB2-39AF-4E0D-B52F-472EA080D6FF}  
TERM.....LACERATION

DATA VALUE: LACKS INTEREST

-----  
ABBREVIATION.....LACKS INT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Lacks Interest  
HELP\_TEXT.....No Help Text  
ID.....{585C5B7D-D229-4C7C-AE3C-5869FF845F0D}  
TERM.....LACKS INTEREST

DATA VALUE: LACTOSE RESTRICTED

-----  
ABBREVIATION.....LACT RESTR  
ACTIVE.....Yes  
DESCRIPTION.....Lactose Restricted  
DISPLAY\_NAME.....Lactose Restricted  
HELP\_TEXT.....No Help Text  
ID.....{78A9066A-48E2-4285-88E3-A180CE771B5A}  
TERM.....LACTOSE RESTRICTED

DATA VALUE: LANGUAGE

-----  
ABBREVIATION.....LANGUAGE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Language  
HELP\_TEXT.....No Help-Text  
ID.....{0666538E-4F44-4D34-A5EA-69CF2F1F86FA}  
TERM.....LANGUAGE

DATA VALUE: LARGE

-----  
ABBREVIATION.....L  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Large  
HELP\_TEXT.....No Help Text  
ID.....{E60D5CF9-920D-44E2-8A83-136DAB1FE251}  
TERM.....LARGE

DATA VALUE: LARYNGECTOMY

-----  
ABBREVIATION.....LARYNGECTO  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Laryngectomy  
HELP\_TEXT.....No Help Text  
ID.....{225CB76A-BAF5-827E-729B-C3D60AD5B50A}  
TERM.....LARYNGECTOMY

DATA VALUE: LEFT

-----  
ABBREVIATION.....LEFT  
ACTIVE.....Yes  
DESCRIPTION.....Left  
DISPLAY\_NAME.....Left  
HELP\_TEXT.....No Help Text  
ID.....{81063CEE-2E83-4C95-AAEE-76B7DCB828D4}  
TERM.....LEFT

DATA VALUE: LESS THAN 30

-----  
ABBREVIATION.....<30  
ACTIVE.....Yes  
DESCRIPTION.....Less Than 30  
DISPLAY\_NAME.....<30  
HELP\_TEXT.....No Help Text  
ID.....{5AA98A0D-9B05-425D-97CA-B8B12425AACA}  
TERM.....LESS THAN 30

DATA VALUE: LETHARGIC

-----  
ABBREVIATION.....L  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Lethargic  
HELP\_TEXT.....No Help Text  
ID.....{4ECA2E44-093A-9354-848E-CDC72B165072}  
TERM.....LETHARGIC

DATA VALUE: LEVELED

-----  
ABBREVIATION.....LEVELED  
ACTIVE.....Yes  
DESCRIPTION.....Leveled  
DISPLAY\_NAME.....Leveled  
HELP\_TEXT.....No Help Text  
ID.....{D72523B2-E0CD-40E1-AC14-2FA43636A456}  
TERM.....LEVELED

DATA VALUE: LIBERAL DIABETIC

-----  
ABBREVIATION.....LIBRL DIAB  
ACTIVE.....Yes  
DESCRIPTION.....Liberal Diabetic  
DISPLAY\_NAME.....Liberal Diabetic  
HELP\_TEXT.....No Help Text  
ID.....{9B66B47B-E77E-4F69-BDFA-5323B28A283D}  
TERM.....LIBERAL DIABETIC

DATA VALUE: LIMITED ASSISTANCE

-----  
ABBREVIATION.....LIM ASSIST  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Limited Assistance  
HELP\_TEXT.....No Help Text  
ID.....{9A8C1310-C840-431E-BA21-065B7EA7B768}  
TERM.....LIMITED ASSISTANCE

DATA VALUE: LIMITED ATTENION SPAN

-----  
ABBREVIATION.....LIM ATT S  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Limited Attention Span  
HELP\_TEXT.....No Help Text  
ID.....{51553020-0CB9-483E-AC1C-15CA6D9380A3}  
TERM.....LIMITED ATTENION SPAN

DATA VALUE: LIMITED TO INCISION

-----  
ABBREVIATION.....ER EXI LIM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Limited to Incision  
HELP\_TEXT.....No Help-Text  
ID.....{C9ED5E2C-8ED1-445A-8A3E-C1207CE3843D}  
TERM.....LIMITED TO INCISION

DATA VALUE: LINE FLUSH

-----  
ABBREVIATION.....LINE FLUSH  
ACTIVE.....Yes  
DESCRIPTION.....Fluid flushed from line  
DISPLAY\_NAME.....Line Flush  
HELP\_TEXT.....No Help Text  
ID.....{8DDD63BE-CF28-4C83-93FF-C81609FFDE51}  
TERM.....LINE FLUSH

DATA VALUE: LIP SALVE

-----  
ABBREVIATION.....LIP SALVE  
ACTIVE.....Yes  
DESCRIPTION.....Lip Salve  
DISPLAY\_NAME.....Lip Salve  
HELP\_TEXT.....No Help Text  
ID.....{612F78AB-7386-4176-9896-02DE5A8E836D}  
TERM.....LIP SALVE

DATA VALUE: LISTED DIET

-----  
ABBREVIATION.....LIST DIET  
ACTIVE.....Yes  
DESCRIPTION.....Listed Diet  
DISPLAY\_NAME.....Listed Diet  
HELP\_TEXT.....No Help Text  
ID.....{24FBBFD9-B426-4DD9-9843-3428A1C91676}  
TERM.....LISTED DIET

DATA VALUE: LOCAL

-----  
ABBREVIATION.....Local  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Local  
HELP\_TEXT.....  
ID.....{64A0FD20-4584-446C-8A2A-57F2FEF3DD10}  
TERM.....LOCAL

DATA VALUE: LOCALIZES TO PAIN

-----  
ABBREVIATION.....LOC TO PAI  
ACTIVE.....Yes  
DESCRIPTION.....Localizes To Pain  
DISPLAY\_NAME.....Localizes To Pain  
HELP\_TEXT.....No Help Text  
ID.....{4072DF9-2D90-423F-B6F0-0F2B924025D9}  
TERM.....LOCALIZES TO PAIN

DATA VALUE: LOCK BOX

-----  
ABBREVIATION.....Lock Box  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Lock Box  
HELP\_TEXT.....  
ID.....{974A2F44-C918-4E85-ABAF-03A21498F98B}  
TERM.....LOCK BOX

DATA VALUE: LOCKER

-----  
ABBREVIATION.....Locker  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Locker  
HELP\_TEXT.....  
ID.....{36AE4F3F-E2BC-4DD3-8DD0-95BD0E38B1B9}  
TERM.....LOCKER

DATA VALUE: LOOSE

-----  
ABBREVIATION.....LOOSE  
ACTIVE.....Yes  
DESCRIPTION.....Loose  
DISPLAY\_NAME.....Loose  
HELP\_TEXT.....No Help Text  
ID.....{1346A74E-3AA1-4D81-876B-91201F9717A5}  
TERM.....LOOSE



DATA VALUE: LOW

-----  
ABBREVIATION.....I  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Low  
HELP\_TEXT.....No Help Text  
ID.....{B6CE52A4-C1DD-FDE3-F87E-8E0C98FCC4AC}  
TERM.....LOW

DATA VALUE: LOW AIRLOSS OVERLAY

-----  
ABBREVIATION.....L AIR OVRL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Low Airloss Overlay  
HELP\_TEXT.....No Help Text  
ID.....{282C439C-A697-43E6-85C1-786352EB56BC}  
TERM.....LOW AIRLOSS OVERLAY

DATA VALUE: LOW AIRLOSS/ALTERNATING PRESSURE COMBO

-----  
ABBREVIATION.....L AIRL ALT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Low Airloss/Alternating Pressure Combo  
HELP\_TEXT.....No Help Text  
ID.....{ADDD1260-C95B-4D7B-AB89-3749D76A3700}  
TERM.....LOW AIRLOSS/ALTERNATING PRESSURE COMBO

DATA VALUE: LOW CHOLESTEROL

-----  
ABBREVIATION.....LOW CHOLE  
ACTIVE.....Yes  
DESCRIPTION.....Low Cholesterol  
DISPLAY\_NAME.....Low Cholesterol  
HELP\_TEXT.....No Help Text  
ID.....{88F0170A-2560-48C6-9604-982453BCFC3B}  
TERM.....LOW CHOLESTEROL

DATA VALUE: LOW FIBER

-----  
ABBREVIATION.....LOW FIBER  
ACTIVE.....Yes  
DESCRIPTION.....Low Fiber  
DISPLAY\_NAME.....Low Fiber  
HELP\_TEXT.....No Help Text  
ID.....{2834D9A0-B6F7-4C21-96ED-1721D32B3B2B}  
TERM.....LOW FIBER

DATA VALUE: LOW PROTEIN 60 GM

-----  
ABBREVIATION.....L PROT 60g  
ACTIVE.....Yes  
DESCRIPTION.....Low Protein 60 GM  
DISPLAY\_NAME.....Low Protein 60 GM  
HELP\_TEXT.....No Help Text  
ID.....{BF5B6D8B-01A4-42D5-8551-6E78A17C43F7}  
TERM.....LOW PROTEIN 60 GM

DATA VALUE: LOW TYRAMINE

-----  
ABBREVIATION.....LOW TYRAMI  
ACTIVE.....Yes  
DESCRIPTION.....Low Tyramine  
DISPLAY\_NAME.....Low Tyramine  
HELP\_TEXT.....No Help Text  
ID.....{DB067A83-9D0C-4FC2-BAC8-6EC5641CAE42}  
TERM.....LOW TYRAMINE

DATA VALUE: LP: LIGHT PERCEPTION

-----  
ABBREVIATION.....LP  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....LP: Light Perception  
HELP\_TEXT.....  
ID.....{3B5C4FA1-6770-4347-B5FB-2B0554D71339}  
TERM.....LP: LIGHT PERCEPTION

DATA VALUE: LPP: LIGHT PERCEPTION WITH PROJECTION

-----  
ABBREVIATION.....LPP  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....LPP: Light Perception with Projection  
HELP\_TEXT.....  
ID.....{29ED35AA-49C5-4835-AF43-17173333A2BB}  
TERM.....LPP: LIGHT PERCEPTION WITH PROJECTION

DATA VALUE: LUMBAR DRAIN

-----  
ABBREVIATION.....LUM DRAIN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Lumbar Drain  
HELP\_TEXT.....No Help Text  
ID.....{FF36A046-1799-0734-286B-2272C34E3056}  
TERM.....LUMBAR DRAIN

DATA VALUE: LUMBAR SPINAL PRECAUTIONS

-----  
ABBREVIATION.....LUMB S P  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Lumbar Spinal Precautions  
HELP\_TEXT.....No Help Text  
ID.....{5AFBBBAB-F919-436A-9C5A-71DCF5D79D18}  
TERM.....LUMBAR SPINAL PRECAUTIONS

DATA VALUE: LVN

-----  
ABBREVIATION.....LVN  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....LVN  
HELP\_TEXT.....  
ID.....{8F8300FD-5A33-4B01-9F0B-901E2A3D96AB}  
TERM.....LVN

DATA VALUE: LVN/LPN

-----  
ABBREVIATION.....LVN/LPN  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....LVN/LPN  
HELP\_TEXT.....  
ID.....{DC89CD99-80C5-4455-A280-3A19E74DCE7D}  
TERM.....LVN/LPN

DATA VALUE: LYING OR SITTING

-----  
ABBREVIATION.....los  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Lying Or Sitting  
HELP\_TEXT.....No Help Text  
ID.....{52A940B1-55AF-AB4E-9E59-F93515A1296C}  
TERM.....LYING OR SITTING

DATA VALUE: MACERATION

-----  
ABBREVIATION.....MACERATION  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....MACERATION  
HELP\_TEXT.....No Help Text  
ID.....{B7C63E05-E193-40AA-A572-F1ACA0162B8B}  
TERM.....MACERATION

DATA VALUE: MAGNACAL

-----  
ABBREVIATION.....MAGNACAL  
ACTIVE.....Yes  
DESCRIPTION.....Enteral Nutrition Type  
DISPLAY\_NAME.....Magnacal  
HELP\_TEXT.....No Help Text  
ID.....{C810A618-FB1A-45BB-A540-EADF1414E343}  
TERM.....MAGNACAL

DATA VALUE: MAINTAINING EFFECTIVE VENTILATION

-----  
ABBREVIATION.....EFF VENT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Maintaining Effective Ventilation  
HELP\_TEXT.....No Help Text  
ID.....{18774550-ECBD-412D-AD95-131FAA4CC549}  
TERM.....MAINTAINING EFFECTIVE VENTILATION

DATA VALUE: MAROON

-----  
ABBREVIATION.....MAROON  
ACTIVE.....Yes  
DESCRIPTION.....Maroon in appearance  
DISPLAY\_NAME.....Maroon  
HELP\_TEXT.....No Help Text  
ID.....{1F518B23-C79F-4F06-9DF3-ADC5C96D3029}  
TERM.....MAROON

DATA VALUE: MD

-----  
ABBREVIATION.....MD  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....MD  
HELP\_TEXT.....  
ID.....{9AE2F01F-A8EF-4423-BD1C-55E9A1CD50A8}  
TERM.....MD

DATA VALUE: MEAL

-----  
ABBREVIATION.....MEAL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Meal  
HELP\_TEXT.....No Help Text  
ID.....{85D38DCD-6A4A-4908-ABE2-4E92ABEF1321}  
TERM.....MEAL

DATA VALUE: MECHANICAL

-----  
ABBREVIATION.....MECHANICAL  
ACTIVE.....Yes  
DESCRIPTION.....Mechanical  
DISPLAY\_NAME.....Mechanical  
HELP\_TEXT.....No Help Text  
ID.....{01EA793C-B524-4685-A3EC-60A317CBACC5}  
TERM.....MECHANICAL

DATA VALUE: MED RECONCILLIATION DOCUMENTED

-----  
ABBREVIATION.....ME RE DO  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Med Reconciliation Documented  
HELP\_TEXT.....No Help Text  
ID.....{9ABF3E2E-4D63-43EF-A80D-9DCF5C820226}  
TERM.....MED RECONCILLIATION DOCUMENTED

DATA VALUE: MEDIAL PORT

-----  
ABBREVIATION.....MEDI PORT  
ACTIVE.....Yes  
DESCRIPTION.....Medial Port  
DISPLAY\_NAME.....Medial Port  
HELP\_TEXT.....No Help Text  
ID.....{D1734D73-2E1A-4538-8884-05D8A78CF556}  
TERM.....MEDIAL PORT

DATA VALUE: MEDICAL DEVICE SAFETY

-----  
ABBREVIATION.....mds  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Medical Device Safety  
HELP\_TEXT.....No Help Text  
ID.....{1985CABE-4F53-9FCF-F444-F9AC53C04DBD}  
TERM.....MEDICAL DEVICE SAFETY

DATA VALUE: MEDICAL NECESSITY CONFIRMED

-----  
ABBREVIATION.....med nec  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Medical Necessity Confirmed  
HELP\_TEXT.....  
ID.....{F380D431-196B-46C9-8A1C-C57E7ECBFC8A}  
TERM.....MEDICAL NECESSITY CONFIRMED

DATA VALUE: MEDICATION

-----  
ABBREVIATION.....MEDICATION  
ACTIVE.....Yes  
DESCRIPTION.....Medication  
DISPLAY\_NAME.....Medication  
HELP\_TEXT.....No Help Text  
ID.....{22E6DEBB-9662-466B-990C-8DA500C0A183}  
TERM.....MEDICATION

DATA VALUE: MEDICATION ROOM

-----  
ABBREVIATION.....mr  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Medication Room  
HELP\_TEXT.....No Help Text  
ID.....{9571AA9F-F538-55EC-37F6-99FE5677B970}  
TERM.....MEDICATION ROOM

DATA VALUE: MEDIUM

-----  
ABBREVIATION.....MEDIUM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Medium  
HELP\_TEXT.....No Help Text  
ID.....{62C94066-A7B6-CAA8-064B-2AF7C0DA6323}  
TERM.....MEDIUM

DATA VALUE: MEMORY

-----  
ABBREVIATION.....MEMORY  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Memory  
HELP\_TEXT.....No Help Text  
ID.....{E30C1930-B98B-4382-B892-34E9C6008D1F}  
TERM.....MEMORY

DATA VALUE: MENTAL HEALTH

-----  
ABBREVIATION.....MEN HEAL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Mental Health  
HELP\_TEXT.....No Help Text  
ID.....{ED55FD2C-E2EC-4458-9AE4-FD7F0E6BD7F6}  
TERM.....MENTAL HEALTH

DATA VALUE: MENTAL HEALTH HOLD

-----  
ABBREVIATION.....MEN HEA HO  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Mental Health Hold  
HELP\_TEXT.....No Help Text  
ID.....{31560285-7BD5-4E53-BB77-93B244367F42}  
TERM.....MENTAL HEALTH HOLD

DATA VALUE: MENTAL STATUS - AWAKE ORIENTED

-----  
ABBREVIATION.....M/S Awake  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Mental Status - Awake Oriented  
HELP\_TEXT.....  
ID.....{37818468-4798-4CBE-82ED-A793BDD10874}  
TERM.....MENTAL STATUS - AWAKE ORIENTED

DATA VALUE: MENTAL STATUS - BASELINE FOR KNOWN DEMENTIA

-----  
ABBREVIATION.....M/S Demen  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Mental Status - Baseline for Dementia  
HELP\_TEXT.....  
ID.....{8E0BBC85-6D9F-44BB-AEAF-41BF561650A9}  
TERM.....MENTAL STATUS - BASELINE FOR KNOWN  
DEMENTIA

DATA VALUE: METAL IMPLANTS

-----  
ABBREVIATION.....Metal Impl  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Metal Implants  
HELP\_TEXT.....  
ID.....{E7424C97-0A81-4C4E-916E-0AE812FFDB96}  
TERM.....METAL IMPLANTS

DATA VALUE: MICU

-----  
ABBREVIATION.....MICU  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....MICU  
HELP\_TEXT.....No Help Text  
ID.....{28244B62-644C-4DA0-A350-EAFE0D2FA7C0}  
TERM.....MICU

DATA VALUE: MILD-MODERATELY REDUCED

-----  
ABBREVIATION.....MI-MOD-RED  
ACTIVE.....Yes  
DESCRIPTION.....Mild-Moderately Reduced  
DISPLAY\_NAME.....Mild-Moderately Reduced  
HELP\_TEXT.....No Help Text  
ID.....{41DA2518-0E49-4E92-A733-3ACF3D2E8199}  
TERM.....MILD-MODERATELY REDUCED

DATA VALUE: MILLIGRAMS PER KILOGRAM PER HOUR

-----  
ABBREVIATION.....mg/kg/hr  
ACTIVE.....Yes  
DESCRIPTION.....Milligrams per Kilogram per Hour  
DISPLAY\_NAME.....mg/kg/hr  
HELP\_TEXT.....Milligrams per Kilogram per Hour  
ID.....{731B270D-C8C5-4910-B4B0-AC101555108A}  
TERM.....MILLIGRAMS PER KILOGRAM PER HOUR

DATA VALUE: MISSING

-----  
ABBREVIATION.....MISSING  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Missing  
HELP\_TEXT.....No Help Text  
ID.....{80019911-11DA-6569-751E-ECAE0D1BC7DF}  
TERM.....MISSING

DATA VALUE: MOANING

-----  
ABBREVIATION.....MOANING  
ACTIVE.....Yes  
DESCRIPTION.....Moaning  
DISPLAY\_NAME.....Moaning  
HELP\_TEXT.....No Help Text  
ID.....{A71E3515-AA20-4B0E-9CB8-1612F76BA406}  
TERM.....MOANING

DATA VALUE: MODERATE

-----  
ABBREVIATION.....MODERATE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Moderate  
HELP\_TEXT.....No Help Text  
ID.....{C3C25064-C9FC-EA1F-3E74-32A25F0FC5E3}  
TERM.....MODERATE



DATA VALUE: MODERATE SEDATION

-----  
ABBREVIATION.....Mod Sed  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Moderate Sedation  
HELP\_TEXT.....  
ID.....{C874EDC7-24BF-4F39-8FDA-E8D1193A8428}  
TERM.....MODERATE SEDATION

DATA VALUE: MODERATE-SEVERELY REDUCED

-----  
ABBREVIATION.....MO-SEV RED  
ACTIVE.....Yes  
DESCRIPTION.....Moderate-Severely Reduced  
DISPLAY\_NAME.....Moderate-Severely Reduced  
HELP\_TEXT.....No Help Text  
ID.....{6D031C70-2F95-4111-A9B8-BC22CEC1B17D}  
TERM.....MODERATE-SEVERELY REDUCED

DATA VALUE: MODERATELY REDUCED

-----  
ABBREVIATION.....MOD REDUCE  
ACTIVE.....Yes  
DESCRIPTION.....Moderately Reduced  
DISPLAY\_NAME.....Moderately Reduced  
HELP\_TEXT.....No Help Text  
ID.....{91F73300-7EE0-428C-A838-5B3BE062D8BC}  
TERM.....MODERATELY REDUCED

DATA VALUE: MOIST

-----  
ABBREVIATION.....m  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Moist  
HELP\_TEXT.....No Help Text  
ID.....{AA73A9D6-6645-5DA7-C602-A537A2DFBFA0}  
TERM.....MOIST

DATA VALUE: MOTORIZED WHEELCHAIR

-----  
ABBREVIATION.....MO WH CH  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Motorized Wheelchair  
HELP\_TEXT.....No Help Text  
ID.....{4CA9CF41-F52C-41A8-B8A7-100A56097FDE}  
TERM.....MOTORIZED WHEELCHAIR

DATA VALUE: MOTTLED

-----  
ABBREVIATION.....m  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Mottled  
HELP\_TEXT.....No Help Text  
ID.....{4B5FFB10-2CFF-2064-1553-4CBB7133B0C5}  
TERM.....MOTTLED

DATA VALUE: MOUTHWASH USED

-----  
ABBREVIATION.....MOUTHWASH  
ACTIVE.....Yes  
DESCRIPTION.....Mouthwash Used  
DISPLAY\_NAME.....Mouthwash Used  
HELP\_TEXT.....No Help Text  
ID.....{9F7EF355-12A4-4EBA-92ED-0F38273E411B}  
TERM.....MOUTHWASH USED

DATA VALUE: MOVE AGAINST GRAVITY AND RESISTANCE, EASILY OVERCOME BY THE EXAMINER

-----  
ABBREVIATION.....magareobte  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Movement Against Gravity And  
Resistance, Easily Overcome By The  
Examiner  
HELP\_TEXT.....No Help Text  
ID.....{2ECCC197-4A81-025F-3EF3-6214004887A9}  
TERM.....MOVE AGAINST GRAVITY AND RESISTANCE,  
EASILY OVERCOME BY THE EXAMINER

DATA VALUE: MOVE AGAINST GRAVITY BUT NO RESISTANCE

-----  
ABBREVIATION.....magbnr  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Movement Against Gravity But No  
Resistance  
HELP\_TEXT.....No Help Text  
ID.....{F995A165-8536-C493-2269-286E8EFF5E43}  
TERM.....MOVE AGAINST GRAVITY BUT NO RESISTANCE

DATA VALUE: MOVEMENT WITH GRAVITY ELIMINATED

-----  
ABBREVIATION.....mwge  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Movement With Gravity Eliminated  
HELP\_TEXT.....No Help Text  
ID.....{B86E1131-1976-FA8E-74FE-B1BC2F0FDE3F}  
TERM.....MOVEMENT WITH GRAVITY ELIMINATED

DATA VALUE: MUCOID

-----  
ABBREVIATION.....MUCOID  
ACTIVE.....Yes  
DESCRIPTION.....Mucoid in appearance  
DISPLAY\_NAME.....Mucoid  
HELP\_TEXT.....No Help Text  
ID.....{073E5514-DACE-4F43-ADE5-91A1A87D6601}  
TERM.....MUCOID

DATA VALUE: MULTIFOCAL PVC

-----  
ABBREVIATION.....Mult PVC  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Multifocal PVC  
HELP\_TEXT.....  
ID.....{34AAE3E0-C0AA-4262-90A2-FCAD14E8023A}  
TERM.....MULTIFOCAL PVC

DATA VALUE: MULTIPODOUS

-----  
ABBREVIATION.....MULTIPODOU  
ACTIVE.....Yes  
DESCRIPTION.....Multipodous  
DISPLAY\_NAME.....Multipodous  
HELP\_TEXT.....No Help Text  
ID.....{ED098F9D-B204-4310-813A-98D91FF9DD25}  
TERM.....MULTIPODOUS

DATA VALUE: MURMUR

-----  
ABBREVIATION.....MURMUR  
ACTIVE.....Yes  
DESCRIPTION.....MURMUR  
DISPLAY\_NAME.....Murmur  
HELP\_TEXT.....MURMUR  
ID.....{591058E0-4B71-E432-B7D7-EBD42A8C4225}  
TERM.....MURMUR

DATA VALUE: MUSHROOM CATHETER

-----  
ABBREVIATION.....MUSH CATHE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Mushroom Catheter  
HELP\_TEXT.....No Help Text  
ID.....{4FFC451E-14B5-46E4-99C9-99A17BF90AE4}  
TERM.....MUSHROOM CATHETER

DATA VALUE: MUTIFOCAL ATRIAL TACHYCARDIA

-----  
ABBREVIATION.....mat  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Multifocal Atrial Tachycardia  
HELP\_TEXT.....  
ID.....{562C6BBD-C210-4467-93B3-8E50640E7CCE}  
TERM.....MUTIFOCAL ATRIAL TACHYCARDIA

DATA VALUE: N/A

-----  
ABBREVIATION.....n  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....N/A  
HELP\_TEXT.....No Help Text  
ID.....{B3DBE3FD-C5CB-3894-607A-BE13319594D0}  
TERM.....N/A

DATA VALUE: NA

-----  
ABBREVIATION.....NA  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....NA  
HELP\_TEXT.....  
ID.....{1993F8D9-D536-4ECA-AFAD-366F79BA0CBD}  
TERM.....NA

DATA VALUE: NAIL CARE

-----  
ABBREVIATION.....NAIL CARE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Nail Care  
HELP\_TEXT.....No Help Text  
ID.....{B783A171-1F0A-4135-8535-85776B81AB68}  
TERM.....NAIL CARE

DATA VALUE: NAPS ONLY

-----  
ABBREVIATION.....NAPS ONLY  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Naps Only  
HELP\_TEXT.....No Help Text  
ID.....{CF749F5C-2E4B-419A-B020-C81B30A7E482}  
TERM.....NAPS ONLY

DATA VALUE: NASAL

-----  
ABBREVIATION.....NASAL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Nasal  
HELP\_TEXT.....No Help Text  
ID.....{0F73F0CB-E25E-4611-9625-62E6F1596A9B}  
TERM.....NASAL

DATA VALUE: NASAL ENDOTRACHEAL

-----  
ABBREVIATION.....Nasal Endo  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Nasal Endotracheal  
HELP\_TEXT.....No Help Text  
ID.....{C631352E-CAB5-0409-6133-E5D11ECC3293}  
TERM.....NASAL ENDOTRACHEAL

DATA VALUE: NASAL INTESTINAL SMALL BORE

-----  
ABBREVIATION.....N INT SM B  
ACTIVE.....Yes  
DESCRIPTION.....Nasal Intestinal Small Bore EF Tube  
DISPLAY\_NAME.....Nasal Intestinal Small Bore  
HELP\_TEXT.....No Help Text  
ID.....{03DCB9EB-9B11-4C3B-869B-B1842BBCC844}  
TERM.....NASAL INTESTINAL SMALL BORE

DATA VALUE: NASAL SMALL BOWEL FEEDING TUBE

-----  
ABBREVIATION.....NAS SM BOW  
ACTIVE.....Yes  
DESCRIPTION.....Nasal Small Bowel Feeding Tube  
DISPLAY\_NAME.....Nasal Small Bowel Feeding Tube  
HELP\_TEXT.....No Help Text  
ID.....{A2968C8A-DDBD-4C3B-AFE6-3E2FED41C413}  
TERM.....NASAL SMALL BOWEL FEEDING TUBE

DATA VALUE: NASAL TRUMPET

-----  
ABBREVIATION.....NASAL TRUM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Nasal Trumpet  
HELP\_TEXT.....No Help Text  
ID.....{2FFB8510-793C-D6A4-838E-20FF7CD26533}  
TERM.....NASAL TRUMPET

DATA VALUE: NASOGASTRIC LARGE BORE TUBE

-----  
ABBREVIATION.....NG LG BORE  
ACTIVE.....Yes  
DESCRIPTION.....Nasogastric Large Bore EF Tube  
DISPLAY\_NAME.....Nasogastric Large Bore  
HELP\_TEXT.....No Help Text  
ID.....{C3FEE16A-9998-455B-BEB8-CC55AF890202}  
TERM.....NASOGASTRIC LARGE BORE TUBE

DATA VALUE: NASOGASTRIC TUBE

-----  
ABBREVIATION.....NG TUBE  
ACTIVE.....Yes  
DESCRIPTION.....Nasogastric Tube  
DISPLAY\_NAME.....Nasogastric Tube  
HELP\_TEXT.....No Help Text  
ID.....{6475E343-A310-48A5-B543-D29573E7684D}  
TERM.....NASOGASTRIC TUBE

DATA VALUE: NASOGASTRIC TUBE SMALL BORE

-----  
ABBREVIATION.....NGT SMALL  
ACTIVE.....Yes  
DESCRIPTION.....NGT Small Bore Tube  
DISPLAY\_NAME.....Nasal Gastric Tube Small Bore  
HELP\_TEXT.....No Help Text  
ID.....{D9FDAD80-9BE2-4B96-87C6-1D25D26D2959}  
TERM.....NASOGASTRIC TUBE SMALL BORE

DATA VALUE: NAUSEA

-----  
ABBREVIATION.....NAUSEA  
ACTIVE.....Yes  
DESCRIPTION.....Nausea  
DISPLAY\_NAME.....Nausea  
HELP\_TEXT.....Nausea  
ID.....{6BF33DC3-0300-4643-BB14-27C35692B806}  
TERM.....NAUSEA

DATA VALUE: NEAT

-----  
ABBREVIATION.....NEAT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Neat  
HELP\_TEXT.....No Help Text  
ID.....{8B01F7B8-B626-4143-A092-1E732EEB4A05}  
TERM.....NEAT

DATA VALUE: NECROTIC

-----  
ABBREVIATION.....NECROTIC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Necrotic  
HELP\_TEXT.....No Help Text  
ID.....{86885A34-E4EC-4BB5-BD94-59DD4B0F1A0E}  
TERM.....NECROTIC

DATA VALUE: NEEDLELESS HUB DISINFECTED AND ACCESSED PER POLICY

-----  
ABBREVIATION.....Needleless  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Needleless Hub Disinfected and  
                  Accessed per Policy  
HELP\_TEXT.....  
ID.....{5D9D8D90-9B39-4DF3-8095-755235687F45}  
TERM.....NEEDLELESS HUB DISINFECTED AND  
                  ACCESSED PER POLICY

DATA VALUE: NEEDS REINFORCEMENT

-----  
ABBREVIATION.....NEED REINF  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Needs Reinforcement  
HELP\_TEXT.....No Help Text  
ID.....{B6171F65-4530-418D-A63B-65DFFF522786}  
TERM.....NEEDS REINFORCEMENT

DATA VALUE: NEGATIVE PRESSURE

-----  
ABBREVIATION.....NEG PRESSU  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Negative Pressure  
HELP\_TEXT.....No Help Text  
ID.....{6D6CB4BB-96B0-4978-97CC-116A03479DAD}  
TERM.....NEGATIVE PRESSURE

DATA VALUE: NEGATIVE PRESSURE THERAPY DEVICE

-----  
ABBREVIATION.....NEG PRE TH  
ACTIVE.....Yes  
DESCRIPTION.....Negative Pressure Therapy Device  
DISPLAY\_NAME.....Negative Pressure Therapy Device  
HELP\_TEXT.....No Help Text  
ID.....{8AF0A069-66AD-4114-99D8-238AD987F0FF}  
TERM.....NEGATIVE PRESSURE THERAPY DEVICE

DATA VALUE: NEPRO

-----  
ABBREVIATION.....Nepro  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Nepro  
HELP\_TEXT.....  
ID.....{00FFC7A1-5CDA-45AF-BD50-D77666A23851}  
TERM.....NEPRO

DATA VALUE: NEPRO CARB STEADY

-----  
ABBREVIATION.....Nepro Car  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Nepro Carb Steady  
HELP\_TEXT.....  
ID.....{6B2B132A-ABCF-4EAC-8B5C-4AFCE885915D}  
TERM.....NEPRO CARB STEADY

DATA VALUE: NEPRO RENAL FORMULA

-----  
ABBREVIATION.....NEP RENAL  
ACTIVE.....Yes  
DESCRIPTION.....Nepro Renal Enteral Nutrition Formula  
DISPLAY\_NAME.....Nepro Renal Enteral Nutrition Formula  
HELP\_TEXT.....No Help Text  
ID.....{B8F28816-BC81-443F-AEC5-DBAE4D435998}  
TERM.....NEPRO RENAL FORMULA

DATA VALUE: NEUROGENIC BLADDER

-----  
ABBREVIATION.....Neuro Blad  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Neurogenic Bladder  
HELP\_TEXT.....  
ID.....{D386E29C-C4B7-4A9D-8092-CA978F9CC4D6}  
TERM.....NEUROGENIC BLADDER

DATA VALUE: NEUTRAL PRESSURE CONNECTOR

-----  
ABBREVIATION.....N Pres Con  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Neutral Pressure Connector  
HELP\_TEXT.....  
ID.....{263F1A2C-4384-4B6A-B930-B9AF85B7F965}  
TERM.....NEUTRAL PRESSURE CONNECTOR



DATA VALUE: NEUTROPENIC

-----  
ABBREVIATION.....NEUTROPENI  
ACTIVE.....Yes  
DESCRIPTION.....Neutropenic Diet  
DISPLAY\_NAME.....Neutropenic  
HELP\_TEXT.....No Help Text  
ID.....{5CB9620C-34D9-48F7-AAF0-C595D06FC8AE}  
TERM.....NEUTROPENIC

DATA VALUE: NEW ADMISSION EDUCATION

-----  
ABBREVIATION.....N ADM EDU  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....New Admission Education  
HELP\_TEXT.....No Help Text  
ID.....{1890BC97-FEEB-4BED-A054-3C754C01817B}  
TERM.....NEW ADMISSION EDUCATION

DATA VALUE: NEW SITE

-----  
ABBREVIATION.....NEW SITE  
ACTIVE.....Yes  
DESCRIPTION.....New Site  
DISPLAY\_NAME.....New Site  
HELP\_TEXT.....No Help Text  
ID.....{10B1FF61-E58E-4608-9CCB-9232BC964418}  
TERM.....NEW SITE

DATA VALUE: NIGHT LIGHT PROVIDED

-----  
ABBREVIATION.....Night Lite  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Night Light Provided  
HELP\_TEXT.....  
ID.....{A6F3E41A-92A1-4AAE-9EFC-84D21E8A32BD}  
TERM.....NIGHT LIGHT PROVIDED

DATA VALUE: NLP: NO LIGHT PERCEPTION

-----  
ABBREVIATION.....NLP  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....NLP: No Light Perception  
HELP\_TEXT.....  
ID.....{F710E3C5-9B5A-4826-876D-B5675F9AA67E}  
TERM.....NLP: NO LIGHT PERCEPTION

DATA VALUE: NO

-----  
ABBREVIATION.....n  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....No  
HELP\_TEXT.....No Help Text  
ID.....{FC275445-1225-9912-18A3-DF36245893D3}  
TERM.....NO

DATA VALUE: NO BREAK IN EPIDERMIS

-----  
ABBREVIATION.....NO BREAK E  
ACTIVE.....Yes  
DESCRIPTION.....No Break In Epidermis  
DISPLAY\_NAME.....No Break In Epidermis  
HELP\_TEXT.....No Help Text  
ID.....{615D96CD-6096-446B-BDDF-3F7242936AF9}  
TERM.....NO BREAK IN EPIDERMIS

DATA VALUE: NO MOVEMENT

-----  
ABBREVIATION.....nm  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....No Movement  
HELP\_TEXT.....No Help Text  
ID.....{2618E301-379B-BAB9-9EDE-98FD4B60D311}  
TERM.....NO MOVEMENT

DATA VALUE: NO SIGNS AND SYMPTOMS

-----  
ABBREVIATION.....NO SIG&SYM  
ACTIVE.....Yes  
DESCRIPTION.....No Signs And Symptoms  
DISPLAY\_NAME.....No Signs And Symptoms  
HELP\_TEXT.....No Help Text  
ID.....{91285942-6C82-4CB1-AAFB-8D3DA15A0B8A}  
TERM.....NO SIGNS AND SYMPTOMS

DATA VALUE: NO WANDERING PREVENTION MEASURES APPLIED

-----  
ABBREVIATION.....No Wan App  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....No Wandering Prevention Measures  
Applied  
HELP\_TEXT.....  
ID.....{6967F640-223E-46CE-82F9-0E4D87745D94}  
TERM.....NO WANDERING PREVENTION MEASURES  
APPLIED

DATA VALUE: NO WANDERING PREVENTION MEASURES NEEDED

-----  
ABBREVIATION.....No Wan Nee  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....No Wandering Prevention Measures Needed  
HELP\_TEXT.....  
ID.....{8D39BA2B-24C6-4355-9AC2-4AAAD2E66D9B}  
TERM.....NO WANDERING PREVENTION MEASURES NEEDED

DATA VALUE: NON-ADHERENT GAUZE

-----  
ABBREVIATION.....NON-AD GAU  
ACTIVE.....Yes  
DESCRIPTION.....Non-Adherent Gauze  
DISPLAY\_NAME.....Non-Adherent gauze  
HELP\_TEXT.....No Help Text  
ID.....{1DBF7C56-F5C2-42B5-AB59-B39A5069A414}  
TERM.....NON-ADHERENT GAUZE

DATA VALUE: NON-PITTING

-----  
ABBREVIATION.....NON-PITTIN  
ACTIVE.....Yes  
DESCRIPTION.....NON-PITTING  
DISPLAY\_NAME.....Non-Pitting  
HELP\_TEXT.....NON-PITTING  
ID.....{786DEC76-36DE-FDC7-885C-F84210F7C7D2}  
TERM.....NON-PITTING

DATA VALUE: NON-PRODUCTIVE

-----  
ABBREVIATION.....np  
ACTIVE.....Yes  
DESCRIPTION.....Does not produce desired result  
DISPLAY\_NAME.....Non-Productive  
HELP\_TEXT.....Does not produce desired result  
ID.....{9AB01716-E395-18B9-6C5D-9BE6AEE41D5B}  
TERM.....NON-PRODUCTIVE

DATA VALUE: NON-SKID SLIPPERS

-----  
ABBREVIATION.....Non-Skid  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Non-Skid Slippers  
HELP\_TEXT.....  
ID.....{D7641B31-2CB9-4284-B505-74B131BF4789}  
TERM.....NON-SKID SLIPPERS

DATA VALUE: NON-TENDER

-----  
ABBREVIATION.....NON-TENDER  
ACTIVE.....Yes  
DESCRIPTION.....Not Tender  
DISPLAY\_NAME.....Non-Tender  
HELP\_TEXT.....No Help Text  
ID.....{5A3D41C0-A83D-4382-A573-D6884BB1670B}  
TERM.....NON-TENDER

DATA VALUE: NONBLANCHABLE ERYTHEMA

-----  
ABBREVIATION.....N BLAN ERY  
ACTIVE.....Yes  
DESCRIPTION.....Nonblanchable Erythema  
DISPLAY\_NAME.....Nonblanchable Erythema  
HELP\_TEXT.....No Help Text  
ID.....{30273ACE-03F3-4818-995B-26586AAB5499}  
TERM.....NONBLANCHABLE ERYTHEMA

DATA VALUE: NONE

-----  
ABBREVIATION.....NONE  
ACTIVE.....Yes  
DESCRIPTION.....None  
DISPLAY\_NAME.....None  
HELP\_TEXT.....No Help Text  
ID.....{D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}  
TERM.....NONE

DATA VALUE: NONE IDENTIFIED AT THIS VISIT

-----  
ABBREVIATION.....None  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....None identified at this Visit  
HELP\_TEXT.....  
ID.....{1B1A62CC-FB5C-4051-95E1-448CA3D0B5F8}  
TERM.....NONE IDENTIFIED AT THIS VISIT

DATA VALUE: NONREACTIVE

-----  
ABBREVIATION.....nonreact  
ACTIVE.....Yes  
DESCRIPTION.....Nonreactive  
DISPLAY\_NAME.....Nonreactive  
HELP\_TEXT.....No Help Text  
ID.....{AB6B7D92-B94C-F3AA-E32F-6035D9CF2EF9}  
TERM.....NONREACTIVE

DATA VALUE: NORMAL

-----  
ABBREVIATION.....n  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Normal  
HELP\_TEXT.....No Help Text  
ID.....{553C9A97-A3D5-C52F-0E4E-C7D168C99A2B}  
TERM.....NORMAL

DATA VALUE: NORMAL FOR ETHNIC GROUP

-----  
ABBREVIATION.....ETH. NORM.  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Normal for ethnic group  
HELP\_TEXT.....Skin color is normal for ethnic group  
ID.....{3B2CD951-58CE-FF53-AD5D-6C65C07658D3}  
TERM.....NORMAL FOR ETHNIC GROUP

DATA VALUE: NORMAL POWER

-----  
ABBREVIATION.....np  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Normal Power  
HELP\_TEXT.....No Help Text  
ID.....{0E45C061-0744-8A21-7256-003702A5005B}  
TERM.....NORMAL POWER

DATA VALUE: NORMAL SINUS

-----  
ABBREVIATION.....ns  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Normal Sinus  
HELP\_TEXT.....No Help Text  
ID.....{735DC86D-08CD-1533-5A36-4775A882C507}  
TERM.....NORMAL SINUS

DATA VALUE: NORMAL STOOL

-----  
ABBREVIATION.....NORM STOOL  
ACTIVE.....Yes  
DESCRIPTION.....Normal Stool  
DISPLAY\_NAME.....Normal Stool  
HELP\_TEXT.....No Help Text  
ID.....{486B0693-CAC5-48EE-B2E6-1F73C651418A}  
TERM.....NORMAL STOOL

DATA VALUE: NOT DONE

-----  
ABBREVIATION.....NOT DONE  
ACTIVE.....Yes  
DESCRIPTION.....Not Done  
DISPLAY\_NAME.....Not Done  
HELP\_TEXT.....No Help Text  
ID.....{C3FE84AC-EEA1-4153-AE6C-69755C2E6BB1}  
TERM.....NOT DONE

DATA VALUE: NOT MONITORED

-----  
ABBREVIATION.....NOT MONITO  
ACTIVE.....Yes  
DESCRIPTION.....Not Monitored  
DISPLAY\_NAME.....Not Monitored  
HELP\_TEXT.....No Help Text  
ID.....{0E6A8A07-48E0-428F-8D86-22572CABCF48}  
TERM.....NOT MONITORED

DATA VALUE: NOVASOURCE 2.0

-----  
ABBREVIATION.....NOVAS 2.0  
ACTIVE.....Yes  
DESCRIPTION.....Novasource 2.0 Enteral Nutrition  
Formula  
DISPLAY\_NAME.....Novasource 2.0 Enteral Nutrition  
Formula  
HELP\_TEXT.....No Help Text  
ID.....{ACE82591-502A-4FBA-9697-B1B330DED78A}  
TERM.....NOVASOURCE 2.0

DATA VALUE: NP

-----  
ABBREVIATION.....NP  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....NP  
HELP\_TEXT.....  
ID.....{E7EA4017-4374-41A1-9D3D-54051CEDF63C}  
TERM.....NP

DATA VALUE: NPO DIET

-----  
ABBREVIATION.....NPO DIET  
ACTIVE.....Yes  
DESCRIPTION.....Nothing By Mouth Diet  
DISPLAY\_NAME.....NPO Diet  
HELP\_TEXT.....Nothing By Mouth Diet  
ID.....{BEC8F771-B082-40A2-9EE8-553222A25FBF}  
TERM.....NPO DIET

DATA VALUE: NSG ASSIST

-----  
ABBREVIATION.....Nsg Assist  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Nsg Assist  
HELP\_TEXT.....  
ID.....{14FA9638-3786-4FEA-BE8F-2345C1F455C3}  
TERM.....NSG ASSIST

DATA VALUE: NUMB

-----  
ABBREVIATION.....NUMB  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....NUMB  
HELP\_TEXT.....No Help Text  
ID.....{652634C0-5334-4D56-BE57-5F58F952FF0D}  
TERM.....NUMB

DATA VALUE: NUMBER: SPECIFY

-----  
ABBREVIATION.....Specify  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Number: Specify  
HELP\_TEXT.....When "Number" is selected, enter  
                  numeric value as a "Comment"  
ID.....{397D773F-78F9-4E69-9A84-B518109A95AB}  
TERM.....NUMBER: SPECIFY

DATA VALUE: NUMBNESS

-----  
ABBREVIATION.....n  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Numbness  
HELP\_TEXT.....No Help Text  
ID.....{037F2F79-36D9-A0B3-C23E-8EE525932AF3}  
TERM.....NUMBNESS

DATA VALUE: NURSING STATION

-----  
ABBREVIATION.....ns  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Nursing Station  
HELP\_TEXT.....No Help Text  
ID.....{B6E48F92-B9BA-488F-3A2E-D74568A6B855}  
TERM.....NURSING STATION

DATA VALUE: NUTRIHEP

-----  
ABBREVIATION.....Nutrihep  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Nutrihep  
HELP\_TEXT.....  
ID.....{ABD98675-EB60-40C7-BCE6-73FD787F6D17}  
TERM.....NUTRIHEP

DATA VALUE: NUTRITION

-----  
ABBREVIATION.....NUTRITION  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Nutrition  
HELP\_TEXT.....No Help Text  
ID.....{429A28EF-8A1F-4247-85E8-68D4EFC189C}  
TERM.....NUTRITION

DATA VALUE: NUTRITIONAL SUPPLEMENT

-----  
ABBREVIATION.....NUTR SUPPL  
ACTIVE.....Yes  
DESCRIPTION.....Nutritional Supplement  
DISPLAY\_NAME.....Nutritional Supplement  
HELP\_TEXT.....No Help Text  
ID.....{4573BD4E-7220-4B88-A2EE-959F8A0E1199}  
TERM.....NUTRITIONAL SUPPLEMENT

DATA VALUE: OBEYS COMMANDS

-----  
ABBREVIATION.....OBEYS COMM  
ACTIVE.....Yes  
DESCRIPTION.....Obeys Commands  
DISPLAY\_NAME.....Obeys Commands  
HELP\_TEXT.....No Help Text  
ID.....{EE88EE08-4C71-4C70-B93E-65F9091544C1}  
TERM.....OBEYS COMMANDS

DATA VALUE: OBSTRUCTIVE SLEEP APNEA

-----  
ABBREVIATION.....Ob Sleep A  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Obstructive Sleep Apnea  
HELP\_TEXT.....  
ID.....{217F9D63-4624-4201-9F98-C0D992B663DE}  
TERM.....OBSTRUCTIVE SLEEP APNEA



DATA VALUE: OCCLUDED

-----  
ABBREVIATION.....Occluded  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Occluded  
HELP\_TEXT.....  
ID.....{25509BE3-F946-4AD7-80E8-5CE8A2E61AD6}  
TERM.....OCCLUDED

DATA VALUE: OFF DURING BATHING

-----  
ABBREVIATION.....OFF BATHIN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Off During Bathing  
HELP\_TEXT.....No Help Text  
ID.....{5D0EFDD9-F170-4F34-9ABE-CE2D376E65F4}  
TERM.....OFF DURING BATHING

DATA VALUE: OFF WHILE AMBULATORY

-----  
ABBREVIATION.....OFF AMBULA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Off While Ambulatory  
HELP\_TEXT.....No Help Text  
ID.....{FF4610BF-8C92-408F-882A-1D6B2C2F146F}  
TERM.....OFF WHILE AMBULATORY

DATA VALUE: OFFER PRN MEDS

-----  
ABBREVIATION.....OFFER PRN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Offer PRN Meds  
HELP\_TEXT.....No Help Text  
ID.....{B7AF4903-913C-4848-A255-0B4767B81AE2}  
TERM.....OFFER PRN MEDS

DATA VALUE: OFFERED

-----  
ABBREVIATION.....OFFERED  
ACTIVE.....Yes  
DESCRIPTION.....Offered  
DISPLAY\_NAME.....Offered  
HELP\_TEXT.....No Help Text  
ID.....{715091C2-9668-4678-900E-6D55D3335A30}  
TERM.....OFFERED

DATA VALUE: ON WHILE SITTING IN CHAIR

-----  
ABBREVIATION.....ON IN CHAI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....On While Sitting In Chair  
HELP\_TEXT.....No Help Text  
ID.....{4139F4C3-BB68-4477-BF63-9BE7813AADA7}  
TERM.....ON WHILE SITTING IN CHAIR

DATA VALUE: ONE PERSON ASSIST

-----  
ABBREVIATION.....1 PERS AS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....One Person Assist  
HELP\_TEXT.....No Help Text  
ID.....{A3D14AA8-D658-412F-BBE8-D8AD195ADA3D}  
TERM.....ONE PERSON ASSIST

DATA VALUE: ONE TO ONE SUPERVISION

-----  
ABBREVIATION.....1 SUPERV  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....One To One Supervision  
HELP\_TEXT.....No Help Text  
ID.....{CFC9C653-BFE9-40B2-BA7A-944116C7FCB7}  
TERM.....ONE TO ONE SUPERVISION

DATA VALUE: OOZING

-----  
ABBREVIATION.....Oozing  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Oozing  
HELP\_TEXT.....  
ID.....{C23549D7-3238-4523-BA31-AB71609F01C2}  
TERM.....OOZING

DATA VALUE: OPEN SPONTANEOUSLY

-----  
ABBREVIATION.....OPEN SPONT  
ACTIVE.....Yes  
DESCRIPTION.....Open Spontaneously  
DISPLAY\_NAME.....Open Spontaneously  
HELP\_TEXT.....No Help Text  
ID.....{013F1473-BDC9-4E7C-9709-C577357709F6}  
TERM.....OPEN SPONTANEOUSLY

DATA VALUE: OPEN SYSTEM

-----  
ABBREVIATION.....Open syste  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Open System  
HELP\_TEXT.....  
ID.....{2047EDE8-5D72-4541-A23A-81F097A6CA2A}  
TERM.....OPEN SYSTEM

DATA VALUE: OPEN TO AIR

-----  
ABBREVIATION.....OPE TO AIR  
ACTIVE.....Yes  
DESCRIPTION.....Open To Air  
DISPLAY\_NAME.....Open To Air  
HELP\_TEXT.....No Help Text  
ID.....{EFB715D0-5666-4281-9034-A0326A14B2BB}  
TERM.....OPEN TO AIR

DATA VALUE: OPERATING ROOM

-----  
ABBREVIATION.....OR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Operating Room  
HELP\_TEXT.....No Help Text  
ID.....{37DACE14-2A28-45F4-9DEA-35A1A433F4E6}  
TERM.....OPERATING ROOM

DATA VALUE: ORAL

-----  
ABBREVIATION.....ORAL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Oral  
HELP\_TEXT.....No Help Text  
ID.....{BBAF0DD0-1F52-49EE-9A47-FA9323D6C529}  
TERM.....ORAL

DATA VALUE: ORAL AIRWAY

-----  
ABBREVIATION.....ORAL AIRWA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Oral Airway  
HELP\_TEXT.....No Help Text  
ID.....{B0CF4270-A68D-1322-0621-1DFF8D677B40}  
TERM.....ORAL AIRWAY

DATA VALUE: ORAL ENDOTRACHEAL

-----  
ABBREVIATION.....ORAL ENDO  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Oral Endotracheal  
HELP\_TEXT.....No Help Text  
ID.....{466CAD16-D7F3-9D73-6500-959C4DA50D9A}  
TERM.....ORAL ENDOTRACHEAL

DATA VALUE: ORAL GASTRIC TUBE

-----  
ABBREVIATION.....OG TUBE  
ACTIVE.....Yes  
DESCRIPTION.....Oral Gastric Enteral Feeding Tube  
DISPLAY\_NAME.....Oral Gastric Tube  
HELP\_TEXT.....No Help Text  
ID.....{C6995D6A-939C-4C31-A3C2-4504CB5597CD}  
TERM.....ORAL GASTRIC TUBE

DATA VALUE: ORANGE

-----  
ABBREVIATION.....ORANGE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Orange  
HELP\_TEXT.....No Help Text  
ID.....{77B1B667-212D-88BA-D8E6-5C22F6ACDC9C}  
TERM.....ORANGE

DATA VALUE: ORIENTATION

-----  
ABBREVIATION.....ORIENTA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Orientation  
HELP\_TEXT.....No Help Text  
ID.....{10FA5B0E-F4A7-44E4-A482-7C80F145C4D2}  
TERM.....ORIENTATION

DATA VALUE: ORIENTED

-----  
ABBREVIATION.....ORIENTED  
ACTIVE.....Yes  
DESCRIPTION.....Oriented  
DISPLAY\_NAME.....Oriented  
HELP\_TEXT.....No Help Text  
ID.....{B8DC9245-6D81-40E4-8B47-B6413E49F152}  
TERM.....ORIENTED

DATA VALUE: ORIENTED TO CALL FOR ASSISTANCE

-----  
ABBREVIATION.....Call/Asst  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Oriented to Call for Assistance  
HELP\_TEXT.....  
ID.....{1F8390CF-2F88-486A-8583-AF7C7A64F6B5}  
TERM.....ORIENTED TO CALL FOR ASSISTANCE

DATA VALUE: ORIENTED TO PERSON

-----  
ABBREVIATION.....otp  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Oriented To Person  
HELP\_TEXT.....No Help Text  
ID.....{9912CFB9-AC55-6445-3EC1-0F1DF98C0561}  
TERM.....ORIENTED TO PERSON

DATA VALUE: ORIENTED TO PLACE

-----  
ABBREVIATION.....otp  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Oriented To Place  
HELP\_TEXT.....No Help Text  
ID.....{8E5933E2-2DF0-C8A4-7A92-20E792F42EA6}  
TERM.....ORIENTED TO PLACE

DATA VALUE: ORIENTED TO SITUATION

-----  
ABBREVIATION.....ots  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Oriented To Situation  
HELP\_TEXT.....No Help Text  
ID.....{955A2F3B-E444-99DF-81DD-5E6BF9D22E49}  
TERM.....ORIENTED TO SITUATION

DATA VALUE: ORIENTED TO SURROUNDINGS/ROUTINE

-----  
ABBREVIATION.....Oriented  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Oriented to Surroundings/Routine  
HELP\_TEXT.....  
ID.....{A8F3F756-E5B5-4B8D-AF21-E45746A64E46}  
TERM.....ORIENTED TO SURROUNDINGS/ROUTINE

DATA VALUE: ORIENTED TO TIME

-----  
ABBREVIATION.....ott  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Oriented To Time  
HELP\_TEXT.....No Help Text  
ID.....{C79EBE25-C51B-478C-CB39-25943792F191}  
TERM.....ORIENTED TO TIME

DATA VALUE: ORTHOPNEIC

-----  
ABBREVIATION.....ORTHOPNIC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Orthopneic  
HELP\_TEXT.....No Help Text  
ID.....{FF373BD7-FD84-DF43-7A44-1B8E62917116}  
TERM.....ORTHOPNEIC

DATA VALUE: ORTHOTIC BOOTS

-----  
ABBREVIATION.....ORTHO BOOT  
ACTIVE.....Yes  
DESCRIPTION.....ORTHOTIC BOOTS  
DISPLAY\_NAME.....Orthotic Boots  
HELP\_TEXT.....ORTHOTIC BOOTS  
ID.....{D4219FE9-7087-B881-6DAC-4E0BB94CDA8E}  
TERM.....ORTHOTIC BOOTS

DATA VALUE: ORTHOTIC COLLAR

-----  
ABBREVIATION.....ORTHO COLL  
ACTIVE.....Yes  
DESCRIPTION.....ORTHOTIC COLLAR  
DISPLAY\_NAME.....Orthotic Collar  
HELP\_TEXT.....ORTHOTIC COLLAR  
ID.....{A74E1C89-2304-0544-9A8B-A23E4E2C6B08}  
TERM.....ORTHOTIC COLLAR

DATA VALUE: ORTHOTIC KNEE IMMOBILIZER

-----  
ABBREVIATION.....KNEE IMMOB  
ACTIVE.....Yes  
DESCRIPTION.....ORTHOTIC KNEE IMMOBILIZER  
DISPLAY\_NAME.....Orthotic Knee Immobilizer  
HELP\_TEXT.....ORTHOTIC KNEE IMMOBILIZER  
ID.....{FF8A8BB4-63FB-3B0D-83E1-5432A0434DAC}  
TERM.....ORTHOTIC KNEE IMMOBILIZER

DATA VALUE: ORTHOTIC SPINE BRACE

-----  
ABBREVIATION.....ORTH SPINE  
ACTIVE.....Yes  
DESCRIPTION.....ORTHOTIC SPINE BRACE  
DISPLAY\_NAME.....Orthotic Spine Brace  
HELP\_TEXT.....ORTHOTIC SPINE BRACE  
ID.....{B47A3536-FEBC-802E-01A1-E909147D18C2}  
TERM.....ORTHOTIC SPINE BRACE

DATA VALUE: ORTHOTIC TED HOSE

-----  
ABBREVIATION.....TED HOSE  
ACTIVE.....Yes  
DESCRIPTION.....TED HOSE  
DISPLAY\_NAME.....Orthotic TED Hose  
HELP\_TEXT.....TED HOSE  
ID.....{FE277269-7BD1-85FB-DC3C-FEDE31C66BB0}  
TERM.....ORTHOTIC TED HOSE

DATA VALUE: ORTHOTIC WEDGE PILLOW

-----  
ABBREVIATION.....WEDG PILLO  
ACTIVE.....Yes  
DESCRIPTION.....WEDGE PILLOW  
DISPLAY\_NAME.....Orthotic Wedge Pillow  
HELP\_TEXT.....WEDGE PILLOW  
ID.....{20E58D21-399B-85E4-BDAB-93B371395BCF}  
TERM.....ORTHOTIC WEDGE PILLOW

DATA VALUE: OSMOLITE

-----  
ABBREVIATION.....OSMOLITE  
ACTIVE.....Yes  
DESCRIPTION.....Enteral Nutrition Type  
DISPLAY\_NAME.....Osmolite  
HELP\_TEXT.....No Help Text  
ID.....{8E41BB55-063E-43C0-BAFE-65150EDD99F8}  
TERM.....OSMOLITE

DATA VALUE: OSMOLITE 1.0

-----  
ABBREVIATION.....Osmo 1.0  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Osmolite 1.0  
HELP\_TEXT.....  
ID.....{D5EE7DA5-707B-4000-A133-0C571824A376}  
TERM.....OSMOLITE 1.0

DATA VALUE: OSMOLITE 1.2

-----  
ABBREVIATION.....Osmo 1.2  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Osmolite 1.2  
HELP\_TEXT.....  
ID.....{CEF57C44-82E5-43F5-8646-F8FBCD4609A0}  
TERM.....OSMOLITE 1.2

DATA VALUE: OSMOLITE 1.5

-----  
ABBREVIATION.....Osmo 1.5  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Osmolite 1.5  
HELP\_TEXT.....  
ID.....{8AF4CD1A-631A-4E95-A631-ACFF775284BA}  
TERM.....OSMOLITE 1.5

DATA VALUE: OSMOLITE HN

-----  
ABBREVIATION.....OSMOLIT HN  
ACTIVE.....Yes  
DESCRIPTION.....Enteral Nutrition Type  
DISPLAY\_NAME.....Osmolite HN  
HELP\_TEXT.....No Help Text  
ID.....{E4186CA6-302F-4513-B125-30B094151423}  
TERM.....OSMOLITE HN

DATA VALUE: OSTOMY CARE

-----  
ABBREVIATION.....Ostomy Car  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Ostomy Care  
HELP\_TEXT.....  
ID.....{BC834BC9-5861-4B94-BA34-FDF5C9FDF979}  
TERM.....OSTOMY CARE

DATA VALUE: OTHER

-----  
ABBREVIATION.....Other  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Other  
HELP\_TEXT.....  
ID.....{4CBA16B3-6BBA-4B8D-AA12-D721AECFE286}  
TERM.....OTHER



DATA VALUE: OTHER (ADD COMMENTS)

-----  
ABBREVIATION.....Other  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Other (Add Comments)  
HELP\_TEXT.....  
ID.....{264AB0B9-64B6-4EC4-BBAC-8A879CCDB23B}  
TERM.....OTHER (ADD COMMENTS)

DATA VALUE: OTHER (SEE COMMENTS)

-----  
ABBREVIATION.....Other  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Other (see comments)  
HELP\_TEXT.....  
ID.....{7908D30B-8BD4-4A02-964D-8E6DBCA24A81}  
TERM.....OTHER (SEE COMMENTS)

DATA VALUE: OVER WIRE

-----  
ABBREVIATION.....OVER WIRE  
ACTIVE.....Yes  
DESCRIPTION.....Over Wire  
DISPLAY\_NAME.....Over Wire  
HELP\_TEXT.....No Help Text  
ID.....{157CD456-1020-43AE-9AAA-1686F1574B6F}  
TERM.....OVER WIRE

DATA VALUE: OXEPA

-----  
ABBREVIATION.....Oxepa  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Oxepa  
HELP\_TEXT.....  
ID.....{1CE54D02-5FBF-43A0-93F4-E5D5DB86FEF6}  
TERM.....OXEPA

DATA VALUE: OXIMETRIC

-----  
ABBREVIATION.....OXIMETRIC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Oximetric Catheter  
HELP\_TEXT.....No Help Text  
ID.....{DA1F5AAA-A349-C17B-74CD-1794492E0A13}  
TERM.....OXIMETRIC

DATA VALUE: PA

-----  
ABBREVIATION.....PA  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....PA  
HELP\_TEXT.....  
ID.....{C057C47D-780A-4230-B577-2C528AD09DEB}  
TERM.....PA

DATA VALUE: PAC

-----  
ABBREVIATION.....p  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....PAC  
HELP\_TEXT.....No Help Text  
ID.....{848AB18A-598D-8E4E-8B1B-ED7478CE340B}  
TERM.....PAC

DATA VALUE: PACED: ATRIAL

-----  
ABBREVIATION.....pa  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Paced: Atrial  
HELP\_TEXT.....No Help Text  
ID.....{88837DD4-098D-B77F-A754-89D5C4F5870F}  
TERM.....PACED: ATRIAL

DATA VALUE: PACED: AV SEQUENTIAL

-----  
ABBREVIATION.....pas  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Paced: Av Sequential  
HELP\_TEXT.....No Help Text  
ID.....{3C485674-224C-47DA-510D-92F1A0CA54A1}  
TERM.....PACED: AV SEQUENTIAL

DATA VALUE: PACED: VENTRICULAR

-----  
ABBREVIATION.....pv  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Paced: Ventricular  
HELP\_TEXT.....No Help Text  
ID.....{A63F6E6E-A7C4-0928-81FB-FE2F70F713C6}  
TERM.....PACED: VENTRICULAR

DATA VALUE: PACER

-----  
ABBREVIATION.....Pacer  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Pacer  
HELP\_TEXT.....  
ID.....{F8C0D639-D7B2-4342-A9FC-8A0F54BEB7DE}  
TERM.....PACER

DATA VALUE: PAIN

-----  
ABBREVIATION.....PAIN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Pain  
HELP\_TEXT.....No Help Text  
ID.....{F2F23EE7-0A06-46B8-A496-B7A4651A551F}  
TERM.....PAIN

DATA VALUE: PAIN CONTROLLED PER PT

-----  
ABBREVIATION.....PAIN CON  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Pain Controlled Per Patient  
HELP\_TEXT.....No Help Text  
ID.....{7EC1AF38-39B1-475E-B58C-E33404A47E96}  
TERM.....PAIN CONTROLLED PER PT

DATA VALUE: PAINAD BL 0

-----  
ABBREVIATION.....PAINAD BLO  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....0 = Relaxed  
HELP\_TEXT.....No Help Text  
ID.....{C77929DB-A59A-4A24-B57B-BF111770C2AB}  
TERM.....PAINAD BL 0

DATA VALUE: PAINAD BL 1

-----  
ABBREVIATION.....PAINAD BL1  
ACTIVE.....Yes  
DESCRIPTION.....1 = Tense, Distressed Pacing, Fidgeting  
DISPLAY\_NAME.....1 = Tense, Distressed Pacing, Fidgeting  
HELP\_TEXT.....1 = Tense, Distressed Pacing, Fidgeting  
ID.....{CC1E1516-086D-4B79-8264-0C3BD67D8552}  
TERM.....PAINAD BL 1

DATA VALUE: PAINAD BL 2

-----  
ABBREVIATION.....PAINAD BL2  
ACTIVE.....Yes  
DESCRIPTION.....2 = Rigid, Fists Clenched, Knees  
                  Pulled Up, Pulling or Pushing Away,  
                  Striking Out  
DISPLAY\_NAME.....2 = Rigid, Fists Clenched, Knees  
                  Pulled Up, Pulling or Pushing  
HELP\_TEXT.....2 = Rigid, Fists Clenched, Knees  
                  Pulled Up, Pulling or Pushing Away,  
                  Striking Out  
ID.....{34C3E304-DF66-40A7-906F-22AEA3F76FAD}  
TERM.....PAINAD BL 2

DATA VALUE: PAINAD BR 0

-----  
ABBREVIATION.....PAINAD B0  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....0 = Normal  
HELP\_TEXT.....No Help Text  
ID.....{58E71B2E-F5C3-419F-8B71-A9FD701E0202}  
TERM.....PAINAD BR 0

DATA VALUE: PAINAD BR 1

-----  
ABBREVIATION.....PAINAD BR1  
ACTIVE.....Yes  
DESCRIPTION.....1 = Occasional Labored Breathing.  
                  Short Period of Hyperventilation  
DISPLAY\_NAME.....1 = Occasional Labored Breathing.  
                  Short Period of Hyperventilation  
HELP\_TEXT.....1 = Occasional Labored Breathing.  
                  Short Period of Hyperventilation  
ID.....{727D1B04-A702-4806-857B-78BCBACB0FE1}  
TERM.....PAINAD BR 1

DATA VALUE: PAINAD BR 2

-----  
ABBREVIATION.....PAINAD BR2  
ACTIVE.....Yes  
DESCRIPTION.....2 = Noisy Labored Breathing.  
                  Cheynes-Stokes Respirations  
DISPLAY\_NAME.....2 = Noisy Labored Breathing.  
                  Cheynes-Stokes Respirations  
HELP\_TEXT.....2 = Noisy Labored Breathing.  
                  Cheynes-Stokes Respirations  
ID.....{37978012-ADAC-416C-8920-48CDAAA5D59D}  
TERM.....PAINAD BR 2

DATA VALUE: PAINAD CO 0

-----  
ABBREVIATION.....PAINAD CO0  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....0 = No Need To Console  
HELP\_TEXT.....No Help Text  
ID.....{06A89B55-5822-42F7-A689-DBD91B62979F}

TERM.....PAINAD CO 0

DATA VALUE: PAINAD CO 1

-----  
ABBREVIATION.....PAINAD CO1  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....1 = Distracted or Reassured By Voice  
                  or Touch  
HELP\_TEXT.....No Help Text  
ID.....{8E5DDB47-1955-40AD-A008-A9CE4646A707}  
TERM.....PAINAD CO 1

DATA VALUE: PAINAD CO 2

-----  
ABBREVIATION.....PAINAD CO2  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....2 = Unable To Console, Distract or  
                  Reassure  
HELP\_TEXT.....No Help Text  
ID.....{1DFFB735-6272-4B60-8F42-CCCAA701C7A7}  
TERM.....PAINAD CO 2

DATA VALUE: PAINAD FE 0

-----  
ABBREVIATION.....PAINAD FE0  
ACTIVE.....Yes  
DESCRIPTION.....0 = Smiling or Inexpressive  
DISPLAY\_NAME.....0 = Smiling or Inexpressive  
HELP\_TEXT.....0 = Smiling or Inexpressive  
ID.....{B3DE1BC6-E7A0-48C5-8C8E-FBE033AD1046}  
TERM.....PAINAD FE 0

DATA VALUE: PAINAD FE 1

-----  
ABBREVIATION.....PAINAD FE1  
ACTIVE.....Yes  
DESCRIPTION.....1 = Sad, Frightened, Frown  
DISPLAY\_NAME.....1 = Sad, Frightened, Frown  
HELP\_TEXT.....1 = Sad, Frightened, Frown  
ID.....{66C01EE9-8F36-40B5-AED2-DBC2E94F860D}  
TERM.....PAINAD FE 1

DATA VALUE: PAINAD FE 2

-----  
ABBREVIATION.....PAINAD FE2  
ACTIVE.....Yes  
DESCRIPTION.....2 = Facial Grimacing  
DISPLAY\_NAME.....2 = Facial Grimacing  
HELP\_TEXT.....2 = Facial Grimacing  
ID.....{35688D06-94F2-4B36-8940-529803D4BC43}  
TERM.....PAINAD FE 2

DATA VALUE: PAINAD NV 0

-----  
ABBREVIATION.....PAINAD NV0  
ACTIVE.....Yes  
DESCRIPTION.....0 = None  
DISPLAY\_NAME.....0 = None  
HELP\_TEXT.....0 = None  
ID.....{63F7F990-2DEA-4B35-AC52-D23B8E13B1FC}  
TERM.....PAINAD NV 0

DATA VALUE: PAINAD NV 1

-----  
ABBREVIATION.....PAINAD NV1  
ACTIVE.....Yes  
DESCRIPTION.....1 = Occasional Moan or Groan. Low  
Level Speech With A Negative or  
Disapproving Quality  
DISPLAY\_NAME.....1 = Occasional Moan or Groan  
HELP\_TEXT.....1 = Occasional Moan or Groan. Low  
Level Speech With A Negative or  
Disapproving Quality  
ID.....{D517E6C9-9D98-4D0F-A8F4-261BED4812D1}  
TERM.....PAINAD NV 1

DATA VALUE: PAINAD NV 2

-----  
ABBREVIATION.....PAINAD NV2  
ACTIVE.....Yes  
DESCRIPTION.....2 = Repeated, Troubled, Calling Out,  
Loud Moaning, Crying  
DISPLAY\_NAME.....2 = Repeated, Troubled, Calling Out,  
Loud Moaning.  
HELP\_TEXT.....2 = Repeated, Troubled, Calling Out,  
Loud Moaning, Crying  
ID.....{3C8D6B1E-53E6-4D2C-BD8A-9C548466E1D4}  
TERM.....PAINAD NV 2

DATA VALUE: PAINFUL EXPRESSION

-----  
ABBREVIATION.....PAINFUL EX  
ACTIVE.....Yes  
DESCRIPTION.....Painful Expression  
DISPLAY\_NAME.....Painful Expression  
HELP\_TEXT.....No Help Text  
ID.....{608554F9-29F9-4B27-BE63-486AF81810BD}  
TERM.....PAINFUL EXPRESSION

DATA VALUE: PAIRS

-----  
ABBREVIATION.....Pairs  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Pairs  
HELP\_TEXT.....  
ID.....{668B0D85-E507-4E12-B07A-A90CB98730DC}  
TERM.....PAIRS

DATA VALUE: PALE

-----  
ABBREVIATION.....p  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Pale  
HELP\_TEXT.....No Help Text  
ID.....{0EB39D28-7ABE-A4AB-D334-3F795E65618D}  
TERM.....PALE

DATA VALUE: PAROXYSMAL ATRIAL TACHYCARDIA

-----  
ABBREVIATION.....pat  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Paroxysmal Atrial Tachycardia  
HELP\_TEXT.....No Help Text  
ID.....{27BC1FAF-A403-6EA8-044C-F1271F12A62A}  
TERM.....PAROXYSMAL ATRIAL TACHYCARDIA

DATA VALUE: PAROXYSMAL SUPRAVENTRICULAR TACHYCARDIA

-----  
ABBREVIATION.....PST  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Paroxysmal Supraventricular Tachycardia  
HELP\_TEXT.....  
ID.....{9C910A04-95D2-4EB7-A4CB-B4A74FE7CCEE}  
TERM.....PAROXYSMAL SUPRAVENTRICULAR TACHYCARDIA

DATA VALUE: PART THICKNESS SKIN LOSS INVOLVE EPIDERMIS/DERMIS

-----  
ABBREVIATION.....PART EPIDE  
ACTIVE.....Yes  
DESCRIPTION.....Partial Thickness Skin Loss Involving  
Epidermis/Dermis  
DISPLAY\_NAME.....Partial Thickness Skin Loss Involving  
Epidermis/Dermis  
HELP\_TEXT.....No Help Text  
ID.....{EAD3A260-CC46-4B47-AA49-91CFE13EEE4B}  
TERM.....PART THICKNESS SKIN LOSS INVOLVE  
EPIDERMIS/DERMIS

DATA VALUE: PARTIAL BATH

-----  
ABBREVIATION.....PART. BATH  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Partial Bath  
HELP\_TEXT.....No Help Text  
ID.....{60B9C7AD-9E6E-4412-9BB3-57AFF956C1FC}  
TERM.....PARTIAL BATH



DATA VALUE: PARTIAL IMPAIRED

-----  
ABBREVIATION.....PART IMPAI  
ACTIVE.....Yes  
DESCRIPTION.....Partial Impaired  
DISPLAY\_NAME.....Partial Impaired  
HELP\_TEXT.....Partial Impaired  
ID.....{2F1B8C52-1C94-4EA7-AB82-117694FD4D52}  
TERM.....PARTIAL IMPAIRED

DATA VALUE: PARTIAL LOWER

-----  
ABBREVIATION.....Part Low  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Partial Lower  
HELP\_TEXT.....  
ID.....{2556CE7E-F211-4992-BCAF-09677C20BEFE}  
TERM.....PARTIAL LOWER

DATA VALUE: PARTIAL SIMPLE

-----  
ABBREVIATION.....PART SIMP  
ACTIVE.....Yes  
DESCRIPTION.....Partial Simple  
DISPLAY\_NAME.....Partial Simple  
HELP\_TEXT.....Partial Simple  
ID.....{F4A5E7ED-92AB-4021-A63C-E1D948EFE6C1}  
TERM.....PARTIAL SIMPLE

DATA VALUE: PARTIAL UPPER

-----  
ABBREVIATION.....Part Up  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Partial Upper  
HELP\_TEXT.....  
ID.....{D7B30583-C22C-4B7E-BE7F-36FA9A577DAB}  
TERM.....PARTIAL UPPER

DATA VALUE: PASSIVE RANGE OF MOTION

-----  
ABBREVIATION.....PASS ROM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Passive Range Of Motion  
HELP\_TEXT.....No Help Text  
ID.....{DD43B9EF-E18F-4203-8564-89EBBAD963A4}  
TERM.....PASSIVE RANGE OF MOTION

DATA VALUE: PATENT

-----  
ABBREVIATION.....Patent  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Patent  
HELP\_TEXT.....  
ID.....{F0C38BCB-1D11-436C-8B8C-4297C27DDE04}  
TERM.....PATENT

DATA VALUE: PATENT CLEAN

-----  
ABBREVIATION.....PAT CLEAN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Patent Clean  
HELP\_TEXT.....No Help Text  
ID.....{68688246-AE58-2EFB-BF9D-5E00D8CFDA8E}  
TERM.....PATENT CLEAN

DATA VALUE: PATIENT

-----  
ABBREVIATION.....PATIENT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Patient  
HELP\_TEXT.....No Help Text  
ID.....{7881F205-BCD1-4948-AE66-AD572F5A02F2}  
TERM.....PATIENT

DATA VALUE: PATIENT INFORMED

-----  
ABBREVIATION.....PAT INFORM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Patient Informed  
HELP\_TEXT.....No Help Text  
ID.....{074A1D19-3920-49F1-8958-2EB6F9CF05BA}  
TERM.....PATIENT INFORMED

DATA VALUE: PATIENT INSTRUCTED

-----  
ABBREVIATION.....pat instru  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Patient Instructed  
HELP\_TEXT.....  
ID.....{D89A50A6-352D-4EB8-ABD7-5D732A060D9A}  
TERM.....PATIENT INSTRUCTED

DATA VALUE: PATIENT LOCATOR DEVICE IN PLACE

-----  
ABBREVIATION.....Loc Device  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Patient Locator Device in Place  
HELP\_TEXT.....  
ID.....{CAF7B45E-0143-472D-BABF-9A92857419E9}  
TERM.....PATIENT LOCATOR DEVICE IN PLACE

DATA VALUE: PATIENT REFUSES MARKING

-----  
ABBREVIATION.....Pat Refuse  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Patient Refuses Marking  
HELP\_TEXT.....  
ID.....{93CE893B-AFE1-4047-9EE9-2AD6388C0086}  
TERM.....PATIENT REFUSES MARKING

DATA VALUE: PATIENT REPORTS NOT PRESENT

-----  
ABBREVIATION.....PT RP N PR  
ACTIVE.....Yes  
DESCRIPTION.....Patient Reports Not Present  
DISPLAY\_NAME.....Patient Reports Not Present  
HELP\_TEXT.....Patient Reports Not Present  
ID.....{AD0EF44D-7231-4E8C-A3AE-0E43DD0C44F2}  
TERM.....PATIENT REPORTS NOT PRESENT

DATA VALUE: PATIENT REPORTS PRESENT

-----  
ABBREVIATION.....PT REP PRE  
ACTIVE.....Yes  
DESCRIPTION.....Patient Reports Present  
DISPLAY\_NAME.....Patient Reports Present  
HELP\_TEXT.....Patient Reports Present  
ID.....{1726BC7E-5A82-41AC-B77A-BF62263E1EF6}  
TERM.....PATIENT REPORTS PRESENT

DATA VALUE: PATIENT RIGHTS AND RESPONSIBILITIES

-----  
ABBREVIATION.....Pat Rights  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Patient Rights and Responsibilities  
HELP\_TEXT.....  
ID.....{DA69AAEE-DBDC-4CAA-849A-7242F0854385}  
TERM.....PATIENT RIGHTS AND RESPONSIBILITIES

DATA VALUE: PATIENT SIGNS SAFETY CONTRACT

-----  
ABBREVIATION.....PT CONTRAC  
ACTIVE.....Yes  
DESCRIPTION.....Patient Signs Safety Contract  
DISPLAY\_NAME.....Patient Signs Safety Contract  
HELP\_TEXT.....No Help Text  
ID.....{C65734EC-2CAF-4547-8FD8-CBCD9C706003}  
TERM.....PATIENT SIGNS SAFETY CONTRACT

DATA VALUE: PATIENT STATED

-----  
ABBREVIATION.....Pat Stated  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Patient Stated  
HELP\_TEXT.....  
ID.....{843654E1-4BC9-4277-9EFC-0B388A70A907}  
TERM.....PATIENT STATED

DATA VALUE: PATIENT STAYS WITHIN LINE OF SITE OF NURSING STATION

-----  
ABBREVIATION.....pswlosons  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Patient Stays Within Line Of Site Of  
Nursing Station  
HELP\_TEXT.....No Help Text  
ID.....{CE3E3531-7D9B-A2B4-00B3-00EBCDD09D3D}  
TERM.....PATIENT STAYS WITHIN LINE OF SITE OF  
NURSING STATION

DATA VALUE: PATIENT'S ROOM

-----  
ABBREVIATION.....pr  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Patient's Room  
HELP\_TEXT.....No Help Text  
ID.....{0F13E11D-9F78-14A3-0568-AF26393F39F1}  
TERM.....PATIENT'S ROOM

DATA VALUE: PCA

-----  
ABBREVIATION.....PCA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Patient Controlled Analgesia  
HELP\_TEXT.....No Help-Text  
ID.....{A1EE27C6-F22A-4109-A8E6-28DCCE599E64}  
TERM.....PCA

DATA VALUE: PCA AND BASAL RATE

-----  
ABBREVIATION.....PCA/BASAL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....PCA and Basal Rate  
HELP\_TEXT.....No Help Text  
ID.....{D711BB51-0662-4DD0-8654-C737F6A6686A}  
TERM.....PCA AND BASAL RATE

DATA VALUE: PCA ONLY

-----  
ABBREVIATION.....PCA ONLY  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....PCA Only  
HELP\_TEXT.....No Help Text  
ID.....{30E0441A-1098-4FE9-B277-BF2BBD6B1667}  
TERM.....PCA ONLY

DATA VALUE: PCA PUMP

-----  
ABBREVIATION.....PCA Pump  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....PCA Pump  
HELP\_TEXT.....  
ID.....{C6F139AF-C42B-4577-83E1-A12675C31A21}  
TERM.....PCA PUMP

DATA VALUE: PEDIATRIC

-----  
ABBREVIATION.....PEDIATRIC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Pediatric  
HELP\_TEXT.....No Help Text  
ID.....{06865E22-3ED9-4EDD-9749-6F8641BF29DE}  
TERM.....PEDIATRIC

DATA VALUE: PEG TUBE

-----  
ABBREVIATION.....PEG TUBE  
ACTIVE.....Yes  
DESCRIPTION.....PEG Tube EF Tube  
DISPLAY\_NAME.....PEG Tube  
HELP\_TEXT.....No Help Text  
ID.....{970C3755-2B3F-462F-AA87-EA33F846843F}  
TERM.....PEG TUBE

DATA VALUE: PENETRATING

-----  
ABBREVIATION.....PENETRATIN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Penetrating  
HELP\_TEXT.....No Help Text  
ID.....{00A2D4C8-EDDE-4593-A1F8-A87A836B3379}  
TERM.....PENETRATING

DATA VALUE: PENROSE DRAIN

-----  
ABBREVIATION.....Penrose D  
ACTIVE.....Yes  
DESCRIPTION.....Penrose Drain  
DISPLAY\_NAME.....Penrose Drain  
HELP\_TEXT.....No Help Text  
ID.....{DE508E95-D044-ACCA-EE69-33B9A8454E6C}  
TERM.....PENROSE DRAIN

DATA VALUE: PERATIVE

-----  
ABBREVIATION.....Perative  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Perative  
HELP\_TEXT.....  
ID.....{9F387502-6FE8-4AFA-AF6E-BE4D0CF60947}  
TERM.....PERATIVE

DATA VALUE: PERCUSSION BED

-----  
ABBREVIATION.....PERC BED  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Percussion Bed  
HELP\_TEXT.....No Help Text  
ID.....{BE0B086B-A2F3-446A-ADE0-9589E92689FE}  
TERM.....PERCUSSION BED

DATA VALUE: PERINEAL CARE

-----  
ABBREVIATION.....PERIN CARE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Perineal Care  
HELP\_TEXT.....No Help Text  
ID.....{D88736E4-E0EB-4FF7-9199-288A4D051A20}  
TERM.....PERINEAL CARE

DATA VALUE: PERMANENT

-----  
ABBREVIATION.....p  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Permanent  
HELP\_TEXT.....No Help Text  
ID.....{A31FF319-F8C3-7AAC-D7FC-765C89732D09}  
TERM.....PERMANENT

DATA VALUE: PERSONAL ARTICLES AT REACH

-----  
ABBREVIATION.....At Reach  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Personal Articles at Reach  
HELP\_TEXT.....  
ID.....{166DE5BA-67F1-4BDB-9378-062862519BE4}  
TERM.....PERSONAL ARTICLES AT REACH

DATA VALUE: PERSONAL HYGIENE AND/OR GROOMING

-----  
ABBREVIATION.....Per Hyg  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Personal Hygiene and/or Grooming  
HELP\_TEXT.....  
ID.....{0DDC2293-90E8-4E0B-A6B1-9C5C946AC831}  
TERM.....PERSONAL HYGIENE AND/OR GROOMING

DATA VALUE: PERSONAL MEDICAL ALERT BRACELET

-----  
ABBREVIATION.....pmab  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Personal Medical Alert Bracelet  
HELP\_TEXT.....No Help Text  
ID.....{E63D5FDE-4509-5751-A11C-C023E296858A}  
TERM.....PERSONAL MEDICAL ALERT BRACELET

DATA VALUE: PETROLATUM

-----  
ABBREVIATION.....PETROLATUM  
ACTIVE.....Yes  
DESCRIPTION.....Petrolatum  
DISPLAY\_NAME.....Petrolatum  
HELP\_TEXT.....No Help Text  
ID.....{DA297C75-4776-4C8D-B0D7-A0FF0E4D1B16}  
TERM.....PETROLATUM

DATA VALUE: PHYSICAL IMPAIRMENT

-----  
ABBREVIATION.....Phys Impai  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Physical Impairment  
HELP\_TEXT.....  
ID.....{4D744081-134E-420D-9712-AF114977E490}  
TERM.....PHYSICAL IMPAIRMENT

DATA VALUE: PHYSICIAN

-----  
ABBREVIATION.....PHYSICIA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Physician  
HELP\_TEXT.....No Help Text  
ID.....{DB0C5AE0-459A-4E04-A202-D1CF54B84A6A}  
TERM.....PHYSICIAN

DATA VALUE: PHYSICIAN NOTIFIED

-----  
ABBREVIATION.....PHYS NOTIF  
ACTIVE.....Yes  
DESCRIPTION.....Physician Notified  
DISPLAY\_NAME.....Physician Notified  
HELP\_TEXT.....No Help Text  
ID.....{4BD68949-74F6-42A6-AEF5-E9A33DC31FFC}  
TERM.....PHYSICIAN NOTIFIED

DATA VALUE: PHYSICIAN ORDER

-----  
ABBREVIATION.....PHY ORDER  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Physician Order  
HELP\_TEXT.....No Help Text  
ID.....{AB89E640-83B7-4E6C-8C8D-2FA8C2564CCB}  
TERM.....PHYSICIAN ORDER

DATA VALUE: PICC

-----  
ABBREVIATION.....PICC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....PICC  
HELP\_TEXT.....Peripherally inserted central catheter  
ID.....{3CFE33D0-5A35-88CC-152C-D7501D1858D9}  
TERM.....PICC



DATA VALUE: PIN

-----  
ABBREVIATION.....PIN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Pin  
HELP\_TEXT.....No Help Text  
ID.....{10D48A4E-46F3-4242-BA56-B9D4B37A4C93}  
TERM.....PIN

DATA VALUE: PIN TRACTION

-----  
ABBREVIATION.....PIN TRACT  
ACTIVE.....Yes  
DESCRIPTION.....Pin Traction  
DISPLAY\_NAME.....Pin Traction  
HELP\_TEXT.....Pin Traction  
ID.....{C75104A5-3B6B-4CFE-8E2E-05877522C811}  
TERM.....PIN TRACTION

DATA VALUE: PINK

-----  
ABBREVIATION.....PINK  
ACTIVE.....Yes  
DESCRIPTION.....The color pink  
DISPLAY\_NAME.....Pink  
HELP\_TEXT.....The color pink  
ID.....{1A7FAF84-D69C-4EEB-7BC5-D21A598A1F4D}  
TERM.....PINK

DATA VALUE: PINK / GRANULATING

-----  
ABBREVIATION.....Pink/Gran  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Pink / Granulating  
HELP\_TEXT.....  
ID.....{05245B27-3292-48B2-9F37-6AB9BB92395B}  
TERM.....PINK / GRANULATING

DATA VALUE: PINPOINT

-----  
ABBREVIATION.....pinpoint  
ACTIVE.....Yes  
DESCRIPTION.....Pinpoint  
DISPLAY\_NAME.....Pinpoint  
HELP\_TEXT.....No Help Text  
ID.....{4E53DDAC-0BF4-66D5-D614-9763BBF1FC8E}  
TERM.....PINPOINT

DATA VALUE: PITTED

-----  
ABBREVIATION.....PITTED  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Pitted  
HELP\_TEXT.....No Help Text  
ID.....{BEA42FD8-91C7-4BF6-936E-726D78064F11}  
TERM.....PITTED

DATA VALUE: PITTING

-----  
ABBREVIATION.....PITTING  
ACTIVE.....Yes  
DESCRIPTION.....PITTING  
DISPLAY\_NAME.....Pitting  
HELP\_TEXT.....PITTING  
ID.....{A9F6ACFF-E427-145A-A168-345ACB87C1DE}  
TERM.....PITTING

DATA VALUE: PIVOT

-----  
ABBREVIATION.....Pivot  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Pivot  
HELP\_TEXT.....  
ID.....{2E68ED54-57F5-490D-A982-668121F0CF5A}  
TERM.....PIVOT

DATA VALUE: PIVOT 1.5

-----  
ABBREVIATION.....Pivot 1.5  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Pivot 1.5  
HELP\_TEXT.....  
ID.....{9126F0C6-9716-4C33-8B7C-178D46899D44}  
TERM.....PIVOT 1.5

DATA VALUE: PJC

-----  
ABBREVIATION.....p  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....PJC  
HELP\_TEXT.....No Help Text  
ID.....{1D2D1FEA-2006-3C6C-FAD4-72384F2033C9}  
TERM.....PJC

DATA VALUE: PLACED NEAR TO NURSE STATION

-----  
ABBREVIATION.....Near N/Sta  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Placed Near to Nurse Station  
HELP\_TEXT.....  
ID.....{7790D901-36F0-4DFD-9BF5-8C3A913F30E1}  
TERM.....PLACED NEAR TO NURSE STATION

DATA VALUE: PLAIN

-----  
ABBREVIATION.....PLAIN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Plain Catheter  
HELP\_TEXT.....No Help Text  
ID.....{483EE48F-F72B-A917-DCD0-E8B69A3FA36F}  
TERM.....PLAIN

DATA VALUE: POOR

-----  
ABBREVIATION.....POOR  
ACTIVE.....Yes  
DESCRIPTION.....Poor  
DISPLAY\_NAME.....Poor  
HELP\_TEXT.....No Help Text  
ID.....{A3D75A12-A791-446D-AF63-9435054DFC4C}  
TERM.....POOR

DATA VALUE: POOR BAL

-----  
ABBREVIATION.....Poor Balan  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Poor Balance  
HELP\_TEXT.....  
ID.....{959AFDB4-6984-452D-9235-B154409DC0ED}  
TERM.....POOR BAL

DATA VALUE: PORCINE SUBMUCOSA

-----  
ABBREVIATION.....PORC SUBMU  
ACTIVE.....Yes  
DESCRIPTION.....Oasis (porcine submucosa)  
DISPLAY\_NAME.....Oasis (porcine submucosa)  
HELP\_TEXT.....No Help Text  
ID.....{6F3917DA-74BD-4790-9D4D-B284E27B7623}  
TERM.....PORCINE SUBMUCOSA

DATA VALUE: POSITIONAL

-----  
ABBREVIATION.....POSITIONAL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Positional  
HELP\_TEXT.....No Help Text  
ID.....{9A74D16B-5E95-4EA8-88AF-E53C80311A5D}  
TERM.....POSITIONAL

DATA VALUE: POSITIVE PRESSURE CONNECTOR

-----  
ABBREVIATION.....P Pres Con  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Positive Pressure Connector  
HELP\_TEXT.....  
ID.....{1A662A2F-4666-4CFD-9E5C-7C5115BB7C2D}  
TERM.....POSITIVE PRESSURE CONNECTOR

DATA VALUE: POST CATH DIET

-----  
ABBREVIATION.....POST CA DT  
ACTIVE.....Yes  
DESCRIPTION.....Post Cath Diet  
DISPLAY\_NAME.....Post Cath Diet  
HELP\_TEXT.....No Help Text  
ID.....{B18FEC50-18E1-4EED-8B37-E7DBAA7D4D3B}  
TERM.....POST CATH DIET

DATA VALUE: POST OPERATION DIAGNOSIS

-----  
ABBREVIATION.....PO OP DIAG  
ACTIVE.....Yes  
DESCRIPTION.....Post Operation Diagnosis  
DISPLAY\_NAME.....Post Operation Diagnosis  
HELP\_TEXT.....Post Operation Diagnosis  
ID.....{A69379C9-C54A-4C55-85B3-C3D92D1705A0}  
TERM.....POST OPERATION DIAGNOSIS

DATA VALUE: POST PROCEDURE DIAGNOSIS

-----  
ABBREVIATION.....PO PRO DIA  
ACTIVE.....Yes  
DESCRIPTION.....Post Procedure Diagnosis  
DISPLAY\_NAME.....Post Procedure Diagnosis  
HELP\_TEXT.....Post Procedure Diagnosis  
ID.....{933C7D2C-88F8-4C18-ACBB-A16A295D1B3A}  
TERM.....POST PROCEDURE DIAGNOSIS

DATA VALUE: POST-FALL

-----  
ABBREVIATION.....POST-FALL  
ACTIVE.....Yes  
DESCRIPTION.....Post-Fall  
DISPLAY\_NAME.....Post-Fall  
HELP\_TEXT.....No Help Text  
ID.....{9BC82B93-CA4A-44E3-A3BA-7AF72BAFD44B}  
TERM.....POST-FALL

DATA VALUE: POST-OPERATIVE

-----  
ABBREVIATION.....POST-OP  
ACTIVE.....Yes  
DESCRIPTION.....Post-Operative  
DISPLAY\_NAME.....Post-Operative  
HELP\_TEXT.....No Help Text  
ID.....{9FA19656-849F-4A65-92FC-755A4C9E8056}  
TERM.....POST-OPERATIVE

DATA VALUE: POST-PROCEDURE

-----  
ABBREVIATION.....POST-PROC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Post-Procedure  
HELP\_TEXT.....No Help Text  
ID.....{0F943BD7-BD13-412E-B4EA-9B5526587340}  
TERM.....POST-PROCEDURE

DATA VALUE: POST-SEDATION

-----  
ABBREVIATION.....POST SED  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Post-Sedation  
HELP\_TEXT.....No Help Text  
ID.....{B829A249-CB14-4919-86F8-4A8C8D6B7753}  
TERM.....POST-SEDATION

DATA VALUE: POTASSIUM MODIFIED 60 MEQ

-----  
ABBREVIATION.....K MOD 60M  
ACTIVE.....Yes  
DESCRIPTION.....Potassium Modified 60 MEQ  
DISPLAY\_NAME.....Potassium Modified 60 MEQ  
HELP\_TEXT.....No Help Text  
ID.....{6EC2F10D-CC44-4872-8473-2445A835ABB0}  
TERM.....POTASSIUM MODIFIED 60 MEQ

DATA VALUE: PRBC

-----  
ABBREVIATION.....PRBC  
ACTIVE.....Yes  
DESCRIPTION.....Packed Red Blood Cells  
DISPLAY\_NAME.....PRBC  
HELP\_TEXT.....No Help Text  
ID.....{F6D9B889-1BCB-4724-B75A-E98743585D55}  
TERM.....PRBC

DATA VALUE: PRE-MEDICATED FOR ACTIVITY

-----  
ABBREVIATION.....PRE-MED AC  
ACTIVE.....Yes  
DESCRIPTION.....Pre-Medicated For Activity  
DISPLAY\_NAME.....Pre-Medicated For Activity  
HELP\_TEXT.....No Help Text  
ID.....{D6DA772A-FA22-4170-B928-855A2F68711F}  
TERM.....PRE-MEDICATED FOR ACTIVITY

DATA VALUE: PRE-OPERATIVE

-----  
ABBREVIATION.....PRE-OP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Pre-Operative  
HELP\_TEXT.....No Help Text  
ID.....{80C418D5-54B1-4CC3-9B19-5653B8D73C45}  
TERM.....PRE-OPERATIVE

DATA VALUE: PRE-PROCEDURE

-----  
ABBREVIATION.....PRE-PROC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Pre-Procedure  
HELP\_TEXT.....No Help Text  
ID.....{D2A7B025-9AB9-4C69-B163-D5EBB78B82BE}  
TERM.....PRE-PROCEDURE

DATA VALUE: PRE-SEDATION

-----  
ABBREVIATION.....PRE SEDA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Pre-Sedation  
HELP\_TEXT.....No Help Text  
ID.....{684C764E-2BE9-4C99-A2B2-C92D75C73751}  
TERM.....PRE-SEDATION

DATA VALUE: PRE-SURGERY

-----  
ABBREVIATION.....PRE-SURG  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Pre-Surgery  
HELP\_TEXT.....No Help Text  
ID.....{0C71761F-3024-4638-A95C-FDF12DC03CF5}  
TERM.....PRE-SURGERY

DATA VALUE: PREGNANT, NO

-----  
ABBREVIATION.....Preg, No  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Pregnant, No  
HELP\_TEXT.....  
ID.....{D08F5D48-965E-4537-83BC-14722E28B0AE}  
TERM.....PREGNANT, NO

DATA VALUE: PREGNANT, YES

-----  
ABBREVIATION.....Preg, Yes  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Pregnant, Yes  
HELP\_TEXT.....  
ID.....{CF4323D7-9C8C-4FA5-9AB5-1AE4FFD1B535}  
TERM.....PREGNANT, YES

DATA VALUE: PREMATURE ATRIAL

-----  
ABBREVIATION.....pa  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Premature Atrial  
HELP\_TEXT.....No Help Text  
ID.....{070FDCB1-2EC4-B4EF-4DFD-81F6CCFFBDED}  
TERM.....PREMATURE ATRIAL

DATA VALUE: PRESENT

-----  
ABBREVIATION.....p  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Present  
HELP\_TEXT.....No Help Text  
ID.....{384DC20F-1625-CE1E-1C85-A3B3924BB2A0}  
TERM.....PRESENT

DATA VALUE: PRESSURE

-----  
ABBREVIATION.....PRESSURE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Pressure  
HELP\_TEXT.....No Help Text  
ID.....{387AFB7F-3973-4D77-B2E6-0724F344D8F6}  
TERM.....PRESSURE

DATA VALUE: PRESSURE CALIBRATION

-----  
ABBREVIATION.....PRESS CALI  
ACTIVE.....Yes  
DESCRIPTION.....Pressure Calibration  
DISPLAY\_NAME.....Pressure Calibration  
HELP\_TEXT.....No Help Text  
ID.....{38C12E69-2CD2-4F61-98D9-06AB857DA51A}  
TERM.....PRESSURE CALIBRATION

DATA VALUE: PRESSURE CONTROL

-----  
ABBREVIATION.....PRESS CONT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Pressure Control  
HELP\_TEXT.....No Help Text  
ID.....{C3CDF756-484F-1483-9883-9B63DEF3B012}  
TERM.....PRESSURE CONTROL

DATA VALUE: PRESSURE DRESSING

-----  
ABBREVIATION.....PRESS DRES  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Pressure Dressing  
HELP\_TEXT.....No Help Text  
ID.....{6D9DB9EF-85EB-8321-03FC-879AA4C5615B}  
TERM.....PRESSURE DRESSING

DATA VALUE: PRESSURE ULCER

-----  
ABBREVIATION.....PRE ULCER  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Pressure Ulcer  
HELP\_TEXT.....No Help Text  
ID.....{7279639E-498D-4EA3-A3FA-4A6D611E6F48}  
TERM.....PRESSURE ULCER



DATA VALUE: PRESYSTOLIC

-----  
ABBREVIATION.....PRESYSTOLI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Presystolic  
HELP\_TEXT.....No Help Text  
ID.....{BB85B8EE-36F8-4AC2-9AB0-0BABCFA9FCBB}  
TERM.....PRESYSTOLIC

DATA VALUE: PREVENT ACCESS TO DEVICE

-----  
ABBREVIATION.....PREV ACCE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Prevent Access To Medical Devices  
HELP\_TEXT.....No Help Text  
ID.....{E7455083-DD2C-4562-A6E2-9626A4CDC4DD}  
TERM.....PREVENT ACCESS TO DEVICE

DATA VALUE: PREVENTIVE HEALTH REHAB TECHNIQUES

-----  
ABBREVIATION.....Rehab Tech  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Preventive Health Rehab Techniques  
HELP\_TEXT.....  
ID.....{2BCABB3F-6688-4BA8-A9E2-D6CA5E38BB7B}  
TERM.....PREVENTIVE HEALTH REHAB TECHNIQUES

DATA VALUE: PREVIOUS PERCUTANEOUS SITE

-----  
ABBREVIATION.....PRE PER SI  
ACTIVE.....Yes  
DESCRIPTION.....Previous Percutaneous Site  
DISPLAY\_NAME.....Previous Percutaneous Site  
HELP\_TEXT.....No Help Text  
ID.....{F2ECCF5C-D91A-40F5-9C09-452C2C81DA66}  
TERM.....PREVIOUS PERCUTANEOUS SITE

DATA VALUE: PRN

-----  
ABBREVIATION.....PRN  
ACTIVE.....Yes  
DESCRIPTION.....As Needed  
DISPLAY\_NAME.....PRN  
HELP\_TEXT.....As Needed  
ID.....{36493D81-954B-4C56-B0A2-9988EA785EBD}  
TERM.....PRN

DATA VALUE: PRO STAT SUGAR FREE 64

-----  
ABBREVIATION.....ProSta 64  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Pro Stat Sugar Free 64  
HELP\_TEXT.....  
ID.....{3622E0A3-8BDF-4ADA-AC4D-20762C10C798}  
TERM.....PRO STAT SUGAR FREE 64

DATA VALUE: PRO-STAT

-----  
ABBREVIATION.....Pro-Stat  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Pro-Stat  
HELP\_TEXT.....  
ID.....{E39E36C7-0A47-459D-BA75-CF6186247B40}  
TERM.....PRO-STAT

DATA VALUE: PROCEDURE ROOM

-----  
ABBREVIATION.....PROC ROO  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Procedure Room  
HELP\_TEXT.....No Help Text  
ID.....{CC9A2AB0-7A73-4DE7-9B7F-DC2C160AE462}  
TERM.....PROCEDURE ROOM

DATA VALUE: PRODUCTIVE

-----  
ABBREVIATION.....p  
ACTIVE.....Yes  
DESCRIPTION.....Produces desired result  
DISPLAY\_NAME.....Productive  
HELP\_TEXT.....No Help Text  
ID.....{22BE0499-F9FE-3A27-2605-61EBEC0ED6AB}  
TERM.....PRODUCTIVE

DATA VALUE: PROJECTILE VOMITING

-----  
ABBREVIATION.....PROJ VOMIT  
ACTIVE.....Yes  
DESCRIPTION.....Projectile Vomiting  
DISPLAY\_NAME.....Projectile Vomiting  
HELP\_TEXT.....Projectile Vomiting  
ID.....{0D2CB23B-B266-40E1-A6FF-20B733655B69}  
TERM.....PROJECTILE VOMITING

DATA VALUE: PROMOTE

-----  
ABBREVIATION.....Promote  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Promote  
HELP\_TEXT.....  
ID.....{6377C8D3-F6DD-4BD2-877F-DB52525DD65D}  
TERM.....PROMOTE

DATA VALUE: PROMOTE WITH FIBER

-----  
ABBREVIATION.....Pro W/Fib  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Promote With Fiber  
HELP\_TEXT.....  
ID.....{33E04B6F-2976-402C-B3B7-B75705430340}  
TERM.....PROMOTE WITH FIBER

DATA VALUE: PRONE

-----  
ABBREVIATION.....PRONE  
ACTIVE.....Yes  
DESCRIPTION.....Prone  
DISPLAY\_NAME.....Prone  
HELP\_TEXT.....No Help Text  
ID.....{41180E09-D9D6-4903-BAE9-28CD5AAC4991}  
TERM.....PRONE

DATA VALUE: PROPASS PROTEIN POWDER

-----  
ABBREVIATION.....Pro Prot  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Propass Protein Powder  
HELP\_TEXT.....  
ID.....{2307C41A-92D2-4DD4-BC1F-2E78F7477D98}  
TERM.....PROPASS PROTEIN POWDER

DATA VALUE: PROSTHESIS

-----  
ABBREVIATION.....PROSTHESIS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Prosthesis  
HELP\_TEXT.....No Help Text  
ID.....{CCD0EF23-FE9C-4B20-942E-9CF1CD7CF6E9}  
TERM.....PROSTHESIS

DATA VALUE: PROTENEX WOUND CARE

-----  
ABBREVIATION.....Protenex  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Protenex Wound Care  
HELP\_TEXT.....  
ID.....{3B5ED9D2-3DBC-42FF-ACC4-3594EFE33627}  
TERM.....PROTENEX WOUND CARE

DATA VALUE: PROVIDED PATIENT EDUCATION

-----  
ABBREVIATION.....Pat Ed  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Provided Patient Education  
HELP\_TEXT.....  
ID.....{25D0998B-3600-40C9-8608-9A0D91CCBEE0}  
TERM.....PROVIDED PATIENT EDUCATION

DATA VALUE: PROXIMAL PORT

-----  
ABBREVIATION.....PROX PORT  
ACTIVE.....Yes  
DESCRIPTION.....Proximal Port  
DISPLAY\_NAME.....Proximal Port  
HELP\_TEXT.....No Help Text  
ID.....{3244EF4F-F18A-4250-96F5-ABCEA7BFA578}  
TERM.....PROXIMAL PORT

DATA VALUE: PT DENIES NAUSEA / NO VOMITING

-----  
ABBREVIATION.....DEN NAUS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Pt. Denies Nausea / No Vomiting  
HELP\_TEXT.....No Help Text  
ID.....{4374D52E-9D4E-41AD-8B52-FCBDF6B5A1E9}  
TERM.....PT DENIES NAUSEA / NO VOMITING

DATA VALUE: PT PAIN AT ACCEPTABLE LEVEL

-----  
ABBREVIATION.....PT PA OK  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Patient's Pain at Acceptable Level  
HELP\_TEXT.....No Help Text  
ID.....{78489498-4A5B-41C5-98C5-95D56C4A0933}  
TERM.....PT PAIN AT ACCEPTABLE LEVEL

DATA VALUE: PULMOCARE

-----  
ABBREVIATION.....Pulmocare  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Pulmocare  
HELP\_TEXT.....  
ID.....{8D683106-B5EE-46E1-B1C6-0C7843EDD15C}  
TERM.....PULMOCARE

DATA VALUE: PULMONARY PROCEDURE ROOM

-----  
ABBREVIATION.....Pul Prc Rm  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Pulmonary Procedure Room  
HELP\_TEXT.....  
ID.....{D082F661-EBA4-4EAA-987C-67E2BDFBE94D}  
TERM.....PULMONARY PROCEDURE ROOM

DATA VALUE: PULSE OXIMETER

-----  
ABBREVIATION.....Pulse Ox  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Pulse Oximeter  
HELP\_TEXT.....  
ID.....{773A956B-0B08-4D72-AAD4-630DF8268694}  
TERM.....PULSE OXIMETER

DATA VALUE: PULSELESS ELECTRICAL ACTIVITY

-----  
ABBREVIATION.....Pulseless  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Pulseless Electrical Activity  
HELP\_TEXT.....  
ID.....{8E21E500-78D4-462D-9483-DD81D3887EE5}  
TERM.....PULSELESS ELECTRICAL ACTIVITY

DATA VALUE: PUMP - CONTINUOUS

-----  
ABBREVIATION.....Pump- cont  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Pump- Continuous  
HELP\_TEXT.....  
ID.....{EEE9625C-7EED-4A31-8C13-5E3F0996EDC8}  
TERM.....PUMP - CONTINUOUS

DATA VALUE: PUMP - INTERMITTENT

-----  
ABBREVIATION.....Pump - int  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Pump - Intermittent  
HELP\_TEXT.....  
ID.....{1D0AB4E0-5DF7-446C-BF2C-8336E7C28FCA}  
TERM.....PUMP - INTERMITTENT

DATA VALUE: PUNCTURE

-----  
ABBREVIATION.....PUNCTURE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Puncture  
HELP\_TEXT.....No Help Text  
ID.....{28CF34AF-BD37-4BA4-A59B-C4AFCE2EB8DD}  
TERM.....PUNCTURE

DATA VALUE: PUREED

-----  
ABBREVIATION.....PUREED  
ACTIVE.....Yes  
DESCRIPTION.....Pureed  
DISPLAY\_NAME.....Pureed  
HELP\_TEXT.....No Help Text  
ID.....{99F8AE5F-E390-4583-8A92-CAC13E21BCF4}  
TERM.....PUREED

DATA VALUE: PURULENT

-----  
ABBREVIATION.....PURULENT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Purulent  
HELP\_TEXT.....No Help Text  
ID.....{86C01494-743A-C6A2-8883-01CA5200C79B}  
TERM.....PURULENT

DATA VALUE: PVC

-----  
ABBREVIATION.....p  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....PVC  
HELP\_TEXT.....No Help Text  
ID.....{33FF92F1-84CC-7CFA-41CF-B271A26C59B8}  
TERM.....PVC

DATA VALUE: QUADRAGEMINY

-----  
ABBREVIATION.....Quadragemi  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Quadrageminy  
HELP\_TEXT.....  
ID.....{899BBEBB-D103-4864-BBE5-C9DF60F1C5BB}  
TERM.....QUADRAGEMINY

DATA VALUE: QUIET OR RESTING

-----  
ABBREVIATION.....qor  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Quiet Or Resting  
HELP\_TEXT.....No Help Text  
ID.....{76E4908F-014C-2088-FF7F-33EBF223DE40}  
TERM.....QUIET OR RESTING

DATA VALUE: QUIET ROOM

-----  
ABBREVIATION.....qr  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Quiet Room  
HELP\_TEXT.....No Help Text  
ID.....{CCD96D41-C273-F6F3-DA18-8818B1CC9980}  
TERM.....QUIET ROOM

DATA VALUE: R BUNDLE BRANCH BLOCK

-----  
ABBREVIATION.....rbbb  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Bundle Branch Block  
HELP\_TEXT.....No Help Text  
ID.....{FDE0D174-5220-F069-D284-20399D9489B0}  
TERM.....R BUNDLE BRANCH BLOCK

DATA VALUE: R FOOT DOWN

-----  
ABBREVIATION.....R FOOT DOW  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Foot Down  
HELP\_TEXT.....No Help Text  
ID.....{DBCC0BA2-3097-42C7-AE9E-DF71FA779162}  
TERM.....R FOOT DOWN

DATA VALUE: R FOOT UP

-----  
ABBREVIATION.....R FOOT UP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Foot Up  
HELP\_TEXT.....No Help Text  
ID.....{5928B710-6EB1-4183-93CB-6C72B37A97C3}  
TERM.....R FOOT UP

DATA VALUE: R HEAD DOWN

-----  
ABBREVIATION.....R HEAD DOW  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Head Down  
HELP\_TEXT.....No Help Text  
ID.....{FDFE0638-444C-4023-B586-985E84EEA1FA}  
TERM.....R HEAD DOWN

DATA VALUE: R HEAD UP

-----  
ABBREVIATION.....R HEAD UP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Head Up  
HELP\_TEXT.....No Help Text  
ID.....{D03DBEC6-F4DD-441A-AAB2-653CC34FCDBF}  
TERM.....R HEAD UP

DATA VALUE: R NECK

-----  
ABBREVIATION.....RIGHT NECK  
ACTIVE.....Yes  
DESCRIPTION.....RIGHT NECK  
DISPLAY\_NAME.....R Neck  
HELP\_TEXT.....RIGHT NECK  
ID.....{907E90CD-87BF-073D-AF76-EA11ABE7FC3C}  
TERM.....R NECK

DATA VALUE: R NEPHROSTOMY

-----  
ABBREVIATION.....R NEPHROST  
ACTIVE.....Yes  
DESCRIPTION.....Right Nephrostomy  
DISPLAY\_NAME.....R Nephrostomy  
HELP\_TEXT.....Right Nephrostomy  
ID.....{33661137-688D-5611-4368-B1DFD8DB1469}  
TERM.....R NEPHROSTOMY



DATA VALUE: R STENT

-----  
ABBREVIATION.....RT. STENT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Stent  
HELP\_TEXT.....No Help Text  
ID.....{6D0EED3B-1F08-11D2-165D-97A9F58327D5}  
TERM.....R STENT

DATA VALUE: RAD TECH

-----  
ABBREVIATION.....Rad Tech  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Rad Tech  
HELP\_TEXT.....  
ID.....{EEA2A91A-DD5A-4062-9DCD-6A909220ABD9}  
TERM.....RAD TECH

DATA VALUE: RADIATING

-----  
ABBREVIATION.....RADIATING  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....RADIATING  
HELP\_TEXT.....No Help Text  
ID.....{B6807939-F809-45C1-912A-ADEF2E0DB0DC}  
TERM.....RADIATING

DATA VALUE: RADIOLOGY

-----  
ABBREVIATION.....Radiology  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Radiology  
HELP\_TEXT.....  
ID.....{AF29AA13-EAA2-4B72-AA9E-994B5CCC4705}  
TERM.....RADIOLOGY

DATA VALUE: RAMSEY ANXIOUS

-----  
ABBREVIATION.....RAM ANXIOU  
ACTIVE.....Yes  
DESCRIPTION.....Anxious, agitated, or restless  
DISPLAY\_NAME.....1 - Anxious  
HELP\_TEXT.....Anxious, agitated, or restless  
ID.....{88205A0C-7065-4616-8FCF-75059462DD15}  
TERM.....RAMSEY ANXIOUS

DATA VALUE: RAMSEY ASLEEP BRISK RESPONSE

-----  
ABBREVIATION.....RAM ASL BR  
ACTIVE.....Yes  
DESCRIPTION.....Asleep; brisk response to light  
                  glabellar tap or loud auditory  
                  stimulus  
DISPLAY\_NAME.....4 - Asleep; Brisk Response  
HELP\_TEXT.....Asleep; brisk response to light  
                  glabellar tap or loud auditory  
                  stimulus  
ID.....{00F4A117-F9E7-4B60-A521-F2E8725AECD3}  
TERM.....RAMSEY ASLEEP BRISK RESPONSE

DATA VALUE: RAMSEY ASLEEP NO RESPONSE

-----  
ABBREVIATION.....RAM ASL NO  
ACTIVE.....Yes  
DESCRIPTION.....Asleep; no response to painful stimuli  
DISPLAY\_NAME.....6 - Asleep No Response  
HELP\_TEXT.....Asleep; no response to painful stimuli  
ID.....{3FF5BA0A-2D3C-4935-A432-23DAA4761AB7}  
TERM.....RAMSEY ASLEEP NO RESPONSE

DATA VALUE: RAMSEY ASLEEP SLUGGISH

-----  
ABBREVIATION.....RAM ASL SL  
ACTIVE.....Yes  
DESCRIPTION.....Asleep; sluggish response to light  
                  glabellar tap or loud auditory  
                  stimulus  
DISPLAY\_NAME.....5 - Asleep Sluggish Response  
HELP\_TEXT.....Asleep; sluggish response to light  
                  glabellar tap or loud auditory  
                  stimulus  
ID.....{AA61604F-00F9-4103-B27F-E817FBC9577B}  
TERM.....RAMSEY ASLEEP SLUGGISH

DATA VALUE: RAMSEY COOPERATIVE

-----  
ABBREVIATION.....RAM COOPER  
ACTIVE.....Yes  
DESCRIPTION.....Cooperative, oriented, and tranquil  
DISPLAY\_NAME.....2 - Cooperative  
HELP\_TEXT.....Cooperative, oriented, and tranquil  
ID.....{56846ADA-0F0E-44AA-9779-A567B9D24C4F}  
TERM.....RAMSEY COOPERATIVE

DATA VALUE: RAMSEY SEDATED

-----  
ABBREVIATION.....RAM SEDATE  
ACTIVE.....Yes  
DESCRIPTION.....Sedated but responds to commands  
DISPLAY\_NAME.....3 - Sedated But Responds  
HELP\_TEXT.....Sedated but responds to commands  
ID.....{CC12405F-60DB-485F-94E5-78C1EB32AB74}  
TERM.....RAMSEY SEDATED



## DATA VALUE: RASS DROWSY

-----  
 ABBREVIATION.....RASS DROWS  
 ACTIVE.....Yes  
 DESCRIPTION.....Not fully alert, but has sustained  
                   awakening (eye-opening/eye contact)  
                   to voice (>10 seconds)  
 DISPLAY\_NAME.....-1: Drowsy  
 HELP\_TEXT.....Not fully alert, but has sustained  
                   awakening (eye-opening/eye contact)  
                   to voice (>10 seconds)  
 ID.....{5E6BDAC6-4F81-4110-A0AF-FED746D9FF0D}  
 TERM.....RASS DROWSY

## DATA VALUE: RASS LIGHT SEDATION

-----  
 ABBREVIATION.....RASS LT SE  
 ACTIVE.....Yes  
 DESCRIPTION.....Briefly awakens with eye contact to  
                   voice (<10 seconds)  
 DISPLAY\_NAME.....-2: Light Sedation  
 HELP\_TEXT.....Briefly awakens with eye contact to  
                   voice (<10 seconds)  
 ID.....{FDEEC084-3922-4065-A863-4893E42F8F93}  
 TERM.....RASS LIGHT SEDATION

## DATA VALUE: RASS MODERATE SEDATION

-----  
 ABBREVIATION.....RASS MD SE  
 ACTIVE.....Yes  
 DESCRIPTION.....Movement or eye opening to voice (but  
                   no eye contact)  
 DISPLAY\_NAME.....-3: Moderate Sedation  
 HELP\_TEXT.....Movement or eye opening to voice (but  
                   no eye contact)  
 ID.....{D12B0A3E-B80E-494F-B89F-A79E8B81B080}  
 TERM.....RASS MODERATE SEDATION

## DATA VALUE: RASS RESTLESS

-----  
 ABBREVIATION.....RASS RESTL  
 ACTIVE.....Yes  
 DESCRIPTION.....Anxious but movements not aggressive  
                   vigorous  
 DISPLAY\_NAME.....+1: Restless  
 HELP\_TEXT.....Anxious but movements not aggressive  
                   vigorous  
 ID.....{2CCFC1ED-D807-45F0-AE69-980E214AD5C3}  
 TERM.....RASS RESTLESS

## DATA VALUE: RASS UNAROUSABLE

-----  
 ABBREVIATION.....RASS UNARO  
 ACTIVE.....Yes  
 DESCRIPTION.....No response to voice or physical  
                   stimulation  
 DISPLAY\_NAME.....-5: Unarousable  
 HELP\_TEXT.....No response to voice or physical  
                   stimulation  
 ID.....{1578E147-7D2C-4469-828E-FCA97109B27C}

TERM.....RASS UNAROUSABLE

DATA VALUE: RASS VERY AGITATED

-----  
ABBREVIATION.....RASS V AGI  
ACTIVE.....Yes  
DESCRIPTION.....Pulls or removes tube(s) or  
catheter(s); aggressive  
DISPLAY\_NAME.....+3: Very Agitated  
HELP\_TEXT.....Pulls or removes tube(s) or  
catheter(s); aggressive  
ID.....{7528CA1D-A9A4-4F4E-8943-E69F3B42E189}  
TERM.....RASS VERY AGITATED

DATA VALUE: REACTIVE

-----  
ABBREVIATION.....Reactive  
ACTIVE.....Yes  
DESCRIPTION.....Reactive  
DISPLAY\_NAME.....Reactive  
HELP\_TEXT.....No Help Text  
ID.....{31B44ADB-0000-A5B0-5F56-59D4A08C05D2}  
TERM.....REACTIVE

DATA VALUE: RECTAL TUBE

-----  
ABBREVIATION.....RECT TUBE  
ACTIVE.....Yes  
DESCRIPTION.....Rectal Tube  
DISPLAY\_NAME.....Rectal Tube  
HELP\_TEXT.....No Help Text  
ID.....{7E492AE5-D3B8-44EC-8841-97770FAD3AFE}  
TERM.....RECTAL TUBE

DATA VALUE: RED

-----  
ABBREVIATION.....RED  
ACTIVE.....Yes  
DESCRIPTION.....The color red  
DISPLAY\_NAME.....Red  
HELP\_TEXT.....The color red  
ID.....{515D7CFB-3527-30E7-EC4B-B607C272DFA0}  
TERM.....RED

DATA VALUE: RED/GRANULATING

-----  
ABBREVIATION.....RED/GRANUL  
ACTIVE.....Yes  
DESCRIPTION.....Red / Granulating  
DISPLAY\_NAME.....Red / Granulating  
HELP\_TEXT.....No Help Text  
ID.....{BB69D2A2-7807-490A-B21D-77B6677796F6}  
TERM.....RED/GRANULATING

DATA VALUE: REDDENING OF THE SKIN

-----  
ABBREVIATION.....RED O SKIN  
ACTIVE.....Yes  
DESCRIPTION.....Reddening Of The Skin  
DISPLAY\_NAME.....Reddening Of The Skin  
HELP\_TEXT.....No Help Text  
ID.....{716AEC57-0835-4D07-A78C-D8453599A3D2}  
TERM.....REDDENING OF THE SKIN

DATA VALUE: REDIRECT PATIENT ACTIVITY

-----  
ABBREVIATION.....REDIR ACT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Redirect Patient Activity  
HELP\_TEXT.....No Help Text  
ID.....{BFF0ABB1-953A-4275-9EF2-08BFF660AA8D}  
TERM.....REDIRECT PATIENT ACTIVITY

DATA VALUE: REFLEX

-----  
ABBREVIATION.....r  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Reflex  
HELP\_TEXT.....No Help Text  
ID.....{F21585F9-5206-A6B3-2BF6-F416E5842366}  
TERM.....REFLEX

DATA VALUE: REFUSE MEDICATION

-----  
ABBREVIATION.....REFUSE MED  
ACTIVE.....Yes  
DESCRIPTION.....Refuse Medication  
DISPLAY\_NAME.....Refuse Medication  
HELP\_TEXT.....No Help Text  
ID.....{17D6956B-4CF6-492C-B113-B5A93D9E215E}  
TERM.....REFUSE MEDICATION

DATA VALUE: REFUSE TO MOVE

-----  
ABBREVIATION.....REF TO MOV  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Refuse To Move  
HELP\_TEXT.....No Help Text  
ID.....{E14588CE-2498-4FE8-AB0B-E954DC1696AC}  
TERM.....REFUSE TO MOVE

DATA VALUE: REFUSED

-----  
ABBREVIATION.....REFUSED  
ACTIVE.....Yes  
DESCRIPTION.....Refused  
DISPLAY\_NAME.....Refused  
HELP\_TEXT.....No Help Text  
ID.....{D5B3E135-4D68-4528-BBC9-6E8A7A997104}  
TERM.....REFUSED

DATA VALUE: REGIONAL BLOCK

-----  
ABBREVIATION.....Region Blk  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Regional Block  
HELP\_TEXT.....  
ID.....{3F829DFC-D3F6-4D08-AACD-A9A29ED82F19}  
TERM.....REGIONAL BLOCK

DATA VALUE: REGISTERED NURSE

-----  
ABBREVIATION.....RN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....RN  
HELP\_TEXT.....No Help Text  
ID.....{524F344F-FDCD-423F-AD34-291D5DB4271B}  
TERM.....REGISTERED NURSE

DATA VALUE: REGULAR

-----  
ABBREVIATION.....REGULAR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Regular  
HELP\_TEXT.....No Help Text  
ID.....{3D3F3D4C-ECD0-2056-A0D1-43CC2F2BA2C8}  
TERM.....REGULAR

DATA VALUE: REINFORCED

-----  
ABBREVIATION.....REINFORCED  
ACTIVE.....Yes  
DESCRIPTION.....Dressing has been reinforced  
DISPLAY\_NAME.....Reinforced  
HELP\_TEXT.....No Help Text  
ID.....{EC524C50-0CDE-4AA1-8059-E4323850575E}  
TERM.....REINFORCED



DATA VALUE: RELATIVE

-----  
ABBREVIATION.....Relative  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Relative  
HELP\_TEXT.....  
ID.....{F4EEB0A5-D25A-45DC-A5AA-BBB7ABA08750}  
TERM.....RELATIVE

DATA VALUE: RELIGIOUS

-----  
ABBREVIATION.....RELIGIOUS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Religious  
HELP\_TEXT.....No Help Text  
ID.....{936C9949-8A30-4BBF-934C-1F194551B8E2}  
TERM.....RELIGIOUS

DATA VALUE: REMNANT P-WAVE

-----  
ABBREVIATION.....Remnant P  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Remnant P-Wave  
HELP\_TEXT.....  
ID.....{3BBAA695-337E-41C4-BB83-8452917C07DF}  
TERM.....REMNANT P-WAVE

DATA VALUE: REMOVE FROM STIM. SETTING

-----  
ABBREVIATION.....REM SET  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Remove From Stimulating Setting  
HELP\_TEXT.....No Help Text  
ID.....{2010BA89-B399-4A30-9ACF-C65AB57B6351}  
TERM.....REMOVE FROM STIM. SETTING

DATA VALUE: REMOVE RECTAL TUBE

-----  
ABBREVIATION.....Rectal Tub  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Remove Rectal Tube  
HELP\_TEXT.....  
ID.....{3D3D75CD-994E-49EE-B51A-8621FF5CA4D8}  
TERM.....REMOVE RECTAL TUBE

DATA VALUE: RENAL DIET

-----  
ABBREVIATION.....REN DIET  
ACTIVE.....Yes  
DESCRIPTION.....Renal Diet  
DISPLAY\_NAME.....Renal Diet  
HELP\_TEXT.....No Help Text  
ID.....{089DCCC2-B12B-4E3B-8D12-FCA5E974A55A}  
TERM.....RENAL DIET

DATA VALUE: REORIENTED TO ENVIRONMENT

-----  
ABBREVIATION.....REORIENT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Reoriented To Environment  
HELP\_TEXT.....No Help Text  
ID.....{4C47B764-4BA7-4D71-9B0B-D42F120E9A4A}  
TERM.....REORIENTED TO ENVIRONMENT

DATA VALUE: REPIRATORY DEPRESSION

-----  
ABBREVIATION.....Resp Depr  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Respiratory Depression  
HELP\_TEXT.....  
ID.....{7FE6D54C-1D7E-4028-B84C-E8EC416F9104}  
TERM.....REPIRATORY DEPRESSION

DATA VALUE: REPORTED

-----  
ABBREVIATION.....REPORTED  
ACTIVE.....Yes  
DESCRIPTION.....Reported as having happened  
DISPLAY\_NAME.....Reported  
HELP\_TEXT.....No Help Text  
ID.....{1D13BAC5-D8F9-455C-B815-8E256943492E}  
TERM.....REPORTED

DATA VALUE: REPOSITION

-----  
ABBREVIATION.....REPOSITION  
ACTIVE.....Yes  
DESCRIPTION.....Reposition  
DISPLAY\_NAME.....Reposition  
HELP\_TEXT.....No Help Text  
ID.....{DF3039CA-5977-4454-81D4-DB0AF966F85D}  
TERM.....REPOSITION

DATA VALUE: REQUEST FOR MORE MEDICATION

-----  
ABBREVIATION.....REQ MOR ME  
ACTIVE.....Yes  
DESCRIPTION.....Request For More Medication  
DISPLAY\_NAME.....Request For More Medication  
HELP\_TEXT.....No Help Text  
ID.....{0639081C-BEEC-4CAD-B420-5F2C4DBA38D0}  
TERM.....REQUEST FOR MORE MEDICATION

DATA VALUE: RES CRIT: AGGRESSIVE

-----  
ABBREVIATION.....R.C. AGGR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Aggressive  
HELP\_TEXT.....No Help Text  
ID.....{7A47B361-85EA-4C88-A1AF-CD8FF7D35D97}  
TERM.....RES CRIT: AGGRESSIVE

DATA VALUE: RES CRIT: AMBULATE

-----  
ABBREVIATION.....R.C. AMBU  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Attempting To Ambulate When  
                  Unable/High Injury Risk  
HELP\_TEXT.....No Help Text  
ID.....{54B72608-2C13-46CD-8034-00179669762C}  
TERM.....RES CRIT: AMBULATE

DATA VALUE: RES CRIT: COMBATIVE

-----  
ABBREVIATION.....R.C. COMB  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Combative  
HELP\_TEXT.....No Help Text  
ID.....{38DBF61A-D364-4215-B4AA-36B660F96462}  
TERM.....RES CRIT: COMBATIVE

DATA VALUE: RES CRIT: EQUIPMENT

-----  
ABBREVIATION.....R.C. EQUI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Attempting To Remove  
                  Equipment/Dressings  
HELP\_TEXT.....No Help Text  
ID.....{38C2C56C-6D57-4750-8E68-F1193B16190A}  
TERM.....RES CRIT: EQUIPMENT

DATA VALUE: RES CRIT: HARM TO SELF

-----  
ABBREVIATION.....R.C. HARM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Risk of Physical Harm To Self  
HELP\_TEXT.....No Help Text  
ID.....{4DC33341-8F53-4D30-BC28-1E32DB54F30D}  
TERM.....RES CRIT: HARM TO SELF

DATA VALUE: RES CRIT: LINES

-----  
ABBREVIATION.....R.C. LINE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Attempting To Pull Lines  
HELP\_TEXT.....No Help Text  
ID.....{27F545F2-5338-4BC8-B5CE-1D321D28F8D6}  
TERM.....RES CRIT: LINES

DATA VALUE: RES CRIT: UNABLE TO FOLLOW DIRECTIONS

-----  
ABBREVIATION.....R.C. DIRE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Unable To Follow Directions While  
Injury Risk  
HELP\_TEXT.....No Help Text  
ID.....{41358EAE-9912-48F8-8931-D3C635406D21}  
TERM.....RES CRIT: UNABLE TO FOLLOW DIRECTIONS

DATA VALUE: RES EVAL: AMBULATION/EXERCISE

-----  
ABBREVIATION.....R.E. EXER  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Ambulation / Exercise  
HELP\_TEXT.....No Help Text  
ID.....{46F4D1EB-3439-4BE4-A78C-453BDBF6C230}  
TERM.....RES EVAL: AMBULATION/EXERCISE

DATA VALUE: RES EVAL: DECREASE STIMULATION

-----  
ABBREVIATION.....R.E. DECR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Decrease Stimulation  
HELP\_TEXT.....No Description  
ID.....{0841626E-994E-4346-B701-1662A00209EB}  
TERM.....RES EVAL: DECREASE STIMULATION

DATA VALUE: RES EVAL: DIVERSIONAL ACTIVITIES

-----  
ABBREVIATION.....R.E. DIVER  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Diversional Activities  
HELP\_TEXT.....No Help Text  
ID.....{569F6D49-C564-4910-835A-B05179086E1F}  
TERM.....RES EVAL: DIVERSIONAL ACTIVITIES

DATA VALUE: RES EVAL: EDUCATION

-----  
ABBREVIATION.....R.E. EDUC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Education  
HELP\_TEXT.....No Help Text  
ID.....{76341E6E-7F0A-429C-B768-9CA83C452FD7}  
TERM.....RES EVAL: EDUCATION

DATA VALUE: RES EVAL: INVOLVEMENT WITH FAMILY

-----  
ABBREVIATION.....R.E. INVO  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Involvement With Family And Others  
HELP\_TEXT.....No Help Text  
ID.....{0EF31047-70A5-4072-953A-B62E1ACDF84B}  
TERM.....RES EVAL: INVOLVEMENT WITH FAMILY

DATA VALUE: RES EVAL: MEDICATION REVIEW

-----  
ABBREVIATION.....R.E. REVI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Medication Review  
HELP\_TEXT.....No Help Text  
ID.....{914A0F82-5863-4C84-B3AB-D5A22BD340D8}  
TERM.....RES EVAL: MEDICATION REVIEW

DATA VALUE: RES EVAL: REALITY ORIENTATION

-----  
ABBREVIATION.....R.E. REALI  
ACTIVE.....Yes  
DESCRIPTION.....None  
DISPLAY\_NAME.....Reality Orientation  
HELP\_TEXT.....None  
ID.....{221C363C-2261-40FC-8400-8B1D2E693F69}  
TERM.....RES EVAL: REALITY ORIENTATION

DATA VALUE: RES EVAL: RELOCATE

-----  
ABBREVIATION.....R.E. RELO  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Relocate Patient  
HELP\_TEXT.....No Help Text  
ID.....{D8165678-ECEF-4F4D-A1C4-C616349E1873}  
TERM.....RES EVAL: RELOCATE

DATA VALUE: RES EVAL: UNABLE TO ATTEMPT

-----  
ABBREVIATION.....R.E. UNAB  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Unable To Attempt Due To Patient  
Safety Concern  
HELP\_TEXT.....No Help Text  
ID.....{979B0477-5789-4319-9562-6752048C4005}  
TERM.....RES EVAL: UNABLE TO ATTEMPT

DATA VALUE: RES LRM: DEVICE ACCESS

-----  
ABBREVIATION.....RLRM ACCE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Prevent Access To Medical Devices  
HELP\_TEXT.....No Help Text  
ID.....{1662F5E1-30E7-4E8A-BCC5-1738868435B8}  
TERM.....RES LRM: DEVICE ACCESS

DATA VALUE: RES LRM: LIMIT SETTING

-----  
ABBREVIATION.....RLRM SETT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Limit Setting  
HELP\_TEXT.....No Help Text  
ID.....{96417DDF-8562-4EFF-B4F6-0933FD91CD36}  
TERM.....RES LRM: LIMIT SETTING

DATA VALUE: RES LRM: PRN MEDS

-----  
ABBREVIATION.....RLRM: MEDS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Offer PRN Meds  
HELP\_TEXT.....No Help Text  
ID.....{CDF02FA3-FC72-4CCC-8197-41710A69A2CB}  
TERM.....RES LRM: PRN MEDS

DATA VALUE: RES LRM: REDIRECT

-----  
ABBREVIATION.....RLRM REDI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Redirect Patient Activity  
HELP\_TEXT.....No Help Text  
ID.....{07DE67AF-E226-4794-ADF5-17E9C56C9FDF}  
TERM.....RES LRM: REDIRECT

DATA VALUE: RES LRM: REMOVE FROM SETTING

-----  
ABBREVIATION.....RLRM REMO  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Remove From Stimulating Setting  
HELP\_TEXT.....No Help Text  
ID.....{857DBF6E-7EFF-4F07-9AC1-7FDB37CC7618}  
TERM.....RES LRM: REMOVE FROM SETTING

DATA VALUE: RES LRM: REORIENT

-----  
ABBREVIATION.....RLRM REOR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Reoriented To Environment  
HELP\_TEXT.....No Help Text  
ID.....{DFCC8245-180D-4F3F-AEFC-5D55AFDDC46D}  
TERM.....RES LRM: REORIENT

DATA VALUE: RES LRM: SUPERVISION

-----  
ABBREVIATION.....RLRM SUPR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....One To One Supervision  
HELP\_TEXT.....No Help Text  
ID.....{18FA8CB0-0155-4D1C-9B14-AE51A0C06280}  
TERM.....RES LRM: SUPERVISION

DATA VALUE: RES LRM: VERBAL DE-ESCALATION

-----  
ABBREVIATION.....RLRM VERB  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Verbal De-escalation  
HELP\_TEXT.....No Help Text  
ID.....{CE8618C8-CAD7-4C74-8566-ACB5E66B082C}  
TERM.....RES LRM: VERBAL DE-ESCALATION

DATA VALUE: RES LRM: VERBAL REMINDERS

-----  
ABBREVIATION.....RLRM REMI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Frequent Verbal Reminders  
HELP\_TEXT.....No Help Text  
ID.....{987B5A32-DE2E-47E5-9B06-1AEE86096AB2}  
TERM.....RES LRM: VERBAL REMINDERS

DATA VALUE: RESOURCE BREEZE

-----  
ABBREVIATION.....Re Breeze  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Resource Breeze  
HELP\_TEXT.....  
ID.....{D3A38687-AD62-4E1D-A945-1201168E1CF6}  
TERM.....RESOURCE BREEZE

DATA VALUE: RESOURCE DM

-----  
ABBREVIATION.....Re Dm  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Resource Dm  
HELP\_TEXT.....  
ID.....{C9D4A461-81F9-4392-98DA-129B00E5807F}  
TERM.....RESOURCE DM

DATA VALUE: RESP TECH

-----  
ABBREVIATION.....Resp Tech  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Resp Tech  
HELP\_TEXT.....  
ID.....{8A7BA130-3EDC-4D07-A852-846B37399FC9}  
TERM.....RESP TECH

DATA VALUE: RESPIRATION

-----  
ABBREVIATION.....Respiratio  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Respiration  
HELP\_TEXT.....  
ID.....{5D6198F0-FE26-4F94-9118-D49A7EC44017}  
TERM.....RESPIRATION



DATA VALUE: RESPIRATORY DEPRESSION

-----  
ABBREVIATION.....Resp Depr  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Respiratory Depression  
HELP\_TEXT.....  
ID.....{A4801B08-0A98-4FDB-B075-C8D84DEAD5A0}  
TERM.....RESPIRATORY DEPRESSION

DATA VALUE: RESTING

-----  
ABBREVIATION.....Resting  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Resting  
HELP\_TEXT.....  
ID.....{D750373F-12E2-4D95-AF0D-F7731905F723}  
TERM.....RESTING

DATA VALUE: RESTING RESPIRATORY RATE BASELINE OR ABOVE

-----  
ABBREVIATION.....Rest Base  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Resting Respiratory Rate Baseline or  
Above  
HELP\_TEXT.....  
ID.....{C3447D0F-AB90-4FA4-A2A2-7681B065AF07}  
TERM.....RESTING RESPIRATORY RATE BASELINE OR  
ABOVE

DATA VALUE: RESTLESS

-----  
ABBREVIATION.....r  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Restless  
HELP\_TEXT.....No Help Text  
ID.....{8FCCAD26-DBA0-7067-DA72-30D342BFA8E8}  
TERM.....RESTLESS

DATA VALUE: RESTRAINT

-----  
ABBREVIATION.....RESTRAINT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Restraint  
HELP\_TEXT.....No Help Text  
ID.....{0429322D-A987-4913-B64D-A6414645B214}  
TERM.....RESTRAINT

DATA VALUE: RESTRAINT BILATERAL ANKLE

-----  
ABBREVIATION.....RES BI ANK  
ACTIVE.....Yes  
DESCRIPTION.....Bilateral Ankle Restraint  
DISPLAY\_NAME.....Bilateral Ankle  
HELP\_TEXT.....No Help Text  
ID.....{8A0B0DE5-E652-40E7-9B67-2F4D27685A68}  
TERM.....RESTRAINT BILATERAL ANKLE

DATA VALUE: RESTRAINT BILATERAL HAND MITT

-----  
ABBREVIATION.....RES BI H M  
ACTIVE.....Yes  
DESCRIPTION.....Bilateral Hand Mitt Restraint  
DISPLAY\_NAME.....Bilateral Hand Mitt  
HELP\_TEXT.....No Help Text  
ID.....{5E293110-1235-4307-80C2-0FFFA007E1DA}  
TERM.....RESTRAINT BILATERAL HAND MITT

DATA VALUE: RESTRAINT GERI CHAIR

-----  
ABBREVIATION.....RES GERI C  
ACTIVE.....Yes  
DESCRIPTION.....Geri Chair Restraint  
DISPLAY\_NAME.....Geri Chair  
HELP\_TEXT.....No Help Text  
ID.....{6B8F87BB-C2AB-4D92-AE04-24F19B26D2F1}  
TERM.....RESTRAINT GERI CHAIR

DATA VALUE: RESTRAINT HARD PLASTIC/LEATHER

-----  
ABBREVIATION.....RES H PL/L  
ACTIVE.....Yes  
DESCRIPTION.....Hard Plastic/Leather Restraint  
DISPLAY\_NAME.....Hard Plastic/Leather  
HELP\_TEXT.....No Help Text  
ID.....{9BB528F4-7016-4030-B73F-4037ADF2621D}  
TERM.....RESTRAINT HARD PLASTIC/LEATHER

DATA VALUE: RESTRAINT HARD PLASTIC/LEATHER 2 POINT

-----  
ABBREVIATION.....RES P/L 2P  
ACTIVE.....Yes  
DESCRIPTION.....Hard Plastic/Leather 2 Point Restraint  
DISPLAY\_NAME.....Hard Plastic/Leather 2 Point  
HELP\_TEXT.....No Help Text  
ID.....{17CB8739-612A-4319-946A-2062C5245C90}  
TERM.....RESTRAINT HARD PLASTIC/LEATHER 2 POINT

DATA VALUE: RESTRAINT HARD PLASTIC/LEATHER 4 POINT

-----  
ABBREVIATION.....RES P/L 4P  
ACTIVE.....Yes  
DESCRIPTION.....Hard Plastic/Leather 4 Point Restraint  
DISPLAY\_NAME.....Hard Plastic/Leather 4 Point  
HELP\_TEXT.....No Help Text  
ID.....{94DFD0BA-E941-4B09-8FD7-87E9AD7C7E6F}  
TERM.....RESTRAINT HARD PLASTIC/LEATHER 4 POINT

DATA VALUE: RESTRAINT ONE ARM

-----  
ABBREVIATION.....RES ONE AR  
ACTIVE.....Yes  
DESCRIPTION.....One Arm Restraint  
DISPLAY\_NAME.....One Arm  
HELP\_TEXT.....No Help Text  
ID.....{839914DA-6D63-43EB-961E-70DA56737A92}  
TERM.....RESTRAINT ONE ARM

DATA VALUE: RESTRAINT ONE LEG

-----  
ABBREVIATION.....RES ONE LE  
ACTIVE.....Yes  
DESCRIPTION.....One Leg Restraint  
DISPLAY\_NAME.....One Leg  
HELP\_TEXT.....No Help Text  
ID.....{2790F0FC-B49A-4C99-9968-BCADD7CC63AC}  
TERM.....RESTRAINT ONE LEG

DATA VALUE: RESTRAINT SINGLE HAND MITT

-----  
ABBREVIATION.....RES S HA M  
ACTIVE.....Yes  
DESCRIPTION.....Single Hand Mitt Restraint  
DISPLAY\_NAME.....Single Hand Mitt  
HELP\_TEXT.....No Help Text  
ID.....{A5F1D037-FB5E-402B-9F25-22DE48404D57}  
TERM.....RESTRAINT SINGLE HAND MITT

DATA VALUE: RESTRAINT SOFT

-----  
ABBREVIATION.....RESTR SOFT  
ACTIVE.....Yes  
DESCRIPTION.....Soft Restraint  
DISPLAY\_NAME.....Soft Restraint  
HELP\_TEXT.....No Help Text  
ID.....{C02CAC4B-5E9C-4060-BECD-0C0820EDA39B}  
TERM.....RESTRAINT SOFT

DATA VALUE: RESTRAINT SOFT 2 POINT

-----  
ABBREVIATION.....RES S2POIN  
ACTIVE.....Yes  
DESCRIPTION.....Soft 2 Point Restraint  
DISPLAY\_NAME.....Soft Restraint 2 Point  
HELP\_TEXT.....No Help Text  
ID.....{66B71CFE-2603-40C9-9625-24448012BAB5}  
TERM.....RESTRAINT SOFT 2 POINT

DATA VALUE: RESTRAINT SOFT 4 POINT

-----  
ABBREVIATION.....RES S4POIN  
ACTIVE.....Yes  
DESCRIPTION.....Soft 4 Point Restraint  
DISPLAY\_NAME.....Soft Restraint 4 Point  
HELP\_TEXT.....No Help Text  
ID.....{6820E997-60A0-4802-BE26-CA29B4F48F32}  
TERM.....RESTRAINT SOFT 4 POINT

DATA VALUE: RESTRAINT SOFT ARMS + LOCKED WAIST

-----  
ABBREVIATION.....RES S A+LW  
ACTIVE.....Yes  
DESCRIPTION.....Soft Arms + Locked Waist Restraint  
DISPLAY\_NAME.....Soft Arms + Locked Waist  
HELP\_TEXT.....No Help Text  
ID.....{C480D67C-958C-42F9-9203-98E2F2BDBEC8}  
TERM.....RESTRAINT SOFT ARMS + LOCKED WAIST

DATA VALUE: RESTRAINT SOFT ARMS + WAIST

-----  
ABBREVIATION.....RES S A+W  
ACTIVE.....Yes  
DESCRIPTION.....Soft Arms + Waist Restraints  
DISPLAY\_NAME.....Soft Arms + Waist  
HELP\_TEXT.....No Help Text  
ID.....{4A322BB6-0F63-4BB7-9290-C0655BD7D772}  
TERM.....RESTRAINT SOFT ARMS + WAIST

DATA VALUE: RESTRAINT STRAP

-----  
ABBREVIATION.....RES STRAP  
ACTIVE.....Yes  
DESCRIPTION.....Restraint Strap  
DISPLAY\_NAME.....Restraint Strap  
HELP\_TEXT.....No Help Text  
ID.....{642FDACE-B320-42BC-9514-F5FCD6583B19}  
TERM.....RESTRAINT STRAP

DATA VALUE: RESTRAINT VEST

-----  
ABBREVIATION.....RES VEST  
ACTIVE.....Yes  
DESCRIPTION.....Restraint Vest  
DISPLAY\_NAME.....Restraint Vest  
HELP\_TEXT.....No Help Text  
ID.....{643378B2-86D7-4DB7-93ED-9DCE16D85E27}  
TERM.....RESTRAINT VEST

DATA VALUE: RESTRAINT WAIST

-----  
ABBREVIATION.....RES WAIST  
ACTIVE.....Yes  
DESCRIPTION.....Waist  
DISPLAY\_NAME.....Waist  
HELP\_TEXT.....No Help Text  
ID.....{253E2757-5A57-4410-9C0E-F4A907A1EC16}  
TERM.....RESTRAINT WAIST

DATA VALUE: RESTRAINT WAIST LOCKED

-----  
ABBREVIATION.....RES WAI LO  
ACTIVE.....Yes  
DESCRIPTION.....Waist Locked Restraint  
DISPLAY\_NAME.....Waist Locked  
HELP\_TEXT.....No Help Text  
ID.....{1EE7AB42-CE17-4D49-849B-CA550966E58C}  
TERM.....RESTRAINT WAIST LOCKED

DATA VALUE: RESTRAINT WRIST

-----  
ABBREVIATION.....RES WRIST  
ACTIVE.....Yes  
DESCRIPTION.....Wrist Restraint  
DISPLAY\_NAME.....Restraint Wrist  
HELP\_TEXT.....No Help Text  
ID.....{03038E52-A66F-4EC0-8DC1-7B0BB596FFC0}  
TERM.....RESTRAINT WRIST

DATA VALUE: RESTRAINTS BROCHURE

-----  
ABBREVIATION.....RES BRO P  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Education Brochure Provided  
HELP\_TEXT.....No Help Text  
ID.....{1E54C2A1-8B8D-4F0D-B247-91E19722CDED}  
TERM.....RESTRAINTS BROCHURE

DATA VALUE: RESTRAINTS ED REMOVAL

-----  
ABBREVIATION.....RES ED REM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Education On Condition Required For  
Removal  
HELP\_TEXT.....No Help Text  
ID.....{B19316C8-687B-48FB-9EBC-9857662B9F85}  
TERM.....RESTRAINTS ED REMOVAL

DATA VALUE: RESTRAINTS ED TYPE

-----  
ABBREVIATION.....RES ED TY  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Education Rational For Type Of  
Restraint  
HELP\_TEXT.....No Help Text  
ID.....{5560A332-C365-4DAF-BA33-610C4AE301F5}  
TERM.....RESTRAINTS ED TYPE

DATA VALUE: RESTRAINTS EXPLAINED

-----  
ABBREVIATION.....RES EXPL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Explain Why Restraints Applied  
HELP\_TEXT.....No Help Text  
ID.....{1C063DD2-C4A3-4027-A245-6DD51DF5A3E9}  
TERM.....RESTRAINTS EXPLAINED

DATA VALUE: RESTRAINTS LIMIT SETTING

-----  
ABBREVIATION.....RES LIM SE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Limit Setting  
HELP\_TEXT.....No Help Text  
ID.....{F97BA57E-ED33-D6FE-5E35-E6EFB0625773}  
TERM.....RESTRAINTS LIMIT SETTING

DATA VALUE: RESTRAINTS TIME LIMIT

-----  
ABBREVIATION.....RES TIME L  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Time Limit On Restraints  
HELP\_TEXT.....No Help Text  
ID.....{513B483E-AD7B-4B1A-BD01-91133C605C3C}  
TERM.....RESTRAINTS TIME LIMIT

DATA VALUE: RHONCHI

-----  
ABBREVIATION.....RHONCHI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Rhonchi  
HELP\_TEXT.....No Help Text  
ID.....{850A9D2A-2F7F-A812-3D56-EE0948B0A0CC}  
TERM.....RHONCHI

DATA VALUE: RIGHT

-----  
ABBREVIATION.....RIGHT  
ACTIVE.....Yes  
DESCRIPTION.....Right  
DISPLAY\_NAME.....Right  
HELP\_TEXT.....No Help Text  
ID.....{AB74C579-DF14-4D7B-A9EF-5F119089F976}  
TERM.....RIGHT

DATA VALUE: RIGID HEEL SPLINT

-----  
ABBREVIATION.....RIG HEEL S  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Rigid Heel Splint  
HELP\_TEXT.....No Help Text  
ID.....{5A231387-9FAC-4B2D-94D6-5DF4312E9AF8}  
TERM.....RIGID HEEL SPLINT

DATA VALUE: RN

-----  
ABBREVIATION.....RN  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....RN  
HELP\_TEXT.....  
ID.....{A103E506-EAA5-4777-8126-268CD723AD87}  
TERM.....RN

DATA VALUE: ROPIVICAINE

-----  
ABBREVIATION.....Ropivicain  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Ropivicaïne  
HELP\_TEXT.....  
ID.....{911BE8BC-6737-4DBF-96AF-A501D3D54493}  
TERM.....ROPIVICAINE

DATA VALUE: ROTATIONAL THERAPY

-----  
ABBREVIATION.....ROT THERAP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Rotational Therapy  
HELP\_TEXT.....No Help Text  
ID.....{BB795594-720C-40D7-9B7C-6792BC849DA0}  
TERM.....ROTATIONAL THERAPY

DATA VALUE: ROUND

-----  
ABBREVIATION.....ROUND  
ACTIVE.....Yes  
DESCRIPTION.....Round  
DISPLAY\_NAME.....Round  
HELP\_TEXT.....No Help Text  
ID.....{E8E741F1-F0C5-4B77-B014-63B545BED755}  
TERM.....ROUND

DATA VALUE: RT

-----  
ABBREVIATION.....RT  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....RT  
HELP\_TEXT.....  
ID.....{6E374EFA-B7CE-44C2-ADCE-DF9C132E951E}  
TERM.....RT

DATA VALUE: RUB

-----  
ABBREVIATION.....r  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Rub  
HELP\_TEXT.....No Help Text  
ID.....{0040835D-1CF0-7D55-EEBB-A94A1E0D1206}  
TERM.....RUB

DATA VALUE: RUNS OF VT

-----  
ABBREVIATION.....Runs of VT  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Runs of VT  
HELP\_TEXT.....  
ID.....{3BA2FF31-1A56-4B4B-AC13-E8608BC3FCC2}  
TERM.....RUNS OF VT



DATA VALUE: S1

-----  
ABBREVIATION.....s  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....S1  
HELP\_TEXT.....No Help Text  
ID.....{1AB9DF0A-8238-E8FB-8520-4749E13DDE1D}  
TERM.....S1

DATA VALUE: S2

-----  
ABBREVIATION.....s  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....S2  
HELP\_TEXT.....No Help Text  
ID.....{57C043D3-B227-3E7A-C1E6-5129A2800718}  
TERM.....S2

DATA VALUE: S3

-----  
ABBREVIATION.....s  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....S3  
HELP\_TEXT.....No Help Text  
ID.....{9374B57F-DE1F-B9C4-26B5-2DD6B194DC26}  
TERM.....S3

DATA VALUE: S4

-----  
ABBREVIATION.....s  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....S4  
HELP\_TEXT.....No Help Text  
ID.....{9C6CA302-85B6-8FFC-EAB4-8FED62C535A5}  
TERM.....S4

DATA VALUE: SAFETY

-----  
ABBREVIATION.....SAFETY  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Safety  
HELP\_TEXT.....No Help Text  
ID.....{AD905A8F-33BD-4EFE-929E-A2B575797402}  
TERM.....SAFETY

DATA VALUE: SALINE LOCK

-----  
ABBREVIATION.....SALIN LOCK  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Saline Lock  
HELP\_TEXT.....No Help Text  
ID.....{A94D7639-2A45-4319-A6C7-5A892AA6B1F4}  
TERM.....SALINE LOCK

DATA VALUE: SALINE MOIST DRESSING

-----  
ABBREVIATION.....SAL MOI DR  
ACTIVE.....Yes  
DESCRIPTION.....Saline Moist Dressing  
DISPLAY\_NAME.....Saline Moist Dressing  
HELP\_TEXT.....No Help Text  
ID.....{19F98514-8ABF-461C-9969-C31177FE7DE7}  
TERM.....SALINE MOIST DRESSING

DATA VALUE: SANDOSOURCE PEPTIDE

-----  
ABBREVIATION.....SAND PEPTI  
ACTIVE.....Yes  
DESCRIPTION.....Sandosource Peptide Enteral Nutrition  
Formula  
DISPLAY\_NAME.....Sandosource Peptide Enteral Nutrition  
Formula  
HELP\_TEXT.....No Help Text  
ID.....{0EF43998-4ED4-47F5-ACC6-5A977B6F9997}  
TERM.....SANDOSOURCE PEPTIDE

DATA VALUE: SAS AGITATED

-----  
ABBREVIATION.....SAS AGITAT  
ACTIVE.....Yes  
DESCRIPTION.....Anxious or physically agitated, calms  
to verbal instructions  
DISPLAY\_NAME.....5 - Agitated  
HELP\_TEXT.....Anxious or physically agitated, calms  
to verbal instructions  
ID.....{107143D4-9527-4A26-8CB7-EEBB656B68B3}  
TERM.....SAS AGITATED

DATA VALUE: SAS CALM AND COOPERATIVE

-----  
ABBREVIATION.....SAS C COOP  
ACTIVE.....Yes  
DESCRIPTION.....Calm, easily arousable, follows  
commands  
DISPLAY\_NAME.....4 - Calm and Cooperative  
HELP\_TEXT.....Calm, easily arousable, follows  
commands  
ID.....{F0633778-280C-43C3-9E6E-9C35C18179F6}  
TERM.....SAS CALM AND COOPERATIVE

DATA VALUE: SAS DANGEROUS AGITATION

-----  
ABBREVIATION.....SAS DA AGI  
ACTIVE.....Yes  
DESCRIPTION.....Pulling at ET tube, trying to remove  
                  catheters, climbing over bedrail,  
                  striking at staff, thrashing  
                  side-to-side  
DISPLAY\_NAME.....7 - Dangerous Agitation  
HELP\_TEXT.....Pulling at ET tube, trying to remove  
                  catheters, climbing over bedrail,  
                  striking at staff, thrashing  
                  side-to-side  
ID.....{0701FB25-1DF0-4771-8C81-AC19A985A3B6}  
TERM.....SAS DANGEROUS AGITATION

DATA VALUE: SAS SEDATED

-----  
ABBREVIATION.....SAS SEDATE  
ACTIVE.....Yes  
DESCRIPTION.....Difficult to arouse but awakens to  
                  verbal stimuli or gentle shaking,  
                  follows simple commands but drifts  
                  off again  
DISPLAY\_NAME.....3 - Sedated  
HELP\_TEXT.....Difficult to arouse but awakens to  
                  verbal stimuli or gentle shaking,  
                  follows simple commands but drifts  
                  off again  
ID.....{AFEC3F29-52D3-4A49-91F8-05AA3891FEC8}  
TERM.....SAS SEDATED

DATA VALUE: SAS UNAROUSABLE

-----  
ABBREVIATION.....SAS UNAROU  
ACTIVE.....Yes  
DESCRIPTION.....Minimal or no response to noxious  
                  stimuli, does not communicate or  
                  follow commands  
DISPLAY\_NAME.....1 - Unarousable  
HELP\_TEXT.....Minimal or no response to noxious  
                  stimuli, does not communicate or  
                  follow commands  
ID.....{3DDFF9FE-A88A-471F-BF36-E5B569DEBDAA}  
TERM.....SAS UNAROUSABLE

DATA VALUE: SAS VERY AGITATED

-----  
ABBREVIATION.....SAS VY AGI  
ACTIVE.....Yes  
DESCRIPTION.....Requiring restraint and frequent  
                  verbal reminding of limits, biting  
                  ETT  
DISPLAY\_NAME.....6 - Very Agitated  
HELP\_TEXT.....Requiring restraint and frequent  
                  verbal reminding of limits, biting  
                  ETT  
ID.....{B501606A-EEA7-470E-81BE-F51D422272F1}  
TERM.....SAS VERY AGITATED

DATA VALUE: SAS VERY SEDATED

-----  
ABBREVIATION.....SAS VY SED

ACTIVE.....Yes

DESCRIPTION.....Arouses to physical stimuli but does  
not communicate or follow commands,  
may move spontaneously

DISPLAY\_NAME.....2 - Very Sedated

HELP\_TEXT.....Arouses to physical stimuli but does  
not communicate or follow commands,  
may move spontaneously

ID.....{4B13529B-61C5-4B68-857A-3A5351DC4BE2}

TERM.....SAS VERY SEDATED

DATA VALUE: SCANT

-----  
ABBREVIATION.....SCANT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Scant  
HELP\_TEXT.....No Help Text  
ID.....{8EF6C62D-0D94-37BB-60AD-86E306E8120D}  
TERM.....SCANT

DATA VALUE: SCAPHOID

-----  
ABBREVIATION.....SCAPHOID  
ACTIVE.....Yes  
DESCRIPTION.....Scaphoid  
DISPLAY\_NAME.....Scaphoid  
HELP\_TEXT.....No Help Text  
ID.....{1BE93489-42FC-44EB-8D07-87F4224979F3}  
TERM.....SCAPHOID

DATA VALUE: SCOOTER

-----  
ABBREVIATION.....Scooter  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Scooter  
HELP\_TEXT.....  
ID.....{663BC475-C20E-4B8B-B8B5-B16CC79A93C6}  
TERM.....SCOOTER

DATA VALUE: SECLUSION

-----  
ABBREVIATION.....SECLUSION  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Seclusion  
HELP\_TEXT.....No Help Text  
ID.....{65B57F6C-C849-4BF7-8705-175361BC900F}  
TERM.....SECLUSION

DATA VALUE: SECOND DEGREE AV BLOCK: TYPE 1

-----  
ABBREVIATION.....sdabt1  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Second Degree Av Block: Type 1  
HELP\_TEXT.....No Help Text  
ID.....{508BF663-96FB-EEBE-4F06-D2F1A02BD099}  
TERM.....SECOND DEGREE AV BLOCK: TYPE 1

DATA VALUE: SECOND DEGREE AV BLOCK: TYPE 2

-----  
ABBREVIATION.....sdabt2  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Second Degree Av Block: Type 2  
HELP\_TEXT.....No Help Text  
ID.....{3367754A-08D1-31A9-E850-1B36358E358A}  
TERM.....SECOND DEGREE AV BLOCK: TYPE 2

DATA VALUE: SECURE TRACK IN PLACE

-----  
ABBREVIATION.....Secure Trk  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Secure Track in Place  
HELP\_TEXT.....  
ID.....{33F72416-BAD4-44AD-825E-93BCBFFB9670}  
TERM.....SECURE TRACK IN PLACE

DATA VALUE: SEDATED

-----  
ABBREVIATION.....s  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Sedated  
HELP\_TEXT.....No Help Text  
ID.....{29D9ADD6-2C1F-F6D5-F4D0-B2BB86416997}  
TERM.....SEDATED

DATA VALUE: SEDIMENT

-----  
ABBREVIATION.....SEDIMENT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Sediment  
HELP\_TEXT.....No Help Text  
ID.....{4D4E4236-55CE-7ECD-FC55-6AD0CFD4D9A5}  
TERM.....SEDIMENT

DATA VALUE: SEIZURES

-----  
ABBREVIATION.....Seizures  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Seizures  
HELP\_TEXT.....  
ID.....{DE974E49-2AB6-4781-9930-468F2B781BA9}  
TERM.....SEIZURES

DATA VALUE: SELECTED SURGERY

-----  
ABBREVIATION.....Sel Surg  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Selected Surgery  
HELP\_TEXT.....  
ID.....{6DBC6356-0F00-4001-AF94-20A479995A6C}  
TERM.....SELECTED SURGERY

DATA VALUE: SELF AM CARE BY PATIENT

-----  
ABBREVIATION.....SELF AM CR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Self AM Care By Patient  
HELP\_TEXT.....No Help Text  
ID.....{C48A46C1-0E68-479C-B803-08C80617ADFF}  
TERM.....SELF AM CARE BY PATIENT

DATA VALUE: SELF CARE

-----  
ABBREVIATION.....SELF CARE  
ACTIVE.....Yes  
DESCRIPTION.....Self Care  
DISPLAY\_NAME.....Self Care  
HELP\_TEXT.....No Help Text  
ID.....{0E068188-75B8-46C5-8947-A6CC5C7849B5}  
TERM.....SELF CARE

DATA VALUE: SELF CARE OF BEDSIDE ENVIRONMENT

-----  
ABBREVIATION.....S CARE BED  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Self Care of Bedside Environment  
HELP\_TEXT.....No Help Text  
ID.....{A1BEC7F4-D404-4824-8BB7-108E0A6C43A2}  
TERM.....SELF CARE OF BEDSIDE ENVIRONMENT

DATA VALUE: SELF FACE WASH

-----  
ABBREVIATION.....S FACE WAS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Self Face Wash  
HELP\_TEXT.....No Help Text  
ID.....{77457D46-5676-46E6-B38F-FE3AF70B5CBC}  
TERM.....SELF FACE WASH

DATA VALUE: SELF FOOT CARE

-----  
ABBREVIATION.....SELF FOOT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Self Foot Care  
HELP\_TEXT.....No Help Text  
ID.....{F2DCAF0E-389D-4987-8C73-66E35266CC78}  
TERM.....SELF FOOT CARE

DATA VALUE: SELF SHAVE

-----  
ABBREVIATION.....SELF SHAVE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Self Shave  
HELP\_TEXT.....No Help Text  
ID.....{90F1760C-4C12-4976-ADCB-DBF456A34023}  
TERM.....SELF SHAVE

DATA VALUE: SELF SHOWER

-----  
ABBREVIATION.....SELF SHOWE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Self Shower  
HELP\_TEXT.....No Help Text  
ID.....{C35CCAE8-27D9-46F2-8971-2C4CF4924CCB}  
TERM.....SELF SHOWER

DATA VALUE: SELF TRACH CARE

-----  
ABBREVIATION.....S TRACH CA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Self Trach Care  
HELP\_TEXT.....No Help Text  
ID.....{70BE79EC-107F-4CD3-B964-7919B3DDE222}  
TERM.....SELF TRACH CARE

DATA VALUE: SELF-CARE SKILLS

-----  
ABBREVIATION.....Self Skill  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Self-Care Skills  
HELP\_TEXT.....  
ID.....{39145E84-85F2-486C-9FF8-27787EC8026A}  
TERM.....SELF-CARE SKILLS



DATA VALUE: SEMIFOWLERS HIGH

-----  
ABBREVIATION.....sh  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Semifowlers High  
HELP\_TEXT.....No Help Text  
ID.....{1BBE7175-2BD2-BFB5-1E1B-E02443433072}  
TERM.....SEMIFOWLERS HIGH

DATA VALUE: SEMIFOWLERS LOW

-----  
ABBREVIATION.....sl  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Semifowlers Low  
HELP\_TEXT.....No Help Text  
ID.....{A265AF02-53FC-B182-8AE3-A06E95718D7E}  
TERM.....SEMIFOWLERS LOW

DATA VALUE: SEQUENTIAL COMPRESSION DEVICE

-----  
ABBREVIATION.....SCD  
ACTIVE.....Yes  
DESCRIPTION.....SEQUENTIAL COMPRESSION DEVICE  
DISPLAY\_NAME.....Sequential Compression Device  
HELP\_TEXT.....SEQUENTIAL COMPRESSION DEVICE  
ID.....{29D657DB-6679-8A55-0E33-235B5DD349F2}  
TERM.....SEQUENTIAL COMPRESSION DEVICE

DATA VALUE: SEROSANGUINOUS

-----  
ABBREVIATION.....SEROSANGUI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Serosanguinous  
HELP\_TEXT.....No Help Text  
ID.....{19A8D332-899D-177A-2247-ABCDAAEE3BFA8}  
TERM.....SEROSANGUINOUS

DATA VALUE: SEROUS

-----  
ABBREVIATION.....SEROUS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Serous  
HELP\_TEXT.....No Help Text  
ID.....{9BDA347E-29E9-80C7-91FA-016192D4257E}  
TERM.....SEROUS

DATA VALUE: SET UP HELP ONLY

-----  
ABBREVIATION.....SETUP ONLY  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Set Up Help only  
HELP\_TEXT.....No Help Text  
ID.....{DACC77AB-5AE7-4F1B-A3E1-C347D3DD7219}  
TERM.....SET UP HELP ONLY

DATA VALUE: SEVERELY REDUCED

-----  
ABBREVIATION.....SEV REDUCE  
ACTIVE.....Yes  
DESCRIPTION.....Severely Reduced  
DISPLAY\_NAME.....Severely Reduced  
HELP\_TEXT.....No Help Text  
ID.....{A895BACE-85AC-45B2-A190-4A6816DEDC1D}  
TERM.....SEVERELY REDUCED

DATA VALUE: SHALLOW

-----  
ABBREVIATION.....SHALLOW  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Shallow  
HELP\_TEXT.....No Help Text  
ID.....{CD205617-BF98-4177-971A-056AFEC823A}  
TERM.....SHALLOW

DATA VALUE: SHAMPOO

-----  
ABBREVIATION.....SHAMPOO  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Shampoo  
HELP\_TEXT.....No Help Text  
ID.....{95CCB23C-E614-40CD-84F0-EE4EFBC95271}  
TERM.....SHAMPOO

DATA VALUE: SHARP

-----  
ABBREVIATION.....SHARP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Sharp  
HELP\_TEXT.....No Help Text  
ID.....{20608864-E5CF-4FE3-8B6D-AE2FEE086503}  
TERM.....SHARP

DATA VALUE: SHAVE

-----  
ABBREVIATION.....SHAVE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Shave  
HELP\_TEXT.....No Help Text  
ID.....{A5CB33DF-C33B-4DD3-B00D-78D78F7966C7}  
TERM.....SHAVE

DATA VALUE: SHOOTING

-----  
ABBREVIATION.....SHOOTING  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Shooting  
HELP\_TEXT.....No Help Text  
ID.....{73C48DB1-6FE6-46B7-B1F6-E93ACBE44C27}  
TERM.....SHOOTING

DATA VALUE: SHOWER

-----  
ABBREVIATION.....SHOWER  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Shower  
HELP\_TEXT.....No Help Text  
ID.....{C259929A-44DB-48FD-A2EC-BC6ECF9D0476}  
TERM.....SHOWER

DATA VALUE: SHUFFLING

-----  
ABBREVIATION.....Shuffling  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Shuffling  
HELP\_TEXT.....  
ID.....{7EA3A898-1D90-472C-BE90-2110330A1ADA}  
TERM.....SHUFFLING

DATA VALUE: SHUNTS

-----  
ABBREVIATION.....Shunts  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Shunts  
HELP\_TEXT.....  
ID.....{8D9CB03B-6D63-4205-86C3-869421BEBD04}  
TERM.....SHUNTS

DATA VALUE: SICU

-----  
ABBREVIATION.....SICU  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....SICU  
HELP\_TEXT.....No Help Text  
ID.....{EAEB2DE6-B057-4143-B44C-99DAB9B7B4AC}  
TERM.....SICU

DATA VALUE: SIGMOIDOSTOMY

-----  
ABBREVIATION.....SIGMOIDOST  
ACTIVE.....Yes  
DESCRIPTION.....Sigmoidostomy  
DISPLAY\_NAME.....Sigmoidostomy  
HELP\_TEXT.....No Help Text  
ID.....{0C84B317-59B3-433B-A9AD-6588C1361E0C}  
TERM.....SIGMOIDOSTOMY

DATA VALUE: SIGNIFICANT OTHER

-----  
ABBREVIATION.....SO  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Significant Other  
HELP\_TEXT.....No Help Text  
ID.....{082C42E8-9E8A-4A96-B3AD-8B0E0A62AF00}  
TERM.....SIGNIFICANT OTHER

DATA VALUE: SILVADENE

-----  
ABBREVIATION.....SILVADENE  
ACTIVE.....Yes  
DESCRIPTION.....Silvadene  
DISPLAY\_NAME.....Silvadene  
HELP\_TEXT.....No Help Text  
ID.....{0889134F-C419-4366-B12E-EAB775719410}  
TERM.....SILVADENE

DATA VALUE: SILVER

-----  
ABBREVIATION.....SILVER  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Silver  
HELP\_TEXT.....No Help Text  
ID.....{AF233570-2BB4-4AE6-8751-056B9717C475}  
TERM.....SILVER

DATA VALUE: SILVER SULFADIAZINE

-----  
ABBREVIATION.....SILVER SUL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Silver Sulfadiazine  
HELP\_TEXT.....No Help Text  
ID.....{E7A99A4E-C9B5-40D6-9820-05F87004F48B}  
TERM.....SILVER SULFADIAZINE

DATA VALUE: SINGLE

-----  
ABBREVIATION.....SINGLE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Single  
HELP\_TEXT.....No Help Text  
ID.....{8D305E0F-9035-E396-E768-40D7E6455BE4}  
TERM.....SINGLE

DATA VALUE: SINGLE USE LUBRICANT PACKET

-----  
ABBREVIATION.....Single Use  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Single Use Lubricant Packet  
HELP\_TEXT.....  
ID.....{7F748EDB-77D4-42D2-BE75-54F73598D4B0}  
TERM.....SINGLE USE LUBRICANT PACKET

DATA VALUE: SINUS ARREST

-----  
ABBREVIATION.....sa  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Sinus Arrest  
HELP\_TEXT.....No Help Text  
ID.....{3CD8F7AC-3FE0-F342-089C-787232D2A5C2}  
TERM.....SINUS ARREST

DATA VALUE: SINUS ARRHYTHMIA

-----  
ABBREVIATION.....sa  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Sinus Arrhythmia  
HELP\_TEXT.....No Help Text  
ID.....{07A24AC9-B486-A4B8-BADB-04D2729542CE}  
TERM.....SINUS ARRHYTHMIA

DATA VALUE: SINUS BLOCK

-----  
ABBREVIATION.....sb  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Sinus Block  
HELP\_TEXT.....No Help Text  
ID.....{94D6EFED-3D08-E724-BC40-9C398919A142}  
TERM.....SINUS BLOCK

DATA VALUE: SINUS BRADYCARDIA

-----  
ABBREVIATION.....sb  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Sinus Bradycardia  
HELP\_TEXT.....No Help Text  
ID.....{D6EF4940-63AE-22AF-1A80-D8E6F0F07839}  
TERM.....SINUS BRADYCARDIA

DATA VALUE: SINUS PAUSE

-----  
ABBREVIATION.....sp  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Sinus Pause  
HELP\_TEXT.....No Help Text  
ID.....{7F8F89A9-6B16-B159-55B3-DEFDF5B4B4F6}  
TERM.....SINUS PAUSE

DATA VALUE: SINUS RHYTHM

-----  
ABBREVIATION.....sr  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Sinus Rhythm  
HELP\_TEXT.....No Help Text  
ID.....{A34F80EA-3D6A-076A-DB2A-E4E5C5AACA11}  
TERM.....SINUS RHYTHM

DATA VALUE: SINUS TACHYCARDIA

-----  
ABBREVIATION.....SINUS T  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Sinus Tachycardia  
HELP\_TEXT.....No Help Text  
ID.....{0AC41AF3-7656-A35D-E397-21D35405812A}  
TERM.....SINUS TACHYCARDIA

DATA VALUE: SITE UNABLE TO BE MARKED

-----  
ABBREVIATION.....Unable/Mar  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Site Unable to be Marked  
HELP\_TEXT.....  
ID.....{E7BA9E0F-1F12-4F74-B48E-B6212D217B30}  
TERM.....SITE UNABLE TO BE MARKED

DATA VALUE: SKIN INCISION WOUND CLOSURE OPEN WOUND

-----  
ABBREVIATION.....INC CL OPE  
ACTIVE.....Yes  
DESCRIPTION.....Skin Incision Closure Open Wound  
DISPLAY\_NAME.....Open Wound  
HELP\_TEXT.....Skin Incision Closure Open Wound  
ID.....{AD018B6F-13E5-41C0-BDBC-306CA0B1B475}  
TERM.....SKIN INCISION WOUND CLOSURE OPEN WOUND

DATA VALUE: SKIN INCISION WOUND CLOSURE STAPLES

-----  
ABBREVIATION.....INC CL STA  
ACTIVE.....Yes  
DESCRIPTION.....Skin Incision Closure Staples  
DISPLAY\_NAME.....Staples  
HELP\_TEXT.....Skin Incision Closure Staples  
ID.....{4F48D42D-C206-423C-970B-02C03B2EE0B6}  
TERM.....SKIN INCISION WOUND CLOSURE STAPLES

DATA VALUE: SKIN INCISION WOUND CLOSURE STERI-STRIPS

-----  
ABBREVIATION.....INC CL STE  
ACTIVE.....Yes  
DESCRIPTION.....Skin Incision Closure Steri-Strips  
DISPLAY\_NAME.....Steri-Strips  
HELP\_TEXT.....Skin Incision Closure Steri-Strips  
ID.....{DF52B6F3-54EE-46B3-AFDF-4A701B065E04}  
TERM.....SKIN INCISION WOUND CLOSURE  
STERI-STRIPS

DATA VALUE: SKIN INCISION WOUND CLOSURE SUTURES

-----  
ABBREVIATION.....INC CL SUT  
ACTIVE.....Yes  
DESCRIPTION.....Skin Incision Closure Sutures  
DISPLAY\_NAME.....Sutures  
HELP\_TEXT.....Skin Incision Closure Sutures  
ID.....{DED8DEDA-82C9-4D3A-A81E-641583314919}  
TERM.....SKIN INCISION WOUND CLOSURE SUTURES

DATA VALUE: SKIN INDURATION AROUND INCISION EXTENDS INTO PERIWOUND LESS THAN 1 cm

-----  
ABBREVIATION.....INC IN <1  
ACTIVE.....Yes  
DESCRIPTION.....Skin Induration Extension Around  
Incision Extends into Periwound Less  
Than 1 cm  
DISPLAY\_NAME.....Extends Into Periwound Less Than 1 cm  
HELP\_TEXT.....Skin Induration Extension Around  
Incision Extends into Periwound Less  
Than 1 cm  
ID.....{2A5097BF-955C-4F79-B483-7EE14AF05786}  
TERM.....SKIN INDURATION AROUND INCISION  
EXTENDS INTO PERIWOUND LESS THAN 1 cm

DATA VALUE: SKIN INDURATION AROUND INCISION EXTENDS INTO PERIWOUND MORE THAN 1 cm

-----  
ABBREVIATION.....INC IN >1  
ACTIVE.....Yes  
DESCRIPTION.....Skin Induration Extension Around  
Incision Extends into Periwound  
Greater Than 1 cm  
DISPLAY\_NAME.....Extends into Periwound Greater Than 1  
cm  
HELP\_TEXT.....Skin Induration Extension Around  
Incision Extends into Periwound  
Greater Than 1 cm  
ID.....{D8EE448D-00EB-412D-95AC-2E2E92C53EEB}  
TERM.....SKIN INDURATION AROUND INCISION  
EXTENDS INTO PERIWOUND MORE THAN 1 cm

DATA VALUE: SKIN INDURATION AROUND INCISION LIMITED TO INCISION

-----  
ABBREVIATION.....INC IN LIM  
ACTIVE.....Yes  
DESCRIPTION.....Skin Induration Extension Around  
Incision Limited to Incision  
DISPLAY\_NAME.....Limited to Incision  
HELP\_TEXT.....Skin Induration Extension Around  
Incision Limited to Incision  
ID.....{88BCE7F9-F21F-427A-A50B-C1BA831639AB}  
TERM.....SKIN INDURATION AROUND INCISION  
LIMITED TO INCISION

DATA VALUE: SKIN INDURATION AROUND INCISION NONE

-----  
ABBREVIATION.....INC IN NON  
ACTIVE.....Yes  
DESCRIPTION.....Skin Induration Extension Around  
Incision None  
DISPLAY\_NAME.....None  
HELP\_TEXT.....Skin Induration Extension Around  
Incision None  
ID.....{A9B95680-F9A0-44C0-929C-70C00B9AA5B7}  
TERM.....SKIN INDURATION AROUND INCISION NONE

DATA VALUE: SKIN INDURATION EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND < 1 cm



-----  
ABBREVIATION.....WOU IN <1  
ACTIVE.....Yes  
DESCRIPTION.....No DescriptionSkin Induration  
                  Extension Around Wound Extends into  
                  Periwound Less Than 1 cm  
DISPLAY\_NAME.....Extends into Periwound Less Than 1 cm  
HELP\_TEXT.....No DescriptionSkin Induration  
                  Extension Around Wound Extends into  
                  Periwound Less Than 1 cm  
ID.....{3594CCCD-FA51-4971-A61E-FFA3B4AC7EE0}  
TERM.....SKIN INDURATION EXTENSION AROUND  
                  WOUND EXTENDS INTO PERIWOUND < 1 cm

DATA VALUE: SKIN INDURATION EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND > 1 cm

-----  
ABBREVIATION.....WOU IN >1  
ACTIVE.....Yes  
DESCRIPTION.....Skin Induration Extension Around  
Wound Extends into Periwound More  
Than 1 cm  
DISPLAY\_NAME.....Extends into Periwound More Than 1 cm  
HELP\_TEXT.....Skin Induration Extension Around  
Wound Extends into Periwound More  
Than 1 cm  
ID.....{7DF7D08D-3E5E-4D18-94E4-A0DFDA5DB72F}  
TERM.....SKIN INDURATION EXTENSION AROUND  
WOUND EXTENDS INTO PERIWOUND > 1 cm

DATA VALUE: SKIN INDURATION EXTENSION AROUND WOUND LIMITED TO WOUND

-----  
ABBREVIATION.....WOU IN LIM  
ACTIVE.....Yes  
DESCRIPTION.....Skin Induration Extension Around  
Wound Limited to Wound  
DISPLAY\_NAME.....Limited to Wound  
HELP\_TEXT.....Skin Induration Extension Around  
Wound Limited to Wound  
ID.....{1828A9EA-FF9D-4A49-9543-15B85FBF211E}  
TERM.....SKIN INDURATION EXTENSION AROUND  
WOUND LIMITED TO WOUND

DATA VALUE: SKIN INDURATION EXTENSION AROUND WOUND NONE

-----  
ABBREVIATION.....WOU IN NON  
ACTIVE.....Yes  
DESCRIPTION.....Skin Induration Extension Around  
Wound None  
DISPLAY\_NAME.....None  
HELP\_TEXT.....Skin Induration Extension Around  
Wound None  
ID.....{2C57CD27-4FDD-4DE8-88A0-D438CBE535D8}  
TERM.....SKIN INDURATION EXTENSION AROUND  
WOUND NONE

DATA VALUE: SKIN SWELLING EXTENSION AROUND INCISION EXTENDS INTO PERIWOUND < 1 cm

-----  
ABBREVIATION.....INC SW <1  
ACTIVE.....Yes  
DESCRIPTION.....Skin Swelling Extension Around  
Incision Extends into Periwound Less  
Than 1 cm  
DISPLAY\_NAME.....Extends into Periwound < 1 cm  
HELP\_TEXT.....Skin Swelling Extension Around  
Incision Extends into Periwound Less  
Than 1 cm  
ID.....{ABA0527C-8FAD-4E34-9F6F-2AA97EA06E20}  
TERM.....SKIN SWELLING EXTENSION AROUND  
INCISION EXTENDS INTO PERIWOUND < 1  
cm

DATA VALUE: SKIN SWELLING EXTENSION AROUND INCISION EXTENDS INTO PERIWOUND

> 1 cm

-----  
ABBREVIATION.....INC SW >1  
ACTIVE.....Yes  
DESCRIPTION.....Skin Swelling Extension Around  
                  Incision Extends into Periwound  
                  Greater Than 1 cm  
DISPLAY\_NAME.....Extends into Periwound More Than 1 cm  
HELP\_TEXT.....Skin Swelling Extension Around  
                  Incision Extends into Periwound  
                  Greater Than 1 cm  
ID.....{106A2AD3-F482-4DEE-967E-6F06C6026944}  
TERM.....SKIN SWELLING EXTENSION AROUND  
                  INCISION EXTENDS INTO PERIWOUND > 1  
                  cm

DATA VALUE: SKIN SWELLING EXTENSION AROUND INCISION FLUCTUANT

-----  
ABBREVIATION.....INC SW FLU  
ACTIVE.....Yes  
DESCRIPTION.....Skin Swelling Extension Around  
                  Incision Fluctuant  
DISPLAY\_NAME.....Fluctuant  
HELP\_TEXT.....Skin Swelling Extension Around  
                  Incision Fluctuant  
ID.....{5D38DC7A-8E40-4D27-A039-2D6A557B4FFD}  
TERM.....SKIN SWELLING EXTENSION AROUND  
                  INCISION FLUCTUANT

DATA VALUE: SKIN SWELLING EXTENSION AROUND INCISION LIMITED TO WOUND

-----  
ABBREVIATION.....INC SW LIM  
ACTIVE.....Yes  
DESCRIPTION.....Skin Swelling Extension Around  
                  Incision Limited to Wound  
DISPLAY\_NAME.....Limited to Wound  
HELP\_TEXT.....Skin Swelling Extension Around  
                  Incision Limited to Wound  
ID.....{B8503D54-6925-430E-8AB6-692412347D5D}  
TERM.....SKIN SWELLING EXTENSION AROUND  
                  INCISION LIMITED TO WOUND

DATA VALUE: SKIN SWELLING EXTENSION AROUND INCISION NONE

-----  
ABBREVIATION.....INC SW NON  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....None  
HELP\_TEXT.....No Help Text  
ID.....{54C9FC46-899C-460B-9524-D6F196B9FF31}  
TERM.....SKIN SWELLING EXTENSION AROUND  
                  INCISION NONE

DATA VALUE: SKIN SWELLING EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND < 1  
cm

-----  
ABBREVIATION.....WOU SW <1  
ACTIVE.....Yes  
DESCRIPTION.....Skin Swelling Extension Around Wound  
                  Extends into Periwound Less Than 1 cm  
DISPLAY\_NAME.....Extends into Periwound Less Than 1 cm  
HELP\_TEXT.....Skin Swelling Extension Around Wound  
                  Extends into Periwound Less Than 1 cm  
ID.....{15F314A0-6B77-4EB5-92C0-FFEC4D82C75C}  
TERM.....SKIN SWELLING EXTENSION AROUND WOUND  
                  EXTENDS INTO PERIWOUND < 1 cm

DATA VALUE: SKIN SWELLING EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND > 1  
cm

-----  
ABBREVIATION.....WOU SW >1  
ACTIVE.....Yes  
DESCRIPTION.....Skin Swelling Extension Around Wound  
                  Extends into Periwound Greater Than 1  
                  cm  
DISPLAY\_NAME.....Extends Into Periwound More Than 1 cm

HELP\_TEXT.....Skin Swelling Extension Around Wound  
Extends into Periwound Greater Than 1  
cm  
ID.....{EA6AC933-0737-4C28-9F31-D36F3452DA96}  
TERM.....SKIN SWELLING EXTENSION AROUND WOUND  
EXTENDS INTO PERIWOUND > 1 cm

DATA VALUE: SKIN SWELLING EXTENSION AROUND WOUND FLUCTANT

-----  
ABBREVIATION.....WOU SW FLU  
ACTIVE.....Yes  
DESCRIPTION.....Skin Swelling Extension Around Wound  
                  Fluctuant  
DISPLAY\_NAME.....Fluctant  
HELP\_TEXT.....Skin Swelling Extension Around Wound  
                  Fluctuant  
ID.....{DB5077E7-D339-4C87-B8CA-D16269DB2935}  
TERM.....SKIN SWELLING EXTENSION AROUND WOUND  
                  FLUCTANT

DATA VALUE: SKIN SWELLING EXTENSION AROUND WOUND LIMITED TO WOUND

-----  
ABBREVIATION.....WOU SW LIM  
ACTIVE.....Yes  
DESCRIPTION.....Skin Swelling Extension Around Wound  
                  Limited to Wound  
DISPLAY\_NAME.....Limited to Wound  
HELP\_TEXT.....Skin Swelling Extension Around Wound  
                  Limited to Wound  
ID.....{413CD1EC-309C-45D0-A1E0-E829EA615E2A}  
TERM.....SKIN SWELLING EXTENSION AROUND WOUND  
                  LIMITED TO WOUND

DATA VALUE: SKIN SWELLING EXTENSION AROUND WOUND NONE

-----  
ABBREVIATION.....WOU SW NON  
ACTIVE.....Yes  
DESCRIPTION.....Skin Swelling Extension Around Wound  
                  None  
DISPLAY\_NAME.....None  
HELP\_TEXT.....Skin Swelling Extension Around Wound  
                  None  
ID.....{550E3B92-2C69-403B-A997-2DBE61374ABF}  
TERM.....SKIN SWELLING EXTENSION AROUND WOUND  
                  NONE

DATA VALUE: SKIN TEAR

-----  
ABBREVIATION.....SKIN TEAR  
ACTIVE.....Yes  
DESCRIPTION.....Skin Tear  
DISPLAY\_NAME.....Skin Tear  
HELP\_TEXT.....No Help Text  
ID.....{2E230C6E-119F-4F9F-8797-2FFA9D76F337}  
TERM.....SKIN TEAR

DATA VALUE: SKIN TEMPERATURE AT INCISION SITE CONSISTENT WITH BODY TEMPERAT  
URE

-----  
ABBREVIATION.....INC TP CON  
ACTIVE.....Yes  
DESCRIPTION.....Skin Temperature at Incision Site  
                  Consistent with Body Temperature  
DISPLAY\_NAME.....Skin Temperature at Incision Site  
                  Consistent with Body Temperature  
HELP\_TEXT.....Skin Temperature at Incision Site  
                  Consistent with Body Temperature

ID.....{D349175A-4B20-4243-BE40-E01E50C5583F}  
TERM.....SKIN TEMPERATURE AT INCISION SITE  
CONSISTENT WITH BODY TEMPERATURE

DATA VALUE: SKIN TEMPERATURE AT INCISION SITE COOLER THAN GENERAL BODY TEMPERATURE

-----  
ABBREVIATION.....INC TP CBT  
ACTIVE.....Yes  
DESCRIPTION.....Skin Temperature at Incision Site  
Cooler than General Body Temperature  
DISPLAY\_NAME.....Cooler Than General Body Temperature  
HELP\_TEXT.....Skin Temperature at Incision Site  
Cooler than General Body Temperature  
ID.....{D357C975-857D-4293-B3EF-EFA7A4586357}  
TERM.....SKIN TEMPERATURE AT INCISION SITE  
COOLER THAN GENERAL BODY TEMPERATURE

T DATA VALUE: SKIN TEMPERATURE AT INCISION SITE COOLER THAN OPPOSITE BODY PART

-----  
ABBREVIATION.....INC TP COP  
ACTIVE.....Yes  
DESCRIPTION.....Skin Temperature at Incision Site  
Cooler than Opposite Body Part  
DISPLAY\_NAME.....Cooler Than Opposite Body Part  
HELP\_TEXT.....Skin Temperature at Incision Site  
Cooler than Opposite Body Part  
ID.....{85307AD8-F730-43F0-BCB3-882594D9DDF8}  
TERM.....SKIN TEMPERATURE AT INCISION SITE  
COOLER THAN OPPOSITE BODY PART

DATA VALUE: SKIN TEMPERATURE AT INCISION SITE WARMER THAN GENERAL BODY TEMPERATURE

-----  
ABBREVIATION.....WARM BD OP  
ACTIVE.....Yes  
DESCRIPTION.....Skin Temperature at Incision Site  
Warmer than General Body Temperature  
DISPLAY\_NAME.....Warmer Than General Body Temp  
HELP\_TEXT.....Skin Temperature at Incision Site  
Warmer than General Body Temperature  
ID.....{5F888E9D-8F65-49ED-BC5D-F6F287A9D478}  
TERM.....SKIN TEMPERATURE AT INCISION SITE  
WARMER THAN GENERAL BODY TEMPERATURE

T DATA VALUE: SKIN TEMPERATURE AT INCISION SITE WARMER THAN OPPOSITE BODY PART

-----  
ABBREVIATION.....WARMER OP  
ACTIVE.....Yes  
DESCRIPTION.....Skin Temperature at Incision Site  
Warmer than Opposite Body Part  
DISPLAY\_NAME.....Warmer Than Opposite Body Part  
HELP\_TEXT.....Skin Temperature at Incision Site  
Warmer than Opposite Body Part  
ID.....{F518787A-325A-45F3-82B4-73C856802073}  
TERM.....SKIN TEMPERATURE AT INCISION SITE  
WARMER THAN OPPOSITE BODY PART

DATA VALUE: SKIN TEMPERATURE AT WOUND SITE CONSISTENT WITH BODY TEMPERATURE

-----  
ABBREVIATION.....WOU TP CBT  
ACTIVE.....Yes



DESCRIPTION.....Skin Temperature at Wound Site  
Consistent with Body Temperature  
DISPLAY\_NAME.....Consistent With Body Temperature  
HELP\_TEXT.....Skin Temperature at Wound Site  
Consistent with Body Temperature  
ID.....{672D56C1-D1A3-4DDF-8815-F05DEE197061}  
TERM.....SKIN TEMPERATURE AT WOUND SITE  
CONSISTENT WITH BODY TEMPERATURE

DATA VALUE: SKIN TEMPERATURE AT WOUND SITE COOLER THAN GENERAL BODY TEMPERATURE

-----  
ABBREVIATION.....WOU TP CGB  
ACTIVE.....Yes  
DESCRIPTION.....Skin Temperature at Wound Site Cooler  
                  than General Body Temp  
DISPLAY\_NAME.....Cooler Than General Body Temperature  
HELP\_TEXT.....Skin Temperature at Wound Site Cooler  
                  than General Body Temp  
ID.....{7203993E-3334-44EC-9FD9-14224253283D}  
TERM.....SKIN TEMPERATURE AT WOUND SITE COOLER  
                  THAN GENERAL BODY TEMPERATURE

DATA VALUE: SKIN TEMPERATURE AT WOUND SITE COOLER THAN OPPOSITE BODY PART

-----  
ABBREVIATION.....WOU TP COB  
ACTIVE.....Yes  
DESCRIPTION.....Skin Temperature at Wound Site Cooler  
                  than Opposite Body Part  
DISPLAY\_NAME.....Cooler Than Opposite Body Part  
HELP\_TEXT.....Skin Temperature at Wound Site Cooler  
                  than Opposite Body Part  
ID.....{A904C666-2810-4211-B53A-2961DCE6894E}  
TERM.....SKIN TEMPERATURE AT WOUND SITE COOLER  
                  THAN OPPOSITE BODY PART

DATA VALUE: SKIN TEMPERATURE AT WOUND SITE WARMER THAN GENERAL BODY TEMPERATURE

-----  
ABBREVIATION.....WOU TP WGB  
ACTIVE.....Yes  
DESCRIPTION.....Skin Temperature at Wound Site Warmer  
                  than General Body Temperature  
DISPLAY\_NAME.....Warmer Than General Body Temperature  
HELP\_TEXT.....Skin Temperature at Wound Site Warmer  
                  than General Body Temperature  
ID.....{03BB8F9A-3DAC-4994-9F09-AA1DEEF48090}  
TERM.....SKIN TEMPERATURE AT WOUND SITE WARMER  
                  THAN GENERAL BODY TEMPERATURE

DATA VALUE: SKIN TEMPERATURE AT WOUND SITE WARMER THAN OPPOSITE BODY PART

-----  
ABBREVIATION.....WOU TP WOB  
ACTIVE.....Yes  
DESCRIPTION.....Skin Temperature at Wound Site Warmer  
                  than Opposite Body Part  
DISPLAY\_NAME.....Warmer Than Opposite Body Part  
HELP\_TEXT.....Skin Temperature at Wound Site Warmer  
                  than Opposite Body Part  
ID.....{6290BD23-24BF-4A3A-8555-DCE52D2274B3}  
TERM.....SKIN TEMPERATURE AT WOUND SITE WARMER  
                  THAN OPPOSITE BODY PART

DATA VALUE: SKIN TRAUMA

-----  
ABBREVIATION.....SKIN TRAUM  
ACTIVE.....Yes  
DESCRIPTION.....Skin Trauma  
DISPLAY\_NAME.....Skin Trauma

HELP\_TEXT.....No Help Text  
ID.....{91C8EE13-2C4A-468C-94EB-C53E0D153EBA}  
TERM.....SKIN TRAUMA

DATA VALUE: SKIN WARM AND DRY COLOR APPROPRIATE

-----  
ABBREVIATION.....S W A DR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Skin Warm and Dry, Color Appropriate  
to Baseline  
HELP\_TEXT.....No Help Text  
ID.....{9637D510-C199-4969-8F2D-0A499151E22B}  
TERM.....SKIN WARM AND DRY COLOR APPROPRIATE

DATA VALUE: SKIN WOUND CLOSURE APPROXIMATED

-----  
ABBREVIATION.....WOU CL APP  
ACTIVE.....Yes  
DESCRIPTION.....Skin Wound Closure Approximated  
DISPLAY\_NAME.....Approximated  
HELP\_TEXT.....Skin Wound Closure Approximated  
ID.....{83D48E0F-8643-436B-9ECC-80795857C322}  
TERM.....SKIN WOUND CLOSURE APPROXIMATED

DATA VALUE: SKIN WOUND CLOSURE DEHISCENCE

-----  
ABBREVIATION.....WOU CL DEH  
ACTIVE.....Yes  
DESCRIPTION.....Skin Wound Closure Dehiscence  
DISPLAY\_NAME.....Dehiscence  
HELP\_TEXT.....Skin Wound Closure Dehiscence  
ID.....{425A460F-19C0-4A4A-8E4E-712C66968C38}  
TERM.....SKIN WOUND CLOSURE DEHISCENCE

DATA VALUE: SKIN WOUND CLOSURE GLUE

-----  
ABBREVIATION.....WOU CL GLU  
ACTIVE.....Yes  
DESCRIPTION.....Skin Wound Closure Glue  
DISPLAY\_NAME.....Glue  
HELP\_TEXT.....Skin Wound Closure Glue  
ID.....{29FE4AF7-39EE-400A-BB3D-E07786F11BFD}  
TERM.....SKIN WOUND CLOSURE GLUE

DATA VALUE: SKIN WOUND CLOSURE NOT CLOSED

-----  
ABBREVIATION.....WOU CL NOT  
ACTIVE.....Yes  
DESCRIPTION.....Skin Wound Closure Not Closed  
DISPLAY\_NAME.....Not Closed  
HELP\_TEXT.....Skin Wound Closure Not Closed  
ID.....{0E69F1FA-35F6-4F34-8547-CC3602670A7B}  
TERM.....SKIN WOUND CLOSURE NOT CLOSED

DATA VALUE: SKIN WOUND CLOSURE OPEN WOUND

-----  
ABBREVIATION.....WOU CL OPE  
ACTIVE.....Yes  
DESCRIPTION.....Skin Wound Closure Open Wound  
DISPLAY\_NAME.....Open Wound  
HELP\_TEXT.....Skin Wound Closure Open Wound  
ID.....{C6F2AD2C-932F-434A-8A7A-513029092FFD}  
TERM.....SKIN WOUND CLOSURE OPEN WOUND

DATA VALUE: SKIN WOUND CLOSURE STAPLES

-----  
ABBREVIATION.....WOU CL STA  
ACTIVE.....Yes  
DESCRIPTION.....Skin Wound Closure Staples  
DISPLAY\_NAME.....Staples  
HELP\_TEXT.....Skin Wound Closure Staples  
ID.....{FE04BFE3-DF81-48E7-9FA1-1CF4ECAD71}  
TERM.....SKIN WOUND CLOSURE STAPLES

DATA VALUE: SKIN WOUND CLOSURE STERI-STRIPS

-----  
ABBREVIATION.....WOU CL STE  
ACTIVE.....Yes  
DESCRIPTION.....Skin Wound Closure Steri - Strips  
DISPLAY\_NAME.....Steri-Strips  
HELP\_TEXT.....Skin Wound Closure Steri-Strips  
ID.....{AD451791-6B49-4A51-ABC3-796D08BEF910}  
TERM.....SKIN WOUND CLOSURE STERI-STRIPS

DATA VALUE: SKIN WOUND CLOSURE SUTURES

-----  
ABBREVIATION.....WOU CL SUT  
ACTIVE.....Yes  
DESCRIPTION.....Skin Wound Closure Sutures  
DISPLAY\_NAME.....Sutures  
HELP\_TEXT.....Skin Wound Closure Sutures  
ID.....{F2840AE4-216B-477C-9FE6-DC91E8CE5364}  
TERM.....SKIN WOUND CLOSURE SUTURES

DATA VALUE: SLEEPING

-----  
ABBREVIATION.....s  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Sleeping  
HELP\_TEXT.....No Help Text  
ID.....{87338420-F351-F98C-7A80-4990B87ACE2B}  
TERM.....SLEEPING

DATA VALUE: SLEEPS POORLY

-----  
ABBREVIATION.....sp  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Sleeps Poorly  
HELP\_TEXT.....No Help Text  
ID.....{AB9424B4-D0B1-148B-C800-FF535F3FAC84}  
TERM.....SLEEPS POORLY

DATA VALUE: SLEPT WELL

-----  
ABBREVIATION.....SLEPT WELL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....States Has Slept Well  
HELP\_TEXT.....No Help Text  
ID.....{18573E82-568F-4D2E-B9E1-6E4589EBEE2D}  
TERM.....SLEPT WELL

DATA VALUE: SLIGHT

-----  
ABBREVIATION.....SLIGHT  
ACTIVE.....Yes  
DESCRIPTION.....Slight  
DISPLAY\_NAME.....Slight  
HELP\_TEXT.....Slight  
ID.....{805549C1-5DC5-4A2D-BC1F-046E44928843}  
TERM.....SLIGHT

DATA VALUE: SLING TRACTION

-----  
ABBREVIATION.....SLING TRAC  
ACTIVE.....Yes  
DESCRIPTION.....Sling Traction  
DISPLAY\_NAME.....Sling Traction  
HELP\_TEXT.....Sling Traction  
ID.....{304E924C-C95A-4DBE-B30B-191BC4D26EF7}  
TERM.....SLING TRACTION

DATA VALUE: SLOUGH

-----  
ABBREVIATION.....SLOUGH  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Slough  
HELP\_TEXT.....No Help Text  
ID.....{548EB536-1AC4-4234-8C4D-DC54E4CAA6F6}  
TERM.....SLOUGH

DATA VALUE: SLUGGISH

-----  
ABBREVIATION.....s  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Sluggish  
HELP\_TEXT.....No Help Text  
ID.....{3215A82E-8DCC-C7C0-8AE1-2D4C9D4A4EFB}  
TERM.....SLUGGISH

DATA VALUE: SLURRED SPEECH

-----  
ABBREVIATION.....Slurred Sp  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Slurred Speech  
HELP\_TEXT.....  
ID.....{6C3D295A-71C9-4DA8-A172-995F431D94C9}  
TERM.....SLURRED SPEECH

DATA VALUE: SMALL

-----  
ABBREVIATION.....s  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Small  
HELP\_TEXT.....No Help Text  
ID.....{5991F892-B28C-3998-E733-D8FFA26D7E21}  
TERM.....SMALL

DATA VALUE: SMALL BORE FEEDING TUBE

-----  
ABBREVIATION.....SM BO TUBE  
ACTIVE.....Yes  
DESCRIPTION.....Small Bore Feeding Tube  
DISPLAY\_NAME.....Small Bore Feeding Tube  
HELP\_TEXT.....No Help Text  
ID.....{7E347B8A-14E7-4093-9A5B-56E6AC0D78DF}  
TERM.....SMALL BORE FEEDING TUBE

DATA VALUE: SMEAR

-----  
ABBREVIATION.....SMEAR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Smear  
HELP\_TEXT.....No Help Text  
ID.....{4B99BBD3-443E-4954-8ADE-5FB3EE50540A}  
TERM.....SMEAR

DATA VALUE: SMOKING AREA

-----  
ABBREVIATION.....sa  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Smoking Area  
HELP\_TEXT.....No Help Text  
ID.....{8AFDBC0-9FC8-523F-5D07-4D8CCBA7036F}  
TERM.....SMOKING AREA

DATA VALUE: SMOKING PORCH

-----  
ABBREVIATION.....sp  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Smoking Porch  
HELP\_TEXT.....No Help Text  
ID.....{5AD3F05B-2A5B-E2B4-FD47-66712CFA44E6}  
TERM.....SMOKING PORCH

DATA VALUE: SMOOTH

-----  
ABBREVIATION.....SMOOTH  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Smooth  
HELP\_TEXT.....No Help Text  
ID.....{5D71A0D5-52ED-4B8D-B359-AC7D500660D8}  
TERM.....SMOOTH

DATA VALUE: SNACK

-----  
ABBREVIATION.....SNACK  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Snack  
HELP\_TEXT.....No Help Text  
ID.....{60B77974-4102-4511-A105-E24850B924D8}  
TERM.....SNACK

DATA VALUE: SNORES

-----  
ABBREVIATION.....Snores  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Snores  
HELP\_TEXT.....  
ID.....{69FE4ECF-E1FE-4235-AD44-5AF9E634E000}  
TERM.....SNORES



DATA VALUE: SODIUM CONTROLLED 2 GM

-----  
ABBREVIATION.....SOD CONT 2  
ACTIVE.....Yes  
DESCRIPTION.....Sodium Controlled 2 GM  
DISPLAY\_NAME.....Sodium Controlled 2 GM  
HELP\_TEXT.....No Help Text  
ID.....{8A1311AF-6C7D-407A-962F-0FAE6C341F9A}  
TERM.....SODIUM CONTROLLED 2 GM

DATA VALUE: SODIUM CONTROLLED 4 GM

-----  
ABBREVIATION.....SOD CON 4G  
ACTIVE.....Yes  
DESCRIPTION.....Sodium Controlled 4 GM Diet  
DISPLAY\_NAME.....Sodium Controlled 4 GM  
HELP\_TEXT.....No Help Text  
ID.....{26538A2F-5502-4D94-A2A7-D9E0BF825C47}  
TERM.....SODIUM CONTROLLED 4 GM

DATA VALUE: SOFT

-----  
ABBREVIATION.....SOFT  
ACTIVE.....Yes  
DESCRIPTION.....Soft  
DISPLAY\_NAME.....Soft  
HELP\_TEXT.....No Help Text  
ID.....{E68B4D4B-B86B-4004-A861-986645267F6E}  
TERM.....SOFT

DATA VALUE: SOMNOLENT

-----  
ABBREVIATION.....Somnolent  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Somnolent  
HELP\_TEXT.....  
ID.....{008A3F20-DD3D-4775-AA18-CCA7B817A77D}  
TERM.....SOMNOLENT

DATA VALUE: SPECIAL PROCEDURE ROOM

-----  
ABBREVIATION.....Special Rm  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Special Procedure Room  
HELP\_TEXT.....  
ID.....{60AA3B10-B3EE-42EF-9CEC-0CAE66DFE164}  
TERM.....SPECIAL PROCEDURE ROOM

DATA VALUE: SPECIALTY MATTRESS/BED

-----  
ABBREVIATION.....Specialty  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Specialty Mattress/Bed  
HELP\_TEXT.....  
ID.....{6B84BE9C-9AAC-430B-8F0D-0001E2F2AAC7}  
TERM.....SPECIALTY MATTRESS/BED

DATA VALUE: SPINAL

-----  
ABBREVIATION.....Spinal  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Spinal  
HELP\_TEXT.....  
ID.....{AC02181C-0A2C-4F7A-B983-B32EF1BCD150}  
TERM.....SPINAL

DATA VALUE: SPIRITUAL

-----  
ABBREVIATION.....SPIRITUAL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Spiritual  
HELP\_TEXT.....No Help Text  
ID.....{9670E2B9-6938-4487-A857-25546775B7A2}  
TERM.....SPIRITUAL

DATA VALUE: SPLINTING

-----  
ABBREVIATION.....SPLINTING  
ACTIVE.....Yes  
DESCRIPTION.....Splinting  
DISPLAY\_NAME.....Splinting  
HELP\_TEXT.....No Help Text  
ID.....{D7FC5AA9-E63A-4A88-BE03-A50D0A8751F9}  
TERM.....SPLINTING

DATA VALUE: SPONTANEOUS

-----  
ABBREVIATION.....SPONTANEOU  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Spontaneous  
HELP\_TEXT.....No Help Text  
ID.....{C137C4A7-08A6-4064-8F95-4A4476E2E3A8}  
TERM.....SPONTANEOUS

DATA VALUE: SPOONING

-----  
ABBREVIATION.....SPOONING  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Spoonng  
HELP\_TEXT.....No Help Text  
ID.....{CA384C9F-D87D-4E2D-B4C0-89BD85692B16}  
TERM.....SPOONING

DATA VALUE: SPOUSE

-----  
ABBREVIATION.....Spouse  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Spouse  
HELP\_TEXT.....  
ID.....{2C601365-23F9-4939-9B63-F73720A17199}  
TERM.....SPOUSE

DATA VALUE: SQ EMPHYSEMA

-----  
ABBREVIATION.....SQ EMPHYS.  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Subcutaneous Emphysema  
HELP\_TEXT.....No Help Text  
ID.....{B8D03538-0168-4CDC-FF41-BA1185EBCC87}  
TERM.....SQ EMPHYSEMA

DATA VALUE: SQUARE

-----  
ABBREVIATION.....SQUARE  
ACTIVE.....Yes  
DESCRIPTION.....Square  
DISPLAY\_NAME.....Square  
HELP\_TEXT.....Square  
ID.....{63F70678-2CAF-48E8-964D-2805D8B9A43B}  
TERM.....SQUARE

DATA VALUE: SQUEEZING

-----  
ABBREVIATION.....SQUEEZING  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....SQUEEZING  
HELP\_TEXT.....No Help Text  
ID.....{A148587A-1B99-4D0D-B11E-87ADC3D2A811}  
TERM.....SQUEEZING

DATA VALUE: STABBING

-----  
ABBREVIATION.....STABBING  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Stabbing  
HELP\_TEXT.....No Help Text  
ID.....{684640D5-B3AB-47DA-98D8-5434ABF22B7E}  
TERM.....STABBING

DATA VALUE: STABLE

-----  
ABBREVIATION.....STABLE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Stable  
HELP\_TEXT.....No Help Text  
ID.....{A474E7BC-47EF-4472-A19A-172C441099EC}  
TERM.....STABLE

DATA VALUE: STABLE BP AND PULSE

-----  
ABBREVIATION.....Stable BP  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Stable BP and Pulse  
HELP\_TEXT.....  
ID.....{F07C2F54-1182-4377-A33A-F529884013E8}  
TERM.....STABLE BP AND PULSE

DATA VALUE: STAFF

-----  
ABBREVIATION.....STAFF  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Staff  
HELP\_TEXT.....No Help Text  
ID.....{BE8D678F-CE19-41FC-8C12-7A5DE7F0C7E2}  
TERM.....STAFF

DATA VALUE: STAGE 3 OR 4 PRESSURE ULCER

-----  
ABBREVIATION.....Press Ulce  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Stage 3 or 4 Pressure Ulcer  
HELP\_TEXT.....  
ID.....{883029EF-BB3A-4146-95DD-3C6F1C0D097D}  
TERM.....STAGE 3 OR 4 PRESSURE ULCER

DATA VALUE: STAGE I

-----  
ABBREVIATION.....STAGE I  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Stage I  
HELP\_TEXT.....No Help Text  
ID.....{B376AAAF-7143-4422-BCA2-CACF1201817C}  
TERM.....STAGE I

DATA VALUE: STAGE II

-----  
ABBREVIATION.....STAGE II  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Stage II  
HELP\_TEXT.....No Help Text  
ID.....{B520BEA3-841E-490A-A3CD-63D1027B561F}  
TERM.....STAGE II

DATA VALUE: STAGE III

-----  
ABBREVIATION.....STAGE III  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Stage III  
HELP\_TEXT.....No Help Text  
ID.....{A2F0866C-ECEE-4FCD-A34B-498232662D7E}  
TERM.....STAGE III

DATA VALUE: STAGE IV

-----  
ABBREVIATION.....STAGE IV  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Stage IV  
HELP\_TEXT.....No Help Text  
ID.....{284334AA-E488-4D45-93EF-9B93B2086A63}  
TERM.....STAGE IV

DATA VALUE: STANDARD

-----  
ABBREVIATION.....STANDARD  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....STANDARD  
HELP\_TEXT.....No Help Text  
ID.....{C4B8C06D-BFFE-4438-A43D-926F039E5EE1}  
TERM.....STANDARD

DATA VALUE: STANDARD PRESSURE REDUCTION

-----  
ABBREVIATION.....STAN P RED  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Standard Pressure Reduction  
HELP\_TEXT.....No Help Text  
ID.....{556537CD-A84D-4493-81F6-51CBD844B5B7}  
TERM.....STANDARD PRESSURE REDUCTION

DATA VALUE: STAPLES

-----  
ABBREVIATION.....STAPLES  
ACTIVE.....Yes  
DESCRIPTION.....Staples to close a wound  
DISPLAY\_NAME.....Staples  
HELP\_TEXT.....No Help Text  
ID.....{344DBBE0-5A94-4E18-BC97-BD7767A83E36}  
TERM.....STAPLES

DATA VALUE: START TIME:

-----  
ABBREVIATION.....Start Time  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Start Time:  
HELP\_TEXT.....When entering data, make sure the  
                  save date/time reflects the start  
                  time.  
ID.....{4D0190FA-FFE4-4D5F-8A15-0517E8AB9B14}  
TERM.....START TIME:

DATA VALUE: STARTED

-----  
ABBREVIATION.....STARTED  
ACTIVE.....Yes  
DESCRIPTION.....Started  
DISPLAY\_NAME.....Started  
HELP\_TEXT.....No Help Text  
ID.....{D6E105C4-5EAD-424E-8E4C-195EB7783A3D}  
TERM.....STARTED

DATA VALUE: STARTED ON LEFT

-----  
ABBREVIATION.....sol  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Started On Left  
HELP\_TEXT.....No Help Text  
ID.....{AA6767EE-57F6-F565-76E9-EB219644492F}  
TERM.....STARTED ON LEFT

DATA VALUE: STARTED ON RIGHT

-----  
ABBREVIATION.....sor  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Started On Right  
HELP\_TEXT.....No Help Text  
ID.....{0C690E36-E710-244D-44FC-5DB5BCC0FF7E}  
TERM.....STARTED ON RIGHT

DATA VALUE: STEADY

-----  
ABBREVIATION.....Steady  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Steady  
HELP\_TEXT.....  
ID.....{42CF2FDA-E774-4784-99D3-24545E30C46F}  
TERM.....STEADY

DATA VALUE: STEADY GATE OR BASELINE

-----  
ABBREVIATION.....ST GA BA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Steady Gate or Baseline  
HELP\_TEXT.....No Help Text  
ID.....{0D47B2CE-88A7-4FE5-964B-7A8107A5DB91}  
TERM.....STEADY GATE OR BASELINE

DATA VALUE: STEP DOWN

-----  
ABBREVIATION.....STEP DOW  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Step Down  
HELP\_TEXT.....No Help Text  
ID.....{F6AC5CA7-0953-4E21-B1A4-6BFE1983D515}  
TERM.....STEP DOWN

DATA VALUE: STERILE GLOVES-DRAPES USED

-----  
ABBREVIATION.....Ster Glove  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Sterile Gloves-Drapes Used  
HELP\_TEXT.....  
ID.....{F26BBB1B-1483-40CD-AC47-828F71B9187C}  
TERM.....STERILE GLOVES-DRAPES USED

DATA VALUE: STINGING

-----  
ABBREVIATION.....STINGING  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Stinging  
HELP\_TEXT.....No Help Text  
ID.....{DF235170-F501-4E8D-89DD-5A3C76F4B13B}  
TERM.....STINGING

DATA VALUE: STOP TIME:

-----  
ABBREVIATION.....Stop Time:  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Stop Time:  
HELP\_TEXT.....When entering data, make sure the  
save date/time reflects the stop  
time.  
ID.....{784CA4D8-988F-4A59-8AB5-48AFC9CA210B}  
TERM.....STOP TIME:

DATA VALUE: STOPCOCK CHANGE

-----  
ABBREVIATION.....Stop Chg  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Stopcock Change  
HELP\_TEXT.....  
ID.....{C4A6E58E-DA4B-478F-823B-D2CFC9860072}  
TERM.....STOPCOCK CHANGE

DATA VALUE: STRAIGHT CATHETER

-----  
ABBREVIATION.....STR. CATH  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Straight Catheter  
HELP\_TEXT.....No Help Text  
ID.....{5F4A8A7F-E720-9458-2360-1DC14E7D05F3}  
TERM.....STRAIGHT CATHETER

DATA VALUE: STRETCHER

-----  
ABBREVIATION.....STRETCHER  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Stretcher  
HELP\_TEXT.....No Help Text  
ID.....{93CB521A-3A67-41FB-A59C-3DD54D810B11}  
TERM.....STRETCHER



DATA VALUE: STRICT

-----  
ABBREVIATION.....STRICT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Strict  
HELP\_TEXT.....No Help Text  
ID.....{5717A3EA-9A1F-44AB-BB0B-E86BBDED3ECF}  
TERM.....STRICT

DATA VALUE: STROKE

-----  
ABBREVIATION.....Stroke  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Stroke  
HELP\_TEXT.....  
ID.....{F0EE98A8-040F-4923-8A21-059DA7B17995}  
TERM.....STROKE

DATA VALUE: STRONG

-----  
ABBREVIATION.....s  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Strong  
HELP\_TEXT.....No Help Text  
ID.....{A12154E1-544E-F372-3839-0D44E5BC476B}  
TERM.....STRONG

DATA VALUE: STRYKER FRAME

-----  
ABBREVIATION.....STRYKER FR  
ACTIVE.....Yes  
DESCRIPTION.....Stryker Fram  
DISPLAY\_NAME.....Stryker Frame  
HELP\_TEXT.....No Help Text  
ID.....{71A1DB72-C860-42D8-BA4F-B71AFA2164B1}  
TERM.....STRYKER FRAME

DATA VALUE: STUDENT

-----  
ABBREVIATION.....STUDENT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Student  
HELP\_TEXT.....No Help Text  
ID.....{613AFF67-9121-4F6C-A44F-1ACFEF10E612}  
TERM.....STUDENT

DATA VALUE: STUPOROUS

-----  
ABBREVIATION.....Stuporous  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Stuporous  
HELP\_TEXT.....  
ID.....{4C25FF58-3346-4B06-ACE5-6ADB942877A8}  
TERM.....STUPOROUS

DATA VALUE: SUBCUTANEOUS TISSUE

-----  
ABBREVIATION.....SUBCUTAN T  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Subcutaneous Tissue  
HELP\_TEXT.....No Help Text  
ID.....{D1345ED8-2CDB-4FD7-84B8-92C9E0872543}  
TERM.....SUBCUTANEOUS TISSUE

DATA VALUE: SUCTIONED ONLY

-----  
ABBREVIATION.....SUCTIONED ONLY  
ACTIVE.....Yes  
DESCRIPTION.....Suctioned Only  
DISPLAY\_NAME.....Suctioned Only  
HELP\_TEXT.....No Help Text  
ID.....{7A8B23FD-69C5-4931-A6D4-E92F8FD70817}  
TERM.....SUCTIONED ONLY

DATA VALUE: SUICIDAL

-----  
ABBREVIATION.....s  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Suicidal  
HELP\_TEXT.....No Help Text  
ID.....{B104E9C7-A78E-82DC-49A4-19B94A639AFE}  
TERM.....SUICIDAL

DATA VALUE: SUMP

-----  
ABBREVIATION.....Sump  
ACTIVE.....Yes  
DESCRIPTION.....Sump Drain  
DISPLAY\_NAME.....Sump Drain  
HELP\_TEXT.....No Help Text  
ID.....{DAD881B9-F1EE-FDAD-AEB0-68C754FDC24E}  
TERM.....SUMP

DATA VALUE: SUPERVISION/SET-UP HELP

-----  
ABBREVIATION.....SUPERVISED  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Supervision/Set-Up Help  
HELP\_TEXT.....No Help Text  
ID.....{596D028B-B60D-4244-A1B0-882177C8D25F}  
TERM.....SUPERVISION/SET-UP HELP

DATA VALUE: SUPINE

-----  
ABBREVIATION.....SUPINE  
ACTIVE.....Yes  
DESCRIPTION.....Supine  
DISPLAY\_NAME.....Supine  
HELP\_TEXT.....No Help Text  
ID.....{A18F9534-17AF-47D7-8D14-2AA457A889BD}  
TERM.....SUPINE

DATA VALUE: SUPLENA

-----  
ABBREVIATION.....Suplena  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Suplena  
HELP\_TEXT.....  
ID.....{F606979A-083E-40C1-9EBB-EA56053EF5FD}  
TERM.....SUPLENA

DATA VALUE: SUPPORT GROUPS

-----  
ABBREVIATION.....Supp Group  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Support Groups  
HELP\_TEXT.....  
ID.....{7D67F542-4210-42FE-B083-7DA9D29055BD}  
TERM.....SUPPORT GROUPS

DATA VALUE: SUPRAPUBIC

-----  
ABBREVIATION.....SUPRAPUBIC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Suprapubic  
HELP\_TEXT.....No Help Text  
ID.....{69E0176F-52FC-A3C4-7A4E-AD57DBD6E572}  
TERM.....SUPRAPUBIC

DATA VALUE: SUPRAVENTRICULAR TACHYCARDIA

-----  
ABBREVIATION.....st  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Supraventricular Tachycardia  
HELP\_TEXT.....No Help Text  
ID.....{F712412F-D39A-58F7-E4F1-6AC8B77445D7}  
TERM.....SUPRAVENTRICULAR TACHYCARDIA

DATA VALUE: SURGICAL

-----  
ABBREVIATION.....SURGICAL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Surgical  
HELP\_TEXT.....No Help-Text  
ID.....{2C6CFEF1-6122-45BA-A9F8-F33B5B3C0E90}  
TERM.....SURGICAL

DATA VALUE: SURGICAL DONOR

-----  
ABBREVIATION.....SURG DONOR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Surgical Donor  
HELP\_TEXT.....No Help Text  
ID.....{F3EE10EA-5776-479B-A58C-F59D6526F90B}  
TERM.....SURGICAL DONOR

DATA VALUE: SURGICAL GRAFT

-----  
ABBREVIATION.....SURG GRAFT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Surgical Graft  
HELP\_TEXT.....No Help Text  
ID.....{328BA3EA-409B-4C5E-92E7-C6ECE8F983E1}  
TERM.....SURGICAL GRAFT

DATA VALUE: SUSPECTED DEEP TISSUE INJURY

-----  
ABBREVIATION.....SUS DEEP T  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Suspected Deep Tissue Injury  
HELP\_TEXT.....No Help Text  
ID.....{E82B4082-646F-4B5A-B68E-BB0872A4F6D2}  
TERM.....SUSPECTED DEEP TISSUE INJURY

DATA VALUE: SUTURES

-----  
ABBREVIATION.....SUTURES  
ACTIVE.....Yes  
DESCRIPTION.....Sutures to close wound  
DISPLAY\_NAME.....Sutures  
HELP\_TEXT.....No Help Text  
ID.....{F6E3535E-2368-4014-8F4E-670EC0CB69A5}  
TERM.....SUTURES

DATA VALUE: SWOLLEN

-----  
ABBREVIATION.....SWOLLEN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Swollen  
HELP\_TEXT.....No Help Text  
ID.....{DCE86B10-C393-8DA9-6518-F2943F0E1A4E}  
TERM.....SWOLLEN

DATA VALUE: SYMPTOMS OF MEDICATION ALLERGY

-----  
ABBREVIATION.....Sym Med Al  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Symptoms of Medication Allergy  
HELP\_TEXT.....  
ID.....{8E73F22D-7583-4F24-B946-091CF8508CED}  
TERM.....SYMPTOMS OF MEDICATION ALLERGY

DATA VALUE: SYNCHRONOUS INTERMITTENT MANDATORY VENTILATION

-----  
ABBREVIATION.....SIMV  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Synchronous Intermittent Mandatory  
Ventilation  
HELP\_TEXT.....No Help Text  
ID.....{A606BB2B-5C00-FEA1-127F-1D3455AE5996}  
TERM.....SYNCHRONOUS INTERMITTENT MANDATORY  
VENTILATION

DATA VALUE: SYRINGE

-----  
ABBREVIATION.....Syringe  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Syringe  
HELP\_TEXT.....  
ID.....{8ED1EF3E-79AC-4507-A35B-473F31510FB7}  
TERM.....SYRINGE

DATA VALUE: SYRINGE CHANGE

-----  
ABBREVIATION.....Syr Chg  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Syringe Change  
HELP\_TEXT.....  
ID.....{1F052EF1-F62C-4442-8B2B-D15280226C69}  
TERM.....SYRINGE CHANGE

DATA VALUE: SYRINGE CHANGED

-----  
ABBREVIATION.....SYR CHANGE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Syringe Changed  
HELP\_TEXT.....No Help Text  
ID.....{BE84F060-7E03-49DD-A6BC-787375C8E1F1}  
TERM.....SYRINGE CHANGED

DATA VALUE: SYSTOLIC

-----  
ABBREVIATION.....SYSTOLIC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Systolic  
HELP\_TEXT.....No Help Text  
ID.....{859B5700-0990-42FF-8129-8DB24EFDA8D1}  
TERM.....SYSTOLIC

DATA VALUE: SaO2% AT ACCEPTABLE LEVEL BY ORDERS

-----  
ABBREVIATION.....SaO2 by Or  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....SaO2% at Acceptable Level by Orders  
HELP\_TEXT.....  
ID.....{408B35FD-0300-4E7E-922A-A3F8CADC3DA9}  
TERM.....SaO2% AT ACCEPTABLE LEVEL BY ORDERS

DATA VALUE: SaO2% AT ACCEPTABLE LEVEL WITH PRESCRIBED O2 LEVEL

-----  
ABBREVIATION.....SaO2 w/Pre  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....SaO2% at Acceptable Level with  
Prescribed O2 Level  
HELP\_TEXT.....  
ID.....{8DE3440F-8049-4846-89D6-2E2DDEEA362C}  
TERM.....SaO2% AT ACCEPTABLE LEVEL WITH  
PRESCRIBED O2 LEVEL

DATA VALUE: T1

-----  
ABBREVIATION.....T1  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....T1  
HELP\_TEXT.....  
ID.....{78DB9C11-B0FE-4A08-A3AE-8BEA2CA4D91F}  
TERM.....T1

DATA VALUE: T10

-----  
ABBREVIATION.....T10  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....T10  
HELP\_TEXT.....  
ID.....{EB6E695A-95AF-4E9D-AEFA-B1484BABFC4E}  
TERM.....T10

DATA VALUE: T11

-----  
ABBREVIATION.....T11  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....T11  
HELP\_TEXT.....  
ID.....{3C40B21A-965B-48C4-A65D-9A01CF2327F1}  
TERM.....T11

DATA VALUE: T12

-----  
ABBREVIATION.....T12  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....T12  
HELP\_TEXT.....  
ID.....{36FA36A9-FBB2-4F16-A3BA-7E3260DA768D}  
TERM.....T12

DATA VALUE: T2

-----  
ABBREVIATION.....T2  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....T2  
HELP\_TEXT.....  
ID.....{7536B154-F3CA-4B21-8CEC-301F204B19A9}  
TERM.....T2

DATA VALUE: T3

-----  
ABBREVIATION.....T3  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....T3  
HELP\_TEXT.....  
ID.....{F8BDC550-D9B9-46CD-9BD1-52E32AA513FD}  
TERM.....T3

DATA VALUE: T4

-----  
ABBREVIATION.....T4  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....T4  
HELP\_TEXT.....  
ID.....{E95970BB-3FDA-4027-AAF0-EE5670642321}  
TERM.....T4

DATA VALUE: T5

-----  
ABBREVIATION.....T5  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....T5  
HELP\_TEXT.....  
ID.....{01B568D4-F9F4-44A7-94A4-EFF830691940}  
TERM.....T5

DATA VALUE: T6

-----  
ABBREVIATION.....T6  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....T6  
HELP\_TEXT.....  
ID.....{0C2C1A00-7CCB-4608-AECB-30BFB52A5342}  
TERM.....T6

DATA VALUE: T7

-----  
ABBREVIATION.....T7  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....T7  
HELP\_TEXT.....  
ID.....{385B254A-F256-4626-9C12-7D8984F4B605}  
TERM.....T7



DATA VALUE: T8

-----  
ABBREVIATION.....T8  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....T8  
HELP\_TEXT.....  
ID.....{68337EB5-902D-4DCB-AFDD-8AC59849F58B}  
TERM.....T8

DATA VALUE: T9

-----  
ABBREVIATION.....T9  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....T9  
HELP\_TEXT.....  
ID.....{06F10D8D-2E9D-454D-AD6B-B750300523DD}  
TERM.....T9

DATA VALUE: TAN

-----  
ABBREVIATION.....TAN  
ACTIVE.....Yes  
DESCRIPTION.....The color tan  
DISPLAY\_NAME.....Tan  
HELP\_TEXT.....The color tan  
ID.....{AB994ED7-8916-FAE0-0E97-CD4E0EF80ED1}  
TERM.....TAN

DATA VALUE: TAPE BURN

-----  
ABBREVIATION.....TAPE BURN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Tape Burn  
HELP\_TEXT.....No Help Text  
ID.....{DDD75AFD-DA8D-4E73-873F-86D4DF4AF47A}  
TERM.....TAPE BURN

DATA VALUE: TARRY

-----  
ABBREVIATION.....TARRY  
ACTIVE.....Yes  
DESCRIPTION.....Tarry in appearance  
DISPLAY\_NAME.....Tarry  
HELP\_TEXT.....No Help Text  
ID.....{09364271-6307-471F-A334-A8F7E051EF46}  
TERM.....TARRY

DATA VALUE: TEA

-----  
ABBREVIATION.....TEA  
ACTIVE.....Yes  
DESCRIPTION.....Tea  
DISPLAY\_NAME.....Tea  
HELP\_TEXT.....Tea  
ID.....{E2C0C287-DFDC-4F7B-8A67-7A339B8E8669}  
TERM.....TEA

DATA VALUE: TECH

-----  
ABBREVIATION.....Tech  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Tech  
HELP\_TEXT.....  
ID.....{13AAF5D5-FEBC-4E09-B972-5D86BA02F307}  
TERM.....TECH

DATA VALUE: TEETH BRUSHED

-----  
ABBREVIATION.....TEETH BRUS  
ACTIVE.....Yes  
DESCRIPTION.....Teeth Brushed  
DISPLAY\_NAME.....Teeth Brushed  
HELP\_TEXT.....No Help Text  
ID.....{20841794-6E6D-4262-B1FE-9E3F20F8CC49}  
TERM.....TEETH BRUSHED

DATA VALUE: TEMPORARY (TRANSVENOUS)

-----  
ABBREVIATION.....tt  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Temporary (Transvenous)  
HELP\_TEXT.....No Help Text  
ID.....{F6A74188-3D79-D749-540F-4EA32201B0E6}  
TERM.....TEMPORARY (TRANSVENOUS)

DATA VALUE: TENACIOUS

-----  
ABBREVIATION.....TENACIOUS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Tenacious  
HELP\_TEXT.....No Help Text  
ID.....{EE90B04C-F74B-43F8-F3F2-8C2D4359D0B3}  
TERM.....TENACIOUS

DATA VALUE: TENDER

-----  
ABBREVIATION.....TENDER  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Tender  
HELP\_TEXT.....No Help Text  
ID.....{30486CC7-3097-4F95-9238-FD6C2126CD10}  
TERM.....TENDER

DATA VALUE: TENDON

-----  
ABBREVIATION.....TENDON  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Tendon  
HELP\_TEXT.....No Help Text  
ID.....{40628D33-E217-4D39-9A82-4C8D7702576F}  
TERM.....TENDON

DATA VALUE: THICK

-----  
ABBREVIATION.....THICK  
ACTIVE.....Yes  
DESCRIPTION.....Thick  
DISPLAY\_NAME.....Thick  
HELP\_TEXT.....No Help Text  
ID.....{1554A49E-1468-46B2-F851-531DDB1D1A98}  
TERM.....THICK

DATA VALUE: THIGH

-----  
ABBREVIATION.....THIGH  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Thigh  
HELP\_TEXT.....No Help Text  
ID.....{65AE788B-52F5-44D3-B9C1-49DAE46FF5AB}  
TERM.....THIGH

DATA VALUE: THIN

-----  
ABBREVIATION.....THIN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Thin  
HELP\_TEXT.....No Help Text  
ID.....{CB8D8B5C-E824-D030-1241-D8EB814671E9}  
TERM.....THIN

DATA VALUE: THIRD DEGREE AV BLOCK

-----  
ABBREVIATION.....tdab  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Third Degree Av Block  
HELP\_TEXT.....No Help Text  
ID.....{1AC7E31B-830E-2309-3AF7-8FB03003DD86}  
TERM.....THIRD DEGREE AV BLOCK

DATA VALUE: THORACIC SPINAL PRECAUTIONS

-----  
ABBREVIATION.....THOR S P  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Thoriacic Spinal Precautions  
HELP\_TEXT.....No Help Text  
ID.....{D55AEF23-7F5E-45DD-A178-F217E6AF3873}  
TERM.....THORACIC SPINAL PRECAUTIONS

DATA VALUE: THREADS

-----  
ABBREVIATION.....Threads  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Threads  
HELP\_TEXT.....  
ID.....{029BB303-B428-41D0-978B-6A8176108D7D}  
TERM.....THREADS

DATA VALUE: THROBBING

-----  
ABBREVIATION.....THROBBING  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Throbbing  
HELP\_TEXT.....No Help Text  
ID.....{3ACE2B10-1418-47F8-A4F3-0C41A699B5D9}  
TERM.....THROBBING

DATA VALUE: TINGLING

-----  
ABBREVIATION.....t  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Tingling  
HELP\_TEXT.....No Help Text  
ID.....{AE8E3101-AC3D-DEC3-0293-58F901F82EB8}  
TERM.....TINGLING

DATA VALUE: TINGLING MOUTH/LIPS

-----  
ABBREVIATION.....Ting Mou/L  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Tingling Mouth/Lips  
HELP\_TEXT.....  
ID.....{E603F979-310D-438B-80A2-1F52E4DE71C6}  
TERM.....TINGLING MOUTH/LIPS

DATA VALUE: TINNITUS

-----  
ABBREVIATION.....Tinnitus  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Tinnitus  
HELP\_TEXT.....  
ID.....{846F5DD2-F8AC-44BA-95C6-BF63A07E7E32}  
TERM.....TINNITUS

DATA VALUE: TO PAIN

-----  
ABBREVIATION.....TO PAIN  
ACTIVE.....Yes  
DESCRIPTION.....To Pain  
DISPLAY\_NAME.....To Pain  
HELP\_TEXT.....No Help Text  
ID.....{AA1F818D-DCFC-4FBC-97D3-BE86D309119E}  
TERM.....TO PAIN

DATA VALUE: TO SPEECH

-----  
ABBREVIATION.....TO SPEECH  
ACTIVE.....Yes  
DESCRIPTION.....To Speech  
DISPLAY\_NAME.....To Speech  
HELP\_TEXT.....No Help Text  
ID.....{2DDD94FA-D4D0-48B0-B4F4-F25C42E4AF96}  
TERM.....TO SPEECH

DATA VALUE: TOENAIL CARE

-----  
ABBREVIATION.....TOEN CARE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Toenail Care  
HELP\_TEXT.....No Help Text  
ID.....{6CF3ABCE-FF74-43E2-81E5-D23339017CA9}  
TERM.....TOENAIL CARE

DATA VALUE: TOILETING BY SCHEDULE

-----  
ABBREVIATION.....Toilet Sch  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Toileting by Schedule  
HELP\_TEXT.....  
ID.....{D47E6118-E5AD-46A0-8DCE-B28CDAFD6154}  
TERM.....TOILETING BY SCHEDULE

DATA VALUE: TOLERATING FLUIDS/NOURISHMENT

-----  
ABBREVIATION.....TO OR FL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Tolerating Fluids/Nourishment  
HELP\_TEXT.....No Help Text  
ID.....{E107F272-CD27-46B8-8C50-F9FF776DDCD9}  
TERM.....TOLERATING FLUIDS/NOURISHMENT

DATA VALUE: TONIC/CLONIC

-----  
ABBREVIATION.....TON/CLON  
ACTIVE.....Yes  
DESCRIPTION.....Tonic/Clonic  
DISPLAY\_NAME.....Tonic/Clonic  
HELP\_TEXT.....Tonic/Clonic  
ID.....{81932BAA-6C79-47A3-85B3-592824B6C09C}  
TERM.....TONIC/CLONIC

DATA VALUE: TOPICAL TREATMENT

-----  
ABBREVIATION.....TOPICAL T  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Topical Treatment  
HELP\_TEXT.....No Help Text  
ID.....{4E995342-8A6F-4B21-8013-F8C6A9FC7658}  
TERM.....TOPICAL TREATMENT

DATA VALUE: TORSADES DE POINTES

-----  
ABBREVIATION.....Tor de Poi  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Torsades de Pointes  
HELP\_TEXT.....  
ID.....{E936B787-7912-4CAE-A465-AD1A1FEA68CA}  
TERM.....TORSADES DE POINTES

DATA VALUE: TOTAL DEPENDENCE

-----  
ABBREVIATION.....TOT DEPEND  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Total Dependence  
HELP\_TEXT.....No Help Text  
ID.....{38702CDD-E440-4EFA-9583-36EEE719F7C3}  
TERM.....TOTAL DEPENDENCE

DATA VALUE: TRACHEOSTOMY

-----  
ABBREVIATION.....TRACHEOSTO  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Tracheostomy  
HELP\_TEXT.....No Help Text  
ID.....{2B0B7599-847E-852F-D16F-091CAAA4B205}  
TERM.....TRACHEOSTOMY

DATA VALUE: TRANSCUTANEOUS

-----  
ABBREVIATION.....t  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Transcutaneous  
HELP\_TEXT.....No Help Text  
ID.....{F4315F92-8804-1CB4-4D03-CCFDC95C6E71}  
TERM.....TRANSCUTANEOUS

DATA VALUE: TRANSFER

-----  
ABBREVIATION.....TRANSFER  
ACTIVE.....Yes  
DESCRIPTION.....Transfer  
DISPLAY\_NAME.....Transfer  
HELP\_TEXT.....No Help Text  
ID.....{DFDB4783-25BE-485B-A9BA-621F1004AAF2}  
TERM.....TRANSFER

DATA VALUE: TRANSFER WITH ASSIST

-----  
ABBREVIATION.....TRANS ASST  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Transfer With Assist  
HELP\_TEXT.....No Help Text  
ID.....{62A5BFF5-9C7C-4412-83A8-F872095DE005}  
TERM.....TRANSFER WITH ASSIST

DATA VALUE: TRANSPARENT

-----  
ABBREVIATION.....TRANSPAREN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Transparent  
HELP\_TEXT.....No Help Text  
ID.....{0E73D24D-0EFC-08F9-196A-D30DD1BB3228}  
TERM.....TRANSPARENT

DATA VALUE: TRANSPARENT DRESSING

-----  
ABBREVIATION.....TRANSP DRE  
ACTIVE.....Yes  
DESCRIPTION.....Transparent Dressing  
DISPLAY\_NAME.....Transparent Dressing  
HELP\_TEXT.....No Help Text  
ID.....{777C7722-AEA8-4A0B-9ADD-1157D37AE539}  
TERM.....TRANSPARENT DRESSING

DATA VALUE: TRANSPORT FACILITATOR PRESENT

-----  
ABBREVIATION.....TR FA PR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Transport Facilitator Present  
HELP\_TEXT.....No Help Text  
ID.....{645BF269-6B23-4377-95ED-71291F075B25}  
TERM.....TRANSPORT FACILITATOR PRESENT

DATA VALUE: TRANSPYLORIC

-----  
ABBREVIATION.....TRANSPYOLO  
ACTIVE.....Yes  
DESCRIPTION.....Transpyloric  
DISPLAY\_NAME.....Transpyloric  
HELP\_TEXT.....No Help Text  
ID.....{78FA5E0A-EA21-40F3-8E8D-3B3FCC1ED060}  
TERM.....TRANSPYLORIC

DATA VALUE: TREATMENTS

-----  
ABBREVIATION.....TREATMENTS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Treatments  
HELP\_TEXT.....No Help Text  
ID.....{7A6879B5-F781-498E-9D71-3FA0243C3F21}  
TERM.....TREATMENTS



DATA VALUE: TRIGEMINY

-----  
ABBREVIATION.....Trigeminy  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Trigeminy  
HELP\_TEXT.....  
ID.....{13D65F9F-B51C-4CEC-8002-037C8323626C}  
TERM.....TRIGEMINY

DATA VALUE: TRIPLE

-----  
ABBREVIATION.....TRIPLE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Triple  
HELP\_TEXT.....No Help Text  
ID.....{077BE8CF-F2AB-56E1-1F25-51DF80F8054D}  
TERM.....TRIPLE

DATA VALUE: TUBED/DYSPHAGIC

-----  
ABBREVIATION.....TUBE/DYSPH  
ACTIVE.....Yes  
DESCRIPTION.....Tubed/Dysphagic  
DISPLAY\_NAME.....Tubed/Dysphagic  
HELP\_TEXT.....No Help Text  
ID.....{82AAFAE3-50B0-4F2A-91CF-7900849A14EB}  
TERM.....TUBED/DYSPHAGIC

DATA VALUE: TUBING CHANGE

-----  
ABBREVIATION.....Tube Chg  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Tubing Change  
HELP\_TEXT.....  
ID.....{6418C2A2-D244-4C5B-A398-904B99AACEA4}  
TERM.....TUBING CHANGE

DATA VALUE: TUBING CHANGED/LABELED

-----  
ABBREVIATION.....TUBE CHANG  
ACTIVE.....Yes  
DESCRIPTION.....Changed or labeled tubing  
DISPLAY\_NAME.....Tubing Changed / Labeled  
HELP\_TEXT.....No Help Text  
ID.....{19BA337C-6E28-411E-9702-48B262C593FA}  
TERM.....TUBING CHANGED/LABELED

DATA VALUE: TURN AND POSITION

-----  
ABBREVIATION.....TURN POSIT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Turn and Position  
HELP\_TEXT.....No Help Text  
ID.....{44A4F470-1760-4073-AB3A-D59CA9EFE7F3}  
TERM.....TURN AND POSITION

DATA VALUE: TURNING SELF INDEPENDENTLY

-----  
ABBREVIATION.....TURN SELF  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Turning Self Independently  
HELP\_TEXT.....No Help Text  
ID.....{DD804CEB-F7CA-4D42-BF6A-447E38B55E98}  
TERM.....TURNING SELF INDEPENDENTLY

DATA VALUE: TWILL TAPE

-----  
ABBREVIATION.....Twill Tape  
ACTIVE.....Yes  
DESCRIPTION.....Twill Tape  
DISPLAY\_NAME.....Twill Tape  
HELP\_TEXT.....No Help Text  
ID.....{C52BE250-732D-03BB-B8F9-72BE4D546807}  
TERM.....TWILL TAPE

DATA VALUE: TWO CAL HN

-----  
ABBREVIATION.....Two Cal HN  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Two Cal HN  
HELP\_TEXT.....  
ID.....{3CAF568F-4844-4B3F-97A2-86EDE50BB2DE}  
TERM.....TWO CAL HN

DATA VALUE: TWO PERSON ASSIST

-----  
ABBREVIATION.....2 PERS AST  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Two Person Assist  
HELP\_TEXT.....No Help Text  
ID.....{0547070B-C9BC-4406-BDC0-BF9C6C23D41A}  
TERM.....TWO PERSON ASSIST

DATA VALUE: ULTRACAL

-----  
ABBREVIATION.....ULTRACAL  
ACTIVE.....Yes  
DESCRIPTION.....Enteral Nutrition Type  
DISPLAY\_NAME.....Ultracal  
HELP\_TEXT.....No Help Text  
ID.....{179AE9A5-154D-4CD2-B70B-66E84F5D0932}  
TERM.....ULTRACAL

DATA VALUE: UNABLE TO AROUSE

-----  
ABBREVIATION.....Un to Ar  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Unable to Arouse  
HELP\_TEXT.....  
ID.....{05B32DD9-E634-4B60-A58B-2E915EA2834D}  
TERM.....UNABLE TO AROUSE

DATA VALUE: UNABLE TO ASSESS

-----  
ABBREVIATION.....UNA TO AS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Unable To Assess  
HELP\_TEXT.....No Help Text  
ID.....{4D11A5E0-1CE8-4390-9EDD-ECEACCAA991B}  
TERM.....UNABLE TO ASSESS

DATA VALUE: UNABLE TO INFORM

-----  
ABBREVIATION.....UNABLE INF  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Unable to Inform  
HELP\_TEXT.....No Help Text  
ID.....{77428C87-B3F9-4BBD-B781-F7D89BCF9492}  
TERM.....UNABLE TO INFORM

DATA VALUE: UNABLE TO READ

-----  
ABBREVIATION.....ILLITERATE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Unable to Read  
HELP\_TEXT.....No Help Text  
ID.....{6D0C2C48-BDD2-489B-94C8-49E6DAC93F2E}  
TERM.....UNABLE TO READ

DATA VALUE: UNABLE TO RESPOND

-----  
ABBREVIATION.....UNA TO RES  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Unable To Respond  
HELP\_TEXT.....No Help Text  
ID.....{1EF29FF5-36AB-4FC0-B5E7-F8C3FC95BA5A}  
TERM.....UNABLE TO RESPOND

DATA VALUE: UNABLE TO VOID

-----  
ABBREVIATION.....Un to Void  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Unable to Void  
HELP\_TEXT.....  
ID.....{1F2A844D-4E14-4759-B8F7-82B58E157720}  
TERM.....UNABLE TO VOID

DATA VALUE: UNABLE TO WEDGE

-----  
ABBREVIATION.....UNABLE WED  
ACTIVE.....Yes  
DESCRIPTION.....Unable to Wedge  
DISPLAY\_NAME.....Unable to Wedge  
HELP\_TEXT.....No Help Text  
ID.....{5D35F546-F1EA-49C6-8B1D-22D24526AF8B}  
TERM.....UNABLE TO WEDGE

DATA VALUE: UNBEARABLE

-----  
ABBREVIATION.....UNBEARABLE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Unbearable  
HELP\_TEXT.....No Help Text  
ID.....{96611993-45F3-4D1D-BF7F-5E5C46E5569E}  
TERM.....UNBEARABLE

DATA VALUE: UNCOOPERATIVE

-----  
ABBREVIATION.....u  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Uncooperative  
HELP\_TEXT.....No Help Text  
ID.....{2D737EFB-F4CA-7AEA-9EB2-6FF133742815}  
TERM.....UNCOOPERATIVE

DATA VALUE: UNDER STRETCHER/BED

-----  
ABBREVIATION.....Stretcher  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Under Stretcher/Bed  
HELP\_TEXT.....  
ID.....{0B945D0B-FE87-46A5-A2A9-A3FC8CC578AC}  
TERM.....UNDER STRETCHER/BED

DATA VALUE: UNEQUAL

-----  
ABBREVIATION.....Unequal  
ACTIVE.....Yes  
DESCRIPTION.....Not the same  
DISPLAY\_NAME.....Unequal  
HELP\_TEXT.....Not the same  
ID.....{73319200-164F-B39F-649C-B7FB1AFBDBB8}  
TERM.....UNEQUAL

DATA VALUE: UNEXPECTED NUMBNESS

-----  
ABBREVIATION.....Unexp Numb  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Unexpected Numbness  
HELP\_TEXT.....  
ID.....{6418B135-56AF-4A76-A207-55D7AE0C5364}  
TERM.....UNEXPECTED NUMBNESS

DATA VALUE: UNIFOCAL PVC

-----  
ABBREVIATION.....Uni PVC  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Unifocal PVC  
HELP\_TEXT.....  
ID.....{5FE7D7FD-048B-4948-A4DC-C9F453AB31C5}  
TERM.....UNIFOCAL PVC

DATA VALUE: UNKNOWN

-----  
ABBREVIATION.....u  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Unknown  
HELP\_TEXT.....No Help Text  
ID.....{1F872397-8635-18A3-E252-AE9B20B88C7D}  
TERM.....UNKNOWN

DATA VALUE: UNSTABLE

-----  
ABBREVIATION.....UNSTABLE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Unstable  
HELP\_TEXT.....No Help Text  
ID.....{5841C521-8070-4034-A69D-27DAB2990EB3}  
TERM.....UNSTABLE

DATA VALUE: UNSTAGEABLE

-----  
ABBREVIATION.....UNSTAGEABL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Unstageable  
HELP\_TEXT.....No Help Text  
ID.....{5C507F64-2EC0-4FD8-A786-94D76D90C0D1}  
TERM.....UNSTAGEABLE

DATA VALUE: UNSTEADY

-----  
ABBREVIATION.....Unsteady  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Unsteady  
HELP\_TEXT.....  
ID.....{7AC26DBA-B7F1-401C-9D4C-5C0F06C87152}  
TERM.....UNSTEADY

DATA VALUE: UNSUCCESSFUL

-----  
ABBREVIATION.....UNSUCCESSF  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Unsuccessful  
HELP\_TEXT.....No Help Text  
ID.....{54DC58A1-7585-4A32-B810-90B586ED9332}  
TERM.....UNSUCCESSFUL

DATA VALUE: UP

-----  
ABBREVIATION.....u  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Up  
HELP\_TEXT.....No Help Text  
ID.....{7BCC74F9-FE2F-FE60-57BA-321BFA316BEC}  
TERM.....UP

DATA VALUE: UP TO CHAIR

-----  
ABBREVIATION.....UP TO CHAI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Up To Chair  
HELP\_TEXT.....No Help Text  
ID.....{42E598E5-82EC-428B-B5FA-1129E3095C78}  
TERM.....UP TO CHAIR

DATA VALUE: URINARY RETENTION

-----  
ABBREVIATION.....Urinary Re  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Urinary Retention  
HELP\_TEXT.....  
ID.....{959EBBDC-7F67-42EC-8CF0-295F709BB433}  
TERM.....URINARY RETENTION

DATA VALUE: UROMETER

-----  
ABBREVIATION.....u  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Urometer  
HELP\_TEXT.....No Help Text  
ID.....{13074A3A-19D7-300C-7BFF-8DC18DDE2951}  
TERM.....UROMETER

DATA VALUE: UROSTOMY

-----  
ABBREVIATION.....UROSTOMY  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Urostomy  
HELP\_TEXT.....No Help Text  
ID.....{1F9E528D-6D78-54EE-94E6-649C6F2056A3}  
TERM.....UROSTOMY

DATA VALUE: V1

-----  
ABBREVIATION.....V1  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....V1  
HELP\_TEXT.....No Help Text  
ID.....{4CAA3A92-F3DB-4847-8523-380EFEE6CFD5}  
TERM.....V1

DATA VALUE: V2

-----  
ABBREVIATION.....V2  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....V2  
HELP\_TEXT.....No Help Text  
ID.....{05A3232E-A30D-4771-A1B5-E7705B717646}  
TERM.....V2

DATA VALUE: V3

-----  
ABBREVIATION.....V3  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....V3  
HELP\_TEXT.....No Help Text  
ID.....{E150126E-4AA7-4F96-961E-FD1169C77158}  
TERM.....V3

DATA VALUE: V4

-----  
ABBREVIATION.....V4  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....V4  
HELP\_TEXT.....No Help Text  
ID.....{255260F9-5CC7-41B1-831C-F9319FA27506}  
TERM.....V4

DATA VALUE: V5

-----  
ABBREVIATION.....V5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....V5  
HELP\_TEXT.....No Help Text  
ID.....{D5C10744-0EB5-46FA-AE6C-CCF571318F8A}  
TERM.....V5

DATA VALUE: V6

-----  
ABBREVIATION.....V6  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....V6  
HELP\_TEXT.....No Help Text  
ID.....{B86F72D7-7884-4470-92A8-EF9277C5B584}  
TERM.....V6



DATA VALUE: VACUUM ASSISTED CLOSURE

-----  
ABBREVIATION.....VAC ASS CL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Vacuum Assisted Closure  
HELP\_TEXT.....No Help Text  
ID.....{9C3DD561-FECF-47FA-BBA5-4C08FF9A4CA1}  
TERM.....VACUUM ASSISTED CLOSURE

DATA VALUE: VASCULAR

-----  
ABBREVIATION.....VASCULAR  
ACTIVE.....Yes  
DESCRIPTION.....Vascular  
DISPLAY\_NAME.....Vascular  
HELP\_TEXT.....No Help Text  
ID.....{60523BF8-6372-4760-881C-A62F85D2A0F6}  
TERM.....VASCULAR

DATA VALUE: VASCULAR LESION ARTERIAL

-----  
ABBREVIATION.....VASLES ART  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Vascular Leasion Arterial  
HELP\_TEXT.....No Help Text  
ID.....{EA565D3F-D984-4553-B597-DFE844554CE9}  
TERM.....VASCULAR LESION ARTERIAL

DATA VALUE: VASCULAR LESION MIXED

-----  
ABBREVIATION.....VASLES MIX  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Vascular Leasion Mixed  
HELP\_TEXT.....No Help Text  
ID.....{8684AB2D-9162-4DC2-A5FA-A14558678689}  
TERM.....VASCULAR LESION MIXED

DATA VALUE: VASCULAR LESION VENOUS

-----  
ABBREVIATION.....VASLES VEN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Vascular Lesion Venous  
HELP\_TEXT.....No Help Text  
ID.....{3F8CEC69-5E13-4F26-AB86-34FE2805E3C6}  
TERM.....VASCULAR LESION VENOUS

DATA VALUE: VASCULAR/NEUROPATHIC

-----  
ABBREVIATION.....VASC/NEURO  
ACTIVE.....Yes  
DESCRIPTION.....Vascular / Neuropathic  
DISPLAY\_NAME.....Vascular / Neuropathic  
HELP\_TEXT.....No Help Text  
ID.....{6B91CEED-4863-46EA-A9AE-3B0D522F1BDC}  
TERM.....VASCULAR/NEUROPATHIC

DATA VALUE: VASOPRESSOR

-----  
ABBREVIATION.....VASOPRESS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Vasopressor  
HELP\_TEXT.....No Help Text  
ID.....{591FCB09-1287-48C6-8007-4CE961296BA9}  
TERM.....VASOPRESSOR

DATA VALUE: VDD

-----  
ABBREVIATION.....VDD  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....VDD  
HELP\_TEXT.....No Help Text  
ID.....{E733E525-EF36-2F93-AECC-BAD1C3E2C633}  
TERM.....VDD

DATA VALUE: VEGETARIAN DIET

-----  
ABBREVIATION.....VEG DIET  
ACTIVE.....Yes  
DESCRIPTION.....Vegetarian Diet  
DISPLAY\_NAME.....Vegetarian Diet  
HELP\_TEXT.....No Help Text  
ID.....{6B8143BE-843F-4944-9ABB-7175AFE2DF17}  
TERM.....VEGETARIAN DIET

DATA VALUE: VENDOR REP

-----  
ABBREVIATION.....Vendor Rep  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Vendor Rep  
HELP\_TEXT.....  
ID.....{4D90C5A3-7395-4BDF-9B5B-BC85F1BE28F7}  
TERM.....VENDOR REP

DATA VALUE: VENTICULAR FIBRILLATION

-----  
ABBREVIATION.....vf  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Ventricular Fibrillation  
HELP\_TEXT.....No Help Text  
ID.....{400B70C2-6E15-458F-70C6-FA2597DE2556}  
TERM.....VENTICULAR FIBRILLATION

DATA VALUE: VENTILATOR

-----  
ABBREVIATION.....v  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Ventilator  
HELP\_TEXT.....No Help Text  
ID.....{A5A647ED-A56D-5012-AD5B-788B3DCB3498}  
TERM.....VENTILATOR

DATA VALUE: VENTILATOR ACHIEVA P5O2

-----  
ABBREVIATION.....ACHIEVA P5  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Achieva P5O2  
HELP\_TEXT.....No Help Text  
ID.....{1690EEA5-3C3B-4D30-AC5C-214C94626699}  
TERM.....VENTILATOR ACHIEVA P5O2

DATA VALUE: VENTILATOR ADULT STAR

-----  
ABBREVIATION.....ADULT STAR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Adult Star Ventilator  
HELP\_TEXT.....No Help Text  
ID.....{879B13DA-4874-42D4-99EC-8E4ADB31B14E}  
TERM.....VENTILATOR ADULT STAR

DATA VALUE: VENTILATOR BEAR 1000

-----  
ABBREVIATION.....BEAR 1000  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bear 1000  
HELP\_TEXT.....No Help Text  
ID.....{03C30F99-A27C-427C-B0ED-872F0DE646BD}  
TERM.....VENTILATOR BEAR 1000

DATA VALUE: VENTILATOR DRAEGER

-----  
ABBREVIATION.....DRAEGER  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Draeger  
HELP\_TEXT.....No Help Text  
ID.....{DD18C42F-EF13-460B-8676-9ED1E429B54B}  
TERM.....VENTILATOR DRAEGER

DATA VALUE: VENTILATOR OTHER

-----  
ABBREVIATION.....VENT OTHER  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Other Ventilator; Use Comments  
HELP\_TEXT.....No Help Text  
ID.....{138B7B3D-9F8B-447A-9785-12AB396AE6DE}  
TERM.....VENTILATOR OTHER

DATA VALUE: VENTILATOR PB 7200ae

-----  
ABBREVIATION.....PB7200ae  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....PB 7200ae  
HELP\_TEXT.....No Help Text  
ID.....{86AE485F-ECD1-4B8E-820F-71CF159F318B}  
TERM.....VENTILATOR PB 7200ae

DATA VALUE: VENTILATOR PB 840

-----  
ABBREVIATION.....PB 840  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....PB 840  
HELP\_TEXT.....No Help Text  
ID.....{25A7B0DA-E437-4201-B8F3-3663B1B91C1A}  
TERM.....VENTILATOR PB 840

DATA VALUE: VENTILATOR SERVO

-----  
ABBREVIATION.....SERVO  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Servo  
HELP\_TEXT.....No Help Text  
ID.....{D6F06715-1ADD-4941-B026-B9A5E463943E}  
TERM.....VENTILATOR SERVO

DATA VALUE: VENTRICULAR

-----  
ABBREVIATION.....v  
ACTIVE.....Yes  
DESCRIPTION.....VENTRICULAR  
DISPLAY\_NAME.....Ventricular  
HELP\_TEXT.....No Help Text  
ID.....{307A4572-3CF0-7875-7313-7FD4AB23F8E7}  
TERM.....VENTRICULAR

DATA VALUE: VENTRICULAR TACHYCARDIA

-----  
ABBREVIATION.....vt  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Ventricular Tachycardia  
HELP\_TEXT.....No Help Text  
ID.....{138258AD-3A2C-2A6C-9DF8-615CCC303DC2}  
TERM.....VENTRICULAR TACHYCARDIA

DATA VALUE: VENTRICULAR WHIP PRESENT

-----  
ABBREVIATION.....VENT WHIP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Ventricular Whip Present  
HELP\_TEXT.....No Help Text  
ID.....{15CB9645-690D-446B-A380-F14CEEC9C6AA}  
TERM.....VENTRICULAR WHIP PRESENT

DATA VALUE: VENTRICULOSTOMY

-----  
ABBREVIATION.....Ventriculo  
ACTIVE.....Yes  
DESCRIPTION.....Ventriculostomy  
DISPLAY\_NAME.....Ventriculostomy  
HELP\_TEXT.....No Help Text  
ID.....{ECF20317-4708-7661-D4C2-9BDB4CAAC4FA}  
TERM.....VENTRICULOSTOMY

DATA VALUE: VERBAL DE-ESCALATION

-----  
ABBREVIATION.....VERB DE-ES  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Verbal De-escalation  
HELP\_TEXT.....No Help Text  
ID.....{8AB46F8C-B4D5-4800-BCF7-52CA3F32E99B}  
TERM.....VERBAL DE-ESCALATION

DATA VALUE: VERBALIZED BY PATIENT

-----  
ABBREVIATION.....Verbal  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Verbalized by Patient  
HELP\_TEXT.....  
ID.....{6AF901B8-A19D-4C89-869F-2ECBF2D105A7}  
TERM.....VERBALIZED BY PATIENT

DATA VALUE: VERY LOW BMI

-----  
ABBREVIATION.....LOW BMI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Very Low BMI  
HELP\_TEXT.....No Help Text  
ID.....{78B3A37B-938B-42E5-A403-75E206722BC3}  
TERM.....VERY LOW BMI

DATA VALUE: VIOLET

-----  
ABBREVIATION.....VIOLET  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Violet  
HELP\_TEXT.....No Help Text  
ID.....{1249AE97-CF1F-4604-A17A-5FC85F6D0185}  
TERM.....VIOLET

DATA VALUE: VIP PORT

-----  
ABBREVIATION.....VIP PORT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....VIP Port Catheter  
HELP\_TEXT.....No Help Text  
ID.....{3D0426B3-F8F3-BAA4-71C7-8A09F911EF08}  
TERM.....VIP PORT

DATA VALUE: VISCOUS

-----  
ABBREVIATION.....VISCOUS  
ACTIVE.....Yes  
DESCRIPTION.....Viscous  
DISPLAY\_NAME.....Viscous  
HELP\_TEXT.....No Help Text  
ID.....{0E8019D5-F2DB-46B8-B750-41568DD0179C}  
TERM.....VISCOUS

DATA VALUE: VISUAL IMPAIRMENT

-----  
ABBREVIATION.....VIS IMPAIR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Visual Impairment  
HELP\_TEXT.....No Help Text  
ID.....{7B2107CE-002D-43BE-9102-FCD5B9FAD617}  
TERM.....VISUAL IMPAIRMENT

DATA VALUE: VITAL 1.0

-----  
ABBREVIATION.....Vital 1.0  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Vital 1.0  
HELP\_TEXT.....  
ID.....{D3CA653C-419D-4186-8704-E171650062DF}  
TERM.....VITAL 1.0

DATA VALUE: VITAL 1.5

-----  
ABBREVIATION.....Vital 1.5  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Vital 1.5  
HELP\_TEXT.....  
ID.....{3665F539-7BBD-4AD9-A4D0-0A33506036D7}  
TERM.....VITAL 1.5

DATA VALUE: VITAL AF

-----  
ABBREVIATION.....Vital AF  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Vital AF  
HELP\_TEXT.....  
ID.....{D7FCB7B2-EA11-4A75-AD23-1BEDEE03DC58}  
TERM.....VITAL AF

DATA VALUE: VITAL ENTERAL NUTRITION TYPE

-----  
ABBREVIATION.....VITAL ENT  
ACTIVE.....Yes  
DESCRIPTION.....Vital Enteral Nutrition Type  
DISPLAY\_NAME.....Vital  
HELP\_TEXT.....No Help Text  
ID.....{2DA036DA-BF27-4E52-86A2-20A312823480}  
TERM.....VITAL ENTERAL NUTRITION TYPE

DATA VALUE: VITAL HN ENTERAL NUTRITION

-----  
ABBREVIATION.....VITAL HN  
ACTIVE.....Yes  
DESCRIPTION.....Vital HN Enteral Nutrition Type  
DISPLAY\_NAME.....Vital HN  
HELP\_TEXT.....No Help Text  
ID.....{E45849D4-05EA-4DC1-BD9C-3FF66DD4D292}  
TERM.....VITAL HN ENTERAL NUTRITION

DATA VALUE: VITAL SIGNS STABLE

-----  
ABBREVIATION.....VI SI ST  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Vital Signs Stable  
HELP\_TEXT.....No Help Text  
ID.....{E6ECD067-54A0-4214-A5F8-43D08C9FEAD6}  
TERM.....VITAL SIGNS STABLE

DATA VALUE: VIVONEX RTF

-----  
ABBREVIATION.....Viv RTF  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Vivonex RTF  
HELP\_TEXT.....  
ID.....{461E7430-CC94-47D1-A90D-9694314C7729}  
TERM.....VIVONEX RTF

DATA VALUE: VOIDED

-----  
ABBREVIATION.....VOIDED  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Voided  
HELP\_TEXT.....No Help Text  
ID.....{7F23E601-A1E5-4016-8A12-06B1099041CF}  
TERM.....VOIDED

DATA VALUE: VOLUNTEER

-----  
ABBREVIATION.....VOLUNTEE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Volunteer  
HELP\_TEXT.....No Help Text  
ID.....{F5406412-DE48-4F84-8A31-0B1D79DB576D}  
TERM.....VOLUNTEER



DATA VALUE: VOMITING

-----  
ABBREVIATION.....VOMITING  
ACTIVE.....Yes  
DESCRIPTION.....Vomiting  
DISPLAY\_NAME.....Vomiting  
HELP\_TEXT.....Vomiting  
ID.....{921F463C-1639-4D13-8E4E-6D3A1B6927F9}  
TERM.....VOMITING

DATA VALUE: VVI

-----  
ABBREVIATION.....VVI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....VVI  
HELP\_TEXT.....No Help Text  
ID.....{2C87F8F3-2B99-80E7-04EB-E6FE88168627}  
TERM.....VVI

DATA VALUE: WAFFLE

-----  
ABBREVIATION.....WAFFLE  
ACTIVE.....Yes  
DESCRIPTION.....Waffle  
DISPLAY\_NAME.....Waffle  
HELP\_TEXT.....No Help Text  
ID.....{29A0FA73-CAA7-4660-9861-79A09F6DD98C}  
TERM.....WAFFLE

DATA VALUE: WAITING ROOM

-----  
ABBREVIATION.....WAIT ROOM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Waiting Room  
HELP\_TEXT.....No Help Text  
ID.....{13476486-F895-4FFC-AB42-B274827B47DA}  
TERM.....WAITING ROOM

DATA VALUE: WALKER

-----  
ABBREVIATION.....WALKER  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Walker  
HELP\_TEXT.....No Help Text  
ID.....{91247001-5E89-46A0-A841-B4FAE9226E63}  
TERM.....WALKER

DATA VALUE: WALKING OR PACING

-----  
ABBREVIATION.....wop  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Walking Or Pacing  
HELP\_TEXT.....No Help Text  
ID.....{CE7AABCB-84B7-A310-3A65-0AF17FB7B0E7}  
TERM.....WALKING OR PACING

DATA VALUE: WANDER GUARD IN PLACE

-----  
ABBREVIATION.....Wander Grd  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Wander Guard in Place  
HELP\_TEXT.....  
ID.....{B49D735F-A6ED-4E8F-9BD0-F8514245B262}  
TERM.....WANDER GUARD IN PLACE

DATA VALUE: WANDERER

-----  
ABBREVIATION.....w  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Wanderer  
HELP\_TEXT.....No Help Text  
ID.....{C5BD737B-F9A6-5800-B45C-2585507C3EDC}  
TERM.....WANDERER

DATA VALUE: WARD

-----  
ABBREVIATION.....WARD  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Ward  
HELP\_TEXT.....No Help Text  
ID.....{A16CB5CB-3401-4F6C-8467-7F78FC929872}  
TERM.....WARD

DATA VALUE: WARM

-----  
ABBREVIATION.....WARM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Warm  
HELP\_TEXT.....No Help Text  
ID.....{E03BBEB4-5643-968F-190F-309A1293FCF9}  
TERM.....WARM

DATA VALUE: WARM HUMIDITY WITH O2

-----  
ABBREVIATION.....WM HUM O2  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Warm Humidity with Oxygen  
HELP\_TEXT.....No Help Text  
ID.....{E62E07C3-4CA2-426B-8EBA-48C406837D9C}  
TERM.....WARM HUMIDITY WITH O2

DATA VALUE: WARMER THAN GENERAL BODY TEMP

-----  
ABBREVIATION.....WARM B TEM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Warmer Than General Body Temperature  
HELP\_TEXT.....No Help Text  
ID.....{AE43F2B4-B585-43F4-9E97-D61E225A803B}  
TERM.....WARMER THAN GENERAL BODY TEMP

DATA VALUE: WARMER THAN OPPOSITE BODY PART

-----  
ABBREVIATION.....WARM OPP P  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Warmer Than Opposite Body Part  
HELP\_TEXT.....No Help Text  
ID.....{5DBF8614-5067-49C5-BD02-C1B21C9DFECC}  
TERM.....WARMER THAN OPPOSITE BODY PART

DATA VALUE: WARMING BLANKET

-----  
ABBREVIATION.....WARM BLANK  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Warming Blanket  
HELP\_TEXT.....Warming Blanket  
ID.....{23D20AE5-DE25-4456-8F77-BF0411FB19A6}  
TERM.....WARMING BLANKET

DATA VALUE: WARMING BOTTLE

-----  
ABBREVIATION.....WARM BOTTL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Warming Bottle  
HELP\_TEXT.....No Help Text  
ID.....{AD1BEFE5-4CE6-4076-8095-9BE7387CFE25}  
TERM.....WARMING BOTTLE

DATA VALUE: WARMING OFF

-----  
ABBREVIATION.....WARMING OF  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Warming Off  
HELP\_TEXT.....No Help Text  
ID.....{98FAC977-8E2E-406A-9217-A77E1DDCDF50}  
TERM.....WARMING OFF

DATA VALUE: WARMING ON

-----  
ABBREVIATION.....WARMING ON  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Warming On  
HELP\_TEXT.....No Help Text  
ID.....{6FFDE928-DE1E-4B25-BA3C-EB29E0CDEB35}  
TERM.....WARMING ON

DATA VALUE: WATERSEAL

-----  
ABBREVIATION.....WATERSEAL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Waterseal  
HELP\_TEXT.....No Help Text  
ID.....{6D91B99D-F199-09ED-78BE-B54F2CAC4E06}  
TERM.....WATERSEAL

DATA VALUE: WAVEFORM REFLECTS ATRIAL PLACEMENT

-----  
ABBREVIATION.....WAVE ATRIA  
ACTIVE.....Yes  
DESCRIPTION.....Waveform Reflects Atrial Placement  
DISPLAY\_NAME.....Waveform Reflects Atrial Placement  
HELP\_TEXT.....No Help Text  
ID.....{A1249745-41C2-4A4F-A86D-7E5E17477B06}  
TERM.....WAVEFORM REFLECTS ATRIAL PLACEMENT

DATA VALUE: WEAK

-----  
ABBREVIATION.....w  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Weak  
HELP\_TEXT.....No Help Text  
ID.....{BE2D9305-E2DA-F492-66CF-164EE58B8179}  
TERM.....WEAK

DATA VALUE: WEAKNESS

-----  
ABBREVIATION.....Weakness  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Weakness  
HELP\_TEXT.....  
ID.....{E4C6D569-7961-4CB4-8E90-A02A31F98094}  
TERM.....WEAKNESS

DATA VALUE: WEDGED

-----  
ABBREVIATION.....WEDGED  
ACTIVE.....Yes  
DESCRIPTION.....Wedged  
DISPLAY\_NAME.....Wedged  
HELP\_TEXT.....No Help Text  
ID.....{03D7F6FD-399D-4A6B-9900-59655D76EE95}  
TERM.....WEDGED

DATA VALUE: WHEELCHAIR

-----  
ABBREVIATION.....WH CHAIR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Wheelchair  
HELP\_TEXT.....No Help Text  
ID.....{EC1B82CE-4037-466F-90AC-52CD5C3B26BE}  
TERM.....WHEELCHAIR

DATA VALUE: WHEELCHAIR EXIT ALARM

-----  
ABBREVIATION.....WHEEL EXIT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Wheelchair Exit Alarm  
HELP\_TEXT.....No Help Text  
ID.....{7D06564E-ECDE-4A11-BEFE-92AB7E448AA3}  
TERM.....WHEELCHAIR EXIT ALARM

DATA VALUE: WHEEZES

-----  
ABBREVIATION.....w  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Wheezes  
HELP\_TEXT.....No Help Text  
ID.....{1B8A9729-554D-59F6-E9E8-40BDFF9C678E}  
TERM.....WHEEZES

DATA VALUE: WHITE

-----  
ABBREVIATION.....WHITE  
ACTIVE.....Yes  
DESCRIPTION.....The color white  
DISPLAY\_NAME.....White  
HELP\_TEXT.....The color white  
ID.....{AB04EF54-5246-EF3A-E1D9-02E9E24369DD}  
TERM.....WHITE

DATA VALUE: WITH PHYSICIAN

-----  
ABBREVIATION.....wp  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....With Physician  
HELP\_TEXT.....No Help Text  
ID.....{1673C53F-DA8C-2CEC-679B-B01D914FE533}  
TERM.....WITH PHYSICIAN

DATA VALUE: WITH THERAPIST

-----  
ABBREVIATION.....wt  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....With Therapist  
HELP\_TEXT.....No Help Text  
ID.....{795A1300-1D60-9D67-692D-A2AB2B190596}  
TERM.....WITH THERAPIST

DATA VALUE: WITHDRAWN

-----  
ABBREVIATION.....w  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Withdrawn  
HELP\_TEXT.....No Help Text  
ID.....{E74B0F04-0E3F-632F-9409-C4953D95B1BA}  
TERM.....WITHDRAWN

DATA VALUE: WITHDRAWS TO PAIN

-----  
ABBREVIATION.....WITHD TO P  
ACTIVE.....Yes  
DESCRIPTION.....Withdraws To Pain  
DISPLAY\_NAME.....Withdraws To Pain  
HELP\_TEXT.....No Help Text  
ID.....{013527C8-FF6B-49CD-BF45-9338F7413E81}  
TERM.....WITHDRAWS TO PAIN

DATA VALUE: WITHOUT ASSISTANCE

-----  
ABBREVIATION.....WO ASSIS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Without Assistance  
HELP\_TEXT.....No Help Text  
ID.....{41CE081D-CB0E-49AF-936A-4FD78E61681B}  
TERM.....WITHOUT ASSISTANCE

DATA VALUE: WOUND

-----  
ABBREVIATION.....WOUND  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Wound  
HELP\_TEXT.....No Help Text  
ID.....{1EC40149-051E-4750-8651-2B2D0DD2DEAC}  
TERM.....WOUND

DATA VALUE: WOUND DISCHARGE CHARACTERISTICS BLOODY

-----  
ABBREVIATION.....WOU DI BLO  
ACTIVE.....Yes  
DESCRIPTION.....Wound Discharge Characteristics Bloody  
DISPLAY\_NAME.....Bloody  
HELP\_TEXT.....Wound Discharge Characteristics Bloody  
ID.....{628DBFBC-9B77-48EA-99FE-0DF000376506}  
TERM.....WOUND DISCHARGE CHARACTERISTICS BLOODY

DATA VALUE: WOUND DISCHARGE CHARACTERISTICS CLOUDY

-----  
ABBREVIATION.....WOU DI CLO  
ACTIVE.....Yes  
DESCRIPTION.....Wound Discharge Characteristics Cloudy  
DISPLAY\_NAME.....Cloudy  
HELP\_TEXT.....Wound Discharge Characteristics Cloudy  
ID.....{159C9993-A2E1-4133-9D5F-315EF9985EC4}  
TERM.....WOUND DISCHARGE CHARACTERISTICS CLOUDY

DATA VALUE: WOUND DISCHARGE CHARACTERISTICS NONE

-----  
ABBREVIATION.....WOU DI NON  
ACTIVE.....Yes  
DESCRIPTION.....Wound Discharge Characteristics None  
DISPLAY\_NAME.....None  
HELP\_TEXT.....Wound Discharge Characteristics None  
ID.....{929F2FA5-923F-4D88-8EDF-E0EBEE9CC7A8}  
TERM.....WOUND DISCHARGE CHARACTERISTICS NONE

DATA VALUE: WOUND DISCHARGE CHARACTERISTICS PURULENT

-----  
ABBREVIATION.....WOU DI PUR  
ACTIVE.....Yes  
DESCRIPTION.....Wound Discharge Characteristics  
Purulent  
DISPLAY\_NAME.....Purulent  
HELP\_TEXT.....Wound Discharge Characteristics  
Purulent  
ID.....{C1E10410-B6D8-4429-BD94-0B179992EAAB}  
TERM.....WOUND DISCHARGE CHARACTERISTICS  
PURULENT

DATA VALUE: WOUND DISCHARGE CHARACTERISTICS SEROSANGUINOUS

-----  
ABBREVIATION.....WOU DI SSA  
ACTIVE.....Yes  
DESCRIPTION.....Wound Discharge Characteristics  
Serosanguinous  
DISPLAY\_NAME.....Serosanguinous  
HELP\_TEXT.....Wound Discharge Characteristics  
Serosanguinous  
ID.....{B431D8EF-5B53-4DCF-A390-431EB5ED3704}  
TERM.....WOUND DISCHARGE CHARACTERISTICS  
SEROSANGUINOUS

DATA VALUE: WOUND DISCHARGE CHARACTERISTICS SEROUS

-----  
ABBREVIATION.....WOU DI SER  
ACTIVE.....Yes  
DESCRIPTION.....Wound Discharge Characteristics Serous  
DISPLAY\_NAME.....Serous  
HELP\_TEXT.....Wound Discharge Characteristics Serous  
ID.....{2124005E-5C36-48AC-AA6B-03E06510D70A}  
TERM.....WOUND DISCHARGE CHARACTERISTICS SEROUS

DATA VALUE: WOUND DISCHARGE CHARACTERISTICS VISCOUS

-----  
ABBREVIATION.....WOU DI VIS  
ACTIVE.....Yes  
DESCRIPTION.....Wound Discharge Characteristics Viscous  
DISPLAY\_NAME.....Viscous  
HELP\_TEXT.....Wound Discharge Characteristics Viscous  
ID.....{382FA6D9-1CB2-4ECC-A458-8598ABA708F8}  
TERM.....WOUND DISCHARGE CHARACTERISTICS VISCOUS

DATA VALUE: WOUND DRAINAGE AMOUNT HEAVY

-----  
ABBREVIATION.....W DR A HEA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Heavy  
HELP\_TEXT.....No Help-Text  
ID.....{BF959172-6672-4FEB-8FB7-0568405F3C83}  
TERM.....WOUND DRAINAGE AMOUNT HEAVY



DATA VALUE: WOUND DRAINAGE AMOUNT MODERATE

-----  
ABBREVIATION.....W DR A MOD  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Moderate  
HELP\_TEXT.....No Help-Text  
ID.....{E5DBFD9C-A2AC-4009-98CB-7246F6AA702D}  
TERM.....WOUND DRAINAGE AMOUNT MODERATE

DATA VALUE: WOUND DRAINAGE AMOUNT NONE

-----  
ABBREVIATION.....WOU DR A N  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....None  
HELP\_TEXT.....No Help-Text  
ID.....{92910398-311F-D4D7-8BD5-AB66C3DC5E72}  
TERM.....WOUND DRAINAGE AMOUNT NONE

DATA VALUE: WOUND DRAINAGE AMOUNT SCANT

-----  
ABBREVIATION.....W DR SCANT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Scant  
HELP\_TEXT.....No Help-Text  
ID.....{9C3CD0BD-F960-4098-B14C-17B20703B242}  
TERM.....WOUND DRAINAGE AMOUNT SCANT

DATA VALUE: WOUND DRAINAGE AMOUNT SMALL

-----  
ABBREVIATION.....W DR A SMA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Small  
HELP\_TEXT.....No Help-Text  
ID.....{3D8EDF1A-EC4B-41EF-A194-73FF60074A57}  
TERM.....WOUND DRAINAGE AMOUNT SMALL

DATA VALUE: WOUND DRESSING TYPE ANTIMICROBIAL

-----  
ABBREVIATION.....WOU DR ANT  
ACTIVE.....No  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Antimicrobial  
HELP\_TEXT.....No Help Text  
ID.....{9DAE650A-3D5C-498B-BC74-F3CE77CD2BA9}  
TERM.....WOUND DRESSING TYPE ANTIMICROBIAL

DATA VALUE: WOUND DRESSING TYPE CALCIUM ALGINATE

-----  
ABBREVIATION.....WOU DR CAL  
ACTIVE.....No  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Calcium Alginate  
HELP\_TEXT.....No Help Text  
ID.....{AA084ED8-07D0-4D8A-84D1-51F025246F40}  
TERM.....WOUND DRESSING TYPE CALCIUM ALGINATE

DATA VALUE: WOUND DRESSING TYPE GAUZE

-----  
ABBREVIATION.....WOU DR GAU  
ACTIVE.....No  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Gauze  
HELP\_TEXT.....No Help Text  
ID.....{A2EE67E2-B006-413F-8F8A-E02C21A39FC1}  
TERM.....WOUND DRESSING TYPE GAUZE

DATA VALUE: WOUND DRESSING TYPE GAUZE PACKING

-----  
ABBREVIATION.....WOU DR GPA  
ACTIVE.....No  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Gauze Packing  
HELP\_TEXT.....No Help Text  
ID.....{ED06B435-B205-4BD2-B3D1-795C40EE016F}  
TERM.....WOUND DRESSING TYPE GAUZE PACKING

DATA VALUE: WOUND DRESSING TYPE GAUZE WRAP

-----  
ABBREVIATION.....WOU DR GWR  
ACTIVE.....No  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Gauze Wrap  
HELP\_TEXT.....No Help Text  
ID.....{4257C446-7627-4632-8507-B5880AD934D5}  
TERM.....WOUND DRESSING TYPE GAUZE WRAP

DATA VALUE: WOUND DRESSING TYPE HYDROCOLLOIDAL

-----  
ABBREVIATION.....WOU DR HYD  
ACTIVE.....No  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Hydrocolloidal  
HELP\_TEXT.....No Help Text  
ID.....{8C20226C-FC6A-49A0-9A72-06876DB1F203}  
TERM.....WOUND DRESSING TYPE HYDROCOLLOIDAL

DATA VALUE: WOUND DRESSING TYPE HYDROGEL

-----  
ABBREVIATION.....WOU DR HYD  
ACTIVE.....No  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Hydrogel  
HELP\_TEXT.....No Help Text  
ID.....{60DDE603-9FF1-44FA-BC4D-45AA8753E254}  
TERM.....WOUND DRESSING TYPE HYDROGEL

DATA VALUE: WOUND DRESSING TYPE NONE

-----  
ABBREVIATION.....WOU DR NON  
ACTIVE.....No  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....None  
HELP\_TEXT.....No Help Text  
ID.....{8A58A493-A02C-4570-9786-9F5AF8D02360}  
TERM.....WOUND DRESSING TYPE NONE

DATA VALUE: WOUND DRESSING TYPE PRESSURE

-----  
ABBREVIATION.....WOU DR PRE  
ACTIVE.....No  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Pressure  
HELP\_TEXT.....No Help Text  
ID.....{520952E4-FB31-4C1C-9FE2-A10E0F1CEE00}  
TERM.....WOUND DRESSING TYPE PRESSURE

DATA VALUE: WOUND DRESSING TYPE TRANSPARENT

-----  
ABBREVIATION.....WOU DR TRA  
ACTIVE.....No  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Transparent  
HELP\_TEXT.....No Help Text  
ID.....{EA57D935-33D8-4F30-8511-706870A27122}  
TERM.....WOUND DRESSING TYPE TRANSPARENT

DATA VALUE: WOUND EDGE INDISTINCT/ATTACHED

-----  
ABBREVIATION.....WOU ED IND  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Indistinct, Attached  
HELP\_TEXT.....No Help-Text  
ID.....{C899235C-4A07-4420-92BC-8F818B8E0A7F}  
TERM.....WOUND EDGE INDISTINCT/ATTACHED

DATA VALUE: WOUND EDGE NON-ATTACHED

-----  
ABBREVIATION.....WOU ED NON  
ACTIVE.....Yes  
DESCRIPTION.....Wound Edge Non Attached  
DISPLAY\_NAME.....Non-Attached  
HELP\_TEXT.....Wound Edge Non Attached  
ID.....{1D5C1317-B4FF-421C-96A0-C20C1089E672}  
TERM.....WOUND EDGE NON-ATTACHED

DATA VALUE: WOUND EDGE ROLLED

-----  
ABBREVIATION.....WOU ED ROL  
ACTIVE.....Yes  
DESCRIPTION.....Wound Edge Rolled  
DISPLAY\_NAME.....Rolled  
HELP\_TEXT.....Wound Edge Rolled  
ID.....{24D9CFEC-1ED6-4069-9715-1F9FCB45BBD3}  
TERM.....WOUND EDGE ROLLED

DATA VALUE: WOUND EDGE TUNNELING

-----  
ABBREVIATION.....WOU ED TUN  
ACTIVE.....Yes  
DESCRIPTION.....Wound Edge Tunneling  
DISPLAY\_NAME.....Tunneling  
HELP\_TEXT.....Wound Edge Tunneling  
ID.....{97425206-5261-46C1-905F-E03E2F071381}  
TERM.....WOUND EDGE TUNNELING

DATA VALUE: WOUND EDGE UNDERMINING

-----  
ABBREVIATION.....WOU ED UND  
ACTIVE.....Yes  
DESCRIPTION.....Wound Edge Undermining  
DISPLAY\_NAME.....Undermining  
HELP\_TEXT.....Wound Edge Undermining  
ID.....{AAA3BB3A-3A3B-4F39-9AC7-3BE1A82BDC02}  
TERM.....WOUND EDGE UNDERMINING

DATA VALUE: WOUND ODOR FOUL

-----  
ABBREVIATION.....WOU OD FOU  
ACTIVE.....Yes  
DESCRIPTION.....Wound Odor Foul  
DISPLAY\_NAME.....Foul  
HELP\_TEXT.....Wound Odor Foul  
ID.....{99A74112-EA71-4896-A4DD-911379E5B453}  
TERM.....WOUND ODOR FOUL

DATA VALUE: WOUND ODOR NONE

-----  
ABBREVIATION.....WOU OD NON  
ACTIVE.....Yes  
DESCRIPTION.....Wound Odor None  
DISPLAY\_NAME.....None  
HELP\_TEXT.....Wound Odor None  
ID.....{D2E388F7-B531-4D60-A2BE-830F11B9AA75}  
TERM.....WOUND ODOR NONE

DATA VALUE: WOUND VAC DEVICE

-----  
ABBREVIATION.....WOUND VAC  
ACTIVE.....Yes  
DESCRIPTION.....Wound Vac Device  
DISPLAY\_NAME.....Wound Vac Device  
HELP\_TEXT.....No Help Text  
ID.....{A03F4978-969D-DCB8-DCE2-26EB33F097B3}  
TERM.....WOUND VAC DEVICE

DATA VALUE: X-RAY

-----  
ABBREVIATION.....X-RAY  
ACTIVE.....Yes  
DESCRIPTION.....X-Ray  
DISPLAY\_NAME.....X-Ray  
HELP\_TEXT.....X-Ray  
ID.....{C82CB7F8-2456-4A08-9FD2-4E498EB1DB20}  
TERM.....X-RAY

DATA VALUE: XENADERM

-----  
ABBREVIATION.....XENADERM  
ACTIVE.....Yes  
DESCRIPTION.....Xenaderm  
DISPLAY\_NAME.....Xenaderm  
HELP\_TEXT.....No Help Text  
ID.....{F181D70E-A231-4D2B-A1E6-DDF4DBED31CD}  
TERM.....XENADERM

DATA VALUE: XEROFORM

-----  
ABBREVIATION.....XEROFORM  
ACTIVE.....Yes  
DESCRIPTION.....Xeroform  
DISPLAY\_NAME.....Xeroform  
HELP\_TEXT.....No Help Text  
ID.....{7D7C3CDC-253E-4262-B8D4-3B70DE9AEBE4}  
TERM.....XEROFORM

DATA VALUE: YELLOW

-----  
ABBREVIATION.....YELLOW  
ACTIVE.....Yes  
DESCRIPTION.....The color yellow  
DISPLAY\_NAME.....Yellow  
HELP\_TEXT.....No Help Text  
ID.....{90DC31AD-6211-9F25-4237-79F7122A2763}  
TERM.....YELLOW

DATA VALUE: YELLOW/FIBRINOUS

-----  
ABBREVIATION.....YELLOW/FIB  
ACTIVE.....Yes  
DESCRIPTION.....Yellow / Fibrinous  
DISPLAY\_NAME.....Yellow / Fibrinous  
HELP\_TEXT.....No Help Text  
ID.....{47CF52FE-0374-4690-ADA1-643E0954E035}  
TERM.....YELLOW/FIBRINOUS

DATA VALUE: YELLOW/GREEN

-----  
ABBREVIATION.....YELLOW/GRN  
ACTIVE.....Yes  
DESCRIPTION.....Yellow/Green in appearance  
DISPLAY\_NAME.....Yellow / Green  
HELP\_TEXT.....No Help Text  
ID.....{1D0E7E81-3331-4D07-9B94-3D4BF72A0F0D}  
TERM.....YELLOW/GREEN

DATA VALUE: YES

-----  
ABBREVIATION.....y  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Yes  
HELP\_TEXT.....No Help Text  
ID.....{483DE806-2EDB-34B0-FE11-C29EEEE649CD}  
TERM.....YES

DATA VALUE: pH

-----  
ABBREVIATION.....pH  
ACTIVE.....Yes  
DESCRIPTION.....pH  
DISPLAY\_NAME.....pH  
HELP\_TEXT.....pH  
ID.....{104FE39C-B0A2-40B5-A523-B11E5F2BA272}  
TERM.....pH

DATA VALUE: {8BC7F294-7126-4843-91A0-14F805EC4578}

-----  
ABBREVIATION.....Morphine  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Morphine  
HELP\_TEXT.....  
ID.....{8BC7F294-7126-4843-91A0-14F805EC4578}  
TERM.....{8BC7F294-7126-4843-91A0-14F805EC4578}

LOCATION: ABDOMEN

-----

ABBREVIATION.....ABDOMEN  
ACTIVE.....Yes  
DESCRIPTION.....ABDOMEN  
DISPLAY\_NAME.....Abdomen  
HELP\_TEXT.....No Help Text  
ID.....{109FE041-2521-6EF2-7F95-1FD1587796AB}  
TERM.....ABDOMEN

LOCATION: ABDOMEN LLQ

-----

ABBREVIATION.....ABD LLQ  
ACTIVE.....Yes  
DESCRIPTION.....Abdomen Left Lower Quarter  
DISPLAY\_NAME.....Abdomen Left Lower Quarter  
HELP\_TEXT.....No Help Text  
ID.....{39831F59-D7B0-4156-AC57-4C622A93CA58}  
TERM.....ABDOMEN LLQ

LOCATION: ABDOMEN LUQ

-----

ABBREVIATION.....ABD LUQ  
ACTIVE.....Yes  
DESCRIPTION.....Abdomen Left Upper Quarter  
DISPLAY\_NAME.....Abdomen Left Upper Quarter  
HELP\_TEXT.....No Help Text  
ID.....{3B31E1EC-CE74-4C76-BB61-10C678431FCA}  
TERM.....ABDOMEN LUQ

LOCATION: ABDOMEN RLQ

-----

ABBREVIATION.....ABD RLQ  
ACTIVE.....Yes  
DESCRIPTION.....Abdomen Right Lower Quarter  
DISPLAY\_NAME.....Abdomen Right Lower Quarter  
HELP\_TEXT.....No Help Text  
ID.....{2C48394F-6877-4738-856C-E933644D8067}  
TERM.....ABDOMEN RLQ

LOCATION: ABDOMEN RUQ

-----

ABBREVIATION.....ABD RUQ  
ACTIVE.....Yes  
DESCRIPTION.....Abdomen Right Upper Quarter  
DISPLAY\_NAME.....Abdomen Right Upper Quarter  
HELP\_TEXT.....No Help Text  
ID.....{F29C9FA1-717A-49B7-8000-95070E9AACA0}  
TERM.....ABDOMEN RUQ



LOCATION: ABDOMINAL

-----  
ABBREVIATION.....Abd  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Abdominal  
HELP\_TEXT.....No Help Text  
ID.....{AE1DA31F-2C28-98C2-501A-F3D5ADB2501F}  
TERM.....ABDOMINAL

LOCATION: ABDOMINAL AORTA

-----  
ABBREVIATION.....ABD AORTA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Abdominal Aorta  
HELP\_TEXT.....No Help Text  
ID.....{C4BD42DD-1D6B-4B2C-B6D5-19E1185E06E3}  
TERM.....ABDOMINAL AORTA

LOCATION: ANKLE

-----  
ABBREVIATION.....Ank  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Ankle  
HELP\_TEXT.....No Help Text  
ID.....{F2A2CF53-BFA5-9279-DB22-EA2627988E5B}  
TERM.....ANKLE

LOCATION: APICAL

-----  
ABBREVIATION.....Ap  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Apical  
HELP\_TEXT.....No Help Text  
ID.....{C5CF5C51-63CC-33FB-9BAE-FB1DD9900911}  
TERM.....APICAL

LOCATION: AV FISTULA

-----  
ABBREVIATION.....AVF  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....AV Fistula  
HELP\_TEXT.....No Help Text  
ID.....{F977EEF6-EA6A-D56D-D1F7-29B91510A02F}  
TERM.....AV FISTULA

LOCATION: AXILLARY

-----  
ABBREVIATION.....Ax  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Axillary  
HELP\_TEXT.....No Help Text  
ID.....{135AA5F4-A1C7-978D-E9AB-6873841AEB57}  
TERM.....AXILLARY

LOCATION: BACK

-----  
ABBREVIATION.....BACK  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Back  
HELP\_TEXT.....No Help Text  
ID.....{22C26239-C692-4E93-B522-E5E8392BA3B1}  
TERM.....BACK

LOCATION: BILATERAL ANKLES

-----  
ABBREVIATION.....BILA ANKLE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bilateral Ankles  
HELP\_TEXT.....No Help-Text  
ID.....{E56A03A5-5CE0-47BF-A14D-F6DF66F4E769}  
TERM.....BILATERAL ANKLES

LOCATION: BILATERAL FOOT

-----  
ABBREVIATION.....BILAT FOOT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bilateral  
HELP\_TEXT.....No Help Text  
ID.....{C758847D-468F-479A-B12D-F4E46F0E5E3A}  
TERM.....BILATERAL FOOT

LOCATION: BILATERAL HAND

-----  
ABBREVIATION.....BILAT HAND  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bilateral  
HELP\_TEXT.....No Help Text  
ID.....{C782FE56-6117-493F-847F-83173407F0B5}  
TERM.....BILATERAL HAND

LOCATION: BILATERAL LEG

-----  
ABBREVIATION.....BILAT LEG  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bilateral  
HELP\_TEXT.....No Help Text  
ID.....{0713761D-E0FC-DA03-02A1-F20EFEB99F45}  
TERM.....BILATERAL LEG

LOCATION: BILATERAL NECK

-----  
ABBREVIATION.....BILAT NECK  
ACTIVE.....Yes  
DESCRIPTION.....Bilateral Neck  
DISPLAY\_NAME.....Bilateral Neck  
HELP\_TEXT.....No Help Text  
ID.....{744BF4A2-C55E-4D55-8711-58700DD1C446}  
TERM.....BILATERAL NECK

LOCATION: BILATERAL PERIPHERALS

-----  
ABBREVIATION.....BILAT PERI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bilateral Peripherals  
HELP\_TEXT.....No Help Text  
ID.....{7F134529-36C5-E9D4-ABD5-64B404F9E706}  
TERM.....BILATERAL PERIPHERALS

LOCATION: BILATERAL WRIST

-----  
ABBREVIATION.....BILAT WRIS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bilateral  
HELP\_TEXT.....No Help Text  
ID.....{89B4D30F-9E46-F8D6-14F6-1CC0843DDC21}  
TERM.....BILATERAL WRIST

LOCATION: BLADDER

-----  
ABBREVIATION.....BLADDER  
ACTIVE.....Yes  
DESCRIPTION.....Bladder  
DISPLAY\_NAME.....Bladder  
HELP\_TEXT.....Bladder  
ID.....{18E083F2-C4BB-407E-B0A3-20C4942FC97F}  
TERM.....BLADDER

LOCATION: BOTH ARMS

-----  
ABBREVIATION.....BOTH ARMS  
ACTIVE.....Yes  
DESCRIPTION.....Both Arms  
DISPLAY\_NAME.....Both Arms  
HELP\_TEXT.....No Help Text  
ID.....{A493AB32-9DE9-40AC-8C46-7EFF395B63B5}  
TERM.....BOTH ARMS

LOCATION: BOTH HANDS

-----  
ABBREVIATION.....BOTH HANDS  
ACTIVE.....Yes  
DESCRIPTION.....Both Hands  
DISPLAY\_NAME.....Both Hands  
HELP\_TEXT.....No Help Text  
ID.....{63AA2685-CB1D-4680-B2F4-2199A22D1EC6}  
TERM.....BOTH HANDS

LOCATION: BOTH LEGS

-----  
ABBREVIATION.....BOTH LEGS  
ACTIVE.....Yes  
DESCRIPTION.....Both Legs  
DISPLAY\_NAME.....Both Legs  
HELP\_TEXT.....No Help Text  
ID.....{5A15A231-4695-4CF3-AD0F-2E7950208967}  
TERM.....BOTH LEGS

LOCATION: BRACHIAL

-----  
ABBREVIATION.....Bra  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Brachial  
HELP\_TEXT.....No Help Text  
ID.....{67426127-09E1-E566-F31E-15714A5B2B7C}  
TERM.....BRACHIAL

LOCATION: CALF

-----  
ABBREVIATION.....Clf  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Calf  
HELP\_TEXT.....No Help Text  
ID.....{7923E496-90D6-B6D5-726A-A5B8DF3436CF}  
TERM.....CALF

LOCATION: CAROTID

-----  
ABBREVIATION.....Car  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Carotid  
HELP\_TEXT.....No Help Text  
ID.....{9536F203-1797-9336-6F69-604249E1B3A0}  
TERM.....CAROTID

LOCATION: CHEST

-----  
ABBREVIATION.....CHEST  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Chest  
HELP\_TEXT.....No Help Text  
ID.....{8C3C30FE-6E06-C35E-8226-CBF0C94014CC}  
TERM.....CHEST

LOCATION: COCCYX

-----  
ABBREVIATION.....COCCYX  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Coccyx  
HELP\_TEXT.....No Help Text  
ID.....{49D65E9D-3D5B-40E3-B7CF-1003F80B852C}  
TERM.....COCCYX

LOCATION: CORE

-----  
ABBREVIATION.....C  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Core  
HELP\_TEXT.....No Help Text  
ID.....{5076C291-4EE7-7F9B-E7A5-3BB0A78639B7}  
TERM.....CORE

LOCATION: DONOR SITE

-----  
ABBREVIATION.....DONOR SITE  
ACTIVE.....Yes  
DESCRIPTION.....Donor Site  
DISPLAY\_NAME.....Donor Site  
HELP\_TEXT.....No Help Text  
ID.....{F0A65974-7650-4A4B-AC88-0945E37A0AD2}  
TERM.....DONOR SITE

LOCATION: DORSALIS PEDIS

-----  
ABBREVIATION.....DP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Dorsalis Pedis  
HELP\_TEXT.....No Help Text  
ID.....{9657D743-1632-2044-B4C5-9EF291A2F6D0}  
TERM.....DORSALIS PEDIS

LOCATION: FACE

-----  
ABBREVIATION.....FACE  
ACTIVE.....Yes  
DESCRIPTION.....PATIENT'S FACE  
DISPLAY\_NAME.....Face  
HELP\_TEXT.....No Help Text  
ID.....{C638EA14-7903-6A94-496A-C6560D51BEE3}  
TERM.....FACE

LOCATION: FEET

-----  
ABBREVIATION.....FEET  
ACTIVE.....Yes  
DESCRIPTION.....BOTH OF PATIENT'S FEET  
DISPLAY\_NAME.....Feet  
HELP\_TEXT.....No Help Text  
ID.....{944CAA0E-5165-A6AF-BE3D-85455440BF40}  
TERM.....FEET

LOCATION: FEMORAL

-----  
ABBREVIATION.....Fem  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Femoral  
HELP\_TEXT.....No Help Text  
ID.....{E3442C9E-8677-7296-D314-7365319C3716}  
TERM.....FEMORAL

LOCATION: FOREHEAD

-----  
ABBREVIATION.....FOREHEAD  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Forehead  
HELP\_TEXT.....No Help Text  
ID.....{5F892722-956E-474C-BA08-4074142D49F6}  
TERM.....FOREHEAD



LOCATION: HEAD

-----  
ABBREVIATION.....Hd  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Head  
HELP\_TEXT.....No Help Text  
ID.....{C3209F44-90C6-FF1B-3169-265B294F276A}  
TERM.....HEAD

LOCATION: HEAD OCCIPITAL

-----  
ABBREVIATION.....HEAD OCCIP  
ACTIVE.....Yes  
DESCRIPTION.....Head Occipital  
DISPLAY\_NAME.....Head Occipital  
HELP\_TEXT.....No Help Text  
ID.....{2066CA47-5A16-447F-9919-24ED435B627F}  
TERM.....HEAD OCCIPITAL

LOCATION: HEART

-----  
ABBREVIATION.....HEART  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Heart  
HELP\_TEXT.....No Help Text  
ID.....{F7D48A36-B512-42DD-B92F-07ACF2B95E08}  
TERM.....HEART

LOCATION: HIP

-----  
ABBREVIATION.....Hip  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Hip  
HELP\_TEXT.....No Help Text  
ID.....{D2C8BD11-7FB2-7217-A663-EBFC8C9428E7}  
TERM.....HIP

LOCATION: INCISIONAL SITE

-----  
ABBREVIATION.....INCISON SI  
ACTIVE.....Yes  
DESCRIPTION.....Incisional Site  
DISPLAY\_NAME.....Incisional Site  
HELP\_TEXT.....No Help Text  
ID.....{EB36A66D-751D-4031-9F68-0BB03E267321}  
TERM.....INCISIONAL SITE



LOCATION: KNEE

-----  
ABBREVIATION.....Kn  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Knee  
HELP\_TEXT.....No Help Text  
ID.....{B577308A-4C37-47BD-AD84-38E6342B6DF7}  
TERM.....KNEE

LOCATION: KNEES

-----  
ABBREVIATION.....KNEES  
ACTIVE.....Yes  
DESCRIPTION.....BOTH OF PATIENT'S KNEES  
DISPLAY\_NAME.....Knees  
HELP\_TEXT.....No Help Text  
ID.....{247D6621-3A85-F1F7-59BF-1FB0361B3C21}  
TERM.....KNEES

LOCATION: L ABDOMEN

-----  
ABBREVIATION.....L ABDOMEN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Abdomen  
HELP\_TEXT.....No Help Text  
ID.....{46B37481-4105-492D-96F0-E2D3895379C1}  
TERM.....L ABDOMEN

LOCATION: L ANKLE

-----  
ABBREVIATION.....L ANKLE  
ACTIVE.....Yes  
DESCRIPTION.....Left Ankle  
DISPLAY\_NAME.....L Ankle  
HELP\_TEXT.....No Help Text  
ID.....{D5A692EF-9D7C-473C-AF96-7EF1C01FA2AF}  
TERM.....L ANKLE

LOCATION: L ANTECUBITAL

-----  
ABBREVIATION.....L ANTECUBI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Antecubital  
HELP\_TEXT.....No Help Text  
ID.....{C21838CC-BAB5-1D7F-E2A4-A7C663FB28C3}  
TERM.....L ANTECUBITAL

LOCATION: L ARM

-----  
ABBREVIATION.....LA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Arm  
HELP\_TEXT.....No Help Text  
ID.....{567F552F-54B5-A8BB-596F-291001746E31}  
TERM.....L ARM

LOCATION: L AXILLARY

-----  
ABBREVIATION.....L AXILLARY  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Axillary  
HELP\_TEXT.....No Help Text  
ID.....{21AF62CA-BEED-2A96-2CD6-F2350A01995E}  
TERM.....L AXILLARY

LOCATION: L BRACHIAL

-----  
ABBREVIATION.....LBRA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Brachial  
HELP\_TEXT.....No Help Text  
ID.....{75F6CE33-1AB3-A1DA-70B8-04825290B577}  
TERM.....L BRACHIAL

LOCATION: L CALF

-----  
ABBREVIATION.....L CALF  
ACTIVE.....Yes  
DESCRIPTION.....Left Calf  
DISPLAY\_NAME.....L Calf  
HELP\_TEXT.....No Help Text  
ID.....{3DAA7817-936B-444F-8A2D-D3AB09144C10}  
TERM.....L CALF

LOCATION: L CAROTID

-----  
ABBREVIATION.....L CAROTID  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Carotid  
HELP\_TEXT.....No Help Text  
ID.....{FCACD1EB-C12B-2515-4A68-4E9E28A72DBC}  
TERM.....L CAROTID

LOCATION: L CHEST

-----  
ABBREVIATION.....L CHEST  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Left Chest  
HELP\_TEXT.....No Help Text  
ID.....{E332E497-FE0C-47CB-A0A7-217FBF52DC97}  
TERM.....L CHEST

LOCATION: L DORSALIS PEDIS

-----  
ABBREVIATION.....L DOR PED  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Dorsalis Pedis  
HELP\_TEXT.....No Help Text  
ID.....{389397B5-2863-B8F2-A779-C38433894A19}  
TERM.....L DORSALIS PEDIS

LOCATION: L EAR

-----  
ABBREVIATION.....L EAR  
ACTIVE.....Yes  
DESCRIPTION.....Left Ear  
DISPLAY\_NAME.....L Ear  
HELP\_TEXT.....No Help Text  
ID.....{DA9C2DD1-D6FF-4899-A62B-1096FAA36C59}  
TERM.....L EAR

LOCATION: L ELBOW

-----  
ABBREVIATION.....L ELBOW  
ACTIVE.....Yes  
DESCRIPTION.....Left Elbow  
DISPLAY\_NAME.....L Elbow  
HELP\_TEXT.....No Help Text  
ID.....{BF4317E7-2024-47EF-9E5F-3D21D0547885}  
TERM.....L ELBOW

LOCATION: L EXTERNAL JUGULAR

-----  
ABBREVIATION.....L EXT JUGU  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L External Jugular  
HELP\_TEXT.....No Help Text  
ID.....{6F7B331B-4F83-0B79-1755-F6F673D647D6}  
TERM.....L EXTERNAL JUGULAR

LOCATION: L EYE

-----  
ABBREVIATION.....L EYE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Eye  
HELP\_TEXT.....No Help Text  
ID.....{D714A0CA-EAD6-EC35-DA54-FF027C23FF45}  
TERM.....L EYE

LOCATION: L FACE

-----  
ABBREVIATION.....L FACE  
ACTIVE.....Yes  
DESCRIPTION.....LEFT SIDE OF PATIENT'S FACE  
DISPLAY\_NAME.....L Face  
HELP\_TEXT.....No Help Text  
ID.....{06F00A70-F475-CD40-665D-AF0569A1E75F}  
TERM.....L FACE

LOCATION: L FEMORAL

-----  
ABBREVIATION.....LFEM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Femoral  
HELP\_TEXT.....No Help Text  
ID.....{50AF6A08-CE37-6859-F4E6-762660A982F7}  
TERM.....L FEMORAL

LOCATION: L FINGER

-----  
ABBREVIATION.....L FINGER  
ACTIVE.....Yes  
DESCRIPTION.....Left Finger  
DISPLAY\_NAME.....L Finger  
HELP\_TEXT.....No Help Text  
ID.....{BB825645-1630-498C-A31B-6728D422EEFC}  
TERM.....L FINGER

LOCATION: L FLANK

-----  
ABBREVIATION.....L FLANK  
ACTIVE.....Yes  
DESCRIPTION.....Left Flank  
DISPLAY\_NAME.....L Flank  
HELP\_TEXT.....No Help Text  
ID.....{6E453E8D-B456-4BFD-A015-F9C013830F30}  
TERM.....L FLANK

LOCATION: L FOOT

-----  
ABBREVIATION.....L FOOT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Foot  
HELP\_TEXT.....No Help Text  
ID.....{2EC3FBEC-4EA2-9E7A-17D4-FDC2F9D28046}  
TERM.....L FOOT

LOCATION: L GLUTEUS

-----  
ABBREVIATION.....L GLUTEUS  
ACTIVE.....Yes  
DESCRIPTION.....Left Gluteus  
DISPLAY\_NAME.....L Gluteus  
HELP\_TEXT.....No Help Text  
ID.....{88C4E547-04CD-4D52-9416-CAEB6D898D21}  
TERM.....L GLUTEUS

LOCATION: L GROIN

-----  
ABBREVIATION.....L GROIN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Groin  
HELP\_TEXT.....No Help Text  
ID.....{5488E4EE-C1A2-4B02-8763-A53C412117A7}  
TERM.....L GROIN

LOCATION: L HAND

-----  
ABBREVIATION.....L HAND  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Hand  
HELP\_TEXT.....No Help Text  
ID.....{639DF25A-19E3-B504-D323-4B3F6ED39A24}  
TERM.....L HAND

LOCATION: L HEAD

-----  
ABBREVIATION.....L HEAD  
ACTIVE.....Yes  
DESCRIPTION.....Left Head  
DISPLAY\_NAME.....L Head  
HELP\_TEXT.....No Help Text  
ID.....{3737DC2D-DC9E-4FC8-839D-0355987B62D1}  
TERM.....L HEAD

LOCATION: L HEEL

-----  
ABBREVIATION.....L HEEL  
ACTIVE.....Yes  
DESCRIPTION.....Left Heel  
DISPLAY\_NAME.....L Heel  
HELP\_TEXT.....No Help Text  
ID.....{C9154881-17C7-4921-A0EE-BFB962C59CEF}  
TERM.....L HEEL

LOCATION: L HIP

-----  
ABBREVIATION.....L HIP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Hip  
HELP\_TEXT.....No Help Text  
ID.....{786CD881-27AD-4136-B9CC-8423CB38A431}  
TERM.....L HIP

LOCATION: L INGUINAL

-----  
ABBREVIATION.....L INGUINAL  
ACTIVE.....Yes  
DESCRIPTION.....Left Inguinal  
DISPLAY\_NAME.....L Inguinal  
HELP\_TEXT.....No Help Text  
ID.....{627D172E-8F33-48C9-BE89-D036E07D46BF}  
TERM.....L INGUINAL

LOCATION: L INNER KNEE

-----  
ABBREVIATION.....L INN KNEE  
ACTIVE.....Yes  
DESCRIPTION.....Left Inner Knee  
DISPLAY\_NAME.....L Inner Knee  
HELP\_TEXT.....No Help Text  
ID.....{79853581-CCF0-49E9-9D91-D16DE02A836D}  
TERM.....L INNER KNEE

LOCATION: L INTERNAL JUGULAR

-----  
ABBREVIATION.....LIJ  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Internal Jugular  
HELP\_TEXT.....No Help Text  
ID.....{AFBB4FE5-8A11-2047-7900-DADF2FC4430D}  
TERM.....L INTERNAL JUGULAR

LOCATION: L ISCHIAL

-----  
ABBREVIATION.....L ISCHIAL  
ACTIVE.....Yes  
DESCRIPTION.....Left Ischial  
DISPLAY\_NAME.....L Ischial  
HELP\_TEXT.....No Help Text  
ID.....{DCDAD741-9599-4479-9D55-3279F800DBE6}  
TERM.....L ISCHIAL

LOCATION: L JAW

-----  
ABBREVIATION.....L JAW  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Jaw  
HELP\_TEXT.....No Help Text  
ID.....{65EBD3CF-76D8-4F45-9B1B-3575D988BE8E}  
TERM.....L JAW

LOCATION: L KNEE

-----  
ABBREVIATION.....L KNEE  
ACTIVE.....Yes  
DESCRIPTION.....PATIENT'S LEFT KNEE  
DISPLAY\_NAME.....L Knee  
HELP\_TEXT.....No Help Text  
ID.....{3A9E30B4-09A9-967A-58EF-69A01EC4F859}  
TERM.....L KNEE

LOCATION: L LATERAL HEEL

-----  
ABBREVIATION.....L LAT HEEL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Lateral Heel  
HELP\_TEXT.....No Help Text  
ID.....{1A10FB20-7BAD-4D2B-9CBD-1A1A01DDCBF1}  
TERM.....L LATERAL HEEL

LOCATION: L LATERAL MALLEOLUS

-----  
ABBREVIATION.....L LAT MALL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Lateral Malleolus  
HELP\_TEXT.....No Help Text  
ID.....{9A523B7C-64DD-4E47-9543-165788108D4C}  
TERM.....L LATERAL MALLEOLUS

LOCATION: L LEG

-----  
ABBREVIATION.....LL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Leg  
HELP\_TEXT.....No Help Text  
ID.....{B088C809-E528-B785-21C2-14E45E1FD8AF}  
TERM.....L LEG

LOCATION: L LOWER ANTERIOR CHEST

-----  
ABBREVIATION.....L L A CHES  
ACTIVE.....Yes  
DESCRIPTION.....Left Lower Anterior Chest  
DISPLAY\_NAME.....L Low. Ant. Chest  
HELP\_TEXT.....No Help Text  
ID.....{6073B3C6-AA4D-4BBA-9FF2-32ED40CF92AD}  
TERM.....L LOWER ANTERIOR CHEST

LOCATION: L LOWER ARM

-----  
ABBREVIATION.....L LOW ARM  
ACTIVE.....Yes  
DESCRIPTION.....Left Lower Arm  
DISPLAY\_NAME.....L Lower Arm  
HELP\_TEXT.....No Help Text  
ID.....{63AA95E8-256D-4390-8CB3-F00A1B3F91CD}  
TERM.....L LOWER ARM

LOCATION: L LOWER EXTREMITY

-----  
ABBREVIATION.....LLEX  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Lower Extremity  
HELP\_TEXT.....No Help Text  
ID.....{FC6F5056-779A-BF76-86D0-A237936C1ECC}  
TERM.....L LOWER EXTREMITY

LOCATION: L LOWER LEG

-----  
ABBREVIATION.....L LOW LEG  
ACTIVE.....Yes  
DESCRIPTION.....Left Lower Leg  
DISPLAY\_NAME.....L Lower Leg  
HELP\_TEXT.....No Help Text  
ID.....{7B8E92FA-E065-4D06-9FB6-4234656CFD61}  
TERM.....L LOWER LEG



LOCATION: L LOWER LIMB

---

ABBREVIATION.....LLLI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Lower Limb  
HELP\_TEXT.....No Help Text  
ID.....{4838CB05-A1DD-F4FB-587C-2D0C17A2484E}  
TERM.....L LOWER LIMB

LOCATION: L LOWER POSTERIOR CHEST

---

ABBREVIATION.....L L P CHES  
ACTIVE.....Yes  
DESCRIPTION.....Left Lower Posterior Chest  
DISPLAY\_NAME.....L Low. Pos. Chest  
HELP\_TEXT.....No Help Text  
ID.....{2C19F10B-7546-4A1F-A5DC-3867359ED971}  
TERM.....L LOWER POSTERIOR CHEST

LOCATION: L MANDIBULAR

---

ABBREVIATION.....L MANDIBUL  
ACTIVE.....Yes  
DESCRIPTION.....Left Mandibular  
DISPLAY\_NAME.....L Mandibular  
HELP\_TEXT.....No Help Text  
ID.....{CB9DB373-5F0A-446B-909F-F76392CB3A04}  
TERM.....L MANDIBULAR

LOCATION: L MEDIAL HEEL

---

ABBREVIATION.....L MED HEEL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Medial Heel  
HELP\_TEXT.....No Help Text  
ID.....{47AD7BA6-4EBB-4A08-B03C-FE73EB240725}  
TERM.....L MEDIAL HEEL

LOCATION: L MEDIAL MALLEOLUS

---

ABBREVIATION.....L MED MALL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Medial Malleolus  
HELP\_TEXT.....No Help Text  
ID.....{330E7682-407F-43CD-84FC-1CDA41F08C68}  
TERM.....L MEDIAL MALLEOLUS

LOCATION: L MIDDLE ANTERIOR CHEST

-----  
ABBREVIATION.....L M A CHES  
ACTIVE.....Yes  
DESCRIPTION.....Left Middle Anterior Chest  
DISPLAY\_NAME.....L Mid. Ant. Chest  
HELP\_TEXT.....No Help Text  
ID.....{36A83F77-117E-474B-AE34-0B2C02E131AE}  
TERM.....L MIDDLE ANTERIOR CHEST

LOCATION: L MIDDLE POSTERIOR CHEST

-----  
ABBREVIATION.....L M P CHES  
ACTIVE.....Yes  
DESCRIPTION.....Left Middle Posterior Chest  
DISPLAY\_NAME.....L Mid. Pos. Chest  
HELP\_TEXT.....No Help Text  
ID.....{5A656939-F82B-4DEC-9F1D-B5FA644E3DA9}  
TERM.....L MIDDLE POSTERIOR CHEST

LOCATION: L NECK

-----  
ABBREVIATION.....L NECK  
ACTIVE.....Yes  
DESCRIPTION.....Left Neck  
DISPLAY\_NAME.....L Neck  
HELP\_TEXT.....No Help Text  
ID.....{C4EA6279-05E2-4195-A638-3794DCBA7B78}  
TERM.....L NECK

LOCATION: L NOSE

-----  
ABBREVIATION.....L NOSE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Nose  
HELP\_TEXT.....No Help Text  
ID.....{C7707847-5135-4749-B196-002E3E59BB54}  
TERM.....L NOSE

LOCATION: L OCCIPITAL

-----  
ABBREVIATION.....L OCCIPTAL  
ACTIVE.....Yes  
DESCRIPTION.....Left Occipital  
DISPLAY\_NAME.....L Occipital  
HELP\_TEXT.....No Help Text  
ID.....{A4234F41-CD6F-4CE1-9210-316DB847BDE2}  
TERM.....L OCCIPITAL

LOCATION: L ORAL

-----  
ABBREVIATION.....L ORAL  
ACTIVE.....Yes  
DESCRIPTION.....Left Oral  
DISPLAY\_NAME.....L Oral  
HELP\_TEXT.....No Help Text  
ID.....{057B787A-8750-4563-80DB-25110CCC1146}  
TERM.....L ORAL

LOCATION: L ORBIT

-----  
ABBREVIATION.....L ORBIT  
ACTIVE.....Yes  
DESCRIPTION.....Left Orbit  
DISPLAY\_NAME.....L Orbit  
HELP\_TEXT.....No Help Text  
ID.....{1D798E35-5095-4BA0-8F9F-FA96668A5239}  
TERM.....L ORBIT

LOCATION: L OUTER KNEE

-----  
ABBREVIATION.....L OUT KNEE  
ACTIVE.....Yes  
DESCRIPTION.....Left Outer Knee  
DISPLAY\_NAME.....L Outer Knee  
HELP\_TEXT.....No Help Text  
ID.....{B143C935-2687-49C3-9F65-8D69C1775ABC}  
TERM.....L OUTER KNEE

LOCATION: L PARIETAL

-----  
ABBREVIATION.....L PARIETAL  
ACTIVE.....Yes  
DESCRIPTION.....Left Parietal  
DISPLAY\_NAME.....L Parietal  
HELP\_TEXT.....No Help Text  
ID.....{2687ED24-F35E-456C-99C9-F6A8EC4364E0}  
TERM.....L PARIETAL

LOCATION: L PEDAL

-----  
ABBREVIATION.....LEFT PEDAL  
ACTIVE.....Yes  
DESCRIPTION.....L Pedal  
DISPLAY\_NAME.....L Pedal  
HELP\_TEXT.....L Pedal  
ID.....{D84E9953-8C17-3130-8F21-A30D2A99FDE4}  
TERM.....L PEDAL

LOCATION: L POST TIBIAL

-----  
ABBREVIATION.....L PO TIB  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Post Tibial  
HELP\_TEXT.....No Help Text  
ID.....{BED130EE-31CC-C7E9-6597-78685BCE91AC}  
TERM.....L POST TIBIAL

LOCATION: L RADIAL

-----  
ABBREVIATION.....L RAD  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Radial  
HELP\_TEXT.....No Help Text  
ID.....{302F78AD-4067-8B6B-BD15-993EC8561AE6}  
TERM.....L RADIAL

LOCATION: L RIBS

-----  
ABBREVIATION.....L RIBS  
ACTIVE.....Yes  
DESCRIPTION.....Left Ribs  
DISPLAY\_NAME.....L Ribs  
HELP\_TEXT.....No Help Text  
ID.....{BDFEE74F-B3CC-4AD9-AE49-78B0BF19A219}  
TERM.....L RIBS

LOCATION: L SCAPULA

-----  
ABBREVIATION.....L SCAPULA  
ACTIVE.....Yes  
DESCRIPTION.....Left Scapula  
DISPLAY\_NAME.....L Scapula  
HELP\_TEXT.....No Help Text  
ID.....{3381EC50-771B-41C8-85C8-6EB0D208D54E}  
TERM.....L SCAPULA

LOCATION: L SHOULDER

-----  
ABBREVIATION.....L SHOULDER  
ACTIVE.....Yes  
DESCRIPTION.....Left Shoulder  
DISPLAY\_NAME.....L Shoulder  
HELP\_TEXT.....No Help Text  
ID.....{91CBD2FA-2743-42EB-946B-6933980D9337}  
TERM.....L SHOULDER

LOCATION: L SUBCLAVIAN

---

ABBREVIATION.....LSCL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Subclavian  
HELP\_TEXT.....No Help Text  
ID.....{71B93594-15F9-AA4A-AEF2-E222EC28C00D}  
TERM.....L SUBCLAVIAN

LOCATION: L TEMPORAL

---

ABBREVIATION.....L TEMPO  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Temporal  
HELP\_TEXT.....No Help Text  
ID.....{375E65CF-A59E-A1E9-EB6D-D53CD818939E}  
TERM.....L TEMPORAL

LOCATION: L TESTICLE

---

ABBREVIATION.....L TESTICLE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Left Testicle  
HELP\_TEXT.....No Help Text  
ID.....{1E376539-E89F-4472-8274-CFD11275E11F}  
TERM.....L TESTICLE

LOCATION: L THIGH

---

ABBREVIATION.....L THIGH  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Thigh  
HELP\_TEXT.....No Help Text  
ID.....{16C2312F-1E0E-45E0-B675-160F40C23F57}  
TERM.....L THIGH

LOCATION: L THORACIC ANTERIOR

---

ABBREVIATION.....L THOR ANT  
ACTIVE.....Yes  
DESCRIPTION.....Left Thoracic Anterior  
DISPLAY\_NAME.....L Thoracic Anterior  
HELP\_TEXT.....No Help Text  
ID.....{C85B5FB9-153B-4A8B-9BE6-5524A572ACAC}  
TERM.....L THORACIC ANTERIOR

LOCATION: L THORACIC LATERAL

-----  
ABBREVIATION.....L THOR LAT  
ACTIVE.....Yes  
DESCRIPTION.....Left Thoracic Lateral  
DISPLAY\_NAME.....L Thoracic Lateral  
HELP\_TEXT.....No Help Text  
ID.....{48D26CEB-2E33-4F57-AB3A-F5A46108FFBF}  
TERM.....L THORACIC LATERAL

LOCATION: L THORACIC POSTERIOR

-----  
ABBREVIATION.....L THOR POS  
ACTIVE.....Yes  
DESCRIPTION.....Left Thoracic Posterior  
DISPLAY\_NAME.....L Thoracic Posterior  
HELP\_TEXT.....No Help Text  
ID.....{68E3381F-7605-4786-85A8-523D215A2B97}  
TERM.....L THORACIC POSTERIOR

LOCATION: L THORACOTOMY

-----  
ABBREVIATION.....L THORACOT  
ACTIVE.....Yes  
DESCRIPTION.....Left Thoracotomy  
DISPLAY\_NAME.....L Thoracotomy  
HELP\_TEXT.....No Help Text  
ID.....{313C11A1-559A-4F68-8718-4012416EC189}  
TERM.....L THORACOTOMY

LOCATION: L THUMB

-----  
ABBREVIATION.....L THUMB  
ACTIVE.....Yes  
DESCRIPTION.....Left Thumb  
DISPLAY\_NAME.....L Thumb  
HELP\_TEXT.....No Help Text  
ID.....{F82F4731-2CC5-4705-91BE-05E40C25B419}  
TERM.....L THUMB

LOCATION: L TOE

-----  
ABBREVIATION.....L TOE  
ACTIVE.....Yes  
DESCRIPTION.....Left Toe  
DISPLAY\_NAME.....L Toe  
HELP\_TEXT.....No Help Text  
ID.....{C666DE86-E7AE-4AB3-8AE1-91496219E63D}  
TERM.....L TOE

LOCATION: L TRAPEZIUS

-----  
ABBREVIATION.....L TRAPEZIU  
ACTIVE.....Yes  
DESCRIPTION.....Left Trapezius  
DISPLAY\_NAME.....L Trapezius  
HELP\_TEXT.....No Help Text  
ID.....{B74CF943-C10F-4A1B-A19A-4975E674A97A}  
TERM.....L TRAPEZIUS

LOCATION: L TROCHANTER

-----  
ABBREVIATION.....L TROCHANT  
ACTIVE.....Yes  
DESCRIPTION.....Left Trochanter  
DISPLAY\_NAME.....L Trochanter  
HELP\_TEXT.....No Help Text  
ID.....{C91647A4-255D-40E4-A24D-01274E5D799F}  
TERM.....L TROCHANTER

LOCATION: L UPPER ANTERIOR CHEST

-----  
ABBREVIATION.....L U A CHES  
ACTIVE.....Yes  
DESCRIPTION.....Left Upper Anterior Chest  
DISPLAY\_NAME.....L. Up. Ant. Chest  
HELP\_TEXT.....No Help Text  
ID.....{FE63B58A-DE39-4C03-AE61-E34A5487C767}  
TERM.....L UPPER ANTERIOR CHEST

LOCATION: L UPPER ARM

-----  
ABBREVIATION.....L UPP ARM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Upper Arm  
HELP\_TEXT.....No Help Text  
ID.....{5F915696-1B08-32E9-6B95-49246117D660}  
TERM.....L UPPER ARM

LOCATION: L UPPER EXTREMITY

-----  
ABBREVIATION.....LUEX  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Upper Extremity  
HELP\_TEXT.....No Help Text  
ID.....{48806E79-B1DA-129B-EB52-453292A3DDC6}  
TERM.....L UPPER EXTREMITY

LOCATION: L UPPER LEG

-----  
ABBREVIATION.....L UPPER LE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Left Upper Leg  
HELP\_TEXT.....No Help Text  
ID.....{1DB35895-F74C-4B88-A3EE-3E556CB6559A}  
TERM.....L UPPER LEG

LOCATION: L UPPER LIMB

-----  
ABBREVIATION.....L Upp Limb  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Upper Limb  
HELP\_TEXT.....No Help Text  
ID.....{E6453716-62D3-DFCB-F84D-F6F2FDFC20B7}  
TERM.....L UPPER LIMB

LOCATION: L UPPER POSTERIOR CHEST

-----  
ABBREVIATION.....L U P CHES  
ACTIVE.....Yes  
DESCRIPTION.....Left Upper Posterior Chest  
DISPLAY\_NAME.....L Up. Pos. Chest  
HELP\_TEXT.....No Help Text  
ID.....{6D5674E9-EF04-4EF1-84A0-3F5ED2FB532A}  
TERM.....L UPPER POSTERIOR CHEST

LOCATION: L WRIST

-----  
ABBREVIATION.....L WRIST  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L Wrist  
HELP\_TEXT.....No Help Text  
ID.....{7BCAFACE-FF3B-46DA-99F9-6437A90B7126}  
TERM.....L WRIST

LOCATION: LEFT

-----  
ABBREVIATION.....Lt  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Left  
HELP\_TEXT.....No Help Text  
ID.....{C9C09768-9EEA-8F2B-86D4-B53458472E4E}  
TERM.....LEFT



LOCATION: LEFT EYE

-----  
ABBREVIATION.....Left Eye  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Left Eye  
HELP\_TEXT.....  
ID.....{D1B20983-6076-47E6-9578-2EFB3CB6A47B}  
TERM.....LEFT EYE

LOCATION: LEFT LOWER LOBE

-----  
ABBREVIATION.....LLL  
ACTIVE.....Yes  
DESCRIPTION.....LEFT LUNG'S LOWER LOBE  
DISPLAY\_NAME.....L Lower Lobe  
HELP\_TEXT.....No Help Text  
ID.....{14155C13-874D-6E26-01D8-2C03B349F913}  
TERM.....LEFT LOWER LOBE

LOCATION: LEFT NECK

-----  
ABBREVIATION.....L NECK  
ACTIVE.....Yes  
DESCRIPTION.....LEFT SIDE OF PATIENT'S NECK  
DISPLAY\_NAME.....L Neck  
HELP\_TEXT.....No Help Text  
ID.....{40522B5F-F23E-EADB-AE64-8EDE3B9E1155}  
TERM.....LEFT NECK

LOCATION: LEFT UPPER LOBE

-----  
ABBREVIATION.....LUL  
ACTIVE.....Yes  
DESCRIPTION.....LEFT LUNG'S UPPER LOBE  
DISPLAY\_NAME.....L Upper Lobe  
HELP\_TEXT.....No Help Text  
ID.....{21BC3BA2-B22E-6A23-A377-5C2D5C1575C7}  
TERM.....LEFT UPPER LOBE

LOCATION: LOWER ABDOMEN

-----  
ABBREVIATION.....LO ABDOM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Lower Abdomen  
HELP\_TEXT.....No Help Text  
ID.....{077C9CF5-DCA7-4D94-A7A3-F3C153CEA47D}  
TERM.....LOWER ABDOMEN

LOCATION: LOWER ARM

-----  
ABBREVIATION.....LrA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Lower Arm  
HELP\_TEXT.....No Help Text  
ID.....{D7627720-D6B0-AC00-8FA3-F7B569ED06FC}  
TERM.....LOWER ARM

LOCATION: LOWER BACK

-----  
ABBREVIATION.....LOWER BACK  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Lower Back  
HELP\_TEXT.....No Help Text  
ID.....{7E761680-61BC-450F-9905-C45CE760EC45}  
TERM.....LOWER BACK

LOCATION: MEDIASTINAL

-----  
ABBREVIATION.....MEDIASTINA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Mediastinal  
HELP\_TEXT.....No Help Text  
ID.....{D99997D2-61EF-CAC2-A854-D906E3A97447}  
TERM.....MEDIASTINAL

LOCATION: MID BACK

-----  
ABBREVIATION.....MID BACK  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Mid Back  
HELP\_TEXT.....No Help Text  
ID.....{F8C515E0-F437-4596-B030-DA1811038CDF}  
TERM.....MID BACK

LOCATION: MID-STERNAL

-----  
ABBREVIATION.....MID-STERN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Mid-Sternal  
HELP\_TEXT.....No Help Text  
ID.....{3F5168F1-9082-4FB2-A364-3CDE33B1214B}  
TERM.....MID-STERNAL

LOCATION: MIDDLE

-----  
ABBREVIATION.....MIDDLE  
ACTIVE.....Yes  
DESCRIPTION.....MIDDLE  
DISPLAY\_NAME.....Middle  
HELP\_TEXT.....No Help Text  
ID.....{654E4D50-B74E-D58B-B3E3-AA3F218A285E}  
TERM.....MIDDLE

LOCATION: MIDLINE ABDOMINAL

-----  
ABBREVIATION.....MID ABDOMI  
ACTIVE.....Yes  
DESCRIPTION.....Midline Abdominal  
DISPLAY\_NAME.....Midline Abdominal  
HELP\_TEXT.....No Help Text  
ID.....{37ED2947-9BDA-4207-94F1-6B8C0FE12ECF}  
TERM.....MIDLINE ABDOMINAL

LOCATION: MIDLINE NECK

-----  
ABBREVIATION.....MID NECK  
ACTIVE.....Yes  
DESCRIPTION.....Midline Neck  
DISPLAY\_NAME.....Midline Neck  
HELP\_TEXT.....No Help Text  
ID.....{FF666713-278D-46D6-8D25-E5D57E78EA08}  
TERM.....MIDLINE NECK

LOCATION: NECK

-----  
ABBREVIATION.....Ne  
ACTIVE.....Yes  
DESCRIPTION.....PATIENTS NECK  
DISPLAY\_NAME.....Neck  
HELP\_TEXT.....No Help Text  
ID.....{D4812C3B-AAC2-1B61-679B-3CD8B5CB947C}  
TERM.....NECK

LOCATION: NOSE

-----  
ABBREVIATION.....NOSE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Nose  
HELP\_TEXT.....No Help Text  
ID.....{85F7003E-6D4E-4A79-A534-8A57FBC2D106}  
TERM.....NOSE

LOCATION: OCCIPITAL

-----  
ABBREVIATION.....OCCIPITAL  
ACTIVE.....Yes  
DESCRIPTION.....Occipital  
DISPLAY\_NAME.....Occipital  
HELP\_TEXT.....No Help Text  
ID.....{5A3E384B-1C87-4A19-A449-DB91AB443FFC}  
TERM.....OCCIPITAL

LOCATION: ORAL

-----  
ABBREVIATION.....O  
ACTIVE.....Yes  
DESCRIPTION.....Oral Description  
DISPLAY\_NAME.....Oral  
HELP\_TEXT.....Oral Help  
ID.....{7A7ACA1B-ACBD-B08D-D40C-54BA8E724B8C}  
TERM.....ORAL

LOCATION: PENIS

-----  
ABBREVIATION.....PENIS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Penis  
HELP\_TEXT.....No Help Text  
ID.....{2D82CE55-83B6-4160-A754-07C71FD388E8}  
TERM.....PENIS

LOCATION: PERINEAL

-----  
ABBREVIATION.....PERINEAL  
ACTIVE.....Yes  
DESCRIPTION.....Perineal  
DISPLAY\_NAME.....Perineal  
HELP\_TEXT.....No Help Text  
ID.....{1801441C-D5EF-4775-9E58-E4985496F7E0}  
TERM.....PERINEAL

LOCATION: PERIORBITAL

-----  
ABBREVIATION.....PERIORBITA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Periorbital  
HELP\_TEXT.....No Help Text  
ID.....{DBAC362B-6278-EDB4-7E28-3180321A38B2}  
TERM.....PERIORBITAL

LOCATION: PERIPHERAL

-----  
ABBREVIATION.....Per  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Peripheral  
HELP\_TEXT.....No Help Text  
ID.....{36B31698-E84A-5792-8EAE-F374DFD1B55A}  
TERM.....PERIPHERAL

LOCATION: PERIRECTAL

-----  
ABBREVIATION.....PERIRECTAL  
ACTIVE.....Yes  
DESCRIPTION.....Perirectal  
DISPLAY\_NAME.....Perirectal  
HELP\_TEXT.....No Help Text  
ID.....{560E1473-7B8F-4A75-8807-29315C908553}  
TERM.....PERIRECTAL

LOCATION: POPLITEAL

-----  
ABBREVIATION.....Pop  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Popliteal  
HELP\_TEXT.....No Help Text  
ID.....{77030DF2-CE06-D9E2-0684-EE70209CBE73}  
TERM.....POPLITEAL

LOCATION: POSTERIOR TIBIAL

-----  
ABBREVIATION.....PT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Posterior Tibial  
HELP\_TEXT.....No Help Text  
ID.....{010F2FFE-DEBC-DF9D-9CF2-AD192B16FFD5}  
TERM.....POSTERIOR TIBIAL

LOCATION: R ABDOMEN

-----  
ABBREVIATION.....R ABDOMEN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Abdomen  
HELP\_TEXT.....No Help Text  
ID.....{5C45EAA8-96F5-4EC7-8AB5-A9543CD6C414}  
TERM.....R ABDOMEN

LOCATION: R ANKLE

-----  
ABBREVIATION.....R ANKLE  
ACTIVE.....Yes  
DESCRIPTION.....Right Ankle  
DISPLAY\_NAME.....R Ankle  
HELP\_TEXT.....No Help Text  
ID.....{4E6DF352-F399-4CE9-908C-48132C3D856B}  
TERM.....R ANKLE

LOCATION: R ANTECUBITAL

-----  
ABBREVIATION.....R ANTECUBI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Antecubital  
HELP\_TEXT.....No Help Text  
ID.....{7FB00A65-B4B2-06A6-BD89-98AC112D8D4E}  
TERM.....R ANTECUBITAL

LOCATION: R ARM

-----  
ABBREVIATION.....RA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Arm  
HELP\_TEXT.....No Help Text  
ID.....{1C09017A-F10A-543D-A7EA-F13790D78692}  
TERM.....R ARM

LOCATION: R ARTERY

-----  
ABBREVIATION.....R ARTERY  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Artery  
HELP\_TEXT.....No Help Text  
ID.....{A2D76F33-A049-B3A0-2CBF-9F83FF52D894}  
TERM.....R ARTERY

LOCATION: R AXILLARY

-----  
ABBREVIATION.....R AXILLARY  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Axillary  
HELP\_TEXT.....No Help Text  
ID.....{3044F28B-FED2-4069-7E3B-2A691C5771F9}  
TERM.....R AXILLARY

LOCATION: R BRACHIAL

-----  
ABBREVIATION.....RBRA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Brachial  
HELP\_TEXT.....No Help Text  
ID.....{A88E69A8-F6AE-80B0-2CD4-687B121F9DDE}  
TERM.....R BRACHIAL

LOCATION: R CALF

-----  
ABBREVIATION.....R CALF  
ACTIVE.....Yes  
DESCRIPTION.....Right Calf  
DISPLAY\_NAME.....R Calf  
HELP\_TEXT.....No Help Text  
ID.....{0298A452-4901-44BF-A7A8-F42DD3BA8B63}  
TERM.....R CALF

LOCATION: R CAROTID

-----  
ABBREVIATION.....R CAROTID  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Carotid  
HELP\_TEXT.....No Help Text  
ID.....{8280ACD2-72D6-202A-89AF-86BCA7736427}  
TERM.....R CAROTID

LOCATION: R CHEST

-----  
ABBREVIATION.....R CHEST  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Right Chest  
HELP\_TEXT.....No Help Text  
ID.....{7E30EC0E-815A-4C2E-A9FC-3C2B2D6E3CD4}  
TERM.....R CHEST

LOCATION: R DORSALIS PEDIS

-----  
ABBREVIATION.....R DOR PED  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Dorsalis Pedis  
HELP\_TEXT.....No Help Text  
ID.....{47D4C64A-FD93-2999-7B46-F25307EA6511}  
TERM.....R DORSALIS PEDIS

LOCATION: R EAR

-----  
ABBREVIATION.....R EAR  
ACTIVE.....Yes  
DESCRIPTION.....Right Ear  
DISPLAY\_NAME.....R Ear  
HELP\_TEXT.....No Help Text  
ID.....{AE31327D-A03C-488C-8A8E-084E7E6E7FDE}  
TERM.....R EAR

LOCATION: R ELBOW

-----  
ABBREVIATION.....R ELBOW  
ACTIVE.....Yes  
DESCRIPTION.....Right Elbow  
DISPLAY\_NAME.....R Elbow  
HELP\_TEXT.....No Help Text  
ID.....{5BAFF87B-6993-4E76-862F-9F386D9B71C5}  
TERM.....R ELBOW

LOCATION: R EXTERNAL JUGULAR

-----  
ABBREVIATION.....R EXT JUGU  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R External Jugular  
HELP\_TEXT.....No Help Text  
ID.....{4D78F704-9E3B-907E-A792-4DC4F4CFF170}  
TERM.....R EXTERNAL JUGULAR

LOCATION: R EYE

-----  
ABBREVIATION.....R EYE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Eye  
HELP\_TEXT.....No Help Text  
ID.....{258F0C1D-67BD-ECF8-3152-EB289DEF4681}  
TERM.....R EYE

LOCATION: R FACE

-----  
ABBREVIATION.....R FACE  
ACTIVE.....Yes  
DESCRIPTION.....RIGHT SIDE OF PATIENT'S FACE  
DISPLAY\_NAME.....R Face  
HELP\_TEXT.....No Help Text  
ID.....{D1B15677-EF6F-0125-3601-68C1FE6386CE}  
TERM.....R FACE



LOCATION: R FEMORAL

-----  
ABBREVIATION.....R FEMORAL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Femoral  
HELP\_TEXT.....No Help Text  
ID.....{B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}  
TERM.....R FEMORAL

LOCATION: R FINGER

-----  
ABBREVIATION.....R FINGER  
ACTIVE.....Yes  
DESCRIPTION.....Right Finger  
DISPLAY\_NAME.....R Finger  
HELP\_TEXT.....No Help Text  
ID.....{E12CA0EF-3267-4B93-8643-230E2226C3D6}  
TERM.....R FINGER

LOCATION: R FLANK

-----  
ABBREVIATION.....R FLANK  
ACTIVE.....Yes  
DESCRIPTION.....Right Flank  
DISPLAY\_NAME.....R Flank  
HELP\_TEXT.....No Help Text  
ID.....{F5700ED2-646F-4B75-A55B-111C27FB939C}  
TERM.....R FLANK

LOCATION: R FOOT

-----  
ABBREVIATION.....R FOOT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Foot  
HELP\_TEXT.....No Help Text  
ID.....{7D9F5FF6-232C-2488-4504-A721FF06928F}  
TERM.....R FOOT

LOCATION: R GLUTEUS

-----  
ABBREVIATION.....R GLUTEUS  
ACTIVE.....Yes  
DESCRIPTION.....Right Gluteus  
DISPLAY\_NAME.....R Gluteus  
HELP\_TEXT.....No Help Text  
ID.....{B4D32045-68D1-4457-9106-7937887FA6F6}  
TERM.....R GLUTEUS

LOCATION: R GROIN

-----  
ABBREVIATION.....R GROIN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Groin  
HELP\_TEXT.....No Help Text  
ID.....{BAE1C889-8CCF-4CEB-9EBD-0243BC6F096A}  
TERM.....R GROIN

LOCATION: R HAND

-----  
ABBREVIATION.....R HAND  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Hand  
HELP\_TEXT.....No Help Text  
ID.....{A54D3A80-8537-D079-9CCE-7DC56346D132}  
TERM.....R HAND

LOCATION: R HEAD

-----  
ABBREVIATION.....R HEAD  
ACTIVE.....Yes  
DESCRIPTION.....Right Head  
DISPLAY\_NAME.....R Head  
HELP\_TEXT.....No Help Text  
ID.....{42A95AAE-3DA9-4C96-A4E7-A15C9840217A}  
TERM.....R HEAD

LOCATION: R HEEL

-----  
ABBREVIATION.....R HEEL  
ACTIVE.....Yes  
DESCRIPTION.....Right Heel  
DISPLAY\_NAME.....R Heel  
HELP\_TEXT.....No Help Text  
ID.....{0020DD20-2EA4-4362-ADF0-07BFA220C60E}  
TERM.....R HEEL

LOCATION: R HIP

-----  
ABBREVIATION.....R HIP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Hip  
HELP\_TEXT.....No Help Text  
ID.....{A48CBCAE-5404-4934-9321-ECEF660A5350}  
TERM.....R HIP

LOCATION: R INGUINAL

-----  
ABBREVIATION.....R INGUINAL  
ACTIVE.....Yes  
DESCRIPTION.....Right Inguinal  
DISPLAY\_NAME.....R Inguinal  
HELP\_TEXT.....No Help Text  
ID.....{395233DF-3B89-49E0-946B-EB5C02E51178}  
TERM.....R INGUINAL

LOCATION: R INNER KNEE

-----  
ABBREVIATION.....R INN KNEE  
ACTIVE.....Yes  
DESCRIPTION.....Right Inner Knee  
DISPLAY\_NAME.....R Inner Knee  
HELP\_TEXT.....No Help Text  
ID.....{2AE24BF1-EA1B-490D-9DB7-EEFC9314C560}  
TERM.....R INNER KNEE

LOCATION: R INTERNAL JUGULAR

-----  
ABBREVIATION.....R INT JUGU  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Internal Jugular  
HELP\_TEXT.....No Help Text  
ID.....{7ABF6916-D379-198D-6454-F0518A74844A}  
TERM.....R INTERNAL JUGULAR

LOCATION: R ISCHIAL

-----  
ABBREVIATION.....R ISCHIAL  
ACTIVE.....Yes  
DESCRIPTION.....Right Ischial  
DISPLAY\_NAME.....R Ischial  
HELP\_TEXT.....No Help Text  
ID.....{97E50F6C-CF75-4E59-AFDE-80FC5CE6D9E3}  
TERM.....R ISCHIAL

LOCATION: R JAW

-----  
ABBREVIATION.....R JAW  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Jaw  
HELP\_TEXT.....No Help Text  
ID.....{7846AEE2-B6F4-4B5A-909F-56E386FD82A8}  
TERM.....R JAW

LOCATION: R KNEE

-----  
ABBREVIATION.....R KNEE  
ACTIVE.....Yes  
DESCRIPTION.....PATIENT'S RIGHT KNEE  
DISPLAY\_NAME.....R Knee  
HELP\_TEXT.....No Help Text  
ID.....{BCD18D53-07B6-2C41-64FC-3A07A7E979E7}  
TERM.....R KNEE

LOCATION: R LATERAL HEEL

-----  
ABBREVIATION.....R LAT HEEL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Lateral Heel  
HELP\_TEXT.....No Help Text  
ID.....{A1896CF9-0246-4323-8DED-E4B5FB12CACE}  
TERM.....R LATERAL HEEL

LOCATION: R LATERAL MALLEOLUS

-----  
ABBREVIATION.....R LAT MALL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Lateral Malleolus  
HELP\_TEXT.....No Help Text  
ID.....{950FC45B-93E8-424F-A5D4-DE646166D7B0}  
TERM.....R LATERAL MALLEOLUS

LOCATION: R LEG

-----  
ABBREVIATION.....RL  
ACTIVE.....Yes  
DESCRIPTION.....Right Leg  
DISPLAY\_NAME.....R Leg  
HELP\_TEXT.....RIGHT LEG  
ID.....{C5E40A7B-6FC4-267D-8CE8-3571D21D81FD}  
TERM.....R LEG

LOCATION: R LOWER ANTERIOR CHEST

-----  
ABBREVIATION.....R L A CHES  
ACTIVE.....Yes  
DESCRIPTION.....Right Lower Anterior Chest  
DISPLAY\_NAME.....R Low. Ant. Chest  
HELP\_TEXT.....No Help Text  
ID.....{96BCF731-B661-4455-B5E7-C4AE2509505F}  
TERM.....R LOWER ANTERIOR CHEST

LOCATION: R LOWER ARM

-----  
ABBREVIATION.....R LOW ARM  
ACTIVE.....Yes  
DESCRIPTION.....Right Lower Arm  
DISPLAY\_NAME.....R Lower Arm  
HELP\_TEXT.....No Help Text  
ID.....{8DE6E277-EFFF-40DD-9E83-D7295AACB6AA}  
TERM.....R LOWER ARM

LOCATION: R LOWER EXTREMITY

-----  
ABBREVIATION.....R LOW EXTR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Lower Extremity  
HELP\_TEXT.....No Help Text  
ID.....{7FC3CECC-1E6E-A686-FF04-2BEF8A9B9DCA}  
TERM.....R LOWER EXTREMITY

LOCATION: R LOWER LEG

-----  
ABBREVIATION.....R LOW LEG  
ACTIVE.....Yes  
DESCRIPTION.....Right Lower Leg  
DISPLAY\_NAME.....R Lower Leg  
HELP\_TEXT.....No Help Text  
ID.....{C8EA4D71-D97E-46BD-AE9D-C10954C525CD}  
TERM.....R LOWER LEG

LOCATION: R LOWER LIMB

-----  
ABBREVIATION.....R LOW LIMB  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Lower Limb  
HELP\_TEXT.....No Help Text  
ID.....{A8E26A05-311A-F196-40B5-B5202D9F8058}  
TERM.....R LOWER LIMB

LOCATION: R LOWER LOBE

-----  
ABBREVIATION.....RLL  
ACTIVE.....Yes  
DESCRIPTION.....RIGHT LUNG'S LOWER LOBE  
DISPLAY\_NAME.....R Lower Lobe  
HELP\_TEXT.....No Help Text  
ID.....{418D350A-3602-1918-22C8-53D093AF49B0}  
TERM.....R LOWER LOBE

LOCATION: R LOWER POSTERIOR CHEST

-----  
ABBREVIATION.....R L P CHES  
ACTIVE.....Yes  
DESCRIPTION.....Right Lower Posterior Chest  
DISPLAY\_NAME.....R Low. Pos. Chest  
HELP\_TEXT.....No Help Text  
ID.....{5163180E-FFD2-4BFB-BF8B-4C6B18EFE5D4}  
TERM.....R LOWER POSTERIOR CHEST

LOCATION: R MANDIBULAR

-----  
ABBREVIATION.....R MANDIBUL  
ACTIVE.....Yes  
DESCRIPTION.....Right Mandibular  
DISPLAY\_NAME.....R Mandibular  
HELP\_TEXT.....No Help Text  
ID.....{303A1DE2-CB9A-498F-9CC5-DFA8FBFF10C1}  
TERM.....R MANDIBULAR

LOCATION: R MEDIAL HEEL

-----  
ABBREVIATION.....R MED HEEL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Medial Heel  
HELP\_TEXT.....No Help Text  
ID.....{C0526527-0AA9-4ECF-B768-99B2A1F5F3BF}  
TERM.....R MEDIAL HEEL

LOCATION: R MEDIAL MALLEOLUS

-----  
ABBREVIATION.....R MED MALL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Medial Malleolus  
HELP\_TEXT.....No Help Text  
ID.....{C70FB4FA-6516-4567-8E26-46BA439F1DF4}  
TERM.....R MEDIAL MALLEOLUS

LOCATION: R MIDDLE ANTERIOR CHEST

-----  
ABBREVIATION.....R M A CHES  
ACTIVE.....Yes  
DESCRIPTION.....Right Middle Anterior Chest  
DISPLAY\_NAME.....R Mid. Ant. Chest  
HELP\_TEXT.....No Help Text  
ID.....{2F751A9C-66CA-4C22-ACC9-6E81FA63FA2B}  
TERM.....R MIDDLE ANTERIOR CHEST

LOCATION: R MIDDLE LOBE

-----  
ABBREVIATION.....RML  
ACTIVE.....Yes  
DESCRIPTION.....RIGHT LUNG'S MIDDLE LOBE  
DISPLAY\_NAME.....R Middle Lobe  
HELP\_TEXT.....No Help Text  
ID.....{57879159-350C-2055-34D3-C6EFA39CA98F}  
TERM.....R MIDDLE LOBE

LOCATION: R MIDDLE POSTERIOR CHEST

-----  
ABBREVIATION.....R M P CHES  
ACTIVE.....Yes  
DESCRIPTION.....Right Middle Posterior Chest  
DISPLAY\_NAME.....R Mid. Pos. Chest  
HELP\_TEXT.....No Help Text  
ID.....{2BC487D6-6CE8-499E-94A0-A79EAA8E199B}  
TERM.....R MIDDLE POSTERIOR CHEST

LOCATION: R NECK

-----  
ABBREVIATION.....R NECK  
ACTIVE.....Yes  
DESCRIPTION.....RIGHT SIDE OF PATIENT'S NECK  
DISPLAY\_NAME.....R Neck  
HELP\_TEXT.....No Help Text  
ID.....{7E4D84D3-535B-6164-C6C1-837C6F997A0D}  
TERM.....R NECK

LOCATION: R NOSE

-----  
ABBREVIATION.....R NOSE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Nose  
HELP\_TEXT.....No Help Text  
ID.....{D80BCB19-E71F-43EA-8990-6C121EE7F7A1}  
TERM.....R NOSE

LOCATION: R OCCIPITAL

-----  
ABBREVIATION.....R OCCIPITA  
ACTIVE.....Yes  
DESCRIPTION.....Right Occipital  
DISPLAY\_NAME.....R Occipital  
HELP\_TEXT.....No Help Text  
ID.....{5DB951DB-CB03-482E-B0DA-CC5EE8FB1A67}  
TERM.....R OCCIPITAL

LOCATION: R ORAL

-----  
ABBREVIATION.....R ORAL  
ACTIVE.....Yes  
DESCRIPTION.....Right Oral  
DISPLAY\_NAME.....R Oral  
HELP\_TEXT.....No Help Text  
ID.....{9CEE934-AA3D-4BDE-BC3B-9E275C51CC9C}  
TERM.....R ORAL

LOCATION: R ORBIT

-----  
ABBREVIATION.....R ORBIT  
ACTIVE.....Yes  
DESCRIPTION.....Right Orbit  
DISPLAY\_NAME.....R Orbit  
HELP\_TEXT.....No Help Text  
ID.....{DFFD974A-1684-49AE-8230-8E12411EA62E}  
TERM.....R ORBIT

LOCATION: R OUTER KNEE

-----  
ABBREVIATION.....R OUT KNEE  
ACTIVE.....Yes  
DESCRIPTION.....Right Outer Knee  
DISPLAY\_NAME.....R Outer Knee  
HELP\_TEXT.....No Help Text  
ID.....{330F1DF1-8C99-4984-9965-7404ABC9562D}  
TERM.....R OUTER KNEE

LOCATION: R PARIETAL

-----  
ABBREVIATION.....R PARIETAL  
ACTIVE.....Yes  
DESCRIPTION.....Right Parietal  
DISPLAY\_NAME.....R Parietal  
HELP\_TEXT.....No Help Text  
ID.....{28C927CD-1342-483B-B366-0B7EA8ACABAD}  
TERM.....R PARIETAL

LOCATION: R PEDAL

-----  
ABBREVIATION.....R PEDAL  
ACTIVE.....Yes  
DESCRIPTION.....R Pedal  
DISPLAY\_NAME.....R Pedal  
HELP\_TEXT.....R Pedal  
ID.....{9532B614-F4E1-316C-1322-905EB58E58C4}  
TERM.....R PEDAL



LOCATION: R POST TIBIAL

-----  
ABBREVIATION.....R PO TIB  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Post Tibial  
HELP\_TEXT.....No Help Text  
ID.....{88CDA53A-5153-4950-701A-3F4F2007CC6B}  
TERM.....R POST TIBIAL

LOCATION: R RADIAL

-----  
ABBREVIATION.....R RADIAL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Radial  
HELP\_TEXT.....No Help Text  
ID.....{BA10E6E5-DDF4-345E-68AF-21D32E9D972D}  
TERM.....R RADIAL

LOCATION: R RIBS

-----  
ABBREVIATION.....R RIBS  
ACTIVE.....Yes  
DESCRIPTION.....Right Ribs  
DISPLAY\_NAME.....R Ribs  
HELP\_TEXT.....No Help Text  
ID.....{260FD564-5E58-419B-9F73-62B6037D1572}  
TERM.....R RIBS

LOCATION: R SCAPULA

-----  
ABBREVIATION.....R SCAPULA  
ACTIVE.....Yes  
DESCRIPTION.....Right Scapula  
DISPLAY\_NAME.....R Scapula  
HELP\_TEXT.....No Help Text  
ID.....{EA657261-4D72-40C2-B70D-2A2E27556A01}  
TERM.....R SCAPULA

LOCATION: R SHOULDER

-----  
ABBREVIATION.....R SHOULDER  
ACTIVE.....Yes  
DESCRIPTION.....Right Shoulder  
DISPLAY\_NAME.....R Shoulder  
HELP\_TEXT.....No Help Text  
ID.....{D27AE130-442E-4DD6-B391-8D7E8AD494A2}  
TERM.....R SHOULDER

LOCATION: R SUBCLAVIAN

-----  
ABBREVIATION.....R SUBCLAVI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Subclavian  
HELP\_TEXT.....No Help Text  
ID.....{53DEAE83-93AE-123B-14E5-4DE38CD4F752}  
TERM.....R SUBCLAVIAN

LOCATION: R TEMPORAL

-----  
ABBREVIATION.....R TEMPOR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Temporal  
HELP\_TEXT.....No Help Text  
ID.....{1F50461D-5ECF-25B2-A29A-E02803CB7814}  
TERM.....R TEMPORAL

LOCATION: R TESTICLE

-----  
ABBREVIATION.....R TESTICLE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Right Testicle  
HELP\_TEXT.....No Help Text  
ID.....{0947D277-50AD-47F8-A76D-0EC34EB3B168}  
TERM.....R TESTICLE

LOCATION: R THIGH

-----  
ABBREVIATION.....R THIGH  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Thigh  
HELP\_TEXT.....No Help Text  
ID.....{B7A2D4A0-801E-4578-967A-F1CFB86FB96B}  
TERM.....R THIGH

LOCATION: R THORACIC ANTERIOR

-----  
ABBREVIATION.....R THOR ANT  
ACTIVE.....Yes  
DESCRIPTION.....Right Thoracic Anterior  
DISPLAY\_NAME.....R Thoracic Anterior  
HELP\_TEXT.....No Help Text  
ID.....{BD3DF757-2C65-4D6A-9C80-5DC95BB41CFB}  
TERM.....R THORACIC ANTERIOR

LOCATION: R THORACIC LATERAL

-----  
ABBREVIATION.....R THOR LAT  
ACTIVE.....Yes  
DESCRIPTION.....Right Thoracic Lateral  
DISPLAY\_NAME.....R Thoracic Lateral  
HELP\_TEXT.....No Help Text  
ID.....{55C48294-D3F3-4013-B552-BA3344D9513B}  
TERM.....R THORACIC LATERAL

LOCATION: R THORACIC POSTERIOR

-----  
ABBREVIATION.....R THOR POS  
ACTIVE.....Yes  
DESCRIPTION.....Right Thoracic Posterior  
DISPLAY\_NAME.....R Thoracic Posterior  
HELP\_TEXT.....No Help Text  
ID.....{E91085D2-E665-4EB6-A086-8CF573DECD8F}  
TERM.....R THORACIC POSTERIOR

LOCATION: R THORACOTOMY

-----  
ABBREVIATION.....R THORACOT  
ACTIVE.....Yes  
DESCRIPTION.....Right Thoracotomy  
DISPLAY\_NAME.....R Thoracotomy  
HELP\_TEXT.....No Help Text  
ID.....{A18E4BB3-0718-48D6-BD49-720A70C0EF60}  
TERM.....R THORACOTOMY

LOCATION: R THUMB

-----  
ABBREVIATION.....R THUMB  
ACTIVE.....Yes  
DESCRIPTION.....Right Thumb  
DISPLAY\_NAME.....R Thumb  
HELP\_TEXT.....No Help Text  
ID.....{B414BD68-4155-45B3-985D-B2F8252A20FE}  
TERM.....R THUMB

LOCATION: R TOE

-----  
ABBREVIATION.....R TOE  
ACTIVE.....Yes  
DESCRIPTION.....Right Toe  
DISPLAY\_NAME.....R Toe  
HELP\_TEXT.....No Help Text  
ID.....{B65727AB-3205-4DCE-833F-6DA8554636F9}  
TERM.....R TOE

LOCATION: R TRAPEZIUS

-----  
ABBREVIATION.....R TRAPEZIU  
ACTIVE.....Yes  
DESCRIPTION.....Right Trapezius  
DISPLAY\_NAME.....R Trapezius  
HELP\_TEXT.....No Help Text  
ID.....{69753B67-AFA1-4C93-BA78-E266F9894E46}  
TERM.....R TRAPEZIUS

LOCATION: R TROCHANTER

-----  
ABBREVIATION.....R TROCHANT  
ACTIVE.....Yes  
DESCRIPTION.....Right Trochanter  
DISPLAY\_NAME.....R Trochanter  
HELP\_TEXT.....No Help Text  
ID.....{1A3E21B2-C391-40D2-8059-12290E667D9E}  
TERM.....R TROCHANTER

LOCATION: R UPPER ANTERIOR CHEST

-----  
ABBREVIATION.....R U A CHES  
ACTIVE.....Yes  
DESCRIPTION.....Right Upper Anterior Chest  
DISPLAY\_NAME.....R Up. Ant. Chest  
HELP\_TEXT.....No Help Text  
ID.....{67890049-A4A8-4241-BDF3-573729BD6EB5}  
TERM.....R UPPER ANTERIOR CHEST

LOCATION: R UPPER ARM

-----  
ABBREVIATION.....R UPP ARM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Upper Arm  
HELP\_TEXT.....No Help Text  
ID.....{3EF20E8E-F68C-E4F2-8F8C-FAD3B44DDAA4}  
TERM.....R UPPER ARM

LOCATION: R UPPER EXTREMITY

-----  
ABBREVIATION.....R UPP EXTR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Upper Extremity  
HELP\_TEXT.....No Help Text  
ID.....{C65BA269-8F05-6FB9-62CB-466A286CF2C0}  
TERM.....R UPPER EXTREMITY

LOCATION: R UPPER LEG

-----  
ABBREVIATION.....R UPPER LE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Right Upper Leg  
HELP\_TEXT.....No Help Text  
ID.....{113BD5ED-5A0A-47CA-B100-91B244BFD04A}  
TERM.....R UPPER LEG

LOCATION: R UPPER LIMB

-----  
ABBREVIATION.....R UPP LIMB  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Upper Limb  
HELP\_TEXT.....No Help Text  
ID.....{5E028BFF-ED9B-014F-74F0-C4137DB074FB}  
TERM.....R UPPER LIMB

LOCATION: R UPPER LOBE

-----  
ABBREVIATION.....RUL  
ACTIVE.....Yes  
DESCRIPTION.....RIGHT LUNG'S UPPER LOBE  
DISPLAY\_NAME.....R Upper Lobe  
HELP\_TEXT.....Enter the Right Upper Lobe of the lung  
ID.....{098C11B4-B556-FBBE-7E21-AEC9D7E4408A}  
TERM.....R UPPER LOBE

LOCATION: R UPPER POSTERIOR CHEST

-----  
ABBREVIATION.....R U P CHES  
ACTIVE.....Yes  
DESCRIPTION.....Right Upper Posterior Chest  
DISPLAY\_NAME.....R Up. Pos. Chest  
HELP\_TEXT.....No Help Text  
ID.....{D8F03FA5-1FD7-46EE-8FA4-17B3B7E62141}  
TERM.....R UPPER POSTERIOR CHEST

LOCATION: R WRIST

-----  
ABBREVIATION.....R WRIST  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....R Wrist  
HELP\_TEXT.....No Help Text  
ID.....{FDD9FE1D-F8A2-09AD-82D7-AA478693DD01}  
TERM.....R WRIST

LOCATION: RADIAL

-----  
ABBREVIATION.....RADIAL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Radial  
HELP\_TEXT.....No Help-Text  
ID.....{D52C4103-2AB3-A0C7-16E3-CF95A1B3B789}  
TERM.....RADIAL

LOCATION: RECTAL

-----  
ABBREVIATION.....RECTAL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Rectal  
HELP\_TEXT.....No Help-Text  
ID.....{0E7069F4-D3E6-67A1-E821-EF1879C3C267}  
TERM.....RECTAL

LOCATION: RIGHT

-----  
ABBREVIATION.....R  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Right  
HELP\_TEXT.....No Help-Text  
ID.....{38528AAC-C0A8-ABB7-A659-7EB54D81EA9F}  
TERM.....RIGHT

LOCATION: RIGHT EYE

-----  
ABBREVIATION.....Right Eye  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Right Eye  
HELP\_TEXT.....  
ID.....{56F81421-F1C7-4FD4-9757-E0EDEF26626E}  
TERM.....RIGHT EYE

LOCATION: SACRAL

-----  
ABBREVIATION.....SACRAL  
ACTIVE.....Yes  
DESCRIPTION.....PATIENT'S SACRAL AREA  
DISPLAY\_NAME.....Sacral  
HELP\_TEXT.....No Help Text  
ID.....{0CF8A8A7-AB35-B803-4544-0C91645D9D1C}  
TERM.....SACRAL

LOCATION: SCROTUM

-----  
ABBREVIATION.....SCROTUM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Scrotum  
HELP\_TEXT.....No Help Text  
ID.....{BDC109F5-0D30-4914-A62C-CEFEB448B547}  
TERM.....SCROTUM

LOCATION: SKIN

-----  
ABBREVIATION.....SKIN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Skin  
HELP\_TEXT.....No Help-Text  
ID.....{0803BC59-1582-AEC8-3DA2-A81363038210}  
TERM.....SKIN

LOCATION: SPINE CERVICAL

-----  
ABBREVIATION.....SPINE CERV  
ACTIVE.....Yes  
DESCRIPTION.....Spine Cervical  
DISPLAY\_NAME.....Spine Cervical  
HELP\_TEXT.....No Help Text  
ID.....{9D28E1BF-D014-4DA4-A832-1F3C25C148E9}  
TERM.....SPINE CERVICAL

LOCATION: SPINE LUMBAR

-----  
ABBREVIATION.....SPINE LUMB  
ACTIVE.....Yes  
DESCRIPTION.....Spine Lumbar  
DISPLAY\_NAME.....Spine Lumbar  
HELP\_TEXT.....No Help Text  
ID.....{E7E411AB-5C29-4906-B0A9-E8658AAFFF91}  
TERM.....SPINE LUMBAR

LOCATION: SPINE THORACIC

-----  
ABBREVIATION.....SPINE THOR  
ACTIVE.....Yes  
DESCRIPTION.....Spine Thoracic  
DISPLAY\_NAME.....Spine Thoracic  
HELP\_TEXT.....No Help Text  
ID.....{2CACBB4F-3EF3-4FBF-AC3E-41D72F0178CB}  
TERM.....SPINE THORACIC

LOCATION: STERNOTOMY

-----  
ABBREVIATION.....STERNOTOMY  
ACTIVE.....Yes  
DESCRIPTION.....Sternotomy  
DISPLAY\_NAME.....Sternotomy  
HELP\_TEXT.....No Help Text  
ID.....{DB97F6A0-9295-42B7-BC29-6A534CE259DD}  
TERM.....STERNOTOMY

LOCATION: STOMACH

-----  
ABBREVIATION.....STOMACH  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Stomach  
HELP\_TEXT.....No Help Text  
ID.....{613C20B9-CC9C-4F5D-B5B1-69B66A663540}  
TERM.....STOMACH

LOCATION: SUBSTERNAL

-----  
ABBREVIATION.....SUBSTERN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Substernal  
HELP\_TEXT.....No Help Text  
ID.....{3DE7E0F1-BCA4-43EC-82CB-A87C6506EF46}  
TERM.....SUBSTERNAL

LOCATION: TEETH

-----  
ABBREVIATION.....TEETH  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Teeth  
HELP\_TEXT.....No Help Text  
ID.....{38B36D8E-A9C3-4EE6-8727-EECD44C6034E}  
TERM.....TEETH

LOCATION: TEMPORAL

-----  
ABBREVIATION.....TEMPORAL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Temporal  
HELP\_TEXT.....No Help-Text  
ID.....{E312E0D7-0FE4-A41A-EAB0-ECCB085C3576}  
TERM.....TEMPORAL



LOCATION: THIGH

-----  
ABBREVIATION.....THIGH  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Thigh  
HELP\_TEXT.....No Help-Text  
ID.....{2DA93C27-12B0-C76C-85D5-B7D6375F3AD3}  
TERM.....THIGH

LOCATION: THYROID

-----  
ABBREVIATION.....THYROID  
ACTIVE.....Yes  
DESCRIPTION.....Thyroid  
DISPLAY\_NAME.....Thyroid  
HELP\_TEXT.....No Help Text  
ID.....{6AC68004-FE58-46D8-9A1E-4EEF72AAD71B}  
TERM.....THYROID

LOCATION: TOOTH

-----  
ABBREVIATION.....TOOTH  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Tooth  
HELP\_TEXT.....No Help Text  
ID.....{8407F780-7C7F-49BF-B1DD-4FDD578CEAB3}  
TERM.....TOOTH

LOCATION: TORSO

-----  
ABBREVIATION.....TORSO  
ACTIVE.....Yes  
DESCRIPTION.....PATIENT'S TORSO  
DISPLAY\_NAME.....Torso  
HELP\_TEXT.....No Help Text  
ID.....{EA1AD477-8BFF-E575-23C1-B9CFAF2B9986}  
TERM.....TORSO

LOCATION: TRACHEOSTOMY

-----  
ABBREVIATION.....TRACHEOSTO  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Tracheostomy  
HELP\_TEXT.....No Help Text  
ID.....{4D12329A-56C8-4175-986F-2F0B68EB7E6D}  
TERM.....TRACHEOSTOMY

LOCATION: TRANSVERSE ABDOMINAL

-----  
ABBREVIATION.....TRANS ABDO  
ACTIVE.....Yes  
DESCRIPTION.....Transverse Abdominal  
DISPLAY\_NAME.....Transverse Abdominal  
HELP\_TEXT.....No Help Text  
ID.....{1C10FF98-7B20-45D0-B420-714981832B3D}  
TERM.....TRANSVERSE ABDOMINAL

LOCATION: TYMPANIC

-----  
ABBREVIATION.....Tym  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Tympanic  
HELP\_TEXT.....No Help Text  
ID.....{0A58C14C-E447-E034-AA93-D41E6E0D2F69}  
TERM.....TYMPANIC

LOCATION: ULNAR

-----  
ABBREVIATION.....Uln  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Ulnar  
HELP\_TEXT.....No Help Text  
ID.....{81FC2EA0-472D-06D7-F5A2-E71B0B7DE091}  
TERM.....ULNAR

LOCATION: UMBILICAL

-----  
ABBREVIATION.....UMBILICAL  
ACTIVE.....Yes  
DESCRIPTION.....UMBILICAL  
DISPLAY\_NAME.....UMBILICAL  
HELP\_TEXT.....No Help Text  
ID.....{8874B300-0751-49AF-8211-CF4C009C68B5}  
TERM.....UMBILICAL

LOCATION: UPPER ABDOMEN

-----  
ABBREVIATION.....UP ABDOM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Upper Abdomen  
HELP\_TEXT.....No Help Text  
ID.....{234C677B-1CE6-498B-9375-BB5DC3958842}  
TERM.....UPPER ABDOMEN

LOCATION: UPPER ARM

-----  
ABBREVIATION.....UpA  
ACTIVE.....Yes  
DESCRIPTION.....Area of arm between elbow and shoulder  
DISPLAY\_NAME.....Upper Arm  
HELP\_TEXT.....UPPER ARM  
ID.....{B43A2B31-1F35-3F4D-5D23-882541ABE477}  
TERM.....UPPER ARM

LOCATION: UPPER BACK

-----  
ABBREVIATION.....UPPER BACK  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Upper Back  
HELP\_TEXT.....No Help Text  
ID.....{C1DA4294-E0DC-438B-810D-7F91E170244A}  
TERM.....UPPER BACK

LOCATION: WAIST

-----  
ABBREVIATION.....WAIST  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Waist  
HELP\_TEXT.....No Help Text  
ID.....{784DD311-8B48-4FE3-96B6-471CCE453CD5}  
TERM.....WAIST

LOCATION: WRIST

-----  
ABBREVIATION.....Wri  
ACTIVE.....Yes  
DESCRIPTION.....Joint between hand and forearm  
DISPLAY\_NAME.....Wrist  
HELP\_TEXT.....WRIST  
ID.....{C4F1453C-A20A-E48C-0BD8-FAF5E5568ADA}  
TERM.....WRIST

METHOD: ADULT

-----

ABBREVIATION.....Ad  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Adult  
HELP\_TEXT.....No Help Text  
ID.....{C2D49445-6556-8E6E-A8E9-CB3FA76BC22E}  
TERM.....ADULT

METHOD: ADULT CUFF

-----

ABBREVIATION.....AC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Adult Cuff  
HELP\_TEXT.....No Help Text  
ID.....{1222302E-6982-C627-8A0B-A651C98138E3}  
TERM.....ADULT CUFF

METHOD: AEROSOL/HUMIDIFIED MASK

-----

ABBREVIATION.....AM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Aerosol/Humidified Mask  
HELP\_TEXT.....No Help Text  
ID.....{CAF54B28-D5D2-8481-FAE3-9C74C2D51551}  
TERM.....AEROSOL/HUMIDIFIED MASK

METHOD: AFTER EXERCISE

-----

ABBREVIATION.....AE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....After Exercise  
HELP\_TEXT.....No Help Text  
ID.....{D129EF4B-3A99-9966-54F1-2C01A36814EE}  
TERM.....AFTER EXERCISE

METHOD: AMBU BAG

-----

ABBREVIATION.....AB  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Ambu Bag  
HELP\_TEXT.....No Help Text  
ID.....{E9D82DB9-D8A3-10E4-CE51-34BFF48F4603}  
TERM.....AMBU BAG

METHOD: ARTERIAL LINE

-----  
ABBREVIATION.....AL  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Arterial Line  
HELP\_TEXT.....No Help Text  
ID.....{D4B3F4ED-8B33-431E-4164-BD165E823776}  
TERM.....ARTERIAL LINE

METHOD: ASSISTED VENTILATOR

-----  
ABBREVIATION.....AV  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Assisted Ventilator  
HELP\_TEXT.....No Help Text  
ID.....{A540F333-5A09-5F21-7B51-2B5DF05E32C0}  
TERM.....ASSISTED VENTILATOR

METHOD: AT REST

-----  
ABBREVIATION.....AR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....At Rest  
HELP\_TEXT.....No Help Text  
ID.....{E8C289A6-0C3D-A7B2-BB85-7F1A23F22B50}  
TERM.....AT REST

METHOD: AUSCULTATE

-----  
ABBREVIATION.....Aus  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Auscultate  
HELP\_TEXT.....No Help Text  
ID.....{0C629C1F-92D3-C641-B02C-9EA43B401DC9}  
TERM.....AUSCULTATE

METHOD: BED

-----  
ABBREVIATION.....B  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bed  
HELP\_TEXT.....No Help Text  
ID.....{20188D9F-3815-1F7F-34EC-BA77F75B310E}  
TERM.....BED

METHOD: BIPAP

-----  
ABBREVIATION.....BIP  
ACTIVE.....Yes  
DESCRIPTION.....Biphasic Positive Airway Pressure  
DISPLAY\_NAME.....BiPAP  
HELP\_TEXT.....No Help Text  
ID.....{3CFD1D11-EF8D-58BA-6381-B24956344D2E}  
TERM.....BIPAP

METHOD: BIPAP-CPAP

-----  
ABBREVIATION.....PAP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....BiPAP-CPAP  
HELP\_TEXT.....No Help Text  
ID.....{2BCBE562-8F44-662F-469B-13AC0E863A3D}  
TERM.....BIPAP-CPAP

METHOD: BLADDER

-----  
ABBREVIATION.....BLAD  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bladder  
HELP\_TEXT.....No Help Text  
ID.....{E3B8D8AC-2846-B84B-75D9-7F50B55B4ADD}  
TERM.....BLADDER

METHOD: CARDIAC MONITOR

-----  
ABBREVIATION.....CARDIAC MO  
ACTIVE.....Yes  
DESCRIPTION.....Cardiac Monitor  
DISPLAY\_NAME.....Cardiac Monitor  
HELP\_TEXT.....Cardiac Monitor  
ID.....{E8B6B5A3-D143-3967-A362-019BD5EB6984}  
TERM.....CARDIAC MONITOR

METHOD: CARDIAC OUTPUT MONITOR

-----  
ABBREVIATION.....COM  
ACTIVE.....Yes  
DESCRIPTION.....A CO Monitor  
DISPLAY\_NAME.....Cardiac Output Monitor  
HELP\_TEXT.....Data from a Cardiac Output Monitor  
ID.....{CCBE5744-7CB9-9DFA-938F-FDDBC73B45B3}  
TERM.....CARDIAC OUTPUT MONITOR

METHOD: CATHETER

-----  
ABBREVIATION.....CATHETER  
ACTIVE.....Yes  
DESCRIPTION.....Catheter used to capture Intake and  
Output  
DISPLAY\_NAME.....Catheter - Method  
HELP\_TEXT.....No Help Text  
ID.....{1AA95835-FB32-4D1A-9C53-271EE38B2F63}  
TERM.....CATHETER

METHOD: CHAIR

-----  
ABBREVIATION.....Ch  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Chair  
HELP\_TEXT.....No Help Text  
ID.....{BBF5C05C-220D-D052-C7E8-236E791A4FC9}  
TERM.....CHAIR

METHOD: CHEST TUBE

-----  
ABBREVIATION.....CH TUBE M  
ACTIVE.....Yes  
DESCRIPTION.....Chest Tube  
DISPLAY\_NAME.....Chest Tube  
HELP\_TEXT.....Chest Tube  
ID.....{F4F467D3-BB57-4E6F-9EA0-C0FBFB66F5B1}  
TERM.....CHEST TUBE

METHOD: COLLOID

-----  
ABBREVIATION.....COLLOID M  
ACTIVE.....Yes  
DESCRIPTION.....Colloid as method for Intake and Output  
DISPLAY\_NAME.....Colloid - Method  
HELP\_TEXT.....No Help Text  
ID.....{2AE33508-5EA2-44DF-A401-74DD7B93AB6E}  
TERM.....COLLOID

METHOD: COMPUTED TOMOGRAPHY

-----  
ABBREVIATION.....COM TOMOGR  
ACTIVE.....Yes  
DESCRIPTION.....Computed Tomography  
DISPLAY\_NAME.....Computed Tomography  
HELP\_TEXT.....No Help Text  
ID.....{B4676AE0-BFAF-4649-AB08-3710E39F9B8D}  
TERM.....COMPUTED TOMOGRAPHY

METHOD: CONTRAST LEFT VENTRICULOGRAPHY

-----  
ABBREVIATION.....CON L VENT  
ACTIVE.....Yes  
DESCRIPTION.....Contrast Left Ventriculography  
DISPLAY\_NAME.....Contrast Left Ventriculography  
HELP\_TEXT.....No Help Text  
ID.....{35BEA63F-9696-4C03-AE1D-D76E4DDBB017}  
TERM.....CONTRAST LEFT VENTRICULOGRAPHY

METHOD: CONTROLLED VENTILATOR

-----  
ABBREVIATION.....CV  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Controlled Ventilator  
HELP\_TEXT.....No Help Text  
ID.....{0CFC2183-5245-E7AB-96ED-93AF4A3FEC64}  
TERM.....CONTROLLED VENTILATOR

METHOD: CPAP

-----  
ABBREVIATION.....CP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Cpap  
HELP\_TEXT.....No Help Text  
ID.....{A041FB71-08EF-F3FE-82FA-051FF74C824F}  
TERM.....CPAP

METHOD: CRRT

-----  
ABBREVIATION.....CRRT  
ACTIVE.....Yes  
DESCRIPTION.....Continuous Renal Replacement Therapy  
DISPLAY\_NAME.....CRRT  
HELP\_TEXT.....No Help Text  
ID.....{AB37FC1D-DD7A-4103-A777-7B8D0F876889}  
TERM.....CRRT

METHOD: CUFF

-----  
ABBREVIATION.....Cu  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Cuff  
HELP\_TEXT.....No Help Text  
ID.....{E59A2543-1860-1E1C-BEF1-6E7F6B6FBC5E}  
TERM.....CUFF



METHOD: CUFF-AUTOMATED

-----  
ABBREVIATION.....CA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Cuff-Automated  
HELP\_TEXT.....No Help Text  
ID.....{088AAD2B-2E7B-F835-376F-3D4340790976}  
TERM.....CUFF-AUTOMATED

METHOD: CUFF-MANUAL

-----  
ABBREVIATION.....CM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Cuff-Manual  
HELP\_TEXT.....No Help Text  
ID.....{5F56C163-2B04-112F-55DE-E31374E0FB58}  
TERM.....CUFF-MANUAL

METHOD: CVVHD

-----  
ABBREVIATION.....CVVHD  
ACTIVE.....Yes  
DESCRIPTION.....Central Venous - Venous Hemodialysis  
DISPLAY\_NAME.....CVVHD  
HELP\_TEXT.....No Help Text  
ID.....{4CC9AB0D-0B6C-43F9-A93E-A8228EC3A77F}  
TERM.....CVVHD

METHOD: DIALYSIS

-----  
ABBREVIATION.....DIALYSIS M  
ACTIVE.....Yes  
DESCRIPTION.....Dialysis used to describe Intake and  
Output  
DISPLAY\_NAME.....Dialysis - Method  
HELP\_TEXT.....No Help Text  
ID.....{1EAE28D4-3AC0-4DA1-9D09-3A65C38267F6}  
TERM.....DIALYSIS

METHOD: DIALYSIS CATHETER

-----  
ABBREVIATION.....DIALYSIS C  
ACTIVE.....Yes  
DESCRIPTION.....Dialysis Catheter  
DISPLAY\_NAME.....Dialysis Catheter  
HELP\_TEXT.....No Help Text  
ID.....{D9FF5CA8-09AF-4318-B784-DEB3A8F819D8}  
TERM.....DIALYSIS CATHETER

METHOD: DOPPLER

-----  
ABBREVIATION.....Dop  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Doppler  
HELP\_TEXT.....No Help Text  
ID.....{52683C19-F9DC-876B-CA3F-53E7C864CF04}  
TERM.....DOPPLER

METHOD: DOUBLE LUMEN CATHETER

-----  
ABBREVIATION.....DOUBLE LUM  
ACTIVE.....Yes  
DESCRIPTION.....Double Lumen Catheter  
DISPLAY\_NAME.....Double Lumen Catheter  
HELP\_TEXT.....No Help Text  
ID.....{ED646D10-3DBE-4049-ADC6-5E3DD2C351DD}  
TERM.....DOUBLE LUMEN CATHETER

METHOD: ECHOCARDIOGRAPHY

-----  
ABBREVIATION.....ECG  
ACTIVE.....Yes  
DESCRIPTION.....Echocardiography  
DISPLAY\_NAME.....Echocardiography  
HELP\_TEXT.....No Help Text  
ID.....{CA44697E-3BBA-4767-A30D-6548A72F9639}  
TERM.....ECHOCARDIOGRAPHY

METHOD: ELECTROCARDIOGRAM

-----  
ABBREVIATION.....ECG  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Electrocardiogram  
HELP\_TEXT.....No Help Text  
ID.....{CF5F1AB9-3B74-7FBA-E38B-C9EF54FE4DC5}  
TERM.....ELECTROCARDIOGRAM

METHOD: ELECTROCARDIOGRAM TRACE 1

-----  
ABBREVIATION.....ECG TR 1  
ACTIVE.....Yes  
DESCRIPTION.....Trace 1 from ECG  
DISPLAY\_NAME.....ECG Trace 1  
HELP\_TEXT.....No Help Text  
ID.....{E0D449BA-9F06-4A27-8FE9-24C1BA0746DB}  
TERM.....ELECTROCARDIOGRAM TRACE 1

METHOD: ELECTROCARDIOGRAM TRACE 2

-----  
ABBREVIATION.....ECG TR 2  
ACTIVE.....Yes  
DESCRIPTION.....Trace 2 from ECG  
DISPLAY\_NAME.....ECG Trace 2  
HELP\_TEXT.....No Help Text  
ID.....{73EACBC1-C1FB-4B2C-96AF-C5C1CA07DA72}  
TERM.....ELECTROCARDIOGRAM TRACE 2

METHOD: ENTERAL FEEDING / SUPPLEMENT

-----  
ABBREVIATION.....EF/SUPPLEM  
ACTIVE.....Yes  
DESCRIPTION.....Enteral Feeding / Supplement - Method  
DISPLAY\_NAME.....Enteral Feeding / Supplement - Method  
HELP\_TEXT.....No Help Text  
ID.....{8389A9F5-0D48-496D-A605-1560B055DCC0}  
TERM.....ENTERAL FEEDING / SUPPLEMENT

METHOD: ENTERAL FEEDING TUBE

-----  
ABBREVIATION.....ENT FEED T  
ACTIVE.....Yes  
DESCRIPTION.....EFT for describing Intakes and Outputs  
DISPLAY\_NAME.....Enteral Feeding Tube - Method  
HELP\_TEXT.....No Help Text  
ID.....{487E3F10-E1F6-4EFF-9C50-9EC5BB7AE3CF}  
TERM.....ENTERAL FEEDING TUBE

METHOD: FACE TENT

-----  
ABBREVIATION.....FT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Face Tent  
HELP\_TEXT.....No Help Text  
ID.....{1329BF1A-DCFD-9AF3-69DB-BA6D4BC3F52A}  
TERM.....FACE TENT

METHOD: GI TUBE

-----  
ABBREVIATION.....GI TUBE  
ACTIVE.....Yes  
DESCRIPTION.....Used to describe Intakes and Outputs  
DISPLAY\_NAME.....GI Tube - Method  
HELP\_TEXT.....No Help Text  
ID.....{AF64F7BE-E72E-4391-BC1F-48163F2F3751}  
TERM.....GI TUBE

METHOD: HEMODIALYSIS

-----  
ABBREVIATION.....HEMODIALYS  
ACTIVE.....Yes  
DESCRIPTION.....Hemodialysis  
DISPLAY\_NAME.....Hemodialysis  
HELP\_TEXT.....No Help Text  
ID.....{DE7CD394-9E69-46B0-A207-5A18F7C5A1A7}  
TERM.....HEMODIALYSIS

METHOD: HEMOVAC

-----  
ABBREVIATION.....HEMOVAC  
ACTIVE.....Yes  
DESCRIPTION.....Hemovac Drain  
DISPLAY\_NAME.....Hemovac Drain  
HELP\_TEXT.....No Help Text  
ID.....{79DA0646-B811-4954-8BE1-C133398CD574}  
TERM.....HEMOVAC

METHOD: IABP

-----  
ABBREVIATION.....IABP  
ACTIVE.....Yes  
DESCRIPTION.....Intra-aortic Balloon Pump  
DISPLAY\_NAME.....IABP  
HELP\_TEXT.....Intra-aortic Balloon Pump  
ID.....{4352A41F-AAF3-4544-9590-BAFD7977DF49}  
TERM.....IABP

METHOD: INTRODUCER CATHETER

-----  
ABBREVIATION.....INTRODUCER  
ACTIVE.....Yes  
DESCRIPTION.....Introducer Catheter  
DISPLAY\_NAME.....Introducer Catheter  
HELP\_TEXT.....No Help Text  
ID.....{60EED7D3-6CB2-4FCE-BFC0-DE0BA8AB0814}  
TERM.....INTRODUCER CATHETER

METHOD: INVASIVE

-----  
ABBREVIATION.....I  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Invasive  
HELP\_TEXT.....No Help Text  
ID.....{5DB8746E-5446-F47A-883C-A44B91275A6E}  
TERM.....INVASIVE

METHOD: IRRIGATION

-----  
ABBREVIATION.....IRRIGATION  
ACTIVE.....Yes  
DESCRIPTION.....Irrigation  
DISPLAY\_NAME.....Irrigation  
HELP\_TEXT.....No Help Text  
ID.....{6C968D0A-400F-448A-93EF-02CB891BA981}  
TERM.....IRRIGATION

METHOD: IV

-----  
ABBREVIATION.....IV  
ACTIVE.....Yes  
DESCRIPTION.....Method to capture Intakes and Outputs  
DISPLAY\_NAME.....IV  
HELP\_TEXT.....No Help Text  
ID.....{52141E2F-7942-4C60-9B08-62F3B4640020}  
TERM.....IV

METHOD: JP DRAIN

-----  
ABBREVIATION.....JP DRAIN  
ACTIVE.....Yes  
DESCRIPTION.....JP Drain  
DISPLAY\_NAME.....JP Drain  
HELP\_TEXT.....No Help Text  
ID.....{94A32247-6190-4047-B581-A212FE4C9453}  
TERM.....JP DRAIN

METHOD: LARYNGEAL MASK

-----  
ABBREVIATION.....LARYNGMASK  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Laryngeal Mask  
HELP\_TEXT.....No Help Text  
ID.....{DE9B1E53-711C-DA9D-59D9-CD0EBB6CD6BE}  
TERM.....LARYNGEAL MASK

METHOD: LEAD

-----  
ABBREVIATION.....L  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Lead  
HELP\_TEXT.....No Help Text  
ID.....{B628FDCE-05E0-DACF-E66E-0F8380AB465A}  
TERM.....LEAD

METHOD: LG ADULT

-----  
ABBREVIATION.....LgA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Lg Adult  
HELP\_TEXT.....No Help Text  
ID.....{0C7B9E22-A426-80D1-80E8-2082F15784BA}  
TERM.....LG ADULT

METHOD: LG ADULT CUFF

-----  
ABBREVIATION.....LAC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Lg Adult Cuff  
HELP\_TEXT.....No Help Text  
ID.....{4FA9374C-AE66-DF6E-E0E8-673D7B9B2E19}  
TERM.....LG ADULT CUFF

METHOD: LIFT SCALE

-----  
ABBREVIATION.....LS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Lift Scale  
HELP\_TEXT.....No Help Text  
ID.....{8D741461-25C5-3955-BE2C-CBD27317CBA1}  
TERM.....LIFT SCALE

METHOD: MAGNETIC RESONANCE IMAGING

-----  
ABBREVIATION.....MRI  
ACTIVE.....Yes  
DESCRIPTION.....Magnetic Resonance Imaging  
DISPLAY\_NAME.....MRI  
HELP\_TEXT.....No Help Text  
ID.....{50EBFFF8-DF0D-41F2-B11B-33A26AA40DBA}  
TERM.....MAGNETIC RESONANCE IMAGING

METHOD: MASK

-----  
ABBREVIATION.....M  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Mask  
HELP\_TEXT.....No Help Text  
ID.....{4FB6B24A-D56E-4EE3-6CBC-E54C9734424A}  
TERM.....MASK

METHOD: MONITOR

-----  
ABBREVIATION.....Mo  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Monitor  
HELP\_TEXT.....No Help Text  
ID.....{63FC24EA-50B8-B302-D948-ABE3517CF6FC}  
TERM.....MONITOR

METHOD: MULTIGATED NUCLEAR VENTRICULOGRAPHY

-----  
ABBREVIATION.....MUGA  
ACTIVE.....Yes  
DESCRIPTION.....Multigated Nuclear Ventriculography  
DISPLAY\_NAME.....MUGA  
HELP\_TEXT.....No Help Text  
ID.....{4D73168B-4DF4-4EA2-83B9-ACE853D5DD42}  
TERM.....MULTIGATED NUCLEAR VENTRICULOGRAPHY

METHOD: MUSTACHE CANNULA

-----  
ABBREVIATION.....MC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Mustache Cannula  
HELP\_TEXT.....No Help Text  
ID.....{DAB3F52F-A8E1-2D9E-C8E5-1D005EC679FC}  
TERM.....MUSTACHE CANNULA

METHOD: NASAL CANNULA

-----  
ABBREVIATION.....NC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Nasal Cannula  
HELP\_TEXT.....No Help Text  
ID.....{EFFF38AE-6AF5-77FD-72E1-E0F7E3DEBA29}  
TERM.....NASAL CANNULA

METHOD: NEBULIZER

-----  
ABBREVIATION.....NEBULIZER  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Nebulizer  
HELP\_TEXT.....No Help Text  
ID.....{3CF184A3-2688-4988-B956-D5B0F303E8E8}  
TERM.....NEBULIZER

METHOD: NON RE-BREATHER

-----  
ABBREVIATION.....NRB  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Non Re-Breather  
HELP\_TEXT.....No Help Text  
ID.....{1B9AA055-FD7D-7DCC-4EDB-BDE3C821AA8E}  
TERM.....NON RE-BREATHER

METHOD: NON-INVASIVE

-----  
ABBREVIATION.....N-I  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Non-Invasive  
HELP\_TEXT.....No Help Text  
ID.....{F50ABFED-5FD7-FC95-20FE-1E3038C9F299}  
TERM.....NON-INVASIVE

METHOD: NONINVASIVE

-----  
ABBREVIATION.....N  
ACTIVE.....No  
DESCRIPTION.....Noninvasive  
DISPLAY\_NAME.....Noninvasive  
HELP\_TEXT.....No Help Text  
ID.....{398E91DF-F425-5083-AA35-C5263792845E}  
TERM.....NONINVASIVE

METHOD: NUCLEAR PERFUSION IMAGING

-----  
ABBREVIATION.....NPI  
ACTIVE.....Yes  
DESCRIPTION.....Nuclear Perfusion Imaging  
DISPLAY\_NAME.....Nuclear Perfusion Imaging  
HELP\_TEXT.....No Help Text  
ID.....{D0F7DB21-0D66-4161-8214-DACCAF41ADAC}  
TERM.....NUCLEAR PERFUSION IMAGING

METHOD: O2 SATS

-----  
ABBREVIATION.....O2  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....O2 Sats  
HELP\_TEXT.....No Help Text  
ID.....{66C2100C-8A52-F36A-F7F1-A55D4E4D7CE6}  
TERM.....O2 SATS



METHOD: OTHER

-----  
ABBREVIATION.....OTHER  
ACTIVE.....Yes  
DESCRIPTION.....Other: Comments Required  
DISPLAY\_NAME.....OTHER: OBSERVATION COMMENTS REQUIRED  
HELP\_TEXT.....No Help Text  
ID.....{AED275BD-8A42-46BE-BCDF-472EF4C31D65}  
TERM.....OTHER

METHOD: P.O.

-----  
ABBREVIATION.....PO  
ACTIVE.....Yes  
DESCRIPTION.....By Mouth, Oral, "Per Os"  
DISPLAY\_NAME.....P.O. - Method  
HELP\_TEXT.....No Help Text  
ID.....{7D50D8E7-DA7B-4605-AB43-1E97DEBEACEE}  
TERM.....P.O.

METHOD: PA CATHETER

-----  
ABBREVIATION.....PA CATHETE  
ACTIVE.....Yes  
DESCRIPTION.....Pulmonary Catheter  
DISPLAY\_NAME.....PA Catheter  
HELP\_TEXT.....No Help Text  
ID.....{5A518696-41B7-4AB8-8FC7-5CF454CE0299}  
TERM.....PA CATHETER

METHOD: PA LINE

-----  
ABBREVIATION.....PALI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....PA Line  
HELP\_TEXT.....No Help Text  
ID.....{1AEA41D7-901B-96FF-BF9B-420C24FF00E4}  
TERM.....PA LINE

METHOD: PALPATED

-----  
ABBREVIATION.....Pal  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Palpated  
HELP\_TEXT.....No Help Text  
ID.....{295D8E72-12A1-08AF-6BE8-0F9CD238F5EF}  
TERM.....PALPATED

METHOD: PARTIAL RE-BREATHER

-----  
ABBREVIATION.....PRB  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Partial Re-Breather  
HELP\_TEXT.....No Help Text  
ID.....{6F1DF170-5128-2E3F-1831-12FB0BB5E95C}  
TERM.....PARTIAL RE-BREATHER

METHOD: PEDIATRIC CUFF

-----  
ABBREVIATION.....PdC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Pediatric Cuff  
HELP\_TEXT.....No Help Text  
ID.....{CA9EA220-EF92-EE77-FBEA-4E49C4730A3E}  
TERM.....PEDIATRIC CUFF

METHOD: PENDANT CANNULA

-----  
ABBREVIATION.....PC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Pendant Cannula  
HELP\_TEXT.....No Help Text  
ID.....{A1451228-4ABB-DBCB-31AD-578BCD93012F}  
TERM.....PENDANT CANNULA

METHOD: PENROSE

-----  
ABBREVIATION.....PENROSE  
ACTIVE.....Yes  
DESCRIPTION.....Penrose Drain  
DISPLAY\_NAME.....Penrose Drain  
HELP\_TEXT.....No Help Text  
ID.....{CBC3C6DE-E5DE-41D5-83AB-30A4F3CFA9E1}  
TERM.....PENROSE

METHOD: PERITONEAL DIALYSIS

-----  
ABBREVIATION.....PERI DIALY  
ACTIVE.....Yes  
DESCRIPTION.....Peritoneal Dialysis  
DISPLAY\_NAME.....Peritoneal Dialysis  
HELP\_TEXT.....No Help Text  
ID.....{6B4AA38B-AFCC-42CE-9F32-59B32E8A6BBB}  
TERM.....PERITONEAL DIALYSIS

METHOD: PICC LINE

-----  
ABBREVIATION.....PICC LINE  
ACTIVE.....Yes  
DESCRIPTION.....Peripherally Inserted Central  
                  Catheter Line  
DISPLAY\_NAME.....PICC Line  
HELP\_TEXT.....No Help Text  
ID.....{BDB5075B-DECC-443A-B774-CFA062DF34F4}  
TERM.....PICC LINE

METHOD: PULSE OX MONITOR

-----  
ABBREVIATION.....POM  
ACTIVE.....Yes  
DESCRIPTION.....A Pulse Ox Monitor  
DISPLAY\_NAME.....Pulse Ox Monitor  
HELP\_TEXT.....Readings that have come from a Pulse  
                  Ox Monitor  
ID.....{A7E3E90C-4C8B-87D7-C2F8-49D9F3846CE8}  
TERM.....PULSE OX MONITOR

METHOD: RECTAL TEMP PROBE

-----  
ABBREVIATION.....RTP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Rectal Temp Probe  
HELP\_TEXT.....No Help Text  
ID.....{0100F644-EF8E-5FE3-0022-804408DCE2DD}  
TERM.....RECTAL TEMP PROBE

METHOD: RESERVOIR CANNULA

-----  
ABBREVIATION.....RC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Reservoir Cannula  
HELP\_TEXT.....No Help Text  
ID.....{B332C44D-29BE-765C-932C-63EBFC9631FC}  
TERM.....RESERVOIR CANNULA

METHOD: RHYTHM

-----  
ABBREVIATION.....RHY  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Rhythm  
HELP\_TEXT.....No Help Text  
ID.....{8A4AC126-12DA-0904-C364-B5D818621E32}  
TERM.....RHYTHM

METHOD: ROOM AIR

-----  
ABBREVIATION.....ROOM AIR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Room Air  
HELP\_TEXT.....No Help Text  
ID.....{B4103F11-4BCE-49ED-AF3E-ACCB526D938C}  
TERM.....ROOM AIR

METHOD: ROOM AIR15

-----  
ABBREVIATION.....RA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Room Air15  
HELP\_TEXT.....No Help Text  
ID.....{55AC16E4-4BAD-EDE4-41E9-853C78D8474B}  
TERM.....ROOM AIR15

METHOD: ROOM AIR20

-----  
ABBREVIATION.....xx  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Room Air20  
HELP\_TEXT.....No Help Text  
ID.....{FBBBD3CC2-53AF-FBCC-EF16-A6300F7F8330}  
TERM.....ROOM AIR20

METHOD: SINGLE LUMEN CATHETER

-----  
ABBREVIATION.....SING LUM C  
ACTIVE.....Yes  
DESCRIPTION.....Single Lumen Catheter  
DISPLAY\_NAME.....Single Lumen Catheter  
HELP\_TEXT.....No Help Text  
ID.....{E5B432BD-C698-4415-98E1-F5F229F54F8A}  
TERM.....SINGLE LUMEN CATHETER

METHOD: SKIN (ADHESIVE INDICATOR)

-----  
ABBREVIATION.....SKIN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Skin (Adhesive Indicator)  
HELP\_TEXT.....No Help Text  
ID.....{C7CE7D84-6FC7-406A-7E56-B080637509DF}  
TERM.....SKIN (ADHESIVE INDICATOR)

METHOD: SM ADULT

-----  
ABBREVIATION.....SmA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Sm Adult  
HELP\_TEXT.....No Help Text  
ID.....{CCB94C2F-CE55-A046-3A21-3ADDB1D23D3C}  
TERM.....SM ADULT

METHOD: SM ADULT CUFF

-----  
ABBREVIATION.....SAC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Sm Adult Cuff  
HELP\_TEXT.....No Help Text  
ID.....{0346341B-6D38-6599-ED9F-24C81A8698A8}  
TERM.....SM ADULT CUFF

METHOD: SPONTANEOUS

-----  
ABBREVIATION.....S  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Spontaneous  
HELP\_TEXT.....No Help Text  
ID.....{5E9463BC-C871-54FD-F800-DCFEDB8A93B}  
TERM.....SPONTANEOUS

METHOD: STANDING WEIGHT

-----  
ABBREVIATION.....SW  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Standing Weight  
HELP\_TEXT.....No Help Text  
ID.....{BE16FA09-16DA-510B-D235-B7F0F93AB105}  
TERM.....STANDING WEIGHT

METHOD: T-PIECE

-----  
ABBREVIATION.....TP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....T-Piece  
HELP\_TEXT.....No Help Text  
ID.....{27FC4BC6-2F25-2154-553C-2B452E28A91C}  
TERM.....T-PIECE

METHOD: THIGH CUFF

-----  
ABBREVIATION.....ThC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Thigh Cuff  
HELP\_TEXT.....No Help Text  
ID.....{C581DFAC-691C-B64A-F058-AF2AA5050147}  
TERM.....THIGH CUFF

METHOD: TPN

-----  
ABBREVIATION.....TPN  
ACTIVE.....Yes  
DESCRIPTION.....Total Parenteral Nutrition Method  
DISPLAY\_NAME.....TPN - Method  
HELP\_TEXT.....No Help Text  
ID.....{B0BE1126-1C5B-4B7E-9B82-5A58AAF63FD2}  
TERM.....TPN

METHOD: TRACHEOSTOMY COLLAR

-----  
ABBREVIATION.....TC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Tracheostomy Collar  
HELP\_TEXT.....No Help Text  
ID.....{CBF8A3B9-D323-E021-FD1B-BF4CE19AD20E}  
TERM.....TRACHEOSTOMY COLLAR

METHOD: TRANSFUSION

-----  
ABBREVIATION.....TRANSFUSIO  
ACTIVE.....Yes  
DESCRIPTION.....Transfusion - Method  
DISPLAY\_NAME.....Transfusion - Method  
HELP\_TEXT.....No Help Text  
ID.....{9CBA9E7F-C0E2-4936-8AD1-04E2BD1D770F}  
TERM.....TRANSFUSION

METHOD: TRANSTRACHEAL

-----  
ABBREVIATION.....Tt  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Transtracheal  
HELP\_TEXT.....No Help Text  
ID.....{48951129-CEEC-10D9-F746-99B4713B9C64}  
TERM.....TRANSTRACHEAL

METHOD: TRIPLE LUMEN CATHETER

-----  
ABBREVIATION.....TRIPLE LUM  
ACTIVE.....Yes  
DESCRIPTION.....Triple Lumen Catheter  
DISPLAY\_NAME.....Triple Lumen Catheter  
HELP\_TEXT.....No Help Text  
ID.....{704E2B3E-5E5D-43FB-B956-CFF90AF7D774}  
TERM.....TRIPLE LUMEN CATHETER

METHOD: VENTILATOR

-----  
ABBREVIATION.....Ven  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Ventilator  
HELP\_TEXT.....No Help Text  
ID.....{416FCBD3-5370-15B9-C480-48855ED486D9}  
TERM.....VENTILATOR

METHOD: VENTRICULAR DRAIN

-----  
ABBREVIATION.....VENT DRAIN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Ventricular Drain  
HELP\_TEXT.....No Help Text  
ID.....{5FF035D1-9C97-415D-BBEB-943506395A0C}  
TERM.....VENTRICULAR DRAIN

METHOD: VENTURI MASK

-----  
ABBREVIATION.....VM  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Venturi Mask  
HELP\_TEXT.....No Help Text  
ID.....{E0401675-9C4E-2AFA-2B93-F352B781E963}  
TERM.....VENTURI MASK

METHOD: WHEELCHAIR SCALE

-----  
ABBREVIATION.....WcS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Wheelchair Scale  
HELP\_TEXT.....No Help Text  
ID.....{8C496753-9F42-2A36-E5CA-750188306ED4}  
TERM.....WHEELCHAIR SCALE

METHOD: WITH ACTIVITY

-----  
ABBREVIATION.....Act  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....With Activity  
HELP\_TEXT.....No Help Text  
ID.....{544A262B-EA34-F82A-46A9-F80740515D65}  
TERM.....WITH ACTIVITY

METHOD: WITH AMBULATION

-----  
ABBREVIATION.....WA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....With Ambulation  
HELP\_TEXT.....No Help Text  
ID.....{4E91DE6A-359E-5AF4-AB49-AFAF735E1B94}  
TERM.....WITH AMBULATION

METHOD: WITH CAST OR BRACE

-----  
ABBREVIATION.....WCB  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....With Cast or Brace  
HELP\_TEXT.....No Help Text  
ID.....{47709C05-DF21-FB7D-DE7F-0699D1B7E457}  
TERM.....WITH CAST OR BRACE

METHOD: WITH PROSTHESIS

-----  
ABBREVIATION.....WP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....With Prosthesis  
HELP\_TEXT.....No Help Text  
ID.....{73C93D2D-934B-CA55-3557-E756A89710F5}  
TERM.....WITH PROSTHESIS

METHOD: WITHOUT PROSTHESIS

-----  
ABBREVIATION.....WOP  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Without Prosthesis  
HELP\_TEXT.....No Help Text  
ID.....{9AACA852-E70D-1048-FCE7-112D46C12619}  
TERM.....WITHOUT PROSTHESIS





NULL FLAVOR: PATIENT ON PASS

-----

ABBREVIATION.....op  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....On Pass  
HELP\_TEXT.....No Help Text  
ID.....{3D2D2E52-D3F2-E02D-B9C7-CEB7C46C807D}  
TERM.....PATIENT ON PASS

NULL FLAVOR: PATIENT REFUSED

-----

ABBREVIATION.....pr  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Refused  
HELP\_TEXT.....No Help Text  
ID.....{E9676179-B942-3851-37AF-75DAEE9E1F8E}  
TERM.....PATIENT REFUSED

POSITION: L MIDDLE POSTERIOR CHEST

-----  
ABBREVIATION.....L M P CHES  
ACTIVE.....No  
DESCRIPTION.....Left Middle Posterior Chest  
DISPLAY\_NAME.....L Mid. Pos. Chest  
HELP\_TEXT.....No Help Text  
ID.....{D6D6483C-90A9-4A32-A101-9DAA13BFD7ED}  
TERM.....L MIDDLE POSTERIOR CHEST

POSITION: LYING

-----  
ABBREVIATION.....Ly  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Lying  
HELP\_TEXT.....No Help Text  
ID.....{DE28D2BF-6F8D-64C3-0092-4CEB5E0902B9}  
TERM.....LYING

POSITION: RES EVAL: BED/WHEELCHAIR ALARM

-----  
ABBREVIATION.....R.E. ALAR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Bed / Wheelchair Alarm  
HELP\_TEXT.....No Help Text  
ID.....{BDF926A9-2BEF-4E7C-A099-D02BD98CA2F2}  
TERM.....RES EVAL: BED/WHEELCHAIR ALARM

POSITION: RES EVAL: REALITY ORIENTATION

-----  
ABBREVIATION.....R.E. REALI  
ACTIVE.....No  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Reality Orientation  
HELP\_TEXT.....No Help Text  
ID.....{04D911E1-AADC-4125-9F93-B298266BA61A}  
TERM.....RES EVAL: REALITY ORIENTATION

POSITION: SEMIFOWLERS

-----  
ABBREVIATION.....Sfl  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Semifowlers  
HELP\_TEXT.....No Help Text  
ID.....{C2F6A945-3829-AABA-3D67-54D1829D3EB9}  
TERM.....SEMIFOWLERS

POSITION: SITTING

-----  
ABBREVIATION.....Sit  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Sitting  
HELP\_TEXT.....No Help Text  
ID.....{DD4DCDF1-0A22-2CC5-404F-0E472C73961C}  
TERM.....SITTING

POSITION: STANDING

-----  
ABBREVIATION.....St  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Standing  
HELP\_TEXT.....No Help Text  
ID.....{0776AE42-C370-B7B3-9B4D-9C914BA8F231}  
TERM.....STANDING

POSITION: SUPINE

-----  
ABBREVIATION.....SUPINE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Supine  
HELP\_TEXT.....No Help Text  
ID.....{DEA8A744-B445-4A9E-875D-485E060F2ACC}  
TERM.....SUPINE

POSITION: TRENDELENBURG

-----  
ABBREVIATION.....Tr  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Trendelenburg  
HELP\_TEXT.....No Help Text  
ID.....{42F3353B-7041-EEE1-A87E-941D3BE1CBCD}  
TERM.....TRENDELENBURG

POSITION: WHEELCHAIR EXIT

-----  
ABBREVIATION.....WHEEL EXIT  
ACTIVE.....No  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Wheelchair Exit  
HELP\_TEXT.....No Help Text  
ID.....{DE13F84A-168F-490D-9AC7-840DAB2D8E22}  
TERM.....WHEELCHAIR EXIT



PRODUCT: 25 PERCENT ALBUMIN

-----

ABBREVIATION.....25 PERC AL  
ACTIVE.....Yes  
DESCRIPTION.....25 Percent Albumin  
DISPLAY\_NAME.....25 Percent Albumin  
HELP\_TEXT.....No Help Text  
ID.....{CBC7E32F-E457-4A7F-8F27-030EBC6ECBD5}  
TERM.....25 PERCENT ALBUMIN

PRODUCT: 5 PERCENT ALBUMIN

-----

ABBREVIATION.....5 PERC ALB  
ACTIVE.....Yes  
DESCRIPTION.....5 Percent Albumin  
DISPLAY\_NAME.....5 Percent Albumin  
HELP\_TEXT.....No Help Text  
ID.....{A2D0F9A0-A9EC-4540-9663-B8B19A189198}  
TERM.....5 PERCENT ALBUMIN

PRODUCT: 6 PERCENT HETASTARCH

-----

ABBREVIATION.....6 PERC HET  
ACTIVE.....Yes  
DESCRIPTION.....6 Percent Hetastarch  
DISPLAY\_NAME.....6 Percent Hetastarch  
HELP\_TEXT.....No Help Text  
ID.....{8A82F108-D2FE-4C21-93BB-344502235AE8}  
TERM.....6 PERCENT HETASTARCH

PRODUCT: ALBUMIN

-----

ABBREVIATION.....ALBUMIN  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Albumin  
HELP\_TEXT.....No Help Text  
ID.....{EF4CC3AE-32B2-4EBC-BBDE-7DE7CCEAF081}  
TERM.....ALBUMIN

PRODUCT: AMPHOTERICIN B

-----

ABBREVIATION.....AMPHOTER B  
ACTIVE.....Yes  
DESCRIPTION.....Amphotericin B  
DISPLAY\_NAME.....Amphotericin B  
HELP\_TEXT.....Amphotericin B  
ID.....{E4ED1488-6A27-4622-A8EE-AB6B53DB86D0}  
TERM.....AMPHOTERICIN B

PRODUCT: ANTIBIOTIC

-----  
ABBREVIATION.....AB  
ACTIVE.....Yes  
DESCRIPTION.....Antibiotic  
DISPLAY\_NAME.....Antibiotic  
HELP\_TEXT.....No Help Text  
ID.....{BC1045A3-FF58-41C0-B48F-C2CFD6E522A7}  
TERM.....ANTIBIOTIC

PRODUCT: BACTERIOSTATIC WATER

-----  
ABBREVIATION.....BACTER H2O  
ACTIVE.....Yes  
DESCRIPTION.....Bacteriostatic Water  
DISPLAY\_NAME.....Bacteriostatic Water  
HELP\_TEXT.....Bacteriostatic Water  
ID.....{EF13BCA6-3AF2-47E9-87A0-82642B529368}  
TERM.....BACTERIOSTATIC WATER

PRODUCT: BLADDER IRRIGATION

-----  
ABBREVIATION.....BLAD IRRIG  
ACTIVE.....Yes  
DESCRIPTION.....Bladder Irrigation Output Product  
DISPLAY\_NAME.....Bladder Irrigation  
HELP\_TEXT.....No Help Text  
ID.....{BCBB106B-ECE4-4232-87C7-D90EDB47F80D}  
TERM.....BLADDER IRRIGATION

PRODUCT: BLOOD PRODUCT

-----  
ABBREVIATION.....bp  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Blood Product  
HELP\_TEXT.....No Help Text  
ID.....{A910DBA4-1CD7-1302-40B2-5BEE95DD3B08}  
TERM.....BLOOD PRODUCT

PRODUCT: CARDIAC OUTPUT INJECTATE

-----  
ABBREVIATION.....CARD OUT I  
ACTIVE.....Yes  
DESCRIPTION.....Cardiac Output Injectate  
DISPLAY\_NAME.....Cardiac Output Injectate  
HELP\_TEXT.....No Help Text  
ID.....{2128B705-DDF4-4776-87C3-AA82F854E313}  
TERM.....CARDIAC OUTPUT INJECTATE

PRODUCT: CHEST TUBE

-----  
ABBREVIATION.....CHEST TUBE  
ACTIVE.....Yes  
DESCRIPTION.....Chest Tube Output  
DISPLAY\_NAME.....Chest Tube  
HELP\_TEXT.....No Help Text  
ID.....{74C0C672-02AD-45A8-99DB-F5256BA14C72}  
TERM.....CHEST TUBE

PRODUCT: CHEST TUBE DRAINAGE

-----  
ABBREVIATION.....CHEST T DR  
ACTIVE.....Yes  
DESCRIPTION.....Chest Tube Drainage  
DISPLAY\_NAME.....Chest Tube Drainage  
HELP\_TEXT.....No Help Text  
ID.....{591D62C9-D173-45AD-8E05-A500F1C4AAA7}  
TERM.....CHEST TUBE DRAINAGE

PRODUCT: CLEAR LIQUIDS

-----  
ABBREVIATION.....Clear Liq  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Clear Liquids  
HELP\_TEXT.....  
ID.....{BBE3A652-ACDC-4797-B140-D4DB4F474B13}  
TERM.....CLEAR LIQUIDS

PRODUCT: COFFEE

-----  
ABBREVIATION.....Coffee  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Coffee  
HELP\_TEXT.....  
ID.....{BEFE282A-F07F-4B3A-A1FC-5FBC0DE75AE1}  
TERM.....COFFEE

PRODUCT: COLLOID

-----  
ABBREVIATION.....COLLOID  
ACTIVE.....Yes  
DESCRIPTION.....Colloid  
DISPLAY\_NAME.....Colloid  
HELP\_TEXT.....No Help Text  
ID.....{5B912A05-9ECA-43C1-A117-BE0843313CDD}  
TERM.....COLLOID



PRODUCT: CRYOPRECIPITATE

-----  
ABBREVIATION.....CRYOPRECIP  
ACTIVE.....Yes  
DESCRIPTION.....Cryoprecipitate Blood Product  
DISPLAY\_NAME.....Cryoprecipitate  
HELP\_TEXT.....No Help Text  
ID.....{5329E737-E059-480D-83F8-1146B091C14E}  
TERM.....CRYOPRECIPITATE

PRODUCT: CRYSTALLOID

-----  
ABBREVIATION.....CRYSTALLOI  
ACTIVE.....Yes  
DESCRIPTION.....Crystalloid  
DISPLAY\_NAME.....Crystalloid  
HELP\_TEXT.....No Help Text  
ID.....{F8B0C9B5-F386-4C88-896A-092C4A95A031}  
TERM.....CRYSTALLOID

PRODUCT: DAKINS SOLUTION

-----  
ABBREVIATION.....DAKINS SOL  
ACTIVE.....Yes  
DESCRIPTION.....Dakins Solution  
DISPLAY\_NAME.....Dakins Solution  
HELP\_TEXT.....Dakins Solution  
ID.....{51651D24-B909-427B-9FF0-3C2859A63D73}  
TERM.....DAKINS SOLUTION

PRODUCT: EF RESIDUAL

-----  
ABBREVIATION.....EF RESID  
ACTIVE.....Yes  
DESCRIPTION.....Enteral Feeding Residual  
DISPLAY\_NAME.....EF Residual  
HELP\_TEXT.....No Help Text  
ID.....{35E269DF-F246-4DF6-9D68-7FD61F8FE088}  
TERM.....EF RESIDUAL

PRODUCT: ELEMENTAL FORMULA

-----  
ABBREVIATION.....ELE FORM  
ACTIVE.....Yes  
DESCRIPTION.....Elemental Enteral Nutrition Formula  
DISPLAY\_NAME.....Elemental Formula  
HELP\_TEXT.....No Help Text  
ID.....{E85DFA33-7F32-4348-9E0E-249AC301EC03}  
TERM.....ELEMENTAL FORMULA

PRODUCT: EMESIS

-----  
ABBREVIATION.....EMESIS  
ACTIVE.....Yes  
DESCRIPTION.....Emesis output  
DISPLAY\_NAME.....Emesis  
HELP\_TEXT.....No Help Text  
ID.....{A7B1605E-1DFE-4FC6-A69B-0E18C1A291FF}  
TERM.....EMESIS

PRODUCT: FECES

-----  
ABBREVIATION.....f  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Feces  
HELP\_TEXT.....No Help Text  
ID.....{C1CED7CD-EF5D-37D8-4EE0-CC9FD6F20DB0}  
TERM.....FECES

PRODUCT: FFP

-----  
ABBREVIATION.....FFP  
ACTIVE.....Yes  
DESCRIPTION.....Fresh Frozen Plasma  
DISPLAY\_NAME.....FFP  
HELP\_TEXT.....No Help Text  
ID.....{9DD0ECF5-634F-49CF-A529-A2657240A80C}  
TERM.....FFP

PRODUCT: FFP - THAWED

-----  
ABBREVIATION.....FFP-THAWED  
ACTIVE.....Yes  
DESCRIPTION.....Fresh Frozen Plasma, Thawed  
DISPLAY\_NAME.....FFP - THAWED  
HELP\_TEXT.....No Help Text  
ID.....{AAF452A3-0DCA-4FF4-BF96-08E60D2D2875}  
TERM.....FFP - THAWED

PRODUCT: FORMULA - HIGH FIBER

-----  
ABBREVIATION.....FORM H FIB  
ACTIVE.....Yes  
DESCRIPTION.....High Fiber Formula  
DISPLAY\_NAME.....High Fiber Formula  
HELP\_TEXT.....No Help Text  
ID.....{249D8C78-B600-4443-A2E2-39F1B67698DD}  
TERM.....FORMULA - HIGH FIBER

PRODUCT: FORMULA - PROBALANCE

-----  
ABBREVIATION.....PROBALANCE  
ACTIVE.....Yes  
DESCRIPTION.....Probalance Formula  
DISPLAY\_NAME.....Probalance Formula  
HELP\_TEXT.....No Help Text  
ID.....{4EA729FE-08EA-490B-822D-696061E5B55F}  
TERM.....FORMULA - PROBALANCE

PRODUCT: FORMULA 1.5 CAL

-----  
ABBREVIATION.....1.5 C FORM  
ACTIVE.....Yes  
DESCRIPTION.....1.5 Cal Enteral Nutrition Formula  
DISPLAY\_NAME.....1.5 Cal Formula  
HELP\_TEXT.....No Help Text  
ID.....{95FAE60B-BEAE-4EB2-88FC-0A199F4534EB}  
TERM.....FORMULA 1.5 CAL

PRODUCT: FORMULA 2.0 CAL

-----  
ABBREVIATION.....2.0 C FORM  
ACTIVE.....Yes  
DESCRIPTION.....2.0 Cal Enteral Nutrition Formula  
DISPLAY\_NAME.....2.0 Cal Formula  
HELP\_TEXT.....No Help Text  
ID.....{E5FB9590-25AA-435A-B9B4-4812EA15C217}  
TERM.....FORMULA 2.0 CAL

PRODUCT: FORMULA D10W ENTERAL NUTRITION

-----  
ABBREVIATION.....D10W ENT N  
ACTIVE.....Yes  
DESCRIPTION.....D10W Enteral Nutrition  
DISPLAY\_NAME.....D10W Enteral Nutrition  
HELP\_TEXT.....10% Dextrose and Water solution  
ID.....{F8FB55C4-6229-427F-987D-3D469E307D0D}  
TERM.....FORMULA D10W ENTERAL NUTRITION

PRODUCT: FORMULA D50W ENTERAL NUTRITION

-----  
ABBREVIATION.....D50W ENT N  
ACTIVE.....Yes  
DESCRIPTION.....D50W Enteral Nutrition  
DISPLAY\_NAME.....D50W Enteral Nutrition  
HELP\_TEXT.....50% Dextrose and Water Solution  
ID.....{19ED9C83-42C1-42AF-AB4D-1A34B621EFB9}  
TERM.....FORMULA D50W ENTERAL NUTRITION

PRODUCT: FORMULA D5W ENTERAL NUTRITION

-----  
ABBREVIATION.....D5W ENT NU  
ACTIVE.....Yes  
DESCRIPTION.....D5W Enteral Nutrition  
DISPLAY\_NAME.....D5W Enteral Nutrition  
HELP\_TEXT.....5% Dextrose and Water Solution  
ID.....{9EA15A1B-DA1C-49EC-B463-82DD5F125F82}  
TERM.....FORMULA D5W ENTERAL NUTRITION

PRODUCT: FORMULA FIBER

-----  
ABBREVIATION.....FIBER FORM  
ACTIVE.....Yes  
DESCRIPTION.....Fiber Enteral Nutrition Formula  
DISPLAY\_NAME.....Fiber Enteral Nutrition Formula  
HELP\_TEXT.....No Help Text  
ID.....{F6D7E75F-917A-4C50-AABB-BE9DE3C1939D}  
TERM.....FORMULA FIBER

PRODUCT: FULL LIQUID

-----  
ABBREVIATION.....FL LIQUID  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Full Liquid  
HELP\_TEXT.....No Help Text  
ID.....{C8405270-CCA3-4BFF-B9D0-398590FBB1D9}  
TERM.....FULL LIQUID

PRODUCT: GASTRIC DRAINAGE

-----  
ABBREVIATION.....GAST DRAIN  
ACTIVE.....Yes  
DESCRIPTION.....Gastric Drainage  
DISPLAY\_NAME.....Gastric Drainage  
HELP\_TEXT.....No Help Text  
ID.....{C01D18F9-191A-4D34-B45D-D8D70ACA67BE}  
TERM.....GASTRIC DRAINAGE

PRODUCT: GLUCERNA DIABETIC

-----  
ABBREVIATION.....GLUCERNA D  
ACTIVE.....Yes  
DESCRIPTION.....Glucerna Diabetic  
DISPLAY\_NAME.....Glucerna Diabetic  
HELP\_TEXT.....No Help Text  
ID.....{C2F2CCFE-6D77-49B4-BF5E-F1A285326440}  
TERM.....GLUCERNA DIABETIC

PRODUCT: GRANULOCYTES PHERESIS

-----  
ABBREVIATION.....GRAN PHERE  
ACTIVE.....Yes  
DESCRIPTION.....Granulocytes Pheresis Blood Product  
DISPLAY\_NAME.....Granulocytes Pheresis  
HELP\_TEXT.....No Help Text  
ID.....{2B90878D-9C29-4731-9525-D7DD7FCB513B}  
TERM.....GRANULOCYTES PHERESIS

PRODUCT: HEME

-----  
ABBREVIATION.....HEME  
ACTIVE.....Yes  
DESCRIPTION.....Heme output  
DISPLAY\_NAME.....Heme  
HELP\_TEXT.....No Help Text  
ID.....{A0ADA229-E779-4562-84B5-67121B88C443}  
TERM.....HEME

PRODUCT: HEMODIALYSIS

-----  
ABBREVIATION.....HEMODIALYS  
ACTIVE.....Yes  
DESCRIPTION.....Hemodialysis Intake and Output Product  
DISPLAY\_NAME.....Hemodialysis  
HELP\_TEXT.....No Help Text  
ID.....{8409B357-C41E-4F4A-BB52-F09AD9704C30}  
TERM.....HEMODIALYSIS

PRODUCT: HEPATIC FORMULA

-----  
ABBREVIATION.....HEPAT FORM  
ACTIVE.....Yes  
DESCRIPTION.....Hepatic Enteral Nutrition Formula  
DISPLAY\_NAME.....Hepatic Enteral Nutrition Formula  
HELP\_TEXT.....No Help Text  
ID.....{41DEECAB-168F-4522-852C-C1ED751BFEDC}  
TERM.....HEPATIC FORMULA

PRODUCT: HIGH PROTEIN FORMULA

-----  
ABBREVIATION.....HIGH PROTE  
ACTIVE.....Yes  
DESCRIPTION.....High Protein  
DISPLAY\_NAME.....High Protein Formula  
HELP\_TEXT.....No Help Text  
ID.....{3618C874-3EAE-41CF-BE2B-3DFB960E034F}  
TERM.....HIGH PROTEIN FORMULA

PRODUCT: JP DRAIN

-----  
ABBREVIATION.....JP DRAIN  
ACTIVE.....Yes  
DESCRIPTION.....Jackson-Pratt Drain  
DISPLAY\_NAME.....JP DRAIN  
HELP\_TEXT.....No Help Text  
ID.....{0889246D-AD40-4CFC-A721-BBCB8C34BE19}  
TERM.....JP DRAIN

PRODUCT: JUICE

-----  
ABBREVIATION.....JUICE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Juice  
HELP\_TEXT.....No Help-Text  
ID.....{91FFF519-1677-BB1D-D959-65B1BFA10285}  
TERM.....JUICE

PRODUCT: LINE FLUSH

-----  
ABBREVIATION.....LINE FLUSH  
ACTIVE.....Yes  
DESCRIPTION.....Fluid flushed from line  
DISPLAY\_NAME.....LINE FLUSH  
HELP\_TEXT.....No Help Text  
ID.....{D3EA8BCD-0EC4-4B39-BA92-BC3FF505C9A4}  
TERM.....LINE FLUSH

PRODUCT: LINE FLUSH HEPARIN

-----  
ABBREVIATION.....LN FL HEPA  
ACTIVE.....Yes  
DESCRIPTION.....Heparin used to flush dialysis lines  
DISPLAY\_NAME.....Line Flush Heparin  
HELP\_TEXT.....No Help Text  
ID.....{E8B68E5A-6B9A-401F-9D1A-E3C907C7B81C}  
TERM.....LINE FLUSH HEPARIN

PRODUCT: LIPIDS

-----  
ABBREVIATION.....LIPIDS  
ACTIVE.....Yes  
DESCRIPTION.....Lipids  
DISPLAY\_NAME.....Lipids  
HELP\_TEXT.....No Help Text  
ID.....{80E456A6-CA3D-4C2D-8DCD-06C63AE1F320}  
TERM.....LIPIDS

PRODUCT: MAGNACAL

-----  
ABBREVIATION.....MAGNACAL  
ACTIVE.....Yes  
DESCRIPTION.....Magnacal Enteral Nutrition Type  
DISPLAY\_NAME.....Maganca  
HELP\_TEXT.....No Help Text  
ID.....{86688633-DFCA-4FBE-99F5-8E0B3FAFE695}  
TERM.....MAGNACAL

PRODUCT: MEDIASTINAL

-----  
ABBREVIATION.....MEDIASTINA  
ACTIVE.....Yes  
DESCRIPTION.....Mediastinal Output Product  
DISPLAY\_NAME.....Mediastinal  
HELP\_TEXT.....No Help Text  
ID.....{A0550B31-B937-4FC7-8F7E-071F75215311}  
TERM.....MEDIASTINAL

PRODUCT: MEDICATION

-----  
ABBREVIATION.....MEDICATION  
ACTIVE.....Yes  
DESCRIPTION.....Medication  
DISPLAY\_NAME.....Medication  
HELP\_TEXT.....Medication  
ID.....{CDD5944E-A0EC-45CE-8EB0-7F9EDF7197E5}  
TERM.....MEDICATION

PRODUCT: MILK

-----  
ABBREVIATION.....Milk  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Milk  
HELP\_TEXT.....  
ID.....{6C4DF9C1-D7BA-429A-BDAF-00FBF35F6B09}  
TERM.....MILK

PRODUCT: MODULAR CUSTOMIZED FORMULA

-----  
ABBREVIATION.....CUST FORMU  
ACTIVE.....Yes  
DESCRIPTION.....Modular Customized Formula  
DISPLAY\_NAME.....Modular Customized Formula  
HELP\_TEXT.....No Help Text  
ID.....{249D095E-9676-4CC5-8A39-097F96A4C032}  
TERM.....MODULAR CUSTOMIZED FORMULA

PRODUCT: NEPRO RENAL FORMULA

-----  
ABBREVIATION.....NEPR RENAL  
ACTIVE.....Yes  
DESCRIPTION.....Nepro Renal Enteral Nutrition Formula  
DISPLAY\_NAME.....Nepro Renal Enteral Nutrition Formula  
HELP\_TEXT.....No Help Text  
ID.....{59C9B396-A705-445B-B566-F030DA4B2B79}  
TERM.....NEPRO RENAL FORMULA

PRODUCT: NG RESIDUAL

-----  
ABBREVIATION.....NG RESID  
ACTIVE.....Yes  
DESCRIPTION.....Naso-Gastric Residual Output  
DISPLAY\_NAME.....NG Residual  
HELP\_TEXT.....No Help Text  
ID.....{B03C9557-84B5-4B05-B784-C87FA90D27A7}  
TERM.....NG RESIDUAL

PRODUCT: NOVASOURCE 2.0

-----  
ABBREVIATION.....NOVAS 2.0  
ACTIVE.....Yes  
DESCRIPTION.....Novasource 2.0 Enteral Nutrition  
Formula  
DISPLAY\_NAME.....Novasource 2.0 Enteral Nutrition  
Formula  
HELP\_TEXT.....No Help Text  
ID.....{67D16273-6271-43E1-80D6-D12C4D55E024}  
TERM.....NOVASOURCE 2.0

PRODUCT: NUTREN 1.5

-----  
ABBREVIATION.....NUTREN 1.5  
ACTIVE.....Yes  
DESCRIPTION.....Nutren 1.5  
DISPLAY\_NAME.....Nutren 1.5  
HELP\_TEXT.....No Help Text  
ID.....{B93EF897-B9A1-45C6-AC9A-B386554A3186}  
TERM.....NUTREN 1.5

PRODUCT: NUTRITIONAL SUPPLEMENT

-----  
ABBREVIATION.....NUTR SUPPL  
ACTIVE.....Yes  
DESCRIPTION.....Nutritional Supplement  
DISPLAY\_NAME.....Nutritional Supplement  
HELP\_TEXT.....No Help Text  
ID.....{17F8F1DA-4FC6-48E5-9E6E-1F67CD06D200}  
TERM.....NUTRITIONAL SUPPLEMENT



PRODUCT: OSMOLITE

-----  
ABBREVIATION.....OSMOLITE  
ACTIVE.....Yes  
DESCRIPTION.....Osmolite Enteral Nutrition  
DISPLAY\_NAME.....Osmolite  
HELP\_TEXT.....No Help Text  
ID.....{18C5843F-BC33-4609-9D58-49F07F12A059}  
TERM.....OSMOLITE

PRODUCT: OSMOLITE HN

-----  
ABBREVIATION.....OSMOLIT HN  
ACTIVE.....Yes  
DESCRIPTION.....Osmolite HN Enteral Nutrition  
DISPLAY\_NAME.....Osmolite HN  
HELP\_TEXT.....No Help Text  
ID.....{3622D5C6-1D7B-4377-94B7-5E33C642C4C5}  
TERM.....OSMOLITE HN

PRODUCT: OTHER PO FLUID

-----  
ABBREVIATION.....PO OTHER  
ACTIVE.....Yes  
DESCRIPTION.....PO Fluid other than those listed.  
DISPLAY\_NAME.....Other - Please Specify  
HELP\_TEXT.....No Help Text  
ID.....{1EA44424-4DE2-42B8-A50F-60DC56204502}  
TERM.....OTHER PO FLUID

PRODUCT: PCA

-----  
ABBREVIATION.....PCA  
ACTIVE.....Yes  
DESCRIPTION.....Patient Controlled Analgesia  
DISPLAY\_NAME.....Patient Controlled Analgesia  
HELP\_TEXT.....No Help Text  
ID.....{D93028EE-9A44-43F5-9C0C-A903BD9B6E6D}  
TERM.....PCA

PRODUCT: PLASMAPHERESIS

-----  
ABBREVIATION.....PLASMAPHER  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Plasmapheresis  
HELP\_TEXT.....No Help Text  
ID.....{E1214367-8BB0-4FAA-A916-EC707569A9C4}  
TERM.....PLASMAPHERESIS

PRODUCT: PLATELETS

-----  
ABBREVIATION.....PLATELETS  
ACTIVE.....Yes  
DESCRIPTION.....Platelets  
DISPLAY\_NAME.....Platelets  
HELP\_TEXT.....No Help Text  
ID.....{3F4FA2EF-DF05-4A6D-839F-F97D3C2AFCF7}  
TERM.....PLATELETS

PRODUCT: PLATELETS PHERESIS

-----  
ABBREVIATION.....PLATELET P  
ACTIVE.....Yes  
DESCRIPTION.....Platelets Pheresis (Apheresis  
Platelets)  
DISPLAY\_NAME.....Platelets Pheresis  
HELP\_TEXT.....No Help Text  
ID.....{65CBDCBD-99DB-44FE-9230-D62886208EBF}  
TERM.....PLATELETS PHERESIS

PRODUCT: PLATELETS PHERESIS IRRADIATED

-----  
ABBREVIATION.....PLATEL P I  
ACTIVE.....Yes  
DESCRIPTION.....Platelets Pheresis (Apheresis  
platelets) Irradiated  
DISPLAY\_NAME.....Platelets Pheresis Irradiated  
HELP\_TEXT.....No Help Text  
ID.....{70C7155A-5AB6-4733-B43F-606F8BBB3B4A}  
TERM.....PLATELETS PHERESIS IRRADIATED

PRODUCT: PPN

-----  
ABBREVIATION.....PPN  
ACTIVE.....Yes  
DESCRIPTION.....Peripheral Parenteral Nutrition  
DISPLAY\_NAME.....PPN  
HELP\_TEXT.....No Help Text  
ID.....{BF402283-3E6F-4DF1-9CC6-4DE35D755C7C}  
TERM.....PPN

PRODUCT: PRBC

-----  
ABBREVIATION.....PRBC  
ACTIVE.....Yes  
DESCRIPTION.....Packed Red Blood Cells  
DISPLAY\_NAME.....PRBC  
HELP\_TEXT.....No Help Text  
ID.....{BCE37EFD-78D6-4A13-A5E4-EA278BF18E02}  
TERM.....PRBC

PRODUCT: PRBC - AUTOLOGOUS

-----  
ABBREVIATION.....PRBC-A  
ACTIVE.....Yes  
DESCRIPTION.....Packed Red Blood Cells, Autologous  
DISPLAY\_NAME.....PRBC - AUTOLOGOUS  
HELP\_TEXT.....No Help Text  
ID.....{1580990C-0B35-4DDD-84AA-682B341C9B76}  
TERM.....PRBC - AUTOLOGOUS

PRODUCT: PRBC - LEUKOREduced

-----  
ABBREVIATION.....PRBC-L  
ACTIVE.....Yes  
DESCRIPTION.....Packed Red Blood Cells, Leukoreduced  
DISPLAY\_NAME.....PRBC - Leukoreduced  
HELP\_TEXT.....No Help Text  
ID.....{13DB3170-6C26-4178-B429-C02AFC541AC7}  
TERM.....PRBC - LEUKOREduced

PRODUCT: PRBC - LEUKOREduced - DIRECTED

-----  
ABBREVIATION.....PRBC-L-D  
ACTIVE.....Yes  
DESCRIPTION.....Packed Red Blood Cells, Leukoreduced,  
Directed  
DISPLAY\_NAME.....PRBC - LEUKOREduced - DIRECTED  
HELP\_TEXT.....No Help Text  
ID.....{9C54B7F7-FBD9-4FA8-B2A1-893A25A9CD2C}  
TERM.....PRBC - LEUKOREduced - DIRECTED

PRODUCT: PRBC - LEUKOREduced - IRRADIATED

-----  
ABBREVIATION.....PRBC-L-I  
ACTIVE.....Yes  
DESCRIPTION.....Packed Red Blood Cells, Leukoreduced,  
Irradiated  
DISPLAY\_NAME.....PRBC - LEUKOREduced - IRRADIATED  
HELP\_TEXT.....No Help Text  
ID.....{B4431519-0102-49ED-94D6-1738E375C5A4}  
TERM.....PRBC - LEUKOREduced - IRRADIATED

PRODUCT: PRBC - LEUKOREduced - WASHED

-----  
ABBREVIATION.....PRBC-L-W  
ACTIVE.....Yes  
DESCRIPTION.....Packed Red Blood Cells, Leukoreduced,  
Washed  
DISPLAY\_NAME.....PRBC - LEUKOREduced - WASHED  
HELP\_TEXT.....No Help Text  
ID.....{BB58F35A-1086-4061-B485-FA35DFB79A0D}  
TERM.....PRBC - LEUKOREduced - WASHED

PRODUCT: Rh IMMUNE GLOBULIN

-----  
ABBREVIATION.....Rh IM GLOB  
ACTIVE.....Yes  
DESCRIPTION.....Rh Immune Globulin  
DISPLAY\_NAME.....Rh Immune Globulin  
HELP\_TEXT.....No Help Text  
ID.....{1FE7BD0C-7252-427A-8216-36FCA4637089}  
TERM.....Rh IMMUNE GLOBULIN

PRODUCT: SALINE

-----  
ABBREVIATION.....SALINE  
ACTIVE.....Yes  
DESCRIPTION.....Saline  
DISPLAY\_NAME.....Saline  
HELP\_TEXT.....No Help Text  
ID.....{19D84001-70DB-47AA-9149-A37B7E0BAEF1}  
TERM.....SALINE

PRODUCT: SANDOSOURCE PEPTIDE

-----  
ABBREVIATION.....SAND PEPTI  
ACTIVE.....Yes  
DESCRIPTION.....Sandosource Peptide Enteral Nutrition  
Formula  
DISPLAY\_NAME.....Sandosource Peptide Enteral Nutrition  
Formula  
HELP\_TEXT.....No Help Text  
ID.....{B543461C-7316-4750-9918-B30E96AB3C95}  
TERM.....SANDOSOURCE PEPTIDE

PRODUCT: SODA: DIET

-----  
ABBREVIATION.....Soda: Diet  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Soda: Diet  
HELP\_TEXT.....  
ID.....{0EFA0EC7-EB52-4F84-B863-5CD8BC5509DA}  
TERM.....SODA: DIET

PRODUCT: SODA: NON-DIET

-----  
ABBREVIATION.....Soda: Non  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Soda: Non-Diet  
HELP\_TEXT.....  
ID.....{B15D88DC-07B6-4980-98B2-A59B5662797A}  
TERM.....SODA: NON-DIET

PRODUCT: SOFT DRINK

-----  
ABBREVIATION.....SOFT DRIN  
ACTIVE.....Yes  
DESCRIPTION.....Soft Drink  
DISPLAY\_NAME.....Soft Drink  
HELP\_TEXT.....No Help Text  
ID.....{9216F0F9-4D5F-4AF7-8D84-04A2CE11BCBE}  
TERM.....SOFT DRINK

PRODUCT: TPN

-----  
ABBREVIATION.....TPN  
ACTIVE.....Yes  
DESCRIPTION.....Total Parenteral Nutrition  
DISPLAY\_NAME.....TPN  
HELP\_TEXT.....No Help Text  
ID.....{55341B62-E356-4C38-A52D-4F196751E9EB}  
TERM.....TPN

PRODUCT: ULTRACAL

-----  
ABBREVIATION.....ULTRACAL  
ACTIVE.....Yes  
DESCRIPTION.....Ultracal Enteral Nutrition  
DISPLAY\_NAME.....Ultracal Enteral Nutrition  
HELP\_TEXT.....No Help Text  
ID.....{B18E6151-9DC6-4380-899A-E84715BFC61B}  
TERM.....ULTRACAL

PRODUCT: URINE

-----  
ABBREVIATION.....u  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Urine  
HELP\_TEXT.....No Help Text  
ID.....{1F61B077-1D34-0CBA-E55B-65A8B6A7CAE7}  
TERM.....URINE

PRODUCT: VITAL ENTERAL NUTRITION TYPE

-----  
ABBREVIATION.....VIT ENT NU  
ACTIVE.....Yes  
DESCRIPTION.....Vital Enteral Nutrition  
DISPLAY\_NAME.....Vital Enteral Nutrition  
HELP\_TEXT.....No Help Text  
ID.....{D77B607C-85C4-428E-AAE6-B319F33F4E9F}  
TERM.....VITAL ENTERAL NUTRITION TYPE

PRODUCT: VITAL HN ENTERAL NUTRITION

-----  
ABBREVIATION.....VIT NUT HN  
ACTIVE.....Yes  
DESCRIPTION.....Vital HN Enteral Nutrition  
DISPLAY\_NAME.....Vital HN Enteral Nutrition  
HELP\_TEXT.....No Help Text  
ID.....{3973F52F-D53F-4D80-928A-BD6154D5C434}  
TERM.....VITAL HN ENTERAL NUTRITION

PRODUCT: WATER

-----  
ABBREVIATION.....WATER  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Water  
HELP\_TEXT.....No Help-Text  
ID.....{ACCC998E-0F64-510D-B843-3D30B2AD67CF}  
TERM.....WATER

PRODUCT: WBC - WASHED

-----  
ABBREVIATION.....WBC WASHED  
ACTIVE.....Yes  
DESCRIPTION.....White Blood Cells, Washed  
DISPLAY\_NAME.....WBC - Washed  
HELP\_TEXT.....No Help Text  
ID.....{F1460F94-203C-473D-AC7C-92E8ADFFC0F5}  
TERM.....WBC - WASHED

PRODUCT: WHOLE BLOOD

-----  
ABBREVIATION.....WHOL BLOOD  
ACTIVE.....Yes  
DESCRIPTION.....Whole Blood  
DISPLAY\_NAME.....Whole Blood  
HELP\_TEXT.....No Help Text  
ID.....{144EE01A-16B5-4F9B-9A1C-D420D0A9A2D0}  
TERM.....WHOLE BLOOD

PRODUCT: WOUND DRAINAGE

-----  
ABBREVIATION.....WOUN DRAIN  
ACTIVE.....Yes  
DESCRIPTION.....Wound Drainage  
DISPLAY\_NAME.....Wound Drainage  
HELP\_TEXT.....Drainage From Wound  
ID.....{73038CBA-1F31-47E2-A78E-50ED1067CBEA}  
TERM.....WOUND DRAINAGE

PRODUCT: WOUND IRRIGATION

-----  
ABBREVIATION.....WOUND IRRI  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Wound Irrigation  
HELP\_TEXT.....No Help Text  
ID.....{192EC6B7-8E2E-418C-9345-BC8A7986CA59}  
TERM.....WOUND IRRIGATION

PRODUCT: WOUND VACUUM

-----  
ABBREVIATION.....WOUND VAC  
ACTIVE.....Yes  
DESCRIPTION.....Wound Vacuum Output Product  
DISPLAY\_NAME.....Wound Vacuum  
HELP\_TEXT.....No Help Text  
ID.....{9F4983A6-B622-4E4D-A7F7-ADA7A4B7CEBC}  
TERM.....WOUND VACUUM

QUALITY: ACTUAL

-----

ABBREVIATION.....A  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Actual  
HELP\_TEXT.....No Help Text  
ID.....{5842A4D9-56A1-7A5A-13C9-F481699C3A09}  
TERM.....ACTUAL

QUALITY: ANTECUBITAL

-----

ABBREVIATION.....ANTECUB  
ACTIVE.....No  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Antecubital  
HELP\_TEXT.....No Help Text  
ID.....{0D56DDB8-EDCD-BF9F-5F77-567BC59F3580}  
TERM.....ANTECUBITAL

QUALITY: BRACHIAL

-----

ABBREVIATION.....BRACHIAL  
ACTIVE.....No  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Brachial  
HELP\_TEXT.....No Help Text  
ID.....{CFE7E2E3-20BC-7ACE-262B-8BAE4FB36B59}  
TERM.....BRACHIAL

QUALITY: DRY

-----

ABBREVIATION.....D  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Dry  
HELP\_TEXT.....No Help Text  
ID.....{B319EFF4-BA4A-85B1-2617-E00E46A54ABE}  
TERM.....DRY

QUALITY: ESTIMATED

-----

ABBREVIATION.....E  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Estimated  
HELP\_TEXT.....No Help Text  
ID.....{9EAFF26B-7D87-DA66-F991-022D095FB5A0}  
TERM.....ESTIMATED



QUALITY: ESTIMATED BY ARM SPAN

-----  
ABBREVIATION.....EAS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Estimated by Arm Span  
HELP\_TEXT.....No Help Text  
ID.....{486CAE5A-BE56-06BA-C65E-1F8131DB1C2D}  
TERM.....ESTIMATED BY ARM SPAN

QUALITY: OTHER

-----  
ABBREVIATION.....Oth  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Other  
HELP\_TEXT.....No Help Text  
ID.....{7C0850DD-931C-65C8-1003-E31FBBADF6E9}  
TERM.....OTHER

QUALITY: PEDIATRIC

-----  
ABBREVIATION.....Ped  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Pediatric  
HELP\_TEXT.....No Help Text  
ID.....{83C8FF79-1919-B001-9F8B-68F014216F0B}  
TERM.....PEDIATRIC

QUALITY: RADIAL

-----  
ABBREVIATION.....RADIAL  
ACTIVE.....No  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Radial  
HELP\_TEXT.....No Help Text  
ID.....{8F407DAB-8C6D-A5B3-524E-444B22613F4D}  
TERM.....RADIAL

QUALITY: STATED

-----  
ABBREVIATION.....Sta  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Stated  
HELP\_TEXT.....No Help Text  
ID.....{6362C194-B7F2-D38C-5B72-3BE5DE09438E}  
TERM.....STATED



SCHEDULE: CATHETER

-----  
ABBREVIATION.....CATHETER  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Catheter  
HELP\_TEXT.....No Help Text  
ID.....{8C788C2B-3C73-7D8D-9032-A1E2188413B8}  
TERM.....CATHETER

SCHEDULE: PRN

-----  
ABBREVIATION.....PRN  
ACTIVE.....Yes  
DESCRIPTION.....AS NEEDED  
DISPLAY\_NAME.....PRN  
HELP\_TEXT.....No Help Text  
ID.....{B9423F6A-166B-CF51-6375-E3928FE0BAD6}  
TERM.....PRN

SCHEDULE: Q2H

-----  
ABBREVIATION.....q2h  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Every 2 hours  
HELP\_TEXT.....No Help Text  
ID.....{8EF08D16-5529-5FF0-5B84-36703D82E01C}  
TERM.....Q2H

SCHEDULE: Q4H

-----  
ABBREVIATION.....q4  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Every 4 Hours  
HELP\_TEXT.....No Help Text  
ID.....{E5638DC5-2772-FD4A-2D81-04B2DA5C8DDC}  
TERM.....Q4H

SCHEDULE: Q7D

-----  
ABBREVIATION.....q7d  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Every 7 Days  
HELP\_TEXT.....No Help Text  
ID.....{E3E171BF-17B3-1C59-112B-042380890C62}  
TERM.....Q7D

SCHEDULE: QD

-----  
ABBREVIATION.....qd  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Everyday  
HELP\_TEXT.....No Help Text  
ID.....{741F0241-5BC5-5AE2-8D8E-543155C46711}  
TERM.....QD

SCHEDULE: RESTRAINT STRAP

-----  
ABBREVIATION.....RES STRAP  
ACTIVE.....No  
DESCRIPTION.....Restraint Strap  
DISPLAY\_NAME.....Restaint Strap  
HELP\_TEXT.....No Help Text  
ID.....{D9162D47-0F50-49F4-837F-70AD65C169BB}  
TERM.....RESTRAINT STRAP

TASK: ACTIVITY - ABD PILLOW

-----  
ABBREVIATION.....ACT ABD PI  
ACTIVE.....Yes  
DESCRIPTION.....ACTIVITY - ABD PILLOW  
DISPLAY\_NAME.....ACTIVITY - ABD PILLOW  
HELP\_TEXT.....ACTIVITY - ABD PILLOW  
ID.....{381AA3BB-8858-4B2B-B607-8FC6C6FF707F}  
TERM.....ACTIVITY - ABD PILLOW

TASK: ACTIVITY - AMBULATE

-----  
ABBREVIATION.....ACT AMBULA  
ACTIVE.....Yes  
DESCRIPTION.....ACTIVITY - AMBULATE  
DISPLAY\_NAME.....ACTIVITY - AMBULATE  
HELP\_TEXT.....ACTIVITY - AMBULATE  
ID.....{E7232201-7FC2-4900-82D0-EF3A10FCB760}  
TERM.....ACTIVITY - AMBULATE

TASK: ACTIVITY - CONT. PASSIVE MOTION DEV.

-----  
ABBREVIATION.....ACT CPM DE  
ACTIVE.....Yes  
DESCRIPTION.....ACTIVITY - CONTINUOUS PASSIVE MOTION  
DEVICE  
DISPLAY\_NAME.....ACTIVITY - CONT. PASS. MOT. DEV.  
HELP\_TEXT.....ACTIVITY - CONTINUOUS PASSIVE MOTION  
DEVICE  
ID.....{7587CDE1-CD6F-4A9C-8E4E-9919485E321E}  
TERM.....ACTIVITY - CONT. PASSIVE MOTION DEV.

TASK: ACTIVITY - ELEVATE HEAD OF BED

-----  
ABBREVIATION.....ACT EL HOB  
ACTIVE.....Yes  
DESCRIPTION.....ACTIVITY - ELEVATE HEAD OF BED  
DISPLAY\_NAME.....ACTIVITY - ELEVATE HEAD OF BED  
HELP\_TEXT.....ACTIVITY - ELEVATE HEAD OF BED  
ID.....{66CDCA34-6B02-4D8B-AD39-FD56BD51EC7E}  
TERM.....ACTIVITY - ELEVATE HEAD OF BED

TASK: ACTIVITY - KNEE IMMOBILIZER

-----  
ABBREVIATION.....ACT KNEE I  
ACTIVE.....Yes  
DESCRIPTION.....ACTIVITY - KNEE IMMOBILIZER  
DISPLAY\_NAME.....ACTIVITY - KNEE IMMOBILIZER  
HELP\_TEXT.....ACTIVITY - KNEE IMMOBILIZER  
ID.....{182692D9-7ECC-45B9-82F6-163DFF937438}  
TERM.....ACTIVITY - KNEE IMMOBILIZER

TASK: ACTIVITY - ORIENT TO WARD

-----  
ABBREVIATION.....ACT ORIENT  
ACTIVE.....Yes  
DESCRIPTION.....Orient patient and family to ward  
                  policies and procedures  
DISPLAY\_NAME.....ACTIVITY - ORIENT TO WARD  
HELP\_TEXT.....ACTIVITY - ORIENT TO WARD  
ID.....{35B092D6-86E0-655B-9B0E-8E819FE6967D}  
TERM.....ACTIVITY - ORIENT TO WARD

TASK: ACTIVITY - RANGE OF MOTION EXERCISES

-----  
ABBREVIATION.....ACT ROM EX  
ACTIVE.....Yes  
DESCRIPTION.....ACTIVITY - RANGE OF MOTION EXERCISES  
DISPLAY\_NAME.....ACTIVITY - RANGE OF MOTION EXERCISES  
HELP\_TEXT.....ACTIVITY - RANGE OF MOTION EXERCISES  
ID.....{D028A8B1-2548-485F-9D87-D9154858E076}  
TERM.....ACTIVITY - RANGE OF MOTION EXERCISES

TASK: ACTIVITY - SIT IN CHAIR

-----  
ABBREVIATION.....ACT SIT CH  
ACTIVE.....Yes  
DESCRIPTION.....ACTIVITY - SIT IN CHAIR  
DISPLAY\_NAME.....ACTIVITY - SIT IN CHAIR  
HELP\_TEXT.....ACTIVITY - SIT IN CHAIR  
ID.....{62765F1D-A170-450C-89C3-5F9EEFAFF71B}  
TERM.....ACTIVITY - SIT IN CHAIR

TASK: ACTIVITY - TURN PATIENT

-----  
ABBREVIATION.....ACT TURN  
ACTIVE.....Yes  
DESCRIPTION.....ACTIVITY - TURN PATIENT  
DISPLAY\_NAME.....ACTIVITY - TURN PATIENT  
HELP\_TEXT.....ACTIVITY - TURN PATIENT  
ID.....{7A545892-7660-4ADE-8017-12DD76762EE8}  
TERM.....ACTIVITY - TURN PATIENT

TASK: ADL - BATH ASSIST

-----  
ABBREVIATION.....ADLBATHASS  
ACTIVE.....Yes  
DESCRIPTION.....ADL - BATH ASSIST  
DISPLAY\_NAME.....ADL - BATH ASSIST  
HELP\_TEXT.....ADL - BATH ASSIST  
ID.....{0DF710C3-7EB8-44EA-8425-C4F5707D1744}  
TERM.....ADL - BATH ASSIST

TASK: ADL - BATH COMPLETE

-----  
ABBREVIATION.....ADLBATHCOM  
ACTIVE.....Yes  
DESCRIPTION.....ADL - BATH COMPLETE  
DISPLAY\_NAME.....ADL - BATH COMPLETE  
HELP\_TEXT.....No Help Text  
ID.....{E48D2F58-C42D-4D61-9E34-19DB98966970}  
TERM.....ADL - BATH COMPLETE

TASK: ADL - BATH PARTIAL

-----  
ABBREVIATION.....ADLBATHPAR  
ACTIVE.....Yes  
DESCRIPTION.....ADL - BATH PARTIAL  
DISPLAY\_NAME.....ADL - BATH PARTIAL  
HELP\_TEXT.....ADL - BATH PARTIAL  
ID.....{6997EBBC-1D7C-46CC-AFC0-0035A777DFBD}  
TERM.....ADL - BATH PARTIAL

TASK: ADL - BATH SELF

-----  
ABBREVIATION.....ADLBATHSEL  
ACTIVE.....Yes  
DESCRIPTION.....ADL - BATH SELF  
DISPLAY\_NAME.....ADL - BATH SELF  
HELP\_TEXT.....ADL - BATH SELF  
ID.....{F79F34A1-F147-4A68-8039-2F80E3BF7CBE}  
TERM.....ADL - BATH SELF

TASK: ADL - CATHETER CARE

-----  
ABBREVIATION.....ADLCATHCAR  
ACTIVE.....Yes  
DESCRIPTION.....ADL - CATHETER CARE  
DISPLAY\_NAME.....ADL - CATHETER CARE  
HELP\_TEXT.....ADL - CATHETER CARE  
ID.....{8AAD5864-D7A5-4886-AD58-BDDDC00D80DC}  
TERM.....ADL - CATHETER CARE

TASK: ADL - FEED PATIENT

-----  
ABBREVIATION.....ADL FEED P  
ACTIVE.....Yes  
DESCRIPTION.....ADL - FEED PATIENT  
DISPLAY\_NAME.....ADL - FEED PATIENT  
HELP\_TEXT.....ADL - FEED PATIENT  
ID.....{D525A017-C515-4FBE-9CC2-FDB5FA6CF010}  
TERM.....ADL - FEED PATIENT

TASK: ADL - ORAL HYGIENE

-----  
ABBREVIATION.....ADLORALHYG  
ACTIVE.....Yes  
DESCRIPTION.....ADL - ORAL HYGIENE  
DISPLAY\_NAME.....ADL - ORAL HYGIENE  
HELP\_TEXT.....ADL - ORAL HYGIENE  
ID.....{27044060-F4C7-4616-9CE4-D2BFE1001EEE}  
TERM.....ADL - ORAL HYGIENE

TASK: ADL - SHAVE

-----  
ABBREVIATION.....ADLSHAVE  
ACTIVE.....Yes  
DESCRIPTION.....ADL - SHAVE  
DISPLAY\_NAME.....ADL - SHAVE  
HELP\_TEXT.....ADL - SHAVE  
ID.....{B5763C28-DEF7-4A6B-86C2-1C1237B6F568}  
TERM.....ADL - SHAVE

TASK: ADL - SKIN CARE

-----  
ABBREVIATION.....ADLSKINCAR  
ACTIVE.....Yes  
DESCRIPTION.....ADL - SKIN CARE  
DISPLAY\_NAME.....ADL - SKIN CARE  
HELP\_TEXT.....ADL - SKIN CARE  
ID.....{7B41CD6E-8F0C-44B5-9415-288ABAF8271A}  
TERM.....ADL - SKIN CARE

TASK: ADL - TURN

-----  
ABBREVIATION.....ADLTURN  
ACTIVE.....Yes  
DESCRIPTION.....ADL - TURN  
DISPLAY\_NAME.....ADL - TURN  
HELP\_TEXT.....ADL - TURN  
ID.....{CCF74486-3ADE-4A67-A6A0-121F195E6FCB}  
TERM.....ADL - TURN

TASK: ASSESSMENT - CO

-----  
ABBREVIATION.....A CO  
ACTIVE.....Yes  
DESCRIPTION.....ASSESSMENT - CO  
DISPLAY\_NAME.....ASSESSMENT - CO  
HELP\_TEXT.....ASSESSMENT - CO  
ID.....{D94C0EFD-CFB8-4249-89D0-1A4BF6C777FC}  
TERM.....ASSESSMENT - CO



TASK: ASSESSMENT - COLLECT VITALS

-----  
ABBREVIATION.....A CO VITAL  
ACTIVE.....Yes  
DESCRIPTION.....ASSESSMENT - COLLECT VITALS  
DISPLAY\_NAME.....ASSESSMENT - COLLECT VITALS  
HELP\_TEXT.....ASSESSMENT - COLLECT VITALS  
ID.....{9A21F145-1DBA-0BA5-4C9D-E9E31B3F3E45}  
TERM.....ASSESSMENT - COLLECT VITALS

TASK: ASSESSMENT - DOPPLERS

-----  
ABBREVIATION.....A DOPPLERS  
ACTIVE.....Yes  
DESCRIPTION.....ASSESSMENT - DOPPLERS  
DISPLAY\_NAME.....ASSESSMENT - DOPPLERS  
HELP\_TEXT.....ASSESSMENT - DOPPLERS  
ID.....{8B1F6E7E-95D5-47A0-A23D-2A8A088930B8}  
TERM.....ASSESSMENT - DOPPLERS

TASK: ASSESSMENT - FINGERSTICK BLOOD SUGAR

-----  
ABBREVIATION.....A BLD SUGR  
ACTIVE.....Yes  
DESCRIPTION.....ASSESSMENT - FINGERSTICK BLOOD SUGAR  
DISPLAY\_NAME.....ASSESSMENT - FINGERSTICK BLOOD SUGAR  
HELP\_TEXT.....ASSESSMENT - FINGERSTICK BLOOD SUGAR  
(FSBS)  
ID.....{50A07B57-D9E6-4F67-B6BC-E21D4A60A56C}  
TERM.....ASSESSMENT - FINGERSTICK BLOOD SUGAR

TASK: ASSESSMENT - HEMODYNAMIC READINGS

-----  
ABBREVIATION.....A HEMO REA  
ACTIVE.....Yes  
DESCRIPTION.....ASSESSMENT - HEMODYNAMIC READINGS  
DISPLAY\_NAME.....ASSESSMENT - HEMODYNAMIC READINGS  
HELP\_TEXT.....ASSESSMENT - HEMODYNAMIC READINGS  
ID.....{98328403-2F7B-7404-34EA-5411633EA55C}  
TERM.....ASSESSMENT - HEMODYNAMIC READINGS

TASK: ASSESSMENT - I&O

-----  
ABBREVIATION.....A I&O  
ACTIVE.....Yes  
DESCRIPTION.....ASSESSMENT - I&O  
DISPLAY\_NAME.....ASSESSMENT - I&O  
HELP\_TEXT.....ASSESSMENT - I&O  
ID.....{53CA663E-40AC-4089-A250-10635CC64F06}  
TERM.....ASSESSMENT - I&O

TASK: ASSESSMENT - NEURO CHECKS

-----  
ABBREVIATION.....A NRO CHKS  
ACTIVE.....Yes  
DESCRIPTION.....ASSESSMENT - NEURO CHECKS  
DISPLAY\_NAME.....ASSESSMENT - NEURO CHECKS  
HELP\_TEXT.....ASSESSMENT - NEURO CHECKS  
ID.....{8E1E8DE5-EFE6-41B9-86CD-484BC9EF4F66}  
TERM.....ASSESSMENT - NEURO CHECKS

TASK: ASSESSMENT - NEUROVASCULAR CHECKS

-----  
ABBREVIATION.....A NROVASCU  
ACTIVE.....Yes  
DESCRIPTION.....ASSESSMENT - NEUROVASCULAR CHECKS  
DISPLAY\_NAME.....ASSESSMENT - NEUROVASCULAR CHECKS  
HELP\_TEXT.....ASSESSMENT - NEUROVASCULAR CHECKS  
ID.....{4C7FBC3F-F88B-426B-9DCB-74EE2FBC40EA}  
TERM.....ASSESSMENT - NEUROVASCULAR CHECKS

TASK: ASSESSMENT - PCA

-----  
ABBREVIATION.....A PCA  
ACTIVE.....Yes  
DESCRIPTION.....ASSESSMENT - PCA  
DISPLAY\_NAME.....ASSESSMENT - PCA  
HELP\_TEXT.....ASSESSMENT - PCA  
ID.....{88464AA6-3DCF-4388-B0EC-6B3D6B098FEC}  
TERM.....ASSESSMENT - PCA

TASK: ASSESSMENT - PULSE OX

-----  
ABBREVIATION.....A PULSE OX  
ACTIVE.....Yes  
DESCRIPTION.....ASSESSMENT - PULSE OX  
DISPLAY\_NAME.....ASSESSMENT - PULSE OX  
HELP\_TEXT.....ASSESSMENT - PULSE OX  
ID.....{FC2A17F3-DC69-427D-946B-C3403B92528B}  
TERM.....ASSESSMENT - PULSE OX

TASK: ASSESSMENT - RESTRAINTS

-----  
ABBREVIATION.....A RESTRAIN  
ACTIVE.....Yes  
DESCRIPTION.....ASSESSMENT - RESTRAINTS  
DISPLAY\_NAME.....ASSESSMENT - RESTRAINTS  
HELP\_TEXT.....ASSESSMENT - RESTRAINTS  
ID.....{3D02BA0C-BFDD-4D82-8B44-40C861E7FD24}  
TERM.....ASSESSMENT - RESTRAINTS

TASK: ASSESSMENT - TELEMETRY

-----  
ABBREVIATION.....A TELEMETR  
ACTIVE.....Yes  
DESCRIPTION.....ASSESSMENT - TELEMETRY  
DISPLAY\_NAME.....ASSESSMENT - TELEMETRY  
HELP\_TEXT.....ASSESSMENT - TELEMETRY  
ID.....{48C7FB6E-74C2-4EA4-984F-19C11C9072D7}  
TERM.....ASSESSMENT - TELEMETRY

TASK: ASSESSMENT - WEIGHT

-----  
ABBREVIATION.....A WEIGHT  
ACTIVE.....Yes  
DESCRIPTION.....ASSESSMENT - GET PATIENTS WEIGHT  
DISPLAY\_NAME.....ASSESSMENT - WEIGHT  
HELP\_TEXT.....ASSESSMENT - GET PATIENTS WEIGHT  
ID.....{CCB2F8DF-9319-4767-AE7D-90A7705C03FD}  
TERM.....ASSESSMENT - WEIGHT

TASK: CALL HOUSE OFFICER - BP DIASTOLIC

-----  
ABBREVIATION.....CHO DBP  
ACTIVE.....Yes  
DESCRIPTION.....CALL HOUSE OFFICER - BP DIASTOLIC  
DISPLAY\_NAME.....CALL HOUSE OFFICER - BP DIASTOLIC  
HELP\_TEXT.....CALL HOUSE OFFICER - BP DIASTOLIC  
ID.....{19A8A7AD-6BAC-4326-806C-BC80E5636FFC}  
TERM.....CALL HOUSE OFFICER - BP DIASTOLIC

TASK: CALL HOUSE OFFICER - BP SYSTOLIC

-----  
ABBREVIATION.....CHO SBP  
ACTIVE.....Yes  
DESCRIPTION.....CALL HOUSE OFFICER - BP SYSTOLIC  
DISPLAY\_NAME.....CALL HOUSE OFFICER - BP SYSTOLIC  
HELP\_TEXT.....CALL HOUSE OFFICER - BP SYSTOLIC  
ID.....{971F71F9-AAD6-46BF-A729-45F0B98FD6D3}  
TERM.....CALL HOUSE OFFICER - BP SYSTOLIC

TASK: CALL HOUSE OFFICER - PULSE

-----  
ABBREVIATION.....CHO PULSE  
ACTIVE.....Yes  
DESCRIPTION.....CALL HOUSE OFFICER - PULSE  
DISPLAY\_NAME.....CALL HOUSE OFFICER - PULSE  
HELP\_TEXT.....CALL HOUSE OFFICER - PULSE  
ID.....{81EED77A-5116-4B2D-B52B-42BAFACC69C9}  
TERM.....CALL HOUSE OFFICER - PULSE

TASK: CALL HOUSE OFFICER - RESPIRATORY RATE

-----  
ABBREVIATION.....CHO RR  
ACTIVE.....Yes  
DESCRIPTION.....CALL HOUSE OFFICER - RESPIRATORY RATE  
DISPLAY\_NAME.....CALL HOUSE OFFICER - RESPIRATORY RATE  
HELP\_TEXT.....CALL HOUSE OFFICER - RESPIRATORY RATE  
ID.....{D2DFE1EC-017F-49B8-B59C-C724F71F2E63}  
TERM.....CALL HOUSE OFFICER - RESPIRATORY RATE

TASK: CALL HOUSE OFFICER - TEMP

-----  
ABBREVIATION.....CHO TEMP  
ACTIVE.....Yes  
DESCRIPTION.....CALL HOUSE OFFICER - TEMP  
DISPLAY\_NAME.....CALL HOUSE OFFICER - TEMP  
HELP\_TEXT.....CALL HOUSE OFFICER - TEMP  
ID.....{E59A3DC0-73CB-4B0F-B031-E3CFF280EA5F}  
TERM.....CALL HOUSE OFFICER - TEMP

TASK: HYDRATION - DIET

-----  
ABBREVIATION.....HY DIET  
ACTIVE.....Yes  
DESCRIPTION.....HYDRATION - DIET  
DISPLAY\_NAME.....HYDRATION - DIET  
HELP\_TEXT.....HYDRATION - DIET  
ID.....{C0F14019-FF1A-4D42-A857-4B1E86E793FA}  
TERM.....HYDRATION - DIET

TASK: HYDRATION - HEPARIN DRIP

-----  
ABBREVIATION.....HY HEPARIN  
ACTIVE.....Yes  
DESCRIPTION.....HYDRATION - HEPARIN DRIP  
DISPLAY\_NAME.....HYDRATION - HEPARIN DRIP  
HELP\_TEXT.....HYDRATION - HEPARIN DRIP  
ID.....{EC49EE9D-C2EB-4F7D-A9D4-9746DDEF84A6}  
TERM.....HYDRATION - HEPARIN DRIP

TASK: HYDRATION - HYPERAL

-----  
ABBREVIATION.....HY HYPERAL  
ACTIVE.....Yes  
DESCRIPTION.....HYDRATION - HYPERAL  
DISPLAY\_NAME.....HYDRATION - HYPERAL  
HELP\_TEXT.....HYDRATION - HYPERAL  
ID.....{5F91CA9E-82F1-4490-8648-896864F1F68C}  
TERM.....HYDRATION - HYPERAL

TASK: HYDRATION - IV FLUIDS

-----  
ABBREVIATION.....HY IV FLU  
ACTIVE.....Yes  
DESCRIPTION.....HYDRATION - IV FLUIDS  
DISPLAY\_NAME.....HYDRATION - IV FLUIDS  
HELP\_TEXT.....HYDRATION - IV FLUIDS  
ID.....{75D6B913-147F-4318-B6A5-985438E8CF66}  
TERM.....HYDRATION - IV FLUIDS

TASK: HYDRATION - NPO

-----  
ABBREVIATION.....HY NPO  
ACTIVE.....Yes  
DESCRIPTION.....HYDRATION - NPO  
DISPLAY\_NAME.....HYDRATION - NPO  
HELP\_TEXT.....HYDRATION - NPO  
ID.....{44028E10-FB0D-4AFD-9B8D-90F157D6329D}  
TERM.....HYDRATION - NPO

TASK: HYDRATION - PERIPHERAL IV

-----  
ABBREVIATION.....HY PIV  
ACTIVE.....Yes  
DESCRIPTION.....HYDRATION - PERIPHERAL IV  
DISPLAY\_NAME.....HYDRATION - PERIPHERAL IV  
HELP\_TEXT.....HYDRATION - PERIPHERAL IV  
ID.....{1ACA57CF-4226-4EF1-AABD-142EE1DE369F}  
TERM.....HYDRATION - PERIPHERAL IV

TASK: HYDRATION - PICC/CENTRAL LINE

-----  
ABBREVIATION.....HY PICC  
ACTIVE.....Yes  
DESCRIPTION.....HYDRATION - PICC/CENTRAL LINE  
DISPLAY\_NAME.....HYDRATION - PICC/CENTRAL LINE  
HELP\_TEXT.....HYDRATION - PICC/CENTRAL LINE  
ID.....{E1D6E811-CE80-4F8A-AB00-570D1C8079CA}  
TERM.....HYDRATION - PICC/CENTRAL LINE

TASK: HYDRATION - SUPPLEMENT SNACK

-----  
ABBREVIATION.....HY SNACK S  
ACTIVE.....Yes  
DESCRIPTION.....HYDRATION - SUPPLEMENT SNACK  
DISPLAY\_NAME.....HYDRATION - SUPPLEMENT SNACK  
HELP\_TEXT.....HYDRATION - SUPPLEMENT SNACK  
ID.....{0916F365-7FBC-4FD3-A55B-CAD4E41A00B1}  
TERM.....HYDRATION - SUPPLEMENT SNACK

TASK: HYDRATION - TUBE FEEDING

-----  
ABBREVIATION.....HY TU FEED  
ACTIVE.....Yes  
DESCRIPTION.....HYDRATION - TUBE FEEDING  
DISPLAY\_NAME.....HYDRATION - TUBE FEEDING  
HELP\_TEXT.....HYDRATION - TUBE FEEDING  
ID.....{6EDD774C-DD5A-4BF6-B06D-725FD7A0B26C}  
TERM.....HYDRATION - TUBE FEEDING

TASK: NURSING INTERVENTIONS - STRAIGHT CATH

-----  
ABBREVIATION.....NI ST CATH  
ACTIVE.....Yes  
DESCRIPTION.....NURSING INTERVENTIONS - STRAIGHT CATH  
DISPLAY\_NAME.....NURSING INTERVENTIONS - STRAIGHT CATH  
HELP\_TEXT.....NURSING INTERVENTIONS - STRAIGHT CATH  
ID.....{A1F54C3F-CFD1-4767-BED8-55BF2BE872DD}  
TERM.....NURSING INTERVENTIONS - STRAIGHT CATH

TASK: NURSING INTERVENTIONS - WOUND CARE

-----  
ABBREVIATION.....NI WOU CAR  
ACTIVE.....Yes  
DESCRIPTION.....NURSING INTERVENTIONS - WOUND CARE  
DISPLAY\_NAME.....NURSING INTERVENTIONS - WOUND CARE  
HELP\_TEXT.....NURSING INTERVENTIONS - WOUND CARE  
ID.....{BFBE2EFF-3644-439F-951C-804B43BFFF30}  
TERM.....NURSING INTERVENTIONS - WOUND CARE

TASK: PRECAUTIONS - 24 HOUR CHART CHECK

-----  
ABBREVIATION.....PREC 24 CH  
ACTIVE.....Yes  
DESCRIPTION.....PRECAUTIONS - 24 HOUR CHART CHECK  
DISPLAY\_NAME.....PRECAUTIONS - 24 HOUR CHART CHECK  
HELP\_TEXT.....PRECAUTIONS - 24 HOUR CHART CHECK  
ID.....{96C39F92-CE4C-4D6B-9837-9C9FECDEC199}  
TERM.....PRECAUTIONS - 24 HOUR CHART CHECK

TASK: PRECAUTIONS - ASPIRATION

-----  
ABBREVIATION.....PREC ASPIR  
ACTIVE.....Yes  
DESCRIPTION.....PRECAUTIONS - ASPIRATION  
DISPLAY\_NAME.....PRECAUTIONS - ASPIRATION  
HELP\_TEXT.....PRECAUTIONS - ASPIRATION  
ID.....{D8D28810-655E-4707-BF6D-B853F18A1E08}  
TERM.....PRECAUTIONS - ASPIRATION

TASK: PRECAUTIONS - ETOH W/D

-----  
ABBREVIATION.....PREC ETOH  
ACTIVE.....Yes  
DESCRIPTION.....PRECAUTIONS - ETOH W/D  
DISPLAY\_NAME.....PRECAUTIONS - ETOH W/D  
HELP\_TEXT.....PRECAUTIONS - ETOH W/D  
ID.....{7800B16E-ECCD-4086-985C-0D12D4671123}  
TERM.....PRECAUTIONS - ETOH W/D

TASK: PRECAUTIONS - FALL

-----  
ABBREVIATION.....PREC FALL  
ACTIVE.....Yes  
DESCRIPTION.....PRECAUTIONS - FALL  
DISPLAY\_NAME.....PRECAUTIONS - FALL  
HELP\_TEXT.....PRECAUTIONS - FALL  
ID.....{CAB2E3CC-18AD-40AB-946B-DB0292E7FA09}  
TERM.....PRECAUTIONS - FALL

TASK: PRECAUTIONS - PATIENT TEACHING

-----  
ABBREVIATION.....PREC TEACH  
ACTIVE.....Yes  
DESCRIPTION.....PRECAUTIONS - PATIENT TEACHING  
DISPLAY\_NAME.....PRECAUTIONS - PATIENT TEACHING  
HELP\_TEXT.....PRECAUTIONS - PATIENT TEACHING  
ID.....{C58B33AB-6442-46A5-B24A-E2B856C72B98}  
TERM.....PRECAUTIONS - PATIENT TEACHING

TASK: PRECAUTIONS - VALIDATE ALARMS

-----  
ABBREVIATION.....PREC ALARM  
ACTIVE.....Yes  
DESCRIPTION.....PRECAUTIONS - VALIDATE ALARMS  
DISPLAY\_NAME.....PRECAUTIONS - VALIDATE ALARMS  
HELP\_TEXT.....PRECAUTIONS - VALIDATE ALARMS  
ID.....{CF386EC1-F5B6-4D15-8E48-8453FE6FC1CC}  
TERM.....PRECAUTIONS - VALIDATE ALARMS

TASK: PULMONARY - CHEST PHYSIO-THERAPY

-----  
ABBREVIATION.....PULMCHPPTH  
ACTIVE.....Yes  
DESCRIPTION.....PULMONARY - CHEST PHYSIO-THERAPY  
DISPLAY\_NAME.....PULMONARY - CHEST PHYSIO-THERAPY  
HELP\_TEXT.....PULMONARY - CHEST PHYSIO-THERAPY  
ID.....{F028FA64-3957-4A15-B30A-277AF594A69D}  
TERM.....PULMONARY - CHEST PHYSIO-THERAPY

TASK: PULMONARY - COUGH AND DEEP BREATH

-----  
ABBREVIATION.....PULMCOUGH  
ACTIVE.....Yes  
DESCRIPTION.....PULMONARY - COUGH AND DEEP BREATH  
DISPLAY\_NAME.....PULMONARY - COUGH AND DEEP BREATH  
HELP\_TEXT.....PULMONARY - COUGH AND DEEP BREATH  
ID.....{17A5D60E-ED50-4FBA-879D-4618563D2B98}  
TERM.....PULMONARY - COUGH AND DEEP BREATH

TASK: PULMONARY - ET TUBE CARE

-----  
ABBREVIATION.....PULM ET TU  
ACTIVE.....Yes  
DESCRIPTION.....PULMONARY - ET TUBE CARE  
DISPLAY\_NAME.....PULMONARY - ET TUBE CARE  
HELP\_TEXT.....PULMONARY - ET TUBE CARE  
ID.....{B2F57E59-33C5-49FB-BEE9-176647329910}  
TERM.....PULMONARY - ET TUBE CARE

TASK: PULMONARY - INCENTIVE SPIROMETER

-----  
ABBREVIATION.....PULMINCSPI  
ACTIVE.....Yes  
DESCRIPTION.....PULMONARY - INCENTIVE SPIROMETER  
DISPLAY\_NAME.....PULMONARY - INCENTIVE SPIROMETER  
HELP\_TEXT.....PULMONARY - INCENTIVE SPIROMETER  
ID.....{F856E287-BAAC-4D8C-AFAD-34912621F3B4}  
TERM.....PULMONARY - INCENTIVE SPIROMETER

TASK: PULMONARY - OXYGEN

-----  
ABBREVIATION.....PULM OXYGE  
ACTIVE.....Yes  
DESCRIPTION.....PULMONARY - OXYGEN  
DISPLAY\_NAME.....PULMONARY - OXYGEN  
HELP\_TEXT.....PULMONARY - OXYGEN  
ID.....{14A986A8-54F1-4D9A-A2EB-43F306704C42}  
TERM.....PULMONARY - OXYGEN

TASK: PULMONARY - SUCTION PATIENT

-----  
ABBREVIATION.....PULM-SUCTI  
ACTIVE.....Yes  
DESCRIPTION.....PULMONARY - SUCTION PATIENT  
DISPLAY\_NAME.....PULMONARY - SUCTION PATIENT  
HELP\_TEXT.....PULMONARY - SUCTION PATIENT  
ID.....{8B74C2F6-8CFB-4A60-9158-5333B980D047}  
TERM.....PULMONARY - SUCTION PATIENT





TASK: RESTRAINT - SEIZURE PRECAUTIONS

-----  
ABBREVIATION.....RES SEI PR  
ACTIVE.....Yes  
DESCRIPTION.....Restraint - Seizure Precautions  
DISPLAY\_NAME.....Restraint - Seizure Precautions  
HELP\_TEXT.....Restraint - Seizure Precautions  
ID.....{71606A0A-05E6-4DFD-9E00-677A07167199}  
TERM.....RESTRAINT - SEIZURE PRECAUTIONS

TASK: TESTS/PROC. - ACCUCHECKS

-----  
ABBREVIATION.....TP ACCUCHE  
ACTIVE.....Yes  
DESCRIPTION.....TESTS/PROC. - ACCUCHECKS  
DISPLAY\_NAME.....TESTS/PROC. - ACCUCHECKS  
HELP\_TEXT.....TESTS/PROC. - ACCUCHECKS  
ID.....{0A5E0495-B990-449B-B5F6-10497842534E}  
TERM.....TESTS/PROC. - ACCUCHECKS

TASK: TESTS/PROC. - BLOOD CULTURES

-----  
ABBREVIATION.....T/P BLOOD  
ACTIVE.....Yes  
DESCRIPTION.....TESTS/PROC. - BLOOD CULTURES  
DISPLAY\_NAME.....TESTS/PROC. - BLOOD CULTURES  
HELP\_TEXT.....TESTS/PROC. - BLOOD CULTURES  
ID.....{BC2F7EFE-3990-47D2-863A-6FB2227746F3}  
TERM.....TESTS/PROC. - BLOOD CULTURES

TASK: TESTS/PROC. - DAILY LABS

-----  
ABBREVIATION.....T/P D LABS  
ACTIVE.....Yes  
DESCRIPTION.....TESTS/PROC. - DAILY LABS  
DISPLAY\_NAME.....TESTS/PROC. - DAILY LABS  
HELP\_TEXT.....TESTS/PROC. - DAILY LABS  
ID.....{03B3A828-D4B0-40CF-B12C-1DD936DE98AC}  
TERM.....TESTS/PROC. - DAILY LABS

TASK: TESTS/PROC. - MORNING LABS

-----  
ABBREVIATION.....T/P M LABS  
ACTIVE.....Yes  
DESCRIPTION.....TESTS/PROC. - MORNING LABS  
DISPLAY\_NAME.....TESTS/PROC. - MORNING LABS  
HELP\_TEXT.....TESTS/PROC. - MORNING LABS  
ID.....{6F81B3D1-CCBC-4BC9-AFF0-6400DF212EB}  
TERM.....TESTS/PROC. - MORNING LABS

TASK: TESTS/PROC. - ZERO BALANCE HEMO. LINES

-----  
ABBREVIATION.....T/P 0 HEMO  
ACTIVE.....Yes  
DESCRIPTION.....TESTS/PROC. - ZERO BALANCE HEMO. LINES  
DISPLAY\_NAME.....TESTS/PROC. - ZERO BALANCE HEMO. LINES  
HELP\_TEXT.....TESTS/PROC. - ZERO BALANCE HEMO. LINES  
ID.....{9CEDB731-A785-48B7-AC92-84B10FCCCE86}  
TERM.....TESTS/PROC. - ZERO BALANCE HEMO. LINES

TASK: TUBES & DRAINS - CONDOM CATH

-----  
ABBREVIATION.....T&D C CATH  
ACTIVE.....Yes  
DESCRIPTION.....TUBES & DRAINS - CONDOM CATH  
DISPLAY\_NAME.....TUBES & DRAINS - CONDOM CATH  
HELP\_TEXT.....TUBES & DRAINS - CONDOM CATH  
ID.....{385EFB2B-55BB-486E-B168-D1FFF71E086C}  
TERM.....TUBES & DRAINS - CONDOM CATH

TASK: TUBES & DRAINS - FOLEY

-----  
ABBREVIATION.....T&D FOLEY  
ACTIVE.....Yes  
DESCRIPTION.....TUBES & DRAINS - FOLEY  
DISPLAY\_NAME.....TUBES & DRAINS - FOLEY  
HELP\_TEXT.....TUBES & DRAINS - FOLEY  
ID.....{14CECC17-A780-41DD-98A4-6E226630BC65}  
TERM.....TUBES & DRAINS - FOLEY

TASK: TUBES & DRAINS - I/O CATH

-----  
ABBREVIATION.....T&D I/O CA  
ACTIVE.....Yes  
DESCRIPTION.....TUBES & DRAINS - I/O CATH  
DISPLAY\_NAME.....TUBES & DRAINS - I/O CATH  
HELP\_TEXT.....TUBES & DRAINS - I/O CATH  
ID.....{5BA34E92-5248-4954-A677-443C4B2183E8}  
TERM.....TUBES & DRAINS - I/O CATH

TASK: TUBES & DRAINS - NG TUBE

-----  
ABBREVIATION.....T&D NG TU  
ACTIVE.....Yes  
DESCRIPTION.....TUBES & DRAINS - NG TUBE  
DISPLAY\_NAME.....TUBES & DRAINS - NG TUBE  
HELP\_TEXT.....TUBES & DRAINS - NG TUBE  
ID.....{BBFCE51A-9CCF-44AC-8C37-4A779232A4B7}  
TERM.....TUBES & DRAINS - NG TUBE

TASK: TUBES & DRAINS - SUPRAPUBIC DRAIN

-----  
ABBREVIATION.....T&D SUP DR  
ACTIVE.....Yes  
DESCRIPTION.....TUBES & DRAINS - SUPRAPUBIC DRAIN  
DISPLAY\_NAME.....TUBES & DRAINS - SUPRAPUBIC DRAIN  
HELP\_TEXT.....TUBES & DRAINS - SUPRAPUBIC DRAIN  
ID.....{6A4A5F28-35D3-494E-BA97-2E0DE4A7B8CC}  
TERM.....TUBES & DRAINS - SUPRAPUBIC DRAIN

UNIT: BEATS PER MINUTE

-----

ABBREVIATION.....bpm  
 ACTIVE.....Yes  
 DESCRIPTION.....No Description  
 DISPLAY\_NAME.....Beats Per Minute  
 HELP\_TEXT.....No Help Text  
 ID.....{8319B3B0-0F31-1393-154D-59A86692D637}  
 TERM.....BEATS PER MINUTE

UNIT: BRADEN SCALE

-----

ABBREVIATION.....BRADEN SCA  
 ACTIVE.....Yes  
 DESCRIPTION.....No Description  
 DISPLAY\_NAME.....Braden Scale  
 HELP\_TEXT.....No Help Text  
 ID.....{B5EE25BB-DA66-4AC9-9A3B-C05497B2647C}  
 TERM.....BRADEN SCALE

UNIT: BREATHS PER MINUTE

-----

ABBREVIATION.....rpm  
 ACTIVE.....Yes  
 DESCRIPTION.....No Description  
 DISPLAY\_NAME.....Resp/Min  
 HELP\_TEXT.....No Help Text  
 ID.....{DCF89527-FFAE-F68E-0871-602F3C32DDB1}  
 TERM.....BREATHS PER MINUTE

UNIT: CALORIE

-----

ABBREVIATION.....CALORIE  
 ACTIVE.....Yes  
 DESCRIPTION.....A unit of food energy  
 DISPLAY\_NAME.....Calorie  
 HELP\_TEXT.....The amount of energy required to  
                   raise the temperature of a liter of  
                   water one degree centigrade at sea  
                   level  
 ID.....{9F7DE371-AB4B-45EE-A8E3-825822419CCA}  
 TERM.....CALORIE

UNIT: CENTIMETER

-----

ABBREVIATION.....cm  
 ACTIVE.....Yes  
 DESCRIPTION.....Measurement of length  
 DISPLAY\_NAME.....Centimeter  
 HELP\_TEXT.....Metric  
 ID.....{EF700458-8C2D-16E8-02DC-9E1711584C53}  
 TERM.....CENTIMETER

UNIT: CENTIMETERS H2O

-----  
ABBREVIATION.....CM H2O  
ACTIVE.....Yes  
DESCRIPTION.....CENTIMETERS OF WATER  
DISPLAY\_NAME.....CM H2O  
HELP\_TEXT.....No Help Text  
ID.....{F7030BB1-AD2C-15F6-FF41-6DA893C851C1}  
TERM.....CENTIMETERS H2O

UNIT: CIWA SCALE

-----  
ABBREVIATION.....CIWA SCALE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....CIWA SCALE  
HELP\_TEXT.....Scale For Addition Of All CIWA Scores  
ID.....{0ED378B9-C1A4-4ABA-82B7-F62BC0D62886}  
TERM.....CIWA SCALE

UNIT: CUBIC CENTIMETER

-----  
ABBREVIATION.....Cubic Cent  
ACTIVE.....Yes  
DESCRIPTION.....Cubic Centimeter  
DISPLAY\_NAME.....Cubic Centimeter  
HELP\_TEXT.....No Help Text  
ID.....{4698E6A0-DA6C-4A0A-800A-395BBCD2871E}  
TERM.....CUBIC CENTIMETER

UNIT: DAY

-----  
ABBREVIATION.....Day  
ACTIVE.....Yes  
DESCRIPTION.....Measurement of passage of time (24  
hours)  
DISPLAY\_NAME.....Day  
HELP\_TEXT.....No Help Text  
ID.....{742619B2-A8C7-1A53-45CD-C3E25A05D58C}  
TERM.....DAY

UNIT: DEGREES (ANGLE)

-----  
ABBREVIATION.....DEG  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Degrees (Angle)  
HELP\_TEXT.....No Help Text  
ID.....{5F407D72-3009-2553-7DDC-82BAF2043FCA}  
TERM.....DEGREES (ANGLE)

UNIT: DEGREES C

-----  
ABBREVIATION.....C  
ACTIVE.....Yes  
DESCRIPTION.....Metric Temperature  
DISPLAY\_NAME.....Celsius  
HELP\_TEXT.....No Help Text  
ID.....{1964CB4C-0584-DA80-C7A9-E7400DCD817D}  
TERM.....DEGREES C

UNIT: DEGREES F

-----  
ABBREVIATION.....F  
ACTIVE.....Yes  
DESCRIPTION.....US measurement of temperature  
DISPLAY\_NAME.....Fahrenheit  
HELP\_TEXT.....No Help Text  
ID.....{59FC6CA4-9389-B770-5A9C-8C5A38C572D8}  
TERM.....DEGREES F

UNIT: DYNE \* SECONDS / CM(-5) \* M2

-----  
ABBREVIATION.....D\*S/CM\*M2  
ACTIVE.....Yes  
DESCRIPTION.....dyne \* seconds / cm(-5) \* m2  
DISPLAY\_NAME.....dyne \* seconds / cm(-5) \* m2  
HELP\_TEXT.....dyne \* seconds / cm(-5) \* m2  
ID.....{1241878F-AECA-B1A8-BF12-D5FA583FC1B2}  
TERM.....DYNE \* SECONDS / CM(-5) \* M2

UNIT: DYNE \* SECONDS / CM(-5) / M2

-----  
ABBREVIATION.....D\*S/CM/M2  
ACTIVE.....Yes  
DESCRIPTION.....dyne \* seconds/cm(-5) / m2  
DISPLAY\_NAME.....dyne \* seconds/cm(-5) / m2  
HELP\_TEXT.....dyne \* seconds/cm(-5) / m2  
ID.....{A642531F-A72E-CF60-D31A-F38AD12A8B30}  
TERM.....DYNE \* SECONDS / CM(-5) / M2

UNIT: DYNE-SEC/CM5

-----  
ABBREVIATION.....dyne-sec/c  
ACTIVE.....Yes  
DESCRIPTION.....dyne-sec/cm5  
DISPLAY\_NAME.....dyne-sec/cm5  
HELP\_TEXT.....dyne-sec/cm5  
ID.....{01A538D8-41AF-4D0B-B827-D557060A618A}  
TERM.....DYNE-SEC/CM5

UNIT: EVENT(S)

-----  
ABBREVIATION.....EVENT(S)  
ACTIVE.....Yes  
DESCRIPTION.....Event(s)  
DISPLAY\_NAME.....Event(s)  
HELP\_TEXT.....No Help Text  
ID.....{723F907D-8B71-4183-AE54-66D0C3C72DCA}  
TERM.....EVENT(S)

UNIT: EVENTS PER HOUR

-----  
ABBREVIATION.....EVENT/H  
ACTIVE.....Yes  
DESCRIPTION.....NUMBER OF TIMES PER HOUR THAT AN  
OBSERVABLE EVENT OCCURS  
DISPLAY\_NAME.....Events Per Hour  
HELP\_TEXT.....No Help Text  
ID.....{3D648D4F-02A9-474A-65BF-70608E75BC29}  
TERM.....EVENTS PER HOUR

UNIT: FOOT

-----  
ABBREVIATION.....FT  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Foot/Feet  
HELP\_TEXT.....No Help Text  
ID.....{66D6F7EC-3884-4790-AC35-10B2AE84257F}  
TERM.....FOOT

UNIT: FRENCH

-----  
ABBREVIATION.....FRENCH  
ACTIVE.....Yes  
DESCRIPTION.....UNIT OF MEASURE FOR ARTIFICIAL AIRWAYS  
DISPLAY\_NAME.....French  
HELP\_TEXT.....No Help Text  
ID.....{3734E0E0-087B-098A-984A-E42264CF6567}  
TERM.....FRENCH

UNIT: GAUGE

-----  
ABBREVIATION.....GAUGE  
ACTIVE.....Yes  
DESCRIPTION.....Gauge For Needles  
DISPLAY\_NAME.....Gauge  
HELP\_TEXT.....No Help Text  
ID.....{5670BBC1-4962-4FA1-A4D0-1E399F2D06B1}  
TERM.....GAUGE



UNIT: GLASGOW COMA SCALE

-----  
ABBREVIATION.....GCS  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Glasgow Coma Scale  
HELP\_TEXT.....No Help Text  
ID.....{CA674718-C98D-669F-8954-F9EC91C2F6F3}  
TERM.....GLASGOW COMA SCALE

UNIT: GM-M/M2/BEAT

-----  
ABBREVIATION.....gm-m/m2/be  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....gm-m/m2/beat  
HELP\_TEXT.....No Help Text  
ID.....{DB5404C8-1880-4B54-B053-DF0875E72C3B}  
TERM.....GM-M/M2/BEAT

UNIT: GRAMS PER HOUR

-----  
ABBREVIATION.....gm/hr  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....gm/hr  
HELP\_TEXT.....No Help Text  
ID.....{A50449B1-DF42-4D18-8F02-59937EA6CDE4}  
TERM.....GRAMS PER HOUR

UNIT: GRAMS PER MINUTE

-----  
ABBREVIATION.....gm/min  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Grams Per Minute  
HELP\_TEXT.....No Help Text  
ID.....{07D005E8-B078-A48F-96B5-8CCAA92EC97F}  
TERM.....GRAMS PER MINUTE

UNIT: HOUR

-----  
ABBREVIATION.....hr  
ACTIVE.....Yes  
DESCRIPTION.....Measurement of time  
DISPLAY\_NAME.....Hour  
HELP\_TEXT.....No Help Text  
ID.....{60034E2A-D088-E9B6-9B98-D0F4357A059F}  
TERM.....HOUR

UNIT: INCH

-----  
ABBREVIATION.....in  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Inch  
HELP\_TEXT.....No Help Text  
ID.....{9CB32FC9-4130-11F8-E750-AF1415556705}  
TERM.....INCH

UNIT: KG/M2

-----  
ABBREVIATION.....KG/M2  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....KG/M2  
HELP\_TEXT.....No Help Text  
ID.....{E3204003-83B3-4D12-872D-C32A50F177D0}  
TERM.....KG/M2

UNIT: KGS/HOUR

-----  
ABBREVIATION.....KGS/HOUR  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....kgs/hour  
HELP\_TEXT.....  
ID.....{704F239F-CAAF-4F32-B705-62906403B874}  
TERM.....KGS/HOUR

UNIT: KILOGRAM

-----  
ABBREVIATION.....kg  
ACTIVE.....Yes  
DESCRIPTION.....Metric measurement of weight  
DISPLAY\_NAME.....Kilogram  
HELP\_TEXT.....No Help Text  
ID.....{0AC4E9AC-41CB-32EA-9674-74FDFAA6B4F4}  
TERM.....KILOGRAM

UNIT: LITER

-----  
ABBREVIATION.....L  
ACTIVE.....Yes  
DESCRIPTION.....Volume, 1000 ml  
DISPLAY\_NAME.....Liter  
HELP\_TEXT.....No Help Text  
ID.....{AFADB383-5867-4D7C-8ADB-E76782E5943D}  
TERM.....LITER

UNIT: LITERS PER MINUTE

-----  
ABBREVIATION.....L/min  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L/min  
HELP\_TEXT.....No Help Text  
ID.....{CFCD64D2-8555-8513-5DB0-F0522C19DCCF}  
TERM.....LITERS PER MINUTE

UNIT: LITERS PER MINUTE PER METER2 (BSA)

-----  
ABBREVIATION.....L/min/m2  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....L/min/m2 (BSA)  
HELP\_TEXT.....No Help Text  
ID.....{F033FDA0-C6D6-0D93-F836-C4A42708FDF2}  
TERM.....LITERS PER MINUTE PER METER2 (BSA)

UNIT: MEWS: SCORE UNIT

-----  
ABBREVIATION.....MEWS-Score  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....MEWS: Score Unit  
HELP\_TEXT.....  
ID.....{C640633E-785A-4ED3-8CD0-356A139ACF36}  
TERM.....MEWS: SCORE UNIT

UNIT: MICROGRAMS

-----  
ABBREVIATION.....mcg  
ACTIVE.....Yes  
DESCRIPTION.....Micrograms  
DISPLAY\_NAME.....mcg  
HELP\_TEXT.....Micrograms  
ID.....{5D2A9B6F-4DD1-48C6-B223-09CAF66A5A50}  
TERM.....MICROGRAMS

UNIT: MICROGRAMS PER HOUR

-----  
ABBREVIATION.....mcg/hr  
ACTIVE.....Yes  
DESCRIPTION.....Micrograms per Hour  
DISPLAY\_NAME.....mcg/hr  
HELP\_TEXT.....Micrograms per Hour  
ID.....{971A5FA5-2FD7-4ABD-9035-229B00DC484E}  
TERM.....MICROGRAMS PER HOUR

UNIT: MICROGRAMS PER KILOGRAM PER HOUR

-----  
ABBREVIATION.....mcg/kg/hr  
ACTIVE.....Yes  
DESCRIPTION.....Micrograms per Kilogram per Hour  
DISPLAY\_NAME.....mcg/kg/hr  
HELP\_TEXT.....Micrograms per Kilogram per Hour  
ID.....{E5C07334-6932-4329-88A9-2F56780C5DA9}  
TERM.....MICROGRAMS PER KILOGRAM PER HOUR

UNIT: MICROGRAMS PER KILOGRAM PER MINUTE

-----  
ABBREVIATION.....mcg/kg/min  
ACTIVE.....Yes  
DESCRIPTION.....Micrograms per Kilogram per Minute  
DISPLAY\_NAME.....mcg/kg/min  
HELP\_TEXT.....Micrograms per Kilogram per Minute  
ID.....{37423DE7-B103-456A-9B36-6C2FBBC55CBE}  
TERM.....MICROGRAMS PER KILOGRAM PER MINUTE

UNIT: MICROGRAMS PER MINUTE

-----  
ABBREVIATION.....mcg/min  
ACTIVE.....Yes  
DESCRIPTION.....Micrograms per Minute  
DISPLAY\_NAME.....mcg/min  
HELP\_TEXT.....Micrograms per Minute  
ID.....{BC684DBB-7DCC-4209-A773-D1243616BD69}  
TERM.....MICROGRAMS PER MINUTE

UNIT: MICROGRAMS PER SHIFT

-----  
ABBREVIATION.....mcg/shift  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Micrograms per Shift  
HELP\_TEXT.....  
ID.....{8FF15E2C-54DF-43C3-B5F1-C4AD323CE043}  
TERM.....MICROGRAMS PER SHIFT

UNIT: MILLIAMPS

-----  
ABBREVIATION.....mA  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....mA  
HELP\_TEXT.....No Help Text  
ID.....{5146A84A-C500-CB21-4718-0E5755691FD5}  
TERM.....MILLIAMPS

UNIT: MILLIBAR

-----  
ABBREVIATION.....mBar  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Millibar  
HELP\_TEXT.....No Help Text  
ID.....{DB8DBC26-6A46-4F2B-AD14-B6A8212F4FCA}  
TERM.....MILLIBAR

UNIT: MILLIGRAMS

-----  
ABBREVIATION.....mg  
ACTIVE.....Yes  
DESCRIPTION.....Milligrams  
DISPLAY\_NAME.....mg  
HELP\_TEXT.....Milligrams  
ID.....{E8CF6B0F-FA19-4F22-BA9D-6EBAA93D7D0D}  
TERM.....MILLIGRAMS

UNIT: MILLIGRAMS PER HOUR

-----  
ABBREVIATION.....mg/hr  
ACTIVE.....Yes  
DESCRIPTION.....Milligrams per Hour  
DISPLAY\_NAME.....mg/hr  
HELP\_TEXT.....Milligrams per Hour  
ID.....{1DCC6F14-E4FF-4BDC-ADBA-99192FB44BE7}  
TERM.....MILLIGRAMS PER HOUR

UNIT: MILLIGRAMS PER KILOGRAM PER HOUR

-----  
ABBREVIATION.....mg/kg/hr  
ACTIVE.....Yes  
DESCRIPTION.....Milligrams per Kilogram per Hour  
DISPLAY\_NAME.....mg/kg/hr  
HELP\_TEXT.....Milligrams per Kilogram per Hour  
ID.....{6900F44E-3581-4558-83E2-76FB376D6711}  
TERM.....MILLIGRAMS PER KILOGRAM PER HOUR

UNIT: MILLIGRAMS PER KILOGRAM PER MINUTE

-----  
ABBREVIATION.....mg/kg/min  
ACTIVE.....Yes  
DESCRIPTION.....Milligrams per Kilogram per Minute  
DISPLAY\_NAME.....mg/kg/min  
HELP\_TEXT.....Milligrams per Kilogram per Minute  
ID.....{7C1C10C2-8B88-4F45-B667-0F988090CD47}  
TERM.....MILLIGRAMS PER KILOGRAM PER MINUTE

UNIT: MILLIGRAMS PER MINUTE

-----  
ABBREVIATION.....mg/min  
ACTIVE.....Yes  
DESCRIPTION.....Milligrams per Minute  
DISPLAY\_NAME.....mg/min  
HELP\_TEXT.....Milligrams per Minute  
ID.....{C5E49AC7-5671-4E84-8047-C7B78607CDA2}  
TERM.....MILLIGRAMS PER MINUTE

UNIT: MILLIGRAMS PER SHIFT

-----  
ABBREVIATION.....mg/shift  
ACTIVE.....Yes  
DESCRIPTION.....  
DISPLAY\_NAME.....Milligrams Per Shift  
HELP\_TEXT.....  
ID.....{518BC06B-84C0-48E3-8241-E123B7B4C642}  
TERM.....MILLIGRAMS PER SHIFT

UNIT: MILLILITER

-----  
ABBREVIATION.....ml  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....ml  
HELP\_TEXT.....No Help Text  
ID.....{EAF55107-F8D5-4F0A-6211-B4AB57A874B4}  
TERM.....MILLILITER

UNIT: MILLILITERS PER BEAT

-----  
ABBREVIATION.....ml/beat  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....ml/beat  
HELP\_TEXT.....No Help Text  
ID.....{548699B0-98CF-2F0F-F289-415228AFAC42}  
TERM.....MILLILITERS PER BEAT

UNIT: MILLILITERS PER BEAT PER METER SQUARE (BSA)

-----  
ABBREVIATION.....ml/beat/m2  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....ml/beat/m2(BSA)  
HELP\_TEXT.....No Help Text  
ID.....{574CDB00-A831-EF35-4F2C-A98C0F58923F}  
TERM.....MILLILITERS PER BEAT PER METER SQUARE  
(BSA)

UNIT: MILLILITERS PER HOUR

-----  
ABBREVIATION.....ml/hr  
ACTIVE.....Yes  
DESCRIPTION.....ml/hr  
DISPLAY\_NAME.....ml/hr  
HELP\_TEXT.....No Help Text  
ID.....{345899A8-9DED-418F-A53F-F8084A5601A1}  
TERM.....MILLILITERS PER HOUR

UNIT: MILLILITERS PER METER SQUARE (BSA)

-----  
ABBREVIATION.....mL/M2  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....mL/M2  
HELP\_TEXT.....No Help Text  
ID.....{42C265B4-76B3-44FC-985F-EA1EFCEDA244}  
TERM.....MILLILITERS PER METER SQUARE (BSA)

UNIT: MILLILITERS PER MINUTE

-----  
ABBREVIATION.....ml/min  
ACTIVE.....Yes  
DESCRIPTION.....Milliliters per Minute  
DISPLAY\_NAME.....ml/min  
HELP\_TEXT.....Milliliters per Minute  
ID.....{25F4A228-D614-457A-9885-8F74D4596897}  
TERM.....MILLILITERS PER MINUTE

UNIT: MILLILITERS PER MINUTE PER METER SQUARE

-----  
ABBREVIATION.....ml/min/m2  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....ml/min/m2(BSA)  
HELP\_TEXT.....No Help Text  
ID.....{6404DFC7-F8B4-3C97-80EF-96627E9C0270}  
TERM.....MILLILITERS PER MINUTE PER METER SQUARE

UNIT: MILLIMETERS

-----  
ABBREVIATION.....mm  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....mm  
HELP\_TEXT.....No Help Text  
ID.....{23F389B2-8308-91DC-2620-EDFFF07B14A5}  
TERM.....MILLIMETERS

UNIT: MILLIMETERS OF MERCURY

-----  
ABBREVIATION.....mm/Hg  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....mm Hg  
HELP\_TEXT.....No Help Text  
ID.....{57CF7550-FD97-351E-DF60-56F8F68BE7C6}  
TERM.....MILLIMETERS OF MERCURY

UNIT: MILLIVOLTS

-----  
ABBREVIATION.....mV  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....mV  
HELP\_TEXT.....No Help Text  
ID.....{40B59357-C9F2-0E5C-DCD7-6C815811F90D}  
TERM.....MILLIVOLTS

UNIT: MILLILITERS PER KILOGRAM PER HOUR

-----  
ABBREVIATION.....ml kg hr  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....ml/kg/hour  
HELP\_TEXT.....No Help Text  
ID.....{2F63C17C-8246-4466-9DC2-C5E5BE1777C2}  
TERM.....MILLILITERS PER KILOGRAM PER HOUR

UNIT: MILLILITERS PER KILOGRAM PER MINUTE

-----  
ABBREVIATION.....ml kg mi  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....ml/kg/min  
HELP\_TEXT.....No Help Text  
ID.....{E564F111-8B23-404A-B610-7707A5B1195D}  
TERM.....MILLILITERS PER KILOGRAM PER MINUTE

UNIT: MINUTE

-----  
ABBREVIATION.....Min  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Minute  
HELP\_TEXT.....No Help Text  
ID.....{0AF0056E-0EB1-3712-6974-A0F5DC60010B}  
TERM.....MINUTE



UNIT: MODIFIED ALDRETE SCORE

-----  
ABBREVIATION.....MOD ALD SC  
ACTIVE.....Yes  
DESCRIPTION.....Modified Aldrete Scoring System Unit  
DISPLAY\_NAME.....Modified Aldrete Score  
HELP\_TEXT.....Unit of measure for Aldrete Scoring  
                  System  
ID.....{8D9A43E8-D4E6-4BBE-A8E6-72677342BC82}  
TERM.....MODIFIED ALDRETE SCORE

UNIT: MONTH

-----  
ABBREVIATION.....Mo  
ACTIVE.....Yes  
DESCRIPTION.....Month of the year  
DISPLAY\_NAME.....Month  
HELP\_TEXT.....1 - 12  
ID.....{8933A469-2D06-D3DA-C2C4-A8FB3DD894D6}  
TERM.....MONTH

UNIT: MORSE SCORE

-----  
ABBREVIATION.....MORSE SCOR  
ACTIVE.....Yes  
DESCRIPTION.....Morse Score  
DISPLAY\_NAME.....Morse Score  
HELP\_TEXT.....Scale for recording Morse Score Values  
ID.....{301282A4-2BC8-430A-853B-C817F9A61878}  
TERM.....MORSE SCORE

UNIT: MULTIVOLTS

-----  
ABBREVIATION.....Multivolts  
ACTIVE.....No  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Multivolts  
HELP\_TEXT.....No Help Text  
ID.....{2882958B-C91D-4B8C-BEE7-382F3D7CCC49}  
TERM.....MULTIVOLTS

UNIT: OUNCE

-----  
ABBREVIATION.....oz  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Ounce  
HELP\_TEXT.....No Help Text  
ID.....{642336FB-7BC9-A1D5-99E5-05447DAD3B75}  
TERM.....OUNCE

UNIT: PAINAD SCALE

-----  
ABBREVIATION.....PAINAD SC  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....PAINAD Scale  
HELP\_TEXT.....No Help Text  
ID.....{D6CE83BB-3DBE-4E94-9E91-E21C07A12CAD}  
TERM.....PAINAD SCALE

UNIT: PERCENTAGE

-----  
ABBREVIATION.....%  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Percentage  
HELP\_TEXT.....No Help Text  
ID.....{0B3163CE-2C5C-B196-DBBE-A1BC6A970B2A}  
TERM.....PERCENTAGE

UNIT: POINTS

-----  
ABBREVIATION.....POINTS  
ACTIVE.....Yes  
DESCRIPTION.....Points  
DISPLAY\_NAME.....Points  
HELP\_TEXT.....No Help Text  
ID.....{A41F9244-F8A9-4CBD-AD62-32165D77414B}  
TERM.....POINTS

UNIT: POUND

-----  
ABBREVIATION.....lb  
ACTIVE.....Yes  
DESCRIPTION.....US measurement of Weight  
DISPLAY\_NAME.....Pound  
HELP\_TEXT.....No Help Text  
ID.....{F9CBB050-E809-52FF-6C2D-51B7E9557F98}  
TERM.....POUND

UNIT: QUAD LUMEN CATHETER

-----  
ABBREVIATION.....QUAD LUMEN  
ACTIVE.....Yes  
DESCRIPTION.....Quad Lumen Catheter  
DISPLAY\_NAME.....Quad Lumen Catheter  
HELP\_TEXT.....No Help Text  
ID.....{C894DDF4-5976-4435-B0E4-70EC7A909565}  
TERM.....QUAD LUMEN CATHETER

UNIT: RANGE

-----  
ABBREVIATION.....range  
ACTIVE.....Yes  
DESCRIPTION.....Defines the minimum and maximum  
                  values in a range  
DISPLAY\_NAME.....Range  
HELP\_TEXT.....No Help Text  
ID.....{DD90689F-3B19-BC1D-1875-E9EAB08C9C9E}  
TERM.....RANGE

UNIT: SCALE

-----  
ABBREVIATION.....SCALE  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Scale  
HELP\_TEXT.....No Help Text  
ID.....{57E423C4-98D9-0891-0982-0B84EFD2D132}  
TERM.....SCALE

UNIT: SECOND

-----  
ABBREVIATION.....sec  
ACTIVE.....Yes  
DESCRIPTION.....Measurement of passage of time (1/60  
                  of a minute)  
DISPLAY\_NAME.....Second  
HELP\_TEXT.....No Help Text  
ID.....{82E7C9B1-38E7-D568-F42D-23F18AA08819}  
TERM.....SECOND

UNIT: SEWS SCALE

-----  
ABBREVIATION.....SEWS SCALE  
ACTIVE.....Yes  
DESCRIPTION.....Severity Of Ethanol Withdrawl Scale  
DISPLAY\_NAME.....SEWS SCALE  
HELP\_TEXT.....Severity Of Ethanol Withdrawl Scale  
ID.....{4C8DFD61-E4FB-4918-93B0-33DAAF4EF1D1}  
TERM.....SEWS SCALE

UNIT: SI UNIT

-----  
ABBREVIATION.....SI UNITS  
ACTIVE.....Yes  
DESCRIPTION.....Unit of measure for Specific Gravity  
DISPLAY\_NAME.....SI Unit  
HELP\_TEXT.....No Help Text  
ID.....{42C760F2-BBFC-4523-82C1-953CB41B4016}  
TERM.....SI UNIT

UNIT: SQUARE METERS

-----  
ABBREVIATION.....M2  
ACTIVE.....Yes  
DESCRIPTION.....Surface area measured in square meters  
DISPLAY\_NAME.....Square Meters  
HELP\_TEXT.....Surface area measured in square meters  
ID.....{5E49A34D-CC24-ABDA-7320-14071BA88781}  
TERM.....SQUARE METERS

UNIT: STROKE WORK PER BEAT

-----  
ABBREVIATION.....gm-M/beat  
ACTIVE.....Yes  
DESCRIPTION.....Stroke Work Per Beat  
DISPLAY\_NAME.....gm-M / beat  
HELP\_TEXT.....Stroke Work Per Beat  
ID.....{61721D2C-02BF-14AA-EB7F-602194924117}  
TERM.....STROKE WORK PER BEAT

UNIT: STROKE WORK PER BEAT PER M2

-----  
ABBREVIATION.....gm-M/B/M2  
ACTIVE.....Yes  
DESCRIPTION.....STROKE WORK PER BEAT PER M2  
DISPLAY\_NAME.....gm-M / Beat / M2  
HELP\_TEXT.....Systemic Vascular Resistance Index  
ID.....{900EBF1F-B706-0419-1697-58FE11257EF3}  
TERM.....STROKE WORK PER BEAT PER M2

UNIT: UNITS

-----  
ABBREVIATION.....UNITS  
ACTIVE.....Yes  
DESCRIPTION.....Units  
DISPLAY\_NAME.....Units  
HELP\_TEXT.....Units  
ID.....{7E9E6B8E-8B47-4B06-81FF-227E46707401}  
TERM.....UNITS

UNIT: UNITS PER HOUR

-----  
ABBREVIATION.....units/hr  
ACTIVE.....Yes  
DESCRIPTION.....Units per Hour  
DISPLAY\_NAME.....units/hr  
HELP\_TEXT.....Units per Hour  
ID.....{7AE61C6C-F8D9-4D6B-8A7C-AE07F93EB3D9}  
TERM.....UNITS PER HOUR

UNIT: UNITS PER MINUTE

-----  
ABBREVIATION.....units/min  
ACTIVE.....Yes  
DESCRIPTION.....Units per Minute  
DISPLAY\_NAME.....Units/min  
HELP\_TEXT.....Units per Minute  
ID.....{3C31F294-AD93-4509-908A-72F3995B01B9}  
TERM.....UNITS PER MINUTE

UNIT: UNITS/kg/HOUR

-----  
ABBREVIATION.....UNkgHR  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Units/kg/Hour  
HELP\_TEXT.....No Help Text  
ID.....{89A25E97-A6B8-49B2-934F-797C80C80711}  
TERM.....UNITS/kg/HOUR

UNIT: YEAR

-----  
ABBREVIATION.....yr  
ACTIVE.....Yes  
DESCRIPTION.....No Description  
DISPLAY\_NAME.....Year  
HELP\_TEXT.....No Help Text  
ID.....{F800BBF4-BC4F-6F5D-BB28-BD3EB5ADC34D}  
TERM.....YEAR

UNIT: mg/dl

-----  
ABBREVIATION.....mg/dl  
ACTIVE.....Yes  
DESCRIPTION.....Milligrams per deciliter  
DISPLAY\_NAME.....mg/dl  
HELP\_TEXT.....No Help Text  
ID.....{796F843B-3EF5-4E9B-AAD6-254088D77FFB}  
TERM.....mg/dl

UNIT: pH

-----  
ABBREVIATION.....pH  
ACTIVE.....Yes  
DESCRIPTION.....pH  
DISPLAY\_NAME.....pH  
HELP\_TEXT.....No Help Text  
ID.....{EF524413-89B4-4ABA-BD68-1FC6713C623D}  
TERM.....pH



- 0 -

0 (data).....	770
0 (q-pair).....	108, 404, 405, 406, 407, 408, 410, 472, 473, 623, 624, 625, 626, 627, 628, 629, 630, 632, 633
0 - ABSENT (data).....	782
0 - ABSENT (q-pair).....	119, 120, 127, 136, 159, 160, 204, 571
0 - ALERT (data).....	782
0 - ALERT (q-pair).....	388
0 - HR 51-100 (data).....	782
0 - HR 51-100 (q-pair).....	389
0 - INDEPENDENT - NO HELP OR STAFF OVERSIGHT (data)....	782
0 - INDEPENDENT - NO HELP OR STAFF OVERSIGHT (q-pair)..	357, 359, 361, 363, 365, 367, 369, 371, 373, 375
0 - INDEPENDENT - NO HELP PROVIDED (data).....	782
0 - INDEPENDENT - NO HELP PROVIDED (q-pair).....	355
0 - NO ANXIETY AT EASE (data).....	783
0 - NO ANXIETY AT EASE (q-pair).....	118
0 - NO NAUSEA/VOMITING (data).....	783
0 - NO NAUSEA/VOMITING (q-pair).....	121
0 - NO SETUP OR PHYSICAL HELP FROM STAFF (data).....	783
0 - NO SETUP OR PHYSICAL HELP FROM STAFF (q-pair).....	356, 358, 360, 362, 364, 366, 368, 370, 372, 374, 376
0 - NO SWEAT VISIBLE (data).....	783
0 - NO SWEAT VISIBLE (q-pair).....	123
0 - NO TREMOR (data).....	783
0 - NO TREMOR (q-pair).....	126
0 - NONE (data).....	784
0 - NONE (q-pair).....	124
0 - NORMAL ACTIVITY (data).....	784
0 - NORMAL ACTIVITY (q-pair).....	117
0 - ORIENTED AND CAN DO SERIAL ADDITIONS (data).....	784
0 - ORIENTED AND CAN DO SERIAL ADDITIONS (q-pair).....	122
0 - RR 9-14 (data).....	784
0 - RR 9-14 (q-pair).....	391
0 - SBP 101-199 (data).....	784
0 - SBP 101-199 (q-pair).....	393
0 - TEMP 36.1-38 C (data).....	785
0 - TEMP 36.1-38 C (q-pair).....	394
0 = 0 EXTREMITIES (data).....	785
0 = 0 EXTREMITIES (q-pair).....	25
0 = <90% (data).....	785
0 = <90% (q-pair).....	28
0 = ABSENT PULSE (data).....	785
0 = APNEIC (data).....	785
0 = APNEIC (q-pair).....	29
0 = BP +/- 50% PREOP (data).....	786
0 = BP +/- 50% PREOP (q-pair).....	26
0 = NOT RESPONDING (data).....	786
0 = NOT RESPONDING (q-pair).....	27

- 1 -

1 (data).....	770
1 (q-pair).....	108, 185, 344, 472,

473, 625, 627, 682,  
715

1 - HR 101-110 (data).....	786
1 - HR 101-110 (q-pair).....	389
1 - BARELY PERCEPTIBLE SWEATING / PALMS MOIST (data)...	786
1 - BARELY PERCEPTIBLE SWEATING / PALMS MOIST (q-pair).....	123
1 - CANNOT DO SERIAL ADDITIONS (data).....	786
1 - CANNOT DO SERIAL ADDITIONS (q-pair).....	122
1 - HR 40-50 (data).....	787
1 - HR 40-50 (q-pair).....	389
1 - LESS THAN 45MLS/HR (data).....	787
1 - LESS THAN 45MLS/HR (q-pair).....	390
1 - MILD NAUSEA / NO VOMITING (data).....	787
1 - MILD NAUSEA / NO VOMITING (q-pair).....	121
1 - MILDLY ANXIOUS (data).....	787
1 - MILDLY ANXIOUS (q-pair).....	118
1 - NEW AGITATION CONFUSION (data).....	787
1 - NEW AGITATION CONFUSION (q-pair).....	388
1 - NO RESPONSE (data).....	788
1 - NO RESPONSE (q-pair).....	204, 205, 207
1 - NOT VISIBLE BUT CAN BE FELT FINGERTIP (data).....	788
1 - NOT VISIBLE BUT CAN BE FELT FINGERTIP (q-pair).....	126
1 - RESPONDS TO VOICE (data).....	788
1 - RESPONDS TO VOICE (q-pair).....	388
1 - RR 15-20 (data).....	788
1 - RR 15-20 (q-pair).....	391
1 - SBP 81-100 (data).....	788
1 - SBP 81-100 (q-pair).....	393
1 - SETUP HELP ONLY (data).....	789
1 - SETUP HELP ONLY (q-pair).....	356, 358, 360, 362, 364, 366, 368, 370, 372, 374, 376
1 - SUPERVISION - OVERSIGHT (data).....	789
1 - SUPERVISION - OVERSIGHT (q-pair).....	357, 359, 361, 363, 365, 367, 369, 371, 373, 375
1 - SUPERVISION - OVERSIGHT HELP ONLY (data).....	789
1 - SUPERVISION - OVERSIGHT HELP ONLY (q-pair).....	355
1 - TEMP 35.1-36 C (data).....	789
1 - TEMP 35.1-36 C (q-pair).....	394
1 - TEMP 38.1-38.5 C (data).....	789
1 - TEMP 38.1-38.5 C (q-pair).....	394
1 - VERY MILD (data).....	790
1 - VERY MILD (q-pair).....	120
1 - VERY MILD HARSHNESS OR ABILITY TO FRIGHTEN (data).....	790
1 - VERY MILD HARSHNESS OR ABILITY TO FRIGHTEN (q-pair).....	119
1 - VERY MILD ITCHING PINS/NEEDLES BURNING (data).....	790
1 - VERY MILD ITCHING PINS/NEEDLES BURNING (q-pair).....	124
1 - VERY MILD SENSITIVITY (data).....	790
1 - VERY MILD SENSITIVITY (q-pair).....	127
1 - WEAK (data).....	790
1 - WEAK (q-pair).....	571
1 = 2 EXTREMITIES (data).....	791
1 = 2 EXTREMITIES (q-pair).....	25
1 = AROUSABLE (data).....	791
1 = AROUSABLE (q-pair).....	27
1 = BP +/- 20-49% PREOP (data).....	791
1 = BP +/- 20-49% PREOP (q-pair).....	26
1 = DYSPNEA (data).....	791
1 = DYSPNEA (q-pair).....	29
1 = NEEDS O2 (data).....	791
1 = NEEDS O2 (q-pair).....	28



1 HOUR (data).....	792
1 HOUR (q-pair).....	492
1+=PALPABLE, THREADY, EASILY OBLITERATED (data).....	792
1- MILD ANXIOUS (data).....	792
1- SOMEWHAT MORE THAN NORMAL ACTIVITY (data).....	792
1- SOMEWHAT MORE THAN NORMAL ACTIVITY (q-pair).....	117
10 (data).....	773
10 (q-pair).....	21, 108, 344, 405, 472, 473, 682, 715
11 (data).....	773
11 (q-pair).....	344
12 (data).....	773
12 (q-pair).....	344
12.5 (data).....	773
13 (data).....	774
13 (q-pair).....	188
13.5 (data).....	774
13.5 (q-pair).....	188
14 (data).....	774
14 (q-pair).....	188
14.5 (data).....	774
14.5 (q-pair).....	188
15 (data).....	774
15 (q-pair).....	188, 404, 408, 410
15.5 (data).....	775
15.5 (q-pair).....	188
16 (data).....	775
16 (q-pair).....	188
16.5 (data).....	775
16.5 (q-pair).....	188
17 (data).....	775
17 (q-pair).....	188
17.5 (data).....	775
17.5 (q-pair).....	188
18 (data).....	776
18 (q-pair).....	188
18.5 (data).....	776
18.5 (q-pair).....	188
1800 CALORIE CONTROLLED (data).....	792
1800 CALORIE CONTROLLED (q-pair).....	150
1800 DIABETIC (data).....	793
1800 DIABETIC (q-pair).....	150
19 (data).....	776
19 (q-pair).....	188
19.5 (data).....	776
19.5 (q-pair).....	188
1:1 (data).....	793
1:1 (q-pair).....	145, 237, 240
1:1 OBSERVATION (data).....	793
1:1 OBSERVATION (q-pair).....	57
1:2 (data).....	793
1:2 (q-pair).....	145, 237, 240
1:3 (data).....	793
1:3 (q-pair).....	145, 237, 240
1:4 (data).....	794
1:4 (q-pair).....	145, 237, 240, 704
1:8 (data).....	794
1:8 (q-pair).....	145, 237, 240
1A + 2V (data).....	794
1A + 2V (q-pair).....	464
1A + 3V (data).....	794
1A + 3V (q-pair).....	464

1A +1V (data).....794  
1A +1V (q-pair).....464

- 2 -

2 (data).....770  
2 (q-pair).....108, 117, 118, 121,  
123, 126, 185, 344,  
472, 473, 682, 715  
2 - ABNORMAL EXTENSION FROM PAIN (data).....795  
2 - ABNORMAL EXTENSION FROM PAIN (q-pair).....205  
2 - BEADS OF SWEAT OBVIOUS ON FOREHEAD (data).....795  
2 - DISORIENTED FOR DATE < 2 CALENDAR DAYS (data).....795  
2 - DISORIENTED FOR DATE < 2 CALENDAR DAYS (q-pair)....122  
2 - HR 0-39 (data).....795  
2 - HR 0-39 (q-pair).....389  
2 - HR 111-129 (data).....795  
2 - HR 111-129 (q-pair).....389  
2 - INCOMPREHENSIBLE SOUNDS (data).....796  
2 - INCOMPREHENSIBLE SOUNDS (q-pair).....207  
2 - INTERMITTENT NAUSEA WITH DRY HEAVES (data).....796  
2 - LESS THAN 30MLS/HR (data).....796  
2 - LESS THAN 30MLS/HR (q-pair).....390  
2 - LIMITED ASSISTANCE (data).....796  
2 - LIMITED ASSISTANCE (q-pair).....357, 359, 361, 363,  
365, 367, 369, 371,  
373, 375  
2 - MILD (data).....796  
2 - MILD (q-pair).....120  
2 - MILD HARSHNESS OR ABILITY TO FRIGHTEN (data).....797  
2 - MILD HARSHNESS OR ABILITY TO FRIGHTEN (q-pair)....119  
2 - MILD ITCHING PINS/NEEDLES BURNING NUMBNESS (data)..797  
2 - MILD ITCHING PINS/NEEDLES BURNING NUMBNESS (q-pair)124  
2 - MILD SENSITIVITY (data).....797  
2 - MILD SENSITIVITY (q-pair).....127  
2 - MODERATE WITH PATIENTS ARMS EXTENDED (data).....797  
2 - NORMAL (data).....797  
2 - NORMAL (q-pair).....571  
2 - ONE PERSON PHYSICAL ASSIST (data).....798  
2 - ONE PERSON PHYSICAL ASSIST (q-pair).....356, 358, 360, 362,  
364, 366, 368, 370,  
372, 374, 376  
2 - PHYSICAL HELP LIMITED TO TRANSFER ONLY (data).....798  
2 - PHYSICAL HELP LIMITED TO TRANSFER ONLY (q-pair)....355  
2 - RESPONDS TO PAIN (data).....798  
2 - RESPONDS TO PAIN (q-pair).....388  
2 - RR 0-8 (data).....798  
2 - RR 0-8 (q-pair).....391  
2 - RR 21-29 (data).....798  
2 - RR 21-29 (q-pair).....391  
2 - SBP 200 OR MORE (data).....799  
2 - SBP 200 OR MORE (q-pair).....393  
2 - SBP 71-80 (data).....799  
2 - SBP 71-80 (q-pair).....393  
2 - TEMP 35.0 C OR LESS (data).....799  
2 - TEMP 35.0 C OR LESS (q-pair).....394  
2 - TEMP MORE THAN 38.6 C (data).....799  
2 - TEMP MORE THAN 38.6 C (q-pair).....394  
2 - TO PAIN (data).....799  
2 - TO PAIN (q-pair).....204  
2 = 4 EXTREMITIES (data).....800  
2 = 4 EXTREMITIES (q-pair).....25

2 = BP +/- 20% PREOP (data).....	800
2 = BP +/- 20% PREOP (q-pair).....	26
2 = BREATHE FREELY (data).....	800
2 = BREATHE FREELY (q-pair).....	29
2 = FULLY AWAKE (data).....	800
2 = FULLY AWAKE (q-pair).....	27
2 = ON ROOM AIR (data).....	800
2 = ON ROOM AIR (q-pair).....	28
2+=NORMAL, NOT EASILY OBLITERATED (data).....	801
2- MILD HARSHNESS OR ABILITY TO FRIGHTEN (data).....	801
2.0 CAL FORMULA (data).....	801
20 (data).....	776
20 (q-pair).....	405, 407
20 MINUTES SINCE LAST IV NARCOTIC/SEDATIVE (data).....	801
20 MINUTES SINCE LAST IV NARCOTIC/SEDATIVE (q-pair)....	465
20.5 (data).....	777
20.5 (q-pair).....	188
20/100 (data).....	801
20/100 (q-pair).....	152
20/15 (data).....	802
20/15 (q-pair).....	152
20/20 (data).....	802
20/20 (q-pair).....	152
20/200 (data).....	802
20/200 (q-pair).....	152
20/25 (data).....	802
20/25 (q-pair).....	152
20/30 (data).....	802
20/30 (q-pair).....	152
20/300 (data).....	803
20/300 (q-pair).....	152
20/40 (data).....	803
20/40 (q-pair).....	152
20/400 (data).....	803
20/400 (q-pair).....	152
20/50 (data).....	803
20/50 (q-pair).....	152
20/60 (data).....	803
20/60 (q-pair).....	152
20/70 (data).....	804
20/70 (q-pair).....	152
20/80 (data).....	804
20/80 (q-pair).....	152
20/800 (data).....	804
20/800 (q-pair).....	152
2000 CALORIE CONTROLLED (data).....	804
2000 CALORIE CONTROLLED (q-pair).....	150
2000 DIABETIC (data).....	804
2000 DIABETIC (q-pair).....	150
21 (data).....	777
21 (q-pair).....	188
21.5 (data).....	777
21.5 (q-pair).....	188
22 (data).....	777
22 (q-pair).....	188
22.5 (data).....	777
22.5 (q-pair).....	188
2200 CALORIE CONTROLLED (data).....	805
2200 CALORIE CONTROLLED (q-pair).....	150
2200 DIABETIC (data).....	805
2200 DIABETIC (q-pair).....	150
23 (data).....	778

23 (q-pair).....	188
23.5 (data).....	778
23.5 (q-pair).....	188
24 (data).....	778
24 (q-pair).....	188
24.5 (data).....	778
24.5 (q-pair).....	188
2400 CALORIE CONTROLLED (data).....	805
2400 CALORIE CONTROLLED (q-pair).....	150
2400 DIABETIC (data).....	805
2400 DIABETIC (q-pair).....	150
25 (data).....	778
25 (q-pair).....	188, 406
25 PERCENT ALBUMIN (data).....	805
25 PERCENT ALBUMIN (prod).....	1205
25 PERCENT ALBUMIN (q-pair).....	132, 309, 317, 320
25.5 (data).....	779
25.5 (q-pair).....	188
26 (data).....	779
26 (q-pair).....	188
26.5 (data).....	779
26.5 (q-pair).....	188
27 (data).....	779
27 (q-pair).....	188
27.5 (data).....	779
27.5 (q-pair).....	188
28 (data).....	780
28 (q-pair).....	188
28.5 (data).....	780
28.5 (q-pair).....	188
29 (data).....	780
2:4 (data).....	806
2:4 (q-pair).....	704
2A (data).....	806
2A (q-pair).....	464
2A + 1V (data).....	806
2A + 1V (q-pair).....	464
2A + 2V (data).....	806
2A + 2V (q-pair).....	464
2V (data).....	806
2V (q-pair).....	464

- 3 -

3 - HR 130 OR MORE (data).....	807
3 - HR 130 OR MORE (q-pair).....	389
3 (data).....	770
3 (q-pair).....	108, 117, 118, 121, 123, 126, 185, 344, 472, 473, 623, 624, 626, 628, 629, 630, 632, 633, 682, 715
3 - ABNORMAL FLEXION FROM PAIN (data).....	807
3 - ABNORMAL FLEXION FROM PAIN (q-pair).....	205
3 - BEADS OF SWEAT OBVIOUS ON FOREHEAD (data).....	807
3 - DISORIENTED FOR DATE > 2 CALENDAR DAYS (data).....	807
3 - DISORIENTED FOR DATE > 2 CALENDAR DAYS (q-pair).....	122
3 - EXTENSIVE ASSISTANCE (data).....	807
3 - EXTENSIVE ASSISTANCE (q-pair).....	357, 359, 361, 363, 365, 367, 369, 371, 373, 375
3 - INAPPROPRIATE WORDS (data).....	808

3 - INAPPROPRIATE WORDS (q-pair).....	207
3 - INTERMITTENT NAUSEA WITH DRY HEAVES (data).....	808
3 - LESS THAN 10ML/HR (data).....	808
3 - LESS THAN 10ML/HR (q-pair).....	390
3 - MODERATE (data).....	808
3 - MODERATE (q-pair).....	120
3 - MODERATE HARSHNESS OR ABILITY TO FRIGHTEN (data)...	808
3 - MODERATE HARSHNESS OR ABILITY TO FRIGHTEN (q-pair).	119
3 - MODERATE ITCHING PINS/NEEDLES BURNING (data).....	809
3 - MODERATE ITCHING PINS/NEEDLES BURNING (q-pair).....	124
3 - MODERATE SENSITIVITY (data).....	809
3 - MODERATE SENSITIVITY (q-pair).....	127
3 - MODERATE WITH PATIENTS ARMS EXTENDED (data).....	809
3 - PHYSICAL HELP IN PART OF BATHING ACTIVITY (data)...	809
3 - PHYSICAL HELP IN PART OF BATHING ACTIVITY (q-pair).	355
3 - RR 30 OR MORE (data).....	809
3 - RR 30 OR MORE (q-pair).....	391
3 - SBP 70 OR LESS (data).....	810
3 - SBP 70 OR LESS (q-pair).....	393
3 - STRONG (data).....	810
3 - STRONG (q-pair).....	571
3 - TO VERBAL COMMAND (data).....	810
3 - TO VERBAL COMMAND (q-pair).....	204
3 - TWO + PERSONS PHYSICAL ASSIST (data).....	810
3 - TWO + PERSONS PHYSICAL ASSIST (q-pair).....	356, 358, 360, 362, 364, 366, 368, 370, 372, 374, 376
3 - UNRESPONSIVE (data).....	810
3 - UNRESPONSIVE (q-pair).....	388
3 WAY FOLEY (data).....	811
3 WAY FOLEY (q-pair).....	219
3+=BOUNDING, EASILY PALPABLE, CANNOT BE OBLITERATED (data)	
.....	811
30 (data).....	780
30 (q-pair).....	404
31 (data).....	780
34 ml (data).....	811
34 ml (q-pair).....	232
3:4 (data).....	811
3:4 (q-pair).....	704

- 4 -

4 (data).....	770
4 (q-pair).....	21, 108, 185, 344, 472, 473, 682, 715
4 - BEADS OF SWEAT OBVIOUS ON FOREHEAD (data).....	811
4 - BEADS OF SWEAT OBVIOUS ON FOREHEAD (q-pair).....	123
4 - DISORIENTED AND CONVERSES (data).....	812
4 - DISORIENTED AND CONVERSES (q-pair).....	207
4 - DISORIENTED FOR PLACE OR PERSON (data).....	812
4 - DISORIENTED FOR PLACE OR PERSON (q-pair).....	122
4 - INTERMITTENT NAUSEA WITH DRY HEAVES (data).....	812
4 - INTERMITTENT NAUSEA WITH DRY HEAVES (q-pair).....	121
4 - MODERATE WITH PATIENTS ARMS EXTENDED (data).....	812
4 - MODERATE WITH PATIENTS ARMS EXTENDED (q-pair).....	126
4 - MODERATELY ANXIOUS (data).....	812
4 - MODERATELY ANXIOUS (q-pair).....	118
4 - MODERATELY FIDGETY AND RESTLESS (data).....	813
4 - MODERATELY FIDGETY AND RESTLESS (q-pair).....	117
4 - MODERATELY SEVERE (data).....	813
4 - MODERATELY SEVERE (q-pair).....	120

4 - MODERATELY SEVERE HALLUCINATIONS (data).....	813
4 - MODERATELY SEVERE HALLUCINATIONS (q-pair).....	119, 124, 127
4 - SPONTANEOUSLY (data).....	813
4 - SPONTANEOUSLY (q-pair).....	204
4 - TOTAL DEPENDENCE (data).....	813
4 - TOTAL DEPENDENCE (q-pair).....	355
4 - TOTAL DEPENDENCE - FULL STAFF PERFORMANCE (data)...	814
4 - TOTAL DEPENDENCE - FULL STAFF PERFORMANCE (q-pair).	357, 359, 361, 363, 365, 367, 369, 371, 373, 375
4 - WITHDRAWS FROM PAIN (data).....	814
4 - WITHDRAWS FROM PAIN (q-pair).....	205
4 HOUR (data).....	814
4 HOUR (q-pair).....	492
4.5 (data).....	771
4.5 (q-pair).....	21
40 ml (data).....	814
40 ml (q-pair).....	232
4:4 (data).....	814
4:4 (q-pair).....	704

- 5 -

5 (data).....	771
5 (q-pair).....	21, 108, 117, 118, 121, 123, 126, 185, 344, 472, 473, 682, 715
5 - CONSTANT NAUSEA W FREQUENT HEAVES AND VOMITING (data) .....	815
5 - DRENCHING SWEAT (data).....	815
5 - LOCALIZES PAIN (data).....	815
5 - LOCALIZES PAIN (q-pair).....	205
5 - ORIENTED AND CONVERSES (data).....	815
5 - ORIENTED AND CONVERSES (q-pair).....	207
5 - SEVERE (data).....	815
5 - SEVERE (q-pair).....	120
5 - SEVERE EVEN WITH ARMS NOT EXTENDED (data).....	816
5 - SEVERE HALLUCINATIONS (data).....	816
5 - SEVERE HALLUCINATIONS (q-pair).....	119, 124, 127
5 PERCENT ALBUMIN (data).....	816
5 PERCENT ALBUMIN (prod).....	1205
5 PERCENT ALBUMIN (q-pair).....	132, 309, 317, 320
5.5 (data).....	771
5.5 (q-pair).....	21

- 6 -

6 (data).....	771
6 (q-pair).....	21, 108, 117, 118, 121, 123, 126, 344, 472, 473, 682, 715
6 - CONSTANT NAUSEA W FREQUENT HEAVES AND VOMITING (data) .....	816
6 - DRENCHING SWEAT (data).....	816
6 - EXTREMELY SEVERE HALLUCINATIONS (data).....	817
6 - EXTREMELY SEVERE HALLUCINATIONS (q-pair).....	119, 124, 127
6 - OBEYS COMMANDS (data).....	817
6 - OBEYS COMMANDS (q-pair).....	205
6 - SEVERE EVEN WITH ARMS NOT EXTENDED (data).....	817
6 - VERY SEVERE (data).....	817
6 - VERY SEVERE (q-pair).....	120

6 HOUR (data).....	817
6 HOUR (q-pair).....	492
6 PERCENT HETASTARCH (data).....	818
6 PERCENT HETASTARCH (prod).....	1205
6 PERCENT HETASTARCH (q-pair).....	132, 309, 320
6.5 (data).....	771
6.5 (q-pair).....	21
60 MINUTES SINCE ADMINISTRATION OF REVERSAL AGENT (data)	
.....	818
60 MINUTES SINCE ADMINISTRATION OF REVERSAL AGENT (q-pair)	
.....	465

- 7 -

7 (data).....	772
7 (q-pair).....	21, 108, 344, 472,
	473, 682, 715
7 - ACUTE PANIC (data).....	818
7 - ACUTE PANIC (q-pair).....	118
7 - CONSTANT NAUSEA W FREQUENT HEAVES AND VOMITING (data)	
.....	818
7 - CONSTANT NAUSEA W FREQUENT HEAVES AND VOMITING (q-pair)	
.....	121
7 - CONTINUOUS HALLUCINATIONS (data).....	818
7 - CONTINUOUS HALLUCINATIONS (q-pair).....	119, 124, 127
7 - DRENCHING SWEAT (data).....	819
7 - DRENCHING SWEAT (q-pair).....	123
7 - EXTREMELY SEVERE (data).....	819
7 - EXTREMELY SEVERE (q-pair).....	120
7 - PACING BACK AND FORTH MOST OF THE INTERVIEW (data).	819
7 - PACING BACK AND FORTH MOST OF THE INTERVIEW (q-pair)	
.....	117
7 - SEVERE EVEN WITH ARMS NOT EXTENDED (data).....	819
7 - SEVERE EVEN WITH ARMS NOT EXTENDED (q-pair).....	126
7.5 (data).....	772
7.5 (q-pair).....	21, 239

- 8 -

8 (data).....	772
8 (q-pair).....	21, 108, 239, 344,
	472, 473, 682, 715
8 - ACTIVITY DID NOT OCCUR (data).....	819
8 - ACTIVITY DID NOT OCCUR (q-pair).....	355, 356, 357, 358,
	359, 360, 361, 362,
	363, 364, 365, 366,
	367, 368, 369, 370,
	371, 372, 373, 374,
	375, 376
8 HOUR (data).....	820
8 HOUR (q-pair).....	492
8.5 (data).....	772
8.5 (q-pair).....	21, 239

- 9 -

9 (data).....	772
9 (q-pair).....	21, 108, 239, 344,
	472, 473, 682, 715
9.5 (data).....	773
9.5 (q-pair).....	21, 239
99 (data).....	781

99 (q-pair).....473

- + -

+1 = SLIGHT (data).....781  
+1 = SLIGHT (q-pair).....159, 160  
+2 = MODERATE (data).....781  
+2 = MODERATE (q-pair).....159, 160  
+3 = MARKED (data).....781  
+3 = MARKED (q-pair).....159, 160  
+4 = VERY MARKED (data).....781  
+4 = VERY MARKED (q-pair).....159, 160

- A -

AAI (data).....820  
AAI (q-pair).....455  
AAIR (data).....820  
AAIR (q-pair).....455  
AAO (data).....820  
AAO (q-pair).....455  
ABDOMEN (loca).....1127  
ABDOMEN (q-pair).....160, 473, 658  
ABDOMEN BOWEL SOUNDS.....1  
ABDOMEN LLQ (loca).....1127  
ABDOMEN LLQ (q-pair).....1, 473, 658  
ABDOMEN LUQ (loca).....1127  
ABDOMEN LUQ (q-pair).....1, 473, 658  
ABDOMEN RLQ (loca).....1127  
ABDOMEN RLQ (q-pair).....1, 473, 658  
ABDOMEN RUQ (loca).....1127  
ABDOMEN RUQ (q-pair).....1, 473, 658  
ABDOMINAL (data).....820  
ABDOMINAL (loca).....1128  
ABDOMINAL (q-pair).....116, 577  
ABDOMINAL AORTA (loca).....1128  
ABDOMINAL AORTA (q-pair).....79, 86  
ABDOMINAL CONDITION.....2  
ABDOMINAL GIRTH.....3  
ABLE TO AMBULATE (data).....821  
ABLE TO AMBULATE (q-pair).....465, 707  
ABRASION (data).....821  
ABRASION (q-pair).....766  
ABSENT (data).....821  
ABSENT (q-pair).....1, 74, 79, 86, 137,  
154, 196, 198, 222,  
350, 352, 618, 619,  
620, 621, 708, 711,  
712  
ABSENT TO MINIMAL BLEEDING (data).....821  
ABSENT TO MINIMAL BLEEDING (q-pair).....465  
ACCELERATING (data).....821  
ACCELERATING (q-pair).....740  
ACHING (data).....822  
ACHING (q-pair).....470  
ACTIVE RANGE OF MOTION (data).....822  
ACTIVE RANGE OF MOTION (q-pair).....12  
ACTIVITIES OF DAILY LIVING (data).....822  
ACTIVITIES OF DAILY LIVING (q-pair).....164  
ACTIVITIES OF DAILY LIVING: WANDERING PREVENTION MEASURES  
.....4  
ACTIVITY - ABD PILLOW (task).....1228



ACTIVITY - AMBULATE (task).....	1228
ACTIVITY - CONT. PASSIVE MOTION DEV. (task).....	1228
ACTIVITY - ELEVATE HEAD OF BED (task).....	1228
ACTIVITY - KNEE IMMOBILIZER (task).....	1228
ACTIVITY - ORIENT TO WARD (task).....	1229
ACTIVITY - RANGE OF MOTION EXERCISES (task).....	1229
ACTIVITY - SIT IN CHAIR (task).....	1229
ACTIVITY - TURN PATIENT (task).....	1229
ACTUAL (q-pair).....	165, 227, 570, 694, 716, 746, 748
ACTUAL (qual).....	1223
ADL - ACTIVITIES OF DAILY LIVING.....	5
ADL - AMBULATION ASSISTANCE.....	6
ADL - AMBULATION DISTANCE.....	7
ADL - BATH ASSIST (task).....	1229
ADL - BATH COMPLETE (task).....	1230
ADL - BATH PARTIAL (task).....	1230
ADL - BATH SELF (task).....	1230
ADL - BED MOBILITY.....	8
ADL - CATHETER CARE (task).....	1230
ADL - DRESSING.....	9
ADL - EATING.....	10
ADL - EMOTIONAL SUPPORT.....	11
ADL - FEED PATIENT (task).....	1230
ADL - ORAL HYGIENE (task).....	1231
ADL - OTHER ACTIVITY.....	12
ADL - REPOSITION.....	13
ADL - SHAVE (task).....	1231
ADL - SKIN CARE (task).....	1231
ADL - TURN (task).....	1231
ADL - VISITORS.....	14
ADL - ORAL CARE.....	15
ADMISSION (data).....	822
ADMISSION (q-pair).....	142, 761
ADMISSION DIAGNOSIS (data).....	822
ADULT (data).....	823
ADULT (meth).....	1179
ADULT (q-pair).....	64
ADULT CUFF (meth).....	1179
ADULT LARGE (data).....	823
ADULT LARGE (q-pair).....	64
ADULT SMALL (data).....	823
ADULT SMALL (q-pair).....	64
ADVANCED CARE PLANNING (data).....	823
ADVANCED CARE PLANNING (q-pair).....	164
ADI - SPIRITUAL CARE.....	16
AEROSOL/HUMIDIFIED MASK (meth).....	1179
AEROSOL/HUMIDIFIED MASK (q-pair).....	686
AFFECT.....	17
AFTER EXERCISE (meth).....	1179
AFTER EXERCISE (q-pair).....	63, 570, 578
AGENT CASHIER (data).....	823
AGENT CASHIER (q-pair).....	516
AGITATED (data).....	824
AGITATED (q-pair).....	56, 351
AGONAL (data).....	824
AGONAL (q-pair).....	84, 577
AGONIZING (data).....	824
AGONIZING (q-pair).....	470
AICD (data).....	824
AICD (q-pair).....	396
AIR FLOW DEVICE (data).....	824

AIR FLOW DEVICE (q-pair).....	744
AIR LEAK (data).....	825
AIR LEAK (q-pair).....	113
AIRBORNE (data).....	825
AIRBORNE (q-pair).....	333
AIRWAY ARTIFICIAL CUFF STATUS.....	18
AIRWAY ARTIFICIAL POSITION AT LIP/NOSE.....	19
AIRWAY ARTIFICIAL SIZE - FRENCH.....	20
AIRWAY ARTIFICIAL SIZE - MILIMETERS.....	21
AIRWAY ARTIFICIAL TYPE.....	22
AIRWAY REFLEXES INTACT / SECRETIONS (data).....	825
AIRWAY REFLEXES INTACT / SECRETIONS (q-pair).....	465
ALARM PARAMETERS VERIFIED.....	23
ALARMS.....	24
ALARMS ON (data).....	825
ALARMS ON (q-pair).....	190
ALBUMIN (data).....	825
ALBUMIN (prod).....	1205
ALBUMIN (q-pair).....	320
ALCOHOL PORT PROTECTOR (data).....	826
ALCOHOL PORT PROTECTOR (q-pair).....	343
ALDRETE - ACTIVITY.....	25
ALDRETE - CIRCULATION.....	26
ALDRETE - CONSCIOUSNESS.....	27
ALDRETE - O2 SATURATION.....	28
ALDRETE - RESPIRATION.....	29
ALDRETE - TOTAL.....	30
ALDRETE SCORE 8 OR GREATER (data).....	826
ALDRETE SCORE 8 OR GREATER (q-pair).....	465
ALDRETE SCORE AT BASELINE (data).....	826
ALDRETE SCORE AT BASELINE (q-pair).....	465
ALERT (data).....	826
ALERT (q-pair).....	56, 351, 490
ALERT AND ORIENTED (data).....	826
ALERT AND ORIENTED (q-pair).....	151, 465
ALL DOWN (data).....	827
ALL DOWN (q-pair).....	634
ALL UP (data).....	827
ALL UP (q-pair).....	634
ALLERGY (data).....	827
ALLERGY (q-pair).....	40
ALLERGY WRISTBAND (data).....	827
ALLERGY WRISTBAND (q-pair).....	509
ALOPECIA PARTIAL (data).....	827
ALOPECIA PARTIAL (q-pair).....	221
ALOPECIA TOTAL (data).....	828
ALOPECIA TOTAL (q-pair).....	221
AMBER (data).....	828
AMBER (q-pair).....	157, 173, 211, 714,
753	
AMBU BAG (meth).....	1179
AMBU BAG (q-pair).....	687
AMBULATE AND ASSIST (data).....	828
AMBULATED (data).....	828
AMBULATED (q-pair).....	5
AMBULATORY CARE UNIT (data).....	828
AMBULATORY CARE UNIT (q-pair).....	467
AMPHOTERICIN B (prod).....	1205
AMPHOTERICIN B (q-pair).....	316
AMPUTEE (data).....	829
AMPUTEE (q-pair).....	571
ANALGESIA (data).....	829

ANALGESIA (q-pair).....	164
ANESTHESIA PLANNED.....	31
ANESTHESIOLOGIST (data).....	829
ANESTHESIOLOGIST (q-pair).....	554
ANGINAL (data).....	829
ANGINAL (q-pair).....	470
ANKLE (loca).....	1128
ANKLE (q-pair).....	116
ANNOYING (data).....	829
ANNOYING (q-pair).....	470
ANTECUBITAL (qual).....	1223
ANTI-EMBOLUS STOCKINGS KNEE HIGH (data).....	830
ANTI-EMBOLUS STOCKINGS KNEE HIGH (q-pair).....	158
ANTI-EMBOLUS STOCKINGS THIGH HIGH (data).....	830
ANTI-EMBOLUS STOCKINGS THIGH HIGH (q-pair).....	158
ANTI-MICROBIAL DRESSING (data).....	830
ANTI-MICROBIAL DRESSING (q-pair).....	96, 758
ANTIBIOTIC (prod).....	1206
ANTIBIOTIC (q-pair).....	316
ANTIBIOTIC OINTMENT (data).....	830
ANTIBIOTIC OINTMENT (q-pair).....	548
ANTIMICROBIAL DISK (data).....	830
ANTIMICROBIAL DISK (q-pair).....	758
ANURIC (data).....	831
ANURIC (q-pair).....	713
ANXIOUS (data).....	831
ANXIOUS (q-pair).....	56, 161
AOO (data).....	831
AOO (q-pair).....	455
APICAL (loca).....	1128
APICAL (q-pair).....	570
APNEA (data).....	831
APNEA (q-pair).....	398
APPEARANCE.....	32
APPEARS ASLEEP (data).....	831
APPEARS ASLEEP (q-pair).....	54, 490
APPLIANCE CHANGED (data).....	832
APPLIANCE CHANGED (q-pair).....	423
APPROPRIATE (data).....	832
APPROPRIATE (q-pair).....	32, 54
ARM BAND - ALLERGIES.....	33
ARM BAND - BLOOD TRANSFUSION IV.....	34
ARM BAND - DNR.....	35
ARM BAND - ELOPMENT.....	36
ARM BAND - HIGH FALL RISK.....	37
ARM BAND - PERSONAL MEDICAL ALERT.....	38
ARM BAND - WANDERER.....	39
ARM BANDS PRESENT.....	40
ARRHYTHMIA (data).....	832
ARRHYTHMIA (q-pair).....	398
ARTERIAL LINE (meth).....	1180
ARTERIAL LINE (q-pair).....	63, 149, 377, 378,
	690
ARTERIAL LINE ALARMS ON.....	41
ARTERIAL LINE CARE (data).....	832
ARTERIAL LINE CARE (q-pair).....	343
ARTERIAL LINE CIRCULATION ASSESSED.....	42
ARTERIAL LINE DRESSING CHANGE.....	43
ARTERIAL LINE DRESSING TYPE.....	44
ARTERIAL LINE SITE CARE.....	45
ARTERIAL LINE SITE CONDITION.....	46
ARTERIAL LINE SIZE.....	47

ARTERIAL LINE WAVEFORM.....	48
ARTERIAL LINE ZEROED.....	49
ARTIFICIAL AIRWAY (data).....	832
ARTIFICIAL AIRWAY (q-pair).....	567
ARTIFICIAL VALVES (data).....	833
ARTIFICIAL VALVES (q-pair).....	396
ASA PHYSICAL STATUS.....	50
ASA PS 1 (data).....	833
ASA PS 1 (q-pair).....	50
ASA PS 2 (data).....	833
ASA PS 2 (q-pair).....	50
ASA PS 3 (data).....	833
ASA PS 3 (q-pair).....	50
ASA PS 4 (data).....	833
ASA PS 4 (q-pair).....	50
ASA PS 5 (data).....	834
ASA PS 5 (q-pair).....	50
ASEPTIC TECHNIQUE FOLLOWED (data).....	834
ASEPTIC TECHNIQUE FOLLOWED (q-pair).....	91
ASHEN (data).....	834
ASHEN (q-pair).....	645
ASLEEP AT ROUNDS (data).....	834
ASLEEP AT ROUNDS (q-pair).....	664
ASPIRATION (data).....	834
ASPIRATION (q-pair).....	212
ASSAULTIVE (data).....	835
ASSAULTIVE (q-pair).....	58, 130
ASSESSMENT - CO (task).....	1231
ASSESSMENT - COLLECT VITALS (task).....	1232
ASSESSMENT - DOPPLERS (task).....	1232
ASSESSMENT - FINGERSTICK BLOOD SUGAR (task).....	1232
ASSESSMENT - HEMODYNAMIC READINGS (task).....	1232
ASSESSMENT - I&O (task).....	1232
ASSESSMENT - NEURO CHECKS (task).....	1233
ASSESSMENT - NEUROVASCULAR CHECKS (task).....	1233
ASSESSMENT - PCA (task).....	1233
ASSESSMENT - PULSE OX (task).....	1233
ASSESSMENT - RESTRAINTS (task).....	1233
ASSESSMENT - TELEMETRY (task).....	1234
ASSESSMENT - WEIGHT (task).....	1234
ASSIGNED ONE TO ONE SITTER (data).....	835
ASSIGNED ONE TO ONE SITTER (q-pair).....	4
ASSIST CONTROL (data).....	835
ASSIST CONTROL (q-pair).....	725
ASSISTED (data).....	835
ASSISTED (q-pair).....	13, 15
ASSISTED VENTILATOR (meth).....	1180
ASSISTED VENTILATOR (q-pair).....	578, 686
ASYSTOLE (data).....	835
ASYSTOLE (q-pair).....	84
AT REST (meth).....	1180
AT REST (q-pair).....	63, 570, 578
ATRIAL (data).....	836
ATRIAL (q-pair).....	461
ATRIAL FIB/FLUTTER (data).....	836
ATRIAL FIB/FLUTTER (q-pair).....	84
ATRIAL FIBRILLATION (data).....	836
ATRIAL FIBRILLATION (q-pair).....	84
ATRIAL FLUTTER (data).....	836
ATRIAL FLUTTER (q-pair).....	84
AUSCULTATE (meth).....	1180
AUSCULTATE (q-pair).....	570

AUSCULTATION (data).....	836
AUSCULTATION (q-pair).....	212
AV FISTULA (loca).....	1128
AV FISTULA (q-pair).....	570
AVF (data).....	837
AVF (q-pair).....	85
AVI (data).....	837
AVI (q-pair).....	455
AVL (data).....	837
AVL (q-pair).....	85
AVO (data).....	837
AVO (q-pair).....	455
AVOIDED LIGHT AND NOISE OVERSTIMULATION (data).....	837
AVOIDED LIGHT AND NOISE OVERSTIMULATION (q-pair).....	190
AVR (data).....	838
AVR (q-pair).....	85
AWAKE (data).....	838
AWAKE (q-pair).....	490
AXILLARY (loca).....	1129
AXILLARY (q-pair).....	694

- B -

BABINSKI.....	51
BACK (loca).....	1129
BACK (q-pair).....	473, 658
BACTERICIDAL (data).....	838
BACTERICIDAL (q-pair).....	548
BACTERIOSTATIC (data).....	838
BACTERIOSTATIC (q-pair).....	548
BACTERIOSTATIC WATER (prod).....	1206
BACTERIOSTATIC WATER (q-pair).....	316
BAG - CONTINUOUS (data).....	838
BAG - CONTINUOUS (q-pair).....	167
BAG - INTERMITTENT (data).....	839
BAG - INTERMITTENT (q-pair).....	167
BAG CHANGE (data).....	839
BAG CHANGE (q-pair).....	170
BAG CHANGED (data).....	839
BAG CHANGED (q-pair).....	343, 493
BAG/MASK & O2 SET UP (data).....	839
BAG/MASK & O2 SET UP (q-pair).....	400
BALSAM OF PERU (data).....	839
BALSAM OF PERU (q-pair).....	758
BARIATRIC FRAME (data).....	840
BARIATRIC FRAME (q-pair).....	665
BARIATRIC LOW AIRLOSS MATTRESS REPLACEMENT (data).....	840
BARIATRIC LOW AIRLOSS MATTRESS REPLACEMENT (q-pair).....	665
BARIATRIC PATIENT (data).....	840
BARIATRIC PATIENT (q-pair).....	545
BARRIERS TO LEARNING.....	52
BASAL RATE ONLY (data).....	840
BASAL RATE ONLY (q-pair).....	495
BATHROOM (data).....	840
BATHROOM (q-pair).....	55, 129
BATHROOM PRIVILEGES (data).....	841
BATHROOM PRIVILEGES (q-pair).....	12
BEATS PER MINUTE (unit).....	1244
BED (data).....	841
BED (meth).....	1180
BED (q-pair).....	746
BED EXIT ALARM (data).....	841

BED EXIT ALARM (q-pair).....	24, 582
BED POSITION.....	53
BEDREST (data).....	841
BEDREST (q-pair).....	12
BEDSIDE (data).....	841
BEDSIDE (q-pair).....	553
BEDSIDE COMMODOE (data).....	842
BEDSIDE COMMODOE (q-pair).....	88
BEHAVIOR.....	54
BEHAVIOR CURRENT LOCATION.....	55
BEHAVIOR OBSERVED.....	56
BEHAVIOR PRECAUTION LEVEL.....	57
BEHAVIOR PRECAUTION REASON.....	58
BIER BLOCK (data).....	842
BIER BLOCK (q-pair).....	31
BIGEMINY (data).....	842
BIGEMINY (q-pair).....	80
BILATERAL (data).....	842
BILATERAL ANKLES (loca).....	1129
BILATERAL ANKLES (q-pair).....	598, 600
BILATERAL FOOT (loca).....	1129
BILATERAL FOOT (q-pair).....	196
BILATERAL HAND (loca).....	1129
BILATERAL HAND (q-pair).....	222
BILATERAL LEG (loca).....	1130
BILATERAL LEG (q-pair).....	352, 598
BILATERAL NECK (loca).....	1130
BILATERAL NECK (q-pair).....	160
BILATERAL PERIPHERALS (loca).....	1130
BILATERAL PERIPHERALS (q-pair).....	570
BILATERAL WRIST (loca).....	1130
BILATERAL WRIST (q-pair).....	598, 599, 712
BIPAP (meth).....	1181
BIPAP (q-pair).....	687
BIPAP - RESPIRONICS (data).....	842
BIPAP - RESPIRONICS (q-pair).....	739
BIPAP - VISION (data).....	843
BIPAP - VISION (q-pair).....	739
BIPAP-CPAP (meth).....	1181
BIPAP-CPAP (q-pair).....	687
BIPAP-EXPIRATORY - DO NOT CHECK TOTAL OR SUBTOTALS....	59
BIPAP-INSPIRATORY - DO NOT CHECK TOTAL OR SUBTOTALS....	60
BIPAP/CPAP (data).....	843
BIPAP/CPAP (q-pair).....	739
BIPHASIC POSITIVE AIRWAY PRESSURE (data).....	843
BIPHASIC POSITIVE AIRWAY PRESSURE (q-pair).....	725, 739
BISMUTH GAUZE (data).....	843
BISMUTH GAUZE (q-pair).....	758
BITE (data).....	843
BITE (q-pair).....	766
BITE BLOCK (data).....	844
BITE BLOCK (q-pair).....	22
BLACK (data).....	844
BLACK (q-pair).....	157, 173, 211, 681, 714, 749, 753
BLACK/NECROTIC (data).....	844
BLACK/NECROTIC (q-pair).....	537
BLADDER (loca).....	1130
BLADDER (meth).....	1181
BLADDER (q-pair).....	316, 473
BLADDER ASSESSMENT.....	61
BLADDER IRRIGATION (prod).....	1206

BLADDER IRRIGATION (q-pair).....	430
BLADDER SCAN (data).....	844
BLADDER SCAN (q-pair).....	88
BLADDER SCAN - DO NOT CHECK TOTAL OR SUBTOTALS.....	62
BLANCHED (data).....	844
BLANCHED (q-pair).....	46
BLAND (data).....	845
BLAND (q-pair).....	150
BLISTERING (data).....	845
BLISTERING (q-pair).....	766
BLOOD GLUCOSE MONITORING (data).....	845
BLOOD GLUCOSE MONITORING (q-pair).....	164
BLOOD PRESSURE.....	63
BLOOD PRESSURE (data).....	845
BLOOD PRESSURE (q-pair).....	402
BLOOD PRESSURE CUFF SIZE.....	64
BLOOD PRODUCT (prod).....	1206
BLOOD PRODUCT (q-pair).....	309
BLOOD STREAKED (data).....	845
BLOOD TRANSFUSION IV (data).....	846
BLOOD TRANSFUSION IV (q-pair).....	40
BLOODY (data).....	846
BLOODY (q-pair).....	109, 296, 539, 680
BLUE (data).....	846
BLUE (q-pair).....	753
BLUNTED (data).....	846
BLUNTED (q-pair).....	17
BLURRED VISION (data).....	846
BLURRED VISION (q-pair).....	183
BMI>35KG/M2 (data).....	847
BMI>35KG/M2 (q-pair).....	396
BODY SURFACE AREA.....	65
BONE (data).....	847
BOOST (data).....	847
BOOST (q-pair).....	172
BOOST PLUS (data).....	847
BOOST PLUS (q-pair).....	172
BORDERLINE REDUCED (data).....	847
BORDERLINE REDUCED (q-pair).....	348
BOTH ARMS (loca).....	1131
BOTH ARMS (q-pair).....	599
BOTH HANDS (loca).....	1131
BOTH HANDS (q-pair).....	160
BOTH LEGS (loca).....	1131
BOTH LEGS (q-pair).....	136, 160, 600
BOWEL MOVEMENT EVENT(S).....	66
BRACHIAL (loca).....	1131
BRACHIAL (q-pair).....	63, 570
BRACHIAL (qual).....	1223
BRADEN FRICTION AND SHEAR 1 PROBLEM (data).....	848
BRADEN FRICTION AND SHEAR 1 PROBLEM (q-pair).....	68
BRADEN SCALE (unit).....	1244
BRADEN SCALE ACTIVITY.....	67
BRADEN SCALE ACTIVITY 1 BEDFAST (data).....	848
BRADEN SCALE ACTIVITY 1 BEDFAST (q-pair).....	67
BRADEN SCALE ACTIVITY 2 CHAIRFAST (data).....	848
BRADEN SCALE ACTIVITY 2 CHAIRFAST (q-pair).....	67
BRADEN SCALE ACTIVITY 3 WALKS OCCASIONALLY (data).....	848
BRADEN SCALE ACTIVITY 3 WALKS OCCASIONALLY (q-pair).....	67
BRADEN SCALE ACTIVITY 4 WALKS FREQUENTLY (data).....	848
BRADEN SCALE ACTIVITY 4 WALKS FREQUENTLY (q-pair).....	67
BRADEN SCALE FRICTION AND SHEAR.....	68

BRADEN SCALE FRICTION AND SHEAR 2 POTENTIAL PROBLEM (data)	849
BRADEN SCALE FRICTION AND SHEAR 2 POTENTIAL PROBLEM (q-pair)	68
BRADEN SCALE FRICTION AND SHEAR 3 NO APPARENT PROBLEM (data)	849
BRADEN SCALE FRICTION AND SHEAR 3 NO APPARENT PROBLEM (q-pair)	68
BRADEN SCALE MOBILITY	69
BRADEN SCALE MOBILITY 1 COMPLETELY IMMOBILE (data)	849
BRADEN SCALE MOBILITY 1 COMPLETELY IMMOBILE (q-pair)	69
BRADEN SCALE MOBILITY 2 VERY LIMITED (data)	849
BRADEN SCALE MOBILITY 2 VERY LIMITED (q-pair)	69
BRADEN SCALE MOBILITY 3 SLIGHTLY LIMITED (data)	849
BRADEN SCALE MOBILITY 3 SLIGHTLY LIMITED (q-pair)	69
BRADEN SCALE MOBILITY 4 NO LIMITATION (data)	850
BRADEN SCALE MOBILITY 4 NO LIMITATION (q-pair)	69
BRADEN SCALE MOISTURE	70
BRADEN SCALE MOISTURE 1 CONSTANTLY MOIST (data)	850
BRADEN SCALE MOISTURE 1 CONSTANTLY MOIST (q-pair)	70
BRADEN SCALE MOISTURE 2 VERY MOIST (data)	850
BRADEN SCALE MOISTURE 2 VERY MOIST (q-pair)	70
BRADEN SCALE MOISTURE 3 OCCASIONALLY MOIST (data)	850
BRADEN SCALE MOISTURE 3 OCCASIONALLY MOIST (q-pair)	70
BRADEN SCALE MOISTURE 4 RARELY MOIST (data)	850
BRADEN SCALE MOISTURE 4 RARELY MOIST (q-pair)	70
BRADEN SCALE NUTRITION	71
BRADEN SCALE NUTRITION 1 VERY POOR (data)	851
BRADEN SCALE NUTRITION 1 VERY POOR (q-pair)	71
BRADEN SCALE NUTRITION 2 PROBABLY INADEQUATE (data)	851
BRADEN SCALE NUTRITION 2 PROBABLY INADEQUATE (q-pair)	71
BRADEN SCALE NUTRITION 3 ADEQUATE (data)	851
BRADEN SCALE NUTRITION 3 ADEQUATE (q-pair)	71
BRADEN SCALE NUTRITION 4 EXCELLENT (data)	851
BRADEN SCALE NUTRITION 4 EXCELLENT (q-pair)	71
BRADEN SCALE SENSORY PERCEPTION	72
BRADEN SCALE SENSORY PERCEPTION 1 COMPLETELY LIMITED (data)	851
BRADEN SCALE SENSORY PERCEPTION 1 COMPLETELY LIMITED (q-pair)	72
BRADEN SCALE SENSORY PERCEPTION 2 VERY LIMITED (data)	852
BRADEN SCALE SENSORY PERCEPTION 2 VERY LIMITED (q-pair)	72
BRADEN SCALE SENSORY PERCEPTION 3 SLIGHTLY LIMITED (data)	852
BRADEN SCALE SENSORY PERCEPTION 3 SLIGHTLY LIMITED (q-pair)	72
BRADEN SCALE SENSORY PERCEPTION 4 NO IMPAIRMENT (data)	852
BRADEN SCALE SENSORY PERCEPTION 4 NO IMPAIRMENT (q-pair)	72
BRADEN SCALE TOTAL	73
BREAKDOWN (data)	852
BREAKDOWN (q-pair)	646
BREATH SOUNDS	74
BREATHS PER MINUTE (unit)	1244
BRIDGE PRESENT (data)	852
BRIDGE PRESENT (q-pair)	423
BRIGHT (data)	853
BRIGHT (q-pair)	17
BRISK (data)	853
BRISK (q-pair)	572
BROWN (data)	853
BROWN (q-pair)	157, 173, 211, 563,



	681, 714, 753
BRUISED (data).....	853
BRUISED (q-pair).....	645, 766
BRUIT (data).....	853
BRUIT (q-pair).....	1
BUCKS TRACTION (data).....	854
BUCKS TRACTION (q-pair).....	702
BULB (data).....	854
BULB (q-pair).....	156
BULB SUCTION (data).....	854
BULB SUCTION (q-pair).....	755
BULLET (data).....	854
BULLET (q-pair).....	766
BUPIVICAINE (data).....	854
BUPIVICAINE (q-pair).....	181
BURN (data).....	855
BURN (q-pair).....	766
BURNING (data).....	855
BURNING (q-pair).....	470
BY PROVIDER ORDER (data).....	855
BY PROVIDER ORDER (q-pair).....	445

- C -

C1 (data).....	855
C1 (q-pair).....	179, 180
C2 (data).....	855
C2 (q-pair).....	179, 180
C3 (data).....	856
C3 (q-pair).....	179, 180
C4 (data).....	856
C4 (q-pair).....	179, 180
C5 (data).....	856
C5 (q-pair).....	179, 180
CADEXOMER IODINE (data).....	856
CADEXOMER IODINE (q-pair).....	548
CALCIUM ALGINATE (data).....	856
CALCIUM ALGINATE (q-pair).....	548, 758
CALF (loca).....	1131
CALF (q-pair).....	116
CALL HOUSE OFFICER - BP DIASTOLIC (task).....	1234
CALL HOUSE OFFICER - BP SYSTOLIC (task).....	1234
CALL HOUSE OFFICER - PULSE (task).....	1234
CALL HOUSE OFFICER - RESPIRATORY RATE (task).....	1235
CALL HOUSE OFFICER - TEMP (task).....	1235
CALL LIGHT WITHIN REACH.....	75
CALM (data).....	857
CALM (q-pair).....	56
CALORIE (unit).....	1244
CANE (data).....	857
CANE (q-pair).....	6, 707
CANES (data).....	857
CANES (q-pair).....	707
CAPILLARY REFILL.....	76
CARDIAC (data).....	857
CARDIAC (q-pair).....	402
CARDIAC ABNORMAL BPM.....	77
CARDIAC ALARM PARAMETERS VERIFIED.....	78
CARDIAC BRUIT.....	79
CARDIAC CHAIR (data).....	857
CARDIAC CHAIR (q-pair).....	53
CARDIAC DIET (data).....	858

CARDIAC DIET (q-pair).....	150
CARDIAC ECTOPY.....	80
CARDIAC HEART SOUNDS.....	81
CARDIAC INDEX.....	82
CARDIAC MONITOR (meth).....	1181
CARDIAC OUTPUT.....	83
CARDIAC OUTPUT INJECTATE (data).....	858
CARDIAC OUTPUT INJECTATE (prod).....	1206
CARDIAC OUTPUT INJECTATE (q-pair).....	430
CARDIAC OUTPUT MONITOR (meth).....	1181
CARDIAC OUTPUT MONITOR (q-pair).....	105, 349, 353, 377, 558, 560, 569, 604, 605, 684, 685, 689
CARDIAC REPOSITION (data).....	858
CARDIAC REPOSITION (q-pair).....	13
CARDIAC RHYTHM.....	84
CARDIAC RHYTHM LEAD.....	85
CARDIAC THRILL.....	86
CARDIAC THRILL TIMING.....	87
CARDIO-RESP MONITOR (data).....	858
CARDIO-RESP MONITOR (q-pair).....	400
CARDIOLOGY PROCEDURE ROOM (data).....	858
CARDIOLOGY PROCEDURE ROOM (q-pair).....	553
CAREGIVER (data).....	859
CAREGIVER (q-pair).....	11
CAREGIVER STRESS (data).....	859
CAREGIVER STRESS (q-pair).....	164
CAROTID (loca).....	1132
CAROTID (q-pair).....	570
CASS ENDOTRACHEAL TUBE (data).....	859
CASS ENDOTRACHEAL TUBE (q-pair).....	22
CATH SECURED WITH DRAIN BAG BELOW PT'S BLADDER LEVEL (data) .....	859
CATH SECURED WITH DRAIN BAG BELOW PT'S BLADDER LEVEL (q-pair) .....	91
CATH SIZE DOCUMENTED (data).....	859
CATH SIZE DOCUMENTED (q-pair).....	91
CATHETER (data).....	860
CATHETER (meth).....	1182
CATHETER (q-pair).....	220
CATHETER (sche).....	1226
CATHETERIZATION TIME DOCUMENTED (data).....	860
CATHETERIZATION TIME DOCUMENTED (q-pair).....	91
CAUTI: ALTERNATIVES ATTEMPTED PRIOR TO CATH.....	88
CAUTI: CATHETER IN PLACE.....	89
CAUTI: INDICATIONS.....	90
CAUTI: INSERTION REQUIREMENTS.....	91
CAUTION-VITALS MAY NOT DISPLAY IN CORRECT SEQUENCE. REVIEW TIMES CAREFULLY. .....	92
CELLULITIS (data).....	860
CELLULITIS (q-pair).....	766
CENTIMETER (unit).....	1244
CENTIMETERS H2O (unit).....	1245
CENTRAL ECG (data).....	860
CENTRAL ECG (q-pair).....	24
CENTRAL LINE ALARMS ON.....	93
CENTRAL LINE CHANGED.....	94
CENTRAL LINE DRESSING CHANGE.....	95
CENTRAL LINE DRESSING TYPE.....	96
CENTRAL LINE INSERTION MARK.....	97
CENTRAL LINE SITE CARE.....	98
CENTRAL LINE SITE CONDITION.....	99

CENTRAL LINE SIZE FRENCH.....	100
CENTRAL LINE SIZE GAUGE.....	101
CENTRAL LINE TYPE.....	102
CENTRAL LINE WAVEFORM.....	103
CENTRAL LINE ZEROED.....	104
CENTRAL VENOUS CATHETER (data).....	860
CENTRAL VENOUS CATHETER (q-pair).....	343, 493
CENTRAL VENOUS PRESSURE.....	105
CEREBRAL PERFUSION PRESSURE - DO NOT CHECK TOTAL OR SUBTOTALS .....	106
CERVICAL SPINAL PRECAUTIONS (data).....	861
CERVICAL SPINAL PRECAUTIONS (q-pair).....	12
CF: COUNTING FINGERS (data).....	861
CF: COUNTING FINGERS (q-pair).....	152
CHAIR (data).....	861
CHAIR (meth).....	1182
CHAIR (q-pair).....	12, 707, 746
CHANGED (data).....	861
CHANGED (q-pair).....	43, 95, 343
CHART VERIFIED CORRECT (data).....	861
CHART VERIFIED CORRECT (q-pair).....	509
CHEST (loca).....	1132
CHEST (q-pair).....	165, 473, 598, 601, 658
CHEST EXPANSION.....	107
CHEST PAIN.....	108
CHEST TUBE (meth).....	1182
CHEST TUBE (prod).....	1207
CHEST TUBE (q-pair).....	431
CHEST TUBE DRAINAGE (prod).....	1207
CHEST TUBE DRAINAGE CHARACTERISTICS.....	109
CHEST TUBE DRAINAGE SUCTION TYPE.....	110
CHEST TUBE DRAINAGE SYSTEM REPLACED.....	111
CHEST TUBE DRESSING CHANGE.....	112
CHEST TUBE SITE ASSESSMENT.....	113
CHEST TUBE SUCTION (CM H2O).....	114
CHEYNE-STOKES (data).....	862
CHEYNE-STOKES (q-pair).....	577
CHIEF CONCERN/COMPLAINT.....	115
CIRCULATION CHECK (data).....	862
CIRCULATION CHECK (q-pair).....	597
CIRCUMFERENCE/GIRTH.....	116
CIWA SCALE (unit).....	1245
CIWA: AGITATION.....	117
CIWA: ANXIETY.....	118
CIWA: AUDITORY DISTURBANCES.....	119
CIWA: HEADACHE FULLNESS IN HEAD.....	120
CIWA: NAUSEA/VOMITING.....	121
CIWA: ORIENTATION AND CLOUDING OF SENSORIUM.....	122
CIWA: PAROXYSMAL SWEATS.....	123
CIWA: TACTILE DISTURBANCES.....	124
CIWA: TOTAL.....	125
CIWA: TREMOR.....	126
CIWA: VISUAL DISTURBANCES.....	127
CLAMPED (data).....	862
CLAMPED (q-pair).....	214
CLASS I (data).....	862
CLASS I (q-pair).....	354
CLASS II (data).....	862
CLASS II (q-pair).....	354
CLASS III (data).....	863
CLASS III (q-pair).....	354

CLASS IV (data).....	863
CLASS IV (q-pair).....	354
CLAY (data).....	863
CLAY (q-pair).....	157, 173, 211, 714, 753
CLEAN (data).....	863
CLEAN (q-pair).....	46, 99, 143, 238, 345, 447
CLEAR (data).....	863
CLEAR (q-pair).....	74, 157, 173, 211, 296, 539, 563, 636, 681, 713, 714
CLEAR LIQUID DIET (data).....	864
CLEAR LIQUID DIET (q-pair).....	150
CLEAR LIQUIDS (prod).....	1207
CLEAR LIQUIDS (q-pair).....	329
CLICK (data).....	864
CLICK (q-pair).....	81
CLINICAL ASSESSMENT (data).....	864
CLINICAL ASSESSMENT (q-pair).....	164
CLOSE OBSERVATION (data).....	864
CLOSE OBSERVATION (q-pair).....	57
CLOSE OBSERVATION BEHAVIOR.....	128
CLOSE OBSERVATION CURRENT LOCATION.....	129
CLOSE OBSERVATION REASON FOR PRECAUTION.....	130
CLOSED (data).....	864
CLOSED (q-pair).....	133
CLOSED SYSTEM (data).....	865
CLOSED SYSTEM (q-pair).....	167
CLOTHED SELF (data).....	865
CLOTHED SELF (q-pair).....	5
CLOTHING ROOM (data).....	865
CLOTHING ROOM (q-pair).....	515, 516
CLOTS (data).....	865
CLOTS (q-pair).....	680, 713
CLOUDY (data).....	865
CLOUDY (q-pair).....	296, 539, 636, 713
CLUBBING (data).....	866
CLUBBING (q-pair).....	415
CMS INTACT DISTAL TO SURGICAL SITE (data).....	866
CMS INTACT DISTAL TO SURGICAL SITE (q-pair).....	465
CNS (data).....	866
CNS (q-pair).....	554
CNS DEPRESSION (data).....	866
CNS DEPRESSION (q-pair).....	183
COCCYX (loca).....	1132
COCCYX (q-pair).....	473, 547, 658, 761
CODE CART IS ASSESSIBLE (data).....	866
CODE CART IS ASSESSIBLE (q-pair).....	400
CODE STATUS.....	131
COFFEE (prod).....	1207
COFFEE (q-pair).....	329
COFFEE GROUND (data).....	867
COGNITIVE IMPAIRMENT (data).....	867
COGNITIVE IMPAIRMENT (q-pair).....	52
COLD (data).....	867
COLD (q-pair).....	475, 658
COLLAGEN (data).....	867
COLLAGEN (q-pair).....	548
COLLOID (meth).....	1182
COLLOID (prod).....	1207
COLLOID (q-pair).....	309, 324

COLLOID TYPE.....	132
COLORLESS (data).....	867
COLORLESS (q-pair).....	753
COLOSTOMY (data).....	868
COLOSTOMY (q-pair).....	424
COMA SCORE EYE.....	133
COMA SCORE MOTOR.....	134
COMA SCORE VERBAL.....	135
COMBATIVE (data).....	868
COMBATIVE (q-pair).....	56
COMFORT CARE (data).....	868
COMFORT CARE (q-pair).....	90
COMFORT MEASURES (data).....	868
COMFORT MEASURES (q-pair).....	190
COMMODO (data).....	868
COMMODO (q-pair).....	12
COMMUNICATION BARRIER (data).....	869
COMMUNICATION BARRIER (q-pair).....	161
COMMUNITY RESOURCES (data).....	869
COMMUNITY RESOURCES (q-pair).....	164
COMPLETE BATH (data).....	869
COMPLETE BATH (q-pair).....	5
COMPLIANCE (data).....	869
COMPLIANCE (q-pair).....	164
COMPR. HOSE SEQ. COMP. DEVICE.....	136
COMPUTED TOMOGRAPHY (meth).....	1182
COMPUTED TOMOGRAPHY (q-pair).....	347, 348
CONFUSED (data).....	869
CONFUSED (q-pair).....	135, 161
CONFUSION (data).....	870
CONFUSION (q-pair).....	183, 498
CONNECTED TO DRAINAGE BAG (data).....	870
CONNECTED TO DRAINAGE BAG (q-pair).....	423
CONSISTENT WITH BODY TEMPERATURE (data).....	870
CONSISTENT WITH BODY TEMPERATURE (q-pair).....	658
CONSTANT (data).....	870
CONSTANT (q-pair).....	470, 471, 767
CONSTIPATION (data).....	870
CONSTIPATION (q-pair).....	183, 209, 680
CONSTRICTED (data).....	871
CONSTRICTED (q-pair).....	17
CONTACT (data).....	871
CONTACT (q-pair).....	333
CONTACTS (data).....	871
CONTACTS (q-pair).....	395
CONTEMPLATIVE (data).....	871
CONTEMPLATIVE (q-pair).....	163
CONTINENT (data).....	871
CONTINENT (q-pair).....	220, 716
CONTINUOUS (data).....	872
CONTINUOUS (q-pair).....	87, 156, 214, 216,
	622
CONTINUOUS POSITIVE AIR PRESSURE (data).....	872
CONTINUOUS POSITIVE AIR PRESSURE (q-pair).....	725, 739
CONTINUOUS SUCTION (data).....	872
CONTINUOUS SUCTION (q-pair).....	110
CONTRAST LEFT VENTRICULOGRAPHY (meth).....	1183
CONTRAST LEFT VENTRICULOGRAPHY (q-pair).....	347, 348
CONTROL MODE VENTILATION (data).....	872
CONTROL MODE VENTILATION (q-pair).....	725
CONTROLLED VENTILATOR (meth).....	1183
CONTROLLED VENTILATOR (q-pair).....	578

CONTUSION (data).....	872
CONTUSION (q-pair).....	766
COOL (data).....	873
COOL (q-pair).....	658
COOL CHEM PACK (data).....	873
COOL CHEM PACK (q-pair).....	744
COOLER THAN GENERAL BODY TEMPERATURE (data).....	873
COOLER THAN GENERAL BODY TEMPERATURE (q-pair).....	658
COOLER THAN OPPOSITE BODY PART (data).....	873
COOLER THAN OPPOSITE BODY PART (q-pair).....	658
COOLING BLANKET (data).....	873
COOLING BLANKET (q-pair).....	744
COOLING OFF (data).....	874
COOLING OFF (q-pair).....	745
COOLING ON (data).....	874
COOLING ON (q-pair).....	745
COOPERATIVE (data).....	874
COOPERATIVE (q-pair).....	56
COPIUS (data).....	874
COPIUS (q-pair).....	435, 562
CORE (loca).....	1132
CORE (q-pair).....	694
COUGH.....	137
CPAP.....	138
CPAP (meth).....	1183
CPAP (q-pair).....	687
CRACKLES (data).....	874
CRACKLES (q-pair).....	74
CRAMPING (data).....	875
CRAMPING (q-pair).....	209, 470
CREPITUS (data).....	875
CREPITUS (q-pair).....	113
CRNA (data).....	875
CRNA (q-pair).....	554
CRRT (meth).....	1183
CRRT (q-pair).....	311, 429
CRUSHING (data).....	875
CRUSHING (q-pair).....	470, 766
CRUTCH (data).....	875
CRUTCH (q-pair).....	707
CRUTCHES (data).....	876
CRUTCHES (q-pair).....	6, 707
CRYING (data).....	876
CRYING (q-pair).....	474
CRYOPRECIPITATE (prod).....	1208
CRYOPRECIPITATE (q-pair).....	309, 317
CRYSTALLOID (data).....	876
CRYSTALLOID (prod).....	1208
CRYSTALLOID (q-pair).....	309, 324
CSM: CENTRAL STEADY MAINTAINED (data).....	876
CSM: CENTRAL STEADY MAINTAINED (q-pair).....	152
CUBIC CENTIMETER (unit).....	1245
CUFF (meth).....	1183
CUFF (q-pair).....	149, 377, 378, 690
CUFF INFLATION CHECKED.....	139
CUFF-AUTOMATED (meth).....	1184
CUFF-AUTOMATED (q-pair).....	63
CUFF-MANUAL (meth).....	1184
CUFF-MANUAL (q-pair).....	63
CULTURAL (data).....	876
CULTURAL (q-pair).....	162
CVVHD (meth).....	1184



DIABETES (data).....	882
DIABETES (q-pair).....	396
DIABETIC/VASCULAR ULCER (data).....	882
DIABETIC/VASCULAR ULCER (q-pair).....	766
DIAGNOSIS.....	141
DIAGNOSIS TYPE.....	142
DIAGNOSIS/HEALTH STATUS (data).....	882
DIAGNOSIS/HEALTH STATUS (q-pair).....	164
DIALYSIS (data).....	882
DIALYSIS (meth).....	1184
DIALYSIS (q-pair).....	102
DIALYSIS CATHETER (meth).....	1184
DIALYSIS CATHETER SITE CONDITION.....	143
DIALYSIS DRESSING TYPE.....	144
DIALYSIS RATIO.....	145
DIALYSIS: ARTERIAL PRESSURE.....	146
DIALYSIS: UFR (ULTRAFILTRATION RATE).....	147
DIALYSIS: VENOUS PRESSURE.....	148
DIAPHORETIC (data).....	883
DIAPHORETIC (q-pair).....	653
DIARRHEA (data).....	883
DIARRHEA (q-pair).....	209, 680
DIASTOLIC (data).....	883
DIASTOLIC (q-pair).....	87
DIASTOLIC PRESSURE.....	149
DIET.....	150
DIFFICULT TO AROUSE (data).....	883
DIFFICULT TO AROUSE (q-pair).....	490
DILATED (data).....	883
DILATED (q-pair).....	572
DILUTE DAKINS SOLUTION (data).....	884
DILUTE DAKINS SOLUTION (q-pair).....	548
DILUTE SODIUM HYPOCHLORITE SOLUTION (data).....	884
DIMINISHED (data).....	884
DIMINISHED (q-pair).....	74
DISCHARGE (data).....	884
DISCHARGE (q-pair).....	142, 164, 761
DISCHARGE ACTIVITIES.....	151
DISCHARGE MEDS PLANNED (data).....	884
DISCHARGE MEDS PLANNED (q-pair).....	151
DISCHARGE ORDER ON CHART (data).....	885
DISCHARGE ORDER ON CHART (q-pair).....	151
DISCOLORATION (data).....	885
DISCOLORATION (q-pair).....	414, 640
DISCOMFORT (data).....	885
DISCOMFORT (q-pair).....	46, 99, 143, 238, 345, 447
DISCONTINUED (data).....	885
DISCONTINUED (q-pair).....	343, 493, 622
DISEASE PROCESS (data).....	885
DISEASE PROCESS (q-pair).....	164
DISHEVELED (data).....	886
DISHEVELED (q-pair).....	32
DISORIENTED (data).....	886
DISORIENTED (q-pair).....	594
DISORIENTED TO PERSON (data).....	886
DISORIENTED TO PERSON (q-pair).....	421
DISORIENTED TO PLACE (data).....	886
DISORIENTED TO PLACE (q-pair).....	421
DISORIENTED TO SITUATION (data).....	886
DISORIENTED TO SITUATION (q-pair).....	421
DISORIENTED TO TIME (data).....	887



DISORIENTED TO TIME (q-pair).....	421
DISTAL PORT (data).....	887
DISTAL PORT (q-pair).....	344
DISTANCE VISION.....	152
DISTENDED (data).....	887
DISTENDED (q-pair).....	2, 417
DIVERSION (data).....	887
DIVERSION (q-pair).....	716
DO (data).....	887
DO (q-pair).....	608
DO NOT PROCESS.....	153
DO NOT RESUSCITATE (data).....	888
DO NOT RESUSCITATE (q-pair).....	40, 131
DONE (data).....	888
DONE (q-pair).....	15
DONOR SITE (loca).....	1132
DONOR SITE (q-pair).....	303
DOO (data).....	888
DOO (q-pair).....	455
DOPPLER (meth).....	1185
DOPPLER (q-pair).....	63, 149, 570, 690
DORSALIS PEDIS (loca).....	1133
DORSALIS PEDIS (q-pair).....	570
DORSI FLEXION.....	154
DOUBLE (data).....	888
DOUBLE (q-pair).....	102
DOUBLE LUMEN CATHETER (meth).....	1185
DOWN (data).....	888
DRAIN SUCTION CM H2O.....	155
DRAIN SUCTION TYPE.....	156
DRAINAGE CHARACTERISTICS.....	157
DRAINAGE FLOW IS UNOBSTRUCTED (data).....	889
DRAINAGE FLOW IS UNOBSTRUCTED (q-pair).....	91
DRESSING (data).....	889
DRESSING (q-pair).....	572, 750
DRESSING CHANGED (data).....	889
DRESSING CHANGED (q-pair).....	343, 493
DRESSING WITH PRESCRIPTION RX (data).....	889
DRESSINGS CHECKED (data).....	889
DRESSINGS CHECKED (q-pair).....	151
DROPLET (data).....	890
DROPLET (q-pair).....	333
DROWSY (data).....	890
DROWSY (q-pair).....	56, 351, 490
DRY (data).....	890
DRY (q-pair).....	43, 95, 638, 653, 746
DRY (qual).....	1223
DRY AND INTACT (data).....	890
DRY SHAMPOO (data).....	890
DRY SHAMPOO (q-pair).....	5
DRY STERILE DRESSING (data).....	891
DULL (data).....	891
DULL (q-pair).....	470
DUMPING SYNDROME (data).....	891
DUMPING SYNDROME (q-pair).....	150
DUODENOSTOMY (data).....	891
DUODENOSTOMY (q-pair).....	166
DUSKY (data).....	891
DUSKY (q-pair).....	645
DVI (data).....	892
DVI (q-pair).....	455

DVT PROPHYLAXIS.....	158
DYE PRESENT (data).....	892
DYE PRESENT (q-pair).....	714
DYNE * SECONDS / CM(-5) * M2 (unit).....	1246
DYNE * SECONDS / CM(-5) / M2 (unit).....	1246
DYNE-SEC/CM5 (unit).....	1246
DYSPHAGIA PHASE I (data).....	892
DYSPHAGIA PHASE I (q-pair).....	150
DYSPHAGIA STAGE II (data).....	892
DYSPHAGIA STAGE II (q-pair).....	150
DYSPHAGIA STAGE III (data).....	892
DYSPHAGIA STAGE III (q-pair).....	150
DYSPHAGIA STAGE IV (data).....	893
DYSPHAGIA STAGE IV (q-pair).....	150
DYSPHAGIA STAGE V (data).....	893
DYSPHAGIA STAGE V (q-pair).....	150
DYSPNEIC (data).....	893
DYSPNEIC (q-pair).....	577
DYSURIA (data).....	893
DYSURIA (q-pair).....	716

- E -

EAGER TO LEARN (data).....	893
EAGER TO LEARN (q-pair).....	163
EASY TO AROUSE (data).....	894
EASY TO AROUSE (q-pair).....	490
ECCHYMOSIS (data).....	894
ECCHYMOSIS (q-pair).....	645
ECHOCARDIOGRAPHY (meth).....	1185
ECHOCARDIOGRAPHY (q-pair).....	347, 348
EDEMA AMOUNT.....	159
EDEMA SEVERITY.....	160
EDUCATION - BARRIERS.....	161
EDUCATION - BELIEFS.....	162
EDUCATION - READINESS TO LEARN.....	163
EDUCATION - TOPIC.....	164
EF RESIDUAL (prod).....	1208
EF RESIDUAL (q-pair).....	430
EKG LEAD.....	165
ELASTIC WRAP (data).....	894
ELASTIC WRAP (q-pair).....	158
ELBOW PAD (data).....	894
ELBOW PAD (q-pair).....	544
ELECTRIC (data).....	894
ELECTRIC (q-pair).....	470
ELECTROCARDIOGRAM (meth).....	1185
ELECTROCARDIOGRAM (q-pair).....	570, 667, 668, 669, 670, 671, 672, 673, 674, 675, 676, 677, 678, 699, 700
ELECTROCARDIOGRAM TRACE 1 (meth).....	1185
ELECTROCARDIOGRAM TRACE 1 (q-pair).....	570
ELECTROCARDIOGRAM TRACE 2 (meth).....	1186
ELECTROCARDIOGRAM TRACE 2 (q-pair).....	570
ELEMENTAL FORMULA (data).....	895
ELEMENTAL FORMULA (prod).....	1208
ELEMENTAL FORMULA (q-pair).....	309
ELIMINATED ENVIRONMENTAL HAZARDS (data).....	895
ELIMINATED ENVIRONMENTAL HAZARDS (q-pair).....	190
ELOPEMENT (data).....	895
ELOPEMENT (q-pair).....	40, 58, 130

EMERGENCY (data).....	895
EMERGENCY (q-pair).....	50
EMERGENCY DEPARTMENT (data).....	895
EMERGENCY DEPARTMENT (q-pair).....	187, 467
EMERGENCY DEPT (data).....	896
EMERGENCY DEPT (q-pair).....	553
EMESIS (prod).....	1209
EMESIS (q-pair).....	432
EMOTIONAL (data).....	896
EMOTIONAL (q-pair).....	161
END OF LIFE CARE (data).....	896
END OF LIFE CARE (q-pair).....	545
END-TIDAL CO2 MONITOR (data).....	896
END-TIDAL CO2 MONITOR (q-pair).....	400
ENSURE PLUS (data).....	896
ENSURE PLUS (q-pair).....	172
ENTERAL FEEDING / SUPPLEMENT (meth).....	1186
ENTERAL FEEDING / SUPPLEMENT (q-pair).....	309
ENTERAL FEEDING FORMULA/WATER ADMINISTRATION (data)....	897
ENTERAL FEEDING FORMULA/WATER ADMINISTRATION (q-pair)..	164
ENTERAL FEEDING PUMP (data).....	897
ENTERAL FEEDING PUMP (q-pair).....	164
ENTERAL FEEDING TUBE (meth).....	1186
ENTERAL FEEDING TUBE (q-pair).....	309
ENTERAL FEEDING TUBE TYPE.....	166
ENTERAL NUTRITION DELIVERY METHOD.....	167
ENTERAL NUTRITION LINE TRACED.....	168
ENTERAL NUTRITION RATE.....	169
ENTERAL NUTRITION SYSTEM MAINTENANCE.....	170
ENTERAL NUTRITION TUBE PLACEMENT VERIFIED.....	171
ENTERAL NUTRITION TYPE.....	172
ENTERAL TUBE ASPIRATE COLOR.....	173
ENTERAL TUBE ASPIRATE PH.....	174
ENTERAL TUBE MEDICATION ADMINISTRATION (data).....	897
ENTERAL TUBE MEDICATION ADMINISTRATION (q-pair).....	164
ENTERAL TUBE RESIDUAL VOLUME.....	175
ENTERAL TUBE SITE CONDITION CHECKED.....	176
ENTERAL TUBE/STOMA CARE (data).....	897
ENTERAL TUBE/STOMA CARE (q-pair).....	164
ENVIRONMENT (data).....	897
ENVIRONMENT (q-pair).....	164
ENZYMATIC OINTMENT (data).....	898
ENZYMATIC OINTMENT (q-pair).....	548
EPICARDIAL LEADS (data).....	898
EPICARDIAL LEADS (q-pair).....	461
EPIDURAL (data).....	898
EPIDURAL (q-pair).....	31
EPIDURAL PUMP (data).....	898
EPIDURAL PUMP (q-pair).....	164
EPIDURAL: BOLUS (ML).....	177
EPIDURAL: CONTINUOUS RATE (ML/HR).....	178
EPIDURAL: DERMATOME SENSORY LEVEL - LEFT.....	179
EPIDURAL: DERMATOME SENSORY LEVEL - RIGHT.....	180
EPIDURAL: MEDICATION.....	181
EPIDURAL: RESERVOIR VOLUME REMAINING.....	182
EPIDURAL: SIDE EFFECTS.....	183
EPIDURAL: TOTAL VOLUME GIVEN.....	184
EPITHELIALIZATION (data).....	898
EPITHELIALIZATION (q-pair).....	760
EQUAL (data).....	899
EQUAL (q-pair).....	352, 712
EQUIPMENT (data).....	899

EQUIPMENT (q-pair).....	164
ER EMERGENCY SEVERITY INDEX.....	185
ER TRIAGE CONDITION.....	186
ER TRIAGE DISPOSITION.....	187
ERYTHEMA (data).....	899
ERYTHEMA (q-pair).....	640
ESCAPE (data).....	899
ESCAPE (q-pair).....	58, 130
ESCHAR (data).....	899
ESCHAR (q-pair).....	760
ESTABLISHED ROUTINE TOILETING SCHEDULE (data).....	900
ESTABLISHED ROUTINE TOILETING SCHEDULE (q-pair).....	190
ESTIMATED (q-pair).....	227, 297, 435, 570, 694, 746, 748
ESTIMATED (qual).....	1223
ESTIMATED BY ARM SPAN (q-pair).....	227
ESTIMATED BY ARM SPAN (qual).....	1224
ET LOCATION.....	188
ETT REPOSITIONED (data).....	900
ETT REPOSITIONED (q-pair).....	15
ETT RETAPED/SECURED (data).....	900
ETT RETAPED/SECURED (q-pair).....	15
EVENT(S) (unit).....	1247
EVENTS PER HOUR (unit).....	1247
EWALD TUBE (data).....	900
EWALD TUBE (q-pair).....	215
EXPIRATORY WHEEZES (data).....	900
EXPIRATORY WHEEZES (q-pair).....	74
EXTENDED PERIOD OF IMMOBILITY (data).....	901
EXTENDED PERIOD OF IMMOBILITY (q-pair).....	183
EXTENDS AROUND INCISION EXTENDS INTO PERIWOUND LESS THAN 1 CM (data) .....	901
EXTENDS AROUND INCISION EXTENDS INTO PERIWOUND LESS THAN 1 CM (q-pair) .....	648
EXTENDS AROUND INCISION MORE THAN 1 CM (data).....	901
EXTENDS AROUND INCISION MORE THAN 1 CM (q-pair).....	648
EXTENDS AROUND INCISION NONE (data).....	901
EXTENDS AROUND INCISION NONE (q-pair).....	648
EXTENDS INTO PERIWOUND <1CM (data).....	901
EXTENDS INTO PERIWOUND <1CM (q-pair).....	301
EXTENDS INTO PERIWOUND>1CM (data).....	902
EXTENDS INTO PERIWOUND>1CM (q-pair).....	301
EXTENDS TO PAIN (data).....	902
EXTENDS TO PAIN (q-pair).....	134
EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND LESS THAN 1 CM (data) .....	902
EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND LESS THAN 1 CM (q-pair) .....	649
EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND MORE THAN 1 CM (data) .....	902
EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND MORE THAN 1 CM (q-pair) .....	649
EXTENSION AROUND WOUND LIMITED TO WOUND (data).....	902
EXTENSION AROUND WOUND LIMITED TO WOUND (q-pair).....	649
EXTENSION AROUND WOUND NONE (data).....	903
EXTENSION AROUND WOUND NONE (q-pair).....	649
EXTENSIVE ASSISTANCE (data).....	903
EXTENSIVE ASSISTANCE (q-pair).....	8, 9, 10
EXTERNAL CATHETER (data).....	903
EXTERNAL CATHETER (q-pair).....	219, 716
EXTERNAL CONDOM CATH (data).....	903
EXTERNAL CONDOM CATH (q-pair).....	88

EXTREMELY DRY (data).....	903
EYE CLOSED (data).....	904
EYE CLOSED (q-pair).....	572
EYEGLASSES (data).....	904
EYEGLASSES (q-pair).....	395

- F -

FACE (loca).....	1133
FACE (q-pair).....	160, 473, 658
FACE TENT (meth).....	1186
FACE TENT (q-pair).....	686
FACE WASH (data).....	904
FACE WASH (q-pair).....	5
FALL PRECAUTIONS.....	189
FALL PREVENTION.....	190
FALL RISK ID (data).....	904
FALL RISK ID (q-pair).....	190
FALL RISK SIGN (data).....	904
FALL RISK SIGN (q-pair).....	190
FAMILY (data).....	905
FAMILY (q-pair).....	11, 515, 516
FAMILY INFORMED (data).....	905
FAMILY INFORMED (q-pair).....	403, 580
FAMILY INSTRUCTED (data).....	905
FAMILY INSTRUCTED (q-pair).....	423
FAMILY NEEDS (data).....	905
FAMILY NEEDS (q-pair).....	164
FAMILY ORIENTED TO UNIT.....	191
FAMILY STATED (data).....	905
FAMILY STATED (q-pair).....	509
FAMILY/FRIEND (data).....	906
FAMILY/FRIEND (q-pair).....	706
FASCIA (data).....	906
FAT RESTRICTED 50 GM DIET (data).....	906
FAT RESTRICTED 50 GM DIET (q-pair).....	150
FECAL FAT TEST DIET (data).....	906
FECAL FAT TEST DIET (q-pair).....	150
FECES (prod).....	1209
FEET (loca).....	1133
FEET (q-pair).....	160
FEMORAL (loca).....	1133
FEMORAL (q-pair).....	63, 570
FENTANYL (data).....	906
FENTANYL (q-pair).....	181
FETAL HEART TONES.....	192
FFP (prod).....	1209
FFP (q-pair).....	309, 317
FFP - THAWED (prod).....	1209
FFP - THAWED (q-pair).....	309, 317
FHWN1 (data).....	907
FIBER FORMULA (data).....	907
FIBER STAT (data).....	907
FIBER STAT (q-pair).....	172
FIBRINOUS (data).....	907
FIBRINOUS (q-pair).....	760
FINANCIAL CONCERNS (data).....	907
FINANCIAL CONCERNS (q-pair).....	164
FINGERSTICK GLUCOSE.....	193
FIRM (data).....	908
FIRM (q-pair).....	2
FIRST DEGREE AV BLOCK (data).....	908

FIRST DEGREE AV BLOCK (q-pair).....	84
FIXED (data).....	908
FIXED (q-pair).....	572
FLAT (data).....	908
FLAT (q-pair).....	2, 17, 417
FLATULENCE (data).....	908
FLATULENCE (q-pair).....	1
FLATUS.....	194
FLEXION TO PAIN (data).....	909
FLEXION TO PAIN (q-pair).....	134
FLICKER OF CONTRACTION (data).....	909
FLICKER OF CONTRACTION (q-pair).....	411
FLOTATION BED (data).....	909
FLOTATION BED (q-pair).....	665
FLOWSHEET ADDENDUM.....	195
FLUCTUANCE (data).....	909
FLUCTUATES (data).....	909
FLUCTUATES (q-pair).....	113
FLUSHED (data).....	910
FLUSHED (q-pair).....	343, 645
FOAM AIR HEEL ELEVATOR (data).....	910
FOAM AIR HEEL ELEVATOR (q-pair).....	544
FOAM DRESSING (data).....	910
FOAM DRESSING (q-pair).....	758
FOLEY (data).....	910
FOLEY (q-pair).....	219
FOOT (unit).....	1247
FOOT CARE (data).....	910
FOOT CARE (q-pair).....	5
FOOT DOWN (data).....	911
FOOT DOWN (q-pair).....	634
FOOT STRENGTH.....	196
FOOT UP (data).....	911
FOOT UP (q-pair).....	634
FOREHEAD (loca).....	1133
FOREHEAD (q-pair).....	547
FORMULA - HIGH FIBER (prod).....	1209
FORMULA - PROBALANCE (prod).....	1210
FORMULA - PROBALANCE (q-pair).....	172
FORMULA 1.5 CAL (prod).....	1210
FORMULA 1.5 CAL (q-pair).....	309
FORMULA 2.0 CAL (prod).....	1210
FORMULA 2.0 CAL (q-pair).....	172, 309
FORMULA D10W ENTERAL NUTRITION (prod).....	1210
FORMULA D10W ENTERAL NUTRITION (q-pair).....	172
FORMULA D50W ENTERAL NUTRITION (prod).....	1210
FORMULA D50W ENTERAL NUTRITION (q-pair).....	172
FORMULA D5W ENTERAL NUTRITION (prod).....	1211
FORMULA D5W ENTERAL NUTRITION (q-pair).....	172
FORMULA FIBER (prod).....	1211
FORMULA FIBER (q-pair).....	309
FOUL SMELLING (data).....	911
FOUL SMELLING (q-pair).....	542, 754
FOWLER'S (data).....	911
FOWLER'S (q-pair).....	13
FRANK BLOOD (data).....	911
FRANK BLOOD (q-pair).....	680
FRENCH (unit).....	1247
FREQUENT VERBAL REMINDERS (data).....	912
FREQUENT VERBAL REMINDERS (q-pair).....	586
FRIEND (data).....	912
FRIEND (q-pair).....	11, 515, 516, 517

FROM A NON-VA OUTSIDE FACILITY (data).....	912
FROM A NON-VA OUTSIDE FACILITY (q-pair).....	89
FROM ANOTHER UNIT (data).....	912
FROM ANOTHER UNIT (q-pair).....	89
FROM ANOTHER VA FACILITY (data).....	912
FROM ANOTHER VA FACILITY (q-pair).....	89
FROM CLC (data).....	913
FROM CLC (q-pair).....	89
FROM EMERGENCY DEPARTMENT (data).....	913
FROM EMERGENCY DEPARTMENT (q-pair).....	89
FROM HOME (data).....	913
FROM HOME (q-pair).....	89
FROM INTERVENTIONAL RADIOLOGY (data).....	913
FROM INTERVENTIONAL RADIOLOGY (q-pair).....	89
FROM MICU (data).....	913
FROM MICU (q-pair).....	89
FROM OPERATING ROOM (data).....	914
FROM OPERATING ROOM (q-pair).....	89
FROM SICU (data).....	914
FROM SICU (q-pair).....	89
FRONTAL (loca).....	1134
FULL CODE (data).....	914
FULL CODE (q-pair).....	131
FULL LIQUID (prod).....	1211
FULL LIQUID (q-pair).....	329
FULL LIQUID DIABETIC (data).....	914
FULL LIQUID DIABETIC (q-pair).....	150
FULL LIQUID DIET (data).....	914
FULL LIQUID DIET (q-pair).....	150
FULL LOWER (data).....	915
FULL LOWER (q-pair).....	140
FULL NAME (data).....	915
FULL NAME (q-pair).....	401
FULL SPINAL PRECAUTIONS (data).....	915
FULL SPINAL PRECAUTIONS (q-pair).....	12
FULL SSN (data).....	915
FULL SSN (q-pair).....	401
FULL THICK SKIN LOSS INVOLVING MUSCLE, BONE, SUPPO (data)	
.....	915
FULL THICK SKIN LOSS INVOLVING MUSCLE, BONE, SUPPO (q-pair)	
.....	547
FULL THICK SKIN LOSS INVOLVING SUBCUTANEOUS TISSUE (data)	
.....	916
FULL THICK SKIN LOSS INVOLVING SUBCUTANEOUS TISSUE (q-pair)	
.....	547
FULL UPPER (data).....	916
FULL UPPER (q-pair).....	140
FUNCTIONAL SECTION SET UP (data).....	916
FUNCTIONAL SECTION SET UP (q-pair).....	400
FUNDAL HEIGHT.....	197

- G -

GAG REFLEX.....	198
GAIT.....	199
GAIT BELT (data).....	916
GAIT BELT (q-pair).....	6
GASTRIC BYPASS CLEAR LIQUID (data).....	916
GASTRIC BYPASS CLEAR LIQUID (q-pair).....	150
GASTRIC BYPASS PUREED (data).....	917
GASTRIC BYPASS PUREED (q-pair).....	150
GASTRIC BYPASS REGULAR (data).....	917

GASTRIC BYPASS REGULAR (q-pair).....	150
GASTRIC DRAINAGE (prod).....	1211
GASTRIC DRAINAGE (q-pair).....	430
GASTRIC RESIDUAL AMOUNT DISCARDED.....	200
GASTRIC RESIDUAL CHECKED.....	201
GASTRIC RESIDUAL NOT DISCARDED.....	202
GASTRIC TUBE (data).....	917
GASTRIC TUBE (q-pair).....	166, 215
GASTRIC pH.....	203
GASTROINTESTINAL (loca).....	1134
GASTROINTESTINAL (q-pair).....	473
GASTROSTOMY FEEDING BUTTON (data).....	917
GASTROSTOMY FEEDING BUTTON (q-pair).....	215
GASTROSTOMY TUBE (data).....	917
GASTROSTOMY TUBE (q-pair).....	215
GAUGE (unit).....	1247
GAUZE (data).....	918
GAUZE (q-pair).....	548, 758
GAUZE AND TAPE (data).....	918
GAUZE AND TAPE (q-pair).....	44, 235
GAUZE PACKING (data).....	918
GAUZE PACKING (q-pair).....	548, 758
GAUZE WRAP (data).....	918
GAUZE WRAP (q-pair).....	548, 758
GCS - EYE OPENING.....	204
GCS - MOTOR RESPONSE.....	205
GCS - TOTAL.....	206
GCS - VERBAL RESPONSE.....	207
GENERAL (data).....	918
GENERAL (q-pair).....	31
GENERAL HEALTH (data).....	919
GENERAL HEALTH (q-pair).....	164
GENERAL UMBILICUS (loca).....	1134
GENERAL UMBILICUS (q-pair).....	1
GENERALIZED (loca).....	1134
GENERALIZED (q-pair).....	160, 761
GENITAL ASSESSMENT.....	208
GENITALS (loca).....	1134
GENITALS (q-pair).....	761
GERI-CHAIR (data).....	919
GERI-CHAIR (q-pair).....	12, 128, 601
GI PROCEDURE ROOM (data).....	919
GI PROCEDURE ROOM (q-pair).....	553
GI SYMPTOMS.....	209
GI TUBE (meth).....	1186
GI TUBE DRAINAGE CHARACTERISTICS.....	210
GI TUBE DRAINAGE COLOR.....	211
GI TUBE PLACEMENT CONFIRMATION.....	212
GI TUBE SIZE.....	213
GI TUBE SUCTION TYPE.....	214
GI TUBE TYPE.....	215
GIVEN DISCHARGE INSTRUCTION AND VERBALIZES UNDERSTANDING (data) .....	919
GIVEN DISCHARGE INSTRUCTION AND VERBALIZES UNDERSTANDING (q-pair) .....	151
GLASGOW COMA SCALE (unit).....	1248
GLUCERNA 1.0 (data).....	919
GLUCERNA 1.0 (q-pair).....	172
GLUCERNA 1.2 (data).....	920
GLUCERNA 1.2 (q-pair).....	172
GLUCERNA DIABETIC (prod).....	1211
GLUCERNA DIABETIC (q-pair).....	172



GLUTEN RESTRICTED (data).....	920
GLUTEN RESTRICTED (q-pair).....	150
GM-M/M2/BEAT (unit).....	1248
GNAWING (data).....	920
GNAWING (q-pair).....	470
GOLD (data).....	920
GOLD (q-pair).....	173
GOOD (data).....	920
GOOD (q-pair).....	661
GRAFTS (data).....	921
GRAFTS (q-pair).....	396
GRAMS PER HOUR (unit).....	1248
GRAMS PER MINUTE (unit).....	1248
GRANULATING (data).....	921
GRANULATING (q-pair).....	760
GRANULOCYTES PHERESIS (prod).....	1212
GRANULOCYTES PHERESIS (q-pair).....	309, 317
GRAVITY (data).....	921
GRAVITY (q-pair).....	110, 214
GRAY (data).....	921
GRAY (q-pair).....	749
GREEN (data).....	921
GREEN (q-pair).....	157, 173, 211, 563, 681, 714, 749, 753
GRIMACING (data).....	922
GRIMACING (q-pair).....	474
GROUP ROOM (data).....	922
GROUP ROOM (q-pair).....	55, 129
GU BLADDER IRRIGATION FREQUENCY.....	216
GU CATHETER CARE PERFORMED.....	217
GU CATHETER SIZE.....	218
GU CATHETER TYPE.....	219
GU PROCEDURE ROOM (data).....	922
GU PROCEDURE ROOM (q-pair).....	553
GU VOIDING.....	220
GUARDING (data).....	922
GUARDING (q-pair).....	2, 474
GURNEY CHAIR (data).....	922
GURNEY CHAIR (q-pair).....	12

- H -

HAIR CONDITION.....	221
HALL (data).....	923
HALL (q-pair).....	55, 129
HAND GRIP.....	222
HAND HYGIENE BEFORE AND AFTER (data).....	923
HAND HYGIENE BEFORE AND AFTER (q-pair).....	91
HANDOFF REPORT.....	223
HARD (data).....	923
HARD (q-pair).....	2, 638
HAS NOT SLEPT (data).....	923
HAS NOT SLEPT (q-pair).....	664
HEAD (loca).....	1135
HEAD (q-pair).....	116, 473, 658
HEAD CIRCUMFERENCE.....	224
HEAD DOWN (data).....	923
HEAD DOWN (q-pair).....	634
HEAD OCCIPITAL (loca).....	1135
HEAD OCCIPITAL (q-pair).....	547, 761
HEAD OF BED ELEVATION.....	225
HEAD UP (data).....	924

HEAD UP (q-pair).....	634
HEADACHE (data).....	924
HEADACHE (q-pair).....	183
HEALTH TECH (data).....	924
HEALTH TECH (q-pair).....	554
HEARING (data).....	924
HEARING (q-pair).....	161
HEARING AIDS (data).....	924
HEARING AIDS (q-pair).....	395
HEART (loca).....	1135
HEART (q-pair).....	86
HEARTBURN (data).....	925
HEARTBURN (q-pair).....	209
HEAT (data).....	925
HEAT (q-pair).....	475
HEAT CHEM PACK (data).....	925
HEAT CHEM PACK (q-pair).....	744
HEAT LAMP (data).....	925
HEAT LAMP (q-pair).....	744
HEAT/MOISTURE EXCHANGER PRESENT.....	226
HEAVY (data).....	925
HEAVY (q-pair).....	470, 647
HEEL LIFT (data).....	926
HEEL LIFT (q-pair).....	544
HEEL PILLOW (data).....	926
HEIGHT.....	227
HELMET (data).....	926
HELMET (q-pair).....	190
HEMATOMA (data).....	926
HEMATOMA (q-pair).....	645, 766
HEMATURIA (data).....	926
HEMATURIA (q-pair).....	713
HEME (prod).....	1212
HEME (q-pair).....	430
HEMODIALYSIS (meth).....	1187
HEMODIALYSIS (prod).....	1212
HEMODIALYSIS (q-pair).....	311, 429, 430
HEMORRHAGE (data).....	927
HEMORRHAGE (q-pair).....	414
HEMOVAC (data).....	927
HEMOVAC (meth).....	1187
HEMOVAC (q-pair).....	431, 755
HEPATIC FORMULA (data).....	927
HEPATIC FORMULA (prod).....	1212
HEPATIC FORMULA (q-pair).....	309
HEPLOCK (data).....	927
HEPLOCK (q-pair).....	343, 493
HIGH (data).....	927
HIGH (q-pair).....	53, 156, 214
HIGH CALORIE (data).....	928
HIGH CALORIE (q-pair).....	150
HIGH FALL RISK (data).....	928
HIGH FALL RISK (q-pair).....	40
HIGH FIBER (data).....	928
HIGH FIBER (q-pair).....	150
HIGH FOWLER'S (data).....	928
HIGH FOWLER'S (q-pair).....	13
HIGH PROTEIN (data).....	928
HIGH PROTEIN (q-pair).....	150
HIGH PROTEIN FORMULA (prod).....	1212
HIGH PROTEIN FORMULA (q-pair).....	172
HIP (loca).....	1135

HIP (q-pair).....	116
HIP-PROTECTORS (data).....	929
HIP-PROTECTORS (q-pair).....	190
HIRSUTISM (data).....	929
HIRSUTISM (q-pair).....	221
HM: HAND MOTIONS (data).....	929
HM: HAND MOTIONS (q-pair).....	152
HOMANS SIGN.....	228
HOME (data).....	929
HOME (q-pair).....	467
HOME HEALTH (data).....	929
HOME HEALTH (q-pair).....	164
HOSTILE (data).....	930
HOSTILE (q-pair).....	56
HOT (data).....	930
HOT (q-pair).....	658
HOUR (unit).....	1248
HUMIDIFIER (data).....	930
HUMIDIFIER (q-pair).....	744
HYDRATION - DIET (task).....	1235
HYDRATION - HEPARIN DRIP (task).....	1235
HYDRATION - HYPERAL (task).....	1235
HYDRATION - IV FLUIDS (task).....	1236
HYDRATION - NPO (task).....	1236
HYDRATION - PERIPHERAL IV (task).....	1236
HYDRATION - PICC/CENTRAL LINE (task).....	1236
HYDRATION - SUPPLEMENT SNACK (task).....	1236
HYDRATION - TUBE FEEDING (task).....	1237
HYDROCOLLOIDAL DRESSING (data).....	930
HYDROCOLLOIDAL DRESSING (q-pair).....	548, 758, 765
HYDROFIBER (data).....	930
HYDROFIBER (q-pair).....	548
HYDROFIBER WITH SILVER (data).....	931
HYDROFIBER WITH SILVER (q-pair).....	548
HYDROGEL DRESSING (data).....	931
HYDROGEL DRESSING (q-pair).....	548, 758
HYDROMORPHONE (data).....	931
HYDROMORPHONE (q-pair).....	181
HYPERACTIVE (data).....	931
HYPERACTIVE (q-pair).....	1
HYPERDYNAMIC (data).....	931
HYPERDYNAMIC (q-pair).....	48, 103, 348, 442,
	450
HYPERGRANULIZATION (data).....	932
HYPERGRANULIZATION (q-pair).....	760
HYPERTONIC HYDROGEL (data).....	932
HYPERTONIC HYDROGEL (q-pair).....	758
HYPOACTIVE (data).....	932
HYPOACTIVE (q-pair).....	1
HYPOALBUMINEMIA (data).....	932
HYPOALBUMINEMIA (q-pair).....	545
HYPOTENSION (data).....	932
HYPOTENSION (q-pair).....	183, 398, 498
HYPOXIA (data).....	933
HYPOXIA (q-pair).....	398

- I -

I (data).....	933
I (q-pair).....	85
I:E RATIO.....	229
IABP (data).....	933

IABP (meth).....	1187
IABP (q-pair).....	24, 378
IABP ASSISTED SYSTOLIC PRESSURE.....	230
IABP BAEDP.....	231
IABP BALLOON SIZE.....	232
IABP CATHETER FLUSH.....	233
IABP CIRCULATION ASSESSED.....	234
IABP DRESSING TYPE.....	235
IABP PEAK DIASTOLIC AUGMENTED PRESSURE.....	236
IABP RATIO.....	237
IABP SITE CONDITION.....	238
IABP SIZE.....	239
IABP TIMING.....	240
IABP TIMING CHECK.....	241
IABP UAEDP.....	242
IABP UNASSISTED SYSTOLIC PRESSURE.....	243
ICE PACK (data).....	933
ICE PACK (q-pair).....	744
ICTERIC (data).....	933
ICTERIC (q-pair).....	645
ICU (data).....	934
ICU (q-pair).....	467, 553
ID BAND (data).....	934
ID BAND (q-pair).....	695
ID CARD (data).....	934
ID CARD (q-pair).....	695
IDENTIFICATION (data).....	934
IDENTIFICATION (q-pair).....	40
IDIOVENTRICULAR (data).....	934
IDIOVENTRICULAR (q-pair).....	84
II (data).....	935
II (q-pair).....	85
III (data).....	935
III (q-pair).....	85
ILEAL CONDUIT (data).....	935
ILEAL CONDUIT (q-pair).....	219, 716
ILEAL-CONDUIT (data).....	935
ILEAL-CONDUIT (q-pair).....	424
ILEOSTOMY (data).....	935
ILEOSTOMY (q-pair).....	424
IMAGING (data).....	936
IMAGING (q-pair).....	467
IMPAIRED LISTENING SKILLS (data).....	936
IMPAIRED LISTENING SKILLS (q-pair).....	161
IMPLANTS (data).....	936
IMPLANTS (q-pair).....	396
INACTIVE - ARTERIAL LINE SITE.....	244
INACTIVE - ARTERIAL LINE TYPE.....	245
INACTIVE - AUDIOMETRY.....	246
INACTIVE - BISPECTRAL INDEX.....	247
INACTIVE - BODY MASS INDEX.....	248
INACTIVE - CARDIAC AV FISTULA BRUIT.....	249
INACTIVE - CATHETER TYPE.....	250
INACTIVE - CENTRAL LINE SITE.....	251
INACTIVE - CHEST TUBE DRESSING INTACT.....	252
INACTIVE - CHEST TUBE SITE.....	253
INACTIVE - COUGH AND DEEP BREATH EVENTS.....	254
INACTIVE - DATE.....	255
INACTIVE - DATE DAY.....	256
INACTIVE - DATE MONTH.....	257
INACTIVE - DATE YEAR.....	258
INACTIVE - DIALYSIS CATHETER SITE.....	259

INACTIVE - EDEMA LOCATION.....	260
INACTIVE - ENTERAL NUTRITION FORMULA.....	261
INACTIVE - GI EMESIS APPEARANCE.....	262
INACTIVE - HEARING.....	263
INACTIVE - IABP INSERTION SITE.....	264
INACTIVE - IABP SITE.....	265
INACTIVE - INCISION LOCATION.....	266
INACTIVE - INCISION SITE.....	267
INACTIVE - IV DOSE RATE.....	268
INACTIVE - IV FLUID.....	269
INACTIVE - IV INSERTION DATE.....	270
INACTIVE - IV SITE.....	271
INACTIVE - OR OUTPUT - CHEST TUBE.....	272
INACTIVE - OR OUTPUT - DRAIN.....	273
INACTIVE - OR OUTPUT - PUMP LOSS.....	274
INACTIVE - OR OUTPUT - URINE.....	275
INACTIVE - OSTOMY TYPE.....	276
INACTIVE - OXYGEN DELIVERED.....	277
INACTIVE - OXYGEN DELIVERY METHOD.....	278
INACTIVE - OXYGEN EXTRACTION RATIO.....	279
INACTIVE - PA CATHETER INSERTION DATE.....	280
INACTIVE - PA CATHETER SITE.....	281
INACTIVE - PACEMAKER DATE INSERTED.....	282
INACTIVE - PACEMAKER SITES.....	283
INACTIVE - PAIN LOCATION.....	284
INACTIVE - POST OPERATION DAY NUMBER.....	285
INACTIVE - PRESSURE ULCER SITE.....	286
INACTIVE - PRESSURE ULCER VISIBLE STRUCTURE.....	287
INACTIVE - PULSE OXIMETRY.....	288
INACTIVE - SKIN DISCHARGE COLOR.....	289
INACTIVE - SKIN DISCHARGE ODOR.....	290
INACTIVE - SKIN DRAINAGE CHARACTERISTICS.....	291
INACTIVE - VENTILATOR IN-LINE NEBULIZER.....	292
INACTIVE - WOUND SITE.....	293
INACTIVE - WOUND TOPICAL TREATMENT TYPE.....	294
INAPPROPRIATE (data).....	936
INAPPROPRIATE (q-pair).....	32, 54, 135
INCENTIVE SPIROMETRY.....	295
INCH (unit).....	1249
INCISION DISCHARGE CHARACTERISTICS.....	296
INCISION DISCHARGE CHARACTERISTICS BLOODY (data).....	936
INCISION DISCHARGE CHARACTERISTICS CLOUDY (data).....	937
INCISION DISCHARGE CHARACTERISTICS NONE (data).....	937
INCISION DISCHARGE CHARACTERISTICS PURULENT (data).....	937
INCISION DISCHARGE CHARACTERISTICS SEROSANGUINOUS (data) .....	937
INCISION DISCHARGE CHARACTERISTICS SEROUS (data).....	937
INCISION DISCHARGE CHARACTERISTICS VISCOUS (data).....	938
INCISION DRAINAGE AMOUNT.....	297
INCISION DRAINAGE AMOUNT HEAVY (data).....	938
INCISION DRAINAGE AMOUNT HEAVY (q-pair).....	297
INCISION DRAINAGE AMOUNT MODERATE (data).....	938
INCISION DRAINAGE AMOUNT MODERATE (q-pair).....	297
INCISION DRAINAGE AMOUNT NONE (data).....	938
INCISION DRAINAGE AMOUNT NONE (q-pair).....	297
INCISION DRAINAGE AMOUNT SCANT (data).....	938
INCISION DRAINAGE AMOUNT SCANT (q-pair).....	297
INCISION DRAINAGE AMOUNT SMALL (data).....	939
INCISION DRAINAGE AMOUNT SMALL (q-pair).....	297
INCISION DRESSING PREVENTS ASSESSMENT.....	298
INCISION DRESSING TYPE.....	299
INCISION DRESSING TYPE ANTIMICROBIAL (data).....	939

INCISION DRESSING TYPE ANTIMICROBIAL (q-pair).....	299
INCISION DRESSING TYPE ANTIMICROBIAL DISK (data).....	939
INCISION DRESSING TYPE ANTIMICROBIAL DISK (q-pair).....	299
INCISION DRESSING TYPE BALSAM OF PERU (data).....	939
INCISION DRESSING TYPE BALSAM OF PERU (q-pair).....	299
INCISION DRESSING TYPE BISMUTH GAUZE (data).....	939
INCISION DRESSING TYPE BISMUTH GAUZE (q-pair).....	299
INCISION DRESSING TYPE CALCIUM ALGINATE (data).....	940
INCISION DRESSING TYPE CALCIUM ALGINATE (q-pair).....	299
INCISION DRESSING TYPE FOAM DRESSING (data).....	940
INCISION DRESSING TYPE FOAM DRESSING (q-pair).....	299
INCISION DRESSING TYPE GAUZE (data).....	940
INCISION DRESSING TYPE GAUZE (q-pair).....	299
INCISION DRESSING TYPE GAUZE PACKING (data).....	940
INCISION DRESSING TYPE GAUZE PACKING (q-pair).....	299
INCISION DRESSING TYPE GAUZE WRAP (data).....	940
INCISION DRESSING TYPE GAUZE WRAP (q-pair).....	299
INCISION DRESSING TYPE HYDROCOLLOIDAL (data).....	941
INCISION DRESSING TYPE HYDROCOLLOIDAL (q-pair).....	299
INCISION DRESSING TYPE HYDROGEL (data).....	941
INCISION DRESSING TYPE HYDROGEL (q-pair).....	299
INCISION DRESSING TYPE HYPERTONIC HYDROGEL (data).....	941
INCISION DRESSING TYPE HYPERTONIC HYDROGEL (q-pair).....	299
INCISION DRESSING TYPE NONE (data).....	941
INCISION DRESSING TYPE NONE (q-pair).....	299
INCISION DRESSING TYPE PRESSURE (data).....	941
INCISION DRESSING TYPE PRESSURE (q-pair).....	299
INCISION DRESSING TYPE TRANSPARENT (data).....	942
INCISION DRESSING TYPE TRANSPARENT (q-pair).....	299
INCISION EDGE.....	300
INCISION EDGE INDISTINCT/ATTACHED (data).....	942
INCISION EDGE INDISTINCT/ATTACHED (q-pair).....	300
INCISION EDGE NON-ATTACHED (data).....	942
INCISION EDGE NON-ATTACHED (q-pair).....	300
INCISION EDGE ROLLED (data).....	942
INCISION EDGE ROLLED (q-pair).....	300
INCISION EDGE TUNNELING (data).....	942
INCISION EDGE TUNNELING (q-pair).....	300
INCISION EDGE UNDERMINING (data).....	943
INCISION EDGE UNDERMINING (q-pair).....	300
INCISION ERYTHEMA.....	301
INCISION ODOR.....	302
INCISION ODOR FOUL (data).....	943
INCISION ODOR FOUL (q-pair).....	302
INCISION ODOR NONE (data).....	943
INCISION ODOR NONE (q-pair).....	302
INCISION PRESENT.....	303
INCISION STAPLES INTACT.....	304
INCISION SUTURES INTACT.....	305
INCISION TENDER TO PALPATION.....	306
INCISION WOUND CLOSURE APPROXIMATED (data).....	943
INCISION WOUND CLOSURE APPROXIMATED (q-pair).....	650
INCISION WOUND CLOSURE DEHISCENCE (data).....	943
INCISION WOUND CLOSURE DEHISCENCE (q-pair).....	650
INCISION WOUND CLOSURE GLUE (data).....	944
INCISION WOUND CLOSURE GLUE (q-pair).....	650
INCISION WOUND CLOSURE NOT CLOSED (data).....	944
INCISION WOUND CLOSURE NOT CLOSED (q-pair).....	650
INCISIONAL SITE (data).....	944
INCISIONAL SITE (loca).....	1135
INCISIONAL SITE (q-pair).....	160, 473
INCOMPREHENSIBLE (data).....	944

INCOMPREHENSIBLE (q-pair).....	135
INCONTENENCE DERMATITIS (data).....	944
INCONTENENCE DERMATITIS (q-pair).....	766
INCONTINENT (data).....	945
INCONTINENT (q-pair).....	220, 716
INCREASE ROUNDS (data).....	945
INCREASE ROUNDS (q-pair).....	190
INDEPENDENT (data).....	945
INDEPENDENT (q-pair).....	6, 8, 9, 10
INDICATES UNDERSTANDING (data).....	945
INDICATES UNDERSTANDING (q-pair).....	591, 596
INDIFFERENT (data).....	945
INDIFFERENT (q-pair).....	163
INDURATION (data).....	946
INDURATION (q-pair).....	640
INDWELLING CATHETER (data).....	946
INDWELLING CATHETER (q-pair).....	219, 716
INEFFECTIVE COPING (data).....	946
INEFFECTIVE COPING (q-pair).....	56
INFECTION (data).....	946
INFECTION (q-pair).....	164
INFECTION PREVENTION AND SCREENING (data).....	946
INFECTION PREVENTION AND SCREENING (q-pair).....	164
INFLATABLE HEEL LIFT (data).....	947
INFLATABLE HEEL LIFT (q-pair).....	544
INFUSION SET CHANGED.....	307
INITIATED (data).....	947
INITIATED (q-pair).....	622
INSERTED THROUGH EXISTING INTRODUCER (data).....	947
INSERTED THROUGH EXISTING INTRODUCER (q-pair).....	94, 444
INSERTED THROUGH NEW INTRODUCER (data).....	947
INSERTED THROUGH NEW INTRODUCER (q-pair).....	94, 444
INSPIRATORY WHEEZES (data).....	947
INSPIRATORY WHEEZES (q-pair).....	74
INSULIN PUMP SQ INFUSION ON?.....	308
INTACT (data).....	948
INTACT (q-pair).....	43, 95, 618, 619, 620, 621, 646
INTAKE.....	309
INTAKE - CBI IRRIGATION TO SUBTRACT (NEGATIVE NUMBER)..	310
INTAKE - DIALYSIS.....	311
INTAKE - EF/SUPPLEMENT CALORIES.....	312
INTAKE - ENTERAL FEEDING TUBE.....	313
INTAKE - ENTERAL TUBE FLUSH.....	314
INTAKE - GI TUBE.....	315
INTAKE - IRRIGATION.....	316
INTAKE - IV BLOOD PRODUCTS.....	317
INTAKE - IV CALORIES.....	318
INTAKE - IV CARDIAC OUTPUT INJECTATE.....	319
INTAKE - IV COLLOID.....	320
INTAKE - IV CRYSTALLOIDS.....	321
INTAKE - IV FLUID.....	322
INTAKE - IV MEDS.....	323
INTAKE - IV NUTRITION.....	324
INTAKE - LINE FLUSH.....	325
INTAKE - PCA VOLUME.....	326
INTAKE - PERCENT MEAL CONSUMED.....	327
INTAKE - PO CALORIES.....	328
INTAKE - PO FLUIDS.....	329
INTAKE - TPN CALORIES.....	330
INTAKE - TRANSFUSION .....	331
INTENTIONAL (data).....	948

INTENTIONAL (q-pair).....	412
INTERMITTENT (data).....	948
INTERMITTENT (q-pair).....	156, 214, 216, 470, 471, 767
INTERMITTENT SLEEP (data).....	948
INTERMITTENT SLEEP (q-pair).....	664
INTERRUPTED (data).....	948
INTERRUPTED (q-pair).....	622
INTERVENTIONAL RADIOLOGY (data).....	949
INTERVENTIONAL RADIOLOGY (q-pair).....	553
INTRACRANIAL PRESSURE - DO NOT CHECK TOTAL OR SUBTOTALS	332
INTRODUCER (data).....	949
INTRODUCER (q-pair).....	102
INTRODUCER CATHETER (meth).....	1187
INVASIVE (meth).....	1187
INVASIVE (q-pair).....	63, 146, 147, 148, 149, 377, 378, 570, 690
IRREGULAR (data).....	949
IRREGULAR (q-pair).....	577
IRRIGATION (meth).....	1188
IRRITABLE (data).....	949
IRRITABLE (q-pair).....	474
ISOLATION PRECAUTIONS.....	333
ISOSOURCE VHN (data).....	949
ISOSOURCE VHN (q-pair).....	172
ITCHING (data).....	950
ITCHING (q-pair).....	183, 498
IV (meth).....	1188
IV (q-pair).....	309
IV DISCONTINUED (data).....	950
IV DISCONTINUED (q-pair).....	151
IV DOSE RATE - UNITS/HOUR - DO NOT CHECK TOTAL OR SUBTOTALS	.....334
IV DOSE RATE - UNITS/kg/HOUR - DO NOT CHECK TOTAL OR SUBTOTALS	.....335
IV DOSE RATE - g/TIME - DO NOT CHECK TOTAL OR SUBTOTALS	336
IV DOSE RATE - mcg/TIME - DO NOT CHECK TOTAL OR SUBTOTALS	.....337
IV DOSE RATE - mcg/kg/HOUR - DO NOT CHECK TOTAL OR SUBTOTALS	.....338
IV DOSE RATE - mg/TIME- DO NOT CHECK TOTAL OR SUBTOTALS	339
IV DOSE RATE - mg/kg/HOUR - DO NOT CHECK TOTAL OR SUBTOTALS	.....340
IV FLOW RATE - ml/TIME - DO NOT CHECK TOTAL OR SUBTOTALS	.....341
IV FLOW RATE - ml/kg/TIME - DO NOT CHECK TOTAL OR SUBTOTALS	.....342
IV LINE CARE.....	343
IV PORT / NUMBER.....	344
IV PUMP (data).....	950
IV PUMP (q-pair).....	24
IV SEDATION (data).....	950
IV SEDATION (q-pair).....	31
IV SITE CONDITION.....	345
IVPB (data).....	950

- J -

JAUNDICE (data).....	951
JAUNDICE (q-pair).....	645
JEJUNOSTOMY (data).....	951



JEJUNOSTOMY (q-pair).....	166, 215, 424
JEVITY 1.0 (data).....	951
JEVITY 1.0 (q-pair).....	172
JEVITY 1.2 (data).....	951
JEVITY 1.2 (q-pair).....	172
JEVITY 1.5 (data).....	951
JEVITY 1.5 (q-pair).....	172
JOINT CAPSULE (data).....	952
JP DRAIN (data).....	952
JP DRAIN (meth).....	1188
JP DRAIN (prod).....	1213
JP DRAIN (q-pair).....	431, 755
JUICE (prod).....	1213
JUICE (q-pair).....	309, 329
JUNCTIONAL (data).....	952
JUNCTIONAL (q-pair).....	84
JUNCTIONAL ESCAPE BEAT (data).....	952
JUNCTIONAL ESCAPE BEAT (q-pair).....	80
JUNCTIONAL TACHYCARDIA (data).....	952
JUNCTIONAL TACHYCARDIA (q-pair).....	84
JUVEN (data).....	953
JUVEN (q-pair).....	172

- K -

K-PAD (data).....	953
K-PAD (q-pair).....	744
KG/M2 (unit).....	1249
KGS/HOUR (unit).....	1249
KILOGRAM (unit).....	1249
KNEE (loca).....	1136
KNEE (q-pair).....	116
KNEES (loca).....	1136
KNEES (q-pair).....	160
KNIFE-LIKE (data).....	953
KNIFE-LIKE (q-pair).....	470
KNOWN/SUSPECTED OBSTRUCTIVE SLEEP APNEA (data).....	953
KNOWN/SUSPECTED OBSTRUCTIVE SLEEP APNEA (q-pair).....	396
KUSSMAUL (data).....	953
KUSSMAUL (q-pair).....	577

- L -

L ABDOMEN (loca).....	1136
L ABDOMEN (q-pair).....	761
L ANKLE (loca).....	1136
L ANKLE (q-pair).....	303, 473, 547, 600, 658, 761
L ANTECUBITAL (loca).....	1136
L ANTECUBITAL (q-pair).....	102, 143, 219, 344, 345, 547, 761
L ARM (loca).....	1137
L ARM (q-pair).....	63, 102, 143, 149, 160, 165, 219, 344, 345, 473, 598, 599, 658, 690
L AXILLARY (loca).....	1137
L AXILLARY (q-pair).....	46, 47, 235, 239
L BRACHIAL (loca).....	1137
L BRACHIAL (q-pair).....	46, 47, 105, 456, 461
L BUNDLE BRANCH BLOCK (data).....	954

L BUNDLE BRANCH BLOCK (q-pair).....	84
L CALF (loca).....	1137
L CALF (q-pair).....	303
L CAROTID (loca).....	1137
L CAROTID (q-pair).....	79, 86
L CHEST (loca).....	1138
L CHEST (q-pair).....	108, 473, 658, 761
L DORSALIS PEDIS (loca).....	1138
L EAR (loca).....	1138
L EAR (q-pair).....	303, 473, 547, 658, 761
L ELBOW (loca).....	1138
L ELBOW (q-pair).....	473, 547, 658, 761
L EXTERNAL JUGULAR (loca).....	1138
L EXTERNAL JUGULAR (q-pair).....	102, 143, 219, 344, 345, 449
L EYE (loca).....	1139
L EYE (q-pair).....	133, 160, 572, 573, 741, 761
L FACE (loca).....	1139
L FACE (q-pair).....	160, 303, 761
L FEMORAL (loca).....	1139
L FEMORAL (q-pair).....	46, 47, 79, 86, 102, 105, 143, 219, 235, 239, 303, 344, 345, 449, 456, 461, 547, 571, 761
L FINGER (loca).....	1139
L FINGER (q-pair).....	303
L FLANK (loca).....	1139
L FLANK (q-pair).....	303
L FOOT (loca).....	1140
L FOOT (q-pair).....	160, 196, 303, 344, 345, 473, 547, 618, 658, 761
L FOOT DOWN (data).....	954
L FOOT DOWN (q-pair).....	634
L FOOT UP (data).....	954
L FOOT UP (q-pair).....	634
L GLUTEUS (loca).....	1140
L GLUTEUS (q-pair).....	303, 547, 761
L GROIN (loca).....	1140
L GROIN (q-pair).....	473, 658
L HAND (loca).....	1140
L HAND (q-pair).....	160, 222, 473, 619, 658, 761
L HEAD (loca).....	1140
L HEAD (q-pair).....	547, 761
L HEAD DOWN (data).....	954
L HEAD DOWN (q-pair).....	634
L HEAD UP (data).....	954
L HEAD UP (q-pair).....	634
L HEEL (loca).....	1141
L HEEL (q-pair).....	547, 761
L HIP (loca).....	1141
L HIP (q-pair).....	303, 473, 658, 761
L INGUINAL (loca).....	1141
L INGUINAL (q-pair).....	303
L INNER KNEE (loca).....	1141
L INNER KNEE (q-pair).....	547, 761
L INTERNAL JUGULAR (loca).....	1141
L INTERNAL JUGULAR (q-pair).....	102, 105, 143, 219,

	344, 345, 449, 456, 461
L ISCHIAL (loca).....	1142
L ISCHIAL (q-pair).....	547, 761
L JAW (loca).....	1142
L JAW (q-pair).....	473, 658
L KNEE (loca).....	1142
L KNEE (q-pair).....	160, 303, 473, 658
L LATERAL HEEL (loca).....	1142
L LATERAL HEEL (q-pair).....	547
L LATERAL MALLEOLUS (loca).....	1142
L LATERAL MALLEOLUS (q-pair).....	547
L LEG (loca).....	1143
L LEG (q-pair).....	63, 136, 149, 160, 165, 344, 345, 473, 598, 600, 658, 690
L LOWER ANTERIOR CHEST (loca).....	1143
L LOWER ANTERIOR CHEST (q-pair).....	110, 303
L LOWER ARM (loca).....	1143
L LOWER ARM (q-pair).....	303, 473, 547, 658, 761
L LOWER EXTREMITY (loca).....	1143
L LOWER EXTREMITY (q-pair).....	411, 620
L LOWER LEG (loca).....	1143
L LOWER LEG (q-pair).....	303, 761
L LOWER LIMB (loca).....	1144
L LOWER LIMB (q-pair).....	352
L LOWER POSTERIOR CHEST (loca).....	1144
L LOWER POSTERIOR CHEST (q-pair).....	110, 303
L MANDIBULAR (loca).....	1144
L MANDIBULAR (q-pair).....	303
L MEDIAL HEEL (loca).....	1144
L MEDIAL HEEL (q-pair).....	547
L MEDIAL MALLEOLUS (loca).....	1144
L MEDIAL MALLEOLUS (q-pair).....	547
L MIDDLE ANTERIOR CHEST (loca).....	1145
L MIDDLE ANTERIOR CHEST (q-pair).....	110, 303
L MIDDLE POSTERIOR CHEST (loca).....	1145
L MIDDLE POSTERIOR CHEST (posi).....	1202
L MIDDLE POSTERIOR CHEST (q-pair).....	110, 303
L NECK (data).....	955
L NECK (loca).....	1145
L NECK (q-pair).....	303
L NEPHROSTOMY (data).....	955
L NEPHROSTOMY (q-pair).....	219, 716
L NOSE (loca).....	1145
L NOSE (q-pair).....	473, 658
L OCCIPITAL (loca).....	1145
L OCCIPITAL (q-pair).....	303
L ORAL (loca).....	1146
L ORAL (q-pair).....	303
L ORBIT (loca).....	1146
L ORBIT (q-pair).....	303
L OUTER KNEE (loca).....	1146
L OUTER KNEE (q-pair).....	547, 761
L PARIETAL (loca).....	1146
L PARIETAL (q-pair).....	303
L PEDAL (loca).....	1146
L PEDAL (q-pair).....	571
L POST TIBIAL (loca).....	1147
L POST TIBIAL (q-pair).....	571
L RADIAL (loca).....	1147

L RADIAL (q-pair).....	46, 47, 571
L RIBS (loca).....	1147
L RIBS (q-pair).....	547, 761
L SCAPULA (loca).....	1147
L SCAPULA (q-pair).....	473, 547, 658, 761
L SHOULDER (loca).....	1147
L SHOULDER (q-pair).....	473, 547, 658, 761
L STENT (data).....	955
L STENT (q-pair).....	219
L SUBCLAVIAN (loca).....	1148
L SUBCLAVIAN (q-pair).....	102, 105, 143, 219, 303, 344, 345, 449, 456, 461
L TEMPORAL (loca).....	1148
L TEMPORAL (q-pair).....	79, 86
L TESTICLE (loca).....	1148
L TESTICLE (q-pair).....	761
L THIGH (loca).....	1148
L THIGH (q-pair).....	303, 473, 658
L THORACIC ANTERIOR (loca).....	1148
L THORACIC ANTERIOR (q-pair).....	303
L THORACIC LATERAL (loca).....	1149
L THORACIC LATERAL (q-pair).....	303
L THORACIC POSTERIOR (loca).....	1149
L THORACIC POSTERIOR (q-pair).....	303
L THORACOTOMY (loca).....	1149
L THORACOTOMY (q-pair).....	303
L THUMB (loca).....	1149
L THUMB (q-pair).....	303
L TOE (loca).....	1149
L TOE (q-pair).....	303
L TRAPEZIUS (loca).....	1150
L TRAPEZIUS (q-pair).....	303
L TROCHANTER (loca).....	1150
L TROCHANTER (q-pair).....	547, 761
L UPPER ANTERIOR CHEST (loca).....	1150
L UPPER ANTERIOR CHEST (q-pair).....	110, 303
L UPPER ARM (loca).....	1150
L UPPER ARM (q-pair).....	102, 143, 219, 303, 344, 345, 473, 547, 658, 761
L UPPER EXTREMITY (loca).....	1150
L UPPER EXTREMITY (q-pair).....	411, 621
L UPPER LEG (loca).....	1151
L UPPER LEG (q-pair).....	761
L UPPER LIMB (loca).....	1151
L UPPER LIMB (q-pair).....	712
L UPPER POSTERIOR CHEST (loca).....	1151
L UPPER POSTERIOR CHEST (q-pair).....	110, 303
L WRIST (loca).....	1151
L WRIST (q-pair).....	102, 143, 219, 344, 345, 598, 599, 761
L1 (data).....	955
L1 (q-pair).....	179, 180
L2 (data).....	955
L2 (q-pair).....	179, 180
L3 (data).....	956
L3 (q-pair).....	179, 180
L4 (data).....	956
L4 (q-pair).....	179, 180
L5 (data).....	956
L5 (q-pair).....	179, 180

LABS PER ORDERS.....	346
LACERATION (data).....	956
LACERATION (q-pair).....	766
LACKS INTEREST (data).....	956
LACKS INTEREST (q-pair).....	163
LACTOSE RESTRICTED (data).....	957
LACTOSE RESTRICTED (q-pair).....	150
LANGUAGE (data).....	957
LANGUAGE (q-pair).....	161
LARGE (data).....	957
LARGE (q-pair).....	428, 435, 437, 562
LARYNGEAL MASK (meth).....	1188
LARYNGEAL MASK (q-pair).....	686
LARYNGECTOMY (data).....	957
LARYNGECTOMY (q-pair).....	22
LEAD (meth).....	1188
LEAD (q-pair).....	570
LEFT (data).....	957
LEFT (loca).....	1151
LEFT (q-pair).....	13
LEFT EYE (loca).....	1152
LEFT EYE (q-pair).....	152
LEFT LOWER LOBE (loca).....	1152
LEFT LOWER LOBE (q-pair).....	74
LEFT NECK (loca).....	1152
LEFT NECK (q-pair).....	160
LEFT UPPER LOBE (loca).....	1152
LEFT UPPER LOBE (q-pair).....	74
LEFT VENTRICULAR EJECTION FRACTION - CATEGORICAL.....	348
LEFT VENTRICULAR EJECTION FRACTION - CONTINUOUS.....	347
LEFT VENTRICULAR STROKE WORK.....	349
LEG MOVEMENTS.....	350
LESS THAN 30 (data).....	958
LESS THAN 30 (q-pair).....	348
LETHARGIC (data).....	958
LETHARGIC (q-pair).....	56, 351
LEVEL OF CONSCIOUSNESS.....	351
LEVELED (data).....	958
LEVELED (q-pair).....	343
LG ADULT (meth).....	1189
LG ADULT CUFF (meth).....	1189
LIBERAL DIABETIC (data).....	958
LIBERAL DIABETIC (q-pair).....	150
LIFT SCALE (meth).....	1189
LIFT SCALE (q-pair).....	746
LIMITED ASSISTANCE (data).....	958
LIMITED ASSISTANCE (q-pair).....	8, 9, 10
LIMITED ATTENTION SPAN (data).....	959
LIMITED ATTENTION SPAN (q-pair).....	161
LIMITED TO INCISION (data).....	959
LIMITED TO INCISION (q-pair).....	301, 648
LINE FLUSH (data).....	959
LINE FLUSH (prod).....	1213
LINE FLUSH HEPARIN (prod).....	1213
LINE FLUSH HEPARIN (q-pair).....	325
LIP SALVE (data).....	959
LIP SALVE (q-pair).....	15
LIPIDS (prod).....	1213
LIPIDS (q-pair).....	324
LISTED DIET (data).....	959
LISTED DIET (q-pair).....	150
LITER (unit).....	1249

LITERS PER MINUTE (unit).....	1250
LITERS PER MINUTE PER METER2 (BSA) (unit).....	1250
LOCAL (data).....	960
LOCAL (q-pair).....	31
LOCALIZES TO PAIN (data).....	960
LOCALIZES TO PAIN (q-pair).....	134
LOCK BOX (data).....	960
LOCK BOX (q-pair).....	516
LOCKER (data).....	960
LOCKER (q-pair).....	515, 516
LOOSE (data).....	960
LOOSE (q-pair).....	680
LOW (data).....	961
LOW (q-pair).....	53, 156, 214
LOW AIRLOSS OVERLAY (data).....	961
LOW AIRLOSS OVERLAY (q-pair).....	665
LOW AIRLOSS/ALTERNATING PRESSURE COMBO (data).....	961
LOW AIRLOSS/ALTERNATING PRESSURE COMBO (q-pair).....	665
LOW CHOLESTEROL (data).....	961
LOW CHOLESTEROL (q-pair).....	150
LOW FIBER (data).....	961
LOW FIBER (q-pair).....	150
LOW PROTEIN 60 GM (data).....	962
LOW PROTEIN 60 GM (q-pair).....	150
LOW TYRAMINE (data).....	962
LOW TYRAMINE (q-pair).....	150
LOWER ABDOMEN (loca).....	1152
LOWER ABDOMEN (q-pair).....	761
LOWER ARM (loca).....	1153
LOWER ARM (q-pair).....	116
LOWER BACK (loca).....	1153
LOWER BACK (q-pair).....	761
LOWER LIMB EXTREMITIES STRENGTH.....	352
LP: LIGHT PERCEPTION (data).....	962
LP: LIGHT PERCEPTION (q-pair).....	152
LPP: LIGHT PERCEPTION WITH PROJECTION (data).....	962
LPP: LIGHT PERCEPTION WITH PROJECTION (q-pair).....	152
LUMBAR DRAIN (data).....	962
LUMBAR DRAIN (q-pair).....	755
LUMBAR SPINAL PRECAUTIONS (data).....	963
LUMBAR SPINAL PRECAUTIONS (q-pair).....	12
LVN (data).....	963
LVN (q-pair).....	608
LVN/LPN (data).....	963
LVN/LPN (q-pair).....	554
LVSWI.....	353
LYING (posi).....	1202
LYING (q-pair).....	63, 149, 570, 578, 690, 746
LYING OR SITTING (data).....	963
LYING OR SITTING (q-pair).....	128

- M -

MACERATION (data).....	963
MACERATION (q-pair).....	640
MAGNACAL (data).....	964
MAGNACAL (prod).....	1214
MAGNACAL (q-pair).....	172
MAGNETIC RESONANCE IMAGING (meth).....	1189
MAGNETIC RESONANCE IMAGING (q-pair).....	347, 348
MAINTAINING EFFECTIVE VENTILATION (data).....	964

MAINTAINING EFFECTIVE VENTILATION (q-pair).....	465
MALLAMPATI SCORE.....	354
MAROON (data).....	964
MAROON (q-pair).....	681, 753
MASK (meth).....	1189
MASK (q-pair).....	686, 687
MD (data).....	964
MD (q-pair).....	608
MDS: BATHING-SELF PERFORMANCE.....	355
MDS: BATHING-SUPPORT.....	356
MDS: BED MOBILITY-SELF PERFORMANCE.....	357
MDS: BED MOBILITY-SUPPORT.....	358
MDS: DRESSING-SELF PERFORMANCE.....	359
MDS: DRESSING-SUPPORT.....	360
MDS: EATING-SELF PERFORMANCE.....	361
MDS: EATING-SUPPORT.....	362
MDS: LOCOMOTION OFF UNIT-SELF PERFORMANCE.....	363
MDS: LOCOMOTION OFF UNIT-SUPPORT.....	364
MDS: LOCOMOTION ON UNIT-SELF PERFORMANCE.....	365
MDS: LOCOMOTION ON UNIT-SUPPORT.....	366
MDS: PERSONAL HYGIENE-SELF PERFORMANCE.....	367
MDS: PERSONAL HYGIENE-SUPPORT.....	368
MDS: TOILET USE-SELF PERFORMANCE.....	369
MDS: TOILET USE-SUPPORT.....	370
MDS: TRANSFER-SELF PERFORMANCE.....	371
MDS: TRANSFER-SUPPORT.....	372
MDS: WALK IN CORRIDOR-SELF PERFORMANCE.....	373
MDS: WALK IN CORRIDOR-SUPPORT.....	374
MDS: WALK IN ROOM-SELF PERFORMANCE.....	375
MDS: WALK IN ROOM-SUPPORT.....	376
MEAL (data).....	964
MEAL (q-pair).....	5
MEAN ARTERIAL PRESSURE.....	377
MEAN BLOOD PRESSURE.....	378
MEAN BLOOD PRESSURE ARTERIAL LINE.....	379
MEAN BLOOD PRESSURE NON-INVASIVE.....	380
MEAN BP CUFF.....	381
MECHANICAL (data).....	965
MECHANICAL (q-pair).....	150
MED LAST DOSE: ANTIBIOTIC.....	382
MED LAST DOSE: ANTICOAGULANT.....	383
MED LAST DOSE: ASPIRIN.....	384
MED LAST DOSE: ESA ERYTHROPOETIN.....	385
MED LAST DOSE: HYPOGLYCEMIC AGENT.....	386
MED LAST DOSE: NON-STEROIDAL ANTI-INFLAMMATORY DRUG....	387
MED RECONCILIATION DOCUMENTED (data).....	965
MED RECONCILIATION DOCUMENTED (q-pair).....	151
MEDIAL PORT (data).....	965
MEDIAL PORT (q-pair).....	344
MEDIASTINAL (loca).....	1153
MEDIASTINAL (prod).....	1214
MEDIASTINAL (q-pair).....	235, 239, 430
MEDICAL DEVICE SAFETY (data).....	965
MEDICAL DEVICE SAFETY (q-pair).....	58, 130
MEDICAL NECESSITY CONFIRMED (data).....	965
MEDICAL NECESSITY CONFIRMED (q-pair).....	343
MEDICATION (data).....	966
MEDICATION (prod).....	1214
MEDICATION (q-pair).....	164, 475
MEDICATION ROOM (data).....	966
MEDICATION ROOM (q-pair).....	55, 129
MEDIUM (data).....	966

MEDIUM (q-pair).....	156
MEMORY (data).....	966
MEMORY (q-pair).....	161
MENTAL HEALTH (data).....	966
MENTAL HEALTH (q-pair).....	164
MENTAL HEALTH HOLD (data).....	967
MENTAL HEALTH HOLD (q-pair).....	164
MENTAL STATUS - AWAKE ORIENTED (data).....	967
MENTAL STATUS - AWAKE ORIENTED (q-pair).....	465
MENTAL STATUS - BASELINE FOR KNOWN DEMENTIA (data).....	967
MENTAL STATUS - BASELINE FOR KNOWN DEMENTIA (q-pair).....	465
METAL IMPLANTS (data).....	967
METAL IMPLANTS (q-pair).....	396
MEWS: CONSCIOUS LEVEL.....	388
MEWS: HEART RATE.....	389
MEWS: HOURLY URINE.....	390
MEWS: RESPIRATORY.....	391
MEWS: SCORE.....	392
MEWS: SCORE UNIT (unit).....	1250
MEWS: SYSTOLIC BLOOD PRESSURE.....	393
MEWS: TEMPERATURE.....	394
MICROGRAMS (unit).....	1250
MICROGRAMS PER HOUR (unit).....	1250
MICROGRAMS PER KILOGRAM PER HOUR (unit).....	1251
MICROGRAMS PER KILOGRAM PER MINUTE (unit).....	1251
MICROGRAMS PER MINUTE (unit).....	1251
MICROGRAMS PER SHIFT (unit).....	1251
MICU (data).....	967
MICU (q-pair).....	467
MID BACK (loca).....	1153
MID BACK (q-pair).....	761
MID-STERNAL (loca).....	1153
MID-STERNAL (q-pair).....	108
MIDDLE (loca).....	1154
MIDLINE ABDOMINAL (loca).....	1154
MIDLINE ABDOMINAL (q-pair).....	303
MIDLINE NECK (loca).....	1154
MIDLINE NECK (q-pair).....	303
MILD-MODERATELY REDUCED (data).....	968
MILD-MODERATELY REDUCED (q-pair).....	348
MILK (prod).....	1214
MILK (q-pair).....	329
MILLIAMPS (unit).....	1251
MILLIBAR (unit).....	1252
MILLIGRAMS (unit).....	1252
MILLIGRAMS PER HOUR (unit).....	1252
MILLIGRAMS PER KILOGRAM PER HOUR (data).....	968
MILLIGRAMS PER KILOGRAM PER HOUR (unit).....	1252
MILLIGRAMS PER KILOGRAM PER MINUTE (unit).....	1252
MILLIGRAMS PER MINUTE (unit).....	1253
MILLIGRAMS PER SHIFT (unit).....	1253
MILLILITER (unit).....	1253
MILLILITERS PER BEAT (unit).....	1253
MILLILITERS PER BEAT PER METER SQUARE (BSA) (unit).....	1253
MILLILITERS PER HOUR (unit).....	1254
MILLILITERS PER METER SQUARE (BSA) (unit).....	1254
MILLILITERS PER MINUTE (unit).....	1254
MILLILITERS PER MINUTE PER METER SQUARE (unit).....	1254
MILLIMETERS (unit).....	1254
MILLIMETERS OF MERCURY (unit).....	1255
MILLIVOLTS (unit).....	1255
MILLILITERS PER KILOGRAM PER HOUR (unit).....	1255



MILLILITERS PER KILOGRAM PER MINUTE (unit).....	1255
MINUTE (unit).....	1255
MISSING (data).....	968
MOANING (data).....	968
MOANING (q-pair).....	474
MODERATE (data).....	968
MODERATE (q-pair).....	428, 435, 437, 562, 647
MODERATE SEDATION (data).....	969
MODERATE SEDATION (q-pair).....	31
MODERATE SEDATION: PRE-ASSESSMENT APPLIANCES REMOVED...	395
MODERATE SEDATION: PRE-ASSESSMENT CPRS HEALTH HX REVIEWED .....	396
MODERATE SEDATION: PRE-ASSESSMENT TIME OUT.....	397
MODERATE SEDATION: PREASSESS: PREVIOUS MOD SED. COMPLICATIONS .....	398
MODERATE SEDATION: PREASSESS: SUBSTANCE ABUSE.....	399
MODERATE SEDATION: PREPARATION.....	400
MODERATE SEDATION: TIME OUT: PATIENT IDENTIFIERS.....	401
MODERATE-SEVERELY REDUCED (data).....	969
MODERATE-SEVERELY REDUCED (q-pair).....	348
MODERATELY REDUCED (data).....	969
MODERATELY REDUCED (q-pair).....	348
MODIFIED ALDRETE SCORE (unit).....	1256
MODULAR CUSTOMIZED FORMULA (prod).....	1214
MODULAR CUSTOMIZED FORMULA (q-pair).....	172
MOIST (data).....	969
MOIST (q-pair).....	653
MONITOR (meth).....	1190
MONITOR (q-pair).....	63, 82, 83, 105, 570, 578, 692
MONITOR PLACED.....	402
MONITORING CAMERA.....	403
MONTH (unit).....	1256
MORSE SCORE (unit).....	1256
MORSE: AMBULATORY AID DURING WALKING.....	404
MORSE: GAIT.....	405
MORSE: HISTORY OF FALLING WITHIN LAST THREE MONTHS.....	406
MORSE: IV THERAPY SALINE LOCK.....	407
MORSE: MENTAL STATUS OF PT OWN ABILITY.....	408
MORSE: SCORE TOTAL.....	409
MORSE: SECONDARY DIAGNOSIS LISTED ON CHART.....	410
MOTOR ACTIVITY.....	411
MOTOR CHECK.....	412
MOTORIZED WHEELCHAIR (data).....	969
MOTORIZED WHEELCHAIR (q-pair).....	707
MOTTLED (data).....	970
MOTTLED (q-pair).....	645
MOUTHWASH USED (data).....	970
MOUTHWASH USED (q-pair).....	15
MOVE AGAINST GRAVITY AND RESISTANCE, EASILY OVERCOME BY THE EXAMINER (data) .....	970
MOVE AGAINST GRAVITY AND RESISTANCE, EASILY OVERCOME BY THE EXAMINER (q-pai r).....	411
MOVE AGAINST GRAVITY BUT NO RESISTANCE (data).....	970
MOVE AGAINST GRAVITY BUT NO RESISTANCE (q-pair).....	411
MOVEMENT WITH GRAVITY ELIMINATED (data).....	970
MOVEMENT WITH GRAVITY ELIMINATED (q-pair).....	411
MUCOID (data).....	971
MUCOID (q-pair).....	680
MUCOUS MEMBRANE CONDITION.....	413
MULTIFOCAL PVC (data).....	971

MULTIFOCAL PVC (q-pair).....	80
MULTIGATED NUCLEAR VENTRICULOGRAPHY (meth).....	1190
MULTIGATED NUCLEAR VENTRICULOGRAPHY (q-pair).....	347, 348
MULTIPODOUS (data).....	971
MULTIPODOUS (q-pair).....	544
MULTIVOLTS (unit).....	1256
MURMUR (data).....	971
MURMUR (q-pair).....	81
MUSHROOM CATHETER (data).....	971
MUSHROOM CATHETER (q-pair).....	219
MUSTACHE CANNULA (meth).....	1190
MUTIFOCAL ATRIAL TACHYCARDIA (data).....	972
MUTIFOCAL ATRIAL TACHYCARDIA (q-pair).....	84
mg/dl (unit).....	1260

- N -

N/A (data).....	972
N/A (q-pair).....	40, 196, 222, 295, 515, 516, 682, 715
NA (data).....	972
NA (q-pair).....	608
NAIL CARE (data).....	972
NAIL CARE (q-pair).....	5
NAIL COLOR.....	414
NAIL SHAPE.....	415
NAIL TEXTURE.....	416
NAPS ONLY (data).....	972
NAPS ONLY (q-pair).....	664
NASAL (data).....	973
NASAL (q-pair).....	567
NASAL CANNULA (meth).....	1190
NASAL CANNULA (q-pair).....	687
NASAL ENDOTRACHEAL (data).....	973
NASAL ENDOTRACHEAL (q-pair).....	22
NASAL INTESTINAL SMALL BORE (data).....	973
NASAL INTESTINAL SMALL BORE (q-pair).....	166
NASAL SMALL BOWEL FEEDING TUBE (data).....	973
NASAL SMALL BOWEL FEEDING TUBE (q-pair).....	215
NASAL TRUMPET (data).....	973
NASAL TRUMPET (q-pair).....	22
NASOGASTRIC LARGE BORE TUBE (data).....	974
NASOGASTRIC LARGE BORE TUBE (q-pair).....	166
NASOGASTRIC TUBE (data).....	974
NASOGASTRIC TUBE (q-pair).....	215
NASOGASTRIC TUBE SMALL BORE (data).....	974
NAUSEA (data).....	974
NAUSEA (q-pair).....	183, 209, 498
NEAT (data).....	974
NEAT (q-pair).....	32
NEBULIZER (meth).....	1190
NEBULIZER (q-pair).....	687
NECK (loca).....	1154
NECK (q-pair).....	116, 473, 658, 761
NECK VEIN CONDITION.....	417
NECROTIC (data).....	975
NECROTIC (q-pair).....	760
NEED FOR FLUIDS/TOILETING ASSESSED.....	418
NEEDLELESS HUB DISINFECTED AND ACCESSED PER POLICY (data) .....	975
NEEDLELESS HUB DISINFECTED AND ACCESSED PER POLICY (q-pair) .....	343

NEEDS REINFORCEMENT (data).....	975
NEEDS REINFORCEMENT (q-pair).....	591, 596
NEGATIVE PRESSURE (data).....	975
NEGATIVE PRESSURE (q-pair).....	333
NEGATIVE PRESSURE THERAPY DEVICE (data).....	975
NEGATIVE PRESSURE THERAPY DEVICE (q-pair).....	548
NEPRO (data).....	976
NEPRO (q-pair).....	172
NEPRO CARB STEADY (data).....	976
NEPRO CARB STEADY (q-pair).....	172
NEPRO RENAL FORMULA (data).....	976
NEPRO RENAL FORMULA (prod).....	1215
NEPRO RENAL FORMULA (q-pair).....	172, 309
NEUROGENIC BLADDER (data).....	976
NEUROGENIC BLADDER (q-pair).....	90
NEUTRAL PRESSURE CONNECTOR (data).....	976
NEUTRAL PRESSURE CONNECTOR (q-pair).....	343
NEUTROPENIC (data).....	977
NEUTROPENIC (q-pair).....	150, 333
NEW ADMISSION EDUCATION (data).....	977
NEW ADMISSION EDUCATION (q-pair).....	164
NEW SITE (data).....	977
NEW SITE (q-pair).....	94, 444
NG RESIDUAL (prod).....	1215
NG RESIDUAL (q-pair).....	430
NIGHT LIGHT PROVIDED (data).....	977
NIGHT LIGHT PROVIDED (q-pair).....	190
NLP: NO LIGHT PERCEPTION (data).....	977
NLP: NO LIGHT PERCEPTION (q-pair).....	152
NO (data).....	978
NO (q-pair).....	14, 16, 18, 23, 33, 34, 35, 36, 37, 38, 39, 41, 42, 45, 49, 51, 61, 75, 76, 78, 93, 98, 104, 107, 111, 112, 139, 189, 191, 201, 208, 217, 226, 228, 233, 234, 241, 295, 298, 303, 304, 305, 306, 307, 418, 440, 446, 451, 454, 456, 458, 466, 468, 482, 496, 505, 506, 531, 540, 543, 546, 549, 551, 561, 576, 583, 584, 585, 587, 588, 589, 590, 606, 616, 617, 637, 639, 641, 642, 643, 654, 655, 701, 709, 710, 718, 757, 764
NO BREAK IN EPIDERMIS (data).....	978
NO BREAK IN EPIDERMIS (q-pair).....	547
NO MOVEMENT (data).....	978
NO MOVEMENT (q-pair).....	411
NO SIGNS AND SYMPTOMS (data).....	978
NO SIGNS AND SYMPTOMS (q-pair).....	474
NO WANDERING PREVENTION MEASURES APPLIED (data).....	978
NO WANDERING PREVENTION MEASURES APPLIED (q-pair).....	4
NO WANDERING PREVENTION MEASURES NEEDED (data).....	979
NO WANDERING PREVENTION MEASURES NEEDED (q-pair).....	4
NON RE-BREATHING (meth).....	1191

NON RE-BREATHING (q-pair).....	687
NON-ADHERENT GAUZE (data).....	979
NON-INVASIVE (meth).....	1191
NON-INVASIVE (q-pair).....	63, 149, 377, 378, 570, 690
NON-PITTING (data).....	979
NON-PRODUCTIVE (data).....	979
NON-PRODUCTIVE (q-pair).....	137
NON-SKID SLIPPERS (data).....	979
NON-SKID SLIPPERS (q-pair).....	190
NON-TENDER (data).....	980
NON-TENDER (q-pair).....	2
NONBLANCHABLE ERYTHEMA (data).....	980
NONBLANCHABLE ERYTHEMA (q-pair).....	547
NONE (data).....	980
NONE (q-pair).....	31, 80, 133, 134, 135, 140, 156, 161, 183, 296, 301, 398, 428, 492, 539, 542, 545, 644, 647, 754, 758
NONE IDENTIFIED AT THIS VISIT (data).....	980
NONE IDENTIFIED AT THIS VISIT (q-pair).....	164
NONINVASIVE (meth).....	1191
NONREACTIVE (data).....	980
NONREACTIVE (q-pair).....	572
NORMAL (data).....	981
NORMAL (q-pair).....	1, 48, 103, 113, 221, 348, 415, 416, 417, 442, 450, 653
NORMAL FOR ETHNIC GROUP (data).....	981
NORMAL FOR ETHNIC GROUP (q-pair).....	414, 645
NORMAL POWER (data).....	981
NORMAL POWER (q-pair).....	411
NORMAL SINUS (data).....	981
NORMAL STOOL (data).....	981
NORMAL STOOL (q-pair).....	680
NOSE (loca).....	1154
NOSE (q-pair).....	473, 658
NOSE SIZE.....	419
NOT DONE (data).....	982
NOT DONE (q-pair).....	8, 9, 10, 15
NOT MONITORED (data).....	982
NOT MONITORED (q-pair).....	48, 103, 450
NOVASOURCE 2.0 (data).....	982
NOVASOURCE 2.0 (prod).....	1215
NOVASOURCE 2.0 (q-pair).....	309
NP (data).....	982
NP (q-pair).....	554, 608
NPO DIET (data).....	982
NPO DIET (q-pair).....	150
NSG ASSIST (data).....	983
NSG ASSIST (q-pair).....	554
NUCLEAR PERFUSION IMAGING (meth).....	1191
NUCLEAR PERFUSION IMAGING (q-pair).....	347, 348
NUMB (data).....	983
NUMB (q-pair).....	470
NUMBER OF DAYS.....	420
NUMBER: SPECIFY (data).....	983
NUMBER: SPECIFY (q-pair).....	156
NUMBNESS (data).....	983
NUMBNESS (q-pair).....	618, 619, 620, 621

NURSING INTERVENTIONS - STRAIGHT CATH (task).....	1237
NURSING INTERVENTIONS - WOUND CARE (task).....	1237
NURSING STATION (data).....	983
NURSING STATION (q-pair).....	55, 129
NUTREN 1.5 (prod).....	1215
NUTREN 1.5 (q-pair).....	172
NUTRIHEP (data).....	984
NUTRIHEP (q-pair).....	172
NUTRITION (data).....	984
NUTRITION (q-pair).....	164
NUTRITIONAL SUPPLEMENT (data).....	984
NUTRITIONAL SUPPLEMENT (prod).....	1215
NUTRITIONAL SUPPLEMENT (q-pair).....	329

- O -

O2 SATS (meth).....	1191
OBEYS COMMANDS (data).....	984
OBEYS COMMANDS (q-pair).....	134
OBSTRUCTIVE SLEEP APNEA (data).....	984
OCCIPITAL (loca).....	1155
OCCIPITAL (q-pair).....	303
OCCLUDED (data).....	985
OCCLUDED (q-pair).....	46, 99, 143, 238, 345, 447
OFF DURING BATHING (data).....	985
OFF DURING BATHING (q-pair).....	622
OFF WHILE AMBULATORY (data).....	985
OFF WHILE AMBULATORY (q-pair).....	622
OFFER PRN MEDS (data).....	985
OFFER PRN MEDS (q-pair).....	586
OFFERED (data).....	985
OFFERED (q-pair).....	15
ON WHILE SITTING IN CHAIR (data).....	986
ON WHILE SITTING IN CHAIR (q-pair).....	622
ONE PERSON ASSIST (data).....	986
ONE PERSON ASSIST (q-pair).....	6, 8, 9, 10
ONE TO ONE SUPERVISION (data).....	986
ONE TO ONE SUPERVISION (q-pair).....	586
OOZING (data).....	986
OOZING (q-pair).....	46, 99, 143, 238, 345, 447
OPEN SPONTANEOUSLY (data).....	986
OPEN SPONTANEOUSLY (q-pair).....	133
OPEN SYSTEM (data).....	987
OPEN SYSTEM (q-pair).....	167
OPEN TO AIR (data).....	987
OPEN TO AIR (q-pair).....	548, 765
OPERATING ROOM (data).....	987
OPERATING ROOM (q-pair).....	467
ORAL (data).....	987
ORAL (loca).....	1155
ORAL (q-pair).....	567, 694, 761
ORAL AIRWAY (data).....	987
ORAL AIRWAY (q-pair).....	22
ORAL ENDOTRACHEAL (data).....	988
ORAL ENDOTRACHEAL (q-pair).....	22
ORAL GASTRIC TUBE (data).....	988
ORAL GASTRIC TUBE (q-pair).....	166, 215
ORANGE (data).....	988
ORANGE (q-pair).....	157, 173, 211, 714, 753

ORIENTATION.....	421
ORIENTATION (data).....	988
ORIENTATION (q-pair).....	164
ORIENTED (data).....	988
ORIENTED (q-pair).....	135, 594
ORIENTED TO CALL FOR ASSISTANCE (data).....	989
ORIENTED TO CALL FOR ASSISTANCE (q-pair).....	190
ORIENTED TO PERSON (data).....	989
ORIENTED TO PERSON (q-pair).....	421
ORIENTED TO PLACE (data).....	989
ORIENTED TO PLACE (q-pair).....	421
ORIENTED TO SITUATION (data).....	989
ORIENTED TO SITUATION (q-pair).....	421
ORIENTED TO SURROUNDINGS/ROUTINE (data).....	989
ORIENTED TO SURROUNDINGS/ROUTINE (q-pair).....	190
ORIENTED TO TIME (data).....	990
ORIENTED TO TIME (q-pair).....	421
ORTHOPNEIC (data).....	990
ORTHOPNEIC (q-pair).....	577
ORTHOTIC BOOTS (data).....	990
ORTHOTIC BOOTS (q-pair).....	422
ORTHOTIC COLLAR (data).....	990
ORTHOTIC COLLAR (q-pair).....	422
ORTHOTIC KNEE IMMOBILIZER (data).....	990
ORTHOTIC KNEE IMMOBILIZER (q-pair).....	422
ORTHOTIC SPINE BRACE (data).....	991
ORTHOTIC SPINE BRACE (q-pair).....	422
ORTHOTIC TED HOSE (data).....	991
ORTHOTIC TED HOSE (q-pair).....	422
ORTHOTIC WEDGE PILLOW (data).....	991
ORTHOTIC WEDGE PILLOW (q-pair).....	422
ORTHOTICS.....	422
OSMOLITE (data).....	991
OSMOLITE (prod).....	1216
OSMOLITE 1.0 (data).....	991
OSMOLITE 1.0 (q-pair).....	172
OSMOLITE 1.2 (data).....	992
OSMOLITE 1.2 (q-pair).....	172
OSMOLITE 1.5 (data).....	992
OSMOLITE 1.5 (q-pair).....	172
OSMOLITE HN (data).....	992
OSMOLITE HN (prod).....	1216
OSMOLITE HN (q-pair).....	172
OSTOMY CARE.....	423
OSTOMY CARE (data).....	992
OSTOMY CARE (q-pair).....	164
OSTOMY TYPE.....	424
OTHER (ADD COMMENTS) (data).....	993
OTHER (ADD COMMENTS) (q-pair).....	544
OTHER (SEE COMMENTS) (data).....	993
OTHER (SEE COMMENTS) (q-pair).....	395, 398, 554
OTHER (data).....	992
OTHER (meth).....	1192
OTHER (q-pair).....	52, 90, 164, 183, 186, 187, 190, 316, 347, 348, 431, 537, 553
OTHER (qual).....	1224
OTHER PO FLUID (prod).....	1216
OTHER PO FLUID (q-pair).....	329
OUNCE (unit).....	1256
OUTPUT - BLOOD.....	425

OUTPUT - CBI DRAINAGE.....	426
OUTPUT - CHEST TUBE.....	427
OUTPUT - CHEST TUBE ESTIMATED.....	428
OUTPUT - DIALYSIS.....	429
OUTPUT - DO NOT CHECK TOTAL OR SUBTOTALS.....	430
OUTPUT - DRAIN.....	431
OUTPUT - EMESIS.....	432
OUTPUT - GI TUBE.....	433
OUTPUT - STOOL.....	434
OUTPUT - STOOL ESTIMATED.....	435
OUTPUT - URINE.....	436
OUTPUT - URINE ESTIMATED.....	437
OVER WIRE (data).....	993
OVER WIRE (q-pair).....	94, 444
OXEPA (data).....	993
OXEPA (q-pair).....	172
OXIMETRIC (data).....	993
OXIMETRIC (q-pair).....	449
OXYGEN CONSUMPTION.....	438
OXYGEN SATURATION.....	439

- P -

P.O. (meth).....	1192
P.O. (q-pair).....	309
PA (data).....	994
PA (q-pair).....	554, 608
PA CATHETER (meth).....	1192
PA CATHETER ALARMS ON.....	440
PA CATHETER DRESSING TYPE.....	441
PA CATHETER DYNAMIC RESPONSE TEST.....	442
PA CATHETER INSERTION MARK.....	443
PA CATHETER INSERTION TECHNIQUE.....	444
PA CATHETER REPOSITIONED/DISCONTINUED REASON.....	445
PA CATHETER SITE CARE.....	446
PA CATHETER SITE CONDITION.....	447
PA CATHETER SIZE.....	448
PA CATHETER TYPE.....	449
PA CATHETER WAVEFORM.....	450
PA CATHETER ZEROED.....	451
PA LINE (meth).....	1192
PAC (data).....	994
PAC (q-pair).....	80
PACED: ATRIAL (data).....	994
PACED: ATRIAL (q-pair).....	84
PACED: AV SEQUENTIAL (data).....	994
PACED: AV SEQUENTIAL (q-pair).....	84
PACED: VENTRICULAR (data).....	994
PACED: VENTRICULAR (q-pair).....	84
PACEMAKER ATRIAL mA.....	452
PACEMAKER AV INTERVAL.....	453
PACEMAKER CAPTURE.....	454
PACEMAKER MODE.....	455
PACEMAKER ON.....	456
PACEMAKER RATE - DO NOT CHECK TOTAL OR SUBTOTALS.....	457
PACEMAKER SENSE.....	458
PACEMAKER SENSITIVITY ATRIAL.....	459
PACEMAKER SENSITIVITY VENTRICULAR.....	460
PACEMAKER TYPE.....	461
PACEMAKER UNIT ID.....	462
PACEMAKER VENTRICULAR mA.....	463
PACEMAKER WIRES.....	464

PACER (data).....	995
PACER (q-pair).....	396, 449
PACU - DISCHARGE CRITERIA.....	465
PACU - DISCHARGE CRITERIA MET.....	466
PACU - DISPOSITION.....	467
PACU - PT RECEIVED REVERSAL AGENT.....	468
PACU - TIME OF REVERSAL AGENT.....	469
PAIN (data).....	995
PAIN (q-pair).....	161, 164
PAIN CHARACTER.....	470
PAIN CONTROLLED PER PT (data).....	995
PAIN CONTROLLED PER PT (q-pair).....	465
PAIN DURATION.....	471
PAIN GOAL.....	472
PAIN LEVEL.....	473
PAIN NON-VERBAL CUES.....	474
PAIN RELIEVING FACTOR INTERVENTION.....	475
PAINAD BL 0 (data).....	995
PAINAD BL 0 (q-pair).....	476
PAINAD BL 1 (data).....	995
PAINAD BL 1 (q-pair).....	476
PAINAD BL 2 (data).....	996
PAINAD BL 2 (q-pair).....	476
PAINAD BODY LANGUAGE.....	476
PAINAD BR 0 (data).....	996
PAINAD BR 0 (q-pair).....	477
PAINAD BR 1 (data).....	996
PAINAD BR 1 (q-pair).....	477
PAINAD BR 2 (data).....	996
PAINAD BR 2 (q-pair).....	477
PAINAD BREATHING.....	477
PAINAD CO 0 (data).....	996
PAINAD CO 0 (q-pair).....	478
PAINAD CO 1 (data).....	997
PAINAD CO 1 (q-pair).....	478
PAINAD CO 2 (data).....	997
PAINAD CO 2 (q-pair).....	478
PAINAD CONSOLABILITY.....	478
PAINAD FACIAL EXPRESSION.....	479
PAINAD FE 0 (data).....	997
PAINAD FE 0 (q-pair).....	479
PAINAD FE 1 (data).....	997
PAINAD FE 1 (q-pair).....	479
PAINAD FE 2 (data).....	997
PAINAD FE 2 (q-pair).....	479
PAINAD NEGATIVE VOCALIZATION.....	480
PAINAD NV 0 (data).....	998
PAINAD NV 0 (q-pair).....	480
PAINAD NV 1 (data).....	998
PAINAD NV 1 (q-pair).....	480
PAINAD NV 2 (data).....	998
PAINAD NV 2 (q-pair).....	480
PAINAD SCALE (unit).....	1257
PAINAD TOTAL.....	481
PAINFUL EXPRESSION (data).....	998
PAINFUL EXPRESSION (q-pair).....	474
PAIRS (data).....	998
PAIRS (q-pair).....	80
PALE (data).....	999
PALE (q-pair).....	413, 414, 645
PALPATED (meth).....	1192
PALPATED (q-pair).....	63, 570



PAROXYSMAL ATRIAL TACHYCARDIA (data).....	999
PAROXYSMAL ATRIAL TACHYCARDIA (q-pair).....	84
PAROXYSMAL SUPRAVENTRICULAR TACHYCARDIA (data).....	999
PAROXYSMAL SUPRAVENTRICULAR TACHYCARDIA (q-pair).....	84
PART THICKNESS SKIN LOSS INVOLVE EPIDERMIS/DERMIS (data)	
.....	999
PART THICKNESS SKIN LOSS INVOLVE EPIDERMIS/DERMIS (q-pair)	
.....	547
PARTIAL BATH (data).....	999
PARTIAL BATH (q-pair).....	5
PARTIAL IMPAIRED (data).....	1000
PARTIAL IMPAIRED (q-pair).....	613
PARTIAL LOWER (data).....	1000
PARTIAL LOWER (q-pair).....	140
PARTIAL RE-BREATHER (meth).....	1193
PARTIAL RE-BREATHER (q-pair).....	687
PARTIAL SIMPLE (data).....	1000
PARTIAL SIMPLE (q-pair).....	613
PARTIAL UPPER (data).....	1000
PARTIAL UPPER (q-pair).....	140
PASSIVE RANGE OF MOTION (data).....	1000
PASSIVE RANGE OF MOTION (q-pair).....	12
PATENT (data).....	1001
PATENT (q-pair).....	46, 99, 143, 238,
	345, 447
PATENT CLEAN (data).....	1001
PATENT CLEAN (q-pair).....	343, 493
PATIENT (data).....	1001
PATIENT (q-pair).....	11
PATIENT ABLE TO RESPOND.....	482
PATIENT EDUCATION PROVIDED PRIOR TO SEDATION.....	483
PATIENT INFORMED (data).....	1001
PATIENT INFORMED (q-pair).....	403, 580
PATIENT INSTRUCTED (data).....	1001
PATIENT INSTRUCTED (q-pair).....	423
PATIENT LOCATOR DEVICE IN PLACE (data).....	1002
PATIENT LOCATOR DEVICE IN PLACE (q-pair).....	4
PATIENT ON PASS (null).....	1201
PATIENT REFUSED (null).....	1201
PATIENT REFUSES MARKING (data).....	1002
PATIENT REFUSES MARKING (q-pair).....	531
PATIENT REPORTS NOT PRESENT (data).....	1002
PATIENT REPORTS NOT PRESENT (q-pair).....	194
PATIENT REPORTS PRESENT (data).....	1002
PATIENT REPORTS PRESENT (q-pair).....	194
PATIENT RIGHTS AND RESPONSIBILITIES (data).....	1002
PATIENT RIGHTS AND RESPONSIBILITIES (q-pair).....	164
PATIENT SIGNS SAFETY CONTRACT (data).....	1003
PATIENT SIGNS SAFETY CONTRACT (q-pair).....	57
PATIENT STATED (data).....	1003
PATIENT STATED (q-pair).....	509
PATIENT STAYS WITHIN LINE OF SITE OF NURSING STATION (data)	
.....	1003
PATIENT STAYS WITHIN LINE OF SITE OF NURSING STATION (q-pair)	
.....	57
PATIENT'S ROOM (data).....	1003
PATIENT'S ROOM (q-pair).....	55, 129
PCA (data).....	1003
PCA (prod).....	1216
PCA (q-pair).....	309
PCA - CLINICAL BOLUS DOSE.....	484
PCA - CONTINUOUS MCG/HR (BASAL).....	485

PCA - CONTINUOUS MG/HR (BASAL).....	486
PCA - DEMAND DOSE - DO NOT CHECK TOTAL OR SUBTOTALS....	487
PCA - DOSAGE DELIVERED.....	488
PCA - DOSES DELIVERED.....	489
PCA - LEVEL OF CONSCIOUSNESS.....	490
PCA - LIMIT DOSE - DO NOT CHECK TOTAL OR SUBTOTALS....	491
PCA - LIMIT TIME.....	492
PCA - LINE CARE.....	493
PCA - LOCK OUT- DO NOT CHECK TOTAL OR SUBTOTALS.....	494
PCA - MODE.....	495
PCA - PUMP CLEARED.....	496
PCA - REQUESTS.....	497
PCA - SIDE EFFECTS.....	498
PCA - TOTAL DOSE 24 HR.....	499
PCA - TOTAL MCG PER SHIFT.....	500
PCA - TOTAL MG PER SHIFT.....	501
PCA - VOLUME REMAINING - DO NOT CHECK TOTAL OR SUBTOTALS	
.....	502
PCA - VOLUME START - DO NOT CHECK TOTAL OR SUBTOTALS...	503
PCA AND BASAL RATE (data).....	1004
PCA AND BASAL RATE (q-pair).....	495
PCA ONLY (data).....	1004
PCA ONLY (q-pair).....	495
PCA PUMP (data).....	1004
PCA PUMP (q-pair).....	164
PEAK INSPIRATORY PRESSURE - DO NOT CHECK TOTAL OR SUBTOTALS	
.....	504
PEDIATRIC (data).....	1004
PEDIATRIC (q-pair).....	64
PEDIATRIC (qual).....	1224
PEDIATRIC CUFF (meth).....	1193
PEG TUBE (data).....	1004
PEG TUBE (q-pair).....	166, 215
PENDANT CANNULA (meth).....	1193
PENDANT CANNULA (q-pair).....	687
PENETRATING (data).....	1005
PENETRATING (q-pair).....	766
PENIS (loca).....	1155
PENIS (q-pair).....	473, 658
PENROSE (meth).....	1193
PENROSE (q-pair).....	431
PENROSE DRAIN (data).....	1005
PENROSE DRAIN (q-pair).....	755
PERATIVE (data).....	1005
PERATIVE (q-pair).....	172
PERCENTAGE (unit).....	1257
PERCUSSION BED (data).....	1005
PERCUSSION BED (q-pair).....	665
PERINEAL (loca).....	1155
PERINEAL (q-pair).....	303, 473, 547, 761
PERINEAL CARE (data).....	1005
PERINEAL CARE (q-pair).....	5
PERIORBITAL (loca).....	1155
PERIORBITAL (q-pair).....	160
PERIPHERAL (loca).....	1156
PERIPHERAL (q-pair).....	570
PERIRECTAL (loca).....	1156
PERITONEAL DIALYSIS (meth).....	1193
PERITONEAL DIALYSIS (q-pair).....	311, 429
PERMANENT (data).....	1006
PERMANENT (q-pair).....	461
PERSONAL ARTICLES AT REACH (data).....	1006

PERSONAL ARTICLES AT REACH (q-pair).....	190
PERSONAL HYGIENE AND/OR GROOMING (data).....	1006
PERSONAL HYGIENE AND/OR GROOMING (q-pair).....	164
PERSONAL MEDICAL ALERT BRACELET (data).....	1006
PERSONAL MEDICAL ALERT BRACELET (q-pair).....	40
PETROLATUM (data).....	1006
PETROLATUM (q-pair).....	548
PHYSICAL IMPAIRMENT (data).....	1007
PHYSICAL IMPAIRMENT (q-pair).....	52
PHYSICIAN (data).....	1007
PHYSICIAN (q-pair).....	554, 706
PHYSICIAN NOTIFIED (data).....	1007
PHYSICIAN NOTIFIED (q-pair).....	475
PHYSICIAN ORDER (data).....	1007
PHYSICIAN ORDER (q-pair).....	592
PICC (data).....	1007
PICC (q-pair).....	102, 343
PICC LINE (meth).....	1194
PIN (data).....	1008
PIN TRACTION (data).....	1008
PIN TRACTION (q-pair).....	702
PINK (data).....	1008
PINK (q-pair).....	157, 173, 211, 413, 563, 645, 714, 749, 753
PINK / GRANULATING (data).....	1008
PINK / GRANULATING (q-pair).....	537
PINPOINT (data).....	1008
PINPOINT (q-pair).....	572
PITTED (data).....	1009
PITTED (q-pair).....	416
PITTING (data).....	1009
PIVOT (data).....	1009
PIVOT (q-pair).....	172
PIVOT 1.5 (data).....	1009
PIVOT 1.5 (q-pair).....	172
PJC (data).....	1009
PJC (q-pair).....	80
PLACED NEAR TO NURSE STATION (data).....	1010
PLACED NEAR TO NURSE STATION (q-pair).....	190
PLAIN (data).....	1010
PLAIN (q-pair).....	449
PLAN OF CARE - REVISION INDICATED.....	505
PLANTAR FLEXION.....	506
PLASMAPHERESIS (prod).....	1216
PLASMAPHERESIS (q-pair).....	317
PLATELETS (prod).....	1217
PLATELETS (q-pair).....	309, 317
PLATELETS PHERESIS (prod).....	1217
PLATELETS PHERESIS (q-pair).....	309, 317
PLATELETS PHERESIS IRRADIATED (prod).....	1217
PLATELETS PHERESIS IRRADIATED (q-pair).....	309, 317
POINTS (unit).....	1257
POOR (data).....	1010
POOR (q-pair).....	661
POOR BAL (data).....	1010
POOR BAL (q-pair).....	183
POPLITEAL (loca).....	1156
POPLITEAL (q-pair).....	570
PORCINE SUBMUCOSA (data).....	1010
PORCINE SUBMUCOSA (q-pair).....	548
POSITIONAL (data).....	1011

POSITIONAL (q-pair).....	48, 103, 450, 470
POSITIVE PRESSURE CONNECTOR (data).....	1011
POSITIVE PRESSURE CONNECTOR (q-pair).....	343
POST CATH DIET (data).....	1011
POST CATH DIET (q-pair).....	150
POST OPERATION DIAGNOSIS (data).....	1011
POST PROCEDURE DIAGNOSIS (data).....	1011
POST PROCEDURE DIAGNOSIS (q-pair).....	142
POST-FALL (data).....	1012
POST-FALL (q-pair).....	761
POST-OPERATIVE (data).....	1012
POST-OPERATIVE (q-pair).....	142, 164, 761
POST-PROCEDURE (data).....	1012
POST-PROCEDURE (q-pair).....	164
POST-SEDATION (data).....	1012
POST-SEDATION (q-pair).....	164
POSTERIOR TIBIAL (loca).....	1156
POSTERIOR TIBIAL (q-pair).....	570
POTASSIUM MODIFIED 60 MEQ (data).....	1012
POTASSIUM MODIFIED 60 MEQ (q-pair).....	150
POUND (unit).....	1257
PPN (prod).....	1217
PPN (q-pair).....	309
PR INTERVAL.....	507
PRBC (data).....	1013
PRBC (prod).....	1217
PRBC (q-pair).....	309, 317
PRBC - AUTOLOGOUS (prod).....	1218
PRBC - AUTOLOGOUS (q-pair).....	309, 317
PRBC - LEUKOREduced (prod).....	1218
PRBC - LEUKOREduced (q-pair).....	309, 317
PRBC - LEUKOREduced - DIRECTED (prod).....	1218
PRBC - LEUKOREduced - DIRECTED (q-pair).....	309, 317
PRBC - LEUKOREduced - IRRADIATED (prod).....	1218
PRBC - LEUKOREduced - IRRADIATED (q-pair).....	309, 317
PRBC - LEUKOREduced - WASHED (prod).....	1218
PRBC - LEUKOREduced - WASHED (q-pair).....	309, 317
PRE-MEDICATED FOR ACTIVITY (data).....	1013
PRE-MEDICATED FOR ACTIVITY (q-pair).....	475
PRE-OPERATIVE (data).....	1013
PRE-OPERATIVE (q-pair).....	164
PRE-PROCEDURE (data).....	1013
PRE-PROCEDURE (q-pair).....	164
PRE-SEDATION (data).....	1013
PRE-SEDATION (q-pair).....	164
PRE-SURGERY (data).....	1014
PRECAUTIONS - 24 HOUR CHART CHECK (task).....	1237
PRECAUTIONS - ASPIRATION (task).....	1237
PRECAUTIONS - ETOH W/D (task).....	1238
PRECAUTIONS - FALL (task).....	1238
PRECAUTIONS - PATIENT TEACHING (task).....	1238
PRECAUTIONS - VALIDATE ALARMS (task).....	1238
PREGNANT, NO (data).....	1014
PREGNANT, NO (q-pair).....	396
PREGNANT, YES (data).....	1014
PREGNANT, YES (q-pair).....	396
PREMATURE ATRIAL (data).....	1014
PREPROCEDURE: 12 LEAD EKG DONE?.....	508
PREPROCEDURE: ALLERGIES REASSESSED AND DOCUMENTED.....	509
PREPROCEDURE: ANTIBIOTIC GIVEN?.....	510
PREPROCEDURE: ANTIBIOTIC ORDERED?.....	511
PREPROCEDURE: BLOOD TRANSFUSION CONSENT CURRENT?.....	512

PREPROCEDURE: CHEST XRAY DONE?.....	513
PREPROCEDURE: DENTURES REMOVED?.....	514
PREPROCEDURE: DISPOSITION OF CLOTHING.....	515
PREPROCEDURE: DISPOSITION OF VALUABLES.....	516
PREPROCEDURE: DRIVER/ESCORT TO HOME CONFIRMED.....	517
PREPROCEDURE: FINGERSTICK GLUCOSE DONE?.....	518
PREPROCEDURE: GLASSES / CONTACTS REMOVED?.....	519
PREPROCEDURE: H AND P CURRENT?.....	520
PREPROCEDURE: HEARING AIDS IN PLACE?.....	521
PREPROCEDURE: HOME MEDICATIONS TAKEN/TIME.....	522
PREPROCEDURE: HOSPITAL GOWN ON?.....	523
PREPROCEDURE: INFORMED CONSENT CURRENT .....	524
PREPROCEDURE: JEWELRY REMOVED?.....	525
PREPROCEDURE: LABS DONE?.....	526
PREPROCEDURE: NPO SINCE (DATE/TIME).....	527
PREPROCEDURE: PLANNED PROCEDURE.....	528
PREPROCEDURE: PREMEDICATION ADMINISTERED?.....	529
PREPROCEDURE: PREMEDICATION ORDERED?.....	530
PREPROCEDURE: PROCEDURE SITE MARKED CORRECTLY.....	531
PREPROCEDURE: PROSTHETICS/IMPLANTS?.....	532
PREPROCEDURE: TYPE AND CROSS MATCH DONE?.....	533
PREPROCEDURE: TYPE AND SCREEN DONE?.....	534
PREPROCEDURE: VITALS SIGNS DONE?.....	535
PREPROCEDURE: VOIDED?.....	536
PRESENT (data).....	1014
PRESENT (q-pair).....	79, 86, 136, 154, 194, 198, 350, 708, 711
PRESSURE (data).....	1015
PRESSURE (q-pair).....	235, 470, 758
PRESSURE CALIBRATION (data).....	1015
PRESSURE CALIBRATION (q-pair).....	343
PRESSURE CONTROL (data).....	1015
PRESSURE CONTROL (q-pair).....	725
PRESSURE DRESSING (data).....	1015
PRESSURE DRESSING (q-pair).....	44, 96, 144, 441
PRESSURE ULCER (data).....	1015
PRESSURE ULCER (q-pair).....	766
PRESSURE ULCER COLOR.....	537
PRESSURE ULCER DEPTH.....	538
PRESSURE ULCER DRAINAGE CHARACTER.....	539
PRESSURE ULCER DRESSING PREVENTS ASSESSMENT.....	540
PRESSURE ULCER LENGTH.....	541
PRESSURE ULCER ODOR.....	542
PRESSURE ULCER PHOTOGRAPHED.....	543
PRESSURE ULCER PREVENTION DEVICE TYPE.....	544
PRESSURE ULCER RISK FACTORS.....	545
PRESSURE ULCER SINUS/UNDERMINING.....	546
PRESSURE ULCER STAGE/SITE.....	547
PRESSURE ULCER TREATMENT.....	548
PRESSURE ULCER TUNNEL.....	549
PRESSURE ULCER TUNNEL LENGTH.....	550
PRESSURE ULCER UNDERMINING.....	551
PRESSURE ULCER WIDTH.....	552
PRESYSTOLIC (data).....	1016
PRESYSTOLIC (q-pair).....	87
PREVENT ACCESS TO DEVICE (data).....	1016
PREVENT ACCESS TO DEVICE (q-pair).....	586
PREVENTIVE HEALTH REHAB TECHNIQUES (data).....	1016
PREVENTIVE HEALTH REHAB TECHNIQUES (q-pair).....	164
PREVIOUS PERCUTANEOUS SITE (data).....	1016
PREVIOUS PERCUTANEOUS SITE (q-pair).....	766

PRN (data).....	1016
PRN (q-pair).....	216
PRN (sche).....	1226
PRO STAT SUGAR FREE 64 (data).....	1017
PRO STAT SUGAR FREE 64 (q-pair).....	172
PRO-STAT (data).....	1017
PRO-STAT (q-pair).....	172
PROCEDURE LOCATION.....	553
PROCEDURE PERSONNEL.....	554
PROCEDURE ROOM (data).....	1017
PROCEDURE ROOM (q-pair).....	467
PROCEDURE TIME.....	555
PROCEDURE, PLANNED.....	556
PRODUCTIVE (data).....	1017
PRODUCTIVE (q-pair).....	137
PROJECTILE VOMITING (data).....	1017
PROJECTILE VOMITING (q-pair).....	209
PROMOTE (data).....	1018
PROMOTE (q-pair).....	172
PROMOTE WITH FIBER (data).....	1018
PROMOTE WITH FIBER (q-pair).....	172
PRONE (data).....	1018
PRONE (q-pair).....	13
PROPASS PROTEIN POWDER (data).....	1018
PROPASS PROTEIN POWDER (q-pair).....	172
PROSTHESIS (data).....	1018
PROSTHESIS (q-pair).....	152, 395
PROTENEX WOUND CARE (data).....	1019
PROTENEX WOUND CARE (q-pair).....	172
PROVIDED PATIENT EDUCATION (data).....	1019
PROVIDED PATIENT EDUCATION (q-pair).....	91
PROXIMAL PORT (data).....	1019
PROXIMAL PORT (q-pair).....	344
PT DENIES NAUSEA / NO VOMITING (data).....	1019
PT DENIES NAUSEA / NO VOMITING (q-pair).....	465
PT PAIN AT ACCEPTABLE LEVEL (data).....	1019
PT PAIN AT ACCEPTABLE LEVEL (q-pair).....	465
PULMOCARE (data).....	1020
PULMOCARE (q-pair).....	172
PULMONARY - CHEST PHYSIO-THERAPY (task).....	1238
PULMONARY - COUGH AND DEEP BREATH (task).....	1239
PULMONARY - ET TUBE CARE (task).....	1239
PULMONARY - INCENTIVE SPIROMETER (task).....	1239
PULMONARY - OXYGEN (task).....	1239
PULMONARY - SUCTION PATIENT (task).....	1239
PULMONARY - TRACH CARE (task).....	1240
PULMONARY ARTERY DIASTOLIC.....	557
PULMONARY ARTERY MEAN PRESSURE.....	558
PULMONARY ARTERY SYSTOLIC.....	559
PULMONARY CAPILLARY WEDGE PRESSURE.....	560
PULMONARY CARE.....	561
PULMONARY PROCEDURE ROOM (data).....	1020
PULMONARY PROCEDURE ROOM (q-pair).....	553
PULMONARY SECRETION AMOUNT.....	562
PULMONARY SECRETION COLOR.....	563
PULMONARY SECRETION CONSISTENCY.....	564
PULMONARY SUCTIONING EVENT.....	565
PULMONARY SUCTIONING FREQUENCY.....	566
PULMONARY SUCTIONING SITE.....	567
PULMONARY VASCULAR RESISTANCE.....	568
PULMONARY VASCULAR RESISTANCE INDEX.....	569
PULSE - DO NOT CHECK TOTAL OR SUBTOTALS.....	570

PULSE OX MONITOR (meth).....1194  
 PULSE OX MONITOR (q-pair).....692  
 PULSE OXIMETER (data).....1020  
 PULSE OXIMETER (q-pair).....400, 402  
 PULSE STRENGTH.....571  
 PULSELESS ELECTRICAL ACTIVITY (data).....1020  
 PULSELESS ELECTRICAL ACTIVITY (q-pair).....84  
 PUMP - CONTINUOUS (data).....1020  
 PUMP - CONTINUOUS (q-pair).....167  
 PUMP - INTERMITTENT (data).....1021  
 PUMP - INTERMITTENT (q-pair).....167  
 PUNCTURE (data).....1021  
 PUNCTURE (q-pair).....766  
 PUPIL REACTION.....572  
 PUPIL SIZE.....573  
 PUREED (data).....1021  
 PUREED (q-pair).....150  
 PURULENT (data).....1021  
 PURULENT (q-pair).....109, 157, 173, 211,  
   296, 539, 713  
 PVC (data).....1021  
 PVC (q-pair).....80  
 pH (data).....1125  
 pH (q-pair).....212  
 pH (unit).....1260

- Q -

Q2H (sche).....1226  
 Q4H (sche).....1226  
 Q7D (sche).....1226  
 QD (sche).....1227  
 QRS.....574  
 QT.....575  
 QUAD LUMEN CATHETER (unit).....1257  
 QUADRAGEMINY (data).....1022  
 QUADRAGEMINY (q-pair).....80  
 QUIET OR RESTING (data).....1022  
 QUIET OR RESTING (q-pair).....128  
 QUIET ROOM (data).....1022  
 QUIET ROOM (q-pair).....55, 129, 475

- R -

R ABDOMEN (loca).....1156  
 R ABDOMEN (q-pair).....761  
 R ANKLE (loca).....1157  
 R ANKLE (q-pair).....303, 473, 547, 600,  
   658, 761  
 R ANTECUBITAL (loca).....1157  
 R ANTECUBITAL (q-pair).....102, 143, 219, 344,  
   345, 547, 761  
 R ARM (loca).....1157  
 R ARM (q-pair).....63, 102, 143, 149,  
   160, 165, 219, 344,  
   345, 473, 598, 599,  
   658, 690  
 R ARTERY (loca).....1157  
 R AXILLARY (loca).....1157  
 R AXILLARY (q-pair).....46, 47, 235, 239  
 R BRACHIAL (loca).....1158  
 R BRACHIAL (q-pair).....46, 47, 105, 456,

R BUNDLE BRANCH BLOCK (data).....	1022
R BUNDLE BRANCH BLOCK (q-pair).....	84
R CALF (loca).....	1158
R CALF (q-pair).....	303
R CAROTID (loca).....	1158
R CAROTID (q-pair).....	79, 86
R CHEST (loca).....	1158
R CHEST (q-pair).....	108, 473, 658, 761
R DORSALIS PEDIS (loca).....	1158
R EAR (loca).....	1159
R EAR (q-pair).....	303, 473, 547, 761
R ELBOW (loca).....	1159
R ELBOW (q-pair).....	473, 547, 658, 761
R EXTERNAL JUGULAR (loca).....	1159
R EXTERNAL JUGULAR (q-pair).....	102, 143, 219, 344, 345, 449
R EYE (loca).....	1159
R EYE (q-pair).....	133, 160, 572, 573, 741, 761
R FACE (loca).....	1159
R FACE (q-pair).....	160, 303, 761
R FEMORAL (loca).....	1160
R FEMORAL (q-pair).....	46, 47, 79, 86, 102, 105, 143, 219, 235, 239, 303, 344, 345, 449, 456, 461, 547, 571, 761
R FINGER (loca).....	1160
R FINGER (q-pair).....	303
R FLANK (loca).....	1160
R FLANK (q-pair).....	303
R FOOT (loca).....	1160
R FOOT (q-pair).....	160, 196, 303, 344, 345, 473, 547, 618, 658, 761
R FOOT DOWN (data).....	1022
R FOOT DOWN (q-pair).....	634
R FOOT UP (data).....	1023
R FOOT UP (q-pair).....	634
R GLUTEUS (loca).....	1160
R GLUTEUS (q-pair).....	303, 547, 761
R GROIN (loca).....	1161
R GROIN (q-pair).....	473, 658
R HAND (loca).....	1161
R HAND (q-pair).....	160, 222, 473, 619, 658, 761
R HEAD (loca).....	1161
R HEAD (q-pair).....	547, 761
R HEAD DOWN (data).....	1023
R HEAD DOWN (q-pair).....	634
R HEAD UP (data).....	1023
R HEAD UP (q-pair).....	634
R HEEL (loca).....	1161
R HEEL (q-pair).....	547, 761
R HIP (loca).....	1161
R HIP (q-pair).....	303, 473, 658, 761
R INGUINAL (loca).....	1162
R INGUINAL (q-pair).....	303
R INNER KNEE (loca).....	1162
R INNER KNEE (q-pair).....	547, 761
R INTERNAL JUGULAR (loca).....	1162



R INTERNAL JUGULAR (q-pair).....102, 105, 143, 219,  
344, 345, 449, 456,  
461  
R ISCHIAL (loca).....1162  
R ISCHIAL (q-pair).....547, 761  
R JAW (loca).....1162  
R JAW (q-pair).....473, 658  
R KNEE (loca).....1163  
R KNEE (q-pair).....160, 303, 473, 658  
R LATERAL HEEL (loca).....1163  
R LATERAL HEEL (q-pair).....547  
R LATERAL MALLEOLUS (loca).....1163  
R LATERAL MALLEOLUS (q-pair).....547  
R LEG (loca).....1163  
R LEG (q-pair).....63, 136, 149, 160,  
165, 344, 345, 473,  
598, 600, 658, 690  
R LOWER ANTERIOR CHEST (loca).....1163  
R LOWER ANTERIOR CHEST (q-pair).....110, 303  
R LOWER ARM (loca).....1164  
R LOWER ARM (q-pair).....303, 473, 547, 658,  
761  
R LOWER EXTREMITY (loca).....1164  
R LOWER EXTREMITY (q-pair).....411, 620  
R LOWER LEG (loca).....1164  
R LOWER LEG (q-pair).....303, 761  
R LOWER LIMB (loca).....1164  
R LOWER LIMB (q-pair).....352  
R LOWER LOBE (loca).....1164  
R LOWER LOBE (q-pair).....74  
R LOWER POSTERIOR CHEST (loca).....1165  
R LOWER POSTERIOR CHEST (q-pair).....110, 303  
R MANDIBULAR (loca).....1165  
R MANDIBULAR (q-pair).....303  
R MEDIAL HEEL (loca).....1165  
R MEDIAL HEEL (q-pair).....547  
R MEDIAL MALLEOLUS (loca).....1165  
R MEDIAL MALLEOLUS (q-pair).....547  
R MIDDLE ANTERIOR CHEST (loca).....1165  
R MIDDLE ANTERIOR CHEST (q-pair).....110, 303  
R MIDDLE LOBE (loca).....1166  
R MIDDLE LOBE (q-pair).....74  
R MIDDLE POSTERIOR CHEST (loca).....1166  
R MIDDLE POSTERIOR CHEST (q-pair).....110, 303  
R NECK (data).....1023  
R NECK (loca).....1166  
R NECK (q-pair).....160, 303  
R NEPHROSTOMY (data).....1023  
R NEPHROSTOMY (q-pair).....219, 716  
R NOSE (loca).....1166  
R NOSE (q-pair).....473, 658  
R OCCIPITAL (loca).....1166  
R OCCIPITAL (q-pair).....303  
R ORAL (loca).....1167  
R ORAL (q-pair).....303  
R ORBIT (loca).....1167  
R ORBIT (q-pair).....303  
R OUTER KNEE (loca).....1167  
R OUTER KNEE (q-pair).....547, 761  
R PARIETAL (loca).....1167  
R PARIETAL (q-pair).....303  
R PEDAL (loca).....1167

R PEDAL (q-pair).....	571
R POST TIBIAL (loca).....	1168
R POST TIBIAL (q-pair).....	571
R RADIAL (loca).....	1168
R RADIAL (q-pair).....	46, 47, 571
R RIBS (loca).....	1168
R RIBS (q-pair).....	547, 761
R SCAPULA (loca).....	1168
R SCAPULA (q-pair).....	473, 547, 658, 761
R SHOULDER (loca).....	1168
R SHOULDER (q-pair).....	473, 547, 658, 761
R STENT (data).....	1024
R STENT (q-pair).....	219
R SUBCLAVIAN (loca).....	1169
R SUBCLAVIAN (q-pair).....	102, 105, 143, 219, 303, 344, 345, 449, 456, 461
R TEMPORAL (loca).....	1169
R TEMPORAL (q-pair).....	79, 86
R TESTICLE (loca).....	1169
R TESTICLE (q-pair).....	761
R THIGH (loca).....	1169
R THIGH (q-pair).....	303, 473, 658
R THORACIC ANTERIOR (loca).....	1169
R THORACIC ANTERIOR (q-pair).....	303
R THORACIC LATERAL (loca).....	1170
R THORACIC LATERAL (q-pair).....	303
R THORACIC POSTERIOR (loca).....	1170
R THORACIC POSTERIOR (q-pair).....	303
R THORACOTOMY (loca).....	1170
R THORACOTOMY (q-pair).....	303
R THUMB (loca).....	1170
R THUMB (q-pair).....	303
R TOE (loca).....	1170
R TOE (q-pair).....	303
R TRAPEZIUS (loca).....	1171
R TRAPEZIUS (q-pair).....	303
R TROCHANTER (loca).....	1171
R TROCHANTER (q-pair).....	547, 761
R UPPER ANTERIOR CHEST (loca).....	1171
R UPPER ANTERIOR CHEST (q-pair).....	110, 303
R UPPER ARM (loca).....	1171
R UPPER ARM (q-pair).....	102, 143, 219, 303, 344, 345, 473, 547, 658, 761
R UPPER EXTREMITY (loca).....	1171
R UPPER EXTREMITY (q-pair).....	411, 621
R UPPER LEG (loca).....	1172
R UPPER LEG (q-pair).....	761
R UPPER LIMB (loca).....	1172
R UPPER LIMB (q-pair).....	712
R UPPER LOBE (loca).....	1172
R UPPER LOBE (q-pair).....	74
R UPPER POSTERIOR CHEST (loca).....	1172
R UPPER POSTERIOR CHEST (q-pair).....	110, 303
R WRIST (loca).....	1172
R WRIST (q-pair).....	102, 143, 219, 344, 345, 598, 599, 761
RAD TECH (data).....	1024
RAD TECH (q-pair).....	554
RADIAL (loca).....	1173
RADIAL (q-pair).....	63, 570

RADIAL (qual).....	1224
RADIATING (data).....	1024
RADIATING (q-pair).....	470
RADIOLOGY (data).....	1024
RADIOLOGY (q-pair).....	553
RAMSEY ANXIOUS (data).....	1024
RAMSEY ANXIOUS (q-pair).....	609
RAMSEY ASLEEP BRISK RESPONSE (data).....	1025
RAMSEY ASLEEP BRISK RESPONSE (q-pair).....	609
RAMSEY ASLEEP NO RESPONSE (data).....	1025
RAMSEY ASLEEP NO RESPONSE (q-pair).....	609
RAMSEY ASLEEP SLUGGISH (data).....	1025
RAMSEY ASLEEP SLUGGISH (q-pair).....	609
RAMSEY COOPERATIVE (data).....	1025
RAMSEY COOPERATIVE (q-pair).....	609
RAMSEY SEDATED (data).....	1025
RAMSEY SEDATED (q-pair).....	609
RANGE (unit).....	1258
RASH/EXFOLIATION (data).....	1026
RASH/EXFOLIATION (q-pair).....	681, 766
RASS AGITATED (data).....	1026
RASS AGITATED (q-pair).....	610
RASS ALERT AND CALM (data).....	1026
RASS ALERT AND CALM (q-pair).....	610
RASS COMBATIVE (data).....	1026
RASS COMBATIVE (q-pair).....	610
RASS DEEP SEDATION (data).....	1026
RASS DEEP SEDATION (q-pair).....	610
RASS DROWSY (data).....	1027
RASS DROWSY (q-pair).....	610
RASS LIGHT SEDATION (data).....	1027
RASS LIGHT SEDATION (q-pair).....	610
RASS MODERATE SEDATION (data).....	1027
RASS MODERATE SEDATION (q-pair).....	610
RASS RESTLESS (data).....	1027
RASS RESTLESS (q-pair).....	610
RASS UNAROUSABLE (data).....	1027
RASS UNAROUSABLE (q-pair).....	610
RASS VERY AGITATED (data).....	1028
RASS VERY AGITATED (q-pair).....	610
REACTIVE (data).....	1028
REACTIVE (q-pair).....	572
RECTAL (loca).....	1173
RECTAL (q-pair).....	694, 761
RECTAL TEMP PROBE (meth).....	1194
RECTAL TUBE (data).....	1028
RECTAL TUBE PRESENT.....	576
RED (data).....	1028
RED (q-pair).....	46, 99, 143, 157, 173, 211, 238, 345, 447, 537, 563, 681, 714, 749, 753
RED/GRANULATING (data).....	1028
REDDENING OF THE SKIN (data).....	1029
REDDENING OF THE SKIN (q-pair).....	547
REDIRECT PATIENT ACTIVITY (data).....	1029
REDIRECT PATIENT ACTIVITY (q-pair).....	586
REFLEX (data).....	1029
REFLEX (q-pair).....	196, 222, 352, 412, 712
REFUSE MEDICATION (data).....	1029
REFUSE MEDICATION (q-pair).....	475

REFUSE TO MOVE (data).....	1029
REFUSE TO MOVE (q-pair).....	545
REFUSED (data).....	1030
REFUSED (q-pair).....	15
REGIONAL BLOCK (data).....	1030
REGIONAL BLOCK (q-pair).....	31
REGISTERED NURSE (data).....	1030
REGISTERED NURSE (q-pair).....	592
REGULAR (data).....	1030
REGULAR (q-pair).....	150, 577
REINFORCED (data).....	1030
REINFORCED (q-pair).....	43, 95
RELATIVE (data).....	1031
RELATIVE (q-pair).....	517
RELIGIOUS (data).....	1031
RELIGIOUS (q-pair).....	162
REMNANT P-WAVE (data).....	1031
REMNANT P-WAVE (q-pair).....	80
REMOVE FROM STIM. SETTING (data).....	1031
REMOVE FROM STIM. SETTING (q-pair).....	586
REMOVE RECTAL TUBE (data).....	1031
RENAL DIET (data).....	1032
RENAL DIET (q-pair).....	150
REORIENTED TO ENVIRONMENT (data).....	1032
REORIENTED TO ENVIRONMENT (q-pair).....	586
REPIRATORY DEPRESSION (data).....	1032
REPIRATORY DEPRESSION (q-pair).....	183
REPORTED (data).....	1032
REPORTED (q-pair).....	716
REPOSITION (data).....	1032
REPOSITION (q-pair).....	475, 586, 593, 622
REQUEST FOR MORE MEDICATION (data).....	1033
REQUEST FOR MORE MEDICATION (q-pair).....	475
RES CRIT: AGGRESSIVE (data).....	1033
RES CRIT: AGGRESSIVE (q-pair).....	579
RES CRIT: AMBULATE (data).....	1033
RES CRIT: AMBULATE (q-pair).....	579
RES CRIT: COMBATIVE (data).....	1033
RES CRIT: COMBATIVE (q-pair).....	579
RES CRIT: EQUIPMENT (data).....	1033
RES CRIT: EQUIPMENT (q-pair).....	579
RES CRIT: HARM TO SELF (data).....	1034
RES CRIT: HARM TO SELF (q-pair).....	579
RES CRIT: LINES (data).....	1034
RES CRIT: LINES (q-pair).....	579
RES CRIT: UNABLE TO FOLLOW DIRECTIONS (data).....	1034
RES CRIT: UNABLE TO FOLLOW DIRECTIONS (q-pair).....	579
RES EVAL: AMBULATION/EXERCISE (data).....	1034
RES EVAL: AMBULATION/EXERCISE (q-pair).....	582
RES EVAL: BED/WHEELCHAIR ALARM (posi).....	1202
RES EVAL: DECREASE STIMULATION (data).....	1034
RES EVAL: DECREASE STIMULATION (q-pair).....	582
RES EVAL: DIVERSIONAL ACTIVITIES (data).....	1035
RES EVAL: DIVERSIONAL ACTIVITIES (q-pair).....	582
RES EVAL: EDUCATION (data).....	1035
RES EVAL: EDUCATION (q-pair).....	582
RES EVAL: INVOLVEMENT WITH FAMILY (data).....	1035
RES EVAL: INVOLVEMENT WITH FAMILY (q-pair).....	582
RES EVAL: MEDICATION REVIEW (data).....	1035
RES EVAL: MEDICATION REVIEW (q-pair).....	582
RES EVAL: REALITY ORIENTATION (data).....	1035
RES EVAL: REALITY ORIENTATION (posi).....	1202

RES EVAL: REALITY ORIENTATION (q-pair).....	582
RES EVAL: RELOCATE (data).....	1036
RES EVAL: RELOCATE (q-pair).....	582
RES EVAL: UNABLE TO ATTEMPT (data).....	1036
RES EVAL: UNABLE TO ATTEMPT (q-pair).....	582
RES LRM: DEVICE ACCESS (data).....	1036
RES LRM: DEVICE ACCESS (q-pair).....	593
RES LRM: LIMIT SETTING (data).....	1036
RES LRM: LIMIT SETTING (q-pair).....	593
RES LRM: PRN MEDS (data).....	1036
RES LRM: PRN MEDS (q-pair).....	593
RES LRM: REDIRECT (data).....	1037
RES LRM: REDIRECT (q-pair).....	593
RES LRM: REMOVE FROM SETTING (data).....	1037
RES LRM: REMOVE FROM SETTING (q-pair).....	593
RES LRM: REORIENT (data).....	1037
RES LRM: REORIENT (q-pair).....	593
RES LRM: SUPERVISION (data).....	1037
RES LRM: SUPERVISION (q-pair).....	593
RES LRM: VERBAL DE-ESCALATION (data).....	1037
RES LRM: VERBAL DE-ESCALATION (q-pair).....	593
RES LRM: VERBAL REMINDERS (data).....	1038
RES LRM: VERBAL REMINDERS (q-pair).....	593
RESERVOIR CANNULA (meth).....	1194
RESERVOIR CANNULA (q-pair).....	687
RESOURCE BREEZE (data).....	1038
RESOURCE BREEZE (q-pair).....	172
RESOURCE DM (data).....	1038
RESOURCE DM (q-pair).....	172
RESP TECH (data).....	1038
RESP TECH (q-pair).....	554
RESPIRATION (data).....	1038
RESPIRATION (q-pair).....	402
RESPIRATORY DEPRESSION (data).....	1039
RESPIRATORY DEPRESSION (q-pair).....	498
RESPIRATORY PATTERN.....	577
RESPIRATORY RATE - DO NOT CHECK TOTAL OR SUBTOTALS.....	578
RESTING (data).....	1039
RESTING (q-pair).....	54
RESTING RESPIRATORY RATE BASELINE OR ABOVE (data).....	1039
RESTING RESPIRATORY RATE BASELINE OR ABOVE (q-pair)....	465
RESTLESS (data).....	1039
RESTLESS (q-pair).....	56, 351, 474, 594
RESTRAINT (data).....	1039
RESTRAINT (q-pair).....	595
RESTRAINT - ASSESS NEED FOR FLUIDS/TOILET (task).....	1240
RESTRAINT - CIRCULATION CHECK (task).....	1240
RESTRAINT - DEVICE RELEASE (task).....	1240
RESTRAINT - FALL PRECAUTIONS (task).....	1240
RESTRAINT - SEIZURE PRECAUTIONS (task).....	1241
RESTRAINT BILATERAL ANKLE (data).....	1040
RESTRAINT BILATERAL ANKLE (q-pair).....	598
RESTRAINT BILATERAL HAND MITT (data).....	1040
RESTRAINT BILATERAL HAND MITT (q-pair).....	598
RESTRAINT CRITERIA.....	579
RESTRAINT EDUCATION.....	580
RESTRAINT EDUCATION TYPE.....	581
RESTRAINT EVAL - ALT. MEASURES ATTEMPTED.....	582
RESTRAINT EVAL - ELIMINATION MONITORED.....	583
RESTRAINT EVAL - FLUID/NUTRITION INTAKE.....	584
RESTRAINT EVAL - HYGIENE/COMFORT.....	585
RESTRAINT EVAL - LEAST RESTRICTIVE MEASURES.....	586

RESTRAINT EVAL - MENTAL.....	587
RESTRAINT EVAL - PRIVACY/DIGNITY.....	588
RESTRAINT EVAL - REPOSITION/RANGE OF MOTION.....	589
RESTRAINT EVAL - SKIN INTEGRETU.....	590
RESTRAINT FAMILY ED. UNDERSTANDING.....	591
RESTRAINT GERI CHAIR (data).....	1040
RESTRAINT GERI CHAIR (q-pair).....	598
RESTRAINT HARD PLASTIC/LEATHER (data).....	1040
RESTRAINT HARD PLASTIC/LEATHER (q-pair).....	598, 599, 600
RESTRAINT HARD PLASTIC/LEATHER 2 POINT (data).....	1040
RESTRAINT HARD PLASTIC/LEATHER 2 POINT (q-pair).....	598
RESTRAINT HARD PLASTIC/LEATHER 4 POINT (data).....	1041
RESTRAINT HARD PLASTIC/LEATHER 4 POINT (q-pair).....	598
RESTRAINT INITIATED BY.....	592
RESTRAINT LEAST RESTRICTIVE MEASURES.....	593
RESTRAINT MENTAL ASSESSMENT.....	594
RESTRAINT ONE ARM (data).....	1041
RESTRAINT ONE ARM (q-pair).....	598
RESTRAINT ONE LEG (data).....	1041
RESTRAINT ONE LEG (q-pair).....	598
RESTRAINT ORDER TYPE.....	595
RESTRAINT PATIENT ED. UNDERSTANDING.....	596
RESTRAINT SAFETY CARE.....	597
RESTRAINT SINGLE HAND MITT (data).....	1041
RESTRAINT SINGLE HAND MITT (q-pair).....	598, 599
RESTRAINT SOFT (data).....	1041
RESTRAINT SOFT (q-pair).....	598, 599, 600
RESTRAINT SOFT 2 POINT (data).....	1042
RESTRAINT SOFT 2 POINT (q-pair).....	598
RESTRAINT SOFT 4 POINT (data).....	1042
RESTRAINT SOFT 4 POINT (q-pair).....	598
RESTRAINT SOFT ARMS + LOCKED WAIST (data).....	1042
RESTRAINT SOFT ARMS + LOCKED WAIST (q-pair).....	598
RESTRAINT SOFT ARMS + WAIST (data).....	1042
RESTRAINT SOFT ARMS + WAIST (q-pair).....	598
RESTRAINT STRAP (data).....	1042
RESTRAINT STRAP (q-pair).....	598, 600
RESTRAINT STRAP (sche).....	1227
RESTRAINT TYPE ALL.....	598
RESTRAINT TYPE ARMS.....	599
RESTRAINT TYPE LEGS.....	600
RESTRAINT TYPE OTHER.....	601
RESTRAINT VEST (data).....	1043
RESTRAINT VEST (q-pair).....	598, 601
RESTRAINT WAIST (data).....	1043
RESTRAINT WAIST (q-pair).....	598, 601
RESTRAINT WAIST LOCKED (data).....	1043
RESTRAINT WAIST LOCKED (q-pair).....	598
RESTRAINT WRIST (data).....	1043
RESTRAINT WRIST (q-pair).....	598
RESTRAINT- ACTIVE ORDER.....	602
RESTRAINTS BROCHURE (data).....	1043
RESTRAINTS BROCHURE (q-pair).....	581
RESTRAINTS ED REMOVAL (data).....	1044
RESTRAINTS ED REMOVAL (q-pair).....	581
RESTRAINTS ED TYPE (data).....	1044
RESTRAINTS ED TYPE (q-pair).....	581
RESTRAINTS EXPLAINED (data).....	1044
RESTRAINTS EXPLAINED (q-pair).....	581
RESTRAINTS LIMIT SETTING (data).....	1044
RESTRAINTS LIMIT SETTING (q-pair).....	586
RESTRAINTS TIME LIMIT (data).....	1044

RESTRAINTS TIME LIMIT (q-pair).....	581
RESTRAINTS: BEHAVIOR CRITERIA FOR DISCONTINUE.....	603
RHONCHI (data).....	1045
RHONCHI (q-pair).....	74
RHYTHM (meth).....	1194
RIGHT (data).....	1045
RIGHT (loca).....	1173
RIGHT (q-pair).....	13
RIGHT EYE (loca).....	1173
RIGHT EYE (q-pair).....	152
RIGHT VENTRICULAR STROKE WORK.....	604
RIGHT VENTRICULAR STROKE WORK INDEX.....	605
RIGID HEEL SPLINT (data).....	1045
RIGID HEEL SPLINT (q-pair).....	544
RN (data).....	1045
RN (q-pair).....	554, 608
ROOM AIR (meth).....	1195
ROOM AIR15 (meth).....	1195
ROOM AIR20 (meth).....	1195
ROPIVICAINE (data).....	1045
ROPIVICAINE (q-pair).....	181
ROTATIONAL THERAPY (data).....	1046
ROTATIONAL THERAPY (q-pair).....	665
ROUND (data).....	1046
ROUND (q-pair).....	2
RT (data).....	1046
RT (q-pair).....	608
RUB (data).....	1046
RUB (q-pair).....	81
RUNS OF VT (data).....	1046
RUNS OF VT (q-pair).....	80
Rh IMMUNE GLOBULIN (prod).....	1219
Rh IMMUNE GLOBULIN (q-pair).....	317

- S -

S1 (data).....	1047
S1 (q-pair).....	81
S2 (data).....	1047
S2 (q-pair).....	81
S3 (data).....	1047
S3 (q-pair).....	81
S4 (data).....	1047
S4 (q-pair).....	81
SACRAL (loca).....	1173
SACRAL (q-pair).....	160, 547, 761
SAFETY (data).....	1047
SAFETY (q-pair).....	164
SAFETY DEVICE WAIST BELT APPLIED.....	606
SALINE (prod).....	1219
SALINE (q-pair).....	316, 325
SALINE LOCK (data).....	1048
SALINE LOCK (q-pair).....	343, 493
SALINE MOIST DRESSING (data).....	1048
SALINE MOIST DRESSING (q-pair).....	548, 765
SANDOSOURCE PEPTIDE (data).....	1048
SANDOSOURCE PEPTIDE (prod).....	1219
SANDOSOURCE PEPTIDE (q-pair).....	309
SAS AGITATED (data).....	1048
SAS AGITATED (q-pair).....	611
SAS CALM AND COOPERATIVE (data).....	1048
SAS CALM AND COOPERATIVE (q-pair).....	611

SAS DANGEROUS AGITATION (data).....	1049
SAS DANGEROUS AGITATION (q-pair).....	611
SAS SEDATED (data).....	1049
SAS SEDATED (q-pair).....	611
SAS UNAROUSABLE (data).....	1049
SAS UNAROUSABLE (q-pair).....	611
SAS VERY AGITATED (data).....	1049
SAS VERY AGITATED (q-pair).....	611
SAS VERY SEDATED (data).....	1049
SAS VERY SEDATED (q-pair).....	611
SCALE (unit).....	1258
SCANT (data).....	1050
SCANT (q-pair).....	428, 562, 647
SCAPHOID (data).....	1050
SCAPHOID (q-pair).....	2
SCOOTER (data).....	1050
SCOOTER (q-pair).....	199
SCROTUM (loca).....	1174
SCROTUM (q-pair).....	473, 547, 658
SECLUSION (data).....	1050
SECLUSION (q-pair).....	595
SECLUSION: BEHAVIOR CRITERIA FOR DISCONTINUATION.....	607
SECOND (unit).....	1258
SECOND DEGREE AV BLOCK: TYPE 1 (data).....	1050
SECOND DEGREE AV BLOCK: TYPE 1 (q-pair).....	84
SECOND DEGREE AV BLOCK: TYPE 2 (data).....	1051
SECOND DEGREE AV BLOCK: TYPE 2 (q-pair).....	84
SECOND VERIFIER.....	608
SECURE TRACK IN PLACE (data).....	1051
SECURE TRACK IN PLACE (q-pair).....	4
SEDATED (data).....	1051
SEDATED (q-pair).....	56, 161, 490, 594
SEDATION SCALE - RAMSEY.....	609
SEDATION SCALE - RASS.....	610
SEDATION SCALE - SAS.....	611
SEDIMENT (data).....	1051
SEDIMENT (q-pair).....	713
SEE CPRS PROGRESS NOTE FOR ADDITIONAL DATA AND ASSESSMENTS .....	612
SEIZURE CHARACTER.....	613
SEIZURE DURATION.....	614
SEIZURE FOCI.....	615
SEIZURE PRECAUTIONS.....	616
SEIZURES.....	617
SEIZURES (data).....	1051
SEIZURES (q-pair).....	396
SELECTED SURGERY (data).....	1052
SELECTED SURGERY (q-pair).....	90
SELF AM CARE BY PATIENT (data).....	1052
SELF AM CARE BY PATIENT (q-pair).....	5
SELF CARE (data).....	1052
SELF CARE (q-pair).....	13, 15
SELF CARE OF BEDSIDE ENVIRONMENT (data).....	1052
SELF CARE OF BEDSIDE ENVIRONMENT (q-pair).....	5
SELF FACE WASH (data).....	1052
SELF FACE WASH (q-pair).....	5
SELF FOOT CARE (data).....	1053
SELF FOOT CARE (q-pair).....	5
SELF SHAVE (data).....	1053
SELF SHAVE (q-pair).....	5
SELF SHOWER (data).....	1053
SELF SHOWER (q-pair).....	5



SELF TRACH CARE (data).....	1053
SELF TRACH CARE (q-pair).....	5
SELF-CARE SKILLS (data).....	1053
SELF-CARE SKILLS (q-pair).....	164
SEMIFOWLERS (posi).....	1202
SEMIFOWLERS (q-pair).....	63, 570, 578
SEMIFOWLERS HIGH (data).....	1054
SEMIFOWLERS HIGH (q-pair).....	53
SEMIFOWLERS LOW (data).....	1054
SEMIFOWLERS LOW (q-pair).....	53
SENSATION FOOT.....	618
SENSATION HAND.....	619
SENSATION LOWER EXTREMITY.....	620
SENSATION UPPER EXTREMITY.....	621
SEQUENTIAL COMPRESSION DEVICE (data).....	1054
SEQUENTIAL COMPRESSION DEVICE (q-pair).....	158
SEQUENTIAL COMPRESSION DEVICE APPLICATION.....	622
SEROSANGUINOUS (data).....	1054
SEROSANGUINOUS (q-pair).....	109, 296, 539, 636
SEROUS (data).....	1054
SEROUS (q-pair).....	109, 296, 539, 636
SET UP HELP ONLY (data).....	1055
SET UP HELP ONLY (q-pair).....	10
SEVERELY REDUCED (data).....	1055
SEVERELY REDUCED (q-pair).....	348
SEWS - AGITATION.....	623
SEWS - ANXIETY.....	624
SEWS - DISORIENTATION NAME PLACE.....	625
SEWS - HALLUCINATION ANY.....	626
SEWS - HALLUCINATION AUDITORY ONLY.....	627
SEWS - NAME ONLY.....	628
SEWS - NAUSEA.....	629
SEWS - SWEATING.....	630
SEWS - TOTAL.....	631
SEWS - TREMOR.....	632
SEWS - VITAL SIGNS.....	633
SEWS SCALE (unit).....	1258
SHALLOW (data).....	1055
SHALLOW (q-pair).....	577
SHAMPOO (data).....	1055
SHAMPOO (q-pair).....	5
SHARP (data).....	1055
SHARP (q-pair).....	470
SHAVE (data).....	1056
SHAVE (q-pair).....	5
SHOOTING (data).....	1056
SHOOTING (q-pair).....	470
SHOWER (data).....	1056
SHOWER (q-pair).....	5
SHUFFLING (data).....	1056
SHUFFLING (q-pair).....	199
SHUNTS (data).....	1056
SHUNTS (q-pair).....	396
SI UNIT (unit).....	1258
SICU (data).....	1057
SICU (q-pair).....	467
SIDE RAILS.....	634
SIGMOIDOSTOMY (data).....	1057
SIGMOIDOSTOMY (q-pair).....	424
SIGNIFICANT EVENT.....	635
SIGNIFICANT OTHER (data).....	1057
SIGNIFICANT OTHER (q-pair).....	11

SILVADENE (data).....	1057
SILVADENE (q-pair).....	548
SILVER (data).....	1057
SILVER SULFADIAZINE (data).....	1058
SINGLE (data).....	1058
SINGLE (q-pair).....	102
SINGLE LUMEN CATHETER (meth).....	1195
SINGLE USE LUBRICANT PACKET (data).....	1058
SINGLE USE LUBRICANT PACKET (q-pair).....	91
SINUS ARREST (data).....	1058
SINUS ARREST (q-pair).....	84
SINUS ARRHYTHMIA (data).....	1058
SINUS ARRHYTHMIA (q-pair).....	84
SINUS BLOCK (data).....	1059
SINUS BLOCK (q-pair).....	84
SINUS BRADYCARDIA (data).....	1059
SINUS BRADYCARDIA (q-pair).....	84
SINUS DRAINAGE CHARACTER.....	636
SINUS DRAINAGE PRESENT.....	637
SINUS NECROTIC TISSUE APPEARANCE.....	638
SINUS NECROTIC TISSUE PRESENT.....	639
SINUS PAUSE (data).....	1059
SINUS PAUSE (q-pair).....	84
SINUS RHYTHM (data).....	1059
SINUS RHYTHM (q-pair).....	84
SINUS TACHYCARDIA (data).....	1059
SINUS TACHYCARDIA (q-pair).....	84
SINUS TRACT SURROUNDING TISSUE CHARACTER.....	640
SINUS TRACT ULCER.....	641
SINUS TRACT ULCER EPITHELIALIZATION.....	642
SINUS TRACT ULCER GRANULIZATION.....	643
SITE UNABLE TO BE MARKED (data).....	1060
SITE UNABLE TO BE MARKED (q-pair).....	531
SITTING (posi).....	1203
SITTING (q-pair).....	63, 149, 570, 578, 690, 746
SKIN (ADHESIVE INDICATOR) (meth).....	1195
SKIN (ADHESIVE INDICATOR) (q-pair).....	165
SKIN (loca).....	1174
SKIN (q-pair).....	694
SKIN ALTERATION.....	644
SKIN COLOR.....	645
SKIN CONDITION.....	646
SKIN DRAINAGE AMOUNT.....	647
SKIN ERYTHEMA EXTENDS AROUND INCISION.....	648
SKIN ERYTHEMA EXTENSION AROUND WOUND.....	649
SKIN INCISION CLOSURE.....	650
SKIN INCISION WOUND CLOSURE OPEN WOUND (data).....	1060
SKIN INCISION WOUND CLOSURE OPEN WOUND (q-pair).....	650
SKIN INCISION WOUND CLOSURE STAPLES (data).....	1060
SKIN INCISION WOUND CLOSURE STAPLES (q-pair).....	650
SKIN INCISION WOUND CLOSURE STERI-STRIPS (data).....	1060
SKIN INCISION WOUND CLOSURE STERI-STRIPS (q-pair).....	650
SKIN INCISION WOUND CLOSURE SUTURES (data).....	1060
SKIN INCISION WOUND CLOSURE SUTURES (q-pair).....	650
SKIN INDURATION AROUND INCISION EXTENDS INTO PERIWOUND LESS THAN 1 cm (data ).....	1061
SKIN INDURATION AROUND INCISION EXTENDS INTO PERIWOUND LESS THAN 1 cm (q-pa ir).....	651
SKIN INDURATION AROUND INCISION EXTENDS INTO PERIWOUND MORE THAN 1 cm (data ).....	1061
SKIN INDURATION AROUND INCISION EXTENDS INTO PERIWOUND MORE THAN 1 cm (q-pa	

ir).....	651
SKIN INDURATION AROUND INCISION LIMITED TO INCISION (data)	.....1061
SKIN INDURATION AROUND INCISION LIMITED TO INCISION (q-pair)	.....651
SKIN INDURATION AROUND INCISION NONE (data).....	1061
SKIN INDURATION AROUND INCISION NONE (q-pair).....	651
SKIN INDURATION EXTENSION AROUND INCISION.....	651
SKIN INDURATION EXTENSION AROUND WOUND.....	652
SKIN INDURATION EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND < 1 cm (data)	.....1061
SKIN INDURATION EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND < 1 cm (q-pai	r).....
SKIN INDURATION EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND > 1 cm (data)	.....1062
SKIN INDURATION EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND > 1 cm (q-pai	r).....
SKIN INDURATION EXTENSION AROUND WOUND LIMITED TO WOUND (data)	.....1062
SKIN INDURATION EXTENSION AROUND WOUND LIMITED TO WOUND (q-pair)	.....652
SKIN INDURATION EXTENSION AROUND WOUND NONE (data).....	1062
SKIN INDURATION EXTENSION AROUND WOUND NONE (q-pair).....	652
SKIN MOISTURE.....	653
SKIN PATCHES PRESENT.....	654
SKIN RISK FACTORS PRESENT.....	655
SKIN SWELLING EXTENSION AROUND INCISION.....	656
SKIN SWELLING EXTENSION AROUND INCISION EXTENDS INTO PERIWOUND < 1 cm (data	).....
SKIN SWELLING EXTENSION AROUND INCISION EXTENDS INTO PERIWOUND < 1 cm (q-pa	ir).....
SKIN SWELLING EXTENSION AROUND INCISION EXTENDS INTO PERIWOUND > 1 cm (data	).....
SKIN SWELLING EXTENSION AROUND INCISION EXTENDS INTO PERIWOUND > 1 cm (q-pa	ir).....
SKIN SWELLING EXTENSION AROUND INCISION FLUCTUANT (data)	.....1063
SKIN SWELLING EXTENSION AROUND INCISION FLUCTUANT (q-pair)	.....656
SKIN SWELLING EXTENSION AROUND INCISION LIMITED TO WOUND (data)	.....1063
SKIN SWELLING EXTENSION AROUND INCISION LIMITED TO WOUND (q-pair)	.....656
SKIN SWELLING EXTENSION AROUND INCISION NONE (data).....	1063
SKIN SWELLING EXTENSION AROUND INCISION NONE (q-pair).....	656
SKIN SWELLING EXTENSION AROUND WOUND.....	657
SKIN SWELLING EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND < 1 cm (data)	.....1063
SKIN SWELLING EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND < 1 cm (q-pair)	.....657
SKIN SWELLING EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND > 1 cm (data)	.....1063
SKIN SWELLING EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND > 1 cm (q-pair)	.....657
SKIN SWELLING EXTENSION AROUND WOUND FLUCTANT (data).....	1064
SKIN SWELLING EXTENSION AROUND WOUND FLUCTANT (q-pair).....	657
SKIN SWELLING EXTENSION AROUND WOUND LIMITED TO WOUND (data)	.....1064
SKIN SWELLING EXTENSION AROUND WOUND LIMITED TO WOUND (q-pair)	.....657
SKIN SWELLING EXTENSION AROUND WOUND NONE (data).....	1064
SKIN SWELLING EXTENSION AROUND WOUND NONE (q-pair).....	657

SKIN TEAR (data).....	1064
SKIN TEAR (q-pair).....	766
SKIN TEMP.....	658
SKIN TEMPERATURE AT INCISION SITE.....	659
SKIN TEMPERATURE AT INCISION SITE CONSISTENT WITH BODY TEMPERATURE (data) .....	1064
SKIN TEMPERATURE AT INCISION SITE CONSISTENT WITH BODY TEMPERATURE (q-pair) .....	659
SKIN TEMPERATURE AT INCISION SITE COOLER THAN GENERAL BODY TEMPERATURE (dat a).....	1065
SKIN TEMPERATURE AT INCISION SITE COOLER THAN GENERAL BODY TEMPERATURE (q-p air).....	659
SKIN TEMPERATURE AT INCISION SITE COOLER THAN OPPOSITE BODY PART (data) .....	1065
SKIN TEMPERATURE AT INCISION SITE COOLER THAN OPPOSITE BODY PART (q-pair) .....	659
SKIN TEMPERATURE AT INCISION SITE WARMER THAN GENERAL BODY TEMPERATURE (dat a).....	1065
SKIN TEMPERATURE AT INCISION SITE WARMER THAN GENERAL BODY TEMPERATURE (q-p air).....	659
SKIN TEMPERATURE AT INCISION SITE WARMER THAN OPPOSITE BODY PART (data) .....	1065
SKIN TEMPERATURE AT INCISION SITE WARMER THAN OPPOSITE BODY PART (q-pair) .....	659
SKIN TEMPERATURE AT WOUND SITE.....	660
SKIN TEMPERATURE AT WOUND SITE CONSISTENT WITH BODY TEMPERATURE (data) .....	1065
SKIN TEMPERATURE AT WOUND SITE CONSISTENT WITH BODY TEMPERATURE (q-pair) .....	660
SKIN TEMPERATURE AT WOUND SITE COOLER THAN GENERAL BODY TEMPERATURE (data) .....	1066
SKIN TEMPERATURE AT WOUND SITE COOLER THAN GENERAL BODY TEMPERATURE (q-pair ).....	660
SKIN TEMPERATURE AT WOUND SITE COOLER THAN OPPOSITE BODY PART (data) .....	1066
SKIN TEMPERATURE AT WOUND SITE COOLER THAN OPPOSITE BODY PART (q-pair) .....	660
SKIN TEMPERATURE AT WOUND SITE WARMER THAN GENERAL BODY TEMPERATURE (data) .....	1066
SKIN TEMPERATURE AT WOUND SITE WARMER THAN GENERAL BODY TEMPERATURE (q-pair ).....	660
SKIN TEMPERATURE AT WOUND SITE WARMER THAN OPPOSITE BODY PART (data) .....	1066
SKIN TEMPERATURE AT WOUND SITE WARMER THAN OPPOSITE BODY PART (q-pair) .....	660
SKIN TRAUMA (data).....	1066
SKIN TRAUMA (q-pair).....	766
SKIN TURGOR.....	661
SKIN WARM AND DRY COLOR APPROPRIATE (data).....	1067
SKIN WARM AND DRY COLOR APPROPRIATE (q-pair).....	465
SKIN WOUND CLOSURE.....	662
SKIN WOUND CLOSURE APPROXIMATED (data).....	1067
SKIN WOUND CLOSURE APPROXIMATED (q-pair).....	662
SKIN WOUND CLOSURE DEHISCENCE (data).....	1067
SKIN WOUND CLOSURE DEHISCENCE (q-pair).....	662
SKIN WOUND CLOSURE GLUE (data).....	1067
SKIN WOUND CLOSURE GLUE (q-pair).....	662
SKIN WOUND CLOSURE NOT CLOSED (data).....	1067
SKIN WOUND CLOSURE NOT CLOSED (q-pair).....	662
SKIN WOUND CLOSURE OPEN WOUND (data).....	1068
SKIN WOUND CLOSURE OPEN WOUND (q-pair).....	662
SKIN WOUND CLOSURE STAPLES (data).....	1068

SKIN WOUND CLOSURE STAPLES (q-pair).....	662
SKIN WOUND CLOSURE STERI-STRIPS (data).....	1068
SKIN WOUND CLOSURE STERI-STRIPS (q-pair).....	662
SKIN WOUND CLOSURE SUTURES (data).....	1068
SKIN WOUND CLOSURE SUTURES (q-pair).....	662
SLEEP AMOUNT.....	663
SLEEP QUALIFIERS.....	664
SLEEPING (data).....	1068
SLEEPING (q-pair).....	128, 474, 594
SLEEPS POORLY (data).....	1069
SLEEPS POORLY (q-pair).....	56
SLEPT WELL (data).....	1069
SLEPT WELL (q-pair).....	664
SLIGHT (data).....	1069
SLIGHT (q-pair).....	437
SLING TRACTION (data).....	1069
SLING TRACTION (q-pair).....	702
SLOUGH (data).....	1069
SLOUGH (q-pair).....	537, 760
SLUGGISH (data).....	1070
SLUGGISH (q-pair).....	572
SLURRED SPEECH (data).....	1070
SLURRED SPEECH (q-pair).....	183
SM ADULT (meth).....	1196
SM ADULT CUFF (meth).....	1196
SMALL (data).....	1070
SMALL (q-pair).....	428, 435, 437, 562,
	647
SMALL BORE FEEDING TUBE (data).....	1070
SMALL BORE FEEDING TUBE (q-pair).....	215
SMEAR (data).....	1070
SMEAR (q-pair).....	435
SMOKING AREA (data).....	1071
SMOKING AREA (q-pair).....	55, 129
SMOKING PORCH (data).....	1071
SMOKING PORCH (q-pair).....	55, 129
SMOOTH (data).....	1071
SMOOTH (q-pair).....	760
SNACK (data).....	1071
SNACK (q-pair).....	5
SNORES (data).....	1071
SNORES (q-pair).....	396
SODA: DIET (prod).....	1219
SODA: DIET (q-pair).....	329
SODA: NON-DIET (prod).....	1219
SODA: NON-DIET (q-pair).....	329
SODIUM CONTROLLED 2 GM (data).....	1072
SODIUM CONTROLLED 2 GM (q-pair).....	150
SODIUM CONTROLLED 4 GM (data).....	1072
SODIUM CONTROLLED 4 GM (q-pair).....	150
SOFT (data).....	1072
SOFT (q-pair).....	2, 638, 680
SOFT DRINK (prod).....	1220
SOMNOLENT (data).....	1072
SOMNOLENT (q-pair).....	498
SPECIAL PROCEDURE ROOM (data).....	1072
SPECIAL PROCEDURE ROOM (q-pair).....	553
SPECIALTY BED/MATTRESS.....	665
SPECIALTY MATTRESS/BED (data).....	1073
SPECIALTY MATTRESS/BED (q-pair).....	544
SPECIFIC GRAVITY.....	666
SPINAL (data).....	1073

SPINAL (q-pair).....	31
SPINE CERVICAL (loca).....	1174
SPINE CERVICAL (q-pair).....	303
SPINE LUMBAR (loca).....	1174
SPINE LUMBAR (q-pair).....	303
SPINE THORACIC (loca).....	1174
SPINE THORACIC (q-pair).....	303
SPIRITUAL (data).....	1073
SPIRITUAL (q-pair).....	162, 164
SPLINTING (data).....	1073
SPLINTING (q-pair).....	474
SPONTANEOUS (data).....	1073
SPONTANEOUS (meth).....	1196
SPONTANEOUS (q-pair).....	578, 725
SPOONING (data).....	1074
SPOONING (q-pair).....	415
SPOUSE (data).....	1074
SPOUSE (q-pair).....	515, 516, 517
SQ EMPHYSEMA (data).....	1074
SQ EMPHYSEMA (q-pair).....	113
SQUARE (data).....	1074
SQUARE (q-pair).....	740
SQUARE METERS (unit).....	1259
SQUEEZING (data).....	1074
SQUEEZING (q-pair).....	470
ST SEGMENT AVF.....	667
ST SEGMENT AVL.....	668
ST SEGMENT AVR.....	669
ST SEGMENT I.....	670
ST SEGMENT II.....	671
ST SEGMENT III.....	672
ST SEGMENT V1.....	673
ST SEGMENT V2.....	674
ST SEGMENT V3.....	675
ST SEGMENT V4.....	676
ST SEGMENT V5.....	677
ST SEGMENT V6.....	678
ST/T WAVE MEASUREMENT.....	679
STABBING (data).....	1075
STABBING (q-pair).....	470
STABLE (data).....	1075
STABLE (q-pair).....	186
STABLE BP AND PULSE (data).....	1075
STABLE BP AND PULSE (q-pair).....	465
STAFF (data).....	1075
STAFF (q-pair).....	706
STAGE 3 OR 4 PRESSURE ULCER (data).....	1075
STAGE 3 OR 4 PRESSURE ULCER (q-pair).....	90
STAGE I (data).....	1076
STAGE I (q-pair).....	547
STAGE II (data).....	1076
STAGE II (q-pair).....	547
STAGE III (data).....	1076
STAGE III (q-pair).....	547
STAGE IV (data).....	1076
STAGE IV (q-pair).....	547
STANDARD (data).....	1076
STANDARD (q-pair).....	333
STANDARD PRESSURE REDUCTION (data).....	1077
STANDARD PRESSURE REDUCTION (q-pair).....	665
STANDING (posi).....	1203
STANDING (q-pair).....	63, 149, 570, 578,

STANDING WEIGHT (meth).....	1196
STANDING WEIGHT (q-pair).....	746
STAPLES (data).....	1077
START TIME: (data).....	1077
START TIME: (q-pair).....	555
STARTED (data).....	1077
STARTED (q-pair).....	343, 493
STARTED ON LEFT (data).....	1077
STARTED ON LEFT (q-pair).....	615
STARTED ON RIGHT (data).....	1078
STARTED ON RIGHT (q-pair).....	615
STATED (q-pair).....	227, 716, 746, 748
STATED (qual).....	1224
STEADY (data).....	1078
STEADY (q-pair).....	199
STEADY GATE OR BASELINE (data).....	1078
STEADY GATE OR BASELINE (q-pair).....	151
STEP DOWN (data).....	1078
STEP DOWN (q-pair).....	467
STERILE GLOVES-DRAPES USED (data).....	1078
STERILE GLOVES-DRAPES USED (q-pair).....	91
STERNOTOMY (loca).....	1175
STINGING (data).....	1079
STINGING (q-pair).....	470
STOMACH (loca).....	1175
STOMACH (q-pair).....	473
STOOL CHARACTERISTICS.....	680
STOOL COLOR.....	681
STOOL NUMBER STOOLS.....	682
STOP TIME: (data).....	1079
STOP TIME: (q-pair).....	555
STOPCOCK CHANGE (data).....	1079
STOPCOCK CHANGE (q-pair).....	170
STRAIGHT CATHETER (data).....	1079
STRAIGHT CATHETER (q-pair).....	219, 716
STRETCHER (data).....	1079
STRETCHER (q-pair).....	707
STRICT (data).....	1080
STRICT (q-pair).....	333
STROKE (data).....	1080
STROKE (q-pair).....	396
STROKE INDEX.....	683
STROKE VOLUME.....	684
STROKE VOLUME INDEX.....	685
STROKE WORK PER BEAT (unit).....	1259
STROKE WORK PER BEAT PER M2 (unit).....	1259
STRONG (data).....	1080
STRONG (q-pair).....	196, 222, 352, 712
STRYKER FRAME (data).....	1080
STRYKER FRAME (q-pair).....	665
STUDENT (data).....	1080
STUDENT (q-pair).....	554, 706
STUPOROUS (data).....	1081
STUPOROUS (q-pair).....	490
SUBCUTANEOUS TISSUE (data).....	1081
SUBSTERNAL (loca).....	1175
SUBSTERNAL (q-pair).....	473
SUCTIONED ONLY (data).....	1081
SUCTIONED ONLY (q-pair).....	15
SUICIDAL (data).....	1081
SUICIDAL (q-pair).....	58, 130

SUMP (data).....	1081
SUMP (q-pair).....	755
SUPERVISION/SET-UP HELP (data).....	1082
SUPERVISION/SET-UP HELP (q-pair).....	8, 9, 10
SUPINE (data).....	1082
SUPINE (posi).....	1203
SUPINE (q-pair).....	13, 63, 578
SUPLENA (data).....	1082
SUPLENA (q-pair).....	172
SUPPLEMENTAL OXYGEN CONCENTRATION.....	686
SUPPLEMENTAL OXYGEN FLOW RATE.....	687
SUPPORT GROUPS (data).....	1082
SUPPORT GROUPS (q-pair).....	164
SUPRAPUBIC (data).....	1082
SUPRAPUBIC (q-pair).....	219, 716
SUPRAVENTRICULAR TACHYCARDIA (data).....	1083
SUPRAVENTRICULAR TACHYCARDIA (q-pair).....	84
SURGICAL (data).....	1083
SURGICAL (q-pair).....	766
SURGICAL DONOR (data).....	1083
SURGICAL DONOR (q-pair).....	766
SURGICAL GRAFT (data).....	1083
SURGICAL GRAFT (q-pair).....	766
SUSPECTED DEEP TISSUE INJURY (data).....	1083
SUSPECTED DEEP TISSUE INJURY (q-pair).....	547
SUTURES (data).....	1084
SWOLLEN (data).....	1084
SWOLLEN (q-pair).....	46, 99, 143, 238, 345, 447
SYMPTOMS OF MEDICATION ALLERGY (data).....	1084
SYMPTOMS OF MEDICATION ALLERGY (q-pair).....	183
SYNCHRONOUS INTERMITTENT MANDATORY VENTILATION (data)..	1084
SYNCHRONOUS INTERMITTENT MANDATORY VENTILATION (q-pair)	725
SYRINGE (data).....	1084
SYRINGE (q-pair).....	167
SYRINGE CHANGE (data).....	1085
SYRINGE CHANGE (q-pair).....	170
SYRINGE CHANGED (data).....	1085
SYRINGE CHANGED (q-pair).....	493
SYSTEMIC VASCULAR RESISTANCE.....	688
SYSTEMIC VASCULAR RESISTANCE INDEX.....	689
SYSTOLIC (data).....	1085
SYSTOLIC (q-pair).....	87
SYSTOLIC PRESSURE.....	690
SaO2%.....	691
SaO2% AT ACCEPTABLE LEVEL BY ORDERS (data).....	1085
SaO2% AT ACCEPTABLE LEVEL BY ORDERS (q-pair).....	465
SaO2% AT ACCEPTABLE LEVEL WITH PRESCRIBED O2 LEVEL (data) .....	1085
SaO2% AT ACCEPTABLE LEVEL WITH PRESCRIBED O2 LEVEL (q-pair) .....	465
SpO2%.....	692
SvO2%.....	693

- T -

T-PIECE (meth).....	1196
T-PIECE (q-pair).....	686
T1 (data).....	1086
T1 (q-pair).....	179, 180
T10 (data).....	1086
T10 (q-pair).....	179, 180



T11 (data).....	1086
T11 (q-pair).....	179, 180
T12 (data).....	1086
T12 (q-pair).....	179, 180
T2 (data).....	1086
T2 (q-pair).....	179, 180
T3 (data).....	1087
T3 (q-pair).....	179, 180
T4 (data).....	1087
T4 (q-pair).....	179, 180
T5 (data).....	1087
T5 (q-pair).....	179, 180
T6 (data).....	1087
T6 (q-pair).....	179, 180
T7 (data).....	1087
T7 (q-pair).....	179, 180
T8 (data).....	1088
T8 (q-pair).....	179, 180
T9 (data).....	1088
T9 (q-pair).....	179, 180
TAN (data).....	1088
TAN (q-pair).....	157, 173, 211, 563, 714, 749, 753
TAPE BURN (data).....	1088
TAPE BURN (q-pair).....	766
TARRY (data).....	1088
TARRY (q-pair).....	680
TEA (data).....	1089
TEA (q-pair).....	714
TECH (data).....	1089
TECH (q-pair).....	554
TEETH (loca).....	1175
TEETH (q-pair).....	473
TEETH BRUSHED (data).....	1089
TEETH BRUSHED (q-pair).....	15
TEMPERATURE.....	694
TEMPORAL (loca).....	1175
TEMPORAL (q-pair).....	694
TEMPORARY (TRANSVENOUS) (data).....	1089
TEMPORARY (TRANSVENOUS) (q-pair).....	461
TENACIOUS (data).....	1089
TENACIOUS (q-pair).....	564
TENDER (data).....	1090
TENDER (q-pair).....	2, 470
TENDON (data).....	1090
TESTS/PROC. - ACCUCHECKS (task).....	1241
TESTS/PROC. - BLOOD CULTURES (task).....	1241
TESTS/PROC. - DAILY LABS (task).....	1241
TESTS/PROC. - MORNING LABS (task).....	1241
TESTS/PROC. - ZERO BALANCE HEMO. LINES (task).....	1242
THICK (data).....	1090
THICK (q-pair).....	157, 210, 416, 564
THIGH (data).....	1090
THIGH (loca).....	1176
THIGH (q-pair).....	64, 116
THIGH CUFF (meth).....	1197
THIN (data).....	1090
THIN (q-pair).....	157, 210, 564
THIRD DEGREE AV BLOCK (data).....	1091
THIRD DEGREE AV BLOCK (q-pair).....	84
THORACIC SPINAL PRECAUTIONS (data).....	1091
THORACIC SPINAL PRECAUTIONS (q-pair).....	12

THREADS (data).....	1091
THREADS (q-pair).....	713
THROBBING (data).....	1091
THROBBING (q-pair).....	470
THYROID (loca).....	1176
THYROID (q-pair).....	303
TIME OUT: PATIENT IDENTIFICATION METHOD.....	695
TIME OUT: VERBAL CONFIRMATION BY TEAM OF PROCEDURE AND LATERALITY .....	697
TIME OUT: VERBAL CONFIRMATION BY TEAM OF PROCEDURE TO BE PERFORMED .....	696
TINGLING (data).....	1091
TINGLING (q-pair).....	470, 618, 619, 620, 621
TINGLING MOUTH/LIPS (data).....	1092
TINGLING MOUTH/LIPS (q-pair).....	183
TINNITUS (data).....	1092
TINNITUS (q-pair).....	498
TO PAIN (data).....	1092
TO PAIN (q-pair).....	133
TO SPEECH (data).....	1092
TO SPEECH (q-pair).....	133
TOENAIL CARE (data).....	1092
TOENAIL CARE (q-pair).....	5
TOILETING BY SCHEDULE (data).....	1093
TOILETING BY SCHEDULE (q-pair).....	88
TOLERATING FLUIDS/NOURISHMENT (data).....	1093
TOLERATING FLUIDS/NOURISHMENT (q-pair).....	151, 465
TONIC/CLONIC (data).....	1093
TONIC/CLONIC (q-pair).....	613
TONOMETRY.....	698
TOOTH (loca).....	1176
TOOTH (q-pair).....	473
TOPICAL TREATMENT (data).....	1093
TOPICAL TREATMENT (q-pair).....	750
TORSADES DE POINTES (data).....	1093
TORSADES DE POINTES (q-pair).....	84
TORSO (loca).....	1176
TORSO (q-pair).....	160, 601
TOTAL DEPENDENCE (data).....	1094
TOTAL DEPENDENCE (q-pair).....	8, 9, 10
TPN (meth).....	1197
TPN (prod).....	1220
TPN (q-pair).....	309
TRACE ST SEGMENT.....	699
TRACE ST SEGMENT 2.....	700
TRACHEA MIDLINE DEVIATED.....	701
TRACHEOSTOMY (data).....	1094
TRACHEOSTOMY (loca).....	1176
TRACHEOSTOMY (q-pair).....	22, 547
TRACHEOSTOMY COLLAR (meth).....	1197
TRACHEOSTOMY COLLAR (q-pair).....	686
TRACTION TYPE.....	702
TRACTION WEIGHT - DO NOT CHECK TOTAL OR SUBTOTALS.....	703
TRAIN OF FOUR RESPONSE.....	704
TRAIN OF FOUR SMS.....	705
TRANSCUTANEOUS (data).....	1094
TRANSCUTANEOUS (q-pair).....	461
TRANSFER (data).....	1094
TRANSFER (q-pair).....	142, 761
TRANSFER WITH ASSIST (data).....	1094
TRANSFER WITH ASSIST (q-pair).....	12

TRANSFUSION (meth).....	1197
TRANSFUSION (q-pair).....	309
TRANSPARENT (data).....	1095
TRANSPARENT (q-pair).....	44, 96, 144, 235, 441, 758
TRANSPARENT DRESSING (data).....	1095
TRANSPARENT DRESSING (q-pair).....	548, 765
TRANSPORT FACILITATOR.....	706
TRANSPORT FACILITATOR PRESENT (data).....	1095
TRANSPORT FACILITATOR PRESENT (q-pair).....	151
TRANSPORTATION MODE/DEVICE.....	707
TRANSPYLORIC (data).....	1095
TRANSPYLORIC (q-pair).....	166
TRANSTRACHEAL (meth).....	1197
TRANSVERSE ABDOMINAL (loca).....	1177
TRANSVERSE ABDOMINAL (q-pair).....	303
TREATMENTS (data).....	1095
TREATMENTS (q-pair).....	164
TREMULOUS.....	708
TRENDELENBURG (posi).....	1203
TRENDELENBURG (q-pair).....	63, 578
TRIGEMINY (data).....	1096
TRIGEMINY (q-pair).....	80
TRIPLE (data).....	1096
TRIPLE (q-pair).....	102
TRIPLE LUMEN CATHETER (meth).....	1198
TUBE POSITION VERIFIED.....	709
TUBED/DYSPHAGIC (data).....	1096
TUBED/DYSPHAGIC (q-pair).....	135
TUBES & DRAINS - CONDOM CATH (task).....	1242
TUBES & DRAINS - FOLEY (task).....	1242
TUBES & DRAINS - I/O CATH (task).....	1242
TUBES & DRAINS - NG TUBE (task).....	1242
TUBES & DRAINS - SUPRAPUBIC DRAIN (task).....	1243
TUBING CHANGE (data).....	1096
TUBING CHANGE (q-pair).....	170
TUBING CHANGED/LABELED (data).....	1096
TUBING CHANGED/LABELED (q-pair).....	343, 493
TURN AND POSITION (data).....	1097
TURN AND POSITION (q-pair).....	5
TURN AND SKIN CARE.....	710
TURNING SELF INDEPENDENTLY (data).....	1097
TURNING SELF INDEPENDENTLY (q-pair).....	5
TWILL TAPE (data).....	1097
TWILL TAPE (q-pair).....	8, 755
TWO CAL HN (data).....	1097
TWO CAL HN (q-pair).....	172
TWO PERSON ASSIST (data).....	1097
TWO PERSON ASSIST (q-pair).....	6, 8, 9, 10
TYMPANIC (loca).....	1177
TYMPANIC (q-pair).....	694

- U -

ULNAR (loca).....	1177
ULNAR (q-pair).....	570
ULNAR DRIFT.....	711
ULTRACAL (data).....	1098
ULTRACAL (prod).....	1220
ULTRACAL (q-pair).....	172
UMBILICAL (loca).....	1177
UMBILICAL (q-pair).....	303

UNABLE TO AROUSE (data).....	1098
UNABLE TO AROUSE (q-pair).....	490
UNABLE TO ASSESS (data).....	1098
UNABLE TO ASSESS (q-pair).....	161, 163
UNABLE TO INFORM (data).....	1098
UNABLE TO INFORM (q-pair).....	403
UNABLE TO READ (data).....	1098
UNABLE TO READ (q-pair).....	161
UNABLE TO RESPOND (data).....	1099
UNABLE TO RESPOND (q-pair).....	163
UNABLE TO VOID (data).....	1099
UNABLE TO VOID (q-pair).....	90
UNABLE TO WEDGE (data).....	1099
UNABLE TO WEDGE (q-pair).....	445
UNBEARABLE (data).....	1099
UNBEARABLE (q-pair).....	470
UNCOOPERATIVE (data).....	1099
UNCOOPERATIVE (q-pair).....	56
UNDER STRETCHER/BED (data).....	1100
UNDER STRETCHER/BED (q-pair).....	515
UNEQUAL (data).....	1100
UNEQUAL (q-pair).....	352, 712
UNEXPECTED NUMBNESS (data).....	1100
UNEXPECTED NUMBNESS (q-pair).....	183
UNIFOCAL PVC (data).....	1100
UNIFOCAL PVC (q-pair).....	80
UNITS (unit).....	1259
UNITS PER HOUR (unit).....	1259
UNITS PER MINUTE (unit).....	1260
UNITS/kg/HOUR (unit).....	1260
UNKNOWN (data).....	1100
UNKNOWN (q-pair).....	194, 615
UNSTABLE (data).....	1101
UNSTABLE (q-pair).....	186, 545
UNSTAGEABLE (data).....	1101
UNSTAGEABLE (q-pair).....	547
UNSTEADY (data).....	1101
UNSTEADY (q-pair).....	199
UNSUCCESSFUL (data).....	1101
UNSUCCESSFUL (q-pair).....	591, 596
UP (data).....	1101
UP TO CHAIR (data).....	1102
UP TO CHAIR (q-pair).....	5
UPPER ABDOMEN (loca).....	1177
UPPER ABDOMEN (q-pair).....	761
UPPER ARM (loca).....	1178
UPPER ARM (q-pair).....	116
UPPER BACK (loca).....	1178
UPPER BACK (q-pair).....	761
UPPER LIMB EXTREMITIES STRENGTH.....	712
URINARY RETENTION (data).....	1102
URINARY RETENTION (q-pair).....	183, 498
URINE (prod).....	1220
URINE (q-pair).....	666
URINE CHARACTERISTICS.....	713
URINE COLOR.....	714
URINE NUMBER VOIDS.....	715
URINE VOIDING METHOD.....	716
UROMETER (data).....	1102
UROMETER (q-pair).....	24
UROSTOMY (data).....	1102
UROSTOMY (q-pair).....	219

V1 (data).....	1102
V1 (q-pair).....	85
V2 (data).....	1103
V2 (q-pair).....	85
V3 (data).....	1103
V3 (q-pair).....	85
V4 (data).....	1103
V4 (q-pair).....	85
V5 (data).....	1103
V5 (q-pair).....	85
V6 (data).....	1103
V6 (q-pair).....	85
VACUUM ASSISTED CLOSURE (data).....	1104
VACUUM ASSISTED CLOSURE (q-pair).....	755
VASCULAR (data).....	1104
VASCULAR (q-pair).....	766
VASCULAR LESION ARTERIAL (data).....	1104
VASCULAR LESION ARTERIAL (q-pair).....	766
VASCULAR LESION MIXED (data).....	1104
VASCULAR LESION MIXED (q-pair).....	766
VASCULAR LESION VENOUS (data).....	1104
VASCULAR LESION VENOUS (q-pair).....	766
VASCULAR/NEUROPATHIC (data).....	1105
VASCULAR/NEUROPATHIC (q-pair).....	766
VASOPRESSOR (data).....	1105
VASOPRESSOR (q-pair).....	545
VDD (data).....	1105
VDD (q-pair).....	455
VEGETARIAN DIET (data).....	1105
VEGETARIAN DIET (q-pair).....	150
VENDOR REP (data).....	1105
VENDOR REP (q-pair).....	554
VENTICULAR FIBRILLATION (data).....	1106
VENTICULAR FIBRILLATION (q-pair).....	84
VENTILATOR (data).....	1106
VENTILATOR (meth).....	1198
VENTILATOR (q-pair).....	24, 504, 577, 578, 686, 717, 718, 719, 720, 721, 722, 723, 724, 725, 726, 727, 728, 729, 730, 731, 732, 733, 734, 735, 736, 737, 738, 739, 740
VENTILATOR ACHIEVA P5O2 (data).....	1106
VENTILATOR ACHIEVA P5O2 (q-pair).....	739
VENTILATOR ADULT STAR (data).....	1106
VENTILATOR ADULT STAR (q-pair).....	739
VENTILATOR AUTO PEEP - DO NOT CHECK TOTAL OR SUBTOTALS.	717
VENTILATOR BEAR 1000 (data).....	1106
VENTILATOR BEAR 1000 (q-pair).....	739
VENTILATOR CIRCUIT CHANGE.....	718
VENTILATOR DRAEGER (data).....	1107
VENTILATOR DRAEGER (q-pair).....	739
VENTILATOR EXHALED VOLUME.....	719
VENTILATOR EXPIRATORY TIME.....	720
VENTILATOR FLOW RATE.....	721
VENTILATOR INSPIRATORY TIME.....	722
VENTILATOR MEAN AIR PRESSURE.....	723

VENTILATOR MINUTE VOLUME.....	724
VENTILATOR MODE.....	725
VENTILATOR OTHER (data).....	1107
VENTILATOR OTHER (q-pair).....	739
VENTILATOR PATIENT EXHALE VOLUME - DO NOT CHECK TOTAL OR SUBTOTALS .....	726
VENTILATOR PATIENT MINUTE VOLUME - DO NOT CHECK TOTAL OR SUBTOTALS .....	727
VENTILATOR PATIENT RATE - DO NOT CHECK TOTAL OR SUBTOTALS .....	728
VENTILATOR PATIENT RATE TOTAL - DO NOT CHECK TOTAL OR SUBTOTALS .....	729
VENTILATOR PB 7200ae (data).....	1107
VENTILATOR PB 7200ae (q-pair).....	739
VENTILATOR PB 840 (data).....	1107
VENTILATOR PB 840 (q-pair).....	739
VENTILATOR PCV SET PRESSURE.....	730
VENTILATOR PEAK PRESSURE - DO NOT CHECK TOTAL OR SUBTOTALS .....	731
VENTILATOR PEEP.....	732
VENTILATOR PRESSURE SUPPORT.....	733
VENTILATOR PSV SET PRESSURE.....	734
VENTILATOR RATE.....	735
VENTILATOR SENSITIVITY.....	736
VENTILATOR SERVO (data).....	1107
VENTILATOR SERVO (q-pair).....	739
VENTILATOR TIDAL VOLUME ACTUAL - DO NOT CHECK TOTAL OR SUBTOTALS .....	737
VENTILATOR TIDAL VOLUME SET.....	738
VENTILATOR TYPE.....	739
VENTILATOR WAVE FORM.....	740
VENTRICULAR (data).....	1108
VENTRICULAR (q-pair).....	461
VENTRICULAR DRAIN (meth).....	1198
VENTRICULAR DRAIN (q-pair).....	332
VENTRICULAR TACHYCARDIA (data).....	1108
VENTRICULAR TACHYCARDIA (q-pair).....	84
VENTRICULAR WHIP PRESENT (data).....	1108
VENTRICULAR WHIP PRESENT (q-pair).....	445
VENTRICULOSTOMY (data).....	1108
VENTRICULOSTOMY (q-pair).....	755
VENTURI MASK (meth).....	1198
VENTURI MASK (q-pair).....	687
VERBAL DE-ESCALATION (data).....	1108
VERBAL DE-ESCALATION (q-pair).....	586
VERBALIZED BY PATIENT (data).....	1109
VERBALIZED BY PATIENT (q-pair).....	695
VERY LOW BMI (data).....	1109
VERY LOW BMI (q-pair).....	545
VIOLET (data).....	1109
VIOLET (q-pair).....	753
VIP PORT (data).....	1109
VIP PORT (q-pair).....	449
VISCOUS (data).....	1109
VISCOUS (q-pair).....	296, 539
VISION CORRECTED.....	741
VISION UNCORRECTED.....	742
VISUAL IMPAIRMENT (data).....	1110
VISUAL IMPAIRMENT (q-pair).....	161
VITAL 1.0 (data).....	1110
VITAL 1.0 (q-pair).....	172
VITAL 1.5 (data).....	1110

VITAL 1.5 (q-pair).....	172
VITAL AF (data).....	1110
VITAL AF (q-pair).....	172
VITAL ENTERAL NUTRITION TYPE (data).....	1110
VITAL ENTERAL NUTRITION TYPE (prod).....	1220
VITAL ENTERAL NUTRITION TYPE (q-pair).....	172
VITAL HN ENTERAL NUTRITION (data).....	1111
VITAL HN ENTERAL NUTRITION (prod).....	1221
VITAL HN ENTERAL NUTRITION (q-pair).....	172
VITAL SIGNS STABLE (data).....	1111
VITAL SIGNS STABLE (q-pair).....	151
VIVONEX RTF (data).....	1111
VIVONEX RTF (q-pair).....	172
VOIDED (data).....	1111
VOIDED (q-pair).....	151, 465
VOLUNTEER (data).....	1111
VOLUNTEER (q-pair).....	517, 706
VOMITING (data).....	1112
VOMITING (q-pair).....	183, 209, 498
VT/VC - DO NOT CHECK TOTAL OR SUBTOTALS.....	743
VVI (data).....	1112
VVI (q-pair).....	455

- W -

WAFFLE (data).....	1112
WAFFLE (q-pair).....	544
WAIST (loca).....	1178
WAIST (q-pair).....	598, 601
WAITING ROOM (data).....	1112
WAITING ROOM (q-pair).....	187, 467
WALKER (data).....	1112
WALKER (q-pair).....	6, 707
WALKING OR PACING (data).....	1113
WALKING OR PACING (q-pair).....	128
WANDER GUARD IN PLACE (data).....	1113
WANDER GUARD IN PLACE (q-pair).....	4
WANDERER (data).....	1113
WANDERER (q-pair).....	40
WARD (data).....	1113
WARD (q-pair).....	467
WARM (data).....	1113
WARM (q-pair).....	658
WARM HUMIDITY WITH O2 (data).....	1114
WARM HUMIDITY WITH O2 (q-pair).....	744
WARMER THAN GENERAL BODY TEMP (data).....	1114
WARMER THAN GENERAL BODY TEMP (q-pair).....	658
WARMER THAN OPPOSITE BODY PART (data).....	1114
WARMER THAN OPPOSITE BODY PART (q-pair).....	658
WARMING BLANKET (data).....	1114
WARMING BLANKET (q-pair).....	744
WARMING BOTTLE (data).....	1114
WARMING BOTTLE (q-pair).....	744
WARMING COOLING DEVICE.....	744
WARMING COOLING DEVICE APPLICATION.....	745
WARMING OFF (data).....	1115
WARMING OFF (q-pair).....	745
WARMING ON (data).....	1115
WARMING ON (q-pair).....	745
WATER (prod).....	1221
WATER (q-pair).....	309, 329
WATERSEAL (data).....	1115

WATERSEAL (q-pair).....	110
WAVEFORM REFLECTS ATRIAL PLACEMENT (data).....	1115
WAVEFORM REFLECTS ATRIAL PLACEMENT (q-pair).....	445
WBC - WASHED (prod).....	1221
WBC - WASHED (q-pair).....	309, 317
WEAK (data).....	1115
WEAK (q-pair).....	196, 222, 352, 712
WEAKNESS (data).....	1116
WEAKNESS (q-pair).....	498
WEDGED (data).....	1116
WEDGED (q-pair).....	445
WEIGHT.....	746
WEIGHT AT ADMISSION.....	747
WEIGHT DRY.....	748
WHEELCHAIR (data).....	1116
WHEELCHAIR (q-pair).....	707
WHEELCHAIR EXIT (posi).....	1203
WHEELCHAIR EXIT ALARM (data).....	1116
WHEELCHAIR EXIT ALARM (q-pair).....	24, 582
WHEELCHAIR SCALE (meth).....	1198
WHEELCHAIR SCALE (q-pair).....	746
WHEEZES (data).....	1116
WHITE (data).....	1117
WHITE (q-pair).....	157, 173, 211, 563, 714, 749, 753
WHOLE BLOOD (prod).....	1221
WHOLE BLOOD (q-pair).....	309
WITH ACTIVITY (meth).....	1199
WITH ACTIVITY (q-pair).....	63, 570, 578
WITH AMBULATION (meth).....	1199
WITH CAST OR BRACE (meth).....	1199
WITH CAST OR BRACE (q-pair).....	746
WITH PHYSICIAN (data).....	1117
WITH PHYSICIAN (q-pair).....	128
WITH PROSTHESIS (meth).....	1199
WITH PROSTHESIS (q-pair).....	746
WITH THERAPIST (data).....	1117
WITH THERAPIST (q-pair).....	128
WITHDRAWN (data).....	1117
WITHDRAWN (q-pair).....	56
WITHDRAWS TO PAIN (data).....	1117
WITHDRAWS TO PAIN (q-pair).....	134
WITHOUT ASSISTANCE (data).....	1118
WITHOUT ASSISTANCE (q-pair).....	6
WITHOUT PROSTHESIS (meth).....	1199
WITHOUT PROSTHESIS (q-pair).....	746
WOUND (data).....	1118
WOUND (q-pair).....	644
WOUND BED COLOR.....	749
WOUND CARE TYPE.....	750
WOUND DEPTH.....	751
WOUND DISCHARGE CHARACTERISTICS.....	752
WOUND DISCHARGE CHARACTERISTICS BLOODY (data).....	1118
WOUND DISCHARGE CHARACTERISTICS BLOODY (q-pair).....	752
WOUND DISCHARGE CHARACTERISTICS CLOUDY (data).....	1118
WOUND DISCHARGE CHARACTERISTICS CLOUDY (q-pair).....	752
WOUND DISCHARGE CHARACTERISTICS NONE (data).....	1118
WOUND DISCHARGE CHARACTERISTICS NONE (q-pair).....	752
WOUND DISCHARGE CHARACTERISTICS PURULENT (data).....	1119
WOUND DISCHARGE CHARACTERISTICS PURULENT (q-pair).....	752
WOUND DISCHARGE CHARACTERISTICS SEROSANGUINOUS (data).....	1119
WOUND DISCHARGE CHARACTERISTICS SEROSANGUINOUS (q-pair).....	752



WOUND DISCHARGE CHARACTERISTICS SEROUS (data).....	1119
WOUND DISCHARGE CHARACTERISTICS SEROUS (q-pair).....	752
WOUND DISCHARGE CHARACTERISTICS VISCOUS (data).....	1119
WOUND DISCHARGE CHARACTERISTICS VISCOUS (q-pair).....	752
WOUND DISCHARGE COLOR.....	753
WOUND DISCHARGE ODOR.....	754
WOUND DRAIN TYPE.....	755
WOUND DRAINAGE (prod).....	1221
WOUND DRAINAGE (q-pair).....	430
WOUND DRAINAGE AMOUNT.....	756
WOUND DRAINAGE AMOUNT HEAVY (data).....	1119
WOUND DRAINAGE AMOUNT HEAVY (q-pair).....	756
WOUND DRAINAGE AMOUNT MODERATE (data).....	1120
WOUND DRAINAGE AMOUNT MODERATE (q-pair).....	756
WOUND DRAINAGE AMOUNT NONE (data).....	1120
WOUND DRAINAGE AMOUNT NONE (q-pair).....	756
WOUND DRAINAGE AMOUNT SCANT (data).....	1120
WOUND DRAINAGE AMOUNT SCANT (q-pair).....	756
WOUND DRAINAGE AMOUNT SMALL (data).....	1120
WOUND DRAINAGE AMOUNT SMALL (q-pair).....	756
WOUND DRESSING CHANGE.....	757
WOUND DRESSING TYPE.....	758
WOUND DRESSING TYPE ANTIMICROBIAL (data).....	1120
WOUND DRESSING TYPE CALCIUM ALGINATE (data).....	1121
WOUND DRESSING TYPE GAUZE (data).....	1121
WOUND DRESSING TYPE GAUZE PACKING (data).....	1121
WOUND DRESSING TYPE GAUZE WRAP (data).....	1121
WOUND DRESSING TYPE HYDROCOLLOIDAL (data).....	1121
WOUND DRESSING TYPE HYDROGEL (data).....	1122
WOUND DRESSING TYPE NONE (data).....	1122
WOUND DRESSING TYPE PRESSURE (data).....	1122
WOUND DRESSING TYPE TRANSPARENT (data).....	1122
WOUND EDGE.....	759
WOUND EDGE INDISTINCT/ATTACHED (data).....	1122
WOUND EDGE INDISTINCT/ATTACHED (q-pair).....	759
WOUND EDGE NON-ATTACHED (data).....	1123
WOUND EDGE NON-ATTACHED (q-pair).....	759
WOUND EDGE ROLLED (data).....	1123
WOUND EDGE ROLLED (q-pair).....	759
WOUND EDGE TUNNELING (data).....	1123
WOUND EDGE TUNNELING (q-pair).....	759
WOUND EDGE UNDERMINING (data).....	1123
WOUND EDGE UNDERMINING (q-pair).....	759
WOUND HEALING STATUS.....	760
WOUND IDENTIFIED/SITE.....	761
WOUND IRRIGATION (prod).....	1222
WOUND LENGTH.....	762
WOUND ODOR.....	763
WOUND ODOR FOUL (data).....	1123
WOUND ODOR FOUL (q-pair).....	763
WOUND ODOR NONE (data).....	1124
WOUND ODOR NONE (q-pair).....	763
WOUND PHOTOGRAPHED.....	764
WOUND TREATMENT.....	765
WOUND TYPE.....	766
WOUND VAC DEVICE (data).....	1124
WOUND VAC DEVICE (q-pair).....	755
WOUND VACUUM (prod).....	1222
WOUND VACUUM (q-pair).....	430
WOUND VACUUM ASSISTED CLOSURE.....	767
WOUND WIDTH.....	768
WRIST (loca).....	1178

WRIST (q-pair).....116

- X -

X-RAY (data).....1124  
X-RAY (q-pair).....212  
XENADERM (data).....1124  
XENADERM (q-pair).....548  
XEROFORM (data).....1124  
XEROFORM (q-pair).....548, 765

- Y -

YEAR (unit).....1260  
YELLOW (data).....1125  
YELLOW (q-pair).....157, 173, 211, 563,  
638, 681, 714, 749,  
753

YELLOW/FIBRINOUS (data).....1125  
YELLOW/FIBRINOUS (q-pair).....537  
YELLOW/GREEN (data).....1125  
YELLOW/GREEN (q-pair).....681  
YES (data).....1125  
YES (q-pair).....14, 16, 18, 23, 33,  
34, 35, 36, 37, 38,  
39, 41, 42, 45, 49,  
51, 61, 75, 76, 78,  
93, 98, 104, 107,  
111, 112, 139, 189,  
191, 201, 208, 217,  
226, 228, 233, 234,  
241, 295, 298, 303,  
304, 305, 306, 307,  
418, 440, 446, 451,  
454, 456, 458, 466,  
468, 482, 496, 505,  
506, 531, 540, 543,  
546, 549, 551, 561,  
576, 583, 584, 585,  
587, 588, 589, 590,  
606, 616, 617, 637,  
639, 641, 642, 643,  
654, 655, 701, 709,  
710, 718, 757, 764

- { -

{8BC7F294-7126-4843-91A0-14F805EC4578} (data).....1126  
{8BC7F294-7126-4843-91A0-14F805EC4578} (q-pair).....181