

=====

TERM: ABDOMEN BOWEL SOUNDS Page: 1
ID: {A7E05486-2BF9-D304-4925-9A32165B19F0}

=====

Abbreviation..... AB BOW S
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Abdominal Bowel Sounds
Display Name..... Abdomen Bowel Sounds
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 6 entries

-
1. NORMAL (n) {553C9A97-A3D5-C52F-0E4E-C7D168C99A2B}
 2. HYPOACTIVE (HYPOACTIVE) {2E0C516D-3858-4A1F-A2F3-BF0AB9E3A7FC}
 3. HYPERACTIVE (HYPERACTIV) {1A548FFE-8645-4C0D-899D-5D91430CC1B1}
 4. FLATULENCE (FLATULENCE) {F6ECE7D8-D9C8-485C-8FAC-03DA9D7E0B66}
 5. ABSENT (a) {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0}
 6. BRUIT (BRUIT) {2627E78A-B8B3-44EE-BC6F-4DAB16394F8C}

LOCATION: 5 entries

-
1. ABDOMEN LUQ (ABD LUQ) {3B31E1EC-CE74-4C76-BB61-10C678431FCA}
 2. ABDOMEN RUQ (ABD RUQ) {F29C9FA1-717A-49B7-8000-95070E9AACAO}
 3. ABDOMEN RLQ (ABD RLQ) {2C48394F-6877-4738-856C-E933644D8067}
 4. ABDOMEN LLQ (ABD LLQ) {39831F59-D7B0-4156-AC57-4C622A93CA58}
 5. GENERAL UMBILICUS (GEN UMBILI)
{34CEEF50-14E8-44B2-9FEA-3BAB36D08DFA}

UNITS:

=====

TERM: ABDOMINAL CONDITION Page: 2
ID: {761ED47D-0598-4B17-EA65-3436FDDB8597}

=====

Abbreviation..... ABDOM COND
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... CURRENT OBSERVATION OF ABDOMINAL

CONDITION

Display Name..... Abdominal Condition
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 10 entries

- 1. SOFT (SOFT) {E68B4D4B-B86B-4004-A861-986645267F6E}
- 2. FIRM (FIRM) {11E23F20-39A0-4A6D-B198-04BB36B0C909}
- 3. HARD (HARD) {7934B543-7028-4CBC-9308-99F5A952A1DE}
- 4. DISTENDED (DISTENDED) {B5C529CA-2BE0-4698-A153-718F4D58E643}
- 5. TENDER (TENDER) {30486CC7-3097-4F95-9238-FD6C2126CD10}
- 6. NON-TENDER (NON-TENDER) {5A3D41C0-A83D-4382-A573-D6884BB1670B}
- 7. GUARDING (GUARDING) {9EFD10A7-B1B7-43DA-B289-646AF5C0A25C}
- 8. FLAT (Flat) {EA4F1980-EDF3-82BC-B677-D747F48F1163}
- 9. ROUND (ROUND) {E8E741F1-F0C5-4B77-B014-63B545BED755}
- 10. SCAPHOID (SCAPHOID) {1BE93489-42FC-44EB-8D07-87F4224979F3}

UNITS:

TERM: ABDOMINAL GIRTH Page: 3
ID: {F70E6642-2719-22BE-BE87-DEF0A884F177}

Abbreviation..... AG
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Girth of patient abdomen
Display Name..... Abdominal Girth
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid..... 4536404

UNITS:

CENTIMETER 0 > 76 * 152 < 300 Dec Prec: 0

TERM: ADL - ACTIVITIES OF DAILY LIVING Page: 4
ID: {B5D8ADB2-0D19-4C07-8A09-6053A4D27003}

Abbreviation..... ADL ADL
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Activities of daily living assessment
Display Name..... Activities of Daily Living

Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 25 entries

- 1. COMPLETE BATH (COMP BATH) {A7CB1F56-F88E-403B-AABC-105BE129FBB0}
- 2. PARTIAL BATH (PART. BATH) {60B9C7AD-9E6E-4412-9BB3-57AFF956C1FC}
- 3. SELF AM CARE BY PATIENT (SELF AM CR)
{C48A46C1-0E68-479C-B803-08C80617ADFF}
- 4. TURN AND POSITION (TURN POSIT)
{44A4F470-1760-4073-AB3A-D59CA9EFE7F3}
- 5. UP TO CHAIR (UP TO CHAI) {42E598E5-82EC-428B-B5FA-1129E3095C78}
- 6. MEAL (MEAL) {85D38DCD-6A4A-4908-ABE2-4E92ABEF1321}
- 7. AMBULATED (AMBULATED) {017D99FE-91D5-49F1-AC15-50620879C466}
- 8. SNACK (SNACK) {60B77974-4102-4511-A105-E24850B924D8}
- 9. SELF CARE OF BEDSIDE ENVIRONMENT (S CARE BED)
{A1BEC7F4-D404-4824-8BB7-108E0A6C43A2}
- 10. SHOWER (SHOWER) {C259929A-44DB-48FD-A2EC-BC6ECF9D0476}
- 11. SELF SHOWER (SELF SHOWE) {C35CCAE8-27D9-46F2-8971-2C4CF4924CCB}
- 12. SHAMPOO (SHAMPOO) {95CCB23C-E614-40CD-84F0-EE4EFBC95271}
- 13. DRY SHAMPOO (DRY SHAMP) {8910A186-2590-4352-A884-7D692A9D2B3C}
- 14. NAIL CARE (NAIL CARE) {B783A171-1F0A-4135-8535-85776B81AB68}
- 15. TOENAIL CARE (TOEN CARE) {6CF3ABCE-FF74-43E2-81E5-D23339017CA9}
- 16. FOOT CARE (FOOT CARE) {332E8338-6DAE-4BC9-B195-63DAD57937D0}
- 17. SELF FOOT CARE (SELF FOOT)
{F2DCAF0E-389D-4987-8C73-66E35266CC78}
- 18. TURNING SELF INDEPENDENTLY (TURN SELF)
{DD804CEB-F7CA-4D42-BF6A-447E38B55E98}
- 19. SHAVE (SHAVE) {A5CB33DF-C33B-4DD3-B00D-78D78F7966C7}
- 20. SELF SHAVE (SELF SHAVE) {90F1760C-4C12-4976-ADCB-DBF456A34023}
- 21. FACE WASH (FACE WASH) {E1162A8B-D26D-4EB3-988F-50583B1849EA}
- 22. SELF FACE WASH (S FACE WAS)
{77457D46-5676-46E6-B38F-FE3AF70B5CBC}
- 23. PERINEAL CARE (PERIN CARE)
{D88736E4-E0EB-4FF7-9199-288A4D051A20}
- 24. SELF TRACH CARE (S TRACH CA)
{70BE79EC-107F-4CD3-B964-7919B3DDE222}
- 25. CLOTHED SELF (CLOTH S) {33802A30-FAA7-426B-ADAC-5297818DA7DE}

UNITS:

=====

TERM: ADL - AMBULATION ASSISTANCE
ID: {BF48D1FD-555C-4714-A624-2A350F74D8A9}

=====

Page: 5

Abbreviation..... ADL AMB AS
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Ambulation Assistance

Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Void.....

DATA VALUE: 8 entries

-
1. INDEPENDENT (INDEPENDEN) {11BC0CF7-1C29-49C7-B7C8-2BCF8FF82250}
 2. WITHOUT ASSISTANCE (WO ASSIS)
 {41CE081D-CB0E-49AF-936A-4FD78E61681B}
 3. ONE PERSON ASSIST (1 PERS AS)
 {A3D14AA8-D658-412F-BBE8-D8AD195ADA3D}
 4. TWO PERSON ASSIST (2 PERS AST)
 {0547070B-C9BC-4406-BDC0-BF9C6C23D41A}
 5. WALKER (WALKER) {91247001-5E89-46A0-A841-B4FAE9226E63}
 6. CANE (CANE) {562B2E0B-0298-44C6-B33B-BC57DE6C3C3C}
 7. CRUTCHES (CRUTCHES) {DDA81269-B25C-4F7D-A059-BB1602B2B179}
 8. GAIT BELT (GAIT BELT) {AD9AB19B-92B6-40BE-B906-FBFC46C39E0D}

UNITS:

=====

TERM: ADL - AMBULATION DISTANCE Page: 6
ID: {49097AB5-ACC4-4FC5-91FB-B62DFF5EBA51}

=====

Abbreviation..... ADL AM D
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Ambulation Distance
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

FOOT -1 > 0 * 500 < 600 Dec Prec: 0

=====

TERM: ADL - BED MOBILITY Page: 7
ID: {ADB040D9-C136-4AF6-809D-CDD868DA3DCA}

=====

Abbreviation..... ADL BED MO
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Bed Mobility
Help Text..... No Help Text

Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 9 entries

- 1. INDEPENDENT (INDEPENDEN) {11BC0CF7-1C29-49C7-B7C8-2BCF8FF82250}
- 2. SUPERVISION/SET-UP HELP (SUPERVISED)
{596D028B-B60D-4244-A1B0-882177C8D25F}
- 3. LIMITED ASSISTANCE (LIM ASSIST)
{9A8C1310-C840-431E-BA21-065B7EA7B768}
- 4. EXTENSIVE ASSISTANCE (EXT ASSIST)
{C402B659-E880-4070-8345-13E5508B4D29}
- 5. TOTAL DEPENDENCE (TOT DEPEND)
{38702CDD-E440-4EFA-9583-36EEE719F7C3}
- 6. ONE PERSON ASSIST (1 PERS AS)
{A3D14AA8-D658-412F-BBE8-D8AD195ADA3D}
- 7. TWILL TAPE (Twill Tape) {C52BE250-732D-03BB-B8F9-72BE4D546807}
- 8. TWO PERSON ASSIST (2 PERS AST)
{0547070B-C9BC-4406-BDC0-BF9C6C23D41A}
- 9. NOT DONE (NOT DONE) {C3FE84AC-EEA1-4153-AE6C-69755C2E6BB1}

UNITS:

TERM: ADL - DRESSING Page: 8
ID: {3FB32B93-096A-4687-B4BE-E3FF103D75F6}

Abbreviation..... ADL DRESS
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Dressing
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 8 entries

- 1. INDEPENDENT (INDEPENDEN) {11BC0CF7-1C29-49C7-B7C8-2BCF8FF82250}
- 2. SUPERVISION/SET-UP HELP (SUPERVISED)
{596D028B-B60D-4244-A1B0-882177C8D25F}
- 3. LIMITED ASSISTANCE (LIM ASSIST)
{9A8C1310-C840-431E-BA21-065B7EA7B768}
- 4. EXTENSIVE ASSISTANCE (EXT ASSIST)
{C402B659-E880-4070-8345-13E5508B4D29}
- 5. TOTAL DEPENDENCE (TOT DEPEND)
{38702CDD-E440-4EFA-9583-36EEE719F7C3}
- 6. ONE PERSON ASSIST (1 PERS AS)
{A3D14AA8-D658-412F-BBE8-D8AD195ADA3D}
- 7. TWO PERSON ASSIST (2 PERS AST)

{0547070B-C9BC-4406-BDC0-BF9C6C23D41A}
8. NOT DONE (NOT DONE) {C3FE84AC-EEA1-4153-AE6C-69755C2E6BB1}

UNITS:

=====

TERM: ADL - EATING Page: 9
ID: {6293D18C-0742-4D1B-B410-E8641798DD2B}

=====

Abbreviation..... ADL EATING
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Eating
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 9 entries

-
1. INDEPENDENT (INDEPENDEN) {11BC0CF7-1C29-49C7-B7C8-2BCF8FF82250}
 2. SUPERVISION/SET-UP HELP (SUPERVISED)
{596D028B-B60D-4244-A1B0-882177C8D25F}
 3. SET UP HELP ONLY (SETUP ONLY)
{DACC77AB-5AE7-4F1B-A3E1-C347D3DD7219}
 4. LIMITED ASSISTANCE (LIM ASSIST)
{9A8C1310-C840-431E-BA21-065B7EA7B768}
 5. EXTENSIVE ASSISTANCE (EXT ASSIST)
{C402B659-E880-4070-8345-13E5508B4D29}
 6. TOTAL DEPENDENCE (TOT DEPEND)
{38702CDD-E440-4EFA-9583-36EEE719F7C3}
 7. ONE PERSON ASSIST (1 PERS AS)
{A3D14AA8-D658-412F-BBE8-D8AD195ADA3D}
 8. TWO PERSON ASSIST (2 PERS AST)
{0547070B-C9BC-4406-BDC0-BF9C6C23D41A}
 9. NOT DONE (NOT DONE) {C3FE84AC-EEA1-4153-AE6C-69755C2E6BB1}

UNITS:

=====

TERM: ADL - EMOTIONAL SUPPORT Page: 10
ID: {8E6573BB-DED8-47D3-AA52-C56404C5D1D9}

=====

Abbreviation..... EMO SUPPT
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description

Display Name..... Emotional Support
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

- 1. PATIENT (PATIENT) {7881F205-BCD1-4948-AE66-AD572F5A02F2}
- 2. FAMILY (FAMILY) {1233C4BB-5AC5-4883-AF2A-F5C93CD8814F}
- 3. FRIEND (FRIEND) {03902E48-AF86-432A-8372-37796A55F625}
- 4. SIGNIFICANT OTHER (SO) {082C42E8-9E8A-4A96-B3AD-8B0E0A62AF00}
- 5. CAREGIVER (CAREGIVER) {F9B8BABE-2F05-437D-85C8-585FC1FFC07F}

UNITS:

TERM: ADL - OTHER ACTIVITY Page: 11
ID: {C0EC0C79-A4AB-4CB6-B026-E193E033A3DE}

Abbreviation..... ADL OTHER
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Other Activity
Help Text..... Other ADL Activities
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 14 entries

- 1. BEDREST (BEDREST) {8CD28057-132E-4CEE-9EA1-F03C8AF50CE8}
- 2. DANGLED (DANGLED) {40665750-112F-401D-AEFA-8E19F45EC0C9}
- 3. CHAIR (CHAIR) {5B8A5DBD-7462-4653-8D1F-CA8B1226DE38}
- 4. COMMODE (COMMODE) {75C896CA-F105-4EF0-B641-3F2B1F984F1B}
- 5. PASSIVE RANGE OF MOTION (PASS ROM)
{DD43B9EF-E18F-4203-8564-89EBBAD963A4}
- 6. ACTIVE RANGE OF MOTION (ACTIVE ROM)
{F55E68F3-6788-47AE-B962-119413B4A089}
- 7. GERI-CHAIR (g) {167F9786-8F24-054F-020A-658A6EE5B0B4}
- 8. GURNEY CHAIR (GURN CHAIR) {B6EA0944-8414-42FD-88C0-50155DE82896}
- 9. BATHROOM PRIVILEGES (BATHR PRIV)
{29C1B6DD-3308-4DED-9970-1322A97C633B}
- 10. TRANSFER WITH ASSIST (TRANS ASST)
{62A5BFF5-9C7C-4412-83A8-F872095DE005}
- 11. FULL SPINAL PRECAUTIONS (FULL S P)
{42ADB74D-12CC-4111-B181-626C0F01BC7C}
- 12. THORACIC SPINAL PRECAUTIONS (THOR S P)
{D55AEF23-7F5E-45DD-A178-F217E6AF3873}
- 13. CERVICAL SPINAL PRECAUTIONS (CERV S P)
{A15F7127-2367-464A-AB7D-1E6A7C09AA21}

14. LUMBAR SPINAL PRECAUTIONS (LUMB S P)
{5AFBBBAB-F919-436A-9C5A-71DCF5D79D18}

UNITS:

=====

TERM: ADL - REPOSITION Page: 12
ID: {AC46FCA8-4D8B-45B9-9938-9CE75B759E76}

=====

Abbreviation..... ADL REPOSI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Reposition
Display Name..... Reposition
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 9 entries

-
1. LEFT (LEFT) {81063CEE-2E83-4C95-AAEE-76B7DCB828D4}
 2. RIGHT (RIGHT) {AB74C579-DF14-4D7B-A9EF-5F119089F976}
 3. SUPINE (SUPINE) {A18F9534-17AF-47D7-8D14-2AA457A889BD}
 4. PRONE (PRONE) {41180E09-D9D6-4903-BAE9-28CD5AAC4991}
 5. SELF CARE (SELF CARE) {0E068188-75B8-46C5-8947-A6CC5C7849B5}
 6. ASSISTED (ASSISTED) {65069804-39BF-41ED-9EA3-BF2C41389F10}
 7. CARDIAC REPOSITION (CARD REPOS)
{E133FE5A-20E1-418F-BD90-D2DFE14908C0}
 8. FOWLER'S (FOWLER'S) {8E691813-A9AC-48BE-9401-528A7D762A32}
 9. HIGH FOWLER'S (HIGH FOWLE)
{EA23C94D-EE91-4552-9F74-E051FAF3E194}

UNITS:

=====

TERM: ADL - VISITORS Page: 13
ID: {C0C8630F-6535-4C89-9556-D97CBCB41952}

=====

Abbreviation..... ADL VISITO
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... Visitors
Display Name..... Visitors
Help Text..... Has the patient had any visitors today.
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
- 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: ADL- ORAL CARE Page: 14
 ID: {0F21C7E0-AD33-46CC-B5BD-40F7A46946D4}

=====

Abbreviation..... ADL ORAL C
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... Oral Care
 Display Name..... Oral Care
 Help Text..... Oral Care
 Multi Select Picklist..... Yes
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 12 entries

- 1. DONE (DONE) {60412C8D-132F-46EA-8DDF-FCC282C4CC8B}
- 2. ASSISTED (ASSISTED) {65069804-39BF-41ED-9EA3-BF2C41389F10}
- 3. SELF CARE (SELF CARE) {0E068188-75B8-46C5-8947-A6CC5C7849B5}
- 4. OFFERED (OFFERED) {715091C2-9668-4678-900E-6D55D3335A30}
- 5. REFUSED (REFUSED) {D5B3E135-4D68-4528-BBC9-6E8A7A997104}
- 6. NOT DONE (NOT DONE) {C3FE84AC-EEA1-4153-AE6C-69755C2E6BB1}
- 7. SUCTIONED ONLY (SUCT ONLY)
 {7A8B23FD-69C5-4931-A6D4-E92F8FD70817}
- 8. LIP SALVE (LIP SALVE) {612F78AB-7386-4176-9896-02DE5A8E836D}
- 9. TEETH BRUSHED (TEETH BRUS)
 {20841794-6E6D-4262-B1FE-9E3F20F8CC49}
- 10. MOUTHWASH USED (MOUTHWASH)
 {9F7EF355-12A4-4EBA-92ED-0F38273E411B}
- 11. ETT REPOSITIONED (ETT REPOS)
 {95B35EDB-76B6-433E-B526-01C543A7CCEA}
- 12. ETT RETAPED/SECURED (ETT RETAPE)
 {6D8A2EDA-CAC9-44D2-BDA9-3B7F20908F1E}

UNITS:

=====

TERM: ADI - SPIRITUAL CARE Page: 15
 ID: {42608477-6B3D-4036-BE05-F0DA3ED218B9}

=====

Abbreviation..... SPIRIT CA
 Active..... Yes

Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... No Description
Display Name..... Spiritual care
Help Text..... Has Patient Had Spiritual Care
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
- 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

TERM: AFFECT Page: 16
ID: {CAE19296-4C9E-0049-DD90-4939FFFAA7F4}

Abbreviation..... AFFECT
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... an expressed or observed emotional response
Display Name..... Affect
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

- 1. FLAT (Flat) {EA4F1980-EDF3-82BC-B677-D747F48F1163}
- 2. BLUNTED (Blunted) {F119F2DA-9182-2480-BA75-02D39F4FF69B}
- 3. BRIGHT (Bright) {D16542CE-0403-2C68-A361-125CB3C75ABA}
- 4. CONSTRICTED (Constrict) {A7F66C5B-89EF-DEEE-A6C0-0DCDFEB3AEB4}

UNITS:

TERM: AIRWAY ARTIFICIAL CUFF STATUS Page: 17
ID: {18958F14-76F1-B4F3-DB01-F5FB4C8A70BC}

Abbreviation..... ART AIR C
Active..... Yes
Boolean Value False..... Deflated
Boolean Value True..... Inflated
Data Type..... Boolean

Description..... AIRWAY CUFF INFLATED OR DEFLATED
Display Name..... Artificial Airway Cuff Status
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: AIRWAY ARTIFICIAL POSITION AT LIP/NOSE
ID: {A9032CCC-5479-87C2-27F8-46D8B4F68272}

=====

Page: 18

Abbreviation..... ART AIR PO
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Artificial Airway Position at Lip/Nose
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

CENTIMETER 11.9 > 12.0 * 28.5 < 28.6 Dec Prec: 1

=====

TERM: AIRWAY ARTIFICIAL SIZE - FRENCH
ID: {48027138-FA36-067B-AC88-61C42902C7CC}

=====

Page: 19

Abbreviation..... ART AIR F
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Artifical Airway Size -French
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

FRENCH 3.9 > 4.0 * 14.0 < 14.1 Dec Prec: 1

=====

TERM: AIRWAY ARTIFICIAL SIZE - MILIMETERS
ID: {C6712EFD-7351-46E2-B31F-3F788CDECA28}

=====

Page: 20

Abbreviation..... ART AIR mm
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Artificial Airway Size - mm
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 13 entries

1. 4 (4)	{D2280C7E-DE50-5A8D-C3FD-5862B69C0E1A}
2. 4.5 (4.5)	{EC6C0C72-8C49-4535-92B2-7E543007EE32}
3. 5 (5)	{188EDCE7-0726-4234-D8D2-92C8052CD221}
4. 5.5 (5.5)	{4878F059-EC9D-44FC-9A61-734C291530E7}
5. 6 (6)	{640CF2F9-3BBF-AB13-185B-CD02A609FB5A}
6. 6.5 (6.5)	{CAFCEBDB1-8533-4C10-8DD1-2BD07A4CF6D3}
7. 7 (7)	{7A010324-D2C2-5D4B-848A-7B8878666D20}
8. 7.5 (7.5)	{7517EE9C-7C8B-4AF0-B60A-88F4E24BFA5E}
9. 8 (8)	{518E7CFC-102E-351E-C7AF-FB575E3C8F4D}
10. 8.5 (8.5)	{96723567-151D-4B5C-9C65-B3EB6D584272}
11. 9 (9)	{58FF2FD4-569F-F454-9F5B-FE6F9E34D8CC}
12. 9.5 (9.5)	{252555A1-ADB2-4948-8764-C952F2587DC5}
13. 10 (10)	{0C827160-18DB-816F-7FAA-A3C175C6D33B}

UNITS:

=====

TERM: AIRWAY ARTIFICIAL TYPE
ID: {2FC0523A-8EBD-7E24-E453-C2DEF54DAE44}

=====

Page: 21

Abbreviation..... ART AIRW T
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... TYPE OF ARTIFICIAL AIRWAY
Display Name..... Artificial Airway Type
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 8 entries

1. NASAL ENDOTRACHEAL (Nasal Endo)
{C631352E-CAB5-0409-6133-E5D11ECC3293}
2. ORAL ENDOTRACHEAL (ORAL ENDO)
{466CAD16-D7F3-9D73-6500-959C4DA50D9A}
3. TRACHEOSTOMY (TRACHEOSTO) {2B0B7599-847E-852F-D16F-091CAAA4B205}
4. LARYNGECTOMY (LARYNGECTO) {225CB76A-BAF5-827E-729B-C3D60AD5B50A}
5. NASAL TRUMPET (NASAL TRUM)
{2FFB8510-793C-D6A4-838E-20FF7CD26533}
6. BITE BLOCK (BITE BLOCK) {E7CC51F3-08A8-5FEE-38B2-16B74E1F542E}
7. ORAL AIRWAY (ORAL AIRWA) {B0CF4270-A68D-1322-0621-1DFF8D677B40}
8. CASS ENDOTRACHEAL TUBE (CASS)
{02BF805C-DF35-4DAB-A99F-2ED673174E49}

UNITS:

=====

TERM: ALARM PARAMETERS VERIFIED
ID: {415CE93B-0510-423D-9CAC-838F59882D94}

=====

Page: 22

Abbreviation..... ALARM PARA
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... Alarm Parameters Verified
Display Name..... Alarm Parameters Verified
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: ALARMS
ID: {97BADE57-C863-F3B4-4507-74BF03D4B0B9}

=====

Page: 23

Abbreviation..... ALARMS
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Alarms
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION

Void.....

DATA VALUE: 7 entries

- 1. CENTRAL ECG (cecg) {26D36175-E4E3-C364-3018-1B0DB4DC4793}
- 2. BED EXIT ALARM (BED EXIT A)
{AC53759E-FFAE-5E0F-2EBA-E56A2712357C}
- 3. WHEELCHAIR EXIT ALARM (WHEEL EXIT)
{7D06564E-ECDE-4A11-BEFE-92AB7E448AA3}
- 4. IV PUMP (IV PUMP) {EFB92C16-CFA3-18BB-9D54-590089DA3962}
- 5. VENTILATOR (v) {A5A647ED-A56D-5012-AD5B-788B3DCB3498}
- 6. UROMETER (u) {13074A3A-19D7-300C-7BFF-8DC18DDE2951}
- 7. IABP (IABP) {763FE975-F5E5-80A6-AB16-E903A829E665}

UNITS:

=====

TERM: ALDRETE - ACTIVITY

Page: 24

ID: {97F85378-357C-4BF5-8198-4323E2656604}

=====

Abbreviation..... ALD ACTIVI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Modified Aldrete Scoring System
Display Name..... Aldrete Score - Activity
Help Text..... Ability to move voluntary: 0=0
extremeties, 1=2 extremeties, 2=4
extremeties
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 3 entries

- 1. 0 = 0 EXTREMITIES (ALD ACT 0)
{6E248A9A-B509-4BD0-A01D-9ED5CE5B8FC9}
- 2. 1 = 2 EXTREMITIES (ALD ACT 1)
{6DCB4762-72CE-48C2-9AFE-B19969ED1DE0}
- 3. 2 = 4 EXTREMITIES (ALD ACT 2)
{3434DF4B-5CB6-417F-B5BE-CEF2D63FE8E0}

UNITS:

=====

TERM: ALDRETE - CIRCULATION

Page: 25

ID: {42D5FBEB-5FBA-4E89-A1CD-D7882740C7B4}

=====

Abbreviation..... ALD CIRCUL
Active..... Yes
Boolean Value False.....

Boolean Value True.....
Data Type..... PickList
Description..... Modified Aldrete Scoring System
Display Name..... Aldrete Score - Circulation
Help Text..... BP +/- 50% Preop=0, 20-49% Preop= 1,
20% Preop = 2
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

-
1. 0 = BP +/- 50% PREOP (ALD CIR 0)
{675F2CD1-DD7B-473B-B8EC-C2BFBF4AC3B2}
 2. 1 = BP +/- 20-49% PREOP (ALD CIR 1)
{D4A6100D-1228-4860-9BD6-B43F30D4FF94}
 3. 2 = BP +/- 20% PREOP (ALD CIR 2)
{04B67413-7A88-409B-8005-471A7C84B879}

UNITS:

=====

TERM: ALDRETE - CONSCIOUSNESS
ID: {8887F08F-099D-4FE7-84EB-FB122BAFF5F5}

=====

Page: 26

Abbreviation..... ALD CONSCI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Modified Aldrete Scoring System
Display Name..... Aldrete Score - Consciousness
Help Text..... Not responding=0, Arousable with
minimal stimulation=1, Fully awake=2
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

-
1. 0 = NOT RESPONDING (ALD CON 0)
{E8765215-B0B0-430C-8C89-AB969A6AD12E}
 2. 1 = AROUSABLE (ALD CON 1) {B2A07804-9829-42D0-9C49-A61FD51EEAA0}
 3. 2 = FULLY AWAKE (ALD CON 2)
{903D0F56-ADEC-4527-BA66-C8A568E3EAC0}

UNITS:

=====

TERM: ALDRETE - O2 SATURATION
ID: {5530278D-7914-4D0B-B3AF-852FFC708808}

=====

Page: 27

Abbreviation..... ALD O2 SAT
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Modified Aldrete Scoring System
Display Name..... Aldrete - O2 Saturation
Help Text..... <90% even with O2=0, Needs O2 to
maintain >90%=1, >92% on room air=2
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

-
1. 0 = <90% (ALD O2 0) {30042730-F600-44FD-A3D0-53172B70082B}
 2. 1 = NEEDS O2 (ALD O2 1) {5D97F6BB-7B02-4644-8CE5-014C89826731}
 3. 2 = ON ROOM AIR (ALD O2 2)
{D58C5B46-CCD3-4E6D-86E8-7A1CFA92BC1C}

UNITS:

=====

TERM: ALDRETE - RESPIRATION Page: 28
ID: {165534D8-3229-4865-B1D5-9B579DBF8E7D}

=====

Abbreviation..... ALD RESPIR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Modified Aldrete Scoring System
Display Name..... Aldrete - Resipration
Help Text..... Apneic=0, Dyspnea, shallow or
limited=1, Breathes deeply, coughs
freely=2
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

-
1. 0 = APNEIC (ALD RESP 0) {649C42B0-6E0A-4E86-AAC5-70A8C20EC511}
 2. 1 = DYSPNEA (ALD RESP 1) {3DAE60A2-57DA-4121-A507-BC43497C7B67}
 3. 2 = BREATHES FREELY (ALD RESP 2)
{278CD44C-538E-4015-9D98-D9312A7039BA}

UNITS:

=====

TERM: ALDRETE - TOTAL Page: 29
ID: {A947CF6C-094E-49B9-A856-E4103162D875}

=====
Abbreviation..... ALD TOTAL
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Modified Aldrete Scoring System Total
Display Name..... Aldrete - Total Score
Help Text..... Total value of Aldrete scores
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

MODIFIED ALDRETE SCORE -1 > 0 * 10 < 11 Dec Prec: 0

=====
TERM: APPEARANCE Page: 30
ID: {5B202547-7AD3-4E86-9FBE-82AEC450510B}

=====
Abbreviation..... APPEARAN
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Appearance
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 4 entries

- 1. NEAT (NEAT) {8B01F7B8-B626-4143-A092-1E732EEB4A05}
2. APPROPRIATE (a) {1CFB285F-34DF-73A2-B89C-F6EE1F40E7E9}
3. INAPPROPRIATE (Inappropri) {7E3B54C2-8A49-28EC-055E-FEA87599421C}
4. DISHEVELED (DISHEVEL) {0455549E-5CCB-470D-9D7F-F3959E6D55FE}

UNITS:

=====
TERM: ARM BAND - ALLERGIES Page: 31
ID: {D72DC694-6192-8D57-FD8B-02E49EF45DE5}

=====
Abbreviation..... AB ALLERGY
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... Is the arm band for allergies present

Display Name..... Arm Band - Allergies
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: ARM BAND - BLOOD TRANSFUSION IV
ID: {66410CB8-1475-12D5-1199-76217CC9B174}

=====

Page: 32

Abbreviation..... AB BL TR
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... No Description
Display Name..... Arm Band - Blood Transfusion IV
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: ARM BAND - DNR
ID: {B7AF531A-26FA-B821-F3B6-F75911A8BBC2}

=====

Page: 33

Abbreviation..... AB DNR
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... Is patient armband for DNR present
Display Name..... Arm Band - DNR
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

-
- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
 - 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: ARM BAND - ELOPMENT

Page: 34

ID: {950320E8-0C68-29C0-E39A-355CEB0119C6}

=====

Abbreviation..... AB ELOPE
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... No Description
Display Name..... Arm Band - Elopment
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

-
- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
 - 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: ARM BAND - HIGH FALL RISK

Page: 35

ID: {48BCE5FD-3703-4F3B-619B-8E5FD2E25D2E}

=====

Abbreviation..... AB FALL RI
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... No Description
Display Name..... Arm Band - High Fall Risk
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

-
- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
 - 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: ARM BAND - PERSONAL MEDICAL ALERT
ID: {4EF7E9C0-4B52-0EA5-B5ED-D29E655089FE}

=====

Page: 36

Abbreviation..... AB PM ALER
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... No Description
Display Name..... Arm Band - Personal Medical Alert
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

-
- | | |
|------------|--|
| 1. NO (n) | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEA649CD} |

UNITS:

=====

TERM: ARM BAND - WANDERER
ID: {F0D2FBB6-0843-239A-356E-D903542E9117}

=====

Page: 37

Abbreviation..... AB WANDER
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... No Description
Display Name..... Arm Band - Wanderer
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

-
- | | |
|------------|--|
| 1. NO (n) | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEA649CD} |

UNITS:

=====

TERM: ARM BANDS PRESENT
ID: {FD6E4F7D-AF64-4208-B6BD-DD965FB234C3}

=====

Page: 38

=====
Abbreviation..... AB PRESENT
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Arm Bands Present
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Void.....

DATA VALUE: 9 entries

-
1. N/A (n) {B3DBE3FD-C5CB-3894-607A-BE13319594D0}
 2. IDENTIFICATION (i) {39274F69-556C-7F05-6D17-CE8A03D1867D}
 3. DO NOT RESUSCITATE (dnr) {0BD8A134-C563-7E00-A6B1-F7C543CEF6FF}
 4. ALLERGY (a) {EBF87F6F-1BDB-6147-C444-2FF787DEAC7E}
 5. HIGH FALL RISK (hfr) {276C15F3-5BED-0FE0-5649-337A23AB1386}
 6. ELOPEMENT (e) {EB030907-FB90-CA9D-CF50-E26C5D9FE671}
 7. PERSONAL MEDICAL ALERT BRACELET (pmab)
{E63D5FDE-4509-5751-A11C-C023E296858A}
 8. BLOOD TRANSFUSION IV (BI TRAN IV)
{27A0BE4C-433C-F920-BAC2-066F304BDF59}
 9. WANDERER (w) {C5BD737B-F9A6-5800-B45C-2585507C3EDC}

UNITS:

=====
TERM: ARTERIAL LINE ALARMS ON
ID: {9267950F-51B3-4C38-B280-C62D7B076F6B}

Page: 39

=====
Abbreviation..... ALINE ALAR
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... Arterial Lines Alarms On
Display Name..... Arterial Line Alarms On
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

-
1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: ARTERIAL LINE CIRCULATION ASSESSED
ID: {CFB536B3-C087-48EF-9417-94CBE1B88911}

=====

Page: 40

Abbreviation..... ART L CR A
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... Has The Arterial Line Circulation
 Been Assessed
Display Name..... Arterial Line Circulation Assessed
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

-
- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
 - 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: ARTERIAL LINE CONDITION
ID: {A87DEFA0-FD22-8AAB-AEEE-61E91FEE2154}

=====

Page: 41

Abbreviation..... ART LI CON
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Arterial Line Condition
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

-
- 1. PATENT CLEAN (PAT CLEAN) {68688246-AE58-2EFB-BF9D-5E00D8CFDA8E}
 - 2. RED (RED) {515D7CFB-3527-30E7-EC4B-B607C272DFA0}
 - 3. SWOLLEN (SWOLLEN) {DCE86B10-C393-8DA9-6518-F2943F0E1A4E}
 - 4. DISCOMFORT (DISCOMFORT) {77D058C9-1EB4-96C0-C9E3-666C840B884C}

LOCATION: 8 entries

-
- 1. L RADIAL (L RAD) {302F78AD-4067-8B6B-BD15-993EC8561AE6}
 - 2. R RADIAL (R RADIAL) {BA10E6E5-DDF4-345E-68AF-21D32E9D972D}
 - 3. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
 - 4. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
 - 5. L AXILLARY (L AXILLARY) {21AF62CA-BEED-2A96-2CD6-F2350A01995E}

- 6. R AXILLARY (R AXILLARY) {3044F28B-FED2-4069-7E3B-2A691C5771F9}
- 7. L BRACHIAL (LBRA) {75F6CE33-1AB3-A1DA-70B8-04825290B577}
- 8. R BRACHIAL (RBRA) {A88E69A8-F6AE-80B0-2CD4-687B121F9DDE}

UNITS:

=====

TERM: ARTERIAL LINE DRESSING CHANGE
ID: {1D4E53DD-2D96-4B9F-9912-920B96FFF942}

=====

Page: 42

Abbreviation..... ART LIN DC
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Arterial Line Dressing Change
Display Name..... Arterial Line Dressing Change
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

- 1. DRY (d) {D8F4F0FD-1B19-2515-3815-EC88D175B305}
- 2. INTACT (i) {935C9CAF-6C88-6610-C6E9-4C1B822923AF}
- 3. CHANGED (CHANGED) {74B14023-D74D-41F6-AA48-B0E5137623CD}
- 4. REINFORCED (REINFORCED) {EC524C50-0CDE-4AA1-8059-E4323850575E}

UNITS:

=====

TERM: ARTERIAL LINE DRESSING TYPE
ID: {A2F0D86E-286A-0D33-2BAA-3DDE78121050}

=====

Page: 43

Abbreviation..... A LINE DRE
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Arterial Line Dressing
Help Text..... Type of dressing used on arterial line
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

- 1. TRANSPARENT (TRANSPAREN) {0E73D24D-0EFC-08F9-196A-D30DD1BB3228}
- 2. PRESSURE DRESSING (PRESS DRES)
{6D9DB9EF-85EB-8321-03FC-879AA4C5615B}

3. GAUZE AND TAPE (GAU & TAPE)
{5BC40405-B5C4-47CF-81C3-520732E100AC}

UNITS:

=====

TERM: ARTERIAL LINE SITE CARE Page: 44
ID: {2CCDDAA0-9563-8B09-11C1-51E1E2823A48}

=====

Abbreviation..... A LINE CAR
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... Has Arterial Line Site Care Been Done?
Display Name..... Arterial Line Site Care
Help Text..... Date arterial line inserted
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

-
- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
 - 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: ARTERIAL LINE SIZE Page: 45
ID: {9D96EBB7-19E0-4832-B61F-D2A42CE15737}

=====

Abbreviation..... ART LIN SI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Arterial Line in size Gauge
Display Name..... Arterial Line Size
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

LOCATION: 8 entries

-
- 1. L RADIAL (L RAD) {302F78AD-4067-8B6B-BD15-993EC8561AE6}
 - 2. R RADIAL (R RADIAL) {BA10E6E5-DDF4-345E-68AF-21D32E9D972D}
 - 3. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
 - 4. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
 - 5. L AXILLARY (L AXILLARY) {21AF62CA-BEED-2A96-2CD6-F2350A01995E}
 - 6. R AXILLARY (R AXILLARY) {3044F28B-FED2-4069-7E3B-2A691C5771F9}

7. L BRACHIAL (LBRA) {75F6CE33-1AB3-A1DA-70B8-04825290B577}
8. R BRACHIAL (RBRA) {A88E69A8-F6AE-80B0-2CD4-687B121F9DDE}

UNITS:

GAUGE 14 > 18 * 26 < 30 Dec Prec: 0

=====

TERM: ARTERIAL LINE WAVEFORM Page: 46
ID: {1B11B2EF-C905-40F6-9F6F-195FEC0CD8E6}

=====

Abbreviation..... ALINE WAVE
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Arterial Line Waveform
Display Name..... Arterial Line Waveform
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

-
1. NORMAL (n) {553C9A97-A3D5-C52F-0E4E-C7D168C99A2B}
 2. DAMPENED (DAMPENED) {77C777B3-7E56-41E3-8E70-DE144E6A3424}
 3. HYPERDYNAMIC (HYPERDYNAM) {818A8A99-04D9-4989-A9BA-ECFACED74605}
 4. POSITIONAL (POSITIONAL) {9A74D16B-5E95-4EA8-88AF-E53C80311A5D}
 5. NOT MONITORED (NOT MONITO)
{0E6A8A07-48E0-428F-8D86-22572CABCF48}

UNITS:

=====

TERM: ARTERIAL LINE ZEROED Page: 47
ID: {A500DBBF-6E69-4CB3-926F-6B9988BB8B68}

=====

Abbreviation..... ALINE ZERO
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... Arterial Line Zeroed
Display Name..... Arterial Line Zeroed
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

-
1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}

2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: BABINSKI Page: 48
ID: {D7BB9690-82D8-F50B-93AA-1937D051AA53}

=====

Abbreviation..... BABINSKI
Active..... Yes
Boolean Value False..... Negative
Boolean Value True..... Positive
Data Type..... Boolean
Description..... Babinski Reflex Test
Display Name..... Babinski
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: BED POSITION Page: 49
ID: {D65573FE-6CAB-3A81-4572-ACC9488D8BAF}

=====

Abbreviation..... BEDPOS
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Bed Position
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

1. HIGH (h) {D88CB5EC-2AD4-2F80-2D43-40ECB0B7D390}
2. LOW (l) {B6CE52A4-C1DD-FDE3-F87E-8E0C98FCC4AC}
3. SEMIFOWLERS HIGH (sh) {1BBE7175-2BD2-BFB5-1E1B-E02443433072}
4. SEMIFOWLERS LOW (sl) {A265AF02-53FC-B182-8AE3-A06E95718D7E}
5. CARDIAC CHAIR (cc) {D4A2CFC1-D3A4-696B-4062-AFB2DAFFD76A}

UNITS:

=====

TERM: BEHAVIOR

Page: 50

ID: {BBF17C69-B1B4-37D9-0469-D124F4ADA673}

=====

Abbreviation..... BEHAVIOR
Active..... Yes
Boolean Value False..... Inappropriate
Boolean Value True..... Appropriate
Data Type..... Boolean
Description..... Is Patients Behavior Appropriate?
Display Name..... Behavior
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: BEHAVIOR CURRENT LOCATION

Page: 51

ID: {A3BE725A-A4D7-6C75-D5CD-453755C5EF21}

=====

Abbreviation..... CUR LOC
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Behavior Current Location
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 10 entries

1. PATIENT'S ROOM (pr) {0F13E11D-9F78-14A3-0568-AF26393F39F1}
2. BATHROOM (b) {D20D91A7-E9F5-D000-5E0B-772BA9F5F72E}
3. HALL (h) {BFCEA10F-2663-199E-88D8-E263122D7D26}
4. DAY ROOM (dr) {269125EF-FA82-B17A-6B56-0CEAB4788899}
5. QUIET ROOM (qr) {CCD96D41-C273-F6F3-DA18-8818B1CC9980}
6. GROUP ROOM (gr) {D41C3D9F-8F10-1A0F-D571-20DE0EA36328}
7. SMOKING PORCH (sp) {5AD3F05B-2A5B-E2B4-FD47-66712CFA44E6}
8. SMOKING AREA (sa) {8AFDDBC0-9FC8-523F-5D07-4D8CCBA7036F}
9. MEDICATION ROOM (mr) {9571AA9F-F538-55EC-37F6-99FE5677B970}
10. NURSING STATION (ns) {B6E48F92-B9BA-488F-3A2E-D74568A6B855}

UNITS:

=====

TERM: BEHAVIOR OBSERVED

Page: 52

ID: {E5D2F4A5-5484-2B9C-EBD0-C0EB0FBC6E8C}

=====

Abbreviation..... BEH OBS
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Behavior
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 15 entries

-
1. COOPERATIVE (c) {C176406E-EBDE-4ADE-53DE-CF191D8AFC85}
 2. CALM (c) {E206C3F5-6977-302C-51FB-AF2155B5FFD9}
 3. ANXIOUS (a) {381BDB41-7E8A-6FDB-841B-2ECC1889B9EF}
 4. HOSTILE (h) {29EF2C04-7DB0-4491-2B5E-D92A6A558721}
 5. COMBATIVE (c) {DF65D42D-6CAA-657E-AEE1-C070C1FDF177}
 6. INEFFECTIVE COPING (ic) {704B1A55-B639-60F4-E3A5-16F942ABC259}
 7. WITHDRAWN (w) {E74B0F04-0E3F-632F-9409-C4953D95B1BA}
 8. UNCOOPERATIVE (u) {2D737EFB-F4CA-7AEA-9EB2-6FF133742815}
 9. SLEEPS POORLY (sp) {AB9424B4-D0B1-148B-C800-FF535F3FAC84}
 10. AGITATED (a) {BAD2CCCE-9083-6F61-1196-C0C9A8DDBFFC}
 11. SEDATED (s) {29D9ADD6-2C1F-F6D5-F4D0-B2BB86416997}
 12. RESTLESS (r) {8FCCAD26-DBA0-7067-DA72-30D342BFA8E8}
 13. ALERT (a) {BC7E2630-179A-5C80-E6C7-598D42F0BFD6}
 14. DROWSY (d) {6B350814-A88A-0F05-C382-9F27830FBC09}
 15. LETHARGIC (l) {4ECA2E44-093A-9354-848E-CDC72B165072}

UNITS:

=====

TERM: BEHAVIOR PRECAUTION LEVEL

Page: 53

ID: {1F322720-9BDF-A3A4-47A6-7DA602FFBBF4}

=====

Abbreviation..... PRE LVL
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Behavior Precaution Level
Help Text..... No Help Text
Multi Select Picklist..... No

Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

- 1. CLOSE OBSERVATION (co) {2679BCAC-8753-FE73-56F8-3AB632161AF4}
- 2. 1:1 OBSERVATION (1o) {B1A7AF4B-9BD2-5F48-79B3-F35861AA5C45}
- 3. PATIENT STAYS WITHIN LINE OF SITE OF NURSING STATION (pswlosons)
{CE3E3531-7D9B-A2B4-00B3-00EBCDD09D3D}
- 4. PATIENT SIGNS SAFETY CONTRACT (PT CONTRAC)
{C65734EC-2CAF-4547-8FD8-CBCD9C706003}

UNITS:

=====

TERM: BEHAVIOR PRECAUTION REASON
ID: {C2808159-80A1-AF7F-7C7F-C3A96096187E}

=====

Page: 54

Abbreviation..... PRE REA
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Behavior Precaution Level
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

- 1. SUICIDAL (s) {B104E9C7-A78E-82DC-49A4-19B94A639AFE}
- 2. ESCAPE (e) {965F3201-9D63-66AB-FA25-99479AE7145A}
- 3. ELOPEMENT (e) {EB030907-FB90-CA9D-CF50-E26C5D9FE671}
- 4. ASSAULTIVE (a) {99ECD8C7-191A-3B9C-2A7A-803A921065FA}
- 5. MEDICAL DEVICE SAFETY (mds)
{1985CABE-4F53-9FCF-F444-F9AC53C04DBD}

UNITS:

=====

TERM: BIPAP-EXPIRATORY
ID: {B52B46E7-A947-447E-A277-71B286B6313A}

=====

Page: 55

Abbreviation..... BIPAP EXP
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Bipap machine setting for expiratory
pressure

UNITS:

=====

TERM: BLOOD PRESSURE

Page: 58

ID: {B15F2DF6-CE99-B847-FE6B-3D5F174A2BCD}

=====

Abbreviation..... BP
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Complex
Description..... Complex observation of Systolic and
Diastolic pressures.
Display Name..... Blood Pressure
Help Text..... Enter in the blood pressure values as
###/###.
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid..... 4500634

LOCATION: 7 entries

1. BRACHIAL (Bra) {67426127-09E1-E566-F31E-15714A5B2B7C}
2. FEMORAL (Fem) {E3442C9E-8677-7296-D314-7365319C3716}
3. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
4. L LEG (LL) {B088C809-E528-B785-21C2-14E45E1FD8AF}
5. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
6. R LEG (RL) {C5E40A7B-6FC4-267D-8CE8-3571D21D81FD}
7. RADIAL (RADIAL) {D52C4103-2AB3-A0C7-16E3-CF95A1B3B789}

METHOD: 11 entries

1. AFTER EXERCISE (AE) {D129EF4B-3A99-9966-54F1-2C01A36814EE}
2. ARTERIAL LINE (AL) {D4B3F4ED-8B33-431E-4164-BD165E823776}
3. AT REST (AR) {E8C289A6-0C3D-A7B2-BB85-7F1A23F22B50}
4. CUFF-AUTOMATED (CA) {088AAD2B-2E7B-F835-376F-3D4340790976}
5. CUFF-MANUAL (CM) {5F56C163-2B04-112F-55DE-E31374E0FB58}
6. DOPPLER (Dop) {52683C19-F9DC-876B-CA3F-53E7C864CF04}
7. INVASIVE (I) {5DB8746E-5446-F47A-883C-A44B91275A6E}
8. MONITOR (Mo) {63FC24EA-50B8-B302-D948-ABE3517CF6FC}
9. NON-INVASIVE (N-I) {F50ABFED-5FD7-FC95-20FE-1E3038C9F299}
10. PALPATED (Pal) {295D8E72-12A1-08AF-6BE8-0F9CD238F5EF}
11. WITH ACTIVITY (Act) {544A262B-EA34-F82A-46A9-F80740515D65}

POSITION: 6 entries

1. LYING (Ly) {DE28D2BF-6F8D-64C3-0092-4CEB5E0902B9}
2. SEMIFOWLERS (Sfl) {C2F6A945-3829-AABA-3D67-54D1829D3EB9}
3. SITTING (Sit) {DD4DCDF1-0A22-2CC5-404F-0E472C73961C}
4. STANDING (St) {0776AE42-C370-B7B3-9B4D-9C914BA8F231}
5. TRENDELENBURG (Tr) {42F3353B-7041-EEE1-A87E-941D3BE1CBCD}
6. SUPINE (SUPINE) {DEA8A744-B445-4A9E-875D-485E060F2ACC}

UNITS:

=====

TERM: BLOOD PRESSURE CUFF SIZE
ID: {46DCF900-6F2C-46D0-ADB8-2EC17E7659FF}

=====

Page: 59

Abbreviation..... BP CUFF SI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Blood Pressure Cuff Size
Help Text..... Blood Pressure Cuff Size
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

-
- 1. ADULT (ADULT) {37B557BE-B932-4757-8E08-5EA315E8AD1F}
 - 2. ADULT LARGE (AD LARGE) {38E64D72-2806-47FA-B82F-51844574F850}
 - 3. ADULT SMALL (AD SMALL) {0161A4E5-713A-406B-B3FD-9ED13BEB5366}
 - 4. THIGH (THIGH) {65AE788B-52F5-44D3-B9C1-49DAE46FF5AB}
 - 5. PEDIATRIC (PEDIATRIC) {06865E22-3ED9-4EDD-9749-6F8641BF29DE}

UNITS:

=====

TERM: BODY SURFACE AREA
ID: {5C9A522B-C5ED-4D59-7BF8-3D88ADEBAB07}

=====

Page: 60

Abbreviation..... M2
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Body Surface Area
Display Name..... BSA
Help Text..... Body Surface Area
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

SQUARE METERS 0.50 > 1.00 * 1.75 < 4.00 Dec Prec: 2

=====

TERM: BOWEL MOVEMENT EVENT(S)
ID: {2EB1DAFB-1E56-40A8-B371-D46EF7AF2435}

=====

Page: 61

=====

TERM: BRADEN SCALE FRICTION AND SHEAR
ID: {03ECD9F4-3C12-A62B-FAAC-0F8A824763AA}

=====

Page: 63

Abbreviation..... BR SC FRIC
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Braden Scale Friction and Shear
Display Name..... Braden Scale Friction and Shear
Help Text..... 1. Problem: Requires moderate to
maximum assistance in moving. 2
Potential Problem: Moves feebly or
requires minimum assistance.3 No
Apparent Problem: Moves in bed and
in chair independently and has
sufficient muscle strength to lift up
completely
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

-
1. BRADEN FRICTION AND SHEAR 1 PROBLEM (BS F/S 1)
{AF8D2952-E0DA-411E-91CC-0114F4EB1C33}
 2. BRADEN SCALE FRICTION AND SHEAR 2 POTENTIAL PROBLEM (BS F/S 2)
{B56338A9-8956-4D27-8BA0-DA50E07359E0}
 3. BRADEN SCALE FRICTION AND SHEAR 3 NO APPARENT PROBLEM (BS F/S 3)
{9850E2BF-6600-4E66-8FCD-1517D903BC0F}

UNITS:

=====

TERM: BRADEN SCALE MOBILITY
ID: {85F9FB55-26A3-6203-C1AE-6EFED5952B3D}

=====

Page: 64

Abbreviation..... BR SC MOBI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Braden Scale Mobility - ability to
change and control body position
Display Name..... Braden Scale Mobility
Help Text..... 1=Completely Immobile, 2=Very
Limited, 3=Slightly Limited, 4=No
Impairment
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

Description..... Braden Scale Nutrition - usual food
intake pattern
Display Name..... Braden Scale Nutrition
Help Text..... 1=Very Poor, 2=Probably Inadequate,
3=Adequate, 4=Excellent
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

-
1. BRADEN SCALE NUTRITION 1 VERY POOR (BS NUT 1)
{28A86074-A6B2-455E-ABA4-19195C2E3679}
 2. BRADEN SCALE NUTRITION 2 PROBABLY INADEQUATE (BS NUT 2)
{7DAB22F6-22E4-4A8B-A6A5-CD1277EB6B87}
 3. BRADEN SCALE NUTRITION 3 ADEQUATE (BS NUT 3)
{3A8A9C59-8DC5-49A9-BAC9-882CB75D8FD5}
 4. BRADEN SCALE NUTRITION 4 EXCELLENT (BS NUT 4)
{F09E4D35-1AA8-49A5-9B38-FBFFAFCFC181}

UNITS:

=====

TERM: BRADEN SCALE SENSORY PERCEPTION
ID: {BC5D5F00-85EE-7FB2-3175-8EDB95485D5E}

=====

Page: 67

Abbreviation..... BR SC SENS
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Sensory Perception - ability to
respond meaningfully
Display Name..... Braden Score - Sensory Perception
Help Text..... 1=Completely Limited, 2=Very Limited,
3=Slightly Limited, 4=No Impairment
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

-
1. BRADEN SCALE SENSORY PERCEPTION 1 COMPLETELY LIMITED (BS SENC 1)
{60F646B3-790F-4C5A-B2A0-845915E120E6}
 2. BRADEN SCALE SENSORY PERCEPTION 2 VERY LIMITED (BS SENC 2)
{98362E63-7245-4E7F-8C30-896A9B923199}
 3. BRADEN SCALE SENSORY PERCEPTION 3 SLIGHTLY LIMITED (BS SENC 3)
{1D01D47D-CDC2-40EE-837B-5417D6EE47D6}
 4. BRADEN SCALE SENSORY PERCEPTION 4 NO IMPAIRMENT (BS SENC 4)
{10304656-E5E8-4F32-8358-A9F80511EFDA}

UNITS:

5. R UPPER LOBE (RUL) {098C11B4-B556-FBBE-7E21-AEC9D7E4408A}

UNITS:

=====

TERM: CALL LIGHT WITHIN REACH Page: 70
ID: {ABAC4D24-C449-117E-BE4C-BFBC1F599A57}

=====

Abbreviation..... CALL LT
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... No Description
Display Name..... Call Light Within Reach
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

-
- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
 - 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: CAPILLARY REFILL Page: 71
ID: {951BBDAA-883A-E4F7-AD05-B086A71A2FEC}

=====

Abbreviation..... CAP REFILL
Active..... Yes
Boolean Value False..... Greater than 3 seconds
Boolean Value True..... Less than 3 seconds
Data Type..... Boolean
Description..... RATE OF CAPILLARY REFILL
Display Name..... Capillary Refill
Help Text..... RATE OF CAPILLARY REFILL
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

-
- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
 - 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

Display Name..... Bruit
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. PRESENT (p) {384DC20F-1625-CE1E-1C85-A3B3924BB2A0}
- 2. ABSENT (a) {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0}

LOCATION: 7 entries

- 1. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
- 2. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
- 3. L CAROTID (L CAROTID) {FCACD1EB-C12B-2515-4A68-4E9E28A72DBC}
- 4. R CAROTID (R CAROTID) {8280ACD2-72D6-202A-89AF-86BCA7736427}
- 5. R TEMPORAL (R TEMPOR) {1F50461D-5ECF-25B2-A29A-E02803CB7814}
- 6. L TEMPORAL (L TEMPO) {375E65CF-A59E-A1E9-EB6D-D53CD818939E}
- 7. ABDOMINAL AORTA (ABD AORTA)
{C4BD42DD-1D6B-4B2C-B6D5-19E1185E06E3}

UNITS:

TERM: CARDIAC ECTOPY Page: 75
ID: {74EC179C-8598-13F7-3D58-442305A3AFDB}

Abbreviation..... CARD E
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Cardiac Ectopy
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

- 1. PJC (p) {1D2D1FEA-2006-3C6C-FAD4-72384F2033C9}
- 2. PAC (p) {848AB18A-598D-8E4E-8B1B-ED7478CE340B}
- 3. PVC (p) {33FF92F1-84CC-7CFA-41CF-B271A26C59B8}
- 4. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}

UNITS:

TERM: CARDIAC HEART SOUNDS Page: 76
ID: {476105AC-C83C-ECEC-6C72-7C7C8EE87853}

=====
Abbreviation..... CARD HS
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Cardiac Heart Sounds
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 7 entries

-
- | | |
|--------------------|--|
| 1. S1 (s) | {1AB9DF0A-8238-E8FB-8520-4749E13DDE1D} |
| 2. S2 (s) | {57C043D3-B227-3E7A-C1E6-5129A2800718} |
| 3. S3 (s) | {9374B57F-DE1F-B9C4-26B5-2DD6B194DC26} |
| 4. S4 (s) | {9C6CA302-85B6-8FFC-EAB4-8FED62C535A5} |
| 5. RUB (r) | {0040835D-1CF0-7D55-EEBB-A94A1E0D1206} |
| 6. CLICK (c) | {F1F3A1A8-35E1-61B4-C783-3314C8247380} |
| 7. MURMUR (MURMUR) | {591058E0-4B71-E432-B7D7-EBD42A8C4225} |

UNITS:

=====
TERM: CARDIAC INDEX Page: 77
ID: {BFC29A9A-C666-68B5-C5CF-36B809CF66C2}

=====
Abbreviation..... CI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Cardiac Index
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

-
- | | |
|-----------------|--|
| 1. MONITOR (Mo) | {63FC24EA-50B8-B302-D948-ABE3517CF6FC} |
|-----------------|--|

UNITS:

LITERS PER MINUTE PER METER2 (BSA)
-0.1 > 0.0 * 8.0 < 8.1 Dec Prec: 1

=====
TERM: CARDIAC OUTPUT Page: 78
ID: {0B7C4991-12E5-10DE-F25D-CB81BFA1ADAD}

Abbreviation..... CO
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Cardiac Output
Display Name..... Cardiac Output
Help Text..... Cardiac Output
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

1. MONITOR (Mo) {63FC24EA-50B8-B302-D948-ABE3517CF6FC}

UNITS:

LITERS PER MINUTE 0.00 > 4.00 * 8.00 < 16.00
Dec Prec: 2

TERM: CARDIAC RHYTHM Page: 79
ID: {D5CF1A63-551A-59D4-E6BB-C851D81BDB58}

Abbreviation..... CARDR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... HEART RHYTHM
Display Name..... Cardiac Rhythm
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 27 entries

1. SINUS RHYTHM (sr) {A34F80EA-3D6A-076A-DB2A-E4E5C5AACA11}
2. SINUS BRADYCARDIA (sb) {D6EF4940-63AE-22AF-1A80-D8E6F0F07839}
3. AGONAL (a) {1621B411-33E3-3A6A-78A6-40B7BDC3E8CB}
4. ASYSTOLE (a) {6E87ECBC-3160-2CEF-D2FF-EEBD1DAA96F2}
5. ATRIAL FIB/FLUTTER (af) {9C1BC2F6-3BF0-BEBB-BE87-50E9F7984544}
6. ATRIAL FIBRILLATION (a fib)
{7B5EBAF2-55B6-2905-3D89-B28CE6A509CB}
7. ATRIAL FLUTTER (af) {9CB8D413-6762-D78E-F40E-25DB66CD23B3}
8. L BUNDLE BRANCH BLOCK (lbbb)
{C3F2E19E-495B-7BF6-99EF-B8236D6B9C3D}
9. R BUNDLE BRANCH BLOCK (rbbb)
{FDE0D174-5220-F069-D284-20399D9489B0}
10. IDIOVENTRICULAR (i) {8BD08966-F536-8E74-0BDA-A95D0D7F8B6C}
11. JUNCTIONAL (j) {4A339D84-7051-1258-877D-85AFC1CA20FE}
12. PACED: AV SEQUENTIAL (pas)

- {3C485674-224C-47DA-510D-92F1A0CA54A1}
- 13. PACED: ATRIAL (pa) {88837DD4-098D-B77F-A754-89D5C4F5870F}
- 14. PACED: VENTRICULAR (pv) {A63F6E6E-A7C4-0928-81FB-FE2F70F713C6}
- 15. PAROXYSMAL ATRIAL TACHYCARDIA (pat)
{27BC1FAF-A403-6EA8-044C-F1271F12A62A}
- 16. SINUS ARREST (sa) {3CD8F7AC-3FE0-F342-089C-787232D2A5C2}
- 17. SINUS ARRHYTHMIA (sa) {07A24AC9-B486-A4B8-BADB-04D2729542CE}
- 18. SUPRAVENTRICULAR TACHYCARDIA (st)
{F712412F-D39A-58F7-E4F1-6AC8B77445D7}
- 19. VENTRICULAR FIBRILLATION (vf)
{400B70C2-6E15-458F-70C6-FA2597DE2556}
- 20. VENTRICULAR TACHYCARDIA (vt)
{138258AD-3A2C-2A6C-9DF8-615CCC303DC2}
- 21. FIRST DEGREE AV BLOCK (fdab)
{4FF046C9-B638-F86C-5D51-006B3D5225E3}
- 22. SECOND DEGREE AV BLOCK: TYPE 1 (sdabt1)
{508BF663-96FB-EEBE-4F06-D2F1A02BD099}
- 23. SECOND DEGREE AV BLOCK: TYPE 2 (sdabt2)
{3367754A-08D1-31A9-E850-1B36358E358A}
- 24. THIRD DEGREE AV BLOCK (tdab)
{1AC7E31B-830E-2309-3AF7-8FB03003DD86}
- 25. SINUS PAUSE (sp) {7F8F89A9-6B16-B159-55B3-DEFDF5B4B4F6}
- 26. SINUS BLOCK (sb) {94D6EFED-3D08-E724-BC40-9C398919A142}
- 27. SINUS TACHYCARDIA (SINUS T)
{0AC41AF3-7656-A35D-E397-21D35405812A}

UNITS:

=====

TERM: CARDIAC RHYTHM LEAD

Page: 80

ID: {3F4A283A-6A9F-345A-1BE0-86980228FA6D}

=====

Abbreviation..... CAR RH LD
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... No Description
 Display Name..... Cardiac Rhythm Lead
 Help Text..... No Help Text
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 12 entries

- 1. I (I) {BEEDFF9A-9401-4D2A-8345-126D8CFFC8D8}
- 2. II (II) {F165C13B-397E-4182-8D72-BD1D74CFBA00}
- 3. III (III) {D8FE8FA9-9BB5-416E-AE2C-60E627F76E5A}
- 4. AVR (AVR) {04AF2551-1185-4BA9-AEED-977D562D188A}
- 5. AVL (AVL) {9124A193-5F50-448D-9E60-C5869EC4029C}
- 6. AVF (AVF) {C8E97206-58AA-47C0-BD2E-7A30E56DB46A}
- 7. V1 (V1) {4CAA3A92-F3DB-4847-8523-380EFEE6CFD5}
- 8. V2 (V2) {05A3232E-A30D-4771-A1B5-E7705B717646}

- 9. V3 (V3) {E150126E-4AA7-4F96-961E-FD1169C77158}
- 10. V4 (V4) {255260F9-5CC7-41B1-831C-F9319FA27506}
- 11. V5 (V5) {D5C10744-0EB5-46FA-AE6C-CCF571318F8A}
- 12. V6 (V6) {B86F72D7-7884-4470-92A8-EF9277C5B584}

UNITS:

=====

TERM: CARDIAC THRILL Page: 81
 ID: {03A59127-CB5C-4DA6-91C4-5AA02DDE2C66}

=====

Abbreviation..... THRILL
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... No Description
 Display Name..... Thrill
 Help Text..... No Help Text
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 2 entries

- 1. PRESENT (p) {384DC20F-1625-CE1E-1C85-A3B3924BB2A0}
- 2. ABSENT (a) {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0}

LOCATION: 8 entries

- 1. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
- 2. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
- 3. R CAROTID (R CAROTID) {8280ACD2-72D6-202A-89AF-86BCA7736427}
- 4. L CAROTID (L CAROTID) {FCACD1EB-C12B-2515-4A68-4E9E28A72DBC}
- 5. R TEMPORAL (R TEMPOR) {1F50461D-5ECF-25B2-A29A-E02803CB7814}
- 6. L TEMPORAL (L TEMPO) {375E65CF-A59E-A1E9-EB6D-D53CD818939E}
- 7. ABDOMINAL AORTA (ABD AORTA)
 {C4BD42DD-1D6B-4B2C-B6D5-19E1185E06E3}
- 8. HEART (HEART) {F7D48A36-B512-42DD-B92F-07ACF2B95E08}

UNITS:

=====

TERM: CARDIAC THRILL TIMING Page: 82
 ID: {D52A1728-1692-4742-886A-B119BB7EE598}

=====

Abbreviation..... THRILL TIM
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... No Description

Display Name..... Thrill Timing
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

- 1. CONTINUOUS (CONTINUOUS) {676CD2DA-F882-4BE5-99BB-510EA7669749}
- 2. SYSTOLIC (SYSTOLIC) {859B5700-0990-42FF-8129-8DB24EFDA8D1}
- 3. DIASTOLIC (DIASTOLIC) {3E9D1DCF-6B3B-482C-AE5F-442814F4B5C0}
- 4. PRESYSTOLIC (PRESYSTOLI) {BB85B8EE-36F8-4AC2-9AB0-0BABCFA9FCBB}

UNITS:

=====

TERM: CENTRAL LINE ALARMS ON
ID: {152A037A-A6CC-4F12-8B85-25CEE78772C4}

=====

Page: 83

Abbreviation..... C-LINE ALA
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... Central Line Alarms On
Display Name..... Central Line Alarms On
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
- 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: CENTRAL LINE CHANGED
ID: {8FA1B26D-CADB-4769-3574-18A6BD86C400}

=====

Page: 84

Abbreviation..... C LINE CHA
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Central Line Changed
Help Text..... Central Line changed to
Multi Select Picklist..... No
Term Type..... OBSERVATION

Vuid.....

DATA VALUE: 4 entries

- 1. NEW SITE (NEW SITE) {10B1FF61-E58E-4608-9CCB-9232BC964418}
- 2. OVER WIRE (OVER WIRE) {157CD456-1020-43AE-9AAA-1686F1574B6F}
- 3. INSERTED THROUGH EXISTING INTRODUCER (INS EXS IN)
{20BC183A-197F-4464-A27C-47F63C5B842A}
- 4. INSERTED THROUGH NEW INTRODUCER (INS NEW IN)
{7FCA403D-41E5-4BAA-B155-B58F35EA49B3}

UNITS:

=====

TERM: CENTRAL LINE CONDITION Page: 85
 ID: {7996453B-0F7C-B4E8-2007-6C26AE26C479}

=====

Abbreviation..... C LINE CON
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... No Description
 Display Name..... Central Line Condition
 Help Text..... Condition of central line
 Multi Select Picklist..... Yes
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 4 entries

- 1. PATENT CLEAN (PAT CLEAN) {68688246-AE58-2EFB-BF9D-5E00D8CFDA8E}
- 2. RED (RED) {515D7CFB-3527-30E7-EC4B-B607C272DFA0}
- 3. SWOLLEN (SWOLLEN) {DCE86B10-C393-8DA9-6518-F2943F0E1A4E}
- 4. DISCOMFORT (DISCOMFORT) {77D058C9-1EB4-96C0-C9E3-666C840B884C}

UNITS:

=====

TERM: CENTRAL LINE DRESSING CHANGE Page: 86
 ID: {DDEA3089-9246-40F3-9A67-56A05B393309}

=====

Abbreviation..... C LI DR CH
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... Central Line Dressing Change
 Display Name..... Central Line Dressing Change
 Help Text..... No Help Text
 Multi Select Picklist..... Yes
 Term Type..... OBSERVATION

Vuid.....

DATA VALUE: 4 entries

- 1. DRY (d) {D8F4F0FD-1B19-2515-3815-EC88D175B305}
- 2. INTACT (i) {935C9CAF-6C88-6610-C6E9-4C1B822923AF}
- 3. CHANGED (CHANGED) {74B14023-D74D-41F6-AA48-B0E5137623CD}
- 4. REINFORCED (REINFORCED) {EC524C50-0CDE-4AA1-8059-E4323850575E}

UNITS:

TERM: CENTRAL LINE DRESSING TYPE

Page: 87

ID: {AF276338-5F2C-1CEF-F009-AC4BAED0DDE3}

Abbreviation..... C LINE DRE
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... No Description
 Display Name..... Central Line Dressing Type
 Help Text..... Type of dressing on central line
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 3 entries

- 1. ANTI-MICROBIAL DRESSING (AB DRESSIN) {9734C5E9-D099-47EB-93B1-D82201F64331}
- 2. PRESSURE DRESSING (PRESS DRES) {6D9DB9EF-85EB-8321-03FC-879AA4C5615B}
- 3. TRANSPARENT (TRANSPAREN) {0E73D24D-0EFC-08F9-196A-D30DD1BB3228}

UNITS:

TERM: CENTRAL LINE INSERTION MARK

Page: 88

ID: {FAF85282-EB13-405A-B78E-36767EE4C33B}

Abbreviation..... C LINE INS
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... Numeric
 Description..... No Description
 Display Name..... Central Line Insertion Mark
 Help Text..... No Help Text
 Multi Select Picklist.....
 Term Type..... OBSERVATION
 Vuid.....

- 9. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
- 10. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
- 11. L ANTECUBITAL (L ANTECUBI)
{C21838CC-BAB5-1D7F-E2A4-A7C663FB28C3}
- 12. R ANTECUBITAL (R ANTECUBI)
{7FB00A65-B4B2-06A6-BD89-98AC112D8D4E}
- 13. L UPPER ARM (L UPP ARM) {5F915696-1B08-32E9-6B95-49246117D660}
- 14. R UPPER ARM (R UPP ARM) {3EF20E8E-F68C-E4F2-8F8C-FAD3B44DDAA4}
- 15. L WRIST (L WRIST) {7BCAFACE-FF3B-46DA-99F9-6437A90B7126}
- 16. R WRIST (R WRIST) {FDD9FE1D-F8A2-09AD-82D7-AA478693DD01}

UNITS:

=====

TERM: CENTRAL LINE WAVEFORM Page: 93
 ID: {4F847833-9A36-4EAE-BC73-B7AB069D39E4}

=====

Abbreviation..... C-LINE WAV
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... Central Line Waveform
 Display Name..... Central Line Waveform
 Help Text..... No Help Text
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 5 entries

- 1. NORMAL (n) {553C9A97-A3D5-C52F-0E4E-C7D168C99A2B}
- 2. DAMPENED (DAMPENED) {77C777B3-7E56-41E3-8E70-DE144E6A3424}
- 3. HYPERDYNAMIC (HYPERDYNAM) {818A8A99-04D9-4989-A9BA-ECFACED74605}
- 4. POSITIONAL (POSITIONAL) {9A74D16B-5E95-4EA8-88AF-E53C80311A5D}
- 5. NOT MONITORED (NOT MONITO)
{0E6A8A07-48E0-428F-8D86-22572CABCF48}

UNITS:

=====

TERM: CENTRAL LINE ZEROED Page: 94
 ID: {7868021A-72D4-451F-A2E3-E3106DEC52B4}

=====

Abbreviation..... C-LINE ZER
 Active..... Yes
 Boolean Value False..... No
 Boolean Value True..... Yes
 Data Type..... Boolean
 Description..... Central Line Zeroed
 Display Name..... Central Line Zeroed
 Help Text..... No Help Text

Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
- 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

TERM: CENTRAL VENOUS PRESSURE
ID: {D30F98A7-4C5D-12E8-AB4D-9C85A4332EC3}

Page: 95

Abbreviation..... CVP
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Central Venous Pressure
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid..... 4688719

LOCATION: 8 entries

- 1. L SUBCLAVIAN (LSCL) {71B93594-15F9-AA4A-AEF2-E222EC28C00D}
- 2. R SUBCLAVIAN (R SUBCLAVI) {53DEAE83-93AE-123B-14E5-4DE38CD4F752}
- 3. L INTERNAL JUGULAR (LIJ) {AFBB4FE5-8A11-2047-7900-DADF2FC4430D}
- 4. R INTERNAL JUGULAR (R INT JUGU)
{7ABF6916-D379-198D-6454-F0518A74844A}
- 5. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
- 6. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
- 7. L BRACHIAL (LBRA) {75F6CE33-1AB3-A1DA-70B8-04825290B577}
- 8. R BRACHIAL (RBRA) {A88E69A8-F6AE-80B0-2CD4-687B121F9DDE}

METHOD: 2 entries

- 1. CARDIAC OUTPUT MONITOR (COM)
{CCBE5744-7CB9-9DFA-938F-FDDBC73B45B3}
- 2. MONITOR (Mo) {63FC24EA-50B8-B302-D948-ABE3517CF6FC}

UNITS:

MILLIMETERS OF MERCURY 0 > 2 * 6 < 19 Dec Prec: 0

TERM: CEREBRAL PERFUSION PRESSURE
ID: {166BD1C4-D983-464A-9A6E-B3A7AC7885EB}

Page: 96

Abbreviation..... CPP
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Cerebral Perfusion Pressure
Display Name..... CPP
Help Text..... Cerebral Perfusion Pressure
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

MILLIMETERS OF MERCURY 0 > 50 * 150 < 300 Dec Prec: 0

=====

TERM: CHEST EXPANSION Page: 97
ID: {BAF272C6-1F0A-40E1-918B-3BABC5F22407}

=====

Abbreviation..... CHEST EXPA
Active..... Yes
Boolean Value False..... Asymmetrical
Boolean Value True..... Symmetrical
Data Type..... Boolean
Description..... No Description
Display Name..... Chest Expansion
Help Text..... Symetrical or Asymetrical Chest
Expansion
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

-
- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
 - 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: CHEST PAIN Page: 98
ID: {57AE16AF-8EC2-5800-131F-0AF254C3FA76}

=====

Abbreviation..... CHEST PN
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... EVALUATION OF PATIENTS CHEST DISCOMFORT
Display Name..... Chest Pain
Help Text..... No Help Text
Multi Select Picklist..... No

Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 11 entries

- 1. 0 (0) {2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
- 2. 1 (1) {8F5D10A8-FE17-47BC-DACA-6972D12CC00A}
- 3. 2 (2) {38E22B79-A2B6-3401-B195-C955F24E1278}
- 4. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}
- 5. 4 (4) {D2280C7E-DE50-5A8D-C3FD-5862B69C0E1A}
- 6. 5 (5) {188EDCE7-0726-4234-D8D2-92C8052CD221}
- 7. 6 (6) {640CF2F9-3BBF-AB13-185B-CD02A609FB5A}
- 8. 7 (7) {7A010324-D2C2-5D4B-848A-7B8878666D20}
- 9. 8 (8) {518E7CFC-102E-351E-C7AF-FB575E3C8F4D}
- 10. 9 (9) {58FF2FD4-569F-F454-9F5B-FE6F9E34D8CC}
- 11. 10 (10) {0C827160-18DB-816F-7FAA-A3C175C6D33B}

LOCATION: 3 entries

- 1. L CHEST (L CHEST) {E332E497-FE0C-47CB-A0A7-217FBF52DC97}
- 2. R CHEST (R CHEST) {7E30EC0E-815A-4C2E-A9FC-3C2B2D6E3CD4}
- 3. MID-STERNAL (MID-STERN) {3F5168F1-9082-4FB2-A364-3CDE33B1214B}

UNITS:

=====

TERM: CHEST TUBE DRAINAGE CHARACTERISTICS
ID: {03144029-2E4A-0794-86A5-DD105F707664}

=====

Page: 99

Abbreviation..... DRAIN CHAR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... DESCRIPTION OF PATIENT'S DRAINAGE
Display Name..... Chest Tube Characteristics
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

- 1. SEROUS (SEROUS) {9BDA347E-29E9-80C7-91FA-016192D4257E}
- 2. SEROSANGUINOUS (SEROSANGUI)
{19A8D332-899D-177A-2247-ABCDAAE3BFA8}
- 3. PURULENT (PURULENT) {86C01494-743A-C6A2-8883-01CA5200C79B}
- 4. BLOODY (BLOODY) {D7904E84-8331-D6DD-F2C0-7F3039C5746A}

UNITS:

=====

Abbreviation..... CH TU SU T
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Chest Tube Drainage Suction Type
Display Name..... Chest Tube Drainage Suction Type
Help Text..... Type of Chest Tube Suction
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

-
1. WATERSEAL (WATERSEAL) {6D91B99D-F199-09ED-78BE-B54F2CAC4E06}
 2. GRAVITY (GRAVITY) {C0233F70-376D-41CE-9F14-714C7B35244F}
 3. CONTINUOUS SUCTION (CONT SUCTI)
{8C10DAC2-9AE5-40B9-99C7-25D47F4F13BA}

LOCATION: 12 entries

-
1. L UPPER ANTERIOR CHEST (L U A CHES)
{FE63B58A-DE39-4C03-AE61-E34A5487C767}
 2. R UPPER ANTERIOR CHEST (R U A CHES)
{67890049-A4A8-4241-BDF3-573729BD6EB5}
 3. L MIDDLE ANTERIOR CHEST (L M A CHES)
{36A83F77-117E-474B-AE34-0B2C02E131AE}
 4. R MIDDLE ANTERIOR CHEST (R M A CHES)
{2F751A9C-66CA-4C22-ACC9-6E81FA63FA2B}
 5. L LOWER ANTERIOR CHEST (L L A CHES)
{6073B3C6-AA4D-4BBA-9FF2-32ED40CF92AD}
 6. R LOWER ANTERIOR CHEST (R L A CHES)
{96BCF731-B661-4455-B5E7-C4AE2509505F}
 7. L UPPER POSTERIOR CHEST (L U P CHES)
{6D5674E9-EF04-4EF1-84A0-3F5ED2FB532A}
 8. R UPPER POSTERIOR CHEST (R U P CHES)
{D8F03FA5-1FD7-46EE-8FA4-17B3B7E62141}
 9. L MIDDLE POSTERIOR CHEST (L M P CHES)
{5A656939-F82B-4DEC-9F1D-B5FA644E3DA9}
 10. R MIDDLE POSTERIOR CHEST (R M P CHES)
{2BC487D6-6CE8-499E-94A0-A79EAA8E199B}
 11. L LOWER POSTERIOR CHEST (L L P CHES)
{2C19F10B-7546-4A1F-A5DC-3867359ED971}
 12. R LOWER POSTERIOR CHEST (R L P CHES)
{5163180E-FFD2-4BFB-BF8B-4C6B18EFE5D4}

UNITS:

Abbreviation..... CH TU S RE
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... No Description
Display Name..... Chest Tube Drainage System Replaced
Help Text..... Chest Tube Drainage System Changed
Replaced
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

-
1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

TERM: CHEST TUBE DRESSING CHANGE
ID: {8DE63ADE-D3AC-49EE-BB24-ED44AB9971F5}

Page: 102

Abbreviation..... CH T DR CH
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... Chest Tube Dressing Change
Display Name..... Chest Tube Dressing Change
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

-
1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

TERM: CHEST TUBE SITE ASSESSMENT
ID: {BA53D6FB-9BE9-18F4-ED0A-28AB94CEB115}

Page: 103

Abbreviation..... CH TU ASSE
Active..... Yes
Boolean Value False.....
Boolean Value True.....

Data Type..... PickList
Description..... No Description
Display Name..... Chest Tube Site Assessment
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

- 1. NORMAL (n) {553C9A97-A3D5-C52F-0E4E-C7D168C99A2B}
- 2. AIR LEAK (AIR LEAK) {52E40E25-6D28-53EC-BC05-D1FEAA4616EF}
- 3. FLUCTUATES (Fluctuates) {C2697650-EF53-5F2E-EEE3-7B2E298F2FA0}
- 4. CREPITUS (CREPITUS) {56A97506-05F0-4F37-8FD9-3AA82191408A}
- 5. SQ EMPHYSEMA (SQ EMPHYS.) {B8D03538-0168-4CDC-FF41-BA1185EBCC87}

UNITS:

=====

TERM: CHEST TUBE SUCTION (CM H2O) Page: 104
ID: {BD9724EC-B5BB-4F95-8C24-B0A1451C33A6}

=====

Abbreviation..... CHST TU SU
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Chest Tube Suction (cm H20)
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

CENTIMETERS H2O 0 > 10 * 20 < 30 Dec Prec: 0

=====

TERM: CIRCUMFERENCE/GIRTH Page: 105
ID: {92A124D4-B75F-9FD9-1A51-605887BCEA79}

=====

Abbreviation..... CG
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... circumference
Display Name..... Circumference/Girth
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid..... 4688720

LOCATION: 11 entries

- 1. ABDOMINAL (Abd) {AE1DA31F-2C28-98C2-501A-F3D5ADB2501F}
- 2. ANKLE (Ank) {F2A2CF53-BFA5-9279-DB22-EA2627988E5B}
- 3. CALF (Clf) {7923E496-90D6-B6D5-726A-A5B8DF3436CF}
- 4. HEAD (Hd) {C3209F44-90C6-FF1B-3169-265B294F276A}
- 5. HIP (Hip) {D2C8BD11-7FB2-7217-A663-EBFC8C9428E7}
- 6. KNEE (Kn) {B577308A-4C37-47BD-AD84-38E6342B6DF7}
- 7. LOWER ARM (LrA) {D7627720-D6B0-AC00-8FA3-F7B569ED06FC}
- 8. NECK (Ne) {D4812C3B-AAC2-1B61-679B-3CD8B5CB947C}
- 9. THIGH (THIGH) {2DA93C27-12B0-C76C-85D5-B7D6375F3AD3}
- 10. UPPER ARM (UpA) {B43A2B31-1F35-3F4D-5D23-882541ABE477}
- 11. WRIST (Wri) {C4F1453C-A20A-E48C-0BD8-FAF5E5568ADA}

UNITS:

INCH 29 > 30 * 60 < 61 Dec Prec: 0
CENTIMETER 75 > 76 * 156 < 157 Dec Prec: 0

TERM: CIWA: AGITATION Page: 106
ID: {D80F9E66-A047-400F-A7B1-D701BF019D2A}

Abbreviation..... CIWA AG
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Agitation
Help Text..... CIWA Agitation Observation
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 8 entries

- 1. 0 - NORMAL ACTIVITY (CIWA AG 0) {9FF32084-A432-467C-B0A2-C8BC9F3D6D13}
- 2. 1- SOMEWHAT MORE THAN NORMAL ACTIVITY (CIWA AG 1) {497F7A3D-72C8-4986-A1E6-74EA0E5704F9}
- 3. 2 (2) {38E22B79-A2B6-3401-B195-C955F24E1278}
- 4. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}
- 5. 4 - MODERATELY FIDGETY AND RESTLESS (CIWA AG 4) {B1515E1F-4305-40D5-84B4-5E37AEF0AD3D}
- 6. 5 (5) {188EDCE7-0726-4234-D8D2-92C8052CD221}
- 7. 6 (6) {640CF2F9-3BBF-AB13-185B-CD02A609FB5A}
- 8. 7 - PACING BACK AND FORTH MOST OF THE INTERVIEW (CIWA AG 7) {9F03A167-CC40-4D65-9FD8-CFB3FEEAF833}

UNITS:

=====
TERM: CIWA: ANXIETY
ID: {9552124C-8650-406F-86FA-AB31D04AC8EF}

Page: 107

=====
Abbreviation..... CIWA AX
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Anxiety
Help Text..... Ask "Do You Feel Nervous?"
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 8 entries

- 1. 0 - NO ANXIETY AT EASE (CIWA AX 0)
 {DA61A13E-6443-4E8A-9C76-762FEACDB4E7}
2. 1 - MILDLY ANXIOUS (CIWA AX 1)
 {401113EE-5CE7-4AE9-80C1-1C3252E706D0}
3. 2 (2) {38E22B79-A2B6-3401-B195-C955F24E1278}
4. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}
5. 4 - MODERATELY ANXIOUS (CIWA AX 4)
 {9A01FEFD-844C-4FAB-BEC3-96333495FA62}
6. 5 (5) {188EDCE7-0726-4234-D8D2-92C8052CD221}
7. 6 (6) {640CF2F9-3BBF-AB13-185B-CD02A609FB5A}
8. 7 - ACUTE PANIC (CIWA AX 7)
 {3F5E6786-FD67-405F-A535-EC1609910F14}

UNITS:

=====
TERM: CIWA: AUDITORY DISTURBANCES
ID: {E6AA330B-96F9-4189-AE97-B93A921055A8}

Page: 108

=====
Abbreviation..... CIWA AU
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Auditory Disturbances
Help Text..... Ask " Are You More Aware Of Sounds
 Around You? Are They Harsh? Do They
 Frighten You? Are You Hearing Things
 You Know Are Not There?"
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 8 entries

- 1. 0 - ABSENT (0-a) {75F14311-F113-963C-4F05-B8AE17F044A9}
- 2. 1 - VERY MILD HARSHNESS OR ABILITY TO FRIGHTEN (CIWA AU 1) {E7D03A9C-B975-48B2-8CA6-2FC6616C6AF6}
- 3. 2 - MILD HARSHNESS OR ABILITY TO FRIGHTEN (CIWA AU 2) {71204969-E614-46A4-ADA5-39FAF81A79B5}
- 4. 3 - MODERATE HARSHNESS OR ABILITY TO FRIGHTEN (CIWA AU 3) {77F59277-6CE4-4323-8531-3CA26E13C9C4}
- 5. 4 - MODERATELY SEVERE HALLUCINATIONS (CIWA 4) {13FCDA97-A967-4993-AED5-13EBA61F3102}
- 6. 5 - SEVERE HALLUCINATIONS (CIWA 5) {14DB90F7-0182-44AC-8B5B-CE86A0E5A9E6}
- 7. 6 - EXTREMELY SEVERE HALLUCINATIONS (CIWA 6) {9C63ED29-2B47-4E59-A026-1FDB0666AF6B}
- 8. 7 - CONTINUOUS HALLUCINATIONS (CIWA 7) {19ECC57C-26A6-4CAE-97C4-2A860F337BE2}

UNITS:

=====

TERM: CIWA: HEADACHE FULLNESS IN HEAD
ID: {BAA1C221-5540-4A46-96AA-32C7BA53317B}

=====

Page: 109

Abbreviation..... CIWA HA
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Headache, Fullness In Head
Help Text..... Ask "Does your head feel different?
Does it feel like there is a band
around your head? Do not rate for
dizziness or lightheadedness."
Otherwise, rate severity.
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 8 entries

- 1. 0 - ABSENT (0-a) {75F14311-F113-963C-4F05-B8AE17F044A9}
- 2. 1 - VERY MILD (CIWA HA 1) {F789C1EA-CC7A-4368-87BA-08AFFE492EE6}
- 3. 2 - MILD (CIWA HA 2) {6B752397-0A4F-4176-90E7-E220DB2BF1C4}
- 4. 3 - MODERATE (CIWA HA 3) {F26EB828-6864-4384-9EBB-35A6CDCACC59}
- 5. 4 - MODERATELY SEVERE (CIWA HA 4) {83B3DE0F-3BA4-4B0A-A222-F372BD403094}
- 6. 5 - SEVERE (CIWA HA 5) {7BD7B91F-DF6F-47B9-BAA6-3EF249340578}
- 7. 6 - VERY SEVERE (CIWA HA 6) {241C1EF4-2E2A-4D65-9991-607BB6B29E01}
- 8. 7 - EXTREMELY SEVERE (CIWA HA 7) {F8DCC910-D36D-480E-867E-B45221170C49}

UNITS:

- 1. 0 - ORIENTED AND CAN DO SERIAL ADDITIONS (CIWA OR 0)
{D8015F54-6E9E-457A-955C-38565BDDDBB6}
- 2. 1 - CANNOT DO SERIAL ADDITIONS (CIWA OR 1)
{F4446E01-E326-420E-8EA8-04451AA9D963}
- 3. 2 - DISORIENTED FOR DATE < 2 CALENDAR DAYS (CIWA OR 2)
{01B894B0-3E70-47FA-820A-83CA48B93EF9}
- 4. 3 - DISORIENTED FOR DATE > 2 CALENDAR DAYS (CIWA OR 3)
{68C351F9-C8DE-49C5-B3BD-D69E97F3E516}
- 5. 4 - DISORIENTED FOR PLACE OR PERSON (CIWA OR 4)
{03BBFFA2-233A-41A1-8212-F988889C5E72}

UNITS:

=====

TERM: CIWA: PAROXYSMAL SWEATS
ID: {5BA369FD-68F0-4FD6-A564-1C50063E69AC}

=====

Page: 112

Abbreviation..... CIWA PS
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Paroxysmal Sweats
Help Text..... Paroxysmal Sweats
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 8 entries

- 1. 0 - NO SWEAT VISIBLE (CIWA PS 0)
{7C83E616-4414-4766-AAF5-5F5442FFBF09}
- 2. 1 - BARELY PERCEPTIBLE SWEATING / PALMS MOIST (CIWA PS 1)
{93B5B833-14B5-4ADD-97E8-7BF14A8BA01A}
- 3. 2 (2)
{38E22B79-A2B6-3401-B195-C955F24E1278}
- 4. 3 (3)
{E714454D-6A67-91F6-2FB7-D8F24F184328}
- 5. 4 - BEADS OF SWEAT OBVIOUS ON FOREHEAD (CIWA PS 4)
{8388B1A0-46EA-4003-A32B-A2C708A6F056}
- 6. 5 (5)
{188EDCE7-0726-4234-D8D2-92C8052CD221}
- 7. 6 (6)
{640CF2F9-3BBF-AB13-185B-CD02A609FB5A}
- 8. 7 - DRENCHING SWEAT (CIWA PS 7)
{15D2E2CC-BAED-403B-A141-B848FFE4A9D0}

UNITS:

=====

TERM: CIWA: TACTILE DISTURBANCES
ID: {306FFC89-BE89-4E74-90BC-AD8EF2273A51}

=====

Page: 113

Abbreviation..... CIWA TC

Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Tactile Disturbances
Help Text..... Ask "Have You Any Itching, Pins and
Needles Sensations, Any Burning, Any
Numbness Or Do You Feel Bugs Crawling
On Or Under Your Skin?"
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 8 entries

- 1. 0 - NONE (CIWA TC 0) {5F9F33C7-EFCE-4AC1-BCA3-1967F68EB597}
- 2. 1 - VERY MILD ITCHING PINS/NEEDLES BURNING (CIWA TC 1)
{00642E0D-E1CF-4EE0-91F4-CBAF209AA54C}
- 3. 2 - MILD ITCHING PINS/NEEDLES BURNING NUMBNESS (CIWA TC 2)
{1A591F17-CFB3-4B24-8AD6-ECF09318E2B6}
- 4. 3 - MODERATE ITCHING PINS/NEEDLES BURNING (CIWA TC 3)
{B7E2BFDD-7D85-4BCB-BE1D-0A220C1CEDAE}
- 5. 4 - MODERATELY SEVERE HALLUCINATIONS (CIWA 4)
{13FCDA97-A967-4993-AED5-13EBA61F3102}
- 6. 5 - SEVERE HALLUCINATIONS (CIWA 5)
{14DB90F7-0182-44AC-8B5B-CE86A0E5A9E6}
- 7. 6 - EXTREMELY SEVERE HALLUCINATIONS (CIWA 6)
{9C63ED29-2B47-4E59-A026-1FDB0666AF6B}
- 8. 7 - CONTINUOUS HALLUCINATIONS (CIWA 7)
{19ECC57C-26A6-4CAE-97C4-2A860F337BE2}

UNITS:

TERM: CIWA: TOTAL Page: 114
ID: {BD47158C-B56B-4B92-AD2E-33B7149EF6C0}

Abbreviation..... CIWA TOTAL
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Total CIWA Score
Help Text..... Add All CIWA Scores Together
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

CIWA SCALE -1 > 0 * 9 < 68 Dec Prec: 0

Abbreviation..... CIWA TR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Tremor
Help Text..... Arms extended and fingers spread apart.
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 8 entries

-
- 1. 0 - NO TREMOR (CIWA TR 0) {2D1A5B39-4C5D-4345-84B4-DDCBC8E25B4D}
 - 2. 1 - NOT VISIBLE BUT CAN BE FELT FINGERTIP (CIWA TR 1)
{31833D55-9018-47EE-8341-A0863744CEF8}
 - 3. 2 (2) {38E22B79-A2B6-3401-B195-C955F24E1278}
 - 4. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}
 - 5. 4 - MODERATE WITH PATIENTS ARMS EXTENDED (CIWA TR 4)
{79CE1486-151C-4492-AFD2-7845DD69DA9E}
 - 6. 5 (5) {188EDCE7-0726-4234-D8D2-92C8052CD221}
 - 7. 6 (6) {640CF2F9-3BBF-AB13-185B-CD02A609FB5A}
 - 8. 7 - SEVERE EVEN WITH ARMS NOT EXTENDED (CIWA TR 7)
{FE0BE72E-913A-48CD-BABF-C5B742AFED93}

UNITS:

Abbreviation..... CIWA VS
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Visual Disturbances
Help Text..... Ask "Does the light appear to be too
bright? Is its color different? Does
it hurt your eyes? Are you seeing
anything that is disturbing to you?
Are you seeing things you know are
not there?"
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 8 entries

- 1. 0 - ABSENT (0-a) {75F14311-F113-963C-4F05-B8AE17F044A9}
- 2. 1 - VERY MILD SENSITIVITY (CIWA VS 1)
{8BF96CF1-3EA7-468C-838A-284430EA850D}
- 3. 2 - MILD SENSITIVITY (CIWA VS 2)
{ED600249-E5E8-4020-955E-A5CDA18B5DE9}
- 4. 3 - MODERATE SENSITIVITY (CIWA VS 3)
{D673A630-1F1E-4734-AADE-17586C48E8B9}
- 5. 4 - MODERATELY SEVERE HALLUCINATIONS (CIWA 4)
{13FCDA97-A967-4993-AED5-13EBA61F3102}
- 6. 5 - SEVERE HALLUCINATIONS (CIWA 5)
{14DB90F7-0182-44AC-8B5B-CE86A0E5A9E6}
- 7. 6 - EXTREMELY SEVERE HALLUCINATIONS (CIWA 6)
{9C63ED29-2B47-4E59-A026-1FDB0666AFCB}
- 8. 7 - CONTINUOUS HALLUCINATIONS (CIWA 7)
{19ECC57C-26A6-4CAE-97C4-2A860F337BE2}

UNITS:

=====

TERM: CLOSE OBSERVATION BEHAVIOR Page: 117
 ID: {42CA9766-A28E-5379-67D7-2BB3BA4C2AE6}

=====

Abbreviation..... CO BEH
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... No Description
 Display Name..... Close Observation Behavior
 Help Text..... No Help Text
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 7 entries

- 1. WALKING OR PACING (wop) {CE7AABCB-84B7-A310-3A65-0AF17FB7B0E7}
- 2. LYING OR SITTING (los) {52A940B1-55AF-AB4E-9E59-F93515A1296C}
- 3. QUIET OR RESTING (qor) {76E4908F-014C-2088-FF7F-33EBF223DE40}
- 4. SLEEPING (s) {87338420-F351-F98C-7A80-4990B87ACE2B}
- 5. WITH PHYSICIAN (wp) {1673C53F-DA8C-2CEC-679B-B01D914FE533}
- 6. WITH THERAPIST (wt) {795A1300-1D60-9D67-692D-A2AB2B190596}
- 7. GERI-CHAIR (g) {167F9786-8F24-054F-020A-658A6EE5B0B4}

UNITS:

=====

TERM: CLOSE OBSERVATION CURRENT LOCATION Page: 118
 ID: {9FFE899C-A69C-4654-9AD7-B0E6F7FF888B}

=====

Abbreviation..... CO CUR LOC
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Close Observation Current Location
Display Name..... Close Observation Current Location
Help Text..... Close Observation Current Location
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 10 entries

-
1. PATIENT'S ROOM (pr) {0F13E11D-9F78-14A3-0568-AF26393F39F1}
 2. BATHROOM (b) {D20D91A7-E9F5-D000-5E0B-772BA9F5F72E}
 3. HALL (h) {BFCEA10F-2663-199E-88D8-E263122D7D26}
 4. DAY ROOM (dr) {269125EF-FA82-B17A-6B56-0CEAB4788899}
 5. QUIET ROOM (qr) {CCD96D41-C273-F6F3-DA18-8818B1CC9980}
 6. GROUP ROOM (gr) {D41C3D9F-8F10-1A0F-D571-20DE0EA36328}
 7. SMOKING AREA (sa) {8AFDBC0-9FC8-523F-5D07-4D8CCBA7036F}
 8. SMOKING PORCH (sp) {5AD3F05B-2A5B-E2B4-FD47-66712CFA44E6}
 9. MEDICATION ROOM (mr) {9571AA9F-F538-55EC-37F6-99FE5677B970}
 10. NURSING STATION (ns) {B6E48F92-B9BA-488F-3A2E-D74568A6B855}

UNITS:

=====

TERM: CLOSE OBSERVATION REASON FOR PRECAUTION
ID: {1B131F5A-CD32-38CA-19A5-1038A07C9867}

=====

Page: 119

Abbreviation..... CO RFP
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Close Observation Reason for Precaution
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

-
1. SUICIDAL (s) {B104E9C7-A78E-82DC-49A4-19B94A639AFE}
 2. ESCAPE (e) {965F3201-9D63-66AB-FA25-99479AE7145A}
 3. ASSAULTIVE (a) {99ECD8C7-191A-3B9C-2A7A-803A921065FA}
 4. ELOPEMENT (e) {EB030907-FB90-CA9D-CF50-E26C5D9FE671}
 5. MEDICAL DEVICE SAFETY (mds)
{1985CABE-4F53-9FCF-F444-F9AC53C04DBD}

UNITS:

=====

TERM: CODE STATUS

Page: 120

ID: {95248A8C-E05A-94F3-CA15-6F113997A621}

=====

Abbreviation..... CS
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Code Status
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

1. FULL CODE (fc) {1DED65E6-2F2E-8C5F-D009-BE6467A9A2AA}
2. DO NOT RESUSCITATE (dnr) {0BD8A134-C563-7E00-A6B1-F7C543CEF6FF}

UNITS:

=====

TERM: COLLOID TYPE

Page: 121

ID: {C22CB099-BC3F-42AE-9A11-67B99093D959}

=====

Abbreviation..... COLLOID TY
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Type of Colloid
Display Name..... Colloid Type
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

1. 5 PERCENT ALBUMIN (5 PERC ALB)
{38CE47CB-736B-475E-8E89-C4C13F2DF184}
2. 25 PERCENT ALBUMIN (25 PERC AL)
{2FCE19CA-B5AB-4F15-8FEA-CF156D1E4CBC}
3. 6 PERCENT HETASTARCH (6 PERC HET)
{D94072C9-87E8-46CA-929F-4904F3BD964C}

UNITS:

=====

TERM: COMA SCORE EYE

Page: 122

ID: {EAAF4E71-3BE6-5058-07FF-1A3D977BC569}

=====

Abbreviation..... CSE
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... 0: CLOSED, 1: NONE, 2: TO PAIN, 3: TO
SPEECH, 4: OPEN SPONTANEOUSLY
Display Name..... Coma Score
Help Text..... 0: CLOSED, 1: NONE, 2: TO PAIN, 3: TO
SPEECH, 4: OPEN SPONTANEOUSLY
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

-
1. CLOSED (CLOSED) {9BA6EE03-F562-4DE2-99A7-337F5A5E34F2}
 2. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
 3. TO PAIN (TO PAIN) {AA1F818D-DCFC-4FBC-97D3-BE86D309119E}
 4. TO SPEECH (TO SPEECH) {2DDD94FA-D4D0-48B0-B4F4-F25C42E4AF96}
 5. OPEN SPONTANEOUSLY (OPEN SPONT)
{013F1473-BDC9-4E7C-9709-C577357709F6}

LOCATION: 2 entries

-
1. L EYE (L EYE) {D714A0CA-EAD6-EC35-DA54-FF027C23FF45}
 2. R EYE (R EYE) {258F0C1D-67BD-ECF8-3152-EB289DEF4681}

UNITS:

=====

TERM: COMA SCORE MOTOR

Page: 123

ID: {73CA2132-12AC-BD84-3847-2E2EA3689D94}

=====

Abbreviation..... CSM
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Coma Score Motor
Help Text..... 1: None, 2: Extends to pain, 3:
Flexion to pain, 4: Withdraws to
pain, 5: Localizes to pain, 6: Obeys
Commands
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 6 entries

=====

TERM: CPAP Page: 127
ID: {11AD0BDB-8E33-C04E-29B2-CCB22FDF9EEC}

=====

Abbreviation..... CPAP
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Continuous Positive Airway Pressure
Display Name..... CPAP
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid..... 4688690

UNITS:

CENTIMETERS H2O -1 > 0 * 15 < 100 Dec Prec: 0

=====

TERM: CUFF INFLATION CHECKED Page: 128
ID: {77D485C1-8A5F-225A-CEA3-1219601C23E9}

=====

Abbreviation..... CUFF INFL
Active..... Yes
Boolean Value False..... Not Checked
Boolean Value True..... Checked
Data Type..... Boolean
Description..... VERIFICATION THAT PATIENT'S CUFF
INFLATION HAS BEEN CHECKED
Display Name..... Cuff Inflation Checked
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: DIAGNOSIS TYPE Page: 129
ID: {3289986D-87C6-4796-AA19-76010682BBAE}

=====

Abbreviation..... DIAG TYPE
Active..... Yes
Boolean Value False.....
Boolean Value True.....

Data Type..... PickList
Description..... Type of Diagnosis
Display Name..... Diagnosis Type
Help Text..... Type of Diagnosis
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

- 1. POST-OPERATIVE (POST-OP) {9FA19656-849F-4A65-92FC-755A4C9E8056}
- 2. ADMISSION (ADMISSION) {483C2830-813A-46C3-A012-CCD70044B2A2}
- 3. POST PROCEDURE DIAGNOSIS (PO PRO DIA)
{933C7D2C-88F8-4C18-ACBB-A16A295D1B3A}
- 4. DISCHARGE (DISCHARGE) {ADF67CC2-E05B-4B8D-A233-4AE1C35B2D20}
- 5. TRANSFER (TRANSFER) {DFDB4783-25BE-485B-A9BA-621F1004AAF2}

UNITS:

TERM: DIALYSIS CATHETER CONDITION Page: 130
ID: {E9D96923-CBA0-67B4-1D5B-3D6134AD8584}

Abbreviation..... D CATH CON
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Dialysis Catheter Condition
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

- 1. PATENT CLEAN (PAT CLEAN) {68688246-AE58-2EFB-BF9D-5E00D8CFDA8E}
- 2. RED (RED) {515D7CFB-3527-30E7-EC4B-B607C272DFA0}
- 3. SWOLLEN (SWOLLEN) {DCE86B10-C393-8DA9-6518-F2943F0E1A4E}
- 4. DISCOMFORT (DISCOMFORT) {77D058C9-1EB4-96C0-C9E3-666C840B884C}

LOCATION: 16 entries

- 1. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
- 2. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
- 3. L INTERNAL JUGULAR (LIJ) {AFBB4FE5-8A11-2047-7900-DADF2FC4430D}
- 4. R INTERNAL JUGULAR (R INT JUGU)
{7ABF6916-D379-198D-6454-F0518A74844A}
- 5. L EXTERNAL JUGULAR (L EXT JUGU)
{6F7B331B-4F83-0B79-1755-F6F673D647D6}
- 6. R EXTERNAL JUGULAR (R EXT JUGU)
{4D78F704-9E3B-907E-A792-4DC4F4CFF170}
- 7. L SUBCLAVIAN (LSCL) {71B93594-15F9-AA4A-AEF2-E222EC28C00D}

- 8. R SUBCLAVIAN (R SUBCLAVI) {53DEAE83-93AE-123B-14E5-4DE38CD4F752}
- 9. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
- 10. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
- 11. L ANTECUBITAL (L ANTECUBI)
{C21838CC-BAB5-1D7F-E2A4-A7C663FB28C3}
- 12. R ANTECUBITAL (R ANTECUBI)
{7FB00A65-B4B2-06A6-BD89-98AC112D8D4E}
- 13. L UPPER ARM (L UPP ARM) {5F915696-1B08-32E9-6B95-49246117D660}
- 14. R UPPER ARM (R UPP ARM) {3EF20E8E-F68C-E4F2-8F8C-FAD3B44DDAA4}
- 15. L WRIST (L WRIST) {7BCAFACE-FF3B-46DA-99F9-6437A90B7126}
- 16. R WRIST (R WRIST) {FDD9FE1D-F8A2-09AD-82D7-AA478693DD01}

UNITS:

=====

TERM: DIALYSIS DRESSING TYPE

Page: 131

ID: {A07D4532-93DF-178D-EB55-070FCF1C570F}

=====

Abbreviation..... DIAL DRESS
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... No Description
 Display Name..... Dialysis Dressing Type
 Help Text..... No Help Text
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 2 entries

- 1. TRANSPARENT (TRANSPAREN) {0E73D24D-0EFC-08F9-196A-D30DD1BB3228}
- 2. PRESSURE DRESSING (PRESS DRES)
{6D9DB9EF-85EB-8321-03FC-879AA4C5615B}

UNITS:

=====

TERM: DIALYSIS RATIO

Page: 132

ID: {F0CDCBBA-E65B-FE01-2E15-640FADF9DF20}

=====

Abbreviation..... DIAL RATIO
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... No Description
 Display Name..... Dialysis Ratio
 Help Text..... No Help Text
 Multi Select Picklist..... No
 Term Type..... OBSERVATION

Vuid.....

DATA VALUE: 5 entries

- 1. 1:1 (1:1) {5C07BBD0-FDB9-0CD7-F2C2-616913E0281B}
- 2. 1:2 (1:2) {287CDCE0-36F3-65F6-E6A5-9EF7CF8976AE}
- 3. 1:3 (1:3) {062B4F93-8F48-7373-F679-11F57CE98BB3}
- 4. 1:4 (1:4) {3F1B75FF-C168-D413-0BDF-6581E7A4EEAD}
- 5. 1:8 (1:8) {94790C4C-3AF2-6F4C-58B8-87021CA97A53}

UNITS:

=====

TERM: DIASTOLIC PRESSURE Page: 133
 ID: {14B24DAF-A0C0-060C-B1C0-ECEC1DF36E8E}

=====

Abbreviation..... DIASTOL PR
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... Numeric
 Description..... No Description
 Display Name..... Diastolic Pressure
 Help Text..... No Help Text
 Multi Select Picklist.....
 Term Type..... OBSERVATION
 Vuid.....

LOCATION: 4 entries

- 1. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
- 2. L LEG (LL) {B088C809-E528-B785-21C2-14E45E1FD8AF}
- 3. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
- 4. R LEG (RL) {C5E40A7B-6FC4-267D-8CE8-3571D21D81FD}

METHOD: 5 entries

- 1. CUFF (Cu) {E59A2543-1860-1E1C-BEF1-6E7F6B6FBC5E}
- 2. DOPPLER (Dop) {52683C19-F9DC-876B-CA3F-53E7C864CF04}
- 3. ARTERIAL LINE (AL) {D4B3F4ED-8B33-431E-4164-BD165E823776}
- 4. NON-INVASIVE (N-I) {F50ABFED-5FD7-FC95-20FE-1E3038C9F299}
- 5. INVASIVE (I) {5DB8746E-5446-F47A-883C-A44B91275A6E}

POSITION: 3 entries

- 1. LYING (Ly) {DE28D2BF-6F8D-64C3-0092-4CEB5E0902B9}
- 2. SITTING (Sit) {DD4DCDF1-0A22-2CC5-404F-0E472C73961C}
- 3. STANDING (St) {0776AE42-C370-B7B3-9B4D-9C914BA8F231}

UNITS:

MILLIMETERS OF MERCURY 0 > 40 * 90 < 300 Dec Prec: 0

Abbreviation..... DIET
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Type of diet ordered
Display Name..... DIET
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 46 entries

1. 1800 CALORIE CONTROLLED (1800 CAL C)
{C564D7FE-EC2C-401E-A77D-79CC5FC2B2F8}
2. 1800 DIABETIC (1800 DIAB) {FC71C7D7-F24E-4316-A455-8674804BE092}
3. 2000 CALORIE CONTROLLED (2000 CAL C)
{DB2FF9DC-298E-4881-94D4-47A7AF0E2BD7}
4. 2000 DIABETIC (2000 DIAB) {DF3056C9-3AEA-4F4E-B758-E35590080AFB}
5. 2200 CALORIE CONTROLLED (2200 CAL C)
{77F7C138-B2CB-4D76-8C0B-D68F317D4DE0}
6. 2200 DIABETIC (2200 DIAB) {1FB3EFB6-B92F-44FB-831E-4641C91C7A8D}
7. 2400 CALORIE CONTROLLED (2400 CAL C)
{F3EA28B5-FC96-4F34-B2BE-A96A196B38A2}
8. 2400 DIABETIC (2400 DIAB) {2018FB65-3D4E-47E0-90E9-205A45EDF130}
9. BLAND (BLAND) {38930AFE-ECF7-4548-857C-854E1B5BE148}
10. CARDIAC DIET (CARD DIET) {218A883D-32A0-438D-B453-7CBCE02693A3}
11. CLEAR LIQUID DIET (CL LI DIET)
{D540348D-3E1D-4E42-9EDC-43C56CF56206}
12. DUMPING SYNDROME (DUMP SYND)
{0173C971-63CB-45A5-9E93-D966EB3B957D}
13. DYSPHAGIA PHASE I (DYSP PH I)
{8011795B-77A5-4FF4-8EB9-EB6E1AA1FD4D}
14. DYSPHAGIA STAGE II (DYSP ST II)
{294BF58A-47FD-43A1-8A1E-9E348987B004}
15. DYSPHAGIA STAGE III (DYS ST III)
{DCF91830-C5CA-47B2-A7C6-4DD971F953EA}
16. DYSPHAGIA STAGE IV (DYSP ST IV)
{E738569A-5104-4C93-8EFC-47D80765F739}
17. DYSPHAGIA STAGE V (DYSPH ST V)
{ABE7D85D-1821-4592-B4E8-84670EE484E4}
18. FAT RESTRICTED 50 GM DIET (F RES 50gm)
{8E97A639-B7A6-4B3B-9BEF-0A46F3B46E43}
19. FECAL FAT TEST DIET (FEC F TEST)
{BA1E8E4A-7698-4EFB-B1DD-18C444EE82D5}
20. FULL LIQUID DIET (FULL LIQ D)
{CAC22234-7B69-44C3-8A0E-9A84A684E429}
21. FULL LIQUID DIABETIC (FULL L DIA)
{B89FAE3C-07EB-4AF6-B123-1153ADADC79F}
22. GASTRIC BYPASS CLEAR LIQUID (GAS BYP LI)
{6A87E820-7BEA-4636-B102-F9CB83BDD4D9}

23. GASTRIC BYPASS PUREED (GAS BY PUR)
{9F9EEB7F-7211-45A2-BCDE-B8398596CBF6}
24. GASTRIC BYPASS REGULAR (GAS BY REG)
{EE831906-EA51-4EE4-AAF2-583DF366A942}
25. GLUTEN RESTRICTED (GLUT RESTR)
{0805A164-D387-4E67-9D42-BAF1F2543337}
26. HIGH FIBER (HIGH FIBER) {3CAF7EE0-6A2D-41A3-BEDC-C57ECBBA72A2}
27. HIGH PROTEIN (HIGH PROTE) {1595A29D-CA91-45AE-893B-856EB909E4DD}
28. HIGH CALORIE (HIGH CALOR) {FC453BF0-1E7E-4733-9605-AC3425F16205}
29. LACTOSE RESTRICTED (LACT RESTR)
{78A9066A-48E2-4285-88E3-A180CE771B5A}
30. LIBERAL DIABETIC (LIBRL DIAB)
{9B66B47B-E77E-4F69-BDFA-5323B28A283D}
31. LISTED DIET (LIST DIET) {24FBBFD9-B426-4DD9-9843-3428A1C91676}
32. LOW CHOLESTEROL (LOW CHOLE)
{88F0170A-2560-48C6-9604-982453BCFC3B}
33. LOW FIBER (LOW FIBER) {2834D9A0-B6F7-4C21-96ED-1721D32B3B2B}
34. LOW PROTEIN 60 GM (L PROT 60g)
{BF5B6D8B-01A4-42D5-8551-6E78A17C43F7}
35. LOW TYRAMINE (LOW TYRAMI) {DB067A83-9D0C-4FC2-BAC8-6EC5641CAE42}
36. MECHANICAL (MECHANICAL) {01EA793C-B524-4685-A3EC-60A317CBACC5}
37. NEUTROPENIC (NEUTROPENI) {5CB9620C-34D9-48F7-AAF0-C595D06FC8AE}
38. NPO DIET (NPO DIET) {BEC8F771-B082-40A2-9EE8-553222A25FBF}
39. POST CATH DIET (POST CA DT)
{B18FEC50-18E1-4EED-8B37-E7DBAA7D4D3B}
40. POTASSIUM MODIFIED 60 MEQ (K MOD 60M)
{6EC2F10D-CC44-4872-8473-2445A835ABB0}
41. PUREED (PUREED) {99F8AE5F-E390-4583-8A92-CAC13E21BCF4}
42. REGULAR (REGULAR) {3D3F3D4C-ECD0-2056-A0D1-43CC2F2BA2C8}
43. RENAL DIET (REN DIET) {089DCCC2-B12B-4E3B-8D12-FCA5E974A55A}
44. SODIUM CONTROLLED 2 GM (SOD CONT 2)
{8A1311AF-6C7D-407A-962F-0FAE6C341F9A}
45. SODIUM CONTROLLED 4 GM (SOD CON 4G)
{26538A2F-5502-4D94-A2A7-D9E0BF825C47}
46. VEGETARIAN DIET (VEG DIET)
{6B8143BE-843F-4944-9ABB-7175AFE2DF17}

UNITS:

=====

TERM: DISCHARGE ACTIVITIES

Page: 135

ID: {22423A29-A1C5-497E-BBB6-109FFBC69D2B}

=====

Abbreviation..... DIS ACTI
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... No Description
 Display Name..... Discharge Activities
 Help Text..... No Help Text
 Multi Select Picklist..... Yes
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 12 entries

- 1. ALERT AND ORIENTED (AL AN OR) {C031F776-36DF-4570-AFAF-A30B8BFA635B}
- 2. VITAL SIGNS STABLE (VI SI ST) {E6ECD067-54A0-4214-A5F8-43D08C9FEAD6}
- 3. STEADY GATE OR BASELINE (ST GA BA) {0D47B2CE-88A7-4FE5-964B-7A8107A5DB91}
- 4. TRANSPORT FACILITATOR PRESENT (TR FA PR) {645BF269-6B23-4377-95ED-71291F075B25}
- 5. TOLERATING FLUIDS/NOURISHMENT (TO OR FL) {E107F272-CD27-46B8-8C50-F9FF776DDCD9}
- 6. IV DISCONTINUED (IV DISC) {FDF951E7-5B72-4166-969B-073D5EC74800}
- 7. VOIDED (VOIDED) {7F23E601-A1E5-4016-8A12-06B1099041CF}
- 8. DISCHARGE MEDS PLANNED (DI ME PL) {04E9DF8B-8807-47A5-AE88-180F431F3BDE}
- 9. MED RECONCILIATION DOCUMENTED (ME RE DO) {9ABF3E2E-4D63-43EF-A80D-9DCF5C820226}
- 10. DISCHARGE ORDER ON CHART (DI OR CH) {F7267E73-07BC-44FB-AB0F-FAA07F6E2C1B}
- 11. GIVEN DISCHARGE INSTRUCTION AND VERBALIZES UNDERSTANDING (DI IN UN) {C9A61418-88B2-42C4-B5F4-FC3CF76AB39D}
- 12. DRESSINGS CHECKED (DRESS CH) {21014F53-29CE-4914-8760-6976D8DBA742}

UNITS:

=====

TERM: DO NOT PROCESS Page: 136
 ID: {4615254C-EC67-46D0-BF70-9A54BEB4B32D}

=====

Abbreviation..... DNP
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... String
 Description..... Do not process
 Display Name..... Do Not Process
 Help Text..... Do Not Process
 Multi Select Picklist.....
 Term Type..... OBSERVATION
 Vuid.....

UNITS:

=====

TERM: DORSI FLEXION Page: 137
 ID: {A1E251EA-9302-8851-A67F-00229BCA7163}

=====

Abbreviation..... DORFL
 Active..... Yes

Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Dorsi Flexion
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

-
- | | |
|----------------|--|
| 1. ABSENT (a) | {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0} |
| 2. PRESENT (p) | {384DC20F-1625-CE1E-1C85-A3B3924BB2A0} |

UNITS:

=====

TERM: DRAINAGE CHARACTERISTICS
ID: {4D556AD1-1850-AB94-0776-B58B1D1713B3}

=====

Page: 138

Abbreviation..... DRAIN CHAR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Drainage Characteristics
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 15 entries

-
- | | |
|-------------------------|--|
| 1. THIN (THIN) | {CB8D8B5C-E824-D030-1241-D8EB814671E9} |
| 2. THICK (THICK) | {1554A49E-1468-46B2-F851-531DDB1D1A98} |
| 3. YELLOW (YELLOW) | {90DC31AD-6211-9F25-4237-79F7122A2763} |
| 4. AMBER (AMBER) | {90A13048-4C35-E01D-3482-64D782001D4A} |
| 5. WHITE (WHITE) | {AB04EF54-5246-EF3A-E1D9-02E9E24369DD} |
| 6. BROWN (BROWN) | {A1C32FBC-7B4D-707C-CA58-92A2E201368E} |
| 7. GREEN (GREEN) | {3D1F1553-56FC-68CE-79A9-E4F6357B6AAF} |
| 8. RED (RED) | {515D7CFB-3527-30E7-EC4B-B607C272DFA0} |
| 9. BLACK (BLACK) | {7A55AFB9-6861-4918-0746-3D67D47FE5C3} |
| 10. CLAY (CLAY) | {A4E4D16C-B45C-0934-ED22-8BF640080108} |
| 11. ORANGE (ORANGE) | {77B1B667-212D-88BA-D8E6-5C22F6ACDC9C} |
| 12. TAN (TAN) | {AB994ED7-8916-FAE0-0E97-CD4E0EF80ED1} |
| 13. PURULENT (PURULENT) | {86C01494-743A-C6A2-8883-01CA5200C79B} |
| 14. PINK (PINK) | {1A7FAF84-D69C-4EEB-7BC5-D21A598A1F4D} |
| 15. CLEAR (c) | {8981C365-2B46-704F-C9F4-975DF1C23C2D} |

UNITS:

ID: {C2EE711E-C076-491F-8063-8DA80BBCA529}

=====

Abbreviation..... DVT PROPHY
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... DVT Prophylaxis devices
Display Name..... DVT Prophylaxis
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

1. ANTI-EMBOLUS STOCKINGS KNEE HIGH (ANT EMB KN)
{4FBFF7A8-D719-98B9-4A46-D76D7BB5213E}
2. ANTI-EMBOLUS STOCKINGS THIGH HIGH (ANT EMB ST)
{9B32D3F3-47F1-886F-A59A-D8AB38A34754}
3. SEQUENTIAL COMPRESSION DEVICE (SCD)
{29D657DB-6679-8A55-0E33-235B5DD349F2}
4. ELASTIC WRAP (ELAST WRAP) {9EB52FA5-32D8-823F-2E3C-83375DA62342}

UNITS:

ID: {64C54E39-9541-4FD2-8F3F-A32053F6889F}

=====

Abbreviation..... EDEMA AMOU
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Edema Amount
Display Name..... Edema Amount
Help Text..... 0=Absent, +1=Slight (0-2 mm),
+2=Moderate (3-5 mm), +3=Marked (6-8
mm), +4=Very Marked (more than 8 mm)
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

1. 0 - ABSENT (0-a) {75F14311-F113-963C-4F05-B8AE17F044A9}
2. +1 = SLIGHT (+1 SLIGHT) {7C0245BA-B261-4B68-A84C-D737CA634230}
3. +2 = MODERATE (+2 = MODER)
{CCF23A60-068C-443F-AFD0-2D56EEE0446C}
4. +3 = MARKED (+3 = MARKE) {41ACC586-0C68-4780-8B2D-F36B21FCEDB2}

5. +4 = VERY MARKED (+4 = VY MA)
{55EF9D42-B4F8-4468-BA1B-98BCD4075E58}

UNITS:

=====

TERM: EDEMA SEVERITY Page: 141
ID: {E60900D9-BDF3-954C-E752-F3D1CAAC8EC1}

=====

Abbreviation..... EDEMA
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Observation Of Patients Edema
Display Name..... Edema Severity
Help Text..... 0=Absent, +1=Slight (0-2 mm),
+2=Moderate (3-5 mm), +3=Marked (6-8
mm), +4=Very Marked (more than 8 mm)
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

-
1. 0 - ABSENT (0-a) {75F14311-F113-963C-4F05-B8AE17F044A9}
 2. +1 = SLIGHT (+1 SLIGHT) {7C0245BA-B261-4B68-A84C-D737CA634230}
 3. +2 = MODERATE (+2 = MODER)
{CCF23A60-068C-443F-AFD0-2D56EEE0446C}
 4. +3 = MARKED (+3 = MARKE) {41ACC586-0C68-4780-8B2D-F36B21FCEDB2}
 5. +4 = VERY MARKED (+4 = VY MA)
{55EF9D42-B4F8-4468-BA1B-98BCD4075E58}

LOCATION: 28 entries

-
1. GENERALIZED (GENERIZED) {17073C6C-8691-CADC-2238-CBAAB0A5C648}
 2. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
 3. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
 4. L HAND (L HAND) {639DF25A-19E3-B504-D323-4B3F6ED39A24}
 5. R HAND (R HAND) {A54D3A80-8537-D079-9CCE-7DC56346D132}
 6. BOTH HANDS (BOTH HANDS) {63AA2685-CB1D-4680-B2F4-2199A22D1EC6}
 7. TORSO (TORSO) {EA1AD477-8BFF-E575-23C1-B9CFAF2B9986}
 8. ABDOMEN (ABDOMEN) {109FE041-2521-6EF2-7F95-1FD1587796AB}
 9. LEFT NECK (L NECK) {40522B5F-F23E-EADB-AE64-8EDE3B9E1155}
 10. R NECK (R NECK) {7E4D84D3-535B-6164-C6C1-837C6F997A0D}
 11. BILATERAL NECK (BILAT NECK)
{744BF4A2-C55E-4D55-8711-58700DD1C446}
 12. L FACE (L FACE) {06F00A70-F475-CD40-665D-AF0569A1E75F}
 13. R FACE (R FACE) {D1B15677-EF6F-0125-3601-68C1FE6386CE}
 14. FACE (FACE) {C638EA14-7903-6A94-496A-C6560D51BEE3}
 15. L FOOT (L FOOT) {2EC3FBEC-4EA2-9E7A-17D4-FDC2F9D28046}
 16. R FOOT (R FOOT) {7D9F5FF6-232C-2488-4504-A721FF06928F}
 17. FEET (FEET) {944CAA0E-5165-A6AF-BE3D-85455440BF40}
 18. L KNEE (L KNEE) {3A9E30B4-09A9-967A-58EF-69A01EC4F859}

- 19. R KNEE (R KNEE) {BCD18D53-07B6-2C41-64FC-3A07A7E979E7}
- 20. KNEES (KNEES) {247D6621-3A85-F1F7-59BF-1FB0361B3C21}
- 21. SACRAL (SACRAL) {0CF8A8A7-AB35-B803-4544-0C91645D9D1C}
- 22. PERIORBITAL (PERIORBITA) {DBAC362B-6278-EDB4-7E28-3180321A38B2}
- 23. L EYE (L EYE) {D714A0CA-EAD6-EC35-DA54-FF027C23FF45}
- 24. R EYE (R EYE) {258F0C1D-67BD-ECF8-3152-EB289DEF4681}
- 25. L LEG (LL) {B088C809-E528-B785-21C2-14E45E1FD8AF}
- 26. R LEG (RL) {C5E40A7B-6FC4-267D-8CE8-3571D21D81FD}
- 27. INCISIONAL SITE (INCISON SI)
{EB36A66D-751D-4031-9F68-0BB03E267321}
- 28. BOTH LEGS (BOTH LEGS) {5A15A231-4695-4CF3-AD0F-2E7950208967}

UNITS:

=====

TERM: EDUCATION - BARRIERS

Page: 142

ID: {6B9914EF-CEA5-4B2A-A9CC-5ED6FBA420F7}

=====

Abbreviation..... EDU BARR
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... No Description
 Display Name..... Barriers
 Help Text..... No Help Text
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 15 entries

- 1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
- 2. HEARING (HEARING) {D1A8AC66-1A5F-4F7D-AA89-529AF7A4AA11}
- 3. LANGUAGE (LANGUAGE) {0666538E-4F44-4D34-A5EA-69CF2F1F86FA}
- 4. LIMITED ATTENTION SPAN (LIM ATT S)
{51553020-0CB9-483E-AC1C-15CA6D9380A3}
- 5. MEMORY (MEMORY) {E30C1930-B98B-4382-B892-34E9C6008D1F}
- 6. PAIN (PAIN) {F2F23EE7-0A06-46B8-A496-B7A4651A551F}
- 7. SEDATED (s) {29D9ADD6-2C1F-F6D5-F4D0-B2BB86416997}
- 8. VISUAL IMPAIRMENT (VIS IMPAIR)
{7B2107CE-002D-43BE-9102-FCD5B9FAD617}
- 9. ANXIOUS (a) {381BDB41-7E8A-6FDB-841B-2ECC1889B9EF}
- 10. CONFUSED (CONFUSED) {31DAAA60-BE89-49A7-A10D-EB0DA82F12C0}
- 11. COMMUNICATION BARRIER (COMM BARRI)
{1A24946D-D966-4FAB-B6EA-940FA9EF60D3}
- 12. IMPAIRED LISTENING SKILLS (IMP LIS SK)
{9CE0E7E7-79E5-4910-918C-2022C9830171}
- 13. EMOTIONAL (EMOTIONAL) {626C22E6-EF2A-489F-AA79-D27541A0BC45}
- 14. UNABLE TO ASSESS (UNA TO AS)
{4D11A5E0-1CE8-4390-9EDD-ECEACCAA991B}
- 15. UNABLE TO READ (ILLITERATE)
{6D0C2C48-BDD2-489B-94C8-49E6DAC93F2E}

- {1EF29FF5-36AB-4FC0-B5E7-F8C3FC95BA5A}
- 6. CONTEMPLATIVE (CONTEMPLAT)
 - {5194DB04-7A68-491E-B117-BD089C912DFB}
- 7. UNABLE TO ASSESS (UNA TO AS)
 - {4D11A5E0-1CE8-4390-9EDD-ECEACCAA991B}

UNITS:

=====

TERM: EDUCATION - TOPIC Page: 145
 ID: {0F8D9E47-CDD7-4652-9E30-0E4355CC2EFA}

=====

Abbreviation..... EDU TOPIC
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... No Description
 Display Name..... Topic
 Help Text..... No Help Text
 Multi Select Picklist..... Yes
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 29 entries

1. ACTIVITIES OF DAILY LIVING (ADL)
 - {12F2CA66-3DCE-443D-883C-119CE4838B4A}
2. ADVANCED CARE PLANNING (ADV CAR PL)
 - {F0ED0F8A-1AD6-4F3E-A4C8-4FB915474FE7}
3. ANALGESIA (ANALGESIA)
 - {5452ED11-25E1-437B-8048-3BB3E03BE3B9}
4. CLINICAL ASSESSMENT (CLIN ASSME)
 - {803BEA4B-E22F-48C2-B20E-D5B3C8C7084C}
5. COMPLIANCE (COMPLIANC)
 - {EF43367E-9898-44C3-8AC8-F5BD28051F1D}
6. DAILY PLAN (DAILY PLAN)
 - {B29357A5-4294-4EA0-A0D3-59182F3C94FB}
7. DISCHARGE (DISCHARGE)
 - {ADF67CC2-E05B-4B8D-A233-4AE1C35B2D20}
8. DISEASE PROCESS (DISEA PROC)
 - {E80807AE-6208-421D-BE4C-445C470D32FA}
9. EQUIPMENT (EQUIPMENT)
 - {04C0F805-C7B3-4A9A-91D8-0F32F3265874}
10. HOME HEALTH (HOME HEALT)
 - {F94FB70A-AFFD-4A5E-A829-0AF933B044BB}
11. INFECTION (INFECTION)
 - {DA00DB58-640D-4294-9468-EEF070DBAE8E}
12. INFECTION PREVENTION AND SCREENING (INF PRE SC)
 - {BA58329F-59EC-49F0-BACE-75AA5EC57F52}
13. MEDICATION (MEDICATION)
 - {22E6DEBB-9662-466B-990C-8DA500C0A183}
14. MENTAL HEALTH (MEN HEAL)
 - {ED55FD2C-E2EC-4458-9AE4-FD7F0E6BD7F6}
15. MENTAL HEALTH HOLD (MEN HEA HO)
 - {31560285-7BD5-4E53-BB77-93B244367F42}
16. NEW ADMISSION EDUCATION (N ADM EDU)
 - {1890BC97-FEEB-4BED-A054-3C754C01817B}
17. NUTRITION (NUTRITION)
 - {429A28EF-8A1F-4247-85E8-68D4EFCDD189C}
18. ORIENTATION (ORIENTA)
 - {10FA5B0E-F4A7-44E4-A482-7C80F145C4D2}
19. PAIN (PAIN)
 - {F2F23EE7-0A06-46B8-A496-B7A4651A551F}
20. POST-OPERATIVE (POST-OP)
 - {9FA19656-849F-4A65-92FC-755A4C9E8056}
21. POST-PROCEDURE (POST-PROC)

- {0F943BD7-BD13-412E-B4EA-9B5526587340}
- 22. PRE-OPERATIVE (PRE-OP) {80C418D5-54B1-4CC3-9B19-5653B8D73C45}
- 23. PRE-PROCEDURE (PRE-PROC) {D2A7B025-9AB9-4C69-B163-D5EBB78B82BE}
- 24. SAFETY (SAFETY) {AD905A8F-33BD-4EFE-929E-A2B575797402}
- 25. SPIRITUAL (SPIRITUAL) {9670E2B9-6938-4487-A857-25546775B7A2}
- 26. TREATMENTS (TREATMENTS) {7A6879B5-F781-498E-9D71-3FA0243C3F21}
- 27. PRE-SEDATION (PRE SEDA) {684C764E-2BE9-4C99-A2B2-C92D75C73751}
- 28. POST-SEDATION (POST SED) {B829A249-CB14-4919-86F8-4A8C8D6B7753}
- 29. ENVIRONMENT (ENVIRO) {C61AA262-F6AF-46AC-A215-0D56FD01B234}

UNITS:

=====

TERM: EKG LEAD Page: 146
 ID: {D7BAE1AC-E3DA-4C9B-1216-7A07AA310A52}

=====

Abbreviation..... EKG
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... Numeric
 Description..... No Description
 Display Name..... EKG Lead
 Help Text..... No Help Text
 Multi Select Picklist.....
 Term Type..... OBSERVATION
 Vuid.....

LOCATION: 5 entries

- 1. CHEST (CHEST) {8C3C30FE-6E06-C35E-8226-CBF0C94014CC}
- 2. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
- 3. L LEG (LL) {B088C809-E528-B785-21C2-14E45E1FD8AF}
- 4. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
- 5. R LEG (RL) {C5E40A7B-6FC4-267D-8CE8-3571D21D81FD}

METHOD: 1 entries

- 1. SKIN (ADHESIVE INDICATOR) (SKIN)
 {C7CE7D84-6FC7-406A-7E56-B080637509DF}

QUALITY: 1 entries

- 1. ACTUAL (A) {5842A4D9-56A1-7A5A-13C9-F481699C3A09}

UNITS:

MILLIMETERS -21.00 > -20.00 * 20.00 < 21.00
Dec Prec: 2

=====

TERM: ENTERAL FEEDING TUBE TYPE Page: 147
 ID: {C0633B19-939C-4213-B7C7-C04DC2ADEF0C}

=====
Abbreviation..... EF TUBE TY
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Type of Enteral Feeding Tube used
Display Name..... E.F. Tube Type
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 8 entries

-
1. NASOGASTRIC LARGE BORE TUBE (NG LG BORE)
{C3FEE16A-9998-455B-BEB8-CC55AF890202}
 2. NASAL INTESTINAL SMALL BORE (N INT SM B)
{03DCB9EB-9B11-4C3B-869B-B1842BBCC844}
 3. PEG TUBE (PEG TUBE) {970C3755-2B3F-462F-AA87-EA33F846843F}
 4. JEJUNOSTOMY (JEJUNOSTOM) {1ED3EC41-3ACF-42C1-B64A-280300EDBCFC}
 5. GASTRIC TUBE (G TUBE) {65F774BA-F977-46D3-AE41-7F7C74CFF47E}
 6. ORAL GASTRIC TUBE (OG TUBE)
{C6995D6A-939C-4C31-A3C2-4504CB5597CD}
 7. DUODENOSTOMY (DUODENOSTO) {FA44135D-6806-42E7-B371-8C3912B3582E}
 8. TRANSPYLORIC (TRANSPYOLO) {78FA5E0A-EA21-40F3-8E8D-3B3FCC1ED060}

UNITS:

=====
TERM: ENTERAL NUTRITION RATE Page: 148
ID: {75A4F240-8D8E-4F7D-84D3-3DBCE1C20F00}

=====
Abbreviation..... ENT NUT RA
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Rate of Enteral Nutrition
Display Name..... Enteral Nutrition Rate
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

MILLILITERS PER HOUR -1 > 0 * 500 < 1000 Dec Prec: 0

=====
TERM: ENTERAL NUTRITION TYPE Page: 149
ID: {F083AE22-811C-432F-951A-3FF9DA485785}

Abbreviation..... ENT NUT TY
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Type of Enteral Nutrition Used
Display Name..... Enteral Nutrition Type
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 9 entries

-
1. MAGNACAL (MAGNACAL) {C810A618-FB1A-45BB-A540-EADF1414E343}
 2. OSMOLITE (OSMOLITE) {8E41BB55-063E-43C0-BAFE-65150EDD99F8}
 3. OSMOLITE HN (OSMOLIT HN) {E4186CA6-302F-4513-B125-30B094151423}
 4. ULTRACAL (ULTRACAL) {179AE9A5-154D-4CD2-B70B-66E84F5D0932}
 5. VITAL ENTERAL NUTRITION TYPE (VITAL ENT)
{2DA036DA-BF27-4E52-86A2-20A312823480}
 6. VITAL HN ENTERAL NUTRITION (VITAL HN)
{E45849D4-05EA-4DC1-BD9C-3FF66DD4D292}
 7. D5W ENTERAL NUTRITION (D5W)
{D6664641-BF43-4157-9A30-8AD8CFE96518}
 8. D10W ENTERAL NUTRITION (D10W ENT N)
{5049BCEF-CAE5-490D-9EB5-4A204A4D85E1}
 9. D50W ENTERAL NUTRITION (D50W ENT N)
{528A4339-0261-4EC5-821C-7CE5592666DC}

PRODUCT: 16 entries

-
1. FORMULA 2.0 CAL (2.0 C FORM)
{E5FB9590-25AA-435A-B9B4-4812EA15C217}
 2. FORMULA D5W ENTERAL NUTRITION (D5W ENT NU)
{9EA15A1B-DA1C-49EC-B463-82DD5F125F82}
 3. FORMULA D10W ENTERAL NUTRITION (D10W ENT N)
{F8FB55C4-6229-427F-987D-3D469E307D0D}
 4. FORMULA D50W ENTERAL NUTRITION (D50W ENT N)
{19ED9C83-42C1-42AF-AB4D-1A34B621EFB9}
 5. FORMULA - PROBALANCE (PROBALANCE)
{4EA729FE-08EA-490B-822D-696061E5B55F}
 6. GLUCERNA DIABETIC (GLUCERNA D)
{C2F2CCFE-6D77-49B4-BF5E-F1A285326440}
 7. HIGH PROTEIN FORMULA (HIGH PROTE)
{3618C874-3EAE-41CF-BE2B-3DFB960E034F}
 8. MAGNACAL (MAGNACAL) {86688633-DFCA-4FBE-99F5-8E0B3FAFE695}
 9. MODULAR CUSTOMIZED FORMULA (CUST FORMU)
{249D095E-9676-4CC5-8A39-097F96A4C032}
 10. NEPRO RENAL FORMULA (NEPR RENAL)
{59C9B396-A705-445B-B566-F030DA4B2B79}
 11. NUTREN 1.5 (NUTREN 1.5) {B93EF897-B9A1-45C6-AC9A-B386554A3186}
 12. OSMOLITE (OSMOLITE) {18C5843F-BC33-4609-9D58-49F07F12A059}
 13. OSMOLITE HN (OSMOLIT HN) {3622D5C6-1D7B-4377-94B7-5E33C642C4C5}
 14. ULTRACAL (ULTRACAL) {B18E6151-9DC6-4380-899A-E84715BFC61B}
 15. VITAL ENTERAL NUTRITION TYPE (VIT ENT NU)
{D77B607C-85C4-428E-AAE6-B319F33F4E9F}

16. VITAL HN ENTERAL NUTRITION (VIT NUT HN)
{3973F52F-D53F-4D80-928A-BD6154D5C434}

UNITS:

=====

TERM: ER EMERGENCY SEVERITY INDEX Page: 150
ID: {62C9FF10-1B3F-48D1-8FE7-ACF8284FD04A}

=====

Abbreviation..... ESI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... EMERGENCY SEVERITY INDEX
Display Name..... ESI
Help Text..... EMERGENCY SEVERITY INDEX
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

-
- 1. 1 (1) {8F5D10A8-FE17-47BC-DACA-6972D12CC00A}
 - 2. 2 (2) {38E22B79-A2B6-3401-B195-C955F24E1278}
 - 3. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}
 - 4. 4 (4) {D2280C7E-DE50-5A8D-C3FD-5862B69C0E1A}
 - 5. 5 (5) {188EDCE7-0726-4234-D8D2-92C8052CD221}

UNITS:

=====

TERM: ER TRIAGE CONDITION Page: 151
ID: {E4848AF5-F7CA-4B23-B436-B902B0F5555A}

=====

Abbreviation..... TRIAGE CON
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Triage Condition
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

-
- 1. STABLE (STABLE) {A474E7BC-47EF-4472-A19A-172C441099EC}
 - 2. UNSTABLE (UNSTABLE) {5841C521-8070-4034-A69D-27DAB2990EB3}
 - 3. OTHER (OTHER) {19094E5B-AEC3-44B0-9BD5-44E79E7EB105}

UNITS:

=====

TERM: ER TRIAGE DISPOSITION

Page: 152

ID: {4D5AB8A3-0883-4D48-AAF7-73933D411038}

=====

Abbreviation..... TRIAGE DIS
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Triage Disposition
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

1. EMERGENCY DEPARTMENT (ER DEPT)
{D2E511F7-8B80-4928-A64B-3F1EE17395F7}
2. WAITING ROOM (WAIT ROOM) {13476486-F895-4FFC-AB42-B274827B47DA}
3. OTHER (OTHER) {19094E5B-AEC3-44B0-9BD5-44E79E7EB105}

UNITS:

=====

TERM: ET LOCATION

Page: 153

ID: {7DA32BAB-D3E3-55F0-CBAC-812B366B15A5}

=====

Abbreviation..... ET LOC
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... LOCATION OF ENDOTRACHEAL TUBE
Display Name..... Endotracheal Tube Location
Help Text..... Measured in centimeters from lip or
nose
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 31 entries

1. 13 (13) {030242DB-3F82-7517-3DEE-4D50E3E8AE23}
2. 13.5 (13.5) {34D24648-B2D9-9C6C-6D71-F80E977CE185}
3. 14 (14) {EDF1C2B1-123B-011F-9B02-7D71088971C1}
4. 14.5 (14.5) {373C8F8C-5D36-4AD6-6A76-513064669CB5}
5. 15 (15) {FE8C049F-B64C-4991-2F50-3C561BA4825C}

- 6. 15.5 (15.5) {1BBF5122-2FFA-6B27-445B-3291C3F83A36}
- 7. 16 (16) {9E4B84FF-430A-39FB-E6AA-E1031814796D}
- 8. 16.5 (16.5) {7D4EA1D4-5233-760F-FD2B-A70FB0969815}
- 9. 17 (17) {2E89B9C8-5E80-F3C9-FE66-4396CB63201A}
- 10. 17.5 (17.5) {FDE6FFD9-00D0-9C55-1A7F-4AA70CE313BA}
- 11. 18 (18) {F2815231-0FCE-4988-95E9-D9AF46C57086}
- 12. 18.5 (18.5) {75ABCB85-4D6E-0BBC-157F-CA69B7737402}
- 13. 19 (19) {478FFB89-4941-1474-04E1-C838B2467439}
- 14. 19.5 (19.5) {58F0F98D-299A-8C4E-F50F-AE21B6E696D6}
- 15. 20.5 (20.5) {DB1F50D2-AF99-7BB3-88F5-B58C0BDE7987}
- 16. 21 (21) {C6188159-9EF1-DA5C-C015-2142A1205C4C}
- 17. 21.5 (21.5) {A259E466-773F-519C-7B5B-1570D1B03B70}
- 18. 22 (22) {DA839156-C7A9-76D7-58CD-7EC2B016EF40}
- 19. 22.5 (22.5) {D402D24F-1EF2-1AA0-DD1C-FA8320AA929A}
- 20. 23 (23) {8FB6A8B1-B692-BDC7-43B9-0332F33D826A}
- 21. 23.5 (23.5) {DDB76295-7E4A-FEAC-324D-37878F62771A}
- 22. 24 (24) {2BE934CF-40B3-2263-9FA7-6AC9DE0C7C8A}
- 23. 24.5 (24.5) {D815830B-5D94-0975-D007-E4FDD966155B}
- 24. 25 (25) {42DD0655-142D-6901-D3FE-66D295BBFCD1}
- 25. 25.5 (25.5) {9D704E9E-FB1D-FB83-2142-EDBF3805E8F7}
- 26. 26 (26) {181395F6-29FD-5A41-5AD7-BFC9CF564F3B}
- 27. 26.5 (26.5) {1D464967-6CFA-E29E-548D-EF68AFCDEC61}
- 28. 27 (27) {64D25A86-6516-0910-B4B5-644B85378876}
- 29. 27.5 (27.5) {8C7BE365-AFD0-3FFD-B8EC-6EC914763CFC}
- 30. 28 (28) {66B5A1F5-5396-10C4-643E-3E8CAB6D6CC8}
- 31. 28.5 (28.5) {C99C9298-D4D2-BC02-36D3-F612B6080A87}

UNITS:

=====

TERM: FALL PRECAUTIONS

Page: 154

ID: {7912B52F-20BC-9128-883B-08AF0EC48045}

=====

Abbreviation..... FALL PRE
 Active..... Yes
 Boolean Value False..... No
 Boolean Value True..... Yes
 Data Type..... Boolean
 Description..... No Description
 Display Name..... Fall Precautions
 Help Text..... Fall Precautions Taken
 Multi Select Picklist.....
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
- 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: FAMILY ORIENTED TO UNIT
ID: {E1CB7972-95FE-FE7B-9E6E-ECF1FF4699EB}

=====

Page: 155

Abbreviation..... FOTU
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... No Description
Display Name..... Family Oriented to Unit
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

1. NO (n)	{FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y)	{483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: FETAL HEART TONES
ID: {A2E22A44-E924-ADDE-2B8E-0251666B4DE6}

=====

Page: 156

Abbreviation..... FT
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... From Vitals
Display Name..... Fetal Heart Tones
Help Text..... From Vitals
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid..... 4688721

UNITS:

BEATS PER MINUTE 50 > 120 * 160 < 225 Dec Prec: 0

=====

TERM: FINGERSTICK GLUCOSE
ID: {805F184D-EF11-446F-A272-B08249A477E9}

=====

Page: 157

Abbreviation..... FINGER GLU
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Fingertstick Glucose Test Result

Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: FOOT STRENGTH

Page: 160

ID: {5B1A78F0-5C93-F5BA-29D5-D9AFC205B3C7}

=====

Abbreviation..... FOOT S
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Foot Strength
Help Text..... Evaluation of foot strength
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 5 entries

1. ABSENT (a) {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0}
2. WEAK (w) {BE2D9305-E2DA-F492-66CF-164EE58B8179}
3. STRONG (s) {A12154E1-544E-F372-3839-0D44E5BC476B}
4. REFLEX (r) {F21585F9-5206-A6B3-2BF6-F416E5842366}
6. N/A (n) {B3DBE3FD-C5CB-3894-607A-BE13319594D0}

LOCATION: 3 entries

1. L FOOT (L FOOT) {2EC3FBEC-4EA2-9E7A-17D4-FDC2F9D28046}
2. R FOOT (R FOOT) {7D9F5FF6-232C-2488-4504-A721FF06928F}
3. BILATERAL FOOT (BILAT FOOT)
{C758847D-468F-479A-B12D-F4E46F0E5E3A}

UNITS:

=====

TERM: FUNDAL HEIGHT

Page: 161

ID: {EEAB8762-624F-7BA3-4001-114FD229BA69}

=====

Abbreviation..... FH
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... From Vitals
Display Name..... Fundal Height
Help Text..... Starting at 24 weeks
Multi Select Picklist.....
Term Type..... OBSERVATION

=====
Abbreviation..... GAS RES CH
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... Was the Gastric Residual Checked?
Display Name..... Gastric Residual Checked
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====
TERM: GASTRIC RESIDUAL NOT DISCARDED Page: 165
ID: {D93EDA7E-706C-4914-A774-5BE46D4A3085}

=====
Abbreviation..... GAS NOT DI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Gastric Residual Amount Not Discarded
Display Name..... Gastric Residual Amount Not Discarded
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

MILLILITER -1 > 0 * 500 < 10000 Dec Prec: 0

=====
TERM: GASTRIC pH Page: 166
ID: {5ECE05DB-47C5-43C9-A3CD-0CD31E285BA8}

=====
Abbreviation..... GASTRIC pH
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Gastric pH
Display Name..... Gastric pH
Help Text..... No Help Text
Multi Select Picklist.....

Term Type..... OBSERVATION
Void.....

UNITS:

pH -1 > 0 * 7 < 15 Dec Prec: 0

=====

TERM: GCS - EYE OPENING Page: 167
ID: {11D9DC0B-FE5C-9374-0904-96896A69E41B}

=====

Abbreviation..... GCS EYE OP
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Eye Opening
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 5 entries

-
- 1. 4 - SPONTANEOUSLY (4-s) {4612FC84-FCA8-8C85-50DC-A9A42C195F51}
 - 2. 3 - TO VERBAL COMMAND (3-tvc)
 {EBF4B748-2861-C81F-C0D0-B9FE99FC7A8F}
 - 3. 2 - TO PAIN (2-tp) {49CE3C40-5FC2-002E-1C89-3E5DFA4CCF0F}
 - 4. 1 - NO RESPONSE (1-nr) {26D3168D-40EA-C918-D321-54BED4837311}
 - 5. 0 - ABSENT (0-a) {75F14311-F113-963C-4F05-B8AE17F044A9}

UNITS:

=====

TERM: GCS - MOTOR RESPONSE Page: 168
ID: {F708D210-5A49-4731-443A-DB3992450B45}

=====

Abbreviation..... GCS MOTOR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Motor Response
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 6 entries

-
- 1. 6 - OBEYS COMMANDS (6-oc) {0D4956F4-7E79-8D17-EDF3-274B96BBC89A}

- 2. 5 - LOCALIZES PAIN (5-lp) {CF13650F-899B-6AE1-5913-DC1F4776E5E4}
- 3. 4 - WITHDRAWS FROM PAIN (4-wfp)
 - {D808ACF1-30CB-29AB-4957-E91D63C4DF18}
- 4. 3 - ABNORMAL FLEXION FROM PAIN (3-affp)
 - {9AED8954-5F5D-D2C2-3CAF-4F8EF92009C2}
- 5. 2 - ABNORMAL EXTENSION FROM PAIN (2-aefp)
 - {9B4951D8-B271-DFC3-0A00-1DC9984039E3}
- 6. 1 - NO RESPONSE (1-nr) {26D3168D-40EA-C918-D321-54BED4837311}

UNITS:

=====

TERM: GCS - TOTAL Page: 169
 ID: {09118FCE-11BC-2679-68FC-67913B3EAAD0}

=====

Abbreviation..... GCS TOTAL
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... Numeric
 Description..... No Description
 Display Name..... Total GCS
 Help Text..... No Help Text
 Multi Select Picklist.....
 Term Type..... OBSERVATION
 Vuid.....

UNITS:

GLASGOW COMA SCALE 2 > 3 * 15 < 16 Dec Prec: 0

=====

TERM: GCS - VERBAL RESPONSE Page: 170
 ID: {B4BA8135-48CB-2930-B5F2-5B93688C07FE}

=====

Abbreviation..... GCS VERBAL
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... No Description
 Display Name..... Verbal Response
 Help Text..... No Help Text
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 5 entries

- 1. 5 - ORIENTED AND CONVERSES (5-oac)
 - {53AD29B1-BA51-4BDE-4F4C-4268619C93A4}
- 2. 4 - DISORIENTED AND CONVERSES (4-dac)
 - {D43F646D-A130-5D06-6DC6-DDB971027231}

- 3. 3 - INAPPROPRIATE WORDS (3-iw) {FE000191-243E-4D32-1E1E-546EE2CF88A5}
- 4. 2 - INCOMPREHENSIBLE SOUNDS (2-is) {F2261C3D-79A3-055B-702D-446B7A6A3113}
- 5. 1 - NO RESPONSE (1-nr) {26D3168D-40EA-C918-D321-54BED4837311}

UNITS:

=====

TERM: GENITAL ASSESSMENT Page: 171
 ID: {139D5B8A-4DF7-4764-A1C9-835C73A13E41}

=====

Abbreviation..... GEN. ASSMT
 Active..... Yes
 Boolean Value False..... Normal
 Boolean Value True..... Abnormal
 Data Type..... Boolean
 Description..... General condition of genitals
 Display Name..... Genital Assessment
 Help Text..... Are genitals normal or abnormal?
 Multi Select Picklist.....
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 2 entries

-
- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
 - 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: GI SYMPTOMS Page: 172
 ID: {2566D442-22EA-4523-8D77-007B858C1EA0}

=====

Abbreviation..... GI SYMPTOM
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... GI Symptoms
 Display Name..... GI Symptoms
 Help Text..... GI Symptoms
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 7 entries

-
- 1. NAUSEA (NAUSEA) {6BF33DC3-0300-4643-BB14-27C35692B806}
 - 2. VOMITING (VOMITING) {921F463C-1639-4D13-8E4E-6D3A1B6927F9}
 - 3. PROJECTILE VOMITING (PROJ VOMIT)

- {0D2CB23B-B266-40E1-A6FF-20B733655B69}
- 4. CRAMPING (CRAMPING) {3EC3E352-5469-4A85-80E4-3AFEB61224C4}
- 5. HEARTBURN (HEARTBURN) {5ECEAF5E-6B18-4177-A18D-2DEE87796CC2}
- 6. DIARRHEA (DIARRHEA) {353472EF-3499-46B1-AD7C-554FE01D8CE2}
- 7. CONSTIPATION (CONSTIPAT.) {4A3B0D47-88BF-4FEB-AA0D-7AD7FED2C5C3}

UNITS:

=====

TERM: GI TUBE DRAINAGE CHARACTERISTICS Page: 173
 ID: {D72728F5-FE40-1697-E896-BFAA1895A4A8}

=====

Abbreviation..... GI DRAIN C
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... No Description
 Display Name..... GI Tube Drainage Characteristics
 Help Text..... No Help Text
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 2 entries

- 1. THICK (THICK) {1554A49E-1468-46B2-F851-531DDB1D1A98}
- 2. THIN (THIN) {CB8D8B5C-E824-D030-1241-D8EB814671E9}

UNITS:

=====

TERM: GI TUBE DRAINAGE COLOR Page: 174
 ID: {4D0CB409-25EA-40CB-ABE5-1BB2E0FDE1C3}

=====

Abbreviation..... GI DRA COL
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... GI Tube Drainage Color
 Display Name..... GI Tube Drainage Color
 Help Text..... No Help Text
 Multi Select Picklist..... Yes
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 13 entries

- 1. YELLOW (YELLOW) {90DC31AD-6211-9F25-4237-79F7122A2763}
- 2. AMBER (AMBER) {90A13048-4C35-E01D-3482-64D782001D4A}
- 3. WHITE (WHITE) {AB04EF54-5246-EF3A-E1D9-02E9E24369DD}

- 4. BROWN (BROWN) {A1C32FBC-7B4D-707C-CA58-92A2E201368E}
- 5. GREEN (GREEN) {3D1F1553-56FC-68CE-79A9-E4F6357B6AAF}
- 6. RED (RED) {515D7CFB-3527-30E7-EC4B-B607C272DFA0}
- 7. BLACK (BLACK) {7A55AFB9-6861-4918-0746-3D67D47FE5C3}
- 8. CLAY (CLAY) {A4E4D16C-B45C-0934-ED22-8BF640080108}
- 9. ORANGE (ORANGE) {77B1B667-212D-88BA-D8E6-5C22F6ACDC9C}
- 10. TAN (TAN) {AB994ED7-8916-FAE0-0E97-CD4E0EF80ED1}
- 11. PURULENT (PURULENT) {86C01494-743A-C6A2-8883-01CA5200C79B}
- 12. PINK (PINK) {1A7FAF84-D69C-4EEB-7BC5-D21A598A1F4D}
- 13. CLEAR (c) {8981C365-2B46-704F-C9F4-975DF1C23C2D}

UNITS:

=====

TERM: GI TUBE PLACEMENT CONFIRMATION

Page: 175

ID: {BFA579B5-CEA7-43EF-84C0-C4893EFD0F3A}

=====

Abbreviation..... GI T PLACE
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... GI Tube Placement Confirmation
 Display Name..... GI Tube Placement Confirmation
 Help Text..... GI Tube Placement Confirmation
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 4 entries

- 1. X-RAY (X-RAY) {C82CB7F8-2456-4A08-9FD2-4E498EB1DB20}
- 2. ASPIRATION (ASPIRATION) {68B91085-8E3C-4B94-9F97-6305F5DFD900}
- 3. AUSCULTATION (AUSCULTATI) {B9D756DF-74F5-4604-B781-D7BC512FC453}
- 4. pH (pH) {104FE39C-B0A2-40B5-A523-B11E5F2BA272}

UNITS:

=====

TERM: GI TUBE SIZE

Page: 176

ID: {9F0E347C-C292-41A6-B213-340B2ECD716E}

=====

Abbreviation..... GI TU SIZE
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... Numeric
 Description..... GI Tube Size
 Display Name..... GI Tube Size
 Help Text..... GI Tube Size
 Multi Select Picklist.....
 Term Type..... OBSERVATION

Void.....

UNITS:

FRENCH 7 > 8 * 18 < 25 Dec Prec: 0

=====

TERM: GI TUBE SUCTION TYPE Page: 177
ID: {884A773C-187A-B910-AEC0-59815DEE579E}

=====

Abbreviation..... TUBE SUC T
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... TYPE OF SUCTION APPLIED TO GI TUBE
Display Name..... GI Tube Suction Type
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Void.....

DATA VALUE: 6 entries

-
- 1. HIGH (h) {D88CB5EC-2AD4-2F80-2D43-40ECB0B7D390}
 - 2. LOW (l) {B6CE52A4-C1DD-FDE3-F87E-8E0C98FCC4AC}
 - 3. CONTINUOUS (CONTINUOUS) {676CD2DA-F882-4BE5-99BB-510EA7669749}
 - 4. INTERMITTENT (INTERMITTE) {DD09ED7C-B977-412F-A678-5130F77D239E}
 - 5. CLAMPED (CLAMPED) {B86865B8-078A-4996-B5EC-46A4A99D5A7D}
 - 6. GRAVITY (GRAVITY) {C0233F70-376D-41CE-9F14-714C7B35244F}

UNITS:

=====

TERM: GI TUBE TYPE Page: 178
ID: {DD63CB7E-0759-F5F2-EFF6-AB6874590B43}

=====

Abbreviation..... GI TUBE TY
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... TYPE OF GI TUBE
Display Name..... GI Tube Type
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 10 entries

-
- 1. NASOGASTRIC TUBE (NG TUBE) {6475E343-A310-48A5-B543-D29573E7684D}

2. GASTRIC TUBE (G TUBE) {65F774BA-F977-46D3-AE41-7F7C74CFF47E}
3. JEJUNOSTOMY (JEJUNOSTOM) {1ED3EC41-3ACF-42C1-B64A-280300EDBCFC}
4. PEG TUBE (PEG TUBE) {970C3755-2B3F-462F-AA87-EA33F846843F}
5. NASAL SMALL BOWEL FEEDING TUBE (NAS SM BOW)
{A2968C8A-DDBD-4C3B-AFE6-3E2FED41C413}
6. ORAL GASTRIC TUBE (OG TUBE)
{C6995D6A-939C-4C31-A3C2-4504CB5597CD}
7. SMALL BORE FEEDING TUBE (SM BO TUBE)
{7E347B8A-14E7-4093-9A5B-56E6AC0D78DF}
8. GASTROSTOMY TUBE (GAST TUBE)
{EF017A39-745E-4C29-821A-5B25F103DF99}
9. GASTROSTOMY FEEDING BUTTON (GAST F BUT)
{1627B53B-374B-42B4-BBFE-463EE37B6110}
10. EWALD TUBE (EWALD TUBE) {B83AB4C1-A098-4FAB-A645-61FE1FEADBBE}

UNITS:

=====

TERM: GU BLADDER IRRIGATION FREQUENCY
ID: {D16099EA-DF97-4481-8B8F-7A25E88D15B1}

=====

Page: 179

Abbreviation..... GU BL IR F
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Bladder Irrigation Frequency
Display Name..... Bladder Irrigation Frequency
Help Text..... Bladder Irrigation Frequency
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

1. CONTINUOUS (CONTINUOUS) {676CD2DA-F882-4BE5-99BB-510EA7669749}
2. INTERMITTENT (INTERMITTE) {DD09ED7C-B977-412F-A678-5130F77D239E}
3. PRN (PRN) {36493D81-954B-4C56-B0A2-9988EA785EBD}

UNITS:

=====

TERM: GU CATHETER CARE PERFORMED
ID: {F99A31A1-6D6D-4532-81F4-35EA957233D7}

=====

Page: 180

Abbreviation..... GU CA CARE
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... Catheter Care Performed
Display Name..... Catheter Care Performed

Help Text..... Catheter Care Performed
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: GU CATHETER SIZE Page: 181
ID: {222FF9DB-7D48-4004-B9B9-49B41D8E510B}

=====

Abbreviation..... GU CA SIZE
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Catheter Size
Display Name..... Catheter Size
Help Text..... Catheter Size
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

FRENCH 11 > 12 * 30 < 31 Dec Prec: 0

=====

TERM: GU CATHETER TYPE Page: 182
ID: {A73423E8-EA6C-007E-1F20-6185BD488AEA}

=====

Abbreviation..... GU CATH TY
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Type of catheter used
Display Name..... Catheter Type
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Void.....

DATA VALUE: 13 entries

- 1. INDWELLING CATHETER (INDW CATH)
{EC16724C-D1EA-5650-58D4-38A89D4456F9}
2. EXTERNAL CATHETER (EXT. CATHE)

- {C8AEE936-DD91-FE91-4DD6-6F210FB59726}
- 3. SUPRAPUBIC (SUPRAPUBIC) {69E0176F-52FC-A3C4-7A4E-AD57DBD6E572}
- 4. L NEPHROSTOMY (L NEPHROST)
 - {D8FC4C18-081E-A7B9-E394-25F194BE9757}
- 5. R NEPHROSTOMY (R NEPHROST)
 - {33661137-688D-5611-4368-B1DFD8DB1469}
- 6. UROSTOMY (UROSTOMY) {1F9E528D-6D78-54EE-94E6-649C6F2056A3}
- 7. STRAIGHT CATHETER (STR. CATH)
 - {5F4A8A7F-E720-9458-2360-1DC14E7D05F3}
- 8. ILEAL CONDUIT (ILEAL COND)
 - {032CB6FE-77D2-702E-A0AE-46176F5CAA57}
- 9. L STENT (LEFT STENT) {FD917F0B-CC96-0CB0-BED3-1D7C51A77B69}
- 10. R STENT (RT. STENT) {6D0EED3B-1F08-11D2-165D-97A9F58327D5}
- 11. FOLEY (FOLEY) {9FE7D3A7-67CF-CF54-F721-2F8035CE2BF1}
- 12. 3 WAY FOLEY (3WAY FOLEY) {32408760-F517-4530-AEC0-54FE0A6A062B}
- 13. MUSHROOM CATHETER (MUSH CATHE)
 - {4FFC451E-14B5-46E4-99C9-99A17BF90AE4}

LOCATION: 16 entries

- 1. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
- 2. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
- 3. L INTERNAL JUGULAR (LIJ) {AFBB4FE5-8A11-2047-7900-DADF2FC4430D}
- 4. R INTERNAL JUGULAR (R INT JUGU)
 - {7ABF6916-D379-198D-6454-F0518A74844A}
- 5. L EXTERNAL JUGULAR (L EXT JUGU)
 - {6F7B331B-4F83-0B79-1755-F6F673D647D6}
- 6. R EXTERNAL JUGULAR (R EXT JUGU)
 - {4D78F704-9E3B-907E-A792-4DC4F4CFF170}
- 7. L SUBCLAVIAN (LSCL) {71B93594-15F9-AA4A-AEF2-E222EC28C00D}
- 8. R SUBCLAVIAN (R SUBCLAVI) {53DEAE83-93AE-123B-14E5-4DE38CD4F752}
- 9. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
- 10. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
- 11. L ANTECUBITAL (L ANTECUBI)
 - {C21838CC-BAB5-1D7F-E2A4-A7C663FB28C3}
- 12. R ANTECUBITAL (R ANTECUBI)
 - {7FB00A65-B4B2-06A6-BD89-98AC112D8D4E}
- 13. L UPPER ARM (L UPP ARM) {5F915696-1B08-32E9-6B95-49246117D660}
- 14. R UPPER ARM (R UPP ARM) {3EF20E8E-F68C-E4F2-8F8C-FAD3B44DDAA4}
- 15. L WRIST (L WRIST) {7BCAFACE-FF3B-46DA-99F9-6437A90B7126}
- 16. R WRIST (R WRIST) {FDD9FE1D-F8A2-09AD-82D7-AA478693DD01}

UNITS:

=====

TERM: GU VOIDING

Page: 183

ID: {BE9D4DD5-CF6E-C973-C1F1-D70D9F554E38}

=====

Abbreviation..... GU VOIDING
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... Method of urine elimination

Display Name..... Voiding
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

-
1. CONTINENT (CONTINENT) {3B658506-F42E-BAF5-04C4-5AD49699E0D2}
 2. INCONTINENT (INCONTINEN) {BA2E9C32-8D03-8FC9-5E02-EAE2026FDC63}
 3. CATHETER (CATHETER) {B0D1AAD7-012D-D3E8-3216-1F71A1EC685E}

UNITS:

=====

TERM: HAIR CONDITION Page: 184
ID: {7FA596FC-81CF-4906-BBA9-96ED54FAC55D}

=====

Abbreviation..... HAIR CONDI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Hair Condition
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

-
1. NORMAL (n) {553C9A97-A3D5-C52F-0E4E-C7D168C99A2B}
 2. ALOPECIA PARTIAL (ALOPE PART) {98EDBD8D-205A-4630-A09D-6B4222F5B2A2}
 3. ALOPECIA TOTAL (ALOPE TOTA) {ECDCE3C4-81F6-4FA1-BADA-97A20BD9A826}
 4. HIRSUTISM (HIRSUTISM) {A909A4FC-080A-4E9A-867A-0B30D4633602}

UNITS:

=====

TERM: HAND GRIP Page: 185
ID: {67289BB5-D6B1-53C5-53E2-6D209C1B0826}

=====

Abbreviation..... HAND GRIP
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Hand Grip

Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

- 1. ABSENT (a) {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0}
- 2. WEAK (w) {BE2D9305-E2DA-F492-66CF-164EE58B8179}
- 3. STRONG (s) {A12154E1-544E-F372-3839-0D44E5BC476B}
- 4. REFLEX (r) {F21585F9-5206-A6B3-2BF6-F416E5842366}
- 7. N/A (n) {B3DBE3FD-C5CB-3894-607A-BE13319594D0}

LOCATION: 3 entries

- 1. L HAND (L HAND) {639DF25A-19E3-B504-D323-4B3F6ED39A24}
- 2. R HAND (R HAND) {A54D3A80-8537-D079-9CCE-7DC56346D132}
- 3. BILATERAL HAND (BILAT HAND) {C782FE56-6117-493F-847F-83173407F0B5}

UNITS:

TERM: HEAD CIRCUMFERENCE Page: 186
ID: {33827E3C-5DBB-083C-D8BE-4DFD7D42071F}

Abbreviation..... HC
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... From Vitals
Display Name..... Head Circumference
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid..... 4536405

UNITS:

CENTIMETER 0 > 1 * 100 < 101 Dec Prec: 0

TERM: HEAD OF BED ELEVATION Page: 187
ID: {97D6AA11-1F1A-1853-23E5-25F2F41E07CF}

Abbreviation..... HOB E
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Head of Bed Elevation

4. STATED (Sta) {486CAE5A-BE56-06BA-C65E-1F8131DB1C2D}
{6362C194-B7F2-D38C-5B72-3BE5DE09438E}

UNITS:

INCH 0.00 > 1.00 * 96.00 < 100.00
Dec Prec: 2
CENTIMETER 0.0 > 1.0 * 2999.0 < 3000.0
Dec Prec: 1

=====

TERM: HOMANS SIGN Page: 190
ID: {996F1DCC-A6EE-4D52-B889-6E0E6185C74C}

=====

Abbreviation..... HOMAN SIGN
Active..... Yes
Boolean Value False..... Negative
Boolean Value True..... Positive
Data Type..... Boolean
Description..... No Description
Display Name..... Homan's Sign
Help Text..... Positive is calf pain with
Dorsiflexion and Plantar Flexion
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: I:E RATIO Page: 191
ID: {7A108E36-2701-4095-B1CE-FC51773D8D10}

=====

Abbreviation..... I:E RATIO
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Complex
Description..... Inspiratory Time to Expiratory Time
Ratio
Display Name..... I:E Ratio
Help Text..... Inspiratory Time to Expiratory Time
Ratio
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

Description..... Size of the IABP
Display Name..... Balloon Size
Help Text..... Size of the IABP
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

-
- | | |
|------------------|--|
| 1. 34 ml (34 ml) | {25952BB5-50B1-4CC0-9483-8AF891767045} |
| 2. 40 ml (40 ml) | {37A2998F-AB1C-4493-9C7A-FFCB05ED2302} |

UNITS:

=====

TERM: IABP CATHETER FLUSH

Page: 195

ID: {D034E8DA-1C70-44D2-B556-F21860827AF9}

=====

Abbreviation..... IABP CAT F
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... Was the Intra-aortic Balloon Pump
Flushed
Display Name..... IABP Catheter Flush
Help Text..... Was the Intra-aortic Balloon Pump
Flushed
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

-
- | | |
|------------|--|
| 1. NO (n) | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEA649CD} |

UNITS:

=====

TERM: IABP CIRCULATION ASSESSED

Page: 196

ID: {59A9ADB9-6630-41B3-AE5E-495AB232E287}

=====

Abbreviation..... IABP CIR A
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... IABP Circulation Assessed
Display Name..... IABP Circulation Assessed
Help Text..... IABP Circulation Assessed
Multi Select Picklist.....

Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: IABP CONDITION Page: 197
ID: {AA81E246-45C3-258F-4FC8-49CDD592B2D4}

=====

Abbreviation..... IABP CONDI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... IABP Condition
Help Text..... Condition of IABP Site
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

- 1. PATENT CLEAN (PAT CLEAN) {68688246-AE58-2EFB-BF9D-5E00D8CFDA8E}
2. RED (RED) {515D7CFB-3527-30E7-EC4B-B607C272DFA0}
3. SWOLLEN (SWOLLEN) {DCE86B10-C393-8DA9-6518-F2943F0E1A4E}
4. DISCOMFORT (DISCOMFORT) {77D058C9-1EB4-96C0-C9E3-666C840B884C}

UNITS:

=====

TERM: IABP DRESSING TYPE Page: 198
ID: {71923854-939D-C730-3A1B-F4A9FA02A6C6}

=====

Abbreviation..... IABP DRESS
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... IABP Dressing Type
Help Text..... Type of Invasive Arterial Blood
Pressure dressing
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

- 1. TRANSPARENT (TRANSPAREN) {0E73D24D-0EFC-08F9-196A-D30DD1BB3228}
- 2. PRESSURE (PRESSURE) {387AFB7F-3973-4D77-B2E6-0724F344D8F6}
- 3. GAUZE AND TAPE (GAU & TAPE)
{5BC40405-B5C4-47CF-81C3-520732E100AC}

LOCATION: 5 entries

- 1. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
- 2. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
- 3. R AXILLARY (R AXILLARY) {3044F28B-FED2-4069-7E3B-2A691C5771F9}
- 4. L AXILLARY (L AXILLARY) {21AF62CA-BEED-2A96-2CD6-F2350A01995E}
- 5. MEDIASTINAL (MEDIASTINA) {D99997D2-61EF-CAC2-A854-D906E3A97447}

UNITS:

=====

TERM: IABP PEAK DIASTOLIC AUGMENTED PRESSURE Page: 199
 ID: {24E18FB7-94D9-4F2D-9D6D-3CAF3C1F9951}

=====

Abbreviation..... IABP PDAP
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... Numeric
 Description..... BP Peak Diastolic Augmented Pressure
 Display Name..... IABP Peak Diastolic Augmented Pressure
 Help Text..... BP Peak Diastolic Augmented Pressure
 Multi Select Picklist.....
 Term Type..... OBSERVATION
 Vuid.....

UNITS:

MILLIMETERS OF MERCURY 0 > 80 * 120 < 200 Dec Prec: 0

=====

TERM: IABP RATIO Page: 200
 ID: {19DC7448-F44F-4E90-9FC0-233738821325}

=====

Abbreviation..... IABP RATIO
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... IABP Ratio
 Display Name..... IABP Ratio
 Help Text..... IABP Ratio
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 5 entries

- 1. 1:1 (1:1) {5C07BBD0-FDB9-0CD7-F2C2-616913E0281B}
- 2. 1:2 (1:2) {287CDCE0-36F3-65F6-E6A5-9EF7CF8976AE}
- 3. 1:3 (1:3) {062B4F93-8F48-7373-F679-11F57CE98BB3}
- 4. 1:4 (1:4) {3F1B75FF-C168-D413-0BDF-6581E7A4EEAD}
- 5. 1:8 (1:8) {94790C4C-3AF2-6F4C-58B8-87021CA97A53}

UNITS:

TERM: IABP SIZE Page: 201
 ID: {3C313604-F012-416F-A000-CCB92D04C30C}

Abbreviation..... IABP SIZE
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... Intra-aortic Balloon Pump Size
 Display Name..... IABP Size
 Help Text..... Intra-aortic Balloon Pump Size
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 5 entries

- 1. 7.5 (7.5) {7517EE9C-7C8B-4AF0-B60A-88F4E24BFA5E}
- 2. 8 (8) {518E7CFC-102E-351E-C7AF-FB575E3C8F4D}
- 3. 8.5 (8.5) {96723567-151D-4B5C-9C65-B3EB6D584272}
- 4. 9 (9) {58FF2FD4-569F-F454-9F5B-FE6F9E34D8CC}
- 5. 9.5 (9.5) {252555A1-ADB2-4948-8764-C952F2587DC5}

LOCATION: 5 entries

- 1. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
- 2. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
- 3. L AXILLARY (L AXILLARY) {21AF62CA-BEED-2A96-2CD6-F2350A01995E}
- 4. R AXILLARY (R AXILLARY) {3044F28B-FED2-4069-7E3B-2A691C5771F9}
- 5. MEDIASTINAL (MEDIASTINA) {D99997D2-61EF-CAC2-A854-D906E3A97447}

UNITS:

TERM: IABP TIMING Page: 202
 ID: {443F4932-9C42-EADA-4A9A-8F8956DF2CE2}

Abbreviation..... IABP TIMIN
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....

Data Type..... PickList
Description..... IABP Timing
Display Name..... IABP Timing
Help Text..... IABP Timing
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

- 1. 1:1 (1:1) {5C07BBD0-FDB9-0CD7-F2C2-616913E0281B}
- 2. 1:2 (1:2) {287CDCE0-36F3-65F6-E6A5-9EF7CF8976AE}
- 3. 1:3 (1:3) {062B4F93-8F48-7373-F679-11F57CE98BB3}
- 4. 1:4 (1:4) {3F1B75FF-C168-D413-0BDF-6581E7A4EEAD}
- 5. 1:8 (1:8) {94790C4C-3AF2-6F4C-58B8-87021CA97A53}

UNITS:

TERM: IABP TIMING CHECK Page: 203
ID: {B2BA8A27-B578-40B6-83A0-77D95C49B59B}

Abbreviation..... IABP TI CH
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... Was The IABP Timing Checked?
Display Name..... IABP Timing Check
Help Text..... Was The IABP Timing Checked?
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
- 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

TERM: IABP UAEDP Page: 204
ID: {6CB7493C-6175-4FFA-AEA9-C60CF2EE0EE4}

Abbreviation..... IABP UAEDP
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... IABP Unassisted Aortic End Diastolic
Pressure

=====
Abbreviation..... ALINE TYPE
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Arterial Line Type
Display Name..... Arterial Line Type
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

UNITS:

=====
TERM: INACTIVE - AUDIOMETRY
ID: {FFD29134-8BB2-248E-0412-93C2C08B076F}

Page: 208

=====
Abbreviation..... AUD
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Audiometry
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void..... 4688718

UNITS:

=====
TERM: INACTIVE - BISPECTRAL INDEX
ID: {4819B466-D0D0-4851-A434-E1691ECA1E65}

Page: 209

=====
Abbreviation..... BIS
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Bispectral Index
Display Name..... Bispectral Index
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - BODY MASS INDEX
ID: {54398789-EC32-4F1D-8E02-913FC96FD198}

=====

Page: 210

Abbreviation..... BMI
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... BMI
Help Text..... Body Mass Index
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - CARDIAC AV FISTULA BRUIT
ID: {A0E99499-8021-47FF-AEA2-EB33B9648B0E}

=====

Page: 211

Abbreviation..... AV FIST BR
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... AV Fistula Bruit
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - CATHETER TYPE
ID: {01941D20-3FD7-A4BF-6108-712C7EBC33A1}

=====

Page: 212

Abbreviation..... CATHETER T
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Catheter Type
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - CENTRAL LINE SITE
ID: {D6FC875B-CC24-2981-B85E-07FF9CD798AB}

=====

Page: 213

Abbreviation..... C LINE SIT
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Central Line Site
Help Text..... Site of Central Line
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - CHEST TUBE DRESSING INTACT
ID: {525B2489-75C2-4AB3-9441-A457BCD7B257}

=====

Page: 214

Abbreviation..... CH TU DR I
Active..... No
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... Chest Tube Dressing Intact
Display Name..... Chest Tube Dressing Intact
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - CHEST TUBE SITE
ID: {C5224B54-65A5-2010-F9E3-1C4570500A11}

=====

Page: 215

Abbreviation..... CHEST TU S
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... SITE OF CHEST TUBE
Display Name..... Chest Tube Site

Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - COUGH AND DEEP BREATH EVENTS Page: 216
ID: {5F7D2087-018E-4DA4-8B36-02743F65A769}

=====

Abbreviation..... C&DB EVENT
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Cough and Deep Breath Events
Display Name..... Cough and Deep Breath Events
Help Text..... Cough and Deep Breath Events
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - DATE Page: 217
ID: {5BB40FB6-173F-9656-8D5B-7A2F5A5F3B82}

=====

Abbreviation..... DATE
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... Complex
Description..... Todays Date
Display Name..... Date
Help Text..... Date: Year, Date: Month, Date: Day
are all available
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - DATE DAY Page: 218
ID: {CD91474C-1EB3-304B-4E3A-B85F6976C1E7}

=====

Abbreviation..... DATE DAY
Active..... No

Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... The day of the month
Display Name..... Date: Day
Help Text..... The day of the month
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - DATE MONTH Page: 219
ID: {D4E1590F-7B12-FD7C-B545-9CF83BA68BAC}

=====

Abbreviation..... DATE MONTH
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Month of the year
Display Name..... Date: Month
Help Text..... Month of the year
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - DATE YEAR Page: 220
ID: {57029BBB-5904-DE1D-F58D-AF13A890C6ED}

=====

Abbreviation..... DATE YEAR
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... The Year
Display Name..... Date: Year
Help Text..... The Year
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - DIALYSIS CATHETER SITE Page: 221

ID: {FC099337-E02E-84E6-7338-E423ECD3DDCB}

Abbreviation..... DIAL CTH S
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Dialysis Catheter Site
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

UNITS:

TERM: INACTIVE - EDEMA LOCATION

Page: 222

ID: {8C4F762B-FD4A-CB97-326F-74C4A2ED5548}

Abbreviation..... EDEMA LOCA
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Location of identified edema
Display Name..... Edema Location
Help Text..... Location of identified edema
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

UNITS:

TERM: INACTIVE - ENTERAL NUTRITION FORMULA

Page: 223

ID: {1C7EF762-4D93-438B-8CBD-2CB98E6C8686}

Abbreviation..... ENT NUT FO
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Enteral Nutrition Formula Used
Display Name..... Enteral Nutrition Formula
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

UNITS:

Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - IABP SITE Page: 227
ID: {EA276ABE-52D0-6D4D-1B0A-7E6443BB264F}

=====

Abbreviation..... IABP SITE
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... IABP Site
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - INCISION LOCATION Page: 228
ID: {8A9D3611-6984-4FA1-9C16-4DC215F9B4D6}

=====

Abbreviation..... IN LOCATIO
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Incision Location
Display Name..... Incision Location
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - INCISION SITE Page: 229
ID: {9DBE601B-4D61-D280-1F13-0FDE3533B0A4}

=====

Abbreviation..... WOUND INCS
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList

Abbreviation..... IV DATE
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... IV Insertion Date
Help Text..... Date IV Inserted
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

TERM: INACTIVE - IV SITE

Page: 233

ID: {7F252F53-B53D-D8BD-097B-78197EE99FBA}

Abbreviation..... IV SITE
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Site of IV
Display Name..... IV Site
Help Text..... Site of IV
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

UNITS:

TERM: INACTIVE - OR OUTPUT - CHEST TUBE

Page: 234

ID: {F6ECA901-354B-4C5B-B03F-C17F7F6448B6}

Abbreviation..... OR OUT CHE
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... OR Output - Chest Tube
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - OR OUTPUT - DRAIN
ID: {1A99C3C7-9FE1-46D3-9A32-57A2C735DFC6}

=====

Page: 235

Abbreviation..... OR OUT DRA
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... OR Output - Drain
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

=====

TERM: INACTIVE - OR OUTPUT - PUMP LOSS
ID: {6A424E02-67FF-49F6-9B33-853BE8545B9B}

=====

Page: 236

Abbreviation..... OR OUT PUM
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... OR Output - Pump Loss
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

=====

TERM: INACTIVE - OR OUTPUT - URINE
ID: {3ECAAD17-1BA5-4D0A-BEDC-EC512D32EFF9}

=====

Page: 237

Abbreviation..... OR OUT URI
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... OR Ouput - Urine
Help Text..... contributes to output total value
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

=====

TERM: INACTIVE - OSTOMY TYPE
ID: {975E1FD8-521C-4818-9975-98A30B12C1EF}

=====

Page: 238

Abbreviation..... OSTOMY TYP
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Ostomy Type
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - OXYGEN DELIVERED
ID: {CFD876D1-AC28-F63F-3B6C-680F5B66A1BD}

=====

Page: 239

Abbreviation..... D02
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Oxygen Delivered
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - OXYGEN DELIVERY METHOD
ID: {F320E0CF-6BAC-3B6D-F5D7-21253FB88A4F}

=====

Page: 240

Abbreviation..... O2 DEL MET
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Type of equipment used to deliver
oxygen to patient
Display Name..... Oxygen Delivery Method

Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - OXYGEN EXTRACTION RATIO Page: 241
ID: {D6C6849F-6134-417F-AD8B-7853F5AE3451}

=====

Abbreviation..... O2ER
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... Complex
Description..... Ratio of Oxygen Consumed (VO2) to
Oxygen delivered (DO2)
Display Name..... O2ER
Help Text..... Ratio of Oxygen Consumed (VO2) to
Oxygen delivered (DO2)
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - PA CATHETER INSERTION DATE Page: 242
ID: {64ADD596-950E-EDB5-D7DE-38F5F8601D9E}

=====

Abbreviation..... PA INS DAT
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... PA Catheter Insertion Date
Help Text..... Date PA Catheter Inserted
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - PA CATHETER SITE Page: 243
ID: {1E2B676E-24D4-28D6-9EC4-2200F3DA31B6}

=====

Abbreviation..... PA SITE

Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... PA Catheter Site
Help Text..... Site of PA Catheter
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - PACEMAKER DATE INSERTED
ID: {B395D1F6-2788-CA17-618C-62719C1CAF4A}

=====

Page: 244

Abbreviation..... PDI
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Pacemaker Date Inserted
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - PACEMAKER SITES
ID: {668C0AE7-061F-C573-A119-8ED8D006F018}

=====

Page: 245

Abbreviation..... PMSI
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Pacemaker Site
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - PAIN LOCATION
ID: {20836431-E2C3-41CA-BE14-A8FB1F6A8637}

Page: 246

Abbreviation..... PAIN LOC
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Location of patients Pain
Display Name..... Pain Location
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Void.....

UNITS:

TERM: INACTIVE - POST OPERATION DAY NUMBER
ID: {82FF0CE1-48CD-4AA8-6292-CB6879B8CD75}

Page: 247

Abbreviation..... POD#
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... NUMBER OF DAYS SINCE OPERATION
Display Name..... Post Op Day #
Help Text..... Number of days since operation.
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

TERM: INACTIVE - PRESSURE ULCER SITE
ID: {580747C6-F7B1-4CCA-AE24-E82C436746B6}

Page: 248

Abbreviation..... PR UL SITE
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Pressure Ulcer Site
Display Name..... Pressure Ulcer Site
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - PRESSURE ULCER VISIBLE STRUCTURE
ID: {C911EB61-F151-42C2-A6EC-1EA357EA6F1D}

=====

Page: 249

Abbreviation..... PR UL VISI
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Pressure Ulcer Visible Anatomical
Structure
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

UNITS:

=====

TERM: INACTIVE - PULSE OXIMETRY
ID: {9B830A66-0278-20E1-DB38-089A952DF381}

=====

Page: 250

Abbreviation..... PO2
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... O2 saturation
Display Name..... Pulse Oximetry
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid..... 4500637

UNITS:

=====

TERM: INACTIVE - SKIN DISCHARGE COLOR
ID: {C4633087-CBC5-4924-ADF6-846D13624D3E}

=====

Page: 251

Abbreviation..... SKI DIS CO
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Skin Discharge Color
Help Text..... No Help Text

Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - SKIN DISCHARGE ODOR
ID: {4E9309BD-6F1E-4D2C-B500-CE724A92A0F6}

=====

Page: 252

Abbreviation..... SKI DIS OD
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Skin Discharge Odor
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - SKIN DRAINAGE CHARACTERISTICS
ID: {9F12C134-BA23-40C9-8CC9-0AC08D8891CA}

=====

Page: 253

Abbreviation..... SKI DRA CH
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Skin Drainage Characteristics
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - VENTILATOR IN-LINE NEBULIZER
ID: {12804E1D-9B65-4027-BCD4-1A28C9C57730}

=====

Page: 254

Abbreviation..... VENT I-L N
Active..... No
Boolean Value False..... Off
Boolean Value True..... On

Data Type..... Boolean
Description..... Ventilator In-Line Nebulizer
Display Name..... Ventilator In-Line Nebulizer
Help Text..... Ventilator In-Line Nebulizer
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - WOUND SITE

Page: 255

ID: {7BFBDAAF-779E-4C7A-82F3-DA9E0E7DD597}

=====

Abbreviation..... WOU SITE
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Wound Site
Display Name..... Wound Site
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INACTIVE - WOUND TOPICAL TREATMENT TYPE

Page: 256

ID: {29E59078-A908-496E-ABC7-B2C48A4E6400}

=====

Abbreviation..... WOU TOP TR
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Topical Treatment
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

UNITS:

=====

TERM: INCENTIVE SPIROMETRY

Page: 257

ID: {C1E6301C-CB9C-4909-9500-9AA01C42583D}

=====

Abbreviation..... INC SPIROM
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Incentive Spirometry
Display Name..... Incentive Spirometry
Help Text..... Incentive Spirometry
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

- 1. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}
- 2. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
- 3. N/A (n) {B3DBE3FD-C5CB-3894-607A-BE13319594D0}

UNITS:

=====

TERM: INCISION DISCHARGE CHARACTERISTICS
ID: {14342EB1-1961-C52E-D00A-DA5E7509C4FD}

=====

Page: 258

Abbreviation..... INC DIS CH
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Incision Discharge Characteristics
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 8 entries

- 1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
- 2. BLOODY (BLOODY) {D7904E84-8331-D6DD-F2C0-7F3039C5746A}
- 3. SEROUS (SEROUS) {9BDA347E-29E9-80C7-91FA-016192D4257E}
- 4. SEROSANGUINOUS (SEROSANGUI)
{19A8D332-899D-177A-2247-ABCDAAEE3BFA8}
- 5. PURULENT (PURULENT) {86C01494-743A-C6A2-8883-01CA5200C79B}
- 6. CLEAR (c) {8981C365-2B46-704F-C9F4-975DF1C23C2D}
- 7. CLOUDY (CLOUDY) {3CCFCCCC-00FC-6B00-1D58-A77C02FDB401}
- 8. VISCOUS (VISCOUS) {0E8019D5-F2DB-46B8-B750-41568DD0179C}

UNITS:

=====

TERM: INCISION DRAINAGE AMOUNT

Page: 259

ID: {77AAAB12-979C-2624-5CFE-7FD5F2310AD5}

Abbreviation..... I DR AMOUN
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Amount of incision drainage
Display Name..... Incision Drainage Amount
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

-
1. INCISION DRAINAGE AMOUNT NONE (I DR A NON)
{5D0D9027-4B20-43A4-B437-EA7B77CB2105}
 2. INCISION DRAINAGE AMOUNT SCANT (I DR A SCA)
{A476E395-46AA-44E7-8A00-CE22C7066BB1}
 3. INCISION DRAINAGE AMOUNT SMALL (I DR A SMA)
{4479479D-8ADA-4FBF-A5A2-426D49FC94D3}
 4. INCISION DRAINAGE AMOUNT MODERATE (I DR A MOD)
{9B0487A9-DCB4-42D3-92E6-805C412B5D34}
 5. INCISION DRAINAGE AMOUNT HEAVY (I DR A HEA)
{D77FCEA0-D247-4DD5-873E-2F95C7FF48A6}

QUALITY: 1 entries

-
1. ESTIMATED (E) {9EAFF26B-7D87-DA66-F991-022D095FB5A0}

UNITS:

TERM: INCISION DRESSING PREVENTS ASSESSMENT
ID: {4BE5369B-FBED-40B9-8950-AF370A31A189}

Page: 260

Abbreviation..... IN DR PREV
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... Incision Dressing Prevents Full
Assessment
Display Name..... Incision Dressing Prevents Full
Assessment
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

-
1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}

2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEE649CD}

UNITS:

=====

TERM: INCISION DRESSING TYPE

Page: 261

ID: {AA1F8811-85BF-187A-7A1D-D55A2B523190}

=====

Abbreviation..... INC DRE TY
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... OBSERVATION OF INCISION WOUND DRESSING
Display Name..... Incision Dressing Type
Help Text..... Incision Dressing Type
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 15 entries

1. INCISION DRESSING TYPE NONE (INC DR NON)
{81F07C8D-A463-48A0-A64A-A143CECB5DD6}
2. INCISION DRESSING TYPE PRESSURE (INC DR PRE)
{B5823373-E5D8-40A5-954C-6A193D7E7511}
3. INCISION DRESSING TYPE ANTIMICROBIAL (INC DR ANT)
{35F5F5C1-B76A-44D2-8576-9D7D3ECA3197}
4. INCISION DRESSING TYPE GAUZE (INC DR GAU)
{E4AC9307-0512-4F9D-A254-FE19E654461F}
5. INCISION DRESSING TYPE TRANSPARENT (INC DR TRA)
{E198EDAD-AA67-49F9-A065-6BB11FE39ED6}
6. INCISION DRESSING TYPE ANTIMICROBIAL DISK (INC DR ADI)
{55A33062-7D6B-4CDF-9131-D4ADC38F6CE3}
7. INCISION DRESSING TYPE GAUZE WRAP (INC DR GWR)
{E1DBEC29-1430-4DA0-99DC-49B9389D2440}
8. INCISION DRESSING TYPE FOAM DRESSING (INC DR FOA)
{5C2C2D24-6022-4BBF-9065-03FB5971CB3F}
9. INCISION DRESSING TYPE GAUZE PACKING (INC DR GPA)
{3561DAEB-8A30-4E16-B920-F7DFE022EBB1}
10. INCISION DRESSING TYPE HYDROCOLLOIDAL (INC DR HYD)
{424CCDCA-5196-4BBD-AE8B-0F68DDCC3FCF}
11. INCISION DRESSING TYPE CALCIUM ALGINATE (INC DR CAL)
{B20558BA-FB9E-499E-82DF-49616ACE9375}
12. INCISION DRESSING TYPE HYPERTONIC HYDROGEL (INC DR HYP)
{25820A63-AD76-4198-BEB1-EF95DDEC5B35}
13. INCISION DRESSING TYPE HYDROGEL (INC DR HYD)
{8CBAC6DF-F51E-456B-8F9A-A68E4452812C}
14. INCISION DRESSING TYPE BALSAM OF PERU (INC DR BAL)
{A1FB45C6-6315-4282-8D89-DCB0C018E357}
15. INCISION DRESSING TYPE BISMUTH GAUZE (INC DR BIS)
{CA88C5DC-3890-49A2-800A-0528CE194BA8}

UNITS:

=====

TERM: INCISION EDGE

Page: 262

ID: {3FEA28A9-1C53-47D7-95F9-BFFFFB25576C6}

=====

Abbreviation..... INC ED
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Incision Edge
Display Name..... Incision Edge
Help Text..... Incision Edge
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

1. INCISION EDGE INDISTINCT/ATTACHED (INC ED ATT)
{AF7E7150-525E-49A4-AC23-042B788E8805}
2. INCISION EDGE NON-ATTACHED (INC ED NON)
{F457F951-F1DC-47AC-9336-AC093E149862}
3. INCISION EDGE TUNNELING (INC ED TUN)
{492CD783-3C32-40F0-AE35-5A89E90409A7}
4. INCISION EDGE UNDERMINING (INC ED UND)
{0F3E968A-B867-4123-8C05-ED6B12D364CD}
5. INCISION EDGE ROLLED (INC ED ROL)
{FFA33CAE-5D8B-4501-8120-A536BE53296F}

UNITS:

=====

TERM: INCISION ERYTHEMA

Page: 263

ID: {D87F3C97-622C-43D3-9C7F-55C2C848941B}

=====

Abbreviation..... INC ERYTHE
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Incision Erythema
Display Name..... Incision Erythema
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
2. LIMITED TO INCISION (LIM TO INC)

- 3. EXTENDS INTO PERIWOUND <1CM (EX IN PER<)
 {81A61254-55B6-47E1-8BE2-F4A47F624BC6}
 {0EBCEAE0-E17A-484F-BFE2-51F17AA761CB}
- 4. EXTENDS INTO PERIWOUND>1CM (EX IN PER>)
 {1F31FEBB-DA19-434C-8C58-67CDC5385231}

UNITS:

=====

TERM: INCISION ODOR Page: 264
 ID: {2584F1F9-CD15-450C-BDDDB-D0B2B1924303}

=====

Abbreviation..... INC OD
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... No Description
 Display Name..... Incision Odor
 Help Text..... No Help Text
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 2 entries

- 1. INCISION ODOR NONE (INC OD NON)
 {CED90890-6FBC-4A4A-B0B8-6F75F0DD567A}
- 2. INCISION ODOR FOUL (INC OD FOU)
 {ECDC0988-BE82-41E3-8104-A66A5F8A1093}

UNITS:

=====

TERM: INCISION PRESENT Page: 265
 ID: {CC0DD035-AAD2-4C9F-B971-E13B2FF0DB9C}

=====

Abbreviation..... INCIS PRES
 Active..... Yes
 Boolean Value False..... No
 Boolean Value True..... Yes
 Data Type..... Boolean
 Description..... Incision Present
 Display Name..... Incision Present
 Help Text..... No Help Text
 Multi Select Picklist.....
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}

2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEE649CD}

LOCATION: 85 entries

-
1. L THORACOTOMY (L THORACOT) {313C11A1-559A-4F68-8718-4012416EC189}
 2. R THORACOTOMY (R THORACOT) {A18E4BB3-0718-48D6-BD49-720A70C0EF60}
 3. L TRAPEZIUS (L TRAPEZIU) {B74CF943-C10F-4A1B-A19A-4975E674A97A}
 4. R TRAPEZIUS (R TRAPEZIU) {69753B67-AFA1-4C93-BA78-E266F9894E46}
 5. DONOR SITE (DONOR SITE) {F0A65974-7650-4A4B-AC88-0945E37A0AD2}
 6. L GLUTEUS (L GLUTEUS) {88C4E547-04CD-4D52-9416-CAEB6D898D21}
 7. R GLUTEUS (R GLUTEUS) {B4D32045-68D1-4457-9106-7937887FA6F6}
 8. L LOWER ANTERIOR CHEST (L L A CHES) {6073B3C6-AA4D-4BBA-9FF2-32ED40CF92AD}
 9. R LOWER ANTERIOR CHEST (R L A CHES) {96BCF731-B661-4455-B5E7-C4AE2509505F}
 10. L MIDDLE ANTERIOR CHEST (L M A CHES) {36A83F77-117E-474B-AE34-0B2C02E131AE}
 11. R MIDDLE ANTERIOR CHEST (R M A CHES) {2F751A9C-66CA-4C22-ACC9-6E81FA63FA2B}
 12. L UPPER ANTERIOR CHEST (L U A CHES) {FE63B58A-DE39-4C03-AE61-E34A5487C767}
 13. R UPPER ANTERIOR CHEST (R U A CHES) {67890049-A4A8-4241-BDF3-573729BD6EB5}
 14. L LOWER POSTERIOR CHEST (L L P CHES) {2C19F10B-7546-4A1F-A5DC-3867359ED971}
 15. R LOWER POSTERIOR CHEST (R L P CHES) {5163180E-FFD2-4BFB-BF8B-4C6B18EFE5D4}
 16. L MIDDLE POSTERIOR CHEST (L M P CHES) {5A656939-F82B-4DEC-9F1D-B5FA644E3DA9}
 17. R MIDDLE POSTERIOR CHEST (R M P CHES) {2BC487D6-6CE8-499E-94A0-A79EAA8E199B}
 18. L UPPER POSTERIOR CHEST (L U P CHES) {6D5674E9-EF04-4EF1-84A0-3F5ED2FB532A}
 19. R UPPER POSTERIOR CHEST (R U P CHES) {D8F03FA5-1FD7-46EE-8FA4-17B3B7E62141}
 20. MIDLINE ABDOMINAL (MID ABDOMI) {37ED2947-9BDA-4207-94F1-6B8C0FE12ECF}
 21. UMBILICAL (UMBILICAL) {8874B300-0751-49AF-8211-CF4C009C68B5}
 22. TRANSVERSE ABDOMINAL (TRANS ABDO) {1C10FF98-7B20-45D0-B420-714981832B3D}
 23. L ANKLE (L ANKLE) {D5A692EF-9D7C-473C-AF96-7EF1C01FA2AF}
 24. R ANKLE (R ANKLE) {4E6DF352-F399-4CE9-908C-48132C3D856B}
 25. L LOWER ARM (L LOW ARM) {63AA95E8-256D-4390-8CB3-F00A1B3F91CD}
 26. R LOWER ARM (R LOW ARM) {8DE6E277-EFFF-40DD-9E83-D7295AACB6AA}
 27. L UPPER ARM (L UPP ARM) {5F915696-1B08-32E9-6B95-49246117D660}
 28. R UPPER ARM (R UPP ARM) {3EF20E8E-F68C-E4F2-8F8C-FAD3B44DDAA4}
 29. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
 30. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
 31. L INGUINAL (L INGUINAL) {627D172E-8F33-48C9-BE89-D036E07D46BF}
 32. R INGUINAL (R INGUINAL) {395233DF-3B89-49E0-946B-EB5C02E51178}
 33. L FLANK (L FLANK) {6E453E8D-B456-4BFD-A015-F9C013830F30}
 34. R FLANK (R FLANK) {F5700ED2-646F-4B75-A55B-111C27FB939C}
 35. L FOOT (L FOOT) {2EC3FBEC-4EA2-9E7A-17D4-FDC2F9D28046}
 36. R FOOT (R FOOT) {7D9F5FF6-232C-2488-4504-A721FF06928F}

37. L HIP (L HIP) {786CD881-27AD-4136-B9CC-8423CB38A431}
38. R HIP (R HIP) {A48CBCAE-5404-4934-9321-ECEF660A5350}
39. L KNEE (L KNEE) {3A9E30B4-09A9-967A-58EF-69A01EC4F859}
40. R KNEE (R KNEE) {BCD18D53-07B6-2C41-64FC-3A07A7E979E7}
41. L THIGH (L THIGH) {16C2312F-1E0E-45E0-B675-160F40C23F57}
42. R THIGH (R THIGH) {B7A2D4A0-801E-4578-967A-F1CFB86FB96B}
43. L CALF (L CALF) {3DAA7817-936B-444F-8A2D-D3AB09144C10}
44. R CALF (R CALF) {0298A452-4901-44BF-A7A8-F42DD3BA8B63}
45. L LOWER LEG (L LOW LEG) {7B8E92FA-E065-4D06-9FB6-4234656CFD61}
46. R LOWER LEG (R LOW LEG) {C8EA4D71-D97E-46BD-AE9D-C10954C525CD}
47. MIDLINE NECK (MID NECK) {FF666713-278D-46D6-8D25-E5D57E78EA08}
48. L NECK (L NECK) {C4EA6279-05E2-4195-A638-3794DCBA7B78}
49. R NECK (R NECK) {7E4D84D3-535B-6164-C6C1-837C6F997A0D}
50. PERINEAL (PERINEAL) {1801441C-D5EF-4775-9E58-E4985496F7E0}
51. SPINE CERVICAL (SPINE CERV)
{9D28E1BF-D014-4DA4-A832-1F3C25C148E9}
52. SPINE THORACIC (SPINE THOR)
{2CACBB4F-3EF3-4FBF-AC3E-41D72F0178CB}
53. SPINE LUMBAR (SPINE LUMB) {E7E411AB-5C29-4906-B0A9-E8658AAFFF91}
54. L SUBCLAVIAN (LSCL) {71B93594-15F9-AA4A-AEF2-E222EC28C00D}
55. R SUBCLAVIAN (R SUBCLAVI) {53DEAE83-93AE-123B-14E5-4DE38CD4F752}
56. L THORACIC ANTERIOR (L THOR ANT)
{C85B5FB9-153B-4A8B-9BE6-5524A572ACAC}
57. R THORACIC ANTERIOR (R THOR ANT)
{BD3DF757-2C65-4D6A-9C80-5DC95BB41CFB}
58. L THORACIC POSTERIOR (L THOR POS)
{68E3381F-7605-4786-85A8-523D215A2B97}
59. R THORACIC POSTERIOR (R THOR POS)
{E91085D2-E665-4EB6-A086-8CF573DECD8F}
60. L THORACIC LATERAL (L THOR LAT)
{48D26CEB-2E33-4F57-AB3A-F5A46108FFBF}
61. R THORACIC LATERAL (R THOR LAT)
{55C48294-D3F3-4013-B552-BA3344D9513B}
62. L THORACOTOMY (L THORACOT)
{313C11A1-559A-4F68-8718-4012416EC189}
63. R THORACOTOMY (R THORACOT)
{A18E4BB3-0718-48D6-BD49-720A70C0EF60}
64. THYROID (THYROID) {6AC68004-FE58-46D8-9A1E-4EEF72AAD71B}
65. L PARIETAL (L PARIETAL) {2687ED24-F35E-456C-99C9-F6A8EC4364E0}
66. R PARIETAL (R PARIETAL) {28C927CD-1342-483B-B366-0B7EA8ACABAD}
67. L ORBIT (L ORBIT) {1D798E35-5095-4BA0-8F9F-FA96668A5239}
68. R ORBIT (R ORBIT) {DFFD974A-1684-49AE-8230-8E12411EA62E}
69. L MANDIBULAR (L MANDIBUL) {CB9DB373-5F0A-446B-909F-F76392CB3A04}
70. R MANDIBULAR (R MANDIBUL) {303A1DE2-CB9A-498F-9CC5-DFA8FBFF10C1}
71. L ORAL (L ORAL) {057B787A-8750-4563-80DB-25110CCC1146}
72. R ORAL (R ORAL) {9CEEF934-AA3D-4BDE-BC3B-9E275C51CC9C}
73. L FACE (L FACE) {06F00A70-F475-CD40-665D-AF0569A1E75F}
74. R FACE (R FACE) {D1B15677-EF6F-0125-3601-68C1FE6386CE}
75. L EAR (L EAR) {DA9C2DD1-D6FF-4899-A62B-1096FAA36C59}
76. R EAR (R EAR) {AE31327D-A03C-488C-8A8E-084E7E6E7FDE}
77. L THUMB (L THUMB) {F82F4731-2CC5-4705-91BE-05E40C25B419}
78. R THUMB (R THUMB) {B414BD68-4155-45B3-985D-B2F8252A20FE}
79. L FINGER (L FINGER) {BB825645-1630-498C-A31B-6728D422EEFC}
80. R FINGER (R FINGER) {E12CA0EF-3267-4B93-8643-230E2226C3D6}
81. L TOE (L TOE) {C666DE86-E7AE-4AB3-8AE1-91496219E63D}
82. R TOE (R TOE) {B65727AB-3205-4DCE-833F-6DA8554636F9}

- 83. L OCCIPITAL (L OCCIPTAL) {A4234F41-CD6F-4CE1-9210-316DB847BDE2}
- 84. R OCCIPITAL (R OCCIPITA) {5DB951DB-CB03-482E-B0DA-CC5EE8FB1A67}
- 85. OCCIPITAL (OCCIPITAL) {5A3E384B-1C87-4A19-A449-DB91AB443FFC}

UNITS:

=====

TERM: INCISION STAPLES INTACT Page: 266
 ID: {FA75C02B-3698-4D85-A293-0ACC8C1D3252}

=====

Abbreviation..... IN STA INT
 Active..... Yes
 Boolean Value False..... No
 Boolean Value True..... Yes
 Data Type..... Boolean
 Description..... Incision Stapes Intact
 Display Name..... Incision Staples Intact
 Help Text..... No Help Text
 Multi Select Picklist.....
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
- 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: INCISION SUTURES INTACT Page: 267
 ID: {A5C84A03-633B-4A9B-891B-292449B8C409}

=====

Abbreviation..... IN SUT INT
 Active..... Yes
 Boolean Value False..... No
 Boolean Value True..... Yes
 Data Type..... Boolean
 Description..... Incision Sutures Intact
 Display Name..... Incision Sutures Intact
 Help Text..... No Help Text
 Multi Select Picklist.....
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
- 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: INCISION TENDER TO PALPATION
ID: {C6ACAC1F-0AEF-4F1B-A487-35A92207AE22}

=====

Page: 268

Abbreviation..... IN TENDER
Active..... Yes
Boolean Value False..... Patient Denies
Boolean Value True..... Appropriate
Data Type..... Boolean
Description..... Incision Tender To Palpation
Display Name..... Incision Tender To Palpation
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

1. NO (n)	{FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y)	{483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: INFUSION SET CHANGED
ID: {BA138A48-8B25-40CF-8345-2AD9E12CC840}

=====

Page: 269

Abbreviation..... INF SET CH
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... Was the infusion set changed?
Display Name..... Infusion Set Changed
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

1. NO (n)	{FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y)	{483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: INTAKE
ID: {062BB93E-EECB-A074-63D7-EB95E176D6A3}

=====

Page: 270

- 17. PCA (PCA) {D93028EE-9A44-43F5-9C0C-A903BD9B6E6D}
- 18. PLATELETS PHERESIS (PLATELET P) {65CBDCBD-99DB-44FE-9230-D62886208EBF}
- 19. PLATELETS PHERESIS IRRADIATED (PLATEL P I) {70C7155A-5AB6-4733-B43F-606F8BBB3B4A}
- 20. PPN (PPN) {BF402283-3E6F-4DF1-9CC6-4DE35D755C7C}
- 21. PRBC (PRBC) {BCE37EFD-78D6-4A13-A5E4-EA278BF18E02}
- 22. PRBC - AUTOLOGOUS (PRBC-A) {1580990C-0B35-4DDD-84AA-682B341C9B76}
- 23. PRBC - LEUKOREduced (PRBC-L) {13DB3170-6C26-4178-B429-C02AFC541AC7}
- 24. PRBC - LEUKOREduced - DIRECTED (PRBC-L-D) {9C54B7F7-FBD9-4FA8-B2A1-893A25A9CD2C}
- 25. PRBC - LEUKOREduced - IRRADIATED (PRBC-L-I) {B4431519-0102-49ED-94D6-1738E375C5A4}
- 26. PRBC - LEUKOREduced - WASHED (PRBC-L-W) {BB58F35A-1086-4061-B485-FA35DFB79A0D}
- 27. SANDOSOURCE PEPTIDE (SAND PEPTI) {B543461C-7316-4750-9918-B30E96AB3C95}
- 28. TPN (TPN) {55341B62-E356-4C38-A52D-4F196751E9EB}
- 29. WATER (WATER) {ACCC998E-0F64-510D-B843-3D30B2AD67CF}
- 30. WBC - WASHED (WBC WASHED) {F1460F94-203C-473D-AC7C-92E8ADFFC0F5}
- 31. CRYSTALLOID (CRYSTALLOI) {F8B0C9B5-F386-4C88-896A-092C4A95A031}
- 32. PLATELETS (PLATELETS) {3F4FA2EF-DF05-4A6D-839F-F97D3C2AFCF7}
- 33. WHOLE BLOOD (WHOL BLOOD) {144EE01A-16B5-4F9B-9A1C-D420D0A9A2D0}

UNITS:

MILLILITER -1 > 0 * 9999 < 10000 Dec Prec: 0
OUNCE -1 > 0 * 64 < 350 Dec Prec: 0
CUBIC CENTIMETER -1 > 0 * 9999 < 10000 Dec Prec: 0

=====

TERM: INTAKE - DIALYSIS Page: 271
ID: {CE3D8DBD-278C-4D3D-AA4F-2E5800C17DBF}

=====

Abbreviation..... INT DIALYS
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Dialysis Intake
Display Name..... Dialysis Intake
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 4 entries

-
- 1. HEMODIALYSIS (HEMODIALYS) {DE7CD394-9E69-46B0-A207-5A18F7C5A1A7}
 - 2. PERITONEAL DIALYSIS (PERI DIALY) {6B4AA38B-AFCC-42CE-9F32-59B32E8A6BBB}
 - 3. CVVHD (CVVHD) {4CC9AB0D-0B6C-43F9-A93E-A8228EC3A77F}
 - 4. CRRT (CRRT) {AB37FC1D-DD7A-4103-A777-7B8D0F876889}

UNITS:

MILLILITER -1 > 0 * 9999 < 10000 Dec Prec: 0
CUBIC CENTIMETER -1 > 0 * 9999 < 10000 Dec Prec: 0

=====

TERM: INTAKE - EF/SUPPLEMENT CALORIES

Page: 272

ID: {58506F9F-9D10-47F6-A195-5CE7D749FFBC}

=====

Abbreviation..... INT EF CAL
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Enteral Feeding and Supplement Calories
Display Name..... Intake - EF/Supplement Calories
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

CALORIE -1 > 0 * 500 < 2001 Dec Prec: 0

=====

TERM: INTAKE - ENTERAL FEEDING TUBE

Page: 273

ID: {769017E8-B98E-45DC-9706-FE6475B7FF8D}

=====

Abbreviation..... INT EF TU
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Intake - EF Tube
Display Name..... Intake - EF Tube
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

MILLILITER -1 > 0 * 9999 < 10000 Dec Prec: 0
CUBIC CENTIMETER -1 > 0 * 9999 < 10000 Dec Prec: 0

=====

TERM: INTAKE - GI TUBE

Page: 274

ID: {38602EA8-600E-4427-A695-029DFE06883F}

=====

Abbreviation..... INT GI TU
Active..... Yes

ID: {80BF8A4E-7536-4325-A803-F66E83C11E3B}

=====
Abbreviation..... INT BLOOD
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Blood Products given to patient
Display Name..... Intake - IV Blood Products
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

PRODUCT: 18 entries

-
1. CRYOPRECIPITATE (CRYOPRECIP)
{5329E737-E059-480D-83F8-1146B091C14E}
 2. FFP (FFP) {9DD0ECF5-634F-49CF-A529-A2657240A80C}
 3. FFP - THAWED (FFP-THAWED) {AAF452A3-0DCA-4FF4-BF96-08E60D2D2875}
 4. GRANULOCYTES PHERESIS (GRAN PHERE)
{2B90878D-9C29-4731-9525-D7DD7FCB513B}
 5. PLATELETS PHERESIS (PLATELET P)
{65CBDCBD-99DB-44FE-9230-D62886208EBF}
 6. PLATELETS PHERESIS IRRADIATED (PLATEL P I)
{70C7155A-5AB6-4733-B43F-606F8BBB3B4A}
 7. PRBC - AUTOLOGOUS (PRBC-A)
{1580990C-0B35-4DDD-84AA-682B341C9B76}
 8. PRBC - LEUKOREduced (PRBC-L)
{13DB3170-6C26-4178-B429-C02AFC541AC7}
 9. PRBC - LEUKOREduced - DIRECTED (PRBC-L-D)
{9C54B7F7-FBD9-4FA8-B2A1-893A25A9CD2C}
 10. PRBC - LEUKOREduced - IRRADIATED (PRBC-L-I)
{B4431519-0102-49ED-94D6-1738E375C5A4}
 11. PRBC - LEUKOREduced - WASHED (PRBC-L-W)
{BB58F35A-1086-4061-B485-FA35DFB79A0D}
 12. WBC - WASHED (WBC WASHED) {F1460F94-203C-473D-AC7C-92E8ADFFC0F5}
 13. PRBC (PRBC) {BCE37EFD-78D6-4A13-A5E4-EA278BF18E02}
 14. PLATELETS (PLATELETS) {3F4FA2EF-DF05-4A6D-839F-F97D3C2AFCF7}
 15. PLASMAPHERESIS (PLASMAPHER)
{E1214367-8BB0-4FAA-A916-EC707569A9C4}
 16. 25 PERCENT ALBUMIN (25 PERC AL)
{CBC7E32F-E457-4A7F-8F27-030EBC6ECBD5}
 17. 5 PERCENT ALBUMIN (5 PERC ALB)
{A2D0F9A0-A9EC-4540-9663-B8B19A189198}
 18. Rh IMMUNE GLOBULIN (Rh IM GLOB)
{1FE7BD0C-7252-427A-8216-36FCA4637089}

UNITS:

MILLILITER -1 > 0 * 9999 < 10000 Dec Prec: 0
CUBIC CENTIMETER -1 > 0 * 9999 < 10000 Dec Prec: 0

ID: {868852E2-3D29-4007-BB36-62F6D90B1685}

=====
Abbreviation..... INT IV CAL
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Calories taken in by IV
Display Name..... Itake - IV Calories
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

CALORIE -1 > 0 * 9999 < 20000 Dec Prec: 0

=====
TERM: INTAKE - IV CARDIAC OUTPUT INJECTATE
ID: {FE023C43-32D0-4937-A100-0F14CDCACEF2}

Page: 278

=====
Abbreviation..... INT CO INJ
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Cardiac Output Injectate
Display Name..... Cardiac Output Injectate
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

MILLILITER -1 > 0 * 9999 < 10000 Dec Prec: 0
CUBIC CENTIMETER -1 > 0 * 9999 < 10000 Dec Prec: 0

=====
TERM: INTAKE - IV COLLOID
ID: {7CD55D76-7D00-414E-BF24-315A5A92A998}

Page: 279

=====
Abbreviation..... INT IV COL
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Intake of IV Colloids
Display Name..... Intake IV Colloids
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

PRODUCT: 4 entries

- 1. ALBUMIN (ALBUMIN) {EF4CC3AE-32B2-4EBC-BBDE-7DE7CCEAF081}
- 2. 5 PERCENT ALBUMIN (5 PERC ALB)
{A2D0F9A0-A9EC-4540-9663-B8B19A189198}
- 3. 25 PERCENT ALBUMIN (25 PERC AL)
{CBC7E32F-E457-4A7F-8F27-030EBC6ECBD5}
- 4. 6 PERCENT HETASTARCH (6 PERC HET)
{8A82F108-D2FE-4C21-93BB-344502235AE8}

UNITS:

MILLILITER -1 > 0 * 9999 < 20000 Dec Prec: 0
 CUBIC CENTIMETER -1 > 0 * 9999 < 10000 Dec Prec: 0

=====

TERM: INTAKE - IV CRYSTALLOIDS Page: 280
 ID: {DD597777-861C-44FD-82E6-1FFF925D0095}

=====

Abbreviation..... INT IV CRY
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... Numeric
 Description..... Intake of IV Crystalloids
 Display Name..... Intake - IV Crystalloids
 Help Text..... No Help Text
 Multi Select Picklist.....
 Term Type..... OBSERVATION
 Vuid.....

UNITS:

MILLILITER -1 > 0 * 9999 < 10000 Dec Prec: 0
 CUBIC CENTIMETER -1 > 0 * 9999 < 10000 Dec Prec: 0

=====

TERM: INTAKE - IV FLUID Page: 281
 ID: {5902ECD7-7C26-481F-AE76-999211321FFD}

=====

Abbreviation..... INT IV FL
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... Numeric
 Description..... Intake - IV Fluids
 Display Name..... Intake - IV Fluids
 Help Text..... Enter Fluid type as 'Comments',
 Continuous or IVPB
 Multi Select Picklist.....
 Term Type..... OBSERVATION
 Vuid.....

=====

TERM: INTAKE - LINE FLUSH

Page: 284

ID: {0DAF17DD-4039-406C-9F53-5AB46E172648}

=====

Abbreviation..... IN LINE FL
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Intake from line flush
Display Name..... Intake - Line Flush
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

PRODUCT: 2 entries

1. SALINE (SALINE) {19D84001-70DB-47AA-9149-A37B7E0BAEF1}
2. LINE FLUSH HEPARIN (LN FL HEPA)
{E8B68E5A-6B9A-401F-9D1A-E3C907C7B81C}

UNITS:

MILLILITER -1 > 0 * 9999 < 10000 Dec Prec: 0
CUBIC CENTIMETER -1 > 0 * 9999 < 10000 Dec Prec: 0

=====

TERM: INTAKE - PCA VOLUME

Page: 285

ID: {86B055A5-04E4-4C85-B20B-9EA8F02E45AB}

=====

Abbreviation..... INT PCA VO
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Intake PCA Volume
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

MILLILITER -1.00 > 0.00 * 5000.00 < 5001.00
Dec Prec: 2

=====

TERM: INTAKE - PERCENT MEAL CONSUMED

Page: 286

ID: {A3C27D46-477C-476E-8D56-A652B8BC3933}

=====

Abbreviation..... INT % MEAL

- {17F8F1DA-4FC6-48E5-9E6E-1F67CD06D200}
- 2. WATER (WATER) {ACCC998E-0F64-510D-B843-3D30B2AD67CF}
- 3. JUICE (JUICE) {91FFF519-1677-BB1D-D959-65B1BFA10285}
- 4. FULL LIQUID (FL LIQUID) {C8405270-CCA3-4BFF-B9D0-398590FBB1D9}
- 5. SOFT DRINK (SOFT DRIN) {9216F0F9-4D5F-4AF7-8D84-04A2CE11BCBE}
- 6. OTHER PO FLUID (PO OTHER) {1EA44424-4DE2-42B8-A50F-60DC56204502}

UNITS:

MILLILITER 0 > 1 * 10000 < 10001 Dec Prec: 0
OUNCE 0 > 1 * 338 < 339 Dec Prec: 0
CUBIC CENTIMETER 0 > 1 * 10000 < 10001 Dec Prec: 0

=====

TERM: INTAKE - TPN CALORIES Page: 289
ID: {6AAC56B6-D4B8-4F1C-A4D0-86FDD3B6C72E}

=====

Abbreviation..... INT TPN CA
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Total parenteral nutrition calories
Display Name..... Intake - TPN Calories
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

CALORIE -1 > 0 * 2000 < 2001 Dec Prec: 0

=====

TERM: INTAKE - TRANSFUSION Page: 290
ID: {DCD2263B-0390-4754-950A-6A410C0C59EB}

=====

Abbreviation..... INT TRANSF
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Transfusion Type
Display Name..... Transfusion
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

CUBIC CENTIMETER -1 > 0 * 9999 < 10000 Dec Prec: 0

=====

TERM: INTRACRANIAL PRESSURE
ID: {E11DD9AE-8220-432E-A855-4697F82AC17A}

=====

Page: 291

Abbreviation..... ICP
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Intracranial Pressure
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

-
1. VENTRICULAR DRAIN (VENT DRAIN)
{5FF035D1-9C97-415D-BBEB-943506395A0C}

UNITS:

MILLIMETERS OF MERCURY 0 > 8 * 20 < 30 Dec Prec: 0

=====

TERM: ISOLATION PRECAUTIONS
ID: {B8AADF9F-BDAC-49F7-BAB8-4B6A1E9865D1}

=====

Page: 292

Abbreviation..... ISOLAT PRE
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Isolation Precautions
Help Text..... Environmental Precautions for
Isolated Patients
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 7 entries

-
1. STRICT (STRICT) {5717A3EA-9A1F-44AB-BB0B-E86BBDED3ECF}
 2. CONTACT (CONTACT) {958F2B5E-F980-47B6-A7F7-EE7D5520B61F}
 3. NEUTROPENIC (NEUTROPENI) {5CB9620C-34D9-48F7-AAF0-C595D06FC8AE}
 4. STANDARD (STANDARD) {C4B8C06D-BFFE-4438-A43D-926F039E5EE1}
 5. AIRBORNE (AIRBORNE) {D2F281F7-DC09-4CA2-B5B3-F92ACBC63A0A}
 6. DROPLET (DROPLET) {CE1DD314-8A30-4785-BD56-1540CDA85146}
 7. NEGATIVE PRESSURE (NEG PRESSU)
{6D6CB4BB-96B0-4978-97CC-116A03479DAD}

UNITS:

=====

TERM: IV DOSE RATE - UNITS/HOUR
ID: {48FBFFFC-1B66-4439-8B74-ECCDDDE37804}

=====

Page: 293

Abbreviation..... DR UN HR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Only to be used with products with
dosage measured in "Units"
Display Name..... IV DOSE RATE - UNITS/HOUR
Help Text..... Only to be used with products with
dosage measured in "Units"
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

UNITS PER HOUR 0 > 1 * 1200 < 2001 Dec Prec: 0

=====

TERM: IV DOSE RATE - UNITS/kg/HOUR
ID: {D1749D90-3EEE-41BE-BBF5-B3B00C18E8D2}

=====

Page: 294

Abbreviation..... IVDR Ukh
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Only to be used with products with
dosage measured in "Units"
Display Name..... IV DOSE RATE - UNITS/kg/HOUR
Help Text..... Only to be used with products with
dosage measured in "Units"
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

UNITS/kg/HOUR 0 > 10 * 51 < 101 Dec Prec: 0

=====

TERM: IV DOSE RATE - g/TIME
ID: {5D83FE3F-45BC-4D82-A7A7-832DA99A751F}

=====

Page: 295

Abbreviation..... IVDR g/t
Active..... Yes
Boolean Value False.....
Boolean Value True.....

MICROGRAMS PER KILOGRAM PER HOUR
0.000 > 1.000 * 2000.000 < 9999.000
Dec Prec: 3

=====

TERM: IV DOSE RATE - mg/TIME Page: 298
ID: {E4E00CBB-AE1C-48BE-A172-CDEFF039A311}

=====

Abbreviation..... nt
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... IV DOSE RATE - mg/TIME
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

MILLIGRAMS PER HOUR	-1.00 > 0.00 * 9999.00 < 10000.00 Dec Prec: 2
MILLIGRAMS PER MINUTE	-1.00 > 0.00 * 9999.00 < 10000.00 Dec Prec: 2

=====

TERM: IV DOSE RATE - mg/kg/HOUR Page: 299
ID: {FE8EB504-6F94-4863-B673-489582F1233C}

=====

Abbreviation..... IV DR MG
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... IV DOSE RATE - mg/kg/HOUR
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

MILLIGRAMS PER KILOGRAM PER HOUR	0 > 1 * 2000 < 9999 Dec Prec: 0
----------------------------------	------------------------------------

=====

TERM: IV FLOW RATE - ml/TIME Page: 300
ID: {0ED682FE-21ED-4D08-9250-46CB6FB85146}

=====

Abbreviation..... IV FLOW R

Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 15 entries

- 1. PATENT CLEAN (PAT CLEAN) {68688246-AE58-2EFB-BF9D-5E00D8CFDA8E}
- 2. DRESSING CHANGED (DRESS CHAN)
{E8CF4B7A-D1AB-49A2-8E7C-3375C3BC436F}
- 3. STARTED (STARTED) {D6E105C4-5EAD-424E-8E4C-195EB7783A3D}
- 4. TUBING CHANGED/LABELED (TUBE CHANG)
{19BA337C-6E28-411E-9702-48B262C593FA}
- 5. HEPLOCK (HEPLOCK) {C827FEDE-59D6-4C45-BDDA-4884FA0882F3}
- 6. BAG CHANGED (BAG CHANGE) {818990DE-4F95-4F8E-92D8-B4D3B975CB36}
- 7. LEVELED (LEVELED) {D72523B2-E0CD-40E1-AC14-2FA43636A456}
- 8. PRESSURE CALIBRATION (PRESS CALI)
{38C12E69-2CD2-4F61-98D9-06AB857DA51A}
- 9. ARTERIAL LINE CARE (ART LINE C)
{787D232D-DFB8-44E2-8A0C-9D7D6B80764E}
- 10. CENTRAL VENOUS CATHETER (CENT VEN C)
{EE91E841-C09B-4432-A121-8EBAB92E9CD2}
- 11. PICC (PICC) {3CFE33D0-5A35-88CC-152C-D7501D1858D9}
- 12. FLUSHED (nt) {C7C54C1A-6E75-4BCF-4C57-0F8622A5BD16}
- 13. CHANGED (CHANGED) {74B14023-D74D-41F6-AA48-B0E5137623CD}
- 14. DISCONTINUED (DISCONTINU) {15E289FB-09DF-4731-9CF3-48A4E0031C31}
- 15. SALINE LOCK (SALIN LOCK) {A94D7639-2A45-4319-A6C7-5A892AA6B1F4}

UNITS:

=====

TERM: IV PORT / NUMBER Page: 303
ID: {84D9EDBF-F38B-496F-90F7-CF99E8D11321}

=====

Abbreviation..... IV PORT #
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... IV Port descriptor or number
Display Name..... IV Port / Number
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 15 entries

- 1. PROXIMAL PORT (PROX PORT) {3244EF4F-F18A-4250-96F5-ABCEA7BFA578}
- 2. MEDIAL PORT (MEDI PORT) {D1734D73-2E1A-4538-8884-05D8A78CF556}
- 3. DISTAL PORT (DIST PORT) {E493DF34-AD0E-4951-9715-5AD4EFA14036}
- 4. 1 (1) {8F5D10A8-FE17-47BC-DACA-6972D12CC00A}
- 5. 2 (2) {38E22B79-A2B6-3401-B195-C955F24E1278}
- 6. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}
- 7. 4 (4) {D2280C7E-DE50-5A8D-C3FD-5862B69C0E1A}

- 8. 5 (5) {188EDCE7-0726-4234-D8D2-92C8052CD221}
- 9. 6 (6) {640CF2F9-3BBF-AB13-185B-CD02A609FB5A}
- 10. 7 (7) {7A010324-D2C2-5D4B-848A-7B8878666D20}
- 11. 8 (8) {518E7CFC-102E-351E-C7AF-FB575E3C8F4D}
- 12. 9 (9) {58FF2FD4-569F-F454-9F5B-FE6F9E34D8CC}
- 13. 10 (10) {0C827160-18DB-816F-7FAA-A3C175C6D33B}
- 14. 11 (11) {BD51A7B3-FACF-890C-472E-E3B28519CC16}
- 15. 12 (12) {ACC5CA06-3ABC-476D-EFB0-D82FEEE129C7}

LOCATION: 20 entries

- 1. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
- 2. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
- 3. L ANTECUBITAL (L ANTECUBI) {C21838CC-BAB5-1D7F-E2A4-A7C663FB28C3}
- 4. R ANTECUBITAL (R ANTECUBI) {7FB00A65-B4B2-06A6-BD89-98AC112D8D4E}
- 5. L UPPER ARM (L UPP ARM) {5F915696-1B08-32E9-6B95-49246117D660}
- 6. R UPPER ARM (R UPP ARM) {3EF20E8E-F68C-E4F2-8F8C-FAD3B44DDAA4}
- 7. L WRIST (L WRIST) {7BCAFACE-FF3B-46DA-99F9-6437A90B7126}
- 8. R WRIST (R WRIST) {FDD9FE1D-F8A2-09AD-82D7-AA478693DD01}
- 9. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
- 10. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
- 11. L LEG (LL) {B088C809-E528-B785-21C2-14E45E1FD8AF}
- 12. R LEG (RL) {C5E40A7B-6FC4-267D-8CE8-3571D21D81FD}
- 13. L FOOT (L FOOT) {2EC3FBEC-4EA2-9E7A-17D4-FDC2F9D28046}
- 14. R FOOT (R FOOT) {7D9F5FF6-232C-2488-4504-A721FF06928F}
- 15. L INTERNAL JUGULAR (LIJ) {AFBB4FE5-8A11-2047-7900-DADF2FC4430D}
- 16. R INTERNAL JUGULAR (R INT JUGU) {7ABF6916-D379-198D-6454-F0518A74844A}
- 17. L EXTERNAL JUGULAR (L EXT JUGU) {6F7B331B-4F83-0B79-1755-F6F673D647D6}
- 18. R EXTERNAL JUGULAR (R EXT JUGU) {4D78F704-9E3B-907E-A792-4DC4F4CFF170}
- 19. L SUBCLAVIAN (LSCL) {71B93594-15F9-AA4A-AEF2-E222EC28C00D}
- 20. R SUBCLAVIAN (R SUBCLAVI) {53DEAE83-93AE-123B-14E5-4DE38CD4F752}

UNITS:

=====

TERM: IV SITE CONDITION

Page: 304

ID: {84C7D095-10E1-AFD2-405D-C44C8278241F}

=====

Abbreviation..... IV ST COND
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... No Description
 Display Name..... IV Site Condition
 Help Text..... Condition of IV Site
 Multi Select Picklist..... Yes
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 4 entries

- 1. PATENT CLEAN (PAT CLEAN) {68688246-AE58-2EFB-BF9D-5E00D8CFDA8E}
- 2. RED (RED) {515D7CFB-3527-30E7-EC4B-B607C272DFA0}
- 3. SWOLLEN (SWOLLEN) {DCE86B10-C393-8DA9-6518-F2943F0E1A4E}
- 4. DISCOMFORT (DISCOMFORT) {77D058C9-1EB4-96C0-C9E3-666C840B884C}

LOCATION: 20 entries

- 1. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
- 2. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
- 3. L ANTECUBITAL (L ANTECUBI) {C21838CC-BAB5-1D7F-E2A4-A7C663FB28C3}
- 4. R ANTECUBITAL (R ANTECUBI) {7FB00A65-B4B2-06A6-BD89-98AC112D8D4E}
- 5. L UPPER ARM (L UPP ARM) {5F915696-1B08-32E9-6B95-49246117D660}
- 6. R UPPER ARM (R UPP ARM) {3EF20E8E-F68C-E4F2-8F8C-FAD3B44DDAA4}
- 7. L WRIST (L WRIST) {7BCAFACE-FF3B-46DA-99F9-6437A90B7126}
- 8. R WRIST (R WRIST) {FDD9FE1D-F8A2-09AD-82D7-AA478693DD01}
- 9. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
- 10. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
- 11. L LEG (LL) {B088C809-E528-B785-21C2-14E45E1FD8AF}
- 12. R LEG (RL) {C5E40A7B-6FC4-267D-8CE8-3571D21D81FD}
- 13. L FOOT (L FOOT) {2EC3FBEC-4EA2-9E7A-17D4-FDC2F9D28046}
- 14. R FOOT (R FOOT) {7D9F5FF6-232C-2488-4504-A721FF06928F}
- 15. L INTERNAL JUGULAR (LIJ) {AFBB4FE5-8A11-2047-7900-DADF2FC4430D}
- 16. R INTERNAL JUGULAR (R INT JUGU) {7ABF6916-D379-198D-6454-F0518A74844A}
- 17. L EXTERNAL JUGULAR (L EXT JUGU) {6F7B331B-4F83-0B79-1755-F6F673D647D6}
- 18. R EXTERNAL JUGULAR (R EXT JUGU) {4D78F704-9E3B-907E-A792-4DC4F4CFF170}
- 19. L SUBCLAVIAN (LSCL) {71B93594-15F9-AA4A-AEF2-E222EC28C00D}
- 20. R SUBCLAVIAN (R SUBCLAVI) {53DEAE83-93AE-123B-14E5-4DE38CD4F752}

UNITS:

=====

TERM: LEFT VENTRICULAR EJECTION FRACTION - CONTINUOUS
ID: {1D980210-D857-457D-A335-DC248B264B10}

=====

Page: 305

Abbreviation..... LVEF - CON
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Continuous Left Ventricular Ejection
Fraction
Display Name..... LVEF - Continuous
Help Text..... If a range is provided use the midpoint
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 7 entries

- 1. ECHOCARDIOGRAPHY (ECG) {CA44697E-3BBA-4767-A30D-6548A72F9639}
- 2. NUCLEAR PERFUSION IMAGING (NPI)
{D0F7DB21-0D66-4161-8214-DACCAF41ADAC}
- 3. MULTIGATED NUCLEAR VENTRICULOGRAPHY (MUGA)
{4D73168B-4DF4-4EA2-83B9-ACE853D5DD42}
- 4. CONTRAST LEFT VENTRICULOGRAPHY (CON L VENT)
{35BEA63F-9696-4C03-AE1D-D76E4DDBB017}
- 5. MAGNETIC RESONANCE IMAGING (MRI)
{50EBFFF8-DF0D-41F2-B11B-33A26AA40DBA}
- 6. COMPUTED TOMOGRAPHY (COM TOMOGR)
{B4676AE0-BFAF-4649-AB08-3710E39F9B8D}
- 7. OTHER (OTHER) {AED275BD-8A42-46BE-BCDF-472EF4C31D65}

UNITS:

PERCENTAGE 0 > 1 * 99 < 100 Dec Prec: 0

TERM: LEFT VENTRICULAR EJECTION FRACTION - CATEGORICAL Page: 306
ID: {8B4D6E3B-1138-4BC0-A303-44B09F05A058}

Abbreviation..... LVEF - CAT
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Categorical Left Ventricular Ejection
Fraction
Display Name..... LVEF - Categorical
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 8 entries

- 1. HYPERDYNAMIC (HYPERDYNAM) {818A8A99-04D9-4989-A9BA-ECFACED74605}
- 2. NORMAL (n) {553C9A97-A3D5-C52F-0E4E-C7D168C99A2B}
- 3. BORDERLINE REDUCED (BORD REDUC)
{63703AB1-6249-403B-BDE8-FECA9FA36CF5}
- 4. MILD-MODERATELY REDUCED (MI-MOD-RED)
{41DA2518-0E49-4E92-A733-3ACF3D2E8199}
- 5. MODERATELY REDUCED (MOD REDUCE)
{91F73300-7EE0-428C-A838-5B3BE062D8BC}
- 6. MODERATE-SEVERELY REDUCED (MO-SEV RED)
{6D031C70-2F95-4111-A9B8-BC22CEC1B17D}
- 7. SEVERELY REDUCED (SEV REDUCE)
{A895BACE-85AC-45B2-A190-4A6816DEDC1D}
- 8. LESS THAN 30 (<30) {5AA98A0D-9B05-425D-97CA-B8B12425AACA}

METHOD: 7 entries

1. ECHOCARDIOGRAPHY (ECG) {CA44697E-3BBA-4767-A30D-6548A72F9639}
2. NUCLEAR PERFUSION IMAGING (NPI)
{D0F7DB21-0D66-4161-8214-DACCAF41ADAC}
3. MULTIGATED NUCLEAR VENTRICULOGRAPHY (MUGA)
{4D73168B-4DF4-4EA2-83B9-ACE853D5DD42}
4. CONTRAST LEFT VENTRICULOGRAPHY (CON L VENT)
{35BEA63F-9696-4C03-AE1D-D76E4DDBB017}
5. MAGNETIC RESONANCE IMAGING (MRI)
{50EBFFF8-DF0D-41F2-B11B-33A26AA40DBA}
6. COMPUTED TOMOGRAPHY (COM TOMOGR)
{B4676AE0-BFAF-4649-AB08-3710E39F9B8D}
7. OTHER (OTHER) {AED275BD-8A42-46BE-BCDF-472EF4C31D65}

UNITS:

=====

TERM: LEFT VENTRICULAR STROKE WORK
ID: {E9D668F0-920C-7E29-A6A1-618FC8ECD19B}

=====

Page: 307

Abbreviation..... LVSW
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Left Ventricular Stroke Work
Display Name..... Left Ventricular Stroke Work
Help Text..... Left Ventricular Stroke Work
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

1. CARDIAC OUTPUT MONITOR (COM)
{CCBE5744-7CB9-9DFA-938F-FDDBC73B45B3}

UNITS:

STROKE WORK PER BEAT 29.9 > 30.0 * 150.0 < 150.1
Dec Prec: 1

=====

TERM: LEG MOVEMENTS
ID: {49E9E897-E19A-8FE9-416C-F98CA66EFBC9}

=====

Page: 308

Abbreviation..... LEGMOV
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Leg Movements
Help Text..... No Help Text

Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

- 1. ABSENT (a) {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0}
- 2. PRESENT (p) {384DC20F-1625-CE1E-1C85-A3B3924BB2A0}

UNITS:

TERM: LEVEL OF CONSCIOUSNESS Page: 309
ID: {48EF6B0D-C435-12D4-EA93-F8F5D344E4C5}

Abbreviation..... LOC
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Level of Consciousness
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 5 entries

- 1. ALERT (a) {BC7E2630-179A-5C80-E6C7-598D42F0BFD6}
- 2. DROWSY (d) {6B350814-A88A-0F05-C382-9F27830FBC09}
- 3. RESTLESS (r) {8FCCAD26-DBA0-7067-DA72-30D342BFA8E8}
- 4. AGITATED (a) {BAD2CCCE-9083-6F61-1196-C0C9A8DDBFFC}
- 5. LETHARGIC (l) {4ECA2E44-093A-9354-848E-CDC72B165072}

UNITS:

TERM: LOWER LIMB EXTREMITIES STRENGTH Page: 310
ID: {9215BD72-FA1A-9146-EA8D-B4A33B455724}

Abbreviation..... LLES
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Lower Limb Extremities Strength
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Void.....

DATA VALUE: 6 entries

- 1. ABSENT (a) {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0}
- 2. WEAK (w) {BE2D9305-E2DA-F492-66CF-164EE58B8179}
- 3. STRONG (s) {A12154E1-544E-F372-3839-0D44E5BC476B}
- 4. REFLEX (r) {F21585F9-5206-A6B3-2BF6-F416E5842366}
- 5. EQUAL (Equal) {AC079584-158D-EBF2-E13A-C2E919375E4D}
- 6. UNEQUAL (Unequal) {73319200-164F-B39F-649C-B7FB1AFBDBB8}

LOCATION: 3 entries

- 1. L LOWER LIMB (LLLI) {4838CB05-A1DD-F4FB-587C-2D0C17A2484E}
- 2. R LOWER LIMB (R LOW LIMB) {A8E26A05-311A-F196-40B5-B5202D9F8058}
- 3. BILATERAL LEG (BILAT LEG) {0713761D-E0FC-DA03-02A1-F20EFEB99F45}

UNITS:

TERM: LVSWI Page: 311
ID: {7945FB00-2FA8-6D4C-BC85-0110EA75AC5C}

Abbreviation..... LVSWI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Left Ventricular Stroke Work Index
Display Name..... LVSWI
Help Text..... Left Ventricular Stroke Work Index
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

- 1. CARDIAC OUTPUT MONITOR (COM)
{CCBE5744-7CB9-9DFA-938F-FDDBC73B45B3}

UNITS:

STROKE WORK PER BEAT PER M2 19 > 20 * 150 < 151 Dec Prec: 0
GM-M/M2/BEAT 19 > 20 * 150 < 151 Dec Prec: 0

TERM: MEAN ARTERIAL PRESSURE Page: 312
ID: {9D9D9FD8-C352-6434-E796-C583407EBD45}

Abbreviation..... MAP
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric

Description..... No Description
Display Name..... Mean Arterial Pressure
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 5 entries

- 1. CUFF (Cu) {E59A2543-1860-1E1C-BEF1-6E7F6B6FBC5E}
- 2. ARTERIAL LINE (AL) {D4B3F4ED-8B33-431E-4164-BD165E823776}
- 3. INVASIVE (I) {5DB8746E-5446-F47A-883C-A44B91275A6E}
- 4. CARDIAC OUTPUT MONITOR (COM)
{CCBE5744-7CB9-9DFA-938F-FDDBC73B45B3}
- 5. NON-INVASIVE (N-I) {F50ABFED-5FD7-FC95-20FE-1E3038C9F299}

UNITS:

MILLIMETERS OF MERCURY 25 > 50 * 100 < 300 Dec Prec: 0

TERM: MEAN BLOOD PRESSURE Page: 313
ID: {9038EF74-932A-42EE-992F-A20BD8784BB7}

Abbreviation..... BP MEAN
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Mean Blood Pressure
Display Name..... Mean Blood Pressure
Help Text..... Mean Blood Pressure
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 5 entries

- 1. ARTERIAL LINE (AL) {D4B3F4ED-8B33-431E-4164-BD165E823776}
- 2. NON-INVASIVE (N-I) {F50ABFED-5FD7-FC95-20FE-1E3038C9F299}
- 3. CUFF (Cu) {E59A2543-1860-1E1C-BEF1-6E7F6B6FBC5E}
- 4. IABP (IABP) {4352A41F-AAF3-4544-9590-BAFD7977DF49}
- 5. INVASIVE (I) {5DB8746E-5446-F47A-883C-A44B91275A6E}

UNITS:

MILLIMETERS OF MERCURY 20 > 50 * 100 < 200 Dec Prec: 0

TERM: MEAN BLOOD PRESSURE ARTERIAL LINE Page: 314
ID: {4493FB14-8C5E-9995-C08F-D0BC9A3B5D0B}

Abbreviation..... MBPAL
Active..... Yes

Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Mean Blood Pressure Arterial Line
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

MILLIMETERS OF MERCURY 24 > 25 * 200 < 201 Dec Prec: 0

=====

TERM: MEAN BLOOD PRESSURE NON-INVASIVE	Page: 315
ID: {A89EC390-60BB-FE79-6F40-F9E2BAFE5BB5}	

=====

Abbreviation..... MBPNI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Mean Blood Pressure Non-Invasive
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

MILLIMETERS OF MERCURY 24 > 25 * 200 < 201 Dec Prec: 0

=====

TERM: MEAN BP CUFF	Page: 316
ID: {1F7DAE88-DE53-4807-F51F-47AB8A672EA1}	

=====

Abbreviation..... MBPC
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Mean BP Cuff
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

MILLIMETERS OF MERCURY 24 > 25 * 200 < 201 Dec Prec: 0

TERM: MORSE: IV THERAPY SALINE LOCK
ID: {5FB26C06-AF8D-ECBE-FC10-442510EC52F1}

Page: 321

Abbreviation..... MORSE IV
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Morse: IV Therapy Saline Lock
Help Text..... 0 if NO, 25 if YES
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

1. 0 (0)	{2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
2. 20 (20)	{532FBDF8-7F7B-9551-300C-F2CFEF656163}

UNITS:

TERM: MORSE: MENTAL STATUS OF PT OWN ABILITY
ID: {E15F905B-5D1C-47CA-A23A-5EC47302D9A3}

Page: 322

Abbreviation..... MORSE PT S
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Morse: Mental Status of Patients Own
Ability
Display Name..... Morse: Mental Status of Patients Own
Ability
Help Text..... 0: Oriented to own ability, 15
Overestimates, forgets limitations
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

1. 0 (0)	{2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
2. 15 (15)	{FE8C049F-B64C-4991-2F50-3C561BA4825C}

UNITS:

TERM: MORSE: SCORE TOTAL
ID: {80F2ACB5-5E42-48B8-8EB8-CADAD149E7BB}

Page: 323

Vuid.....

DATA VALUE: 6 entries

- 1. NORMAL POWER (np) {0E45C061-0744-8A21-7256-003702A5005B}
- 2. NO MOVEMENT (nm) {2618E301-379B-BAB9-9EDE-98FD4B60D311}
- 3. FLICKER OF CONTRACTION (foc)
{218B322B-84AA-F6F6-5207-D86F9A4311F0}
- 4. MOVEMENT WITH GRAVITY ELIMINATED (mwge)
{B86E1131-1976-FA8E-74FE-B1BC2F0FDE3F}
- 5. MOVE AGAINST GRAVITY BUT NO RESISTANCE (magbnr)
{F995A165-8536-C493-2269-286E8EFF5E43}
- 6. MOVE AGAINST GRAVITY AND RESISTANCE, EASILY OVERCOME BY THE EXAMINER (magareobte)
{2ECCC197-4A81-025F-3EF3-6214004887A9}

LOCATION: 4 entries

- 1. L LOWER EXTREMITY (LLEX) {FC6F5056-779A-BF76-86D0-A237936C1ECC}
- 2. L UPPER EXTREMITY (LUEX) {48806E79-B1DA-129B-EB52-453292A3DDC6}
- 3. R LOWER EXTREMITY (R LOW EXTR)
{7FC3CECC-1E6E-A686-FF04-2BEF8A9B9DCA}
- 4. R UPPER EXTREMITY (R UPP EXTR)
{C65BA269-8F05-6FB9-62CB-466A286CF2C0}

UNITS:

=====

TERM: MOTOR CHECK Page: 326
 ID: {521171CF-BD9D-2F28-4042-EA885C40E556}

=====

Abbreviation..... MOCH
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... No Description
 Display Name..... Motor Check
 Help Text..... No Help Text
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 2 entries

- 1. INTENTIONAL (i) {D5BBBFFC-F739-1E69-FED4-288D8E803D04}
- 2. REFLEX (r) {F21585F9-5206-A6B3-2BF6-F416E5842366}

UNITS:

=====

TERM: MUCOUS MEMBRANE CONDITION Page: 327
 ID: {4D9B0E42-4A15-4699-A80E-31588D8C13FD}

=====

Abbreviation..... MUC MEMBRA
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Mucous Membrane Condition
Display Name..... Mucous Membrane Condition
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

-
1. PINK (PINK) {1A7FAF84-D69C-4EEB-7BC5-D21A598A1F4D}
 2. CYANOTIC (c) {06BD77DB-1849-D731-0F27-F8CE48755EA9}
 3. PALE (p) {0EB39D28-7ABE-A4AB-D334-3F795E65618D}

UNITS:

TERM: NAIL COLOR Page: 328
ID: {63FF0B1B-B57D-43A9-A05F-3F2CF6FE8926}

Abbreviation..... NAIL COLOR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Nail Color
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

-
1. NORMAL FOR ETHNIC GROUP (ETH. NORM.)
{3B2CD951-58CE-FF53-AD5D-6C65C07658D3}
 2. CYANOTIC (c) {06BD77DB-1849-D731-0F27-F8CE48755EA9}
 3. DISCOLORATION (DISCOLORAT)
{BBE3E910-CC49-430C-BD90-FB0A97F1B183}
 4. PALE (p) {0EB39D28-7ABE-A4AB-D334-3F795E65618D}
 5. HEMORRHAGE (HEMORRHAGE) {E00AC828-83F3-4F69-B97B-D6AD1D87A440}

UNITS:

TERM: NAIL SHAPE Page: 329
ID: {6351A78B-67BC-46A4-B7E4-72D63C8C2917}

Abbreviation..... NAIL SHAPE
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Nail Shape
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

-
1. NORMAL (n) {553C9A97-A3D5-C52F-0E4E-C7D168C99A2B}
 2. CLUBBING (CLUBBING) {AE77D344-6F0F-452B-B1AA-EC8D9F092CFC}
 3. SPOONING (SPOONING) {CA384C9F-D87D-4E2D-B4C0-89BD85692B16}

UNITS:

TERM: NAIL TEXTURE Page: 330
ID: {AD3F008F-6AA8-483F-AAA6-5F797109F4A6}

Abbreviation..... NAIL TEXTU
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Nail Texture
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

-
1. NORMAL (n) {553C9A97-A3D5-C52F-0E4E-C7D168C99A2B}
 2. THICK (THICK) {1554A49E-1468-46B2-F851-531DDB1D1A98}
 3. PITTED (PITTED) {BEA42FD8-91C7-4BF6-936E-726D78064F11}

UNITS:

TERM: NECK VEIN CONDITION Page: 331
ID: {D76459ED-8572-409E-85A4-F64228AE9A36}

Abbreviation..... NECK VEIN
Active..... Yes
Boolean Value False.....

Boolean Value True.....
Data Type..... PickList
Description..... Neck Vein Condition
Display Name..... Neck Vein Condition
Help Text..... Condition of Neck Veins
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

-
- 1. NORMAL (n) {553C9A97-A3D5-C52F-0E4E-C7D168C99A2B}
 - 2. FLAT (Flat) {EA4F1980-EDF3-82BC-B677-D747F48F1163}
 - 3. DISTENDED (DISTENDED) {B5C529CA-2BE0-4698-A153-718F4D58E643}

UNITS:

=====

TERM: NEED FOR FLUIDS/TOILETING ASSESSED
ID: {816D5DB3-7AEA-E73F-08D8-68C58D4FB4FE}

=====

Page: 332

Abbreviation..... NFFTA
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... No Description
Display Name..... Need for Fluids/Toileting Assessed
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

-
- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
 - 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: NOSE SIZE
ID: {0A7CD784-279C-9F56-67CC-1C0624B26F4A}

=====

Page: 333

Abbreviation..... ns
Active..... No
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Describe the size of the honker
Display Name..... Nose Size
Help Text..... Honker etiquette

human function
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 6 entries

-
1. ORTHOTIC BOOTS (ORTHO BOOT)
{D4219FE9-7087-B881-6DAC-4E0BB94CDA8E}
 2. ORTHOTIC COLLAR (ORTHO COLL)
{A74E1C89-2304-0544-9A8B-A23E4E2C6B08}
 3. ORTHOTIC SPINE BRACE (ORTH SPINE)
{B47A3536-FEBC-802E-01A1-E909147D18C2}
 4. ORTHOTIC KNEE IMMOBILIZER (KNEE IMMOB)
{FF8A8BB4-63FB-3B0D-83E1-5432A0434DAC}
 5. ORTHOTIC WEDGE PILLOW (WEDG PILLO)
{20E58D21-399B-85E4-BDAB-93B371395BCF}
 6. ORTHOTIC TED HOSE (TED HOSE)
{FE277269-7BD1-85FB-DC3C-FEDE31C66BB0}

UNITS:

=====

TERM: OSTOMY TYPE Page: 336
ID: {CA5CC473-EE95-4E81-8E24-DC23B47A7F42}

=====

Abbreviation..... OSTOMY TYP
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Type of Ostomy
Display Name..... OSTOMY TYPE
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

-
1. COLOSTOMY (COLOSTOMY) {53F939AB-6AF4-4DD5-8D64-17AC1D1736B7}
 2. JEJUNOSTOMY (JEJUNOSTOM) {1ED3EC41-3ACF-42C1-B64A-280300EDBCFC}
 3. SIGMOIDOSTOMY (SIGMOIDOST)
{0C84B317-59B3-433B-A9AD-6588C1361E0C}
 4. ILEOSTOMY (ILEOSTOMY) {D5DF49C8-C6BA-44C1-9840-7765907730A3}
 5. RECTAL TUBE (RECT TUBE) {7E492AE5-D3B8-44EC-8841-97770FAD3AFE}

UNITS:

=====

TERM: OUTPUT Page: 337
ID: {733293B4-0FEC-9D4B-191C-85D490209348}

Abbreviation..... OUTPUT
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Output
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

PRODUCT: 10 entries

1. CARDIAC OUTPUT INJECTATE (CARD OUT I)
{2128B705-DDF4-4776-87C3-AA82F854E313}
2. EF RESIDUAL (EF RESID) {35E269DF-F246-4DF6-9D68-7FD61F8FE088}
3. NG RESIDUAL (NG RESID) {B03C9557-84B5-4B05-B784-C87FA90D27A7}
4. HEME (HEME) {A0ADA229-E779-4562-84B5-67121B88C443}
5. HEMODIALYSIS (HEMODIALYS) {8409B357-C41E-4F4A-BB52-F09AD9704C30}
6. MEDIASTINAL (MEDIASTINA) {A0550B31-B937-4FC7-8F7E-071F75215311}
7. BLADDER IRRIGATION (BLAD IRRIG)
{BCBB106B-ECE4-4232-87C7-D90EDB47F80D}
8. WOUND VACUUM (WOUND VAC) {9F4983A6-B622-4E4D-A7F7-ADA7A4B7CEBC}
9. GASTRIC DRAINAGE (GAST DRAIN)
{C01D18F9-191A-4D34-B45D-D8D70ACA67BE}
10. WOUND DRAINAGE (WOUN DRAIN)
{73038CBA-1F31-47E2-A78E-50ED1067CBEA}

UNITS:

MILLILITER -1 > 0 * 9999 < 10000 Dec Prec: 0
OUNCE -1 > 0 * 338 < 339 Dec Prec: 0

TERM: OUTPUT - BLADDER SCAN Page: 338
ID: {403C285B-40ED-4ECA-BFD0-4404BE4FB59E}

Abbreviation..... BLAD SCAN
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Bladder Scan
Display Name..... Bladder Scan
Help Text..... Bladder Scan Urine Volume Results
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

MILLILITER -1 > 0 * 100 < 1000 Dec Prec: 0

Description..... AMOUNT OF CHEST TUBE DRAINAGE
Display Name..... Chest Tube Drainage Amount Estimated
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

- 1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
- 2. SCANT (SCANT) {8EF6C62D-0D94-37BB-60AD-86E306E8120D}
- 3. SMALL (s) {5991F892-B28C-3998-E733-D8FFA26D7E21}
- 4. MODERATE (MODERATE) {C3C25064-C9FC-EA1F-3E74-32A25F0FC5E3}
- 5. LARGE (l) {E60D5CF9-920D-44E2-8A83-136DAB1FE251}

UNITS:

=====

TERM: OUTPUT - DIALYSIS Page: 342
ID: {C7954A6C-4214-4CB2-B80E-8F4128D05847}

=====

Abbreviation..... OUT DIALYS
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Dialysis Output
Display Name..... Dialysis Output
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 4 entries

- 1. HEMODIALYSIS (HEMODIALYS) {DE7CD394-9E69-46B0-A207-5A18F7C5A1A7}
- 2. PERITONEAL DIALYSIS (PERI DIALY) {6B4AA38B-AFCC-42CE-9F32-59B32E8A6BBB}
- 3. CVVHD (CVVHD) {4CC9AB0D-0B6C-43F9-A93E-A8228EC3A77F}
- 4. CRRT (CRRT) {AB37FC1D-DD7A-4103-A777-7B8D0F876889}

UNITS:

MILLILITER -1 > 0 * 9999 < 10000 Dec Prec: 0
CUBIC CENTIMETER -1 > 0 * 9999 < 10000 Dec Prec: 0

=====

TERM: OUTPUT - DRAIN Page: 343
ID: {B8FD2F72-FCDD-484A-98E1-14B3637B797A}

=====

Abbreviation..... OUT DRAIN
Active..... Yes
Boolean Value False.....

Boolean Value True.....
Data Type..... Numeric
Description..... Drain Output
Display Name..... Drain Output
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 5 entries

- 1. JP DRAIN (JP DRAIN) {94A32247-6190-4047-B581-A212FE4C9453}
- 2. HEMOVAC (HEMOVAC) {79DA0646-B811-4954-8BE1-C133398CD574}
- 3. PENROSE (PENROSE) {CBC3C6DE-E5DE-41D5-83AB-30A4F3CFA9E1}
- 4. CHEST TUBE (CH TUBE M) {F4F467D3-BB57-4E6F-9EA0-C0FBFB66F5B1}
- 5. OTHER (OTHER) {AED275BD-8A42-46BE-BCDF-472EF4C31D65}

UNITS:

MILLILITER -1 > 0 * 9999 < 10000 Dec Prec: 0
CUBIC CENTIMETER -1 > 0 * 9999 < 10000 Dec Prec: 0

TERM: OUTPUT - EMESIS Page: 344
ID: {56BD11AB-FE61-4785-B3CA-C528A5F4EBD6}

Abbreviation..... OUT EMESIS
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Emesis Output
Display Name..... Emesis Output
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

PRODUCT: 1 entries

- 1. EMESIS (EMESIS) {A7B1605E-1DFE-4FC6-A69B-0E18C1A291FF}

UNITS:

MILLILITER -1 > 0 * 9999 < 10000 Dec Prec: 0
CUBIC CENTIMETER -1 > 0 * 9999 < 10000 Dec Prec: 0

TERM: OUTPUT - GI TUBE Page: 345
ID: {CB54872E-CF17-4BA5-BF3B-80B0E6C745F8}

Abbreviation..... OUT GI TU
Active..... Yes
Boolean Value False.....

- 2. SMALL (s) {5991F892-B28C-3998-E733-D8FFA26D7E21}
- 3. MODERATE (MODERATE) {C3C25064-C9FC-EA1F-3E74-32A25F0FC5E3}
- 4. LARGE (l) {E60D5CF9-920D-44E2-8A83-136DAB1FE251}
- 5. COPIUS (COPIUS) {C74A6979-0E35-AA12-F992-69ACA3B0FA23}

QUALITY: 1 entries

- 1. ESTIMATED (E) {9EAFF26B-7D87-DA66-F991-022D095FB5A0}

UNITS:

TERM: OUTPUT - URINE

Page: 348

ID: {F4998CFC-19CA-402A-AB0F-C32FA194DAA5}

Abbreviation..... OUT URINE
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... Numeric
 Description..... Urine Output
 Display Name..... Urine Output
 Help Text..... No Help Text
 Multi Select Picklist.....
 Term Type..... OBSERVATION
 Vuid.....

UNITS:

MILLILITER -1 > 0 * 9999 < 10000 Dec Prec: 0
 CUBIC CENTIMETER -1 > 0 * 9999 < 10000 Dec Prec: 0

TERM: OUTPUT - URINE ESTIMATED

Page: 349

ID: {C5116BFA-332E-4C05-8FF7-4C828D8999BF}

Abbreviation..... OUT UR EST
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... Output - Urine Estimated
 Display Name..... Output - Urine Estimated
 Help Text..... Output - Urine Estimated
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 4 entries

- 1. SLIGHT (SLIGHT) {805549C1-5DC5-4A2D-BC1F-046E44928843}
- 2. SMALL (s) {5991F892-B28C-3998-E733-D8FFA26D7E21}
- 3. MODERATE (MODERATE) {C3C25064-C9FC-EA1F-3E74-32A25F0FC5E3}

Boolean Value False..... Off
Boolean Value True..... On
Data Type..... Boolean
Description..... PA Catheter Alarms On
Display Name..... PA Catheter Alarms On
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

-
- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
 - 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: PA CATHETER CONDITION
ID: {00A7358A-DBE8-CB1D-FFC0-046F3C5E266F}

=====

Page: 353

Abbreviation..... PA CONDITI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... PA Catheter Condition
Help Text..... Condition of PA catheter
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

-
- 1. PATENT CLEAN (PAT CLEAN) {68688246-AE58-2EFB-BF9D-5E00D8CFDA8E}
 - 2. RED (RED) {515D7CFB-3527-30E7-EC4B-B607C272DFA0}
 - 3. SWOLLEN (SWOLLEN) {DCE86B10-C393-8DA9-6518-F2943F0E1A4E}
 - 4. DISCOMFORT (DISCOMFORT) {77D058C9-1EB4-96C0-C9E3-666C840B884C}

UNITS:

=====

TERM: PA CATHETER DRESSING TYPE
ID: {16F7EF6A-E0B6-9931-1F74-039E9D18AC57}

=====

Page: 354

Abbreviation..... PA DRESSIN
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description

Display Name..... PA Catheter Dressing Type
Help Text..... Type of dressing on PA Catheter site
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. TRANSPARENT (TRANSPAREN) {0E73D24D-0EFC-08F9-196A-D30DD1BB3228}
- 2. PRESSURE DRESSING (PRESS DRES)
{6D9DB9EF-85EB-8321-03FC-879AA4C5615B}

UNITS:

TERM: PA CATHETER DYNAMIC RESPONSE TEST
ID: {7ED7B728-5683-4660-B301-79E72C3FD1C9}

Page: 355

Abbreviation..... PA CATH DR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... PA Catheter Dynamic Response Test
Display Name..... PA Catheter Dynamic Response Test
Help Text..... PA Catheter Dynamic Response Test
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

- 1. NORMAL (n) {553C9A97-A3D5-C52F-0E4E-C7D168C99A2B}
- 2. DAMPENED (DAMPENED) {77C777B3-7E56-41E3-8E70-DE144E6A3424}
- 3. HYPERDYNAMIC (HYPERDYNAM) {818A8A99-04D9-4989-A9BA-ECFACED74605}

UNITS:

TERM: PA CATHETER INSERTION MARK
ID: {7E62F0A7-9233-4BD3-870B-189BC86A179E}

Page: 356

Abbreviation..... PA INS MAR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... PA Catheter Insertion Mark
Display Name..... PA Catheter Insertion Mark (cm)
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION

Description..... PA Catheter Waveform
Display Name..... PA Catheter Waveform
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

- 1. NORMAL (n) {553C9A97-A3D5-C52F-0E4E-C7D168C99A2B}
- 2. DAMPENED (DAMPENED) {77C777B3-7E56-41E3-8E70-DE144E6A3424}
- 3. HYPERDYNAMIC (HYPERDYNAM) {818A8A99-04D9-4989-A9BA-ECFACED74605}
- 4. POSITIONAL (POSITIONAL) {9A74D16B-5E95-4EA8-88AF-E53C80311A5D}
- 5. NOT MONITORED (NOT MONITO) {0E6A8A07-48E0-428F-8D86-22572CABCF48}

UNITS:

TERM: PA CATHETER ZEROED Page: 363
ID: {E8B0247F-BD12-4902-87F7-012AA4C1D73D}

Abbreviation..... PA ZEROED
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... PA Catheter Zeroed
Display Name..... PA Catheter Zeroed
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
- 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

TERM: PACEMAKER ATRIAL mA Page: 364
ID: {2551E3D7-CD46-909E-1D79-082B92AEDF7A}

Abbreviation..... PAmA
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Pacemaker Atrial mA

ID: {AC25DEC8-680C-EDD4-1141-43F841651D2A}

=====

Abbreviation..... PMM
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Pacemaker Mode
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 13 entries

- 1. AOO (AOO) {86F25E48-AADB-5B30-4B89-E8BE63DAE7D9}
- 2. AAO (AAO) {E87E0670-2C0F-6E6A-DEC9-F52E5C530A2A}
- 3. AAI (AAI) {2DD18ACA-B98E-24D3-E760-6D082A76074B}
- 4. AAIR (AAIR) {3B4B6F37-B631-BE32-F1C1-65B99D0C5018}
- 5. AVO (AVO) {02EA752A-D156-96E6-3D28-3D389FFC7B58}
- 6. AVI (AVI) {8BB49456-FB80-4916-6EBE-06837BF12B0C}
- 7. VVI (VVI) {2C87F8F3-2B99-80E7-04EB-E6FE88168627}
- 8. VDD (VDD) {E733E525-EF36-2F93-AECC-BAD1C3E2C633}
- 9. DVI (DVI) {4D6E3081-7BB5-45DB-E8E6-0D3D7869E6BD}
- 10. DOO (DOO) {749BCE3D-FF09-B796-FF2A-ECD2D116C96F}
- 11. DDI (DDI) {0C1E359F-B1F4-7994-D337-37622CB14BC8}
- 12. DDD (DDD) {7AA7152E-48F0-8550-AD0D-9F298E90B765}
- 13. DDDR (DDDR) {A3CC0D9D-D329-1A52-2DD1-911A1E04B871}

UNITS:

ID: {A7D28B32-845F-0E16-FDAF-F654FA6EE8A4}

=====

Abbreviation..... PMO
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... No Description
Display Name..... Pacemaker On
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}

2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

LOCATION: 8 entries

- 1. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
- 2. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
- 3. L BRACHIAL (LBRA) {75F6CE33-1AB3-A1DA-70B8-04825290B577}
- 4. R BRACHIAL (RBRA) {A88E69A8-F6AE-80B0-2CD4-687B121F9DDE}
- 5. L INTERNAL JUGULAR (LIJ) {AFBB4FE5-8A11-2047-7900-DADF2FC4430D}
- 6. R INTERNAL JUGULAR (R INT JUGU)
{7ABF6916-D379-198D-6454-F0518A74844A}
- 7. L SUBCLAVIAN (LSCL) {71B93594-15F9-AA4A-AEF2-E222EC28C00D}
- 8. R SUBCLAVIAN (R SUBCLAVI) {53DEAE83-93AE-123B-14E5-4DE38CD4F752}

UNITS:

=====

TERM: PACEMAKER RATE Page: 369
ID: {FCA825C5-12A2-132D-8842-23D56ED6660B}

=====

Abbreviation..... PMR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Pacemaker Rate
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

BEATS PER MINUTE 0 > 45 * 100 < 200 Dec Prec: 0

=====

TERM: PACEMAKER SENSE Page: 370
ID: {7F716D0A-ACA4-7E33-62D7-988BA22F0D7F}

=====

Abbreviation..... PMS
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... No Description
Display Name..... Pacemaker Sense
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: PACEMAKER SENSITIVITY ATRIAL Page: 371
ID: {6CBB234D-694E-CB88-3844-6B9E9E262320}

=====

Abbreviation..... PSA
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Pacemaker Sensitivity Atrial
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

MILLIVOLTS 0.01 > 0.05 * 2.00 < 6.00 Dec Prec: 2

=====

TERM: PACEMAKER SENSITIVITY VENTRICULAR Page: 372
ID: {305A8823-01D3-F3FE-D66E-3CACE736CD30}

=====

Abbreviation..... PSV
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Pacemaker Sensitivity Ventricular
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

MILLIVOLTS 0.01 > 0.05 * 2.00 < 6.00 Dec Prec: 2

=====

TERM: PACEMAKER TYPE Page: 373
ID: {3B5B64D7-9091-8511-64CE-848456CAD2CD}

=====

Abbreviation..... PACEM TYPE
Active..... Yes

Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Pacemaker Type
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Void.....

DATA VALUE: 6 entries

-
1. TRANSCUTANEOUS (t) {F4315F92-8804-1CB4-4D03-CCFDC95C6E71}
 2. PERMANENT (p) {A31FF319-F8C3-7AAC-D7FC-765C89732D09}
 3. TEMPORARY (TRANSVENOUS) (tt)
{F6A74188-3D79-D749-540F-4EA32201B0E6}
 4. EPICARDIAL LEADS (el) {2B551575-3D50-96BD-7BAA-50CC9796361E}
 5. ATRIAL (a) {266D2FDA-5D6A-203D-6FF0-8282D49BFD67}
 6. VENTRICULAR (v) {307A4572-3CF0-7875-7313-7FD4AB23F8E7}

LOCATION: 8 entries

-
1. L BRACHIAL (LBRA) {75F6CE33-1AB3-A1DA-70B8-04825290B577}
 2. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
 3. L INTERNAL JUGULAR (LIJ) {AFBB4FE5-8A11-2047-7900-DADF2FC4430D}
 4. L SUBCLAVIAN (LSCL) {71B93594-15F9-AA4A-AEF2-E222EC28C00D}
 5. R BRACHIAL (RBRA) {A88E69A8-F6AE-80B0-2CD4-687B121F9DDE}
 6. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
 7. R INTERNAL JUGULAR (R INT JUGU)
{7ABF6916-D379-198D-6454-F0518A74844A}
 8. R SUBCLAVIAN (R SUBCLAVI) {53DEAE83-93AE-123B-14E5-4DE38CD4F752}

UNITS:

=====

TERM: PACEMAKER UNIT ID

Page: 374

ID: {CDA2992F-349C-9904-45E3-20D01DC87917}

=====

Abbreviation..... PUID
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... String
Description..... No Description
Display Name..... Pacemaker Unit ID
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

Abbreviation..... P DIS C
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... PACU - Discharge Criteria
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 14 entries

- 1. TOLERATING FLUIDS/NOURISHMENT (TO OR FL)
{E107F272-CD27-46B8-8C50-F9FF776DDCD9}
- 2. ABLE TO AMBULATE (AMBULAT)
{4E1D40CA-F5C4-442E-8BBA-C6E4C4781F57}
- 3. ALERT AND ORIENTED (AL AN OR)
{C031F776-36DF-4570-AFAF-A30B8BFA635B}
- 4. PT DENIES NAUSEA / NO VOMITING (DEN NAUS)
{4374D52E-9D4E-41AD-8B52-FCBDF6B5A1E9}
- 5. ALDRETE SCORE 8 OR GREATER (ALD 8+)
{1FACE6D7-BC92-4579-8665-C5E1D3475461}
- 6. AIRWAY REFLEXES INTACT / SECRETIONS (AIR REFL)
{88A7E7D3-6283-4640-A5FD-0F1E37FB8249}
- 7. PAIN CONTROLLED PER PT (PAIN CON)
{7EC1AF38-39B1-475E-B58C-E33404A47E96}
- 8. PT PAIN AT ACCEPTABLE LEVEL (PT PA OK)
{78489498-4A5B-41C5-98C5-95D56C4A0933}
- 9. VOIDED (VOIDED) {7F23E601-A1E5-4016-8A12-06B1099041CF}
- 10. ABSENT TO MINIMAL BLEEDING (ABS BLEE)
{123460FB-1A8D-45AB-B0A2-489EAB482F67}
- 11. CMS INTACT DISTAL TO SURGICAL SITE (CMS INTA)
{CE2D2CDE-71C5-4993-8E46-C4799301A851}
- 12. SKIN WARM AND DRY COLOR APPROPRIATE (S W A DR)
{9637D510-C199-4969-8F2D-0A499151E22B}
- 13. MAINTAINING EFFECTIVE VENTILATION (EFF VENT)
{18774550-ECBD-412D-AD95-131FAA4CC549}
- 14. ALDRETE SCORE AT BASELINE (ALD BASE)
{CC86D48F-9473-4AD8-AA5B-3BA0D1CCA8FD}

UNITS:

=====

TERM: PACU - DISCHARGE CRITERIA MET
ID: {3B2F659B-ABA4-41E4-9861-B4B0FE3C873D}

=====

Abbreviation..... P DIS ME
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description

Display Name..... Discharge Criteria Met
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

-
- 1. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}
 - 2. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}

UNITS:

=====

TERM: PACU - DISPOSITION Page: 379
ID: {9EADBFE2-2072-4B3F-A077-354B76E5F190}

=====

Abbreviation..... PAC DISP
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... PACU - Disposition
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 12 entries

-
- 1. ICU (ICU) {E7FA69FE-1088-4853-AFBD-2D815E3C2EB0}
 - 2. MICU (MICU) {28244B62-644C-4DA0-A350-EAFE0D2FA7C0}
 - 3. SICU (SICU) {EAEB2DE6-B057-4143-B44C-99DAB9B7B4AC}
 - 4. STEP DOWN (STEP DOW) {F6AC5CA7-0953-4E21-B1A4-6BFE1983D515}
 - 5. WARD (WARD) {A16CB5CB-3401-4F6C-8467-7F78FC929872}
 - 6. AMBULATORY CARE UNIT (ACU) {7FB7FBC5-42A6-4BA4-9CA6-494D77C4C549}
 - 7. PROCEDURE ROOM (PROC ROO) {CC9A2AB0-7A73-4DE7-9B7F-DC2C160AE462}
 - 8. IMAGING (IMAGING) {791A322E-F45E-40EB-B95D-E5B62E1D0A75}
 - 9. OPERATING ROOM (OR) {37DACE14-2A28-45F4-9DEA-35A1A433F4E6}
 - 10. EMERGENCY DEPARTMENT (ER DEPT) {D2E511F7-8B80-4928-A64B-3F1EE17395F7}
 - 11. WAITING ROOM (WAIT ROOM) {13476486-F895-4FFC-AB42-B274827B47DA}
 - 12. HOME (HOME) {C8E5E320-361B-421C-A6A9-13DC16125440}

UNITS:

=====

TERM: PACU - PT RECEIVED REVERSAL AGENT Page: 380
ID: {08A85E74-4365-41AA-8C26-4D712C9364DE}

=====

- 17. POSITIONAL (POSITIONAL) {9A74D16B-5E95-4EA8-88AF-E53C80311A5D}
- 18. PRESSURE (PRESSURE) {387AFB7F-3973-4D77-B2E6-0724F344D8F6}
- 19. RADIATING (RADIATING) {B6807939-F809-45C1-912A-ADEF2E0DB0DC}
- 20. SHARP (SHARP) {20608864-E5CF-4FE3-8B6D-AE2FEE086503}
- 21. SHOOTING (SHOOTING) {73C48DB1-6FE6-46B7-B1F6-E93ACBE44C27}
- 22. SQUEEZING (SQUEEZING) {A148587A-1B99-4D0D-B11E-87ADC3D2A811}
- 23. STABBING (STABBING) {684640D5-B3AB-47DA-98D8-5434ABF22B7E}
- 24. STINGING (STINGING) {DF235170-F501-4E8D-89DD-5A3C76F4B13B}
- 25. TENDER (TENDER) {30486CC7-3097-4F95-9238-FD6C2126CD10}
- 26. THROBBING (THROBBING) {3ACE2B10-1418-47F8-A4F3-0C41A699B5D9}
- 27. TINGLING (t) {AE8E3101-AC3D-DEC3-0293-58F901F82EB8}
- 28. UNBEARABLE (UNBEARABLE) {96611993-45F3-4D1D-BF7F-5E5C46E5569E}

UNITS:

=====

TERM: PAIN DURATION Page: 382
 ID: {DF714694-E547-4D39-9355-ADF3B419F9CD}

=====

Abbreviation..... PAIN DURAT
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... Pain Duration
 Display Name..... Pain Duration
 Help Text..... No Help Text
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 2 entries

- 1. CONSTANT (CONSTANT) {0FB8E43C-2106-4C37-9430-38675953A85B}
- 2. INTERMITTENT (INTERMITTE) {DD09ED7C-B977-412F-A678-5130F77D239E}

UNITS:

=====

TERM: PAIN GOAL Page: 383
 ID: {709DFACE-B211-4160-B291-359777ED8D0F}

=====

Abbreviation..... PAIN GOAL
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... No Description
 Display Name..... Pain Goal
 Help Text..... No Help Text
 Multi Select Picklist..... No
 Term Type..... OBSERVATION

Vuid.....

DATA VALUE: 11 entries

- 1. 0 (0) {2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
- 2. 1 (1) {8F5D10A8-FE17-47BC-DACA-6972D12CC00A}
- 3. 2 (2) {38E22B79-A2B6-3401-B195-C955F24E1278}
- 4. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}
- 5. 4 (4) {D2280C7E-DE50-5A8D-C3FD-5862B69C0E1A}
- 6. 5 (5) {188EDCE7-0726-4234-D8D2-92C8052CD221}
- 7. 6 (6) {640CF2F9-3BBF-AB13-185B-CD02A609FB5A}
- 8. 7 (7) {7A010324-D2C2-5D4B-848A-7B8878666D20}
- 9. 8 (8) {518E7CFC-102E-351E-C7AF-FB575E3C8F4D}
- 10. 9 (9) {58FF2FD4-569F-F454-9F5B-FE6F9E34D8CC}
- 11. 10 (10) {0C827160-18DB-816F-7FAA-A3C175C6D33B}

UNITS:

TERM: PAIN LEVEL Page: 384
ID: {47A83DEA-BA95-38AD-DF2E-1F209122E684}

Abbreviation..... PAIN LEVEL
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Standard pain scale
Display Name..... Pain Level
Help Text..... Enter 0-10, 99 or Unable to respond
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid..... 4500635

DATA VALUE: 12 entries

- 1. 0 (0) {2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
- 2. 1 (1) {8F5D10A8-FE17-47BC-DACA-6972D12CC00A}
- 3. 2 (2) {38E22B79-A2B6-3401-B195-C955F24E1278}
- 4. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}
- 5. 4 (4) {D2280C7E-DE50-5A8D-C3FD-5862B69C0E1A}
- 6. 5 (5) {188EDCE7-0726-4234-D8D2-92C8052CD221}
- 7. 6 (6) {640CF2F9-3BBF-AB13-185B-CD02A609FB5A}
- 8. 7 (7) {7A010324-D2C2-5D4B-848A-7B8878666D20}
- 9. 8 (8) {518E7CFC-102E-351E-C7AF-FB575E3C8F4D}
- 10. 9 (9) {58FF2FD4-569F-F454-9F5B-FE6F9E34D8CC}
- 11. 10 (10) {0C827160-18DB-816F-7FAA-A3C175C6D33B}
- 12. 99 (9) {B013C05B-1AF0-1547-4F03-378883F887C5}

LOCATION: 58 entries

- 1. HEAD (Hd) {C3209F44-90C6-FF1B-3169-265B294F276A}
- 2. FACE (FACE) {C638EA14-7903-6A94-496A-C6560D51BEE3}
- 3. NECK (Ne) {D4812C3B-AAC2-1B61-679B-3CD8B5CB947C}

4. STOMACH (STOMACH) {613C20B9-CC9C-4F5D-B5B1-69B66A663540}
 5. ABDOMEN (ABDOMEN) {109FE041-2521-6EF2-7F95-1FD1587796AB}
 6. BACK (BACK) {22C26239-C692-4E93-B522-E5E8392BA3B1}
 7. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
 8. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
 9. L HAND (L HAND) {639DF25A-19E3-B504-D323-4B3F6ED39A24}
 10. R HAND (R HAND) {A54D3A80-8537-D079-9CCE-7DC56346D132}
 11. L HIP (L HIP) {786CD881-27AD-4136-B9CC-8423CB38A431}
 12. R HIP (R HIP) {A48CBCAE-5404-4934-9321-ECEF660A5350}
 13. L GROIN (L GROIN) {5488E4EE-C1A2-4B02-8763-A53C412117A7}
 14. R GROIN (R GROIN) {BAE1C889-8CCF-4CEB-9EBD-0243BC6F096A}
 15. L THIGH (L THIGH) {16C2312F-1E0E-45E0-B675-160F40C23F57}
 16. R THIGH (R THIGH) {B7A2D4A0-801E-4578-967A-F1CFB86FB96B}
 17. L KNEE (L KNEE) {3A9E30B4-09A9-967A-58EF-69A01EC4F859}
 18. R KNEE (R KNEE) {BCD18D53-07B6-2C41-64FC-3A07A7E979E7}
 19. L FOOT (L FOOT) {2EC3FBEC-4EA2-9E7A-17D4-FDC2F9D28046}
 20. R FOOT (R FOOT) {7D9F5FF6-232C-2488-4504-A721FF06928F}
 21. L LEG (LL) {B088C809-E528-B785-21C2-14E45E1FD8AF}
 22. R LEG (RL) {C5E40A7B-6FC4-267D-8CE8-3571D21D81FD}
 23. COCCYX (COCCYX) {49D65E9D-3D5B-40E3-B7CF-1003F80B852C}
 24. GASTROINTESTINAL (GI) {C27C66DB-3A6A-487E-9981-FA7FD323C0C3}
 25. L EAR (L EAR) {DA9C2DD1-D6FF-4899-A62B-1096FAA36C59}
 26. R EAR (R EAR) {AE31327D-A03C-488C-8A8E-084E7E6E7FDE}
 27. L SHOULDER (L SHOULDER) {91CBD2FA-2743-42EB-946B-6933980D9337}
 28. R SHOULDER (R SHOULDER) {D27AE130-442E-4DD6-B391-8D7E8AD494A2}
 29. L ELBOW (L ELBOW) {BF4317E7-2024-47EF-9E5F-3D21D0547885}
 30. R ELBOW (R ELBOW) {5BAFF87B-6993-4E76-862F-9F386D9B71C5}
 31. L ANKLE (L ANKLE) {D5A692EF-9D7C-473C-AF96-7EF1C01FA2AF}
 32. R ANKLE (R ANKLE) {4E6DF352-F399-4CE9-908C-48132C3D856B}
 33. PERINEAL (PERINEAL) {1801441C-D5EF-4775-9E58-E4985496F7E0}
 34. BLADDER (BLADDER) {18E083F2-C4BB-407E-B0A3-20C4942FC97F}
 35. L JAW (L JAW) {65EBD3CF-76D8-4F45-9B1B-3575D988BE8E}
 36. R JAW (R JAW) {7846AEE2-B6F4-4B5A-909F-56E386FD82A8}
 37. TOOTH (TOOTH) {8407F780-7C7F-49BF-B1DD-4FDD578CEAB3}
 38. TEETH (TEETH) {38B36D8E-A9C3-4EE6-8727-EECDA4C6034E}
 39. L NOSE (L NOSE) {C7707847-5135-4749-B196-002E3E59BB54}
 40. R NOSE (R NOSE) {D80BCB19-E71F-43EA-8990-6C121EE7F7A1}
 41. NOSE (NOSE) {85F7003E-6D4E-4A79-A534-8A57FBC2D106}
 42. PENIS (PENIS) {2D82CE55-83B6-4160-A754-07C71FD388E8}
 43. ABDOMEN RUQ (ABD RUQ) {F29C9FA1-717A-49B7-8000-95070E9AACAO}
 44. ABDOMEN LUQ (ABD LUQ) {3B31E1EC-CE74-4C76-BB61-10C678431FCA}
 45. ABDOMEN RLQ (ABD RLQ) {2C48394F-6877-4738-856C-E933644D8067}
 46. ABDOMEN LLQ (ABD LLQ) {39831F59-D7B0-4156-AC57-4C622A93CA58}
 47. INCISIONAL SITE (INCISON SI) {EB36A66D-751D-4031-9F68-0BB03E267321}
 48. CHEST (CHEST) {8C3C30FE-6E06-C35E-8226-CBF0C94014CC}
 49. L CHEST (L CHEST) {E332E497-FE0C-47CB-A0A7-217FBF52DC97}
 50. R CHEST (R CHEST) {7E30EC0E-815A-4C2E-A9FC-3C2B2D6E3CD4}
 51. L SCAPULA (L SCAPULA) {3381EC50-771B-41C8-85C8-6EB0D208D54E}
 52. R SCAPULA (R SCAPULA) {EA657261-4D72-40C2-B70D-2A2E27556A01}
 53. L UPPER ARM (L UPP ARM) {5F915696-1B08-32E9-6B95-49246117D660}
 54. R UPPER ARM (R UPP ARM) {3EF20E8E-F68C-E4F2-8F8C-FAD3B44DDAA4}
 55. SUBSTERNAL (SUBSTERN) {3DE7E0F1-BCA4-43EC-82CB-A87C6506EF46}
 56. L LOWER ARM (L LOW ARM) {63AA95E8-256D-4390-8CB3-F00A1B3F91CD}
 57. R LOWER ARM (R LOW ARM) {8DE6E277-EFFF-40DD-9E83-D7295AACB6AA}
 58. SCROTUM (SCROTUM) {BDC109F5-0D30-4914-A62C-CEFEB448B547}

UNITS:

=====

TERM: PAIN NON-VERBAL CUES

Page: 385

ID: {F875330C-D06C-41A7-877F-FBF07596EFF6}

=====

Abbreviation..... PAI NONVER
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Pain Non-Verbal Cues
Display Name..... Pain Non-Verbal Cues
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 10 entries

1. CRYING (CRYING) {D87AB50C-0B18-40DD-884C-6AAEE4DD97FC}
2. GRIMACING (GRIMACING) {69EFDB2D-A0C4-48D4-A7E0-A862F34545C8}
3. MOANING (MOANING) {A71E3515-AA20-4B0E-9CB8-1612F76BA406}
4. RESTLESS (r) {8FCCAD26-DBA0-7067-DA72-30D342BFA8E8}
5. SPLINTING (SPLINTING) {D7FC5AA9-E63A-4A88-BE03-A50D0A8751F9}
6. GUARDING (GUARDING) {9EFD10A7-B1B7-43DA-B289-646AF5C0A25C}
7. PAINFUL EXPRESSION (PAINFUL EX)
{608554F9-29F9-4B27-BE63-486AF81810BD}
8. IRRITABLE (Irritable) {10651C11-2736-4056-BA84-E544358BA877}
9. SLEEPING (s) {87338420-F351-F98C-7A80-4990B87ACE2B}
10. NO SIGNS AND SYMPTOMS (NO SIG&SYM)
{91285942-6C82-4CB1-AAFB-8D3DA15A0B8A}

UNITS:

=====

TERM: PAIN RELIEVING FACTOR INTERVENTION

Page: 386

ID: {325DB732-E996-4DF5-9F5F-EF3B55D80935}

=====

Abbreviation..... PN FACT IN
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Pain Relieving Factor Intervention
Display Name..... Pain Relieving Factor Intervention
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 9 entries

- 1. COLD (COLD) {B7B422BA-22EC-117D-0A0C-8419E9447D46}
- 2. HEAT (HEAT) {E0FCA383-3AB8-4E03-9ABC-00DE3FE01431}
- 3. QUIET ROOM (qr) {CCD96D41-C273-F6F3-DA18-8818B1CC9980}
- 4. REPOSITION (REPOSITION) {DF3039CA-5977-4454-81D4-DB0AF966F85D}
- 5. MEDICATION (MEDICATION) {22E6DEBB-9662-466B-990C-8DA500C0A183}
- 6. REFUSE MEDICATION (REFUSE MED)
{17D6956B-4CF6-492C-B113-B5A93D9E215E}
- 7. REQUEST FOR MORE MEDICATION (REQ MOR ME)
{0639081C-BEEC-4CAD-B420-5F2C4DBA38D0}
- 8. PRE-MEDICATED FOR ACTIVITY (PRE-MED AC)
{D6DA772A-FA22-4170-B928-855A2F68711F}
- 9. PHYSICIAN NOTIFIED (PHYS NOTIF)
{4BD68949-74F6-42A6-AEF5-E9A33DC31FFC}

UNITS:

=====

TERM: PAINAD BODY LANGUAGE Page: 387
 ID: {6DF1DD0D-C735-4D34-ABF6-A274621EC177}

=====

Abbreviation..... PAINAD BL
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... No Description
 Display Name..... PAINAD Body Language
 Help Text..... Pain Assessment in Advanced Dementia
 (PAINAD)
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 3 entries

- 1. PAINAD BL 0 (PAINAD BL0) {C77929DB-A59A-4A24-B57B-BF111770C2AB}
- 2. PAINAD BL 1 (PAINAD BL1) {CC1E1516-086D-4B79-8264-0C3BD67D8552}
- 3. PAINAD BL 2 (PAINAD BL2) {34C3E304-DF66-40A7-906F-22AEA3F76FAD}

UNITS:

=====

TERM: PAINAD BREATHING Page: 388
 ID: {EBF3C2FD-5906-4FC8-8E75-CF998FA797C3}

=====

Abbreviation..... PAINAD BR
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList

Description..... No Description
Display Name..... PAINAD Breathing
Help Text..... Pain Assessment in Advanced Dementia
(PAINAD)
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

- 1. PAINAD BR 0 (PAINAD B0) {58E71B2E-F5C3-419F-8B71-A9FD701E0202}
- 2. PAINAD BR 1 (PAINAD BR1) {727D1B04-A702-4806-857B-78BCBACB0FE1}
- 3. PAINAD BR 2 (PAINAD BR2) {37978012-ADAC-416C-8920-48CDAAA5D59D}

UNITS:

TERM: PAINAD CONSOLABILITY Page: 389
ID: {869E7E2E-04C5-46A5-8811-C7F94E066A2C}

Abbreviation..... PAINAD CO
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... PAINAD Consolability
Help Text..... Pain Assessment in Advanced Dementia
(PAINAD)
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

- 1. PAINAD CO 0 (PAINAD CO0) {06A89B55-5822-42F7-A689-DBD91B62979F}
- 2. PAINAD CO 1 (PAINAD CO1) {8E5DDB47-1955-40AD-A008-A9CE4646A707}
- 3. PAINAD CO 2 (PAINAD CO2) {1DFFB735-6272-4B60-8F42-CCCAA701C7A7}

UNITS:

TERM: PAINAD FACIAL EXPRESSION Page: 390
ID: {CB17C4F3-2E4B-4355-9E63-DE165CE47A25}

Abbreviation..... PAINAD FE
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Pain Assessment in Advanced Dementia
(PAINAD)

Display Name..... Painad Facial Expression
Help Text..... Pain Assessment in Advanced Dementia
(PAINAD)
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

- 1. PAINAD FE 0 (PAINAD FE0) {B3DE1BC6-E7A0-48C5-8C8E-FBE033AD1046}
- 2. PAINAD FE 1 (PAINAD FE1) {66C01EE9-8F36-40B5-AED2-DBC2E94F860D}
- 3. PAINAD FE 2 (PAINAD FE2) {35688D06-94F2-4B36-8940-529803D4BC43}

UNITS:

TERM: PAINAD NEGATIVE VOCALIZATION Page: 391
ID: {815392EF-49E0-48D9-BCC0-39FD3333B255}

Abbreviation..... PAINAD NV
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Pain Assessment in Advanced Dementia
(PAINAD)
Display Name..... PAINAD Negative Vocalization
Help Text..... Pain Assessment in Advanced Dementia
(PAINAD)
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

- 1. PAINAD NV 0 (PAINAD NV0) {63F7F990-2DEA-4B35-AC52-D23B8E13B1FC}
- 2. PAINAD NV 1 (PAINAD NV1) {D517E6C9-9D98-4D0F-A8F4-261BED4812D1}
- 3. PAINAD NV 2 (PAINAD NV2) {3C8D6B1E-53E6-4D2C-BD8A-9C548466E1D4}

UNITS:

TERM: PAINAD TOTAL Page: 392
ID: {4E35BBC8-DD3C-4EEA-A0BD-DADBD1252F28}

Abbreviation..... PAINAD TO
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... PAINAD Total Score

passed
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

- 1. 1 HOUR (1 HOUR) {319BA174-1624-483C-B666-4998127B0800}
- 2. 4 HOUR (4 HOUR) {91B2A49D-0308-47C9-B220-63A0E318A201}
- 3. 6 HOUR (6 HOUR) {0D1F000A-01F0-422C-9CB8-2A04CD479444}
- 4. 8 HOUR (8 HOUR) {05E46D8F-C165-4552-8063-B735DECD4AF3}
- 5. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}

UNITS:

TERM: PCA - LINE CARE Page: 400
ID: {B84F9A0A-4353-4480-9DCB-F44437E41F49}

Abbreviation..... PCA LINE C
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... PCA Line Care
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 10 entries

- 1. BAG CHANGED (BAG CHANGE) {818990DE-4F95-4F8E-92D8-B4D3B975CB36}
- 2. CENTRAL VENOUS CATHETER (CENT VEN C)
{EE91E841-C09B-4432-A121-8EBAB92E9CD2}
- 3. DISCONTINUED (DISCONTINU) {15E289FB-09DF-4731-9CF3-48A4E0031C31}
- 4. DRESSING CHANGED (DRESS CHAN)
{E8CF4B7A-D1AB-49A2-8E7C-3375C3BC436F}
- 5. HEPLOCK (HEPLOCK) {C827FEDE-59D6-4C45-BDDA-4884FA0882F3}
- 6. PATENT CLEAN (PAT CLEAN) {68688246-AE58-2EFB-BF9D-5E00D8CFDA8E}
- 7. SALINE LOCK (SALIN LOCK) {A94D7639-2A45-4319-A6C7-5A892AA6B1F4}
- 8. STARTED (STARTED) {D6E105C4-5EAD-424E-8E4C-195EB7783A3D}
- 9. SYRINGE CHANGED (SYR CHANGE)
{BE84F060-7E03-49DD-A6BC-787375C8E1F1}
- 10. TUBING CHANGED/LABELED (TUBE CHANG)
{19BA337C-6E28-411E-9702-48B262C593FA}

UNITS:

TERM: PCA - LOCK OUT Page: 401

Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... PCA - Pump Cleared
Help Text..... PCA Pump Usage History Cleared
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

1. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}
2. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}

UNITS:

=====

TERM: PCA - REQUESTS Page: 404
ID: {851F9560-C9AB-4F6C-AB24-283F64B5943A}

=====

Abbreviation..... PCA REQUES
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... PCA - Requests
Help Text..... Number of times patient requested
analgesia
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

EVENT(S) -1 > 0 * 1000 < 1001 Dec Prec: 0

=====

TERM: PCA - VOLUME REMAINING Page: 405
ID: {7A6CB341-8E80-4F5C-B399-857409777BB1}

=====

Abbreviation..... PCA VOL RE
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... PCA - Volume Remaining
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

MILLILITER -1.00 > 0.00 * 8000.00 < 8001.00
Dec Prec: 2

=====

TERM: PCA - VOLUME START Page: 406
ID: {B5644934-F363-4BE1-BC2D-473BDECBA9AC}

=====

Abbreviation..... PCA VOL ST
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... PCA - Volume Start
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

MILLILITER -1.00 > 0.00 * 8000.00 < 8001.00
Dec Prec: 2

=====

TERM: PEAK INSPIRATORY PRESSURE Page: 407
ID: {05BD86DC-10C3-C4EB-1359-C2406381F357}

=====

Abbreviation..... PK INSP PR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Peak Inspiratory Pressure
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

METHOD: 1 entries

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

CENTIMETERS H2O 0 > 5 * 100 < 200 Dec Prec: 0

=====

TERM: PLAN OF CARE - REVISION INDICATED Page: 408

ID: {DE8BBB70-21B4-4955-BEC4-C1E71BBA3E31}

Abbreviation..... POC REV I
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... No Description
Display Name..... Plan Of Care Revision Indicated
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
- 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

TERM: PLANTAR FLEXION Page: 409
ID: {4BBDBB76-053B-CA20-908E-017F5956DF5E}

Abbreviation..... PLANFL
Active..... Yes
Boolean Value False..... Absent
Boolean Value True..... Present
Data Type..... Boolean
Description..... No Description
Display Name..... Plantar Flexion
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
- 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

TERM: PR INTERVAL Page: 410
ID: {3657545A-7138-94C0-03A5-3136C51B0369}

Abbreviation..... PR INTERVA
Active..... Yes
Boolean Value False.....
Boolean Value True.....

Assessment
Help Text..... Does the pressure ulcer dressing
prevent a full assessment?
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: PRESSURE ULCER LENGTH Page: 415
ID: {B31F418F-99E1-4273-BB5B-3C56745B38C5}

=====

Abbreviation..... PR UL LENG
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Pressure Ulcer Length
Display Name..... Pressure Ulcer Length
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

CENTIMETER -1.00 > 0.00 * 50.00 < 100.00
Dec Prec: 2

=====

TERM: PRESSURE ULCER ODOR Page: 416
ID: {48FE49CA-D1FE-4812-ACB2-D4540B0EC2B9}

=====

Abbreviation..... PR UL ODOR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Pressure Ulcer Odor
Display Name..... Pressure Ulcer Odor
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
- 2. FOUL SMELLING (FOUL SMELL) {761346A9-329D-4F42-8D80-0E9D8D40653C}

UNITS:

=====

TERM: PRESSURE ULCER PHOTOGRAPHED Page: 417
 ID: {453C53C4-516D-4FE6-9032-496517880D3F}

=====

Abbreviation..... PR UL PHOT
 Active..... Yes
 Boolean Value False..... No
 Boolean Value True..... Yes
 Data Type..... Boolean
 Description..... Was the pressure ulcer photographed?
 Display Name..... Pressure Ulcer Photographed
 Help Text..... Was the pressure ulcer photographed?
 Multi Select Picklist.....
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
- 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: PRESSURE ULCER PREVENTION DEVICE TYPE Page: 418
 ID: {4C16776B-B53E-441C-AA73-BA1083D9F626}

=====

Abbreviation..... PR ULC BOO
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... Pressure Ulcer Prevention Boots
 Display Name..... Pressure Ulcer Prevention Devices
 Help Text..... No Help Text
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 8 entries

- 1. ELBOW PAD (ELBOW PAD) {26942F81-1E33-4962-B772-A96FFBCED814}
- 2. FOAM AIR HEEL ELEVATOR (FO AI HEEL) {4EC26AA3-4577-4808-8728-F360C5C5A0D4}
- 3. HEEL LIFT (HEEL LIFT) {974F90E6-382D-421A-A4CE-94A74CBEFC32}
- 4. HEEL PILLOW (HEEL PILLO) {F0F1D997-C780-40BF-AA88-78262E81BFAA}

- 5. INFLATABLE HEEL LIFT (INF HEEL L) {793B34A4-2F21-42DA-881D-1E6E0B4EAD6B}
- 6. RIGID HEEL SPLINT (RIG HEEL S) {5A231387-9FAC-4B2D-94D6-5DF4312E9AF8}
- 7. MULTIPODUS (MULTIPODUS) {ED098F9D-B204-4310-813A-98D91FF9DD25}
- 8. WAFFLE (WAFFLE) {29A0FA73-CAA7-4660-9861-79A09F6DD98C}

UNITS:

=====

TERM: PRESSURE ULCER RISK FACTORS
ID: {A6658726-4361-4B3C-932D-AF2C0F9D6E54}

=====

Page: 419

Abbreviation..... PRE ULC RF
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Pressure Ulcer Risk Factors
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 9 entries

- 1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
- 2. BARIATRIC PATIENT (BARIAC PAT) {6F614334-685E-4533-A475-BDF5440855EA}
- 3. DEVICE RELATED PRESSURE (DEV REL PR) {3AED6107-A173-4B53-8EB5-2FB2B3336FC5}
- 4. END OF LIFE CARE (END LIFE) {047D8977-0567-4B90-AD6A-3AAF8DB35238}
- 5. HYPOALBUMINEMIA (HYPOALBUM) {A3397F38-2491-48E7-8366-A6FA2D91E866}
- 6. VASOPRESSOR (VASOPRESS) {591FCB09-1287-48C6-8007-4CE961296BA9}
- 7. REFUSE TO MOVE (REF TO MOV) {E14588CE-2498-4FE8-AB0B-E954DC1696AC}
- 8. UNSTABLE (UNSTABLE) {5841C521-8070-4034-A69D-27DAB2990EB3}
- 9. VERY LOW BMI (LOW BMI) {78B3A37B-938B-42E5-A403-75E206722BC3}

UNITS:

=====

TERM: PRESSURE ULCER SINUS/UNDERMINING
ID: {CAE9260D-3B9E-433E-8FCD-F5D98A0EB02D}

=====

Page: 420

Abbreviation..... PR UL SINU
Active..... Yes
Boolean Value False..... Not Present
Boolean Value True..... Present

Data Type..... Boolean
Description..... Pressure Ulcer Sinus/Undermining
Display Name..... Pressure Ulcer Sinus/Undermining
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

-
- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
 - 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: PRESSURE ULCER STAGE/SITE
ID: {4DA7EC17-AF79-4FCC-A75F-F1CEB4FE1172}

=====

Page: 421

Abbreviation..... PR UL STAG
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Pressure Ulcer Stage
Display Name..... Pressure Ulcer Stage
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 12 entries

-
- 1. NONBLANCHABLE ERYTHEMA (N BLAN ERY)
{30273ACE-03F3-4818-995B-26586AAB5499}
 - 2. REDDENING OF THE SKIN (RED O SKIN)
{716AEC57-0835-4D07-A78C-D8453599A3D2}
 - 3. NO BREAK IN EPIDERMIS (NO BREAK E)
{615D96CD-6096-446B-BDDF-3F7242936AF9}
 - 4. PART THICKNESS SKIN LOSS INVOLVE EPIDERMIS/DERMIS (PART EPIDE)
{EAD3A260-CC46-4B47-AA49-91CFE13EEE4B}
 - 5. FULL THICK SKIN LOSS INVOLVING SUBCUTANEOUS TISSUE (FULL SKINL)
{14C31B14-1D35-4CCB-81A3-C6E30801950D}
 - 6. FULL THICK SKIN LOSS INVOLVING MUSCLE, BONE, SUPPO (F SK L MBS)
{E027FBEF-1C83-41B7-B848-26D9D62F4266}
 - 7. SUSPECTED DEEP TISSUE INJURY (SUS DEEP T)
{E82B4082-646F-4B5A-B68E-BB0872A4F6D2}
 - 8. STAGE I (STAGE I) {B376AAAF-7143-4422-BCA2-CACF1201817C}
 - 9. STAGE II (STAGE II) {B520BEA3-841E-490A-A3CD-63D1027B561F}
 - 10. STAGE III (STAGE III) {A2F0866C-ECEE-4FCD-A34B-498232662D7E}
 - 11. STAGE IV (STAGE IV) {284334AA-E488-4D45-93EF-9B93B2086A63}
 - 12. UNSTAGEABLE (UNSTAGEABL) {5C507F64-2EC0-4FD8-A786-94D76D90C0D1}

LOCATION: 51 entries

-
1. COCCYX (COCCYX) {49D65E9D-3D5B-40E3-B7CF-1003F80B852C}
 2. HEAD OCCIPITAL (HEAD OCCIP)
{2066CA47-5A16-447F-9919-24ED435B627F}
 3. L ANKLE (L ANKLE) {D5A692EF-9D7C-473C-AF96-7EF1C01FA2AF}
 4. L ANTECUBITAL (L ANTECUBI)
{C21838CC-BAB5-1D7F-E2A4-A7C663FB28C3}
 5. L EAR (L EAR) {DA9C2DD1-D6FF-4899-A62B-1096FAA36C59}
 6. L ELBOW (L ELBOW) {BF4317E7-2024-47EF-9E5F-3D21D0547885}
 7. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
 8. L FOOT (L FOOT) {2EC3FBEC-4EA2-9E7A-17D4-FDC2F9D28046}
 9. L GLUTEUS (L GLUTEUS) {88C4E547-04CD-4D52-9416-CAEB6D898D21}
 10. L HEAD (L HEAD) {3737DC2D-DC9E-4FC8-839D-0355987B62D1}
 11. L HEEL (L HEEL) {C9154881-17C7-4921-A0EE-BFB962C59CEF}
 12. L INNER KNEE (L INN KNEE) {79853581-CCF0-49E9-9D91-D16DE02A836D}
 13. L ISCHIAL (L ISCHIAL) {DCDAD741-9599-4479-9D55-3279F800DBE6}
 14. L LOWER ARM (L LOW ARM) {63AA95E8-256D-4390-8CB3-F00A1B3F91CD}
 15. L OUTER KNEE (L OUT KNEE) {B143C935-2687-49C3-9F65-8D69C1775ABC}
 16. L RIBS (L RIBS) {BDFEE74F-B3CC-4AD9-AE49-78B0BF19A219}
 17. L SCAPULA (L SCAPULA) {3381EC50-771B-41C8-85C8-6EB0D208D54E}
 18. L SHOULDER (L SHOULDER) {91CBD2FA-2743-42EB-946B-6933980D9337}
 19. L TROCHANTER (L TROCHANT) {C91647A4-255D-40E4-A24D-01274E5D799F}
 20. L UPPER ARM (L UPP ARM) {5F915696-1B08-32E9-6B95-49246117D660}
 21. PERINEAL (PERINEAL) {1801441C-D5EF-4775-9E58-E4985496F7E0}
 22. R ANKLE (R ANKLE) {4E6DF352-F399-4CE9-908C-48132C3D856B}
 23. R ANTECUBITAL (R ANTECUBI)
{7FB00A65-B4B2-06A6-BD89-98AC112D8D4E}
 24. R EAR (R EAR) {AE31327D-A03C-488C-8A8E-084E7E6E7FDE}
 25. R ELBOW (R ELBOW) {5BAFF87B-6993-4E76-862F-9F386D9B71C5}
 26. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
 27. R FOOT (R FOOT) {7D9F5FF6-232C-2488-4504-A721FF06928F}
 28. R GLUTEUS (R GLUTEUS) {B4D32045-68D1-4457-9106-7937887FA6F6}
 29. R HEAD (R HEAD) {42A95AAE-3DA9-4C96-A4E7-A15C9840217A}
 30. R HEEL (R HEEL) {0020DD20-2EA4-4362-ADF0-07BFA220C60E}
 31. R INNER KNEE (R INN KNEE) {2AE24BF1-EA1B-490D-9DB7-EEFC9314C560}
 32. R ISCHIAL (R ISCHIAL) {97E50F6C-CF75-4E59-AFDE-80FC5CE6D9E3}
 33. R LOWER ARM (R LOW ARM) {8DE6E277-EFFF-40DD-9E83-D7295AACB6AA}
 34. R OUTER KNEE (R OUT KNEE) {330F1DF1-8C99-4984-9965-7404ABC9562D}
 35. R RIBS (R RIBS) {260FD564-5E58-419B-9F73-62B6037D1572}
 36. R SCAPULA (R SCAPULA) {EA657261-4D72-40C2-B70D-2A2E27556A01}
 37. R SHOULDER (R SHOULDER) {D27AE130-442E-4DD6-B391-8D7E8AD494A2}
 38. R TROCHANTER (R TROCHANT) {1A3E21B2-C391-40D2-8059-12290E667D9E}
 39. R UPPER ARM (R UPP ARM) {3EF20E8E-F68C-E4F2-8F8C-FAD3B44DDAA4}
 40. SACRAL (SACRAL) {0CF8A8A7-AB35-B803-4544-0C91645D9D1C}
 41. FOREHEAD (FOREHEAD) {5F892722-956E-474C-BA08-4074142D49F6}
 42. L LATERAL HEEL (L LAT HEEL)
{1A10FB20-7BAD-4D2B-9CBD-1A1A01DDCBF1}
 43. R LATERAL HEEL (R LAT HEEL)
{A1896CF9-0246-4323-8DED-E4B5FB12CACE}
 44. L MEDIAL HEEL (L MED HEEL)
{47AD7BA6-4EBB-4A08-B03C-FE73EB240725}
 45. R MEDIAL HEEL (R MED HEEL)
{C0526527-0AA9-4ECF-B768-99B2A1F5F3BF}
 46. L LATERAL MALLEOLUS (L LAT MALL)
{9A523B7C-64DD-4E47-9543-165788108D4C}
 47. R LATERAL MALLEOLUS (R LAT MALL)

- {F6192D10-D6EA-42C2-AB15-270659D9925C}
- 16. OPEN TO AIR (OPE TO AIR) {EFB715D0-5666-4281-9034-A0326A14B2BB}
- 17. NEGATIVE PRESSURE THERAPY DEVICE (NEG PRE TH)
{8AF0A069-66AD-4114-99D8-238AD987F0FF}
- 18. SALINE MOIST DRESSING (SAL MOI DR)
{19F98514-8ABF-461C-9969-C31177FE7DE7}
- 19. TRANSPARENT DRESSING (TRANSP DRE)
{777C7722-AEA8-4A0B-9ADD-1157D37AE539}
- 20. PETROLATUM (PETROLATUM) {DA297C75-4776-4C8D-B0D7-A0FF0E4D1B16}
- 21. PORCINE SUBMUCOSA (PORC SUBMU)
{6F3917DA-74BD-4790-9D4D-B284E27B7623}
- 22. SILVADENE (SILVADENE) {0889134F-C419-4366-B12E-EAB775719410}
- 23. XENADERM (XENADERM) {F181D70E-A231-4D2B-A1E6-DDF4DBED31CD}
- 24. XEROFORM (XEROFORM) {7D7C3CDC-253E-4262-B8D4-3B70DE9AEBE4}

UNITS:

=====

TERM: PRESSURE ULCER TUNNEL Page: 423
 ID: {584775CB-EB35-47BB-A7D4-62966A23D0DF}

=====

Abbreviation..... PRE ULC TU
 Active..... Yes
 Boolean Value False..... No
 Boolean Value True..... Yes
 Data Type..... Boolean
 Description..... No Description
 Display Name..... Pressure Ulcer Tunnelling
 Help Text..... Is the Pressure Ulcer Tunnelling?
 Multi Select Picklist.....
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
- 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: PRESSURE ULCER TUNNEL LENGTH Page: 424
 ID: {D7FC23D9-8C2B-42F0-A545-5102AFC22335}

=====

Abbreviation..... PR UL TU L
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... Numeric
 Description..... No Description
 Display Name..... Pressure Ulcer Tunnel Length
 Help Text..... No Help Text

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
- 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: PULMONARY SECRETION AMOUNT Page: 432
 ID: {FEED2CE3-F697-D1AB-AB1B-72A7968F0EB0}

=====

Abbreviation..... PUL SEC AM
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... AMOUNT OF SPUTUM FROM PATIENT
 Display Name..... Pulmonary Secretion Amount
 Help Text..... No Help Text
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 5 entries

- 1. SCANT (SCANT) {8EF6C62D-0D94-37BB-60AD-86E306E8120D}
- 2. SMALL (s) {5991F892-B28C-3998-E733-D8FFA26D7E21}
- 3. MODERATE (MODERATE) {C3C25064-C9FC-EA1F-3E74-32A25F0FC5E3}
- 4. LARGE (l) {E60D5CF9-920D-44E2-8A83-136DAB1FE251}
- 5. COPIUS (COPIUS) {C74A6979-0E35-AA12-F992-69ACA3B0FA23}

UNITS:

=====

TERM: PULMONARY SECRETION COLOR Page: 433
 ID: {52529F3D-B461-1C5E-892A-1E46DFB8C003}

=====

Abbreviation..... SPUTUM COL
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... Color of sputum
 Display Name..... Pulmonary Secretion Color
 Help Text..... No Help Text
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 8 entries

- 1. CLEAR (c) {8981C365-2B46-704F-C9F4-975DF1C23C2D}
- 2. WHITE (WHITE) {AB04EF54-5246-EF3A-E1D9-02E9E24369DD}
- 3. YELLOW (YELLOW) {90DC31AD-6211-9F25-4237-79F7122A2763}

Vuid..... 4500636

LOCATION: 12 entries

- 1. APICAL (Ap) {C5CF5C51-63CC-33FB-9BAE-FB1DD9900911}
- 2. AV FISTULA (AVF) {F977EEF6-EA6A-D56D-D1F7-29B91510A02F}
- 3. BILATERAL PERIPHERALS (BILAT PERI)
{7F134529-36C5-E9D4-ABD5-64B404F9E706}
- 4. BRACHIAL (Bra) {67426127-09E1-E566-F31E-15714A5B2B7C}
- 5. CAROTID (Car) {9536F203-1797-9336-6F69-604249E1B3A0}
- 6. DORSALIS PEDIS (DP) {9657D743-1632-2044-B4C5-9EF291A2F6D0}
- 7. FEMORAL (Fem) {E3442C9E-8677-7296-D314-7365319C3716}
- 8. PERIPHERAL (Per) {36B31698-E84A-5792-8EAE-F374DFD1B55A}
- 9. POPLITEAL (Pop) {77030DF2-CE06-D9E2-0684-EE70209CBE73}
- 10. POSTERIOR TIBIAL (PT) {010F2FFE-DEBC-DF9D-9CF2-AD192B16FFD5}
- 11. RADIAL (RADIAL) {D52C4103-2AB3-A0C7-16E3-CF95A1B3B789}
- 12. ULNAR (Uln) {81FC2EA0-472D-06D7-F5A2-E71B0B7DE091}

METHOD: 13 entries

- 1. AFTER EXERCISE (AE) {D129EF4B-3A99-9966-54F1-2C01A36814EE}
- 2. AT REST (AR) {E8C289A6-0C3D-A7B2-BB85-7F1A23F22B50}
- 3. AUSCULTATE (Aus) {0C629C1F-92D3-C641-B02C-9EA43B401DC9}
- 4. DOPPLER (Dop) {52683C19-F9DC-876B-CA3F-53E7C864CF04}
- 5. INVASIVE (I) {5DB8746E-5446-F47A-883C-A44B91275A6E}
- 6. LEAD (L) {B628FDCE-05E0-DACF-E66E-0F8380AB465A}
- 7. MONITOR (Mo) {63FC24EA-50B8-B302-D948-ABE3517CF6FC}
- 8. NON-INVASIVE (N-I) {F50ABFED-5FD7-FC95-20FE-1E3038C9F299}
- 9. PALPATED (Pal) {295D8E72-12A1-08AF-6BE8-0F9CD238F5EF}
- 10. WITH ACTIVITY (Act) {544A262B-EA34-F82A-46A9-F80740515D65}
- 11. ELECTROCARDIOGRAM (ECG) {CF5F1AB9-3B74-7FBA-E38B-C9EF54FE4DC5}
- 12. ELECTROCARDIOGRAM TRACE 1 (ECG TR 1)
{E0D449BA-9F06-4A27-8FE9-24C1BA0746DB}
- 13. ELECTROCARDIOGRAM TRACE 2 (ECG TR 2)
{73EACBC1-C1FB-4B2C-96AF-C5C1CA07DA72}

POSITION: 4 entries

- 1. LYING (Ly) {DE28D2BF-6F8D-64C3-0092-4CEB5E0902B9}
- 2. SEMIFOWLERS (Sfl) {C2F6A945-3829-AABA-3D67-54D1829D3EB9}
- 3. SITTING (Sit) {DD4DCDF1-0A22-2CC5-404F-0E472C73961C}
- 4. STANDING (St) {0776AE42-C370-B7B3-9B4D-9C914BA8F231}

QUALITY: 2 entries

- 1. ACTUAL (A) {5842A4D9-56A1-7A5A-13C9-F481699C3A09}
- 2. ESTIMATED (E) {9EAF26B-7D87-DA66-F991-022D095FB5A0}

UNITS:

BEATS PER MINUTE 0 > 60 * 100 < 300 Dec Prec: 0

Abbreviation..... PULSE STR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Quality of pulse
Display Name..... Pulse Strength
Help Text..... Select one of the measurements for
pulse quality
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 6 entries

-
- 1. 3 - STRONG (3-s) {C55EAB72-7816-5440-84E9-6E7E190ED197}
 - 2. 2 - NORMAL (2-n) {25FAADFC-D280-78AC-4F01-CFCF34ECD3FF}
 - 3. 1 - WEAK (1 - WEAK) {A41570DC-2F35-DDF6-977B-7C262150DFA2}
 - 4. 0 - ABSENT (0-a) {75F14311-F113-963C-4F05-B8AE17F044A9}
 - 5. D - DOPPLER (d-d) {2B77F099-0520-474D-A4BC-22BF824C8364}
 - 6. AMPUTEE (AMPUTEE) {5804286F-2B23-91BD-FCB5-6AD18C918203}

LOCATION: 8 entries

-
- 1. R RADIAL (R RADIAL) {BA10E6E5-DDF4-345E-68AF-21D32E9D972D}
 - 2. L RADIAL (L RAD) {302F78AD-4067-8B6B-BD15-993EC8561AE6}
 - 3. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
 - 4. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
 - 5. R POST TIBIAL (R PO TIB) {88CDA53A-5153-4950-701A-3F4F2007CC6B}
 - 6. L POST TIBIAL (L PO TIB) {BED130EE-31CC-C7E9-6597-78685BCE91AC}
 - 7. R PEDAL (R PEDAL) {9532B614-F4E1-316C-1322-905EB58E58C4}
 - 8. L PEDAL (LEFT PEDAL) {D84E9953-8C17-3130-8F21-A30D2A99FDE4}

UNITS:

TERM: PUPIL REACTION
ID: {9FCBF7CF-39B8-9710-6F41-9B7DDDC3BDD2}

Page: 442

Abbreviation..... PREA
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Pupil Reaction
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 9 entries

1. FIXED (f) {697D8C5A-976F-F423-30EC-A04102CF85F0}
2. SLUGGISH (s) {3215A82E-8DCC-C7C0-8AE1-2D4C9D4A4EFB}
3. EYE CLOSED (ec) {140B7D2F-F9C2-0B4C-3CDC-7355EDDB5030}
4. DRESSING (d) {FFEF0714-3B5B-E2B4-D760-3C03D45E59EA}
5. REACTIVE (Reactive) {31B44ADB-0000-A5B0-5F56-59D4A08C05D2}
6. NONREACTIVE (nonreact) {AB6B7D92-B94C-F3AA-E32F-6035D9CF2EF9}
7. BRISK (brisk) {980AE687-0A78-458B-BA7B-5D646C38A5F5}
8. PINPOINT (pinpoint) {4E53DDAC-0BF4-66D5-D614-9763BBF1FC8E}
9. DILATED (Dilated) {BF98ADC4-C69E-B17A-9BDE-D426D7F4C395}

LOCATION: 2 entries

-
1. L EYE (L EYE) {D714A0CA-EAD6-EC35-DA54-FF027C23FF45}
 2. R EYE (R EYE) {258F0C1D-67BD-ECF8-3152-EB289DEF4681}

UNITS:

TERM: PUPIL SIZE Page: 443
 ID: {80D7541F-0EE5-D362-E781-F825A9848E8C}

Abbreviation..... PS
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... Numeric
 Description..... No Description
 Display Name..... Pupil Size
 Help Text..... No Help Text
 Multi Select Picklist.....
 Term Type..... OBSERVATION
 Vuid.....

LOCATION: 2 entries

-
1. L EYE (L EYE) {D714A0CA-EAD6-EC35-DA54-FF027C23FF45}
 2. R EYE (R EYE) {258F0C1D-67BD-ECF8-3152-EB289DEF4681}

UNITS:

MILLIMETERS 0 > 1 * 9 < 11 Dec Prec: 0

TERM: QRS Page: 444
 ID: {7E76309C-0486-40CC-B14F-0D9B042C1D83}

Abbreviation..... QRS
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... Numeric
 Description..... The QRS complex is a structure on the
 ECG that corresponds to the

depolarization of the ventricles.
Display Name..... QRS Complex
Help Text..... The QRS complex is a structure on the
ECG that corresponds to the
depolarization of the ventricles.
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

SECOND 0.02 > 0.04 * 0.12 < 0.24 Dec Prec: 2

TERM: QT Page: 445
ID: {5226FEC7-8AD3-AE64-4122-6F9D2925448A}

Abbreviation..... QT
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... The QT interval is measured from the
beginning of the QRS complex to the
end of the T wave.
Display Name..... QT Interval
Help Text..... The QT interval is measured from the
beginning of the QRS complex to the
end of the T wave.
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

SECOND 0.10 > 0.30 * 0.46 < 1.00 Dec Prec: 2

TERM: RECTAL TUBE PRESENT Page: 446
ID: {D05DB527-8EDD-4324-B7B5-989AC004486E}

Abbreviation..... REC TU PRE
Active..... Yes
Boolean Value False..... Not Present
Boolean Value True..... Present
Data Type..... Boolean
Description..... Rectal Tube Present
Display Name..... Rectal Tube Present
Help Text..... Rectal Tube Present
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

-
- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
 - 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: RESPIRATORY PATTERN Page: 447
 ID: {5E7469A0-F426-2342-DA03-7054F44D7614}

=====

Abbreviation..... RESP PATT
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... No Description
 Display Name..... Respiratory Pattern
 Help Text..... No Help Text
 Multi Select Picklist..... Yes
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 11 entries

-
- 1. REGULAR (REGULAR) {3D3F3D4C-ECD0-2056-A0D1-43CC2F2BA2C8}
 - 2. IRREGULAR (IRREGULAR) {6CCE27A2-F1EB-3EC8-B140-CDA241D159F2}
 - 3. DEEP (DEEP) {78AD569C-16D4-CFC3-6B6F-41FB9FFAB9D7}
 - 4. SHALLOW (SHALLOW) {CD205617-BF98-4177-971A-056AFECC823A}
 - 5. DYSPNEIC (DYSPENIC) {4C6B42D9-694F-5EB2-6396-7ED85D5B154E}
 - 6. ORTHOPNEIC (ORTHOPNIC) {FF373BD7-FD84-DF43-7A44-1B8E62917116}
 - 7. ABDOMINAL (ABDOMINAL) {29E279F4-5D2E-4644-A2FA-CF85BEDFD480}
 - 8. CHEYNE-STOKES (CHEYNE-STO)
 {58DD2B80-B522-499C-2C60-C48A654BAD84}
 - 9. KUSSMAUL (KUSSMAUL) {0C5ED066-F445-ABFE-F7B1-D158EC7C2976}
 - 10. AGONAL (a) {1621B411-33E3-3A6A-78A6-40B7BDC3E8CB}
 - 11. VENTILATOR (v) {A5A647ED-A56D-5012-AD5B-788B3DCB3498}

UNITS:

=====

TERM: RESPIRATORY RATE Page: 448
 ID: {973ED2C0-0625-7DF9-17DC-8FFF7E376F23}

=====

Abbreviation..... RR
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... Numeric
 Description..... RESPIRATIONS PER MINUTE
 Display Name..... Respiratory Rate
 Help Text..... Breaths per minute
 Multi Select Picklist.....

Term Type..... OBSERVATION
Vuid..... 4688725

METHOD: 8 entries

- 1. AFTER EXERCISE (AE) {D129EF4B-3A99-9966-54F1-2C01A36814EE}
- 2. ASSISTED VENTILATOR (AV) {A540F333-5A09-5F21-7B51-2B5DF05E32C0}
- 3. AT REST (AR) {E8C289A6-0C3D-A7B2-BB85-7F1A23F22B50}
- 4. CONTROLLED VENTILATOR (CV) {0CFC2183-5245-E7AB-96ED-93AF4A3FEC64}
- 5. MONITOR (Mo) {63FC24EA-50B8-B302-D948-ABE3517CF6FC}
- 6. SPONTANEOUS (S) {5E9463BC-C871-54FD-F800-DCFDEDB8A93B}
- 7. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}
- 8. WITH ACTIVITY (Act) {544A262B-EA34-F82A-46A9-F80740515D65}

POSITION: 6 entries

- 1. LYING (Ly) {DE28D2BF-6F8D-64C3-0092-4CEB5E0902B9}
- 2. SITTING (Sit) {DD4DCDF1-0A22-2CC5-404F-0E472C73961C}
- 3. STANDING (St) {0776AE42-C370-B7B3-9B4D-9C914BA8F231}
- 4. SEMIFOWLERS (Sfl) {C2F6A945-3829-AABA-3D67-54D1829D3EB9}
- 5. SUPINE (SUPINE) {DEA8A744-B445-4A9E-875D-485E060F2ACC}
- 6. TRENDELENBURG (Tr) {42F3353B-7041-EEE1-A87E-941D3BE1CBCD}

UNITS:

BREATHS PER MINUTE -1 > 12 * 20 < 201 Dec Prec: 0

=====

TERM: RESTRAINT CRITERIA Page: 449
ID: {4B4C42F3-0F9E-40BE-90AB-C744E91228B0}

=====

Abbreviation..... RES CRITER
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Restraint Criteria
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 7 entries

- 1. RES CRIT: LINES (R.C. LINE) {27F545F2-5338-4BC8-B5CE-1D321D28F8D6}
- 2. RES CRIT: EQUIPMENT (R.C. EQUI) {38C2C56C-6D57-4750-8E68-F1193B16190A}
- 3. RES CRIT: AMBULATE (R.C. AMBU) {54B72608-2C13-46CD-8034-00179669762C}
- 4. RES CRIT: UNABLE TO FOLLOW DIRECTIONS (R.C. DIRE) {41358EAE-9912-48F8-8931-D3C635406D21}
- 5. RES CRIT: COMBATIVE (R.C. COMB)

- {38DBF61A-D364-4215-B4AA-36B660F96462}
- 6. RES CRIT: AGGRESSIVE (R.C. AGGR)
{7A47B361-85EA-4C88-A1AF-CD8FF7D35D97}
- 7. RES CRIT: HARM TO SELF (R.C. HARM)
{4DC33341-8F53-4D30-BC28-1E32DB54F30D}

UNITS:

=====

TERM: RESTRAINT EDUCATION Page: 450
 ID: {08B0A919-200B-4874-993B-E9924E9A439F}

=====

Abbreviation..... RES EDUCAT
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... No Description
 Display Name..... Restraint Education
 Help Text..... No Help Text
 Multi Select Picklist..... Yes
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 2 entries

- 1. FAMILY INFORMED (FAM INFORM)
{77135D32-9CDD-4D3F-9B62-6E8189B77844}
- 2. PATIENT INFORMED (PAT INFORM)
{074A1D19-3920-49F1-8958-2EB6F9CF05BA}

UNITS:

=====

TERM: RESTRAINT EDUCATION TYPE Page: 451
 ID: {05693D33-977E-4089-988F-4FE029D64BE4}

=====

Abbreviation..... RES EDU
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... No Description
 Display Name..... Restraint Education Type
 Help Text..... No Help Text
 Multi Select Picklist..... Yes
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 5 entries

- 1. RESTRAINTS EXPLAINED (RES EXPL)

- {1C063DD2-C4A3-4027-A245-6DD51DF5A3E9}
- 2. RESTRAINTS ED TYPE (RES ED TY)
{5560A332-C365-4DAF-BA33-610C4AE301F5}
- 3. RESTRAINTS ED REMOVAL (RES ED REM)
{B19316C8-687B-48FB-9EBC-9857662B9F85}
- 4. RESTRAINTS BROCHURE (RES BRO P)
{1E54C2A1-8B8D-4F0D-B247-91E19722CDED}
- 5. RESTRAINTS TIME LIMIT (RES TIME L)
{513B483E-AD7B-4B1A-BD01-91133C605C3C}

UNITS:

=====

TERM: RESTRAINT EVAL - ALT. MEASURES ATTEMPTED
ID: {B3B065C3-2612-45C0-9127-3CB42BDA99B6}

=====

Page: 452

Abbreviation..... R.E. ALT
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Alternative Measures Attempted
Help Text..... Restraint Evaluation: Alternative
Measures Attempted
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 11 entries

- 1. RES EVAL: DECREASE STIMULATION (R.E. DECR)
{0841626E-994E-4346-B701-1662A00209EB}
- 2. RES EVAL: REALITY ORIENTATION (R.E. REALI)
{221C363C-2261-40FC-8400-8B1D2E693F69}
- 3. RES EVAL: RELOCATE (R.E. RELO)
{D8165678-ECEF-4F4D-A1C4-C616349E1873}
- 4. RES EVAL: INVOLVEMENT WITH FAMILY (R.E. INVO)
{0EF31047-70A5-4072-953A-B62E1ACDF84B}
- 5. RES EVAL: AMBULATION/EXERCISE (R.E. EXER)
{46F4D1EB-3439-4BE4-A78C-453BDBF6C230}
- 6. RES EVAL: EDUCATION (R.E. EDUC)
{76341E6E-7F0A-429C-B768-9CA83C452FD7}
- 7. RES EVAL: DIVERSIONAL ACTIVITIES (R.E. DIVER)
{569F6D49-C564-4910-835A-B05179086E1F}
- 8. RES EVAL: MEDICATION REVIEW (R.E. REVI)
{914A0F82-5863-4C84-B3AB-D5A22BD340D8}
- 9. RES EVAL: UNABLE TO ATTEMPT (R.E. UNAB)
{979B0477-5789-4319-9562-6752048C4005}
- 10. BED EXIT ALARM (BED EXIT A)
{AC53759E-FFAE-5E0F-2EBA-E56A2712357C}
- 11. WHEELCHAIR EXIT ALARM (WHEEL EXIT)
{7D06564E-ECDE-4A11-BEFE-92AB7E448AA3}

1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: RESTRAINT EVAL - REPOSITION/RANGE OF MOTION Page: 459
ID: {54C82D71-B32F-48A2-87A9-683D6023D905}

=====

Abbreviation..... R.E. REPOS
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... No Description
Display Name..... Reposition / Range of Motion Evaluated
Help Text..... Restraint Evaluation: Reposition /
Range of Motion Evaluated
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: RESTRAINT EVAL - SKIN INTEGRITY Page: 460
ID: {1EE22D46-855D-4705-858B-D7623CC79F83}

=====

Abbreviation..... R.E. SKIN
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... No Description
Display Name..... Skin Integrity Evaluated
Help Text..... Restraint Evaluation: Assess Skin
Integrity Evaluated
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

=====

TERM: RESTRAINT LEAST RESTRICTIVE MEASURES
ID: {1132E26A-19E0-40D0-9602-68AB41F415FA}

=====

Page: 463

Abbreviation..... RES LRME
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Least Restrictive Measures Evaluated
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 10 entries

1. RES LRM: LIMIT SETTING (RLRM SETT)
{96417DDF-8562-4EFF-B4F6-0933FD91CD36}
2. RES LRM: REDIRECT (RLRM REDI)
{07DE67AF-E226-4794-ADF5-17E9C56C9FDF}
3. RES LRM: REMOVE FROM SETTING (RLRM REMO)
{857DBF6E-7EFF-4F07-9AC1-7FDB37CC7618}
4. RES LRM: PRN MEDS (RLRM: MEDS)
{CDF02FA3-FC72-4CCC-8197-41710A69A2CB}
5. RES LRM: VERBAL DE-ESCALATION (RLRM VERB)
{CE8618C8-CAD7-4C74-8566-ACB5E66B082C}
6. RES LRM: SUPERVISION (RLRM SUPR)
{18FA8CB0-0155-4D1C-9B14-AE51A0C06280}
7. RES LRM: VERBAL REMINDERS (RLRM REMI)
{987B5A32-DE2E-47E5-9B06-1AEE86096AB2}
8. RES LRM: DEVICE ACCESS (RLRM ACCE)
{1662F5E1-30E7-4E8A-BCC5-1738868435B8}
9. RES LRM: REORIENT (RLRM REOR)
{DFCC8245-180D-4F3F-AEFC-5D55AFDDC46D}
10. REPOSITION (REPOSITION) {DF3039CA-5977-4454-81D4-DB0AF966F85D}

UNITS:

=====

TERM: RESTRAINT MENTAL ASSESSMENT
ID: {78C70CF3-3253-46B1-AD61-4A6E7C72442E}

=====

Page: 464

Abbreviation..... RES MENT A
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Mental Assessment Patient Evaluation
Display Name..... Mental Assessment

Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

- 1. ORIENTED (ORIENTED) {B8DC9245-6D81-40E4-8B47-B6413E49F152}
- 2. DISORIENTED (DISORIENTE) {A7D1F2EA-CD7C-4100-8552-44834684C6DE}
- 3. RESTLESS (r) {8FCCAD26-DBA0-7067-DA72-30D342BFA8E8}
- 4. SLEEPING (s) {87338420-F351-F98C-7A80-4990B87ACE2B}
- 5. SEDATED (s) {29D9ADD6-2C1F-F6D5-F4D0-B2BB86416997}

UNITS:

=====

TERM: RESTRAINT ORDER TYPE Page: 465
ID: {843976A3-7C4A-4E09-A273-4589923DDDC2}

=====

Abbreviation..... RES ORD TY
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Order Type
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. SECLUSION (SECLUSION) {65B57F6C-C849-4BF7-8705-175361BC900F}
- 2. RESTRAINT (RESTRAINT) {0429322D-A987-4913-B64D-A6414645B214}

UNITS:

=====

TERM: RESTRAINT PATIENT ED. UNDERSTANDING Page: 466
ID: {5002371C-E683-459E-A4B0-3CC914081E78}

=====

Abbreviation..... RES PT U
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Patient Understanding Of Restraints
Education
Help Text..... No Help Text
Multi Select Picklist..... No

Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

- 1. INDICATES UNDERSTANDING (UNDERSTAND)
{C4AEFE5F-C388-4E00-9663-6B308934284B}
- 2. NEEDS REINFORCEMENT (NEED REINF)
{B6171F65-4530-418D-A63B-65DFFF522786}
- 3. UNSUCCESSFUL (UNSUCCESSF) {54DC58A1-7585-4A32-B810-90B586ED9332}

UNITS:

=====

TERM: RESTRAINT SAFETY CARE Page: 467
ID: {1450E52B-CBF3-BA02-F4E0-4AE793AE9614}

=====

Abbreviation..... RESSC
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Restraint Safety Care
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. CIRCULATION CHECK (ccq) {BA32C2A8-F15C-84EF-4859-6277EFA6463A}
- 2. DEVICE RELEASE (drq) {9F97947F-CFFE-3285-F4C3-A8DF294DA46C}

UNITS:

=====

TERM: RESTRAINT TYPE ALL Page: 468
ID: {696245B7-717E-4CF4-3C92-7E0F4DFDEDDE}

=====

Abbreviation..... RES TYPE A
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Restraint Type
Help Text..... Location and type of restraint used
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE:

19 entries

-
1. RESTRAINT SOFT (RESTR SOFT)
{C02CAC4B-5E9C-4060-BECD-0C0820EDA39B}
 2. RESTRAINT SOFT 2 POINT (RES S2POIN)
{66B71CFE-2603-40C9-9625-24448012BAB5}
 3. RESTRAINT SOFT 4 POINT (RES S4POIN)
{6820E997-60A0-4802-BE26-CA29B4F48F32}
 4. RESTRAINT SOFT ARMS + WAIST (RES S A+W)
{4A322BB6-0F63-4BB7-9290-C0655BD7D772}
 5. RESTRAINT SOFT ARMS + LOCKED WAIST (RES S A+LW)
{C480D67C-958C-42F9-9203-98E2F2BDBEC8}
 6. RESTRAINT SINGLE HAND MITT (RES S HA M)
{A5F1D037-FB5E-402B-9F25-22DE48404D57}
 7. RESTRAINT BILATERAL HAND MITT (RES BI H M)
{5E293110-1235-4307-80C2-0FFFA007E1DA}
 8. RESTRAINT VEST (RES VEST) {643378B2-86D7-4DB7-93ED-9DCE16D85E27}
 9. RESTRAINT BILATERAL ANKLE (RES BI ANK)
{8A0B0DE5-E652-40E7-9B67-2F4D27685A68}
 10. RESTRAINT WRIST (RES WRIST)
{03038E52-A66F-4EC0-8DC1-7B0BB596FFC0}
 11. RESTRAINT GERI CHAIR (RES GERI C)
{6B8F87BB-C2AB-4D92-AE04-24F19B26D2F1}
 12. RESTRAINT HARD PLASTIC/LEATHER 2 POINT (RES P/L 2P)
{17CB8739-612A-4319-946A-2062C5245C90}
 13. RESTRAINT HARD PLASTIC/LEATHER 4 POINT (RES P/L 4P)
{94DFD0BA-E941-4B09-8FD7-87E9AD7C7E6F}
 14. RESTRAINT ONE ARM (RES ONE AR)
{839914DA-6D63-43EB-961E-70DA56737A92}
 15. RESTRAINT ONE LEG (RES ONE LE)
{2790F0FC-B49A-4C99-9968-BCADD7CC63AC}
 16. RESTRAINT WAIST LOCKED (RES WAI LO)
{1EE7AB42-CE17-4D49-849B-CA550966E58C}
 17. RESTRAINT WAIST (RES WAIST)
{253E2757-5A57-4410-9C0E-F4A907A1EC16}
 18. RESTRAINT HARD PLASTIC/LEATHER (RES H PL/L)
{9BB528F4-7016-4030-B73F-4037ADF2621D}
 19. RESTRAINT STRAP (RES STRAP)
{642FDACE-B320-42BC-9514-F5FCD6583B19}

LOCATION:

11 entries

-
1. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
 2. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
 3. L WRIST (L WRIST) {7BCAFACE-FF3B-46DA-99F9-6437A90B7126}
 4. BILATERAL WRIST (BILAT WRIS)
{89B4D30F-9E46-F8D6-14F6-1CC0843DDC21}
 5. BILATERAL LEG (BILAT LEG) {0713761D-E0FC-DA03-02A1-F20EFEB99F45}
 6. BILATERAL ANKLES (BILA ANKLE)
{E56A03A5-5CE0-47BF-A14D-F6DF66F4E769}
 7. CHEST (CHEST) {8C3C30FE-6E06-C35E-8226-CBF0C94014CC}
 8. L LEG (LL) {B088C809-E528-B785-21C2-14E45E1FD8AF}
 9. R LEG (RL) {C5E40A7B-6FC4-267D-8CE8-3571D21D81FD}
 10. R WRIST (R WRIST) {FDD9FE1D-F8A2-09AD-82D7-AA478693DD01}
 11. WAIST (WAIST) {784DD311-8B48-4FE3-96B6-471CCE453CD5}

UNITS:

=====

TERM: RESTRAINT TYPE ARMS

Page: 469

ID: {0A2CB443-AA1A-4526-9BF0-9C63B49D7603}

=====

Abbreviation..... RES TY ARM
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Arm Restraints
Display Name..... Arm Restraints
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

1. RESTRAINT SINGLE HAND MITT (RES S HA M)
{A5F1D037-FB5E-402B-9F25-22DE48404D57}
2. RESTRAINT HARD PLASTIC/LEATHER (RES H PL/L)
{9BB528F4-7016-4030-B73F-4037ADF2621D}
3. RESTRAINT SOFT (RESTR SOFT)
{C02CAC4B-5E9C-4060-BECD-0C0820EDA39B}

LOCATION: 6 entries

1. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
2. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
3. BOTH ARMS (BOTH ARMS) {A493AB32-9DE9-40AC-8C46-7EFF395B63B5}
4. L WRIST (L WRIST) {7BCAFACE-FF3B-46DA-99F9-6437A90B7126}
5. R WRIST (R WRIST) {FDD9FE1D-F8A2-09AD-82D7-AA478693DD01}
6. BILATERAL WRIST (BILAT WRIS)
{89B4D30F-9E46-F8D6-14F6-1CC0843DDC21}

UNITS:

=====

TERM: RESTRAINT TYPE LEGS

Page: 470

ID: {C3B1D18A-FF73-4362-8FCF-7AFB5304A4CF}

=====

Abbreviation..... RES TY LEG
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Leg Restraints
Display Name..... Leg Restraints
Help Text..... No Help Text
Multi Select Picklist..... No

Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

- 1. RESTRAINT HARD PLASTIC/LEATHER (RES H PL/L)
{9BB528F4-7016-4030-B73F-4037ADF2621D}
- 2. RESTRAINT SOFT (RESTR SOFT)
{C02CAC4B-5E9C-4060-BECD-0C0820EDA39B}
- 3. RESTRAINT STRAP (RES STRAP)
{642FDACE-B320-42BC-9514-F5FCD6583B19}

LOCATION: 6 entries

- 1. L LEG (LL) {B088C809-E528-B785-21C2-14E45E1FD8AF}
- 2. R LEG (RL) {C5E40A7B-6FC4-267D-8CE8-3571D21D81FD}
- 3. BOTH LEGS (BOTH LEGS) {5A15A231-4695-4CF3-AD0F-2E7950208967}
- 4. L ANKLE (L ANKLE) {D5A692EF-9D7C-473C-AF96-7EF1C01FA2AF}
- 5. R ANKLE (R ANKLE) {4E6DF352-F399-4CE9-908C-48132C3D856B}
- 6. BILATERAL ANKLES (BILA ANKLE)
{E56A03A5-5CE0-47BF-A14D-F6DF66F4E769}

UNITS:

TERM: RESTRAINT TYPE OTHER

Page: 471

ID: {40448F75-66C2-4FA1-AA02-6B827C9A06F8}

Abbreviation..... RES TYP OT
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Restraint Type Other
Display Name..... Restraint Type Other
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

- 1. GERI-CHAIR (g) {167F9786-8F24-054F-020A-658A6EE5B0B4}
- 2. RESTRAINT VEST (RES VEST) {643378B2-86D7-4DB7-93ED-9DCE16D85E27}
- 3. RESTRAINT WAIST (RES WAIST)
{253E2757-5A57-4410-9C0E-F4A907A1EC16}

LOCATION: 3 entries

- 1. TORSO (TORSO) {EA1AD477-8BFF-E575-23C1-B9CFAF2B9986}
- 2. CHEST (CHEST) {8C3C30FE-6E06-C35E-8226-CBF0C94014CC}
- 3. WAIST (WAIST) {784DD311-8B48-4FE3-96B6-471CCE453CD5}

UNITS:

=====

TERM: RIGHT VENTRICULAR STROKE WORK
ID: {3A72DA75-5112-7980-4F4F-E1E34FFCD1A6}

=====

Page: 472

Abbreviation..... RVSW
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Right Ventricular Stroke Work
Display Name..... Right Ventricular Stroke Work
Help Text..... Right Ventricular Stroke Work
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

1. CARDIAC OUTPUT MONITOR (COM)
 {CCBE5744-7CB9-9DFA-938F-FDDBC73B45B3}

UNITS:

STROKE WORK PER BEAT 0.0 > 5.0 * 30.0 < 50.0 Dec Prec: 1

=====

TERM: RIGHT VENTRICULAR STROKE WORK INDEX
ID: {EE450A40-E87C-1F5C-FDEF-BAAD030B8149}

=====

Page: 473

Abbreviation..... RVSWI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Right Ventricular Stroke Work Index
Display Name..... Right Ventricular Stroke Work Index
Help Text..... Right Ventricular Stroke Work Index
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

1. CARDIAC OUTPUT MONITOR (COM)
 {CCBE5744-7CB9-9DFA-938F-FDDBC73B45B3}

UNITS:

STROKE WORK PER BEAT PER M2 1.9 > 2.0 * 18.0 < 18.1 Dec Prec: 1
GM-M/M2/BEAT 3.9 > 4.0 * 18.0 < 18.1 Dec Prec: 1

=====

TERM: SAFETY DEVICE WAIST BELT APPLIED
ID: {27E75399-989F-4598-A972-646FBFF51A20}

=====

Page: 474

Abbreviation..... WAIST BE
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Safety Device Waist Belt Applied
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

-
- 1. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}
 - 2. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}

UNITS:

=====

TERM: SEDATION SCALE - RAMSEY
ID: {CF2E7E29-8C96-461B-9074-CA94455F62D6}

=====

Page: 475

Abbreviation..... SED RAMSEY
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Ramsey Sedation Scale
Display Name..... RAMSEY
Help Text..... Ramsey Sedation Scale
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 6 entries

-
- 1. RAMSEY ANXIOUS (RAM ANXIOU)
{88205A0C-7065-4616-8FCF-75059462DD15}
 - 2. RAMSEY COOPERATIVE (RAM COOPER)
{56846ADA-0F0E-44AA-9779-A567B9D24C4F}
 - 3. RAMSEY SEDATED (RAM SEDATE)
{CC12405F-60DB-485F-94E5-78C1EB32AB74}
 - 4. RAMSEY ASLEEP BRISK RESPONSE (RAM ASL BR)
{00F4A117-F9E7-4B60-A521-F2E8725AECD3}
 - 5. RAMSEY ASLEEP SLUGGISH (RAM ASL SL)
{AA61604F-00F9-4103-B27F-E817FBC9577B}
 - 6. RAMSEY ASLEEP NO RESPONSE (RAM ASL NO)
{3FF5BA0A-2D3C-4935-A432-23DAA4761AB7}

UNITS:

=====

TERM: SEDATION SCALE - RASS
ID: {F4900650-2584-4A5B-9F97-171A9552C7BA}

=====

Page: 476

Abbreviation..... RASS
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Richmond Agitation Sedation Scale
Display Name..... RASS
Help Text..... Richmond Agitation Sedation Scale
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 10 entries

-
1. RASS COMBATIVE (RASS COMBA)
{A97E064A-7709-42F9-9A88-4AAFF1020883}
 2. RASS VERY AGITATED (RASS V AGI)
{7528CA1D-A9A4-4F4E-8943-E69F3B42E189}
 3. RASS AGITATED (RASS AGITA)
{D81A5D5D-E7FC-445A-A814-495731ABA0F8}
 4. RASS RESTLESS (RASS RESTL)
{2CCFC1ED-D807-45F0-AE69-980E214AD5C3}
 5. RASS ALERT AND CALM (RASS AL CA)
{F1834BF6-F8F1-4E13-BC25-B13E756E0FF0}
 6. RASS DROWSY (RASS DROWS) {5E6BDAC6-4F81-4110-A0AF-FED746D9FF0D}
 7. RASS LIGHT SEDATION (RASS LT SE)
{FDEEC084-3922-4065-A863-4893E42F8F93}
 8. RASS MODERATE SEDATION (RASS MD SE)
{D12B0A3E-B80E-494F-B89F-A79E8B81B080}
 9. RASS DEEP SEDATION (RASS DP SE)
{83CB4E64-3581-452F-987A-D31CB3AB38D2}
 10. RASS UNAROUSABLE (RASS UNARO)
{1578E147-7D2C-4469-828E-FCA97109B27C}

UNITS:

=====

TERM: SEDATION SCALE - SAS
ID: {8281E806-CEC9-4236-B984-495B8A3973E5}

=====

Page: 477

Abbreviation..... SED SC SAS
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Riker Sedation-Agitation Scale

Display Name..... SAS
Help Text..... Riker Sedation-Agitation Scale
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 7 entries

- 1. SAS DANGEROUS AGITATION (SAS DA AGI)
{0701FB25-1DF0-4771-8C81-AC19A985A3B6}
- 2. SAS VERY AGITATED (SAS VY AGI)
{B501606A-EEA7-470E-81BE-F51D422272F1}
- 3. SAS AGITATED (SAS AGITAT) {107143D4-9527-4A26-8CB7-EEBB656B68B3}
- 4. SAS CALM AND COOPERATIVE (SAS C COOP)
{F0633778-280C-43C3-9E6E-9C35C18179F6}
- 5. SAS SEDATED (SAS SEDATE) {AFEC3F29-52D3-4A49-91F8-05AA3891FEC8}
- 6. SAS VERY SEDATED (SAS VY SED)
{4B13529B-61C5-4B68-857A-3A5351DC4BE2}
- 7. SAS UNAROUSABLE (SAS UNAROU)
{3DDFF9FE-A88A-471F-BF36-E5B569DEBDAA}

UNITS:

TERM: SEIZURE CHARACTER Page: 478
ID: {15B23878-6E4C-4FCE-9F49-1C728A4A65D2}

Abbreviation..... SEIZU CHAR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Seizure Character
Display Name..... Seizure Character
Help Text..... Character of Seizure
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

- 1. TONIC/CLONIC (TON/CLON) {81932BAA-6C79-47A3-85B3-592824B6C09C}
- 2. PARTIAL IMPAIRED (PART IMPAI)
{2F1B8C52-1C94-4EA7-AB82-117694FD4D52}
- 3. PARTIAL SIMPLE (PART SIMP)
{F4A5E7ED-92AB-4021-A63C-E1D948EFE6C1}

UNITS:

TERM: SEIZURE DURATION Page: 479
ID: {DAABDF54-F64C-6306-87EA-B8F5DE034DA0}

Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: SEIZURES Page: 482
ID: {BB006255-D69A-E29C-1286-D71E6FF15830}

=====

Abbreviation..... SEI
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... No Description
Display Name..... Seizures
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: SENSATION FOOT Page: 483
ID: {CCF961DC-8590-2498-3836-0F010461F69A}

=====

Abbreviation..... SENSF
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Sensation Foot
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 4 entries

- 1. INTACT (i) {935C9CAF-6C88-6610-C6E9-4C1B822923AF}
- 2. NUMBNESS (n) {037F2F79-36D9-A0B3-C23E-8EE525932AF3}
- 3. TINGLING (t) {AE8E3101-AC3D-DEC3-0293-58F901F82EB8}
- 4. ABSENT (a) {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0}

LOCATION: 2 entries

- 1. L FOOT (L FOOT) {2EC3FBEC-4EA2-9E7A-17D4-FDC2F9D28046}
- 2. R FOOT (R FOOT) {7D9F5FF6-232C-2488-4504-A721FF06928F}

UNITS:

TERM: SENSATION HAND Page: 484
ID: {7B880F9C-83F8-3395-9983-F10A1EF2F97A}

Abbreviation..... SENSH
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Sensation Hand
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

- 1. INTACT (i) {935C9CAF-6C88-6610-C6E9-4C1B822923AF}
- 2. NUMBNESS (n) {037F2F79-36D9-A0B3-C23E-8EE525932AF3}
- 3. TINGLING (t) {AE8E3101-AC3D-DEC3-0293-58F901F82EB8}
- 4. ABSENT (a) {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0}

LOCATION: 2 entries

- 1. L HAND (L HAND) {639DF25A-19E3-B504-D323-4B3F6ED39A24}
- 2. R HAND (R HAND) {A54D3A80-8537-D079-9CCE-7DC56346D132}

UNITS:

TERM: SENSATION LOWER EXTREMITY Page: 485
ID: {EB69E23D-9C3B-F2AB-B7F6-CAE402C87DE9}

Abbreviation..... SLEX
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description

Display Name..... Sensation Lower Extremity
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

- 1. INTACT (i) {935C9CAF-6C88-6610-C6E9-4C1B822923AF}
- 2. NUMBNESS (n) {037F2F79-36D9-A0B3-C23E-8EE525932AF3}
- 3. TINGLING (t) {AE8E3101-AC3D-DEC3-0293-58F901F82EB8}
- 4. ABSENT (a) {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0}

LOCATION: 2 entries

- 1. L LOWER EXTREMITY (LLEX) {FC6F5056-779A-BF76-86D0-A237936C1ECC}
- 2. R LOWER EXTREMITY (R LOW EXTR)
{7FC3CECC-1E6E-A686-FF04-2BEF8A9B9DCA}

UNITS:

TERM: SENSATION UPPER EXTREMITY
ID: {25313956-8B3E-FA21-0AA7-8E9DD4ABFF8D}

Page: 486

Abbreviation..... SUEX
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Sensation Upper Extremity
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

- 1. INTACT (i) {935C9CAF-6C88-6610-C6E9-4C1B822923AF}
- 2. NUMBNESS (n) {037F2F79-36D9-A0B3-C23E-8EE525932AF3}
- 3. TINGLING (t) {AE8E3101-AC3D-DEC3-0293-58F901F82EB8}
- 4. ABSENT (a) {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0}

LOCATION: 2 entries

- 1. L UPPER EXTREMITY (LUEX) {48806E79-B1DA-129B-EB52-453292A3DDC6}
- 2. R UPPER EXTREMITY (R UPP EXTR)
{C65BA269-8F05-6FB9-62CB-466A286CF2C0}

UNITS:

UNITS:

=====

TERM: SEWS - ANXIETY Page: 489
ID: {1F7184F3-45F0-4D93-808C-FB797589814F}

=====

Abbreviation..... SEWS ANXIE
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Severity of Ethanol Withdrawl Scale
Display Name..... Anxiety
Help Text..... Do you feel that something bad is
about to happen to you right now?
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

1. 0 (0)	{2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
2. 3 (3)	{E714454D-6A67-91F6-2FB7-D8F24F184328}

UNITS:

=====

TERM: SEWS - DISORIENTATION NAME PLACE Page: 490
ID: {A6E459CA-ADEF-4EA5-8B24-11300836DE9A}

=====

Abbreviation..... SEWS N PL
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Severity Of Ethanol Withdrawl Scale
Display Name..... Disorientation: Knows Name and Place
Help Text..... Disorientation: Knows name and place
but not date
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

1. 0 (0)	{2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
2. 1 (1)	{8F5D10A8-FE17-47BC-DACA-6972D12CC00A}

UNITS:

=====
Abbreviation..... SEWS N ONL
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Severity Of Ethanol Withdrawl Scale
Display Name..... Disorientation: Knows Name Only
Help Text..... Disorientation: Knows Name Only
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

1. 0 (0) {2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
2. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}

UNITS:

=====
TERM: SEWS - NAUSEA Page: 494
ID: {168E97DD-48FA-45BF-A3C0-8FCB5B316B6D}

=====
Abbreviation..... SEWS NAUSE
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Severity Of Ethanol Withdrawl Scale
Display Name..... Nausea and Dry Heaves or Vomitting?
Help Text..... Nausea and Dry Heaves or Vomitting?
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

1. 0 (0) {2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
2. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}

UNITS:

=====
TERM: SEWS - SWEATING Page: 495
ID: {4B7D9035-79F8-4736-9A20-87671D76E12B}

=====
Abbreviation..... SEWS SWEAT
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList

Description..... Severity Of Ethanol Withdrawl Scale
Display Name..... Sweating
Help Text..... Sweating: Includes moist palms,
sweating now?
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. 0 (0) {2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
2. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}

UNITS:

=====

TERM: SEWS - TOTAL Page: 496
ID: {6220504D-6961-4138-B1EC-1257142F5E49}

=====

Abbreviation..... SEWS TOTAL
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Severity Of Ethanol Withdrawl Scale
Total
Display Name..... SEWS TOTAL
Help Text..... 5 Or Less = Lower Risk For Withdrawl.
6 Or More = Higher Risk For Withdrawl
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

SEWS SCALE -1 > 0 * 5 < 24 Dec Prec: 0

=====

TERM: SEWS - TREMOR Page: 497
ID: {8BC41692-A9A8-4900-BB2D-5E60E5625798}

=====

Abbreviation..... SEWS TREMO
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Severity Of Ethanol Withdrawl Scale
Display Name..... Tremor
Help Text..... Tremor: with arms extended, eyes
closed?
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

1. 0 (0) {2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
2. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}

UNITS:

TERM: SEWS - VITAL SIGNS Page: 498
ID: {657A9E0B-12CC-4D33-B4A8-99827CAE5CFF}

Abbreviation..... SEWS VITAL
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Severity Of Ethanol Withdrawl Scale
Display Name..... Vitals: Pulse > 110, Diastolic BP >
90, Temp > 99.6
Help Text..... Vitals: Pulse > 110, Diastolic BP >
90, Temp > 99.6
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

1. 0 (0) {2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
2. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}

UNITS:

TERM: SIDE RAILS Page: 499
ID: {1BE3FF8D-BC95-CE9E-ABFB-9C94CF591F24}

Abbreviation..... SD RAIL
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Side Rails
Help Text..... Side Rail Positions
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Void.....

DATA VALUE: 14 entries

1. ALL UP (ALL UP) {48199083-ACC2-408C-9242-7A6313CD7A7F}

- 2. ALL DOWN (ALL DOWN) {76C14E79-95E2-49B7-B847-3B114653E98C}
- 3. HEAD UP (hu) {BB1C79C0-80F1-7CDE-CC10-7CA0C704AE35}
- 4. R HEAD UP (R HEAD UP) {D03DBEC6-F4DD-441A-AAB2-653CC34FCDBF}
- 5. L HEAD UP (L HEAD UP) {BB39792D-B572-4A10-BA81-8685809E7437}
- 6. HEAD DOWN (hd) {80BF3425-D252-F917-052E-DC54D39E464E}
- 7. R HEAD DOWN (R HEAD DOW) {DFFE0638-444C-4023-B586-985E84EEA1FA}
- 8. L HEAD DOWN (L HEAD DOW) {3C919CB8-1C0A-498D-95A6-C46D25744CEC}
- 9. FOOT UP (fu) {49CCC59E-002E-8387-E8D6-D57AAC089EF3}
- 10. R FOOT UP (R FOOT UP) {5928B710-6EB1-4183-93CB-6C72B37A97C3}
- 11. L FOOT UP (L FOOT UP) {A15CD8C5-5CD1-4B90-B79C-5D3D331C7466}
- 12. FOOT DOWN (fd) {13F73BE9-2FA8-FE88-443B-9BA3FE6EDFD7}
- 13. R FOOT DOWN (R FOOT DOW) {DBCC0BA2-3097-42C7-AE9E-DF71FA779162}
- 14. L FOOT DOWN (L FOOT DOW) {5B18895B-DAD5-4874-87DB-DACFFF88B170}

UNITS:

=====

TERM: SIGNIFICANT EVENT

Page: 500

ID: {0DBEF78C-8CFC-0173-2A39-02C24EB10A5F}

=====

Abbreviation..... se
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... String
 Description..... No Description
 Display Name..... Significant Event
 Help Text..... No Help Text
 Multi Select Picklist.....
 Term Type..... OBSERVATION
 Vuid.....

UNITS:

=====

TERM: SINUS DRAINAGE CHARACTER

Page: 501

ID: {05B77F59-77DE-424F-ACD1-413E70D52429}

=====

Abbreviation..... SIN DRA CH
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... No Description
 Display Name..... Sinus Drainage Character
 Help Text..... Describe The Sinus Drainage
 Multi Select Picklist..... Yes
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE:

4 entries

1. SEROUS (SEROUS) {9BDA347E-29E9-80C7-91FA-016192D4257E}
2. SEROSANGUINOUS (SEROSANGUI) {19A8D332-899D-177A-2247-ABCDAAEE3BFA8}
3. CLEAR (c) {8981C365-2B46-704F-C9F4-975DF1C23C2D}
4. CLOUDY (CLOUDY) {3CCFCCCC-00FC-6B00-1D58-A77C02FDB401}

UNITS:

=====

TERM: SINUS DRAINAGE PRESENT Page: 502
 ID: {0F9CDF0D-A7A8-408D-BBFD-A6F74E5B8D6F}

=====

Abbreviation..... SIN DRAI P
 Active..... Yes
 Boolean Value False..... No
 Boolean Value True..... Yes
 Data Type..... Boolean
 Description..... No Description
 Display Name..... Sinus Drainage Present
 Help Text..... Is Sinus Drainage Present?
 Multi Select Picklist.....
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 2 entries

-
1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
 2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: SINUS NECROTIC TISSUE APPEARANCE Page: 503
 ID: {0578E21E-CFBE-440A-A532-8135AD414085}

=====

Abbreviation..... SNEC TI AP
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... No Description
 Display Name..... Sinus Necrotic Tissue Appearance
 Help Text..... Sinus Necrotic Tissue Appearance
 Multi Select Picklist..... Yes
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 4 entries

-
1. SOFT (SOFT) {E68B4D4B-B86B-4004-A861-986645267F6E}
 2. YELLOW (YELLOW) {90DC31AD-6211-9F25-4237-79F7122A2763}
 3. HARD (HARD) {7934B543-7028-4CBC-9308-99F5A952A1DE}

4. DRY (d) {D8F4F0FD-1B19-2515-3815-EC88D175B305}

UNITS:

=====

TERM: SINUS NECROTIC TISSUE PRESENT

Page: 504

ID: {299921C9-5A13-4D9E-BC0D-8F777ED490D4}

=====

Abbreviation..... SIN NEC TI
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... No Description
Display Name..... Sinus Necrotic Tissue Present
Help Text..... Is Sinus Necrotic Tissue Present?
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: SINUS TRACT SURROUNDING TISSUE CHARACTER

Page: 505

ID: {59D298AC-1BE0-4C13-B75B-1CE331454151}

=====

Abbreviation..... SIN TI CHA
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Sinus Tract Surrounding Tissue
Character
Help Text..... Sinus Tract Surrounding Tissue
Character
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

1. ERYTHEMA (ERYTHEMA) {0083E1B3-AE2B-4CCF-8794-86E612B9DCEA}
2. DISCOLORATION (DISCOLORAT)
{BBE3E910-CC49-430C-BD90-FB0A97F1B183}
3. INDURATION (INDURAITON) {A1669576-45D9-435C-87B9-488E9CF6BBA4}
4. MACERATION (MACERATION) {B7C63E05-E193-40AA-A572-F1ACA0162B8B}

UNITS:

=====

TERM: SINUS TRACT ULCER

Page: 506

ID: {D89C6750-E1F0-4817-9EDD-205EE310A09A}

=====

Abbreviation..... SIN TRA UL
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... No Description
Display Name..... Sinus Tract Ulcer
Help Text..... Is a Sinus Tract Ulcer Present?
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

- | | |
|------------|--|
| 1. NO (n) | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEA649CD} |

UNITS:

=====

TERM: SINUS TRACT ULCER EPITHELIALIZATION

Page: 507

ID: {63375C18-43BF-4969-989F-8239440E4040}

=====

Abbreviation..... SIN ULC EP
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... No Description
Display Name..... Sinus Tract Ulcer Epithelialization
Help Text..... Is Sinus Tract Ulcer
Epithelialization Present?
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

- | | |
|------------|--|
| 1. NO (n) | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEA649CD} |

UNITS:

=====

TERM: SINUS TRACT ULCER GRANULIZATION
ID: {64BE9437-2EAB-40F2-93BC-2D0C9DBDC1EB}

=====

Page: 508

Abbreviation..... SIN ULC GR
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... No Description
Display Name..... Sinus Tract Ulcer Granulization
Help Text..... Is Sinus Tract Ulcer Granulization
Present?
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

-
- | | |
|------------|--|
| 1. NO (n) | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEA649CD} |

UNITS:

=====

TERM: SKIN ALTERATION
ID: {6B4325F4-F776-48C1-902D-2CFFC73A22BB}

=====

Page: 509

Abbreviation..... SKIN ALTER
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Skin Alteration
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

-
- | | |
|------------------|--|
| 1. NONE (NONE) | {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65} |
| 2. WOUND (WOUND) | {1EC40149-051E-4750-8651-2B2D0DD2DEAC} |

UNITS:

=====

TERM: SKIN COLOR
ID: {8D3A8036-39B5-AA2B-BF62-CC1EB134BB19}

=====

Page: 510

Abbreviation..... SKIN COLOR

Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Observation of Patients Skin Color
Display Name..... Skin Color
Help Text..... Observation of Patients Skin Color
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 13 entries

- 1. PALE (p) {0EB39D28-7ABE-A4AB-D334-3F795E65618D}
- 2. MOTTLED (m) {4B5FFB10-2CFF-2064-1553-4CBB7133B0C5}
- 3. ICTERIC (ICTERIC) {274ED5BB-5A34-20EE-A9B0-81AB7803780E}
- 4. JAUNDICE (j) {B17F4603-7C68-660D-96A7-119140A998D3}
- 5. FLUSHED (nt) {C7C54C1A-6E75-4BCF-4C57-0F8622A5BD16}
- 6. DUSKY (d) {2ABF845E-4E2F-883F-3FCC-F72B8507FF85}
- 7. CYANOTIC (c) {06BD77DB-1849-D731-0F27-F8CE48755EA9}
- 8. ASHEN (ASHEN) {7431E4EC-DCF1-4966-85CB-29BF65311855}
- 9. NORMAL FOR ETHNIC GROUP (ETH. NORM.)
{3B2CD951-58CE-FF53-AD5D-6C65C07658D3}
- 10. PINK (PINK) {1A7FAF84-D69C-4EEB-7BC5-D21A598A1F4D}
- 11. BRUISED (BRUISED) {FA5B38A7-7744-186F-055A-8458D696AF33}
- 12. HEMATOMA (HEMATOMA) {6A24B7EE-62F3-C952-C1A3-086CDBB1FE68}
- 13. ECCHYMOSIS (ECCHYMOSIS) {DF076906-CC17-6AD1-76E7-B07678A126DE}

UNITS:

=====

TERM: SKIN CONDITION Page: 511
ID: {831BEB82-92CE-07A0-F36B-09C94030D088}

=====

Abbreviation..... SKIN COND.
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Skin Condition
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. INTACT (i) {935C9CAF-6C88-6610-C6E9-4C1B822923AF}
- 2. BREAKDOWN (BREAKDOWN) {9CFD21C0-A805-47C4-E559-8033C754DE80}

UNITS:

=====

TERM: SKIN DRAINAGE AMOUNT
ID: {B973CF56-27AD-447B-9B12-A4326997E027}

=====

Page: 512

Abbreviation..... SKI DRA AM
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Skin Device Drainage Amount
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

-
1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
 2. SCANT (SCANT) {8EF6C62D-0D94-37BB-60AD-86E306E8120D}
 3. SMALL (s) {5991F892-B28C-3998-E733-D8FFA26D7E21}
 4. MODERATE (MODERATE) {C3C25064-C9FC-EA1F-3E74-32A25F0FC5E3}
 5. HEAVY (HEAVY) {1F10CD98-03FA-4705-BF7D-941E881C4E3D}

UNITS:

=====

TERM: SKIN ERYTHEMA EXTENDS AROUND INCISION
ID: {A88D0DBE-7CC6-4983-B9AC-AF9CCF9653AB}

=====

Page: 513

Abbreviation..... ER EXI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Skin Erythema Extension Around Incision
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

-
1. EXTENDS AROUND INCISION NONE (ER EXI NON)
{91F1A7C5-885C-4B13-9F7F-189C93E505CC}
 - 1) 2. EXTENDS AROUND INCISION EXTENDS INTO PERIWOUND LESS THAN 1 CM (ER EXI <
{095709D0-6190-49A3-AB08-BF4575DC8396}
 3. EXTENDS AROUND INCISION MORE THAN 1 CM (ER EXI >1)
{11F13A1F-7AC2-4629-A960-CD92655EE101}
 4. LIMITED TO INCISION (ER EXI LIM)
{C9ED5E2C-8ED1-445A-8A3E-C1207CE3843D}

UNITS:

=====
TERM: SKIN ERYTHEMA EXTENSION AROUND WOUND
ID: {6BABC7F1-2C62-41A7-83E2-F5F4B5ACE087}

Page: 514

Abbreviation..... ERY EX
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Skin Erythema Extension Around Wound
Display Name..... Erythema Extension Around Wound
Help Text..... Skin Erythema Extension Around Wound
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

- 1. EXTENSION AROUND WOUND NONE (ERY EX NON)
2. EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND MORE THAN 1 CM (ERY EX >1)
3. EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND LESS THAN 1 CM (ERY EX <1)
4. EXTENSION AROUND WOUND LIMITED TO WOUND (ERY EX LIM)

UNITS:

=====
TERM: SKIN INCISION CLOSURE
ID: {210023AD-813A-05F9-29E5-9F4339FD8BFE}

Page: 515

Abbreviation..... INCI CLOSU
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Skin Incision Closure
Help Text..... Type of material used to close incision
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 8 entries

- 1. SKIN INCISION WOUND CLOSURE SUTURES (INC CL SUT)
2. INCISION WOUND CLOSURE GLUE (IN CL GLU)

- {9EEB436E-2991-4E9A-9F8A-BFEC196784CB}
- 3. SKIN INCISION WOUND CLOSURE STAPLES (INC CL STA)
{4F48D42D-C206-423C-970B-02C03B2EE0B6}
- 4. SKIN INCISION WOUND CLOSURE STERI-STRIPS (INC CL STE)
{DF52B6F3-54EE-46B3-AFDF-4A701B065E04}
- 5. INCISION WOUND CLOSURE NOT CLOSED (INC CL)
{D63A05BF-1205-45F6-A78D-15D198BCB1C6}
- 6. SKIN INCISION WOUND CLOSURE OPEN WOUND (INC CL OPE)
{AD018B6F-13E5-41C0-BDBC-306CA0B1B475}
- 7. INCISION WOUND CLOSURE APPROXIMATED (INC CL APP)
{2DD5713E-D81C-484B-A0DF-6A5730B396DF}
- 8. INCISION WOUND CLOSURE DEHISCENCE (INC CL DEH)
{DCA1D137-0D9A-42D2-9270-7B3A82DCB6EB}

UNITS:

=====

TERM: SKIN INDURATION EXTENSION AROUND INCISION
ID: {36580CA5-4FB0-4562-9652-0B3BD7D04307}

=====

Page: 516

Abbreviation..... INC IN
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Induration Extension Around Incision
Display Name..... Induration Extension Around Incision
Help Text..... Induration Extension Around Incision
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

- 1. SKIN INDURATION AROUND INCISION NONE (INC IN NON)
{A9B95680-F9A0-44C0-929C-70C00B9AA5B7}
- 2. SKIN INDURATION AROUND INCISION EXTENDS INTO PERIWOUND LESS THAN 1 cm (INC IN <1)
{2A5097BF-955C-4F79-B483-7EE14AF05786}
- 3. SKIN INDURATION AROUND INCISION EXTENDS INTO PERIWOUND MORE THAN 1 cm (INC IN >1)
{D8EE448D-00EB-412D-95AC-2E2E92C53EEB}
- 4. SKIN INDURATION AROUND INCISION LIMITED TO INCISION (INC IN LIM)
{88BCE7F9-F21F-427A-A50B-C1BA831639AB}

UNITS:

=====

TERM: SKIN INDURATION EXTENSION AROUND WOUND
ID: {20391515-609B-4D18-AE79-533E47AFCAA8}

=====

Page: 517

Abbreviation..... WOU IN
Active..... Yes
Boolean Value False.....

Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Skin Induration Extension Around Wound
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 4 entries

-
1. SKIN INDURATION EXTENSION AROUND WOUND NONE (WOU IN NON)
{2C57CD27-4FDD-4DE8-88A0-D438CBE535D8}
 2. SKIN INDURATION EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND < 1 cm (WOU IN <1)
{3594CCCD-FA51-4971-A61E-FFA3B4AC7EE0}
 3. SKIN INDURATION EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND > 1 cm (WOU IN >1)
{7DF7D08D-3E5E-4D18-94E4-A0DFDA5DB72F}
 4. SKIN INDURATION EXTENSION AROUND WOUND LIMITED TO WOUND (WOU IN LIM)
{1828A9EA-FF9D-4A49-9543-15B85FBF211E}

UNITS:

=====

TERM: SKIN MOISTURE Page: 518
ID: {BE934424-A7CF-EEF0-8098-26B3AFC4B2FE}

=====

Abbreviation..... SKIN MOIST
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Skin Moisture
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 4 entries

-
1. DRY (d) {D8F4F0FD-1B19-2515-3815-EC88D175B305}
 2. MOIST (m) {AA73A9D6-6645-5DA7-C602-A537A2DFBFA0}
 3. DIAPHORETIC (d) {F5F1E644-A1BD-0137-54B5-1C02678C1D5A}
 4. NORMAL (n) {553C9A97-A3D5-C52F-0E4E-C7D168C99A2B}

UNITS:

=====

TERM: SKIN PATCHES PRESENT Page: 519
ID: {D48AF8EF-A439-4D1B-9987-A04D9E030B45}

=====

Abbreviation..... SKIN PAT P

Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... Are Skin Patches Present
Display Name..... Skin Patches Present
Help Text..... Are Skin Patches Present? (EKG,
Medication...)
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: SKIN RISK FACTORS PRESENT Page: 520
ID: {8258C84C-8C09-417D-8FB2-FC1999456654}

=====

Abbreviation..... SK RISK PR
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... Are Skin Risk Factors Present?
Display Name..... Skin Risk Factors Present
Help Text..... Does The Patient Have Spinal Cord
Injury, Paralysis, Or Neurological
Disease?
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: SKIN SWELLING EXTENSION AROUND INCISION Page: 521
ID: {AB9716EF-235C-4434-A96D-5BC289C9458A}

=====

Abbreviation..... INC SW
Active..... Yes
Boolean Value False.....
Boolean Value True.....

Data Type..... PickList
Description..... Skin Swelling Extension Around Incision
Display Name..... Skin Swelling Extension Around Incision
Help Text..... Skin Swelling Extension Around Incision
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

-
1. SKIN SWELLING EXTENSION AROUND INCISION NONE (INC SW NON)
{54C9FC46-899C-460B-9524-D6F196B9FF31}
 2. SKIN SWELLING EXTENSION AROUND INCISION EXTENDS INTO PERIWOUND < 1 cm (INC SW <1)
{ABA0527C-8FAD-4E34-9F6F-2AA97EA06E20}
 3. SKIN SWELLING EXTENSION AROUND INCISION EXTENDS INTO PERIWOUND > 1 cm (INC SW >1)
{106A2AD3-F482-4DEE-967E-6F06C6026944}
 4. SKIN SWELLING EXTENSION AROUND INCISION LIMITED TO WOUND (INC SW LIM)
{B8503D54-6925-430E-8AB6-692412347D5D}
 5. SKIN SWELLING EXTENSION AROUND INCISION FLUCTUANT (INC SW FLU)
{5D38DC7A-8E40-4D27-A039-2D6A557B4FFD}

UNITS:

=====

TERM: SKIN SWELLING EXTENSION AROUND WOUND
ID: {FE9F8966-B828-473C-9285-E189A448E304}

=====

Page: 522

Abbreviation..... WOU SW
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Skin Swelling Extension Around Wound
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

-
1. SKIN SWELLING EXTENSION AROUND WOUND NONE (WOU SW NON)
{550E3B92-2C69-403B-A997-2DBE61374ABF}
 2. SKIN SWELLING EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND < 1 cm (WOU SW <1)
{15F314A0-6B77-4EB5-92C0-FFEC4D82C75C}
 3. SKIN SWELLING EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND > 1 cm (WOU SW >1)
{EA6AC933-0737-4C28-9F31-D36F3452DA96}
 4. SKIN SWELLING EXTENSION AROUND WOUND LIMITED TO WOUND (WOU SW LIM)
{413CD1EC-309C-45D0-A1E0-E829EA615E2A}
 5. SKIN SWELLING EXTENSION AROUND WOUND FLUCTANT (WOU SW FLU)
{DB5077E7-D339-4C87-B8CA-D16269DB2935}

UNITS:

ID: {7E791776-2C73-E8F3-64A6-D965319D79BA}

=====

Abbreviation..... SKIN TEMP
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Skin Temperature
Display Name..... Skin Temperature
Help Text..... Skin Temperature
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 9 entries

1. COOL (COOL) {48A1FEA8-FD75-DAF8-174A-9C6AB2044439}
2. WARM (WARM) {E03BBEB4-5643-968F-190F-309A1293FCF9}
3. COLD (COLD) {B7B422BA-22EC-117D-0A0C-8419E9447D46}
4. HOT (HOT) {401EA6CD-12DE-AB0E-2B52-8B2E62F7C34E}
5. CONSISTENT WITH BODY TEMPERATURE (CON W TEMP)
{941FC5EF-6AF8-4DA1-A43C-95575AC74B57}
6. WARMER THAN GENERAL BODY TEMP (WARM B TEM)
{AE43F2B4-B585-43F4-9E97-D61E225A803B}
7. WARMER THAN OPPOSITE BODY PART (WARM OPP P)
{5DBF8614-5067-49C5-BD02-C1B21C9DFECC}
8. COOLER THAN OPPOSITE BODY PART (COOL OPP P)
{AFA2635C-8E80-4C71-A1B4-4AB117096FD2}
9. COOLER THAN GENERAL BODY TEMPERATURE (CoolGenT)
{03D7D4CF-9B3C-49CA-9650-213E9BABE064}

LOCATION: 49 entries

1. HEAD (Hd) {C3209F44-90C6-FF1B-3169-265B294F276A}
2. FACE (FACE) {C638EA14-7903-6A94-496A-C6560D51BEE3}
3. NECK (Ne) {D4812C3B-AAC2-1B61-679B-3CD8B5CB947C}
4. ABDOMEN (ABDOMEN) {109FE041-2521-6EF2-7F95-1FD1587796AB}
5. BACK (BACK) {22C26239-C692-4E93-B522-E5E8392BA3B1}
6. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
7. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
8. L HAND (L HAND) {639DF25A-19E3-B504-D323-4B3F6ED39A24}
9. R HAND (R HAND) {A54D3A80-8537-D079-9CCE-7DC56346D132}
10. L HIP (L HIP) {786CD881-27AD-4136-B9CC-8423CB38A431}
11. R HIP (R HIP) {A48CBCAE-5404-4934-9321-ECEF660A5350}
12. L GROIN (L GROIN) {5488E4EE-C1A2-4B02-8763-A53C412117A7}
13. R GROIN (R GROIN) {BAE1C889-8CCF-4CEB-9EBD-0243BC6F096A}
14. L THIGH (L THIGH) {16C2312F-1E0E-45E0-B675-160F40C23F57}
15. R THIGH (R THIGH) {B7A2D4A0-801E-4578-967A-F1CFB86FB96B}
16. L KNEE (L KNEE) {3A9E30B4-09A9-967A-58EF-69A01EC4F859}
17. R KNEE (R KNEE) {BCD18D53-07B6-2C41-64FC-3A07A7E979E7}
18. L FOOT (L FOOT) {2EC3FBEC-4EA2-9E7A-17D4-FDC2F9D28046}
19. R FOOT (R FOOT) {7D9F5FF6-232C-2488-4504-A721FF06928F}

- 20. L LEG (LL) {B088C809-E528-B785-21C2-14E45E1FD8AF}
- 21. R LEG (RL) {C5E40A7B-6FC4-267D-8CE8-3571D21D81FD}
- 22. L EAR (L EAR) {DA9C2DD1-D6FF-4899-A62B-1096FAA36C59}
- 23. L SHOULDER (L SHOULDER) {91CBD2FA-2743-42EB-946B-6933980D9337}
- 24. R SHOULDER (R SHOULDER) {D27AE130-442E-4DD6-B391-8D7E8AD494A2}
- 25. L ELBOW (L ELBOW) {BF4317E7-2024-47EF-9E5F-3D21D0547885}
- 26. R ELBOW (R ELBOW) {5BAFF87B-6993-4E76-862F-9F386D9B71C5}
- 27. L ANKLE (L ANKLE) {D5A692EF-9D7C-473C-AF96-7EF1C01FA2AF}
- 28. R ANKLE (R ANKLE) {4E6DF352-F399-4CE9-908C-48132C3D856B}
- 29. COCCYX (COCCYX) {49D65E9D-3D5B-40E3-B7CF-1003F80B852C}
- 30. L JAW (L JAW) {65EBD3CF-76D8-4F45-9B1B-3575D988BE8E}
- 31. R JAW (R JAW) {7846AEE2-B6F4-4B5A-909F-56E386FD82A8}
- 32. L NOSE (L NOSE) {C7707847-5135-4749-B196-002E3E59BB54}
- 33. R NOSE (R NOSE) {D80BCB19-E71F-43EA-8990-6C121EE7F7A1}
- 34. NOSE (NOSE) {85F7003E-6D4E-4A79-A534-8A57FBC2D106}
- 35. PENIS (PENIS) {2D82CE55-83B6-4160-A754-07C71FD388E8}
- 36. ABDOMEN LUQ (ABD LUQ) {3B31E1EC-CE74-4C76-BB61-10C678431FCA}
- 37. ABDOMEN RUQ (ABD RUQ) {F29C9FA1-717A-49B7-8000-95070E9AACA0}
- 38. ABDOMEN LLQ (ABD LLQ) {39831F59-D7B0-4156-AC57-4C622A93CA58}
- 39. ABDOMEN RLQ (ABD RLQ) {2C48394F-6877-4738-856C-E933644D8067}
- 40. L CHEST (L CHEST) {E332E497-FE0C-47CB-A0A7-217FBF52DC97}
- 41. R CHEST (R CHEST) {7E30EC0E-815A-4C2E-A9FC-3C2B2D6E3CD4}
- 42. CHEST (CHEST) {8C3C30FE-6E06-C35E-8226-CBF0C94014CC}
- 43. L SCAPULA (L SCAPULA) {3381EC50-771B-41C8-85C8-6EB0D208D54E}
- 44. R SCAPULA (R SCAPULA) {EA657261-4D72-40C2-B70D-2A2E27556A01}
- 45. L UPPER ARM (L UPP ARM) {5F915696-1B08-32E9-6B95-49246117D660}
- 46. R UPPER ARM (R UPP ARM) {3EF20E8E-F68C-E4F2-8F8C-FAD3B44DDAA4}
- 47. L LOWER ARM (L LOW ARM) {63AA95E8-256D-4390-8CB3-F00A1B3F91CD}
- 48. R LOWER ARM (R LOW ARM) {8DE6E277-EFFF-40DD-9E83-D7295AACB6AA}
- 49. SCROTUM (SCROTUM) {BDC109F5-0D30-4914-A62C-CEFEB448B547}

UNITS:

=====

TERM: SKIN TEMPERATURE AT INCISION SITE
ID: {832CCFB8-A5D4-49C0-A97D-2079E698D0EA}

=====

Page: 524

Abbreviation..... INC TP
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Skin Temperature at Incision Site
Display Name..... Skin Temperature at Incision Site
Help Text..... Skin Temperature at Incision Site
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

- 1. SKIN TEMPERATURE AT INCISION SITE CONSISTENT WITH BODY TEMPERATURE (INC TP CON) {D349175A-4B20-4243-BE40-E01E50C5583F}
- 2. SKIN TEMPERATURE AT INCISION SITE WARMER THAN GENERAL BODY TEMPERATURE

- (WARM BD OP) {5F888E9D-8F65-49ED-BC5D-F6F287A9D478}
3. SKIN TEMPERATURE AT INCISION SITE WARMER THAN OPPOSITE BODY PART (WARME R OP) {F518787A-325A-45F3-82B4-73C856802073}
4. SKIN TEMPERATURE AT INCISION SITE COOLER THAN GENERAL BODY TEMPERATURE (INC TP CBT) {D357C975-857D-4293-B3EF-EFA7A4586357}
5. SKIN TEMPERATURE AT INCISION SITE COOLER THAN OPPOSITE BODY PART (INC T P COP) {85307AD8-F730-43F0-BCB3-882594D9DDF8}

UNITS:

=====

TERM: SKIN TEMPERATURE AT WOUND SITE Page: 525
ID: {52AE1F41-553B-4D21-8A65-33B4CEDFA783}

=====

Abbreviation..... WOU TP
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Skin Temperature at Wound Site
Display Name..... Skin Temperature at Wound Site
Help Text..... Skin Temperature at Wound Site
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

1. SKIN TEMPERATURE AT WOUND SITE CONSISTENT WITH BODY TEMPERATURE (WOU TP CBT) {672D56C1-D1A3-4DDF-8815-F05DEE197061}
2. SKIN TEMPERATURE AT WOUND SITE WARMER THAN GENERAL BODY TEMPERATURE (WOU TP WGB) {03BB8F9A-3DAC-4994-9F09-AA1DEEF48090}
3. SKIN TEMPERATURE AT WOUND SITE WARMER THAN OPPOSITE BODY PART (WOU TP WOB) {6290BD23-24BF-4A3A-8555-DCE52D2274B3}
4. SKIN TEMPERATURE AT WOUND SITE COOLER THAN OPPOSITE BODY PART (WOU TP COB) {A904C666-2810-4211-B53A-2961DCE6894E}
5. SKIN TEMPERATURE AT WOUND SITE COOLER THAN GENERAL BODY TEMPERATURE (WOU TP CGB) {7203993E-3334-44EC-9FD9-14224253283D}

UNITS:

=====

TERM: SKIN TURGOR Page: 526
ID: {4D90E9F2-5E61-67AA-1ACA-874AED4E6C79}

=====

Abbreviation..... SKIN TURG
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList

Description..... No Description
Display Name..... Skin Turgor
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

- 1. GOOD (GOOD) {39ECA304-49F9-44B7-9735-B33CEC9D8119}
- 2. DECREASED (DECREASED) {D910FE54-A6AC-4BE3-AE35-C9E6679A93BD}
- 3. POOR (POOR) {A3D75A12-A791-446D-AF63-9435054DFC4C}

UNITS:

TERM: SKIN WOUND CLOSURE Page: 527
ID: {975D025B-F852-4212-A182-1D9469BCA5B1}

Abbreviation..... WOU CLOSUR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Wound Closure
Help Text..... Type of material used to close wound
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 8 entries

- 1. SKIN WOUND CLOSURE APPROXIMATED (WOU CL APP)
{83D48E0F-8643-436B-9ECC-80795857C322}
- 2. SKIN WOUND CLOSURE DEHISCENCE (WOU CL DEH)
{425A460F-19C0-4A4A-8E4E-712C66968C38}
- 3. SKIN WOUND CLOSURE GLUE (WOU CL GLU)
{29FE4AF7-39EE-400A-BB3D-E07786F11BFD}
- 4. SKIN WOUND CLOSURE NOT CLOSED (WOU CL NOT)
{0E69F1FA-35F6-4F34-8547-CC3602670A7B}
- 5. SKIN WOUND CLOSURE OPEN WOUND (WOU CL OPE)
{C6F2AD2C-932F-434A-8A7A-513029092FFD}
- 6. SKIN WOUND CLOSURE STAPLES (WOU CL STA)
{FE04BFE3-DF81-48E7-9FA1-1CF4ECAD71}
- 7. SKIN WOUND CLOSURE STERI-STRIPS (WOU CL STE)
{AD451791-6B49-4A51-ABC3-796D08BEF910}
- 8. SKIN WOUND CLOSURE SUTURES (WOU CL SUT)
{F2840AE4-216B-477C-9FE6-DC91E8CE5364}

UNITS:

=====

TERM: SLEEP AMOUNT Page: 528
ID: {E68F85F4-35A1-4FAC-BA86-90BF1FD5D63F}

=====

Abbreviation..... SLEEP AMOU
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Sleep Amount
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

HOUR -0.1 > 0.0 * 24.0 < 24.1 Dec Prec: 1

=====

TERM: SLEEP QUALIFIERS Page: 529
ID: {7727D9C8-C3F2-4676-A2B8-23F4DAC72B90}

=====

Abbreviation..... SLEEP QUAL
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Sleep Qualifiers
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 6 entries

-
1. NAPS ONLY (NAPS ONLY) {CF749F5C-2E4B-419A-B020-C81B30A7E482}
 2. HAS NOT SLEPT (NO SLEEP) {3BF9DB54-5D1E-422D-90AF-918F81C9EA85}
 3. ASLEEP AT ROUNDS (SLEEP ROUN)
 {24E1D9F0-D8F1-4DC8-911F-6170F343C087}
 4. DENIES HAVING SLEPT (DENIES SL)
 {1AC4B0CF-EC41-426B-89AB-68D5A9BAFADA}
 5. SLEPT WELL (SLEPT WELL) {18573E82-568F-4D2E-B9E1-6E4589EBEE2D}
 6. INTERMITTENT SLEEP (INT SLEEP)
 {5005C433-5892-4F4D-A64D-4C94AD1853E2}

UNITS:

=====

TERM: SPECIALTY BED/MATTRESS Page: 530
ID: {274B7CFD-D091-47F2-9F35-602F14D75340}

=====

Abbreviation..... SPEC BED/M
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Type of specialty Bed or Mattress
Display Name..... Specialty Bed/Mattress
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Void.....

DATA VALUE: 9 entries

-
1. BARIATRIC FRAME (BARI FRAME)
{CF95B426-F280-451B-86F2-CEF7A1A2B967}
 2. BARIATRIC LOW AIRLOSS MATTRESS REPLACEMENT (BAR L AIR)
{220C6E92-1D0A-4513-B3A9-8065A6EC7DB1}
 3. LOW AIRLOSS/ALTERNATING PRESSURE COMBO (L AIRL ALT)
{ADDD1260-C95B-4D7B-AB89-3749D76A3700}
 4. LOW AIRLOSS OVERLAY (L AIR OVRL)
{282C439C-A697-43E6-85C1-786352EB56BC}
 5. PERCUSSION BED (PERC BED) {BE0B086B-A2F3-446A-ADE0-9589E92689FE}
 6. ROTATIONAL THERAPY (ROT THERAP)
{BB795594-720C-40D7-9B7C-6792BC849DA0}
 7. STANDARD PRESSURE REDUCTION (STAN P RED)
{556537CD-A84D-4493-81F6-51CBD844B5B7}
 8. STRYKER FRAME (STRYKER FR)
{71A1DB72-C860-42D8-BA4F-B71AFA2164B1}
 9. FLOTATION BED (FLOT BED) {FA6329AB-85A6-4301-BC14-D0D089249356}

UNITS:

TERM: SPECIFIC GRAVITY
ID: {CFB61154-54BD-4DE2-B83B-BC497AEC0F33}

Page: 531

Abbreviation..... SPEC GRAVI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Specific Gravity
Display Name..... Specific Gravity
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

PRODUCT: 1 entries

-
1. URINE (u) {1F61B077-1D34-0CBA-E55B-65A8B6A7CAE7}

Void.....

METHOD: 1 entries

1. ELECTROCARDIOGRAM (ECG) {CF5F1AB9-3B74-7FBA-E38B-C9EF54FE4DC5}

UNITS:

MILLIVOLTS -0.50 > -0.05 * 0.10 < 1.00
Dec Prec: 2
MILLIMETERS -5.00 > -0.50 * 1.00 < 10.00
Dec Prec: 2

TERM: ST SEGMENT II Page: 536
ID: {6F42CE65-5097-918C-93B6-7D3000CA402E}

Abbreviation..... ST2
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... ST Segment 2
Display Name..... ST Segment 2
Help Text..... The part of an electrocardiogram
immediately following the QRS complex
and merging into the T wave
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

METHOD: 1 entries

1. ELECTROCARDIOGRAM (ECG) {CF5F1AB9-3B74-7FBA-E38B-C9EF54FE4DC5}

UNITS:

MILLIVOLTS -0.50 > -0.05 * 0.10 < 1.00
Dec Prec: 2
MILLIMETERS -5.00 > -0.50 * 1.00 < 10.00
Dec Prec: 2

TERM: ST SEGMENT III Page: 537
ID: {C2D177EB-A435-9EAA-E518-2039E8592E69}

Abbreviation..... ST 3
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... ST Segment 3
Display Name..... ST Segment 3
Help Text..... The part of an electrocardiogram

immediately following the QRS complex
and merging into the T wave

Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

1. ELECTROCARDIOGRAM (ECG) {CF5F1AB9-3B74-7FBA-E38B-C9EF54FE4DC5}

UNITS:

MILLIVOLTS -0.50 > -0.05 * 0.10 < 1.00
Dec Prec: 2
MILLIMETERS -5.00 > -0.50 * 1.00 < 10.00
Dec Prec: 2

TERM: ST SEGMENT V1 Page: 538
ID: {6D9794A3-5374-CE1D-6AD4-D4601683F6F0}

Abbreviation..... STV1
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... ST Segment from Lead V1
Display Name..... ST Segment V1
Help Text..... The part of an electrocardiogram
immediately following the QRS complex
and merging into the T wave from Lead
V1

Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

1. ELECTROCARDIOGRAM (ECG) {CF5F1AB9-3B74-7FBA-E38B-C9EF54FE4DC5}

UNITS:

MILLIVOLTS -0.50 > -0.05 * 0.10 < 1.00
Dec Prec: 2
MILLIMETERS -5.00 > -0.50 * 1.00 < 10.00
Dec Prec: 2

TERM: ST SEGMENT V2 Page: 539
ID: {BA2FE2E9-D587-5F6C-66D4-B046A3746ADE}

Abbreviation..... STV2
Active..... Yes
Boolean Value False.....

Boolean Value True.....
Data Type..... Numeric
Description..... ST Segment from Lead V2
Display Name..... ST Segment V2
Help Text..... The part of an electrocardiogram
immediately following the QRS complex
and merging into the T wave from Lead
V2
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

1. ELECTROCARDIOGRAM (ECG) {CF5F1AB9-3B74-7FBA-E38B-C9EF54FE4DC5}

UNITS:

MILLIVOLTS -0.50 > -0.05 * 0.10 < 1.00
Dec Prec: 2
MILLIMETERS -5.00 > -0.50 * 1.00 < 10.00
Dec Prec: 2

=====

TERM: ST SEGMENT V3 Page: 540
ID: {ACEEBB4D-FDFD-31BA-7151-446EAB274549}

=====

Abbreviation..... STV3
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... ST Segment from Lead v3
Display Name..... ST Segment V3
Help Text..... The part of an electrocardiogram
immediately following the QRS complex
and merging into the T wave from Lead
V3
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

1. ELECTROCARDIOGRAM (ECG) {CF5F1AB9-3B74-7FBA-E38B-C9EF54FE4DC5}

UNITS:

MILLIVOLTS -0.50 > -0.05 * 0.10 < 1.00
Dec Prec: 2
MILLIMETERS -5.00 > -0.50 * 1.00 < 10.00
Dec Prec: 2

=====

Abbreviation..... STV4
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... ST Segment from Lead 4
Display Name..... ST Segment V4
Help Text..... The part of an electrocardiogram
immediately following the QRS complex
and merging into the T wave from Lead
V4
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

METHOD: 1 entries

1. ELECTROCARDIOGRAM (ECG) {CF5F1AB9-3B74-7FBA-E38B-C9EF54FE4DC5}

UNITS:

MILLIVOLTS -0.50 > -0.05 * 0.10 < 1.00
Dec Prec: 2
MILLIMETERS -5.00 > -0.50 * 1.00 < 10.00
Dec Prec: 2

Abbreviation..... STV5
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... ST Segment from Lead 5
Display Name..... ST Segment V5
Help Text..... The part of an electrocardiogram
immediately following the QRS complex
and merging into the T wave from Lead
V5
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

METHOD: 1 entries

1. ELECTROCARDIOGRAM (ECG) {CF5F1AB9-3B74-7FBA-E38B-C9EF54FE4DC5}

UNITS:

MILLIVOLTS -0.50 > -0.05 * 0.10 < 1.00

- 3. SOFT (SOFT) {E68B4D4B-B86B-4004-A861-986645267F6E}
- 4. DIARRHEA (DIARRHEA) {353472EF-3499-46B1-AD7C-554FE01D8CE2}
- 5. CONSTIPATION (CONSTIPAT.) {4A3B0D47-88BF-4FEB-AA0D-7AD7FED2C5C3}
- 6. TARRY (TARRY) {09364271-6307-471F-A334-A8F7E051EF46}
- 7. MUCOID (MUCOID) {073E5514-DACE-4F43-ADE5-91A1A87D6601}
- 8. BLOODY (BLOODY) {D7904E84-8331-D6DD-F2C0-7F3039C5746A}
- 9. FRANK BLOOD (FRANK BLD) {CFAE6E6B-87B2-4CBA-A280-C5233B7995ED}
- 10. CLOTS (CLOTS) {B52B88E0-6D23-4657-97CD-81E3EF64991B}

UNITS:

=====

TERM: STOOL COLOR Page: 545
 ID: {85012C5D-09F3-4078-9A4D-82564BF57235}

=====

Abbreviation..... STOOL COLO
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... Stool Color
 Display Name..... Stool Color
 Help Text..... No Help Text
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 9 entries

- 1. BROWN (BROWN) {A1C32FBC-7B4D-707C-CA58-92A2E201368E}
- 2. YELLOW/GREEN (YELLOW/GRN) {1D0E7E81-3331-4D07-9B94-3D4BF72A0F0D}
- 3. BLACK (BLACK) {7A55AFB9-6861-4918-0746-3D67D47FE5C3}
- 4. RED (RED) {515D7CFB-3527-30E7-EC4B-B607C272DFA0}
- 5. MAROON (MAROON) {1F518B23-C79F-4F06-9DF3-ADC5C96D3029}
- 6. CLEAR (c) {8981C365-2B46-704F-C9F4-975DF1C23C2D}
- 7. GREEN (GREEN) {3D1F1553-56FC-68CE-79A9-E4F6357B6AAF}
- 8. RASH/EXFOLIATION (RASH/EXFOL)
 {AB83ABA9-28A7-6484-67E7-4E66B5526942}
- 9. YELLOW (YELLOW) {90DC31AD-6211-9F25-4237-79F7122A2763}

UNITS:

=====

TERM: STOOL NUMBER STOOLS Page: 546
 ID: {85348D14-F821-4043-A1C6-07D5BA407198}

=====

Abbreviation..... NUM STOOLS
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... Number of stools

Display Name..... Number of stools
Help Text..... Number of stools
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 11 entries

1. 1 (1) {8F5D10A8-FE17-47BC-DACA-6972D12CC00A}
2. 2 (2) {38E22B79-A2B6-3401-B195-C955F24E1278}
3. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}
4. 4 (4) {D2280C7E-DE50-5A8D-C3FD-5862B69C0E1A}
5. 5 (5) {188EDCE7-0726-4234-D8D2-92C8052CD221}
6. 6 (6) {640CF2F9-3BBF-AB13-185B-CD02A609FB5A}
7. 7 (7) {7A010324-D2C2-5D4B-848A-7B8878666D20}
8. 8 (8) {518E7CFC-102E-351E-C7AF-FB575E3C8F4D}
9. 9 (9) {58FF2FD4-569F-F454-9F5B-FE6F9E34D8CC}
10. 10 (10) {0C827160-18DB-816F-7FAA-A3C175C6D33B}
11. N/A (n) {B3DBE3FD-C5CB-3894-607A-BE13319594D0}

UNITS:

=====

TERM: STROKE INDEX Page: 547
ID: {14D70939-EDAB-9318-200B-D71EF2F38570}

=====

Abbreviation..... SI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Stroke Index
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

MILLILITERS PER BEAT 10 > 25 * 45 < 100 Dec Prec: 0

=====

TERM: STROKE VOLUME Page: 548
ID: {97938D03-0837-552B-5BC2-3F8D1A38E386}

=====

Abbreviation..... SV
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Stroke Volume

Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

METHOD: 8 entries

- 1. AEROSOL/HUMIDIFIED MASK (AM) {CAF54B28-D5D2-8481-FAE3-9C74C2D51551}
- 2. ASSISTED VENTILATOR (AV) {A540F333-5A09-5F21-7B51-2B5DF05E32C0}
- 3. FACE TENT (FT) {1329BF1A-DCFD-9AF3-69DB-BA6D4BC3F52A}
- 4. LARYNGEAL MASK (LARYNGMASK) {DE9B1E53-711C-DA9D-59D9-CD0EBB6CD6BE}
- 5. MASK (M) {4FB6B24A-D56E-4EE3-6CBC-E54C9734424A}
- 6. T-PIECE (TP) {27FC4BC6-2F25-2154-553C-2B452E28A91C}
- 7. TRACHEOSTOMY COLLAR (TC) {CBF8A3B9-D323-E021-FD1B-BF4CE19AD20E}
- 8. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

PERCENTAGE 0 > 21 * 100 < 101 Dec Prec: 0

=====
TERM: SUPPLEMENTAL OXYGEN FLOW RATE Page: 551
ID: {3BB314E8-9BBB-480E-B34E-B56EDE43BAC4}
=====

Abbreviation..... O2FLOWRATE
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Supplemental Oxygen Flow Rate
Display Name..... Supplemental Oxygen Flow Rate
Help Text..... Supplemental Oxygen Flow Rate
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

METHOD: 12 entries

- 1. AMBU BAG (AB) {E9D82DB9-D8A3-10E4-CE51-34BFF48F4603}
- 2. BIPAP (BIP) {3CFD1D11-EF8D-58BA-6381-B24956344D2E}
- 3. BIPAP-CPAP (PAP) {2BCBE562-8F44-662F-469B-13AC0E863A3D}
- 4. CPAP (CP) {A041FB71-08EF-F3FE-82FA-051FF74C824F}
- 5. MASK (M) {4FB6B24A-D56E-4EE3-6CBC-E54C9734424A}
- 6. NASAL CANNULA (NC) {EFFD38AE-6AF5-77FD-72E1-E0F7E3DEBA29}
- 7. NEBULIZER (NEBULIZER) {3CF184A3-2688-4988-B956-D5B0F303E8E8}
- 8. NON RE-BREATHING (NRB) {1B9AA055-FD7D-7DCC-4EDB-BDE3C821AA8E}
- 9. PARTIAL RE-BREATHING (PRB) {6F1DF170-5128-2E3F-1831-12FB0BB5E95C}
- 10. PENDANT CANNULA (PC) {A1451228-4ABB-DBCB-31AD-578BCD93012F}
- 11. RESERVOIR CANNULA (RC) {B332C44D-29BE-765C-932C-63EBFC9631FC}
- 12. VENTURI MASK (VM) {E0401675-9C4E-2AFA-2B93-F352B781E963}

UNITS:

LITERS PER MINUTE -1 > 0 * 20 < 21 Dec Prec: 0

=====

TERM: SYSTEMIC VASCULAR RESISTANCE
ID: {BD2C3F43-186B-A312-2514-186EFE223B65}

=====

Page: 552

Abbreviation..... SVR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Systemic Vascular Resistance
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

DYNE * SECONDS / CM(-5) * M2 400 > 800 * 1200 < 3000 Dec Prec: 0

=====

TERM: SYSTEMIC VASCULAR RESISTANCE INDEX
ID: {62C5A9F1-C0D2-DD0C-5A01-459C833375EE}

=====

Page: 553

Abbreviation..... SVRI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Systemic Vascular Resistance Index
Display Name..... Systemic Vascular Resistance Index
Help Text..... Systemic Vascular Resistance Index
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

1. CARDIAC OUTPUT MONITOR (COM)
 {CCBE5744-7CB9-9DFA-938F-FDDBC73B45B3}

UNITS:

DYNE * SECONDS / CM(-5) / M2 1500 > 1760 * 2600 < 4000 Dec Prec: 0

=====

TERM: SYSTOLIC PRESSURE
ID: {D66590AE-5D02-A206-A7CD-7EFB22AE591D}

=====

Page: 554

Abbreviation..... SYSTOL PRE
Active..... Yes
Boolean Value False.....

Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Systolic Pressure
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

LOCATION: 4 entries

1. L ARM (LA) {567F552F-54B5-A8BB-596F-291001746E31}
2. L LEG (LL) {B088C809-E528-B785-21C2-14E45E1FD8AF}
3. R ARM (RA) {1C09017A-F10A-543D-A7EA-F13790D78692}
4. R LEG (RL) {C5E40A7B-6FC4-267D-8CE8-3571D21D81FD}

METHOD: 5 entries

1. CUFF (Cu) {E59A2543-1860-1E1C-BEF1-6E7F6B6FBC5E}
2. DOPPLER (Dop) {52683C19-F9DC-876B-CA3F-53E7C864CF04}
3. ARTERIAL LINE (AL) {D4B3F4ED-8B33-431E-4164-BD165E823776}
4. INVASIVE (I) {5DB8746E-5446-F47A-883C-A44B91275A6E}
5. NON-INVASIVE (N-I) {F50ABFED-5FD7-FC95-20FE-1E3038C9F299}

POSITION: 3 entries

1. LYING (Ly) {DE28D2BF-6F8D-64C3-0092-4CEB5E0902B9}
2. SITTING (Sit) {DD4DCDF1-0A22-2CC5-404F-0E472C73961C}
3. STANDING (St) {0776AE42-C370-B7B3-9B4D-9C914BA8F231}

UNITS:

MILLIMETERS OF MERCURY 0 > 80 * 150 < 300 Dec Prec: 0

=====

TERM: SaO2% Page: 555
ID: {D2D98D5C-552B-9913-31FB-E4FE5BB32912}

=====

Abbreviation..... SaO2
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Arterial Oxygen Saturation
Display Name..... SaO2%
Help Text..... Arterial Oxygen Saturation
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

PERCENTAGE 50 > 92 * 100 < 101 Dec Prec: 0

=====

TERM: SpO2% Page: 556
ID: {5F84DD55-3CCF-094C-2536-B51EB7FAD999}

=====

Abbreviation..... SpO2
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Spot Oxygen Saturation
Display Name..... SpO2%
Help Text..... Oxygen Saturation
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid..... 4500637

METHOD: 2 entries

13. MONITOR (Mo) {63FC24EA-50B8-B302-D948-ABE3517CF6FC}
19. PULSE OX MONITOR (POM) {A7E3E90C-4C8B-87D7-C2F8-49D9F3846CE8}

UNITS:

PERCENTAGE 50 > 92 * 100 < 101 Dec Prec: 0

=====

TERM: SvO2% Page: 557
ID: {A7862F66-E515-838A-6F71-2099DCC02C9E}

=====

Abbreviation..... SvO2%
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Venous Oxygen Saturation
Display Name..... SvO2%
Help Text..... Venous Oxygen Saturation
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

PERCENTAGE 30 > 60 * 75 < 100 Dec Prec: 0

=====

TERM: TEMPERATURE Page: 558
ID: {0F33223E-DF2C-6B8B-5201-5E091C5F9065}

=====

Abbreviation..... T
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric

Description..... Internal body temperature of a patient.
Display Name..... Temperature
Help Text..... Enter in the numeric value of the
patients temperature with up to one
decimal.

Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid..... 4500638

LOCATION: 7 entries

-
- 1. AXILLARY (Ax) {135AA5F4-A1C7-978D-E9AB-6873841AEB57}
 - 2. CORE (C) {5076C291-4EE7-7F9B-E7A5-3BB0A78639B7}
 - 3. ORAL (O) {7A7ACA1B-ACBD-B08D-D40C-54BA8E724B8C}
 - 4. RECTAL (RECTAL) {0E7069F4-D3E6-67A1-E821-EF1879C3C267}
 - 5. SKIN (SKIN) {0803BC59-1582-AEC8-3DA2-A81363038210}
 - 6. TEMPORAL (TEMPORAL) {E312E0D7-0FE4-A41A-EAB0-ECCB085C3576}
 - 7. TYMPANIC (Tym) {0A58C14C-E447-E034-AA93-D41E6E0D2F69}

QUALITY: 2 entries

-
- 1. ACTUAL (A) {5842A4D9-56A1-7A5A-13C9-F481699C3A09}
 - 2. ESTIMATED (E) {9EAFF26B-7D87-DA66-F991-022D095FB5A0}

UNITS:

DEGREES F 70.0 > 95.0 * 100.9 < 108.0
Dec Prec: 1

DEGREES C 21.1 > 35.0 * 38.3 < 42.2 Dec Prec: 1

=====

TERM: TONOMETRY Page: 559
ID: {C06989EF-4B0F-4941-B1A7-FA9D81A480FF}

=====

Abbreviation..... TON
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Tonometry
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid..... 4688726

UNITS:

MILLIMETERS OF MERCURY 5 > 10 * 21 < 35 Dec Prec: 0

=====

TERM: TRACE ST SEGMENT Page: 560
ID: {8A29B49B-7338-7943-B6C0-9027BE094607}

=====

Abbreviation..... TSS
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Electrocardiogram Reading
Display Name..... Trace ST Segment
Help Text..... The part of an electrocardiogram
immediately following the QRS complex
and merging into the T wave
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

METHOD: 1 entries

1. ELECTROCARDIOGRAM (ECG) {CF5F1AB9-3B74-7FBA-E38B-C9EF54FE4DC5}

UNITS:

MILLIVOLTS -0.50 > -0.05 * 0.10 < 1.00
Dec Prec: 2
MILLIMETERS -5.00 > -0.50 * 1.00 < 10.00
Dec Prec: 2

TERM: TRACE ST SEGMENT 2

Page: 561

ID: {ED106F43-5AB7-8B5C-D282-8458B35169C2}

Abbreviation..... TSS2
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... ECG Reading
Display Name..... Trace ST Segment 2
Help Text..... The part of an electrocardiogram
immediately following the QRS complex
and merging into the T wave
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

METHOD: 1 entries

1. ELECTROCARDIOGRAM (ECG) {CF5F1AB9-3B74-7FBA-E38B-C9EF54FE4DC5}

UNITS:

MILLIVOLTS -0.50 > -0.05 * 0.10 < 1.00
Dec Prec: 2
MILLIMETERS -5.00 > -0.50 * 1.00 < 10.00
Dec Prec: 2

=====

TERM: TRACHEA MIDLINE DEVIATED
ID: {E952DF23-0FC7-5809-2135-C8AA60FE159B}

=====

Page: 562

Abbreviation..... TRACH MIDL
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... No Description
Display Name..... Trachea Midline Deviated
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

-
- | | |
|------------|--|
| 1. NO (n) | {FC275445-1225-9912-18A3-DF36245893D3} |
| 2. YES (y) | {483DE806-2EDB-34B0-FE11-C29EEEA649CD} |

UNITS:

=====

TERM: TRACTION TYPE
ID: {5A9CB3E0-7141-45E9-BE82-AA026A7889AB}

=====

Page: 563

Abbreviation..... TRACT TYPE
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Traction Type
Display Name..... Traction Type
Help Text..... Type of traction used
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

-
- | | |
|--------------------------------|--|
| 1. BUCKS TRACTION (BUCK TRACT) | {5E63F920-02D7-4729-8F30-A53A0A6E1EA0} |
| 2. SLING TRACTION (SLING TRAC) | {304E924C-C95A-4DBE-B30B-191BC4D26EF7} |
| 3. PIN TRACTION (PIN TRACT) | {C75104A5-3B6B-4CFE-8E2E-05877522C811} |

UNITS:

=====

TERM: TRACTION WEIGHT
ID: {0C913807-C678-4E94-BD20-B8B2EBE697BE}

=====

Page: 564

=====
Abbreviation..... TRACT WEIG
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Weight used in traction
Display Name..... Traction Weight
Help Text..... Weight used in traction
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

UNITS:

POUND -1 > 0 * 30 < 31 Dec Prec: 0

=====
TERM: TRAIN OF FOUR RESPONSE Page: 565
ID: {86E4352A-A72E-4287-8C17-ABCF97DCAF94}

=====
Abbreviation..... TOF RESP
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Train of Four response
Display Name..... Train of Four Response
Help Text..... Record the number of responses
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 4 entries

- 1. 1:4 (1:4) {3F1B75FF-C168-D413-0BDF-6581E7A4EEAD}
2. 2:4 (2:4) {7D423C9A-42B3-4FDD-999B-889193492F5C}
3. 3:4 (3:4) {E4EBD0E2-FE8C-40FA-9DCA-475DCE43B1A2}
4. 4:4 (4:4) {E9FCD5A0-995A-401D-8E85-6F0C753609B6}

UNITS:

=====
TERM: TRAIN OF FOUR SMS Page: 566
ID: {63E1D230-3238-4D03-95C5-24BF51FD0CEF}

=====
Abbreviation..... TOF SMS
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Train of Four sms
Display Name..... Train of Four SMS

Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: TURN AND SKIN CARE Page: 571
ID: {42B46B96-C43B-4F36-B7F9-A3413F02EEA7}

Abbreviation..... TURN SKIN
Active..... Yes
Boolean Value False..... Not Done
Boolean Value True..... Done
Data Type..... Boolean
Description..... Turn and Skin Care
Display Name..... Turn and Skin Care
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: ULNAR DRIFT Page: 572
ID: {70284DFB-2F51-403D-54BE-DA51C83F2DDD}

Abbreviation..... ULDT
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Ulnar Drift
Help Text..... position of the metacarpophalangeal
joints
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

-
- 1. ABSENT (a) {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0}
 - 2. PRESENT (p) {384DC20F-1625-CE1E-1C85-A3B3924BB2A0}

UNITS:

=====

TERM: UPPER LIMB EXTREMITIES STRENGTH
ID: {C713F418-605C-1E77-7834-610B6394D5F2}

=====

Page: 573

Abbreviation..... UP. LIM. S
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Upper Limb Extremities Strength
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 6 entries

-
- 1. ABSENT (a) {6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0}
 - 2. WEAK (w) {BE2D9305-E2DA-F492-66CF-164EE58B8179}
 - 3. STRONG (s) {A12154E1-544E-F372-3839-0D44E5BC476B}
 - 4. REFLEX (r) {F21585F9-5206-A6B3-2BF6-F416E5842366}
 - 5. EQUAL (Equal) {AC079584-158D-EBF2-E13A-C2E919375E4D}
 - 6. UNEQUAL (Unequal) {73319200-164F-B39F-649C-B7FB1AFBDBB8}

LOCATION: 3 entries

-
- 1. L UPPER LIMB (L Upp Limb) {E6453716-62D3-DFCB-F84D-F6F2FDFC20B7}
 - 2. R UPPER LIMB (R UPP LIMB) {5E028BFF-ED9B-014F-74F0-C4137DB074FB}
 - 3. BILATERAL WRIST (BILAT WRIS)
{89B4D30F-9E46-F8D6-14F6-1CC0843DDC21}

UNITS:

=====

TERM: URINE CHARACTERISTICS
ID: {3112FF47-AF91-5D4B-AD9B-1A2D95AA9CA7}

=====

Page: 574

Abbreviation..... URINE CHAR
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Urine Characteristics
Help Text..... No Help Text

Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 6 entries

- 1. ANURIC (ANURIC) {0CCA4B44-7A59-4042-BAB9-EB2395A63C51}
- 2. CLEAR (c) {8981C365-2B46-704F-C9F4-975DF1C23C2D}
- 3. CLOUDY (CLOUDY) {3CCFCCCC-00FC-6B00-1D58-A77C02FDB401}
- 4. HEMATURIA (HEMATURIA) {ABD21F4E-BC48-D43A-38CF-75DB7704C4B8}
- 5. PURULENT (PURULENT) {86C01494-743A-C6A2-8883-01CA5200C79B}
- 6. SEDIMENT (SEDIMENT) {4D4E4236-55CE-7ECD-FC55-6AD0CFD4D9A5}

UNITS:

TERM: URINE COLOR Page: 575
ID: {80D67904-3FE2-98BA-6810-F15D88A637B3}

Abbreviation..... URINE COLO
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Color of urine
Display Name..... Urine Color
Help Text..... Color of urine
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 14 entries

- 1. YELLOW (YELLOW) {90DC31AD-6211-9F25-4237-79F7122A2763}
- 2. WHITE (WHITE) {AB04EF54-5246-EF3A-E1D9-02E9E24369DD}
- 3. BROWN (BROWN) {A1C32FBC-7B4D-707C-CA58-92A2E201368E}
- 4. GREEN (GREEN) {3D1F1553-56FC-68CE-79A9-E4F6357B6AAF}
- 5. RED (RED) {515D7CFB-3527-30E7-EC4B-B607C272DFA0}
- 6. TAN (TAN) {AB994ED7-8916-FAE0-0E97-CD4E0EF80ED1}
- 7. PINK (PINK) {1A7FAF84-D69C-4EEB-7BC5-D21A598A1F4D}
- 8. CLEAR (c) {8981C365-2B46-704F-C9F4-975DF1C23C2D}
- 9. AMBER (AMBER) {90A13048-4C35-E01D-3482-64D782001D4A}
- 10. BLACK (BLACK) {7A55AFB9-6861-4918-0746-3D67D47FE5C3}
- 11. CLAY (CLAY) {A4E4D16C-B45C-0934-ED22-8BF640080108}
- 12. ORANGE (ORANGE) {77B1B667-212D-88BA-D8E6-5C22F6ACDC9C}
- 13. DYE PRESENT (DYE PRESEN) {E865AD3E-6DAD-C3F9-6F93-7AD6F71234A4}
- 14. TEA (TEA) {E2C0C287-DFDC-4F7B-8A67-7A339B8E8669}

UNITS:

TERM: URINE NUMBER VOIDS Page: 576

ID: {298781A7-DC46-4065-AA8E-1567F72FBCB8}

Abbreviation..... URINE VOID
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Urine Voids
Help Text..... Number of urine voids since last entry
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 11 entries

1. 1 (1) {8F5D10A8-FE17-47BC-DACA-6972D12CC00A}
2. 2 (2) {38E22B79-A2B6-3401-B195-C955F24E1278}
3. 3 (3) {E714454D-6A67-91F6-2FB7-D8F24F184328}
4. 4 (4) {D2280C7E-DE50-5A8D-C3FD-5862B69C0E1A}
5. 5 (5) {188EDCE7-0726-4234-D8D2-92C8052CD221}
6. 6 (6) {640CF2F9-3BBF-AB13-185B-CD02A609FB5A}
7. 7 (7) {7A010324-D2C2-5D4B-848A-7B8878666D20}
8. 8 (8) {518E7CFC-102E-351E-C7AF-FB575E3C8F4D}
9. 9 (9) {58FF2FD4-569F-F454-9F5B-FE6F9E34D8CC}
10. 10 (10) {0C827160-18DB-816F-7FAA-A3C175C6D33B}
11. N/A (n) {B3DBE3FD-C5CB-3894-607A-BE13319594D0}

UNITS:

TERM: URINE VOIDING METHOD

Page: 577

ID: {CF4BE586-DB30-CD14-BEE1-704461E412C5}

Abbreviation..... VOIDING
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Urine Voiding Method
Help Text..... No Help Text
Multi Select Picklist..... Yes
Term Type..... OBSERVATION
Void.....

DATA VALUE: 12 entries

1. CONTINENT (CONTINENT) {3B658506-F42E-BAF5-04C4-5AD49699E0D2}
2. DYSURIA (DYSURIA) {75AC5ECD-8975-406A-30C7-E2EB23F10F84}
3. INCONTINENT (INCONTINEN) {BA2E9C32-8D03-8FC9-5E02-EAE2026FDC63}
4. DIVERSION (DIVERSION) {087AAE9F-7FE1-5717-6423-07574C637875}
5. INDWELLING CATHETER (INDW CATH)
{EC16724C-D1EA-5650-58D4-38A89D4456F9}

- 6. EXTERNAL CATHETER (EXT. CATHE)
{C8AEE936-DD91-FE91-4DD6-6F210FB59726}
- 7. STRAIGHT CATHETER (STR. CATH)
{5F4A8A7F-E720-9458-2360-1DC14E7D05F3}
- 8. ILEAL CONDUIT (ILEAL COND)
{032CB6FE-77D2-702E-A0AE-46176F5CAA57}
- 9. SUPRAPUBIC (SUPRAPUBIC) {69E0176F-52FC-A3C4-7A4E-AD57DBD6E572}
- 10. L NEPHROSTOMY (L NEPHROST)
{D8FC4C18-081E-A7B9-E394-25F194BE9757}
- 11. R NEPHROSTOMY (R NEPHROST)
{33661137-688D-5611-4368-B1DFD8DB1469}
- 12. REPORTED (REPORTED) {1D13BAC5-D8F9-455C-B815-8E256943492E}

QUALITY: 2 entries

- 1. ACTUAL (A) {5842A4D9-56A1-7A5A-13C9-F481699C3A09}
- 2. STATED (Sta) {6362C194-B7F2-D38C-5B72-3BE5DE09438E}

UNITS:

=====
 TERM: VENTILATOR AUTO PEEP Page: 578
 ID: {EA42E7E0-23CA-43D3-984E-CC421E0A90DF}
 =====

Abbreviation..... VENT AU PE
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... Numeric
 Description..... Ventilator Auto Peep
 Display Name..... Ventilator Auto Peep
 Help Text..... Ventilator Auto Peep
 Multi Select Picklist.....
 Term Type..... OBSERVATION
 Vuid.....

METHOD: 1 entries

- 1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

 CENTIMETERS H2O -1.00 > 0.00 * 10.00 < 100.00
 Dec Prec: 2

=====
 TERM: VENTILATOR CIRCUIT CHANGE Page: 579
 ID: {C0494BB3-5896-41BD-89B4-86C684078604}
 =====

Abbreviation..... VENT CIR C
 Active..... Yes
 Boolean Value False..... No
 Boolean Value True..... Yes

Data Type..... Boolean
Description..... Ventilator Circuit Change
Display Name..... Ventilator Circuit Change
Help Text..... Ventilator Circuit Change
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

METHOD: 1 entries

- 1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

=====

TERM: VENTILATOR EXHALED VOLUME
ID: {404A3850-382D-4E86-B7F3-29DDAF1616BC}

=====

Page: 580

Abbreviation..... VENT EXHAL
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Ventilator Exhaled Volume
Display Name..... Ventilator Exhaled Volume
Help Text..... Ventilator Exhaled Volume
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

- 1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

LITER 0.0 > 1.0 * 99.0 < 100.0 Dec Prec: 1

=====

TERM: VENTILATOR EXPIRATORY TIME
ID: {004C6736-FFB2-4002-B92C-E2F17D03AF4E}

=====

Page: 581

Abbreviation..... VENT EXP T
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Ventilator Expiratory Time

Display Name..... Ventilator Expiratory Time
Help Text..... Ventilator Expiratory Time
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

SECOND 0.00 > 0.01 * 9.99 < 10.00
Dec Prec: 2

=====

TERM: VENTILATOR FLOW RATE

Page: 582

ID: {548E1F72-B501-49DD-B22B-1366CE06C2BE}

=====

Abbreviation..... VENT FL RA
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Ventilator Flow Rate
Display Name..... Ventilator Flow Rate
Help Text..... Ventilator Flow Rate
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

LITERS PER MINUTE -1 > 0 * 150 < 151 Dec Prec: 0

=====

TERM: VENTILATOR INSPIRATORY TIME

Page: 583

ID: {35948B65-6D8F-4BD4-8F4E-CE785D8D829D}

=====

Abbreviation..... VENT INS T
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Ventilator Inspiratory Time
Display Name..... Ventilator Inspiratory Time
Help Text..... Ventilator Inspiratory Time
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

SECOND 0.00 > 0.01 * 9.99 < 10.00
Dec Prec: 2

=====

TERM: VENTILATOR MEAN AIR PRESSURE
ID: {BD28B791-7E79-4116-8DF1-6296B60FCC57}

=====

Page: 584

Abbreviation..... VENT MAP
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Ventilator Mean Air Pressure
Display Name..... Ventilator Mean Air Pressure
Help Text..... Ventilator Mean Air Pressure
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

METHOD: 1 entries

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

CENTIMETERS H2O 0 > 1 * 99 < 100 Dec Prec: 0
MILLIBAR 0 > 20 * 80 < 100 Dec Prec: 0

=====

TERM: VENTILATOR MINUTE VOLUME
ID: {1C718F1D-EE59-4743-888D-7767864C4EB0}

=====

Page: 585

Abbreviation..... VENT MIN V
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Ventilator Minute Volume
Display Name..... Ventilator Minute Volume
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

METHOD: 1 entries

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

minute minus the breaths set to be delivered by the machine

Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

BREATHS PER MINUTE 0 > 1 * 119 < 120 Dec Prec: 0

=====

TERM: VENTILATOR PATIENT RATE TOTAL Page: 590
ID: {B2681753-B7DE-4E0B-8901-B7CDF22D0E45}

=====

Abbreviation..... VT PT R TO
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Total rate that the patient actually achieved
Display Name..... Ventilator Patient Rate Total
Help Text..... Value of breaths delivered by the machine plus any additional breaths that the patient takes on their own
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

BREATHS PER MINUTE 0 > 1 * 119 < 120 Dec Prec: 0

=====

TERM: VENTILATOR PCV SET PRESSURE Page: 591
ID: {BC521E68-B7AA-4B6E-AAC1-0594EE27D94E}

=====

Abbreviation..... VENT PCV S
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Pressure Control Ventilator Pressure
Display Name..... Ventilator PCV Set Pressure
Help Text..... Pressure Control Ventilator Pressure
Multi Select Picklist.....

Term Type..... OBSERVATION
Void.....

METHOD: 1 entries

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

CENTIMETERS H2O 0 > 1 * 99 < 100 Dec Prec: 0

TERM: VENTILATOR PEAK PRESSURE Page: 592
ID: {F955312E-79E0-A72D-1F93-533B1B6EAD04}

Abbreviation..... PEAK PRESS
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Ventilator Peak Pressure
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

METHOD: 1 entries

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

CENTIMETERS H2O 0 > 1 * 99 < 100 Dec Prec: 0

TERM: VENTILATOR PEEP Page: 593
ID: {4847847A-7EB7-9A8F-D565-1DF26438A616}

Abbreviation..... VENT PEEP
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... PEEP: Positive End Expiratory Pressure
Display Name..... Ventilator PEEP
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

METHOD: 1 entries

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

=====

TERM: VENTILATOR RATE

Page: 596

ID: {001D591B-6CE7-F51F-098C-79AC6B3E0E09}

=====

Abbreviation..... VENT RATE
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Ventilator Rate Set
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

BREATHS PER MINUTE 0 > 8 * 40 < 120 Dec Prec: 0

=====

TERM: VENTILATOR SENSITIVITY

Page: 597

ID: {52F9AF8D-CF31-430B-8A11-A8115BFF4B2E}

=====

Abbreviation..... VENT SENSI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Ventilator Sensitivity
Display Name..... Ventilator Sensitivity
Help Text..... Ventilator Sensitivity
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

CENTIMETERS H2O -10.0 > -9.9 * -0.1 < 0.0 Dec Prec: 1

=====

TERM: VENTILATOR TIDAL VOLUME ACTUAL

Page: 598

ID: {0592309F-98B7-DBF1-F9F7-94051584BF9F}

=====

Abbreviation..... TIDL VOL A
Active..... Yes

Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Tidal Volume Actual
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

MILLILITER 0 > 400 * 1200 < 5000 Dec Prec: 0
LITER 0.000 > 0.400 * 1.200 < 5.000
Dec Prec: 3

=====

TERM: VENTILATOR TIDAL VOLUME SET Page: 599
ID: {90D4A5CA-E21B-06F6-7AAC-FC6D4B11E7AC}

=====

Abbreviation..... TIDL VOL S
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Tidal Volume Set
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

METHOD: 1 entries

1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

MILLILITER 0 > 800 * 1200 < 5000 Dec Prec: 0
LITER 0.000 > 0.800 * 1.200 < 5.000
Dec Prec: 3

=====

TERM: VENTILATOR TYPE Page: 600
ID: {58151FCB-646B-4CC5-A89F-6D34DFB7E622}

=====

Abbreviation..... VENT TYPE
Active..... Yes
Boolean Value False.....
Boolean Value True.....

Data Type..... PickList
Description..... Type of Ventilator
Display Name..... Ventilator Type
Help Text..... Type of Ventilator
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 13 entries

-
1. VENTILATOR PB 840 (PB 840)
 {25A7B0DA-E437-4201-B8F3-3663B1B91C1A}
 2. VENTILATOR ACHIEVA P5O2 (ACHIEVA P5)
 {1690EEA5-3C3B-4D30-AC5C-214C94626699}
 3. VENTILATOR ADULT STAR (ADULT STAR)
 {879B13DA-4874-42D4-99EC-8E4ADB31B14E}
 4. VENTILATOR BEAR 1000 (BEAR 1000)
 {03C30F99-A27C-427C-B0ED-872F0DE646BD}
 5. VENTILATOR DRAEGER (DRAEGER)
 {DD18C42F-EF13-460B-8676-9ED1E429B54B}
 6. VENTILATOR PB 7200ae (PB7200ae)
 {86AE485F-ECD1-4B8E-820F-71CF159F318B}
 7. VENTILATOR SERVO (SERVO) {D6F06715-1ADD-4941-B026-B9A5E463943E}
 8. BIPHASIC POSITIVE AIRWAY PRESSURE (BiPAP)
 {4802AB1E-A0AD-BB9F-E8BA-33863BCD3466}
 9. CONTINUOUS POSITIVE AIR PRESSURE (CPAP)
 {387B42F0-F587-AAFC-ECA0-66067F18BEC6}
 10. VENTILATOR OTHER (VENT OTHER)
 {138B7B3D-9F8B-447A-9785-12AB396AE6DE}
 11. BIPAP - RESPIRONICS (BIPAP RESP)
 {AA94A723-4C68-489F-88AC-2C6F276EAB76}
 12. BIPAP - VISION (BIPAP VISI)
 {AA39C765-BA1F-4289-B459-200A57C07FFF}
 13. BIPAP/CPAP (BIPAP/CPAP) {8D767C67-3712-410F-A96B-69C731E03C5D}

METHOD: 1 entries

-
1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

=====

TERM: VENTILATOR WAVE FORM Page: 601
ID: {17F52B3E-B4C9-483A-B02F-0F4DD809E673}

=====

Abbreviation..... VENT WA FO
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Ventilator Wave Form Setting
Display Name..... Ventilator Wave Form
Help Text..... Ventilator Wave Form Setting
Multi Select Picklist..... No

Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 3 entries

- 1. SQUARE (SQUARE) {63F70678-2CAF-48E8-964D-2805D8B9A43B}
2. ACCELERATING (ACCELERATI) {514FE300-87EE-4C2C-AC70-5E57DBECAFF7}
3. DECELERATING (DECELERATI) {B259D4D4-82B8-454A-B53D-F12F4952BB61}

METHOD: 1 entries

- 1. VENTILATOR (Ven) {416FCBD3-5370-15B9-C480-48855ED486D9}

UNITS:

=====

TERM: VISION CORRECTED Page: 602
ID: {ED02AC1-EBE4-E708-684D-63D00628A94C}

=====

Abbreviation..... VC
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Vision Corrected
Help Text..... Enter Acuity value XX/20
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid..... 4688727

LOCATION: 2 entries

- 1. L EYE (L EYE) {D714A0CA-EAD6-EC35-DA54-FF027C23FF45}
2. R EYE (R EYE) {258F0C1D-67BD-ECF8-3152-EB289DEF4681}

UNITS:

FOOT -1 > 0 * 200 < 201 Dec Prec: 0

=====

TERM: VISION UNCORRECTED Page: 603
ID: {BEA5E565-D728-F5B3-0A3A-0528C42A45C4}

=====

Abbreviation..... VU
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Vision Uncorrected
Help Text..... Enter Acuity value XX/20
Multi Select Picklist.....

Term Type..... OBSERVATION
Vuid..... 4688728

UNITS:

FOOT -1 > 0 * 200 < 201 Dec Prec: 0

TERM: VT/VC Page: 604
ID: {FF8D8ABA-E89E-429A-AE65-7C25E22EF199}

Abbreviation..... VT/VC
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Tidal Volume / Vital Capacity
Display Name..... Vt/VC
Help Text..... Tidal Volume / Vital Capacity
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

LITER 0.00 > 1.00 * 99.00 < 100.00
Dec Prec: 2

TERM: WARMING COOLING DEVICE Page: 605
ID: {34F4BAF5-8388-4431-8AE1-F04670A25FE4}

Abbreviation..... WM CL DEV
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Warming, Cooling Device
Help Text..... Warming or Cooling Device
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 11 entries

-
1. AIR FLOW DEVICE (AIR FLOW D)
{60AF6215-CE2A-4EEE-9596-909E3E50BE71}
 2. K-PAD (K-PAD) {10A0C893-F279-48ED-ACAA-1DCA4AD9DB24}
 3. WARMING BLANKET (WARM BLANK)
{23D20AE5-DE25-4456-8F77-BF0411FB19A6}
 4. HUMIDIFIER (HUMIDIFIER) {A4CE67D1-F885-423B-B914-392D8CAC697E}
 5. HEAT LAMP (HEAT LAMP) {4E9DE229-F810-4AB0-9FD7-4596759CEA11}
 6. HEAT CHEM PACK (HT CHEM PA)

- {D759F891-2535-40C2-A68C-678170E56F7A}
- 7. COOL CHEM PACK (CO CHEM PA) {6DE5E644-28C2-4190-A2EF-1A766E014A65}
- 8. ICE PACK (ICE PACK) {E46ED309-D60D-45DF-BFA3-F74DD13FC4AE}
- 9. WARMING BOTTLE (WARM BOTTL) {AD1BEFE5-4CE6-4076-8095-9BE7387CFE25}
- 10. WARM HUMIDITY WITH O2 (WM HUM O2) {E62E07C3-4CA2-426B-8EBA-48C406837D9C}
- 11. COOLING BLANKET (COOL BLANK) {13A42B12-4959-4674-AF97-2D1242CFD934}

UNITS:

=====

TERM: WARMING COOLING DEVICE APPLICATION
ID: {17908A29-AF84-4217-A035-1CE6CA0F7316}

=====

Page: 606

Abbreviation..... WM CL APPL
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Warming/Cooling Device Application
Help Text..... Current State of Warming/Cooling Device
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 4 entries

- 1. WARMING ON (WARMING ON) {6FFDE928-DE1E-4B25-BA3C-EB29E0CDEB35}
- 2. WARMING OFF (WARMING OF) {98FAC977-8E2E-406A-9217-A77E1DDCDF50}
- 3. COOLING ON (COOLING ON) {911CA17D-614E-445A-969E-4F5220E469EB}
- 4. COOLING OFF (COOLING OF) {14F573ED-0D10-4E0F-A1F3-A11F29D6DB11}

UNITS:

=====

TERM: WEIGHT Page: 607
ID: {CD2D8263-6B71-0E1C-0AFE-87B4B2C12632}

=====

Abbreviation..... WT
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Weight of patient
Display Name..... Weight
Help Text..... Patients weight in lb or kg.
Multi Select Picklist.....
Term Type..... OBSERVATION

Vuid..... 4500639

METHOD: 8 entries

- 1. BED (B) {20188D9F-3815-1F7F-34EC-BA77F75B310E}
- 2. CHAIR (Ch) {BBF5C05C-220D-D052-C7E8-236E791A4FC9}
- 3. LIFT SCALE (LS) {8D741461-25C5-3955-BE2C-CBD27317CBA1}
- 4. STANDING WEIGHT (SW) {BE16FA09-16DA-510B-D235-B7F0F93AB105}
- 5. WHEELCHAIR SCALE (WcS) {8C496753-9F42-2A36-E5CA-750188306ED4}
- 6. WITH CAST OR BRACE (WCB) {47709C05-DF21-FB7D-DE7F-0699D1B7E457}
- 7. WITH PROSTHESIS (WP) {73C93D2D-934B-CA55-3557-E756A89710F5}
- 8. WITHOUT PROSTHESIS (WOP) {9AACA852-E70D-1048-FCE7-112D46C12619}

POSITION: 3 entries

- 1. STANDING (St) {0776AE42-C370-B7B3-9B4D-9C914BA8F231}
- 2. SITTING (Sit) {DD4DCDF1-0A22-2CC5-404F-0E472C73961C}
- 3. LYING (Ly) {DE28D2BF-6F8D-64C3-0092-4CEB5E0902B9}

QUALITY: 4 entries

- 1. ACTUAL (A) {5842A4D9-56A1-7A5A-13C9-F481699C3A09}
- 2. DRY (D) {B319EFF4-BA4A-85B1-2617-E00E46A54ABE}
- 3. ESTIMATED (E) {9EAF26B-7D87-DA66-F991-022D095FB5A0}
- 4. STATED (Sta) {6362C194-B7F2-D38C-5B72-3BE5DE09438E}

UNITS:

POUND 0.0 > 75.0 * 200.0 < 1500.0
Dec Prec: 1

KILOGRAM 0.0 > 34.0 * 90.0 < 680.4 Dec Prec: 1

=====

TERM: WEIGHT AT ADMISSION Page: 608
ID: {F9C859E2-ACD7-44F5-8E0D-3166C3D0A61F}

=====

Abbreviation..... WT AT ADMI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Weight at Admission
Help Text..... Weight of patient at admission time
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

POUND 0 > 1 * 500 < 501 Dec Prec: 0

KILOGRAM 0.0 > 1.0 * 226.7 < 226.8 Dec Prec: 1

=====

Abbreviation..... WEIGHT DRY
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... No Description
Display Name..... Dry Weight
Help Text..... Patient Dry Weight
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

QUALITY: 3 entries

-
- | | |
|------------------|--|
| 1. ACTUAL (A) | {5842A4D9-56A1-7A5A-13C9-F481699C3A09} |
| 2. ESTIMATED (E) | {9EAF26B-7D87-DA66-F991-022D095FB5A0} |
| 3. STATED (Sta) | {6362C194-B7F2-D38C-5B72-3BE5DE09438E} |

UNITS:

POUND	0 > 1 * 500 < 501	Dec Prec: 0
KILOGRAM	0.0 > 1.0 * 226.7 < 226.8	Dec Prec: 1

Abbreviation..... WOU BED CO
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Wound Bed Color
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 8 entries

-
- | | |
|--------------------|--|
| 1. BLACK (BLACK) | {7A55AFB9-6861-4918-0746-3D67D47FE5C3} |
| 2. GRAY (GRAY) | {D38ABEE0-52D7-4316-A713-40ACF4250599} |
| 3. GREEN (GREEN) | {3D1F1553-56FC-68CE-79A9-E4F6357B6AAF} |
| 4. PINK (PINK) | {1A7FAF84-D69C-4EEB-7BC5-D21A598A1F4D} |
| 5. RED (RED) | {515D7CFB-3527-30E7-EC4B-B607C272DFA0} |
| 6. TAN (TAN) | {AB994ED7-8916-FAE0-0E97-CD4E0EF80ED1} |
| 7. WHITE (WHITE) | {AB04EF54-5246-EF3A-E1D9-02E9E24369DD} |
| 8. YELLOW (YELLOW) | {90DC31AD-6211-9F25-4237-79F7122A2763} |

UNITS:

=====

TERM: WOUND CARE TYPE
ID: {09114B46-DF13-431B-B9F3-C15667811A2C}

=====

Page: 611

Abbreviation..... WOU CARE T
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Wound Care Type
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

-
1. TOPICAL TREATMENT (TOPICAL T)
{4E995342-8A6F-4B21-8013-F8C6A9FC7658}
 2. DRESSING (d) {FFEF0714-3B5B-E2B4-D760-3C03D45E59EA}

UNITS:

=====

TERM: WOUND DEPTH
ID: {B8FDF2BE-E6ED-49E3-8F42-4CEF94E29C47}

=====

Page: 612

Abbreviation..... WOU DEPTH
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Wound Depth
Display Name..... Wound Depth
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

CENTIMETER -1.00 > 0.00 * 200.00 < 400.00
Dec Prec: 2

=====

TERM: WOUND DISCHARGE CHARACTERISTICS
ID: {02DD62D6-E771-43DB-8CDF-F45A27263CFC}

=====

Page: 613

Abbreviation..... WOU DR CHA
Active..... Yes

Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Wound Drainage
Display Name..... Wound Discharge Characteristics
Help Text..... Wound Discharge Characteristics Pick
One
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 7 entries

- 1. WOUND DISCHARGE CHARACTERISTICS NONE (WOU DI NON)
{929F2FA5-923F-4D88-8EDF-E0EBEE9CC7A8}
- 2. WOUND DISCHARGE CHARACTERISTICS BLOODY (WOU DI BLO)
{628DBFBC-9B77-48EA-99FE-0DF000376506}
- 3. WOUND DISCHARGE CHARACTERISTICS CLOUDY (WOU DI CLO)
{159C9993-A2E1-4133-9D5F-315EF9985EC4}
- 4. WOUND DISCHARGE CHARACTERISTICS PURULENT (WOU DI PUR)
{C1E10410-B6D8-4429-BD94-0B179992EAAB}
- 5. WOUND DISCHARGE CHARACTERISTICS SEROSANGUINOUS (WOU DI SSA)
{B431D8EF-5B53-4DCF-A390-431EB5ED3704}
- 6. WOUND DISCHARGE CHARACTERISTICS SEROUS (WOU DI SER)
{2124005E-5C36-48AC-AA6B-03E06510D70A}
- 7. WOUND DISCHARGE CHARACTERISTICS VISCOUS (WOU DI VIS)
{382FA6D9-1CB2-4ECC-A458-8598ABA708F8}

UNITS:

=====

TERM: WOUND DISCHARGE COLOR Page: 614
ID: {6BF2AAEF-0D62-4D2D-9E41-F2EDCA714150}

=====

Abbreviation..... WOU DIS CO
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Wound Discharge Color
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 15 entries

- 1. AMBER (AMBER) {90A13048-4C35-E01D-3482-64D782001D4A}
- 2. BLACK (BLACK) {7A55AFB9-6861-4918-0746-3D67D47FE5C3}
- 3. BLUE (BLUE) {8FA78575-E1DD-4E6E-BC43-9E5C92E2B7E5}
- 4. BROWN (BROWN) {A1C32FBC-7B4D-707C-CA58-92A2E201368E}
- 5. COLORLESS (COLORLESS) {D093D8C5-E284-47D4-8FA6-DCD5B030B10E}
- 6. CLAY (CLAY) {A4E4D16C-B45C-0934-ED22-8BF640080108}

- 7. GREEN (GREEN) {3D1F1553-56FC-68CE-79A9-E4F6357B6AAF}
- 8. MAROON (MAROON) {1F518B23-C79F-4F06-9DF3-ADC5C96D3029}
- 9. ORANGE (ORANGE) {77B1B667-212D-88BA-D8E6-5C22F6ACDC9C}
- 10. PINK (PINK) {1A7FAF84-D69C-4EEB-7BC5-D21A598A1F4D}
- 11. RED (RED) {515D7CFB-3527-30E7-EC4B-B607C272DFA0}
- 12. TAN (TAN) {AB994ED7-8916-FAE0-0E97-CD4E0EF80ED1}
- 13. VIOLET (VIOLET) {1249AE97-CF1F-4604-A17A-5FC85F6D0185}
- 14. WHITE (WHITE) {AB04EF54-5246-EF3A-E1D9-02E9E24369DD}
- 15. YELLOW (YELLOW) {90DC31AD-6211-9F25-4237-79F7122A2763}

UNITS:

=====

TERM: WOUND DISCHARGE ODOR
ID: {67D752F7-CC83-4D76-96EB-ECE6DBA40210}

=====

Page: 615

Abbreviation..... WOU DIS OD
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Wound Discharge Odor
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
- 2. FOUL SMELLING (FOUL SMELL) {761346A9-329D-4F42-8D80-0E9D8D40653C}

UNITS:

=====

TERM: WOUND DRAIN TYPE
ID: {8EC56666-A337-AA7F-5109-A10155D4FD58}

=====

Page: 616

Abbreviation..... WOU DR TY
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Type of drain used
Display Name..... Wound Drain Type
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 10 entries

- 1. BULB SUCTION (Bulb Sucti) {B0839A82-BB79-BB9F-5E4E-0ED67F716A50}
- 2. HEMOVAC (HEMOVAC) {82DA2DA7-E20B-E571-659C-457DF42472C2}
- 3. JP DRAIN (JP Drain) {20F3B3D5-D3A4-96FD-5446-253FD031DF2B}
- 4. LUMBAR DRAIN (LUM DRAIN) {FF36A046-1799-0734-286B-2272C34E3056}
- 5. PENROSE DRAIN (Penrose D) {DE508E95-D044-ACCA-EE69-33B9A8454E6C}
- 6. SUMP (Sump) {DAD881B9-F1EE-FDAD-AEB0-68C754FDC24E}
- 7. TWILL TAPE (Twill Tape) {C52BE250-732D-03BB-B8F9-72BE4D546807}
- 8. VENTRICULOSTOMY (Ventriculo) {ECF20317-4708-7661-D4C2-9BDB4CAAC4FA}
- 9. WOUND VAC DEVICE (WOUND VAC) {A03F4978-969D-DCB8-DCE2-26EB33F097B3}
- 10. VACUUM ASSISTED CLOSURE (VAC ASS CL) {9C3DD561-FECF-47FA-BBA5-4C08FF9A4CA1}

UNITS:

=====

TERM: WOUND DRAINAGE AMOUNT

Page: 617

ID: {3BABF5D4-BD30-4814-A397-EEEEAE0EE1D13}

=====

Abbreviation..... WOU DRA AM
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Amount of drainage from wound
Display Name..... Wound Drainage Amount
Help Text..... Amount of drainage from wound
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

- 1. WOUND DRAINAGE AMOUNT NONE (WOU DR A N) {92910398-311F-D4D7-8BD5-AB66C3DC5E72}
- 2. WOUND DRAINAGE AMOUNT SCANT (W DR SCANT) {9C3CD0BD-F960-4098-B14C-17B20703B242}
- 3. WOUND DRAINAGE AMOUNT SMALL (W DR A SMA) {3D8EDF1A-EC4B-41EF-A194-73FF60074A57}
- 4. WOUND DRAINAGE AMOUNT MODERATE (W DR A MOD) {E5DBFD9C-A2AC-4009-98CB-7246F6AA702D}
- 5. WOUND DRAINAGE AMOUNT HEAVY (W DR A HEA) {BF959172-6672-4FEB-8FB7-0568405F3C83}

UNITS:

=====

TERM: WOUND DRESSING CHANGE

Page: 618

ID: {01A45F49-9EBB-44E7-8B9A-6920886F51BB}

=====

=====
Abbreviation..... WOU DR CHA
Active..... Yes
Boolean Value False..... No
Boolean Value True..... Yes
Data Type..... Boolean
Description..... Has Wound Dressing Been Changed?
Display Name..... Wound Dressing Changed
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Void.....

DATA VALUE: 2 entries

- 1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====
TERM: WOUND DRESSING TYPE

Page: 619

ID: {4C47F5C1-BA38-4FDC-A334-32EA297D8A16}

=====
Abbreviation..... WOU DR TY
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Wound Dressing Type
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Void.....

DATA VALUE: 15 entries

- 1. NONE (NONE) {D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
2. ANTI-MICROBIAL DRESSING (AB DRESSIN)
{9734C5E9-D099-47EB-93B1-D82201F64331}
3. ANTIMICROBIAL DISK (WOU DR ADI)
{1F1745EA-D763-4AFD-9CED-82E974E57692}
4. BALSAM OF PERU (WOU DR PER)
{DD412BB1-B1C1-4A48-A38F-832629448E06}
5. BISMUTH GAUZE (WOU DR BIS)
{010DED59-C808-4D92-A7F8-87AEB86B4439}
6. CALCIUM ALGINATE (CA ALGINAT)
{2B2ACD25-A77B-4D39-8598-B26FD9BAFFC2}
7. FOAM DRESSING (WOU DR FOA)
{D18A5EDD-0EF6-470E-8370-4EDB6D322477}
8. GAUZE (GAUZE) {CF10C658-4334-4F3B-A38E-E98894458333}
9. GAUZE PACKING (GAUZE PACK)
{6C093717-82CB-4966-929E-C909310649C1}

- 10. GAUZE WRAP (GAUZE WRAP) {040A9648-E6E1-4AC5-93C7-53A49537B15C}
- 11. HYDROCOLLOIDAL DRESSING (HCOLL DRES)
{D85741EF-0DB3-48CF-BC9A-88C834A72F36}
- 12. HYDROGEL DRESSING (H-GEL DRES)
{F6192D10-D6EA-42C2-AB15-270659D9925C}
- 13. HYPERTONIC HYDROGEL (WOU DR HYP)
{533D6FC3-5183-4417-AB89-9FAD9BAD0F44}
- 14. PRESSURE (PRESSURE) {387AFB7F-3973-4D77-B2E6-0724F344D8F6}
- 15. TRANSPARENT (TRANSPAREN) {0E73D24D-0EFC-08F9-196A-D30DD1BB3228}

UNITS:

=====

TERM: WOUND EDGE Page: 620
 ID: {BC1B996C-35E6-4A05-BDF6-00DF84EA90EF}

=====

Abbreviation..... WOU ED
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... No Description
 Display Name..... Wound Edge
 Help Text..... No Help Text
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 5 entries

- 1. WOUND EDGE INDISTINCT/ATTACHED (WOU ED IND)
{C899235C-4A07-4420-92BC-8F818B8E0A7F}
- 2. WOUND EDGE NON-ATTACHED (WOU ED NON)
{1D5C1317-B4FF-421C-96A0-C20C1089E672}
- 3. WOUND EDGE TUNNELING (WOU ED TUN)
{97425206-5261-46C1-905F-E03E2F071381}
- 4. WOUND EDGE UNDERMINING (WOU ED UND)
{AAA3BB3A-3A3B-4F39-9AC7-3BE1A82BDC02}
- 5. WOUND EDGE ROLLED (WOU ED ROL)
{24D9CFEC-1ED6-4069-9715-1F9FCB45BBD3}

UNITS:

=====

TERM: WOUND HEALING STATUS Page: 621
 ID: {687C2787-FDDF-4796-8195-23152E8A9FD9}

=====

Abbreviation..... WOU HEA ST
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList

Description..... No Description
Display Name..... Wound Healing Status
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 8 entries

-
1. EPITHELIALIZATION (EPITHELIAL) {5FC43A6D-921D-4736-A234-9993DB6499B7}
 2. ESCHAR (ESCHAR) {555E03C5-FBFE-4EB9-816F-25CE34476B3C}
 3. FIBRINOUS (FIBRINOUS) {D286EAE2-6902-4B7E-845D-BBB316941B2C}
 4. GRANULATING (GRANULATI) {F8930276-3954-4BFE-82B0-6C44A8E9CAE0}
 5. HYPERGRANULIZATION (HYPERGRANU) {C14B9FB0-6F33-4ACB-88B3-4B799594F845}
 6. NECROTIC (NECROTIC) {86885A34-E4EC-4BB5-BD94-59DD4B0F1A0E}
 7. SLOUGH (SLOUGH) {548EB536-1AC4-4234-8C4D-DC54E4CAA6F6}
 8. SMOOTH (SMOOTH) {5D71A0D5-52ED-4B8D-B359-AC7D500660D8}

UNITS:

=====

TERM: WOUND IDENTIFIED/SITE

Page: 622

ID: {7DDD37DD-4323-4704-BC23-7722732D9B32}

=====

Abbreviation..... WOU IDENTI
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Wound Identified
Display Name..... Wound Identified
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 5 entries

-
1. ADMISSION (ADMISSION) {483C2830-813A-46C3-A012-CCD70044B2A2}
 2. DISCHARGE (DISCHARGE) {ADF67CC2-E05B-4B8D-A233-4AE1C35B2D20}
 3. TRANSFER (TRANSFER) {DFDB4783-25BE-485B-A9BA-621F1004AAF2}
 4. POST-OPERATIVE (POST-OP) {9FA19656-849F-4A65-92FC-755A4C9E8056}
 5. POST-FALL (POST-FALL) {9BC82B93-CA4A-44E3-A3BA-7AF72BAFD44B}

LOCATION: 70 entries

-
1. R WRIST (R WRIST) {FDD9FE1D-F8A2-09AD-82D7-AA478693DD01}
 2. L WRIST (L WRIST) {7BCAFACE-FF3B-46DA-99F9-6437A90B7126}
 3. R TESTICLE (R TESTICLE) {0947D277-50AD-47F8-A76D-0EC34EB3B168}
 4. L TESTICLE (L TESTICLE) {1E376539-E89F-4472-8274-CFD11275E11F}
 5. R SHOULDER (R SHOULDER) {D27AE130-442E-4DD6-B391-8D7E8AD494A2}
 6. L SHOULDER (L SHOULDER) {91CBD2FA-2743-42EB-946B-6933980D9337}

7. RECTAL (RECTAL) {0E7069F4-D3E6-67A1-E821-EF1879C3C267}
8. NECK (Ne) {D4812C3B-AAC2-1B61-679B-3CD8B5CB947C}
9. ORAL (O) {7A7ACA1B-ACBD-B08D-D40C-54BA8E724B8C}
10. R UPPER LEG (R UPPER LE) {113BD5ED-5A0A-47CA-B100-91B244BFD04A}
11. L UPPER LEG (L UPPER LE) {1DB35895-F74C-4B88-A3EE-3E556CB6559A}
12. R LOWER LEG (R LOW LEG) {C8EA4D71-D97E-46BD-AE9D-C10954C525CD}
13. L LOWER LEG (L LOW LEG) {7B8E92FA-E065-4D06-9FB6-4234656CFD61}
14. R INNER KNEE (R INN KNEE) {2AE24BF1-EA1B-490D-9DB7-EEFC9314C560}
15. R OUTER KNEE (R OUT KNEE) {330F1DF1-8C99-4984-9965-7404ABC9562D}
16. L INNER KNEE (L INN KNEE) {79853581-CCF0-49E9-9D91-D16DE02A836D}
17. L OUTER KNEE (L OUT KNEE) {B143C935-2687-49C3-9F65-8D69C1775ABC}
18. R HIP (R HIP) {A48CBCAE-5404-4934-9321-ECEF660A5350}
19. L HIP (L HIP) {786CD881-27AD-4136-B9CC-8423CB38A431}
20. L HEAD (L HEAD) {3737DC2D-DC9E-4FC8-839D-0355987B62D1}
21. R HEAD (R HEAD) {42A95AAE-3DA9-4C96-A4E7-A15C9840217A}
22. R HAND (R HAND) {A54D3A80-8537-D079-9CCE-7DC56346D132}
23. L HAND (L HAND) {639DF25A-19E3-B504-D323-4B3F6ED39A24}
24. GENITALS (GENITALS) {83437106-3336-4925-92C8-F959561687FC}
25. GENERALIZED (GENERIZED) {17073C6C-8691-CADC-2238-CBAAB0A5C648}
26. R FOOT (R FOOT) {7D9F5FF6-232C-2488-4504-A721FF06928F}
27. L FOOT (L FOOT) {2EC3FBEC-4EA2-9E7A-17D4-FDC2F9D28046}
28. R FACE (R FACE) {D1B15677-EF6F-0125-3601-68C1FE6386CE}
29. L FACE (L FACE) {06F00A70-F475-CD40-665D-AF0569A1E75F}
30. R EYE (R EYE) {258F0C1D-67BD-ECF8-3152-EB289DEF4681}
31. L EYE (L EYE) {D714A0CA-EAD6-EC35-DA54-FF027C23FF45}
32. R ELBOW (R ELBOW) {5BAFF87B-6993-4E76-862F-9F386D9B71C5}
33. L ELBOW (L ELBOW) {BF4317E7-2024-47EF-9E5F-3D21D0547885}
34. R EAR (R EAR) {AE31327D-A03C-488C-8A8E-084E7E6E7FDE}
35. L EAR (L EAR) {DA9C2DD1-D6FF-4899-A62B-1096FAA36C59}
36. R CHEST (R CHEST) {7E30EC0E-815A-4C2E-A9FC-3C2B2D6E3CD4}
37. L CHEST (L CHEST) {E332E497-FE0C-47CB-A0A7-217FBF52DC97}
38. R GLUTEUS (R GLUTEUS) {B4D32045-68D1-4457-9106-7937887FA6F6}
39. L GLUTEUS (L GLUTEUS) {88C4E547-04CD-4D52-9416-CAEB6D898D21}
40. UPPER BACK (UPPER BACK) {C1DA4294-E0DC-438B-810D-7F91E170244A}
41. MID BACK (MID BACK) {F8C515E0-F437-4596-B030-DA1811038CDF}
42. LOWER BACK (LOWER BACK) {7E761680-61BC-450F-9905-C45CE760EC45}
43. R UPPER ARM (R UPP ARM) {3EF20E8E-F68C-E4F2-8F8C-FAD3B44DDAA4}
44. L UPPER ARM (L UPP ARM) {5F915696-1B08-32E9-6B95-49246117D660}
45. R LOWER ARM (R LOW ARM) {8DE6E277-EFFF-40DD-9E83-D7295AACB6AA}
46. L LOWER ARM (L LOW ARM) {63AA95E8-256D-4390-8CB3-F00A1B3F91CD}
47. R ANKLE (R ANKLE) {4E6DF352-F399-4CE9-908C-48132C3D856B}
48. L ANKLE (L ANKLE) {D5A692EF-9D7C-473C-AF96-7EF1C01FA2AF}
49. UPPER ABDOMEN (UP ABDOM) {234C677B-1CE6-498B-9375-BB5DC3958842}
50. R ABDOMEN (R ABDOMEN) {5C45EAA8-96F5-4EC7-8AB5-A9543CD6C414}
51. LOWER ABDOMEN (LO ABDOM) {077C9CF5-DCA7-4D94-A7A3-F3C153CEA47D}
52. L ABDOMEN (L ABDOMEN) {46B37481-4105-492D-96F0-E2D3895379C1}
53. COCCYX (COCCYX) {49D65E9D-3D5B-40E3-B7CF-1003F80B852C}
54. HEAD OCCIPITAL (HEAD OCCIP)
{2066CA47-5A16-447F-9919-24ED435B627F}
55. L ANTECUBITAL (L ANTECUBI)
{C21838CC-BAB5-1D7F-E2A4-A7C663FB28C3}
56. R ANTECUBITAL (R ANTECUBI)
{7FB00A65-B4B2-06A6-BD89-98AC112D8D4E}
57. L FEMORAL (LFEM) {50AF6A08-CE37-6859-F4E6-762660A982F7}
58. R FEMORAL (R FEMORAL) {B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
59. L HEEL (L HEEL) {C9154881-17C7-4921-A0EE-BFB962C59CEF}

- 60. R HEEL (R HEEL) {0020DD20-2EA4-4362-ADF0-07BFA220C60E}
- 61. L ISCHIAL (L ISCHIAL) {DCDAD741-9599-4479-9D55-3279F800DBE6}
- 62. R ISCHIAL (R ISCHIAL) {97E50F6C-CF75-4E59-AFDE-80FC5CE6D9E3}
- 63. L RIBS (L RIBS) {BDFEE74F-B3CC-4AD9-AE49-78B0BF19A219}
- 64. R RIBS (R RIBS) {260FD564-5E58-419B-9F73-62B6037D1572}
- 65. L SCAPULA (L SCAPULA) {3381EC50-771B-41C8-85C8-6EB0D208D54E}
- 66. R SCAPULA (R SCAPULA) {EA657261-4D72-40C2-B70D-2A2E27556A01}
- 67. L TROCHANTER (L TROCHANT) {C91647A4-255D-40E4-A24D-01274E5D799F}
- 68. R TROCHANTER (R TROCHANT) {1A3E21B2-C391-40D2-8059-12290E667D9E}
- 69. PERINEAL (PERINEAL) {1801441C-D5EF-4775-9E58-E4985496F7E0}
- 70. SACRAL (SACRAL) {0CF8A8A7-AB35-B803-4544-0C91645D9D1C}

UNITS:

=====

TERM: WOUND LENGTH Page: 623
 ID: {B43EAD3E-A282-4351-B83A-1A4F9ADD0B46}

=====

Abbreviation..... WOU LEN
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... Numeric
 Description..... Wound Length
 Display Name..... Wound Length
 Help Text..... No Help Text
 Multi Select Picklist.....
 Term Type..... OBSERVATION
 Vuid.....

UNITS:

CENTIMETER -1.00 > 0.00 * 200.00 < 400.00
Dec Prec: 2

=====

TERM: WOUND ODOR Page: 624
 ID: {F71814E1-2D36-4074-B8F4-253AC6BDA1AA}

=====

Abbreviation..... WOU ODOR
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... Wound Odor
 Display Name..... Wound Odor
 Help Text..... Wound Odor
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 2 entries

1. WOUND ODOR NONE (WOU OD NON)
{D2E388F7-B531-4D60-A2BE-830F11B9AA75}
2. WOUND ODOR FOUL (WOU OD FOU)
{99A74112-EA71-4896-A4DD-911379E5B453}

UNITS:

=====

TERM: WOUND PHOTOGRAPHED Page: 625
 ID: {398E171E-93A6-4DB3-82E7-5C0D118A80EE}

=====

Abbreviation..... WOU PHOTOG
 Active..... Yes
 Boolean Value False..... No
 Boolean Value True..... Yes
 Data Type..... Boolean
 Description..... Has Wound Been Photographed
 Display Name..... Wound Photographed
 Help Text..... No Help Text
 Multi Select Picklist.....
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 2 entries

1. NO (n) {FC275445-1225-9912-18A3-DF36245893D3}
2. YES (y) {483DE806-2EDB-34B0-FE11-C29EEEA649CD}

UNITS:

=====

TERM: WOUND TREATMENT Page: 626
 ID: {4883E272-E8D4-48D4-A7C5-B00F5315B6DC}

=====

Abbreviation..... WOU TREATM
 Active..... Yes
 Boolean Value False.....
 Boolean Value True.....
 Data Type..... PickList
 Description..... Wound Treatment
 Display Name..... Wound Treatment
 Help Text..... No Help Text
 Multi Select Picklist..... No
 Term Type..... OBSERVATION
 Vuid.....

DATA VALUE: 5 entries

1. HYDROCOLLOIDAL DRESSING (HCOLL DRES)
{D85741EF-0DB3-48CF-BC9A-88C834A72F36}
2. TRANSPARENT DRESSING (TRANSP DRE)
{777C7722-AEA8-4A0B-9ADD-1157D37AE539}

- 3. SALINE MOIST DRESSING (SAL MOI DR) {19F98514-8ABF-461C-9969-C31177FE7DE7}
- 4. OPEN TO AIR (OPE TO AIR) {EFB715D0-5666-4281-9034-A0326A14B2BB}
- 5. XEROFORM (XEROFORM) {7D7C3CDC-253E-4262-B8D4-3B70DE9AEBE4}

UNITS:

=====

TERM: WOUND TYPE Page: 627
ID: {1DCC28AC-700F-4E61-85BC-71D00758F9EB}

=====

Abbreviation..... WOU TYPE
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... Type of wound
Display Name..... Wound Type
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 29 entries

- 1. ABRASION (ABRASION) {7E962B43-24B1-46D2-A79A-54811117271A}
- 2. BLISTERING (BLISTERING) {3A4212D5-F2A8-4636-BB7F-F36FA8CED192}
- 3. BRUISED (BRUISED) {FA5B38A7-7744-186F-055A-8458D696AF33}
- 4. BURN (BURN) {97B37D7A-28C6-4ABF-8D09-1FC49D5F160F}
- 5. CELLULITIS (CELLULITIS) {709F8795-0D4A-4E98-AB82-DD3E2AA7EA7B}
- 6. CRUSHING (CRUSHING) {92D4BAF1-63E3-4FC0-B096-B5CBC4736245}
- 7. DIABETIC/VASCULAR ULCER (DIA/VAS UL) {E209E53E-EFF1-48F5-9581-75236C89D9B9}
- 8. HEMATOMA (HEMATOMA) {6A24B7EE-62F3-C952-C1A3-086CDBB1FE68}
- 9. PREVIOUS PERCUTANEOUS SITE (PRE PER SI) {F2ECCF5C-D91A-40F5-9C09-452C2C81DA66}
- 10. RASH/EXFOLIATION (RASH/EXFOL) {AB83ABA9-28A7-6484-67E7-4E66B5526942}
- 11. SKIN TEAR (SKIN TEAR) {2E230C6E-119F-4F9F-8797-2FFA9D76F337}
- 12. SKIN TRAUMA (SKIN TRAUM) {91C8EE13-2C4A-468C-94EB-C53E0D153EBA}
- 13. SURGICAL DONOR (SURG DONOR) {F3EE10EA-5776-479B-A58C-F59D6526F90B}
- 14. SURGICAL GRAFT (SURG GRAFT) {328BA3EA-409B-4C5E-92E7-C6ECE8F983E1}
- 15. SURGICAL (SURGICAL) {2C6CFEF1-6122-45BA-A9F8-F33B5B3C0E90}
- 16. VASCULAR (VASCULAR) {60523BF8-6372-4760-881C-A62F85D2A0F6}
- 17. VASCULAR/NEUROPATHIC (VASC/NEURO) {6B91CEED-4863-46EA-A9AE-3B0D522F1BDC}
- 18. BITE (Bite) {FC33E0AE-3130-4881-B9EF-9B81FBEC6AC0}
- 19. BULLET (BULLET) {1570D8E2-D936-494D-AF6E-13838EA34EFD}
- 20. CONTUSION (Contusion) {5285197A-CC0A-4D12-A1F8-966765C687F1}
- 21. INCONTENENCE DERMATITIS (INCON DERM) {9F6C5BA2-79C7-4B09-B19A-6DDC444936E5}
- 22. LACERATION (LACERATION) {F4548BB2-39AF-4E0D-B52F-472EA080D6FF}

- 23. PENETRATING (PENETRATIN) {00A2D4C8-EDDE-4593-A1F8-A87A836B3379}
- 24. PRESSURE ULCER (PRE ULCER)
{7279639E-498D-4EA3-A3FA-4A6D611E6F48}
- 25. PUNCTURE (PUNCTURE) {28CF34AF-BD37-4BA4-A59B-C4AFCE2EB8DD}
- 26. TAPE BURN (TAPE BURN) {DDD75AFD-DA8D-4E73-873F-86D4DF4AF47A}
- 27. VASCULAR LESION ARTERIAL (VASLES ART)
{EA565D3F-D984-4553-B597-DFE844554CE9}
- 28. VASCULAR LESION MIXED (VASLES MIX)
{8684AB2D-9162-4DC2-A5FA-A14558678689}
- 29. VASCULAR LESION VENOUS (VASLES VEN)
{3F8CEC69-5E13-4F26-AB86-34FE2805E3C6}

UNITS:

=====

TERM: WOUND VACUUM ASSISTED CLOSURE
ID: {04748983-B61B-4121-AF49-FA1FDA1D9646}

=====

Page: 628

Abbreviation..... WOU VAC AS
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... PickList
Description..... No Description
Display Name..... Vacuum Assisted Closure Type
Help Text..... No Help Text
Multi Select Picklist..... No
Term Type..... OBSERVATION
Vuid.....

DATA VALUE: 2 entries

- 1. INTERMITTENT (INTERMITTE) {DD09ED7C-B977-412F-A678-5130F77D239E}
- 2. CONSTANT (CONSTANT) {0FB8E43C-2106-4C37-9430-38675953A85B}

UNITS:

=====

TERM: WOUND WIDTH
ID: {7EED882C-3CAE-4A54-88FA-D0B380D81658}

=====

Page: 629

Abbreviation..... WOU WIDTH
Active..... Yes
Boolean Value False.....
Boolean Value True.....
Data Type..... Numeric
Description..... Wound Width
Display Name..... Wound Width
Help Text..... No Help Text
Multi Select Picklist.....
Term Type..... OBSERVATION
Vuid.....

UNITS:

CENTIMETER -1.00 > 0.00 * 200.00 < 400.00
 Dec Prec: 2

=====

TYPE: DATA VALUE

Page: 631

=====

DATA VALUE: 0

ABBREVIATION.....0
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....0
HELP_TEXT.....No Help Text
ID.....{2B57EA26-2ACE-1CD8-EC5A-73E999C4F3FF}
TERM.....0

DATA VALUE: 1

ABBREVIATION.....1
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....1
HELP_TEXT.....No Help Text
ID.....{8F5D10A8-FE17-47BC-DACA-6972D12CC00A}
TERM.....1

DATA VALUE: 2

ABBREVIATION.....2
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....2
HELP_TEXT.....No Help Text
ID.....{38E22B79-A2B6-3401-B195-C955F24E1278}
TERM.....2

DATA VALUE: 3

ABBREVIATION.....3
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....3
HELP_TEXT.....No Help Text
ID.....{E714454D-6A67-91F6-2FB7-D8F24F184328}
TERM.....3

DATA VALUE: 4

ABBREVIATION.....4

ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....4
HELP_TEXT.....No Help Text
ID.....{D2280C7E-DE50-5A8D-C3FD-5862B69C0E1A}
TERM.....4

Page: 632

DATA VALUE: 4.5

ABBREVIATION.....4.5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....4.5
HELP_TEXT.....No Help Text
ID.....{EC6C0C72-8C49-4535-92B2-7E543007EE32}
TERM.....4.5

DATA VALUE: 5

ABBREVIATION.....5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....5
HELP_TEXT.....No Help Text
ID.....{188EDCE7-0726-4234-D8D2-92C8052CD221}
TERM.....5

DATA VALUE: 5.5

ABBREVIATION.....5.5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....5.5
HELP_TEXT.....No Help Text
ID.....{4878F059-EC9D-44FC-9A61-734C291530E7}
TERM.....5.5

DATA VALUE: 6

ABBREVIATION.....6
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....6
HELP_TEXT.....No Help Text
ID.....{640CF2F9-3BBF-AB13-185B-CD02A609FB5A}
TERM.....6

DATA VALUE: 6.5

ABBREVIATION.....6.5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....6.5

HELP_TEXT.....No Help Text
ID.....{CAFCBDB1-8533-4C10-8DD1-2BD07A4CF6D3}
TERM.....6.5

Page: 633

DATA VALUE: 7

ABBREVIATION.....7
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....7
HELP_TEXT.....No Help Text
ID.....{7A010324-D2C2-5D4B-848A-7B8878666D20}
TERM.....7

DATA VALUE: 7.5

ABBREVIATION.....7.5
ACTIVE.....Yes
DESCRIPTION.....7.5
DISPLAY_NAME.....7.5
HELP_TEXT.....No Help Text
ID.....{7517EE9C-7C8B-4AF0-B60A-88F4E24BFA5E}
TERM.....7.5

DATA VALUE: 8

ABBREVIATION.....8
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....8
HELP_TEXT.....No Help Text
ID.....{518E7CFC-102E-351E-C7AF-FB575E3C8F4D}
TERM.....8

DATA VALUE: 8.5

ABBREVIATION.....8.5
ACTIVE.....Yes
DESCRIPTION.....8.5
DISPLAY_NAME.....8.5
HELP_TEXT.....8.5
ID.....{96723567-151D-4B5C-9C65-B3EB6D584272}
TERM.....8.5

DATA VALUE: 9

ABBREVIATION.....9
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....9
HELP_TEXT.....No Help Text
ID.....{58FF2FD4-569F-F454-9F5B-FE6F9E34D8CC}
TERM.....9

DATA VALUE: 9.5

ABBREVIATION.....9.5
ACTIVE.....Yes
DESCRIPTION.....9.5
DISPLAY_NAME.....9.5
HELP_TEXT.....9.5
ID.....{252555A1-ADB2-4948-8764-C952F2587DC5}
TERM.....9.5

DATA VALUE: 10

ABBREVIATION.....10
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....10
HELP_TEXT.....No Help-Text
ID.....{0C827160-18DB-816F-7FAA-A3C175C6D33B}
TERM.....10

DATA VALUE: 11

ABBREVIATION.....11
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....11
HELP_TEXT.....11
ID.....{BD51A7B3-FACF-890C-472E-E3B28519CC16}
TERM.....11

DATA VALUE: 12

ABBREVIATION.....12
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....12
HELP_TEXT.....No Help Text
ID.....{ACC5CA06-3ABC-476D-EFB0-D82FEEEE129C7}
TERM.....12

DATA VALUE: 12.5

ABBREVIATION.....12.5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....12.5
HELP_TEXT.....No Help Text
ID.....{97A947AE-2BCC-F75D-1794-6AF9BAC1497B}
TERM.....12.5

DATA VALUE: 13

ABBREVIATION.....13
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....13
HELP_TEXT.....No Help Text
ID.....{030242DB-3F82-7517-3DEE-4D50E3E8AE23}
TERM.....13

DATA VALUE: 13.5

ABBREVIATION.....13.5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....13.5
HELP_TEXT.....No Help Text
ID.....{34D24648-B2D9-9C6C-6D71-F80E977CE185}
TERM.....13.5

DATA VALUE: 14

ABBREVIATION.....14
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....14
HELP_TEXT.....No Help Text
ID.....{EDF1C2B1-123B-011F-9B02-7D71088971C1}
TERM.....14

DATA VALUE: 14.5

ABBREVIATION.....14.5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....14.5
HELP_TEXT.....No Help Text
ID.....{373C8F8C-5D36-4AD6-6A76-513064669CB5}
TERM.....14.5

DATA VALUE: 15

ABBREVIATION.....15
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....15
HELP_TEXT.....No Help Text
ID.....{FE8C049F-B64C-4991-2F50-3C561BA4825C}
TERM.....15

DATA VALUE: 15.5

ABBREVIATION.....15.5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....15.5
HELP_TEXT.....No Help Text
ID.....{1BBF5122-2FFA-6B27-445B-3291C3F83A36}
TERM.....15.5

DATA VALUE: 16

ABBREVIATION.....16
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....16
HELP_TEXT.....No Help Text
ID.....{9E4B84FF-430A-39FB-E6AA-E1031814796D}
TERM.....16

DATA VALUE: 16.5

ABBREVIATION.....16.5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....16.5
HELP_TEXT.....No Help Text
ID.....{7D4EA1D4-5233-760F-FD2B-A70FB0969815}
TERM.....16.5

DATA VALUE: 17

ABBREVIATION.....17
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....17
HELP_TEXT.....No Help Text
ID.....{2E89B9C8-5E80-F3C9-FE66-4396CB63201A}
TERM.....17

DATA VALUE: 17.5

ABBREVIATION.....17.5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....17.5
HELP_TEXT.....No Help Text
ID.....{FDE6FFD9-00D0-9C55-1A7F-4AA70CE313BA}
TERM.....17.5

Page: 637

DATA VALUE: 18

ABBREVIATION.....18
ACTIVE.....Yes
DESCRIPTION.....No Description

DISPLAY_NAME.....18
HELP_TEXT.....No Help Text
ID.....{F2815231-0FCE-4988-95E9-D9AF46C57086}
TERM.....18

DATA VALUE: 18.5

ABBREVIATION.....18.5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....18.5
HELP_TEXT.....No Help Text
ID.....{75ABCB85-4D6E-0BBC-157F-CA69B7737402}
TERM.....18.5

DATA VALUE: 19

ABBREVIATION.....19
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....19
HELP_TEXT.....No Help Text
ID.....{478FFB89-4941-1474-04E1-C838B2467439}
TERM.....19

DATA VALUE: 19.5

ABBREVIATION.....19.5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....19.5
HELP_TEXT.....No Help Text
ID.....{58F0F98D-299A-8C4E-F50F-AE21B6E696D6}
TERM.....19.5

DATA VALUE: 20

ABBREVIATION.....20
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....20
HELP_TEXT.....No Help Text
ID.....{532FBDF8-7F7B-9551-300C-F2CFEF656163}
TERM.....20

Page: 638

DATA VALUE: 20.5

ABBREVIATION.....20.5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....20.5
HELP_TEXT.....No Help Text
ID.....{DB1F50D2-AF99-7BB3-88F5-B58C0BDE7987}

TERM.....20.5

DATA VALUE: 21

ABBREVIATION.....21
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....21
HELP_TEXT.....No Help Text
ID.....{C6188159-9EF1-DA5C-C015-2142A1205C4C}
TERM.....21

DATA VALUE: 21.5

ABBREVIATION.....21.5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....21.5
HELP_TEXT.....No Help Text
ID.....{A259E466-773F-519C-7B5B-1570D1B03B70}
TERM.....21.5

DATA VALUE: 22

ABBREVIATION.....22
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....22
HELP_TEXT.....No Help Text
ID.....{DA839156-C7A9-76D7-58CD-7EC2B016EF40}
TERM.....22

DATA VALUE: 22.5

ABBREVIATION.....22.5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....22.5
HELP_TEXT.....No Help Text
ID.....{D402D24F-1EF2-1AA0-DD1C-FA8320AA929A}
TERM.....22.5

Page: 639

DATA VALUE: 23

ABBREVIATION.....23
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....23
HELP_TEXT.....No Help Text
ID.....{8FB6A8B1-B692-BDC7-43B9-0332F33D826A}
TERM.....23

DATA VALUE: 23.5

ABBREVIATION.....23.5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....23.5
HELP_TEXT.....No Help Text
ID.....{DDB76295-7E4A-FEAC-324D-37878F62771A}
TERM.....23.5

DATA VALUE: 24

ABBREVIATION.....24
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....24
HELP_TEXT.....No Help Text
ID.....{DBE934CF-40B3-2263-9FA7-6AC9DE0C7C8A}
TERM.....24

DATA VALUE: 24.5

ABBREVIATION.....24.5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....24.5
HELP_TEXT.....No Help Text
ID.....{D815830B-5D94-0975-D007-E4FDD966155B}
TERM.....24.5

DATA VALUE: 25

ABBREVIATION.....25
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....25
HELP_TEXT.....No Help Text
ID.....{42DD0655-142D-6901-D3FE-66D295BBFCD1}
TERM.....25

Page: 640

DATA VALUE: 25.5

ABBREVIATION.....25.5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....25.5
HELP_TEXT.....No Help Text
ID.....{9D704E9E-FB1D-FB83-2142-EDBF3805E8F7}
TERM.....25.5

DATA VALUE: 26

ABBREVIATION.....26
ACTIVE.....Yes

DESCRIPTION.....No Description
DISPLAY_NAME.....26
HELP_TEXT.....No Help Text
ID.....{181395F6-29FD-5A41-5AD7-BFC9CF564F3B}
TERM.....26

DATA VALUE: 26.5

ABBREVIATION.....26.5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....26.6
HELP_TEXT.....No Help Text
ID.....{1D464967-6CFA-E29E-548D-EF68AFCDEC61}
TERM.....26.5

DATA VALUE: 27

ABBREVIATION.....27
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....27
HELP_TEXT.....No Help Text
ID.....{64D25A86-6516-0910-B4B5-644B85378876}
TERM.....27

DATA VALUE: 27.5

ABBREVIATION.....27.5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....27.5
HELP_TEXT.....No Help Text
ID.....{8C7BE365-AFD0-3FFD-B8EC-6EC914763CFC}
TERM.....27.5

Page: 641

DATA VALUE: 28

ABBREVIATION.....28
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....28
HELP_TEXT.....No Help Text
ID.....{66B5A1F5-5396-10C4-643E-3E8CAB6D6CC8}
TERM.....28

DATA VALUE: 28.5

ABBREVIATION.....28.5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....28.5
HELP_TEXT.....No Help Text

ID.....{C99C9298-D4D2-BC02-36D3-F612B6080A87}
TERM.....28.5

DATA VALUE: 29

ABBREVIATION.....29
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....29
HELP_TEXT.....29
ID.....{93F62F6F-5EFE-C451-5908-47C895B7177C}
TERM.....29

DATA VALUE: 30

ABBREVIATION.....30
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....30
HELP_TEXT.....30
ID.....{1CF40F33-4B6B-FE03-98FB-55DFA8FB4145}
TERM.....30

DATA VALUE: 31

ABBREVIATION.....31
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....31
HELP_TEXT.....31
ID.....{5EEDDCF5-7D59-1526-2522-FA59099C8C21}
TERM.....31

Page: 642

DATA VALUE: 99

ABBREVIATION.....9
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....99
HELP_TEXT.....No Help Text
ID.....{B013C05B-1AF0-1547-4F03-378883F887C5}
TERM.....99

DATA VALUE: +1 = SLIGHT

ABBREVIATION.....+1 SLIGHT
ACTIVE.....Yes
DESCRIPTION.....+1 = Slight
DISPLAY_NAME.....+1 = Slight
HELP_TEXT.....No Help Text
ID.....{7C0245BA-B261-4B68-A84C-D737CA634230}
TERM.....+1 = SLIGHT

DATA VALUE: +2 = MODERATE

ABBREVIATION.....+2 = MODER
ACTIVE.....Yes
DESCRIPTION.....+2 = Moderate
DISPLAY_NAME.....+2 = Moderate
HELP_TEXT.....No Help Text
ID.....{CCF23A60-068C-443F-AFD0-2D56EEE0446C}
TERM.....+2 = MODERATE

DATA VALUE: +3 = MARKED

ABBREVIATION.....+3 = MARKE
ACTIVE.....Yes
DESCRIPTION.....+3 = Marked
DISPLAY_NAME.....+3 = Marked
HELP_TEXT.....No Help Text
ID.....{41ACC586-0C68-4780-8B2D-F36B21FCEDB2}
TERM.....+3 = MARKED

DATA VALUE: +4 = VERY MARKED

ABBREVIATION.....+4 = VY MA
ACTIVE.....Yes
DESCRIPTION.....+4 = Very Marked
DISPLAY_NAME.....+4 = Very Marked
HELP_TEXT.....No Help Text
ID.....{55EF9D42-B4F8-4468-BA1B-98BCD4075E58}
TERM.....+4 = VERY MARKED

Page: 643

DATA VALUE: 0 - ABSENT

ABBREVIATION.....0-a
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....0 - Absent
HELP_TEXT.....No Help Text
ID.....{75F14311-F113-963C-4F05-B8AE17F044A9}
TERM.....0 - ABSENT

DATA VALUE: 0 - NO ANXIETY AT EASE

ABBREVIATION.....CIWA AX 0
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....0 - No Anxiety, At Ease
HELP_TEXT.....No Help Text
ID.....{DA61A13E-6443-4E8A-9C76-762FEACDB4E7}
TERM.....0 - NO ANXIETY AT EASE

DATA VALUE: 0 - NO NAUSEA/VOMITING

ABBREVIATION.....CIWA NV 0

ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....0 - No Nausea / Vomiting
HELP_TEXT.....No Help Text
ID.....{0569555E-E7B3-49E3-9A18-68D23CFB2B8D}
TERM.....0 - NO NAUSEA/VOMITING

DATA VALUE: 0 - NO SWEAT VISIBLE

ABBREVIATION.....CIWA PS 0
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....0 - No Sweat Visible
HELP_TEXT.....No Help Text
ID.....{7C83E616-4414-4766-AAF5-5F5442FFBF09}
TERM.....0 - NO SWEAT VISIBLE

DATA VALUE: 0 - NO TREMOR

ABBREVIATION.....CIWA TR 0
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....0 - No Tremor
HELP_TEXT.....No Help Text
ID.....{2D1A5B39-4C5D-4345-84B4-DDCBC8E25B4D}
TERM.....0 - NO TREMOR

Page: 644

DATA VALUE: 0 - NONE

ABBREVIATION.....CIWA TC 0
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....0 - None
HELP_TEXT.....No Help Text
ID.....{5F9F33C7-EFCE-4AC1-BCA3-1967F68EB597}
TERM.....0 - NONE

DATA VALUE: 0 - NORMAL ACTIVITY

ABBREVIATION.....CIWA AG 0
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....0 - Normal Activity
HELP_TEXT.....No Help Text
ID.....{9FF32084-A432-467C-B0A2-C8BC9F3D6D13}
TERM.....0 - NORMAL ACTIVITY

DATA VALUE: 0 - ORIENTED AND CAN DO SERIAL ADDITIONS

ABBREVIATION.....CIWA OR 0
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....0 - Oriented And Can Do Serial

Additions

HELP_TEXT.....No Help Text
ID.....{D8015F54-6E9E-457A-955C-38565BDDDBB6}
TERM.....0 - ORIENTED AND CAN DO SERIAL
ADDITIONS

DATA VALUE: 0 = 0 EXTREMITIES

ABBREVIATION.....ALD ACT 0
ACTIVE.....Yes
DESCRIPTION.....0 = 0 Extremities
DISPLAY_NAME.....0 = Moves 0 Extremities
HELP_TEXT.....0 = 0 Extremities
ID.....{6E248A9A-B509-4BD0-A01D-9ED5CE5B8FC9}
TERM.....0 = 0 EXTREMITIES

DATA VALUE: 0 = <90%

ABBREVIATION.....ALD O2 0
ACTIVE.....Yes
DESCRIPTION.....0 = <90% with O2
DISPLAY_NAME.....0 = <90% with O2
HELP_TEXT.....0 = <90% with O2
ID.....{30042730-F600-44FD-A3D0-53172B70082B}
TERM.....0 = <90%

Page: 645

DATA VALUE: 0 = ABSENT PULSE

ABBREVIATION.....0p
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....0 = Absent Pulse
HELP_TEXT.....No Help Text
ID.....{4253CB74-51AE-BBFA-214F-937762F7638C}
TERM.....0 = ABSENT PULSE

DATA VALUE: 0 = APNEIC

ABBREVIATION.....ALD RESP 0
ACTIVE.....Yes
DESCRIPTION.....0 = Apneic
DISPLAY_NAME.....0 = Apneic
HELP_TEXT.....0 = Apneic
ID.....{649C42B0-6E0A-4E86-AAC5-70A8C20EC511}
TERM.....0 = APNEIC

DATA VALUE: 0 = BP +/- 50% PREOP

ABBREVIATION.....ALD CIR 0
ACTIVE.....Yes
DESCRIPTION.....0 = BP +/- 50% Preop
DISPLAY_NAME.....0 = BP +/- 50% Preop
HELP_TEXT.....0 = BP +/- 50% Preop

ID.....{675F2CD1-DD7B-473B-B8EC-C2BFBF4AC3B2}
TERM.....0 = BP +/- 50% PREOP

DATA VALUE: 0 = NOT RESPONDING

ABBREVIATION.....ALD CON 0
ACTIVE.....Yes
DESCRIPTION.....0 = Not Responding
DISPLAY_NAME.....0 = Not Responding
HELP_TEXT.....0 = Not Responding
ID.....{E8765215-B0B0-430C-8C89-AB969A6AD12E}
TERM.....0 = NOT RESPONDING

DATA VALUE: 1 - BARELY PERCEPTIBLE SWEATING / PALMS MOIST

ABBREVIATION.....CIWA PS 1
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....1 - Barely Perceptible Sweating /
Palms Moist
HELP_TEXT.....No Help Text
ID.....{93B5B833-14B5-4ADD-97E8-7BF14A8BA01A}
TERM.....1 - BARELY PERCEPTIBLE SWEATING /
PALMS MOIST

Page: 646

DATA VALUE: 1 - CANNOT DO SERIAL ADDITIONS

ABBREVIATION.....CIWA OR 1
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....1 - Cannot Do Serial Additions Or Is
Uncertain About Date
HELP_TEXT.....No Help Text
ID.....{F4446E01-E326-420E-8EA8-04451AA9D963}
TERM.....1 - CANNOT DO SERIAL ADDITIONS

DATA VALUE: 1 - MILD NAUSEA / NO VOMITING

ABBREVIATION.....CIWA NV 1
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....1 - Mild Nausea / No Vomiting
HELP_TEXT.....No Help Text
ID.....{582CBFB6-60A0-4F6D-936A-9D2D6F871B4A}
TERM.....1 - MILD NAUSEA / NO VOMITING

DATA VALUE: 1 - MILDLY ANXIOUS

ABBREVIATION.....CIWA AX 1
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....1 - Mildly Anxious
HELP_TEXT.....No Help Text

ID.....{401113EE-5CE7-4AE9-80C1-1C3252E706D0}
TERM.....1 - MILDLY ANXIOUS

DATA VALUE: 1 - NO RESPONSE

ABBREVIATION.....1-nr
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....1 - No Response
HELP_TEXT.....No Help Text
ID.....{26D3168D-40EA-C918-D321-54BED4837311}
TERM.....1 - NO RESPONSE

DATA VALUE: 1 - NOT VISIBLE BUT CAN BE FELT FINGERTIP

ABBREVIATION.....CIWA TR 1
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....1 - Not Visible But Can Be Felt
Fingertip To Fingertip
HELP_TEXT.....No Help Text
ID.....{31833D55-9018-47EE-8341-A0863744CEF8}
TERM.....1 - NOT VISIBLE BUT CAN BE FELT
FINGERTIP

Page: 647

DATA VALUE: 1 - VERY MILD

ABBREVIATION.....CIWA HA 1
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....1 - Very Mild
HELP_TEXT.....No Help Text
ID.....{F789C1EA-CC7A-4368-87BA-08AFFE492EE6}
TERM.....1 - VERY MILD

DATA VALUE: 1 - VERY MILD HARSHNESS OR ABILITY TO FRIGHTEN

ABBREVIATION.....CIWA AU 1
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....1 - Very Mild Harshness Or Ability To
Frighten
HELP_TEXT.....No Help Text
ID.....{E7D03A9C-B975-48B2-8CA6-2FC6616C6AF6}
TERM.....1 - VERY MILD HARSHNESS OR ABILITY TO
FRIGHTEN

DATA VALUE: 1 - VERY MILD ITCHING PINS/NEEDLES BURNING

ABBREVIATION.....CIWA TC 1
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....1 - Very Mild Itching, Pins/Needles,

Burning, Numbness
HELP_TEXT.....No Help Text
ID.....{00642E0D-E1CF-4EE0-91F4-CBAF209AA54C}
TERM.....1 - VERY MILD ITCHING PINS/NEEDLES
BURNING

DATA VALUE: 1 - VERY MILD SENSITIVITY

ABBREVIATION.....CIWA VS 1
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....1 - Very Mild Sensitivity
HELP_TEXT.....No Help Text
ID.....{8BF96CF1-3EA7-468C-838A-284430EA850D}
TERM.....1 - VERY MILD SENSITIVITY

DATA VALUE: 1 - WEAK

ABBREVIATION.....1 - WEAK
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....1 - Weak
HELP_TEXT.....No Help Text
ID.....{A41570DC-2F35-DDF6-977B-7C262150DFA2}
TERM.....1 - WEAK

Page: 648

DATA VALUE: 1 = 2 EXTREMITIES

ABBREVIATION.....ALD ACT 1
ACTIVE.....Yes
DESCRIPTION.....1 = 2 Extremities
DISPLAY_NAME.....1 = Moves 2 Extremities
HELP_TEXT.....1 = 2 Extremities
ID.....{6DCB4762-72CE-48C2-9AFE-B19969ED1DE0}
TERM.....1 = 2 EXTREMITIES

DATA VALUE: 1 = AROUSABLE

ABBREVIATION.....ALD CON 1
ACTIVE.....Yes
DESCRIPTION.....1 = Arousable with minimal stimulation
DISPLAY_NAME.....1 = Arousable with minimal stimulation
HELP_TEXT.....1 = Arousable with minimal stimulation
ID.....{B2A07804-9829-42D0-9C49-A61FD51EEAA0}
TERM.....1 = AROUSABLE

DATA VALUE: 1 = BP +/- 20-49% PREOP

ABBREVIATION.....ALD CIR 1
ACTIVE.....Yes
DESCRIPTION.....1 = BP +/- 20-49% Preop
DISPLAY_NAME.....1 = BP +/- 20-49% Preop
HELP_TEXT.....1 = BP +/- 20-49% Preop

ID.....{D4A6100D-1228-4860-9BD6-B43F30D4FF94}
TERM.....1 = BP +/- 20-49% PREOP

DATA VALUE: 1 = DYSPNEA

ABBREVIATION.....ALD RESP 1
ACTIVE.....Yes
DESCRIPTION.....1 = Dyspnea, shallow or limited
DISPLAY_NAME.....1 = Dyspnea, shallow or limited
HELP_TEXT.....1 = Dyspnea, shallow or limited
ID.....{3DAE60A2-57DA-4121-A507-BC43497C7B67}
TERM.....1 = DYSPNEA

DATA VALUE: 1 = NEEDS O2

ABBREVIATION.....ALD O2 1
ACTIVE.....Yes
DESCRIPTION.....1 = Needs O2 to maintain >90%
DISPLAY_NAME.....1 = Needs O2 to maintain >90%
HELP_TEXT.....1 = Needs O2 to maintain >90%
ID.....{5D97F6BB-7B02-4644-8CE5-014C89826731}
TERM.....1 = NEEDS O2

Page: 649

DATA VALUE: 1 HOUR

ABBREVIATION.....1 HOUR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....1 Hour
HELP_TEXT.....No Help Text
ID.....{319BA174-1624-483C-B666-4998127B0800}
TERM.....1 HOUR

DATA VALUE: 1+=PALPABLE, THREADY, EASILY OBLITERATED

ABBREVIATION.....1teo
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....1+=Palpable, Thready, Easily
Obliterated
HELP_TEXT.....No Help Text
ID.....{3AE92556-1337-BFE2-C3EE-9676327935C3}
TERM.....1+=PALPABLE, THREADY, EASILY
OBLITERATED

DATA VALUE: 1- MILD ANXIOUS

ABBREVIATION.....CIWA AX 1
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....1 - Mild Anxious
HELP_TEXT.....No Help Text
ID.....{CAB6FB3A-4088-4D2C-AE1E-623352511D73}

TERM.....1- MILD ANXIOUS

DATA VALUE: 1- SOMEWHAT MORE THAN NORMAL ACTIVITY

ABBREVIATION.....CIWA AG 1
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....1 - Somewhat More Than Normal Activity
HELP_TEXT.....No Help Text
ID.....{497F7A3D-72C8-4986-A1E6-74EA0E5704F9}
TERM.....1- SOMEWHAT MORE THAN NORMAL ACTIVITY

DATA VALUE: 1800 CALORIE CONTROLLED

ABBREVIATION.....1800 CAL C
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....1800 Calorie Controlled Diet
HELP_TEXT.....No Help Text
ID.....{C564D7FE-EC2C-401E-A77D-79CC5FC2B2F8}
TERM.....1800 CALORIE CONTROLLED

Page: 650

DATA VALUE: 1800 DIABETIC

ABBREVIATION.....1800 DIAB
ACTIVE.....Yes
DESCRIPTION.....1800 Calorie Diabetic Diet
DISPLAY_NAME.....1800 Diabetic
HELP_TEXT.....No Help Text
ID.....{FC71C7D7-F24E-4316-A455-8674804BE092}
TERM.....1800 DIABETIC

DATA VALUE: 1:1

ABBREVIATION.....1:1
ACTIVE.....Yes
DESCRIPTION.....1 to 1 ratio
DISPLAY_NAME.....1:1
HELP_TEXT.....No Help Text
ID.....{5C07BBD0-FDB9-0CD7-F2C2-616913E0281B}
TERM.....1:1

DATA VALUE: 1:1 OBSERVATION

ABBREVIATION.....1o
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....1:1 Observation
HELP_TEXT.....No Help Text
ID.....{B1A7AF4B-9BD2-5F48-79B3-F35861AA5C45}
TERM.....1:1 OBSERVATION

DATA VALUE: 1:2

ABBREVIATION.....1:2
ACTIVE.....Yes
DESCRIPTION.....1 to 2 ratio
DISPLAY_NAME.....1:2
HELP_TEXT.....No Help Text
ID.....{287CDCE0-36F3-65F6-E6A5-9EF7CF8976AE}
TERM.....1:2

DATA VALUE: 1:3

ABBREVIATION.....1:3
ACTIVE.....Yes
DESCRIPTION.....1 to 3 ratio
DISPLAY_NAME.....1:3
HELP_TEXT.....No Help Text
ID.....{062B4F93-8F48-7373-F679-11F57CE98BB3}
TERM.....1:3

Page: 651

DATA VALUE: 1:4

ABBREVIATION.....1:4
ACTIVE.....Yes
DESCRIPTION.....1 to 4 ratio
DISPLAY_NAME.....1:4
HELP_TEXT.....No Help Text
ID.....{3F1B75FF-C168-D413-0BDF-6581E7A4EEAD}
TERM.....1:4

DATA VALUE: 1:8

ABBREVIATION.....1:8
ACTIVE.....Yes
DESCRIPTION.....1 to 8 ratio
DISPLAY_NAME.....1:8
HELP_TEXT.....No Help Text
ID.....{94790C4C-3AF2-6F4C-58B8-87021CA97A53}
TERM.....1:8

DATA VALUE: 1A + 2V

ABBREVIATION.....1+2
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....1A + 2V
HELP_TEXT.....No Help Text
ID.....{FE47CEC3-C3A5-0F5C-3EE8-50F30FD4E0C4}
TERM.....1A + 2V

DATA VALUE: 1A + 3V

ABBREVIATION.....1+3
ACTIVE.....Yes

DESCRIPTION.....No Description
DISPLAY_NAME.....1A + 3V
HELP_TEXT.....No Help Text
ID.....{CEDF220D-E593-B8F9-DB08-D32800AF3F5E}
TERM.....1A + 3V

DATA VALUE: 1A +1V

ABBREVIATION.....1+
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....1A +1V
HELP_TEXT.....No Help Text
ID.....{31711946-9F1B-0AD9-BDBE-846EDBC86AE2}
TERM.....1A +1V

Page: 652

DATA VALUE: 2 - ABNORMAL EXTENSION FROM PAIN

ABBREVIATION.....2-aeftp
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....2 - Abnormal Extension From Pain
HELP_TEXT.....No Help Text
ID.....{9B4951D8-B271-DFC3-0A00-1DC9984039E3}
TERM.....2 - ABNORMAL EXTENSION FROM PAIN

DATA VALUE: 2 - BEADS OF SWEAT OBVIOUS ON FOREHEAD

ABBREVIATION.....CIWA PS 2
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....2 - Beads Of Sweat Obvious On Forehead
HELP_TEXT.....No Help Text
ID.....{38EA4194-F5F4-42CB-8DE2-715167888F86}
TERM.....2 - BEADS OF SWEAT OBVIOUS ON FOREHEAD

DATA VALUE: 2 - DISORIENTED FOR DATE < 2 CALENDAR DAYS

ABBREVIATION.....CIWA OR 2
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....2 - Disoriented For Date By No More
Than 2 Calendar Days
HELP_TEXT.....No Help Text
ID.....{01B894B0-3E70-47FA-820A-83CA48B93EF9}
TERM.....2 - DISORIENTED FOR DATE < 2 CALENDAR
DAYS

DATA VALUE: 2 - INCOMPREHENSIBLE SOUNDS

ABBREVIATION.....2-is
ACTIVE.....Yes
DESCRIPTION.....No Description

DISPLAY_NAME.....2 - Incomprehensible Sounds
HELP_TEXT.....No Help Text
ID.....{F2261C3D-79A3-055B-702D-446B7A6A3113}
TERM.....2 - INCOMPREHENSIBLE SOUNDS

DATA VALUE: 2 - INTERMITTENT NAUSEA WITH DRY HEAVES

ABBREVIATION.....CIWA NV 2
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....2 - Intermittent Nausea With Dry Heaves
HELP_TEXT.....No Help Text
ID.....{E92221EE-1298-4D9D-9EFD-35F615118C96}
TERM.....2 - INTERMITTENT NAUSEA WITH DRY HEAVES

Page: 653

DATA VALUE: 2 - MILD

ABBREVIATION.....CIWA HA 2
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....2 - Mild
HELP_TEXT.....No Help Text
ID.....{6B752397-0A4F-4176-90E7-E220DB2BF1C4}
TERM.....2 - MILD

DATA VALUE: 2 - MILD HARSHNESS OR ABILITY TO FRIGHTEN

ABBREVIATION.....CIWA AU 2
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....2 - Mild Harshness Or Ability To
Frighten
HELP_TEXT.....No Help-Text
ID.....{71204969-E614-46A4-ADA5-39FAF81A79B5}
TERM.....2 - MILD HARSHNESS OR ABILITY TO
FRIGHTEN

DATA VALUE: 2 - MILD ITCHING PINS/NEEDLES BURNING NUMBNESS

ABBREVIATION.....CIWA TC 2
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....2 - Mild Itching, Pins/Needles,
Burning, Numbness
HELP_TEXT.....No Help Text
ID.....{1A591F17-CFB3-4B24-8AD6-ECF09318E2B6}
TERM.....2 - MILD ITCHING PINS/NEEDLES BURNING
NUMBNESS

DATA VALUE: 2 - MILD SENSITIVITY

ABBREVIATION.....CIWA VS 2
ACTIVE.....Yes

DESCRIPTION.....No Description
DISPLAY_NAME.....2 - Mild Sensitivity
HELP_TEXT.....No Help Text
ID.....{ED600249-E5E8-4020-955E-A5CDA18B5DE9}
TERM.....2 - MILD SENSITIVITY

DATA VALUE: 2 - MODERATE WITH PATIENTS ARMS EXTENDED

ABBREVIATION.....CIWA TR 2
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....2 - Moderate With Patient's Arms
Extended
HELP_TEXT.....No Help Text
ID.....{B4E7B7FB-B88E-435A-904E-C613AC937323}
TERM.....2 - MODERATE WITH PATIENTS ARMS
EXTENDED

Page: 654

DATA VALUE: 2 - NORMAL

ABBREVIATION.....2-n
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....2 - Normal
HELP_TEXT.....No Help Text
ID.....{25FAADFC-D280-78AC-4F01-CFCF34ECD3FF}
TERM.....2 - NORMAL

DATA VALUE: 2 - TO PAIN

ABBREVIATION.....2-tp
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....2 - To Pain
HELP_TEXT.....No Help Text
ID.....{49CE3C40-5FC2-002E-1C89-3E5DFA4CCF0F}
TERM.....2 - TO PAIN

DATA VALUE: 2 = 4 EXTREMITIES

ABBREVIATION.....ALD ACT 2
ACTIVE.....Yes
DESCRIPTION.....2 = 4 Extremities
DISPLAY_NAME.....2 = Moves 4 Extremities
HELP_TEXT.....2 = 4 Extremities
ID.....{3434DF4B-5CB6-417F-B5BE-CEF2D63FE8E0}
TERM.....2 = 4 EXTREMITIES

DATA VALUE: 2 = BP +/- 20% PREOP

ABBREVIATION.....ALD CIR 2
ACTIVE.....Yes
DESCRIPTION.....2 = BP +/- 20% Preop

DISPLAY_NAME.....2 = BP +/- 20% Preop
HELP_TEXT.....2 = BP +/- 20% Preop
ID.....{04B67413-7A88-409B-8005-471A7C84B879}
TERM.....2 = BP +/- 20% PREOP

DATA VALUE: 2 = BREATHES FREELY

ABBREVIATION.....ALD RESP 2
ACTIVE.....Yes
DESCRIPTION.....2 = Breathes freely, coughs freely
DISPLAY_NAME.....2 = Breathes freely, coughs freely
HELP_TEXT.....2 = Breathes freely, coughs freely
ID.....{278CD44C-538E-4015-9D98-D9312A7039BA}
TERM.....2 = BREATHES FREELY

Page: 655

DATA VALUE: 2 = FULLY AWAKE

ABBREVIATION.....ALD CON 2
ACTIVE.....Yes
DESCRIPTION.....2 = Fully awake
DISPLAY_NAME.....2 = Fully awake
HELP_TEXT.....2 = Fully awake
ID.....{903D0F56-ADEC-4527-BA66-C8A568E3EAC0}
TERM.....2 = FULLY AWAKE

DATA VALUE: 2 = ON ROOM AIR

ABBREVIATION.....ALD O2 2
ACTIVE.....Yes
DESCRIPTION.....2 = On room air
DISPLAY_NAME.....2 = >92% on room air
HELP_TEXT.....2 = On room air
ID.....{D58C5B46-CCD3-4E6D-86E8-7A1CFA92BC1C}
TERM.....2 = ON ROOM AIR

DATA VALUE: 2+=NORMAL, NOT EASILY OBLITERATED

ABBREVIATION.....2neo
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....2+=Normal, Not Easily Obliterated
HELP_TEXT.....No Help Text
ID.....{50484DA4-EC15-47E5-3923-268C4B5D6428}
TERM.....2+=NORMAL, NOT EASILY OBLITERATED

DATA VALUE: 2- MILD HARSHNESS OR ABILITY TO FRIGHTEN

ABBREVIATION.....CIWA AU 2
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....2 - Mild Harshness Or Ability To
Frighten
HELP_TEXT.....No Help Text

ID.....{42DAFB3B-E9AF-4558-95D4-2FFF7BE779F1}
TERM.....2- MILD HARSHNESS OR ABILITY TO
FRIGHTEN

DATA VALUE: 2.0 CAL FORMULA

ABBREVIATION.....2.0 C FORM
ACTIVE.....Yes
DESCRIPTION.....2.0 Cal Enteral Nutrition Formula
DISPLAY_NAME.....2.0 Cal Formula
HELP_TEXT.....No Help Text
ID.....{4A595230-45D1-40F0-AEF3-11E75DF418CB}
TERM.....2.0 CAL FORMULA

Page: 656

DATA VALUE: 2000 CALORIE CONTROLLED

ABBREVIATION.....2000 CAL C
ACTIVE.....Yes
DESCRIPTION.....2000 Calorie Controlled Diet
DISPLAY_NAME.....2000 Calorie Controlled Diet
HELP_TEXT.....No Help Text
ID.....{DB2FF9DC-298E-4881-94D4-47A7AF0E2BD7}
TERM.....2000 CALORIE CONTROLLED

DATA VALUE: 2000 DIABETIC

ABBREVIATION.....2000 DIAB
ACTIVE.....Yes
DESCRIPTION.....2000 Diabetic Diet
DISPLAY_NAME.....2000 Diabetic Diet
HELP_TEXT.....No Help Text
ID.....{DF3056C9-3AEA-4F4E-B758-E35590080AFB}
TERM.....2000 DIABETIC

DATA VALUE: 2200 CALORIE CONTROLLED

ABBREVIATION.....2200 CAL C
ACTIVE.....Yes
DESCRIPTION.....2200 Calorie Controlled Diet
DISPLAY_NAME.....2200 Calorie Controlled Diet
HELP_TEXT.....No Help Text
ID.....{77F7C138-B2CB-4D76-8C0B-D68F317D4DE0}
TERM.....2200 CALORIE CONTROLLED

DATA VALUE: 2200 DIABETIC

ABBREVIATION.....2200 DIAB
ACTIVE.....Yes
DESCRIPTION.....2200 Diabetic Diet
DISPLAY_NAME.....2200 Diabetic Diet
HELP_TEXT.....No Help Text
ID.....{1FB3EFB6-B92F-44FB-831E-4641C91C7A8D}
TERM.....2200 DIABETIC

DATA VALUE: 2400 CALORIE CONTROLLED

ABBREVIATION.....2400 CAL C
ACTIVE.....Yes
DESCRIPTION.....2400 Calorie Controlled Diet
DISPLAY_NAME.....2400 Calorie Controlled Diet
HELP_TEXT.....No Help Text
ID.....{F3EA28B5-FC96-4F34-B2BE-A96A196B38A2}
TERM.....2400 CALORIE CONTROLLED

Page: 657

DATA VALUE: 2400 DIABETIC

ABBREVIATION.....2400 DIAB
ACTIVE.....Yes
DESCRIPTION.....2400 Diabetic Diet
DISPLAY_NAME.....2400 Diabetic Diet
HELP_TEXT.....No Help Text
ID.....{2018FB65-3D4E-47E0-90E9-205A45EDF130}
TERM.....2400 DIABETIC

DATA VALUE: 25 PERCENT ALBUMIN

ABBREVIATION.....25 PERC AL
ACTIVE.....Yes
DESCRIPTION.....25 Percent Albumin
DISPLAY_NAME.....25 Percent Albumin
HELP_TEXT.....No Help Text
ID.....{2FCE19CA-B5AB-4F15-8FEA-CF156D1E4CBC}
TERM.....25 PERCENT ALBUMIN

DATA VALUE: 2:4

ABBREVIATION.....2:4
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....2:4
HELP_TEXT.....No Help Text
ID.....{7D423C9A-42B3-4FDD-999B-889193492F5C}
TERM.....2:4

DATA VALUE: 2A

ABBREVIATION.....2
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....2A
HELP_TEXT.....No Help Text
ID.....{C75BDC9F-C756-9051-CAA3-A27DDD94B35E}
TERM.....2A

DATA VALUE: 2A + 1V

ABBREVIATION.....2+1
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....2A + 1V
HELP_TEXT.....No Help Text
ID.....{1849AA0D-1572-55E0-A780-ACE602E2C72D}
TERM.....2A + 1V

Page: 658

DATA VALUE: 2A + 2V

ABBREVIATION.....2+2
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....2A + 2V
HELP_TEXT.....No Help Text
ID.....{17203251-68E8-A00A-CF54-E2CC55B8EEDD}
TERM.....2A + 2V

DATA VALUE: 2V

ABBREVIATION.....2
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....2V
HELP_TEXT.....No Help Text
ID.....{D24B5FA4-1DA6-36DE-49DB-627A9B82E968}
TERM.....2V

DATA VALUE: 3 - ABNORMAL FLEXION FROM PAIN

ABBREVIATION.....3-affp
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....3 - Abnormal Flexion From Pain
HELP_TEXT.....No Help Text
ID.....{9AED8954-5F5D-D2C2-3CAF-4F8EF92009C2}
TERM.....3 - ABNORMAL FLEXION FROM PAIN

DATA VALUE: 3 - BEADS OF SWEAT OBVIOUS ON FOREHEAD

ABBREVIATION.....CIWA PS 3
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....3 - Beads Of Sweat Obvious On Forehead
HELP_TEXT.....No Help Text
ID.....{663591B6-3FB4-4060-8A49-18747F857B16}
TERM.....3 - BEADS OF SWEAT OBVIOUS ON FOREHEAD

DATA VALUE: 3 - DISORIENTED FOR DATE > 2 CALENDAR DAYS

ABBREVIATION.....CIWA OR 3
ACTIVE.....Yes
DESCRIPTION.....No Description

DISPLAY_NAME.....3 - Disoriented For Date By More Than
2 Calendar Days
HELP_TEXT.....No Help Text
ID.....{68C351F9-C8DE-49C5-B3BD-D69E97F3E516}
TERM.....3 - DISORIENTED FOR DATE > 2 CALENDAR
DAYS

Page: 659

DATA VALUE: 3 - INAPPROPRIATE WORDS

ABBREVIATION.....3-iw
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....3 - Inappropriate Words
HELP_TEXT.....No Help Text
ID.....{FE00191-243E-4D32-1E1E-546EE2CF88A5}
TERM.....3 - INAPPROPRIATE WORDS

DATA VALUE: 3 - INTERMITTENT NAUSEA WITH DRY HEAVES

ABBREVIATION.....CIWA NV 3
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....3 - Intermittent Nausea With Dry Heaves
HELP_TEXT.....No Help Text
ID.....{D299B4E1-0FF4-474D-BECD-5537391EEE7A}
TERM.....3 - INTERMITTENT NAUSEA WITH DRY HEAVES

DATA VALUE: 3 - MODERATE

ABBREVIATION.....CIWA HA 3
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....3 - Moderate
HELP_TEXT.....No Help Text
ID.....{F26EB828-6864-4384-9EBB-35A6CDCACC59}
TERM.....3 - MODERATE

DATA VALUE: 3 - MODERATE HARSHNESS OR ABILITY TO FRIGHTEN

ABBREVIATION.....CIWA AU 3
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....3 - Moderate Harshness Or Ability To
Frighten
HELP_TEXT.....No Help Text
ID.....{77F59277-6CE4-4323-8531-3CA26E13C9C4}
TERM.....3 - MODERATE HARSHNESS OR ABILITY TO
FRIGHTEN

DATA VALUE: 3 - MODERATE ITCHING PINS/NEEDLES BURNING

ABBREVIATION.....CIWA TC 3
ACTIVE.....Yes

DESCRIPTION.....No Description
DISPLAY_NAME.....3 - Moderate Itching, Pins/Needles,
Burning, Numbness
HELP_TEXT.....No Help Text
ID.....{B7E2BFDD-7D85-4BCB-BE1D-0A220C1CEDAE}
TERM.....3 - MODERATE ITCHING PINS/NEEDLES
BURNING

Page: 660

DATA VALUE: 3 - MODERATE SENSITIVITY

ABBREVIATION.....CIWA VS 3
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....3 - Moderate Sensitivity
HELP_TEXT.....No Help Text
ID.....{D673A630-1F1E-4734-AADE-17586C48E8B9}
TERM.....3 - MODERATE SENSITIVITY

DATA VALUE: 3 - MODERATE WITH PATIENTS ARMS EXTENDED

ABBREVIATION.....CIWA TR 3
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....3 - Moderate With Patient's Arms
Extended
HELP_TEXT.....No Help Text
ID.....{36FF0AB3-F5ED-4CBA-A6F3-110DFCB65244}
TERM.....3 - MODERATE WITH PATIENTS ARMS
EXTENDED

DATA VALUE: 3 - STRONG

ABBREVIATION.....3-s
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....3 - Strong
HELP_TEXT.....No Help Text
ID.....{C55EAB72-7816-5440-84E9-6E7E190ED197}
TERM.....3 - STRONG

DATA VALUE: 3 - TO VERBAL COMMAND

ABBREVIATION.....3-tvc
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....3 - To Verbal Command
HELP_TEXT.....No Help Text
ID.....{EBF4B748-2861-C81F-C0D0-B9FE99FC7A8F}
TERM.....3 - TO VERBAL COMMAND

DATA VALUE: 3 WAY FOLEY

ABBREVIATION.....3WAY FOLEY

ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....3 Way Foley
HELP_TEXT.....No Help Text
ID.....{32408760-F517-4530-AEC0-54FE0A6A062B}
TERM.....3 WAY FOLEY

Page: 661

DATA VALUE: 3+=BOUNDING, EASILY PALPABLE, CANNOT BE OBLITERATED

ABBREVIATION.....3epcbo
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....3+=Bounding, Easily Palpable, Cannot
Be Obliterated
HELP_TEXT.....No Help Text
ID.....{7A7A07ED-26B6-A3BF-C564-94BD4F6D4E1C}
TERM.....3+=BOUNDING, EASILY PALPABLE, CANNOT
BE OBLITERATED

DATA VALUE: 34 ml

ABBREVIATION.....34 ml
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....34 ml
HELP_TEXT.....No Help Text
ID.....{25952BB5-50B1-4CC0-9483-8AF891767045}
TERM.....34 ml

DATA VALUE: 3:4

ABBREVIATION.....3:4
ACTIVE.....Yes
DESCRIPTION.....3:4
DISPLAY_NAME.....3:4
HELP_TEXT.....No Help Text
ID.....{E4EBD0E2-FE8C-40FA-9DCA-475DCE43B1A2}
TERM.....3:4

DATA VALUE: 4 - BEADS OF SWEAT OBVIOUS ON FOREHEAD

ABBREVIATION.....CIWA PS 4
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....4 - Beads Of Sweat Obvious On Forehead
HELP_TEXT.....No Help Text
ID.....{8388B1A0-46EA-4003-A32B-A2C708A6F056}
TERM.....4 - BEADS OF SWEAT OBVIOUS ON FOREHEAD

DATA VALUE: 4 - DISORIENTED AND CONVERSES

ABBREVIATION.....4-dac
ACTIVE.....Yes

DESCRIPTION.....No Description
DISPLAY_NAME.....4 - Disoriented And Converses
HELP_TEXT.....No Help Text
ID.....{D43F646D-A130-5D06-6DC6-DDB971027231}
TERM.....4 - DISORIENTED AND CONVERSES

Page: 662

DATA VALUE: 4 - DISORIENTED FOR PLACE OR PERSON

ABBREVIATION.....CIWA OR 4
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....4 - Disoriented For Place/Or Person
HELP_TEXT.....No Help Text
ID.....{03BBFFA2-233A-41A1-8212-F988889C5E72}
TERM.....4 - DISORIENTED FOR PLACE OR PERSON

DATA VALUE: 4 - INTERMITTENT NAUSEA WITH DRY HEAVES

ABBREVIATION.....CIWA NV 4
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....4 - Intermittent Nausea With Dry Heaves
HELP_TEXT.....No Help Text
ID.....{7FBE10F0-472F-4F21-BE55-0892259141B5}
TERM.....4 - INTERMITTENT NAUSEA WITH DRY HEAVES

DATA VALUE: 4 - MODERATE WITH PATIENTS ARMS EXTENDED

ABBREVIATION.....CIWA TR 4
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....4 - Moderate With Patient's Arms
Extended
HELP_TEXT.....No Help Text
ID.....{79CE1486-151C-4492-AFD2-7845DD69DA9E}
TERM.....4 - MODERATE WITH PATIENTS ARMS
EXTENDED

DATA VALUE: 4 - MODERATELY ANXIOUS

ABBREVIATION.....CIWA AX 4
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....4 - Moderately Anxious, Or Guarded,
Anxiety Is Inferred
HELP_TEXT.....No Help Text
ID.....{9A01FEFD-844C-4FAB-BEC3-96333495FA62}
TERM.....4 - MODERATELY ANXIOUS

DATA VALUE: 4 - MODERATELY FIDGETY AND RESTLESS

ABBREVIATION.....CIWA AG 4
ACTIVE.....Yes

DESCRIPTION.....No Description
DISPLAY_NAME.....4 - Moderately Fidgety And Restless
HELP_TEXT.....No Help Text
ID.....{B1515E1F-4305-40D5-84B4-5E37AEF0AD3D}
TERM.....4 - MODERATELY FIDGETY AND RESTLESS

Page: 663

DATA VALUE: 4 - MODERATELY SEVERE

ABBREVIATION.....CIWA HA 4
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....4 - Moderately Severe
HELP_TEXT.....No Help Text
ID.....{83B3DE0F-3BA4-4B0A-A222-F372BD403094}
TERM.....4 - MODERATELY SEVERE

DATA VALUE: 4 - MODERATELY SEVERE HALLUCINATIONS

ABBREVIATION.....CIWA 4
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....4 - Moderately Severe Hallucinations
HELP_TEXT.....No Help Text
ID.....{13FCDA97-A967-4993-AED5-13EBA61F3102}
TERM.....4 - MODERATELY SEVERE HALLUCINATIONS

DATA VALUE: 4 - SPONTANEOUSLY

ABBREVIATION.....4-s
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....4 - Spontaneously
HELP_TEXT.....No Help Text
ID.....{4612FC84-FCA8-8C85-50DC-A9A42C195F51}
TERM.....4 - SPONTANEOUSLY

DATA VALUE: 4 - WITHDRAWS FROM PAIN

ABBREVIATION.....4-wfp
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....4 - Withdraws From Pain
HELP_TEXT.....No Help Text
ID.....{D808ACF1-30CB-29AB-4957-E91D63C4DF18}
TERM.....4 - WITHDRAWS FROM PAIN

DATA VALUE: 4 HOUR

ABBREVIATION.....4 HOUR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....4 Hour
HELP_TEXT.....No Help Text

ID.....{91B2A49D-0308-47C9-B220-63A0E318A201}
TERM.....4 HOUR

Page: 664

DATA VALUE: 40 ml

ABBREVIATION.....40 ml
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....40 ml
HELP_TEXT.....No Help Text
ID.....{37A2998F-AB1C-4493-9C7A-FFCB05ED2302}
TERM.....40 ml

DATA VALUE: 4:4

ABBREVIATION.....4:4
ACTIVE.....Yes
DESCRIPTION.....4:4
DISPLAY_NAME.....4:4
HELP_TEXT.....No Help Text
ID.....{E9FCD5A0-995A-401D-8E85-6F0C753609B6}
TERM.....4:4

DATA VALUE: 5 - CONSTANT NAUSEA W FREQUENT HEAVES AND VOMITING

ABBREVIATION.....CIWA NV 5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....5 - Constant Nausea W Frequent Heaves
and Vomiting
HELP_TEXT.....No Help Text
ID.....{3A90064F-0B58-48AB-ABEE-188F215F4E33}
TERM.....5 - CONSTANT NAUSEA W FREQUENT HEAVES
AND VOMITING

DATA VALUE: 5 - DRENCHING SWEAT

ABBREVIATION.....CIWA PS 5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....5 - Drenching Sweat
HELP_TEXT.....No Help Text
ID.....{4CAA2CFF-6DA5-4238-A25B-4D3A960808A9}
TERM.....5 - DRENCHING SWEAT

DATA VALUE: 5 - LOCALIZES PAIN

ABBREVIATION.....5-lp
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....5 - Localizes Pain
HELP_TEXT.....No Help Text
ID.....{CF13650F-899B-6AE1-5913-DC1F4776E5E4}

TERM.....5 - LOCALIZES PAIN

Page: 665

DATA VALUE: 5 - ORIENTED AND CONVERSES

ABBREVIATION.....5-oac
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....5 - Oriented And Converses
HELP_TEXT.....No Help Text
ID.....{53AD29B1-BA51-4BDE-4F4C-4268619C93A4}
TERM.....5 - ORIENTED AND CONVERSES

DATA VALUE: 5 - SEVERE

ABBREVIATION.....CIWA HA 5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....5 - Severe
HELP_TEXT.....No Help Text
ID.....{7BD7B91F-DF6F-47B9-BAA6-3EF249340578}
TERM.....5 - SEVERE

DATA VALUE: 5 - SEVERE EVEN WITH ARMS NOT EXTENDED

ABBREVIATION.....CIWA TR 5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....5 - Severe Even With Arms Not Extended
HELP_TEXT.....No Help Text
ID.....{611A344A-C540-4861-81D7-681DAF1AA473}
TERM.....5 - SEVERE EVEN WITH ARMS NOT EXTENDED

DATA VALUE: 5 - SEVERE HALLUCINATIONS

ABBREVIATION.....CIWA 5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....5 - Severe Hallucinations
HELP_TEXT.....No Help Text
ID.....{14DB90F7-0182-44AC-8B5B-CE86A0E5A9E6}
TERM.....5 - SEVERE HALLUCINATIONS

DATA VALUE: 5 PERCENT ALBUMIN

ABBREVIATION.....5 PERC ALB
ACTIVE.....Yes
DESCRIPTION.....5 Percent Albumin
DISPLAY_NAME.....5 Percent Albumin
HELP_TEXT.....No Help Text
ID.....{38CE47CB-736B-475E-8E89-C4C13F2DF184}
TERM.....5 PERCENT ALBUMIN

DATA VALUE: 6 - CONSTANT NAUSEA W FREQUENT HEAVES AND VOMITING

ABBREVIATION.....CIWA NV 6
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....6 - Constant Nausea W Frequent Heaves
and Vomiting
HELP_TEXT.....No Help Text
ID.....{446FE39A-EB53-4D6B-B5B8-1370A5D55EEB}
TERM.....6 - CONSTANT NAUSEA W FREQUENT HEAVES
AND VOMITING

DATA VALUE: 6 - DRENCHING SWEAT

ABBREVIATION.....CIWA PS 6
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....6 - Drenching Sweat
HELP_TEXT.....No Help Text
ID.....{7294B311-9032-4FE1-9779-363D19FD1FD3}
TERM.....6 - DRENCHING SWEAT

DATA VALUE: 6 - EXTREMELY SEVERE HALLUCINATIONS

ABBREVIATION.....CIWA 6
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....6 - Extremely Severe Hallucinations
HELP_TEXT.....No Help Text
ID.....{9C63ED29-2B47-4E59-A026-1FDB0666AFCB}
TERM.....6 - EXTREMELY SEVERE HALLUCINATIONS

DATA VALUE: 6 - OBEYS COMMANDS

ABBREVIATION.....6-oc
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....6 - Obeys Commands
HELP_TEXT.....No Help Text
ID.....{0D4956F4-7E79-8D17-EDF3-274B96BBC89A}
TERM.....6 - OBEYS COMMANDS

DATA VALUE: 6 - SEVERE EVEN WITH ARMS NOT EXTENDED

ABBREVIATION.....CIWA TR 6
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....6 - Severe Even With Arms Not Extended
HELP_TEXT.....No Help Text
ID.....{9BCF893C-76D0-4A3E-80D0-281A917F1A08}
TERM.....6 - SEVERE EVEN WITH ARMS NOT EXTENDED

DATA VALUE: 6 - VERY SEVERE

ABBREVIATION.....CIWA HA 6
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....6 - Very Severe
HELP_TEXT.....No Help Text
ID.....{241C1EF4-2E2A-4D65-9991-607BB6B29E01}
TERM.....6 - VERY SEVERE

DATA VALUE: 6 HOUR

ABBREVIATION.....6 HOUR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....6 Hour
HELP_TEXT.....No Help Text
ID.....{0D1F000A-01F0-422C-9CB8-2A04CD479444}
TERM.....6 HOUR

DATA VALUE: 6 PERCENT HETASTARCH

ABBREVIATION.....6 PERC HET
ACTIVE.....Yes
DESCRIPTION.....6 Percent Hetastarch
DISPLAY_NAME.....6 Percent Hetastarch
HELP_TEXT.....No Help Text
ID.....{D94072C9-87E8-46CA-929F-4904F3BD964C}
TERM.....6 PERCENT HETASTARCH

DATA VALUE: 7 - ACUTE PANIC

ABBREVIATION.....CIWA AX 7
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....7 - Equivalent To Acute Panic States
HELP_TEXT.....No Help Text
ID.....{3F5E6786-FD67-405F-A535-EC1609910F14}
TERM.....7 - ACUTE PANIC

DATA VALUE: 7 - CONSTANT NAUSEA W FREQUENT HEAVES AND VOMITING

ABBREVIATION.....CIWA NV 7
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....7 - Constant Nausea W Frequent Heaves
and Vomiting
HELP_TEXT.....No Help Text
ID.....{7F1076C4-C655-4BF2-8BB1-036F5169B619}
TERM.....7 - CONSTANT NAUSEA W FREQUENT HEAVES
AND VOMITING

DATA VALUE: 7 - CONTINUOUS HALLUCINATIONS

ABBREVIATION.....CIWA 7
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....7 - Continuous Hallucinations
HELP_TEXT.....No Help Text
ID.....{19ECC57C-26A6-4CAE-97C4-2A860F337BE2}
TERM.....7 - CONTINUOUS HALLUCINATIONS

DATA VALUE: 7 - DRENCHING SWEAT

ABBREVIATION.....CIWA PS 7
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....7 - Drenching Sweat
HELP_TEXT.....No Help Text
ID.....{15D2E2CC-BAED-403B-A141-B848FFE4A9D0}
TERM.....7 - DRENCHING SWEAT

DATA VALUE: 7 - EXTREMELY SEVERE

ABBREVIATION.....CIWA HA 7
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....7 - Extremely Severe
HELP_TEXT.....No Help Text
ID.....{F8DCC910-D36D-480E-867E-B45221170C49}
TERM.....7 - EXTREMELY SEVERE

DATA VALUE: 7 - PACING BACK AND FORTH MOST OF THE INTERVIEW

ABBREVIATION.....CIWA AG 7
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....7 - Pacing Back And Forth The
Interview/Thrashing
HELP_TEXT.....No Help Text
ID.....{9F03A167-CC40-4D65-9FD8-CFB3FEEAF833}
TERM.....7 - PACING BACK AND FORTH MOST OF THE
INTERVIEW

DATA VALUE: 7 - SEVERE EVEN WITH ARMS NOT EXTENDED

ABBREVIATION.....CIWA TR 7
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....7 - Severe Even With Arms Not Extended
HELP_TEXT.....No Help Text
ID.....{FE0BE72E-913A-48CD-BABF-C5B742AFED93}
TERM.....7 - SEVERE EVEN WITH ARMS NOT EXTENDED

DATA VALUE: 8 HOUR

ABBREVIATION.....8 HOUR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....8 Hour
HELP_TEXT.....No Help Text
ID.....{05E46D8F-C165-4552-8063-B735DECD4AF3}
TERM.....8 HOUR

DATA VALUE: AAI

ABBREVIATION.....AAI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....AAI
HELP_TEXT.....No Help Text
ID.....{DD18ACA-B98E-24D3-E760-6D082A76074B}
TERM.....AAI

DATA VALUE: AAIR

ABBREVIATION.....AAIR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....AAIR
HELP_TEXT.....No Help Text
ID.....{3B4B6F37-B631-BE32-F1C1-65B99D0C5018}
TERM.....AAIR

DATA VALUE: AAO

ABBREVIATION.....AAO
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....AAO
HELP_TEXT.....No Help Text
ID.....{E87E0670-2C0F-6E6A-DEC9-F52E5C530A2A}
TERM.....AAO

DATA VALUE: ABDOMINAL

ABBREVIATION.....ABDOMINAL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Abdominal
HELP_TEXT.....No Help Text
ID.....{29E279F4-5D2E-4644-A2FA-CF85BEDFD480}
TERM.....ABDOMINAL

Page: 670

DATA VALUE: ABLE TO AMBULATE

ABBREVIATION.....AMBULAT
ACTIVE.....Yes

DESCRIPTION.....No Description
DISPLAY_NAME.....Able to Ambulate
HELP_TEXT.....No Help Text
ID.....{4E1D40CA-F5C4-442E-8BBA-C6E4C4781F57}
TERM.....ABLE TO AMBULATE

DATA VALUE: ABRASION

ABBREVIATION.....ABRASION
ACTIVE.....Yes
DESCRIPTION.....Abrasion
DISPLAY_NAME.....Abrasion
HELP_TEXT.....No Help Text
ID.....{7E962B43-24B1-46D2-A79A-54811117271A}
TERM.....ABRASION

DATA VALUE: ABSENT

ABBREVIATION.....a
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Absent
HELP_TEXT.....No Help Text
ID.....{6DB90FB5-3CC2-B19D-1265-48EAC64F1BC0}
TERM.....ABSENT

DATA VALUE: ABSENT TO MINIMAL BLEEDING

ABBREVIATION.....ABS BLEE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Absent to Minimal Bleeding
HELP_TEXT.....No Help Text
ID.....{123460FB-1A8D-45AB-B0A2-489EAB482F67}
TERM.....ABSENT TO MINIMAL BLEEDING

DATA VALUE: ACCELERATING

ABBREVIATION.....ACCELERATI
ACTIVE.....Yes
DESCRIPTION.....Accelerating
DISPLAY_NAME.....Accelerating
HELP_TEXT.....Accelerating
ID.....{514FE300-87EE-4C2C-AC70-5E57DBECAFF7}
TERM.....ACCELERATING

DATA VALUE: ACHING

ABBREVIATION.....ACHING
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Aching
HELP_TEXT.....No Help Text

ID.....{736BC1A4-2DA7-4A56-ACC0-28086EA103E0}
TERM.....ACHING

DATA VALUE: ACTIVE RANGE OF MOTION

ABBREVIATION.....ACTIVE ROM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Active Range Of Motion
HELP_TEXT.....No Help Text
ID.....{F55E68F3-6788-47AE-B962-119413B4A089}
TERM.....ACTIVE RANGE OF MOTION

DATA VALUE: ACTIVITIES OF DAILY LIVING

ABBREVIATION.....ADL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Activities of Daily Living
HELP_TEXT.....No Help Text
ID.....{12F2CA66-3DCE-443D-883C-119CE4838B4A}
TERM.....ACTIVITIES OF DAILY LIVING

DATA VALUE: ADMISSION

ABBREVIATION.....ADMISSION
ACTIVE.....Yes
DESCRIPTION.....Admission
DISPLAY_NAME.....Admission
HELP_TEXT.....No Help Text
ID.....{483C2830-813A-46C3-A012-CCD70044B2A2}
TERM.....ADMISSION

DATA VALUE: ADMISSION DIAGNOSIS

ABBREVIATION.....ADM DIAG
ACTIVE.....Yes
DESCRIPTION.....Admission Diagnosis
DISPLAY_NAME.....Admission Diagnosis
HELP_TEXT.....Admission Diagnosis
ID.....{2EEC3F2F-9802-46C5-AB47-008AFA093A7F}
TERM.....ADMISSION DIAGNOSIS

Page: 672

DATA VALUE: ADULT

ABBREVIATION.....ADULT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Adult
HELP_TEXT.....No Help Text
ID.....{37B557BE-B932-4757-8E08-5EA315E8AD1F}
TERM.....ADULT

DATA VALUE: ADULT LARGE

ABBREVIATION.....AD LARGE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Adult Large
HELP_TEXT.....No Help Text
ID.....{38E64D72-2806-47FA-B82F-51844574F850}
TERM.....ADULT LARGE

DATA VALUE: ADULT SMALL

ABBREVIATION.....AD SMALL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Adult Small
HELP_TEXT.....No Help Text
ID.....{0161A4E5-713A-406B-B3FD-9ED13BEB5366}
TERM.....ADULT SMALL

DATA VALUE: ADVANCED CARE PLANNING

ABBREVIATION.....ADV CAR PL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Advanced Care Planning
HELP_TEXT.....No Help Text
ID.....{F0ED0F8A-1AD6-4F3E-A4C8-4FB915474FE7}
TERM.....ADVANCED CARE PLANNING

DATA VALUE: AGITATED

ABBREVIATION.....a
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Agitated
HELP_TEXT.....No Help Text
ID.....{BAD2CCCE-9083-6F61-1196-C0C9A8DDBFFC}
TERM.....AGITATED

Page: 673

DATA VALUE: AGONAL

ABBREVIATION.....a
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Agonal
HELP_TEXT.....No Help Text
ID.....{1621B411-33E3-3A6A-78A6-40B7BDC3E8CB}
TERM.....AGONAL

DATA VALUE: AGONIZING

ABBREVIATION.....AGONIZING

ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Agonizing
HELP_TEXT.....No Help Text
ID.....{F43DF9E7-555B-49FF-987C-35F083D547AA}
TERM.....AGONIZING

DATA VALUE: AIR FLOW DEVICE

ABBREVIATION.....AIR FLOW D
ACTIVE.....Yes
DESCRIPTION.....Air Flow Device
DISPLAY_NAME.....Air Flow Device
HELP_TEXT.....Air Flow Device
ID.....{60AF6215-CE2A-4EEE-9596-909E3E50BE71}
TERM.....AIR FLOW DEVICE

DATA VALUE: AIR LEAK

ABBREVIATION.....AIR LEAK
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Air Leak
HELP_TEXT.....No Help Text
ID.....{52E40E25-6D28-53EC-BC05-D1FEAA4616EF}
TERM.....AIR LEAK

DATA VALUE: AIRBORNE

ABBREVIATION.....AIRBORNE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Airborne
HELP_TEXT.....No Help Text
ID.....{D2F281F7-DC09-4CA2-B5B3-F92ACBC63A0A}
TERM.....AIRBORNE

Page: 674

DATA VALUE: AIRWAY REFLEXES INTACT / SECRETIONS

ABBREVIATION.....AIR REFL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Airway Reflexes Intact Manages
Secretions
HELP_TEXT.....No Help Text
ID.....{88A7E7D3-6283-4640-A5FD-0F1E37FB8249}
TERM.....AIRWAY REFLEXES INTACT / SECRETIONS

DATA VALUE: ALBUMIN

ABBREVIATION.....ALBUMIN
ACTIVE.....Yes
DESCRIPTION.....Albumin

DISPLAY_NAME.....Albumin
HELP_TEXT.....No Help Text
ID.....{4CC265EE-F506-4752-8309-E569A25B40E5}
TERM.....ALBUMIN

DATA VALUE: ALDRETE SCORE 8 OR GREATER

ABBREVIATION.....ALD 8+
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Aldrete Score 8 or Greater
HELP_TEXT.....No Help Text
ID.....{1FACE6D7-BC92-4579-8665-C5E1D3475461}
TERM.....ALDRETE SCORE 8 OR GREATER

DATA VALUE: ALDRETE SCORE AT BASELINE

ABBREVIATION.....ALD BASE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Aldrete Score at Baseline
HELP_TEXT.....No Help Text
ID.....{CC86D48F-9473-4AD8-AA5B-3BA0D1CCA8FD}
TERM.....ALDRETE SCORE AT BASELINE

DATA VALUE: ALERT

ABBREVIATION.....a
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Alert
HELP_TEXT.....No Help Text
ID.....{BC7E2630-179A-5C80-E6C7-598D42F0BFD6}
TERM.....ALERT

Page: 675

DATA VALUE: ALERT AND ORIENTED

ABBREVIATION.....AL AN OR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Alert and Oriented
HELP_TEXT.....No Help Text
ID.....{C031F776-36DF-4570-AFAF-A30B8BFA635B}
TERM.....ALERT AND ORIENTED

DATA VALUE: ALL DOWN

ABBREVIATION.....ALL DOWN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....All Down
HELP_TEXT.....No Help Text
ID.....{76C14E79-95E2-49B7-B847-3B114653E98C}

TERM.....ALL DOWN

DATA VALUE: ALL UP

ABBREVIATION.....ALL UP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....All Up
HELP_TEXT.....No Help Text
ID.....{48199083-ACC2-408C-9242-7A6313CD7A7F}
TERM.....ALL UP

DATA VALUE: ALLERGY

ABBREVIATION.....a
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Allergy
HELP_TEXT.....No Help Text
ID.....{EBF87F6F-1BDB-6147-C444-2FF787DEAC7E}
TERM.....ALLERGY

DATA VALUE: ALOPECIA PARTIAL

ABBREVIATION.....ALOPE PART
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Alopecia Partial
HELP_TEXT.....No Help Text
ID.....{98EDBD8D-205A-4630-A09D-6B4222F5B2A2}
TERM.....ALOPECIA PARTIAL

Page: 676

DATA VALUE: ALOPECIA TOTAL

ABBREVIATION.....ALOPE TOTA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Alopecia Total
HELP_TEXT.....No Help Text
ID.....{ECDCE3C4-81F6-4FA1-BADA-97A20BD9A826}
TERM.....ALOPECIA TOTAL

DATA VALUE: AMBER

ABBREVIATION.....AMBER
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Amber
HELP_TEXT.....No Help Text
ID.....{90A13048-4C35-E01D-3482-64D782001D4A}
TERM.....AMBER

DATA VALUE: AMBULATE AND ASSIST

ABBREVIATION.....AMB ASSIS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Ambulate and Assist
HELP_TEXT.....No Help Text
ID.....{EB53ADD9-B5B0-4174-95C7-C4B33E80DBDC}
TERM.....AMBULATE AND ASSIST

DATA VALUE: AMBULATED

ABBREVIATION.....AMBULATED
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Ambulated
HELP_TEXT.....No Help Text
ID.....{017D99FE-91D5-49F1-AC15-50620879C466}
TERM.....AMBULATED

DATA VALUE: AMBULATORY CARE UNIT

ABBREVIATION.....ACU
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Ambulatory Care Unit
HELP_TEXT.....No Help Text
ID.....{7FB7FBC5-42A6-4BA4-9CA6-494D77C4C549}
TERM.....AMBULATORY CARE UNIT

Page: 677

DATA VALUE: AMPUTEE

ABBREVIATION.....AMPUTEE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Amputee
HELP_TEXT.....No Help-Text
ID.....{5804286F-2B23-91BD-FCB5-6AD18C918203}
TERM.....AMPUTEE

DATA VALUE: ANALGESIA

ABBREVIATION.....ANALGESIA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Analgesia
HELP_TEXT.....No Help Text
ID.....{5452ED11-25E1-437B-8048-3BB3E03BE3B9}
TERM.....ANALGESIA

DATA VALUE: ANGINAL

ABBREVIATION.....ANGINAL
ACTIVE.....Yes

DESCRIPTION.....No Description
DISPLAY_NAME.....ANGINAL
HELP_TEXT.....No Help Text
ID.....{ECF02415-4124-46D3-AA01-F21AC40FEFF7}
TERM.....ANGINAL

DATA VALUE: ANNOYING

ABBREVIATION.....ANNOYING
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Annoying
HELP_TEXT.....No Help Text
ID.....{53437530-57C9-4EA8-AC0D-36B5F2BD051D}
TERM.....ANNOYING

DATA VALUE: ANTI-EMBOLUS STOCKINGS KNEE HIGH

ABBREVIATION.....ANT EMB KN
ACTIVE.....Yes
DESCRIPTION.....ANTI-EMBOLUS STOCKINGS KNEE HIGH
DISPLAY_NAME.....Anti-Embolus Stockings Knee High
HELP_TEXT.....ANTI-EMBOLUS STOCKINGS KNEE HIGH
ID.....{4FBFF7A8-D719-98B9-4A46-D76D7BB5213E}
TERM.....ANTI-EMBOLUS STOCKINGS KNEE HIGH

Page: 678

DATA VALUE: ANTI-EMBOLUS STOCKINGS THIGH HIGH

ABBREVIATION.....ANT EMB ST
ACTIVE.....Yes
DESCRIPTION.....ANTI-EMBOLUS STOCKINGS THIGH HIGH
DISPLAY_NAME.....Anti-Embolus Stockings Thigh High
HELP_TEXT.....ANTI-EMBOLUS STOCKINGS THIGH HIGH
ID.....{9B32D3F3-47F1-886F-A59A-D8AB38A34754}
TERM.....ANTI-EMBOLUS STOCKINGS THIGH HIGH

DATA VALUE: ANTI-MICROBIAL DRESSING

ABBREVIATION.....AB DRESSIN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Anti-Microbial Dressing
HELP_TEXT.....No Help Text
ID.....{9734C5E9-D099-47EB-93B1-D82201F64331}
TERM.....ANTI-MICROBIAL DRESSING

DATA VALUE: ANTIBIOTIC OINTMENT

ABBREVIATION.....AB OINTMEN
ACTIVE.....Yes
DESCRIPTION.....Antibiotic Ointment
DISPLAY_NAME.....Antibiotic Ointment
HELP_TEXT.....No Help Text

ID.....{352D9690-77BB-44DC-919E-5357876DEAC4}
TERM.....ANTIBIOTIC OINTMENT

DATA VALUE: ANTIMICROBIAL DISK

ABBREVIATION.....WOU DR ADI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Antimicrobial Disk
HELP_TEXT.....No Help-Text
ID.....{1F1745EA-D763-4AFD-9CED-82E974E57692}
TERM.....ANTIMICROBIAL DISK

DATA VALUE: ANURIC

ABBREVIATION.....ANURIC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Anuric
HELP_TEXT.....No Help Text
ID.....{0CCA4B44-7A59-4042-BAB9-EB2395A63C51}
TERM.....ANURIC

Page: 679

DATA VALUE: ANXIOUS

ABBREVIATION.....a
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Anxious
HELP_TEXT.....No Help Text
ID.....{381BDB41-7E8A-6FDB-841B-2ECC1889B9EF}
TERM.....ANXIOUS

DATA VALUE: AOO

ABBREVIATION.....AOO
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....AOO
HELP_TEXT.....No Help Text
ID.....{86F25E48-AADB-5B30-4B89-E8BE63DAE7D9}
TERM.....AOO

DATA VALUE: APPROPRIATE

ABBREVIATION.....a
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Appropriate
HELP_TEXT.....No Help Text
ID.....{1CFB285F-34DF-73A2-B89C-F6EE1F40E7E9}
TERM.....APPROPRIATE

DATA VALUE: ARTERIAL LINE CARE

ABBREVIATION.....ART LINE C
ACTIVE.....Yes
DESCRIPTION.....Arterial Line Care
DISPLAY_NAME.....Arterial Line Care
HELP_TEXT.....No Help Text
ID.....{787D232D-DFB8-44E2-8A0C-9D7D6B80764E}
TERM.....ARTERIAL LINE CARE

DATA VALUE: ARTIFICIAL AIRWAY

ABBREVIATION.....ART AIRWAY
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Artificial Airway
HELP_TEXT.....No Help Text
ID.....{8DEE2358-8892-4BFB-AFAE-6C134D1CCAD4}
TERM.....ARTIFICIAL AIRWAY

Page: 680

DATA VALUE: ASHEN

ABBREVIATION.....ASHEN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Ashen
HELP_TEXT.....No Help-Text
ID.....{7431E4EC-DCF1-4966-85CB-29BF65311855}
TERM.....ASHEN

DATA VALUE: ASLEEP AT ROUNDS

ABBREVIATION.....SLEEP ROUN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Appeared Asleep At All Rounds
HELP_TEXT.....No Help Text
ID.....{24E1D9F0-D8F1-4DC8-911F-6170F343C087}
TERM.....ASLEEP AT ROUNDS

DATA VALUE: ASPIRATION

ABBREVIATION.....ASPIRATION
ACTIVE.....Yes
DESCRIPTION.....Aspiration
DISPLAY_NAME.....Aspiration
HELP_TEXT.....No Help Text
ID.....{68B91085-8E3C-4B94-9F97-6305F5DFD900}
TERM.....ASPIRATION

DATA VALUE: ASSAULTIVE

ABBREVIATION.....a

ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Assaultive
HELP_TEXT.....No Help Text
ID.....{99ECD8C7-191A-3B9C-2A7A-803A921065FA}
TERM.....ASSAULTIVE

DATA VALUE: ASSIST CONTROL

ABBREVIATION.....AC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Assist Control
HELP_TEXT.....No Help Text
ID.....{862ED94C-9DC3-DFAF-F824-0F7DE348C2CB}
TERM.....ASSIST CONTROL

Page: 681

DATA VALUE: ASSISTED

ABBREVIATION.....ASSISTED
ACTIVE.....Yes
DESCRIPTION.....Assisted
DISPLAY_NAME.....Assisted
HELP_TEXT.....No Help Text
ID.....{65069804-39BF-41ED-9EA3-BF2C41389F10}
TERM.....ASSISTED

DATA VALUE: ASYSTOLE

ABBREVIATION.....a
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Asystole
HELP_TEXT.....No Help Text
ID.....{6E87ECBC-3160-2CEF-D2FF-EEBD1DAA96F2}
TERM.....ASYSTOLE

DATA VALUE: ATRIAL

ABBREVIATION.....a
ACTIVE.....Yes
DESCRIPTION.....ATRIAL
DISPLAY_NAME.....Atrial
HELP_TEXT.....No Help Text
ID.....{266D2FDA-5D6A-203D-6FF0-8282D49BFD67}
TERM.....ATRIAL

DATA VALUE: ATRIAL FIB/FLUTTER

ABBREVIATION.....af
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Atrial Fib/Flutter

HELP_TEXT.....No Help Text
ID.....{9C1BC2F6-3BF0-BEBB-BE87-50E9F7984544}
TERM.....ATRIAL FIB/FLUTTER

DATA VALUE: ATRIAL FIBRILLATION

ABBREVIATION.....a fib
ACTIVE.....Yes
DESCRIPTION.....Atrial Fibrillation
DISPLAY_NAME.....Atrial Fibrillation
HELP_TEXT.....Atrial Fibrillation
ID.....{7B5EBAF2-55B6-2905-3D89-B28CE6A509CB}
TERM.....ATRIAL FIBRILLATION

Page: 682

DATA VALUE: ATRIAL FLUTTER

ABBREVIATION.....af
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Atrial Flutter
HELP_TEXT.....No Help Text
ID.....{9CB8D413-6762-D78E-F40E-25DB66CD23B3}
TERM.....ATRIAL FLUTTER

DATA VALUE: AUSCULTATION

ABBREVIATION.....AUSCULTATI
ACTIVE.....Yes
DESCRIPTION.....Auscultation
DISPLAY_NAME.....Auscultation
HELP_TEXT.....Auscultation
ID.....{B9D756DF-74F5-4604-B781-D7BC512FC453}
TERM.....AUSCULTATION

DATA VALUE: AVF

ABBREVIATION.....AVF
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....AVF
HELP_TEXT.....No Help Text
ID.....{C8E97206-58AA-47C0-BD2E-7A30E56DB46A}
TERM.....AVF

DATA VALUE: AVI

ABBREVIATION.....AVI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....AVI
HELP_TEXT.....No Help Text
ID.....{8BB49456-FB80-4916-6EBE-06837BF12B0C}
TERM.....AVI

DATA VALUE: AVL

ABBREVIATION.....AVL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....AVL
HELP_TEXT.....No Help Text
ID.....{9124A193-5F50-448D-9E60-C5869EC4029C}
TERM.....AVL

Page: 683

DATA VALUE: AVO

ABBREVIATION.....AVO
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....AVO
HELP_TEXT.....No Help Text
ID.....{02EA752A-D156-96E6-3D28-3D389FFC7B58}
TERM.....AVO

DATA VALUE: AVR

ABBREVIATION.....AVR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....AVR
HELP_TEXT.....No Help Text
ID.....{04AF2551-1185-4BA9-AEED-977D562D188A}
TERM.....AVR

DATA VALUE: BACTERICIDAL

ABBREVIATION.....BACTERIOCI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bactericidal
HELP_TEXT.....No Help-Text
ID.....{680E4A3C-A40A-461E-818D-003AE8A1EFC8}
TERM.....BACTERICIDAL

DATA VALUE: BACTERIOSTATIC

ABBREVIATION.....BACTERIOST
ACTIVE.....Yes
DESCRIPTION.....Bacteriostatic
DISPLAY_NAME.....Bacteriostatic
HELP_TEXT.....No Help Text
ID.....{AB74D109-8ED7-4A6C-B028-7A041781DF67}
TERM.....BACTERIOSTATIC

DATA VALUE: BAG CHANGED

ABBREVIATION.....BAG CHANGE
ACTIVE.....Yes
DESCRIPTION.....Changed Bag
DISPLAY_NAME.....Bag Changed
HELP_TEXT.....No Help Text
ID.....{818990DE-4F95-4F8E-92D8-B4D3B975CB36}
TERM.....BAG CHANGED

Page: 684

DATA VALUE: BALSAM OF PERU

ABBREVIATION.....WOU DR PER
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Balsam Of Peru
HELP_TEXT.....No Help-Text
ID.....{DD412BB1-B1C1-4A48-A38F-832629448E06}
TERM.....BALSAM OF PERU

DATA VALUE: BARIATRIC FRAME

ABBREVIATION.....BARI FRAME
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bariatric Frame
HELP_TEXT.....No Help Text
ID.....{CF95B426-F280-451B-86F2-CEF7A1A2B967}
TERM.....BARIATRIC FRAME

DATA VALUE: BARIATRIC LOW AIRLOSS MATTRESS REPLACEMENT

ABBREVIATION.....BAR L AIR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bariatric Low Airloss Mattress
Replacement
HELP_TEXT.....No Help Text
ID.....{220C6E92-1D0A-4513-B3A9-8065A6EC7DB1}
TERM.....BARIATRIC LOW AIRLOSS MATTRESS
REPLACEMENT

DATA VALUE: BARIATRIC PATIENT

ABBREVIATION.....BARIAC PAT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bariatric Patient
HELP_TEXT.....No Help Text
ID.....{6F614334-685E-4533-A475-BDF5440855EA}
TERM.....BARIATRIC PATIENT

DATA VALUE: BASAL RATE ONLY

ABBREVIATION.....BASAL ONL

ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Basal Rate Only
HELP_TEXT.....No Help Text
ID.....{4BFFF2BC-837C-4D62-B811-0E547247D4D0}
TERM.....BASAL RATE ONLY

Page: 685

DATA VALUE: BATHROOM

ABBREVIATION.....b
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bathroom
HELP_TEXT.....No Help Text
ID.....{D20D91A7-E9F5-D000-5E0B-772BA9F5F72E}
TERM.....BATHROOM

DATA VALUE: BATHROOM PRIVILEGES

ABBREVIATION.....BATHR PRIV
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bathroom Privileges
HELP_TEXT.....No Help Text
ID.....{29C1B6DD-3308-4DED-9970-1322A97C633B}
TERM.....BATHROOM PRIVILEGES

DATA VALUE: BED

ABBREVIATION.....BED
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bed
HELP_TEXT.....No Help Text
ID.....{7E447C9C-41A0-43C9-AF63-38D6EA68279C}
TERM.....BED

DATA VALUE: BED EXIT ALARM

ABBREVIATION.....BED EXIT A
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bed Exit Alarm
HELP_TEXT.....No Help Text
ID.....{AC53759E-FFAE-5E0F-2EBA-E56A2712357C}
TERM.....BED EXIT ALARM

DATA VALUE: BEDREST

ABBREVIATION.....BEDREST
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bedrest

HELP_TEXT.....No Help Text
ID.....{8CD28057-132E-4CEE-9EA1-F03C8AF50CE8}
TERM.....BEDREST

Page: 686

DATA VALUE: BILATERAL

ABBREVIATION.....Bilateral
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bilateral
HELP_TEXT.....No Help Text
ID.....{FCD5361B-C674-473C-8F29-7B08650BA338}
TERM.....BILATERAL

DATA VALUE: BIPAP - RESPIRONICS

ABBREVIATION.....BIPAP RESP
ACTIVE.....Yes
DESCRIPTION.....BiPap machine by Respironics
DISPLAY_NAME.....Bipap - Respironics
HELP_TEXT.....No Help Text
ID.....{AA94A723-4C68-489F-88AC-2C6F276EAB76}
TERM.....BIPAP - RESPIRONICS

DATA VALUE: BIPAP - VISION

ABBREVIATION.....BIPAP VISI
ACTIVE.....Yes
DESCRIPTION.....BiPap ventilator made by Vision
DISPLAY_NAME.....Bipap - Vision
HELP_TEXT.....No Help Text
ID.....{AA39C765-BA1F-4289-B459-200A57C07FFF}
TERM.....BIPAP - VISION

DATA VALUE: BIPAP/CPAP

ABBREVIATION.....BIPAP/CPAP
ACTIVE.....Yes
DESCRIPTION.....BiPap / CPAP machine
DISPLAY_NAME.....BiPap / CPAP machine
HELP_TEXT.....No Help Text
ID.....{8D767C67-3712-410F-A96B-69C731E03C5D}
TERM.....BIPAP/CPAP

DATA VALUE: BIPHASIC POSITIVE AIRWAY PRESSURE

ABBREVIATION.....BiPAP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....BiPap
HELP_TEXT.....No Help Text
ID.....{4802AB1E-A0AD-BB9F-E8BA-33863BCD3466}
TERM.....BIPHASIC POSITIVE AIRWAY PRESSURE

DATA VALUE: BISMUTH GAUZE

ABBREVIATION.....WOU DR BIS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bismuth Gauze
HELP_TEXT.....No Help-Text
ID.....{010DED59-C808-4D92-A7F8-87AEB86B4439}
TERM.....BISMUTH GAUZE

DATA VALUE: BITE

ABBREVIATION.....Bite
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bite
HELP_TEXT.....No Help Text
ID.....{FC33E0AE-3130-4881-B9EF-9B81FBEC6AC0}
TERM.....BITE

DATA VALUE: BITE BLOCK

ABBREVIATION.....BITE BLOCK
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bite Block
HELP_TEXT.....No Help Text
ID.....{E7CC51F3-08A8-5FEE-38B2-16B74E1F542E}
TERM.....BITE BLOCK

DATA VALUE: BLACK

ABBREVIATION.....BLACK
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Black
HELP_TEXT.....No Help Text
ID.....{7A55AFB9-6861-4918-0746-3D67D47FE5C3}
TERM.....BLACK

DATA VALUE: BLACK/NECROTIC

ABBREVIATION.....BLACK/NECR
ACTIVE.....Yes
DESCRIPTION.....Black / Necrotic
DISPLAY_NAME.....Black / Necrotic
HELP_TEXT.....No Help Text
ID.....{53519F9A-B844-41C7-A013-F9D2EA8E4BC0}
TERM.....BLACK/NECROTIC

DATA VALUE: BLAND

ABBREVIATION.....BLAND
ACTIVE.....Yes
DESCRIPTION.....BLAND DIET
DISPLAY_NAME.....BLAND
HELP_TEXT.....No Help Text
ID.....{38930AFE-ECF7-4548-857C-854E1B5BE148}
TERM.....BLAND

DATA VALUE: BLISTERING

ABBREVIATION.....BLISTERING
ACTIVE.....Yes
DESCRIPTION.....Blistering
DISPLAY_NAME.....Blistering
HELP_TEXT.....No Help Text
ID.....{3A4212D5-F2A8-4636-BB7F-F36FA8CED192}
TERM.....BLISTERING

DATA VALUE: BLOOD STREAKED

ABBREVIATION.....BLOOD STRE
ACTIVE.....Yes
DESCRIPTION.....Blood Streaked
DISPLAY_NAME.....Blood Streaked
HELP_TEXT.....Blood Streaked
ID.....{AEC5607D-806F-49AE-8436-DC47F49F39BC}
TERM.....BLOOD STREAKED

DATA VALUE: BLOOD TRANSFUSION IV

ABBREVIATION.....BI TRAN IV
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Blood Transfusion IV
HELP_TEXT.....No Help Text
ID.....{27A0BE4C-433C-F920-BAC2-066F304BDF59}
TERM.....BLOOD TRANSFUSION IV

DATA VALUE: BLOODY

ABBREVIATION.....BLOODY
ACTIVE.....Yes
DESCRIPTION.....Containing blood
DISPLAY_NAME.....Bloody
HELP_TEXT.....No Help Text
ID.....{D7904E84-8331-D6DD-F2C0-7F3039C5746A}
TERM.....BLOODY

DATA VALUE: BLUE

ABBREVIATION.....BLUE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Blue
HELP_TEXT.....No Help Text
ID.....{8FA78575-E1DD-4E6E-BC43-9E5C92E2B7E5}
TERM.....BLUE

DATA VALUE: BLUNTED

ABBREVIATION.....Blunted
ACTIVE.....Yes
DESCRIPTION.....Blunted
DISPLAY_NAME.....Blunted
HELP_TEXT.....No Help Text
ID.....{F119F2DA-9182-2480-BA75-02D39F4FF69B}
TERM.....BLUNTED

DATA VALUE: BONE

ABBREVIATION.....BONE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bone
HELP_TEXT.....No Help Text
ID.....{3729DEB9-0967-4B32-86B7-409CF27C8625}
TERM.....BONE

DATA VALUE: BORDERLINE REDUCED

ABBREVIATION.....BORD REDUC
ACTIVE.....Yes
DESCRIPTION.....Borderline Reduced
DISPLAY_NAME.....Borderline Reduced
HELP_TEXT.....No Help Text
ID.....{63703AB1-6249-403B-BDE8-FECA9FA36CF5}
TERM.....BORDERLINE REDUCED

DATA VALUE: BRADEN FRICTION AND SHEAR 1 PROBLEM

ABBREVIATION.....BS F/S 1
ACTIVE.....Yes
DESCRIPTION.....Braden Scale Friction and Shear 1
Problem
DISPLAY_NAME.....1 - Problem
HELP_TEXT.....Requires moderate to maximum
assistance in moving. Complete
lifting without sliding against
sheets is impossible. Frequently
slides down in bed or chair,
requiring frequent repositioning with
maximum assistance.
ID.....{AF8D2952-E0DA-411E-91CC-0114F4EB1C33}
TERM.....BRADEN FRICTION AND SHEAR 1 PROBLEM

DATA VALUE: BRADEN SCALE ACTIVITY 1 BEDFAST

ABBREVIATION.....BS ACT 1
ACTIVE.....Yes
DESCRIPTION.....Braden Scale Activity 1 Bedfast
DISPLAY_NAME.....1 - Bedfast
HELP_TEXT.....Confined to bed
ID.....{B9382966-5DE9-4D1C-BEC7-4DEBC5D6DCF5}
TERM.....BRADEN SCALE ACTIVITY 1 BEDFAST

DATA VALUE: BRADEN SCALE ACTIVITY 2 CHAIRFAST

ABBREVIATION.....BS ACT 2
ACTIVE.....Yes
DESCRIPTION.....Braden Scale Activity 2 Chairfast
DISPLAY_NAME.....2 - Chairfast
HELP_TEXT.....Ability to walk severely limited or
non-existent. Cannot bear own weight
and/or must be assisted into chair or
wheelchair.
ID.....{CDBB8C3B-B53A-4A0D-91ED-8A53CFE1485F}
TERM.....BRADEN SCALE ACTIVITY 2 CHAIRFAST

DATA VALUE: BRADEN SCALE ACTIVITY 3 WALKS OCCASIONALLY

ABBREVIATION.....BS ACT 3
ACTIVE.....Yes
DESCRIPTION.....Braden Scale Activity 3 Walks
Occasionally
DISPLAY_NAME.....3 - Walks Occasionally
HELP_TEXT.....Walks occasionally during day, but
for very short distances, with or
without assistance. Spends majority
of each shift in bed or chair.
ID.....{E47B47BC-6097-4602-9184-BF5A4F6D686A}
TERM.....BRADEN SCALE ACTIVITY 3 WALKS
OCCASIONALLY

DATA VALUE: BRADEN SCALE ACTIVITY 4 WALKS FREQUENTLY

ABBREVIATION.....BS ACT 4
ACTIVE.....Yes
DESCRIPTION.....Braden Scale Activity 4 Walks
Frequently
DISPLAY_NAME.....4 - Walks Frequently
HELP_TEXT.....Walks outside room at least twice a
day and inside room at least once
every two hours during waking hours.
ID.....{C5A8C2EA-E042-4804-9812-77603C47FC6A}
TERM.....BRADEN SCALE ACTIVITY 4 WALKS
FREQUENTLY

DATA VALUE: BRADEN SCALE FRICTION AND SHEAR 2 POTENTIAL PROBLEM

ABBREVIATION.....BS F/S 2
ACTIVE.....Yes
DESCRIPTION.....Braden Scale Friction and Shear 2
Potential Problem
DISPLAY_NAME.....2 - Potential Problem
HELP_TEXT.....Moves feebly or requires minimum
assistance. During a move skin
probably slides to some extent
against sheets, chair, restraints or
other devices. Maintains relatively
good position in chair or bed most of
the time but occasionally slides down
ID.....{B56338A9-8956-4D27-8BA0-DA50E07359E0}
TERM.....BRADEN SCALE FRICTION AND SHEAR 2
POTENTIAL PROBLEM

Page: 691

DATA VALUE: BRADEN SCALE FRICTION AND SHEAR 3 NO APPARENT PROBLEM

ABBREVIATION.....BS F/S 3
ACTIVE.....Yes
DESCRIPTION.....Braden Scale Friction and Shear 3 No
Apparent Problem
DISPLAY_NAME.....3 - No Apparent Problem
HELP_TEXT.....No Apparent Problem
ID.....{9850E2BF-6600-4E66-8FCD-1517D903BC0F}
TERM.....BRADEN SCALE FRICTION AND SHEAR 3 NO
APPARENT PROBLEM

DATA VALUE: BRADEN SCALE MOBILITY 1 COMPLETELY IMMOBILE

ABBREVIATION.....BS MOB 1
ACTIVE.....Yes
DESCRIPTION.....Braden Scale Mobility 1 Completely
Immobile
DISPLAY_NAME.....1 - Completely Immobile
HELP_TEXT.....Does not make even slight changes in
body or extremity position without
assistance
ID.....{7607574E-D4C2-4E64-8DC5-ABE7D00917A6}
TERM.....BRADEN SCALE MOBILITY 1 COMPLETELY
IMMOBILE

DATA VALUE: BRADEN SCALE MOBILITY 2 VERY LIMITED

ABBREVIATION.....BS MOB 2
ACTIVE.....Yes
DESCRIPTION.....Braden Scale Mobility 2 Very Limited
DISPLAY_NAME.....2 - Very Limited
HELP_TEXT.....Makes occasional slight changes in
body or extremity position but unable
to make frequent or significant
changes independently
ID.....{868EFA33-877E-4B91-B4FD-9655F66A7160}

TERM.....BRADEN SCALE MOBILITY 2 VERY LIMITED

DATA VALUE: BRADEN SCALE MOBILITY 3 SLIGHTLY LIMITED

ABBREVIATION.....BS MOB 3
ACTIVE.....Yes
DESCRIPTION.....Braden Scale Mobility 3 Slightly
 Limited
DISPLAY_NAME.....3 - Slightly Limited
HELP_TEXT.....Makes frequent though slight changes
 in body or extremity position
 independently
ID.....{C54B0A8F-46F1-4DFA-9C24-AB6C2A29CA82}
TERM.....BRADEN SCALE MOBILITY 3 SLIGHTLY
 LIMITED

DATA VALUE: BRADEN SCALE MOBILITY 4 NO LIMITATION

ABBREVIATION.....BS MOB 4
ACTIVE.....Yes
DESCRIPTION.....Braden Scale Mobility 4 No Limitation
DISPLAY_NAME.....4 - No Limitation
HELP_TEXT.....Makes major and frequent changes in
 position without assistance
ID.....{84BFB8D0-D6E4-4384-88EE-28B7C4EDB275}
TERM.....BRADEN SCALE MOBILITY 4 NO LIMITATION

Page: 692

DATA VALUE: BRADEN SCALE MOISTURE 1 CONSTANTLY MOIST

ABBREVIATION.....BS MOI 1
ACTIVE.....Yes
DESCRIPTION.....Braden Scale Moisture 1 Constantly
 Moist
DISPLAY_NAME.....1 - Constantly Moist
HELP_TEXT.....Skin is kept moist almost constantly
 by perspiration, urine, etc. Dampness
 is detected every time patient is
 moved or turned
ID.....{3E2569BA-ECF0-4622-8F4A-26378A2CB466}
TERM.....BRADEN SCALE MOISTURE 1 CONSTANTLY
 MOIST

DATA VALUE: BRADEN SCALE MOISTURE 2 VERY MOIST

ABBREVIATION.....BS MOI 2
ACTIVE.....Yes
DESCRIPTION.....Braden Scale Moisture 2 Very Moist
DISPLAY_NAME.....2 - Very Moist
HELP_TEXT.....Skin is often, but not always moist.
 Linen must be changed at least once a
 shift
ID.....{2C6D15A3-6166-4D7D-B200-B58EC9C2CABE}
TERM.....BRADEN SCALE MOISTURE 2 VERY MOIST

DATA VALUE: BRADEN SCALE MOISTURE 3 OCCASIONALLY MOIST

ABBREVIATION.....BS MOI 3
ACTIVE.....Yes
DESCRIPTION.....Braden Scale Moisture 3 Occasionally
Moist
DISPLAY_NAME.....3 - Occasionally Moist
HELP_TEXT.....Skin is occasionally moist, requiring
an extra linen change approximately
once a day
ID.....{D7F2714B-B0C7-4000-8588-8582C1B32ECE}
TERM.....BRADEN SCALE MOISTURE 3 OCCASIONALLY
MOIST

DATA VALUE: BRADEN SCALE MOISTURE 4 RARELY MOIST

ABBREVIATION.....BS MOI 4
ACTIVE.....Yes
DESCRIPTION.....Braden Scale Moisture 4 Rarely Moist
DISPLAY_NAME.....4 - Rarely Moist
HELP_TEXT.....Skin is usually dry, linen only
requires changing at routine
intervals
ID.....{745EA3B9-18C9-4CF3-9DB8-9806A02656E8}
TERM.....BRADEN SCALE MOISTURE 4 RARELY MOIST

DATA VALUE: BRADEN SCALE NUTRITION 1 VERY POOR

ABBREVIATION.....BS NUT 1
ACTIVE.....Yes
DESCRIPTION.....Braden Scale Nutrition 1 Very Poor
DISPLAY_NAME.....1 - Very Poor
HELP_TEXT.....Never eats a complete meal. Rarely
eats more than 1/2 of any food
offered. Eats 2 servings or less of
protein per day. Takes fluids poorly,
does not take a liquid supplement. OR
is NPO maintained on clear liquids or
IV's for more than 5 days
ID.....{28A86074-A6B2-455E-ABA4-19195C2E3679}
TERM.....BRADEN SCALE NUTRITION 1 VERY POOR

Page: 693

DATA VALUE: BRADEN SCALE NUTRITION 2 PROBABLY INADEQUATE

ABBREVIATION.....BS NUT 2
ACTIVE.....Yes
DESCRIPTION.....Braden Scale Nutrition 2 Probably
Inadequate
DISPLAY_NAME.....2 - Probably Inadequate
HELP_TEXT.....Rarely eats a complete meal and
generally eats only about 2 of any
food offered. Protein intake includes

only 3 servings of meat or dairy products per day. Occasionally will take a dietary supplement. OR receives less than optimum amount of liquid.

ID.....{7DAB22F6-22E4-4A8B-A6A5-CD1277EB6B87}
TERM.....BRADEN SCALE NUTRITION 2 PROBABLY INADEQUATE

DATA VALUE: BRADEN SCALE NUTRITION 3 ADEQUATE

ABBREVIATION.....BS NUT 3
ACTIVE.....Yes
DESCRIPTION.....Braden Scale Nutrition 3 Adequate
DISPLAY_NAME.....3 - Adequate
HELP_TEXT.....Eats over half of most meals. Eats a total of 4 servings of protein (meat, dairy products) per day. Occasionally will refuse a meal, but will usually take a supplement when offered. OR is on a tube feeding or TPN regimen.

ID.....{3A8A9C59-8DC5-49A9-BAC9-882CB75D8FD5}
TERM.....BRADEN SCALE NUTRITION 3 ADEQUATE

DATA VALUE: BRADEN SCALE NUTRITION 4 EXCELLENT

ABBREVIATION.....BS NUT 4
ACTIVE.....Yes
DESCRIPTION.....Braden Scale Nutrition 4 Excellent
DISPLAY_NAME.....4 - Excellent
HELP_TEXT.....Eats most of every meal. Never refuses a meal. Usually eats a total of 4 or more servings of meat and dairy products. Occasionally eats between meals. Does NOT require supplements

ID.....{F09E4D35-1AA8-49A5-9B38-FBFFAFCFC181}
TERM.....BRADEN SCALE NUTRITION 4 EXCELLENT

DATA VALUE: BRADEN SCALE SENSORY PERCEPTION 1 COMPLETELY LIMITED

ABBREVIATION.....BS SENC 1
ACTIVE.....Yes
DESCRIPTION.....Braden Scale Sensory Perception 1 Completely Limited
DISPLAY_NAME.....1 - Completely Limited
HELP_TEXT.....Unresponsive (does not moan, flinch or grasp) to painful stimuli, due to diminished level of consciousness or sedation. OR limited ability to feel pain over most of body

ID.....{60F646B3-790F-4C5A-B2A0-845915E120E6}
TERM.....BRADEN SCALE SENSORY PERCEPTION 1 COMPLETELY LIMITED

DATA VALUE: BRADEN SCALE SENSORY PERCEPTION 2 VERY LIMITED

ABBREVIATION.....BS SENC 2
ACTIVE.....Yes
DESCRIPTION.....Braden Scale Sensory Perception 2
Very Limited
DISPLAY_NAME.....2 - Very Limited
HELP_TEXT.....Responds only to painful stimuli.
Cannot Communicate discomfort except
by moaning or restlessness. OR has a
sensory impairment which limits the
ability to feel pain or discomfort
over 1/2 of the body
ID.....{98362E63-7245-4E7F-8C30-896A9B923199}
TERM.....BRADEN SCALE SENSORY PERCEPTION 2
VERY LIMITED

Page: 694

DATA VALUE: BRADEN SCALE SENSORY PERCEPTION 3 SLIGHTLY LIMITED

ABBREVIATION.....BS SENC 3
ACTIVE.....Yes
DESCRIPTION.....Braden Scale Sensory Perception 3
Slightly Limited
DISPLAY_NAME.....3 - Slightly Limited
HELP_TEXT.....Responds to verbal commands, but
cannot always communicate discomfort
or the need to be turned. OR has some
sensory impairment which limits
ability to feel pain or discomfort in
one or two extremities.
ID.....{1D01D47D-CDC2-40EE-837B-5417D6EE47D6}
TERM.....BRADEN SCALE SENSORY PERCEPTION 3
SLIGHTLY LIMITED

DATA VALUE: BRADEN SCALE SENSORY PERCEPTION 4 NO IMPAIRMENT

ABBREVIATION.....BS SENC 4
ACTIVE.....Yes
DESCRIPTION.....Braden Scale Sensory Perception 4 No
Impairment
DISPLAY_NAME.....4 - No Impairment
HELP_TEXT.....Responds to verbal commands. Has no
sensory deficit which would limit
ability to feel or voice pain or
discomfort.
ID.....{10304656-E5E8-4F32-8358-A9F80511EFDA}
TERM.....BRADEN SCALE SENSORY PERCEPTION 4 NO
IMPAIRMENT

DATA VALUE: BREAKDOWN

ABBREVIATION.....BREAKDOWN
ACTIVE.....Yes
DESCRIPTION.....No Description

DISPLAY_NAME.....Breakdown
HELP_TEXT.....Use comment to describe breakdown
ID.....{9CFD21C0-A805-47C4-E559-8033C754DE80}
TERM.....BREAKDOWN

DATA VALUE: BRIGHT

ABBREVIATION.....Bright
ACTIVE.....Yes
DESCRIPTION.....Bright
DISPLAY_NAME.....Bright
HELP_TEXT.....No Help Text
ID.....{D16542CE-0403-2C68-A361-125CB3C75ABA}
TERM.....BRIGHT

DATA VALUE: BRISK

ABBREVIATION.....brisk
ACTIVE.....Yes
DESCRIPTION.....Brisk
DISPLAY_NAME.....Brisk
HELP_TEXT.....No Help Text
ID.....{980AE687-0A78-458B-BA7B-5D646C38A5F5}
TERM.....BRISK

Page: 695

DATA VALUE: BROWN

ABBREVIATION.....BROWN
ACTIVE.....Yes
DESCRIPTION.....The color brown
DISPLAY_NAME.....Brown
HELP_TEXT.....The color brown
ID.....{A1C32FBC-7B4D-707C-CA58-92A2E201368E}
TERM.....BROWN

DATA VALUE: BRUISED

ABBREVIATION.....BRUISED
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bruised
HELP_TEXT.....No Help Text
ID.....{FA5B38A7-7744-186F-055A-8458D696AF33}
TERM.....BRUISED

DATA VALUE: BRUIT

ABBREVIATION.....BRUIT
ACTIVE.....Yes
DESCRIPTION.....Bruit
DISPLAY_NAME.....Bruit
HELP_TEXT.....No Help Text
ID.....{2627E78A-B8B3-44EE-BC6F-4DAB16394F8C}

TERM.....BRUIT

DATA VALUE: BUCKS TRACTION

ABBREVIATION.....BUCK TRACT
ACTIVE.....Yes
DESCRIPTION.....Buck's Traction
DISPLAY_NAME.....Buck's Traction
HELP_TEXT.....Buck's Traction
ID.....{5E63F920-02D7-4729-8F30-A53A0A6E1EA0}
TERM.....BUCKS TRACTION

DATA VALUE: BULB SUCTION

ABBREVIATION.....Bulb Sucti
ACTIVE.....Yes
DESCRIPTION.....Bulb Suction Drain
DISPLAY_NAME.....Bulb Suction Drain
HELP_TEXT.....No Help Text
ID.....{B0839A82-BB79-BB9F-5E4E-0ED67F716A50}
TERM.....BULB SUCTION

Page: 696

DATA VALUE: BULLET

ABBREVIATION.....BULLET
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bullet
HELP_TEXT.....No Help Text
ID.....{1570D8E2-D936-494D-AF6E-13838EA34EFD}
TERM.....BULLET

DATA VALUE: BURN

ABBREVIATION.....BURN
ACTIVE.....Yes
DESCRIPTION.....Burn
DISPLAY_NAME.....Burn
HELP_TEXT.....No Help Text
ID.....{97B37D7A-28C6-4ABF-8D09-1FC49D5F160F}
TERM.....BURN

DATA VALUE: BURNING

ABBREVIATION.....BURNING
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Burning
HELP_TEXT.....No Help Text
ID.....{D75E2D5D-9CF9-4E70-A2EB-0632EBDB82C9}
TERM.....BURNING

DATA VALUE: BY PROVIDER ORDER

ABBREVIATION.....BY ORDER
ACTIVE.....Yes
DESCRIPTION.....By Provider Order
DISPLAY_NAME.....By Provider Order
HELP_TEXT.....No Help Text
ID.....{C038BC16-44CF-4933-BECC-B0F1A1B12603}
TERM.....BY PROVIDER ORDER

DATA VALUE: CADEXOMER IODINE

ABBREVIATION.....CADEXO IOD
ACTIVE.....Yes
DESCRIPTION.....Cadexomer Iodine
DISPLAY_NAME.....Cadexomer Iodine
HELP_TEXT.....No Help Text
ID.....{C83E41AA-6A1F-43CC-AB2E-7F6555FECD92}
TERM.....CADEXOMER IODINE

Page: 697

DATA VALUE: CALCIUM ALGINATE

ABBREVIATION.....CA ALGINAT
ACTIVE.....Yes
DESCRIPTION.....Calcium Alginate
DISPLAY_NAME.....Calcium Alginate
HELP_TEXT.....No Help Text
ID.....{2B2ACD25-A77B-4D39-8598-B26FD9BAFFC2}
TERM.....CALCIUM ALGINATE

DATA VALUE: CALM

ABBREVIATION.....c
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Calm
HELP_TEXT.....No Help Text
ID.....{E206C3F5-6977-302C-51FB-AF2155B5FFD9}
TERM.....CALM

DATA VALUE: CANE

ABBREVIATION.....CANE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Cane
HELP_TEXT.....No Help Text
ID.....{562B2E0B-0298-44C6-B33B-BC57DE6CBC3C}
TERM.....CANE

DATA VALUE: CANES

ABBREVIATION.....CANES
ACTIVE.....Yes

DESCRIPTION.....No Description
DISPLAY_NAME.....Canes
HELP_TEXT.....No Help Text
ID.....{DF37C541-1C7D-4527-A4D2-31867B5ED438}
TERM.....CANES

DATA VALUE: CARDIAC CHAIR

ABBREVIATION.....cc
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Cardiac Chair
HELP_TEXT.....No Help Text
ID.....{D4A2CFC1-D3A4-696B-4062-AFB2DAFFD76A}
TERM.....CARDIAC CHAIR

Page: 698

DATA VALUE: CARDIAC DIET

ABBREVIATION.....CARD DIET
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....CARDIAC DIET
HELP_TEXT.....No Help Text
ID.....{218A883D-32A0-438D-B453-7CBCE02693A3}
TERM.....CARDIAC DIET

DATA VALUE: CARDIAC OUTPUT INJECTATE

ABBREVIATION.....CARD OUT I
ACTIVE.....Yes
DESCRIPTION.....Cardiac Output Injectate
DISPLAY_NAME.....Cardiac Output Injectate
HELP_TEXT.....No Help Text
ID.....{19605CAB-8D54-4704-819B-20ABDAB6C604}
TERM.....CARDIAC OUTPUT INJECTATE

DATA VALUE: CARDIAC REPOSITION

ABBREVIATION.....CARD REPOS
ACTIVE.....Yes
DESCRIPTION.....Cardiac Reposition
DISPLAY_NAME.....Cardiac Reposition
HELP_TEXT.....No Help Text
ID.....{E133FE5A-20E1-418F-BD90-D2DFE14908C0}
TERM.....CARDIAC REPOSITION

DATA VALUE: CAREGIVER

ABBREVIATION.....CAREGIVER
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Caregiver
HELP_TEXT.....No Help Text

ID.....{F9B8BABE-2F05-437D-85C8-585FC1FFC07F}
TERM.....CAREGIVER

DATA VALUE: CASS ENDOTRACHEAL TUBE

ABBREVIATION.....CASS
ACTIVE.....Yes
DESCRIPTION.....CASS Endotracheal Tube
DISPLAY_NAME.....CASS Endotracheal Tube
HELP_TEXT.....Mucous Slurping Endotracheal Tube
ID.....{02BF805C-DF35-4DAB-A99F-2ED673174E49}
TERM.....CASS ENDOTRACHEAL TUBE

Page: 699

DATA VALUE: CATHETER

ABBREVIATION.....CATHETER
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Catheter
HELP_TEXT.....No Help Text
ID.....{B0D1AAD7-012D-D3E8-3216-1F71A1EC685E}
TERM.....CATHETER

DATA VALUE: CELLULITIS

ABBREVIATION.....CELLULITIS
ACTIVE.....Yes
DESCRIPTION.....Cellulitis
DISPLAY_NAME.....Cellulitis
HELP_TEXT.....No Help Text
ID.....{709F8795-0D4A-4E98-AB82-DD3E2AA7EA7B}
TERM.....CELLULITIS

DATA VALUE: CENTRAL ECG

ABBREVIATION.....cecg
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Central ECG
HELP_TEXT.....No Help Text
ID.....{26D36175-E4E3-C364-3018-1B0DB4DC4793}
TERM.....CENTRAL ECG

DATA VALUE: CENTRAL VENOUS CATHETER

ABBREVIATION.....CENT VEN C
ACTIVE.....Yes
DESCRIPTION.....Central Venous Catheter
DISPLAY_NAME.....Central Venous Catheter
HELP_TEXT.....No Help Text
ID.....{EE91E841-C09B-4432-A121-8EBAB92E9CD2}
TERM.....CENTRAL VENOUS CATHETER

DATA VALUE: CERVICAL SPINAL PRECAUTIONS

ABBREVIATION.....CERV S P
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Cervical Spinal Precautions
HELP_TEXT.....No Help Text
ID.....{A15F7127-2367-464A-AB7D-1E6A7C09AA21}
TERM.....CERVICAL SPINAL PRECAUTIONS

Page: 700

DATA VALUE: CHAIR

ABBREVIATION.....CHAIR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Chair
HELP_TEXT.....No Help Text
ID.....{5B8A5DBD-7462-4653-8D1F-CA8B1226DE38}
TERM.....CHAIR

DATA VALUE: CHANGED

ABBREVIATION.....CHANGED
ACTIVE.....Yes
DESCRIPTION.....Changed
DISPLAY_NAME.....Changed
HELP_TEXT.....No Help Text
ID.....{74B14023-D74D-41F6-AA48-B0E5137623CD}
TERM.....CHANGED

DATA VALUE: CHEYNE-STOKES

ABBREVIATION.....CHEYNE-STO
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Cheyne-Stokes
HELP_TEXT.....No Help Text
ID.....{58DD2B80-B522-499C-2C60-C48A654BAD84}
TERM.....CHEYNE-STOKES

DATA VALUE: CIRCULATION CHECK

ABBREVIATION.....ccq
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Circulation Check
HELP_TEXT.....No Help Text
ID.....{BA32C2A8-F15C-84EF-4859-6277EFA6463A}
TERM.....CIRCULATION CHECK

DATA VALUE: CLAMPED

ABBREVIATION.....CLAMPED

ACTIVE.....Yes
DESCRIPTION.....Clamped
DISPLAY_NAME.....Clamped
HELP_TEXT.....No Help Text
ID.....{B86865B8-078A-4996-B5EC-46A4A99D5A7D}
TERM.....CLAMPED

Page: 701

DATA VALUE: CLAY

ABBREVIATION.....CLAY
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Clay
HELP_TEXT.....No Help Text
ID.....{A4E4D16C-B45C-0934-ED22-8BF640080108}
TERM.....CLAY

DATA VALUE: CLEAR

ABBREVIATION.....c
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Clear
HELP_TEXT.....No Help Text
ID.....{8981C365-2B46-704F-C9F4-975DF1C23C2D}
TERM.....CLEAR

DATA VALUE: CLEAR LIQUID DIET

ABBREVIATION.....CL LI DIET
ACTIVE.....Yes
DESCRIPTION.....Clear Liquid Diet
DISPLAY_NAME.....Clear Liquid Diet
HELP_TEXT.....No Help Text
ID.....{D540348D-3E1D-4E42-9EDC-43C56CF56206}
TERM.....CLEAR LIQUID DIET

DATA VALUE: CLICK

ABBREVIATION.....c
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Click
HELP_TEXT.....No Help Text
ID.....{F1F3A1A8-35E1-61B4-C783-3314C8247380}
TERM.....CLICK

DATA VALUE: CLINICAL ASSESSMENT

ABBREVIATION.....CLIN ASSME
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Clinical Assessment

HELP_TEXT.....No Help Text
ID.....{803BEA4B-E22F-48C2-B20E-D5B3C8C7084C}
TERM.....CLINICAL ASSESSMENT

Page: 702

DATA VALUE: CLOSE OBSERVATION

ABBREVIATION.....co
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Close Observation
HELP_TEXT.....No Help Text
ID.....{2679BCAC-8753-FE73-56F8-3AB632161AF4}
TERM.....CLOSE OBSERVATION

DATA VALUE: CLOSED

ABBREVIATION.....CLOSED
ACTIVE.....Yes
DESCRIPTION.....CLOSED
DISPLAY_NAME.....CLOSED
HELP_TEXT.....No Help Text
ID.....{9BA6EE03-F562-4DE2-99A7-337F5A5E34F2}
TERM.....CLOSED

DATA VALUE: CLOTHED SELF

ABBREVIATION.....CLOTH S
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Clothed Self
HELP_TEXT.....No Help Text
ID.....{33802A30-FAA7-426B-ADAC-5297818DA7DE}
TERM.....CLOTHED SELF

DATA VALUE: CLOTS

ABBREVIATION.....CLOTS
ACTIVE.....Yes
DESCRIPTION.....Clots Present
DISPLAY_NAME.....Clots
HELP_TEXT.....No Help Text
ID.....{B52B88E0-6D23-4657-97CD-81E3EF64991B}
TERM.....CLOTS

DATA VALUE: CLOUDY

ABBREVIATION.....CLOUDY
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Cloudy
HELP_TEXT.....No Help Text
ID.....{3CCFCCCC-00FC-6B00-1D58-A77C02FDB401}
TERM.....CLOUDY

DATA VALUE: CLUBBING

ABBREVIATION.....CLUBBING
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Clubbing
HELP_TEXT.....No Help Text
ID.....{AE77D344-6F0F-452B-B1AA-EC8D9F092CFC}
TERM.....CLUBBING

DATA VALUE: CMS INTACT DISTAL TO SURGICAL SITE

ABBREVIATION.....CMS INTA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....CMS Intact Distal to Surgical Site
HELP_TEXT.....No Help Text
ID.....{CE2D2CDE-71C5-4993-8E46-C4799301A851}
TERM.....CMS INTACT DISTAL TO SURGICAL SITE

DATA VALUE: COFFEE GROUND

ABBREVIATION.....COFFEE GRO
ACTIVE.....Yes
DESCRIPTION.....Coffee Ground
DISPLAY_NAME.....Coffee Ground
HELP_TEXT.....Coffee Ground
ID.....{A5B3119A-C739-4B8D-87C2-7D858B94F380}
TERM.....COFFEE GROUND

DATA VALUE: COLD

ABBREVIATION.....COLD
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Cold
HELP_TEXT.....No Help Text
ID.....{B7B422BA-22EC-117D-0A0C-8419E9447D46}
TERM.....COLD

DATA VALUE: COLLAGEN

ABBREVIATION.....COLLAGEN
ACTIVE.....Yes
DESCRIPTION.....Collagen
DISPLAY_NAME.....Collagen
HELP_TEXT.....No Help Text
ID.....{8D734495-2FAF-41C0-9146-55B22583A830}
TERM.....COLLAGEN

DATA VALUE: COLORLESS

ABBREVIATION.....COLORLESS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Colorless
HELP_TEXT.....No Help Text
ID.....{D093D8C5-E284-47D4-8FA6-DCD5B030B10E}
TERM.....COLORLESS

DATA VALUE: COLOSTOMY

ABBREVIATION.....COLOSTOMY
ACTIVE.....Yes
DESCRIPTION.....Colostomy
DISPLAY_NAME.....Colostomy
HELP_TEXT.....No Help Text
ID.....{53F939AB-6AF4-4DD5-8D64-17AC1D1736B7}
TERM.....COLOSTOMY

DATA VALUE: COMBATIVE

ABBREVIATION.....c
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Combative
HELP_TEXT.....No Help Text
ID.....{DF65D42D-6CAA-657E-AEE1-C070C1FDF177}
TERM.....COMBATIVE

DATA VALUE: COMM MODE

ABBREVIATION.....COMM MODE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Commode
HELP_TEXT.....No Help-Text
ID.....{75C896CA-F105-4EF0-B641-3F2B1F984F1B}
TERM.....COMM MODE

DATA VALUE: COMMUNICATION BARRIER

ABBREVIATION.....COMM BARRI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Communication Barrier
HELP_TEXT.....No Help Text
ID.....{1A24946D-D966-4FAB-B6EA-940FA9EF60D3}
TERM.....COMMUNICATION BARRIER

DATA VALUE: COMPLETE BATH

ABBREVIATION.....COMP BATH
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Complete Bath
HELP_TEXT.....No Help Text
ID.....{A7CB1F56-F88E-403B-AABC-105BE129FBB0}
TERM.....COMPLETE BATH

DATA VALUE: COMPLIANCE

ABBREVIATION.....COMPLIANC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Compliance
HELP_TEXT.....No Help Text
ID.....{EF43367E-9898-44C3-8AC8-F5BD28051F1D}
TERM.....COMPLIANCE

DATA VALUE: CONFUSED

ABBREVIATION.....CONFUSED
ACTIVE.....Yes
DESCRIPTION.....Confused
DISPLAY_NAME.....Confused
HELP_TEXT.....No Help Text
ID.....{31DAAA60-BE89-49A7-A10D-EB0DA82F12C0}
TERM.....CONFUSED

DATA VALUE: CONSISTENT WITH BODY TEMPERATURE

ABBREVIATION.....CON W TEMP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Consistent With Body Temperature
HELP_TEXT.....No Help Text
ID.....{941FC5EF-6AF8-4DA1-A43C-95575AC74B57}
TERM.....CONSISTENT WITH BODY TEMPERATURE

DATA VALUE: CONSTANT

ABBREVIATION.....CONSTANT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Constant
HELP_TEXT.....No Help Text
ID.....{0FB8E43C-2106-4C37-9430-38675953A85B}
TERM.....CONSTANT

DATA VALUE: CONSTIPATION

ABBREVIATION.....CONSTIPAT.
ACTIVE.....Yes
DESCRIPTION.....Constipation

DISPLAY_NAME.....Constipation
HELP_TEXT.....No Help Text
ID.....{4A3B0D47-88BF-4FEB-AA0D-7AD7FED2C5C3}
TERM.....CONSTIPATION

DATA VALUE: CONSTRICTED

ABBREVIATION.....Constrict
ACTIVE.....Yes
DESCRIPTION.....Constricted
DISPLAY_NAME.....Constricted
HELP_TEXT.....No Help Text
ID.....{A7F66C5B-89EF-DEEE-A6C0-0DCDFEB3AEB4}
TERM.....CONSTRICTED

DATA VALUE: CONTACT

ABBREVIATION.....CONTACT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Contact
HELP_TEXT.....No Help Text
ID.....{958F2B5E-F980-47B6-A7F7-EE7D5520B61F}
TERM.....CONTACT

DATA VALUE: CONTEMPLATIVE

ABBREVIATION.....CONTEMPLAT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Contemplative
HELP_TEXT.....No Help Text
ID.....{5194DB04-7A68-491E-B117-BD089C912DFB}
TERM.....CONTEMPLATIVE

DATA VALUE: CONTINENT

ABBREVIATION.....CONTINENT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Continent
HELP_TEXT.....No Help Text
ID.....{3B658506-F42E-BAF5-04C4-5AD49699E0D2}
TERM.....CONTINENT

Page: 707

DATA VALUE: CONTINUOUS

ABBREVIATION.....CONTINUOUS
ACTIVE.....Yes
DESCRIPTION.....Continuous
DISPLAY_NAME.....Continuous
HELP_TEXT.....No Help Text
ID.....{676CD2DA-F882-4BE5-99BB-510EA7669749}

TERM.....CONTINUOUS

DATA VALUE: CONTINUOUS POSITIVE AIR PRESSURE

ABBREVIATION.....CPAP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....CPAP
HELP_TEXT.....No Help Text
ID.....{387B42F0-F587-AAFC-ECA0-66067F18BEC6}
TERM.....CONTINUOUS POSITIVE AIR PRESSURE

DATA VALUE: CONTINUOUS SUCTION

ABBREVIATION.....CONT SUCTI
ACTIVE.....Yes
DESCRIPTION.....Continuous Suction
DISPLAY_NAME.....Continuous Suction
HELP_TEXT.....No Help Text
ID.....{8C10DAC2-9AE5-40B9-99C7-25D47F4F13BA}
TERM.....CONTINUOUS SUCTION

DATA VALUE: CONTROL MODE VENTILATION

ABBREVIATION.....CMV
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Control Mode Ventilation
HELP_TEXT.....No Help Text
ID.....{2E79351C-BA31-FC04-B636-F39728B23EA6}
TERM.....CONTROL MODE VENTILATION

DATA VALUE: CONTUSION

ABBREVIATION.....Contusion
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Contusion
HELP_TEXT.....No Help Text
ID.....{5285197A-CC0A-4D12-A1F8-966765C687F1}
TERM.....CONTUSION

Page: 708

DATA VALUE: COOL

ABBREVIATION.....COOL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Cool
HELP_TEXT.....No Help Text
ID.....{48A1FEA8-FD75-DAF8-174A-9C6AB2044439}
TERM.....COOL

DATA VALUE: COOL CHEM PACK

ABBREVIATION.....CO CHEM PA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Cool Chem Pack
HELP_TEXT.....No Help Text
ID.....{6DE5E644-28C2-4190-A2EF-1A766E014A65}
TERM.....COOL CHEM PACK

DATA VALUE: COOLER THAN GENERAL BODY TEMPERATURE

ABBREVIATION.....CoolGenT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Cooler Than General Body Temperature
HELP_TEXT.....No Help Text
ID.....{03D7D4CF-9B3C-49CA-9650-213E9BABE064}
TERM.....COOLER THAN GENERAL BODY TEMPERATURE

DATA VALUE: COOLER THAN OPPOSITE BODY PART

ABBREVIATION.....COOL OPP P
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Cooler Than Opposite Body Part
HELP_TEXT.....No Help Text
ID.....{AFA2635C-8E80-4C71-A1B4-4AB117096FD2}
TERM.....COOLER THAN OPPOSITE BODY PART

DATA VALUE: COOLING BLANKET

ABBREVIATION.....COOL BLANK
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Cooling Blanket
HELP_TEXT.....No Help Text
ID.....{13A42B12-4959-4674-AF97-2D1242CFD934}
TERM.....COOLING BLANKET

Page: 709

DATA VALUE: COOLING OFF

ABBREVIATION.....COOLING OF
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Cooling Off
HELP_TEXT.....No Help Text
ID.....{14F573ED-0D10-4E0F-A1F3-A11F29D6DB11}
TERM.....COOLING OFF

DATA VALUE: COOLING ON

ABBREVIATION.....COOLING ON
ACTIVE.....Yes

DESCRIPTION.....No Description
DISPLAY_NAME.....Cooling On
HELP_TEXT.....No Help Text
ID.....{911CA17D-614E-445A-969E-4F5220E469EB}
TERM.....COOLING ON

DATA VALUE: COOPERATIVE

ABBREVIATION.....c
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Cooperative
HELP_TEXT.....No Help Text
ID.....{C174A06E-EBDE-4ADE-53DE-CF191D8AFC85}
TERM.....COOPERATIVE

DATA VALUE: COPIUS

ABBREVIATION.....COPIUS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Copius
HELP_TEXT.....No Help Text
ID.....{C74A6979-0E35-AA12-F992-69ACA3B0FA23}
TERM.....COPIUS

DATA VALUE: CRACKLES

ABBREVIATION.....c
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Crackles
HELP_TEXT.....No Help Text
ID.....{6B64C108-4424-00D9-525F-D89D549476D6}
TERM.....CRACKLES

Page: 710

DATA VALUE: CRAMPING

ABBREVIATION.....CRAMPING
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Cramping
HELP_TEXT.....No Help Text
ID.....{3EC3E352-5469-4A85-80E4-3AFEB61224C4}
TERM.....CRAMPING

DATA VALUE: CREPITUS

ABBREVIATION.....CREPITUS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Crepitus
HELP_TEXT.....No Help Text

ID.....{56A97506-05F0-4F37-8FD9-3AA82191408A}
TERM.....CREPITUS

DATA VALUE: CRUSHING

ABBREVIATION.....CRUSHING
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Crushing
HELP_TEXT.....No Help Text
ID.....{92D4BAF1-63E3-4FC0-B096-B5CBC4736245}
TERM.....CRUSHING

DATA VALUE: CRUTCH

ABBREVIATION.....CRUTCH
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Crutch
HELP_TEXT.....No Help Text
ID.....{2F040C2C-7A95-4FFA-8DAC-8B2EA9A90259}
TERM.....CRUTCH

DATA VALUE: CRUTCHES

ABBREVIATION.....CRUTCHES
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Crutches
HELP_TEXT.....No Help Text
ID.....{DDA81269-B25C-4F7D-A059-BB1602B2B179}
TERM.....CRUTCHES

Page: 711

DATA VALUE: CRYING

ABBREVIATION.....CRYING
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Crying
HELP_TEXT.....No Help Text
ID.....{D87AB50C-0B18-40DD-884C-6AAEE4DD97FC}
TERM.....CRYING

DATA VALUE: CRYSTALLOID

ABBREVIATION.....CRYSTALLOI
ACTIVE.....Yes
DESCRIPTION.....Crystalloid
DISPLAY_NAME.....Crystalloid
HELP_TEXT.....No Help Text
ID.....{4EC07E19-88B2-4C13-99B6-8233DF612997}
TERM.....CRYSTALLOID

DATA VALUE: CULTURAL

ABBREVIATION.....CULTURAL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Cultural
HELP_TEXT.....No Help Text
ID.....{530D166D-D7C1-4204-94DC-83F4E140F11C}
TERM.....CULTURAL

DATA VALUE: CYANOTIC

ABBREVIATION.....c
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Cyanotic
HELP_TEXT.....No Help Text
ID.....{06BD77DB-1849-D731-0F27-F8CE48755EA9}
TERM.....CYANOTIC

DATA VALUE: D - DOPPLER

ABBREVIATION.....d-d
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....D - Doppler
HELP_TEXT.....No Help Text
ID.....{2B77F099-0520-474D-A4BC-22BF824C8364}
TERM.....D - DOPPLER

Page: 712

DATA VALUE: D10W ENTERAL NUTRITION

ABBREVIATION.....D10W ENT N
ACTIVE.....Yes
DESCRIPTION.....D10W Enteral Nutrition Type
DISPLAY_NAME.....D10W Enteral Nutrition
HELP_TEXT.....No Help Text
ID.....{5049BCEF-CAE5-490D-9EB5-4A204A4D85E1}
TERM.....D10W ENTERAL NUTRITION

DATA VALUE: D50W ENTERAL NUTRITION

ABBREVIATION.....D50W ENT N
ACTIVE.....Yes
DESCRIPTION.....D50W Enteral Nutrition Type
DISPLAY_NAME.....D50W Enteral Nutrition
HELP_TEXT.....No Help Text
ID.....{528A4339-0261-4EC5-821C-7CE5592666DC}
TERM.....D50W ENTERAL NUTRITION

DATA VALUE: D5W ENTERAL NUTRITION

ABBREVIATION.....D5W

ACTIVE.....Yes
DESCRIPTION.....D5W Enteral Nutrition Type
DISPLAY_NAME.....D5W Enteral Nutrition
HELP_TEXT.....No Help Text
ID.....{D6664641-BF43-4157-9A30-8AD8CFE96518}
TERM.....D5W ENTERAL NUTRITION

DATA VALUE: D=DOPPLER

ABBREVIATION.....d
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....D=Doppler
HELP_TEXT.....No Help Text
ID.....{BCF37855-1E6E-9544-B316-D80A651697E4}
TERM.....D=DOPPLER

DATA VALUE: DAILY PLAN

ABBREVIATION.....DAILY PLAN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Daily Plan
HELP_TEXT.....No Help Text
ID.....{B29357A5-4294-4EA0-A0D3-59182F3C94FB}
TERM.....DAILY PLAN

Page: 713

DATA VALUE: DAMPENED

ABBREVIATION.....DAMPENED
ACTIVE.....Yes
DESCRIPTION.....Dampened
DISPLAY_NAME.....Dampened
HELP_TEXT.....No Help Text
ID.....{77C777B3-7E56-41E3-8E70-DE144E6A3424}
TERM.....DAMPENED

DATA VALUE: DANGLED

ABBREVIATION.....DANGLED
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Dangled
HELP_TEXT.....No Help Text
ID.....{40665750-112F-401D-AEFA-8E19F45EC0C9}
TERM.....DANGLED

DATA VALUE: DAY ROOM

ABBREVIATION.....dr
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Day Room

HELP_TEXT.....No Help Text
ID.....{269125EF-FA82-B17A-6B56-0CEAB4788899}
TERM.....DAY ROOM

DATA VALUE: DDD

ABBREVIATION.....DDD
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....DDD
HELP_TEXT.....No Help Text
ID.....{7AA7152E-48F0-8550-AD0D-9F298E90B765}
TERM.....DDD

DATA VALUE: DDDR

ABBREVIATION.....DDDR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....DDDR
HELP_TEXT.....No Help Text
ID.....{A3CC0D9D-D329-1A52-2DD1-911A1E04B871}
TERM.....DDDR

Page: 714

DATA VALUE: DDI

ABBREVIATION.....DDI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....DDI
HELP_TEXT.....No Help Text
ID.....{0C1E359F-B1F4-7994-D337-37622CB14BC8}
TERM.....DDI

DATA VALUE: DECELERATING

ABBREVIATION.....DECELERATI
ACTIVE.....Yes
DESCRIPTION.....Decelerating
DISPLAY_NAME.....Decelerating
HELP_TEXT.....Decelerating
ID.....{B259D4D4-82B8-454A-B53D-F12F4952BB61}
TERM.....DECELERATING

DATA VALUE: DECREASED

ABBREVIATION.....DECREASED
ACTIVE.....Yes
DESCRIPTION.....Decreased
DISPLAY_NAME.....Decreased
HELP_TEXT.....No Help Text
ID.....{D910FE54-A6AC-4BE3-AE35-C9E6679A93BD}
TERM.....DECREASED

DATA VALUE: DEEP

ABBREVIATION.....DEEP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Deep
HELP_TEXT.....No Help Text
ID.....{78AD569C-16D4-CFC3-6B6F-41FB9FFAB9D7}
TERM.....DEEP

DATA VALUE: DENIES HAVING SLEPT

ABBREVIATION.....DENIES SL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Denies Having Slept
HELP_TEXT.....No Help Text
ID.....{1AC4B0CF-EC41-426B-89AB-68D5A9BAFADA}
TERM.....DENIES HAVING SLEPT

Page: 715

DATA VALUE: DENIES NEED

ABBREVIATION.....DENIES NE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Denies Need
HELP_TEXT.....No Help Text
ID.....{28C8C733-5A28-4DDC-8D2E-F9D0E6CD9B28}
TERM.....DENIES NEED

DATA VALUE: DESCRIP

ABBREVIATION.....N
ACTIVE.....Yes
DESCRIPTION.....1.5 Cal Enteral Nutrition Formula
DISPLAY_NAME.....N
HELP_TEXT.....No Help Text
ID.....{47FD2D75-2C45-476F-82AF-7B40108E9C07}
TERM.....DESCRIP

DATA VALUE: DEVICE RELATED PRESSURE

ABBREVIATION.....DEV REL PR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Device Related Pressure
HELP_TEXT.....No Help Text
ID.....{3AED6107-A173-4B53-8EB5-2FB2B3336FC5}
TERM.....DEVICE RELATED PRESSURE

DATA VALUE: DEVICE RELEASE

ABBREVIATION.....drq
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Device Release
HELP_TEXT.....No Help Text
ID.....{9F97947F-CFFE-3285-F4C3-A8DF294DA46C}
TERM.....DEVICE RELEASE

DATA VALUE: DIABETIC/VASCULAR ULCER

ABBREVIATION.....DIA/VAS UL
ACTIVE.....Yes
DESCRIPTION.....Type of wound
DISPLAY_NAME.....Diabetic / Vascular Ulcer
HELP_TEXT.....No Help Text
ID.....{E209E53E-EFF1-48F5-9581-75236C89D9B9}
TERM.....DIABETIC/VASCULAR ULCER

Page: 716

DATA VALUE: DIALYSIS

ABBREVIATION.....DIALYSIS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Dialysis Line
HELP_TEXT.....No Help Text
ID.....{D218E23C-5A34-09B0-DE43-DC792DCD0F1C}
TERM.....DIALYSIS

DATA VALUE: DIAPHORETIC

ABBREVIATION.....d
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Diaphoretic
HELP_TEXT.....No Help Text
ID.....{F5F1E644-A1BD-0137-54B5-1C02678C1D5A}
TERM.....DIAPHORETIC

DATA VALUE: DIARRHEA

ABBREVIATION.....DIARRHEA
ACTIVE.....Yes
DESCRIPTION.....Diarrhea
DISPLAY_NAME.....Diarrhea
HELP_TEXT.....No Help Text
ID.....{353472EF-3499-46B1-AD7C-554FE01D8CE2}
TERM.....DIARRHEA

DATA VALUE: DIASTOLIC

ABBREVIATION.....DIASTOLIC
ACTIVE.....Yes
DESCRIPTION.....No Description

DISPLAY_NAME.....Diastolic
HELP_TEXT.....No Help Text
ID.....{3E9D1DCF-6B3B-482C-AE5F-442814F4B5C0}
TERM.....DIASTOLIC

DATA VALUE: DILATED

ABBREVIATION.....Dilated
ACTIVE.....Yes
DESCRIPTION.....Dilated
DISPLAY_NAME.....Dilated
HELP_TEXT.....No Help Text
ID.....{BF98ADC4-C69E-B17A-9BDE-D426D7F4C395}
TERM.....DILATED

Page: 717

DATA VALUE: DILUTE DAKINS SOLUTION

ABBREVIATION.....DL DAKINS
ACTIVE.....Yes
DESCRIPTION.....Dilute Dakins Solution
DISPLAY_NAME.....Dilute Dakins Solution
HELP_TEXT.....No Help Text
ID.....{66DEF4D0-54C0-4ADF-8FFD-F4E022C30F43}
TERM.....DILUTE DAKINS SOLUTION

DATA VALUE: DILUTE SODIUM HYPOCHLORITE SOLUTION

ABBREVIATION.....D NA HYPO
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Dilute Sodium Hypochlorite Solution
HELP_TEXT.....No Help Text
ID.....{14D2B1FB-1A5D-42AA-ACCF-A5033741DB61}
TERM.....DILUTE SODIUM HYPOCHLORITE SOLUTION

DATA VALUE: DIMINISHED

ABBREVIATION.....d
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Diminished
HELP_TEXT.....No Help Text
ID.....{49DE419E-DC3C-F7CD-771D-742031E0DEBC}
TERM.....DIMINISHED

DATA VALUE: DISCHARGE

ABBREVIATION.....DISCHARGE
ACTIVE.....Yes
DESCRIPTION.....Discharge
DISPLAY_NAME.....Discharge
HELP_TEXT.....No Help Text
ID.....{ADF67CC2-E05B-4B8D-A233-4AE1C35B2D20}

TERM.....DISCHARGE

DATA VALUE: DISCHARGE MEDS PLANNED

ABBREVIATION.....DI ME PL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Discharge Meds Planned
HELP_TEXT.....No Help Text
ID.....{04E9DF8B-8807-47A5-AE88-180F431F3BDE}
TERM.....DISCHARGE MEDS PLANNED

Page: 718

DATA VALUE: DISCHARGE ORDER ON CHART

ABBREVIATION.....DI OR CH
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Discharge Order On Chart
HELP_TEXT.....No Help Text
ID.....{F7267E73-07BC-44FB-AB0F-FAA07F6E2C1B}
TERM.....DISCHARGE ORDER ON CHART

DATA VALUE: DISCOLORATION

ABBREVIATION.....DISCOLORAT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Discoloration
HELP_TEXT.....No Help Text
ID.....{BBE3E910-CC49-430C-BD90-FB0A97F1B183}
TERM.....DISCOLORATION

DATA VALUE: DISCOMFORT

ABBREVIATION.....DISCOMFORT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Discomfort
HELP_TEXT.....No Help Text
ID.....{77D058C9-1EB4-96C0-C9E3-666C840B884C}
TERM.....DISCOMFORT

DATA VALUE: DISCONTINUED

ABBREVIATION.....DISCONTINU
ACTIVE.....Yes
DESCRIPTION.....Discontinued
DISPLAY_NAME.....Discontinued
HELP_TEXT.....No Help Text
ID.....{15E289FB-09DF-4731-9CF3-48A4E0031C31}
TERM.....DISCONTINUED

DATA VALUE: DISEASE PROCESS

ABBREVIATION.....DISEA PROC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Disease Process
HELP_TEXT.....No Help Text
ID.....{E80807AE-6208-421D-BE4C-445C470D32FA}
TERM.....DISEASE PROCESS

Page: 719

DATA VALUE: DISHEVELED

ABBREVIATION.....DISHEVEL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Disheveled
HELP_TEXT.....No Help Text
ID.....{0455549E-5CCB-470D-9D7F-F3959E6D55FE}
TERM.....DISHEVELED

DATA VALUE: DISORIENTED

ABBREVIATION.....DISORIENTE
ACTIVE.....Yes
DESCRIPTION.....Disoriented
DISPLAY_NAME.....Disoriented
HELP_TEXT.....No Help Text
ID.....{A7D1F2EA-CD7C-4100-8552-44834684C6DE}
TERM.....DISORIENTED

DATA VALUE: DISORIENTED TO PERSON

ABBREVIATION.....DISOR TO P
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Disoriented To Person
HELP_TEXT.....No Help Text
ID.....{E5D83C03-7EF5-770B-930F-7B3B3345048C}
TERM.....DISORIENTED TO PERSON

DATA VALUE: DISORIENTED TO PLACE

ABBREVIATION.....dtp
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Disoriented To Place
HELP_TEXT.....No Help Text
ID.....{00C98809-C0D0-C7B4-0615-8F7D31017FE8}
TERM.....DISORIENTED TO PLACE

DATA VALUE: DISORIENTED TO SITUATION

ABBREVIATION.....dts
ACTIVE.....Yes

DESCRIPTION.....No Description
DISPLAY_NAME.....Disoriented To Situation
HELP_TEXT.....No Help Text
ID.....{7143B3D7-5569-4378-D3AA-10C4675F0BA5}
TERM.....DISORIENTED TO SITUATION

Page: 720

DATA VALUE: DISORIENTED TO TIME

ABBREVIATION.....dtt
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Disoriented To Time
HELP_TEXT.....No Help Text
ID.....{2AD1A28A-4121-4B82-CFB8-E5A8D3095704}
TERM.....DISORIENTED TO TIME

DATA VALUE: DISTAL PORT

ABBREVIATION.....DIST PORT
ACTIVE.....Yes
DESCRIPTION.....Distal Port
DISPLAY_NAME.....Distal Port
HELP_TEXT.....No Help Text
ID.....{E493DF34-AD0E-4951-9715-5AD4EFA14036}
TERM.....DISTAL PORT

DATA VALUE: DISTENDED

ABBREVIATION.....DISTENDED
ACTIVE.....Yes
DESCRIPTION.....Distended
DISPLAY_NAME.....Distended
HELP_TEXT.....No Help Text
ID.....{B5C529CA-2BE0-4698-A153-718F4D58E643}
TERM.....DISTENDED

DATA VALUE: DIVERSION

ABBREVIATION.....DIVERSION
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Diversion
HELP_TEXT.....No Help Text
ID.....{087AAE9F-7FE1-5717-6423-07574C637875}
TERM.....DIVERSION

DATA VALUE: DO NOT RESUSCITATE

ABBREVIATION.....dnr
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Do Not Resuscitate
HELP_TEXT.....No Help Text

ID.....{0BD8A134-C563-7E00-A6B1-F7C543CEF6FF}
TERM.....DO NOT RESUSCITATE

Page: 721

DATA VALUE: DONE

ABBREVIATION.....DONE
ACTIVE.....Yes
DESCRIPTION.....Done
DISPLAY_NAME.....Done
HELP_TEXT.....No Help Text
ID.....{60412C8D-132F-46EA-8DDF-FCC282C4CC8B}
TERM.....DONE

DATA VALUE: DOO

ABBREVIATION.....DOO
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....DOO
HELP_TEXT.....No Help Text
ID.....{749BCE3D-FF09-B796-FF2A-ECD2D116C96F}
TERM.....DOO

DATA VALUE: DOUBLE

ABBREVIATION.....DOUBLE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Double
HELP_TEXT.....No Help Text
ID.....{A486852F-150B-24E3-8EB3-8CAABB3D0520}
TERM.....DOUBLE

DATA VALUE: DOWN

ABBREVIATION.....d
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Down
HELP_TEXT.....No Help Text
ID.....{36ECC86E-2D73-D4A9-5E13-9972E2ECF2E2}
TERM.....DOWN

DATA VALUE: DRESSING

ABBREVIATION.....d
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Dressing
HELP_TEXT.....No Help Text
ID.....{FFEF0714-3B5B-E2B4-D760-3C03D45E59EA}
TERM.....DRESSING

DATA VALUE: DRESSING CHANGED

ABBREVIATION.....DRESS CHAN
ACTIVE.....Yes
DESCRIPTION.....Changed Dressing
DISPLAY_NAME.....Dressing Changed
HELP_TEXT.....No Help Text
ID.....{E8CF4B7A-D1AB-49A2-8E7C-3375C3BC436F}
TERM.....DRESSING CHANGED

DATA VALUE: DRESSING WITH PRESCRIPTION RX

ABBREVIATION.....DR W RX
ACTIVE.....Yes
DESCRIPTION.....Dressing With Prescription Rx
DISPLAY_NAME.....Dressing With Prescription Rx
HELP_TEXT.....No Help Text
ID.....{AD5F9FDC-F280-4F1C-99E9-2BE0F9FFB15C}
TERM.....DRESSING WITH PRESCRIPTION RX

DATA VALUE: DRESSINGS CHECKED

ABBREVIATION.....DRESS CH
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Dressings Checked
HELP_TEXT.....No Help Text
ID.....{21014F53-29CE-4914-8760-6976D8DBA742}
TERM.....DRESSINGS CHECKED

DATA VALUE: DROPLET

ABBREVIATION.....DROPLET
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Droplet
HELP_TEXT.....No Help Text
ID.....{CE1DD314-8A30-4785-BD56-1540CDA85146}
TERM.....DROPLET

DATA VALUE: DROWSY

ABBREVIATION.....d
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Drowsy
HELP_TEXT.....No Help Text
ID.....{6B350814-A88A-0F05-C382-9F27830FBC09}
TERM.....DROWSY

DATA VALUE: DRY

ABBREVIATION.....d
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Dry
HELP_TEXT.....No Help Text
ID.....{D8F4F0FD-1B19-2515-3815-EC88D175B305}
TERM.....DRY

DATA VALUE: DRY AND INTACT

ABBREVIATION.....DRY INTACT
ACTIVE.....Yes
DESCRIPTION.....Dry and Intact
DISPLAY_NAME.....Dry and Intact
HELP_TEXT.....No Help Text
ID.....{77B62A18-F098-4764-86FA-66373E08E8C3}
TERM.....DRY AND INTACT

DATA VALUE: DRY SHAMPOO

ABBREVIATION.....DRY SHAMP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Dry Shampoo
HELP_TEXT.....No Help Text
ID.....{8910A186-2590-4352-A884-7D692A9D2B3C}
TERM.....DRY SHAMPOO

DATA VALUE: DRY STERILE DRESSING

ABBREVIATION.....DRY ST DRE
ACTIVE.....Yes
DESCRIPTION.....Dry Sterile Dressing
DISPLAY_NAME.....Dry Sterile Dressing
HELP_TEXT.....No Help Text
ID.....{3069DFCE-1B59-467C-97CE-9DBAF9AF42CC}
TERM.....DRY STERILE DRESSING

DATA VALUE: DULL

ABBREVIATION.....DULL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Dull
HELP_TEXT.....No Help Text
ID.....{F0F80976-F6A0-4363-9C75-AAAE9937DB49}
TERM.....DULL

Page: 724

DATA VALUE: DUMPING SYNDROME

ABBREVIATION.....DUMP SYND

ACTIVE.....Yes
DESCRIPTION.....Dumping Syndrome Diet
DISPLAY_NAME.....Dumping Syndrome Diet
HELP_TEXT.....No Help Text
ID.....{0173C971-63CB-45A5-9E93-D966EB3B957D}
TERM.....DUMPING SYNDROME

DATA VALUE: DUODENOSTOMY

ABBREVIATION.....DUODENOSTO
ACTIVE.....Yes
DESCRIPTION.....Duodenostomy
DISPLAY_NAME.....Duodenostomy
HELP_TEXT.....No Help Text
ID.....{FA44135D-6806-42E7-B371-8C3912B3582E}
TERM.....DUODENOSTOMY

DATA VALUE: DUSKY

ABBREVIATION.....d
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Dusky
HELP_TEXT.....No Help Text
ID.....{2ABF845E-4E2F-883F-3FCC-F72B8507FF85}
TERM.....DUSKY

DATA VALUE: DVI

ABBREVIATION.....DVI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....DVI
HELP_TEXT.....No Help Text
ID.....{4D6E3081-7BB5-45DB-E8E6-0D3D7869E6BD}
TERM.....DVI

DATA VALUE: DYE PRESENT

ABBREVIATION.....DYE PRESEN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Dye Present
HELP_TEXT.....No Help Text
ID.....{E865AD3E-6DAD-C3F9-6F93-7AD6F71234A4}
TERM.....DYE PRESENT

Page: 725

DATA VALUE: DYSPHAGIA PHASE I

ABBREVIATION.....DYSP PH I
ACTIVE.....Yes
DESCRIPTION.....Dysphagia Stage I Diet
DISPLAY_NAME.....Dysphagia Stage I Diet

HELP_TEXT.....No Help Text
ID.....{8011795B-77A5-4FF4-8EB9-EB6E1AA1FD4D}
TERM.....DYSPHAGIA PHASE I

DATA VALUE: DYSPHAGIA STAGE II

ABBREVIATION.....DYSP ST II
ACTIVE.....Yes
DESCRIPTION.....Dysphagia Stage II Diet
DISPLAY_NAME.....Dysphagia Stage II Diet
HELP_TEXT.....No Help Text
ID.....{294BF58A-47FD-43A1-8A1E-9E348987B004}
TERM.....DYSPHAGIA STAGE II

DATA VALUE: DYSPHAGIA STAGE III

ABBREVIATION.....DYS ST III
ACTIVE.....Yes
DESCRIPTION.....Dysphagia Stage III Diet
DISPLAY_NAME.....Dysphagia Stage III Diet
HELP_TEXT.....No Help Text
ID.....{DCF91830-C5CA-47B2-A7C6-4DD971F953EA}
TERM.....DYSPHAGIA STAGE III

DATA VALUE: DYSPHAGIA STAGE IV

ABBREVIATION.....DYSP ST IV
ACTIVE.....Yes
DESCRIPTION.....Dysphagia Stage IV Diet
DISPLAY_NAME.....Dysphagia Stage IV Diet
HELP_TEXT.....No Help Text
ID.....{E738569A-5104-4C93-8EFC-47D80765F739}
TERM.....DYSPHAGIA STAGE IV

DATA VALUE: DYSPHAGIA STAGE V

ABBREVIATION.....DYS PH ST V
ACTIVE.....Yes
DESCRIPTION.....Dysphagia Stage V Diet
DISPLAY_NAME.....Dysphagia Stage V Diet
HELP_TEXT.....No Help Text
ID.....{ABE7D85D-1821-4592-B4E8-84670EE484E4}
TERM.....DYSPHAGIA STAGE V

Page: 726

DATA VALUE: DYSPNEIC

ABBREVIATION.....DYSPENIC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Dyspenic
HELP_TEXT.....No Help Text
ID.....{4C6B42D9-694F-5EB2-6396-7ED85D5B154E}
TERM.....DYSPNEIC

DATA VALUE: DYSURIA

ABBREVIATION.....DYSURIA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Dysuria
HELP_TEXT.....No Help Text
ID.....{75AC5ECD-8975-406A-30C7-E2EB23F10F84}
TERM.....DYSURIA

DATA VALUE: EAGER TO LEARN

ABBREVIATION.....EAGER
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Eager To Learn
HELP_TEXT.....No Help Text
ID.....{A20AF51B-C89C-4308-AE17-C66144FE68BE}
TERM.....EAGER TO LEARN

DATA VALUE: ECCHYMOSIS

ABBREVIATION.....ECCHYMOSIS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Ecchymosis
HELP_TEXT.....No Help Text
ID.....{DF076906-CC17-6AD1-76E7-B07678A126DE}
TERM.....ECCHYMOSIS

DATA VALUE: ELASTIC WRAP

ABBREVIATION.....ELAST WRAP
ACTIVE.....Yes
DESCRIPTION.....ELASTIC WRAP
DISPLAY_NAME.....Elastic Wrap
HELP_TEXT.....ELASTIC WRAP
ID.....{9EB52FA5-32D8-823F-2E3C-83375DA62342}
TERM.....ELASTIC WRAP

Page: 727

DATA VALUE: ELBOW PAD

ABBREVIATION.....ELBOW PAD
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Elbow Pad
HELP_TEXT.....No Help Text
ID.....{26942F81-1E33-4962-B772-A96FFBCED814}
TERM.....ELBOW PAD

DATA VALUE: ELECTRIC

ABBREVIATION.....ELECTRIC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....ELECTRIC
HELP_TEXT.....No Help Text
ID.....{0C221385-61AF-4660-A129-35E610EA71AB}
TERM.....ELECTRIC

DATA VALUE: ELEMENTAL FORMULA

ABBREVIATION.....ELE FORM
ACTIVE.....Yes
DESCRIPTION.....Elemental Enteral Nutrition Formula
DISPLAY_NAME.....Elemental Enteral Nutrition Formula
HELP_TEXT.....No Help Text
ID.....{DCE1E81C-1AC7-48F4-97AA-F1FC4D0601D8}
TERM.....ELEMENTAL FORMULA

DATA VALUE: ELOPEMENT

ABBREVIATION.....e
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Elopement
HELP_TEXT.....No Help Text
ID.....{EB030907-FB90-CA9D-CF50-E26C5D9FE671}
TERM.....ELOPEMENT

DATA VALUE: EMERGENCY DEPARTMENT

ABBREVIATION.....ER DEPT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Emergency Department
HELP_TEXT.....No Help Text
ID.....{D2E511F7-8B80-4928-A64B-3F1EE17395F7}
TERM.....EMERGENCY DEPARTMENT

Page: 728

DATA VALUE: EMOTIONAL

ABBREVIATION.....EMOTIONAL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Emotional
HELP_TEXT.....No Help Text
ID.....{626C22E6-EF2A-489F-AA79-D27541A0BC45}
TERM.....EMOTIONAL

DATA VALUE: END OF LIFE CARE

ABBREVIATION.....END LIFE
ACTIVE.....Yes
DESCRIPTION.....No Description

DISPLAY_NAME.....End Of Life Care
HELP_TEXT.....No Help Text
ID.....{047D8977-0567-4B90-AD6A-3AAF8DB35238}
TERM.....END OF LIFE CARE

DATA VALUE: ENVIRONMENT

ABBREVIATION.....ENVIRO
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Environment
HELP_TEXT.....No Help Text
ID.....{C61AA262-F6AF-46AC-A215-0D56FD01B234}
TERM.....ENVIRONMENT

DATA VALUE: ENZYMATIC OINTMENT

ABBREVIATION.....ENZY OINTM
ACTIVE.....Yes
DESCRIPTION.....Enzymatic Ointment
DISPLAY_NAME.....Enzymatic Ointment
HELP_TEXT.....No Help Text
ID.....{E781DA75-8777-4884-A093-FB623774DB8F}
TERM.....ENZYMATIC OINTMENT

DATA VALUE: EPICARDIAL LEADS

ABBREVIATION.....el
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Epicardial Leads
HELP_TEXT.....No Help Text
ID.....{2B551575-3D50-96BD-7BAA-50CC9796361E}
TERM.....EPICARDIAL LEADS

Page: 729

DATA VALUE: EPITHELIALIZATION

ABBREVIATION.....EPITHELIAL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Epithelialization
HELP_TEXT.....No Help Text
ID.....{5FC43A6D-921D-4736-A234-9993DB6499B7}
TERM.....EPITHELIALIZATION

DATA VALUE: EQUAL

ABBREVIATION.....Equal
ACTIVE.....Yes
DESCRIPTION.....Equal
DISPLAY_NAME.....Equal
HELP_TEXT.....No Help Text
ID.....{AC079584-158D-EBF2-E13A-C2E919375E4D}

TERM.....EQUAL

DATA VALUE: EQUIPMENT

ABBREVIATION.....EQUIPMENT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Equipment
HELP_TEXT.....No Help Text
ID.....{04C0F805-C7B3-4A9A-91D8-0F32F3265874}
TERM.....EQUIPMENT

DATA VALUE: ERYTHEMA

ABBREVIATION.....ERYTHEMA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Erythema
HELP_TEXT.....No Help Text
ID.....{0083E1B3-AE2B-4CCF-8794-86E612B9DCEA}
TERM.....ERYTHEMA

DATA VALUE: ESCAPE

ABBREVIATION.....e
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Escape
HELP_TEXT.....No Help Text
ID.....{965F3201-9D63-66AB-FA25-99479AE7145A}
TERM.....ESCAPE

Page: 730

DATA VALUE: ESCHAR

ABBREVIATION.....ESCHAR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Eschar
HELP_TEXT.....No Help Text
ID.....{555E03C5-FBFE-4EB9-816F-25CE34476B3C}
TERM.....ESCHAR

DATA VALUE: ETT REPOSITIONED

ABBREVIATION.....ETT REPOSI
ACTIVE.....Yes
DESCRIPTION.....ETT Repositioned
DISPLAY_NAME.....ETT Repositioned
HELP_TEXT.....No Help Text
ID.....{95B35EDB-76B6-433E-B526-01C543A7CCEA}
TERM.....ETT REPOSITIONED

DATA VALUE: ETT RETAPED/SECURED

ABBREVIATION.....ETT RETAPE
ACTIVE.....Yes
DESCRIPTION.....ETT Retaped / Secured
DISPLAY_NAME.....ETT Retaped / Secured
HELP_TEXT.....No Help Text
ID.....{6D8A2EDA-CAC9-44D2-BDA9-3B7F20908F1E}
TERM.....ETT RETAPED/SECURED

DATA VALUE: EWALD TUBE

ABBREVIATION.....EWALD TUBE
ACTIVE.....Yes
DESCRIPTION.....Ewald Tube
DISPLAY_NAME.....Ewald Tube
HELP_TEXT.....No Help Text
ID.....{B83AB4C1-A098-4FAB-A645-61FE1FEADBBE}
TERM.....EWALD TUBE

DATA VALUE: EXPIRATORY WHEEZES

ABBREVIATION.....EXP WHEEZE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Expiratory Wheezes
HELP_TEXT.....No Help Text
ID.....{CFC103A8-8B02-4942-BBDD-8D24C930883E}
TERM.....EXPIRATORY WHEEZES

Page: 731

DATA VALUE: EXTENDS AROUND INCISION EXTENDS INTO PERIWOUND LESS THAN 1 CM

ABBREVIATION.....ER EXI <1
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Less Than 1 cm
HELP_TEXT.....No Help-Text
ID.....{095709D0-6190-49A3-AB08-BF4575DC8396}
TERM.....EXTENDS AROUND INCISION EXTENDS INTO
PERIWOUND LESS THAN 1 CM

DATA VALUE: EXTENDS AROUND INCISION MORE THAN 1 CM

ABBREVIATION.....ER EXI >1
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Extends Into Periwound More Than 1 cm
HELP_TEXT.....No Help-Text
ID.....{11F13A1F-7AC2-4629-A960-CD92655EE101}
TERM.....EXTENDS AROUND INCISION MORE THAN 1 CM

DATA VALUE: EXTENDS AROUND INCISION NONE

ABBREVIATION.....ER EXI NON

ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....None
HELP_TEXT.....No Help-Text
ID.....{91F1A7C5-885C-4B13-9F7F-189C93E505CC}
TERM.....EXTENDS AROUND INCISION NONE

DATA VALUE: EXTENDS INTO PERIWOUND <1CM

ABBREVIATION.....EX IN PER<
ACTIVE.....Yes
DESCRIPTION.....Extends Into Periwound <1cm
DISPLAY_NAME.....Extends Into Periwound <1cm
HELP_TEXT.....No Help Text
ID.....{0EBCEAE0-E17A-484F-BFE2-51F17AA761CB}
TERM.....EXTENDS INTO PERIWOUND <1CM

DATA VALUE: EXTENDS INTO PERIWOUND>1CM

ABBREVIATION.....EX IN PER>
ACTIVE.....Yes
DESCRIPTION.....Extends Into Periwound > 1cm
DISPLAY_NAME.....Extends Into Periwound > 1cm
HELP_TEXT.....No Help Text
ID.....{1F31FE9B-DA19-434C-8C58-67CDC5385231}
TERM.....EXTENDS INTO PERIWOUND>1CM

Page: 732

DATA VALUE: EXTENDS TO PAIN

ABBREVIATION.....EXT TO PAI
ACTIVE.....Yes
DESCRIPTION.....Extends To Pain
DISPLAY_NAME.....Extends To Pain
HELP_TEXT.....No Help Text
ID.....{A62B2C4D-BDB7-49FA-9669-DF6732063DA6}
TERM.....EXTENDS TO PAIN

DATA VALUE: EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND LESS THAN 1 CM

ABBREVIATION.....ERY EX <1
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Extends Into Periwound Less Than 1 cm
HELP_TEXT.....No Help-Text
ID.....{48C7FA69-D31A-4AAA-84F3-79AC863B71FD}
TERM.....EXTENSION AROUND WOUND EXTENDS INTO
PERIWOUND LESS THAN 1 CM

DATA VALUE: EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND MORE THAN 1 CM

ABBREVIATION.....ERY EX >1
ACTIVE.....Yes
DESCRIPTION.....No Description

DISPLAY_NAME.....Extends Into Periwound Greater Than 1
cm
HELP_TEXT.....No Help-Text
ID.....{175B2505-C596-4036-8FC3-DFEA8EFB76F9}
TERM.....EXTENSION AROUND WOUND EXTENDS INTO
PERIWOUND MORE THAN 1 CM

DATA VALUE: EXTENSION AROUND WOUND LIMITED TO WOUND

ABBREVIATION.....ERY EX LIM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Limited To Wound
HELP_TEXT.....No Help-Text
ID.....{2DC83B85-B3B6-4679-9F48-7DD26EFFF4D4}
TERM.....EXTENSION AROUND WOUND LIMITED TO WOUND

DATA VALUE: EXTENSION AROUND WOUND NONE

ABBREVIATION.....ERY EX NON
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....None
HELP_TEXT.....No Help-Text
ID.....{815C7BDF-C9A8-4AFE-B2BE-2C2FF24AA495}
TERM.....EXTENSION AROUND WOUND NONE

Page: 733

DATA VALUE: EXTENSIVE ASSISTANCE

ABBREVIATION.....EXT ASSIST
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Extensive Assistance
HELP_TEXT.....No Help Text
ID.....{C402B659-E880-4070-8345-13E5508B4D29}
TERM.....EXTENSIVE ASSISTANCE

DATA VALUE: EXTERNAL CATHETER

ABBREVIATION.....EXT. CATHE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....External Catheter
HELP_TEXT.....No Help Text
ID.....{C8AEE936-DD91-FE91-4DD6-6F210FB59726}
TERM.....EXTERNAL CATHETER

DATA VALUE: EXTREMELY DRY

ABBREVIATION.....EXTR DRY
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Extremely Dry

HELP_TEXT.....No Help Text
ID.....{51677B19-80B9-4787-B6E2-E0BCF6601AA1}
TERM.....EXTREMELY DRY

DATA VALUE: EYE CLOSED

ABBREVIATION.....ec
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Eye Closed
HELP_TEXT.....No Help Text
ID.....{140B7D2F-F9C2-0B4C-3CDC-7355EDDB5030}
TERM.....EYE CLOSED

DATA VALUE: FACE WASH

ABBREVIATION.....FACE WASH
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Face Wash
HELP_TEXT.....No Help Text
ID.....{E1162A8B-D26D-4EB3-988F-50583B1849EA}
TERM.....FACE WASH

Page: 734

DATA VALUE: FAMILY

ABBREVIATION.....FAMILY
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Family
HELP_TEXT.....No Help Text
ID.....{1233C4BB-5AC5-4883-AF2A-F5C93CD8814F}
TERM.....FAMILY

DATA VALUE: FAMILY INFORMED

ABBREVIATION.....FAM INFORM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Family Informed
HELP_TEXT.....No Help Text
ID.....{77135D32-9CDD-4D3F-9B62-6E8189B77844}
TERM.....FAMILY INFORMED

DATA VALUE: FAMILY/FRIEND

ABBREVIATION.....FAM/FRIE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Family/Friend
HELP_TEXT.....No Help Text
ID.....{FA8FE272-E659-41EF-A394-F907A5F55ED8}
TERM.....FAMILY/FRIEND

DATA VALUE: FASCIA

ABBREVIATION.....FASCIA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Fascia
HELP_TEXT.....No Help Text
ID.....{EEE8535C-B887-4A34-A93F-C1C4AFC1296F}
TERM.....FASCIA

DATA VALUE: FAT RESTRICTED 50 GM DIET

ABBREVIATION.....F RES 50gm
ACTIVE.....Yes
DESCRIPTION.....Fat Restricted 50 gm Diet
DISPLAY_NAME.....Fat Restricted 50 gm Diet
HELP_TEXT.....No Help Text
ID.....{8E97A639-B7A6-4B3B-9BEF-0A46F3B46E43}
TERM.....FAT RESTRICTED 50 GM DIET

Page: 735

DATA VALUE: FECAL FAT TEST DIET

ABBREVIATION.....FEC F TEST
ACTIVE.....Yes
DESCRIPTION.....Fecal Fat Test Diet
DISPLAY_NAME.....Fecal Fat Test Diet
HELP_TEXT.....No Help Text
ID.....{BA1E8E4A-7698-4EFB-B1DD-18C444EE82D5}
TERM.....FECAL FAT TEST DIET

DATA VALUE: FHWN1

ABBREVIATION.....FHWN1
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....FHWN1
HELP_TEXT.....No Help Text
ID.....{93B28875-51BF-4294-B2F9-58BE9365DC81}
TERM.....FHWN1

DATA VALUE: FIBER FORMULA

ABBREVIATION.....FIBER FORM
ACTIVE.....Yes
DESCRIPTION.....Fiber Enteral Nutrition Formula
DISPLAY_NAME.....Fiber Enteral Nutrition Formula
HELP_TEXT.....No Help Text
ID.....{1AF14DA4-8D59-420B-9895-4256E7A818BE}
TERM.....FIBER FORMULA

DATA VALUE: FIBRINOUS

ABBREVIATION.....FIBRINOUS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Fibrinous
HELP_TEXT.....No Help Text
ID.....{D286EAE2-6902-4B7E-845D-BBB316941B2C}
TERM.....FIBRINOUS

DATA VALUE: FIRM

ABBREVIATION.....FIRM
ACTIVE.....Yes
DESCRIPTION.....Firm
DISPLAY_NAME.....Firm
HELP_TEXT.....No Help Text
ID.....{11E23F20-39A0-4A6D-B198-04BB36B0C909}
TERM.....FIRM

Page: 736

DATA VALUE: FIRST DEGREE AV BLOCK

ABBREVIATION.....fdab
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....First Degree Av Block
HELP_TEXT.....No Help Text
ID.....{4FF046C9-B638-F86C-5D51-006B3D5225E3}
TERM.....FIRST DEGREE AV BLOCK

DATA VALUE: FIXED

ABBREVIATION.....f
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Fixed
HELP_TEXT.....No Help Text
ID.....{697D8C5A-976F-F423-30EC-A04102CF85F0}
TERM.....FIXED

DATA VALUE: FLAT

ABBREVIATION.....Flat
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Flat
HELP_TEXT.....No Help Text
ID.....{EA4F1980-EDF3-82BC-B677-D747F48F1163}
TERM.....FLAT

DATA VALUE: FLATULENCE

ABBREVIATION.....FLATULENCE
ACTIVE.....Yes
DESCRIPTION.....Flatulence

DISPLAY_NAME.....Flatulence
HELP_TEXT.....No Help Text
ID.....{F6ECE7D8-D9C8-485C-8FAC-03DA9D7E0B66}
TERM.....FLATULENCE

DATA VALUE: FLEXION TO PAIN

ABBREVIATION.....FLEX TO PA
ACTIVE.....Yes
DESCRIPTION.....Flexion To Pain
DISPLAY_NAME.....Flexion To Pain
HELP_TEXT.....No Help Text
ID.....{B9F5046F-4F8A-4A20-9D47-90432DC8FEC2}
TERM.....FLEXION TO PAIN

Page: 737

DATA VALUE: FLICKER OF CONTRACTION

ABBREVIATION.....foc
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Flicker Of Contraction
HELP_TEXT.....No Help Text
ID.....{218B322B-84AA-F6F6-5207-D86F9A4311F0}
TERM.....FLICKER OF CONTRACTION

DATA VALUE: FLOTATION BED

ABBREVIATION.....FLOT BED
ACTIVE.....Yes
DESCRIPTION.....Flotation Bed
DISPLAY_NAME.....Flotation Bed
HELP_TEXT.....No Help Text
ID.....{FA6329AB-85A6-4301-BC14-D0D089249356}
TERM.....FLOTATION BED

DATA VALUE: FLUCTUANCE

ABBREVIATION.....FLUCTUANCE
ACTIVE.....Yes
DESCRIPTION.....Fluctuance
DISPLAY_NAME.....Fluctuance
HELP_TEXT.....No Help Text
ID.....{A6AD2660-621D-4F39-8485-A013D7491566}
TERM.....FLUCTUANCE

DATA VALUE: FLUCTUATES

ABBREVIATION.....Fluctuates
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Fluctuates
HELP_TEXT.....No Help Text
ID.....{C2697650-EF53-5F2E-EEE3-7B2E298F2FA0}

TERM.....FLUCTUATES

DATA VALUE: FLUSHED

ABBREVIATION.....nt
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Flushed
HELP_TEXT.....No Help Text
ID.....{C7C54C1A-6E75-4BCF-4C57-0F8622A5BD16}
TERM.....FLUSHED

Page: 738

DATA VALUE: FOAM AIR HEEL ELEVATOR

ABBREVIATION.....FO AI HEEL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Foam Air Heel Elevator
HELP_TEXT.....No Help Text
ID.....{4EC26AA3-4577-4808-8728-F360C5C5A0D4}
TERM.....FOAM AIR HEEL ELEVATOR

DATA VALUE: FOAM DRESSING

ABBREVIATION.....WOU DR FOA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Foam Dressing
HELP_TEXT.....No Help-Text
ID.....{D18A5EDD-0EF6-470E-8370-4EDB6D322477}
TERM.....FOAM DRESSING

DATA VALUE: FOLEY

ABBREVIATION.....FOLEY
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Foley
HELP_TEXT.....No Help Text
ID.....{9FE7D3A7-67CF-CF54-F721-2F8035CE2BF1}
TERM.....FOLEY

DATA VALUE: FOOT CARE

ABBREVIATION.....FOOT CARE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Foot Care
HELP_TEXT.....No Help Text
ID.....{332E8338-6DAE-4BC9-B195-63DAD57937D0}
TERM.....FOOT CARE

DATA VALUE: FOOT DOWN

ABBREVIATION.....fd
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Foot Down
HELP_TEXT.....No Help Text
ID.....{13F73BE9-2FA8-FE88-443B-9BA3FE6EDFD7}
TERM.....FOOT DOWN

Page: 739

DATA VALUE: FOOT UP

ABBREVIATION.....fu
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Foot Up
HELP_TEXT.....No Help Text
ID.....{49CCC59E-002E-8387-E8D6-D57AAC089EF3}
TERM.....FOOT UP

DATA VALUE: FOUL SMELLING

ABBREVIATION.....FOUL SMELL
ACTIVE.....Yes
DESCRIPTION.....Foul Smelling
DISPLAY_NAME.....Foul Smelling
HELP_TEXT.....No Help Text
ID.....{761346A9-329D-4F42-8D80-0E9D8D40653C}
TERM.....FOUL SMELLING

DATA VALUE: FOWLER'S

ABBREVIATION.....FOWLER'S
ACTIVE.....Yes
DESCRIPTION.....Fowler's
DISPLAY_NAME.....Fowler's
HELP_TEXT.....No Help Text
ID.....{8E691813-A9AC-48BE-9401-528A7D762A32}
TERM.....FOWLER'S

DATA VALUE: FRANK BLOOD

ABBREVIATION.....FRANK BLD
ACTIVE.....Yes
DESCRIPTION.....Frank Blood
DISPLAY_NAME.....Frank Blood
HELP_TEXT.....No Help Text
ID.....{CFAE6E6B-87B2-4CBA-A280-C5233B7995ED}
TERM.....FRANK BLOOD

DATA VALUE: FREQUENT VERBAL REMINDERS

ABBREVIATION.....FR VER RE
ACTIVE.....Yes

DESCRIPTION.....No Description
DISPLAY_NAME.....Frequent Verbal Reminders
HELP_TEXT.....No Help Text
ID.....{F30D791D-BB80-4AB3-B37F-22B6B547EA65}
TERM.....FREQUENT VERBAL REMINDERS

Page: 740

DATA VALUE: FRIEND

ABBREVIATION.....FRIEND
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Friend
HELP_TEXT.....No Help Text
ID.....{03902E48-AF86-432A-8372-37796A55F625}
TERM.....FRIEND

DATA VALUE: FULL CODE

ABBREVIATION.....fc
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Full Code
HELP_TEXT.....No Help Text
ID.....{1DED65E6-2F2E-8C5F-D009-BE6467A9A2AA}
TERM.....FULL CODE

DATA VALUE: FULL LIQUID DIABETIC

ABBREVIATION.....FULL L DIA
ACTIVE.....Yes
DESCRIPTION.....Full Liquid Diabetic
DISPLAY_NAME.....Full Liquid Diabetic
HELP_TEXT.....No Help Text
ID.....{B89FAE3C-07EB-4AF6-B123-1153ADADC79F}
TERM.....FULL LIQUID DIABETIC

DATA VALUE: FULL LIQUID DIET

ABBREVIATION.....FULL LIQ D
ACTIVE.....Yes
DESCRIPTION.....Full Liquid Diet
DISPLAY_NAME.....Full Liquid Diet
HELP_TEXT.....No Help Text
ID.....{CAC22234-7B69-44C3-8A0E-9A84A684E429}
TERM.....FULL LIQUID DIET

DATA VALUE: FULL SPINAL PRECAUTIONS

ABBREVIATION.....FULL S P
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Full Spinal Precautions
HELP_TEXT.....No Help Text

ID.....{42ADB74D-12CC-4111-B181-626C0F01BC7C}
TERM.....FULL SPINAL PRECAUTIONS

Page: 741

DATA VALUE: FULL THICK SKIN LOSS INVOLVING MUSCLE, BONE, SUPPO

ABBREVIATION.....F SK L MBS
ACTIVE.....Yes
DESCRIPTION.....Full Thickness Skin Loss Involving
Muscle, Bone, or Supporting Tissue
DISPLAY_NAME.....Full Thickness Skin Loss Involving
Muscle, Bone, or Supporting Tissue
HELP_TEXT.....No Help Text
ID.....{E027FBEF-1C83-41B7-B848-26D9D62F4266}
TERM.....FULL THICK SKIN LOSS INVOLVING
MUSCLE, BONE, SUPPO

DATA VALUE: FULL THICK SKIN LOSS INVOLVING SUBCUTANEOUS TISSUE

ABBREVIATION.....FULL SKINL
ACTIVE.....Yes
DESCRIPTION.....Full Thickness Skin Loss Involving
Subcutaneous Tissue
DISPLAY_NAME.....Full Thickness Skin Loss Involving
Subcutaneous Tissue
HELP_TEXT.....No Help Text
ID.....{14C31B14-1D35-4CCB-81A3-C6E30801950D}
TERM.....FULL THICK SKIN LOSS INVOLVING
SUBCUTANEOUS TISSUE

DATA VALUE: GAIT BELT

ABBREVIATION.....GAIT BELT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Gait Belt
HELP_TEXT.....No Help Text
ID.....{AD9AB19B-92B6-40BE-B906-FBFC46C39E0D}
TERM.....GAIT BELT

DATA VALUE: GASTRIC BYPASS CLEAR LIQUID

ABBREVIATION.....GAS BYP LI
ACTIVE.....Yes
DESCRIPTION.....Gastric Bypass Clear Liquid
DISPLAY_NAME.....Gastric Bypass Clear Liquid
HELP_TEXT.....No Help Text
ID.....{6A87E820-7BEA-4636-B102-F9CB83BDD4D9}
TERM.....GASTRIC BYPASS CLEAR LIQUID

DATA VALUE: GASTRIC BYPASS PUREED

ABBREVIATION.....GAS BY PUR
ACTIVE.....Yes

DESCRIPTION.....Gastric Bypass Pureed
DISPLAY_NAME.....Gastric Bypass Pureed
HELP_TEXT.....No Help Text
ID.....{9F9EEB7F-7211-45A2-BCDE-B8398596CBF6}
TERM.....GASTRIC BYPASS PUREED

Page: 742

DATA VALUE: GASTRIC BYPASS REGULAR

ABBREVIATION.....GAS BY REG
ACTIVE.....Yes
DESCRIPTION.....Gastric Bypass Regular
DISPLAY_NAME.....Gastric Bypass Regular
HELP_TEXT.....No Help Text
ID.....{EE831906-EA51-4EE4-AAF2-583DF366A942}
TERM.....GASTRIC BYPASS REGULAR

DATA VALUE: GASTRIC TUBE

ABBREVIATION.....G TUBE
ACTIVE.....Yes
DESCRIPTION.....Gastric Feeding Tube
DISPLAY_NAME.....G Tube
HELP_TEXT.....No Help Text
ID.....{65F774BA-F977-46D3-AE41-7F7C74CFF47E}
TERM.....GASTRIC TUBE

DATA VALUE: GASTROSTOMY FEEDING BUTTON

ABBREVIATION.....GAST F BUT
ACTIVE.....Yes
DESCRIPTION.....Gastrostomy Feeding Button
DISPLAY_NAME.....Gastrostomy Feeding Button
HELP_TEXT.....No Help Text
ID.....{1627B53B-374B-42B4-BBFE-463EE37B6110}
TERM.....GASTROSTOMY FEEDING BUTTON

DATA VALUE: GASTROSTOMY TUBE

ABBREVIATION.....GAST TUBE
ACTIVE.....Yes
DESCRIPTION.....Gastrostomy Tube
DISPLAY_NAME.....Gastrostomy Tube
HELP_TEXT.....No Help Text
ID.....{EF017A39-745E-4C29-821A-5B25F103DF99}
TERM.....GASTROSTOMY TUBE

DATA VALUE: GAUZE

ABBREVIATION.....GAUZE
ACTIVE.....Yes
DESCRIPTION.....Gauze
DISPLAY_NAME.....Gauze
HELP_TEXT.....No Help Text

ID.....{CF10C658-4334-4F3B-A38E-E98894458333}
TERM.....GAUZE

Page: 743

DATA VALUE: GAUZE AND TAPE

ABBREVIATION.....GAU & TAPE
ACTIVE.....Yes
DESCRIPTION.....Gauze And Tape Dressing
DISPLAY_NAME.....Gauze and Tape
HELP_TEXT.....No Help Text
ID.....{5BC40405-B5C4-47CF-81C3-520732E100AC}
TERM.....GAUZE AND TAPE

DATA VALUE: GAUZE PACKING

ABBREVIATION.....GAUZE PACK
ACTIVE.....Yes
DESCRIPTION.....Gauze Packing
DISPLAY_NAME.....Gauze Packing
HELP_TEXT.....No Help Text
ID.....{6C093717-82CB-4966-929E-C909310649C1}
TERM.....GAUZE PACKING

DATA VALUE: GAUZE WRAP

ABBREVIATION.....GAUZE WRAP
ACTIVE.....Yes
DESCRIPTION.....Gauze Wrap
DISPLAY_NAME.....Gauze Wrap
HELP_TEXT.....No Help Text
ID.....{040A9648-E6E1-4AC5-93C7-53A49537B15C}
TERM.....GAUZE WRAP

DATA VALUE: GERI-CHAIR

ABBREVIATION.....g
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Geri-Chair
HELP_TEXT.....No Help Text
ID.....{167F9786-8F24-054F-020A-658A6EE5B0B4}
TERM.....GERI-CHAIR

DATA VALUE: GIVEN DISCHARGE INSTRUCTION AND VERBALIZES UNDERSTANDING

ABBREVIATION.....DI IN UN
ACTIVE.....No
DESCRIPTION.....No Description
DISPLAY_NAME.....Given Discharge Instructions and
Verbalizes Understanding
HELP_TEXT.....No Help Text
ID.....{C9A61418-88B2-42C4-B5F4-FC3CF76AB39D}
TERM.....GIVEN DISCHARGE INSTRUCTION AND

VERBALIZES UNDERSTANDING

Page: 744

DATA VALUE: GLUTEN RESTRICTED

ABBREVIATION.....GLUT RESTR
ACTIVE.....Yes
DESCRIPTION.....Gluten Restricted
DISPLAY_NAME.....Gluten Restricted
HELP_TEXT.....No Help Text
ID.....{0805A164-D387-4E67-9D42-BAF1F2543337}
TERM.....GLUTEN RESTRICTED

DATA VALUE: GNAWING

ABBREVIATION.....GNAWING
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Gnawing
HELP_TEXT.....No Help Text
ID.....{8AC8EDB0-D7B7-4368-BE95-C35CBADD43E0}
TERM.....GNAWING

DATA VALUE: GOOD

ABBREVIATION.....GOOD
ACTIVE.....Yes
DESCRIPTION.....Good
DISPLAY_NAME.....Good
HELP_TEXT.....No Help Text
ID.....{39ECA304-49F9-44B7-9735-B33CEC9D8119}
TERM.....GOOD

DATA VALUE: GRANULATING

ABBREVIATION.....GRANULATI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Granulating
HELP_TEXT.....No Help Text
ID.....{F8930276-3954-4BFE-82B0-6C44A8E9CAE0}
TERM.....GRANULATING

DATA VALUE: GRAVITY

ABBREVIATION.....GRAVITY
ACTIVE.....Yes
DESCRIPTION.....Gravity
DISPLAY_NAME.....Gravity
HELP_TEXT.....No Help Text
ID.....{C0233F70-376D-41CE-9F14-714C7B35244F}
TERM.....GRAVITY

DATA VALUE: GRAY

ABBREVIATION.....GRAY
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Gray
HELP_TEXT.....No Help Text
ID.....{D38ABEE0-52D7-4316-A713-40ACF4250599}
TERM.....GRAY

DATA VALUE: GREEN

ABBREVIATION.....GREEN
ACTIVE.....Yes
DESCRIPTION.....The color green
DISPLAY_NAME.....Green
HELP_TEXT.....The color green
ID.....{3D1F1553-56FC-68CE-79A9-E4F6357B6AAF}
TERM.....GREEN

DATA VALUE: GRIMACING

ABBREVIATION.....GRIMACING
ACTIVE.....Yes
DESCRIPTION.....Grimacing
DISPLAY_NAME.....Grimacing
HELP_TEXT.....No Help Text
ID.....{69EFDDB2D-A0C4-48D4-A7E0-A862F34545C8}
TERM.....GRIMACING

DATA VALUE: GROUP ROOM

ABBREVIATION.....gr
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Group Room
HELP_TEXT.....No Help Text
ID.....{D41C3D9F-8F10-1A0F-D571-20DE0EA36328}
TERM.....GROUP ROOM

DATA VALUE: GUARDING

ABBREVIATION.....GUARDING
ACTIVE.....Yes
DESCRIPTION.....Guarding
DISPLAY_NAME.....Guarding
HELP_TEXT.....No Help Text
ID.....{9EFD10A7-B1B7-43DA-B289-646AF5C0A25C}
TERM.....GUARDING

DATA VALUE: GURNEY CHAIR

ABBREVIATION.....GURN CHAIR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Gurney Chair
HELP_TEXT.....No Help Text
ID.....{B6EA0944-8414-42FD-88C0-50155DE82896}
TERM.....GURNEY CHAIR

DATA VALUE: HALL

ABBREVIATION.....h
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Hall
HELP_TEXT.....No Help Text
ID.....{BFCEA10F-2663-199E-88D8-E263122D7D26}
TERM.....HALL

DATA VALUE: HARD

ABBREVIATION.....HARD
ACTIVE.....Yes
DESCRIPTION.....Hard
DISPLAY_NAME.....Hard
HELP_TEXT.....No Help Text
ID.....{7934B543-7028-4CBC-9308-99F5A952A1DE}
TERM.....HARD

DATA VALUE: HAS NOT SLEPT

ABBREVIATION.....NO SLEEP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Has Not Slept
HELP_TEXT.....No Help Text
ID.....{3BF9DB54-5D1E-422D-90AF-918F81C9EA85}
TERM.....HAS NOT SLEPT

DATA VALUE: HEAD DOWN

ABBREVIATION.....hd
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Head Down
HELP_TEXT.....No Help Text
ID.....{80BF3425-D252-F917-052E-DC54D39E464E}
TERM.....HEAD DOWN

Page: 747

DATA VALUE: HEAD UP

ABBREVIATION.....hu
ACTIVE.....Yes

DESCRIPTION.....No Description
DISPLAY_NAME.....Head Up
HELP_TEXT.....No Help Text
ID.....{BB1C79C0-80F1-7CDE-CC10-7CA0C704AE35}
TERM.....HEAD UP

DATA VALUE: HEARING

ABBREVIATION.....HEARING
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Hearing
HELP_TEXT.....No Help Text
ID.....{D1A8AC66-1A5F-4F7D-AA89-529AF7A4AA11}
TERM.....HEARING

DATA VALUE: HEARTBURN

ABBREVIATION.....HEARTBURN
ACTIVE.....Yes
DESCRIPTION.....Heartburn
DISPLAY_NAME.....Heartburn
HELP_TEXT.....Heartburn
ID.....{5ECEAF5E-6B18-4177-A18D-2DEE87796CC2}
TERM.....HEARTBURN

DATA VALUE: HEAT

ABBREVIATION.....HEAT
ACTIVE.....Yes
DESCRIPTION.....Heat
DISPLAY_NAME.....Heat
HELP_TEXT.....No Help Text
ID.....{E0FCA383-3AB8-4E03-9ABC-00DE3FE01431}
TERM.....HEAT

DATA VALUE: HEAT CHEM PACK

ABBREVIATION.....HT CHEM PA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Heat Chemi Pack
HELP_TEXT.....No Help Text
ID.....{D759F891-2535-40C2-A68C-678170E56F7A}
TERM.....HEAT CHEM PACK

DATA VALUE: HEAT LAMP

ABBREVIATION.....HEAT LAMP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Heat Lamp
HELP_TEXT.....Heat Lamp

ID.....{4E9DE229-F810-4AB0-9FD7-4596759CEA11}
TERM.....HEAT LAMP

DATA VALUE: HEAVY

ABBREVIATION.....HEAVY
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Heavy
HELP_TEXT.....No Help Text
ID.....{1F10CD98-03FA-4705-BF7D-941E881C4E3D}
TERM.....HEAVY

DATA VALUE: HEEL LIFT

ABBREVIATION.....HEEL LIFT
ACTIVE.....Yes
DESCRIPTION.....Heel Lift
DISPLAY_NAME.....Heel Lift
HELP_TEXT.....No Help Text
ID.....{974F90E6-382D-421A-A4CE-94A74CBEFC32}
TERM.....HEEL LIFT

DATA VALUE: HEEL PILLOW

ABBREVIATION.....HEEL PILLO
ACTIVE.....Yes
DESCRIPTION.....Heel Pillow
DISPLAY_NAME.....Heel Pillow
HELP_TEXT.....No Help Text
ID.....{F0F1D997-C780-40BF-AA88-78262E81BFAA}
TERM.....HEEL PILLOW

DATA VALUE: HEMATOMA

ABBREVIATION.....HEMATOMA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Hematoma
HELP_TEXT.....No Help Text
ID.....{6A24B7EE-62F3-C952-C1A3-086CDBB1FE68}
TERM.....HEMATOMA

Page: 749

DATA VALUE: HEMATURIA

ABBREVIATION.....HEMATURIA
ACTIVE.....Yes
DESCRIPTION.....Blood present in urine
DISPLAY_NAME.....Hematuria
HELP_TEXT.....Blood present in urine
ID.....{ABD21F4E-BC48-D43A-38CF-75DB7704C4B8}
TERM.....HEMATURIA

DATA VALUE: HEMORRHAGE

ABBREVIATION.....HEMORRHAGE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Hemorrhage
HELP_TEXT.....No Help Text
ID.....{E00AC828-83F3-4F69-B97B-D6AD1D87A440}
TERM.....HEMORRHAGE

DATA VALUE: HEMOVAC

ABBREVIATION.....HEMOVAC
ACTIVE.....Yes
DESCRIPTION.....Hemovac
DISPLAY_NAME.....Hemovac
HELP_TEXT.....No Help Text
ID.....{82DA2DA7-E20B-E571-659C-457DF42472C2}
TERM.....HEMOVAC

DATA VALUE: HEPATIC FORMULA

ABBREVIATION.....HEPAT FORM
ACTIVE.....Yes
DESCRIPTION.....Hepatic Enteral Nutrition Formula
DISPLAY_NAME.....Hepatic Enteral Nutrition Formula
HELP_TEXT.....No Help Text
ID.....{1BC0A09A-CCD7-457F-BC08-9B6FC5144BF3}
TERM.....HEPATIC FORMULA

DATA VALUE: HEPLOCK

ABBREVIATION.....HEPLOCK
ACTIVE.....Yes
DESCRIPTION.....Heplock
DISPLAY_NAME.....Heplock
HELP_TEXT.....No Help Text
ID.....{C827FEDE-59D6-4C45-BDDA-4884FA0882F3}
TERM.....HEPLOCK

Page: 750

DATA VALUE: HIGH

ABBREVIATION.....h
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....High
HELP_TEXT.....No Help Text
ID.....{D88CB5EC-2AD4-2F80-2D43-40ECB0B7D390}
TERM.....HIGH

DATA VALUE: HIGH CALORIE

ABBREVIATION.....HIGH CALOR

ACTIVE.....Yes
DESCRIPTION.....High Calorie
DISPLAY_NAME.....High Calorie
HELP_TEXT.....No Help Text
ID.....{FC453BF0-1E7E-4733-9605-AC3425F16205}
TERM.....HIGH CALORIE

DATA VALUE: HIGH FALL RISK

ABBREVIATION.....hfr
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....High Fall Risk
HELP_TEXT.....No Help Text
ID.....{276C15F3-5BED-0FE0-5649-337A23AB1386}
TERM.....HIGH FALL RISK

DATA VALUE: HIGH FIBER

ABBREVIATION.....HIGH FIBER
ACTIVE.....Yes
DESCRIPTION.....High Fiber
DISPLAY_NAME.....High Fiber
HELP_TEXT.....No Help Text
ID.....{3CAF7EE0-6A2D-41A3-BEDC-C57ECBBA72A2}
TERM.....HIGH FIBER

DATA VALUE: HIGH FOWLER'S

ABBREVIATION.....HIGH FOWLE
ACTIVE.....Yes
DESCRIPTION.....High Fowler's
DISPLAY_NAME.....High Fowler's
HELP_TEXT.....No Help Text
ID.....{EA23C94D-EE91-4552-9F74-E051FAF3E194}
TERM.....HIGH FOWLER'S

Page: 751

DATA VALUE: HIGH PROTEIN

ABBREVIATION.....HIGH PROTE
ACTIVE.....Yes
DESCRIPTION.....High Protein
DISPLAY_NAME.....High Protein
HELP_TEXT.....No Help Text
ID.....{1595A29D-CA91-45AE-893B-856EB909E4DD}
TERM.....HIGH PROTEIN

DATA VALUE: HIRSUTISM

ABBREVIATION.....HIRSUTISM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Hirsutism

HELP_TEXT.....No Help Text
ID.....{A909A4FC-080A-4E9A-867A-0B30D4633602}
TERM.....HIRSUTISM

DATA VALUE: HOME

ABBREVIATION.....HOME
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Home
HELP_TEXT.....No Help Text
ID.....{C8E5E320-361B-421C-A6A9-13DC16125440}
TERM.....HOME

DATA VALUE: HOME HEALTH

ABBREVIATION.....HOME HEALT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Home Health
HELP_TEXT.....No Help Text
ID.....{F94FB70A-AFFD-4A5E-A829-0AF933B044BB}
TERM.....HOME HEALTH

DATA VALUE: HOSTILE

ABBREVIATION.....h
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Hostile
HELP_TEXT.....No Help Text
ID.....{29EF2C04-7DB0-4491-2B5E-D92A6A558721}
TERM.....HOSTILE

Page: 752

DATA VALUE: HOT

ABBREVIATION.....HOT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Hot
HELP_TEXT.....No Help Text
ID.....{401EA6CD-12DE-AB0E-2B52-8B2E62F7C34E}
TERM.....HOT

DATA VALUE: HUMIDIFIER

ABBREVIATION.....HUMIDIFIER
ACTIVE.....Yes
DESCRIPTION.....Humidifier
DISPLAY_NAME.....Humidifier
HELP_TEXT.....Humidifier
ID.....{A4CE67D1-F885-423B-B914-392D8CAC697E}
TERM.....HUMIDIFIER

DATA VALUE: HYDROCOLLOIDAL DRESSING

ABBREVIATION.....HCOLL DRES
ACTIVE.....Yes
DESCRIPTION.....Hydrocolloidal Dressing
DISPLAY_NAME.....Hydrocolloidal Dressing
HELP_TEXT.....No Help Text
ID.....{D85741EF-0DB3-48CF-BC9A-88C834A72F36}
TERM.....HYDROCOLLOIDAL DRESSING

DATA VALUE: HYDROFIBER

ABBREVIATION.....H-FIBER
ACTIVE.....Yes
DESCRIPTION.....Hydrofiber
DISPLAY_NAME.....Hydrofiber
HELP_TEXT.....No Help Text
ID.....{42740A7F-A831-4FBD-9E59-88F315176674}
TERM.....HYDROFIBER

DATA VALUE: HYDROFIBER WITH SILVER

ABBREVIATION.....HYFIB W SI
ACTIVE.....Yes
DESCRIPTION.....Hydrofiber With Silver
DISPLAY_NAME.....Hydrofiber With Silver
HELP_TEXT.....No Help Text
ID.....{04BD9512-DE5D-421A-9818-8A33D1F0693D}
TERM.....HYDROFIBER WITH SILVER

Page: 753

DATA VALUE: HYDROGEL DRESSING

ABBREVIATION.....H-GEL DRES
ACTIVE.....Yes
DESCRIPTION.....Hydrogel Dressing
DISPLAY_NAME.....Hydrogel Dressing
HELP_TEXT.....No Help Text
ID.....{F6192D10-D6EA-42C2-AB15-270659D9925C}
TERM.....HYDROGEL DRESSING

DATA VALUE: HYPERACTIVE

ABBREVIATION.....HYPERACTIV
ACTIVE.....Yes
DESCRIPTION.....Hyperactive
DISPLAY_NAME.....Hyperactive
HELP_TEXT.....No Help Text
ID.....{1A548FFE-8645-4C0D-899D-5D91430CC1B1}
TERM.....HYPERACTIVE

DATA VALUE: HYPERDYNAMIC

ABBREVIATION.....HYPERDYNAM
ACTIVE.....Yes
DESCRIPTION.....Hyperdynamic
DISPLAY_NAME.....Hyperdynamic
HELP_TEXT.....No Help Text
ID.....{818A8A99-04D9-4989-A9BA-ECFACED74605}
TERM.....HYPERDYNAMIC

DATA VALUE: HYPERGRANULIZATION

ABBREVIATION.....HYPERGRANU
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Hypergranulization
HELP_TEXT.....No Help Text
ID.....{C14B9FB0-6F33-4ACB-88B3-4B799594F845}
TERM.....HYPERGRANULIZATION

DATA VALUE: HYPERTONIC HYDROGEL

ABBREVIATION.....WOU DR HYP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Hypertonic Hydrogel
HELP_TEXT.....No Help-Text
ID.....{533D6FC3-5183-4417-AB89-9FAD9BAD0F44}
TERM.....HYPERTONIC HYDROGEL

Page: 754

DATA VALUE: HYPOACTIVE

ABBREVIATION.....HYPOACTIVE
ACTIVE.....Yes
DESCRIPTION.....Hypoactive
DISPLAY_NAME.....Hypoactive
HELP_TEXT.....No Help Text
ID.....{2E0C516D-3858-4A1F-A2F3-BF0AB9E3A7FC}
TERM.....HYPOACTIVE

DATA VALUE: HYPOALBUMINEMIA

ABBREVIATION.....HYPOALBUM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Hypoalbuminemia
HELP_TEXT.....No Help Text
ID.....{A3397F38-2491-48E7-8366-A6FA2D91E866}
TERM.....HYPOALBUMINEMIA

DATA VALUE: I

ABBREVIATION.....I
ACTIVE.....Yes
DESCRIPTION.....No Description

DISPLAY_NAME.....I
HELP_TEXT.....No Help Text
ID.....{BEEDFF9A-9401-4D2A-8345-126D8CFFC8D8}
TERM.....I

DATA VALUE: IABP

ABBREVIATION.....IABP
ACTIVE.....Yes
DESCRIPTION.....Intra-Aortic Balloon Pump
DISPLAY_NAME.....IABP
HELP_TEXT.....No Help Text
ID.....{763FE975-F5E5-80A6-AB16-E903A829E665}
TERM.....IABP

DATA VALUE: ICE PACK

ABBREVIATION.....ICE PACK
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Ice Pack
HELP_TEXT.....No Help Text
ID.....{E46ED309-D60D-45DF-BFA3-F74DD13FC4AE}
TERM.....ICE PACK

Page: 755

DATA VALUE: ICTERIC

ABBREVIATION.....ICTERIC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Icteric
HELP_TEXT.....Relating to or affected with jaundice
ID.....{274ED5BB-5A34-20EE-A9B0-81AB7803780E}
TERM.....ICTERIC

DATA VALUE: ICU

ABBREVIATION.....ICU
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....ICU
HELP_TEXT.....No Help Text
ID.....{E7FA69FE-1088-4853-AFBD-2D815E3C2EB0}
TERM.....ICU

DATA VALUE: IDENTIFICATION

ABBREVIATION.....i
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Identification
HELP_TEXT.....No Help Text
ID.....{39274F69-556C-7F05-6D17-CE8A03D1867D}

TERM.....IDENTIFICATION

DATA VALUE: IDIOVENTRICULAR

ABBREVIATION.....i
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Idioventricular
HELP_TEXT.....No Help Text
ID.....{8BD08966-F536-8E74-0BDA-A95D0D7F8B6C}
TERM.....IDIOVENTRICULAR

DATA VALUE: II

ABBREVIATION.....II
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....II
HELP_TEXT.....No Help Text
ID.....{F165C13B-397E-4182-8D72-BD1D74CFBA00}
TERM.....II

Page: 756

DATA VALUE: III

ABBREVIATION.....III
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....III
HELP_TEXT.....No Help Text
ID.....{D8FE8FA9-9BB5-416E-AE2C-60E627F76E5A}
TERM.....III

DATA VALUE: ILEAL CONDUIT

ABBREVIATION.....ILEAL COND
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Ileal Conduit
HELP_TEXT.....No Help Text
ID.....{032CB6FE-77D2-702E-A0AE-46176F5CAA57}
TERM.....ILEAL CONDUIT

DATA VALUE: ILEOSTOMY

ABBREVIATION.....ILEOSTOMY
ACTIVE.....Yes
DESCRIPTION.....Ileostomy
DISPLAY_NAME.....Ileostomy
HELP_TEXT.....No Help Text
ID.....{D5DF49C8-C6BA-44C1-9840-7765907730A3}
TERM.....ILEOSTOMY

DATA VALUE: IMAGING

ABBREVIATION.....IMAGING
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Imaging
HELP_TEXT.....No Help Text
ID.....{791A322E-F45E-40EB-B95D-E5B62E1D0A75}
TERM.....IMAGING

DATA VALUE: IMPAIRED LISTENING SKILLS

ABBREVIATION.....IMP LIS SK
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Impaired Listening Skills
HELP_TEXT.....No Help Text
ID.....{9CE0E7E7-79E5-4910-918C-2022C9830171}
TERM.....IMPAIRED LISTENING SKILLS

Page: 757

DATA VALUE: INAPPROPRIATE

ABBREVIATION.....Inappropri
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Inappropriate
HELP_TEXT.....No Help Text
ID.....{7E3B54C2-8A49-28EC-055E-FEA87599421C}
TERM.....INAPPROPRIATE

DATA VALUE: INCISION DISCHARGE CHARACTERISTICS BLOODY

ABBREVIATION.....INC DI BLO
ACTIVE.....Yes
DESCRIPTION.....Incision Discharge Characteristics
Bloody
DISPLAY_NAME.....Bloody
HELP_TEXT.....Incision Discharge Characteristics
Bloody
ID.....{D185E55B-EC21-49B6-A040-5422382C3B90}
TERM.....INCISION DISCHARGE CHARACTERISTICS
BLOODY

DATA VALUE: INCISION DISCHARGE CHARACTERISTICS CLOUDY

ABBREVIATION.....INC DI CLO
ACTIVE.....Yes
DESCRIPTION.....Incision Discharge Characteristics
Cloudy
DISPLAY_NAME.....Cloudy
HELP_TEXT.....Incision Discharge Characteristics
Cloudy
ID.....{3EBBB25E-1B90-415B-BCB5-3C3DAC33C118}
TERM.....INCISION DISCHARGE CHARACTERISTICS

CLOUDY

DATA VALUE: INCISION DISCHARGE CHARACTERISTICS NONE

ABBREVIATION.....INC DI NON
ACTIVE.....Yes
DESCRIPTION.....Incision Discharge Characteristics None
DISPLAY_NAME.....None
HELP_TEXT.....Incision Discharge Characteristics None
ID.....{4F118F41-8C09-400E-A303-0955DFDF54E2}
TERM.....INCISION DISCHARGE CHARACTERISTICS NONE

DATA VALUE: INCISION DISCHARGE CHARACTERISTICS PURULENT

ABBREVIATION.....INC DI PUR
ACTIVE.....Yes
DESCRIPTION.....Incision Discharge Characteristics
Purulent
DISPLAY_NAME.....Purulent
HELP_TEXT.....Incision Discharge Characteristics
Purulent
ID.....{B2D71BDF-41AB-4F10-ABE3-A9D9F338F6F1}
TERM.....INCISION DISCHARGE CHARACTERISTICS
PURULENT

Page: 758

DATA VALUE: INCISION DISCHARGE CHARACTERISTICS SEROSANGUINOUS

ABBREVIATION.....INC DI SSA
ACTIVE.....Yes
DESCRIPTION.....Incision Discharge Characteristics
Serosanguinous
DISPLAY_NAME.....Serosanguinous
HELP_TEXT.....Incision Discharge Characteristics
Serosanguinous
ID.....{102F19A9-CC28-4156-88DE-1B4C8A682AC2}
TERM.....INCISION DISCHARGE CHARACTERISTICS
SEROSANGUINOUS

DATA VALUE: INCISION DISCHARGE CHARACTERISTICS SEROUS

ABBREVIATION.....INC DI SER
ACTIVE.....Yes
DESCRIPTION.....Incision Discharge Characteristics
Serous
DISPLAY_NAME.....Serous
HELP_TEXT.....Incision Discharge Characteristics
Serous
ID.....{CC4F3BA1-AD0D-42BA-9D76-1FE66F864FC8}
TERM.....INCISION DISCHARGE CHARACTERISTICS
SEROUS

DATA VALUE: INCISION DISCHARGE CHARACTERISTICS VISCOUS

ABBREVIATION.....INC DI VIS
ACTIVE.....Yes
DESCRIPTION.....Incision Discharge Characteristics
Viscous
DISPLAY_NAME.....Viscous
HELP_TEXT.....Incision Discharge Characteristics
Viscous
ID.....{A38C2D8C-25B1-41E5-A69F-7A5C439C6EB9}
TERM.....INCISION DISCHARGE CHARACTERISTICS
VISCOUS

DATA VALUE: INCISION DRAINAGE AMOUNT HEAVY

ABBREVIATION.....I DR A HEA
ACTIVE.....Yes
DESCRIPTION.....Incision Drainage Amount Heavy
DISPLAY_NAME.....Heavy
HELP_TEXT.....Incision Drainage Amount Heavy
ID.....{D77FCEA0-D247-4DD5-873E-2F95C7FF48A6}
TERM.....INCISION DRAINAGE AMOUNT HEAVY

DATA VALUE: INCISION DRAINAGE AMOUNT MODERATE

ABBREVIATION.....I DR A MOD
ACTIVE.....Yes
DESCRIPTION.....Incision Drainage Amount Moderate
DISPLAY_NAME.....Moderate
HELP_TEXT.....Incision Drainage Amount Moderate
ID.....{9B0487A9-DCB4-42D3-92E6-805C412B5D34}
TERM.....INCISION DRAINAGE AMOUNT MODERATE

Page: 759

DATA VALUE: INCISION DRAINAGE AMOUNT NONE

ABBREVIATION.....I DR A NON
ACTIVE.....Yes
DESCRIPTION.....Incision Drainage Amount None
DISPLAY_NAME.....None
HELP_TEXT.....Incision Drainage Amount None
ID.....{5D0D9027-4B20-43A4-B437-EA7B77CB2105}
TERM.....INCISION DRAINAGE AMOUNT NONE

DATA VALUE: INCISION DRAINAGE AMOUNT SCANT

ABBREVIATION.....I DR A SCA
ACTIVE.....Yes
DESCRIPTION.....Incision Drainage Amount Scant
DISPLAY_NAME.....Scant
HELP_TEXT.....Incision Drainage Amount Scant
ID.....{A476E395-46AA-44E7-8A00-CE22C7066BB1}
TERM.....INCISION DRAINAGE AMOUNT SCANT

DATA VALUE: INCISION DRAINAGE AMOUNT SMALL

ABBREVIATION.....I DR A SMA
ACTIVE.....Yes
DESCRIPTION.....Incision Drainage Amount Small
DISPLAY_NAME.....Small
HELP_TEXT.....Incision Drainage Amount Small
ID.....{4479479D-8ADA-4FBF-A5A2-426D49FC94D3}
TERM.....INCISION DRAINAGE AMOUNT SMALL

DATA VALUE: INCISION DRESSING TYPE ANTIMICROBIAL

ABBREVIATION.....INC DR ANT
ACTIVE.....Yes
DESCRIPTION.....Incision Dressing Type Antimicrobial
DISPLAY_NAME.....Antimicrobial
HELP_TEXT.....Incision Dressing Type Antimicrobial
ID.....{35F5F5C1-B76A-44D2-8576-9D7D3ECA3197}
TERM.....INCISION DRESSING TYPE ANTIMICROBIAL

DATA VALUE: INCISION DRESSING TYPE ANTIMICROBIAL DISK

ABBREVIATION.....INC DR ADI
ACTIVE.....Yes
DESCRIPTION.....Incision Dressing Type Antimicrobial
Disk
DISPLAY_NAME.....Antimicrobial Disk
HELP_TEXT.....Incision Dressing Type Antimicrobial
Disk
ID.....{55A33062-7D6B-4CDF-9131-D4ADC38F6CE3}
TERM.....INCISION DRESSING TYPE ANTIMICROBIAL
DISK

Page: 760

DATA VALUE: INCISION DRESSING TYPE BALSAM OF PERU

ABBREVIATION.....INC DR BAL
ACTIVE.....Yes
DESCRIPTION.....Incision Dressing Type Balsam of Peru
DISPLAY_NAME.....Balsam Of Peru
HELP_TEXT.....Incision Dressing Type Balsam of Peru
ID.....{A1FB45C6-6315-4282-8D89-DCB0C018E357}
TERM.....INCISION DRESSING TYPE BALSAM OF PERU

DATA VALUE: INCISION DRESSING TYPE BISMUTH GAUZE

ABBREVIATION.....INC DR BIS
ACTIVE.....Yes
DESCRIPTION.....Incision Dressing Type Bismuth Gauze
DISPLAY_NAME.....Bismuth Gauze
HELP_TEXT.....Incision Dressing Type Bismuth Gauze
ID.....{CA88C5DC-3890-49A2-800A-0528CE194BA8}
TERM.....INCISION DRESSING TYPE BISMUTH GAUZE

DATA VALUE: INCISION DRESSING TYPE CALCIUM ALGINATE

ABBREVIATION.....INC DR CAL
ACTIVE.....Yes
DESCRIPTION.....Incision Dressing Type Calcium Alginate
DISPLAY_NAME.....Calcium Alginate
HELP_TEXT.....Incision Dressing Type Calcium Alginate
ID.....{B20558BA-FB9E-499E-82DF-49616ACE9375}
TERM.....INCISION DRESSING TYPE CALCIUM ALGINATE

DATA VALUE: INCISION DRESSING TYPE FOAM DRESSING

ABBREVIATION.....INC DR FOA
ACTIVE.....Yes
DESCRIPTION.....Incision Dressing Type Foam dressing
DISPLAY_NAME.....Foam Dressing
HELP_TEXT.....Incision Dressing Type Foam dressing
ID.....{5C2C2D24-6022-4BBF-9065-03FB5971CB3F}
TERM.....INCISION DRESSING TYPE FOAM DRESSING

DATA VALUE: INCISION DRESSING TYPE GAUZE

ABBREVIATION.....INC DR GAU
ACTIVE.....Yes
DESCRIPTION.....Incision Dressing Type Gauze
DISPLAY_NAME.....Gauze
HELP_TEXT.....Incision Dressing Type Gauze
ID.....{E4AC9307-0512-4F9D-A254-FE19E654461F}
TERM.....INCISION DRESSING TYPE GAUZE

Page: 761

DATA VALUE: INCISION DRESSING TYPE GAUZE PACKING

ABBREVIATION.....INC DR GPA
ACTIVE.....Yes
DESCRIPTION.....Incision Dressing Type Gauze Packing
DISPLAY_NAME.....Gauze Packing
HELP_TEXT.....Incision Dressing Type Gauze Packing
ID.....{3561DAEB-8A30-4E16-B920-F7DFE022EBB1}
TERM.....INCISION DRESSING TYPE GAUZE PACKING

DATA VALUE: INCISION DRESSING TYPE GAUZE WRAP

ABBREVIATION.....INC DR GWR
ACTIVE.....Yes
DESCRIPTION.....Incision Dressing Type Gauze Wrap
DISPLAY_NAME.....Gauze Wrap
HELP_TEXT.....Incision Dressing Type Gauze Wrap
ID.....{E1DBEC29-1430-4DA0-99DC-49B9389D2440}
TERM.....INCISION DRESSING TYPE GAUZE WRAP

DATA VALUE: INCISION DRESSING TYPE HYDROCOLLOIDAL

ABBREVIATION.....INC DR HYD
ACTIVE.....Yes
DESCRIPTION.....Incision Dressing Type Hydrocolloidal

DISPLAY_NAME.....Hydrocolloidal
HELP_TEXT.....Incision Dressing Type Hydrocolloidal
ID.....{424CCDCA-5196-4BBD-AE8B-0F68DDCC3FCF}
TERM.....INCISION DRESSING TYPE HYDROCOLLOIDAL

DATA VALUE: INCISION DRESSING TYPE HYDROGEL

ABBREVIATION.....INC DR HYD
ACTIVE.....Yes
DESCRIPTION.....Incision Dressing Type Hydrogel
DISPLAY_NAME.....Hydrogel
HELP_TEXT.....Incision Dressing Type Hydrogel
ID.....{8CBAC6DF-F51E-456B-8F9A-A68E4452812C}
TERM.....INCISION DRESSING TYPE HYDROGEL

DATA VALUE: INCISION DRESSING TYPE HYPERTONIC HYDROGEL

ABBREVIATION.....INC DR HYP
ACTIVE.....Yes
DESCRIPTION.....Incision Dressing Type Hypertonic
Hydrogel
DISPLAY_NAME.....Hypertonic Hydrogel
HELP_TEXT.....Incision Dressing Type Hypertonic
Hydrogel
ID.....{25820A63-AD76-4198-BEB1-EF95DDEC5B35}
TERM.....INCISION DRESSING TYPE HYPERTONIC
HYDROGEL

Page: 762

DATA VALUE: INCISION DRESSING TYPE NONE

ABBREVIATION.....INC DR NON
ACTIVE.....Yes
DESCRIPTION.....Incision Dressing Type None
DISPLAY_NAME.....None
HELP_TEXT.....Incision Dressing Type None
ID.....{81F07C8D-A463-48A0-A64A-A143CECB5DD6}
TERM.....INCISION DRESSING TYPE NONE

DATA VALUE: INCISION DRESSING TYPE PRESSURE

ABBREVIATION.....INC DR PRE
ACTIVE.....Yes
DESCRIPTION.....Incision Dressing Type Pressure
DISPLAY_NAME.....Pressure
HELP_TEXT.....Incision Dressing Type Pressure
ID.....{B5823373-E5D8-40A5-954C-6A193D7E7511}
TERM.....INCISION DRESSING TYPE PRESSURE

DATA VALUE: INCISION DRESSING TYPE TRANSPARENT

ABBREVIATION.....INC DR TRA
ACTIVE.....Yes
DESCRIPTION.....Incision Dressing Type Gauze

DISPLAY_NAME.....Transparent
HELP_TEXT.....Incision Dressing Type Gauze
ID.....{E198EDAD-AA67-49F9-A065-6BB11FE39ED6}
TERM.....INCISION DRESSING TYPE TRANSPARENT

DATA VALUE: INCISION EDGE INDISTINCT/ATTACHED

ABBREVIATION.....INC ED ATT
ACTIVE.....Yes
DESCRIPTION.....Incision Edge Indistinct, Attached
DISPLAY_NAME.....Indistinct, Attached
HELP_TEXT.....Incision Edge Indistinct, Attached
ID.....{AF7E7150-525E-49A4-AC23-042B788E8805}
TERM.....INCISION EDGE INDISTINCT/ATTACHED

DATA VALUE: INCISION EDGE NON-ATTACHED

ABBREVIATION.....INC ED NON
ACTIVE.....Yes
DESCRIPTION.....Incision Edge Non Attached
DISPLAY_NAME.....Non-Attached
HELP_TEXT.....Incision Edge Non Attached
ID.....{F457F951-F1DC-47AC-9336-AC093E149862}
TERM.....INCISION EDGE NON-ATTACHED

Page: 763

DATA VALUE: INCISION EDGE ROLLED

ABBREVIATION.....INC ED ROL
ACTIVE.....Yes
DESCRIPTION.....Incision Edge Rolled
DISPLAY_NAME.....Rolled
HELP_TEXT.....Incision Edge Rolled
ID.....{FFA33CAE-5D8B-4501-8120-A536BE53296F}
TERM.....INCISION EDGE ROLLED

DATA VALUE: INCISION EDGE TUNNELING

ABBREVIATION.....INC ED TUN
ACTIVE.....Yes
DESCRIPTION.....Incision Edge Tunneling
DISPLAY_NAME.....Tunneling
HELP_TEXT.....Incision Edge Tunneling
ID.....{492CD783-3C32-40F0-AE35-5A89E90409A7}
TERM.....INCISION EDGE TUNNELING

DATA VALUE: INCISION EDGE UNDERMINING

ABBREVIATION.....INC ED UND
ACTIVE.....Yes
DESCRIPTION.....Incision Edge Undermining
DISPLAY_NAME.....Undermining
HELP_TEXT.....Incision Edge Undermining
ID.....{0F3E968A-B867-4123-8C05-ED6B12D364CD}

TERM.....INCISION EDGE UNDERMINING

DATA VALUE: INCISION ODOR FOUL

ABBREVIATION.....INC OD FOU
ACTIVE.....Yes
DESCRIPTION.....Incision Odor Foul
DISPLAY_NAME.....Foul
HELP_TEXT.....Incision Odor Foul
ID.....{ECDC0988-BE82-41E3-8104-A66A5F8A1093}
TERM.....INCISION ODOR FOUL

DATA VALUE: INCISION ODOR NONE

ABBREVIATION.....INC OD NON
ACTIVE.....Yes
DESCRIPTION.....Incision Odor None
DISPLAY_NAME.....None
HELP_TEXT.....Incision Odor None
ID.....{CED90890-6FBC-4A4A-B0B8-6F75F0DD567A}
TERM.....INCISION ODOR NONE

Page: 764

DATA VALUE: INCISION WOUND CLOSURE APPROXIMATED

ABBREVIATION.....INC CL APP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Approximated
HELP_TEXT.....No Help-Text
ID.....{2DD5713E-D81C-484B-A0DF-6A5730B396DF}
TERM.....INCISION WOUND CLOSURE APPROXIMATED

DATA VALUE: INCISION WOUND CLOSURE DEHISCENCE

ABBREVIATION.....INC CL DEH
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Dehiscence
HELP_TEXT.....No Help-Text
ID.....{DCA1D137-0D9A-42D2-9270-7B3A82DCB6EB}
TERM.....INCISION WOUND CLOSURE DEHISCENCE

DATA VALUE: INCISION WOUND CLOSURE GLUE

ABBREVIATION.....IN CL GLU
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Glue
HELP_TEXT.....No Help-Text
ID.....{9EEB436E-2991-4E9A-9F8A-BFEC196784CB}
TERM.....INCISION WOUND CLOSURE GLUE

DATA VALUE: INCISION WOUND CLOSURE NOT CLOSED

ABBREVIATION.....INC CL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Not Closed
HELP_TEXT.....No Help-Text
ID.....{D63A05BF-1205-45F6-A78D-15D198BCB1C6}
TERM.....INCISION WOUND CLOSURE NOT CLOSED

DATA VALUE: INCISIONAL SITE

ABBREVIATION.....INCIS SITE
ACTIVE.....Yes
DESCRIPTION.....INCISIONAL SITE
DISPLAY_NAME.....Incisional Site
HELP_TEXT.....INCISIONAL SITE
ID.....{19E69DE2-097B-0FA5-449D-B2E9161A3E05}
TERM.....INCISIONAL SITE

Page: 765

DATA VALUE: INCOMPREHENSIBLE

ABBREVIATION.....INCOMPREHE
ACTIVE.....Yes
DESCRIPTION.....Incomprehensible
DISPLAY_NAME.....Incomprehensible
HELP_TEXT.....No Help Text
ID.....{916CD781-563D-4F54-8752-800D37BFAA3B}
TERM.....INCOMPREHENSIBLE

DATA VALUE: INCONTENENCE DERMATITIS

ABBREVIATION.....INCON DERM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Incontenance Dermatitis
HELP_TEXT.....No Help Text
ID.....{9F6C5BA2-79C7-4B09-B19A-6DDC444936E5}
TERM.....INCONTENENCE DERMATITIS

DATA VALUE: INCONTINENT

ABBREVIATION.....INCONTINEN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Incontinent
HELP_TEXT.....No Help Text
ID.....{BA2E9C32-8D03-8FC9-5E02-EAE2026FDC63}
TERM.....INCONTINENT

DATA VALUE: INDEPENDENT

ABBREVIATION.....INDEPENDEN
ACTIVE.....Yes

DESCRIPTION.....No Description
DISPLAY_NAME.....Independent
HELP_TEXT.....No Help-Text
ID.....{11BC0CF7-1C29-49C7-B7C8-2BCF8FF82250}
TERM.....INDEPENDENT

DATA VALUE: INDICATES UNDERSTANDING

ABBREVIATION.....UNDERSTAND
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Indicates Understanding
HELP_TEXT.....No Help Text
ID.....{C4AEFE5F-C388-4E00-9663-6B308934284B}
TERM.....INDICATES UNDERSTANDING

Page: 766

DATA VALUE: INDIFFERENT

ABBREVIATION.....INDIFF
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Indifferent
HELP_TEXT.....No Help Text
ID.....{7C3B96CF-C3D0-4117-A952-94BD9FBE05A5}
TERM.....INDIFFERENT

DATA VALUE: INDURATION

ABBREVIATION.....INDURAITON
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Induration (Hardness)
HELP_TEXT.....No Help Text
ID.....{A1669576-45D9-435C-87B9-488E9CF6BBA4}
TERM.....INDURATION

DATA VALUE: INDWELLING CATHETER

ABBREVIATION.....INDW CATH
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Indwelling Catheter
HELP_TEXT.....No Help Text
ID.....{EC16724C-D1EA-5650-58D4-38A89D4456F9}
TERM.....INDWELLING CATHETER

DATA VALUE: INEFFECTIVE COPING

ABBREVIATION.....ic
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Ineffective Coping
HELP_TEXT.....No Help Text

ID.....{704B1A55-B639-60F4-E3A5-16F942ABC259}
TERM.....INEFFECTIVE COPING

DATA VALUE: INFECTION

ABBREVIATION.....INFECTION
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Infection
HELP_TEXT.....No Help Text
ID.....{DA00DB58-640D-4294-9468-EEF070DBAE8E}
TERM.....INFECTION

Page: 767

DATA VALUE: INFECTION PREVENTION AND SCREENING

ABBREVIATION.....INF PRE SC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Infection Prevention and Screening
HELP_TEXT.....No Help Text
ID.....{BA58329F-59EC-49F0-BACE-75AA5EC57F52}
TERM.....INFECTION PREVENTION AND SCREENING

DATA VALUE: INFLATABLE HEEL LIFT

ABBREVIATION.....INF HEEL L
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Inflatable Heel Lift
HELP_TEXT.....No Help Text
ID.....{793B34A4-2F21-42DA-881D-1E6E0B4EAD6B}
TERM.....INFLATABLE HEEL LIFT

DATA VALUE: INITIATED

ABBREVIATION.....INITIATED
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Initiated
HELP_TEXT.....No Help Text
ID.....{562E5FFB-4E53-4614-ACAC-E7D9B6E1B781}
TERM.....INITIATED

DATA VALUE: INSERTED THROUGH EXISTING INTRODUCER

ABBREVIATION.....INS EXS IN
ACTIVE.....Yes
DESCRIPTION.....Inserted Through Existing Introducer
DISPLAY_NAME.....Inserted Through Existing Introducer
HELP_TEXT.....No Help Text
ID.....{20BC183A-197F-4464-A27C-47F63C5B842A}
TERM.....INSERTED THROUGH EXISTING INTRODUCER

DATA VALUE: INSERTED THROUGH NEW INTRODUCER

ABBREVIATION.....INS NEW IN
ACTIVE.....Yes
DESCRIPTION.....Inserted Through New Introducer
DISPLAY_NAME.....Inserted Through New Introducer
HELP_TEXT.....No Help Text
ID.....{7FCA403D-41E5-4BAA-B155-B58F35EA49B3}
TERM.....INSERTED THROUGH NEW INTRODUCER

Page: 768

DATA VALUE: INSPIRATORY WHEEZES

ABBREVIATION.....INSP WHEEZ
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Inspiratory Wheezes
HELP_TEXT.....No Help Text
ID.....{3FA7F54C-DE17-4DFF-B5F6-C2CD05EC0B24}
TERM.....INSPIRATORY WHEEZES

DATA VALUE: INTACT

ABBREVIATION.....i
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Intact
HELP_TEXT.....No apparent breakdown
ID.....{935C9CAF-6C88-6610-C6E9-4C1B822923AF}
TERM.....INTACT

DATA VALUE: INTENTIONAL

ABBREVIATION.....i
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Intentional
HELP_TEXT.....No Help Text
ID.....{D5BBBFFC-F739-1E69-FED4-288D8E803D04}
TERM.....INTENTIONAL

DATA VALUE: INTERMITTENT

ABBREVIATION.....INTERMITTE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Intermittent
HELP_TEXT.....No Help Text
ID.....{DD09ED7C-B977-412F-A678-5130F77D239E}
TERM.....INTERMITTENT

DATA VALUE: INTERMITTENT SLEEP

ABBREVIATION.....INT SLEEP

ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....States Has Slept Intermittently
HELP_TEXT.....No Help Text
ID.....{5005C433-5892-4F4D-A64D-4C94AD1853E2}
TERM.....INTERMITTENT SLEEP

Page: 769

DATA VALUE: INTERRUPTED

ABBREVIATION.....INTERRUPTED
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Interrupted
HELP_TEXT.....No Help Text
ID.....{AB53B9AB-DA3F-42D3-BEB2-E02A452ED777}
TERM.....INTERRUPTED

DATA VALUE: INTRODUCER

ABBREVIATION.....INTRODUCER
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Introducer
HELP_TEXT.....No Help Text
ID.....{B46F7CA4-C94E-68FD-2892-52866763B5A0}
TERM.....INTRODUCER

DATA VALUE: IRREGULAR

ABBREVIATION.....IRREGULAR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Irregular
HELP_TEXT.....No Help Text
ID.....{6CCE27A2-F1EB-3EC8-B140-CDA241D159F2}
TERM.....IRREGULAR

DATA VALUE: IRRITABLE

ABBREVIATION.....Irritable
ACTIVE.....Yes
DESCRIPTION.....Irritable
DISPLAY_NAME.....Irritable
HELP_TEXT.....No Help Text
ID.....{10651C11-2736-4056-BA84-E544358BA877}
TERM.....IRRITABLE

DATA VALUE: IV DISCONTINUED

ABBREVIATION.....IV DISC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....IV Discontinued

HELP_TEXT.....No Help Text
ID.....{FDF951E7-5B72-4166-969B-073D5EC74800}
TERM.....IV DISCONTINUED

Page: 770

DATA VALUE: IV PUMP

ABBREVIATION.....IV PUMP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....IV Pump
HELP_TEXT.....No Help Text
ID.....{EFB92C16-CFA3-18BB-9D54-590089DA3962}
TERM.....IV PUMP

DATA VALUE: IVPB

ABBREVIATION.....IVPB
ACTIVE.....Yes
DESCRIPTION.....IV Piggy Back
DISPLAY_NAME.....IV Piggy Back
HELP_TEXT.....No Help Text
ID.....{50EE76CF-5A9C-4EEB-9201-37CD90509355}
TERM.....IVPB

DATA VALUE: JAUNDICE

ABBREVIATION.....j
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Jaundice
HELP_TEXT.....No Help Text
ID.....{B17F4603-7C68-660D-96A7-119140A998D3}
TERM.....JAUNDICE

DATA VALUE: JEJUNOSTOMY

ABBREVIATION.....JEJUNOSTOM
ACTIVE.....Yes
DESCRIPTION.....Jejunostomy EF Tube
DISPLAY_NAME.....Jejunostomy
HELP_TEXT.....No Help Text
ID.....{1ED3EC41-3ACF-42C1-B64A-280300EDBCFC}
TERM.....JEJUNOSTOMY

DATA VALUE: JOINT CAPSULE

ABBREVIATION.....JOINT CAPS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Joint Capsule
HELP_TEXT.....No Help Text
ID.....{F10E4DA8-DF84-4625-A317-724E7C9E028C}
TERM.....JOINT CAPSULE

DATA VALUE: JP DRAIN

ABBREVIATION.....JP Drain
ACTIVE.....Yes
DESCRIPTION.....Jackson Pratt Drain
DISPLAY_NAME.....JP Drain
HELP_TEXT.....No Help Text
ID.....{20F3B3D5-D3A4-96FD-5446-253FD031DF2B}
TERM.....JP DRAIN

DATA VALUE: JUNCTIONAL

ABBREVIATION.....j
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Junctional
HELP_TEXT.....No Help Text
ID.....{4A339D84-7051-1258-877D-85AFC1CA20FE}
TERM.....JUNCTIONAL

DATA VALUE: K-PAD

ABBREVIATION.....K-PAD
ACTIVE.....Yes
DESCRIPTION.....K-pad
DISPLAY_NAME.....K-pad
HELP_TEXT.....K-pad
ID.....{10A0C893-F279-48ED-ACAA-1DCA4AD9DB24}
TERM.....K-PAD

DATA VALUE: KNIFE-LIKE

ABBREVIATION.....KNIFE-LIKE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Knife-Like
HELP_TEXT.....No Help Text
ID.....{FBA0EE3F-134B-431B-879A-9B17D65CD464}
TERM.....KNIFE-LIKE

DATA VALUE: KUSSMAUL

ABBREVIATION.....KUSSMAUL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Kussmaul
HELP_TEXT.....No Help Text
ID.....{0C5ED066-F445-ABFE-F7B1-D158EC7C2976}
TERM.....KUSSMAUL

DATA VALUE: L BUNDLE BRANCH BLOCK

ABBREVIATION.....lbbb
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Bundle Branch Block
HELP_TEXT.....No Help Text
ID.....{C3F2E19E-495B-7BF6-99EF-B8236D6B9C3D}
TERM.....L BUNDLE BRANCH BLOCK

DATA VALUE: L FOOT DOWN

ABBREVIATION.....L FOOT DOW
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Foot Down
HELP_TEXT.....No Help Text
ID.....{5B18895B-DAD5-4874-87DB-DACFFF88B170}
TERM.....L FOOT DOWN

DATA VALUE: L FOOT UP

ABBREVIATION.....L FOOT UP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Foot Up
HELP_TEXT.....No Help Text
ID.....{A15CD8C5-5CD1-4B90-B79C-5D3D331C7466}
TERM.....L FOOT UP

DATA VALUE: L HEAD DOWN

ABBREVIATION.....L HEAD DOW
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Head Down
HELP_TEXT.....No Help Text
ID.....{3C919CB8-1C0A-498D-95A6-C46D25744CEC}
TERM.....L HEAD DOWN

DATA VALUE: L HEAD UP

ABBREVIATION.....L HEAD UP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Head Up
HELP_TEXT.....No Help Text
ID.....{BB39792D-B572-4A10-BA81-8685809E7437}
TERM.....L HEAD UP

DATA VALUE: L NECK

ABBREVIATION.....LEFT NECK
ACTIVE.....Yes
DESCRIPTION.....Left Neck
DISPLAY_NAME.....L Neck
HELP_TEXT.....Left Neck
ID.....{6A53228A-5BD1-6AAF-029F-701C7EFE7C48}
TERM.....L NECK

DATA VALUE: L NEPHROSTOMY

ABBREVIATION.....L NEPHROST
ACTIVE.....Yes
DESCRIPTION.....Left Nephrostomy
DISPLAY_NAME.....L Nephrostomy
HELP_TEXT.....Left Nephrostomy
ID.....{D8FC4C18-081E-A7B9-E394-25F194BE9757}
TERM.....L NEPHROSTOMY

DATA VALUE: L STENT

ABBREVIATION.....LEFT STENT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Stent
HELP_TEXT.....No Help Text
ID.....{FD917F0B-CC96-0CB0-BED3-1D7C51A77B69}
TERM.....L STENT

DATA VALUE: LACERATION

ABBREVIATION.....LACERATION
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Laceration
HELP_TEXT.....No Help Text
ID.....{F4548BB2-39AF-4E0D-B52F-472EA080D6FF}
TERM.....LACERATION

DATA VALUE: LACKS INTEREST

ABBREVIATION.....LACKS INT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Lacks Interest
HELP_TEXT.....No Help Text
ID.....{585C5B7D-D229-4C7C-AE3C-5869FF845F0D}
TERM.....LACKS INTEREST

Page: 774

DATA VALUE: LACTOSE RESTRICTED

ABBREVIATION.....LACT RESTR
ACTIVE.....Yes
DESCRIPTION.....Lactose Restricted

DISPLAY_NAME.....Lactose Restricted
HELP_TEXT.....No Help Text
ID.....{78A9066A-48E2-4285-88E3-A180CE771B5A}
TERM.....LACTOSE RESTRICTED

DATA VALUE: LANGUAGE

ABBREVIATION.....LANGUAGE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Language
HELP_TEXT.....No Help-Text
ID.....{0666538E-4F44-4D34-A5EA-69CF2F1F86FA}
TERM.....LANGUAGE

DATA VALUE: LARGE

ABBREVIATION.....L
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Large
HELP_TEXT.....No Help Text
ID.....{E60D5CF9-920D-44E2-8A83-136DAB1FE251}
TERM.....LARGE

DATA VALUE: LARYNGECTOMY

ABBREVIATION.....LARYNGECTO
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Laryngectomy
HELP_TEXT.....No Help Text
ID.....{225CB76A-BAF5-827E-729B-C3D60AD5B50A}
TERM.....LARYNGECTOMY

DATA VALUE: LEFT

ABBREVIATION.....LEFT
ACTIVE.....Yes
DESCRIPTION.....Left
DISPLAY_NAME.....Left
HELP_TEXT.....No Help Text
ID.....{81063CEE-2E83-4C95-AAEE-76B7DCB828D4}
TERM.....LEFT

Page: 775

DATA VALUE: LESS THAN 30

ABBREVIATION.....<30
ACTIVE.....Yes
DESCRIPTION.....Less Than 30
DISPLAY_NAME.....<30
HELP_TEXT.....No Help Text
ID.....{5AA98A0D-9B05-425D-97CA-B8B12425AACA}

TERM.....LESS THAN 30

DATA VALUE: LETHARGIC

ABBREVIATION.....L
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Lethargic
HELP_TEXT.....No Help Text
ID.....{4ECA2E44-093A-9354-848E-CDC72B165072}
TERM.....LETHARGIC

DATA VALUE: LEVELED

ABBREVIATION.....LEVELED
ACTIVE.....Yes
DESCRIPTION.....Leveled
DISPLAY_NAME.....Leveled
HELP_TEXT.....No Help Text
ID.....{D72523B2-E0CD-40E1-AC14-2FA43636A456}
TERM.....LEVELED

DATA VALUE: LIBERAL DIABETIC

ABBREVIATION.....LIBRL DIAB
ACTIVE.....Yes
DESCRIPTION.....Liberal Diabetic
DISPLAY_NAME.....Liberal Diabetic
HELP_TEXT.....No Help Text
ID.....{9B66B47B-E77E-4F69-BDFA-5323B28A283D}
TERM.....LIBERAL DIABETIC

DATA VALUE: LIMITED ASSISTANCE

ABBREVIATION.....LIM ASSIST
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Limited Assistance
HELP_TEXT.....No Help Text
ID.....{9A8C1310-C840-431E-BA21-065B7EA7B768}
TERM.....LIMITED ASSISTANCE

Page: 776

DATA VALUE: LIMITED ATTENTION SPAN

ABBREVIATION.....LIM ATT S
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Limited Attention Span
HELP_TEXT.....No Help Text
ID.....{51553020-0CB9-483E-AC1C-15CA6D9380A3}
TERM.....LIMITED ATTENTION SPAN

DATA VALUE: LIMITED TO INCISION

ABBREVIATION.....ER EXI LIM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Limited to Incision
HELP_TEXT.....No Help-Text
ID.....{C9ED5E2C-8ED1-445A-8A3E-C1207CE3843D}
TERM.....LIMITED TO INCISION

DATA VALUE: LINE FLUSH

ABBREVIATION.....LINE FLUSH
ACTIVE.....Yes
DESCRIPTION.....Fluid flushed from line
DISPLAY_NAME.....Line Flush
HELP_TEXT.....No Help Text
ID.....{8DDD63BE-CF28-4C83-93FF-C81609FFDE51}
TERM.....LINE FLUSH

DATA VALUE: LIP SALVE

ABBREVIATION.....LIP SALVE
ACTIVE.....Yes
DESCRIPTION.....Lip Salve
DISPLAY_NAME.....Lip Salve
HELP_TEXT.....No Help Text
ID.....{612F78AB-7386-4176-9896-02DE5A8E836D}
TERM.....LIP SALVE

DATA VALUE: LISTED DIET

ABBREVIATION.....LIST DIET
ACTIVE.....Yes
DESCRIPTION.....Listed Diet
DISPLAY_NAME.....Listed Diet
HELP_TEXT.....No Help Text
ID.....{24FBBFD9-B426-4DD9-9843-3428A1C91676}
TERM.....LISTED DIET

Page: 777

DATA VALUE: LOCALIZES TO PAIN

ABBREVIATION.....LOC TO PAI
ACTIVE.....Yes
DESCRIPTION.....Localizes To Pain
DISPLAY_NAME.....Localizes To Pain
HELP_TEXT.....No Help Text
ID.....{4072DFF9-2D90-423F-B6F0-0F2B924025D9}
TERM.....LOCALIZES TO PAIN

DATA VALUE: LOOSE

ABBREVIATION.....LOOSE
ACTIVE.....Yes

DESCRIPTION.....Loose
DISPLAY_NAME.....Loose
HELP_TEXT.....No Help Text
ID.....{1346A74E-3AA1-4D81-876B-91201F9717A5}
TERM.....LOOSE

DATA VALUE: LOW

ABBREVIATION.....I
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Low
HELP_TEXT.....No Help Text
ID.....{B6CE52A4-C1DD-FDE3-F87E-8E0C98FCC4AC}
TERM.....LOW

DATA VALUE: LOW AIRLOSS OVERLAY

ABBREVIATION.....L AIR OVRL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Low Airloss Overlay
HELP_TEXT.....No Help Text
ID.....{282C439C-A697-43E6-85C1-786352EB56BC}
TERM.....LOW AIRLOSS OVERLAY

DATA VALUE: LOW AIRLOSS/ALTERNATING PRESSURE COMBO

ABBREVIATION.....L AIRL ALT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Low Airloss/Alternating Pressure Combo
HELP_TEXT.....No Help Text
ID.....{ADDD1260-C95B-4D7B-AB89-3749D76A3700}
TERM.....LOW AIRLOSS/ALTERNATING PRESSURE COMBO

Page: 778

DATA VALUE: LOW CHOLESTEROL

ABBREVIATION.....LOW CHOLES
ACTIVE.....Yes
DESCRIPTION.....Low Cholesterol
DISPLAY_NAME.....Low Cholesterol
HELP_TEXT.....No Help Text
ID.....{88F0170A-2560-48C6-9604-982453BCFC3B}
TERM.....LOW CHOLESTEROL

DATA VALUE: LOW FIBER

ABBREVIATION.....LOW FIBER
ACTIVE.....Yes
DESCRIPTION.....Low Fiber
DISPLAY_NAME.....Low Fiber
HELP_TEXT.....No Help Text

ID.....{2834D9A0-B6F7-4C21-96ED-1721D32B3B2B}
TERM.....LOW FIBER

DATA VALUE: LOW PROTEIN 60 GM

ABBREVIATION.....L PROT 60g
ACTIVE.....Yes
DESCRIPTION.....Low Protein 60 GM
DISPLAY_NAME.....Low Protein 60 GM
HELP_TEXT.....No Help Text
ID.....{BF5B6D8B-01A4-42D5-8551-6E78A17C43F7}
TERM.....LOW PROTEIN 60 GM

DATA VALUE: LOW TYRAMINE

ABBREVIATION.....LOW TYRAMI
ACTIVE.....Yes
DESCRIPTION.....Low Tyramine
DISPLAY_NAME.....Low Tyramine
HELP_TEXT.....No Help Text
ID.....{DB067A83-9D0C-4FC2-BAC8-6EC5641CAE42}
TERM.....LOW TYRAMINE

DATA VALUE: LUMBAR DRAIN

ABBREVIATION.....LUM DRAIN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Lumbar Drain
HELP_TEXT.....No Help Text
ID.....{FF36A046-1799-0734-286B-2272C34E3056}
TERM.....LUMBAR DRAIN

Page: 779

DATA VALUE: LUMBAR SPINAL PRECAUTIONS

ABBREVIATION.....LUMB S P
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Lumbar Spinal Precautions
HELP_TEXT.....No Help Text
ID.....{5AFBBBAB-F919-436A-9C5A-71DCF5D79D18}
TERM.....LUMBAR SPINAL PRECAUTIONS

DATA VALUE: LYING OR SITTING

ABBREVIATION.....los
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Lying Or Sitting
HELP_TEXT.....No Help Text
ID.....{52A940B1-55AF-AB4E-9E59-F93515A1296C}
TERM.....LYING OR SITTING

DATA VALUE: MACERATION

ABBREVIATION.....MACERATION
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....MACERATION
HELP_TEXT.....No Help Text
ID.....{B7C63E05-E193-40AA-A572-F1ACA0162B8B}
TERM.....MACERATION

DATA VALUE: MAGNACAL

ABBREVIATION.....MAGNACAL
ACTIVE.....Yes
DESCRIPTION.....Enteral Nutrition Type
DISPLAY_NAME.....Magnacal
HELP_TEXT.....No Help Text
ID.....{C810A618-FB1A-45BB-A540-EADF1414E343}
TERM.....MAGNACAL

DATA VALUE: MAINTAINING EFFECTIVE VENTILATION

ABBREVIATION.....EFF VENT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Maintaining Effective Ventilation
HELP_TEXT.....No Help Text
ID.....{18774550-ECBD-412D-AD95-131FAA4CC549}
TERM.....MAINTAINING EFFECTIVE VENTILATION

Page: 780

DATA VALUE: MAROON

ABBREVIATION.....MAROON
ACTIVE.....Yes
DESCRIPTION.....Maroon in appearance
DISPLAY_NAME.....Maroon
HELP_TEXT.....No Help Text
ID.....{1F518B23-C79F-4F06-9DF3-ADC5C96D3029}
TERM.....MAROON

DATA VALUE: MEAL

ABBREVIATION.....MEAL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Meal
HELP_TEXT.....No Help Text
ID.....{85D38DCD-6A4A-4908-ABE2-4E92ABEF1321}
TERM.....MEAL

DATA VALUE: MECHANICAL

ABBREVIATION.....MECHANICAL

ACTIVE.....Yes
DESCRIPTION.....Mechanical
DISPLAY_NAME.....Mechanical
HELP_TEXT.....No Help Text
ID.....{01EA793C-B524-4685-A3EC-60A317CBACC5}
TERM.....MECHANICAL

DATA VALUE: MED RECONCILLIATION DOCUMENTED

ABBREVIATION.....ME RE DO
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Med Reconciliation Documented
HELP_TEXT.....No Help Text
ID.....{9ABF3E2E-4D63-43EF-A80D-9DCF5C820226}
TERM.....MED RECONCILLIATION DOCUMENTED

DATA VALUE: MEDIAL PORT

ABBREVIATION.....MEDI PORT
ACTIVE.....Yes
DESCRIPTION.....Medial Port
DISPLAY_NAME.....Medial Port
HELP_TEXT.....No Help Text
ID.....{D1734D73-2E1A-4538-8884-05D8A78CF556}
TERM.....MEDIAL PORT

Page: 781

DATA VALUE: MEDICAL DEVICE SAFETY

ABBREVIATION.....mds
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Medical Device Safety
HELP_TEXT.....No Help Text
ID.....{1985CABE-4F53-9FCF-F444-F9AC53C04DBD}
TERM.....MEDICAL DEVICE SAFETY

DATA VALUE: MEDICATION

ABBREVIATION.....MEDICATION
ACTIVE.....Yes
DESCRIPTION.....Medication
DISPLAY_NAME.....Medication
HELP_TEXT.....No Help Text
ID.....{22E6DEBB-9662-466B-990C-8DA500C0A183}
TERM.....MEDICATION

DATA VALUE: MEDICATION ROOM

ABBREVIATION.....mr
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Medication Room

HELP_TEXT.....No Help Text
ID.....{9571AA9F-F538-55EC-37F6-99FE5677B970}
TERM.....MEDICATION ROOM

DATA VALUE: MEDIUM

ABBREVIATION.....MEDIUM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Medium
HELP_TEXT.....No Help Text
ID.....{62C94066-A7B6-CAA8-064B-2AF7C0DA6323}
TERM.....MEDIUM

DATA VALUE: MEMORY

ABBREVIATION.....MEMORY
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Memory
HELP_TEXT.....No Help Text
ID.....{E30C1930-B98B-4382-B892-34E9C6008D1F}
TERM.....MEMORY

Page: 782

DATA VALUE: MENTAL HEALTH

ABBREVIATION.....MEN HEAL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Mental Health
HELP_TEXT.....No Help Text
ID.....{ED55FD2C-E2EC-4458-9AE4-FD7F0E6BD7F6}
TERM.....MENTAL HEALTH

DATA VALUE: MENTAL HEALTH HOLD

ABBREVIATION.....MEN HEA HO
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Mental Health Hold
HELP_TEXT.....No Help Text
ID.....{31560285-7BD5-4E53-BB77-93B244367F42}
TERM.....MENTAL HEALTH HOLD

DATA VALUE: MICU

ABBREVIATION.....MICU
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....MICU
HELP_TEXT.....No Help Text
ID.....{28244B62-644C-4DA0-A350-EAFE0D2FA7C0}
TERM.....MICU

DATA VALUE: MILD-MODERATELY REDUCED

ABBREVIATION.....MI-MOD-RED
ACTIVE.....Yes
DESCRIPTION.....Mild-Moderately Reduced
DISPLAY_NAME.....Mild-Moderately Reduced
HELP_TEXT.....No Help Text
ID.....{41DA2518-0E49-4E92-A733-3ACF3D2E8199}
TERM.....MILD-MODERATELY REDUCED

DATA VALUE: MILLIGRAMS PER KILOGRAM PER HOUR

ABBREVIATION.....mg/kg/hr
ACTIVE.....Yes
DESCRIPTION.....Milligrams per Kilogram per Hour
DISPLAY_NAME.....mg/kg/hr
HELP_TEXT.....Milligrams per Kilogram per Hour
ID.....{731B270D-C8C5-4910-B4B0-AC101555108A}
TERM.....MILLIGRAMS PER KILOGRAM PER HOUR

Page: 783

DATA VALUE: MISSING

ABBREVIATION.....MISSING
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Missing
HELP_TEXT.....No Help Text
ID.....{80019911-11DA-6569-751E-ECAE0D1BC7DF}
TERM.....MISSING

DATA VALUE: MOANING

ABBREVIATION.....MOANING
ACTIVE.....Yes
DESCRIPTION.....Moaning
DISPLAY_NAME.....Moaning
HELP_TEXT.....No Help Text
ID.....{A71E3515-AA20-4B0E-9CB8-1612F76BA406}
TERM.....MOANING

DATA VALUE: MODERATE

ABBREVIATION.....MODERATE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Moderate
HELP_TEXT.....No Help Text
ID.....{C3C25064-C9FC-EA1F-3E74-32A25F0FC5E3}
TERM.....MODERATE

DATA VALUE: MODERATE-SEVERELY REDUCED

ABBREVIATION.....MO-SEV RED
ACTIVE.....Yes
DESCRIPTION.....Moderate-Severely Reduced
DISPLAY_NAME.....Moderate-Severely Reduced
HELP_TEXT.....No Help Text
ID.....{6D031C70-2F95-4111-A9B8-BC22CEC1B17D}
TERM.....MODERATE-SEVERELY REDUCED

DATA VALUE: MODERATELY REDUCED

ABBREVIATION.....MOD REDUCE
ACTIVE.....Yes
DESCRIPTION.....Moderately Reduced
DISPLAY_NAME.....Moderately Reduced
HELP_TEXT.....No Help Text
ID.....{91F73300-7EE0-428C-A838-5B3BE062D8BC}
TERM.....MODERATELY REDUCED

Page: 784

DATA VALUE: MOIST

ABBREVIATION.....m
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Moist
HELP_TEXT.....No Help Text
ID.....{AA73A9D6-6645-5DA7-C602-A537A2DFBFA0}
TERM.....MOIST

DATA VALUE: MOTORIZED WHEELCHAIR

ABBREVIATION.....MO WH CH
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Motorized Wheelchair
HELP_TEXT.....No Help Text
ID.....{4CA9CF41-F52C-41A8-B8A7-100A56097FDE}
TERM.....MOTORIZED WHEELCHAIR

DATA VALUE: MOTTLED

ABBREVIATION.....m
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Mottled
HELP_TEXT.....No Help Text
ID.....{4B5FFB10-2CFF-2064-1553-4CBB7133B0C5}
TERM.....MOTTLED

DATA VALUE: MOUTHWASH USED

ABBREVIATION.....MOUTHWASH
ACTIVE.....Yes
DESCRIPTION.....Mouthwash Used

DISPLAY_NAME.....Mouthwash Used
HELP_TEXT.....No Help Text
ID.....{9F7EF355-12A4-4EBA-92ED-0F38273E411B}
TERM.....MOUTHWASH USED

DATA VALUE: MOVE AGAINST GRAVITY AND RESISTANCE, EASILY OVERCOME BY THE EXA
MINER

ABBREVIATION.....magareobte
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Movement Against Gravity And
Resistance, Easily Overcome By The
Examiner
HELP_TEXT.....No Help Text
ID.....{2ECCC197-4A81-025F-3EF3-6214004887A9}
TERM.....MOVE AGAINST GRAVITY AND RESISTANCE,
EASILY OVERCOME BY THE EXAMINER

Page: 785

DATA VALUE: MOVE AGAINST GRAVITY BUT NO RESISTANCE

ABBREVIATION.....magbnr
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Movement Against Gravity But No
Resistance
HELP_TEXT.....No Help Text
ID.....{F995A165-8536-C493-2269-286E8EFF5E43}
TERM.....MOVE AGAINST GRAVITY BUT NO RESISTANCE

DATA VALUE: MOVEMENT WITH GRAVITY ELIMINATED

ABBREVIATION.....mwge
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Movement With Gravity Eliminated
HELP_TEXT.....No Help Text
ID.....{B86E1131-1976-FA8E-74FE-B1BC2F0FDE3F}
TERM.....MOVEMENT WITH GRAVITY ELIMINATED

DATA VALUE: MUCOID

ABBREVIATION.....MUCOID
ACTIVE.....Yes
DESCRIPTION.....Muroid in appearance
DISPLAY_NAME.....Muroid
HELP_TEXT.....No Help Text
ID.....{073E5514-DACE-4F43-ADE5-91A1A87D6601}
TERM.....MUCOID

DATA VALUE: MULTIPODUS

ABBREVIATION.....MULTIPODUS

ACTIVE.....Yes
DESCRIPTION.....Multipodus
DISPLAY_NAME.....Multipodus
HELP_TEXT.....No Help Text
ID.....{ED098F9D-B204-4310-813A-98D91FF9DD25}
TERM.....MULTIPODUS

DATA VALUE: MURMUR

ABBREVIATION.....MURMUR
ACTIVE.....Yes
DESCRIPTION.....MURMUR
DISPLAY_NAME.....Murmur
HELP_TEXT.....MURMUR
ID.....{591058E0-4B71-E432-B7D7-EBD42A8C4225}
TERM.....MURMUR

Page: 786

DATA VALUE: MUSHROOM CATHETER

ABBREVIATION.....MUSH CATHE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Mushroom Catheter
HELP_TEXT.....No Help Text
ID.....{4FFC451E-14B5-46E4-99C9-99A17BF90AE4}
TERM.....MUSHROOM CATHETER

DATA VALUE: N/A

ABBREVIATION.....n
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....N/A
HELP_TEXT.....No Help Text
ID.....{B3DBE3FD-C5CB-3894-607A-BE13319594D0}
TERM.....N/A

DATA VALUE: NAIL CARE

ABBREVIATION.....NAIL CARE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Nail Care
HELP_TEXT.....No Help Text
ID.....{B783A171-1F0A-4135-8535-85776B81AB68}
TERM.....NAIL CARE

DATA VALUE: NAPS ONLY

ABBREVIATION.....NAPS ONLY
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Naps Only

HELP_TEXT.....No Help Text
ID.....{CF749F5C-2E4B-419A-B020-C81B30A7E482}
TERM.....NAPS ONLY

DATA VALUE: NASAL

ABBREVIATION.....NASAL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Nasal
HELP_TEXT.....No Help Text
ID.....{0F73F0CB-E25E-4611-9625-62E6F1596A9B}
TERM.....NASAL

Page: 787

DATA VALUE: NASAL ENDOTRACHEAL

ABBREVIATION.....Nasal Endo
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Nasal Endotracheal
HELP_TEXT.....No Help Text
ID.....{C631352E-CAB5-0409-6133-E5D11ECC3293}
TERM.....NASAL ENDOTRACHEAL

DATA VALUE: NASAL INTESTINAL SMALL BORE

ABBREVIATION.....N INT SM B
ACTIVE.....Yes
DESCRIPTION.....Nasal Intestinal Small Bore EF Tube
DISPLAY_NAME.....Nasal Intestinal Small Bore
HELP_TEXT.....No Help Text
ID.....{03DCB9EB-9B11-4C3B-869B-B1842BBCC844}
TERM.....NASAL INTESTINAL SMALL BORE

DATA VALUE: NASAL SMALL BOWEL FEEDING TUBE

ABBREVIATION.....NAS SM BOW
ACTIVE.....Yes
DESCRIPTION.....Nasal Small Bowel Feeding Tube
DISPLAY_NAME.....Nasal Small Bowel Feeding Tube
HELP_TEXT.....No Help Text
ID.....{A2968C8A-DDBD-4C3B-AFE6-3E2FED41C413}
TERM.....NASAL SMALL BOWEL FEEDING TUBE

DATA VALUE: NASAL TRUMPET

ABBREVIATION.....NASAL TRUM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Nasal Trumpet
HELP_TEXT.....No Help Text
ID.....{2FFB8510-793C-D6A4-838E-20FF7CD26533}
TERM.....NASAL TRUMPET

DATA VALUE: NASOGASTRIC LARGE BORE TUBE

ABBREVIATION.....NG LG BORE
ACTIVE.....Yes
DESCRIPTION.....Nasogastric Large Bore EF Tube
DISPLAY_NAME.....Nasogastric Large Bore
HELP_TEXT.....No Help Text
ID.....{C3FEE16A-9998-455B-BEB8-CC55AF890202}
TERM.....NASOGASTRIC LARGE BORE TUBE

Page: 788

DATA VALUE: NASOGASTRIC TUBE

ABBREVIATION.....NG TUBE
ACTIVE.....Yes
DESCRIPTION.....Nasogastric Tube
DISPLAY_NAME.....Nasogastric Tube
HELP_TEXT.....No Help Text
ID.....{6475E343-A310-48A5-B543-D29573E7684D}
TERM.....NASOGASTRIC TUBE

DATA VALUE: NASOGASTRIC TUBE SMALL BORE

ABBREVIATION.....NGT SMALL
ACTIVE.....Yes
DESCRIPTION.....NGT Small Bore Tube
DISPLAY_NAME.....Nasal Gastric Tube Small Bore
HELP_TEXT.....No Help Text
ID.....{D9FDAD80-9BE2-4B96-87C6-1D25D26D2959}
TERM.....NASOGASTRIC TUBE SMALL BORE

DATA VALUE: NAUSEA

ABBREVIATION.....NAUSEA
ACTIVE.....Yes
DESCRIPTION.....Nausea
DISPLAY_NAME.....Nausea
HELP_TEXT.....Nausea
ID.....{8BF33DC3-0300-4643-BB14-27C35692B806}
TERM.....NAUSEA

DATA VALUE: NEAT

ABBREVIATION.....NEAT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Neat
HELP_TEXT.....No Help Text
ID.....{8B01F7B8-B626-4143-A092-1E732EEB4A05}
TERM.....NEAT

DATA VALUE: NECROTIC

ABBREVIATION.....NECROTIC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Necrotic
HELP_TEXT.....No Help Text
ID.....{86885A34-E4EC-4BB5-BD94-59DD4B0F1A0E}
TERM.....NECROTIC

Page: 789

DATA VALUE: NEEDS REINFORCEMENT

ABBREVIATION.....NEED REINF
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Needs Reinforcement
HELP_TEXT.....No Help Text
ID.....{B6171F65-4530-418D-A63B-65DFFF522786}
TERM.....NEEDS REINFORCEMENT

DATA VALUE: NEGATIVE PRESSURE

ABBREVIATION.....NEG PRESSU
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Negative Pressure
HELP_TEXT.....No Help Text
ID.....{6D6CB4BB-96B0-4978-97CC-116A03479DAD}
TERM.....NEGATIVE PRESSURE

DATA VALUE: NEGATIVE PRESSURE THERAPY DEVICE

ABBREVIATION.....NEG PRE TH
ACTIVE.....Yes
DESCRIPTION.....Negative Pressure Therapy Device
DISPLAY_NAME.....Negative Pressure Therapy Device
HELP_TEXT.....No Help Text
ID.....{8AF0A069-66AD-4114-99D8-238AD987F0FF}
TERM.....NEGATIVE PRESSURE THERAPY DEVICE

DATA VALUE: NEPRO RENAL FORMULA

ABBREVIATION.....NEP RENAL
ACTIVE.....Yes
DESCRIPTION.....Nepro Renal Enteral Nutrition Formula
DISPLAY_NAME.....Nepro Renal Enteral Nutrition Formula
HELP_TEXT.....No Help Text
ID.....{B8F28816-BC81-443F-AEC5-DBAE4D435998}
TERM.....NEPRO RENAL FORMULA

DATA VALUE: NEUTROPENIC

ABBREVIATION.....NEUTROPENI
ACTIVE.....Yes
DESCRIPTION.....Neutropenic Diet

DISPLAY_NAME.....Neutropenic
HELP_TEXT.....No Help Text
ID.....{5CB9620C-34D9-48F7-AAF0-C595D06FC8AE}
TERM.....NEUTROPENIC

Page: 790

DATA VALUE: NEW ADMISSION EDUCATION

ABBREVIATION.....N ADM EDU
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....New Admission Education
HELP_TEXT.....No Help Text
ID.....{1890BC97-FEEB-4BED-A054-3C754C01817B}
TERM.....NEW ADMISSION EDUCATION

DATA VALUE: NEW SITE

ABBREVIATION.....NEW SITE
ACTIVE.....Yes
DESCRIPTION.....New Site
DISPLAY_NAME.....New Site
HELP_TEXT.....No Help Text
ID.....{10B1FF61-E58E-4608-9CCB-9232BC964418}
TERM.....NEW SITE

DATA VALUE: NO

ABBREVIATION.....n
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....No
HELP_TEXT.....No Help Text
ID.....{FC275445-1225-9912-18A3-DF36245893D3}
TERM.....NO

DATA VALUE: NO BREAK IN EPIDERMIS

ABBREVIATION.....NO BREAK E
ACTIVE.....Yes
DESCRIPTION.....No Break In Epidermis
DISPLAY_NAME.....No Break In Epidermis
HELP_TEXT.....No Help Text
ID.....{615D96CD-6096-446B-BDDF-3F7242936AF9}
TERM.....NO BREAK IN EPIDERMIS

DATA VALUE: NO MOVEMENT

ABBREVIATION.....nm
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....No Movement
HELP_TEXT.....No Help Text
ID.....{2618E301-379B-BAB9-9EDE-98FD4B60D311}

TERM.....NO MOVEMENT

Page: 791

DATA VALUE: NO SIGNS AND SYMPTOMS

ABBREVIATION.....NO SIG&SYM
ACTIVE.....Yes
DESCRIPTION.....No Signs And Symptoms
DISPLAY_NAME.....No Signs And Symptoms
HELP_TEXT.....No Help Text
ID.....{91285942-6C82-4CB1-AAFB-8D3DA15A0B8A}
TERM.....NO SIGNS AND SYMPTOMS

DATA VALUE: NON-ADHERENT GAUZE

ABBREVIATION.....NON-AD GAU
ACTIVE.....Yes
DESCRIPTION.....Non-Adherent Gauze
DISPLAY_NAME.....Non-Adherent gauze
HELP_TEXT.....No Help Text
ID.....{1DBF7C56-F5C2-42B5-AB59-B39A5069A414}
TERM.....NON-ADHERENT GAUZE

DATA VALUE: NON-PITTING

ABBREVIATION.....NON-PITTIN
ACTIVE.....Yes
DESCRIPTION.....NON-PITTING
DISPLAY_NAME.....Non-Pitting
HELP_TEXT.....NON-PITTING
ID.....{786DEC76-36DE-FDC7-885C-F84210F7C7D2}
TERM.....NON-PITTING

DATA VALUE: NON-PRODUCTIVE

ABBREVIATION.....np
ACTIVE.....Yes
DESCRIPTION.....Does not produce desired result
DISPLAY_NAME.....Non-Productive
HELP_TEXT.....Does not produce desired result
ID.....{9AB01716-E395-18B9-6C5D-9BE6AEE41D5B}
TERM.....NON-PRODUCTIVE

DATA VALUE: NON-TENDER

ABBREVIATION.....NON-TENDER
ACTIVE.....Yes
DESCRIPTION.....Not Tender
DISPLAY_NAME.....Non-Tender
HELP_TEXT.....No Help Text
ID.....{5A3D41C0-A83D-4382-A573-D6884BB1670B}
TERM.....NON-TENDER

DATA VALUE: NONBLANCHABLE ERYTHEMA

ABBREVIATION.....N BLAN ERY
ACTIVE.....Yes
DESCRIPTION.....Nonblanchable Erythema
DISPLAY_NAME.....Nonblanchable Erythema
HELP_TEXT.....No Help Text
ID.....{30273ACE-03F3-4818-995B-26586AAB5499}
TERM.....NONBLANCHABLE ERYTHEMA

DATA VALUE: NONE

ABBREVIATION.....NONE
ACTIVE.....Yes
DESCRIPTION.....None
DISPLAY_NAME.....None
HELP_TEXT.....No Help Text
ID.....{D3118C4E-FEB4-4E9E-B191-9A31A7A0BA65}
TERM.....NONE

DATA VALUE: NONREACTIVE

ABBREVIATION.....nonreact
ACTIVE.....Yes
DESCRIPTION.....Nonreactive
DISPLAY_NAME.....Nonreactive
HELP_TEXT.....No Help Text
ID.....{AB6B7D92-B94C-F3AA-E32F-6035D9CF2EF9}
TERM.....NONREACTIVE

DATA VALUE: NORMAL

ABBREVIATION.....n
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Normal
HELP_TEXT.....No Help Text
ID.....{553C9A97-A3D5-C52F-0E4E-C7D168C99A2B}
TERM.....NORMAL

DATA VALUE: NORMAL FOR ETHNIC GROUP

ABBREVIATION.....ETH. NORM.
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Normal for ethnic group
HELP_TEXT.....Skin color is normal for ethnic group
ID.....{3B2CD951-58CE-FF53-AD5D-6C65C07658D3}
TERM.....NORMAL FOR ETHNIC GROUP

DATA VALUE: NORMAL POWER

ABBREVIATION.....np
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Normal Power
HELP_TEXT.....No Help Text
ID.....{0E45C061-0744-8A21-7256-003702A5005B}
TERM.....NORMAL POWER

DATA VALUE: NORMAL SINUS

ABBREVIATION.....ns
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Normal Sinus
HELP_TEXT.....No Help Text
ID.....{735DC86D-08CD-1533-5A36-4775A882C507}
TERM.....NORMAL SINUS

DATA VALUE: NORMAL STOOL

ABBREVIATION.....NORM STOOL
ACTIVE.....Yes
DESCRIPTION.....Normal Stool
DISPLAY_NAME.....Normal Stool
HELP_TEXT.....No Help Text
ID.....{486B0693-CAC5-48EE-B2E6-1F73C651418A}
TERM.....NORMAL STOOL

DATA VALUE: NOT DONE

ABBREVIATION.....NOT DONE
ACTIVE.....Yes
DESCRIPTION.....Not Done
DISPLAY_NAME.....Not Done
HELP_TEXT.....No Help Text
ID.....{C3FE84AC-EEA1-4153-AE6C-69755C2E6BB1}
TERM.....NOT DONE

DATA VALUE: NOT MONITORED

ABBREVIATION.....NOT MONITO
ACTIVE.....Yes
DESCRIPTION.....Not Monitored
DISPLAY_NAME.....Not Monitored
HELP_TEXT.....No Help Text
ID.....{0E6A8A07-48E0-428F-8D86-22572CABCF48}
TERM.....NOT MONITORED

Page: 794

DATA VALUE: NOVASOURCE 2.0

ABBREVIATION.....NOVAS 2.0
ACTIVE.....Yes

DESCRIPTION.....Novasource 2.0 Enteral Nutrition
Formula
DISPLAY_NAME.....Novasource 2.0 Enteral Nutrition
Formula
HELP_TEXT.....No Help Text
ID.....{ACE82591-502A-4FBA-9697-B1B330DED78A}
TERM.....NOVASOURCE 2.0

DATA VALUE: NPO DIET

ABBREVIATION.....NPO DIET
ACTIVE.....Yes
DESCRIPTION.....Nothing By Mouth Diet
DISPLAY_NAME.....NPO Diet
HELP_TEXT.....Nothing By Mouth Diet
ID.....{BEC8F771-B082-40A2-9EE8-553222A25FBF}
TERM.....NPO DIET

DATA VALUE: NUMB

ABBREVIATION.....NUMB
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....NUMB
HELP_TEXT.....No Help Text
ID.....{652634C0-5334-4D56-BE57-5F58F952FF0D}
TERM.....NUMB

DATA VALUE: NUMBNESS

ABBREVIATION.....n
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Numbness
HELP_TEXT.....No Help Text
ID.....{037F2F79-36D9-A0B3-C23E-8EE525932AF3}
TERM.....NUMBNESS

DATA VALUE: NURSING STATION

ABBREVIATION.....ns
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Nursing Station
HELP_TEXT.....No Help Text
ID.....{B6E48F92-B9BA-488F-3A2E-D74568A6B855}
TERM.....NURSING STATION

Page: 795

DATA VALUE: NUTRITION

ABBREVIATION.....NUTRITION
ACTIVE.....Yes
DESCRIPTION.....No Description

DISPLAY_NAME.....Nutrition
HELP_TEXT.....No Help Text
ID.....{429A28EF-8A1F-4247-85E8-68D4EFCD189C}
TERM.....NUTRITION

DATA VALUE: NUTRITIONAL SUPPLEMENT

ABBREVIATION.....NUTR SUPPL
ACTIVE.....Yes
DESCRIPTION.....Nutritional Supplement
DISPLAY_NAME.....Nutritional Supplement
HELP_TEXT.....No Help Text
ID.....{4573BD4E-7220-4B88-A2EE-959F8A0E1199}
TERM.....NUTRITIONAL SUPPLEMENT

DATA VALUE: OBEYS COMMANDS

ABBREVIATION.....OBEYS COMM
ACTIVE.....Yes
DESCRIPTION.....Obeys Commands
DISPLAY_NAME.....Obeys Commands
HELP_TEXT.....No Help Text
ID.....{EE88EE08-4C71-4C70-B93E-65F9091544C1}
TERM.....OBEYS COMMANDS

DATA VALUE: OFF DURING BATHING

ABBREVIATION.....OFF BATHIN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Off During Bathing
HELP_TEXT.....No Help Text
ID.....{5D0EFDD9-F170-4F34-9ABE-CE2D376E65F4}
TERM.....OFF DURING BATHING

DATA VALUE: OFF WHILE AMBULATORY

ABBREVIATION.....OFF AMBULA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Off While Ambulatory
HELP_TEXT.....No Help Text
ID.....{FF4610BF-8C92-408F-882A-1D6B2C2F146F}
TERM.....OFF WHILE AMBULATORY

Page: 796

DATA VALUE: OFFER PRN MEDS

ABBREVIATION.....OFFER PRN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Offer PRN Meds
HELP_TEXT.....No Help Text
ID.....{B7AF4903-913C-4848-A255-0B4767B81AE2}

TERM.....OFFER PRN MEDS

DATA VALUE: OFFERED

ABBREVIATION.....OFFERED
ACTIVE.....Yes
DESCRIPTION.....Offered
DISPLAY_NAME.....Offered
HELP_TEXT.....No Help Text
ID.....{715091C2-9668-4678-900E-6D55D3335A30}
TERM.....OFFERED

DATA VALUE: ON WHILE SITTING IN CHAIR

ABBREVIATION.....ON IN CHAI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....On While Sitting In Chair
HELP_TEXT.....No Help Text
ID.....{4139F4C3-BB68-4477-BF63-9BE7813AADA7}
TERM.....ON WHILE SITTING IN CHAIR

DATA VALUE: ONE PERSON ASSIST

ABBREVIATION.....1 PERS AS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....One Person Assist
HELP_TEXT.....No Help Text
ID.....{A3D14AA8-D658-412F-BBE8-D8AD195ADA3D}
TERM.....ONE PERSON ASSIST

DATA VALUE: ONE TO ONE SUPERVISION

ABBREVIATION.....1 SUPERV
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....One To One Supervision
HELP_TEXT.....No Help Text
ID.....{CFC9C653-BFE9-40B2-BA7A-944116C7FCB7}
TERM.....ONE TO ONE SUPERVISION

Page: 797

DATA VALUE: OPEN SPONTANEOUSLY

ABBREVIATION.....OPEN SPONT
ACTIVE.....Yes
DESCRIPTION.....Open Spontaneously
DISPLAY_NAME.....Open Spontaneously
HELP_TEXT.....No Help Text
ID.....{013F1473-BDC9-4E7C-9709-C577357709F6}
TERM.....OPEN SPONTANEOUSLY

DATA VALUE: OPEN TO AIR

ABBREVIATION.....OPE TO AIR
ACTIVE.....Yes
DESCRIPTION.....Open To Air
DISPLAY_NAME.....Open To Air
HELP_TEXT.....No Help Text
ID.....{EFB715D0-5666-4281-9034-A0326A14B2BB}
TERM.....OPEN TO AIR

DATA VALUE: OPERATING ROOM

ABBREVIATION.....OR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Operating Room
HELP_TEXT.....No Help Text
ID.....{37DACE14-2A28-45F4-9DEA-35A1A433F4E6}
TERM.....OPERATING ROOM

DATA VALUE: ORAL

ABBREVIATION.....ORAL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Oral
HELP_TEXT.....No Help Text
ID.....{BBAF0DD0-1F52-49EE-9A47-FA9323D6C529}
TERM.....ORAL

DATA VALUE: ORAL AIRWAY

ABBREVIATION.....ORAL AIRWA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Oral Airway
HELP_TEXT.....No Help Text
ID.....{B0CF4270-A68D-1322-0621-1DFF8D677B40}
TERM.....ORAL AIRWAY

Page: 798

DATA VALUE: ORAL ENDOTRACHEAL

ABBREVIATION.....ORAL ENDO
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Oral Endotracheal
HELP_TEXT.....No Help Text
ID.....{466CAD16-D7F3-9D73-6500-959C4DA50D9A}
TERM.....ORAL ENDOTRACHEAL

DATA VALUE: ORAL GASTRIC TUBE

ABBREVIATION.....OG TUBE
ACTIVE.....Yes

DESCRIPTION.....Oral Gastric Enteral Feeding Tube
DISPLAY_NAME.....Oral Gastric Tube
HELP_TEXT.....No Help Text
ID.....{C6995D6A-939C-4C31-A3C2-4504CB5597CD}
TERM.....ORAL GASTRIC TUBE

DATA VALUE: ORANGE

ABBREVIATION.....ORANGE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Orange
HELP_TEXT.....No Help Text
ID.....{77B1B667-212D-88BA-D8E6-5C22F6ACDC9C}
TERM.....ORANGE

DATA VALUE: ORIENTATION

ABBREVIATION.....ORIENTA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Orientation
HELP_TEXT.....No Help Text
ID.....{10FA5B0E-F4A7-44E4-A482-7C80F145C4D2}
TERM.....ORIENTATION

DATA VALUE: ORIENTED

ABBREVIATION.....ORIENTED
ACTIVE.....Yes
DESCRIPTION.....Oriented
DISPLAY_NAME.....Oriented
HELP_TEXT.....No Help Text
ID.....{B8DC9245-6D81-40E4-8B47-B6413E49F152}
TERM.....ORIENTED

Page: 799

DATA VALUE: ORIENTED TO PERSON

ABBREVIATION.....otp
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Oriented To Person
HELP_TEXT.....No Help Text
ID.....{9912CFB9-AC55-6445-3EC1-0F1DF98C0561}
TERM.....ORIENTED TO PERSON

DATA VALUE: ORIENTED TO PLACE

ABBREVIATION.....otp
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Oriented To Place
HELP_TEXT.....No Help Text

ID.....{8E5933E2-2DF0-C8A4-7A92-20E792F42EA6}
TERM.....ORIENTED TO PLACE

DATA VALUE: ORIENTED TO SITUATION

ABBREVIATION.....ots
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Oriented To Situation
HELP_TEXT.....No Help Text
ID.....{955A2F3B-E444-99DF-81DD-5E6BF9D22E49}
TERM.....ORIENTED TO SITUATION

DATA VALUE: ORIENTED TO TIME

ABBREVIATION.....ott
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Oriented To Time
HELP_TEXT.....No Help Text
ID.....{C79EBE25-C51B-478C-CB39-25943792F191}
TERM.....ORIENTED TO TIME

DATA VALUE: ORTHOPNEIC

ABBREVIATION.....ORTHOPNIC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Orthopneic
HELP_TEXT.....No Help Text
ID.....{FF373BD7-FD84-DF43-7A44-1B8E62917116}
TERM.....ORTHOPNEIC

Page: 800

DATA VALUE: ORTHOTIC BOOTS

ABBREVIATION.....ORTHO BOOT
ACTIVE.....Yes
DESCRIPTION.....ORTHOTIC BOOTS
DISPLAY_NAME.....Orthotic Boots
HELP_TEXT.....ORTHOTIC BOOTS
ID.....{D4219FE9-7087-B881-6DAC-4E0BB94CDA8E}
TERM.....ORTHOTIC BOOTS

DATA VALUE: ORTHOTIC COLLAR

ABBREVIATION.....ORTHO COLL
ACTIVE.....Yes
DESCRIPTION.....ORTHOTIC COLLAR
DISPLAY_NAME.....Orthotic Collar
HELP_TEXT.....ORTHOTIC COLLAR
ID.....{A74E1C89-2304-0544-9A8B-A23E4E2C6B08}
TERM.....ORTHOTIC COLLAR

DATA VALUE: ORTHOTIC KNEE IMMOBILIZER

ABBREVIATION.....KNEE IMMOB
ACTIVE.....Yes
DESCRIPTION.....ORTHOTIC KNEE IMMOBILIZER
DISPLAY_NAME.....Orthotic Knee Immobilizer
HELP_TEXT.....ORTHOTIC KNEE IMMOBILIZER
ID.....{FF8A8BB4-63FB-3B0D-83E1-5432A0434DAC}
TERM.....ORTHOTIC KNEE IMMOBILIZER

DATA VALUE: ORTHOTIC SPINE BRACE

ABBREVIATION.....ORTH SPINE
ACTIVE.....Yes
DESCRIPTION.....ORTHOTIC SPINE BRACE
DISPLAY_NAME.....Orthotic Spine Brace
HELP_TEXT.....ORTHOTIC SPINE BRACE
ID.....{B47A3536-FEBC-802E-01A1-E909147D18C2}
TERM.....ORTHOTIC SPINE BRACE

DATA VALUE: ORTHOTIC TED HOSE

ABBREVIATION.....TED HOSE
ACTIVE.....Yes
DESCRIPTION.....TED HOSE
DISPLAY_NAME.....Orthotic TED Hose
HELP_TEXT.....TED HOSE
ID.....{FE277269-7BD1-85FB-DC3C-FEDE31C66BB0}
TERM.....ORTHOTIC TED HOSE

Page: 801

DATA VALUE: ORTHOTIC WEDGE PILLOW

ABBREVIATION.....WEDG PILLO
ACTIVE.....Yes
DESCRIPTION.....WEDGE PILLOW
DISPLAY_NAME.....Orthotic Wedge Pillow
HELP_TEXT.....WEDGE PILLOW
ID.....{20E58D21-399B-85E4-BDAB-93B371395BCF}
TERM.....ORTHOTIC WEDGE PILLOW

DATA VALUE: OSMOLITE

ABBREVIATION.....OSMOLITE
ACTIVE.....Yes
DESCRIPTION.....Enteral Nutrition Type
DISPLAY_NAME.....Osmolite
HELP_TEXT.....No Help Text
ID.....{8E41BB55-063E-43C0-BAFE-65150EDD99F8}
TERM.....OSMOLITE

DATA VALUE: OSMOLITE HN

ABBREVIATION.....OSMOLIT HN

ACTIVE.....Yes
DESCRIPTION.....Enteral Nutrition Type
DISPLAY_NAME.....Osmolite HN
HELP_TEXT.....No Help Text
ID.....{E4186CA6-302F-4513-B125-30B094151423}
TERM.....OSMOLITE HN

DATA VALUE: OTHER

ABBREVIATION.....OTHER
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Other
HELP_TEXT.....No Help Text
ID.....{19094E5B-AEC3-44B0-9BD5-44E79E7EB105}
TERM.....OTHER

DATA VALUE: OVER WIRE

ABBREVIATION.....OVER WIRE
ACTIVE.....Yes
DESCRIPTION.....Over Wire
DISPLAY_NAME.....Over Wire
HELP_TEXT.....No Help Text
ID.....{157CD456-1020-43AE-9AAA-1686F1574B6F}
TERM.....OVER WIRE

Page: 802

DATA VALUE: OXIMETRIC

ABBREVIATION.....OXIMETRIC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Oximetric Catheter
HELP_TEXT.....No Help Text
ID.....{DA1F5AAA-A349-C17B-74CD-1794492E0A13}
TERM.....OXIMETRIC

DATA VALUE: PAC

ABBREVIATION.....p
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....PAC
HELP_TEXT.....No Help Text
ID.....{848AB18A-598D-8E4E-8B1B-ED7478CE340B}
TERM.....PAC

DATA VALUE: PACED: ATRIAL

ABBREVIATION.....pa
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Paced: Atrial

HELP_TEXT.....No Help Text
ID.....{88837DD4-098D-B77F-A754-89D5C4F5870F}
TERM.....PACED: ATRIAL

DATA VALUE: PACED: AV SEQUENTIAL

ABBREVIATION.....pas
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Paced: Av Sequential
HELP_TEXT.....No Help Text
ID.....{3C485674-224C-47DA-510D-92F1A0CA54A1}
TERM.....PACED: AV SEQUENTIAL

DATA VALUE: PACED: VENTRICULAR

ABBREVIATION.....pv
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Paced: Ventricular
HELP_TEXT.....No Help Text
ID.....{A63F6E6E-A7C4-0928-81FB-FE2F70F713C6}
TERM.....PACED: VENTRICULAR

Page: 803

DATA VALUE: PACER

ABBREVIATION.....PACER
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Pacer Catheter
HELP_TEXT.....No Help Text
ID.....{79499460-ACAB-33B3-52CF-D2D2E5BCD0EE}
TERM.....PACER

DATA VALUE: PAIN

ABBREVIATION.....PAIN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Pain
HELP_TEXT.....No Help Text
ID.....{F2F23EE7-0A06-46B8-A496-B7A4651A551F}
TERM.....PAIN

DATA VALUE: PAIN CONTROLLED PER PT

ABBREVIATION.....PAIN CON
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Pain Controlled Per Patient
HELP_TEXT.....No Help Text
ID.....{7EC1AF38-39B1-475E-B58C-E33404A47E96}
TERM.....PAIN CONTROLLED PER PT

DATA VALUE: PAINAD BL 0

ABBREVIATION.....PAINAD BL0
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....0 = Relaxed
HELP_TEXT.....No Help Text
ID.....{C77929DB-A59A-4A24-B57B-BF111770C2AB}
TERM.....PAINAD BL 0

DATA VALUE: PAINAD BL 1

ABBREVIATION.....PAINAD BL1
ACTIVE.....Yes
DESCRIPTION.....1 = Tense, Distressed Pacing, Fidgeting
DISPLAY_NAME.....1 = Tense, Distressed Pacing, Fidgeting
HELP_TEXT.....1 = Tense, Distressed Pacing, Fidgeting
ID.....{CC1E1516-086D-4B79-8264-0C3BD67D8552}
TERM.....PAINAD BL 1

Page: 804

DATA VALUE: PAINAD BL 2

ABBREVIATION.....PAINAD BL2
ACTIVE.....Yes
DESCRIPTION.....2 = Rigid, Fists Clenched, Knees
Pulled Up, Pulling or Pushing Away,
Striking Out
DISPLAY_NAME.....2 = Rigid, Fists Clenched, Knees
Pulled Up, Pulling or Pushing
HELP_TEXT.....2 = Rigid, Fists Clenched, Knees
Pulled Up, Pulling or Pushing Away,
Striking Out
ID.....{34C3E304-DF66-40A7-906F-22AEA3F76FAD}
TERM.....PAINAD BL 2

DATA VALUE: PAINAD BR 0

ABBREVIATION.....PAINAD B0
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....0 = Normal
HELP_TEXT.....No Help Text
ID.....{58E71B2E-F5C3-419F-8B71-A9FD701E0202}
TERM.....PAINAD BR 0

DATA VALUE: PAINAD BR 1

ABBREVIATION.....PAINAD BR1
ACTIVE.....Yes
DESCRIPTION.....1 = Occasional Labored Breathing.
Short Period of Hyperventilation
DISPLAY_NAME.....1 = Occasional Labored Breathing.

Short Period of Hyperventilation
HELP_TEXT.....1 = Occasional Labored Breathing.
Short Period of Hyperventilation
ID.....{727D1B04-A702-4806-857B-78BCBACB0FE1}
TERM.....PAINAD BR 1

DATA VALUE: PAINAD BR 2

ABBREVIATION.....PAINAD BR2
ACTIVE.....Yes
DESCRIPTION.....2 = Noisy Labored Breathing.
Cheynes-Stokes Respirations
DISPLAY_NAME.....2 = Noisy Labored Breathing.
Cheynes-Stokes Respirations
HELP_TEXT.....2 = Noisy Labored Breathing.
Cheynes-Stokes Respirations
ID.....{37978012-ADAC-416C-8920-48CDAAA5D59D}
TERM.....PAINAD BR 2

DATA VALUE: PAINAD CO 0

ABBREVIATION.....PAINAD CO0
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....0 = No Need To Console
HELP_TEXT.....No Help Text
ID.....{06A89B55-5822-42F7-A689-DBD91B62979F}
TERM.....PAINAD CO 0

Page: 805

DATA VALUE: PAINAD CO 1

ABBREVIATION.....PAINAD CO1
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....1 = Distracted or Reassured By Voice
or Touch
HELP_TEXT.....No Help Text
ID.....{8E5DDB47-1955-40AD-A008-A9CE4646A707}
TERM.....PAINAD CO 1

DATA VALUE: PAINAD CO 2

ABBREVIATION.....PAINAD CO2
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....2 = Unable To Console, Distract or
Reassure
HELP_TEXT.....No Help Text
ID.....{1DFFB735-6272-4B60-8F42-CCCAA701C7A7}
TERM.....PAINAD CO 2

DATA VALUE: PAINAD FE 0

ABBREVIATION.....PAINAD FE0
ACTIVE.....Yes
DESCRIPTION.....0 = Smiling or Inexpressive
DISPLAY_NAME.....0 = Smiling or Inexpressive
HELP_TEXT.....0 = Smiling or Inexpressive
ID.....{B3DE1BC6-E7A0-48C5-8C8E-FBE033AD1046}
TERM.....PAINAD FE 0

DATA VALUE: PAINAD FE 1

ABBREVIATION.....PAINAD FE1
ACTIVE.....Yes
DESCRIPTION.....1 = Sad, Frightened, Frown
DISPLAY_NAME.....1 = Sad, Frightened, Frown
HELP_TEXT.....1 = Sad, Frightened, Frown
ID.....{66C01EE9-8F36-40B5-AED2-DBC2E94F860D}
TERM.....PAINAD FE 1

DATA VALUE: PAINAD FE 2

ABBREVIATION.....PAINAD FE2
ACTIVE.....Yes
DESCRIPTION.....2 = Facial Grimacing
DISPLAY_NAME.....2 = Facial Grimacing
HELP_TEXT.....2 = Facial Grimacing
ID.....{35688D06-94F2-4B36-8940-529803D4BC43}
TERM.....PAINAD FE 2

Page: 806

DATA VALUE: PAINAD NV 0

ABBREVIATION.....PAINAD NV0
ACTIVE.....Yes
DESCRIPTION.....0 = None
DISPLAY_NAME.....0 = None
HELP_TEXT.....0 = None
ID.....{63F7F990-2DEA-4B35-AC52-D23B8E13B1FC}
TERM.....PAINAD NV 0

DATA VALUE: PAINAD NV 1

ABBREVIATION.....PAINAD NV1
ACTIVE.....Yes
DESCRIPTION.....1 = Occasional Moan or Groan. Low
Level Speech With A Negative or
Disapproving Quality
DISPLAY_NAME.....1 = Occasional Moan or Groan
HELP_TEXT.....1 = Occasional Moan or Groan. Low
Level Speech With A Negative or
Disapproving Quality
ID.....{D517E6C9-9D98-4D0F-A8F4-261BED4812D1}
TERM.....PAINAD NV 1

DATA VALUE: PAINAD NV 2

ABBREVIATION.....PAINAD NV2
ACTIVE.....Yes
DESCRIPTION.....2 = Repeated, Troubled, Calling Out,
Loud Moaning, Crying
DISPLAY_NAME.....2 = Repeated, Troubled, Calling Out,
Loud Moaning.
HELP_TEXT.....2 = Repeated, Troubled, Calling Out,
Loud Moaning, Crying
ID.....{3C8D6B1E-53E6-4D2C-BD8A-9C548466E1D4}
TERM.....PAINAD NV 2

DATA VALUE: PAINFUL EXPRESSION

ABBREVIATION.....PAINFUL EX
ACTIVE.....Yes
DESCRIPTION.....Painful Expression
DISPLAY_NAME.....Painful Expression
HELP_TEXT.....No Help Text
ID.....{608554F9-29F9-4B27-BE63-486AF81810BD}
TERM.....PAINFUL EXPRESSION

DATA VALUE: PALE

ABBREVIATION.....p
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Pale
HELP_TEXT.....No Help Text
ID.....{0EB39D28-7ABE-A4AB-D334-3F795E65618D}
TERM.....PALE

Page: 807

DATA VALUE: PAROXYSMAL ATRIAL TACHYCARDIA

ABBREVIATION.....pat
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Paroxysmal Atrial Tachycardia
HELP_TEXT.....No Help Text
ID.....{27BC1FAF-A403-6EA8-044C-F1271F12A62A}
TERM.....PAROXYSMAL ATRIAL TACHYCARDIA

DATA VALUE: PART THICKNESS SKIN LOSS INVOLVE EPIDERMIS/DERMIS

ABBREVIATION.....PART EPIDE
ACTIVE.....Yes
DESCRIPTION.....Partial Thickness Skin Loss Involving
Epidermis/Dermis
DISPLAY_NAME.....Partial Thickness Skin Loss Involving
Epidermis/Dermis
HELP_TEXT.....No Help Text
ID.....{EAD3A260-CC46-4B47-AA49-91CFE13EEE4B}
TERM.....PART THICKNESS SKIN LOSS INVOLVE

EPIDERMIS/DERMIS

DATA VALUE: PARTIAL BATH

ABBREVIATION.....PART. BATH
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Partial Bath
HELP_TEXT.....No Help Text
ID.....{60B9C7AD-9E6E-4412-9BB3-57AFF956C1FC}
TERM.....PARTIAL BATH

DATA VALUE: PARTIAL IMPAIRED

ABBREVIATION.....PART IMPAI
ACTIVE.....Yes
DESCRIPTION.....Partial Impaired
DISPLAY_NAME.....Partial Impaired
HELP_TEXT.....Partial Impaired
ID.....{F1B8C52-1C94-4EA7-AB82-117694FD4D52}
TERM.....PARTIAL IMPAIRED

DATA VALUE: PARTIAL SIMPLE

ABBREVIATION.....PART SIMP
ACTIVE.....Yes
DESCRIPTION.....Partial Simple
DISPLAY_NAME.....Partial Simple
HELP_TEXT.....Partial Simple
ID.....{F4A5E7ED-92AB-4021-A63C-E1D948EFE6C1}
TERM.....PARTIAL SIMPLE

Page: 808

DATA VALUE: PASSIVE RANGE OF MOTION

ABBREVIATION.....PASS ROM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Passive Range Of Motion
HELP_TEXT.....No Help Text
ID.....{DD43B9EF-E18F-4203-8564-89EBBAD963A4}
TERM.....PASSIVE RANGE OF MOTION

DATA VALUE: PATENT CLEAN

ABBREVIATION.....PAT CLEAN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Patent Clean
HELP_TEXT.....No Help Text
ID.....{68688246-AE58-2EFB-BF9D-5E00D8CFDA8E}
TERM.....PATENT CLEAN

DATA VALUE: PATIENT

ABBREVIATION.....PATIENT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Patient
HELP_TEXT.....No Help Text
ID.....{7881F205-BCD1-4948-AE66-AD572F5A02F2}
TERM.....PATIENT

DATA VALUE: PATIENT INFORMED

ABBREVIATION.....PAT INFORM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Patient Informed
HELP_TEXT.....No Help Text
ID.....{074A1D19-3920-49F1-8958-2EB6F9CF05BA}
TERM.....PATIENT INFORMED

DATA VALUE: PATIENT REPORTS NOT PRESENT

ABBREVIATION.....PT RP N PR
ACTIVE.....Yes
DESCRIPTION.....Patient Reports Not Present
DISPLAY_NAME.....Patient Reports Not Present
HELP_TEXT.....Patient Reports Not Present
ID.....{AD0EF44D-7231-4E8C-A3AE-0E43DD0C44F2}
TERM.....PATIENT REPORTS NOT PRESENT

Page: 809

DATA VALUE: PATIENT REPORTS PRESENT

ABBREVIATION.....PT REP PRE
ACTIVE.....Yes
DESCRIPTION.....Patient Reports Present
DISPLAY_NAME.....Patient Reports Present
HELP_TEXT.....Patient Reports Present
ID.....{1726BC7E-5A82-41AC-B77A-BF62263E1EF6}
TERM.....PATIENT REPORTS PRESENT

DATA VALUE: PATIENT SIGNS SAFETY CONTRACT

ABBREVIATION.....PT CONTRAC
ACTIVE.....Yes
DESCRIPTION.....Patient Signs Safety Contract
DISPLAY_NAME.....Patient Signs Safety Contract
HELP_TEXT.....No Help Text
ID.....{C65734EC-2CAF-4547-8FD8-CBCD9C706003}
TERM.....PATIENT SIGNS SAFETY CONTRACT

DATA VALUE: PATIENT STAYS WITHIN LINE OF SITE OF NURSING STATION

ABBREVIATION.....pswlosons
ACTIVE.....Yes

DESCRIPTION.....No Description
DISPLAY_NAME.....Patient Stays Within Line Of Site Of
Nursing Station
HELP_TEXT.....No Help Text
ID.....{CE3E3531-7D9B-A2B4-00B3-00EBCDD09D3D}
TERM.....PATIENT STAYS WITHIN LINE OF SITE OF
NURSING STATION

DATA VALUE: PATIENT'S ROOM

ABBREVIATION.....pr
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Patient's Room
HELP_TEXT.....No Help Text
ID.....{0F13E11D-9F78-14A3-0568-AF26393F39F1}
TERM.....PATIENT'S ROOM

DATA VALUE: PCA

ABBREVIATION.....PCA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Patient Controlled Analgesia
HELP_TEXT.....No Help-Text
ID.....{A1EE27C6-F22A-4109-A8E6-28DCCE599E64}
TERM.....PCA

Page: 810

DATA VALUE: PCA AND BASAL RATE

ABBREVIATION.....PCA/BASAL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....PCA and Basal Rate
HELP_TEXT.....No Help Text
ID.....{D711BB51-0662-4DD0-8654-C737F6A6686A}
TERM.....PCA AND BASAL RATE

DATA VALUE: PCA ONLY

ABBREVIATION.....PCA ONLY
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....PCA Only
HELP_TEXT.....No Help Text
ID.....{30E0441A-1098-4FE9-B277-BF2BBD6B1667}
TERM.....PCA ONLY

DATA VALUE: PEDIATRIC

ABBREVIATION.....PEDIATRIC
ACTIVE.....Yes
DESCRIPTION.....No Description

DISPLAY_NAME.....Pediatric
HELP_TEXT.....No Help Text
ID.....{06865E22-3ED9-4EDD-9749-6F8641BF29DE}
TERM.....PEDIATRIC

DATA VALUE: PEG TUBE

ABBREVIATION.....PEG TUBE
ACTIVE.....Yes
DESCRIPTION.....PEG Tube EF Tube
DISPLAY_NAME.....PEG Tube
HELP_TEXT.....No Help Text
ID.....{970C3755-2B3F-462F-AA87-EA33F846843F}
TERM.....PEG TUBE

DATA VALUE: PENETRATING

ABBREVIATION.....PENETRATIN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Penetrating
HELP_TEXT.....No Help Text
ID.....{00A2D4C8-EDDE-4593-A1F8-A87A836B3379}
TERM.....PENETRATING

Page: 811

DATA VALUE: PENROSE DRAIN

ABBREVIATION.....Penrose D
ACTIVE.....Yes
DESCRIPTION.....Penrose Drain
DISPLAY_NAME.....Penrose Drain
HELP_TEXT.....No Help Text
ID.....{DE508E95-D044-ACCA-EE69-33B9A8454E6C}
TERM.....PENROSE DRAIN

DATA VALUE: PERCUSSION BED

ABBREVIATION.....PERC BED
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Percussion Bed
HELP_TEXT.....No Help Text
ID.....{BE0B086B-A2F3-446A-ADE0-9589E92689FE}
TERM.....PERCUSSION BED

DATA VALUE: PERINEAL CARE

ABBREVIATION.....PERIN CARE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Perineal Care
HELP_TEXT.....No Help Text
ID.....{D88736E4-E0EB-4FF7-9199-288A4D051A20}

TERM.....PERINEAL CARE

DATA VALUE: PERMANENT

ABBREVIATION.....p
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Permanent
HELP_TEXT.....No Help Text
ID.....{A31FF319-F8C3-7AAC-D7FC-765C89732D09}
TERM.....PERMANENT

DATA VALUE: PERSONAL MEDICAL ALERT BRACELET

ABBREVIATION.....pmab
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Personal Medical Alert Bracelet
HELP_TEXT.....No Help Text
ID.....{E63D5FDE-4509-5751-A11C-C023E296858A}
TERM.....PERSONAL MEDICAL ALERT BRACELET

Page: 812

DATA VALUE: PETROLATUM

ABBREVIATION.....PETROLATUM
ACTIVE.....Yes
DESCRIPTION.....Petrolatum
DISPLAY_NAME.....Petrolatum
HELP_TEXT.....No Help Text
ID.....{DA297C75-4776-4C8D-B0D7-A0FF0E4D1B16}
TERM.....PETROLATUM

DATA VALUE: PHYSICIAN

ABBREVIATION.....PHYSICIA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Physician
HELP_TEXT.....No Help Text
ID.....{DB0C5AE0-459A-4E04-A202-D1CF54B84A6A}
TERM.....PHYSICIAN

DATA VALUE: PHYSICIAN NOTIFIED

ABBREVIATION.....PHYS NOTIF
ACTIVE.....Yes
DESCRIPTION.....Physician Notified
DISPLAY_NAME.....Physician Notified
HELP_TEXT.....No Help Text
ID.....{4BD68949-74F6-42A6-AEF5-E9A33DC31FFC}
TERM.....PHYSICIAN NOTIFIED

DATA VALUE: PHYSICIAN ORDER

ABBREVIATION.....PHY ORDER
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Physician Order
HELP_TEXT.....No Help Text
ID.....{AB89E640-83B7-4E6C-8C8D-2FA8C2564CCB}
TERM.....PHYSICIAN ORDER

DATA VALUE: PICC

ABBREVIATION.....PICC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....PICC
HELP_TEXT.....Peripherally inserted central catheter
ID.....{3CFE33D0-5A35-88CC-152C-D7501D1858D9}
TERM.....PICC

Page: 813

DATA VALUE: PIN

ABBREVIATION.....PIN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Pin
HELP_TEXT.....No Help Text
ID.....{10D48A4E-46F3-4242-BA56-B9D4B37A4C93}
TERM.....PIN

DATA VALUE: PIN TRACTION

ABBREVIATION.....PIN TRACT
ACTIVE.....Yes
DESCRIPTION.....Pin Traction
DISPLAY_NAME.....Pin Traction
HELP_TEXT.....Pin Traction
ID.....{C75104A5-3B6B-4CFE-8E2E-05877522C811}
TERM.....PIN TRACTION

DATA VALUE: PINK

ABBREVIATION.....PINK
ACTIVE.....Yes
DESCRIPTION.....The color pink
DISPLAY_NAME.....Pink
HELP_TEXT.....The color pink
ID.....{1A7FAF84-D69C-4EEB-7BC5-D21A598A1F4D}
TERM.....PINK

DATA VALUE: PINPOINT

ABBREVIATION.....pinpoint
ACTIVE.....Yes

DESCRIPTION.....Pinpoint
DISPLAY_NAME.....Pinpoint
HELP_TEXT.....No Help Text
ID.....{4E53DDAC-0BF4-66D5-D614-9763BBF1FC8E}
TERM.....PINPOINT

DATA VALUE: PITTED

ABBREVIATION.....PITTED
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Pitted
HELP_TEXT.....No Help Text
ID.....{BEA42FD8-91C7-4BF6-936E-726D78064F11}
TERM.....PITTED

Page: 814

DATA VALUE: PITTING

ABBREVIATION.....PITTING
ACTIVE.....Yes
DESCRIPTION.....PITTING
DISPLAY_NAME.....Pitting
HELP_TEXT.....PITTING
ID.....{A9F6ACFF-E427-145A-A168-345ACB87C1DE}
TERM.....PITTING

DATA VALUE: PJC

ABBREVIATION.....p
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....PJC
HELP_TEXT.....No Help Text
ID.....{1D2D1FEA-2006-3C6C-FAD4-72384F2033C9}
TERM.....PJC

DATA VALUE: PLAIN

ABBREVIATION.....PLAIN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Plain Catheter
HELP_TEXT.....No Help Text
ID.....{483EE48F-F72B-A917-DCD0-E8B69A3FA36F}
TERM.....PLAIN

DATA VALUE: POOR

ABBREVIATION.....POOR
ACTIVE.....Yes
DESCRIPTION.....Poor
DISPLAY_NAME.....Poor
HELP_TEXT.....No Help Text

ID.....{A3D75A12-A791-446D-AF63-9435054DFC4C}
TERM.....POOR

DATA VALUE: PORCINE SUBMUCOSA

ABBREVIATION.....PORC SUBMU
ACTIVE.....Yes
DESCRIPTION.....Oasis (porcine submucosa)
DISPLAY_NAME.....Oasis (porcine submucosa)
HELP_TEXT.....No Help Text
ID.....{6F3917DA-74BD-4790-9D4D-B284E27B7623}
TERM.....PORCINE SUBMUCOSA

Page: 815

DATA VALUE: POSITIONAL

ABBREVIATION.....POSITIONAL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Positional
HELP_TEXT.....No Help Text
ID.....{9A74D16B-5E95-4EA8-88AF-E53C80311A5D}
TERM.....POSITIONAL

DATA VALUE: POST CATH DIET

ABBREVIATION.....POST CA DT
ACTIVE.....Yes
DESCRIPTION.....Post Cath Diet
DISPLAY_NAME.....Post Cath Diet
HELP_TEXT.....No Help Text
ID.....{B18FEC50-18E1-4EED-8B37-E7DBAA7D4D3B}
TERM.....POST CATH DIET

DATA VALUE: POST OPERATION DIAGNOSIS

ABBREVIATION.....PO OP DIAG
ACTIVE.....Yes
DESCRIPTION.....Post Operation Diagnosis
DISPLAY_NAME.....Post Operation Diagnosis
HELP_TEXT.....Post Operation Diagnosis
ID.....{A69379C9-C54A-4C55-85B3-C3D92D1705A0}
TERM.....POST OPERATION DIAGNOSIS

DATA VALUE: POST PROCEDURE DIAGNOSIS

ABBREVIATION.....PO PRO DIA
ACTIVE.....Yes
DESCRIPTION.....Post Procedure Diagnosis
DISPLAY_NAME.....Post Procedure Diagnosis
HELP_TEXT.....Post Procedure Diagnosis
ID.....{933C7D2C-88F8-4C18-ACBB-A16A295D1B3A}
TERM.....POST PROCEDURE DIAGNOSIS

DATA VALUE: POST-FALL

ABBREVIATION.....POST-FALL
ACTIVE.....Yes
DESCRIPTION.....Post-Fall
DISPLAY_NAME.....Post-Fall
HELP_TEXT.....No Help Text
ID.....{9BC82B93-CA4A-44E3-A3BA-7AF72BAFD44B}
TERM.....POST-FALL

Page: 816

DATA VALUE: POST-OPERATIVE

ABBREVIATION.....POST-OP
ACTIVE.....Yes
DESCRIPTION.....Post-Operative
DISPLAY_NAME.....Post-Operative
HELP_TEXT.....No Help Text
ID.....{9FA19656-849F-4A65-92FC-755A4C9E8056}
TERM.....POST-OPERATIVE

DATA VALUE: POST-PROCEDURE

ABBREVIATION.....POST-PROC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Post-Procedure
HELP_TEXT.....No Help Text
ID.....{0F943BD7-BD13-412E-B4EA-9B5526587340}
TERM.....POST-PROCEDURE

DATA VALUE: POST-SEDATION

ABBREVIATION.....POST SED
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Post-Sedation
HELP_TEXT.....No Help Text
ID.....{B829A249-CB14-4919-86F8-4A8C8D6B7753}
TERM.....POST-SEDATION

DATA VALUE: POTASSIUM MODIFIED 60 MEQ

ABBREVIATION.....K MOD 60M
ACTIVE.....Yes
DESCRIPTION.....Potassium Modified 60 MEQ
DISPLAY_NAME.....Potassium Modified 60 MEQ
HELP_TEXT.....No Help Text
ID.....{6EC2F10D-CC44-4872-8473-2445A835ABB0}
TERM.....POTASSIUM MODIFIED 60 MEQ

DATA VALUE: PRBC

ABBREVIATION.....PRBC

ACTIVE.....Yes
DESCRIPTION.....Packed Red Blood Cells
DISPLAY_NAME.....PRBC
HELP_TEXT.....No Help Text
ID.....{F6D9B889-1BCB-4724-B75A-E98743585D55}
TERM.....PRBC

Page: 817

DATA VALUE: PRE-MEDICATED FOR ACTIVITY

ABBREVIATION.....PRE-MED AC
ACTIVE.....Yes
DESCRIPTION.....Pre-Medicated For Activity
DISPLAY_NAME.....Pre-Medicated For Activity
HELP_TEXT.....No Help Text
ID.....{D6DA772A-FA22-4170-B928-855A2F68711F}
TERM.....PRE-MEDICATED FOR ACTIVITY

DATA VALUE: PRE-OPERATIVE

ABBREVIATION.....PRE-OP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Pre-Operative
HELP_TEXT.....No Help Text
ID.....{80C418D5-54B1-4CC3-9B19-5653B8D73C45}
TERM.....PRE-OPERATIVE

DATA VALUE: PRE-PROCEDURE

ABBREVIATION.....PRE-PROC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Pre-Procedure
HELP_TEXT.....No Help Text
ID.....{D2A7B025-9AB9-4C69-B163-D5EBB78B82BE}
TERM.....PRE-PROCEDURE

DATA VALUE: PRE-SEDATION

ABBREVIATION.....PRE SEDA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Pre-Sedation
HELP_TEXT.....No Help Text
ID.....{684C764E-2BE9-4C99-A2B2-C92D75C73751}
TERM.....PRE-SEDATION

DATA VALUE: PRE-SURGERY

ABBREVIATION.....PRE-SURG
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Pre-Surgery

HELP_TEXT.....No Help Text
ID.....{0C71761F-3024-4638-A95C-FDF12DC03CF5}
TERM.....PRE-SURGERY

Page: 818

DATA VALUE: PREMATURE ATRIAL

ABBREVIATION.....pa
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Premature Atrial
HELP_TEXT.....No Help Text
ID.....{070FDCB1-2EC4-B4EF-4DFD-81F6CCFFBDED}
TERM.....PREMATURE ATRIAL

DATA VALUE: PRESENT

ABBREVIATION.....p
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Present
HELP_TEXT.....No Help Text
ID.....{384DC20F-1625-CE1E-1C85-A3B3924BB2A0}
TERM.....PRESENT

DATA VALUE: PRESSURE

ABBREVIATION.....PRESSURE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Pressure
HELP_TEXT.....No Help Text
ID.....{387AFB7F-3973-4D77-B2E6-0724F344D8F6}
TERM.....PRESSURE

DATA VALUE: PRESSURE CALIBRATION

ABBREVIATION.....PRESS CALI
ACTIVE.....Yes
DESCRIPTION.....Pressure Calibration
DISPLAY_NAME.....Pressure Calibration
HELP_TEXT.....No Help Text
ID.....{38C12E69-2CD2-4F61-98D9-06AB857DA51A}
TERM.....PRESSURE CALIBRATION

DATA VALUE: PRESSURE CONTROL

ABBREVIATION.....PRESS CONT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Pressure Control
HELP_TEXT.....No Help Text
ID.....{C3CDF756-484F-1483-9883-9B63DEF3B012}
TERM.....PRESSURE CONTROL

DATA VALUE: PRESSURE DRESSING

ABBREVIATION.....PRESS DRES
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Pressure Dressing
HELP_TEXT.....No Help Text
ID.....{6D9DB9EF-85EB-8321-03FC-879AA4C5615B}
TERM.....PRESSURE DRESSING

DATA VALUE: PRESSURE ULCER

ABBREVIATION.....PRE ULCER
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Pressure Ulcer
HELP_TEXT.....No Help Text
ID.....{7279639E-498D-4EA3-A3FA-4A6D611E6F48}
TERM.....PRESSURE ULCER

DATA VALUE: PRESYSTOLIC

ABBREVIATION.....PRESYSTOLI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Presystolic
HELP_TEXT.....No Help Text
ID.....{BB85B8EE-36F8-4AC2-9AB0-0BABCFA9FCBB}
TERM.....PRESYSTOLIC

DATA VALUE: PREVENT ACCESS TO DEVICE

ABBREVIATION.....PREV ACCE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Prevent Access To Medical Devices
HELP_TEXT.....No Help Text
ID.....{E7455083-DD2C-4562-A6E2-9626A4CDC4DD}
TERM.....PREVENT ACCESS TO DEVICE

DATA VALUE: PREVIOUS PERCUTANEOUS SITE

ABBREVIATION.....PRE PER SI
ACTIVE.....Yes
DESCRIPTION.....Previous Percutaneous Site
DISPLAY_NAME.....Previous Percutaneous Site
HELP_TEXT.....No Help Text
ID.....{F2ECCF5C-D91A-40F5-9C09-452C2C81DA66}
TERM.....PREVIOUS PERCUTANEOUS SITE

DATA VALUE: PRN

ABBREVIATION.....PRN
ACTIVE.....Yes
DESCRIPTION.....As Needed
DISPLAY_NAME.....PRN
HELP_TEXT.....As Needed
ID.....{36493D81-954B-4C56-B0A2-9988EA785EBD}
TERM.....PRN

DATA VALUE: PROCEDURE ROOM

ABBREVIATION.....PROC ROO
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Procedure Room
HELP_TEXT.....No Help Text
ID.....{CC9A2AB0-7A73-4DE7-9B7F-DC2C160AE462}
TERM.....PROCEDURE ROOM

DATA VALUE: PRODUCTIVE

ABBREVIATION.....p
ACTIVE.....Yes
DESCRIPTION.....Produces desired result
DISPLAY_NAME.....Productive
HELP_TEXT.....No Help Text
ID.....{22BE0499-F9FE-3A27-2605-61EBEC0ED6AB}
TERM.....PRODUCTIVE

DATA VALUE: PROJECTILE VOMITING

ABBREVIATION.....PROJ VOMIT
ACTIVE.....Yes
DESCRIPTION.....Projectile Vomiting
DISPLAY_NAME.....Projectile Vomiting
HELP_TEXT.....Projectile Vomiting
ID.....{0D2CB23B-B266-40E1-A6FF-20B733655B69}
TERM.....PROJECTILE VOMITING

DATA VALUE: PRONE

ABBREVIATION.....PRONE
ACTIVE.....Yes
DESCRIPTION.....Prone
DISPLAY_NAME.....Prone
HELP_TEXT.....No Help Text
ID.....{41180E09-D9D6-4903-BAE9-28CD5AAC4991}
TERM.....PRONE

DATA VALUE: PROSTHESIS

ABBREVIATION.....PROSTHESIS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Prosthesis
HELP_TEXT.....No Help Text
ID.....{CCD0EF23-FE9C-4B20-942E-9CF1CD7CF6E9}
TERM.....PROSTHESIS

DATA VALUE: PROXIMAL PORT

ABBREVIATION.....PROX PORT
ACTIVE.....Yes
DESCRIPTION.....Proximal Port
DISPLAY_NAME.....Proximal Port
HELP_TEXT.....No Help Text
ID.....{3244EF4F-F18A-4250-96F5-ABCEA7BFA578}
TERM.....PROXIMAL PORT

DATA VALUE: PT DENIES NAUSEA / NO VOMITING

ABBREVIATION.....DEN NAUS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Pt. Denies Nausea / No Vomiting
HELP_TEXT.....No Help Text
ID.....{4374D52E-9D4E-41AD-8B52-FCBDF6B5A1E9}
TERM.....PT DENIES NAUSEA / NO VOMITING

DATA VALUE: PT PAIN AT ACCEPTABLE LEVEL

ABBREVIATION.....PT PA OK
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Patient's Pain at Acceptable Level
HELP_TEXT.....No Help Text
ID.....{78489498-4A5B-41C5-98C5-95D56C4A0933}
TERM.....PT PAIN AT ACCEPTABLE LEVEL

DATA VALUE: PUNCTURE

ABBREVIATION.....PUNCTURE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Puncture
HELP_TEXT.....No Help Text
ID.....{28CF34AF-BD37-4BA4-A59B-C4AFCE2EB8DD}
TERM.....PUNCTURE

Page: 822

DATA VALUE: PUREED

ABBREVIATION.....PUREED
ACTIVE.....Yes
DESCRIPTION.....Pureed

DISPLAY_NAME.....Pureed
HELP_TEXT.....No Help Text
ID.....{99F8AE5F-E390-4583-8A92-CAC13E21BCF4}
TERM.....PUREED

DATA VALUE: PURULENT

ABBREVIATION.....PURULENT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Purulent
HELP_TEXT.....No Help Text
ID.....{86C01494-743A-C6A2-8883-01CA5200C79B}
TERM.....PURULENT

DATA VALUE: PVC

ABBREVIATION.....p
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....PVC
HELP_TEXT.....No Help Text
ID.....{33FF92F1-84CC-7CFA-41CF-B271A26C59B8}
TERM.....PVC

DATA VALUE: QUIET OR RESTING

ABBREVIATION.....qor
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Quiet Or Resting
HELP_TEXT.....No Help Text
ID.....{76E4908F-014C-2088-FF7F-33EBF223DE40}
TERM.....QUIET OR RESTING

DATA VALUE: QUIET ROOM

ABBREVIATION.....qr
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Quiet Room
HELP_TEXT.....No Help Text
ID.....{CCD96D41-C273-F6F3-DA18-8818B1CC9980}
TERM.....QUIET ROOM

Page: 823

DATA VALUE: R BUNDLE BRANCH BLOCK

ABBREVIATION.....rbbb
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Bundle Branch Block
HELP_TEXT.....No Help Text
ID.....{FDE0D174-5220-F069-D284-20399D9489B0}

TERM.....R BUNDLE BRANCH BLOCK

DATA VALUE: R FOOT DOWN

ABBREVIATION.....R FOOT DOW
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Foot Down
HELP_TEXT.....No Help Text
ID.....{DBCC0BA2-3097-42C7-AE9E-DF71FA779162}
TERM.....R FOOT DOWN

DATA VALUE: R FOOT UP

ABBREVIATION.....R FOOT UP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Foot Up
HELP_TEXT.....No Help Text
ID.....{5928B710-6EB1-4183-93CB-6C72B37A97C3}
TERM.....R FOOT UP

DATA VALUE: R HEAD DOWN

ABBREVIATION.....R HEAD DOW
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Head Down
HELP_TEXT.....No Help Text
ID.....{DFDE0638-444C-4023-B586-985E84EEA1FA}
TERM.....R HEAD DOWN

DATA VALUE: R HEAD UP

ABBREVIATION.....R HEAD UP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Head Up
HELP_TEXT.....No Help Text
ID.....{D03DBEC6-F4DD-441A-AAB2-653CC34FCDBF}
TERM.....R HEAD UP

Page: 824

DATA VALUE: R NECK

ABBREVIATION.....RIGHT NECK
ACTIVE.....Yes
DESCRIPTION.....RIGHT NECK
DISPLAY_NAME.....R Neck
HELP_TEXT.....RIGHT NECK
ID.....{907E90CD-87BF-073D-AF76-EA11ABE7FC3C}
TERM.....R NECK

DATA VALUE: R NEPHROSTOMY

ABBREVIATION.....R NEPHROST
ACTIVE.....Yes
DESCRIPTION.....Right Nephrostomy
DISPLAY_NAME.....R Nephrostomy
HELP_TEXT.....Right Nephrostomy
ID.....{33661137-688D-5611-4368-B1DFD8DB1469}
TERM.....R NEPHROSTOMY

DATA VALUE: R STENT

ABBREVIATION.....RT. STENT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Stent
HELP_TEXT.....No Help Text
ID.....{B6D0EED3B-1F08-11D2-165D-97A9F58327D5}
TERM.....R STENT

DATA VALUE: RADIATING

ABBREVIATION.....RADIATING
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....RADIATING
HELP_TEXT.....No Help Text
ID.....{B6807939-F809-45C1-912A-ADEF2E0DB0DC}
TERM.....RADIATING

DATA VALUE: RAMSEY ANXIOUS

ABBREVIATION.....RAM ANXIOU
ACTIVE.....Yes
DESCRIPTION.....Anxious, agitated, or restless
DISPLAY_NAME.....1 - Anxious
HELP_TEXT.....Anxious, agitated, or restless
ID.....{88205A0C-7065-4616-8FCF-75059462DD15}
TERM.....RAMSEY ANXIOUS

Page: 825

DATA VALUE: RAMSEY ASLEEP BRISK RESPONSE

ABBREVIATION.....RAM ASL BR
ACTIVE.....Yes
DESCRIPTION.....Asleep; brisk response to light
glabellar tap or loud auditory
stimulus
DISPLAY_NAME.....4 - Asleep; Brisk Response
HELP_TEXT.....Asleep; brisk response to light
glabellar tap or loud auditory
stimulus
ID.....{00F4A117-F9E7-4B60-A521-F2E8725AEC3}
TERM.....RAMSEY ASLEEP BRISK RESPONSE

DATA VALUE: RAMSEY ASLEEP NO RESPONSE

ABBREVIATION.....RAM ASL NO
ACTIVE.....Yes
DESCRIPTION.....Asleep; no response to painful stimuli
DISPLAY_NAME.....6 - Asleep No Response
HELP_TEXT.....Asleep; no response to painful stimuli
ID.....{3FF5BA0A-2D3C-4935-A432-23DAA4761AB7}
TERM.....RAMSEY ASLEEP NO RESPONSE

DATA VALUE: RAMSEY ASLEEP SLUGGISH

ABBREVIATION.....RAM ASL SL
ACTIVE.....Yes
DESCRIPTION.....Asleep; sluggish response to light
glabellar tap or loud auditory
stimulus
DISPLAY_NAME.....5 - Asleep Sluggish Response
HELP_TEXT.....Asleep; sluggish response to light
glabellar tap or loud auditory
stimulus
ID.....{AA61604F-00F9-4103-B27F-E817FBC9577B}
TERM.....RAMSEY ASLEEP SLUGGISH

DATA VALUE: RAMSEY COOPERATIVE

ABBREVIATION.....RAM COOPER
ACTIVE.....Yes
DESCRIPTION.....Cooperative, oriented, and tranquil
DISPLAY_NAME.....2 - Cooperative
HELP_TEXT.....Cooperative, oriented, and tranquil
ID.....{56846ADA-0F0E-44AA-9779-A567B9D24C4F}
TERM.....RAMSEY COOPERATIVE

DATA VALUE: RAMSEY SEDATED

ABBREVIATION.....RAM SEDATE
ACTIVE.....Yes
DESCRIPTION.....Sedated but responds to commands
DISPLAY_NAME.....3 - Sedated But Responds
HELP_TEXT.....Sedated but responds to commands
ID.....{CC12405F-60DB-485F-94E5-78C1EB32AB74}
TERM.....RAMSEY SEDATED

Page: 826

DATA VALUE: RASH/EXFOLIATION

ABBREVIATION.....RASH/EXFOL
ACTIVE.....Yes
DESCRIPTION.....Rash / Exfoliation
DISPLAY_NAME.....Rash / Exfoliation
HELP_TEXT.....No Help Text
ID.....{AB83ABA9-28A7-6484-67E7-4E66B5526942}
TERM.....RASH/EXFOLIATION

DATA VALUE: RASS AGITATED

ABBREVIATION.....RASS AGITA
ACTIVE.....Yes
DESCRIPTION.....Frequent non-purposeful movement,
 fights ventilator
DISPLAY_NAME.....+2: Agitated
HELP_TEXT.....Frequent non-purposeful movement,
 fights ventilator
ID.....{D81A5D5D-E7FC-445A-A814-495731ABA0F8}
TERM.....RASS AGITATED

DATA VALUE: RASS ALERT AND CALM

ABBREVIATION.....RASS AL CA
ACTIVE.....Yes
DESCRIPTION.....Alert and Calm
DISPLAY_NAME.....0: Alert and Calm
HELP_TEXT.....Alert and Calm
ID.....{F1834BF6-F8F1-4E13-BC25-B13E756E0FF0}
TERM.....RASS ALERT AND CALM

DATA VALUE: RASS COMBATIVE

ABBREVIATION.....RASS COMBA
ACTIVE.....Yes
DESCRIPTION.....Overtly combative, violent, immediate
 danger to staff
DISPLAY_NAME.....+4: Combative
HELP_TEXT.....Overtly combative, violent, immediate
 danger to staff
ID.....{A97E064A-7709-42F9-9A88-4AAFF1020883}
TERM.....RASS COMBATIVE

DATA VALUE: RASS DEEP SEDATION

ABBREVIATION.....RASS DP SE
ACTIVE.....Yes
DESCRIPTION.....No response to voice, but movement or
 eye opening to physical stimulation
DISPLAY_NAME.....-4: Deep Sedation
HELP_TEXT.....No response to voice, but movement or
 eye opening to physical stimulation
ID.....{83CB4E64-3581-452F-987A-D31CB3AB38D2}
TERM.....RASS DEEP SEDATION

Page: 827

DATA VALUE: RASS DROWSY

ABBREVIATION.....RASS DROWS
ACTIVE.....Yes
DESCRIPTION.....Not fully alert, but has sustained
 awakening (eye-opening/eye contact)

to voice (>10 seconds)
DISPLAY_NAME.....-1: Drowsy
HELP_TEXT.....Not fully alert, but has sustained
awakening (eye-opening/eye contact)
to voice (>10 seconds)
ID.....{5E6BDAC6-4F81-4110-A0AF-FED746D9FF0D}
TERM.....RASS DROWSY

DATA VALUE: RASS LIGHT SEDATION

ABBREVIATION.....RASS LT SE
ACTIVE.....Yes
DESCRIPTION.....Briefly awakens with eye contact to
voice (<10 seconds)
DISPLAY_NAME.....-2: Light Sedation
HELP_TEXT.....Briefly awakens with eye contact to
voice (<10 seconds)
ID.....{FDEEC084-3922-4065-A863-4893E42F8F93}
TERM.....RASS LIGHT SEDATION

DATA VALUE: RASS MODERATE SEDATION

ABBREVIATION.....RASS MD SE
ACTIVE.....Yes
DESCRIPTION.....Movement or eye opening to voice (but
no eye contact)
DISPLAY_NAME.....-3: Moderate Sedation
HELP_TEXT.....Movement or eye opening to voice (but
no eye contact)
ID.....{D12B0A3E-B80E-494F-B89F-A79E8B81B080}
TERM.....RASS MODERATE SEDATION

DATA VALUE: RASS RESTLESS

ABBREVIATION.....RASS RESTL
ACTIVE.....Yes
DESCRIPTION.....Anxious but movements not aggressive
vigorous
DISPLAY_NAME.....+1: Restless
HELP_TEXT.....Anxious but movements not aggressive
vigorous
ID.....{2CCFC1ED-D807-45F0-AE69-980E214AD5C3}
TERM.....RASS RESTLESS

DATA VALUE: RASS UNAROUSABLE

ABBREVIATION.....RASS UNARO
ACTIVE.....Yes
DESCRIPTION.....No response to voice or physical
stimulation
DISPLAY_NAME.....-5: Unarousable
HELP_TEXT.....No response to voice or physical
stimulation
ID.....{1578E147-7D2C-4469-828E-FCA97109B27C}
TERM.....RASS UNAROUSABLE

DATA VALUE: RASS VERY AGITATED

ABBREVIATION.....RASS V AGI
ACTIVE.....Yes
DESCRIPTION.....Pulls or removes tube(s) or
 catheter(s); aggressive
DISPLAY_NAME.....+3: Very Agitated
HELP_TEXT.....Pulls or removes tube(s) or
 catheter(s); aggressive
ID.....{7528CA1D-A9A4-4F4E-8943-E69F3B42E189}
TERM.....RASS VERY AGITATED

DATA VALUE: REACTIVE

ABBREVIATION.....Reactive
ACTIVE.....Yes
DESCRIPTION.....Reactive
DISPLAY_NAME.....Reactive
HELP_TEXT.....No Help Text
ID.....{31B44ADB-0000-A5B0-5F56-59D4A08C05D2}
TERM.....REACTIVE

DATA VALUE: RECTAL TUBE

ABBREVIATION.....RECT TUBE
ACTIVE.....Yes
DESCRIPTION.....Rectal Tube
DISPLAY_NAME.....Rectal Tube
HELP_TEXT.....No Help Text
ID.....{7E492AE5-D3B8-44EC-8841-97770FAD3AFE}
TERM.....RECTAL TUBE

DATA VALUE: RED

ABBREVIATION.....RED
ACTIVE.....Yes
DESCRIPTION.....The color red
DISPLAY_NAME.....Red
HELP_TEXT.....The color red
ID.....{515D7CFB-3527-30E7-EC4B-B607C272DFA0}
TERM.....RED

DATA VALUE: RED/GRANULATING

ABBREVIATION.....RED/GRANUL
ACTIVE.....Yes
DESCRIPTION.....Red / Granulating
DISPLAY_NAME.....Red / Granulating
HELP_TEXT.....No Help Text
ID.....{BB69D2A2-7807-490A-B21D-77B6677796F6}
TERM.....RED/GRANULATING

DATA VALUE: REDDENING OF THE SKIN

ABBREVIATION.....RED O SKIN
ACTIVE.....Yes
DESCRIPTION.....Reddening Of The Skin
DISPLAY_NAME.....Reddening Of The Skin
HELP_TEXT.....No Help Text
ID.....{716AEC57-0835-4D07-A78C-D8453599A3D2}
TERM.....REDDENING OF THE SKIN

DATA VALUE: REDIRECT PATIENT ACTIVITY

ABBREVIATION.....REDIR ACT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Redirect Patient Activity
HELP_TEXT.....No Help Text
ID.....{BFF0ABB1-953A-4275-9EF2-08BFF660AA8D}
TERM.....REDIRECT PATIENT ACTIVITY

DATA VALUE: REFLEX

ABBREVIATION.....r
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Reflex
HELP_TEXT.....No Help Text
ID.....{F21585F9-5206-A6B3-2BF6-F416E5842366}
TERM.....REFLEX

DATA VALUE: REFUSE MEDICATION

ABBREVIATION.....REFUSE MED
ACTIVE.....Yes
DESCRIPTION.....Refuse Medication
DISPLAY_NAME.....Refuse Medication
HELP_TEXT.....No Help Text
ID.....{17D6956B-4CF6-492C-B113-B5A93D9E215E}
TERM.....REFUSE MEDICATION

DATA VALUE: REFUSE TO MOVE

ABBREVIATION.....REF TO MOV
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Refuse To Move
HELP_TEXT.....No Help Text
ID.....{E14588CE-2498-4FE8-AB0B-E954DC1696AC}
TERM.....REFUSE TO MOVE

DATA VALUE: REFUSED

ABBREVIATION.....REFUSED
ACTIVE.....Yes
DESCRIPTION.....Refused
DISPLAY_NAME.....Refused
HELP_TEXT.....No Help Text
ID.....{D5B3E135-4D68-4528-BBC9-6E8A7A997104}
TERM.....REFUSED

DATA VALUE: REGISTERED NURSE

ABBREVIATION.....RN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....RN
HELP_TEXT.....No Help Text
ID.....{524F344F-FDCD-423F-AD34-291D5DB4271B}
TERM.....REGISTERED NURSE

DATA VALUE: REGULAR

ABBREVIATION.....REGULAR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Regular
HELP_TEXT.....No Help Text
ID.....{3D3F3D4C-ECD0-2056-A0D1-43CC2F2BA2C8}
TERM.....REGULAR

DATA VALUE: REINFORCED

ABBREVIATION.....REINFORCED
ACTIVE.....Yes
DESCRIPTION.....Dressing has been reinforced
DISPLAY_NAME.....Reinforced
HELP_TEXT.....No Help Text
ID.....{EC524C50-0CDE-4AA1-8059-E4323850575E}
TERM.....REINFORCED

DATA VALUE: RELIGIOUS

ABBREVIATION.....RELIGIOUS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Religious
HELP_TEXT.....No Help Text
ID.....{936C9949-8A30-4BBF-934C-1F194551B8E2}
TERM.....RELIGIOUS

DATA VALUE: REMOVE FROM STIM. SETTING

ABBREVIATION.....REM SET
ACTIVE.....Yes

DESCRIPTION.....No Description
DISPLAY_NAME.....Remove From Stimulating Setting
HELP_TEXT.....No Help Text
ID.....{2010BA89-B399-4A30-9ACF-C65AB57B6351}
TERM.....REMOVE FROM STIM. SETTING

DATA VALUE: RENAL DIET

ABBREVIATION.....REN DIET
ACTIVE.....Yes
DESCRIPTION.....Renal Diet
DISPLAY_NAME.....Renal Diet
HELP_TEXT.....No Help Text
ID.....{089DCCC2-B12B-4E3B-8D12-FCA5E974A55A}
TERM.....RENAL DIET

DATA VALUE: REORIENTED TO ENVIRONMENT

ABBREVIATION.....REORIENT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Reoriented To Environment
HELP_TEXT.....No Help Text
ID.....{4C47B764-4BA7-4D71-9B0B-D42F120E9A4A}
TERM.....REORIENTED TO ENVIRONMENT

DATA VALUE: REPORTED

ABBREVIATION.....REPORTED
ACTIVE.....Yes
DESCRIPTION.....Reported as having happened
DISPLAY_NAME.....Reported
HELP_TEXT.....No Help Text
ID.....{1D13BAC5-D8F9-455C-B815-8E256943492E}
TERM.....REPORTED

DATA VALUE: REPOSITION

ABBREVIATION.....REPOSITION
ACTIVE.....Yes
DESCRIPTION.....Reposition
DISPLAY_NAME.....Reposition
HELP_TEXT.....No Help Text
ID.....{DF3039CA-5977-4454-81D4-DB0AF966F85D}
TERM.....REPOSITION

Page: 832

DATA VALUE: REQUEST FOR MORE MEDICATION

ABBREVIATION.....REQ MOR ME
ACTIVE.....Yes
DESCRIPTION.....Request For More Medication
DISPLAY_NAME.....Request For More Medication
HELP_TEXT.....No Help Text

ID.....{0639081C-BEEC-4CAD-B420-5F2C4DBA38D0}
TERM.....REQUEST FOR MORE MEDICATION

DATA VALUE: RES CRIT: AGGRESSIVE

ABBREVIATION.....R.C. AGGR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Aggressive
HELP_TEXT.....No Help Text
ID.....{7A47B361-85EA-4C88-A1AF-CD8FF7D35D97}
TERM.....RES CRIT: AGGRESSIVE

DATA VALUE: RES CRIT: AMBULATE

ABBREVIATION.....R.C. AMBU
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Attempting To Ambulate When
 Unable/High Injury Risk
HELP_TEXT.....No Help Text
ID.....{54B72608-2C13-46CD-8034-00179669762C}
TERM.....RES CRIT: AMBULATE

DATA VALUE: RES CRIT: COMBATIVE

ABBREVIATION.....R.C. COMB
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Combative
HELP_TEXT.....No Help Text
ID.....{38DBF61A-D364-4215-B4AA-36B660F96462}
TERM.....RES CRIT: COMBATIVE

DATA VALUE: RES CRIT: EQUIPMENT

ABBREVIATION.....R.C. EQUI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Attempting To Remove
 Equipment/Dressings
HELP_TEXT.....No Help Text
ID.....{38C2C56C-6D57-4750-8E68-F1193B16190A}
TERM.....RES CRIT: EQUIPMENT

Page: 833

DATA VALUE: RES CRIT: HARM TO SELF

ABBREVIATION.....R.C. HARM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Risk of Physical Harm To Self
HELP_TEXT.....No Help Text
ID.....{4DC33341-8F53-4D30-BC28-1E32DB54F30D}

TERM.....RES CRIT: HARM TO SELF

DATA VALUE: RES CRIT: LINES

ABBREVIATION.....R.C. LINE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Attempting To Pull Lines
HELP_TEXT.....No Help Text
ID.....{27F545F2-5338-4BC8-B5CE-1D321D28F8D6}
TERM.....RES CRIT: LINES

DATA VALUE: RES CRIT: UNABLE TO FOLLOW DIRECTIONS

ABBREVIATION.....R.C. DIRE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Unable To Follow Directions While
Injury Risk
HELP_TEXT.....No Help Text
ID.....{41358EAE-9912-48F8-8931-D3C635406D21}
TERM.....RES CRIT: UNABLE TO FOLLOW DIRECTIONS

DATA VALUE: RES EVAL: AMBULATION/EXERCISE

ABBREVIATION.....R.E. EXER
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Ambulation / Exercise
HELP_TEXT.....No Help Text
ID.....{46F4D1EB-3439-4BE4-A78C-453BDBF6C230}
TERM.....RES EVAL: AMBULATION/EXERCISE

DATA VALUE: RES EVAL: DECREASE STIMULATION

ABBREVIATION.....R.E. DECR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Decrease Stimulation
HELP_TEXT.....No Description
ID.....{0841626E-994E-4346-B701-1662A00209EB}
TERM.....RES EVAL: DECREASE STIMULATION

Page: 834

DATA VALUE: RES EVAL: DIVERSIONAL ACTIVITIES

ABBREVIATION.....R.E. DIVER
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Diversional Activities
HELP_TEXT.....No Help Text
ID.....{569F6D49-C564-4910-835A-B05179086E1F}
TERM.....RES EVAL: DIVERSIONAL ACTIVITIES

DATA VALUE: RES EVAL: EDUCATION

ABBREVIATION.....R.E. EDUC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Education
HELP_TEXT.....No Help Text
ID.....{76F341E6E-7F0A-429C-B768-9CA83C452FD7}
TERM.....RES EVAL: EDUCATION

DATA VALUE: RES EVAL: INVOLVEMENT WITH FAMILY

ABBREVIATION.....R.E. INVO
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Involvement With Family And Others
HELP_TEXT.....No Help Text
ID.....{0EF31047-70A5-4072-953A-B62E1ACDF84B}
TERM.....RES EVAL: INVOLVEMENT WITH FAMILY

DATA VALUE: RES EVAL: MEDICATION REVIEW

ABBREVIATION.....R.E. REVI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Medication Review
HELP_TEXT.....No Help Text
ID.....{914A0F82-5863-4C84-B3AB-D5A22BD340D8}
TERM.....RES EVAL: MEDICATION REVIEW

DATA VALUE: RES EVAL: REALITY ORIENTATION

ABBREVIATION.....R.E. REALI
ACTIVE.....Yes
DESCRIPTION.....None
DISPLAY_NAME.....Reality Orientation
HELP_TEXT.....None
ID.....{221C363C-2261-40FC-8400-8B1D2E693F69}
TERM.....RES EVAL: REALITY ORIENTATION

Page: 835

DATA VALUE: RES EVAL: RELOCATE

ABBREVIATION.....R.E. RELO
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Relocate Patient
HELP_TEXT.....No Help Text
ID.....{D8165678-ECEF-4F4D-A1C4-C616349E1873}
TERM.....RES EVAL: RELOCATE

DATA VALUE: RES EVAL: UNABLE TO ATTEMPT

ABBREVIATION.....R.E. UNAB

ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Unable To Attempt Due To Patient
Safety Concern
HELP_TEXT.....No Help Text
ID.....{979B0477-5789-4319-9562-6752048C4005}
TERM.....RES EVAL: UNABLE TO ATTEMPT

DATA VALUE: RES LRM: DEVICE ACCESS

ABBREVIATION.....RLRM ACCE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Prevent Access To Medical Devices
HELP_TEXT.....No Help Text
ID.....{1662F5E1-30E7-4E8A-BCC5-1738868435B8}
TERM.....RES LRM: DEVICE ACCESS

DATA VALUE: RES LRM: LIMIT SETTING

ABBREVIATION.....RLRM SETT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Limit Setting
HELP_TEXT.....No Help Text
ID.....{96417DDF-8562-4EFF-B4F6-0933FD91CD36}
TERM.....RES LRM: LIMIT SETTING

DATA VALUE: RES LRM: PRN MEDS

ABBREVIATION.....RLRM: MEDS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Offer PRN Meds
HELP_TEXT.....No Help Text
ID.....{CDF02FA3-FC72-4CCC-8197-41710A69A2CB}
TERM.....RES LRM: PRN MEDS

Page: 836

DATA VALUE: RES LRM: REDIRECT

ABBREVIATION.....RLRM REDI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Redirect Patient Activity
HELP_TEXT.....No Help Text
ID.....{07DE67AF-E226-4794-ADF5-17E9C56C9FDF}
TERM.....RES LRM: REDIRECT

DATA VALUE: RES LRM: REMOVE FROM SETTING

ABBREVIATION.....RLRM REMO
ACTIVE.....Yes
DESCRIPTION.....No Description

DISPLAY_NAME.....Remove From Stimulating Setting
HELP_TEXT.....No Help Text
ID.....{857DBF6E-7EFF-4F07-9AC1-7FDB37CC7618}
TERM.....RES LRM: REMOVE FROM SETTING

DATA VALUE: RES LRM: REORIENT

ABBREVIATION.....RLRM REOR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Reoriented To Environment
HELP_TEXT.....No Help Text
ID.....{DFCC8245-180D-4F3F-AEFC-5D55AFDDC46D}
TERM.....RES LRM: REORIENT

DATA VALUE: RES LRM: SUPERVISION

ABBREVIATION.....RLRM SUPR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....One To One Supervision
HELP_TEXT.....No Help Text
ID.....{18FA8CB0-0155-4D1C-9B14-AE51A0C06280}
TERM.....RES LRM: SUPERVISION

DATA VALUE: RES LRM: VERBAL DE-ESCALATION

ABBREVIATION.....RLRM VERB
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Verbal De-escalation
HELP_TEXT.....No Help Text
ID.....{CE8618C8-CAD7-4C74-8566-ACB5E66B082C}
TERM.....RES LRM: VERBAL DE-ESCALATION

Page: 837

DATA VALUE: RES LRM: VERBAL REMINDERS

ABBREVIATION.....RLRM REMI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Frequent Verbal Reminders
HELP_TEXT.....No Help Text
ID.....{987B5A32-DE2E-47E5-9B06-1AEE86096AB2}
TERM.....RES LRM: VERBAL REMINDERS

DATA VALUE: RESTLESS

ABBREVIATION.....r
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Restless
HELP_TEXT.....No Help Text
ID.....{8FCCAD26-DBA0-7067-DA72-30D342BFA8E8}

TERM.....RESTLESS

DATA VALUE: RESTRAINT

ABBREVIATION.....RESTRAINT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Restraint
HELP_TEXT.....No Help Text
ID.....{0429322D-A987-4913-B64D-A6414645B214}
TERM.....RESTRAINT

DATA VALUE: RESTRAINT BILATERAL ANKLE

ABBREVIATION.....RES BI ANK
ACTIVE.....Yes
DESCRIPTION.....Bilateral Ankle Restraint
DISPLAY_NAME.....Bilateral Ankle
HELP_TEXT.....No Help Text
ID.....{8A0B0DE5-E652-40E7-9B67-2F4D27685A68}
TERM.....RESTRAINT BILATERAL ANKLE

DATA VALUE: RESTRAINT BILATERAL HAND MITT

ABBREVIATION.....RES BI H M
ACTIVE.....Yes
DESCRIPTION.....Bilateral Hand Mitt Restraint
DISPLAY_NAME.....Bilateral Hand Mitt
HELP_TEXT.....No Help Text
ID.....{5E293110-1235-4307-80C2-0FFFA007E1DA}
TERM.....RESTRAINT BILATERAL HAND MITT

Page: 838

DATA VALUE: RESTRAINT GERI CHAIR

ABBREVIATION.....RES GERI C
ACTIVE.....Yes
DESCRIPTION.....Geri Chair Restraint
DISPLAY_NAME.....Geri Chair
HELP_TEXT.....No Help Text
ID.....{6B8F87BB-C2AB-4D92-AE04-24F19B26D2F1}
TERM.....RESTRAINT GERI CHAIR

DATA VALUE: RESTRAINT HARD PLASTIC/LEATHER

ABBREVIATION.....RES H PL/L
ACTIVE.....Yes
DESCRIPTION.....Hard Plastic/Leather Restraint
DISPLAY_NAME.....Hard Plastic/Leather
HELP_TEXT.....No Help Text
ID.....{9BB528F4-7016-4030-B73F-4037ADF2621D}
TERM.....RESTRAINT HARD PLASTIC/LEATHER

DATA VALUE: RESTRAINT HARD PLASTIC/LEATHER 2 POINT

ABBREVIATION.....RES P/L 2P
ACTIVE.....Yes
DESCRIPTION.....Hard Plastic/Leather 2 Point Restraint
DISPLAY_NAME.....Hard Plastic/Leather 2 Point
HELP_TEXT.....No Help Text
ID.....{17CB8739-612A-4319-946A-2062C5245C90}
TERM.....RESTRAINT HARD PLASTIC/LEATHER 2 POINT

DATA VALUE: RESTRAINT HARD PLASTIC/LEATHER 4 POINT

ABBREVIATION.....RES P/L 4P
ACTIVE.....Yes
DESCRIPTION.....Hard Plastic/Leather 4 Point Restraint
DISPLAY_NAME.....Hard Plastic/Leather 4 Point
HELP_TEXT.....No Help Text
ID.....{94DFD0BA-E941-4B09-8FD7-87E9AD7C7E6F}
TERM.....RESTRAINT HARD PLASTIC/LEATHER 4 POINT

DATA VALUE: RESTRAINT ONE ARM

ABBREVIATION.....RES ONE AR
ACTIVE.....Yes
DESCRIPTION.....One Arm Restraint
DISPLAY_NAME.....One Arm
HELP_TEXT.....No Help Text
ID.....{839914DA-6D63-43EB-961E-70DA56737A92}
TERM.....RESTRAINT ONE ARM

Page: 839

DATA VALUE: RESTRAINT ONE LEG

ABBREVIATION.....RES ONE LE
ACTIVE.....Yes
DESCRIPTION.....One Leg Restraint
DISPLAY_NAME.....One Leg
HELP_TEXT.....No Help Text
ID.....{2790F0FC-B49A-4C99-9968-BCADD7CC63AC}
TERM.....RESTRAINT ONE LEG

DATA VALUE: RESTRAINT SINGLE HAND MITT

ABBREVIATION.....RES S HA M
ACTIVE.....Yes
DESCRIPTION.....Single Hand Mitt Restraint
DISPLAY_NAME.....Single Hand Mitt
HELP_TEXT.....No Help Text
ID.....{A5F1D037-FB5E-402B-9F25-22DE48404D57}
TERM.....RESTRAINT SINGLE HAND MITT

DATA VALUE: RESTRAINT SOFT

ABBREVIATION.....RESTR SOFT
ACTIVE.....Yes

DESCRIPTION.....Soft Restraint
DISPLAY_NAME.....Soft Restraint
HELP_TEXT.....No Help Text
ID.....{C02CAC4B-5E9C-4060-BECD-0C0820EDA39B}
TERM.....RESTRAINT SOFT

DATA VALUE: RESTRAINT SOFT 2 POINT

ABBREVIATION.....RES S2POIN
ACTIVE.....Yes
DESCRIPTION.....Soft 2 Point Restraint
DISPLAY_NAME.....Soft Restraint 2 Point
HELP_TEXT.....No Help Text
ID.....{66B71CFE-2603-40C9-9625-24448012BAB5}
TERM.....RESTRAINT SOFT 2 POINT

DATA VALUE: RESTRAINT SOFT 4 POINT

ABBREVIATION.....RES S4POIN
ACTIVE.....Yes
DESCRIPTION.....Soft 4 Point Restraint
DISPLAY_NAME.....Soft Restraint 4 Point
HELP_TEXT.....No Help Text
ID.....{6820E997-60A0-4802-BE26-CA29B4F48F32}
TERM.....RESTRAINT SOFT 4 POINT

Page: 840

DATA VALUE: RESTRAINT SOFT ARMS + LOCKED WAIST

ABBREVIATION.....RES S A+LW
ACTIVE.....Yes
DESCRIPTION.....Soft Arms + Locked Waist Restraint
DISPLAY_NAME.....Soft Arms + Locked Waist
HELP_TEXT.....No Help Text
ID.....{C480D67C-958C-42F9-9203-98E2F2BDBEC8}
TERM.....RESTRAINT SOFT ARMS + LOCKED WAIST

DATA VALUE: RESTRAINT SOFT ARMS + WAIST

ABBREVIATION.....RES S A+W
ACTIVE.....Yes
DESCRIPTION.....Soft Arms + Waist Restraints
DISPLAY_NAME.....Soft Arms + Waist
HELP_TEXT.....No Help Text
ID.....{4A322BB6-0F63-4BB7-9290-C0655BD7D772}
TERM.....RESTRAINT SOFT ARMS + WAIST

DATA VALUE: RESTRAINT STRAP

ABBREVIATION.....RES STRAP
ACTIVE.....Yes
DESCRIPTION.....Restraint Strap
DISPLAY_NAME.....Restraint Strap
HELP_TEXT.....No Help Text

ID.....{642FDACE-B320-42BC-9514-F5FCD6583B19}
TERM.....RESTRAINT STRAP

DATA VALUE: RESTRAINT VEST

ABBREVIATION.....RES VEST
ACTIVE.....Yes
DESCRIPTION.....Restraint Vest
DISPLAY_NAME.....Restraint Vest
HELP_TEXT.....No Help Text
ID.....{643378B2-86D7-4DB7-93ED-9DCE16D85E27}
TERM.....RESTRAINT VEST

DATA VALUE: RESTRAINT WAIST

ABBREVIATION.....RES WAIST
ACTIVE.....Yes
DESCRIPTION.....Waist
DISPLAY_NAME.....Waist
HELP_TEXT.....No Help Text
ID.....{253E2757-5A57-4410-9C0E-F4A907A1EC16}
TERM.....RESTRAINT WAIST

Page: 841

DATA VALUE: RESTRAINT WAIST LOCKED

ABBREVIATION.....RES WAI LO
ACTIVE.....Yes
DESCRIPTION.....Waist Locked Restraint
DISPLAY_NAME.....Waist Locked
HELP_TEXT.....No Help Text
ID.....{1EE7AB42-CE17-4D49-849B-CA550966E58C}
TERM.....RESTRAINT WAIST LOCKED

DATA VALUE: RESTRAINT WRIST

ABBREVIATION.....RES WRIST
ACTIVE.....Yes
DESCRIPTION.....Wrist Restraint
DISPLAY_NAME.....Restraint Wrist
HELP_TEXT.....No Help Text
ID.....{03038E52-A66F-4EC0-8DC1-7B0BB596FFC0}
TERM.....RESTRAINT WRIST

DATA VALUE: RESTRAINTS BROCHURE

ABBREVIATION.....RES BRO P
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Education Brochure Provided
HELP_TEXT.....No Help Text
ID.....{1E54C2A1-8B8D-4F0D-B247-91E19722CDED}
TERM.....RESTRAINTS BROCHURE

DATA VALUE: RESTRAINTS ED REMOVAL

ABBREVIATION.....RES ED REM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Education On Condition Required For
Removal
HELP_TEXT.....No Help Text
ID.....{B19316C8-687B-48FB-9EBC-9857662B9F85}
TERM.....RESTRAINTS ED REMOVAL

DATA VALUE: RESTRAINTS ED TYPE

ABBREVIATION.....RES ED TY
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Education Rational For Type Of
Restraint
HELP_TEXT.....No Help Text
ID.....{5560A332-C365-4DAF-BA33-610C4AE301F5}
TERM.....RESTRAINTS ED TYPE

Page: 842

DATA VALUE: RESTRAINTS EXPLAINED

ABBREVIATION.....RES EXPL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Explain Why Restraints Applied
HELP_TEXT.....No Help Text
ID.....{1C063DD2-C4A3-4027-A245-6DD51DF5A3E9}
TERM.....RESTRAINTS EXPLAINED

DATA VALUE: RESTRAINTS LIMIT SETTING

ABBREVIATION.....RES LIM SE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Limit Setting
HELP_TEXT.....No Help Text
ID.....{F97BA57E-ED33-D6FE-5E35-E6EFB0625773}
TERM.....RESTRAINTS LIMIT SETTING

DATA VALUE: RESTRAINTS TIME LIMIT

ABBREVIATION.....RES TIME L
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Time Limit On Restraints
HELP_TEXT.....No Help Text
ID.....{513B483E-AD7B-4B1A-BD01-91133C605C3C}
TERM.....RESTRAINTS TIME LIMIT

DATA VALUE: RHONCHI

ABBREVIATION.....RHONCHI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Rhonchi
HELP_TEXT.....No Help Text
ID.....{850A9D2A-2F7F-A812-3D56-EE0948B0A0CC}
TERM.....RHONCHI

DATA VALUE: RIGHT

ABBREVIATION.....RIGHT
ACTIVE.....Yes
DESCRIPTION.....Right
DISPLAY_NAME.....Right
HELP_TEXT.....No Help Text
ID.....{AB74C579-DF14-4D7B-A9EF-5F119089F976}
TERM.....RIGHT

Page: 843

DATA VALUE: RIGID HEEL SPLINT

ABBREVIATION.....RIG HEEL S
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Rigid Heel Splint
HELP_TEXT.....No Help Text
ID.....{5A231387-9FAC-4B2D-94D6-5DF4312E9AF8}
TERM.....RIGID HEEL SPLINT

DATA VALUE: ROTATIONAL THERAPY

ABBREVIATION.....ROT THERAP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Rotational Therapy
HELP_TEXT.....No Help Text
ID.....{BB795594-720C-40D7-9B7C-6792BC849DA0}
TERM.....ROTATIONAL THERAPY

DATA VALUE: ROUND

ABBREVIATION.....ROUND
ACTIVE.....Yes
DESCRIPTION.....Round
DISPLAY_NAME.....Round
HELP_TEXT.....No Help Text
ID.....{E8E741F1-F0C5-4B77-B014-63B545BED755}
TERM.....ROUND

DATA VALUE: RUB

ABBREVIATION.....r
ACTIVE.....Yes

DESCRIPTION.....No Description
DISPLAY_NAME.....Rub
HELP_TEXT.....No Help Text
ID.....{0040835D-1CF0-7D55-EEBB-A94A1E0D1206}
TERM.....RUB

DATA VALUE: S1

ABBREVIATION.....s
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....S1
HELP_TEXT.....No Help Text
ID.....{1AB9DF0A-8238-E8FB-8520-4749E13DDE1D}
TERM.....S1

Page: 844

DATA VALUE: S2

ABBREVIATION.....s
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....S2
HELP_TEXT.....No Help Text
ID.....{57C043D3-B227-3E7A-C1E6-5129A2800718}
TERM.....S2

DATA VALUE: S3

ABBREVIATION.....s
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....S3
HELP_TEXT.....No Help Text
ID.....{9374B57F-DE1F-B9C4-26B5-2DD6B194DC26}
TERM.....S3

DATA VALUE: S4

ABBREVIATION.....s
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....S4
HELP_TEXT.....No Help Text
ID.....{9C6CA302-85B6-8FFC-EAB4-8FED62C535A5}
TERM.....S4

DATA VALUE: SAFETY

ABBREVIATION.....SAFETY
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Safety
HELP_TEXT.....No Help Text

ID.....{AD905A8F-33BD-4EFE-929E-A2B575797402}
TERM.....SAFETY

DATA VALUE: SALINE LOCK

ABBREVIATION.....SALIN LOCK
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Saline Lock
HELP_TEXT.....No Help Text
ID.....{A94D7639-2A45-4319-A6C7-5A892AA6B1F4}
TERM.....SALINE LOCK

Page: 845

DATA VALUE: SALINE MOIST DRESSING

ABBREVIATION.....SAL MOI DR
ACTIVE.....Yes
DESCRIPTION.....Saline Moist Dressing
DISPLAY_NAME.....Saline Moist Dressing
HELP_TEXT.....No Help Text
ID.....{19F98514-8ABF-461C-9969-C31177FE7DE7}
TERM.....SALINE MOIST DRESSING

DATA VALUE: SANDOSOURCE PEPTIDE

ABBREVIATION.....SAND PEPTI
ACTIVE.....Yes
DESCRIPTION.....Sandosource Peptide Enteral Nutrition
Formula
DISPLAY_NAME.....Sandosource Peptide Enteral Nutrition
Formula
HELP_TEXT.....No Help Text
ID.....{0EF43998-4ED4-47F5-ACC6-5A977B6F9997}
TERM.....SANDOSOURCE PEPTIDE

DATA VALUE: SAS AGITATED

ABBREVIATION.....SAS AGITAT
ACTIVE.....Yes
DESCRIPTION.....Anxious or physically agitated, calms
to verbal instructions
DISPLAY_NAME.....5 - Agitated
HELP_TEXT.....Anxious or physically agitated, calms
to verbal instructions
ID.....{107143D4-9527-4A26-8CB7-EEBB656B68B3}
TERM.....SAS AGITATED

DATA VALUE: SAS CALM AND COOPERATIVE

ABBREVIATION.....SAS C COOP
ACTIVE.....Yes
DESCRIPTION.....Calm, easily arousable, follows
commands

DISPLAY_NAME.....4 - Calm and Cooperative
HELP_TEXT.....Calm, easily arousable, follows
commands
ID.....{F0633778-280C-43C3-9E6E-9C35C18179F6}
TERM.....SAS CALM AND COOPERATIVE

DATA VALUE: SAS DANGEROUS AGITATION

ABBREVIATION.....SAS DA AGI
ACTIVE.....Yes
DESCRIPTION.....Pulling at ET tube, trying to remove
catheters, climbing over bedrail,
striking at staff, thrashing
side-to-side
DISPLAY_NAME.....7 - Dangerous Agitation
HELP_TEXT.....Pulling at ET tube, trying to remove
catheters, climbing over bedrail,
striking at staff, thrashing
side-to-side
ID.....{0701FB25-1DF0-4771-8C81-AC19A985A3B6}
TERM.....SAS DANGEROUS AGITATION

Page: 846

DATA VALUE: SAS SEDATED

ABBREVIATION.....SAS SEDATE
ACTIVE.....Yes
DESCRIPTION.....Difficult to arouse but awakens to
verbal stimuli or gentle shaking,
follows simple commands but drifts
off again
DISPLAY_NAME.....3 - Sedated
HELP_TEXT.....Difficult to arouse but awakens to
verbal stimuli or gentle shaking,
follows simple commands but drifts
off again
ID.....{AFEC3F29-52D3-4A49-91F8-05AA3891FEC8}
TERM.....SAS SEDATED

DATA VALUE: SAS UNAROUSABLE

ABBREVIATION.....SAS UNAROU
ACTIVE.....Yes
DESCRIPTION.....Minimal or no response to noxious
stimuli, does not communicate or
follow commands
DISPLAY_NAME.....1 - Unarousable
HELP_TEXT.....Minimal or no response to noxious
stimuli, does not communicate or
follow commands
ID.....{3DDFF9FE-A88A-471F-BF36-E5B569DEBDAA}
TERM.....SAS UNAROUSABLE

DATA VALUE: SAS VERY AGITATED

ABBREVIATION.....SAS VY AGI
ACTIVE.....Yes
DESCRIPTION.....Requiring restraint and frequent
verbal reminding of limits, biting
ETT
DISPLAY_NAME.....6 - Very Agitated
HELP_TEXT.....Requiring restraint and frequent
verbal reminding of limits, biting
ETT
ID.....{B501606A-EEA7-470E-81BE-F51D422272F1}
TERM.....SAS VERY AGITATED

DATA VALUE: SAS VERY SEDATED

ABBREVIATION.....SAS VY SED
ACTIVE.....Yes
DESCRIPTION.....Arouses to physical stimuli but does
not communicate or follow commands,
may move spontaneously
DISPLAY_NAME.....2 - Very Sedated
HELP_TEXT.....Arouses to physical stimuli but does
not communicate or follow commands,
may move spontaneously
ID.....{4B13529B-61C5-4B68-857A-3A5351DC4BE2}
TERM.....SAS VERY SEDATED

DATA VALUE: SCANT

ABBREVIATION.....SCANT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Scant
HELP_TEXT.....No Help Text
ID.....{8EF6C62D-0D94-37BB-60AD-86E306E8120D}
TERM.....SCANT

Page: 847

DATA VALUE: SCAPHOID

ABBREVIATION.....SCAPHOID
ACTIVE.....Yes
DESCRIPTION.....Scaphoid
DISPLAY_NAME.....Scaphoid
HELP_TEXT.....No Help Text
ID.....{1BE93489-42FC-44EB-8D07-87F4224979F3}
TERM.....SCAPHOID

DATA VALUE: SECLUSION

ABBREVIATION.....SECLUSION
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Seclusion

HELP_TEXT.....No Help Text
ID.....{65B57F6C-C849-4BF7-8705-175361BC900F}
TERM.....SECLUSION

DATA VALUE: SECOND DEGREE AV BLOCK: TYPE 1

ABBREVIATION.....sdabt1
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Second Degree Av Block: Type 1
HELP_TEXT.....No Help Text
ID.....{508BF663-96FB-EEBE-4F06-D2F1A02BD099}
TERM.....SECOND DEGREE AV BLOCK: TYPE 1

DATA VALUE: SECOND DEGREE AV BLOCK: TYPE 2

ABBREVIATION.....sdabt2
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Second Degree Av Block: Type 2
HELP_TEXT.....No Help Text
ID.....{3367754A-08D1-31A9-E850-1B36358E358A}
TERM.....SECOND DEGREE AV BLOCK: TYPE 2

DATA VALUE: SEDATED

ABBREVIATION.....s
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Sedated
HELP_TEXT.....No Help Text
ID.....{29D9ADD6-2C1F-F6D5-F4D0-B2BB86416997}
TERM.....SEDATED

Page: 848

DATA VALUE: SEDIMENT

ABBREVIATION.....SEDIMENT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Sediment
HELP_TEXT.....No Help Text
ID.....{4D4E4236-55CE-7ECD-FC55-6AD0CFD4D9A5}
TERM.....SEDIMENT

DATA VALUE: SELF AM CARE BY PATIENT

ABBREVIATION.....SELF AM CR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Self AM Care By Patient
HELP_TEXT.....No Help Text
ID.....{C48A46C1-0E68-479C-B803-08C80617ADFF}
TERM.....SELF AM CARE BY PATIENT

DATA VALUE: SELF CARE

ABBREVIATION.....SELF CARE
ACTIVE.....Yes
DESCRIPTION.....Self Care
DISPLAY_NAME.....Self Care
HELP_TEXT.....No Help Text
ID.....{0E068188-75B8-46C5-8947-A6CC5C7849B5}
TERM.....SELF CARE

DATA VALUE: SELF CARE OF BEDSIDE ENVIRONMENT

ABBREVIATION.....S CARE BED
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Self Care of Bedside Environment
HELP_TEXT.....No Help Text
ID.....{A1BEC7F4-D404-4824-8BB7-108E0A6C43A2}
TERM.....SELF CARE OF BEDSIDE ENVIRONMENT

DATA VALUE: SELF FACE WASH

ABBREVIATION.....S FACE WAS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Self Face Wash
HELP_TEXT.....No Help Text
ID.....{77457D46-5676-46E6-B38F-FE3AF70B5CBC}
TERM.....SELF FACE WASH

Page: 849

DATA VALUE: SELF FOOT CARE

ABBREVIATION.....SELF FOOT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Self Foot Care
HELP_TEXT.....No Help Text
ID.....{F2DCAF0E-389D-4987-8C73-66E35266CC78}
TERM.....SELF FOOT CARE

DATA VALUE: SELF SHAVE

ABBREVIATION.....SELF SHAVE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Self Shave
HELP_TEXT.....No Help Text
ID.....{90F1760C-4C12-4976-ADCB-DBF456A34023}
TERM.....SELF SHAVE

DATA VALUE: SELF SHOWER

ABBREVIATION.....SELF SHOWER
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Self Shower
HELP_TEXT.....No Help Text
ID.....{C35CCAE8-27D9-46F2-8971-2C4CF4924CCB}
TERM.....SELF SHOWER

DATA VALUE: SELF TRACH CARE

ABBREVIATION.....S TRACH CA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Self Trach Care
HELP_TEXT.....No Help Text
ID.....{70BE79EC-107F-4CD3-B964-7919B3DDE222}
TERM.....SELF TRACH CARE

DATA VALUE: SEMIFOWLERS HIGH

ABBREVIATION.....sh
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Semifowlers High
HELP_TEXT.....No Help Text
ID.....{1BBE7175-2BD2-BFB5-1E1B-E02443433072}
TERM.....SEMIFOWLERS HIGH

Page: 850

DATA VALUE: SEMIFOWLERS LOW

ABBREVIATION.....sl
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Semifowlers Low
HELP_TEXT.....No Help Text
ID.....{A265AF02-53FC-B182-8AE3-A06E95718D7E}
TERM.....SEMIFOWLERS LOW

DATA VALUE: SEQUENTIAL COMPRESSION DEVICE

ABBREVIATION.....SCD
ACTIVE.....Yes
DESCRIPTION.....SEQUENTIAL COMPRESSION DEVICE
DISPLAY_NAME.....Sequential Compression Device
HELP_TEXT.....SEQUENTIAL COMPRESSION DEVICE
ID.....{29D657DB-6679-8A55-0E33-235B5DD349F2}
TERM.....SEQUENTIAL COMPRESSION DEVICE

DATA VALUE: SEROSANGUINOUS

ABBREVIATION.....SEROSANGUI
ACTIVE.....Yes
DESCRIPTION.....No Description

DISPLAY_NAME.....Serosanguinous
HELP_TEXT.....No Help Text
ID.....{19A8D332-899D-177A-2247-ABCDAAEE3BFA8}
TERM.....SEROSANGUINOUS

DATA VALUE: SEROUS

ABBREVIATION.....SEROUS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Serous
HELP_TEXT.....No Help Text
ID.....{9BDA347E-29E9-80C7-91FA-016192D4257E}
TERM.....SEROUS

DATA VALUE: SET UP HELP ONLY

ABBREVIATION.....SETUP ONLY
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Set Up Help only
HELP_TEXT.....No Help Text
ID.....{DACC77AB-5AE7-4F1B-A3E1-C347D3DD7219}
TERM.....SET UP HELP ONLY

Page: 851

DATA VALUE: SEVERELY REDUCED

ABBREVIATION.....SEV REDUCE
ACTIVE.....Yes
DESCRIPTION.....Severely Reduced
DISPLAY_NAME.....Severely Reduced
HELP_TEXT.....No Help Text
ID.....{A895BACE-85AC-45B2-A190-4A6816DEDC1D}
TERM.....SEVERELY REDUCED

DATA VALUE: SHALLOW

ABBREVIATION.....SHALLOW
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Shallow
HELP_TEXT.....No Help Text
ID.....{CD205617-BF98-4177-971A-056AFECC823A}
TERM.....SHALLOW

DATA VALUE: SHAMPOO

ABBREVIATION.....SHAMPOO
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Shampoo
HELP_TEXT.....No Help Text
ID.....{95CCB23C-E614-40CD-84F0-EE4EFBC95271}

TERM.....SHAMPOO

DATA VALUE: SHARP

ABBREVIATION.....SHARP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Sharp
HELP_TEXT.....No Help Text
ID.....{20608864-E5CF-4FE3-8B6D-AE2FEE086503}
TERM.....SHARP

DATA VALUE: SHAVE

ABBREVIATION.....SHAVE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Shave
HELP_TEXT.....No Help Text
ID.....{A5CB33DF-C33B-4DD3-B00D-78D78F7966C7}
TERM.....SHAVE

Page: 852

DATA VALUE: SHOOTING

ABBREVIATION.....SHOOTING
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Shooting
HELP_TEXT.....No Help Text
ID.....{73C48DB1-6FE6-46B7-B1F6-E93ACBE44C27}
TERM.....SHOOTING

DATA VALUE: SHOWER

ABBREVIATION.....SHOWER
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Shower
HELP_TEXT.....No Help Text
ID.....{C259929A-44DB-48FD-A2EC-BC6ECF9D0476}
TERM.....SHOWER

DATA VALUE: SICU

ABBREVIATION.....SICU
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....SICU
HELP_TEXT.....No Help Text
ID.....{EAEB2DE6-B057-4143-B44C-99DAB9B7B4AC}
TERM.....SICU

DATA VALUE: SIGMOIDOSTOMY

ABBREVIATION.....SIGMOIDOST
ACTIVE.....Yes
DESCRIPTION.....Sigmoidostomy
DISPLAY_NAME.....Sigmoidostomy
HELP_TEXT.....No Help Text
ID.....{0C84B317-59B3-433B-A9AD-6588C1361E0C}
TERM.....SIGMOIDOSTOMY

DATA VALUE: SIGNIFICANT OTHER

ABBREVIATION.....SO
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Significant Other
HELP_TEXT.....No Help Text
ID.....{082C42E8-9E8A-4A96-B3AD-8B0E0A62AF00}
TERM.....SIGNIFICANT OTHER

Page: 853

DATA VALUE: SILVADENE

ABBREVIATION.....SILVADENE
ACTIVE.....Yes
DESCRIPTION.....Silvadene
DISPLAY_NAME.....Silvadene
HELP_TEXT.....No Help Text
ID.....{0889134F-C419-4366-B12E-EAB775719410}
TERM.....SILVADENE

DATA VALUE: SILVER

ABBREVIATION.....SILVER
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Silver
HELP_TEXT.....No Help Text
ID.....{AF233570-2BB4-4AE6-8751-056B9717C475}
TERM.....SILVER

DATA VALUE: SILVER SULFADIAZINE

ABBREVIATION.....SILVER SUL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Silver Sulfadiazine
HELP_TEXT.....No Help Text
ID.....{E7A99A4E-C9B5-40D6-9820-05F87004F48B}
TERM.....SILVER SULFADIAZINE

DATA VALUE: SINGLE

ABBREVIATION.....SINGLE
ACTIVE.....Yes

DESCRIPTION.....No Description
DISPLAY_NAME.....Single
HELP_TEXT.....No Help Text
ID.....{8D305E0F-9035-E396-E768-40D7E6455BE4}
TERM.....SINGLE

DATA VALUE: SINUS ARREST

ABBREVIATION.....sa
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Sinus Arrest
HELP_TEXT.....No Help Text
ID.....{3CD8F7AC-3FE0-F342-089C-787232D2A5C2}
TERM.....SINUS ARREST

Page: 854

DATA VALUE: SINUS ARRHYTHMIA

ABBREVIATION.....sa
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Sinus Arrhythmia
HELP_TEXT.....No Help Text
ID.....{07A24AC9-B486-A4B8-BADB-04D2729542CE}
TERM.....SINUS ARRHYTHMIA

DATA VALUE: SINUS BLOCK

ABBREVIATION.....sb
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Sinus Block
HELP_TEXT.....No Help Text
ID.....{94D6EFED-3D08-E724-BC40-9C398919A142}
TERM.....SINUS BLOCK

DATA VALUE: SINUS BRADYCARDIA

ABBREVIATION.....sb
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Sinus Bradycardia
HELP_TEXT.....No Help Text
ID.....{D6EF4940-63AE-22AF-1A80-D8E6F0F07839}
TERM.....SINUS BRADYCARDIA

DATA VALUE: SINUS PAUSE

ABBREVIATION.....sp
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Sinus Pause
HELP_TEXT.....No Help Text

ID.....{7F8F89A9-6B16-B159-55B3-DEFDF5B4B4F6}
TERM.....SINUS PAUSE

DATA VALUE: SINUS RHYTHM

ABBREVIATION.....sr
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Sinus Rhythm
HELP_TEXT.....No Help Text
ID.....{A34F80EA-3D6A-076A-DB2A-E4E5C5AACA11}
TERM.....SINUS RHYTHM

Page: 855

DATA VALUE: SINUS TACHYCARDIA

ABBREVIATION.....SINUS T
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Sinus Tachycardia
HELP_TEXT.....No Help Text
ID.....{0AC41AF3-7656-A35D-E397-21D35405812A}
TERM.....SINUS TACHYCARDIA

DATA VALUE: SKIN INCISION WOUND CLOSURE OPEN WOUND

ABBREVIATION.....INC CL OPE
ACTIVE.....Yes
DESCRIPTION.....Skin Incision Closure Open Wound
DISPLAY_NAME.....Open Wound
HELP_TEXT.....Skin Incision Closure Open Wound
ID.....{AD018B6F-13E5-41C0-BDBC-306CA0B1B475}
TERM.....SKIN INCISION WOUND CLOSURE OPEN WOUND

DATA VALUE: SKIN INCISION WOUND CLOSURE STAPLES

ABBREVIATION.....INC CL STA
ACTIVE.....Yes
DESCRIPTION.....Skin Incision Closure Staples
DISPLAY_NAME.....Staples
HELP_TEXT.....Skin Incision Closure Staples
ID.....{4F48D42D-C206-423C-970B-02C03B2EE0B6}
TERM.....SKIN INCISION WOUND CLOSURE STAPLES

DATA VALUE: SKIN INCISION WOUND CLOSURE STERI-STRIPS

ABBREVIATION.....INC CL STE
ACTIVE.....Yes
DESCRIPTION.....Skin Incision Closure Steri-Strips
DISPLAY_NAME.....Steri-Strips
HELP_TEXT.....Skin Incision Closure Steri-Strips
ID.....{DF52B6F3-54EE-46B3-AFDF-4A701B065E04}
TERM.....SKIN INCISION WOUND CLOSURE
STERI-STRIPS

DATA VALUE: SKIN INCISION WOUND CLOSURE SUTURES

ABBREVIATION.....INC CL SUT
ACTIVE.....Yes
DESCRIPTION.....Skin Incision Closure Sutures
DISPLAY_NAME.....Sutures
HELP_TEXT.....Skin Incision Closure Sutures
ID.....{DED8DEDA-82C9-4D3A-A81E-641583314919}
TERM.....SKIN INCISION WOUND CLOSURE SUTURES

Page: 856

DATA VALUE: SKIN INDURATION AROUND INCISION EXTENDS INTO PERIWOUND LESS THAN 1 cm

ABBREVIATION.....INC IN <1
ACTIVE.....Yes
DESCRIPTION.....Skin Induration Extension Around
Incision Extends into Periwound Less
Than1 cm
DISPLAY_NAME.....Extends Into Periwound Less Than 1 cm
HELP_TEXT.....Skin Induration Extension Around
Incision Extends into Periwound Less
Than1 cm
ID.....{2A5097BF-955C-4F79-B483-7EE14AF05786}
TERM.....SKIN INDURATION AROUND INCISION
EXTENDS INTO PERIWOUND LESS THAN 1 cm

DATA VALUE: SKIN INDURATION AROUND INCISION EXTENDS INTO PERIWOUND MORE THAN 1 cm

ABBREVIATION.....INC IN >1
ACTIVE.....Yes
DESCRIPTION.....Skin Induration Extension Around
Incision Extends into Periwound
Greater Than 1 cm
DISPLAY_NAME.....Extends into Periwound Greater Than 1
cm
HELP_TEXT.....Skin Induration Extension Around
Incision Extends into Periwound
Greater Than 1 cm
ID.....{D8EE448D-00EB-412D-95AC-2E2E92C53EEB}
TERM.....SKIN INDURATION AROUND INCISION
EXTENDS INTO PERIWOUND MORE THAN 1 cm

DATA VALUE: SKIN INDURATION AROUND INCISION LIMITED TO INCISION

ABBREVIATION.....INC IN LIM
ACTIVE.....Yes
DESCRIPTION.....Skin Induration Extension Around
Incision Limited to Incision
DISPLAY_NAME.....Limited to Incision
HELP_TEXT.....Skin Induration Extension Around
Incision Limited to Incision

ID.....{88BCE7F9-F21F-427A-A50B-C1BA831639AB}
TERM.....SKIN INDURATION AROUND INCISION
LIMITED TO INCISION

DATA VALUE: SKIN INDURATION AROUND INCISION NONE

ABBREVIATION.....INC IN NON
ACTIVE.....Yes
DESCRIPTION.....Skin Induration Extension Around
Incision None
DISPLAY_NAME.....None
HELP_TEXT.....Skin Induration Extension Around
Incision None

ID.....{A9B95680-F9A0-44C0-929C-70C00B9AA5B7}
TERM.....SKIN INDURATION AROUND INCISION NONE

DATA VALUE: SKIN INDURATION EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND

<
1 cm

ABBREVIATION.....WOU IN <1
ACTIVE.....Yes
DESCRIPTION.....No DescriptionSkin Induration
Extension Around Wound Extends into
Periwound Less Than 1 cm
DISPLAY_NAME.....Extends into Periwound Less Than 1 cm
HELP_TEXT.....No DescriptionSkin Induration
Extension Around Wound Extends into
Periwound Less Than 1 cm
ID.....{3594CCCD-FA51-4971-A61E-FFA3B4AC7EE0}
TERM.....SKIN INDURATION EXTENSION AROUND
WOUND EXTENDS INTO PERIWOUND < 1 cm

Page: 857

DATA VALUE: SKIN INDURATION EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND

>
1 cm

ABBREVIATION.....WOU IN >1
ACTIVE.....Yes
DESCRIPTION.....Skin Induration Extension Around
Wound Extends into Periwound More
Than 1 cm
DISPLAY_NAME.....Extends into Periwound More Than 1 cm
HELP_TEXT.....Skin Induration Extension Around
Wound Extends into Periwound More
Than 1 cm
ID.....{7DF7D08D-3E5E-4D18-94E4-A0DFDA5DB72F}
TERM.....SKIN INDURATION EXTENSION AROUND
WOUND EXTENDS INTO PERIWOUND > 1 cm

DATA VALUE: SKIN INDURATION EXTENSION AROUND WOUND LIMITED TO WOUND

ABBREVIATION.....WOU IN LIM

ACTIVE.....Yes
DESCRIPTION.....Skin Induration Extension Around
Wound Limited to Wound
DISPLAY_NAME.....Limited to Wound
HELP_TEXT.....Skin Induration Extension Around
Wound Limited to Wound
ID.....{1828A9EA-FF9D-4A49-9543-15B85FBF211E}
TERM.....SKIN INDURATION EXTENSION AROUND
WOUND LIMITED TO WOUND

DATA VALUE: SKIN INDURATION EXTENSION AROUND WOUND NONE

ABBREVIATION.....WOU IN NON
ACTIVE.....Yes
DESCRIPTION.....Skin Induration Extension Around
Wound None
DISPLAY_NAME.....None
HELP_TEXT.....Skin Induration Extension Around
Wound None
ID.....{2C57CD27-4FDD-4DE8-88A0-D438CBE535D8}
TERM.....SKIN INDURATION EXTENSION AROUND
WOUND NONE

DATA VALUE: SKIN SWELLING EXTENSION AROUND INCISION EXTENDS INTO PERIWOUND
< 1 cm

ABBREVIATION.....INC SW <1
ACTIVE.....Yes
DESCRIPTION.....Skin Swelling Extension Around
Incision Extends into Periwound Less
Than 1 cm
DISPLAY_NAME.....Extends into Periwound < 1 cm
HELP_TEXT.....Skin Swelling Extension Around
Incision Extends into Periwound Less
Than 1 cm
ID.....{ABA0527C-8FAD-4E34-9F6F-2AA97EA06E20}
TERM.....SKIN SWELLING EXTENSION AROUND
INCISION EXTENDS INTO PERIWOUND < 1
cm

DATA VALUE: SKIN SWELLING EXTENSION AROUND INCISION EXTENDS INTO PERIWOUND
> 1 cm

ABBREVIATION.....INC SW >1
ACTIVE.....Yes
DESCRIPTION.....Skin Swelling Extension Around
Incision Extends into Periwound
Greater Than 1 cm
DISPLAY_NAME.....Extends into Periwound More Than 1 cm
HELP_TEXT.....Skin Swelling Extension Around
Incision Extends into Periwound
Greater Than 1 cm
ID.....{106A2AD3-F482-4DEE-967E-6F06C6026944}
TERM.....SKIN SWELLING EXTENSION AROUND
INCISION EXTENDS INTO PERIWOUND > 1
cm

DATA VALUE: SKIN SWELLING EXTENSION AROUND INCISION FLUCTUANT

ABBREVIATION.....INC SW FLU
ACTIVE.....Yes
DESCRIPTION.....Skin Swelling Extension Around
 Incision Fluctuant
DISPLAY_NAME.....Fluctuant
HELP_TEXT.....Skin Swelling Extension Around
 Incision Fluctuant
ID.....{5D38DC7A-8E40-4D27-A039-2D6A557B4FFD}
TERM.....SKIN SWELLING EXTENSION AROUND
 INCISION FLUCTUANT

DATA VALUE: SKIN SWELLING EXTENSION AROUND INCISION LIMITED TO WOUND

ABBREVIATION.....INC SW LIM
ACTIVE.....Yes
DESCRIPTION.....Skin Swelling Extension Around
 Incision Limited to Wound
DISPLAY_NAME.....Limited to Wound
HELP_TEXT.....Skin Swelling Extension Around
 Incision Limited to Wound
ID.....{B8503D54-6925-430E-8AB6-692412347D5D}
TERM.....SKIN SWELLING EXTENSION AROUND
 INCISION LIMITED TO WOUND

DATA VALUE: SKIN SWELLING EXTENSION AROUND INCISION NONE

ABBREVIATION.....INC SW NON
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....None
HELP_TEXT.....No Help Text
ID.....{54C9FC46-899C-460B-9524-D6F196B9FF31}
TERM.....SKIN SWELLING EXTENSION AROUND
 INCISION NONE

DATA VALUE: SKIN SWELLING EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND < 1
cm

ABBREVIATION.....WOU SW <1
ACTIVE.....Yes
DESCRIPTION.....Skin Swelling Extension Around Wound
 Extends into Periwound Less Than 1 cm
DISPLAY_NAME.....Extends into Periwound Less Than 1 cm
HELP_TEXT.....Skin Swelling Extension Around Wound
 Extends into Periwound Less Than 1 cm
ID.....{15F314A0-6B77-4EB5-92C0-FFEC4D82C75C}
TERM.....SKIN SWELLING EXTENSION AROUND WOUND
 EXTENDS INTO PERIWOUND < 1 cm

DATA VALUE: SKIN SWELLING EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND > 1

cm

ABBREVIATION.....WOU SW >1
ACTIVE.....Yes
DESCRIPTION.....Skin Swelling Extension Around Wound
 Extends into Periwound Greater Than 1
 cm
DISPLAY_NAME.....Extends Into Periwound More Than 1 cm
HELP_TEXT.....Skin Swelling Extension Around Wound
 Extends into Periwound Greater Than 1
 cm
ID.....{EA6AC933-0737-4C28-9F31-D36F3452DA96}
TERM.....SKIN SWELLING EXTENSION AROUND WOUND
 EXTENDS INTO PERIWOUND > 1 cm

Page: 859

DATA VALUE: SKIN SWELLING EXTENSION AROUND WOUND FLUCTANT

ABBREVIATION.....WOU SW FLU
ACTIVE.....Yes
DESCRIPTION.....Skin Swelling Extension Around Wound
 Fluctuant
DISPLAY_NAME.....Fluctant
HELP_TEXT.....Skin Swelling Extension Around Wound
 Fluctuant
ID.....{DB5077E7-D339-4C87-B8CA-D16269DB2935}
TERM.....SKIN SWELLING EXTENSION AROUND WOUND
 FLUCTANT

DATA VALUE: SKIN SWELLING EXTENSION AROUND WOUND LIMITED TO WOUND

ABBREVIATION.....WOU SW LIM
ACTIVE.....Yes
DESCRIPTION.....Skin Swelling Extension Around Wound
 Limited to Wound
DISPLAY_NAME.....Limited to Wound
HELP_TEXT.....Skin Swelling Extension Around Wound
 Limited to Wound
ID.....{413CD1EC-309C-45D0-A1E0-E829EA615E2A}
TERM.....SKIN SWELLING EXTENSION AROUND WOUND
 LIMITED TO WOUND

DATA VALUE: SKIN SWELLING EXTENSION AROUND WOUND NONE

ABBREVIATION.....WOU SW NON
ACTIVE.....Yes
DESCRIPTION.....Skin Swelling Extension Around Wound
 None
DISPLAY_NAME.....None
HELP_TEXT.....Skin Swelling Extension Around Wound
 None
ID.....{550E3B92-2C69-403B-A997-2DBE61374ABF}
TERM.....SKIN SWELLING EXTENSION AROUND WOUND
 NONE

DATA VALUE: SKIN TEAR

ABBREVIATION.....SKIN TEAR
ACTIVE.....Yes
DESCRIPTION.....Skin Tear
DISPLAY_NAME.....Skin Tear
HELP_TEXT.....No Help Text
ID.....{2E230C6E-119F-4F9F-8797-2FFA9D76F337}
TERM.....SKIN TEAR

DATA VALUE: SKIN TEMPERATURE AT INCISION SITE CONSISTENT WITH BODY TEMPERATURE

ABBREVIATION.....INC TP CON
ACTIVE.....Yes
DESCRIPTION.....Skin Temperature at Incision Site
Consistent with Body Temperature
DISPLAY_NAME.....Skin Temperature at Incision Site
Consistent with Body Temperature
HELP_TEXT.....Skin Temperature at Incision Site
Consistent with Body Temperature
ID.....{D349175A-4B20-4243-BE40-E01E50C5583F}
TERM.....SKIN TEMPERATURE AT INCISION SITE
CONSISTENT WITH BODY TEMPERATURE

Page: 860

DATA VALUE: SKIN TEMPERATURE AT INCISION SITE COOLER THAN GENERAL BODY TEMPERATURE

ABBREVIATION.....INC TP CBT
ACTIVE.....Yes
DESCRIPTION.....Skin Temperature at Incision Site
Cooler than General Body Temperature
DISPLAY_NAME.....Cooler Than General Body Temperature
HELP_TEXT.....Skin Temperature at Incision Site
Cooler than General Body Temperature
ID.....{D357C975-857D-4293-B3EF-EFA7A4586357}
TERM.....SKIN TEMPERATURE AT INCISION SITE
COOLER THAN GENERAL BODY TEMPERATURE

DATA VALUE: SKIN TEMPERATURE AT INCISION SITE COOLER THAN OPPOSITE BODY PART

ABBREVIATION.....INC TP COP
ACTIVE.....Yes
DESCRIPTION.....Skin Temperature at Incision Site
Cooler than Opposite Body Part
DISPLAY_NAME.....Cooler Than Opposite Body Part
HELP_TEXT.....Skin Temperature at Incision Site
Cooler than Opposite Body Part
ID.....{85307AD8-F730-43F0-BCB3-882594D9DDF8}
TERM.....SKIN TEMPERATURE AT INCISION SITE
COOLER THAN OPPOSITE BODY PART

DATA VALUE: SKIN TEMPERATURE AT INCISION SITE WARMER THAN GENERAL BODY TEMPERATURE

ABBREVIATION.....WARM BD OP
ACTIVE.....Yes
DESCRIPTION.....Skin Temperature at Incision Site
Warmer than General Body Temperture
DISPLAY_NAME.....Warmer Than General Body Temp
HELP_TEXT.....Skin Temperature at Incision Site
Warmer than General Body Temperture
ID.....{5F888E9D-8F65-49ED-BC5D-F6F287A9D478}
TERM.....SKIN TEMPERATURE AT INCISION SITE
WARMER THAN GENERAL BODY TEMPERATURE

T DATA VALUE: SKIN TEMPERATURE AT INCISION SITE WARMER THAN OPPOSITE BODY PART

ABBREVIATION.....WARMER OP
ACTIVE.....Yes
DESCRIPTION.....Skin Temperature at Incision Site
Warmer than Opposite Body Part
DISPLAY_NAME.....Warmer Than Opposite Body Part
HELP_TEXT.....Skin Temperature at Incision Site
Warmer than Opposite Body Part
ID.....{F518787A-325A-45F3-82B4-73C856802073}
TERM.....SKIN TEMPERATURE AT INCISION SITE
WARMER THAN OPPOSITE BODY PART

DATA VALUE: SKIN TEMPERATURE AT WOUND SITE CONSISTENT WITH BODY TEMPERATURE

ABBREVIATION.....WOU TP CBT
ACTIVE.....Yes
DESCRIPTION.....Skin Temperature at Wound Site
Consistent with Body Temperature
DISPLAY_NAME.....Consistent With Body Temperature
HELP_TEXT.....Skin Temperature at Wound Site
Consistent with Body Temperature
ID.....{672D56C1-D1A3-4DDF-8815-F05DEE197061}
TERM.....SKIN TEMPERATURE AT WOUND SITE
CONSISTENT WITH BODY TEMPERATURE

Page: 861

DATA VALUE: SKIN TEMPERATURE AT WOUND SITE COOLER THAN GENERAL BODY TEMPERATURE

ABBREVIATION.....WOU TP CGB
ACTIVE.....Yes
DESCRIPTION.....Skin Temperature at Wound Site Cooler
than General Body Temp
DISPLAY_NAME.....Cooler Than General Body Temperature
HELP_TEXT.....Skin Temperature at Wound Site Cooler

than General Body Temp
ID.....{7203993E-3334-44EC-9FD9-14224253283D}
TERM.....SKIN TEMPERATURE AT WOUND SITE COOLER
THAN GENERAL BODY TEMPERATURE

DATA VALUE: SKIN TEMPERATURE AT WOUND SITE COOLER THAN OPPOSITE BODY PART

ABBREVIATION.....WOU TP COB
ACTIVE.....Yes
DESCRIPTION.....Skin Temperature at Wound Site Cooler
than Opposite Body Part
DISPLAY_NAME.....Cooler Than Opposite Body Part
HELP_TEXT.....Skin Temperature at Wound Site Cooler
than Opposite Body Part
ID.....{A904C666-2810-4211-B53A-2961DCE6894E}
TERM.....SKIN TEMPERATURE AT WOUND SITE COOLER
THAN OPPOSITE BODY PART

DATA VALUE: SKIN TEMPERATURE AT WOUND SITE WARMER THAN GENERAL BODY
TEMPERA
TURE

ABBREVIATION.....WOU TP WGB
ACTIVE.....Yes
DESCRIPTION.....Skin Temperature at Wound Site Warmer
than General Body Temperature
DISPLAY_NAME.....Warmer Than General Body Temperature
HELP_TEXT.....Skin Temperature at Wound Site Warmer
than General Body Temperature
ID.....{03BB8F9A-3DAC-4994-9F09-AA1DEEF48090}
TERM.....SKIN TEMPERATURE AT WOUND SITE WARMER
THAN GENERAL BODY TEMPERATURE

DATA VALUE: SKIN TEMPERATURE AT WOUND SITE WARMER THAN OPPOSITE BODY PART

ABBREVIATION.....WOU TP WOB
ACTIVE.....Yes
DESCRIPTION.....Skin Temperature at Wound Site Warmer
than Opposite Body Part
DISPLAY_NAME.....Warmer Than Opposite Body Part
HELP_TEXT.....Skin Temperature at Wound Site Warmer
than Opposite Body Part
ID.....{6290BD23-24BF-4A3A-8555-DCE52D2274B3}
TERM.....SKIN TEMPERATURE AT WOUND SITE WARMER
THAN OPPOSITE BODY PART

DATA VALUE: SKIN TRAUMA

ABBREVIATION.....SKIN TRAUM
ACTIVE.....Yes
DESCRIPTION.....Skin Trauma
DISPLAY_NAME.....Skin Trauma
HELP_TEXT.....No Help Text
ID.....{91C8EE13-2C4A-468C-94EB-C53E0D153EBA}
TERM.....SKIN TRAUMA

DATA VALUE: SKIN WARM AND DRY COLOR APPROPRIATE

ABBREVIATION.....S W A DR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Skin Warm and Dry, Color Appropriate
to Baseline
HELP_TEXT.....No Help Text
ID.....{9637D510-C199-4969-8F2D-0A499151E22B}
TERM.....SKIN WARM AND DRY COLOR APPROPRIATE

DATA VALUE: SKIN WOUND CLOSURE APPROXIMATED

ABBREVIATION.....WOU CL APP
ACTIVE.....Yes
DESCRIPTION.....Skin Wound Closure Approximated
DISPLAY_NAME.....Approximated
HELP_TEXT.....Skin Wound Closure Approximated
ID.....{83D48E0F-8643-436B-9ECC-80795857C322}
TERM.....SKIN WOUND CLOSURE APPROXIMATED

DATA VALUE: SKIN WOUND CLOSURE DEHISCENCE

ABBREVIATION.....WOU CL DEH
ACTIVE.....Yes
DESCRIPTION.....Skin Wound Closure Dehiscence
DISPLAY_NAME.....Dehiscence
HELP_TEXT.....Skin Wound Closure Dehiscence
ID.....{425A460F-19C0-4A4A-8E4E-712C66968C38}
TERM.....SKIN WOUND CLOSURE DEHISCENCE

DATA VALUE: SKIN WOUND CLOSURE GLUE

ABBREVIATION.....WOU CL GLU
ACTIVE.....Yes
DESCRIPTION.....Skin Wound Closure Glue
DISPLAY_NAME.....Glue
HELP_TEXT.....Skin Wound Closure Glue
ID.....{29FE4AF7-39EE-400A-BB3D-E07786F11BFD}
TERM.....SKIN WOUND CLOSURE GLUE

DATA VALUE: SKIN WOUND CLOSURE NOT CLOSED

ABBREVIATION.....WOU CL NOT
ACTIVE.....Yes
DESCRIPTION.....Skin Wound Closure Not Closed
DISPLAY_NAME.....Not Closed
HELP_TEXT.....Skin Wound Closure Not Closed
ID.....{0E69F1FA-35F6-4F34-8547-CC3602670A7B}
TERM.....SKIN WOUND CLOSURE NOT CLOSED

DATA VALUE: SKIN WOUND CLOSURE OPEN WOUND

ABBREVIATION.....WOU CL OPE
ACTIVE.....Yes
DESCRIPTION.....Skin Wound Closure Open Wound
DISPLAY_NAME.....Open Wound
HELP_TEXT.....Skin Wound Closure Open Wound
ID.....{C6F2AD2C-932F-434A-8A7A-513029092FFD}
TERM.....SKIN WOUND CLOSURE OPEN WOUND

DATA VALUE: SKIN WOUND CLOSURE STAPLES

ABBREVIATION.....WOU CL STA
ACTIVE.....Yes
DESCRIPTION.....Skin Wound Closure Staples
DISPLAY_NAME.....Staples
HELP_TEXT.....Skin Wound Closure Staples
ID.....{FE04BFE3-DF81-48E7-9FA1-1CFCF4ECAD71}
TERM.....SKIN WOUND CLOSURE STAPLES

DATA VALUE: SKIN WOUND CLOSURE STERI-STRIPS

ABBREVIATION.....WOU CL STE
ACTIVE.....Yes
DESCRIPTION.....Skin Wound Closure Steri - Strips
DISPLAY_NAME.....Steri-Strips
HELP_TEXT.....Skin Wound Closure Steri-Strips
ID.....{AD451791-6B49-4A51-ABC3-796D08BEF910}
TERM.....SKIN WOUND CLOSURE STERI-STRIPS

DATA VALUE: SKIN WOUND CLOSURE SUTURES

ABBREVIATION.....WOU CL SUT
ACTIVE.....Yes
DESCRIPTION.....Skin Wound Closure Sutures
DISPLAY_NAME.....Sutures
HELP_TEXT.....Skin Wound Closure Sutures
ID.....{F2840AE4-216B-477C-9FE6-DC91E8CE5364}
TERM.....SKIN WOUND CLOSURE SUTURES

DATA VALUE: SLEEPING

ABBREVIATION.....s
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Sleeping
HELP_TEXT.....No Help Text
ID.....{87338420-F351-F98C-7A80-4990B87ACE2B}
TERM.....SLEEPING

DATA VALUE: SLEEPS POORLY

ABBREVIATION.....sp
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Sleeps Poorly
HELP_TEXT.....No Help Text
ID.....{AB9424B4-D0B1-148B-C800-FF535F3FAC84}
TERM.....SLEEPS POORLY

DATA VALUE: SLEPT WELL

ABBREVIATION.....SLEPT WELL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....States Has Slept Well
HELP_TEXT.....No Help Text
ID.....{18573E82-568F-4D2E-B9E1-6E4589EBEE2D}
TERM.....SLEPT WELL

DATA VALUE: SLIGHT

ABBREVIATION.....SLIGHT
ACTIVE.....Yes
DESCRIPTION.....Slight
DISPLAY_NAME.....Slight
HELP_TEXT.....Slight
ID.....{805549C1-5DC5-4A2D-BC1F-046E44928843}
TERM.....SLIGHT

DATA VALUE: SLING TRACTION

ABBREVIATION.....SLING TRAC
ACTIVE.....Yes
DESCRIPTION.....Sling Traction
DISPLAY_NAME.....Sling Traction
HELP_TEXT.....Sling Traction
ID.....{304E924C-C95A-4DBE-B30B-191BC4D26EF7}
TERM.....SLING TRACTION

DATA VALUE: SLOUGH

ABBREVIATION.....SLOUGH
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Slough
HELP_TEXT.....No Help Text
ID.....{548EB536-1AC4-4234-8C4D-DC54E4CAA6F6}
TERM.....SLOUGH

Page: 865

DATA VALUE: SLUGGISH

ABBREVIATION.....s
ACTIVE.....Yes
DESCRIPTION.....No Description

DISPLAY_NAME.....Sluggish
HELP_TEXT.....No Help Text
ID.....{3215A82E-8DCC-C7C0-8AE1-2D4C9D4A4EFB}
TERM.....SLUGGISH

DATA VALUE: SMALL

ABBREVIATION.....s
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Small
HELP_TEXT.....No Help Text
ID.....{5991F892-B28C-3998-E733-D8FFA26D7E21}
TERM.....SMALL

DATA VALUE: SMALL BORE FEEDING TUBE

ABBREVIATION.....SM BO TUBE
ACTIVE.....Yes
DESCRIPTION.....Small Bore Feeding Tube
DISPLAY_NAME.....Small Bore Feeding Tube
HELP_TEXT.....No Help Text
ID.....{7E347B8A-14E7-4093-9A5B-56E6AC0D78DF}
TERM.....SMALL BORE FEEDING TUBE

DATA VALUE: SMEAR

ABBREVIATION.....SMEAR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Smear
HELP_TEXT.....No Help Text
ID.....{4B99BBD3-443E-4954-8ADE-5FB3EE50540A}
TERM.....SMEAR

DATA VALUE: SMOKING AREA

ABBREVIATION.....sa
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Smoking Area
HELP_TEXT.....No Help Text
ID.....{8AFDBC0-9FC8-523F-5D07-4D8CCBA7036F}
TERM.....SMOKING AREA

Page: 866

DATA VALUE: SMOKING PORCH

ABBREVIATION.....sp
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Smoking Porch
HELP_TEXT.....No Help Text
ID.....{5AD3F05B-2A5B-E2B4-FD47-66712CFA44E6}

TERM.....SMOKING PORCH

DATA VALUE: SMOOTH

ABBREVIATION.....SMOOTH
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Smooth
HELP_TEXT.....No Help Text
ID.....{5D71A0D5-52ED-4B8D-B359-AC7D500660D8}
TERM.....SMOOTH

DATA VALUE: SNACK

ABBREVIATION.....SNACK
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Snack
HELP_TEXT.....No Help Text
ID.....{60B77974-4102-4511-A105-E24850B924D8}
TERM.....SNACK

DATA VALUE: SODIUM CONTROLLED 2 GM

ABBREVIATION.....SOD CONT 2
ACTIVE.....Yes
DESCRIPTION.....Sodium Controlled 2 GM
DISPLAY_NAME.....Sodium Controlled 2 GM
HELP_TEXT.....No Help Text
ID.....{8A1311AF-6C7D-407A-962F-0FAE6C341F9A}
TERM.....SODIUM CONTROLLED 2 GM

DATA VALUE: SODIUM CONTROLLED 4 GM

ABBREVIATION.....SOD CON 4G
ACTIVE.....Yes
DESCRIPTION.....Sodium Controlled 4 GM Diet
DISPLAY_NAME.....Sodium Controlled 4 GM
HELP_TEXT.....No Help Text
ID.....{26538A2F-5502-4D94-A2A7-D9E0BF825C47}
TERM.....SODIUM CONTROLLED 4 GM

Page: 867

DATA VALUE: SOFT

ABBREVIATION.....SOFT
ACTIVE.....Yes
DESCRIPTION.....Soft
DISPLAY_NAME.....Soft
HELP_TEXT.....No Help Text
ID.....{E68B4D4B-B86B-4004-A861-986645267F6E}
TERM.....SOFT

DATA VALUE: SPIRITUAL

ABBREVIATION.....SPIRITUAL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Spiritual
HELP_TEXT.....No Help Text
ID.....{9670E2B9-6938-4487-A857-25546775B7A2}
TERM.....SPIRITUAL

DATA VALUE: SPLINTING

ABBREVIATION.....SPLINTING
ACTIVE.....Yes
DESCRIPTION.....Splinting
DISPLAY_NAME.....Splinting
HELP_TEXT.....No Help Text
ID.....{D7FC5AA9-E63A-4A88-BE03-A50D0A8751F9}
TERM.....SPLINTING

DATA VALUE: SPONTANEOUS

ABBREVIATION.....SPONTANEOU
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Spontaneous
HELP_TEXT.....No Help Text
ID.....{C137C4A7-08A6-4064-8F95-4A4476E2E3A8}
TERM.....SPONTANEOUS

DATA VALUE: SPOONING

ABBREVIATION.....SPOONING
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Spoonning
HELP_TEXT.....No Help Text
ID.....{CA384C9F-D87D-4E2D-B4C0-89BD85692B16}
TERM.....SPOONING

Page: 868

DATA VALUE: SQ EMPHYSEMA

ABBREVIATION.....SQ EMPHYS.
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Subcutaneous Emphysema
HELP_TEXT.....No Help Text
ID.....{B8D03538-0168-4CDC-FF41-BA1185EBCC87}
TERM.....SQ EMPHYSEMA

DATA VALUE: SQUARE

ABBREVIATION.....SQUARE
ACTIVE.....Yes

DESCRIPTION.....Square
DISPLAY_NAME.....Square
HELP_TEXT.....Square
ID.....{63F70678-2CAF-48E8-964D-2805D8B9A43B}
TERM.....SQUARE

DATA VALUE: SQUEEZING

ABBREVIATION.....SQUEEZING
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....SQUEEZING
HELP_TEXT.....No Help Text
ID.....{A148587A-1B99-4D0D-B11E-87ADC3D2A811}
TERM.....SQUEEZING

DATA VALUE: STABBING

ABBREVIATION.....STABBING
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Stabbing
HELP_TEXT.....No Help Text
ID.....{684640D5-B3AB-47DA-98D8-5434ABF22B7E}
TERM.....STABBING

DATA VALUE: STABLE

ABBREVIATION.....STABLE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Stable
HELP_TEXT.....No Help Text
ID.....{A474E7BC-47EF-4472-A19A-172C441099EC}
TERM.....STABLE

Page: 869

DATA VALUE: STAFF

ABBREVIATION.....STAFF
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Staff
HELP_TEXT.....No Help Text
ID.....{BE8D678F-CE19-41FC-8C12-7A5DE7F0C7E2}
TERM.....STAFF

DATA VALUE: STAGE I

ABBREVIATION.....STAGE I
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Stage I
HELP_TEXT.....No Help Text

ID.....{B376AAAF-7143-4422-BCA2-CACF1201817C}
TERM.....STAGE I

DATA VALUE: STAGE II

ABBREVIATION.....STAGE II
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Stage II
HELP_TEXT.....No Help Text
ID.....{B520BEA3-841E-490A-A3CD-63D1027B561F}
TERM.....STAGE II

DATA VALUE: STAGE III

ABBREVIATION.....STAGE III
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Stage III
HELP_TEXT.....No Help Text
ID.....{A2F0866C-ECEE-4FCD-A34B-498232662D7E}
TERM.....STAGE III

DATA VALUE: STAGE IV

ABBREVIATION.....STAGE IV
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Stage IV
HELP_TEXT.....No Help Text
ID.....{284334AA-E488-4D45-93EF-9B93B2086A63}
TERM.....STAGE IV

Page: 870

DATA VALUE: STANDARD

ABBREVIATION.....STANDARD
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....STANDARD
HELP_TEXT.....No Help Text
ID.....{C4B8C06D-BFFE-4438-A43D-926F039E5EE1}
TERM.....STANDARD

DATA VALUE: STANDARD PRESSURE REDUCTION

ABBREVIATION.....STAN P RED
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Standard Pressure Reduction
HELP_TEXT.....No Help Text
ID.....{556537CD-A84D-4493-81F6-51CBD844B5B7}
TERM.....STANDARD PRESSURE REDUCTION

DATA VALUE: STAPLES

ABBREVIATION.....STAPLES
ACTIVE.....Yes
DESCRIPTION.....Staples to close a wound
DISPLAY_NAME.....Staples
HELP_TEXT.....No Help Text
ID.....{344DBBE0-5A94-4E18-BC97-BD7767A83E36}
TERM.....STAPLES

DATA VALUE: STARTED

ABBREVIATION.....STARTED
ACTIVE.....Yes
DESCRIPTION.....Started
DISPLAY_NAME.....Started
HELP_TEXT.....No Help Text
ID.....{D6E105C4-5EAD-424E-8E4C-195EB7783A3D}
TERM.....STARTED

DATA VALUE: STARTED ON LEFT

ABBREVIATION.....sol
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Started On Left
HELP_TEXT.....No Help Text
ID.....{AA6767EE-57F6-F565-76E9-EB219644492F}
TERM.....STARTED ON LEFT

Page: 871

DATA VALUE: STARTED ON RIGHT

ABBREVIATION.....sor
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Started On Right
HELP_TEXT.....No Help Text
ID.....{0C690E36-E710-244D-44FC-5DB5BCC0FF7E}
TERM.....STARTED ON RIGHT

DATA VALUE: STEADY GATE OR BASELINE

ABBREVIATION.....ST GA BA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Steady Gate or Baseline
HELP_TEXT.....No Help Text
ID.....{0D47B2CE-88A7-4FE5-964B-7A8107A5DB91}
TERM.....STEADY GATE OR BASELINE

DATA VALUE: STEP DOWN

ABBREVIATION.....STEP DOW

ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Step Down
HELP_TEXT.....No Help Text
ID.....{F6AC5CA7-0953-4E21-B1A4-6BFE1983D515}
TERM.....STEP DOWN

DATA VALUE: STINGING

ABBREVIATION.....STINGING
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Stinging
HELP_TEXT.....No Help Text
ID.....{DF235170-F501-4E8D-89DD-5A3C76F4B13B}
TERM.....STINGING

DATA VALUE: STRAIGHT CATHETER

ABBREVIATION.....STR. CATH
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Straight Catheter
HELP_TEXT.....No Help Text
ID.....{5F4A8A7F-E720-9458-2360-1DC14E7D05F3}
TERM.....STRAIGHT CATHETER

Page: 872

DATA VALUE: STRETCHER

ABBREVIATION.....STRETCHER
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Stretcher
HELP_TEXT.....No Help Text
ID.....{93CB521A-3A67-41FB-A59C-3DD54D810B11}
TERM.....STRETCHER

DATA VALUE: STRICT

ABBREVIATION.....STRICT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Strict
HELP_TEXT.....No Help Text
ID.....{5717A3EA-9A1F-44AB-BB0B-E86BBDED3ECF}
TERM.....STRICT

DATA VALUE: STRONG

ABBREVIATION.....s
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Strong

HELP_TEXT.....No Help Text
ID.....{A12154E1-544E-F372-3839-0D44E5BC476B}
TERM.....STRONG

DATA VALUE: STRYKER FRAME

ABBREVIATION.....STRYKER FR
ACTIVE.....Yes
DESCRIPTION.....Stryker Fram
DISPLAY_NAME.....Stryker Frame
HELP_TEXT.....No Help Text
ID.....{71A1DB72-C860-42D8-BA4F-B71AFA2164B1}
TERM.....STRYKER FRAME

DATA VALUE: STUDENT

ABBREVIATION.....STUDENT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Student
HELP_TEXT.....No Help Text
ID.....{613AFF67-9121-4F6C-A44F-1ACFEF10E612}
TERM.....STUDENT

Page: 873

DATA VALUE: SUBCUTANEOUS TISSUE

ABBREVIATION.....SUBCUTAN T
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Subcutaneous Tissue
HELP_TEXT.....No Help Text
ID.....{D1345ED8-2CDB-4FD7-84B8-92C9E0872543}
TERM.....SUBCUTANEOUS TISSUE

DATA VALUE: SUCTIONED ONLY

ABBREVIATION.....SUCT ONLY
ACTIVE.....Yes
DESCRIPTION.....Suctioned Only
DISPLAY_NAME.....Suctioned Only
HELP_TEXT.....No Help Text
ID.....{7A8B23FD-69C5-4931-A6D4-E92F8FD70817}
TERM.....SUCTIONED ONLY

DATA VALUE: SUICIDAL

ABBREVIATION.....s
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Suicidal
HELP_TEXT.....No Help Text
ID.....{B104E9C7-A78E-82DC-49A4-19B94A639AFE}
TERM.....SUICIDAL

DATA VALUE: SUMP

ABBREVIATION.....Sump
ACTIVE.....Yes
DESCRIPTION.....Sump Drain
DISPLAY_NAME.....Sump Drain
HELP_TEXT.....No Help Text
ID.....{DAD881B9-F1EE-FDAD-AEB0-68C754FDC24E}
TERM.....SUMP

DATA VALUE: SUPERVISION/SET-UP HELP

ABBREVIATION.....SUPERVISED
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Supervision/Set-Up Help
HELP_TEXT.....No Help Text
ID.....{596D028B-B60D-4244-A1B0-882177C8D25F}
TERM.....SUPERVISION/SET-UP HELP

Page: 874

DATA VALUE: SUPINE

ABBREVIATION.....SUPINE
ACTIVE.....Yes
DESCRIPTION.....Supine
DISPLAY_NAME.....Supine
HELP_TEXT.....No Help Text
ID.....{A18F9534-17AF-47D7-8D14-2AA457A889BD}
TERM.....SUPINE

DATA VALUE: SUPRAPUBIC

ABBREVIATION.....SUPRAPUBIC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Suprapubic
HELP_TEXT.....No Help Text
ID.....{69E0176F-52FC-A3C4-7A4E-AD57DBD6E572}
TERM.....SUPRAPUBIC

DATA VALUE: SUPRAVENTRICULAR TACHYCARDIA

ABBREVIATION.....st
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Supraventricular Tachycardia
HELP_TEXT.....No Help Text
ID.....{F712412F-D39A-58F7-E4F1-6AC8B77445D7}
TERM.....SUPRAVENTRICULAR TACHYCARDIA

DATA VALUE: SURGICAL

ABBREVIATION.....SURGICAL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Surgical
HELP_TEXT.....No Help-Text
ID.....{2C6CFEF1-6122-45BA-A9F8-F33B5B3C0E90}
TERM.....SURGICAL

DATA VALUE: SURGICAL DONOR

ABBREVIATION.....SURG DONOR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Surgical Donor
HELP_TEXT.....No Help Text
ID.....{F3EE10EA-5776-479B-A58C-F59D6526F90B}
TERM.....SURGICAL DONOR

Page: 875

DATA VALUE: SURGICAL GRAFT

ABBREVIATION.....SURG GRAFT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Surgical Graft
HELP_TEXT.....No Help Text
ID.....{328BA3EA-409B-4C5E-92E7-C6ECE8F983E1}
TERM.....SURGICAL GRAFT

DATA VALUE: SUSPECTED DEEP TISSUE INJURY

ABBREVIATION.....SUS DEEP T
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Suspected Deep Tissue Injury
HELP_TEXT.....No Help Text
ID.....{E82B4082-646F-4B5A-B68E-BB0872A4F6D2}
TERM.....SUSPECTED DEEP TISSUE INJURY

DATA VALUE: SUTURES

ABBREVIATION.....SUTURES
ACTIVE.....Yes
DESCRIPTION.....Sutures to close wound
DISPLAY_NAME.....Sutures
HELP_TEXT.....No Help Text
ID.....{F6E3535E-2368-4014-8F4E-670EC0CB69A5}
TERM.....SUTURES

DATA VALUE: SWOLLEN

ABBREVIATION.....SWOLLEN
ACTIVE.....Yes
DESCRIPTION.....No Description

DISPLAY_NAME.....Swollen
HELP_TEXT.....No Help Text
ID.....{DCE86B10-C393-8DA9-6518-F2943F0E1A4E}
TERM.....SWOLLEN

DATA VALUE: SYNCHRONOUS INTERMITTENT MANDATORY VENTILATION

ABBREVIATION.....SIMV
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Synchronous Intermittent Mandatory
Ventilation
HELP_TEXT.....No Help Text
ID.....{A606BB2B-5C00-FEA1-127F-1D3455AE5996}
TERM.....SYNCHRONOUS INTERMITTENT MANDATORY
VENTILATION

Page: 876

DATA VALUE: SYRINGE CHANGED

ABBREVIATION.....SYR CHANGE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Syringe Changed
HELP_TEXT.....No Help Text
ID.....{BE84F060-7E03-49DD-A6BC-787375C8E1F1}
TERM.....SYRINGE CHANGED

DATA VALUE: SYSTOLIC

ABBREVIATION.....SYSTOLIC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Systolic
HELP_TEXT.....No Help Text
ID.....{859B5700-0990-42FF-8129-8DB24EFDA8D1}
TERM.....SYSTOLIC

DATA VALUE: TAN

ABBREVIATION.....TAN
ACTIVE.....Yes
DESCRIPTION.....The color tan
DISPLAY_NAME.....Tan
HELP_TEXT.....The color tan
ID.....{AB994ED7-8916-FAE0-0E97-CD4E0EF80ED1}
TERM.....TAN

DATA VALUE: TAPE BURN

ABBREVIATION.....TAPE BURN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Tape Burn

HELP_TEXT.....No Help Text
ID.....{DDD75AFD-DA8D-4E73-873F-86D4DF4AF47A}
TERM.....TAPE BURN

DATA VALUE: TARRY

ABBREVIATION.....TARRY
ACTIVE.....Yes
DESCRIPTION.....Tarry in appearance
DISPLAY_NAME.....Tarry
HELP_TEXT.....No Help Text
ID.....{09364271-6307-471F-A334-A8F7E051EF46}
TERM.....TARRY

Page: 877

DATA VALUE: TEA

ABBREVIATION.....TEA
ACTIVE.....Yes
DESCRIPTION.....Tea
DISPLAY_NAME.....Tea
HELP_TEXT.....Tea
ID.....{E2C0C287-DFDC-4F7B-8A67-7A339B8E8669}
TERM.....TEA

DATA VALUE: TEETH BRUSHED

ABBREVIATION.....TEETH BRUS
ACTIVE.....Yes
DESCRIPTION.....Teeth Brushed
DISPLAY_NAME.....Teeth Brushed
HELP_TEXT.....No Help Text
ID.....{20841794-6E6D-4262-B1FE-9E3F20F8CC49}
TERM.....TEETH BRUSHED

DATA VALUE: TEMPORARY (TRANSVENOUS)

ABBREVIATION.....tt
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Temporary (Transvenous)
HELP_TEXT.....No Help Text
ID.....{F6A74188-3D79-D749-540F-4EA32201B0E6}
TERM.....TEMPORARY (TRANSVENOUS)

DATA VALUE: TENACIOUS

ABBREVIATION.....TENACIOUS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Tenacious
HELP_TEXT.....No Help Text
ID.....{EE90B04C-F74B-43F8-F3F2-8C2D4359D0B3}
TERM.....TENACIOUS

DATA VALUE: TENDER

ABBREVIATION.....TENDER
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Tender
HELP_TEXT.....No Help Text
ID.....{30486CC7-3097-4F95-9238-FD6C2126CD10}
TERM.....TENDER

Page: 878

DATA VALUE: TENDON

ABBREVIATION.....TENDON
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Tendon
HELP_TEXT.....No Help Text
ID.....{40628D33-E217-4D39-9A82-4C8D7702576F}
TERM.....TENDON

DATA VALUE: THICK

ABBREVIATION.....THICK
ACTIVE.....Yes
DESCRIPTION.....Thick
DISPLAY_NAME.....Thick
HELP_TEXT.....No Help Text
ID.....{1554A49E-1468-46B2-F851-531DDB1D1A98}
TERM.....THICK

DATA VALUE: THIGH

ABBREVIATION.....THIGH
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Thigh
HELP_TEXT.....No Help Text
ID.....{65AE788B-52F5-44D3-B9C1-49DAE46FF5AB}
TERM.....THIGH

DATA VALUE: THIN

ABBREVIATION.....THIN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Thin
HELP_TEXT.....No Help Text
ID.....{CB8D8B5C-E824-D030-1241-D8EB814671E9}
TERM.....THIN

DATA VALUE: THIRD DEGREE AV BLOCK

ABBREVIATION.....tdab
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Third Degree Av Block
HELP_TEXT.....No Help Text
ID.....{1AC7E31B-830E-2309-3AF7-8FB03003DD86}
TERM.....THIRD DEGREE AV BLOCK

Page: 879

DATA VALUE: THORACIC SPINAL PRECAUTIONS

ABBREVIATION.....THOR S P
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Thoriacic Spinal Precautions
HELP_TEXT.....No Help Text
ID.....{D55AEF23-7F5E-45DD-A178-F217E6AF3873}
TERM.....THORACIC SPINAL PRECAUTIONS

DATA VALUE: THROBBING

ABBREVIATION.....THROBBING
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Throbbing
HELP_TEXT.....No Help Text
ID.....{3ACE2B10-1418-47F8-A4F3-0C41A699B5D9}
TERM.....THROBBING

DATA VALUE: TINGLING

ABBREVIATION.....t
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Tingling
HELP_TEXT.....No Help Text
ID.....{AE8E3101-AC3D-DEC3-0293-58F901F82EB8}
TERM.....TINGLING

DATA VALUE: TO PAIN

ABBREVIATION.....TO PAIN
ACTIVE.....Yes
DESCRIPTION.....To Pain
DISPLAY_NAME.....To Pain
HELP_TEXT.....No Help Text
ID.....{AA1F818D-DCFC-4FBC-97D3-BE86D309119E}
TERM.....TO PAIN

DATA VALUE: TO SPEECH

ABBREVIATION.....TO SPEECH
ACTIVE.....Yes
DESCRIPTION.....To Speech

DISPLAY_NAME.....To Speech
HELP_TEXT.....No Help Text
ID.....{2DDD94FA-D4D0-48B0-B4F4-F25C42E4AF96}
TERM.....TO SPEECH

Page: 880

DATA VALUE: TOENAIL CARE

ABBREVIATION.....TOEN CARE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Toenail Care
HELP_TEXT.....No Help Text
ID.....{6CF3ABCE-FF74-43E2-81E5-D23339017CA9}
TERM.....TOENAIL CARE

DATA VALUE: TOLERATING FLUIDS/NOURISHMENT

ABBREVIATION.....TO OR FL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Tolerating Fluids/Nourishment
HELP_TEXT.....No Help Text
ID.....{E107F272-CD27-46B8-8C50-F9FF776DDCD9}
TERM.....TOLERATING FLUIDS/NOURISHMENT

DATA VALUE: TONIC/CLONIC

ABBREVIATION.....TON/CLON
ACTIVE.....Yes
DESCRIPTION.....Tonic/Clonic
DISPLAY_NAME.....Tonic/Clonic
HELP_TEXT.....Tonic/Clonic
ID.....{81932BAA-6C79-47A3-85B3-592824B6C09C}
TERM.....TONIC/CLONIC

DATA VALUE: TOPICAL TREATMENT

ABBREVIATION.....TOPICAL T
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Topical Treatment
HELP_TEXT.....No Help Text
ID.....{4E995342-8A6F-4B21-8013-F8C6A9FC7658}
TERM.....TOPICAL TREATMENT

DATA VALUE: TOTAL DEPENDENCE

ABBREVIATION.....TOT DEPEND
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Total Dependence
HELP_TEXT.....No Help Text
ID.....{38702CDD-E440-4EFA-9583-36EEE719F7C3}

TERM.....TOTAL DEPENDENCE

Page: 881

DATA VALUE: TRACHEOSTOMY

ABBREVIATION.....TRACHEOSTO
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Tracheostomy
HELP_TEXT.....No Help Text
ID.....{2B0B7599-847E-852F-D16F-091CAA4B205}
TERM.....TRACHEOSTOMY

DATA VALUE: TRANSCUTANEOUS

ABBREVIATION.....t
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Transcutaneous
HELP_TEXT.....No Help Text
ID.....{F4315F92-8804-1CB4-4D03-CCFDC95C6E71}
TERM.....TRANSCUTANEOUS

DATA VALUE: TRANSFER

ABBREVIATION.....TRANSFER
ACTIVE.....Yes
DESCRIPTION.....Transfer
DISPLAY_NAME.....Transfer
HELP_TEXT.....No Help Text
ID.....{DFDB4783-25BE-485B-A9BA-621F1004AAF2}
TERM.....TRANSFER

DATA VALUE: TRANSFER WITH ASSIST

ABBREVIATION.....TRANS ASST
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Transfer With Assist
HELP_TEXT.....No Help Text
ID.....{62A5BFF5-9C7C-4412-83A8-F872095DE005}
TERM.....TRANSFER WITH ASSIST

DATA VALUE: TRANSPARENT

ABBREVIATION.....TRANSPAREN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Transparent
HELP_TEXT.....No Help Text
ID.....{0E73D24D-0EFC-08F9-196A-D30DD1BB3228}
TERM.....TRANSPARENT

DATA VALUE: TRANSPARENT DRESSING

ABBREVIATION.....TRANSP DRE
ACTIVE.....Yes
DESCRIPTION.....Transparent Dressing
DISPLAY_NAME.....Transparent Dressing
HELP_TEXT.....No Help Text
ID.....{777C7722-AEA8-4A0B-9ADD-1157D37AE539}
TERM.....TRANSPARENT DRESSING

DATA VALUE: TRANSPORT FACILITATOR PRESENT

ABBREVIATION.....TR FA PR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Transport Facilitator Present
HELP_TEXT.....No Help Text
ID.....{645BF269-6B23-4377-95ED-71291F075B25}
TERM.....TRANSPORT FACILITATOR PRESENT

DATA VALUE: TRANSPYLORIC

ABBREVIATION.....TRANSPYOLO
ACTIVE.....Yes
DESCRIPTION.....Transpyloric
DISPLAY_NAME.....Transpyloric
HELP_TEXT.....No Help Text
ID.....{78FA5E0A-EA21-40F3-8E8D-3B3FCC1ED060}
TERM.....TRANSPYLORIC

DATA VALUE: TREATMENTS

ABBREVIATION.....TREATMENTS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Treatments
HELP_TEXT.....No Help Text
ID.....{7A6879B5-F781-498E-9D71-3FA0243C3F21}
TERM.....TREATMENTS

DATA VALUE: TRIPLE

ABBREVIATION.....TRIPLE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Triple
HELP_TEXT.....No Help Text
ID.....{077BE8CF-F2AB-56E1-1F25-51DF80F8054D}
TERM.....TRIPLE

DATA VALUE: TUBED/DYSPHAGIC

ABBREVIATION.....TUBE/DYSPH
ACTIVE.....Yes
DESCRIPTION.....Tubed/Dysphagic
DISPLAY_NAME.....Tubed/Dysphagic
HELP_TEXT.....No Help Text
ID.....{82AAFAE3-50B0-4F2A-91CF-7900849A14EB}
TERM.....TUBED/DYSPHAGIC

DATA VALUE: TUBING CHANGED/LABELED

ABBREVIATION.....TUBE CHANG
ACTIVE.....Yes
DESCRIPTION.....Changed or labeled tubing
DISPLAY_NAME.....Tubing Changed / Labeled
HELP_TEXT.....No Help Text
ID.....{19BA337C-6E28-411E-9702-48B262C593FA}
TERM.....TUBING CHANGED/LABELED

DATA VALUE: TURN AND POSITION

ABBREVIATION.....TURN POSIT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Turn and Position
HELP_TEXT.....No Help Text
ID.....{44A4F470-1760-4073-AB3A-D59CA9EFE7F3}
TERM.....TURN AND POSITION

DATA VALUE: TURNING SELF INDEPENDENTLY

ABBREVIATION.....TURN SELF
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Turning Self Independently
HELP_TEXT.....No Help Text
ID.....{DD804CEB-F7CA-4D42-BF6A-447E38B55E98}
TERM.....TURNING SELF INDEPENDENTLY

DATA VALUE: TWILL TAPE

ABBREVIATION.....Twill Tape
ACTIVE.....Yes
DESCRIPTION.....Twill Tape
DISPLAY_NAME.....Twill Tape
HELP_TEXT.....No Help Text
ID.....{C52BE250-732D-03BB-B8F9-72BE4D546807}
TERM.....TWILL TAPE

Page: 884

DATA VALUE: TWO PERSON ASSIST

ABBREVIATION.....2 PERS AST
ACTIVE.....Yes

DESCRIPTION.....No Description
DISPLAY_NAME.....Two Person Assist
HELP_TEXT.....No Help Text
ID.....{0547070B-C9BC-4406-BDC0-BF9C6C23D41A}
TERM.....TWO PERSON ASSIST

DATA VALUE: ULTRACAL

ABBREVIATION.....ULTRACAL
ACTIVE.....Yes
DESCRIPTION.....Enteral Nutrition Type
DISPLAY_NAME.....Ultracal
HELP_TEXT.....No Help Text
ID.....{179AE9A5-154D-4CD2-B70B-66E84F5D0932}
TERM.....ULTRACAL

DATA VALUE: UNABLE TO ASSESS

ABBREVIATION.....UNA TO AS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Unable To Assess
HELP_TEXT.....No Help Text
ID.....{4D11A5E0-1CE8-4390-9EDD-ECEACCAA991B}
TERM.....UNABLE TO ASSESS

DATA VALUE: UNABLE TO INFORM

ABBREVIATION.....UNABLE INF
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Unable to Inform
HELP_TEXT.....No Help Text
ID.....{77428C87-B3F9-4BBD-B781-F7D89BCF9492}
TERM.....UNABLE TO INFORM

DATA VALUE: UNABLE TO READ

ABBREVIATION.....ILLITERATE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Unable to Read
HELP_TEXT.....No Help Text
ID.....{6D0C2C48-BDD2-489B-94C8-49E6DAC93F2E}
TERM.....UNABLE TO READ

Page: 885

DATA VALUE: UNABLE TO RESPOND

ABBREVIATION.....UNA TO RES
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Unable To Respond
HELP_TEXT.....No Help Text

ID.....{1EF29FF5-36AB-4FC0-B5E7-F8C3FC95BA5A}
TERM.....UNABLE TO RESPOND

DATA VALUE: UNABLE TO WEDGE

ABBREVIATION.....UNABLE WED
ACTIVE.....Yes
DESCRIPTION.....Unable to Wedge
DISPLAY_NAME.....Unable to Wedge
HELP_TEXT.....No Help Text
ID.....{5D35F546-F1EA-49C6-8B1D-22D24526AF8B}
TERM.....UNABLE TO WEDGE

DATA VALUE: UNBEARABLE

ABBREVIATION.....UNBEARABLE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Unbearable
HELP_TEXT.....No Help Text
ID.....{96611993-45F3-4D1D-BF7F-5E5C46E5569E}
TERM.....UNBEARABLE

DATA VALUE: UNCOOPERATIVE

ABBREVIATION.....u
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Uncooperative
HELP_TEXT.....No Help Text
ID.....{2D737EFB-F4CA-7AEA-9EB2-6FF133742815}
TERM.....UNCOOPERATIVE

DATA VALUE: UNEQUAL

ABBREVIATION.....Unequal
ACTIVE.....Yes
DESCRIPTION.....Not the same
DISPLAY_NAME.....Unequal
HELP_TEXT.....Not the same
ID.....{73319200-164F-B39F-649C-B7FB1AFBDBB8}
TERM.....UNEQUAL

DATA VALUE: UNKNOWN

ABBREVIATION.....u
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Unknown
HELP_TEXT.....No Help Text
ID.....{1F872397-8635-18A3-E252-AE9B20B88C7D}
TERM.....UNKNOWN

DATA VALUE: UNSTABLE

ABBREVIATION.....UNSTABLE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Unstable
HELP_TEXT.....No Help Text
ID.....{5841C521-8070-4034-A69D-27DAB2990EB3}
TERM.....UNSTABLE

DATA VALUE: UNSTAGEABLE

ABBREVIATION.....UNSTAGEABL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Unstageable
HELP_TEXT.....No Help Text
ID.....{5C507F64-2EC0-4FD8-A786-94D76D90C0D1}
TERM.....UNSTAGEABLE

DATA VALUE: UNSUCCESSFUL

ABBREVIATION.....UNSUCCESSF
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Unsuccessful
HELP_TEXT.....No Help Text
ID.....{54DC58A1-7585-4A32-B810-90B586ED9332}
TERM.....UNSUCCESSFUL

DATA VALUE: UP

ABBREVIATION.....u
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Up
HELP_TEXT.....No Help Text
ID.....{7BCC74F9-FE2F-FE60-57BA-321BFA316BEC}
TERM.....UP

Page: 887

DATA VALUE: UP TO CHAIR

ABBREVIATION.....UP TO CHAI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Up To Chair
HELP_TEXT.....No Help Text
ID.....{42E598E5-82EC-428B-B5FA-1129E3095C78}
TERM.....UP TO CHAIR

DATA VALUE: UROMETER

ABBREVIATION.....u

ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Urometer
HELP_TEXT.....No Help Text
ID.....{13074A3A-19D7-300C-7BFF-8DC18DDE2951}
TERM.....UROMETER

DATA VALUE: UROSTOMY

ABBREVIATION.....UROSTOMY
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Urostomy
HELP_TEXT.....No Help Text
ID.....{1F9E528D-6D78-54EE-94E6-649C6F2056A3}
TERM.....UROSTOMY

DATA VALUE: V1

ABBREVIATION.....V1
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....V1
HELP_TEXT.....No Help Text
ID.....{4CAA3A92-F3DB-4847-8523-380EFEE6CFD5}
TERM.....V1

DATA VALUE: V2

ABBREVIATION.....V2
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....V2
HELP_TEXT.....No Help Text
ID.....{05A3232E-A30D-4771-A1B5-E7705B717646}
TERM.....V2

Page: 888

DATA VALUE: V3

ABBREVIATION.....V3
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....V3
HELP_TEXT.....No Help Text
ID.....{E150126E-4AA7-4F96-961E-FD1169C77158}
TERM.....V3

DATA VALUE: V4

ABBREVIATION.....V4
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....V4

HELP_TEXT.....No Help Text
ID.....{255260F9-5CC7-41B1-831C-F9319FA27506}
TERM.....V4

DATA VALUE: V5

ABBREVIATION.....V5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....V5
HELP_TEXT.....No Help Text
ID.....{D5C10744-0EB5-46FA-AE6C-CCF571318F8A}
TERM.....V5

DATA VALUE: V6

ABBREVIATION.....V6
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....V6
HELP_TEXT.....No Help Text
ID.....{B86F72D7-7884-4470-92A8-EF9277C5B584}
TERM.....V6

DATA VALUE: VACUUM ASSISTED CLOSURE

ABBREVIATION.....VAC ASS CL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Vacuum Assisted Closure
HELP_TEXT.....No Help Text
ID.....{9C3DD561-FECF-47FA-BBA5-4C08FF9A4CA1}
TERM.....VACUUM ASSISTED CLOSURE

Page: 889

DATA VALUE: VASCULAR

ABBREVIATION.....VASCULAR
ACTIVE.....Yes
DESCRIPTION.....Vascular
DISPLAY_NAME.....Vascular
HELP_TEXT.....No Help Text
ID.....{60523BF8-6372-4760-881C-A62F85D2A0F6}
TERM.....VASCULAR

DATA VALUE: VASCULAR LESION ARTERIAL

ABBREVIATION.....VASLES ART
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Vascular Leasion Arterial
HELP_TEXT.....No Help Text
ID.....{EA565D3F-D984-4553-B597-DFE844554CE9}
TERM.....VASCULAR LESION ARTERIAL

DATA VALUE: VASCULAR LESION MIXED

ABBREVIATION.....VASLES MIX
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Vascular Leasion Mixed
HELP_TEXT.....No Help Text
ID.....{8684AB2D-9162-4DC2-A5FA-A14558678689}
TERM.....VASCULAR LESION MIXED

DATA VALUE: VASCULAR LESION VENOUS

ABBREVIATION.....VASLES VEN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Vascular Lesion Venous
HELP_TEXT.....No Help Text
ID.....{3F8CEC69-5E13-4F26-AB86-34FE2805E3C6}
TERM.....VASCULAR LESION VENOUS

DATA VALUE: VASCULAR/NEUROPATHIC

ABBREVIATION.....VASC/NEURO
ACTIVE.....Yes
DESCRIPTION.....Vascular / Neuropathic
DISPLAY_NAME.....Vascular / Neuropathic
HELP_TEXT.....No Help Text
ID.....{6B91CEED-4863-46EA-A9AE-3B0D522F1BDC}
TERM.....VASCULAR/NEUROPATHIC

Page: 890

DATA VALUE: VASOPRESSOR

ABBREVIATION.....VASOPRESS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Vasopressor
HELP_TEXT.....No Help Text
ID.....{591FCB09-1287-48C6-8007-4CE961296BA9}
TERM.....VASOPRESSOR

DATA VALUE: VDD

ABBREVIATION.....VDD
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....VDD
HELP_TEXT.....No Help Text
ID.....{E733E525-EF36-2F93-AECC-BAD1C3E2C633}
TERM.....VDD

DATA VALUE: VEGETARIAN DIET

ABBREVIATION.....VEG DIET
ACTIVE.....Yes
DESCRIPTION.....Vegetarian Diet
DISPLAY_NAME.....Vegetarian Diet
HELP_TEXT.....No Help Text
ID.....{6B8143BE-843F-4944-9ABB-7175AFE2DF17}
TERM.....VEGETARIAN DIET

DATA VALUE: VENTICULAR FIBRILLATION

ABBREVIATION.....vf
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Ventricular Fibrillation
HELP_TEXT.....No Help Text
ID.....{400B70C2-6E15-458F-70C6-FA2597DE2556}
TERM.....VENTICULAR FIBRILLATION

DATA VALUE: VENTILATOR

ABBREVIATION.....v
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Ventilator
HELP_TEXT.....No Help Text
ID.....{A5A647ED-A56D-5012-AD5B-788B3DCB3498}
TERM.....VENTILATOR

Page: 891

DATA VALUE: VENTILATOR ACHIEVA P5O2

ABBREVIATION.....ACHIEVA P5
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Achieva P5O2
HELP_TEXT.....No Help Text
ID.....{1690EEA5-3C3B-4D30-AC5C-214C94626699}
TERM.....VENTILATOR ACHIEVA P5O2

DATA VALUE: VENTILATOR ADULT STAR

ABBREVIATION.....ADULT STAR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Adult Star Ventilator
HELP_TEXT.....No Help Text
ID.....{879B13DA-4874-42D4-99EC-8E4ADB31B14E}
TERM.....VENTILATOR ADULT STAR

DATA VALUE: VENTILATOR BEAR 1000

ABBREVIATION.....BEAR 1000
ACTIVE.....Yes
DESCRIPTION.....No Description

DISPLAY_NAME.....Bear 1000
HELP_TEXT.....No Help Text
ID.....{03C30F99-A27C-427C-B0ED-872F0DE646BD}
TERM.....VENTILATOR BEAR 1000

DATA VALUE: VENTILATOR DRAEGER

ABBREVIATION.....DRAEGER
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Draeger
HELP_TEXT.....No Help Text
ID.....{DD18C42F-EF13-460B-8676-9ED1E429B54B}
TERM.....VENTILATOR DRAEGER

DATA VALUE: VENTILATOR OTHER

ABBREVIATION.....VENT OTHER
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Other Ventilator; Use Comments
HELP_TEXT.....No Help Text
ID.....{138B7B3D-9F8B-447A-9785-12AB396AE6DE}
TERM.....VENTILATOR OTHER

Page: 892

DATA VALUE: VENTILATOR PB 7200ae

ABBREVIATION.....PB7200ae
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....PB 7200ae
HELP_TEXT.....No Help Text
ID.....{86AE485F-ECD1-4B8E-820F-71CF159F318B}
TERM.....VENTILATOR PB 7200ae

DATA VALUE: VENTILATOR PB 840

ABBREVIATION.....PB 840
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....PB 840
HELP_TEXT.....No Help Text
ID.....{25A7B0DA-E437-4201-B8F3-3663B1B91C1A}
TERM.....VENTILATOR PB 840

DATA VALUE: VENTILATOR SERVO

ABBREVIATION.....SERVO
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Servo
HELP_TEXT.....No Help Text
ID.....{D6F06715-1ADD-4941-B026-B9A5E463943E}

TERM.....VENTILATOR SERVO

DATA VALUE: VENTRICULAR

ABBREVIATION.....v
ACTIVE.....Yes
DESCRIPTION.....VENTRICULAR
DISPLAY_NAME.....Ventricular
HELP_TEXT.....No Help Text
ID.....{307A4572-3CF0-7875-7313-7FD4AB23F8E7}
TERM.....VENTRICULAR

DATA VALUE: VENTRICULAR TACHYCARDIA

ABBREVIATION.....vt
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Ventricular Tachycardia
HELP_TEXT.....No Help Text
ID.....{138258AD-3A2C-2A6C-9DF8-615CCC303DC2}
TERM.....VENTRICULAR TACHYCARDIA

Page: 893

DATA VALUE: VENTRICULAR WHIP PRESENT

ABBREVIATION.....VENT WHIP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Ventricular Whip Present
HELP_TEXT.....No Help Text
ID.....{15CB9645-690D-446B-A380-F14CEEC9C6AA}
TERM.....VENTRICULAR WHIP PRESENT

DATA VALUE: VENTRICULOSTOMY

ABBREVIATION.....Ventriculo
ACTIVE.....Yes
DESCRIPTION.....Ventriculostomy
DISPLAY_NAME.....Ventriculostomy
HELP_TEXT.....No Help Text
ID.....{ECF20317-4708-7661-D4C2-9BDB4CAAC4FA}
TERM.....VENTRICULOSTOMY

DATA VALUE: VERBAL DE-ESCALATION

ABBREVIATION.....VERB DE-ES
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Verbal De-escalation
HELP_TEXT.....No Help Text
ID.....{8AB46F8C-B4D5-4800-BCF7-52CA3F32E99B}
TERM.....VERBAL DE-ESCALATION

DATA VALUE: VERY LOW BMI

ABBREVIATION.....LOW BMI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Very Low BMI
HELP_TEXT.....No Help Text
ID.....{78B3A37B-938B-42E5-A403-75E206722BC3}
TERM.....VERY LOW BMI

DATA VALUE: VIOLET

ABBREVIATION.....VIOLET
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Violet
HELP_TEXT.....No Help Text
ID.....{1249AE97-CF1F-4604-A17A-5FC85F6D0185}
TERM.....VIOLET

Page: 894

DATA VALUE: VIP PORT

ABBREVIATION.....VIP PORT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....VIP Port Catheter
HELP_TEXT.....No Help Text
ID.....{3D0426B3-F8F3-BAA4-71C7-8A09F911EF08}
TERM.....VIP PORT

DATA VALUE: VISCOUS

ABBREVIATION.....VISCOUS
ACTIVE.....Yes
DESCRIPTION.....Viscous
DISPLAY_NAME.....Viscous
HELP_TEXT.....No Help Text
ID.....{0E8019D5-F2DB-46B8-B750-41568DD0179C}
TERM.....VISCOUS

DATA VALUE: VISUAL IMPAIRMENT

ABBREVIATION.....VIS IMPAIR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Visual Impairment
HELP_TEXT.....No Help Text
ID.....{7B2107CE-002D-43BE-9102-FCD5B9FAD617}
TERM.....VISUAL IMPAIRMENT

DATA VALUE: VITAL ENTERAL NUTRITION TYPE

ABBREVIATION.....VITAL ENT
ACTIVE.....Yes

DESCRIPTION.....Vital Enteral Nutrition Type
DISPLAY_NAME.....Vital
HELP_TEXT.....No Help Text
ID.....{2DA036DA-BF27-4E52-86A2-20A312823480}
TERM.....VITAL ENTERAL NUTRITION TYPE

DATA VALUE: VITAL HN ENTERAL NUTRITION

ABBREVIATION.....VITAL HN
ACTIVE.....Yes
DESCRIPTION.....Vital HN Enteral Nutrition Type
DISPLAY_NAME.....Vital HN Enteral Nutrition Type
HELP_TEXT.....No Help Text
ID.....{E45849D4-05EA-4DC1-BD9C-3FF66DD4D292}
TERM.....VITAL HN ENTERAL NUTRITION

Page: 895

DATA VALUE: VITAL SIGNS STABLE

ABBREVIATION.....VI SI ST
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Vital Signs Stable
HELP_TEXT.....No Help Text
ID.....{E6ECD067-54A0-4214-A5F8-43D08C9FEAD6}
TERM.....VITAL SIGNS STABLE

DATA VALUE: VOIDED

ABBREVIATION.....VOIDED
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Voided
HELP_TEXT.....No Help Text
ID.....{7F23E601-A1E5-4016-8A12-06B1099041CF}
TERM.....VOIDED

DATA VALUE: VOLUNTEER

ABBREVIATION.....VOLUNTEE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Volunteer
HELP_TEXT.....No Help Text
ID.....{F5406412-DE48-4F84-8A31-0B1D79DB576D}
TERM.....VOLUNTEER

DATA VALUE: VOMITING

ABBREVIATION.....VOMITING
ACTIVE.....Yes
DESCRIPTION.....Vomiting
DISPLAY_NAME.....Vomiting
HELP_TEXT.....Vomiting

ID.....{921F463C-1639-4D13-8E4E-6D3A1B6927F9}
TERM.....VOMITING

DATA VALUE: VVI

ABBREVIATION.....VVI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....VVI
HELP_TEXT.....No Help Text
ID.....{2C87F8F3-2B99-80E7-04EB-E6FE88168627}
TERM.....VVI

Page: 896

DATA VALUE: WAFFLE

ABBREVIATION.....WAFFLE
ACTIVE.....Yes
DESCRIPTION.....Waffle
DISPLAY_NAME.....Waffle
HELP_TEXT.....No Help Text
ID.....{29A0FA73-CAA7-4660-9861-79A09F6DD98C}
TERM.....WAFFLE

DATA VALUE: WAITING ROOM

ABBREVIATION.....WAIT ROOM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Waiting Room
HELP_TEXT.....No Help Text
ID.....{13476486-F895-4FFC-AB42-B274827B47DA}
TERM.....WAITING ROOM

DATA VALUE: WALKER

ABBREVIATION.....WALKER
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Walker
HELP_TEXT.....No Help Text
ID.....{91247001-5E89-46A0-A841-B4FAE9226E63}
TERM.....WALKER

DATA VALUE: WALKING OR PACING

ABBREVIATION.....wop
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Walking Or Pacing
HELP_TEXT.....No Help Text
ID.....{CE7AABCB-84B7-A310-3A65-0AF17FB7B0E7}
TERM.....WALKING OR PACING

DATA VALUE: WANDERER

ABBREVIATION.....w
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Wanderer
HELP_TEXT.....No Help Text
ID.....{C5BD737B-F9A6-5800-B45C-2585507C3EDC}
TERM.....WANDERER

Page: 897

DATA VALUE: WARD

ABBREVIATION.....WARD
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Ward
HELP_TEXT.....No Help Text
ID.....{A16CB5CB-3401-4F6C-8467-7F78FC929872}
TERM.....WARD

DATA VALUE: WARM

ABBREVIATION.....WARM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Warm
HELP_TEXT.....No Help Text
ID.....{E03BBEB4-5643-968F-190F-309A1293FCF9}
TERM.....WARM

DATA VALUE: WARM HUMIDITY WITH O2

ABBREVIATION.....WM HUM O2
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Warm Humidity with Oxygen
HELP_TEXT.....No Help Text
ID.....{E62E07C3-4CA2-426B-8EBA-48C406837D9C}
TERM.....WARM HUMIDITY WITH O2

DATA VALUE: WARMER THAN GENERAL BODY TEMP

ABBREVIATION.....WARM B TEM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Warmer Than General Body Temperature
HELP_TEXT.....No Help Text
ID.....{AE43F2B4-B585-43F4-9E97-D61E225A803B}
TERM.....WARMER THAN GENERAL BODY TEMP

DATA VALUE: WARMER THAN OPPOSITE BODY PART

ABBREVIATION.....WARM OPP P

ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Warmer Than Opposite Body Part
HELP_TEXT.....No Help Text
ID.....{5DBF8614-5067-49C5-BD02-C1B21C9DFECC}
TERM.....WARMER THAN OPPOSITE BODY PART

Page: 898

DATA VALUE: WARMING BLANKET

ABBREVIATION.....WARM BLANK
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Warming Blanket
HELP_TEXT.....Warming Blanket
ID.....{23D20AE5-DE25-4456-8F77-BF0411FB19A6}
TERM.....WARMING BLANKET

DATA VALUE: WARMING BOTTLE

ABBREVIATION.....WARM BOTTL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Warming Bottle
HELP_TEXT.....No Help Text
ID.....{AD1BEFE5-4CE6-4076-8095-9BE7387CFE25}
TERM.....WARMING BOTTLE

DATA VALUE: WARMING OFF

ABBREVIATION.....WARMING OF
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Warming Off
HELP_TEXT.....No Help Text
ID.....{98FAC977-8E2E-406A-9217-A77E1DDCDF50}
TERM.....WARMING OFF

DATA VALUE: WARMING ON

ABBREVIATION.....WARMING ON
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Warming On
HELP_TEXT.....No Help Text
ID.....{6FFDE928-DE1E-4B25-BA3C-EB29E0CDEB35}
TERM.....WARMING ON

DATA VALUE: WATERSEAL

ABBREVIATION.....WATERSEAL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Waterseal

HELP_TEXT.....No Help Text
ID.....{6D91B99D-F199-09ED-78BE-B54F2CAC4E06}
TERM.....WATERSEAL

Page: 899

DATA VALUE: WAVEFORM REFLECTS ATRIAL PLACEMENT

ABBREVIATION.....WAVE ATRIA
ACTIVE.....Yes
DESCRIPTION.....Waveform Reflects Atrial Placement
DISPLAY_NAME.....Waveform Reflects Atrial Placement
HELP_TEXT.....No Help Text
ID.....{A1249745-41C2-4A4F-A86D-7E5E17477B06}
TERM.....WAVEFORM REFLECTS ATRIAL PLACEMENT

DATA VALUE: WEAK

ABBREVIATION.....w
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Weak
HELP_TEXT.....No Help Text
ID.....{BE2D9305-E2DA-F492-66CF-164EE58B8179}
TERM.....WEAK

DATA VALUE: WEDGED

ABBREVIATION.....WEDGED
ACTIVE.....Yes
DESCRIPTION.....Wedged
DISPLAY_NAME.....Wedged
HELP_TEXT.....No Help Text
ID.....{03D7F6FD-399D-4A6B-9900-59655D76EE95}
TERM.....WEDGED

DATA VALUE: WHEELCHAIR

ABBREVIATION.....WH CHAIR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Wheelchair
HELP_TEXT.....No Help Text
ID.....{EC1B82CE-4037-466F-90AC-52CD5C3B26BE}
TERM.....WHEELCHAIR

DATA VALUE: WHEELCHAIR EXIT ALARM

ABBREVIATION.....WHEEL EXIT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Wheelchair Exit Alarm
HELP_TEXT.....No Help Text
ID.....{7D06564E-ECDE-4A11-BEFE-92AB7E448AA3}
TERM.....WHEELCHAIR EXIT ALARM

DATA VALUE: WHEEZES

ABBREVIATION.....w
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Wheezes
HELP_TEXT.....No Help Text
ID.....{1B8A9729-554D-59F6-E9E8-40BDFF9C678E}
TERM.....WHEEZES

DATA VALUE: WHITE

ABBREVIATION.....WHITE
ACTIVE.....Yes
DESCRIPTION.....The color white
DISPLAY_NAME.....White
HELP_TEXT.....The color white
ID.....{AB04EF54-5246-EF3A-E1D9-02E9E24369DD}
TERM.....WHITE

DATA VALUE: WITH PHYSICIAN

ABBREVIATION.....wp
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....With Physician
HELP_TEXT.....No Help Text
ID.....{1673C53F-DA8C-2CEC-679B-B01D914FE533}
TERM.....WITH PHYSICIAN

DATA VALUE: WITH THERAPIST

ABBREVIATION.....wt
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....With Therapist
HELP_TEXT.....No Help Text
ID.....{795A1300-1D60-9D67-692D-A2AB2B190596}
TERM.....WITH THERAPIST

DATA VALUE: WITHDRAWN

ABBREVIATION.....w
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Withdrawn
HELP_TEXT.....No Help Text
ID.....{E74B0F04-0E3F-632F-9409-C4953D95B1BA}
TERM.....WITHDRAWN

DATA VALUE: WITHDRAWS TO PAIN

ABBREVIATION.....WITHD TO P
ACTIVE.....Yes
DESCRIPTION.....Withdraws To Pain
DISPLAY_NAME.....Withdraws To Pain
HELP_TEXT.....No Help Text
ID.....{013527C8-FF6B-49CD-BF45-9338F7413E81}
TERM.....WITHDRAWS TO PAIN

DATA VALUE: WITHOUT ASSISTANCE

ABBREVIATION.....WO ASSIS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Without Assistance
HELP_TEXT.....No Help Text
ID.....{41CE081D-CB0E-49AF-936A-4FD78E61681B}
TERM.....WITHOUT ASSISTANCE

DATA VALUE: WOUND

ABBREVIATION.....WOUND
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Wound
HELP_TEXT.....No Help Text
ID.....{1EC40149-051E-4750-8651-2B2D0DD2DEAC}
TERM.....WOUND

DATA VALUE: WOUND DISCHARGE CHARACTERISTICS BLOODY

ABBREVIATION.....WOU DI BLO
ACTIVE.....Yes
DESCRIPTION.....Wound Discharge Characteristics Bloody
DISPLAY_NAME.....Bloody
HELP_TEXT.....Wound Discharge Characteristics Bloody
ID.....{628DBFBC-9B77-48EA-99FE-0DF000376506}
TERM.....WOUND DISCHARGE CHARACTERISTICS BLOODY

DATA VALUE: WOUND DISCHARGE CHARACTERISTICS CLOUDY

ABBREVIATION.....WOU DI CLO
ACTIVE.....Yes
DESCRIPTION.....Wound Discharge Characteristics Cloudy
DISPLAY_NAME.....Cloudy
HELP_TEXT.....Wound Discharge Characteristics Cloudy
ID.....{159C9993-A2E1-4133-9D5F-315EF9985EC4}
TERM.....WOUND DISCHARGE CHARACTERISTICS CLOUDY

DATA VALUE: WOUND DISCHARGE CHARACTERISTICS NONE

ABBREVIATION.....WOU DI NON
ACTIVE.....Yes
DESCRIPTION.....Wound Discharge Characteristics None
DISPLAY_NAME.....None
HELP_TEXT.....Wound Discharge Characteristics None
ID.....{929F2FA5-923F-4D88-8EDF-E0EBEE9CC7A8}
TERM.....WOUND DISCHARGE CHARACTERISTICS NONE

DATA VALUE: WOUND DISCHARGE CHARACTERISTICS PURULENT

ABBREVIATION.....WOU DI PUR
ACTIVE.....Yes
DESCRIPTION.....Wound Discharge Characteristics
Purulent
DISPLAY_NAME.....Purulent
HELP_TEXT.....Wound Discharge Characteristics
Purulent
ID.....{C1E10410-B6D8-4429-BD94-0B179992EAAB}
TERM.....WOUND DISCHARGE CHARACTERISTICS
PURULENT

DATA VALUE: WOUND DISCHARGE CHARACTERISTICS SEROSANGUINOUS

ABBREVIATION.....WOU DI SSA
ACTIVE.....Yes
DESCRIPTION.....Wound Discharge Characteristics
Serosanguinous
DISPLAY_NAME.....Serosanguinous
HELP_TEXT.....Wound Discharge Characteristics
Serosanguinous
ID.....{B431D8EF-5B53-4DCF-A390-431EB5ED3704}
TERM.....WOUND DISCHARGE CHARACTERISTICS
SEROSANGUINOUS

DATA VALUE: WOUND DISCHARGE CHARACTERISTICS SEROUS

ABBREVIATION.....WOU DI SER
ACTIVE.....Yes
DESCRIPTION.....Wound Discharge Characteristics Serous
DISPLAY_NAME.....Serous
HELP_TEXT.....Wound Discharge Characteristics Serous
ID.....{2124005E-5C36-48AC-AA6B-03E06510D70A}
TERM.....WOUND DISCHARGE CHARACTERISTICS SEROUS

DATA VALUE: WOUND DISCHARGE CHARACTERISTICS VISCOUS

ABBREVIATION.....WOU DI VIS
ACTIVE.....Yes
DESCRIPTION.....Wound Discharge Characteristics Viscous
DISPLAY_NAME.....Viscous
HELP_TEXT.....Wound Discharge Characteristics Viscous
ID.....{382FA6D9-1CB2-4ECC-A458-8598ABA708F8}
TERM.....WOUND DISCHARGE CHARACTERISTICS VISCOUS

DATA VALUE: WOUND DRAINAGE AMOUNT HEAVY

ABBREVIATION.....W DR A HEA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Heavy
HELP_TEXT.....No Help-Text
ID.....{BF959172-6672-4FEB-8FB7-0568405F3C83}
TERM.....WOUND DRAINAGE AMOUNT HEAVY

DATA VALUE: WOUND DRAINAGE AMOUNT MODERATE

ABBREVIATION.....W DR A MOD
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Moderate
HELP_TEXT.....No Help-Text
ID.....{E5DBFD9C-A2AC-4009-98CB-7246F6AA702D}
TERM.....WOUND DRAINAGE AMOUNT MODERATE

DATA VALUE: WOUND DRAINAGE AMOUNT NONE

ABBREVIATION.....WOU DR A N
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....None
HELP_TEXT.....No Help-Text
ID.....{92910398-311F-D4D7-8BD5-AB66C3DC5E72}
TERM.....WOUND DRAINAGE AMOUNT NONE

DATA VALUE: WOUND DRAINAGE AMOUNT SCANT

ABBREVIATION.....W DR SCANT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Scant
HELP_TEXT.....No Help-Text
ID.....{9C3CD0BD-F960-4098-B14C-17B20703B242}
TERM.....WOUND DRAINAGE AMOUNT SCANT

DATA VALUE: WOUND DRAINAGE AMOUNT SMALL

ABBREVIATION.....W DR A SMA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Small
HELP_TEXT.....No Help-Text
ID.....{3D8EDF1A-EC4B-41EF-A194-73FF60074A57}
TERM.....WOUND DRAINAGE AMOUNT SMALL

DATA VALUE: WOUND DRESSING TYPE ANTIMICROBIAL

ABBREVIATION.....WOU DR ANT
ACTIVE.....No
DESCRIPTION.....No Description
DISPLAY_NAME.....Antimicrobial
HELP_TEXT.....No Help Text
ID.....{9DAE650A-3D5C-498B-BC74-F3CE77CD2BA9}
TERM.....WOUND DRESSING TYPE ANTIMICROBIAL

DATA VALUE: WOUND DRESSING TYPE CALCIUM ALGINATE

ABBREVIATION.....WOU DR CAL
ACTIVE.....No
DESCRIPTION.....No Description
DISPLAY_NAME.....Calcium Alginate
HELP_TEXT.....No Help Text
ID.....{AA084ED8-07D0-4D8A-84D1-51F025246F40}
TERM.....WOUND DRESSING TYPE CALCIUM ALGINATE

DATA VALUE: WOUND DRESSING TYPE GAUZE

ABBREVIATION.....WOU DR GAU
ACTIVE.....No
DESCRIPTION.....No Description
DISPLAY_NAME.....Gauze
HELP_TEXT.....No Help Text
ID.....{A2EE67E2-B006-413F-8F8A-E02C21A39FC1}
TERM.....WOUND DRESSING TYPE GAUZE

DATA VALUE: WOUND DRESSING TYPE GAUZE PACKING

ABBREVIATION.....WOU DR GPA
ACTIVE.....No
DESCRIPTION.....No Description
DISPLAY_NAME.....Gauze Packing
HELP_TEXT.....No Help Text
ID.....{ED06B435-B205-4BD2-B3D1-795C40EE016F}
TERM.....WOUND DRESSING TYPE GAUZE PACKING

DATA VALUE: WOUND DRESSING TYPE GAUZE WRAP

ABBREVIATION.....WOU DR GWR
ACTIVE.....No
DESCRIPTION.....No Description
DISPLAY_NAME.....Gauze Wrap
HELP_TEXT.....No Help Text
ID.....{4257C446-7627-4632-8507-B5880AD934D5}
TERM.....WOUND DRESSING TYPE GAUZE WRAP

DATA VALUE: WOUND DRESSING TYPE HYDROCOLLOIDAL

ABBREVIATION.....WOU DR HYD
ACTIVE.....No
DESCRIPTION.....No Description

DISPLAY_NAME.....Hydrocolloidal
HELP_TEXT.....No Help Text
ID.....{8C20226C-FC6A-49A0-9A72-06876DB1F203}
TERM.....WOUND DRESSING TYPE HYDROCOLLOIDAL

DATA VALUE: WOUND DRESSING TYPE HYDROGEL

ABBREVIATION.....WOU DR HYD
ACTIVE.....No
DESCRIPTION.....No Description
DISPLAY_NAME.....Hydrogel
HELP_TEXT.....No Help Text
ID.....{60DDE603-9FF1-44FA-BC4D-45AA8753E254}
TERM.....WOUND DRESSING TYPE HYDROGEL

DATA VALUE: WOUND DRESSING TYPE NONE

ABBREVIATION.....WOU DR NON
ACTIVE.....No
DESCRIPTION.....No Description
DISPLAY_NAME.....None
HELP_TEXT.....No Help Text
ID.....{8A58A493-A02C-4570-9786-9F5AF8D02360}
TERM.....WOUND DRESSING TYPE NONE

DATA VALUE: WOUND DRESSING TYPE PRESSURE

ABBREVIATION.....WOU DR PRE
ACTIVE.....No
DESCRIPTION.....No Description
DISPLAY_NAME.....Pressure
HELP_TEXT.....No Help Text
ID.....{520952E4-FB31-4C1C-9FE2-A10E0F1CEE00}
TERM.....WOUND DRESSING TYPE PRESSURE

DATA VALUE: WOUND DRESSING TYPE TRANSPARENT

ABBREVIATION.....WOU DR TRA
ACTIVE.....No
DESCRIPTION.....No Description
DISPLAY_NAME.....Transparent
HELP_TEXT.....No Help Text
ID.....{EA57D935-33D8-4F30-8511-706870A27122}
TERM.....WOUND DRESSING TYPE TRANSPARENT

Page: 906

DATA VALUE: WOUND EDGE INDISTINCT/ATTACHED

ABBREVIATION.....WOU ED IND
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Indistinct, Attached
HELP_TEXT.....No Help-Text
ID.....{C899235C-4A07-4420-92BC-8F818B8E0A7F}

TERM.....WOUND EDGE INDISTINCT/ATTACHED

DATA VALUE: WOUND EDGE NON-ATTACHED

ABBREVIATION.....WOU ED NON
ACTIVE.....Yes
DESCRIPTION.....Wound Edge Non Attached
DISPLAY_NAME.....Non-Attached
HELP_TEXT.....Wound Edge Non Attached
ID.....{1D5C1317-B4FF-421C-96A0-C20C1089E672}
TERM.....WOUND EDGE NON-ATTACHED

DATA VALUE: WOUND EDGE ROLLED

ABBREVIATION.....WOU ED ROL
ACTIVE.....Yes
DESCRIPTION.....Wound Edge Rolled
DISPLAY_NAME.....Rolled
HELP_TEXT.....Wound Edge Rolled
ID.....{24D9CFEC-1ED6-4069-9715-1F9FCB45BBD3}
TERM.....WOUND EDGE ROLLED

DATA VALUE: WOUND EDGE TUNNELING

ABBREVIATION.....WOU ED TUN
ACTIVE.....Yes
DESCRIPTION.....Wound Edge Tunneling
DISPLAY_NAME.....Tunneling
HELP_TEXT.....Wound Edge Tunneling
ID.....{97425206-5261-46C1-905F-E03E2F071381}
TERM.....WOUND EDGE TUNNELING

DATA VALUE: WOUND EDGE UNDERMINING

ABBREVIATION.....WOU ED UND
ACTIVE.....Yes
DESCRIPTION.....Wound Edge Undermining
DISPLAY_NAME.....Undermining
HELP_TEXT.....Wound Edge Undermining
ID.....{AAA3BB3A-3A3B-4F39-9AC7-3BE1A82BDC02}
TERM.....WOUND EDGE UNDERMINING

Page: 907

DATA VALUE: WOUND ODOR FOUL

ABBREVIATION.....WOU OD FOU
ACTIVE.....Yes
DESCRIPTION.....Wound Odor Foul
DISPLAY_NAME.....Foul
HELP_TEXT.....Wound Odor Foul
ID.....{99A74112-EA71-4896-A4DD-911379E5B453}
TERM.....WOUND ODOR FOUL

DATA VALUE: WOUND ODOR NONE

ABBREVIATION.....WOU OD NON
ACTIVE.....Yes
DESCRIPTION.....Wound Odor None
DISPLAY_NAME.....None
HELP_TEXT.....Wound Odor None
ID.....{D2E388F7-B531-4D60-A2BE-830F11B9AA75}
TERM.....WOUND ODOR NONE

DATA VALUE: WOUND VAC DEVICE

ABBREVIATION.....WOUND VAC
ACTIVE.....Yes
DESCRIPTION.....Wound Vac Device
DISPLAY_NAME.....Wound Vac Device
HELP_TEXT.....No Help Text
ID.....{A03F4978-969D-DCB8-DCE2-26EB33F097B3}
TERM.....WOUND VAC DEVICE

DATA VALUE: X-RAY

ABBREVIATION.....X-RAY
ACTIVE.....Yes
DESCRIPTION.....X-Ray
DISPLAY_NAME.....X-Ray
HELP_TEXT.....X-Ray
ID.....{C82CB7F8-2456-4A08-9FD2-4E498EB1DB20}
TERM.....X-RAY

DATA VALUE: XENADERM

ABBREVIATION.....XENADERM
ACTIVE.....Yes
DESCRIPTION.....Xenaderm
DISPLAY_NAME.....Xenaderm
HELP_TEXT.....No Help Text
ID.....{F181D70E-A231-4D2B-A1E6-DDF4DBED31CD}
TERM.....XENADERM

Page: 908

DATA VALUE: XEROFORM

ABBREVIATION.....XEROFORM
ACTIVE.....Yes
DESCRIPTION.....Xeroform
DISPLAY_NAME.....Xeroform
HELP_TEXT.....No Help Text
ID.....{7D7C3CDC-253E-4262-B8D4-3B70DE9AEBE4}
TERM.....XEROFORM

DATA VALUE: YELLOW

ABBREVIATION.....YELLOW
ACTIVE.....Yes

DESCRIPTION.....The color yellow
DISPLAY_NAME.....Yellow
HELP_TEXT.....No Help Text
ID.....{90DC31AD-6211-9F25-4237-79F7122A2763}
TERM.....YELLOW

DATA VALUE: YELLOW/FIBRINOUS

ABBREVIATION.....YELLOW/FIB
ACTIVE.....Yes
DESCRIPTION.....Yellow / Fibrinous
DISPLAY_NAME.....Yellow / Fibrinous
HELP_TEXT.....No Help Text
ID.....{47CF52FE-0374-4690-ADA1-643E0954E035}
TERM.....YELLOW/FIBRINOUS

DATA VALUE: YELLOW/GREEN

ABBREVIATION.....YELLOW/GRN
ACTIVE.....Yes
DESCRIPTION.....Yellow/Green in appearance
DISPLAY_NAME.....Yellow / Green
HELP_TEXT.....No Help Text
ID.....{1D0E7E81-3331-4D07-9B94-3D4BF72A0F0D}
TERM.....YELLOW/GREEN

DATA VALUE: YES

ABBREVIATION.....y
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Yes
HELP_TEXT.....No Help Text
ID.....{483DE806-2EDB-34B0-FE11-C29EEEE649CD}
TERM.....YES

Page: 909

DATA VALUE: pH

ABBREVIATION.....pH
ACTIVE.....Yes
DESCRIPTION.....pH
DISPLAY_NAME.....pH
HELP_TEXT.....pH
ID.....{104FE39C-B0A2-40B5-A523-B11E5F2BA272}
TERM.....pH

=====

TYPE: LOCATION

Page: 910

=====

LOCATION: ABDOMEN

ABBREVIATION.....ABDOMEN
ACTIVE.....Yes
DESCRIPTION.....ABDOMEN
DISPLAY_NAME.....Abdomen
HELP_TEXT.....No Help Text
ID.....{109FE041-2521-6EF2-7F95-1FD1587796AB}
TERM.....ABDOMEN

LOCATION: ABDOMEN LLQ

ABBREVIATION.....ABD LLQ
ACTIVE.....Yes
DESCRIPTION.....Abdomen Left Lower Quarter
DISPLAY_NAME.....Abdomen Left Lower Quarter
HELP_TEXT.....No Help Text
ID.....{39831F59-D7B0-4156-AC57-4C622A93CA58}
TERM.....ABDOMEN LLQ

LOCATION: ABDOMEN LUQ

ABBREVIATION.....ABD LUQ
ACTIVE.....Yes
DESCRIPTION.....Abdomen Left Upper Quarter
DISPLAY_NAME.....Abdomen Left Upper Quarter
HELP_TEXT.....No Help Text
ID.....{3B31E1EC-CE74-4C76-BB61-10C678431FCA}
TERM.....ABDOMEN LUQ

LOCATION: ABDOMEN RLQ

ABBREVIATION.....ABD RLQ
ACTIVE.....Yes
DESCRIPTION.....Abdomen Right Lower Quarter
DISPLAY_NAME.....Abdomen Right Lower Quarter
HELP_TEXT.....No Help Text
ID.....{2C48394F-6877-4738-856C-E933644D8067}
TERM.....ABDOMEN RLQ

LOCATION: ABDOMEN RUQ

ABBREVIATION.....ABD RUQ
ACTIVE.....Yes
DESCRIPTION.....Abdomen Right Upper Quarter
DISPLAY_NAME.....Abdomen Right Upper Quarter
HELP_TEXT.....No Help Text
ID.....{F29C9FA1-717A-49B7-8000-95070E9AACA0}
TERM.....ABDOMEN RUQ

Page: 911

LOCATION: ABDOMINAL

ABBREVIATION.....Abd
ACTIVE.....Yes
DESCRIPTION.....No Description

DISPLAY_NAME.....Abdominal
HELP_TEXT.....No Help Text
ID.....{AE1DA31F-2C28-98C2-501A-F3D5ADB2501F}
TERM.....ABDOMINAL

LOCATION: ABDOMINAL AORTA

ABBREVIATION.....ABD AORTA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Abdominal Aorta
HELP_TEXT.....No Help Text
ID.....{C4BD42DD-1D6B-4B2C-B6D5-19E1185E06E3}
TERM.....ABDOMINAL AORTA

LOCATION: ANKLE

ABBREVIATION.....Ank
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Ankle
HELP_TEXT.....No Help Text
ID.....{F2A2CF53-BFA5-9279-DB22-EA2627988E5B}
TERM.....ANKLE

LOCATION: APICAL

ABBREVIATION.....Ap
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Apical
HELP_TEXT.....No Help Text
ID.....{C5CF5C51-63CC-33FB-9BAE-FB1DD9900911}
TERM.....APICAL

LOCATION: AV FISTULA

ABBREVIATION.....AVF
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....AV Fistula
HELP_TEXT.....No Help Text
ID.....{F977EEF6-EA6A-D56D-D1F7-29B91510A02F}
TERM.....AV FISTULA

Page: 912

LOCATION: AXILLARY

ABBREVIATION.....Ax
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Axillary
HELP_TEXT.....No Help Text
ID.....{135AA5F4-A1C7-978D-E9AB-6873841AEB57}

TERM.....AXILLARY

LOCATION: BACK

ABBREVIATION.....BACK
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Back
HELP_TEXT.....No Help Text
ID.....{22C26239-C692-4E93-B522-E5E8392BA3B1}
TERM.....BACK

LOCATION: BILATERAL ANKLES

ABBREVIATION.....BILA ANKLE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bilateral Ankles
HELP_TEXT.....No Help-Text
ID.....{E56A03A5-5CE0-47BF-A14D-F6DF66F4E769}
TERM.....BILATERAL ANKLES

LOCATION: BILATERAL FOOT

ABBREVIATION.....BILAT FOOT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bilateral
HELP_TEXT.....No Help Text
ID.....{C758847D-468F-479A-B12D-F4E46F0E5E3A}
TERM.....BILATERAL FOOT

LOCATION: BILATERAL HAND

ABBREVIATION.....BILAT HAND
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bilateral
HELP_TEXT.....No Help Text
ID.....{C782FE56-6117-493F-847F-83173407F0B5}
TERM.....BILATERAL HAND

Page: 913

LOCATION: BILATERAL LEG

ABBREVIATION.....BILAT LEG
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bilateral
HELP_TEXT.....No Help Text
ID.....{0713761D-E0FC-DA03-02A1-F20EFEB99F45}
TERM.....BILATERAL LEG

LOCATION: BILATERAL NECK

ABBREVIATION.....BILAT NECK
ACTIVE.....Yes
DESCRIPTION.....Bilateral Neck
DISPLAY_NAME.....Bilateral Neck
HELP_TEXT.....No Help Text
ID.....{744BF4A2-C55E-4D55-8711-58700DD1C446}
TERM.....BILATERAL NECK

LOCATION: BILATERAL PERIPHERALS

ABBREVIATION.....BILAT PERI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bilateral Peripherals
HELP_TEXT.....No Help Text
ID.....{7F134529-36C5-E9D4-ABD5-64B404F9E706}
TERM.....BILATERAL PERIPHERALS

LOCATION: BILATERAL WRIST

ABBREVIATION.....BILAT WRIS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bilateral
HELP_TEXT.....No Help Text
ID.....{89B4D30F-9E46-F8D6-14F6-1CC0843DDC21}
TERM.....BILATERAL WRIST

LOCATION: BLADDER

ABBREVIATION.....BLADDER
ACTIVE.....Yes
DESCRIPTION.....Bladder
DISPLAY_NAME.....Bladder
HELP_TEXT.....Bladder
ID.....{18E083F2-C4BB-407E-B0A3-20C4942FC97F}
TERM.....BLADDER

Page: 914

LOCATION: BOTH ARMS

ABBREVIATION.....BOTH ARMS
ACTIVE.....Yes
DESCRIPTION.....Both Arms
DISPLAY_NAME.....Both Arms
HELP_TEXT.....No Help Text
ID.....{A493AB32-9DE9-40AC-8C46-7EFF395B63B5}
TERM.....BOTH ARMS

LOCATION: BOTH HANDS

ABBREVIATION.....BOTH HANDS
ACTIVE.....Yes

DESCRIPTION.....Both Hands
DISPLAY_NAME.....Both Hands
HELP_TEXT.....No Help Text
ID.....{63AA2685-CB1D-4680-B2F4-2199A22D1EC6}
TERM.....BOTH HANDS

LOCATION: BOTH LEGS

ABBREVIATION.....BOTH LEGS
ACTIVE.....Yes
DESCRIPTION.....Both Legs
DISPLAY_NAME.....Both Legs
HELP_TEXT.....No Help Text
ID.....{5A15A231-4695-4CF3-AD0F-2E7950208967}
TERM.....BOTH LEGS

LOCATION: BRACHIAL

ABBREVIATION.....Bra
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Brachial
HELP_TEXT.....No Help Text
ID.....{67426127-09E1-E566-F31E-15714A5B2B7C}
TERM.....BRACHIAL

LOCATION: CALF

ABBREVIATION.....Clf
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Calf
HELP_TEXT.....No Help Text
ID.....{7923E496-90D6-B6D5-726A-A5B8DF3436CF}
TERM.....CALF

Page: 915

LOCATION: CAROTID

ABBREVIATION.....Car
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Carotid
HELP_TEXT.....No Help Text
ID.....{9536F203-1797-9336-6F69-604249E1B3A0}
TERM.....CAROTID

LOCATION: CHEST

ABBREVIATION.....CHEST
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Chest
HELP_TEXT.....No Help Text

ID.....{8C3C30FE-6E06-C35E-8226-CBF0C94014CC}
TERM.....CHEST

LOCATION: COCCYX

ABBREVIATION.....COCCYX
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Coccyx
HELP_TEXT.....No Help Text
ID.....{49D65E9D-3D5B-40E3-B7CF-1003F80B852C}
TERM.....COCCYX

LOCATION: CORE

ABBREVIATION.....C
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Core
HELP_TEXT.....No Help Text
ID.....{5076C291-4EE7-7F9B-E7A5-3BB0A78639B7}
TERM.....CORE

LOCATION: DONOR SITE

ABBREVIATION.....DONOR SITE
ACTIVE.....Yes
DESCRIPTION.....Donor Site
DISPLAY_NAME.....Donor Site
HELP_TEXT.....No Help Text
ID.....{F0A65974-7650-4A4B-AC88-0945E37A0AD2}
TERM.....DONOR SITE

Page: 916

LOCATION: DORSALIS PEDIS

ABBREVIATION.....DP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Dorsalis Pedis
HELP_TEXT.....No Help Text
ID.....{9657D743-1632-2044-B4C5-9EF291A2F6D0}
TERM.....DORSALIS PEDIS

LOCATION: FACE

ABBREVIATION.....FACE
ACTIVE.....Yes
DESCRIPTION.....PATIENT'S FACE
DISPLAY_NAME.....Face
HELP_TEXT.....No Help Text
ID.....{C638EA14-7903-6A94-496A-C6560D51BEE3}
TERM.....FACE

LOCATION: FEET

ABBREVIATION.....FEET
ACTIVE.....Yes
DESCRIPTION.....BOTH OF PATIENT'S FEET
DISPLAY_NAME.....Feet
HELP_TEXT.....No Help Text
ID.....{944CAA0E-5165-A6AF-BE3D-85455440BF40}
TERM.....FEET

LOCATION: FEMORAL

ABBREVIATION.....Fem
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Femoral
HELP_TEXT.....No Help Text
ID.....{E3442C9E-8677-7296-D314-7365319C3716}
TERM.....FEMORAL

LOCATION: FOREHEAD

ABBREVIATION.....FOREHEAD
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Forehead
HELP_TEXT.....No Help Text
ID.....{5F892722-956E-474C-BA08-4074142D49F6}
TERM.....FOREHEAD

Page: 917

LOCATION: FRONTAL

ABBREVIATION.....FRONTAL
ACTIVE.....Yes
DESCRIPTION.....Frontal
DISPLAY_NAME.....Frontal
HELP_TEXT.....No Help Text
ID.....{CF8F7565-13BF-4180-A222-44FFFA06992A}
TERM.....FRONTAL

LOCATION: GASTROINTESTINAL

ABBREVIATION.....GI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....GI
HELP_TEXT.....No Help Text
ID.....{C27C66DB-3A6A-487E-9981-FA7FD323C0C3}
TERM.....GASTROINTESTINAL

LOCATION: GENERAL UMBILICUS

ABBREVIATION.....GEN UMBILI

ACTIVE.....Yes
DESCRIPTION.....General Umbilicus
DISPLAY_NAME.....General Umbilicus
HELP_TEXT.....General Umbilicus
ID.....{34CEEF50-14E8-44B2-9FEA-3BAB36D08DFA}
TERM.....GENERAL UMBILICUS

LOCATION: GENERALIZED

ABBREVIATION.....GENERIZED
ACTIVE.....Yes
DESCRIPTION.....TERM DESCRIBES GENERAL OBSERVATIONS
 OF PATIENT'S ENTIRE BODY
DISPLAY_NAME.....Generalized
HELP_TEXT.....No Help Text
ID.....{17073C6C-8691-CADC-2238-CBAAB0A5C648}
TERM.....GENERALIZED

LOCATION: GENITALS

ABBREVIATION.....GENITALS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Genitals
HELP_TEXT.....No Help Text
ID.....{83437106-3336-4925-92C8-F959561687FC}
TERM.....GENITALS

Page: 918

LOCATION: HEAD

ABBREVIATION.....Hd
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Head
HELP_TEXT.....No Help Text
ID.....{C3209F44-90C6-FF1B-3169-265B294F276A}
TERM.....HEAD

LOCATION: HEAD OCCIPITAL

ABBREVIATION.....HEAD OCCIP
ACTIVE.....Yes
DESCRIPTION.....Head Occipital
DISPLAY_NAME.....Head Occipital
HELP_TEXT.....No Help Text
ID.....{2066CA47-5A16-447F-9919-24ED435B627F}
TERM.....HEAD OCCIPITAL

LOCATION: HEART

ABBREVIATION.....HEART
ACTIVE.....Yes
DESCRIPTION.....No Description

DISPLAY_NAME.....Heart
HELP_TEXT.....No Help Text
ID.....{F7D48A36-B512-42DD-B92F-07ACF2B95E08}
TERM.....HEART

LOCATION: HIP

ABBREVIATION.....Hip
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Hip
HELP_TEXT.....No Help Text
ID.....{D2C8BD11-7FB2-7217-A663-EBFC8C9428E7}
TERM.....HIP

LOCATION: INCISIONAL SITE

ABBREVIATION.....INCISON SI
ACTIVE.....Yes
DESCRIPTION.....Incisional Site
DISPLAY_NAME.....Incisional Site
HELP_TEXT.....No Help Text
ID.....{EB36A66D-751D-4031-9F68-0BB03E267321}
TERM.....INCISIONAL SITE

Page: 919

LOCATION: KNEE

ABBREVIATION.....Kn
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Knee
HELP_TEXT.....No Help Text
ID.....{B577308A-4C37-47BD-AD84-38E6342B6DF7}
TERM.....KNEE

LOCATION: KNEES

ABBREVIATION.....KNEES
ACTIVE.....Yes
DESCRIPTION.....BOTH OF PATIENT'S KNEES
DISPLAY_NAME.....Knees
HELP_TEXT.....No Help Text
ID.....{247D6621-3A85-F1F7-59BF-1FB0361B3C21}
TERM.....KNEES

LOCATION: L ABDOMEN

ABBREVIATION.....L ABDOMEN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Abdomen
HELP_TEXT.....No Help Text
ID.....{46B37481-4105-492D-96F0-E2D3895379C1}

TERM.....L ABDOMEN

LOCATION: L ANKLE

ABBREVIATION.....L ANKLE
ACTIVE.....Yes
DESCRIPTION.....Left Ankle
DISPLAY_NAME.....L Ankle
HELP_TEXT.....No Help Text
ID.....{D5A692EF-9D7C-473C-AF96-7EF1C01FA2AF}
TERM.....L ANKLE

LOCATION: L ANTECUBITAL

ABBREVIATION.....L ANTECUBI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Antecubital
HELP_TEXT.....No Help Text
ID.....{C21838CC-BAB5-1D7F-E2A4-A7C663FB28C3}
TERM.....L ANTECUBITAL

Page: 920

LOCATION: L ARM

ABBREVIATION.....LA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Arm
HELP_TEXT.....No Help Text
ID.....{567F552F-54B5-A8BB-596F-291001746E31}
TERM.....L ARM

LOCATION: L AXILLARY

ABBREVIATION.....L AXILLARY
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Axillary
HELP_TEXT.....No Help Text
ID.....{21AF62CA-BEED-2A96-2CD6-F2350A01995E}
TERM.....L AXILLARY

LOCATION: L BRACHIAL

ABBREVIATION.....LBRA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Brachial
HELP_TEXT.....No Help Text
ID.....{75F6CE33-1AB3-A1DA-70B8-04825290B577}
TERM.....L BRACHIAL

LOCATION: L CALF

ABBREVIATION.....L CALF
ACTIVE.....Yes
DESCRIPTION.....Left Calf
DISPLAY_NAME.....L Calf
HELP_TEXT.....No Help Text
ID.....{3DAA7817-936B-444F-8A2D-D3AB09144C10}
TERM.....L CALF

LOCATION: L CAROTID

ABBREVIATION.....L CAROTID
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Carotid
HELP_TEXT.....No Help Text
ID.....{FCACD1EB-C12B-2515-4A68-4E9E28A72DBC}
TERM.....L CAROTID

Page: 921

LOCATION: L CHEST

ABBREVIATION.....L CHEST
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Left Chest
HELP_TEXT.....No Help Text
ID.....{E332E497-FE0C-47CB-A0A7-217FBF52DC97}
TERM.....L CHEST

LOCATION: L DORSALIS PEDIS

ABBREVIATION.....L DOR PED
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Dorsalis Pedis
HELP_TEXT.....No Help Text
ID.....{389397B5-2863-B8F2-A779-C38433894A19}
TERM.....L DORSALIS PEDIS

LOCATION: L EAR

ABBREVIATION.....L EAR
ACTIVE.....Yes
DESCRIPTION.....Left Ear
DISPLAY_NAME.....L Ear
HELP_TEXT.....No Help Text
ID.....{DA9C2DD1-D6FF-4899-A62B-1096FAA36C59}
TERM.....L EAR

LOCATION: L ELBOW

ABBREVIATION.....L ELBOW
ACTIVE.....Yes

DESCRIPTION.....Left Elbow
DISPLAY_NAME.....L Elbow
HELP_TEXT.....No Help Text
ID.....{BF4317E7-2024-47EF-9E5F-3D21D0547885}
TERM.....L ELBOW

LOCATION: L EXTERNAL JUGULAR

ABBREVIATION.....L EXT JUGU
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L External Jugular
HELP_TEXT.....No Help Text
ID.....{6F7B331B-4F83-0B79-1755-F6F673D647D6}
TERM.....L EXTERNAL JUGULAR

Page: 922

LOCATION: L EYE

ABBREVIATION.....L EYE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Eye
HELP_TEXT.....No Help Text
ID.....{D714A0CA-EAD6-EC35-DA54-FF027C23FF45}
TERM.....L EYE

LOCATION: L FACE

ABBREVIATION.....L FACE
ACTIVE.....Yes
DESCRIPTION.....LEFT SIDE OF PATIENT'S FACE
DISPLAY_NAME.....L Face
HELP_TEXT.....No Help Text
ID.....{06F00A70-F475-CD40-665D-AF0569A1E75F}
TERM.....L FACE

LOCATION: L FEMORAL

ABBREVIATION.....LFEM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Femoral
HELP_TEXT.....No Help Text
ID.....{50AF6A08-CE37-6859-F4E6-762660A982F7}
TERM.....L FEMORAL

LOCATION: L FINGER

ABBREVIATION.....L FINGER
ACTIVE.....Yes
DESCRIPTION.....Left Finger
DISPLAY_NAME.....L Finger
HELP_TEXT.....No Help Text

ID.....{BB825645-1630-498C-A31B-6728D422EEFC}
TERM.....L FINGER

LOCATION: L FLANK

ABBREVIATION.....L FLANK
ACTIVE.....Yes
DESCRIPTION.....Left Flank
DISPLAY_NAME.....L Flank
HELP_TEXT.....No Help Text
ID.....{6E453E8D-B456-4BFD-A015-F9C013830F30}
TERM.....L FLANK

Page: 923

LOCATION: L FOOT

ABBREVIATION.....L FOOT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Foot
HELP_TEXT.....No Help Text
ID.....{2EC3FBEC-4EA2-9E7A-17D4-FDC2F9D28046}
TERM.....L FOOT

LOCATION: L GLUTEUS

ABBREVIATION.....L GLUTEUS
ACTIVE.....Yes
DESCRIPTION.....Left Gluteus
DISPLAY_NAME.....L Gluteus
HELP_TEXT.....No Help Text
ID.....{88C4E547-04CD-4D52-9416-CAEB6D898D21}
TERM.....L GLUTEUS

LOCATION: L GROIN

ABBREVIATION.....L GROIN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Groin
HELP_TEXT.....No Help Text
ID.....{5488E4EE-C1A2-4B02-8763-A53C412117A7}
TERM.....L GROIN

LOCATION: L HAND

ABBREVIATION.....L HAND
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Hand
HELP_TEXT.....No Help Text
ID.....{639DF25A-19E3-B504-D323-4B3F6ED39A24}
TERM.....L HAND

LOCATION: L HEAD

ABBREVIATION.....L HEAD
ACTIVE.....Yes
DESCRIPTION.....Left Head
DISPLAY_NAME.....L Head
HELP_TEXT.....No Help Text
ID.....{3737DC2D-DC9E-4FC8-839D-0355987B62D1}
TERM.....L HEAD

Page: 924

LOCATION: L HEEL

ABBREVIATION.....L HEEL
ACTIVE.....Yes
DESCRIPTION.....Left Heel
DISPLAY_NAME.....L Heel
HELP_TEXT.....No Help Text
ID.....{C9154881-17C7-4921-A0EE-BFB962C59CEF}
TERM.....L HEEL

LOCATION: L HIP

ABBREVIATION.....L HIP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Hip
HELP_TEXT.....No Help Text
ID.....{786CD881-27AD-4136-B9CC-8423CB38A431}
TERM.....L HIP

LOCATION: L INGUINAL

ABBREVIATION.....L INGUINAL
ACTIVE.....Yes
DESCRIPTION.....Left Inguinal
DISPLAY_NAME.....L Inguinal
HELP_TEXT.....No Help Text
ID.....{627D172E-8F33-48C9-BE89-D036E07D46BF}
TERM.....L INGUINAL

LOCATION: L INNER KNEE

ABBREVIATION.....L INN KNEE
ACTIVE.....Yes
DESCRIPTION.....Left Inner Knee
DISPLAY_NAME.....L Inner Knee
HELP_TEXT.....No Help Text
ID.....{79853581-CCF0-49E9-9D91-D16DE02A836D}
TERM.....L INNER KNEE

LOCATION: L INTERNAL JUGULAR

ABBREVIATION.....LIJ

ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Internal Jugular
HELP_TEXT.....No Help Text
ID.....{AFBB4FE5-8A11-2047-7900-DADF2FC4430D}
TERM.....L INTERNAL JUGULAR

Page: 925

LOCATION: L ISCHIAL

ABBREVIATION.....L ISCHIAL
ACTIVE.....Yes
DESCRIPTION.....Left Ischial
DISPLAY_NAME.....L Ischial
HELP_TEXT.....No Help Text
ID.....{DCDAD741-9599-4479-9D55-3279F800DBE6}
TERM.....L ISCHIAL

LOCATION: L JAW

ABBREVIATION.....L JAW
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Jaw
HELP_TEXT.....No Help Text
ID.....{65EBD3CF-76D8-4F45-9B1B-3575D988BE8E}
TERM.....L JAW

LOCATION: L KNEE

ABBREVIATION.....L KNEE
ACTIVE.....Yes
DESCRIPTION.....PATIENT'S LEFT KNEE
DISPLAY_NAME.....L Knee
HELP_TEXT.....No Help Text
ID.....{3A9E30B4-09A9-967A-58EF-69A01EC4F859}
TERM.....L KNEE

LOCATION: L LATERAL HEEL

ABBREVIATION.....L LAT HEEL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Lateral Heel
HELP_TEXT.....No Help Text
ID.....{1A10FB20-7BAD-4D2B-9CBD-1A1A01DDCBF1}
TERM.....L LATERAL HEEL

LOCATION: L LATERAL MALLEOLUS

ABBREVIATION.....L LAT MALL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Lateral Malleolus

HELP_TEXT.....No Help Text
ID.....{9A523B7C-64DD-4E47-9543-165788108D4C}
TERM.....L LATERAL MALLEOLUS

Page: 926

LOCATION: L LEG

ABBREVIATION.....LL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Leg
HELP_TEXT.....No Help Text
ID.....{B088C809-E528-B785-21C2-14E45E1FD8AF}
TERM.....L LEG

LOCATION: L LOWER ANTERIOR CHEST

ABBREVIATION.....L L A CHES
ACTIVE.....Yes
DESCRIPTION.....Left Lower Anterior Chest
DISPLAY_NAME.....L Low. Ant. Chest
HELP_TEXT.....No Help Text
ID.....{6073B3C6-AA4D-4BBA-9FF2-32ED40CF92AD}
TERM.....L LOWER ANTERIOR CHEST

LOCATION: L LOWER ARM

ABBREVIATION.....L LOW ARM
ACTIVE.....Yes
DESCRIPTION.....Left Lower Arm
DISPLAY_NAME.....L Lower Arm
HELP_TEXT.....No Help Text
ID.....{63AA95E8-256D-4390-8CB3-F00A1B3F91CD}
TERM.....L LOWER ARM

LOCATION: L LOWER EXTREMITY

ABBREVIATION.....LLEX
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Lower Extremity
HELP_TEXT.....No Help Text
ID.....{FC6F5056-779A-BF76-86D0-A237936C1ECC}
TERM.....L LOWER EXTREMITY

LOCATION: L LOWER LEG

ABBREVIATION.....L LOW LEG
ACTIVE.....Yes
DESCRIPTION.....Left Lower Leg
DISPLAY_NAME.....L Lower Leg
HELP_TEXT.....No Help Text
ID.....{7B8E92FA-E065-4D06-9FB6-4234656CFD61}
TERM.....L LOWER LEG

LOCATION: L LOWER LIMB

ABBREVIATION.....LLLI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Lower Limb
HELP_TEXT.....No Help Text
ID.....{4838CB05-A1DD-F4FB-587C-2D0C17A2484E}
TERM.....L LOWER LIMB

LOCATION: L LOWER POSTERIOR CHEST

ABBREVIATION.....L L P CHES
ACTIVE.....Yes
DESCRIPTION.....Left Lower Posterior Chest
DISPLAY_NAME.....L Low. Pos. Chest
HELP_TEXT.....No Help Text
ID.....{2C19F10B-7546-4A1F-A5DC-3867359ED971}
TERM.....L LOWER POSTERIOR CHEST

LOCATION: L MANDIBULAR

ABBREVIATION.....L MANDIBUL
ACTIVE.....Yes
DESCRIPTION.....Left Mandibular
DISPLAY_NAME.....L Mandibular
HELP_TEXT.....No Help Text
ID.....{CB9DB373-5F0A-446B-909F-F76392CB3A04}
TERM.....L MANDIBULAR

LOCATION: L MEDIAL HEEL

ABBREVIATION.....L MED HEEL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Medial Heel
HELP_TEXT.....No Help Text
ID.....{47AD7BA6-4EBB-4A08-B03C-FE73EB240725}
TERM.....L MEDIAL HEEL

LOCATION: L MEDIAL MALLEOLUS

ABBREVIATION.....L MED MALL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Medial Malleolus
HELP_TEXT.....No Help Text
ID.....{330E7682-407F-43CD-84FC-1CDA41F08C68}
TERM.....L MEDIAL MALLEOLUS

LOCATION: L MIDDLE ANTERIOR CHEST

ABBREVIATION.....L M A CHES
ACTIVE.....Yes
DESCRIPTION.....Left Middle Anterior Chest
DISPLAY_NAME.....L Mid. Ant. Chest
HELP_TEXT.....No Help Text
ID.....{36A83F77-117E-474B-AE34-0B2C02E131AE}
TERM.....L MIDDLE ANTERIOR CHEST

LOCATION: L MIDDLE POSTERIOR CHEST

ABBREVIATION.....L M P CHES
ACTIVE.....Yes
DESCRIPTION.....Left Middle Posterior Chest
DISPLAY_NAME.....L Mid. Pos. Chest
HELP_TEXT.....No Help Text
ID.....{5A656939-F82B-4DEC-9F1D-B5FA644E3DA9}
TERM.....L MIDDLE POSTERIOR CHEST

LOCATION: L NECK

ABBREVIATION.....L NECK
ACTIVE.....Yes
DESCRIPTION.....Left Neck
DISPLAY_NAME.....L Neck
HELP_TEXT.....No Help Text
ID.....{C4EA6279-05E2-4195-A638-3794DCBA7B78}
TERM.....L NECK

LOCATION: L NOSE

ABBREVIATION.....L NOSE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Nose
HELP_TEXT.....No Help Text
ID.....{C7707847-5135-4749-B196-002E3E59BB54}
TERM.....L NOSE

LOCATION: L OCCIPITAL

ABBREVIATION.....L OCCIPTAL
ACTIVE.....Yes
DESCRIPTION.....Left Occipital
DISPLAY_NAME.....L Occipital
HELP_TEXT.....No Help Text
ID.....{A4234F41-CD6F-4CE1-9210-316DB847BDE2}
TERM.....L OCCIPITAL

LOCATION: L ORAL

ABBREVIATION.....L ORAL
ACTIVE.....Yes
DESCRIPTION.....Left Oral
DISPLAY_NAME.....L Oral
HELP_TEXT.....No Help Text
ID.....{057B787A-8750-4563-80DB-25110CCC1146}
TERM.....L ORAL

LOCATION: L ORBIT

ABBREVIATION.....L ORBIT
ACTIVE.....Yes
DESCRIPTION.....Left Orbit
DISPLAY_NAME.....L Orbit
HELP_TEXT.....No Help Text
ID.....{1D798E35-5095-4BA0-8F9F-FA96668A5239}
TERM.....L ORBIT

LOCATION: L OUTER KNEE

ABBREVIATION.....L OUT KNEE
ACTIVE.....Yes
DESCRIPTION.....Left Outer Knee
DISPLAY_NAME.....L Outer Knee
HELP_TEXT.....No Help Text
ID.....{B143C935-2687-49C3-9F65-8D69C1775ABC}
TERM.....L OUTER KNEE

LOCATION: L PARIETAL

ABBREVIATION.....L PARIETAL
ACTIVE.....Yes
DESCRIPTION.....Left Parietal
DISPLAY_NAME.....L Parietal
HELP_TEXT.....No Help Text
ID.....{2687ED24-F35E-456C-99C9-F6A8EC4364E0}
TERM.....L PARIETAL

LOCATION: L PEDAL

ABBREVIATION.....LEFT PEDAL
ACTIVE.....Yes
DESCRIPTION.....L Pedal
DISPLAY_NAME.....L Pedal
HELP_TEXT.....L Pedal
ID.....{D84E9953-8C17-3130-8F21-A30D2A99FDE4}
TERM.....L PEDAL

Page: 930

LOCATION: L POST TIBIAL

ABBREVIATION.....L PO TIB
ACTIVE.....Yes
DESCRIPTION.....No Description

DISPLAY_NAME.....L Post Tibial
HELP_TEXT.....No Help Text
ID.....{BED130EE-31CC-C7E9-6597-78685BCE91AC}
TERM.....L POST TIBIAL

LOCATION: L RADIAL

ABBREVIATION.....L RAD
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Radial
HELP_TEXT.....No Help Text
ID.....{302F78AD-4067-8B6B-BD15-993EC8561AE6}
TERM.....L RADIAL

LOCATION: L RIBS

ABBREVIATION.....L RIBS
ACTIVE.....Yes
DESCRIPTION.....Left Ribs
DISPLAY_NAME.....L Ribs
HELP_TEXT.....No Help Text
ID.....{BDFEE74F-B3CC-4AD9-AE49-78B0BF19A219}
TERM.....L RIBS

LOCATION: L SCAPULA

ABBREVIATION.....L SCAPULA
ACTIVE.....Yes
DESCRIPTION.....Left Scapula
DISPLAY_NAME.....L Scapula
HELP_TEXT.....No Help Text
ID.....{3381EC50-771B-41C8-85C8-6EB0D208D54E}
TERM.....L SCAPULA

LOCATION: L SHOULDER

ABBREVIATION.....L SHOULDER
ACTIVE.....Yes
DESCRIPTION.....Left Shoulder
DISPLAY_NAME.....L Shoulder
HELP_TEXT.....No Help Text
ID.....{91CBD2FA-2743-42EB-946B-6933980D9337}
TERM.....L SHOULDER

LOCATION: L SUBCLAVIAN

ABBREVIATION.....LSCL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Subclavian
HELP_TEXT.....No Help Text
ID.....{71B93594-15F9-AA4A-AEF2-E222EC28C00D}

TERM.....L SUBCLAVIAN

LOCATION: L TEMPORAL

ABBREVIATION.....L TEMPO
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Temporal
HELP_TEXT.....No Help Text
ID.....{375E65CF-A59E-A1E9-EB6D-D53CD818939E}
TERM.....L TEMPORAL

LOCATION: L TESTICLE

ABBREVIATION.....L TESTICLE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Left Testicle
HELP_TEXT.....No Help Text
ID.....{1E376539-E89F-4472-8274-CFD11275E11F}
TERM.....L TESTICLE

LOCATION: L THIGH

ABBREVIATION.....L THIGH
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Thigh
HELP_TEXT.....No Help Text
ID.....{16C2312F-1E0E-45E0-B675-160F40C23F57}
TERM.....L THIGH

LOCATION: L THORACIC ANTERIOR

ABBREVIATION.....L THOR ANT
ACTIVE.....Yes
DESCRIPTION.....Left Thoracic Anterior
DISPLAY_NAME.....L Thoracic Anterior
HELP_TEXT.....No Help Text
ID.....{C85B5FB9-153B-4A8B-9BE6-5524A572ACAC}
TERM.....L THORACIC ANTERIOR

Page: 932

LOCATION: L THORACIC LATERAL

ABBREVIATION.....L THOR LAT
ACTIVE.....Yes
DESCRIPTION.....Left Thoracic Lateral
DISPLAY_NAME.....L Thoracic Lateral
HELP_TEXT.....No Help Text
ID.....{48D26CEB-2E33-4F57-AB3A-F5A46108FFBF}
TERM.....L THORACIC LATERAL

LOCATION: L THORACIC POSTERIOR

ABBREVIATION.....L THOR POS
ACTIVE.....Yes
DESCRIPTION.....Left Thoracic Posterior
DISPLAY_NAME.....L Thoracic Posterior
HELP_TEXT.....No Help Text
ID.....{68E3381F-7605-4786-85A8-523D215A2B97}
TERM.....L THORACIC POSTERIOR

LOCATION: L THORACOTOMY

ABBREVIATION.....L THORACOT
ACTIVE.....Yes
DESCRIPTION.....Left Thoracotomy
DISPLAY_NAME.....L Thoracotomy
HELP_TEXT.....No Help Text
ID.....{313C11A1-559A-4F68-8718-4012416EC189}
TERM.....L THORACOTOMY

LOCATION: L THUMB

ABBREVIATION.....L THUMB
ACTIVE.....Yes
DESCRIPTION.....Left Thumb
DISPLAY_NAME.....L Thumb
HELP_TEXT.....No Help Text
ID.....{F82F4731-2CC5-4705-91BE-05E40C25B419}
TERM.....L THUMB

LOCATION: L TOE

ABBREVIATION.....L TOE
ACTIVE.....Yes
DESCRIPTION.....Left Toe
DISPLAY_NAME.....L Toe
HELP_TEXT.....No Help Text
ID.....{C666DE86-E7AE-4AB3-8AE1-91496219E63D}
TERM.....L TOE

Page: 933

LOCATION: L TRAPEZIUS

ABBREVIATION.....L TRAPEZIU
ACTIVE.....Yes
DESCRIPTION.....Left Trapezius
DISPLAY_NAME.....L Trapezius
HELP_TEXT.....No Help Text
ID.....{B74CF943-C10F-4A1B-A19A-4975E674A97A}
TERM.....L TRAPEZIUS

LOCATION: L TROCHANTER

ABBREVIATION.....L TROCHANT
ACTIVE.....Yes

DESCRIPTION.....Left Trochanter
DISPLAY_NAME.....L Trochanter
HELP_TEXT.....No Help Text
ID.....{C91647A4-255D-40E4-A24D-01274E5D799F}
TERM.....L TROCHANTER

LOCATION: L UPPER ANTERIOR CHEST

ABBREVIATION.....L U A CHES
ACTIVE.....Yes
DESCRIPTION.....Left Upper Anterior Chest
DISPLAY_NAME.....L. Up. Ant. Chest
HELP_TEXT.....No Help Text
ID.....{FE63B58A-DE39-4C03-AE61-E34A5487C767}
TERM.....L UPPER ANTERIOR CHEST

LOCATION: L UPPER ARM

ABBREVIATION.....L UPP ARM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Upper Arm
HELP_TEXT.....No Help Text
ID.....{5F915696-1B08-32E9-6B95-49246117D660}
TERM.....L UPPER ARM

LOCATION: L UPPER EXTREMITY

ABBREVIATION.....LUEX
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Upper Extremity
HELP_TEXT.....No Help Text
ID.....{48806E79-B1DA-129B-EB52-453292A3DDC6}
TERM.....L UPPER EXTREMITY

Page: 934

LOCATION: L UPPER LEG

ABBREVIATION.....L UPPER LE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Left Upper Leg
HELP_TEXT.....No Help Text
ID.....{1DB35895-F74C-4B88-A3EE-3E556CB6559A}
TERM.....L UPPER LEG

LOCATION: L UPPER LIMB

ABBREVIATION.....L Upp Limb
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Upper Limb
HELP_TEXT.....No Help Text

ID.....{E6453716-62D3-DFCB-F84D-F6F2FDFC20B7}
TERM.....L UPPER LIMB

LOCATION: L UPPER POSTERIOR CHEST

ABBREVIATION.....L U P CHES
ACTIVE.....Yes
DESCRIPTION.....Left Upper Posterior Chest
DISPLAY_NAME.....L Up. Pos. Chest
HELP_TEXT.....No Help Text
ID.....{6D5674E9-EF04-4EF1-84A0-3F5ED2FB532A}
TERM.....L UPPER POSTERIOR CHEST

LOCATION: L WRIST

ABBREVIATION.....L WRIST
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L Wrist
HELP_TEXT.....No Help Text
ID.....{7BCAFACE-FF3B-46DA-99F9-6437A90B7126}
TERM.....L WRIST

LOCATION: LEFT

ABBREVIATION.....Lt
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Left
HELP_TEXT.....No Help Text
ID.....{C9C09768-9EEA-8F2B-86D4-B53458472E4E}
TERM.....LEFT

Page: 935

LOCATION: LEFT LOWER LOBE

ABBREVIATION.....LLL
ACTIVE.....Yes
DESCRIPTION.....LEFT LUNG'S LOWER LOBE
DISPLAY_NAME.....L Lower Lobe
HELP_TEXT.....No Help Text
ID.....{14155C13-874D-6E26-01D8-2C03B349F913}
TERM.....LEFT LOWER LOBE

LOCATION: LEFT NECK

ABBREVIATION.....L NECK
ACTIVE.....Yes
DESCRIPTION.....LEFT SIDE OF PATIENT'S NECK
DISPLAY_NAME.....L Neck
HELP_TEXT.....No Help Text
ID.....{40522B5F-F23E-EADB-AE64-8EDE3B9E1155}
TERM.....LEFT NECK

LOCATION: LEFT UPPER LOBE

ABBREVIATION.....LUL
ACTIVE.....Yes
DESCRIPTION.....LEFT LUNG'S UPPER LOBE
DISPLAY_NAME.....L Upper Lobe
HELP_TEXT.....No Help Text
ID.....{21BC3BA2-B22E-6A23-A377-5C2D5C1575C7}
TERM.....LEFT UPPER LOBE

LOCATION: LOWER ABDOMEN

ABBREVIATION.....LO ABDOM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Lower Abdomen
HELP_TEXT.....No Help Text
ID.....{077C9CF5-DCA7-4D94-A7A3-F3C153CEA47D}
TERM.....LOWER ABDOMEN

LOCATION: LOWER ARM

ABBREVIATION.....LrA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Lower Arm
HELP_TEXT.....No Help Text
ID.....{D7627720-D6B0-AC00-8FA3-F7B569ED06FC}
TERM.....LOWER ARM

Page: 936

LOCATION: LOWER BACK

ABBREVIATION.....LOWER BACK
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Lower Back
HELP_TEXT.....No Help Text
ID.....{7E761680-61BC-450F-9905-C45CE760EC45}
TERM.....LOWER BACK

LOCATION: MEDIASTINAL

ABBREVIATION.....MEDIASTINA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Mediastinal
HELP_TEXT.....No Help Text
ID.....{D99997D2-61EF-CAC2-A854-D906E3A97447}
TERM.....MEDIASTINAL

LOCATION: MID BACK

ABBREVIATION.....MID BACK

ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Mid Back
HELP_TEXT.....No Help Text
ID.....{F8C515E0-F437-4596-B030-DA1811038CDF}
TERM.....MID BACK

LOCATION: MID-STERNAL

ABBREVIATION.....MID-STERN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Mid-Sternal
HELP_TEXT.....No Help Text
ID.....{3F5168F1-9082-4FB2-A364-3CDE33B1214B}
TERM.....MID-STERNAL

LOCATION: MIDDLE

ABBREVIATION.....MIDDLE
ACTIVE.....Yes
DESCRIPTION.....MIDDLE
DISPLAY_NAME.....Middle
HELP_TEXT.....No Help Text
ID.....{654E4D50-B74E-D58B-B3E3-AA3F218A285E}
TERM.....MIDDLE

Page: 937

LOCATION: MIDLINE ABDOMINAL

ABBREVIATION.....MID ABDOMI
ACTIVE.....Yes
DESCRIPTION.....Midline Abdominal
DISPLAY_NAME.....Midline Abdominal
HELP_TEXT.....No Help Text
ID.....{37ED2947-9BDA-4207-94F1-6B8C0FE12ECF}
TERM.....MIDLINE ABDOMINAL

LOCATION: MIDLINE NECK

ABBREVIATION.....MID NECK
ACTIVE.....Yes
DESCRIPTION.....Midline Neck
DISPLAY_NAME.....Midline Neck
HELP_TEXT.....No Help Text
ID.....{FF666713-278D-46D6-8D25-E5D57E78EA08}
TERM.....MIDLINE NECK

LOCATION: NECK

ABBREVIATION.....Ne
ACTIVE.....Yes
DESCRIPTION.....PATIENTS NECK
DISPLAY_NAME.....Neck

HELP_TEXT.....No Help Text
ID.....{D4812C3B-AAC2-1B61-679B-3CD8B5CB947C}
TERM.....NECK

LOCATION: NOSE

ABBREVIATION.....NOSE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Nose
HELP_TEXT.....No Help Text
ID.....{85F7003E-6D4E-4A79-A534-8A57FBC2D106}
TERM.....NOSE

LOCATION: OCCIPITAL

ABBREVIATION.....OCCIPITAL
ACTIVE.....Yes
DESCRIPTION.....Occipital
DISPLAY_NAME.....Occipital
HELP_TEXT.....No Help Text
ID.....{5A3E384B-1C87-4A19-A449-DB91AB443FFC}
TERM.....OCCIPITAL

Page: 938

LOCATION: ORAL

ABBREVIATION.....O
ACTIVE.....Yes
DESCRIPTION.....Oral Description
DISPLAY_NAME.....Oral
HELP_TEXT.....Oral Help
ID.....{7A7ACA1B-ACBD-B08D-D40C-54BA8E724B8C}
TERM.....ORAL

LOCATION: PENIS

ABBREVIATION.....PENIS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Penis
HELP_TEXT.....No Help Text
ID.....{2D82CE55-83B6-4160-A754-07C71FD388E8}
TERM.....PENIS

LOCATION: PERINEAL

ABBREVIATION.....PERINEAL
ACTIVE.....Yes
DESCRIPTION.....Perineal
DISPLAY_NAME.....Perineal
HELP_TEXT.....No Help Text
ID.....{1801441C-D5EF-4775-9E58-E4985496F7E0}
TERM.....PERINEAL

LOCATION: PERIORBITAL

ABBREVIATION.....PERIORBITA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Periorbital
HELP_TEXT.....No Help Text
ID.....{DBAC362B-6278-EDB4-7E28-3180321A38B2}
TERM.....PERIORBITAL

LOCATION: PERIPHERAL

ABBREVIATION.....Per
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Peripheral
HELP_TEXT.....No Help Text
ID.....{36B31698-E84A-5792-8EAE-F374DFD1B55A}
TERM.....PERIPHERAL

Page: 939

LOCATION: PERIRECTAL

ABBREVIATION.....PERIRECTAL
ACTIVE.....Yes
DESCRIPTION.....Perirectal
DISPLAY_NAME.....Perirectal
HELP_TEXT.....No Help Text
ID.....{560E1473-7B8F-4A75-8807-29315C908553}
TERM.....PERIRECTAL

LOCATION: POPLITEAL

ABBREVIATION.....Pop
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Popliteal
HELP_TEXT.....No Help Text
ID.....{77030DF2-CE06-D9E2-0684-EE70209CBE73}
TERM.....POPLITEAL

LOCATION: POSTERIOR TIBIAL

ABBREVIATION.....PT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Posterior Tibial
HELP_TEXT.....No Help Text
ID.....{010F2FFE-DEBC-DF9D-9CF2-AD192B16FFD5}
TERM.....POSTERIOR TIBIAL

LOCATION: R ABDOMEN

ABBREVIATION.....R ABDOMEN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Abdomen
HELP_TEXT.....No Help Text
ID.....{5C45EAA8-96F5-4EC7-8AB5-A9543CD6C414}
TERM.....R ABDOMEN

LOCATION: R ANKLE

ABBREVIATION.....R ANKLE
ACTIVE.....Yes
DESCRIPTION.....Right Ankle
DISPLAY_NAME.....R Ankle
HELP_TEXT.....No Help Text
ID.....{4E6DF352-F399-4CE9-908C-48132C3D856B}
TERM.....R ANKLE

Page: 940

LOCATION: R ANTECUBITAL

ABBREVIATION.....R ANTECUBI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Antecubital
HELP_TEXT.....No Help Text
ID.....{7FB00A65-B4B2-06A6-BD89-98AC112D8D4E}
TERM.....R ANTECUBITAL

LOCATION: R ARM

ABBREVIATION.....RA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Arm
HELP_TEXT.....No Help Text
ID.....{1C09017A-F10A-543D-A7EA-F13790D78692}
TERM.....R ARM

LOCATION: R ARTERY

ABBREVIATION.....R ARTERY
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Artery
HELP_TEXT.....No Help Text
ID.....{A2D76F33-A049-B3A0-2CBF-9F83FF52D894}
TERM.....R ARTERY

LOCATION: R AXILLARY

ABBREVIATION.....R AXILLARY
ACTIVE.....Yes
DESCRIPTION.....No Description

DISPLAY_NAME.....R Axillary
HELP_TEXT.....No Help Text
ID.....{3044F28B-FED2-4069-7E3B-2A691C5771F9}
TERM.....R AXILLARY

LOCATION: R BRACHIAL

ABBREVIATION.....RBRA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Brachial
HELP_TEXT.....No Help Text
ID.....{A88E69A8-F6AE-80B0-2CD4-687B121F9DDE}
TERM.....R BRACHIAL

Page: 941

LOCATION: R CALF

ABBREVIATION.....R CALF
ACTIVE.....Yes
DESCRIPTION.....Right Calf
DISPLAY_NAME.....R Calf
HELP_TEXT.....No Help Text
ID.....{0298A452-4901-44BF-A7A8-F42DD3BA8B63}
TERM.....R CALF

LOCATION: R CAROTID

ABBREVIATION.....R CAROTID
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Carotid
HELP_TEXT.....No Help Text
ID.....{8280ACD2-72D6-202A-89AF-86BCA7736427}
TERM.....R CAROTID

LOCATION: R CHEST

ABBREVIATION.....R CHEST
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Right Chest
HELP_TEXT.....No Help Text
ID.....{7E30EC0E-815A-4C2E-A9FC-3C2B2D6E3CD4}
TERM.....R CHEST

LOCATION: R DORSALIS PEDIS

ABBREVIATION.....R DOR PED
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Dorsalis Pedis
HELP_TEXT.....No Help Text
ID.....{47D4C64A-FD93-2999-7B46-F25307EA6511}

TERM.....R DORSALIS PEDIS

LOCATION: R EAR

ABBREVIATION.....R EAR
ACTIVE.....Yes
DESCRIPTION.....Right Ear
DISPLAY_NAME.....R Ear
HELP_TEXT.....No Help Text
ID.....{AE31327D-A03C-488C-8A8E-084E7E6E7FDE}
TERM.....R EAR

Page: 942

LOCATION: R ELBOW

ABBREVIATION.....R ELBOW
ACTIVE.....Yes
DESCRIPTION.....Right Elbow
DISPLAY_NAME.....R Elbow
HELP_TEXT.....No Help Text
ID.....{5BAFF87B-6993-4E76-862F-9F386D9B71C5}
TERM.....R ELBOW

LOCATION: R EXTERNAL JUGULAR

ABBREVIATION.....R EXT JUGU
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R External Jugular
HELP_TEXT.....No Help Text
ID.....{4D78F704-9E3B-907E-A792-4DC4F4CFF170}
TERM.....R EXTERNAL JUGULAR

LOCATION: R EYE

ABBREVIATION.....R EYE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Eye
HELP_TEXT.....No Help Text
ID.....{258F0C1D-67BD-ECF8-3152-EB289DEF4681}
TERM.....R EYE

LOCATION: R FACE

ABBREVIATION.....R FACE
ACTIVE.....Yes
DESCRIPTION.....RIGHT SIDE OF PATIENT'S FACE
DISPLAY_NAME.....R Face
HELP_TEXT.....No Help Text
ID.....{D1B15677-EF6F-0125-3601-68C1FE6386CE}
TERM.....R FACE

LOCATION: R FEMORAL

ABBREVIATION.....R FEMORAL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Femoral
HELP_TEXT.....No Help Text
ID.....{B1F9EF50-B022-2BA9-CDA9-374E6A1557CC}
TERM.....R FEMORAL

Page: 943

LOCATION: R FINGER

ABBREVIATION.....R FINGER
ACTIVE.....Yes
DESCRIPTION.....Right Finger
DISPLAY_NAME.....R Finger
HELP_TEXT.....No Help Text
ID.....{E12CA0EF-3267-4B93-8643-230E2226C3D6}
TERM.....R FINGER

LOCATION: R FLANK

ABBREVIATION.....R FLANK
ACTIVE.....Yes
DESCRIPTION.....Right Flank
DISPLAY_NAME.....R Flank
HELP_TEXT.....No Help Text
ID.....{F5700ED2-646F-4B75-A55B-111C27FB939C}
TERM.....R FLANK

LOCATION: R FOOT

ABBREVIATION.....R FOOT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Foot
HELP_TEXT.....No Help Text
ID.....{7D9F5FF6-232C-2488-4504-A721FF06928F}
TERM.....R FOOT

LOCATION: R GLUTEUS

ABBREVIATION.....R GLUTEUS
ACTIVE.....Yes
DESCRIPTION.....Right Gluteus
DISPLAY_NAME.....R Gluteus
HELP_TEXT.....No Help Text
ID.....{B4D32045-68D1-4457-9106-7937887FA6F6}
TERM.....R GLUTEUS

LOCATION: R GROIN

ABBREVIATION.....R GROIN
ACTIVE.....Yes

DESCRIPTION.....No Description
DISPLAY_NAME.....R Groin
HELP_TEXT.....No Help Text
ID.....{BAE1C889-8CCF-4CEB-9EBD-0243BC6F096A}
TERM.....R GROIN

Page: 944

LOCATION: R HAND

ABBREVIATION.....R HAND
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Hand
HELP_TEXT.....No Help Text
ID.....{A54D3A80-8537-D079-9CCE-7DC56346D132}
TERM.....R HAND

LOCATION: R HEAD

ABBREVIATION.....R HEAD
ACTIVE.....Yes
DESCRIPTION.....Right Head
DISPLAY_NAME.....R Head
HELP_TEXT.....No Help Text
ID.....{42A95AAE-3DA9-4C96-A4E7-A15C9840217A}
TERM.....R HEAD

LOCATION: R HEEL

ABBREVIATION.....R HEEL
ACTIVE.....Yes
DESCRIPTION.....Right Heel
DISPLAY_NAME.....R Heel
HELP_TEXT.....No Help Text
ID.....{0020DD20-2EA4-4362-ADF0-07BFA220C60E}
TERM.....R HEEL

LOCATION: R HIP

ABBREVIATION.....R HIP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Hip
HELP_TEXT.....No Help Text
ID.....{A48CBCAE-5404-4934-9321-ECEF660A5350}
TERM.....R HIP

LOCATION: R INGUINAL

ABBREVIATION.....R INGUINAL
ACTIVE.....Yes
DESCRIPTION.....Right Inguinal
DISPLAY_NAME.....R Inguinal
HELP_TEXT.....No Help Text

ID.....{395233DF-3B89-49E0-946B-EB5C02E51178}
TERM.....R INGUINAL

Page: 945

LOCATION: R INNER KNEE

ABBREVIATION.....R INN KNEE
ACTIVE.....Yes
DESCRIPTION.....Right Inner Knee
DISPLAY_NAME.....R Inner Knee
HELP_TEXT.....No Help Text
ID.....{2AE24BF1-EA1B-490D-9DB7-EEFC9314C560}
TERM.....R INNER KNEE

LOCATION: R INTERNAL JUGULAR

ABBREVIATION.....R INT JUGU
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Internal Jugular
HELP_TEXT.....No Help Text
ID.....{7ABF6916-D379-198D-6454-F0518A74844A}
TERM.....R INTERNAL JUGULAR

LOCATION: R ISCHIAL

ABBREVIATION.....R ISCHIAL
ACTIVE.....Yes
DESCRIPTION.....Right Ischial
DISPLAY_NAME.....R Ischial
HELP_TEXT.....No Help Text
ID.....{97E50F6C-CF75-4E59-AFDE-80FC5CE6D9E3}
TERM.....R ISCHIAL

LOCATION: R JAW

ABBREVIATION.....R JAW
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Jaw
HELP_TEXT.....No Help Text
ID.....{7846AEE2-B6F4-4B5A-909F-56E386FD82A8}
TERM.....R JAW

LOCATION: R KNEE

ABBREVIATION.....R KNEE
ACTIVE.....Yes
DESCRIPTION.....PATIENT'S RIGHT KNEE
DISPLAY_NAME.....R Knee
HELP_TEXT.....No Help Text
ID.....{BCD18D53-07B6-2C41-64FC-3A07A7E979E7}
TERM.....R KNEE

LOCATION: R LATERAL HEEL

ABBREVIATION.....R LAT HEEL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Lateral Heel
HELP_TEXT.....No Help Text
ID.....{A1896CF9-0246-4323-8DED-E4B5FB12CACE}
TERM.....R LATERAL HEEL

LOCATION: R LATERAL MALLEOLUS

ABBREVIATION.....R LAT MALL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Lateral Malleolus
HELP_TEXT.....No Help Text
ID.....{950FC45B-93E8-424F-A5D4-DE646166D7B0}
TERM.....R LATERAL MALLEOLUS

LOCATION: R LEG

ABBREVIATION.....RL
ACTIVE.....Yes
DESCRIPTION.....Right Leg
DISPLAY_NAME.....R Leg
HELP_TEXT.....RIGHT LEG
ID.....{C5E40A7B-6FC4-267D-8CE8-3571D21D81FD}
TERM.....R LEG

LOCATION: R LOWER ANTERIOR CHEST

ABBREVIATION.....R L A CHES
ACTIVE.....Yes
DESCRIPTION.....Right Lower Anterior Chest
DISPLAY_NAME.....R Low. Ant. Chest
HELP_TEXT.....No Help Text
ID.....{96BCF731-B661-4455-B5E7-C4AE2509505F}
TERM.....R LOWER ANTERIOR CHEST

LOCATION: R LOWER ARM

ABBREVIATION.....R LOW ARM
ACTIVE.....Yes
DESCRIPTION.....Right Lower Arm
DISPLAY_NAME.....R Lower Arm
HELP_TEXT.....No Help Text
ID.....{8DE6E277-EFFF-40DD-9E83-D7295AACB6AA}
TERM.....R LOWER ARM

LOCATION: R LOWER EXTREMITY

ABBREVIATION.....R LOW EXTR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Lower Extremity
HELP_TEXT.....No Help Text
ID.....{7FC3CECC-1E6E-A686-FF04-2BEF8A9B9DCA}
TERM.....R LOWER EXTREMITY

LOCATION: R LOWER LEG

ABBREVIATION.....R LOW LEG
ACTIVE.....Yes
DESCRIPTION.....Right Lower Leg
DISPLAY_NAME.....R Lower Leg
HELP_TEXT.....No Help Text
ID.....{C8EA4D71-D97E-46BD-AE9D-C10954C525CD}
TERM.....R LOWER LEG

LOCATION: R LOWER LIMB

ABBREVIATION.....R LOW LIMB
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Lower Limb
HELP_TEXT.....No Help Text
ID.....{A8E26A05-311A-F196-40B5-B5202D9F8058}
TERM.....R LOWER LIMB

LOCATION: R LOWER LOBE

ABBREVIATION.....RLL
ACTIVE.....Yes
DESCRIPTION.....RIGHT LUNG'S LOWER LOBE
DISPLAY_NAME.....R Lower Lobe
HELP_TEXT.....No Help Text
ID.....{418D350A-3602-1918-22C8-53D093AF49B0}
TERM.....R LOWER LOBE

LOCATION: R LOWER POSTERIOR CHEST

ABBREVIATION.....R L P CHES
ACTIVE.....Yes
DESCRIPTION.....Right Lower Posterior Chest
DISPLAY_NAME.....R Low. Pos. Chest
HELP_TEXT.....No Help Text
ID.....{5163180E-FFD2-4BFB-BF8B-4C6B18EFE5D4}
TERM.....R LOWER POSTERIOR CHEST

LOCATION: R MANDIBULAR

ABBREVIATION.....R MANDIBUL

ACTIVE.....Yes
DESCRIPTION.....Right Mandibular
DISPLAY_NAME.....R Mandibular
HELP_TEXT.....No Help Text
ID.....{303A1DE2-CB9A-498F-9CC5-DFA8FBFF10C1}
TERM.....R MANDIBULAR

LOCATION: R MEDIAL HEEL

ABBREVIATION.....R MED HEEL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Medial Heel
HELP_TEXT.....No Help Text
ID.....{C0526527-0AA9-4ECF-B768-99B2A1F5F3BF}
TERM.....R MEDIAL HEEL

LOCATION: R MEDIAL MALLEOLUS

ABBREVIATION.....R MED MALL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Medial Malleolus
HELP_TEXT.....No Help Text
ID.....{C70FB4FA-6516-4567-8E26-46BA439F1DF4}
TERM.....R MEDIAL MALLEOLUS

LOCATION: R MIDDLE ANTERIOR CHEST

ABBREVIATION.....R M A CHES
ACTIVE.....Yes
DESCRIPTION.....Right Middle Anterior Chest
DISPLAY_NAME.....R Mid. Ant. Chest
HELP_TEXT.....No Help Text
ID.....{2F751A9C-66CA-4C22-ACC9-6E81FA63FA2B}
TERM.....R MIDDLE ANTERIOR CHEST

LOCATION: R MIDDLE LOBE

ABBREVIATION.....RML
ACTIVE.....Yes
DESCRIPTION.....RIGHT LUNG'S MIDDLE LOBE
DISPLAY_NAME.....R Middle Lobe
HELP_TEXT.....No Help Text
ID.....{57879159-350C-2055-34D3-C6EFA39CA98F}
TERM.....R MIDDLE LOBE

Page: 949

LOCATION: R MIDDLE POSTERIOR CHEST

ABBREVIATION.....R M P CHES
ACTIVE.....Yes
DESCRIPTION.....Right Middle Posterior Chest
DISPLAY_NAME.....R Mid. Pos. Chest

HELP_TEXT.....No Help Text
ID.....{2BC487D6-6CE8-499E-94A0-A79EAA8E199B}
TERM.....R MIDDLE POSTERIOR CHEST

LOCATION: R NECK

ABBREVIATION.....R NECK
ACTIVE.....Yes
DESCRIPTION.....RIGHT SIDE OF PATIENT'S NECK
DISPLAY_NAME.....R Neck
HELP_TEXT.....No Help Text
ID.....{7E4D84D3-535B-6164-C6C1-837C6F997A0D}
TERM.....R NECK

LOCATION: R NOSE

ABBREVIATION.....R NOSE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Nose
HELP_TEXT.....No Help Text
ID.....{D80BCB19-E71F-43EA-8990-6C121EE7F7A1}
TERM.....R NOSE

LOCATION: R OCCIPITAL

ABBREVIATION.....R OCCIPITA
ACTIVE.....Yes
DESCRIPTION.....Right Occipital
DISPLAY_NAME.....R Occipital
HELP_TEXT.....No Help Text
ID.....{5DB951DB-CB03-482E-B0DA-CC5EE8FB1A67}
TERM.....R OCCIPITAL

LOCATION: R ORAL

ABBREVIATION.....R ORAL
ACTIVE.....Yes
DESCRIPTION.....Right Oral
DISPLAY_NAME.....R Oral
HELP_TEXT.....No Help Text
ID.....{9CEEF934-AA3D-4BDE-BC3B-9E275C51CC9C}
TERM.....R ORAL

Page: 950

LOCATION: R ORBIT

ABBREVIATION.....R ORBIT
ACTIVE.....Yes
DESCRIPTION.....Right Orbit
DISPLAY_NAME.....R Orbit
HELP_TEXT.....No Help Text
ID.....{FFFD974A-1684-49AE-8230-8E12411EA62E}
TERM.....R ORBIT

LOCATION: R OUTER KNEE

ABBREVIATION.....R OUT KNEE
ACTIVE.....Yes
DESCRIPTION.....Right Outer Knee
DISPLAY_NAME.....R Outer Knee
HELP_TEXT.....No Help Text
ID.....{330F1DF1-8C99-4984-9965-7404ABC9562D}
TERM.....R OUTER KNEE

LOCATION: R PARIETAL

ABBREVIATION.....R PARIETAL
ACTIVE.....Yes
DESCRIPTION.....Right Parietal
DISPLAY_NAME.....R Parietal
HELP_TEXT.....No Help Text
ID.....{28C927CD-1342-483B-B366-0B7EA8ACABAD}
TERM.....R PARIETAL

LOCATION: R PEDAL

ABBREVIATION.....R PEDAL
ACTIVE.....Yes
DESCRIPTION.....R Pedal
DISPLAY_NAME.....R Pedal
HELP_TEXT.....R Pedal
ID.....{9532B614-F4E1-316C-1322-905EB58E58C4}
TERM.....R PEDAL

LOCATION: R POST TIBIAL

ABBREVIATION.....R PO TIB
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Post Tibial
HELP_TEXT.....No Help Text
ID.....{88CDA53A-5153-4950-701A-3F4F2007CC6B}
TERM.....R POST TIBIAL

Page: 951

LOCATION: R RADIAL

ABBREVIATION.....R RADIAL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Radial
HELP_TEXT.....No Help Text
ID.....{BA10E6E5-DDF4-345E-68AF-21D32E9D972D}
TERM.....R RADIAL

LOCATION: R RIBS

ABBREVIATION.....R RIBS
ACTIVE.....Yes
DESCRIPTION.....Right Ribs
DISPLAY_NAME.....R Ribs
HELP_TEXT.....No Help Text
ID.....{260FD564-5E58-419B-9F73-62B6037D1572}
TERM.....R RIBS

LOCATION: R SCAPULA

ABBREVIATION.....R SCAPULA
ACTIVE.....Yes
DESCRIPTION.....Right Scapula
DISPLAY_NAME.....R Scapula
HELP_TEXT.....No Help Text
ID.....{EA657261-4D72-40C2-B70D-2A2E27556A01}
TERM.....R SCAPULA

LOCATION: R SHOULDER

ABBREVIATION.....R SHOULDER
ACTIVE.....Yes
DESCRIPTION.....Right Shoulder
DISPLAY_NAME.....R Shoulder
HELP_TEXT.....No Help Text
ID.....{D27AE130-442E-4DD6-B391-8D7E8AD494A2}
TERM.....R SHOULDER

LOCATION: R SUBCLAVIAN

ABBREVIATION.....R SUBCLAVI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Subclavian
HELP_TEXT.....No Help Text
ID.....{53DEAE83-93AE-123B-14E5-4DE38CD4F752}
TERM.....R SUBCLAVIAN

Page: 952

LOCATION: R TEMPORAL

ABBREVIATION.....R TEMPOR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Temporal
HELP_TEXT.....No Help Text
ID.....{1F50461D-5ECF-25B2-A29A-E02803CB7814}
TERM.....R TEMPORAL

LOCATION: R TESTICLE

ABBREVIATION.....R TESTICLE
ACTIVE.....Yes
DESCRIPTION.....No Description

DISPLAY_NAME.....Right Testicle
HELP_TEXT.....No Help Text
ID.....{0947D277-50AD-47F8-A76D-0EC34EB3B168}
TERM.....R TESTICLE

LOCATION: R THIGH

ABBREVIATION.....R THIGH
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Thigh
HELP_TEXT.....No Help Text
ID.....{B7A2D4A0-801E-4578-967A-F1CFB86FB96B}
TERM.....R THIGH

LOCATION: R THORACIC ANTERIOR

ABBREVIATION.....R THOR ANT
ACTIVE.....Yes
DESCRIPTION.....Right Thoracic Anterior
DISPLAY_NAME.....R Thoracic Anterior
HELP_TEXT.....No Help Text
ID.....{BD3DF757-2C65-4D6A-9C80-5DC95BB41CFB}
TERM.....R THORACIC ANTERIOR

LOCATION: R THORACIC LATERAL

ABBREVIATION.....R THOR LAT
ACTIVE.....Yes
DESCRIPTION.....Right Thoracic Lateral
DISPLAY_NAME.....R Thoracic Lateral
HELP_TEXT.....No Help Text
ID.....{55C48294-D3F3-4013-B552-BA3344D9513B}
TERM.....R THORACIC LATERAL

Page: 953

LOCATION: R THORACIC POSTERIOR

ABBREVIATION.....R THOR POS
ACTIVE.....Yes
DESCRIPTION.....Right Thoracic Posterior
DISPLAY_NAME.....R Thoracic Posterior
HELP_TEXT.....No Help Text
ID.....{E91085D2-E665-4EB6-A086-8CF573DECD8F}
TERM.....R THORACIC POSTERIOR

LOCATION: R THORACOTOMY

ABBREVIATION.....R THORACOT
ACTIVE.....Yes
DESCRIPTION.....Right Thoracotomy
DISPLAY_NAME.....R Thoracotomy
HELP_TEXT.....No Help Text
ID.....{A18E4BB3-0718-48D6-BD49-720A70C0EF60}

TERM.....R THORACOTOMY

LOCATION: R THUMB

ABBREVIATION.....R THUMB
ACTIVE.....Yes
DESCRIPTION.....Right Thumb
DISPLAY_NAME.....R Thumb
HELP_TEXT.....No Help Text
ID.....{B414BD68-4155-45B3-985D-B2F8252A20FE}
TERM.....R THUMB

LOCATION: R TOE

ABBREVIATION.....R TOE
ACTIVE.....Yes
DESCRIPTION.....Right Toe
DISPLAY_NAME.....R Toe
HELP_TEXT.....No Help Text
ID.....{B65727AB-3205-4DCE-833F-6DA8554636F9}
TERM.....R TOE

LOCATION: R TRAPEZIUS

ABBREVIATION.....R TRAPEZIU
ACTIVE.....Yes
DESCRIPTION.....Right Trapezius
DISPLAY_NAME.....R Trapezius
HELP_TEXT.....No Help Text
ID.....{69753B67-AFA1-4C93-BA78-E266F9894E46}
TERM.....R TRAPEZIUS

Page: 954

LOCATION: R TROCHANTER

ABBREVIATION.....R TROCHANT
ACTIVE.....Yes
DESCRIPTION.....Right Trochanter
DISPLAY_NAME.....R Trochanter
HELP_TEXT.....No Help Text
ID.....{1A3E21B2-C391-40D2-8059-12290E667D9E}
TERM.....R TROCHANTER

LOCATION: R UPPER ANTERIOR CHEST

ABBREVIATION.....R U A CHES
ACTIVE.....Yes
DESCRIPTION.....Right Upper Anterior Chest
DISPLAY_NAME.....R Up. Ant. Chest
HELP_TEXT.....No Help Text
ID.....{67890049-A4A8-4241-BDF3-573729BD6EB5}
TERM.....R UPPER ANTERIOR CHEST

LOCATION: R UPPER ARM

ABBREVIATION.....R UPP ARM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Upper Arm
HELP_TEXT.....No Help Text
ID.....{3EF20E8E-F68C-E4F2-8F8C-FAD3B44DDAA4}
TERM.....R UPPER ARM

LOCATION: R UPPER EXTREMITY

ABBREVIATION.....R UPP EXTR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Upper Extremity
HELP_TEXT.....No Help Text
ID.....{C65BA269-8F05-6FB9-62CB-466A286CF2C0}
TERM.....R UPPER EXTREMITY

LOCATION: R UPPER LEG

ABBREVIATION.....R UPPER LE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Right Upper Leg
HELP_TEXT.....No Help Text
ID.....{113BD5ED-5A0A-47CA-B100-91B244BFD04A}
TERM.....R UPPER LEG

Page: 955

LOCATION: R UPPER LIMB

ABBREVIATION.....R UPP LIMB
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Upper Limb
HELP_TEXT.....No Help Text
ID.....{5E028BFF-ED9B-014F-74F0-C4137DB074FB}
TERM.....R UPPER LIMB

LOCATION: R UPPER LOBE

ABBREVIATION.....RUL
ACTIVE.....Yes
DESCRIPTION.....RIGHT LUNG'S UPPER LOBE
DISPLAY_NAME.....R Upper Lobe
HELP_TEXT.....Enter the Right Upper Lobe of the lung
ID.....{098C11B4-B556-FBBE-7E21-AEC9D7E4408A}
TERM.....R UPPER LOBE

LOCATION: R UPPER POSTERIOR CHEST

ABBREVIATION.....R U P CHES
ACTIVE.....Yes

DESCRIPTION.....Right Upper Posterior Chest
DISPLAY_NAME.....R Up. Pos. Chest
HELP_TEXT.....No Help Text
ID.....{D8F03FA5-1FD7-46EE-8FA4-17B3B7E62141}
TERM.....R UPPER POSTERIOR CHEST

LOCATION: R WRIST

ABBREVIATION.....R WRIST
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....R Wrist
HELP_TEXT.....No Help Text
ID.....{FDD9FE1D-F8A2-09AD-82D7-AA478693DD01}
TERM.....R WRIST

LOCATION: RADIAL

ABBREVIATION.....RADIAL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Radial
HELP_TEXT.....No Help-Text
ID.....{D52C4103-2AB3-A0C7-16E3-CF95A1B3B789}
TERM.....RADIAL

Page: 956

LOCATION: RECTAL

ABBREVIATION.....RECTAL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Rectal
HELP_TEXT.....No Help-Text
ID.....{0E7069F4-D3E6-67A1-E821-EF1879C3C267}
TERM.....RECTAL

LOCATION: RIGHT

ABBREVIATION.....R
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Right
HELP_TEXT.....No Help-Text
ID.....{38528AAC-C0A8-ABB7-A659-7EB54D81EA9F}
TERM.....RIGHT

LOCATION: SACRAL

ABBREVIATION.....SACRAL
ACTIVE.....Yes
DESCRIPTION.....PATIENT'S SACRAL AREA
DISPLAY_NAME.....Sacral
HELP_TEXT.....No Help Text

ID.....{0CF8A8A7-AB35-B803-4544-0C91645D9D1C}
TERM.....SACRAL

LOCATION: SCROTUM

ABBREVIATION.....SCROTUM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Scrotum
HELP_TEXT.....No Help Text
ID.....{BDC109F5-0D30-4914-A62C-CEFE448B547}
TERM.....SCROTUM

LOCATION: SKIN

ABBREVIATION.....SKIN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Skin
HELP_TEXT.....No Help-Text
ID.....{0803BC59-1582-AEC8-3DA2-A81363038210}
TERM.....SKIN

Page: 957

LOCATION: SPINE CERVICAL

ABBREVIATION.....SPINE CERV
ACTIVE.....Yes
DESCRIPTION.....Spine Cervical
DISPLAY_NAME.....Spine Cervical
HELP_TEXT.....No Help Text
ID.....{9D28E1BF-D014-4DA4-A832-1F3C25C148E9}
TERM.....SPINE CERVICAL

LOCATION: SPINE LUMBAR

ABBREVIATION.....SPINE LUMB
ACTIVE.....Yes
DESCRIPTION.....Spine Lumbar
DISPLAY_NAME.....Spine Lumbar
HELP_TEXT.....No Help Text
ID.....{E7E411AB-5C29-4906-B0A9-E8658AAFFF91}
TERM.....SPINE LUMBAR

LOCATION: SPINE THORACIC

ABBREVIATION.....SPINE THOR
ACTIVE.....Yes
DESCRIPTION.....Spine Thoracic
DISPLAY_NAME.....Spine Thoracic
HELP_TEXT.....No Help Text
ID.....{2CACBB4F-3EF3-4FBF-AC3E-41D72F0178CB}
TERM.....SPINE THORACIC

LOCATION: STERNOTOMY

ABBREVIATION.....STERNOTOMY
ACTIVE.....Yes
DESCRIPTION.....Sternotomy
DISPLAY_NAME.....Sternotomy
HELP_TEXT.....No Help Text
ID.....{DB97F6A0-9295-42B7-BC29-6A534CE259DD}
TERM.....STERNOTOMY

LOCATION: STOMACH

ABBREVIATION.....STOMACH
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Stomach
HELP_TEXT.....No Help Text
ID.....{613C20B9-CC9C-4F5D-B5B1-69B66A663540}
TERM.....STOMACH

Page: 958

LOCATION: SUBSTERNAL

ABBREVIATION.....SUBSTERN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Substernal
HELP_TEXT.....No Help Text
ID.....{3DE7E0F1-BCA4-43EC-82CB-A87C6506EF46}
TERM.....SUBSTERNAL

LOCATION: TEETH

ABBREVIATION.....TEETH
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Teeth
HELP_TEXT.....No Help Text
ID.....{38B36D8E-A9C3-4EE6-8727-EECDA4C6034E}
TERM.....TEETH

LOCATION: TEMPORAL

ABBREVIATION.....TEMPORAL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Temporal
HELP_TEXT.....No Help-Text
ID.....{E312E0D7-0FE4-A41A-EAB0-ECCB085C3576}
TERM.....TEMPORAL

LOCATION: THIGH

ABBREVIATION.....THIGH

ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Thigh
HELP_TEXT.....No Help-Text
ID.....{2DA93C27-12B0-C76C-85D5-B7D6375F3AD3}
TERM.....THIGH

LOCATION: THYROID

ABBREVIATION.....THYROID
ACTIVE.....Yes
DESCRIPTION.....Thyroid
DISPLAY_NAME.....Thyroid
HELP_TEXT.....No Help Text
ID.....{6AC68004-FE58-46D8-9A1E-4EEF72AAD71B}
TERM.....THYROID

Page: 959

LOCATION: TOOTH

ABBREVIATION.....TOOTH
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Tooth
HELP_TEXT.....No Help Text
ID.....{8407F780-7C7F-49BF-B1DD-4FDD578CEAB3}
TERM.....TOOTH

LOCATION: TORSO

ABBREVIATION.....TORSO
ACTIVE.....Yes
DESCRIPTION.....PATIENT'S TORSO
DISPLAY_NAME.....Torso
HELP_TEXT.....No Help Text
ID.....{EA1AD477-8BFF-E575-23C1-B9CFAF2B9986}
TERM.....TORSO

LOCATION: TRACHEOSTOMY

ABBREVIATION.....TRACHEOSTO
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Tracheostomy
HELP_TEXT.....No Help Text
ID.....{4D12329A-56C8-4175-986F-2F0B68EB7E6D}
TERM.....TRACHEOSTOMY

LOCATION: TRANSVERSE ABDOMINAL

ABBREVIATION.....TRANS ABDO
ACTIVE.....Yes
DESCRIPTION.....Transverse Abdominal
DISPLAY_NAME.....Transverse Abdominal

HELP_TEXT.....No Help Text
ID.....{1C10FF98-7B20-45D0-B420-714981832B3D}
TERM.....TRANSVERSE ABDOMINAL

LOCATION: TYMPANIC

ABBREVIATION.....Tym
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Tympanic
HELP_TEXT.....No Help Text
ID.....{0A58C14C-E447-E034-AA93-D41E6E0D2F69}
TERM.....TYMPANIC

Page: 960

LOCATION: ULNAR

ABBREVIATION.....Uln
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Ulnar
HELP_TEXT.....No Help Text
ID.....{81FC2EA0-472D-06D7-F5A2-E71B0B7DE091}
TERM.....ULNAR

LOCATION: UMBILICAL

ABBREVIATION.....UMBILICAL
ACTIVE.....Yes
DESCRIPTION.....UMBILICAL
DISPLAY_NAME.....UMBILICAL
HELP_TEXT.....No Help Text
ID.....{874B300-0751-49AF-8211-CF4C009C68B5}
TERM.....UMBILICAL

LOCATION: UPPER ABDOMEN

ABBREVIATION.....UP ABDOM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Upper Abdomen
HELP_TEXT.....No Help Text
ID.....{234C677B-1CE6-498B-9375-BB5DC3958842}
TERM.....UPPER ABDOMEN

LOCATION: UPPER ARM

ABBREVIATION.....UpA
ACTIVE.....Yes
DESCRIPTION.....Area of arm between elbow and shoulder
DISPLAY_NAME.....Upper Arm
HELP_TEXT.....UPPER ARM
ID.....{B43A2B31-1F35-3F4D-5D23-882541ABE477}
TERM.....UPPER ARM

LOCATION: UPPER BACK

ABBREVIATION.....UPPER BACK
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Upper Back
HELP_TEXT.....No Help Text
ID.....{C1DA4294-E0DC-438B-810D-7F91E170244A}
TERM.....UPPER BACK

Page: 961

LOCATION: WAIST

ABBREVIATION.....WAIST
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Waist
HELP_TEXT.....No Help Text
ID.....{784DD311-8B48-4FE3-96B6-471CCE453CD5}
TERM.....WAIST

LOCATION: WRIST

ABBREVIATION.....Wri
ACTIVE.....Yes
DESCRIPTION.....Joint between hand and forearm
DISPLAY_NAME.....Wrist
HELP_TEXT.....WRIST
ID.....{C4F1453C-A20A-E48C-0BD8-FAF5E5568ADA}
TERM.....WRIST

=====

TYPE: METHOD

Page: 962

=====

METHOD: ADULT

ABBREVIATION.....Ad
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Adult
HELP_TEXT.....No Help Text
ID.....{C2D49445-6556-8E6E-A8E9-CB3FA76BC22E}
TERM.....ADULT

METHOD: ADULT CUFF

ABBREVIATION.....AC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Adult Cuff
HELP_TEXT.....No Help Text

ID.....{1222302E-6982-C627-8A0B-A651C98138E3}
TERM.....ADULT CUFF

METHOD: AEROSOL/HUMIDIFIED MASK

ABBREVIATION.....AM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Aerosol/Humidified Mask
HELP_TEXT.....No Help Text
ID.....{CAF54B28-D5D2-8481-FAE3-9C74C2D51551}
TERM.....AEROSOL/HUMIDIFIED MASK

METHOD: AFTER EXERCISE

ABBREVIATION.....AE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....After Exercise
HELP_TEXT.....No Help Text
ID.....{D129EF4B-3A99-9966-54F1-2C01A36814EE}
TERM.....AFTER EXERCISE

METHOD: AMBU BAG

ABBREVIATION.....AB
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Ambu Bag
HELP_TEXT.....No Help Text
ID.....{E9D82DB9-D8A3-10E4-CE51-34BFF48F4603}
TERM.....AMBU BAG

Page: 963

METHOD: ARTERIAL LINE

ABBREVIATION.....AL
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Arterial Line
HELP_TEXT.....No Help Text
ID.....{D4B3F4ED-8B33-431E-4164-BD165E823776}
TERM.....ARTERIAL LINE

METHOD: ASSISTED VENTILATOR

ABBREVIATION.....AV
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Assisted Ventilator
HELP_TEXT.....No Help Text
ID.....{A540F333-5A09-5F21-7B51-2B5DF05E32C0}
TERM.....ASSISTED VENTILATOR

METHOD: AT REST

ABBREVIATION.....AR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....At Rest
HELP_TEXT.....No Help Text
ID.....{E8C289A6-0C3D-A7B2-BB85-7F1A23F22B50}
TERM.....AT REST

METHOD: AUSCULTATE

ABBREVIATION.....Aus
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Auscultate
HELP_TEXT.....No Help Text
ID.....{0C629C1F-92D3-C641-B02C-9EA43B401DC9}
TERM.....AUSCULTATE

METHOD: BED

ABBREVIATION.....B
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bed
HELP_TEXT.....No Help Text
ID.....{20188D9F-3815-1F7F-34EC-BA77F75B310E}
TERM.....BED

Page: 964

METHOD: BIPAP

ABBREVIATION.....BIP
ACTIVE.....Yes
DESCRIPTION.....Biphasic Positive Airway Pressure
DISPLAY_NAME.....BiPAP
HELP_TEXT.....No Help Text
ID.....{3CFD1D11-EF8D-58BA-6381-B24956344D2E}
TERM.....BIPAP

METHOD: BIPAP-CPAP

ABBREVIATION.....PAP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....BiPAP-CPAP
HELP_TEXT.....No Help Text
ID.....{2BCBE562-8F44-662F-469B-13AC0E863A3D}
TERM.....BIPAP-CPAP

METHOD: BLADDER

ABBREVIATION.....BLAD

ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bladder
HELP_TEXT.....No Help Text
ID.....{E3B8D8AC-2846-B84B-75D9-7F50B55B4ADD}
TERM.....BLADDER

METHOD: CARDIAC MONITOR

ABBREVIATION.....CARDIAC MO
ACTIVE.....Yes
DESCRIPTION.....Cardiac Monitor
DISPLAY_NAME.....Cardiac Monitor
HELP_TEXT.....Cardiac Monitor
ID.....{E8B6B5A3-D143-3967-A362-019BD5EB6984}
TERM.....CARDIAC MONITOR

METHOD: CARDIAC OUTPUT MONITOR

ABBREVIATION.....COM
ACTIVE.....Yes
DESCRIPTION.....A CO Monitor
DISPLAY_NAME.....Cardiac Output Monitor
HELP_TEXT.....Data from a Cardiac Output Monitor
ID.....{CCBE5744-7CB9-9DFA-938F-FDDBC73B45B3}
TERM.....CARDIAC OUTPUT MONITOR

Page: 965

METHOD: CATHETER

ABBREVIATION.....CATHETER
ACTIVE.....Yes
DESCRIPTION.....Catheter used to capture Intake and
Output
DISPLAY_NAME.....Catheter - Method
HELP_TEXT.....No Help Text
ID.....{1AA95835-FB32-4D1A-9C53-271EE38B2F63}
TERM.....CATHETER

METHOD: CHAIR

ABBREVIATION.....Ch
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Chair
HELP_TEXT.....No Help Text
ID.....{BBF5C05C-220D-D052-C7E8-236E791A4FC9}
TERM.....CHAIR

METHOD: CHEST TUBE

ABBREVIATION.....CH TUBE M
ACTIVE.....Yes
DESCRIPTION.....Chest Tube

DISPLAY_NAME.....Chest Tube
HELP_TEXT.....Chest Tube
ID.....{F4F467D3-BB57-4E6F-9EA0-C0FBFB66F5B1}
TERM.....CHEST TUBE

METHOD: COLLOID

ABBREVIATION.....COLLOID M
ACTIVE.....Yes
DESCRIPTION.....Colloid as method for Intake and Output
DISPLAY_NAME.....Colloid - Method
HELP_TEXT.....No Help Text
ID.....{2AE33508-5EA2-44DF-A401-74DD7B93AB6E}
TERM.....COLLOID

METHOD: COMPUTED TOMOGRAPHY

ABBREVIATION.....COM TOMOGR
ACTIVE.....Yes
DESCRIPTION.....Computed Tomography
DISPLAY_NAME.....Computed Tomography
HELP_TEXT.....No Help Text
ID.....{B4676AE0-BFAF-4649-AB08-3710E39F9B8D}
TERM.....COMPUTED TOMOGRAPHY

Page: 966

METHOD: CONTRAST LEFT VENTRICULOGRAPHY

ABBREVIATION.....CON L VENT
ACTIVE.....Yes
DESCRIPTION.....Contrast Left Ventriculography
DISPLAY_NAME.....Contrast Left Ventriculography
HELP_TEXT.....No Help Text
ID.....{35BEA63F-9696-4C03-AE1D-D76E4DDBB017}
TERM.....CONTRAST LEFT VENTRICULOGRAPHY

METHOD: CONTROLLED VENTILATOR

ABBREVIATION.....CV
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Controlled Ventilator
HELP_TEXT.....No Help Text
ID.....{0CFC2183-5245-E7AB-96ED-93AF4A3FEC64}
TERM.....CONTROLLED VENTILATOR

METHOD: CPAP

ABBREVIATION.....CP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Cpap
HELP_TEXT.....No Help Text
ID.....{A041FB71-08EF-F3FE-82FA-051FF74C824F}

TERM.....CPAP

METHOD: CRRT

ABBREVIATION.....CRRT
ACTIVE.....Yes
DESCRIPTION.....Continuous Renal Replacement Therapy
DISPLAY_NAME.....CRRT
HELP_TEXT.....No Help Text
ID.....{AB37FC1D-DD7A-4103-A777-7B8D0F876889}
TERM.....CRRT

METHOD: CUFF

ABBREVIATION.....Cu
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Cuff
HELP_TEXT.....No Help Text
ID.....{E59A2543-1860-1E1C-BEF1-6E7F6B6FBC5E}
TERM.....CUFF

Page: 967

METHOD: CUFF-AUTOMATED

ABBREVIATION.....CA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Cuff-Automated
HELP_TEXT.....No Help Text
ID.....{088AAD2B-2E7B-F835-376F-3D4340790976}
TERM.....CUFF-AUTOMATED

METHOD: CUFF-MANUAL

ABBREVIATION.....CM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Cuff-Manual
HELP_TEXT.....No Help Text
ID.....{5F56C163-2B04-112F-55DE-E31374E0FB58}
TERM.....CUFF-MANUAL

METHOD: CVVHD

ABBREVIATION.....CVVHD
ACTIVE.....Yes
DESCRIPTION.....Central Venous - Venous Hemodialysis
DISPLAY_NAME.....CVVHD
HELP_TEXT.....No Help Text
ID.....{4CC9AB0D-0B6C-43F9-A93E-A8228EC3A77F}
TERM.....CVVHD

METHOD: DIALYSIS

ABBREVIATION.....DIALYSIS M
ACTIVE.....Yes
DESCRIPTION.....Dialysis used to describe Intake and
Output
DISPLAY_NAME.....Dialysis - Method
HELP_TEXT.....No Help Text
ID.....{1EAE28D4-3AC0-4DA1-9D09-3A65C38267F6}
TERM.....DIALYSIS

METHOD: DIALYSIS CATHETER

ABBREVIATION.....DIALYSIS C
ACTIVE.....Yes
DESCRIPTION.....Dialysis Catheter
DISPLAY_NAME.....Dialysis Catheter
HELP_TEXT.....No Help Text
ID.....{D9FF5CA8-09AF-4318-B784-DEB3A8F819D8}
TERM.....DIALYSIS CATHETER

Page: 968

METHOD: DOPPLER

ABBREVIATION.....Dop
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Doppler
HELP_TEXT.....No Help Text
ID.....{52683C19-F9DC-876B-CA3F-53E7C864CF04}
TERM.....DOPPLER

METHOD: DOUBLE LUMEN CATHETER

ABBREVIATION.....DOUBLE LUM
ACTIVE.....Yes
DESCRIPTION.....Double Lumen Catheter
DISPLAY_NAME.....Double Lumen Catheter
HELP_TEXT.....No Help Text
ID.....{ED646D10-3DBE-4049-ADC6-5E3DD2C351DD}
TERM.....DOUBLE LUMEN CATHETER

METHOD: ECHOCARDIOGRAPHY

ABBREVIATION.....ECG
ACTIVE.....Yes
DESCRIPTION.....Echocardiography
DISPLAY_NAME.....Echocardiography
HELP_TEXT.....No Help Text
ID.....{CA44697E-3BBA-4767-A30D-6548A72F9639}
TERM.....ECHOCARDIOGRAPHY

METHOD: ELECTROCARDIOGRAM

ABBREVIATION.....ECG

ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Electrocardiogram
HELP_TEXT.....No Help Text
ID.....{CF5F1AB9-3B74-7FBA-E38B-C9EF54FE4DC5}
TERM.....ELECTROCARDIOGRAM

METHOD: ELECTROCARDIOGRAM TRACE 1

ABBREVIATION.....ECG TR 1
ACTIVE.....Yes
DESCRIPTION.....Trace 1 from ECG
DISPLAY_NAME.....ECG Trace 1
HELP_TEXT.....No Help Text
ID.....{E0D449BA-9F06-4A27-8FE9-24C1BA0746DB}
TERM.....ELECTROCARDIOGRAM TRACE 1

Page: 969

METHOD: ELECTROCARDIOGRAM TRACE 2

ABBREVIATION.....ECG TR 2
ACTIVE.....Yes
DESCRIPTION.....Trace 2 from ECG
DISPLAY_NAME.....ECG Trace 2
HELP_TEXT.....No Help Text
ID.....{73EACBC1-C1FB-4B2C-96AF-C5C1CA07DA72}
TERM.....ELECTROCARDIOGRAM TRACE 2

METHOD: ENTERAL FEEDING / SUPPLEMENT

ABBREVIATION.....EF/SUPPLEM
ACTIVE.....Yes
DESCRIPTION.....Enteral Feeding / Supplement - Method
DISPLAY_NAME.....Enteral Feeding / Supplement - Method
HELP_TEXT.....No Help Text
ID.....{8389A9F5-0D48-496D-A605-1560B055DCC0}
TERM.....ENTERAL FEEDING / SUPPLEMENT

METHOD: ENTERAL FEEDING TUBE

ABBREVIATION.....ENT FEED T
ACTIVE.....Yes
DESCRIPTION.....EFT for describing Intakes and Outputs
DISPLAY_NAME.....Enteral Feeding Tube - Method
HELP_TEXT.....No Help Text
ID.....{487E3F10-E1F6-4EFF-9C50-9EC5BB7AE3CF}
TERM.....ENTERAL FEEDING TUBE

METHOD: FACE TENT

ABBREVIATION.....FT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Face Tent

HELP_TEXT.....No Help Text
ID.....{1329BF1A-DCFD-9AF3-69DB-BA6D4BC3F52A}
TERM.....FACE TENT

METHOD: GI TUBE

ABBREVIATION.....GI TUBE
ACTIVE.....Yes
DESCRIPTION.....Used to describe Intakes and Outputs
DISPLAY_NAME.....GI Tube - Method
HELP_TEXT.....No Help Text
ID.....{AF64F7BE-E72E-4391-BC1F-48163F2F3751}
TERM.....GI TUBE

Page: 970

METHOD: HEMODIALYSIS

ABBREVIATION.....HEMODIALYS
ACTIVE.....Yes
DESCRIPTION.....Hemodialysis
DISPLAY_NAME.....Hemodialysis
HELP_TEXT.....No Help Text
ID.....{DE7CD394-9E69-46B0-A207-5A18F7C5A1A7}
TERM.....HEMODIALYSIS

METHOD: HEMOVAC

ABBREVIATION.....HEMOVAC
ACTIVE.....Yes
DESCRIPTION.....Hemovac Drain
DISPLAY_NAME.....Hemovac Drain
HELP_TEXT.....No Help Text
ID.....{79DA0646-B811-4954-8BE1-C133398CD574}
TERM.....HEMOVAC

METHOD: IABP

ABBREVIATION.....IABP
ACTIVE.....Yes
DESCRIPTION.....Intra-aortic Balloon Pump
DISPLAY_NAME.....IABP
HELP_TEXT.....Intra-aortic Balloon Pump
ID.....{4352A41F-AAF3-4544-9590-BAFD7977DF49}
TERM.....IABP

METHOD: INTRODUCER CATHETER

ABBREVIATION.....INTRODUCER
ACTIVE.....Yes
DESCRIPTION.....Introducer Catheter
DISPLAY_NAME.....Introducer Catheter
HELP_TEXT.....No Help Text
ID.....{60EED7D3-6CB2-4FCE-BFC0-DE0BA8AB0814}
TERM.....INTRODUCER CATHETER

METHOD: INVASIVE

ABBREVIATION.....I
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Invasive
HELP_TEXT.....No Help Text
ID.....{5DB8746E-5446-F47A-883C-A44B91275A6E}
TERM.....INVASIVE

Page: 971

METHOD: IRRIGATION

ABBREVIATION.....IRRIGATION
ACTIVE.....Yes
DESCRIPTION.....Irrigation
DISPLAY_NAME.....Irrigation
HELP_TEXT.....No Help Text
ID.....{6C968D0A-400F-448A-93EF-02CB891BA981}
TERM.....IRRIGATION

METHOD: IV

ABBREVIATION.....IV
ACTIVE.....Yes
DESCRIPTION.....Method to capture Intakes and Outputs
DISPLAY_NAME.....IV
HELP_TEXT.....No Help Text
ID.....{52141E2F-7942-4C60-9B08-62F3B4640020}
TERM.....IV

METHOD: JP DRAIN

ABBREVIATION.....JP DRAIN
ACTIVE.....Yes
DESCRIPTION.....JP Drain
DISPLAY_NAME.....JP Drain
HELP_TEXT.....No Help Text
ID.....{94A32247-6190-4047-B581-A212FE4C9453}
TERM.....JP DRAIN

METHOD: LARYNGEAL MASK

ABBREVIATION.....LARYNGMASK
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Laryngeal Mask
HELP_TEXT.....No Help Text
ID.....{DE9B1E53-711C-DA9D-59D9-CD0EBB6CD6BE}
TERM.....LARYNGEAL MASK

METHOD: LEAD

ABBREVIATION.....L
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Lead
HELP_TEXT.....No Help Text
ID.....{B628FDCE-05E0-DACF-E66E-0F8380AB465A}
TERM.....LEAD

Page: 972

METHOD: LG ADULT

ABBREVIATION.....LgA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Lg Adult
HELP_TEXT.....No Help Text
ID.....{0C7B9E22-A426-80D1-80E8-2082F15784BA}
TERM.....LG ADULT

METHOD: LG ADULT CUFF

ABBREVIATION.....LAC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Lg Adult Cuff
HELP_TEXT.....No Help Text
ID.....{4FA9374C-AE66-DF6E-E0E8-673D7B9B2E19}
TERM.....LG ADULT CUFF

METHOD: LIFT SCALE

ABBREVIATION.....LS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Lift Scale
HELP_TEXT.....No Help Text
ID.....{8D741461-25C5-3955-BE2C-CBD27317CBA1}
TERM.....LIFT SCALE

METHOD: MAGNETIC RESONANCE IMAGING

ABBREVIATION.....MRI
ACTIVE.....Yes
DESCRIPTION.....Magnetic Resonance Imaging
DISPLAY_NAME.....MRI
HELP_TEXT.....No Help Text
ID.....{50EBFFF8-DF0D-41F2-B11B-33A26AA40DBA}
TERM.....MAGNETIC RESONANCE IMAGING

METHOD: MASK

ABBREVIATION.....M
ACTIVE.....Yes
DESCRIPTION.....No Description

DISPLAY_NAME.....Mask
HELP_TEXT.....No Help Text
ID.....{4FB6B24A-D56E-4EE3-6CBC-E54C9734424A}
TERM.....MASK

Page: 973

METHOD: MONITOR

ABBREVIATION.....Mo
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Monitor
HELP_TEXT.....No Help Text
ID.....{63FC24EA-50B8-B302-D948-ABE3517CF6FC}
TERM.....MONITOR

METHOD: MULTIGATED NUCLEAR VENTRICULOGRAPHY

ABBREVIATION.....MUGA
ACTIVE.....Yes
DESCRIPTION.....Multigated Nuclear Ventriculography
DISPLAY_NAME.....MUGA
HELP_TEXT.....No Help Text
ID.....{4D73168B-4DF4-4EA2-83B9-ACE853D5DD42}
TERM.....MULTIGATED NUCLEAR VENTRICULOGRAPHY

METHOD: MUSTACHE CANNULA

ABBREVIATION.....MC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Mustache Cannula
HELP_TEXT.....No Help Text
ID.....{DAB3F52F-A8E1-2D9E-C8E5-1D005EC679FC}
TERM.....MUSTACHE CANNULA

METHOD: NASAL CANNULA

ABBREVIATION.....NC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Nasal Cannula
HELP_TEXT.....No Help Text
ID.....{EFFD38AE-6AF5-77FD-72E1-E0F7E3DEBA29}
TERM.....NASAL CANNULA

METHOD: NEBULIZER

ABBREVIATION.....NEBULIZER
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Nebulizer
HELP_TEXT.....No Help Text
ID.....{3CF184A3-2688-4988-B956-D5B0F303E8E8}

TERM.....NEBULIZER

Page: 974

METHOD: NON RE-BREATHER

ABBREVIATION.....NRB
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Non Re-Breather
HELP_TEXT.....No Help Text
ID.....{1B9AA055-FD7D-7DCC-4EDB-BDE3C821AA8E}
TERM.....NON RE-BREATHER

METHOD: NON-INVASIVE

ABBREVIATION.....N-I
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Non-Invasive
HELP_TEXT.....No Help Text
ID.....{F50ABFED-5FD7-FC95-20FE-1E3038C9F299}
TERM.....NON-INVASIVE

METHOD: NONINVASIVE

ABBREVIATION.....N
ACTIVE.....No
DESCRIPTION.....Noninvasive
DISPLAY_NAME.....Noninvasive
HELP_TEXT.....No Help Text
ID.....{398E91DF-F425-5083-AA35-C5263792845E}
TERM.....NONINVASIVE

METHOD: NUCLEAR PERFUSION IMAGING

ABBREVIATION.....NPI
ACTIVE.....Yes
DESCRIPTION.....Nuclear Perfusion Imaging
DISPLAY_NAME.....Nuclear Perfusion Imaging
HELP_TEXT.....No Help Text
ID.....{D0F7DB21-0D66-4161-8214-DACCAF41ADAC}
TERM.....NUCLEAR PERFUSION IMAGING

METHOD: O2 SATS

ABBREVIATION.....O2
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....O2 Sats
HELP_TEXT.....No Help Text
ID.....{66C2100C-8A52-F36A-F7F1-A55D4E4D7CE6}
TERM.....O2 SATS

METHOD: OTHER

ABBREVIATION.....OTHER
ACTIVE.....Yes
DESCRIPTION.....Other: Comments Required
DISPLAY_NAME.....OTHER: OBSERVATION COMMENTS REQUIRED
HELP_TEXT.....No Help Text
ID.....{AED275BD-8A42-46BE-BCDF-472EF4C31D65}
TERM.....OTHER

METHOD: P.O.

ABBREVIATION.....PO
ACTIVE.....Yes
DESCRIPTION.....By Mouth, Oral, "Per Os"
DISPLAY_NAME.....P.O. - Method
HELP_TEXT.....No Help Text
ID.....{7D50D8E7-DA7B-4605-AB43-1E97DEBEACEE}
TERM.....P.O.

METHOD: PA CATHETER

ABBREVIATION.....PA CATHETE
ACTIVE.....Yes
DESCRIPTION.....Pulmonary Catheter
DISPLAY_NAME.....PA Catheter
HELP_TEXT.....No Help Text
ID.....{5A518696-41B7-4AB8-8FC7-5CF454CE0299}
TERM.....PA CATHETER

METHOD: PA LINE

ABBREVIATION.....PALI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....PA Line
HELP_TEXT.....No Help Text
ID.....{1AEA41D7-901B-96FF-BF9B-420C24FF00E4}
TERM.....PA LINE

METHOD: PALPATED

ABBREVIATION.....Pal
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Palpated
HELP_TEXT.....No Help Text
ID.....{295D8E72-12A1-08AF-6BE8-0F9CD238F5EF}
TERM.....PALPATED

METHOD: PARTIAL RE-BREATH

ABBREVIATION.....PRB
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Partial Re-Breather
HELP_TEXT.....No Help Text
ID.....{6F1DF170-5128-2E3F-1831-12FB0BB5E95C}
TERM.....PARTIAL RE-BREATHER

METHOD: PEDIATRIC CUFF

ABBREVIATION.....PdC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Pediatric Cuff
HELP_TEXT.....No Help Text
ID.....{CA9EA220-EF92-EE77-FBEA-4E49C4730A3E}
TERM.....PEDIATRIC CUFF

METHOD: PENDANT CANNULA

ABBREVIATION.....PC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Pendant Cannula
HELP_TEXT.....No Help Text
ID.....{A1451228-4ABB-DBCB-31AD-578BCD93012F}
TERM.....PENDANT CANNULA

METHOD: PENROSE

ABBREVIATION.....PENROSE
ACTIVE.....Yes
DESCRIPTION.....Penrose Drain
DISPLAY_NAME.....Penrose Drain
HELP_TEXT.....No Help Text
ID.....{CBC3C6DE-E5DE-41D5-83AB-30A4F3CFA9E1}
TERM.....PENROSE

METHOD: PERITONEAL DIALYSIS

ABBREVIATION.....PERI DIALY
ACTIVE.....Yes
DESCRIPTION.....Peritoneal Dialysis
DISPLAY_NAME.....Peritoneal Dialysis
HELP_TEXT.....No Help Text
ID.....{6B4AA38B-AFCC-42CE-9F32-59B32E8A6BBB}
TERM.....PERITONEAL DIALYSIS

METHOD: PICC LINE

ABBREVIATION.....PICC LINE
ACTIVE.....Yes

DESCRIPTION.....Peripherally Inserted Central
Catheter Line
DISPLAY_NAME.....PICC Line
HELP_TEXT.....No Help Text
ID.....{BDB5075B-DECC-443A-B774-CFA062DF34F4}
TERM.....PICC LINE

METHOD: PULSE OX MONITOR

ABBREVIATION.....POM
ACTIVE.....Yes
DESCRIPTION.....A Pulse Ox Monitor
DISPLAY_NAME.....Pulse Ox Monitor
HELP_TEXT.....Readings that have come from a Pulse
Ox Monitor
ID.....{A7E3E90C-4C8B-87D7-C2F8-49D9F3846CE8}
TERM.....PULSE OX MONITOR

METHOD: RECTAL TEMP PROBE

ABBREVIATION.....RTP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Rectal Temp Probe
HELP_TEXT.....No Help Text
ID.....{0100F644-EF8E-5FE3-0022-804408DCE2DD}
TERM.....RECTAL TEMP PROBE

METHOD: RESERVOIR CANNULA

ABBREVIATION.....RC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Reservoir Cannula
HELP_TEXT.....No Help Text
ID.....{B332C44D-29BE-765C-932C-63EBFC9631FC}
TERM.....RESERVOIR CANNULA

METHOD: RHYTHM

ABBREVIATION.....RHY
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Rhythm
HELP_TEXT.....No Help Text
ID.....{8A4AC126-12DA-0904-C364-B5D818621E32}
TERM.....RHYTHM

Page: 978

METHOD: ROOM AIR

ABBREVIATION.....ROOM AIR
ACTIVE.....Yes
DESCRIPTION.....No Description

DISPLAY_NAME.....Room Air
HELP_TEXT.....No Help Text
ID.....{B4103F11-4BCE-49ED-AF3E-ACCB526D938C}
TERM.....ROOM AIR

METHOD: ROOM AIR15

ABBREVIATION.....RA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Room Air15
HELP_TEXT.....No Help Text
ID.....{55AC16E4-4BAD-EDE4-41E9-853C78D8474B}
TERM.....ROOM AIR15

METHOD: ROOM AIR20

ABBREVIATION.....xx
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Room Air20
HELP_TEXT.....No Help Text
ID.....{FBBD3CC2-53AF-FBCC-EF16-A6300F7F8330}
TERM.....ROOM AIR20

METHOD: SINGLE LUMEN CATHETER

ABBREVIATION.....SING LUM C
ACTIVE.....Yes
DESCRIPTION.....Single Lumen Catheter
DISPLAY_NAME.....Single Lumen Catheter
HELP_TEXT.....No Help Text
ID.....{E5B432BD-C698-4415-98E1-F5F229F54F8A}
TERM.....SINGLE LUMEN CATHETER

METHOD: SKIN (ADHESIVE INDICATOR)

ABBREVIATION.....SKIN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Skin (Adhesive Indicator)
HELP_TEXT.....No Help Text
ID.....{C7CE7D84-6FC7-406A-7E56-B080637509DF}
TERM.....SKIN (ADHESIVE INDICATOR)

Page: 979

METHOD: SM ADULT

ABBREVIATION.....SmA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Sm Adult
HELP_TEXT.....No Help Text
ID.....{CCB94C2F-CE55-A046-3A21-3ADDB1D23D3C}

TERM.....SM ADULT

METHOD: SM ADULT CUFF

ABBREVIATION.....SAC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Sm Adult Cuff
HELP_TEXT.....No Help Text
ID.....{0346341B-6D38-6599-ED9F-24C81A8698A8}
TERM.....SM ADULT CUFF

METHOD: SPONTANEOUS

ABBREVIATION.....S
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Spontaneous
HELP_TEXT.....No Help Text
ID.....{5E9463BC-C871-54FD-F800-DCFDEDB8A93B}
TERM.....SPONTANEOUS

METHOD: STANDING WEIGHT

ABBREVIATION.....SW
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Standing Weight
HELP_TEXT.....No Help Text
ID.....{BE16FA09-16DA-510B-D235-B7F0F93AB105}
TERM.....STANDING WEIGHT

METHOD: T-PIECE

ABBREVIATION.....TP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....T-Piece
HELP_TEXT.....No Help Text
ID.....{27FC4BC6-2F25-2154-553C-2B452E28A91C}
TERM.....T-PIECE

Page: 980

METHOD: THIGH CUFF

ABBREVIATION.....ThC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Thigh Cuff
HELP_TEXT.....No Help Text
ID.....{C581DFAC-691C-B64A-F058-AF2AA5050147}
TERM.....THIGH CUFF

METHOD: TPN

ABBREVIATION.....TPN
ACTIVE.....Yes
DESCRIPTION.....Total Parenteral Nutrition Method
DISPLAY_NAME.....TPN - Method
HELP_TEXT.....No Help Text
ID.....{B0BE1126-1C5B-4B7E-9B82-5A58AAF63FD2}
TERM.....TPN

METHOD: TRACHEOSTOMY COLLAR

ABBREVIATION.....TC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Tracheostomy Collar
HELP_TEXT.....No Help Text
ID.....{CBF8A3B9-D323-E021-FD1B-BF4CE19AD20E}
TERM.....TRACHEOSTOMY COLLAR

METHOD: TRANSFUSION

ABBREVIATION.....TRANSFUSIO
ACTIVE.....Yes
DESCRIPTION.....Transfusion - Method
DISPLAY_NAME.....Transfusion - Method
HELP_TEXT.....No Help Text
ID.....{9CBA9E7F-C0E2-4936-8AD1-04E2BD1D770F}
TERM.....TRANSFUSION

METHOD: TRANSTRACHEAL

ABBREVIATION.....Tt
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Transtracheal
HELP_TEXT.....No Help Text
ID.....{48951129-CEEC-10D9-F746-99B4713B9C64}
TERM.....TRANSTRACHEAL

Page: 981

METHOD: TRIPLE LUMEN CATHETER

ABBREVIATION.....TRIPLE LUM
ACTIVE.....Yes
DESCRIPTION.....Triple Lumen Catheter
DISPLAY_NAME.....Triple Lumen Catheter
HELP_TEXT.....No Help Text
ID.....{704E2B3E-5E5D-43FB-B956-CFF90AF7D774}
TERM.....TRIPLE LUMEN CATHETER

METHOD: VENTILATOR

ABBREVIATION.....Ven
ACTIVE.....Yes

DESCRIPTION.....No Description
DISPLAY_NAME.....Ventilator
HELP_TEXT.....No Help Text
ID.....{416FCBD3-5370-15B9-C480-48855ED486D9}
TERM.....VENTILATOR

METHOD: VENTRICULAR DRAIN

ABBREVIATION.....VENT DRAIN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Ventricular Drain
HELP_TEXT.....No Help Text
ID.....{5FF035D1-9C97-415D-BBEB-943506395A0C}
TERM.....VENTRICULAR DRAIN

METHOD: VENTURI MASK

ABBREVIATION.....VM
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Venturi Mask
HELP_TEXT.....No Help Text
ID.....{E0401675-9C4E-2AFA-2B93-F352B781E963}
TERM.....VENTURI MASK

METHOD: WHEELCHAIR SCALE

ABBREVIATION.....WcS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Wheelchair Scale
HELP_TEXT.....No Help Text
ID.....{8C496753-9F42-2A36-E5CA-750188306ED4}
TERM.....WHEELCHAIR SCALE

Page: 982

METHOD: WITH ACTIVITY

ABBREVIATION.....Act
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....With Activity
HELP_TEXT.....No Help Text
ID.....{544A262B-EA34-F82A-46A9-F80740515D65}
TERM.....WITH ACTIVITY

METHOD: WITH AMBULATION

ABBREVIATION.....WA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....With Ambulation
HELP_TEXT.....No Help Text

ID.....{4E91DE6A-359E-5AF4-AB49-AFAF735E1B94}
TERM.....WITH AMBULATION

METHOD: WITH CAST OR BRACE

ABBREVIATION.....WCB
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....With Cast or Brace
HELP_TEXT.....No Help Text
ID.....{47709C05-DF21-FB7D-DE7F-0699D1B7E457}
TERM.....WITH CAST OR BRACE

METHOD: WITH PROSTHESIS

ABBREVIATION.....WP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....With Prosthesis
HELP_TEXT.....No Help Text
ID.....{73C93D2D-934B-CA55-3557-E756A89710F5}
TERM.....WITH PROSTHESIS

METHOD: WITHOUT PROSTHESIS

ABBREVIATION.....WOP
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Without Prosthesis
HELP_TEXT.....No Help Text
ID.....{9AACA852-E70D-1048-FCE7-112D46C12619}
TERM.....WITHOUT PROSTHESIS

Page: 983

=====

TYPE: NULL FLAVOR

Page: 984

=====

NULL FLAVOR: PATIENT ON PASS

ABBREVIATION.....op
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....On Pass
HELP_TEXT.....No Help Text
ID.....{3D2D2E52-D3F2-E02D-B9C7-CEB7C46C807D}
TERM.....PATIENT ON PASS

NULL FLAVOR: PATIENT REFUSED

ABBREVIATION.....pr
ACTIVE.....Yes
DESCRIPTION.....No Description

DISPLAY_NAME.....Refused
HELP_TEXT.....No Help Text
ID.....{E9676179-B942-3851-37AF-75DAEE9E1F8E}
TERM.....PATIENT REFUSED

=====

TYPE: POSITION

Page: 985

=====

POSITION: L MIDDLE POSTERIOR CHEST

ABBREVIATION.....L M P CHES
ACTIVE.....No
DESCRIPTION.....Left Middle Posterior Chest
DISPLAY_NAME.....L Mid. Pos. Chest
HELP_TEXT.....No Help Text
ID.....{D6D6483C-90A9-4A32-A101-9DAA13BFD7ED}
TERM.....L MIDDLE POSTERIOR CHEST

POSITION: LYING

ABBREVIATION.....Ly
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Lying
HELP_TEXT.....No Help Text
ID.....{DE28D2BF-6F8D-64C3-0092-4CEB5E0902B9}
TERM.....LYING

POSITION: RES EVAL: BED/WHEELCHAIR ALARM

ABBREVIATION.....R.E. ALAR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Bed / Wheelchair Alarm
HELP_TEXT.....No Help Text
ID.....{BDF926A9-2BEF-4E7C-A099-D02BD98CA2F2}
TERM.....RES EVAL: BED/WHEELCHAIR ALARM

POSITION: RES EVAL: REALITY ORIENTATION

ABBREVIATION.....R.E. REALI
ACTIVE.....No
DESCRIPTION.....No Description
DISPLAY_NAME.....Reality Orientation
HELP_TEXT.....No Help Text
ID.....{04D911E1-AA DC-4125-9F93-B298266BA61A}
TERM.....RES EVAL: REALITY ORIENTATION

POSITION: SEMIFOWLERS

ABBREVIATION.....Sfl
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Semifowlers

HELP_TEXT.....No Help Text
ID.....{C2F6A945-3829-AABA-3D67-54D1829D3EB9}
TERM.....SEMIFOWLERS

Page: 986

POSITION: SITTING

ABBREVIATION.....Sit
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Sitting
HELP_TEXT.....No Help Text
ID.....{DD4DCDF1-0A22-2CC5-404F-0E472C73961C}
TERM.....SITTING

POSITION: STANDING

ABBREVIATION.....St
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Standing
HELP_TEXT.....No Help Text
ID.....{0776AE42-C370-B7B3-9B4D-9C914BA8F231}
TERM.....STANDING

POSITION: SUPINE

ABBREVIATION.....SUPINE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Supine
HELP_TEXT.....No Help Text
ID.....{DEA8A744-B445-4A9E-875D-485E060F2ACC}
TERM.....SUPINE

POSITION: TRENDELENBURG

ABBREVIATION.....Tr
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Trendelenburg
HELP_TEXT.....No Help Text
ID.....{42F3353B-7041-EEE1-A87E-941D3BE1CBCD}
TERM.....TRENDELENBURG

POSITION: WHEELCHAIR EXIT

ABBREVIATION.....WHEEL EXIT
ACTIVE.....No
DESCRIPTION.....No Description
DISPLAY_NAME.....Wheelchair Exit
HELP_TEXT.....No Help Text
ID.....{DE13F84A-168F-490D-9AC7-840DAB2D8E22}
TERM.....WHEELCHAIR EXIT

=====

TYPE: PRODUCT

Page: 988

=====

PRODUCT: 25 PERCENT ALBUMIN

ABBREVIATION.....25 PERC AL
ACTIVE.....Yes
DESCRIPTION.....25 Percent Albumin
DISPLAY_NAME.....25 Percent Albumin
HELP_TEXT.....No Help Text
ID.....{CBC7E32F-E457-4A7F-8F27-030EBC6ECBD5}
TERM.....25 PERCENT ALBUMIN

PRODUCT: 5 PERCENT ALBUMIN

ABBREVIATION.....5 PERC ALB
ACTIVE.....Yes
DESCRIPTION.....5 Percent Albumin
DISPLAY_NAME.....5 Percent Albumin
HELP_TEXT.....No Help Text
ID.....{A2D0F9A0-A9EC-4540-9663-B8B19A189198}
TERM.....5 PERCENT ALBUMIN

PRODUCT: 6 PERCENT HETASTARCH

ABBREVIATION.....6 PERC HET
ACTIVE.....Yes
DESCRIPTION.....6 Percent Hetastarch
DISPLAY_NAME.....6 Percent Hetastarch
HELP_TEXT.....No Help Text
ID.....{8A82F108-D2FE-4C21-93BB-344502235AE8}
TERM.....6 PERCENT HETASTARCH

PRODUCT: ALBUMIN

ABBREVIATION.....ALBUMIN
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Albumin
HELP_TEXT.....No Help Text
ID.....{EF4CC3AE-32B2-4EBC-BBDE-7DE7CCEAF081}
TERM.....ALBUMIN

PRODUCT: AMPHOTERICIN B

ABBREVIATION.....AMPHOTER B
ACTIVE.....Yes
DESCRIPTION.....Amphotericin B
DISPLAY_NAME.....Amphotericin B
HELP_TEXT.....Amphotericin B

ID.....{E4ED1488-6A27-4622-A8EE-AB6B53DB86D0}
TERM.....AMPHOTERICIN B

Page: 989

PRODUCT: ANTIBIOTIC

ABBREVIATION.....AB
ACTIVE.....Yes
DESCRIPTION.....Antibiotic
DISPLAY_NAME.....Antibiotic
HELP_TEXT.....No Help Text
ID.....{BC1045A3-FF58-41C0-B48F-C2CFD6E522A7}
TERM.....ANTIBIOTIC

PRODUCT: BACTERIOSTATIC WATER

ABBREVIATION.....BACTER H2O
ACTIVE.....Yes
DESCRIPTION.....Bacteriostatic Water
DISPLAY_NAME.....Bacteriostatic Water
HELP_TEXT.....Bacteriostatic Water
ID.....{EF13BCA6-3AF2-47E9-87A0-82642B529368}
TERM.....BACTERIOSTATIC WATER

PRODUCT: BLADDER IRRIGATION

ABBREVIATION.....BLAD IRRIG
ACTIVE.....Yes
DESCRIPTION.....Bladder Irrigation Output Product
DISPLAY_NAME.....Bladder Irrigation
HELP_TEXT.....No Help Text
ID.....{BCBB106B-ECE4-4232-87C7-D90EDB47F80D}
TERM.....BLADDER IRRIGATION

PRODUCT: BLOOD PRODUCT

ABBREVIATION.....bp
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Blood Product
HELP_TEXT.....No Help Text
ID.....{A910DBA4-1CD7-1302-40B2-5BEE95DD3B08}
TERM.....BLOOD PRODUCT

PRODUCT: CARDIAC OUTPUT INJECTATE

ABBREVIATION.....CARD OUT I
ACTIVE.....Yes
DESCRIPTION.....Cardiac Output Injectate
DISPLAY_NAME.....Cardiac Output Injectate
HELP_TEXT.....No Help Text
ID.....{2128B705-DDF4-4776-87C3-AA82F854E313}
TERM.....CARDIAC OUTPUT INJECTATE

PRODUCT: CHEST TUBE

ABBREVIATION.....CHEST TUBE
ACTIVE.....Yes
DESCRIPTION.....Chest Tube Output
DISPLAY_NAME.....Chest Tube
HELP_TEXT.....No Help Text
ID.....{74C0C672-02AD-45A8-99DB-F5256BA14C72}
TERM.....CHEST TUBE

PRODUCT: CHEST TUBE DRAINAGE

ABBREVIATION.....CHEST T DR
ACTIVE.....Yes
DESCRIPTION.....Chest Tube Drainage
DISPLAY_NAME.....Chest Tube Drainage
HELP_TEXT.....No Help Text
ID.....{591D62C9-D173-45AD-8E05-A500F1C4AAA7}
TERM.....CHEST TUBE DRAINAGE

PRODUCT: COLLOID

ABBREVIATION.....COLLOID
ACTIVE.....Yes
DESCRIPTION.....Colloid
DISPLAY_NAME.....Colloid
HELP_TEXT.....No Help Text
ID.....{5B912A05-9ECA-43C1-A117-BE0843313CDD}
TERM.....COLLOID

PRODUCT: CRYOPRECIPITATE

ABBREVIATION.....CRYOPRECIP
ACTIVE.....Yes
DESCRIPTION.....Cryoprecipitate Blood Product
DISPLAY_NAME.....Cryoprecipitate
HELP_TEXT.....No Help Text
ID.....{5329E737-E059-480D-83F8-1146B091C14E}
TERM.....CRYOPRECIPITATE

PRODUCT: CRYSTALLOID

ABBREVIATION.....CRYSTALLOI
ACTIVE.....Yes
DESCRIPTION.....Crystalloid
DISPLAY_NAME.....Crystalloid
HELP_TEXT.....No Help Text
ID.....{F8B0C9B5-F386-4C88-896A-092C4A95A031}
TERM.....CRYSTALLOID

PRODUCT: DAKINS SOLUTION

ABBREVIATION.....DAKINS SOL
ACTIVE.....Yes
DESCRIPTION.....Dakins Solution
DISPLAY_NAME.....Dakins Solution
HELP_TEXT.....Dakins Solution
ID.....{51651D24-B909-427B-9FF0-3C2859A63D73}
TERM.....DAKINS SOLUTION

PRODUCT: EF RESIDUAL

ABBREVIATION.....EF RESID
ACTIVE.....Yes
DESCRIPTION.....Enteral Feeding Residual
DISPLAY_NAME.....EF Residual
HELP_TEXT.....No Help Text
ID.....{35E269DF-F246-4DF6-9D68-7FD61F8FE088}
TERM.....EF RESIDUAL

PRODUCT: ELEMENTAL FORMULA

ABBREVIATION.....ELE FORM
ACTIVE.....Yes
DESCRIPTION.....Elemental Enteral Nutrition Formula
DISPLAY_NAME.....Elemental Formula
HELP_TEXT.....No Help Text
ID.....{E85DFA33-7F32-4348-9E0E-249AC301EC03}
TERM.....ELEMENTAL FORMULA

PRODUCT: EMESIS

ABBREVIATION.....EMESIS
ACTIVE.....Yes
DESCRIPTION.....Emesis output
DISPLAY_NAME.....Emesis
HELP_TEXT.....No Help Text
ID.....{A7B1605E-1DFE-4FC6-A69B-0E18C1A291FF}
TERM.....EMESIS

PRODUCT: FECES

ABBREVIATION.....f
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Feces
HELP_TEXT.....No Help Text
ID.....{C1CED7CD-EF5D-37D8-4EE0-CC9FD6F20DB0}
TERM.....FECES

PRODUCT: FFP

ABBREVIATION.....FFP

ACTIVE.....Yes
DESCRIPTION.....Fresh Frozen Plasma
DISPLAY_NAME.....FFP
HELP_TEXT.....No Help Text
ID.....{9DD0ECF5-634F-49CF-A529-A2657240A80C}
TERM.....FFP

PRODUCT: FFP - THAWED

ABBREVIATION.....FFP-THAWED
ACTIVE.....Yes
DESCRIPTION.....Fresh Frozen Plasma, Thawed
DISPLAY_NAME.....FFP - THAWED
HELP_TEXT.....No Help Text
ID.....{AAF452A3-0DCA-4FF4-BF96-08E60D2D2875}
TERM.....FFP - THAWED

PRODUCT: FORMULA - HIGH FIBER

ABBREVIATION.....FORM H FIB
ACTIVE.....Yes
DESCRIPTION.....High Fiber Formula
DISPLAY_NAME.....High Fiber Formula
HELP_TEXT.....No Help Text
ID.....{249D8C78-B600-4443-A2E2-39F1B67698DD}
TERM.....FORMULA - HIGH FIBER

PRODUCT: FORMULA - PROBALANCE

ABBREVIATION.....PROBALANCE
ACTIVE.....Yes
DESCRIPTION.....Probalance Formula
DISPLAY_NAME.....Probalance Formula
HELP_TEXT.....No Help Text
ID.....{4EA729FE-08EA-490B-822D-696061E5B55F}
TERM.....FORMULA - PROBALANCE

PRODUCT: FORMULA 1.5 CAL

ABBREVIATION.....1.5 C FORM
ACTIVE.....Yes
DESCRIPTION.....1.5 Cal Enteral Nutrition Formula
DISPLAY_NAME.....1.5 Cal Formula
HELP_TEXT.....No Help Text
ID.....{95FAE60B-BEAE-4EB2-88FC-0A199F4534EB}
TERM.....FORMULA 1.5 CAL

Page: 993

PRODUCT: FORMULA 2.0 CAL

ABBREVIATION.....2.0 C FORM
ACTIVE.....Yes
DESCRIPTION.....2.0 Cal Enteral Nutrition Formula
DISPLAY_NAME.....2.0 Cal Formula

HELP_TEXT.....No Help Text
ID.....{E5FB9590-25AA-435A-B9B4-4812EA15C217}
TERM.....FORMULA 2.0 CAL

PRODUCT: FORMULA D10W ENTERAL NUTRITION

ABBREVIATION.....D10W ENT N
ACTIVE.....Yes
DESCRIPTION.....D10W Enteral Nutrition
DISPLAY_NAME.....D10W Enteral Nutrition
HELP_TEXT.....10% Dextrose and Water solution
ID.....{F8FB55C4-6229-427F-987D-3D469E307D0D}
TERM.....FORMULA D10W ENTERAL NUTRITION

PRODUCT: FORMULA D50W ENTERAL NUTRITION

ABBREVIATION.....D50W ENT N
ACTIVE.....Yes
DESCRIPTION.....D50W Enteral Nutrition
DISPLAY_NAME.....D50W Enteral Nutrition
HELP_TEXT.....50% Dextrose and Water Solution
ID.....{19ED9C83-42C1-42AF-AB4D-1A34B621EFB9}
TERM.....FORMULA D50W ENTERAL NUTRITION

PRODUCT: FORMULA D5W ENTERAL NUTRITION

ABBREVIATION.....D5W ENT NU
ACTIVE.....Yes
DESCRIPTION.....D5W Enteral Nutrition
DISPLAY_NAME.....D5W Enteral Nutrition
HELP_TEXT.....5% Dextrose and Water Solution
ID.....{9EA15A1B-DA1C-49EC-B463-82DD5F125F82}
TERM.....FORMULA D5W ENTERAL NUTRITION

PRODUCT: FORMULA FIBER

ABBREVIATION.....FIBER FORM
ACTIVE.....Yes
DESCRIPTION.....Fiber Enteral Nutrition Formula
DISPLAY_NAME.....Fiber Enteral Nutrition Formula
HELP_TEXT.....No Help Text
ID.....{F6D7E75F-917A-4C50-AABB-BE9DE3C1939D}
TERM.....FORMULA FIBER

PRODUCT: FULL LIQUID

ABBREVIATION.....FL LIQUID
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Full Liquid
HELP_TEXT.....No Help Text
ID.....{C8405270-CCA3-4BFF-B9D0-398590FBB1D9}
TERM.....FULL LIQUID

PRODUCT: GASTRIC DRAINAGE

ABBREVIATION.....GAST DRAIN
ACTIVE.....Yes
DESCRIPTION.....Gastric Drainage
DISPLAY_NAME.....Gastric Drainage
HELP_TEXT.....No Help Text
ID.....{C01D18F9-191A-4D34-B45D-D8D70ACA67BE}
TERM.....GASTRIC DRAINAGE

PRODUCT: GLUCERNA DIABETIC

ABBREVIATION.....GLUCERNA D
ACTIVE.....Yes
DESCRIPTION.....Glucerna Diabetic
DISPLAY_NAME.....Glucerna Diabetic
HELP_TEXT.....No Help Text
ID.....{C2F2CCFE-6D77-49B4-BF5E-F1A285326440}
TERM.....GLUCERNA DIABETIC

PRODUCT: GRANULOCYTES PHERESIS

ABBREVIATION.....GRAN PHERE
ACTIVE.....Yes
DESCRIPTION.....Granulocytes Pheresis Blood Product
DISPLAY_NAME.....Granulocytes Pheresis
HELP_TEXT.....No Help Text
ID.....{2B90878D-9C29-4731-9525-D7DD7FCB513B}
TERM.....GRANULOCYTES PHERESIS

PRODUCT: HEME

ABBREVIATION.....HEME
ACTIVE.....Yes
DESCRIPTION.....Heme output
DISPLAY_NAME.....Heme
HELP_TEXT.....No Help Text
ID.....{A0ADA229-E779-4562-84B5-67121B88C443}
TERM.....HEME

Page: 995

PRODUCT: HEMODIALYSIS

ABBREVIATION.....HEMODIALYS
ACTIVE.....Yes
DESCRIPTION.....Hemodialysis Intake and Output Product
DISPLAY_NAME.....Hemodialysis
HELP_TEXT.....No Help Text
ID.....{8409B357-C41E-4F4A-BB52-F09AD9704C30}
TERM.....HEMODIALYSIS

PRODUCT: HEPATIC FORMULA

ABBREVIATION.....HEPAT FORM
ACTIVE.....Yes
DESCRIPTION.....Hepatic Enteral Nutrition Formula
DISPLAY_NAME.....Hepatic Enteral Nutrition Formula
HELP_TEXT.....No Help Text
ID.....{41DEECAB-168F-4522-852C-C1ED751BFEDC}
TERM.....HEPATIC FORMULA

PRODUCT: HIGH PROTEIN FORMULA

ABBREVIATION.....HIGH PROTE
ACTIVE.....Yes
DESCRIPTION.....High Protein
DISPLAY_NAME.....High Protein Formula
HELP_TEXT.....No Help Text
ID.....{3618C874-3EAE-41CF-BE2B-3DFB960E034F}
TERM.....HIGH PROTEIN FORMULA

PRODUCT: JP DRAIN

ABBREVIATION.....JP DRAIN
ACTIVE.....Yes
DESCRIPTION.....Jackson-Pratt Drain
DISPLAY_NAME.....JP DRAIN
HELP_TEXT.....No Help Text
ID.....{0889246D-AD40-4CFC-A721-BBCB8C34BE19}
TERM.....JP DRAIN

PRODUCT: JUICE

ABBREVIATION.....JUICE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Juice
HELP_TEXT.....No Help-Text
ID.....{91FFF519-1677-BB1D-D959-65B1BFA10285}
TERM.....JUICE

Page: 996

PRODUCT: LINE FLUSH

ABBREVIATION.....LINE FLUSH
ACTIVE.....Yes
DESCRIPTION.....Fluid flushed from line
DISPLAY_NAME.....LINE FLUSH
HELP_TEXT.....No Help Text
ID.....{D3EA8BCD-0EC4-4B39-BA92-BC3FF505C9A4}
TERM.....LINE FLUSH

PRODUCT: LINE FLUSH HEPARIN

ABBREVIATION.....LN FL HEPA
ACTIVE.....Yes
DESCRIPTION.....Heparin used to flush dialysis lines

DISPLAY_NAME.....Line Flush Heparin
HELP_TEXT.....No Help Text
ID.....{E8B68E5A-6B9A-401F-9D1A-E3C907C7B81C}
TERM.....LINE FLUSH HEPARIN

PRODUCT: LIPIDS

ABBREVIATION.....LIPIDS
ACTIVE.....Yes
DESCRIPTION.....Lipids
DISPLAY_NAME.....Lipids
HELP_TEXT.....No Help Text
ID.....{80E456A6-CA3D-4C2D-8DCD-06C63AE1F320}
TERM.....LIPIDS

PRODUCT: MAGNACAL

ABBREVIATION.....MAGNACAL
ACTIVE.....Yes
DESCRIPTION.....Magnacal Enteral Nutrition Type
DISPLAY_NAME.....Magancal
HELP_TEXT.....No Help Text
ID.....{86688633-DFCA-4FBE-99F5-8E0B3FAFE695}
TERM.....MAGNACAL

PRODUCT: MEDIASTINAL

ABBREVIATION.....MEDIASTINA
ACTIVE.....Yes
DESCRIPTION.....Mediastinal Output Product
DISPLAY_NAME.....Mediastinal
HELP_TEXT.....No Help Text
ID.....{A0550B31-B937-4FC7-8F7E-071F75215311}
TERM.....MEDIASTINAL

Page: 997

PRODUCT: MEDICATION

ABBREVIATION.....MEDICATION
ACTIVE.....Yes
DESCRIPTION.....Medication
DISPLAY_NAME.....Medication
HELP_TEXT.....Medication
ID.....{CDD5944E-A0EC-45CE-8EB0-7F9EDF7197E5}
TERM.....MEDICATION

PRODUCT: MODULAR CUSTOMIZED FORMULA

ABBREVIATION.....CUST FORMU
ACTIVE.....Yes
DESCRIPTION.....Modular Customized Formula
DISPLAY_NAME.....Modular Customized Formula
HELP_TEXT.....No Help Text
ID.....{249D095E-9676-4CC5-8A39-097F96A4C032}

TERM.....MODULAR CUSTOMIZED FORMULA

PRODUCT: NEPRO RENAL FORMULA

ABBREVIATION.....NEPR RENAL
ACTIVE.....Yes
DESCRIPTION.....Nepro Renal Enteral Nutrition Formula
DISPLAY_NAME.....Nepro Renal Enteral Nutrition Formula
HELP_TEXT.....No Help Text
ID.....{59C9B396-A705-445B-B566-F030DA4B2B79}
TERM.....NEPRO RENAL FORMULA

PRODUCT: NG RESIDUAL

ABBREVIATION.....NG RESID
ACTIVE.....Yes
DESCRIPTION.....Naso-Gastric Residual Output
DISPLAY_NAME.....NG Residual
HELP_TEXT.....No Help Text
ID.....{B03C9557-84B5-4B05-B784-C87FA90D27A7}
TERM.....NG RESIDUAL

PRODUCT: NOVASOURCE 2.0

ABBREVIATION.....NOVAS 2.0
ACTIVE.....Yes
DESCRIPTION.....Novasource 2.0 Enteral Nutrition
Formula
DISPLAY_NAME.....Novasource 2.0 Enteral Nutrition
Formula
HELP_TEXT.....No Help Text
ID.....{67D16273-6271-43E1-80D6-D12C4D55E024}
TERM.....NOVASOURCE 2.0

Page: 998

PRODUCT: NUTREN 1.5

ABBREVIATION.....NUTREN 1.5
ACTIVE.....Yes
DESCRIPTION.....Nutren 1.5
DISPLAY_NAME.....Nutren 1.5
HELP_TEXT.....No Help Text
ID.....{B93EF897-B9A1-45C6-AC9A-B386554A3186}
TERM.....NUTREN 1.5

PRODUCT: NUTRITIONAL SUPPLEMENT

ABBREVIATION.....NUTR SUPPL
ACTIVE.....Yes
DESCRIPTION.....Nutritional Supplement
DISPLAY_NAME.....Nutritional Supplement
HELP_TEXT.....No Help Text
ID.....{17F8F1DA-4FC6-48E5-9E6E-1F67CD06D200}
TERM.....NUTRITIONAL SUPPLEMENT

PRODUCT: OSMOLITE

ABBREVIATION.....OSMOLITE
ACTIVE.....Yes
DESCRIPTION.....Osmolite Enteral Nutrition
DISPLAY_NAME.....Osmolite
HELP_TEXT.....No Help Text
ID.....{18C5843F-BC33-4609-9D58-49F07F12A059}
TERM.....OSMOLITE

PRODUCT: OSMOLITE HN

ABBREVIATION.....OSMOLIT HN
ACTIVE.....Yes
DESCRIPTION.....Osmolite HN Enteral Nutrition
DISPLAY_NAME.....Osmolite HN
HELP_TEXT.....No Help Text
ID.....{3622D5C6-1D7B-4377-94B7-5E33C642C4C5}
TERM.....OSMOLITE HN

PRODUCT: OTHER PO FLUID

ABBREVIATION.....PO OTHER
ACTIVE.....Yes
DESCRIPTION.....PO Fluid other than those listed.
DISPLAY_NAME.....Other - Please Specify
HELP_TEXT.....No Help Text
ID.....{1EA44424-4DE2-42B8-A50F-60DC56204502}
TERM.....OTHER PO FLUID

Page: 999

PRODUCT: PCA

ABBREVIATION.....PCA
ACTIVE.....Yes
DESCRIPTION.....Patient Controlled Analgesia
DISPLAY_NAME.....Patient Controlled Analgesia
HELP_TEXT.....No Help Text
ID.....{D93028EE-9A44-43F5-9C0C-A903BD9B6E6D}
TERM.....PCA

PRODUCT: PLASMAPHERESIS

ABBREVIATION.....PLASMAPHER
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Plasmapheresis
HELP_TEXT.....No Help Text
ID.....{E1214367-8BB0-4FAA-A916-EC707569A9C4}
TERM.....PLASMAPHERESIS

PRODUCT: PLATELETS

ABBREVIATION.....PLATELETS
ACTIVE.....Yes
DESCRIPTION.....Platelets
DISPLAY_NAME.....Platelets
HELP_TEXT.....No Help Text
ID.....{3F4FA2EF-DF05-4A6D-839F-F97D3C2AFCF7}
TERM.....PLATELETS

PRODUCT: PLATELETS PHERESIS

ABBREVIATION.....PLATELET P
ACTIVE.....Yes
DESCRIPTION.....Platelets Pheresis (Apheresis
Platelets)
DISPLAY_NAME.....Platelets Pheresis
HELP_TEXT.....No Help Text
ID.....{65CBDCBD-99DB-44FE-9230-D62886208EBF}
TERM.....PLATELETS PHERESIS

PRODUCT: PLATELETS PHERESIS IRRADIATED

ABBREVIATION.....PLATEL P I
ACTIVE.....Yes
DESCRIPTION.....Platelets Pheresis (Apheresis
platelets) Irradiated
DISPLAY_NAME.....Platelets Pheresis Irradiated
HELP_TEXT.....No Help Text
ID.....{70C7155A-5AB6-4733-B43F-606F8BBB3B4A}
TERM.....PLATELETS PHERESIS IRRADIATED

Page: 1000

PRODUCT: PPN

ABBREVIATION.....PPN
ACTIVE.....Yes
DESCRIPTION.....Peripheral Parenteral Nutrition
DISPLAY_NAME.....PPN
HELP_TEXT.....No Help Text
ID.....{BF402283-3E6F-4DF1-9CC6-4DE35D755C7C}
TERM.....PPN

PRODUCT: PRBC

ABBREVIATION.....PRBC
ACTIVE.....Yes
DESCRIPTION.....Packed Red Blood Cells
DISPLAY_NAME.....PRBC
HELP_TEXT.....No Help Text
ID.....{BCE37EFD-78D6-4A13-A5E4-EA278BF18E02}
TERM.....PRBC

PRODUCT: PRBC - AUTOLOGOUS

ABBREVIATION.....PRBC-A

ACTIVE.....Yes
DESCRIPTION.....Packed Red Blood Cells, Autologous
DISPLAY_NAME.....PRBC - AUTOLOGOUS
HELP_TEXT.....No Help Text
ID.....{1580990C-0B35-4DDD-84AA-682B341C9B76}
TERM.....PRBC - AUTOLOGOUS

PRODUCT: PRBC - LEUKOREDUCED

ABBREVIATION.....PRBC-L
ACTIVE.....Yes
DESCRIPTION.....Packed Red Blood Cells, Leukoreduced
DISPLAY_NAME.....PRBC - Leukoreduced
HELP_TEXT.....No Help Text
ID.....{13DB3170-6C26-4178-B429-C02AFC541AC7}
TERM.....PRBC - LEUKOREDUCED

PRODUCT: PRBC - LEUKOREDUCED - DIRECTED

ABBREVIATION.....PRBC-L-D
ACTIVE.....Yes
DESCRIPTION.....Packed Red Blood Cells, Leukoreduced,
Directed
DISPLAY_NAME.....PRBC - LEUKOREDUCED - DIRECTED
HELP_TEXT.....No Help Text
ID.....{9C54B7F7-FBD9-4FA8-B2A1-893A25A9CD2C}
TERM.....PRBC - LEUKOREDUCED - DIRECTED

Page: 1001

PRODUCT: PRBC - LEUKOREDUCED - IRRADIATED

ABBREVIATION.....PRBC-L-I
ACTIVE.....Yes
DESCRIPTION.....Packed Red Blood Cells, Leukoreduced,
Irradiated
DISPLAY_NAME.....PRBC - LEUKOREDUCED - IRRADIATED
HELP_TEXT.....No Help Text
ID.....{B4431519-0102-49ED-94D6-1738E375C5A4}
TERM.....PRBC - LEUKOREDUCED - IRRADIATED

PRODUCT: PRBC - LEUKOREDUCED - WASHED

ABBREVIATION.....PRBC-L-W
ACTIVE.....Yes
DESCRIPTION.....Packed Red Blood Cells, Leukoreduced,
Washed
DISPLAY_NAME.....PRBC - LEUKOREDUCED - WASHED
HELP_TEXT.....No Help Text
ID.....{BB58F35A-1086-4061-B485-FA35DFB79A0D}
TERM.....PRBC - LEUKOREDUCED - WASHED

PRODUCT: Rh IMMUNE GLOBULIN

ABBREVIATION.....Rh IM GLOB

ACTIVE.....Yes
DESCRIPTION.....Rh Immune Globulin
DISPLAY_NAME.....Rh Immune Globulin
HELP_TEXT.....No Help Text
ID.....{1FE7BD0C-7252-427A-8216-36FCA4637089}
TERM.....Rh IMMUNE GLOBULIN

PRODUCT: SALINE

ABBREVIATION.....SALINE
ACTIVE.....Yes
DESCRIPTION.....Saline
DISPLAY_NAME.....Saline
HELP_TEXT.....No Help Text
ID.....{19D84001-70DB-47AA-9149-A37B7E0BAEF1}
TERM.....SALINE

PRODUCT: SANDOSOURCE PEPTIDE

ABBREVIATION.....SAND PEPTI
ACTIVE.....Yes
DESCRIPTION.....Sandosource Peptide Enteral Nutrition
Formula
DISPLAY_NAME.....Sandosource Peptide Enteral Nutrition
Formula
HELP_TEXT.....No Help Text
ID.....{B543461C-7316-4750-9918-B30E96AB3C95}
TERM.....SANDOSOURCE PEPTIDE

Page: 1002

PRODUCT: SOFT DRINK

ABBREVIATION.....SOFT DRIN
ACTIVE.....Yes
DESCRIPTION.....Soft Drink
DISPLAY_NAME.....Soft Drink
HELP_TEXT.....No Help Text
ID.....{9216F0F9-4D5F-4AF7-8D84-04A2CE11BCBE}
TERM.....SOFT DRINK

PRODUCT: TPN

ABBREVIATION.....TPN
ACTIVE.....Yes
DESCRIPTION.....Total Parenteral Nutrition
DISPLAY_NAME.....TPN
HELP_TEXT.....No Help Text
ID.....{55341B62-E356-4C38-A52D-4F196751E9EB}
TERM.....TPN

PRODUCT: ULTRACAL

ABBREVIATION.....ULTRACAL
ACTIVE.....Yes

DESCRIPTION.....Ultracal Enteral Nutrition
DISPLAY_NAME.....Ultracal Enteral Nutrition
HELP_TEXT.....No Help Text
ID.....{B18E6151-9DC6-4380-899A-E84715BFC61B}
TERM.....ULTRACAL

PRODUCT: URINE

ABBREVIATION.....u
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Urine
HELP_TEXT.....No Help Text
ID.....{1F61B077-1D34-0CBA-E55B-65A8B6A7CAE7}
TERM.....URINE

PRODUCT: VITAL ENTERAL NUTRITION TYPE

ABBREVIATION.....VIT ENT NU
ACTIVE.....Yes
DESCRIPTION.....Vital Enteral Nutrition
DISPLAY_NAME.....Vital Enteral Nutrition
HELP_TEXT.....No Help Text
ID.....{D77B607C-85C4-428E-AAE6-B319F33F4E9F}
TERM.....VITAL ENTERAL NUTRITION TYPE

Page: 1003

PRODUCT: VITAL HN ENTERAL NUTRITION

ABBREVIATION.....VIT NUT HN
ACTIVE.....Yes
DESCRIPTION.....Vital HN Enteral Nutrition
DISPLAY_NAME.....Vital HN Enteral Nutrition
HELP_TEXT.....No Help Text
ID.....{3973F52F-D53F-4D80-928A-BD6154D5C434}
TERM.....VITAL HN ENTERAL NUTRITION

PRODUCT: WATER

ABBREVIATION.....WATER
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Water
HELP_TEXT.....No Help-Text
ID.....{ACCC998E-0F64-510D-B843-3D30B2AD67CF}
TERM.....WATER

PRODUCT: WBC - WASHED

ABBREVIATION.....WBC WASHED
ACTIVE.....Yes
DESCRIPTION.....White Blood Cells, Washed
DISPLAY_NAME.....WBC - Washed
HELP_TEXT.....No Help Text

ID.....{F1460F94-203C-473D-AC7C-92E8ADFFC0F5}
TERM.....WBC - WASHED

PRODUCT: WHOLE BLOOD

ABBREVIATION.....WHOL BLOOD
ACTIVE.....Yes
DESCRIPTION.....Whole Blood
DISPLAY_NAME.....Whole Blood
HELP_TEXT.....No Help Text
ID.....{144EE01A-16B5-4F9B-9A1C-D420D0A9A2D0}
TERM.....WHOLE BLOOD

PRODUCT: WOUND DRAINAGE

ABBREVIATION.....WOUN DRAIN
ACTIVE.....Yes
DESCRIPTION.....Wound Drainage
DISPLAY_NAME.....Wound Drainage
HELP_TEXT.....Drainage From Wound
ID.....{73038CBA-1F31-47E2-A78E-50ED1067CBEA}
TERM.....WOUND DRAINAGE

Page: 1004

PRODUCT: WOUND IRRIGATION

ABBREVIATION.....WOUND IRRI
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Wound Irrigation
HELP_TEXT.....No Help Text
ID.....{192EC6B7-8E2E-418C-9345-BC8A7986CA59}
TERM.....WOUND IRRIGATION

PRODUCT: WOUND VACUUM

ABBREVIATION.....WOUND VAC
ACTIVE.....Yes
DESCRIPTION.....Wound Vacuum Output Product
DISPLAY_NAME.....Wound Vacuum
HELP_TEXT.....No Help Text
ID.....{9F4983A6-B622-4E4D-A7F7-ADA7A4B7CEBC}
TERM.....WOUND VACUUM

=====

TYPE: QUALITY

Page: 1005

QUALITY: ACTUAL

ABBREVIATION.....A
ACTIVE.....Yes
DESCRIPTION.....No Description

DISPLAY_NAME.....Actual
HELP_TEXT.....No Help Text
ID.....{5842A4D9-56A1-7A5A-13C9-F481699C3A09}
TERM.....ACTUAL

QUALITY: ANTECUBITAL

ABBREVIATION.....ANTECUB
ACTIVE.....No
DESCRIPTION.....No Description
DISPLAY_NAME.....Antecubital
HELP_TEXT.....No Help Text
ID.....{0D56DDB8-EDCD-BF9F-5F77-567BC59F3580}
TERM.....ANTECUBITAL

QUALITY: BRACHIAL

ABBREVIATION.....BRACHIAL
ACTIVE.....No
DESCRIPTION.....No Description
DISPLAY_NAME.....Brachial
HELP_TEXT.....No Help Text
ID.....{CFE7E2E3-20BC-7ACE-262B-8BAE4FB36B59}
TERM.....BRACHIAL

QUALITY: DRY

ABBREVIATION.....D
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Dry
HELP_TEXT.....No Help Text
ID.....{B319EFF4-BA4A-85B1-2617-E00E46A54ABE}
TERM.....DRY

QUALITY: ESTIMATED

ABBREVIATION.....E
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Estimated
HELP_TEXT.....No Help Text
ID.....{9EAFF26B-7D87-DA66-F991-022D095FB5A0}
TERM.....ESTIMATED

Page: 1006

QUALITY: ESTIMATED BY ARM SPAN

ABBREVIATION.....EAS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Estimated by Arm Span
HELP_TEXT.....No Help Text
ID.....{486CAE5A-BE56-06BA-C65E-1F8131DB1C2D}

TERM.....ESTIMATED BY ARM SPAN

QUALITY: OTHER

ABBREVIATION.....Oth
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Other
HELP_TEXT.....No Help Text
ID.....{7C0850DD-931C-65C8-1003-E31FBBADF6E9}
TERM.....OTHER

QUALITY: PEDIATRIC

ABBREVIATION.....Ped
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Pediatric
HELP_TEXT.....No Help Text
ID.....{83C8FF79-1919-B001-9F8B-68F014216F0B}
TERM.....PEDIATRIC

QUALITY: RADIAL

ABBREVIATION.....RADIAL
ACTIVE.....No
DESCRIPTION.....No Description
DISPLAY_NAME.....Radial
HELP_TEXT.....No Help Text
ID.....{8F407DAB-8C6D-A5B3-524E-444B22613F4D}
TERM.....RADIAL

QUALITY: STATED

ABBREVIATION.....Sta
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Stated
HELP_TEXT.....No Help Text
ID.....{6362C194-B7F2-D38C-5B72-3BE5DE09438E}
TERM.....STATED

=====

SCHEDULE: CATHETER

ABBREVIATION.....CATHETER
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Catheter

HELP_TEXT.....No Help Text
ID.....{8C788C2B-3C73-7D8D-9032-A1E2188413B8}
TERM.....CATHETER

SCHEDULE: PRN

ABBREVIATION.....PRN
ACTIVE.....Yes
DESCRIPTION.....AS NEEDED
DISPLAY_NAME.....PRN
HELP_TEXT.....No Help Text
ID.....{B9423F6A-166B-CF51-6375-E3928FE0BAD6}
TERM.....PRN

SCHEDULE: Q2H

ABBREVIATION.....q2h
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Every 2 hours
HELP_TEXT.....No Help Text
ID.....{8EF08D16-5529-5FF0-5B84-36703D82E01C}
TERM.....Q2H

SCHEDULE: Q4H

ABBREVIATION.....q4
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Every 4 Hours
HELP_TEXT.....No Help Text
ID.....{E5638DC5-2772-FD4A-2D81-04B2DA5C8DDC}
TERM.....Q4H

SCHEDULE: Q7D

ABBREVIATION.....q7d
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Every 7 Days
HELP_TEXT.....No Help Text
ID.....{E3E171BF-17B3-1C59-112B-042380890C62}
TERM.....Q7D

SCHEDULE: QD

ABBREVIATION.....qd
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Everyday
HELP_TEXT.....No Help Text
ID.....{741F0241-5BC5-5AE2-8D8E-543155C46711}
TERM.....QD

SCHEDULE: RESTRAINT STRAP

ABBREVIATION.....RES STRAP
ACTIVE.....No
DESCRIPTION.....Restraint Strap
DISPLAY_NAME.....Restaint Strap
HELP_TEXT.....No Help Text
ID.....{D9162D47-0F50-49F4-837F-70AD65C169BB}
TERM.....RESTRAINT STRAP

=====

TYPE: TASK

Page: 1010

=====

TASK: ACTIVITY - ABD PILLOW

ABBREVIATION.....ACT ABD PI
ACTIVE.....Yes
DESCRIPTION.....ACTIVITY - ABD PILLOW
DISPLAY_NAME.....ACTIVITY - ABD PILLOW
HELP_TEXT.....ACTIVITY - ABD PILLOW
ID.....{381AA3BB-8858-4B2B-B607-8FC6C6FF707F}
TERM.....ACTIVITY - ABD PILLOW

TASK: ACTIVITY - AMBULATE

ABBREVIATION.....ACT AMBULA
ACTIVE.....Yes
DESCRIPTION.....ACTIVITY - AMBULATE
DISPLAY_NAME.....ACTIVITY - AMBULATE
HELP_TEXT.....ACTIVITY - AMBULATE
ID.....{E7232201-7FC2-4900-82D0-EF3A10FCB760}
TERM.....ACTIVITY - AMBULATE

TASK: ACTIVITY - CONT. PASSIVE MOTION DEV.

ABBREVIATION.....ACT CPM DE
ACTIVE.....Yes
DESCRIPTION.....ACTIVITY - CONTINUOUS PASSIVE MOTION
DEVICE
DISPLAY_NAME.....ACTIVITY - CONT. PASS. MOT. DEV.
HELP_TEXT.....ACTIVITY - CONTINUOUS PASSIVE MOTION
DEVICE
ID.....{7587CDE1-CD6F-4A9C-8E4E-9919485E321E}
TERM.....ACTIVITY - CONT. PASSIVE MOTION DEV.

TASK: ACTIVITY - ELEVATE HEAD OF BED

ABBREVIATION.....ACT EL HOB
ACTIVE.....Yes
DESCRIPTION.....ACTIVITY - ELEVATE HEAD OF BED
DISPLAY_NAME.....ACTIVITY - ELEVATE HEAD OF BED
HELP_TEXT.....ACTIVITY - ELEVATE HEAD OF BED
ID.....{66CDCA34-6B02-4D8B-AD39-FD56BD51EC7E}

TERM.....ACTIVITY - ELEVATE HEAD OF BED

TASK: ACTIVITY - KNEE IMMOBILIZER

ABBREVIATION.....ACT KNEE I
ACTIVE.....Yes
DESCRIPTION.....ACTIVITY - KNEE IMMOBILIZER
DISPLAY_NAME.....ACTIVITY - KNEE IMMOBILIZER
HELP_TEXT.....ACTIVITY - KNEE IMMOBILIZER
ID.....{182692D9-7ECC-45B9-82F6-163DFF937438}
TERM.....ACTIVITY - KNEE IMMOBILIZER

Page: 1011

TASK: ACTIVITY - ORIENT TO WARD

ABBREVIATION.....ACT ORIENT
ACTIVE.....Yes
DESCRIPTION.....Orient patient and family to ward
policies and procedures
DISPLAY_NAME.....ACTIVITY - ORIENT TO WARD
HELP_TEXT.....ACTIVITY - ORIENT TO WARD
ID.....{35B092D6-86E0-655B-9B0E-8E819FE6967D}
TERM.....ACTIVITY - ORIENT TO WARD

TASK: ACTIVITY - RANGE OF MOTION EXERCISES

ABBREVIATION.....ACT ROM EX
ACTIVE.....Yes
DESCRIPTION.....ACTIVITY - RANGE OF MOTION EXERCISES
DISPLAY_NAME.....ACTIVITY - RANGE OF MOTION EXERCISES
HELP_TEXT.....ACTIVITY - RANGE OF MOTION EXERCISES
ID.....{D028A8B1-2548-485F-9D87-D9154858E076}
TERM.....ACTIVITY - RANGE OF MOTION EXERCISES

TASK: ACTIVITY - SIT IN CHAIR

ABBREVIATION.....ACT SIT CH
ACTIVE.....Yes
DESCRIPTION.....ACTIVITY - SIT IN CHAIR
DISPLAY_NAME.....ACTIVITY - SIT IN CHAIR
HELP_TEXT.....ACTIVITY - SIT IN CHAIR
ID.....{62765F1D-A170-450C-89C3-5F9EEFAFF71B}
TERM.....ACTIVITY - SIT IN CHAIR

TASK: ACTIVITY - TURN PATIENT

ABBREVIATION.....ACT TURN
ACTIVE.....Yes
DESCRIPTION.....ACTIVITY - TURN PATIENT
DISPLAY_NAME.....ACTIVITY - TURN PATIENT
HELP_TEXT.....ACTIVITY - TURN PATIENT
ID.....{7A545892-7660-4ADE-8017-12DD76762EE8}
TERM.....ACTIVITY - TURN PATIENT

TASK: ADL - BATH ASSIST

ABBREVIATION.....ADLBATHASS
ACTIVE.....Yes
DESCRIPTION.....ADL - BATH ASSIST
DISPLAY_NAME.....ADL - BATH ASSIST
HELP_TEXT.....ADL - BATH ASSIST
ID.....{0DF710C3-7EB8-44EA-8425-C4F5707D1744}
TERM.....ADL - BATH ASSIST

Page: 1012

TASK: ADL - BATH COMPLETE

ABBREVIATION.....ADLBATHCOM
ACTIVE.....Yes
DESCRIPTION.....ADL - BATH COMPLETE
DISPLAY_NAME.....ADL - BATH COMPLETE
HELP_TEXT.....No Help Text
ID.....{E48D2F58-C42D-4D61-9E34-19DB98966970}
TERM.....ADL - BATH COMPLETE

TASK: ADL - BATH PARTIAL

ABBREVIATION.....ADLBATHPAR
ACTIVE.....Yes
DESCRIPTION.....ADL - BATH PARTIAL
DISPLAY_NAME.....ADL - BATH PARTIAL
HELP_TEXT.....ADL - BATH PARTIAL
ID.....{6997EBBC-1D7C-46CC-AFC0-0035A777DFBD}
TERM.....ADL - BATH PARTIAL

TASK: ADL - BATH SELF

ABBREVIATION.....ADLBATHSEL
ACTIVE.....Yes
DESCRIPTION.....ADL - BATH SELF
DISPLAY_NAME.....ADL - BATH SELF
HELP_TEXT.....ADL - BATH SELF
ID.....{F79F34A1-F147-4A68-8039-2F80E3BF7CBE}
TERM.....ADL - BATH SELF

TASK: ADL - CATHETER CARE

ABBREVIATION.....ADLCATHCAR
ACTIVE.....Yes
DESCRIPTION.....ADL - CATHETER CARE
DISPLAY_NAME.....ADL - CATHETER CARE
HELP_TEXT.....ADL - CATHETER CARE
ID.....{8AAD5864-D7A5-4886-AD58-BDDDC00D80DC}
TERM.....ADL - CATHETER CARE

TASK: ADL - FEED PATIENT

ABBREVIATION.....ADL FEED P

ACTIVE.....Yes
DESCRIPTION.....ADL - FEED PATIENT
DISPLAY_NAME.....ADL - FEED PATIENT
HELP_TEXT.....ADL - FEED PATIENT
ID.....{D525A017-C515-4FBE-9CC2-FDB5FA6CF010}
TERM.....ADL - FEED PATIENT

Page: 1013

TASK: ADL - ORAL HYGIENE

ABBREVIATION.....ADLORALHYG
ACTIVE.....Yes
DESCRIPTION.....ADL - ORAL HYGIENE
DISPLAY_NAME.....ADL - ORAL HYGIENE
HELP_TEXT.....ADL - ORAL HYGIENE
ID.....{27044060-F4C7-4616-9CE4-D2BFE1001EEE}
TERM.....ADL - ORAL HYGIENE

TASK: ADL - SHAVE

ABBREVIATION.....ADLSHAVE
ACTIVE.....Yes
DESCRIPTION.....ADL - SHAVE
DISPLAY_NAME.....ADL - SHAVE
HELP_TEXT.....ADL - SHAVE
ID.....{B5763C28-DEF7-4A6B-86C2-1C1237B6F568}
TERM.....ADL - SHAVE

TASK: ADL - SKIN CARE

ABBREVIATION.....ADLSKINCAR
ACTIVE.....Yes
DESCRIPTION.....ADL - SKIN CARE
DISPLAY_NAME.....ADL - SKIN CARE
HELP_TEXT.....ADL - SKIN CARE
ID.....{7B41CD6E-8F0C-44B5-9415-288ABAF8271A}
TERM.....ADL - SKIN CARE

TASK: ADL - TURN

ABBREVIATION.....ADLTURN
ACTIVE.....Yes
DESCRIPTION.....ADL - TURN
DISPLAY_NAME.....ADL - TURN
HELP_TEXT.....ADL - TURN
ID.....{CCF74486-3ADE-4A67-A6A0-121F195E6FCB}
TERM.....ADL - TURN

TASK: ASSESSMENT - CO

ABBREVIATION.....A CO
ACTIVE.....Yes
DESCRIPTION.....ASSESSMENT - CO
DISPLAY_NAME.....ASSESSMENT - CO

HELP_TEXT.....ASSESSMENT - CO
ID.....{D94C0EFD-CFB8-4249-89D0-1A4BF6C777FC}
TERM.....ASSESSMENT - CO

Page: 1014

TASK: ASSESSMENT - COLLECT VITALS

ABBREVIATION.....A CO VITAL
ACTIVE.....Yes
DESCRIPTION.....ASSESSMENT - COLLECT VITALS
DISPLAY_NAME.....ASSESSMENT - COLLECT VITALS
HELP_TEXT.....ASSESSMENT - COLLECT VITALS
ID.....{9A21F145-1DBA-0BA5-4C9D-E9E31B3F3E45}
TERM.....ASSESSMENT - COLLECT VITALS

TASK: ASSESSMENT - DOPPLERS

ABBREVIATION.....A DOPPLERS
ACTIVE.....Yes
DESCRIPTION.....ASSESSMENT - DOPPLERS
DISPLAY_NAME.....ASSESSMENT - DOPPLERS
HELP_TEXT.....ASSESSMENT - DOPPLERS
ID.....{8B1F6E7E-95D5-47A0-A23D-2A8A088930B8}
TERM.....ASSESSMENT - DOPPLERS

TASK: ASSESSMENT - FINGERSTICK BLOOD SUGAR

ABBREVIATION.....A BLD SUGR
ACTIVE.....Yes
DESCRIPTION.....ASSESSMENT - FINGERSTICK BLOOD SUGAR
DISPLAY_NAME.....ASSESSMENT - FINGERSTICK BLOOD SUGAR
HELP_TEXT.....ASSESSMENT - FINGERSTICK BLOOD SUGAR
(FSBS)
ID.....{50A07B57-D9E6-4F67-B6BC-E21D4A60A56C}
TERM.....ASSESSMENT - FINGERSTICK BLOOD SUGAR

TASK: ASSESSMENT - HEMODYNAMIC READINGS

ABBREVIATION.....A HEMO REA
ACTIVE.....Yes
DESCRIPTION.....ASSESSMENT - HEMODYNAMIC READINGS
DISPLAY_NAME.....ASSESSMENT - HEMODYNAMIC READINGS
HELP_TEXT.....ASSESSMENT - HEMODYNAMIC READINGS
ID.....{98328403-2F7B-7404-34EA-5411633EA55C}
TERM.....ASSESSMENT - HEMODYNAMIC READINGS

TASK: ASSESSMENT - I&O

ABBREVIATION.....A I&O
ACTIVE.....Yes
DESCRIPTION.....ASSESSMENT - I&O
DISPLAY_NAME.....ASSESSMENT - I&O
HELP_TEXT.....ASSESSMENT - I&O
ID.....{53CA663E-40AC-4089-A250-10635CC64F06}

TERM.....ASSESSMENT - I&O

Page: 1015

TASK: ASSESSMENT - NEURO CHECKS

ABBREVIATION.....A NRO CHKS
ACTIVE.....Yes
DESCRIPTION.....ASSESSMENT - NEURO CHECKS
DISPLAY_NAME.....ASSESSMENT - NEURO CHECKS
HELP_TEXT.....ASSESSMENT - NEURO CHECKS
ID.....{8E1E8DE5-EFE6-41B9-86CD-484BC9EF4F66}
TERM.....ASSESSMENT - NEURO CHECKS

TASK: ASSESSMENT - NEUROVASCULAR CHECKS

ABBREVIATION.....A NROVASCU
ACTIVE.....Yes
DESCRIPTION.....ASSESSMENT - NEUROVASCULAR CHECKS
DISPLAY_NAME.....ASSESSMENT - NEUROVASCULAR CHECKS
HELP_TEXT.....ASSESSMENT - NEUROVASCULAR CHECKS
ID.....{4C7FBC3F-F88B-426B-9DCB-74EE2FBC40EA}
TERM.....ASSESSMENT - NEUROVASCULAR CHECKS

TASK: ASSESSMENT - PCA

ABBREVIATION.....A PCA
ACTIVE.....Yes
DESCRIPTION.....ASSESSMENT - PCA
DISPLAY_NAME.....ASSESSMENT - PCA
HELP_TEXT.....ASSESSMENT - PCA
ID.....{88464AA6-3DCF-4388-B0EC-6B3D6B098FEC}
TERM.....ASSESSMENT - PCA

TASK: ASSESSMENT - PULSE OX

ABBREVIATION.....A PULSE OX
ACTIVE.....Yes
DESCRIPTION.....ASSESSMENT - PULSE OX
DISPLAY_NAME.....ASSESSMENT - PULSE OX
HELP_TEXT.....ASSESSMENT - PULSE OX
ID.....{FC2A17F3-DC69-427D-946B-C3403B92528B}
TERM.....ASSESSMENT - PULSE OX

TASK: ASSESSMENT - RESTRAINTS

ABBREVIATION.....A RESTRAIN
ACTIVE.....Yes
DESCRIPTION.....ASSESSMENT - RESTRAINTS
DISPLAY_NAME.....ASSESSMENT - RESTRAINTS
HELP_TEXT.....ASSESSMENT - RESTRAINTS
ID.....{3D02BA0C-BFDD-4D82-8B44-40C861E7FD24}
TERM.....ASSESSMENT - RESTRAINTS

TASK: ASSESSMENT - TELEMETRY

ABBREVIATION.....A TELEMETR
ACTIVE.....Yes
DESCRIPTION.....ASSESSMENT - TELEMETRY
DISPLAY_NAME.....ASSESSMENT - TELEMETRY
HELP_TEXT.....ASSESSMENT - TELEMETRY
ID.....{48C7FB6E-74C2-4EA4-984F-19C11C9072D7}
TERM.....ASSESSMENT - TELEMETRY

TASK: ASSESSMENT - WEIGHT

ABBREVIATION.....A WEIGHT
ACTIVE.....Yes
DESCRIPTION.....ASSESSMENT - GET PATIENTS WEIGHT
DISPLAY_NAME.....ASSESSMENT - WEIGHT
HELP_TEXT.....ASSESSMENT - GET PATIENTS WEIGHT
ID.....{CCB2F8DF-9319-4767-AE7D-90A7705C03FD}
TERM.....ASSESSMENT - WEIGHT

TASK: CALL HOUSE OFFICER - BP DIASTOLIC

ABBREVIATION.....CHO DBP
ACTIVE.....Yes
DESCRIPTION.....CALL HOUSE OFFICER - BP DIASTOLIC
DISPLAY_NAME.....CALL HOUSE OFFICER - BP DIASTOLIC
HELP_TEXT.....CALL HOUSE OFFICER - BP DIASTOLIC
ID.....{19A8A7AD-6BAC-4326-806C-BC80E5636FFC}
TERM.....CALL HOUSE OFFICER - BP DIASTOLIC

TASK: CALL HOUSE OFFICER - BP SYSTOLIC

ABBREVIATION.....CHO SBP
ACTIVE.....Yes
DESCRIPTION.....CALL HOUSE OFFICER - BP SYSTOLIC
DISPLAY_NAME.....CALL HOUSE OFFICER - BP SYSTOLIC
HELP_TEXT.....CALL HOUSE OFFICER - BP SYSTOLIC
ID.....{971F71F9-AAD6-46BF-A729-45F0B98FD6D3}
TERM.....CALL HOUSE OFFICER - BP SYSTOLIC

TASK: CALL HOUSE OFFICER - PULSE

ABBREVIATION.....CHO PULSE
ACTIVE.....Yes
DESCRIPTION.....CALL HOUSE OFFICER - PULSE
DISPLAY_NAME.....CALL HOUSE OFFICER - PULSE
HELP_TEXT.....CALL HOUSE OFFICER - PULSE
ID.....{81EED77A-5116-4B2D-B52B-42BAFACC69C9}
TERM.....CALL HOUSE OFFICER - PULSE

TASK: CALL HOUSE OFFICER - RESPIRATORY RATE

ABBREVIATION.....CHO RR
ACTIVE.....Yes
DESCRIPTION.....CALL HOUSE OFFICER - RESPIRATORY RATE
DISPLAY_NAME.....CALL HOUSE OFFICER - RESPIRATORY RATE
HELP_TEXT.....CALL HOUSE OFFICER - RESPIRATORY RATE
ID.....{D2DFE1EC-017F-49B8-B59C-C724F71F2E63}
TERM.....CALL HOUSE OFFICER - RESPIRATORY RATE

TASK: CALL HOUSE OFFICER - TEMP

ABBREVIATION.....CHO TEMP
ACTIVE.....Yes
DESCRIPTION.....CALL HOUSE OFFICER - TEMP
DISPLAY_NAME.....CALL HOUSE OFFICER - TEMP
HELP_TEXT.....CALL HOUSE OFFICER - TEMP
ID.....{E59A3DC0-73CB-4B0F-B031-E3CFF280EA5F}
TERM.....CALL HOUSE OFFICER - TEMP

TASK: HYDRATION - DIET

ABBREVIATION.....HY DIET
ACTIVE.....Yes
DESCRIPTION.....HYDRATION - DIET
DISPLAY_NAME.....HYDRATION - DIET
HELP_TEXT.....HYDRATION - DIET
ID.....{C0F14019-FF1A-4D42-A857-4B1E86E793FA}
TERM.....HYDRATION - DIET

TASK: HYDRATION - HEPARIN DRIP

ABBREVIATION.....HY HEPARIN
ACTIVE.....Yes
DESCRIPTION.....HYDRATION - HEPARIN DRIP
DISPLAY_NAME.....HYDRATION - HEPARIN DRIP
HELP_TEXT.....HYDRATION - HEPARIN DRIP
ID.....{EC49EE9D-C2EB-4F7D-A9D4-9746DDEF84A6}
TERM.....HYDRATION - HEPARIN DRIP

TASK: HYDRATION - HYPERAL

ABBREVIATION.....HY HYPERAL
ACTIVE.....Yes
DESCRIPTION.....HYDRATION - HYPERAL
DISPLAY_NAME.....HYDRATION - HYPERAL
HELP_TEXT.....HYDRATION - HYPERAL
ID.....{5F91CA9E-82F1-4490-8648-896864F1F68C}
TERM.....HYDRATION - HYPERAL

Page: 1018

TASK: HYDRATION - IV FLUIDS

ABBREVIATION.....HY IV FLU
ACTIVE.....Yes

DESCRIPTION.....HYDRATION - IV FLUIDS
DISPLAY_NAME.....HYDRATION - IV FLUIDS
HELP_TEXT.....HYDRATION - IV FLUIDS
ID.....{75D6B913-147F-4318-B6A5-985438E8CF66}
TERM.....HYDRATION - IV FLUIDS

TASK: HYDRATION - NPO

ABBREVIATION.....HY NPO
ACTIVE.....Yes
DESCRIPTION.....HYDRATION - NPO
DISPLAY_NAME.....HYDRATION - NPO
HELP_TEXT.....HYDRATION - NPO
ID.....{44028E10-FB0D-4AFD-9B8D-90F157D6329D}
TERM.....HYDRATION - NPO

TASK: HYDRATION - PERIPHERAL IV

ABBREVIATION.....HY PIV
ACTIVE.....Yes
DESCRIPTION.....HYDRATION - PERIPHERAL IV
DISPLAY_NAME.....HYDRATION - PERIPHERAL IV
HELP_TEXT.....HYDRATION - PERIPHERAL IV
ID.....{1ACA57CF-4226-4EF1-AABD-142EE1DE369F}
TERM.....HYDRATION - PERIPHERAL IV

TASK: HYDRATION - PICC/CENTRAL LINE

ABBREVIATION.....HY PICC
ACTIVE.....Yes
DESCRIPTION.....HYDRATION - PICC/CENTRAL LINE
DISPLAY_NAME.....HYDRATION - PICC/CENTRAL LINE
HELP_TEXT.....HYDRATION - PICC/CENTRAL LINE
ID.....{E1D6E811-CE80-4F8A-AB00-570D1C8079CA}
TERM.....HYDRATION - PICC/CENTRAL LINE

TASK: HYDRATION - SUPPLEMENT SNACK

ABBREVIATION.....HY SNACK S
ACTIVE.....Yes
DESCRIPTION.....HYDRATION - SUPPLEMENT SNACK
DISPLAY_NAME.....HYDRATION - SUPPLEMENT SNACK
HELP_TEXT.....HYDRATION - SUPPLEMENT SNACK
ID.....{0916F365-7FBC-4FD3-A55B-CAD4E41A00B1}
TERM.....HYDRATION - SUPPLEMENT SNACK

Page: 1019

TASK: HYDRATION - TUBE FEEDING

ABBREVIATION.....HY TU FEED
ACTIVE.....Yes
DESCRIPTION.....HYDRATION - TUBE FEEDING
DISPLAY_NAME.....HYDRATION - TUBE FEEDING
HELP_TEXT.....HYDRATION - TUBE FEEDING

ID.....{6EDD774C-DD5A-4BF6-B06D-725FD7A0B26C}
TERM.....HYDRATION - TUBE FEEDING

TASK: NURSING INTERVENTIONS - STRAIGHT CATH

ABBREVIATION.....NI ST CATH
ACTIVE.....Yes
DESCRIPTION.....NURSING INTERVENTIONS - STRAIGHT CATH
DISPLAY_NAME.....NURSING INTERVENTIONS - STRAIGHT CATH
HELP_TEXT.....NURSING INTERVENTIONS - STRAIGHT CATH
ID.....{A1F54C3F-CFD1-4767-BED8-55BF2BE872DD}
TERM.....NURSING INTERVENTIONS - STRAIGHT CATH

TASK: NURSING INTERVENTIONS - WOUND CARE

ABBREVIATION.....NI WOU CAR
ACTIVE.....Yes
DESCRIPTION.....NURSING INTERVENTIONS - WOUND CARE
DISPLAY_NAME.....NURSING INTERVENTIONS - WOUND CARE
HELP_TEXT.....NURSING INTERVENTIONS - WOUND CARE
ID.....{BFB2EFF-3644-439F-951C-804B43BFFF30}
TERM.....NURSING INTERVENTIONS - WOUND CARE

TASK: PRECAUTIONS - 24 HOUR CHART CHECK

ABBREVIATION.....PREC 24 CH
ACTIVE.....Yes
DESCRIPTION.....PRECAUTIONS - 24 HOUR CHART CHECK
DISPLAY_NAME.....PRECAUTIONS - 24 HOUR CHART CHECK
HELP_TEXT.....PRECAUTIONS - 24 HOUR CHART CHECK
ID.....{96C39F92-CE4C-4D6B-9837-9C9FECDEC199}
TERM.....PRECAUTIONS - 24 HOUR CHART CHECK

TASK: PRECAUTIONS - ASPIRATION

ABBREVIATION.....PREC ASPIR
ACTIVE.....Yes
DESCRIPTION.....PRECAUTIONS - ASPIRATION
DISPLAY_NAME.....PRECAUTIONS - ASPIRATION
HELP_TEXT.....PRECAUTIONS - ASPIRATION
ID.....{D8D28810-655E-4707-BF6D-B853F18A1E08}
TERM.....PRECAUTIONS - ASPIRATION

TASK: PRECAUTIONS - ETOH W/D

ABBREVIATION.....PREC ETOH
ACTIVE.....Yes
DESCRIPTION.....PRECAUTIONS - ETOH W/D
DISPLAY_NAME.....PRECAUTIONS - ETOH W/D
HELP_TEXT.....PRECAUTIONS - ETOH W/D
ID.....{7800B16E-ECCD-4086-985C-0D12D4671123}
TERM.....PRECAUTIONS - ETOH W/D

TASK: PRECAUTIONS - FALL

ABBREVIATION.....PREC FALL
ACTIVE.....Yes
DESCRIPTION.....PRECAUTIONS - FALL
DISPLAY_NAME.....PRECAUTIONS - FALL
HELP_TEXT.....PRECAUTIONS - FALL
ID.....{CAB2E3CC-18AD-40AB-946B-DB0292E7FA09}
TERM.....PRECAUTIONS - FALL

TASK: PRECAUTIONS - PATIENT TEACHING

ABBREVIATION.....PREC TEACH
ACTIVE.....Yes
DESCRIPTION.....PRECAUTIONS - PATIENT TEACHING
DISPLAY_NAME.....PRECAUTIONS - PATIENT TEACHING
HELP_TEXT.....PRECAUTIONS - PATIENT TEACHING
ID.....{C58B33AB-6442-46A5-B24A-E2B856C72B98}
TERM.....PRECAUTIONS - PATIENT TEACHING

TASK: PRECAUTIONS - VALIDATE ALARMS

ABBREVIATION.....PREC ALARM
ACTIVE.....Yes
DESCRIPTION.....PRECAUTIONS - VALIDATE ALARMS
DISPLAY_NAME.....PRECAUTIONS - VALIDATE ALARMS
HELP_TEXT.....PRECAUTIONS - VALIDATE ALARMS
ID.....{CF386EC1-F5B6-4D15-8E48-8453FE6FC1CC}
TERM.....PRECAUTIONS - VALIDATE ALARMS

TASK: PULMONARY - CHEST PHYSIO-THERAPY

ABBREVIATION.....PULMCHPPTH
ACTIVE.....Yes
DESCRIPTION.....PULMONARY - CHEST PHYSIO-THERAPY
DISPLAY_NAME.....PULMONARY - CHEST PHYSIO-THERAPY
HELP_TEXT.....PULMONARY - CHEST PHYSIO-THERAPY
ID.....{F028FA64-3957-4A15-B30A-277AF594A69D}
TERM.....PULMONARY - CHEST PHYSIO-THERAPY

Page: 1021

TASK: PULMONARY - COUGH AND DEEP BREATH

ABBREVIATION.....PULMCOUGH
ACTIVE.....Yes
DESCRIPTION.....PULMONARY - COUGH AND DEEP BREATH
DISPLAY_NAME.....PULMONARY - COUGH AND DEEP BREATH
HELP_TEXT.....PULMONARY - COUGH AND DEEP BREATH
ID.....{17A5D60E-ED50-4FBA-879D-4618563D2B98}
TERM.....PULMONARY - COUGH AND DEEP BREATH

TASK: PULMONARY - ET TUBE CARE

ABBREVIATION.....PULM ET TU

ACTIVE.....Yes
DESCRIPTION.....PULMONARY - ET TUBE CARE
DISPLAY_NAME.....PULMONARY - ET TUBE CARE
HELP_TEXT.....PULMONARY - ET TUBE CARE
ID.....{B2F57E59-33C5-49FB-BEE9-176647329910}
TERM.....PULMONARY - ET TUBE CARE

TASK: PULMONARY - INCENTIVE SPIROMETER

ABBREVIATION.....PULMINCSPI
ACTIVE.....Yes
DESCRIPTION.....PULMONARY - INCENTIVE SPIROMETER
DISPLAY_NAME.....PULMONARY - INCENTIVE SPIROMETER
HELP_TEXT.....PULMONARY - INCENTIVE SPIROMETER
ID.....{F856E287-BAAC-4D8C-AFAD-34912621F3B4}
TERM.....PULMONARY - INCENTIVE SPIROMETER

TASK: PULMONARY - OXYGEN

ABBREVIATION.....PULM OXYGE
ACTIVE.....Yes
DESCRIPTION.....PULMONARY - OXYGEN
DISPLAY_NAME.....PULMONARY - OXYGEN
HELP_TEXT.....PULMONARY - OXYGEN
ID.....{14A986A8-54F1-4D9A-A2EB-43F306704C42}
TERM.....PULMONARY - OXYGEN

TASK: PULMONARY - SUCTION PATIENT

ABBREVIATION.....PULM-SUCTI
ACTIVE.....Yes
DESCRIPTION.....PULMONARY - SUCTION PATIENT
DISPLAY_NAME.....PULMONARY - SUCTION PATIENT
HELP_TEXT.....PULMONARY - SUCTION PATIENT
ID.....{8B74C2F6-8CFB-4A60-9158-5333B980D047}
TERM.....PULMONARY - SUCTION PATIENT

Page: 1022

TASK: PULMONARY - TRACH CARE

ABBREVIATION.....PULMTRACHC
ACTIVE.....Yes
DESCRIPTION.....PULMONARY - TRACH CARE
DISPLAY_NAME.....PULMONARY - TRACH CARE
HELP_TEXT.....PULMONARY - TRACH CARE
ID.....{EAAFB44D-359B-4FB0-9683-BC44F7D59C23}
TERM.....PULMONARY - TRACH CARE

TASK: RESTRAINT - ASSESS NEED FOR FLUIDS/TOILET

ABBREVIATION.....RES TOILET
ACTIVE.....Yes
DESCRIPTION.....Restraint - Assess Need for
Fluids/Toileting

DISPLAY_NAME.....Restraint - Assess Need for
Fluids/Toileting
HELP_TEXT.....Restraint - Assess Need for
Fluids/Toileting
ID.....{58A09EA6-3336-4985-83A9-B765E3CB9C9E}
TERM.....RESTRAINT - ASSESS NEED FOR
FLUIDS/TOILET

TASK: RESTRAINT - CIRCULATION CHECK

ABBREVIATION.....RES CIR CH
ACTIVE.....Yes
DESCRIPTION.....Restraint - Circulation Check
DISPLAY_NAME.....Restraint - Circulation Check
HELP_TEXT.....Restraint - Circulation Check
ID.....{0A1DA91C-BC0E-4B07-967D-2BAB7E64B028}
TERM.....RESTRAINT - CIRCULATION CHECK

TASK: RESTRAINT - DEVICE RELEASE

ABBREVIATION.....RES DEV RE
ACTIVE.....Yes
DESCRIPTION.....Restraint - Device Release
DISPLAY_NAME.....Restraint - Device Release
HELP_TEXT.....Restraint - Device Release
ID.....{7B0C29DA-6593-472C-9716-24F7E2B7ABD7}
TERM.....RESTRAINT - DEVICE RELEASE

TASK: RESTRAINT - FALL PRECAUTIONS

ABBREVIATION.....RES FALL P
ACTIVE.....Yes
DESCRIPTION.....Restraint - Fall Precautions
DISPLAY_NAME.....Restraint - Fall Precautions
HELP_TEXT.....Restraint - Fall Precautions
ID.....{D812036A-7B6B-4C9C-9EE2-616C5EFC4789}
TERM.....RESTRAINT - FALL PRECAUTIONS

Page: 1023

TASK: RESTRAINT - SEIZURE PRECAUTIONS

ABBREVIATION.....RES SEI PR
ACTIVE.....Yes
DESCRIPTION.....Restraint - Seizure Precautions
DISPLAY_NAME.....Restraint - Seizure Precautions
HELP_TEXT.....Restraint - Seizure Precautions
ID.....{71606A0A-05E6-4DFD-9E00-677A07167199}
TERM.....RESTRAINT - SEIZURE PRECAUTIONS

TASK: TESTS/PROC. - ACCUCHECKS

ABBREVIATION.....TP ACCUCHE
ACTIVE.....Yes
DESCRIPTION.....TESTS/PROC. - ACCUCHECKS

DISPLAY_NAME.....TESTS/PROC. - ACCUCHECKS
HELP_TEXT.....TESTS/PROC. - ACCUCHECKS
ID.....{0A5E0495-B990-449B-B5F6-10497842534E}
TERM.....TESTS/PROC. - ACCUCHECKS

TASK: TESTS/PROC. - BLOOD CULTURES

ABBREVIATION.....T/P BLOOD
ACTIVE.....Yes
DESCRIPTION.....TESTS/PROC. - BLOOD CULTURES
DISPLAY_NAME.....TESTS/PROC. - BLOOD CULTURES
HELP_TEXT.....TESTS/PROC. - BLOOD CULTURES
ID.....{BC2F7EFE-3990-47D2-863A-6FB2227746F3}
TERM.....TESTS/PROC. - BLOOD CULTURES

TASK: TESTS/PROC. - DAILY LABS

ABBREVIATION.....T/P D LABS
ACTIVE.....Yes
DESCRIPTION.....TESTS/PROC. - DAILY LABS
DISPLAY_NAME.....TESTS/PROC. - DAILY LABS
HELP_TEXT.....TESTS/PROC. - DAILY LABS
ID.....{03B3A828-D4B0-40CF-B12C-1DD936DE98AC}
TERM.....TESTS/PROC. - DAILY LABS

TASK: TESTS/PROC. - MORNING LABS

ABBREVIATION.....T/P M LABS
ACTIVE.....Yes
DESCRIPTION.....TESTS/PROC. - MORNING LABS
DISPLAY_NAME.....TESTS/PROC. - MORNING LABS
HELP_TEXT.....TESTS/PROC. - MORNING LABS
ID.....{6F81B3D1-CCBC-4BC9-AFF0-64000DF212EB}
TERM.....TESTS/PROC. - MORNING LABS

Page: 1024

TASK: TESTS/PROC. - ZERO BALANCE HEMO. LINES

ABBREVIATION.....T/P 0 HEMO
ACTIVE.....Yes
DESCRIPTION.....TESTS/PROC. - ZERO BALANCE HEMO. LINES
DISPLAY_NAME.....TESTS/PROC. - ZERO BALANCE HEMO. LINES
HELP_TEXT.....TESTS/PROC. - ZERO BALANCE HEMO. LINES
ID.....{9CEDB731-A785-48B7-AC92-84B10FCCCE86}
TERM.....TESTS/PROC. - ZERO BALANCE HEMO. LINES

TASK: TUBES & DRAINS - CONDOM CATH

ABBREVIATION.....T&D C CATH
ACTIVE.....Yes
DESCRIPTION.....TUBES & DRAINS - CONDOM CATH
DISPLAY_NAME.....TUBES & DRAINS - CONDOM CATH
HELP_TEXT.....TUBES & DRAINS - CONDOM CATH
ID.....{385EFB2B-55BB-486E-B168-D1FFF71E086C}

TERM.....TUBES & DRAINS - CONDOM CATH

TASK: TUBES & DRAINS - FOLEY

ABBREVIATION.....T&D FOLEY
ACTIVE.....Yes
DESCRIPTION.....TUBES & DRAINS - FOLEY
DISPLAY_NAME.....TUBES & DRAINS - FOLEY
HELP_TEXT.....TUBES & DRAINS - FOLEY
ID.....{14CECC17-A780-41DD-98A4-6E226630BC65}
TERM.....TUBES & DRAINS - FOLEY

TASK: TUBES & DRAINS - I/O CATH

ABBREVIATION.....T&D I/O CA
ACTIVE.....Yes
DESCRIPTION.....TUBES & DRAINS - I/O CATH
DISPLAY_NAME.....TUBES & DRAINS - I/O CATH
HELP_TEXT.....TUBES & DRAINS - I/O CATH
ID.....{5BA34E92-5248-4954-A677-443C4B2183E8}
TERM.....TUBES & DRAINS - I/O CATH

TASK: TUBES & DRAINS - NG TUBE

ABBREVIATION.....T&D NG TU
ACTIVE.....Yes
DESCRIPTION.....TUBES & DRAINS - NG TUBE
DISPLAY_NAME.....TUBES & DRAINS - NG TUBE
HELP_TEXT.....TUBES & DRAINS - NG TUBE
ID.....{BBFCE51A-9CCF-44AC-8C37-4A779232A4B7}
TERM.....TUBES & DRAINS - NG TUBE

Page: 1025

TASK: TUBES & DRAINS - SUPRAPUBIC DRAIN

ABBREVIATION.....T&D SUP DR
ACTIVE.....Yes
DESCRIPTION.....TUBES & DRAINS - SUPRAPUBIC DRAIN
DISPLAY_NAME.....TUBES & DRAINS - SUPRAPUBIC DRAIN
HELP_TEXT.....TUBES & DRAINS - SUPRAPUBIC DRAIN
ID.....{6A4A5F28-35D3-494E-BA97-2E0DE4A7B8CC}
TERM.....TUBES & DRAINS - SUPRAPUBIC DRAIN

=====

TYPE: UNIT

Page: 1026

=====

UNIT: BEATS PER MINUTE

ABBREVIATION.....bpm
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Beats Per Minute

HELP_TEXT.....No Help Text
ID.....{8319B3B0-0F31-1393-154D-59A86692D637}
TERM.....BEATS PER MINUTE

UNIT: BRADEN SCALE

ABBREVIATION.....BRADEN SCA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Braden Scale
HELP_TEXT.....No Help Text
ID.....{B5EE25BB-DA66-4AC9-9A3B-C05497B2647C}
TERM.....BRADEN SCALE

UNIT: BREATHS PER MINUTE

ABBREVIATION.....rpm
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Resp/Min
HELP_TEXT.....No Help Text
ID.....{DCF89527-FFAE-F68E-0871-602F3C32DDB1}
TERM.....BREATHS PER MINUTE

UNIT: CALORIE

ABBREVIATION.....CALORIE
ACTIVE.....Yes
DESCRIPTION.....A unit of food energy
DISPLAY_NAME.....Calorie
HELP_TEXT.....The amount of energy required to
raise the temperature of a liter of
water one degree centigrade at sea
level
ID.....{9F7DE371-AB4B-45EE-A8E3-825822419CCA}
TERM.....CALORIE

UNIT: CENTIMETER

ABBREVIATION.....cm
ACTIVE.....Yes
DESCRIPTION.....Measurment of length
DISPLAY_NAME.....Centimeter
HELP_TEXT.....Metric
ID.....{EF700458-8C2D-16E8-02DC-9E1711584C53}
TERM.....CENTIMETER

UNIT: CENTIMETERS H2O

ABBREVIATION.....CM H2O
ACTIVE.....Yes
DESCRIPTION.....CENTIMETERS OF WATER
DISPLAY_NAME.....CM H2O

HELP_TEXT.....No Help Text
ID.....{F7030BB1-AD2C-15F6-FF41-6DA893C851C1}
TERM.....CENTIMETERS H2O

UNIT: CIWA SCALE

ABBREVIATION.....CIWA SCALE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....CIWA SCALE
HELP_TEXT.....Scale For Addition Of All CIWA Scores
ID.....{0ED378B9-C1A4-4ABA-82B7-F62BC0D62886}
TERM.....CIWA SCALE

UNIT: CUBIC CENTIMETER

ABBREVIATION.....Cubic Cent
ACTIVE.....Yes
DESCRIPTION.....Cubic Centimeter
DISPLAY_NAME.....Cubic Centimeter
HELP_TEXT.....No Help Text
ID.....{4698E6A0-DA6C-4A0A-800A-395BBCD2871E}
TERM.....CUBIC CENTIMETER

UNIT: DAY

ABBREVIATION.....Day
ACTIVE.....Yes
DESCRIPTION.....Measurement of passage of time (24
hours)
DISPLAY_NAME.....Day
HELP_TEXT.....No Help Text
ID.....{742619B2-A8C7-1A53-45CD-C3E25A05D58C}
TERM.....DAY

UNIT: DEGREES (ANGLE)

ABBREVIATION.....DEG
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Degrees (Angle)
HELP_TEXT.....No Help Text
ID.....{5F407D72-3009-2553-7DDC-82BAF2043FCA}
TERM.....DEGREES (ANGLE)

Page: 1028

UNIT: DEGREES C

ABBREVIATION.....c
ACTIVE.....Yes
DESCRIPTION.....Metric Temperature
DISPLAY_NAME.....Celsius
HELP_TEXT.....No Help Text
ID.....{1964CB4C-0584-DA80-C7A9-E7400DCD817D}

TERM.....DEGREES C

UNIT: DEGREES F

ABBREVIATION.....F
ACTIVE.....Yes
DESCRIPTION.....US measurement of temperature
DISPLAY_NAME.....Fahrenheit
HELP_TEXT.....No Help Text
ID.....{59FC6CA4-9389-B770-5A9C-8C5A38C572D8}
TERM.....DEGREES F

UNIT: DYNE * SECONDS / CM(-5) * M2

ABBREVIATION.....D*S/CM*M2
ACTIVE.....Yes
DESCRIPTION.....dyne * seconds / cm(-5) * m2
DISPLAY_NAME.....dyne * seconds / cm(-5) * m2
HELP_TEXT.....dyne * seconds / cm(-5) * m2
ID.....{1241878F-AECA-B1A8-BF12-D5FA583FC1B2}
TERM.....DYNE * SECONDS / CM(-5) * M2

UNIT: DYNE * SECONDS / CM(-5) / M2

ABBREVIATION.....D*S/CM/M2
ACTIVE.....Yes
DESCRIPTION.....dyne * seconds/cm(-5) / m2
DISPLAY_NAME.....dyne * seconds/cm(-5) / m2
HELP_TEXT.....dyne * seconds/cm(-5) / m2
ID.....{A642531F-A72E-CF60-D31A-F38AD12A8B30}
TERM.....DYNE * SECONDS / CM(-5) / M2

UNIT: DYNE-SEC/CM5

ABBREVIATION.....dyne-sec/c
ACTIVE.....Yes
DESCRIPTION.....dyne-sec/cm5
DISPLAY_NAME.....dyne-sec/cm5
HELP_TEXT.....dyne-sec/cm5
ID.....{01A538D8-41AF-4D0B-B827-D557060A618A}
TERM.....DYNE-SEC/CM5

Page: 1029

UNIT: EVENT(S)

ABBREVIATION.....EVENT(S)
ACTIVE.....Yes
DESCRIPTION.....Event(s)
DISPLAY_NAME.....Event(s)
HELP_TEXT.....No Help Text
ID.....{723F907D-8B71-4183-AE54-66D0C3C72DCA}
TERM.....EVENT(S)

UNIT: EVENTS PER HOUR

ABBREVIATION.....EVENT/H
ACTIVE.....Yes
DESCRIPTION.....NUMBER OF TIMES PER HOUR THAT AN
OBSERVABLE EVENT OCCURS
DISPLAY_NAME.....Events Per Hour
HELP_TEXT.....No Help Text
ID.....{3D648D4F-02A9-474A-65BF-70608E75BC29}
TERM.....EVENTS PER HOUR

UNIT: FOOT

ABBREVIATION.....FT
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Foot/Feet
HELP_TEXT.....No Help Text
ID.....{66D6F7EC-3884-4790-AC35-10B2AE84257F}
TERM.....FOOT

UNIT: FRENCH

ABBREVIATION.....FRENCH
ACTIVE.....Yes
DESCRIPTION.....UNIT OF MEASURE FOR ARTIFICIAL AIRWAYS
DISPLAY_NAME.....French
HELP_TEXT.....No Help Text
ID.....{3734E0E0-087B-098A-984A-E42264CF6567}
TERM.....FRENCH

UNIT: GAUGE

ABBREVIATION.....GAUGE
ACTIVE.....Yes
DESCRIPTION.....Gauge For Needles
DISPLAY_NAME.....Gauge
HELP_TEXT.....No Help Text
ID.....{5670BBC1-4962-4FA1-A4D0-1E399F2D06B1}
TERM.....GAUGE

Page: 1030

UNIT: GLASGOW COMA SCALE

ABBREVIATION.....GCS
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Glasgow Coma Scale
HELP_TEXT.....No Help Text
ID.....{CA674718-C98D-669F-8954-F9EC91C2F6F3}
TERM.....GLASGOW COMA SCALE

UNIT: GM-M/M2/BEAT

ABBREVIATION.....gm-m/m2/be

ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....gm-m/m2/beat
HELP_TEXT.....No Help Text
ID.....{DB5404C8-1880-4B54-B053-DF0875E72C3B}
TERM.....GM-M/M2/BEAT

UNIT: GRAMS PER HOUR

ABBREVIATION.....gm/hr
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....gm/hr
HELP_TEXT.....No Help Text
ID.....{A50449B1-DF42-4D18-8F02-59937EA6CDE4}
TERM.....GRAMS PER HOUR

UNIT: GRAMS PER MINUTE

ABBREVIATION.....gm/min
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Grams Per Minute
HELP_TEXT.....No Help Text
ID.....{07D005E8-B078-A48F-96B5-8CCAA92EC97F}
TERM.....GRAMS PER MINUTE

UNIT: HOUR

ABBREVIATION.....hr
ACTIVE.....Yes
DESCRIPTION.....Measurement of time
DISPLAY_NAME.....Hour
HELP_TEXT.....No Help Text
ID.....{60034E2A-D088-E9B6-9B98-D0F4357A059F}
TERM.....HOUR

Page: 1031

UNIT: INCH

ABBREVIATION.....in
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Inch
HELP_TEXT.....No Help Text
ID.....{9CB32FC9-4130-11F8-E750-AF1415556705}
TERM.....INCH

UNIT: KG/M2

ABBREVIATION.....KG/M2
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....KG/M2

HELP_TEXT.....No Help Text
ID.....{E3204003-83B3-4D12-872D-C32A50F177D0}
TERM.....KG/M2

UNIT: KILOGRAM

ABBREVIATION.....kg
ACTIVE.....Yes
DESCRIPTION.....Metric measurement of weight
DISPLAY_NAME.....Kilogram
HELP_TEXT.....No Help Text
ID.....{0AC4E9AC-41CB-32EA-9674-74FDFAA6B4F4}
TERM.....KILOGRAM

UNIT: LITER

ABBREVIATION.....L
ACTIVE.....Yes
DESCRIPTION.....Volume, 1000 ml
DISPLAY_NAME.....Liter
HELP_TEXT.....No Help Text
ID.....{AFADB383-5867-4D7C-8ADB-E76782E5943D}
TERM.....LITER

UNIT: LITERS PER MINUTE

ABBREVIATION.....L/min
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L/min
HELP_TEXT.....No Help Text
ID.....{CFCD64D2-8555-8513-5DB0-F0522C19DCCF}
TERM.....LITERS PER MINUTE

Page: 1032

UNIT: LITERS PER MINUTE PER METER2 (BSA)

ABBREVIATION.....L/min/m2
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....L/min/m2 (BSA)
HELP_TEXT.....No Help Text
ID.....{F033FDA0-C6D6-0D93-F836-C4A42708FDF2}
TERM.....LITERS PER MINUTE PER METER2 (BSA)

UNIT: MICROGRAMS

ABBREVIATION.....mcg
ACTIVE.....Yes
DESCRIPTION.....Micrograms
DISPLAY_NAME.....mcg
HELP_TEXT.....Micrograms
ID.....{5D2A9B6F-4DD1-48C6-B223-09CAF66A5A50}
TERM.....MICROGRAMS

UNIT: MICROGRAMS PER HOUR

ABBREVIATION.....mcg/hr
ACTIVE.....Yes
DESCRIPTION.....Micrograms per Hour
DISPLAY_NAME.....mcg/hr
HELP_TEXT.....Micrograms per Hour
ID.....{971A5FA5-2FD7-4ABD-9035-229B00DC484E}
TERM.....MICROGRAMS PER HOUR

UNIT: MICROGRAMS PER KILOGRAM PER HOUR

ABBREVIATION.....mcg/kg/hr
ACTIVE.....Yes
DESCRIPTION.....Micrograms per Kilogram per Hour
DISPLAY_NAME.....mcg/kg/hr
HELP_TEXT.....Micrograms per Kilogram per Hour
ID.....{E5C07334-6932-4329-88A9-2F56780C5DA9}
TERM.....MICROGRAMS PER KILOGRAM PER HOUR

UNIT: MICROGRAMS PER KILOGRAM PER MINUTE

ABBREVIATION.....mcg/kg/min
ACTIVE.....Yes
DESCRIPTION.....Micrograms per Kilogram per Minute
DISPLAY_NAME.....mcg/kg/min
HELP_TEXT.....Micrograms per Kilogram per Minute
ID.....{37423DE7-B103-456A-9B36-6C2FBBC55CBE}
TERM.....MICROGRAMS PER KILOGRAM PER MINUTE

Page: 1033

UNIT: MICROGRAMS PER MINUTE

ABBREVIATION.....mcg/min
ACTIVE.....Yes
DESCRIPTION.....Micrograms per Minute
DISPLAY_NAME.....mcg/min
HELP_TEXT.....Micrograms per Minute
ID.....{BC684DBB-7DCC-4209-A773-D1243616BD69}
TERM.....MICROGRAMS PER MINUTE

UNIT: MILLIAMPS

ABBREVIATION.....mA
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....mA
HELP_TEXT.....No Help Text
ID.....{5146A84A-C500-CB21-4718-0E5755691FD5}
TERM.....MILLIAMPS

UNIT: MILLIBAR

ABBREVIATION.....mBar
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Millibar
HELP_TEXT.....No Help Text
ID.....{DB8DBC26-6A46-4F2B-AD14-B6A8212F4FCA}
TERM.....MILLIBAR

UNIT: MILLIGRAMS

ABBREVIATION.....mg
ACTIVE.....Yes
DESCRIPTION.....Milligrams
DISPLAY_NAME.....mg
HELP_TEXT.....Milligrams
ID.....{E8CF6B0F-FA19-4F22-BA9D-6EBAA93D7D0D}
TERM.....MILLIGRAMS

UNIT: MILLIGRAMS PER HOUR

ABBREVIATION.....mg/hr
ACTIVE.....Yes
DESCRIPTION.....Milligrams per Hour
DISPLAY_NAME.....mg/hr
HELP_TEXT.....Milligrams per Hour
ID.....{1DCC6F14-E4FF-4BDC-ADBA-99192FB44BE7}
TERM.....MILLIGRAMS PER HOUR

Page: 1034

UNIT: MILLIGRAMS PER KILOGRAM PER HOUR

ABBREVIATION.....mg/kg/hr
ACTIVE.....Yes
DESCRIPTION.....Milligrams per Kilogram per Hour
DISPLAY_NAME.....mg/kg/hr
HELP_TEXT.....Milligrams per Kilogram per Hour
ID.....{6900F44E-3581-4558-83E2-76FB376D6711}
TERM.....MILLIGRAMS PER KILOGRAM PER HOUR

UNIT: MILLIGRAMS PER KILOGRAM PER MINUTE

ABBREVIATION.....mg/kg/min
ACTIVE.....Yes
DESCRIPTION.....Milligrams per Kilogram per Minute
DISPLAY_NAME.....mg/kg/min
HELP_TEXT.....Milligrams per Kilogram per Minute
ID.....{7C1C10C2-8B88-4F45-B667-0F988090CD47}
TERM.....MILLIGRAMS PER KILOGRAM PER MINUTE

UNIT: MILLIGRAMS PER MINUTE

ABBREVIATION.....mg/min
ACTIVE.....Yes
DESCRIPTION.....Milligrams per Minute

DISPLAY_NAME.....mg/min
HELP_TEXT.....Milligrams per Minute
ID.....{C5E49AC7-5671-4E84-8047-C7B78607CDA2}
TERM.....MILLIGRAMS PER MINUTE

UNIT: MILLILITER

ABBREVIATION.....ml
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....ml
HELP_TEXT.....No Help Text
ID.....{EAF55107-F8D5-4F0A-6211-B4AB57A874B4}
TERM.....MILLILITER

UNIT: MILLILITERS PER BEAT

ABBREVIATION.....ml/beat
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....ml/beat
HELP_TEXT.....No Help Text
ID.....{548699B0-98CF-2F0F-F289-415228AFAC42}
TERM.....MILLILITERS PER BEAT

Page: 1035

UNIT: MILLILITERS PER BEAT PER METER SQUARE (BSA)

ABBREVIATION.....ml/beat/m2
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....ml/beat/m2(BSA)
HELP_TEXT.....No Help Text
ID.....{574CDB00-A831-EF35-4F2C-A98C0F58923F}
TERM.....MILLILITERS PER BEAT PER METER SQUARE
(BSA)

UNIT: MILLILITERS PER HOUR

ABBREVIATION.....ml/hr
ACTIVE.....Yes
DESCRIPTION.....ml/hr
DISPLAY_NAME.....ml/hr
HELP_TEXT.....No Help Text
ID.....{345899A8-9DED-418F-A53F-F8084A5601A1}
TERM.....MILLILITERS PER HOUR

UNIT: MILLILITERS PER METER SQUARE (BSA)

ABBREVIATION.....mL/M2
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....mL/M2
HELP_TEXT.....No Help Text

ID.....{42C265B4-76B3-44FC-985F-EA1EFCEDA244}
TERM.....MILLILITERS PER METER SQUARE (BSA)

UNIT: MILLILITERS PER MINUTE

ABBREVIATION.....ml/min
ACTIVE.....Yes
DESCRIPTION.....Milliliters per Minute
DISPLAY_NAME.....ml/min
HELP_TEXT.....Milliliters per Minute
ID.....{25F4A228-D614-457A-9885-8F74D4596897}
TERM.....MILLILITERS PER MINUTE

UNIT: MILLILITERS PER MINUTE PER METER SQUARE

ABBREVIATION.....ml/min/m2
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....ml/min/m2(BSA)
HELP_TEXT.....No Help Text
ID.....{6404DFC7-F8B4-3C97-80EF-96627E9C0270}
TERM.....MILLILITERS PER MINUTE PER METER SQUARE

Page: 1036

UNIT: MILLIMETERS

ABBREVIATION.....mm
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....mm
HELP_TEXT.....No Help Text
ID.....{23F389B2-8308-91DC-2620-EDFFF07B14A5}
TERM.....MILLIMETERS

UNIT: MILLIMETERS OF MERCURY

ABBREVIATION.....mm/Hg
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....mm Hg
HELP_TEXT.....No Help Text
ID.....{57CF7550-FD97-351E-DF60-56F8F68BE7C6}
TERM.....MILLIMETERS OF MERCURY

UNIT: MILLIVOLTS

ABBREVIATION.....mV
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....mV
HELP_TEXT.....No Help Text
ID.....{40B59357-C9F2-0E5C-DCD7-6C815811F90D}
TERM.....MILLIVOLTS

UNIT: MILLILITERS PER KILOGRAM PER HOUR

ABBREVIATION.....ml kg hr
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....ml/kg/hour
HELP_TEXT.....No Help Text
ID.....{2F63C17C-8246-4466-9DC2-C5E5BE1777C2}
TERM.....MILLILITERS PER KILOGRAM PER HOUR

UNIT: MILLILITERS PER KILOGRAM PER MINUTE

ABBREVIATION.....ml kg mi
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....ml/kg/min
HELP_TEXT.....No Help Text
ID.....{E564F111-8B23-404A-B610-7707A5B1195D}
TERM.....MILLILITERS PER KILOGRAM PER MINUTE

Page: 1037

UNIT: MINUTE

ABBREVIATION.....Min
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Minute
HELP_TEXT.....No Help Text
ID.....{0AF0056E-0EB1-3712-6974-A0F5DC60010B}
TERM.....MINUTE

UNIT: MODIFIED ALDRETE SCORE

ABBREVIATION.....MOD ALD SC
ACTIVE.....Yes
DESCRIPTION.....Modified Aldrete Scoring System Unit
DISPLAY_NAME.....Modified Aldrete Score
HELP_TEXT.....Unit of measure for Aldrete Scoring
System
ID.....{8D9A43E8-D4E6-4BBE-A8E6-72677342BC82}
TERM.....MODIFIED ALDRETE SCORE

UNIT: MONTH

ABBREVIATION.....Mo
ACTIVE.....Yes
DESCRIPTION.....Month of the year
DISPLAY_NAME.....Month
HELP_TEXT.....1 - 12
ID.....{8933A469-2D06-D3DA-C2C4-A8FB3DD894D6}
TERM.....MONTH

UNIT: MORSE SCORE

ABBREVIATION.....MORSE SCOR
ACTIVE.....Yes
DESCRIPTION.....Morse Score
DISPLAY_NAME.....Morse Score
HELP_TEXT.....Scale for recording Morse Score Values
ID.....{301282A4-2BC8-430A-853B-C817F9A61878}
TERM.....MORSE SCORE

UNIT: MULTIVOLTS

ABBREVIATION.....Multivolts
ACTIVE.....No
DESCRIPTION.....No Description
DISPLAY_NAME.....Multivolts
HELP_TEXT.....No Help Text
ID.....{2882958B-C91D-4B8C-BEE7-382F3D7CCC49}
TERM.....MULTIVOLTS

Page: 1038

UNIT: OUNCE

ABBREVIATION.....oz
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Ounce
HELP_TEXT.....No Help Text
ID.....{642336FB-7BC9-A1D5-99E5-05447DAD3B75}
TERM.....OUNCE

UNIT: PAINAD SCALE

ABBREVIATION.....PAINAD SC
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....PAINAD Scale
HELP_TEXT.....No Help Text
ID.....{D6CE83BB-3DBE-4E94-9E91-E21C07A12CAD}
TERM.....PAINAD SCALE

UNIT: PERCENTAGE

ABBREVIATION.....%
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Percentage
HELP_TEXT.....No Help Text
ID.....{0B3163CE-2C5C-B196-DBBE-A1BC6A970B2A}
TERM.....PERCENTAGE

UNIT: POINTS

ABBREVIATION.....POINTS
ACTIVE.....Yes
DESCRIPTION.....Points

DISPLAY_NAME.....Points
HELP_TEXT.....No Help Text
ID.....{A41F9244-F8A9-4CBD-AD62-32165D77414B}
TERM.....POINTS

UNIT: POUND

ABBREVIATION.....lb
ACTIVE.....Yes
DESCRIPTION.....US measurement of Weight
DISPLAY_NAME.....Pound
HELP_TEXT.....No Help Text
ID.....{F9CBB050-E809-52FF-6C2D-51B7E9557F98}
TERM.....POUND

Page: 1039

UNIT: QUAD LUMEN CATHETER

ABBREVIATION.....QUAD LUMEN
ACTIVE.....Yes
DESCRIPTION.....Quad Lumen Catheter
DISPLAY_NAME.....Quad Lumen Catheter
HELP_TEXT.....No Help Text
ID.....{C894DDF4-5976-4435-B0E4-70EC7A909565}
TERM.....QUAD LUMEN CATHETER

UNIT: RANGE

ABBREVIATION.....range
ACTIVE.....Yes
DESCRIPTION.....Defines the minimum and maximum
values in a range
DISPLAY_NAME.....Range
HELP_TEXT.....No Help Text
ID.....{DD90689F-3B19-BC1D-1875-E9EAB08C9C9E}
TERM.....RANGE

UNIT: SCALE

ABBREVIATION.....SCALE
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Scale
HELP_TEXT.....No Help Text
ID.....{57E423C4-98D9-0891-0982-0B84EFD2D132}
TERM.....SCALE

UNIT: SECOND

ABBREVIATION.....sec
ACTIVE.....Yes
DESCRIPTION.....Measurement of passage of time (1/60
of a minute)
DISPLAY_NAME.....Second

HELP_TEXT.....No Help Text
ID.....{82E7C9B1-38E7-D568-F42D-23F18AA08819}
TERM.....SECOND

UNIT: SEWS SCALE

ABBREVIATION.....SEWS SCALE
ACTIVE.....Yes
DESCRIPTION.....Severity Of Ethanol Withdrawl Scale
DISPLAY_NAME.....SEWS SCALE
HELP_TEXT.....Severity Of Ethanol Withdrawl Scale
ID.....{4C8DFD61-E4FB-4918-93B0-33DAAF4EF1D1}
TERM.....SEWS SCALE

Page: 1040

UNIT: SI UNIT

ABBREVIATION.....SI UNITS
ACTIVE.....Yes
DESCRIPTION.....Unit of measure for Specific Gravity
DISPLAY_NAME.....SI Unit
HELP_TEXT.....No Help Text
ID.....{42C760F2-BBFC-4523-82C1-953CB41B4016}
TERM.....SI UNIT

UNIT: SQUARE METERS

ABBREVIATION.....M2
ACTIVE.....Yes
DESCRIPTION.....Surface area measured in square meters
DISPLAY_NAME.....Square Meters
HELP_TEXT.....Surface area measured in square meters
ID.....{5E49A34D-CC24-ABDA-7320-14071BA88781}
TERM.....SQUARE METERS

UNIT: STROKE WORK PER BEAT

ABBREVIATION.....gm-M/beat
ACTIVE.....Yes
DESCRIPTION.....Stroke Work Per Beat
DISPLAY_NAME.....gm-M / beat
HELP_TEXT.....Stroke Work Per Beat
ID.....{61721D2C-02BF-14AA-EB7F-602194924117}
TERM.....STROKE WORK PER BEAT

UNIT: STROKE WORK PER BEAT PER M2

ABBREVIATION.....gm-M/B/M2
ACTIVE.....Yes
DESCRIPTION.....STROKE WORK PER BEAT PER M2
DISPLAY_NAME.....gm-M / Beat / M2
HELP_TEXT.....Systemic Vascular Resistance Index
ID.....{900EBF1F-B706-0419-1697-58FE11257EF3}
TERM.....STROKE WORK PER BEAT PER M2

UNIT: UNITS

ABBREVIATION.....UNITS
ACTIVE.....Yes
DESCRIPTION.....Units
DISPLAY_NAME.....Units
HELP_TEXT.....Units
ID.....{7E9E6B8E-8B47-4B06-81FF-227E46707401}
TERM.....UNITS

Page: 1041

UNIT: UNITS PER HOUR

ABBREVIATION.....units/hr
ACTIVE.....Yes
DESCRIPTION.....Units per Hour
DISPLAY_NAME.....units/hr
HELP_TEXT.....Units per Hour
ID.....{7AE61C6C-F8D9-4D6B-8A7C-AE07F93EB3D9}
TERM.....UNITS PER HOUR

UNIT: UNITS PER MINUTE

ABBREVIATION.....units/min
ACTIVE.....Yes
DESCRIPTION.....Units per Minute
DISPLAY_NAME.....Units/min
HELP_TEXT.....Units per Minute
ID.....{3C31F294-AD93-4509-908A-72F3995B01B9}
TERM.....UNITS PER MINUTE

UNIT: UNITS/kg/HOUR

ABBREVIATION.....UNkgHR
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Units/kg/Hour
HELP_TEXT.....No Help Text
ID.....{89A25E97-A6B8-49B2-934F-797C80C80711}
TERM.....UNITS/kg/HOUR

UNIT: YEAR

ABBREVIATION.....yr
ACTIVE.....Yes
DESCRIPTION.....No Description
DISPLAY_NAME.....Year
HELP_TEXT.....No Help Text
ID.....{F800BBF4-BC4F-6F5D-BB28-BD3EB5ADC34D}
TERM.....YEAR

UNIT: mg/dl

ABBREVIATION.....mg/dl
ACTIVE.....Yes
DESCRIPTION.....Milligrams per deciliter
DISPLAY_NAME.....mg/dl
HELP_TEXT.....No Help Text
ID.....{796F843B-3EF5-4E9B-AAD6-254088D77FFB}
TERM.....mg/dl

Page: 1042

UNIT: pH

ABBREVIATION.....pH
ACTIVE.....Yes
DESCRIPTION.....pH
DISPLAY_NAME.....pH
HELP_TEXT.....No Help Text
ID.....{EF524413-89B4-4ABA-BD68-1FC6713C623D}
TERM.....pH

- 0 -

0 (data).....631
0 (q-pair).....98, 318, 319, 320,
321, 322, 324, 383,
384, 488, 489, 490,
491, 492, 493, 494,
495, 497, 498
0 - ABSENT (data).....643
0 - ABSENT (q-pair).....108, 109, 116, 125,
140, 141, 167, 441
0 - NO ANXIETY AT EASE (data).....643
0 - NO ANXIETY AT EASE (q-pair).....107
0 - NO NAUSEA/VOMITING (data).....643
0 - NO NAUSEA/VOMITING (q-pair).....110
0 - NO SWEAT VISIBLE (data).....643
0 - NO SWEAT VISIBLE (q-pair).....112
0 - NO TREMOR (data).....643
0 - NO TREMOR (q-pair).....115
0 - NONE (data).....644
0 - NONE (q-pair).....113
0 - NORMAL ACTIVITY (data).....644
0 - NORMAL ACTIVITY (q-pair).....106
0 - ORIENTED AND CAN DO SERIAL ADDITIONS (data).....644
0 - ORIENTED AND CAN DO SERIAL ADDITIONS (q-pair).....111
0 = 0 EXTREMITIES (data).....644
0 = 0 EXTREMITIES (q-pair).....24
0 = <90% (data).....644
0 = <90% (q-pair).....27
0 = ABSENT PULSE (data).....645
0 = APNEIC (data).....645
0 = APNEIC (q-pair).....28
0 = BP +/- 50% PREOP (data).....645

0 = BP +/- 50% PREOP (q-pair).....25
 0 = NOT RESPONDING (data).....645
 0 = NOT RESPONDING (q-pair).....26

- 1 -

1 (data).....631
 1 (q-pair).....98, 150, 303, 383,
 384, 490, 492, 546,
 576
 1 - BARELY PERCEPTIBLE SWEATING / PALMS MOIST (data)..645
 1 - BARELY PERCEPTIBLE SWEATING / PALMS MOIST (q-pair).112
 1 - CANNOT DO SERIAL ADDITIONS (data).....646
 1 - CANNOT DO SERIAL ADDITIONS (q-pair).....111
 1 - MILD NAUSEA / NO VOMITING (data).....646
 1 - MILD NAUSEA / NO VOMITING (q-pair).....110
 1 - MILDLY ANXIOUS (data).....646
 1 - MILDLY ANXIOUS (q-pair).....107
 1 - NO RESPONSE (data).....646
 1 - NO RESPONSE (q-pair).....167, 168, 170
 1 - NOT VISIBLE BUT CAN BE FELT FINGERTIP (data).....646
 1 - NOT VISIBLE BUT CAN BE FELT FINGERTIP (q-pair).....115
 1 - VERY MILD (data).....647
 1 - VERY MILD (q-pair).....109
 1 - VERY MILD HARSHNESS OR ABILITY TO FRIGHTEN (data)..647
 1 - VERY MILD HARSHNESS OR ABILITY TO FRIGHTEN (q-pair)108
 1 - VERY MILD ITCHING PINS/NEEDLES BURNING (data).....647
 1 - VERY MILD ITCHING PINS/NEEDLES BURNING (q-pair)....113
 1 - VERY MILD SENSITIVITY (data).....647
 1 - VERY MILD SENSITIVITY (q-pair).....116
 1 - WEAK (data).....647
 1 - WEAK (q-pair).....441
 1 = 2 EXTREMITIES (data).....648
 1 = 2 EXTREMITIES (q-pair).....24
 1 = AROUSABLE (data).....648
 1 = AROUSABLE (q-pair).....26
 1 = BP +/- 20-49% PREOP (data).....648
 1 = BP +/- 20-49% PREOP (q-pair).....25
 1 = DYSPNEA (data).....648
 1 = DYSPNEA (q-pair).....28
 1 = NEEDS O2 (data).....648
 1 = NEEDS O2 (q-pair).....27
 1 HOUR (data).....649
 1 HOUR (q-pair).....399
 1+=PALPABLE, THRADY, EASILY OBLITERATED (data).....649
 1- MILD ANXIOUS (data).....649
 1- SOMEWHAT MORE THAN NORMAL ACTIVITY (data).....649
 1- SOMEWHAT MORE THAN NORMAL ACTIVITY (q-pair).....106
 10 (data).....634
 10 (q-pair).....20, 98, 303, 319,
 383, 384, 546, 576
 11 (data).....634
 11 (q-pair).....303
 12 (data).....634
 12 (q-pair).....303
 12.5 (data).....634

13 (data).....	635
13 (q-pair).....	153
13.5 (data).....	635
13.5 (q-pair).....	153
14 (data).....	635
14 (q-pair).....	153
14.5 (data).....	635
14.5 (q-pair).....	153
15 (data).....	635
15 (q-pair).....	153, 318, 322, 324
15.5 (data).....	636
15.5 (q-pair).....	153
16 (data).....	636
16 (q-pair).....	153
16.5 (data).....	636
16.5 (q-pair).....	153
17 (data).....	636
17 (q-pair).....	153
17.5 (data).....	636
17.5 (q-pair).....	153
18 (data).....	637
18 (q-pair).....	153
18.5 (data).....	637
18.5 (q-pair).....	153
1800 CALORIE CONTROLLED (data).....	649
1800 CALORIE CONTROLLED (q-pair).....	134
1800 DIABETIC (data).....	650
1800 DIABETIC (q-pair).....	134
19 (data).....	637
19 (q-pair).....	153
19.5 (data).....	637
19.5 (q-pair).....	153
1:1 (data).....	650
1:1 (q-pair).....	132, 200, 202
1:1 OBSERVATION (data).....	650
1:1 OBSERVATION (q-pair).....	53
1:2 (data).....	650
1:2 (q-pair).....	132, 200, 202
1:3 (data).....	650
1:3 (q-pair).....	132, 200, 202
1:4 (data).....	651
1:4 (q-pair).....	132, 200, 202, 565
1:8 (data).....	651
1:8 (q-pair).....	132, 200, 202
1A + 2V (data).....	651
1A + 2V (q-pair).....	376
1A + 3V (data).....	651
1A + 3V (q-pair).....	376
1A +1V (data).....	651
1A +1V (q-pair).....	376

- 2 -

2 (data).....	631
2 (q-pair).....	98, 106, 107, 110, 112, 115, 150, 303,

383, 384, 546, 576

2 - ABNORMAL EXTENSION FROM PAIN (data).....	652
2 - ABNORMAL EXTENSION FROM PAIN (q-pair).....	168
2 - BEADS OF SWEAT OBVIOUS ON FOREHEAD (data).....	652
2 - DISORIENTED FOR DATE < 2 CALENDAR DAYS (data).....	652
2 - DISORIENTED FOR DATE < 2 CALENDAR DAYS (q-pair)....	111
2 - INCOMPREHENSIBLE SOUNDS (data).....	652
2 - INCOMPREHENSIBLE SOUNDS (q-pair).....	170
2 - INTERMITTENT NAUSEA WITH DRY HEAVES (data).....	652
2 - MILD (data).....	653
2 - MILD (q-pair).....	109
2 - MILD HARSHNESS OR ABILITY TO FRIGHTEN (data).....	653
2 - MILD HARSHNESS OR ABILITY TO FRIGHTEN (q-pair).....	108
2 - MILD ITCHING PINS/NEEDLES BURNING NUMBNESS (data)..	653
2 - MILD ITCHING PINS/NEEDLES BURNING NUMBNESS (q-pair)	113
2 - MILD SENSITIVITY (data).....	653
2 - MILD SENSITIVITY (q-pair).....	116
2 - MODERATE WITH PATIENTS ARMS EXTENDED (data).....	653
2 - NORMAL (data).....	654
2 - NORMAL (q-pair).....	441
2 - TO PAIN (data).....	654
2 - TO PAIN (q-pair).....	167
2 = 4 EXTREMITIES (data).....	654
2 = 4 EXTREMITIES (q-pair).....	24
2 = BP +/- 20% PREOP (data).....	654
2 = BP +/- 20% PREOP (q-pair).....	25
2 = BREATHE FREELY (data).....	654
2 = BREATHE FREELY (q-pair).....	28
2 = FULLY AWAKE (data).....	655
2 = FULLY AWAKE (q-pair).....	26
2 = ON ROOM AIR (data).....	655
2 = ON ROOM AIR (q-pair).....	27
2+=NORMAL, NOT EASILY OBLITERATED (data).....	655
2- MILD HARSHNESS OR ABILITY TO FRIGHTEN (data).....	655
2.0 CAL FORMULA (data).....	655
20 (data).....	637
20 (q-pair).....	319, 321
20.5 (data).....	638
20.5 (q-pair).....	153
2000 CALORIE CONTROLLED (data).....	656
2000 CALORIE CONTROLLED (q-pair).....	134
2000 DIABETIC (data).....	656
2000 DIABETIC (q-pair).....	134
21 (data).....	638
21 (q-pair).....	153
21.5 (data).....	638
21.5 (q-pair).....	153
22 (data).....	638
22 (q-pair).....	153
22.5 (data).....	638
22.5 (q-pair).....	153
2200 CALORIE CONTROLLED (data).....	656
2200 CALORIE CONTROLLED (q-pair).....	134
2200 DIABETIC (data).....	656
2200 DIABETIC (q-pair).....	134
23 (data).....	639

23 (q-pair).....	153
23.5 (data).....	639
23.5 (q-pair).....	153
24 (data).....	639
24 (q-pair).....	153
24.5 (data).....	639
24.5 (q-pair).....	153
2400 CALORIE CONTROLLED (data).....	656
2400 CALORIE CONTROLLED (q-pair).....	134
2400 DIABETIC (data).....	657
2400 DIABETIC (q-pair).....	134
25 (data).....	639
25 (q-pair).....	153, 320
25 PERCENT ALBUMIN (data).....	657
25 PERCENT ALBUMIN (prod).....	988
25 PERCENT ALBUMIN (q-pair).....	121, 270, 276, 279
25.5 (data).....	640
25.5 (q-pair).....	153
26 (data).....	640
26 (q-pair).....	153
26.5 (data).....	640
26.5 (q-pair).....	153
27 (data).....	640
27 (q-pair).....	153
27.5 (data).....	640
27.5 (q-pair).....	153
28 (data).....	641
28 (q-pair).....	153
28.5 (data).....	641
28.5 (q-pair).....	153
29 (data).....	641
2:4 (data).....	657
2:4 (q-pair).....	565
2A (data).....	657
2A (q-pair).....	376
2A + 1V (data).....	657
2A + 1V (q-pair).....	376
2A + 2V (data).....	658
2A + 2V (q-pair).....	376
2V (data).....	658
2V (q-pair).....	376

- 3 -

3 (data).....	631
3 (q-pair).....	98, 106, 107, 110, 112, 115, 150, 303, 383, 384, 488, 489, 491, 493, 494, 495, 497, 498, 546, 576
3 - ABNORMAL FLEXION FROM PAIN (data).....	658
3 - ABNORMAL FLEXION FROM PAIN (q-pair).....	168
3 - BEADS OF SWEAT OBVIOUS ON FOREHEAD (data).....	658
3 - DISORIENTED FOR DATE > 2 CALENDAR DAYS (data).....	658
3 - DISORIENTED FOR DATE > 2 CALENDAR DAYS (q-pair)....	111
3 - INAPPROPRIATE WORDS (data).....	659

3 - INAPPROPRIATE WORDS (q-pair).....	170
3 - INTERMITTENT NAUSEA WITH DRY HEAVES (data).....	659
3 - MODERATE (data).....	659
3 - MODERATE (q-pair).....	109
3 - MODERATE HARSHNESS OR ABILITY TO FRIGHTEN (data)...	659
3 - MODERATE HARSHNESS OR ABILITY TO FRIGHTEN (q-pair).....	108
3 - MODERATE ITCHING PINS/NEEDLES BURNING (data).....	659
3 - MODERATE ITCHING PINS/NEEDLES BURNING (q-pair).....	113
3 - MODERATE SENSITIVITY (data).....	660
3 - MODERATE SENSITIVITY (q-pair).....	116
3 - MODERATE WITH PATIENTS ARMS EXTENDED (data).....	660
3 - STRONG (data).....	660
3 - STRONG (q-pair).....	441
3 - TO VERBAL COMMAND (data).....	660
3 - TO VERBAL COMMAND (q-pair).....	167
3 WAY FOLEY (data).....	660
3 WAY FOLEY (q-pair).....	182
3+=BOUNDING, EASILY PALPABLE, CANNOT BE OBLITERATED (data)	661
30 (data).....	641
30 (q-pair).....	318
31 (data).....	641
34 ml (data).....	661
34 ml (q-pair).....	194
3:4 (data).....	661
3:4 (q-pair).....	565

- 4 -

4 (data).....	631
4 (q-pair).....	20, 98, 150, 303, 383, 384, 546, 576
4 - BEADS OF SWEAT OBVIOUS ON FOREHEAD (data).....	661
4 - BEADS OF SWEAT OBVIOUS ON FOREHEAD (q-pair).....	112
4 - DISORIENTED AND CONVERSES (data).....	661
4 - DISORIENTED AND CONVERSES (q-pair).....	170
4 - DISORIENTED FOR PLACE OR PERSON (data).....	662
4 - DISORIENTED FOR PLACE OR PERSON (q-pair).....	111
4 - INTERMITTENT NAUSEA WITH DRY HEAVES (data).....	662
4 - INTERMITTENT NAUSEA WITH DRY HEAVES (q-pair).....	110
4 - MODERATE WITH PATIENTS ARMS EXTENDED (data).....	662
4 - MODERATE WITH PATIENTS ARMS EXTENDED (q-pair).....	115
4 - MODERATELY ANXIOUS (data).....	662
4 - MODERATELY ANXIOUS (q-pair).....	107
4 - MODERATELY FIDGETY AND RESTLESS (data).....	662
4 - MODERATELY FIDGETY AND RESTLESS (q-pair).....	106
4 - MODERATELY SEVERE (data).....	663
4 - MODERATELY SEVERE (q-pair).....	109
4 - MODERATELY SEVERE HALLUCINATIONS (data).....	663
4 - MODERATELY SEVERE HALLUCINATIONS (q-pair).....	108, 113, 116
4 - SPONTANEOUSLY (data).....	663
4 - SPONTANEOUSLY (q-pair).....	167
4 - WITHDRAWS FROM PAIN (data).....	663
4 - WITHDRAWS FROM PAIN (q-pair).....	168
4 HOUR (data).....	663
4 HOUR (q-pair).....	399

4.5 (data).....	632
4.5 (q-pair).....	20
40 ml (data).....	664
40 ml (q-pair).....	194
4:4 (data).....	664
4:4 (q-pair).....	565

- 5 -

5 (data).....	632
5 (q-pair).....	20, 98, 106, 107, 110, 112, 115, 150, 303, 383, 384, 546, 576
5 - CONSTANT NAUSEA W FREQUENT HEAVES AND VOMITING (data)	664
5 - DRENCHING SWEAT (data).....	664
5 - LOCALIZES PAIN (data).....	664
5 - LOCALIZES PAIN (q-pair).....	168
5 - ORIENTED AND CONVERSES (data).....	665
5 - ORIENTED AND CONVERSES (q-pair).....	170
5 - SEVERE (data).....	665
5 - SEVERE (q-pair).....	109
5 - SEVERE EVEN WITH ARMS NOT EXTENDED (data).....	665
5 - SEVERE HALLUCINATIONS (data).....	665
5 - SEVERE HALLUCINATIONS (q-pair).....	108, 113, 116
5 PERCENT ALBUMIN (data).....	665
5 PERCENT ALBUMIN (prod).....	988
5 PERCENT ALBUMIN (q-pair).....	121, 270, 276, 279
5.5 (data).....	632
5.5 (q-pair).....	20

- 6 -

6 (data).....	632
6 (q-pair).....	20, 98, 106, 107, 110, 112, 115, 303, 383, 384, 546, 576
6 - CONSTANT NAUSEA W FREQUENT HEAVES AND VOMITING (data)	666
6 - DRENCHING SWEAT (data).....	666
6 - EXTREMELY SEVERE HALLUCINATIONS (data).....	666
6 - EXTREMELY SEVERE HALLUCINATIONS (q-pair).....	108, 113, 116
6 - OBEYS COMMANDS (data).....	666
6 - OBEYS COMMANDS (q-pair).....	168
6 - SEVERE EVEN WITH ARMS NOT EXTENDED (data).....	666
6 - VERY SEVERE (data).....	667
6 - VERY SEVERE (q-pair).....	109
6 HOUR (data).....	667
6 HOUR (q-pair).....	399
6 PERCENT HETASTARCH (data).....	667
6 PERCENT HETASTARCH (prod).....	988
6 PERCENT HETASTARCH (q-pair).....	121, 270, 279
6.5 (data).....	632
6.5 (q-pair).....	20

- 7 -

7 (data).....	633
7 (q-pair).....	20, 98, 303, 383, 384, 546, 576
7 - ACUTE PANIC (data).....	667
7 - ACUTE PANIC (q-pair).....	107
7 - CONSTANT NAUSEA W FREQUENT HEAVES AND VOMITING (data)	667
7 - CONSTANT NAUSEA W FREQUENT HEAVES AND VOMITING (q-pair)	110
7 - CONTINUOUS HALLUCINATIONS (data).....	668
7 - CONTINUOUS HALLUCINATIONS (q-pair).....	108, 113, 116
7 - DRENCHING SWEAT (data).....	668
7 - DRENCHING SWEAT (q-pair).....	112
7 - EXTREMELY SEVERE (data).....	668
7 - EXTREMELY SEVERE (q-pair).....	109
7 - PACING BACK AND FORTH MOST OF THE INTERVIEW (data).....	668
7 - PACING BACK AND FORTH MOST OF THE INTERVIEW (q-pair)	106
7 - SEVERE EVEN WITH ARMS NOT EXTENDED (data).....	668
7 - SEVERE EVEN WITH ARMS NOT EXTENDED (q-pair).....	115
7.5 (data).....	633
7.5 (q-pair).....	20, 201

- 8 -

8 (data).....	633
8 (q-pair).....	20, 98, 201, 303, 383, 384, 546, 576
8 HOUR (data).....	669
8 HOUR (q-pair).....	399
8.5 (data).....	633
8.5 (q-pair).....	20, 201

- 9 -

9 (data).....	633
9 (q-pair).....	20, 98, 201, 303, 383, 384, 546, 576
9.5 (data).....	634
9.5 (q-pair).....	20, 201
99 (data).....	642
99 (q-pair).....	384

- + -

+1 = SLIGHT (data).....	642
+1 = SLIGHT (q-pair).....	140, 141
+2 = MODERATE (data).....	642
+2 = MODERATE (q-pair).....	140, 141
+3 = MARKED (data).....	642
+3 = MARKED (q-pair).....	140, 141
+4 = VERY MARKED (data).....	642
+4 = VERY MARKED (q-pair).....	140, 141

AAI (data).....	669
AAI (q-pair).....	367
AAIR (data).....	669
AAIR (q-pair).....	367
AAO (data).....	669
AAO (q-pair).....	367
ABDOMEN (loca).....	910
ABDOMEN (q-pair).....	141, 384, 523
ABDOMEN BOWEL SOUNDS.....	1
ABDOMEN LLQ (loca).....	910
ABDOMEN LLQ (q-pair).....	1, 384, 523
ABDOMEN LUQ (loca).....	910
ABDOMEN LUQ (q-pair).....	1, 384, 523
ABDOMEN RLQ (loca).....	910
ABDOMEN RLQ (q-pair).....	1, 384, 523
ABDOMEN RUQ (loca).....	910
ABDOMEN RUQ (q-pair).....	1, 384, 523
ABDOMINAL (data).....	669
ABDOMINAL (loca).....	911
ABDOMINAL (q-pair).....	105, 447
ABDOMINAL AORTA (loca).....	911
ABDOMINAL AORTA (q-pair).....	74, 81
ABDOMINAL CONDITION.....	2
ABDOMINAL GIRTH.....	3
ABLE TO AMBULATE (data).....	670
ABLE TO AMBULATE (q-pair).....	377, 568
ABRASION (data).....	670
ABRASION (q-pair).....	627
ABSENT (data).....	670
ABSENT (q-pair).....	1, 69, 74, 81, 126, 137, 160, 162, 185, 308, 310, 483, 484, 485, 486, 569, 572, 573
ABSENT TO MINIMAL BLEEDING (data).....	670
ABSENT TO MINIMAL BLEEDING (q-pair).....	377
ACCELERATING (data).....	670
ACCELERATING (q-pair).....	601
ACHING (data).....	671
ACHING (q-pair).....	381
ACTIVE RANGE OF MOTION (data).....	671
ACTIVE RANGE OF MOTION (q-pair).....	11
ACTIVITIES OF DAILY LIVING (data).....	671
ACTIVITIES OF DAILY LIVING (q-pair).....	145
ACTIVITY - ABD PILLOW (task).....	1010
ACTIVITY - AMBULATE (task).....	1010
ACTIVITY - CONT. PASSIVE MOTION DEV. (task).....	1010
ACTIVITY - ELEVATE HEAD OF BED (task).....	1010
ACTIVITY - KNEE IMMOBILIZER (task).....	1010
ACTIVITY - ORIENT TO WARD (task).....	1011
ACTIVITY - RANGE OF MOTION EXERCISES (task).....	1011
ACTIVITY - SIT IN CHAIR (task).....	1011
ACTIVITY - TURN PATIENT (task).....	1011
ACTUAL (q-pair).....	146, 189, 440, 558,

577, 607, 609

ACTUAL (qual).....	1005
ADL - ACTIVITIES OF DAILY LIVING.....	4
ADL - AMBULATION ASSISTANCE.....	5
ADL - AMBULATION DISTANCE.....	6
ADL - BATH ASSIST (task).....	1011
ADL - BATH COMPLETE (task).....	1012
ADL - BATH PARTIAL (task).....	1012
ADL - BATH SELF (task).....	1012
ADL - BED MOBILITY.....	7
ADL - CATHETER CARE (task).....	1012
ADL - DRESSING.....	8
ADL - EATING.....	9
ADL - EMOTIONAL SUPPORT.....	10
ADL - FEED PATIENT (task).....	1012
ADL - ORAL HYGIENE (task).....	1013
ADL - OTHER ACTIVITY.....	11
ADL - REPOSITION.....	12
ADL - SHAVE (task).....	1013
ADL - SKIN CARE (task).....	1013
ADL - TURN (task).....	1013
ADL - VISITORS.....	13
ADL- ORAL CARE.....	14
ADMISSION (data).....	671
ADMISSION (q-pair).....	129, 622
ADMISSION DIAGNOSIS (data).....	671
ADULT (data).....	672
ADULT (meth).....	962
ADULT (q-pair).....	59
ADULT CUFF (meth).....	962
ADULT LARGE (data).....	672
ADULT LARGE (q-pair).....	59
ADULT SMALL (data).....	672
ADULT SMALL (q-pair).....	59
ADVANCED CARE PLANNING (data).....	672
ADVANCED CARE PLANNING (q-pair).....	145
ADI - SPIRITUAL CARE.....	15
AEROSOL/HUMIDIFIED MASK (meth).....	962
AEROSOL/HUMIDIFIED MASK (q-pair).....	550
AFFECT.....	16
AFTER EXERCISE (meth).....	962
AFTER EXERCISE (q-pair).....	58, 440, 448
AGITATED (data).....	672
AGITATED (q-pair).....	52, 309
AGONAL (data).....	673
AGONAL (q-pair).....	79, 447
AGONIZING (data).....	673
AGONIZING (q-pair).....	381
AIR FLOW DEVICE (data).....	673
AIR FLOW DEVICE (q-pair).....	605
AIR LEAK (data).....	673
AIR LEAK (q-pair).....	103
AIRBORNE (data).....	673
AIRBORNE (q-pair).....	292
AIRWAY ARTIFICIAL CUFF STATUS.....	17
AIRWAY ARTIFICIAL POSITION AT LIP/NOSE.....	18

AIRWAY ARTIFICIAL SIZE - FRENCH.....	19
AIRWAY ARTIFICIAL SIZE - MILIMETERS.....	20
AIRWAY ARTIFICIAL TYPE.....	21
AIRWAY REFLEXES INTACT / SECRETIONS (data).....	674
AIRWAY REFLEXES INTACT / SECRETIONS (q-pair).....	377
ALARM PARAMETERS VERIFIED.....	22
ALARMS.....	23
ALBUMIN (data).....	674
ALBUMIN (prod).....	988
ALBUMIN (q-pair).....	279
ALDRETE - ACTIVITY.....	24
ALDRETE - CIRCULATION.....	25
ALDRETE - CONSCIOUSNESS.....	26
ALDRETE - O2 SATURATION.....	27
ALDRETE - RESPIRATION.....	28
ALDRETE - TOTAL.....	29
ALDRETE SCORE 8 OR GREATER (data).....	674
ALDRETE SCORE 8 OR GREATER (q-pair).....	377
ALDRETE SCORE AT BASELINE (data).....	674
ALDRETE SCORE AT BASELINE (q-pair).....	377
ALERT (data).....	674
ALERT (q-pair).....	52, 309
ALERT AND ORIENTED (data).....	675
ALERT AND ORIENTED (q-pair).....	135, 377
ALL DOWN (data).....	675
ALL DOWN (q-pair).....	499
ALL UP (data).....	675
ALL UP (q-pair).....	499
ALLERGY (data).....	675
ALLERGY (q-pair).....	38
ALOPECIA PARTIAL (data).....	675
ALOPECIA PARTIAL (q-pair).....	184
ALOPECIA TOTAL (data).....	676
ALOPECIA TOTAL (q-pair).....	184
AMBER (data).....	676
AMBER (q-pair).....	138, 174, 575, 614
AMBU BAG (meth).....	962
AMBU BAG (q-pair).....	551
AMBULATE AND ASSIST (data).....	676
AMBULATED (data).....	676
AMBULATED (q-pair).....	4
AMBULATORY CARE UNIT (data).....	676
AMBULATORY CARE UNIT (q-pair).....	379
AMPHOTERICIN B (prod).....	988
AMPHOTERICIN B (q-pair).....	275
AMPUTEE (data).....	677
AMPUTEE (q-pair).....	441
ANALGESIA (data).....	677
ANALGESIA (q-pair).....	145
ANGINAL (data).....	677
ANGINAL (q-pair).....	381
ANKLE (loca).....	911
ANKLE (q-pair).....	105
ANNOYING (data).....	677
ANNOYING (q-pair).....	381
ANTECUBITAL (qual).....	1005

ANTI-EMBOLUS STOCKINGS KNEE HIGH (data).....	677
ANTI-EMBOLUS STOCKINGS KNEE HIGH (q-pair).....	139
ANTI-EMBOLUS STOCKINGS THIGH HIGH (data).....	678
ANTI-EMBOLUS STOCKINGS THIGH HIGH (q-pair).....	139
ANTI-MICROBIAL DRESSING (data).....	678
ANTI-MICROBIAL DRESSING (q-pair).....	87, 619
ANTIBIOTIC (prod).....	989
ANTIBIOTIC (q-pair).....	275
ANTIBIOTIC OINTMENT (data).....	678
ANTIBIOTIC OINTMENT (q-pair).....	422
ANTIMICROBIAL DISK (data).....	678
ANTIMICROBIAL DISK (q-pair).....	619
ANURIC (data).....	678
ANURIC (q-pair).....	574
ANXIOUS (data).....	679
ANXIOUS (q-pair).....	52, 142
AOO (data).....	679
AOO (q-pair).....	367
APICAL (loca).....	911
APICAL (q-pair).....	440
APPEARANCE.....	30
APPROPRIATE (data).....	679
APPROPRIATE (q-pair).....	30
ARM BAND - ALLERGIES.....	31
ARM BAND - BLOOD TRANSFUSION IV.....	32
ARM BAND - DNR.....	33
ARM BAND - ELOPMENT.....	34
ARM BAND - HIGH FALL RISK.....	35
ARM BAND - PERSONAL MEDICAL ALERT.....	36
ARM BAND - WANDERER.....	37
ARM BANDS PRESENT.....	38
ARTERIAL LINE (meth).....	963
ARTERIAL LINE (q-pair).....	58, 133, 312, 313, 554
ARTERIAL LINE ALARMS ON.....	39
ARTERIAL LINE CARE (data).....	679
ARTERIAL LINE CARE (q-pair).....	302
ARTERIAL LINE CIRCULATION ASSESSED.....	40
ARTERIAL LINE CONDITION.....	41
ARTERIAL LINE DRESSING CHANGE.....	42
ARTERIAL LINE DRESSING TYPE.....	43
ARTERIAL LINE SITE CARE.....	44
ARTERIAL LINE SIZE.....	45
ARTERIAL LINE WAVEFORM.....	46
ARTERIAL LINE ZEROED.....	47
ARTIFICIAL AIRWAY (data).....	679
ARTIFICIAL AIRWAY (q-pair).....	437
ASHEN (data).....	680
ASHEN (q-pair).....	510
ASLEEP AT ROUNDS (data).....	680
ASLEEP AT ROUNDS (q-pair).....	529
ASPIRATION (data).....	680
ASPIRATION (q-pair).....	175
ASSAULTIVE (data).....	680
ASSAULTIVE (q-pair).....	54, 119
ASSESSMENT - CO (task).....	1013

ASSESSMENT - COLLECT VITALS (task).....	1014
ASSESSMENT - DOPPLERS (task).....	1014
ASSESSMENT - FINGERSTICK BLOOD SUGAR (task).....	1014
ASSESSMENT - HEMODYNAMIC READINGS (task).....	1014
ASSESSMENT - I&O (task).....	1014
ASSESSMENT - NEURO CHECKS (task).....	1015
ASSESSMENT - NEUROVASCULAR CHECKS (task).....	1015
ASSESSMENT - PCA (task).....	1015
ASSESSMENT - PULSE OX (task).....	1015
ASSESSMENT - RESTRAINTS (task).....	1015
ASSESSMENT - TELEMETRY (task).....	1016
ASSESSMENT - WEIGHT (task).....	1016
ASSIST CONTROL (data).....	680
ASSIST CONTROL (q-pair).....	586
ASSISTED (data).....	681
ASSISTED (q-pair).....	12, 14
ASSISTED VENTILATOR (meth).....	963
ASSISTED VENTILATOR (q-pair).....	448, 550
ASYSTOLE (data).....	681
ASYSTOLE (q-pair).....	79
AT REST (meth).....	963
AT REST (q-pair).....	58, 440, 448
ATRIAL (data).....	681
ATRIAL (q-pair).....	373
ATRIAL FIB/FLUTTER (data).....	681
ATRIAL FIB/FLUTTER (q-pair).....	79
ATRIAL FIBRILLATION (data).....	681
ATRIAL FIBRILLATION (q-pair).....	79
ATRIAL FLUTTER (data).....	682
ATRIAL FLUTTER (q-pair).....	79
AUSCULTATE (meth).....	963
AUSCULTATE (q-pair).....	440
AUSCULTATION (data).....	682
AUSCULTATION (q-pair).....	175
AV FISTULA (loca).....	911
AV FISTULA (q-pair).....	440
AVF (data).....	682
AVF (q-pair).....	80
AVI (data).....	682
AVI (q-pair).....	367
AVL (data).....	682
AVL (q-pair).....	80
AVO (data).....	683
AVO (q-pair).....	367
AVR (data).....	683
AVR (q-pair).....	80
AXILLARY (loca).....	912
AXILLARY (q-pair).....	558

- B -

BABINSKI.....	48
BACK (loca).....	912
BACK (q-pair).....	384, 523
BACTERICIDAL (data).....	683
BACTERICIDAL (q-pair).....	422

BACTERIOSTATIC (data).....	683
BACTERIOSTATIC (q-pair).....	422
BACTERIOSTATIC WATER (prod).....	989
BACTERIOSTATIC WATER (q-pair).....	275
BAG CHANGED (data).....	683
BAG CHANGED (q-pair).....	302, 400
BALSAM OF PERU (data).....	684
BALSAM OF PERU (q-pair).....	619
BARIATRIC FRAME (data).....	684
BARIATRIC FRAME (q-pair).....	530
BARIATRIC LOW AIRLOSS MATTRESS REPLACEMENT (data).....	684
BARIATRIC LOW AIRLOSS MATTRESS REPLACEMENT (q-pair)....	530
BARIATRIC PATIENT (data).....	684
BARIATRIC PATIENT (q-pair).....	419
BASAL RATE ONLY (data).....	684
BASAL RATE ONLY (q-pair).....	402
BATHROOM (data).....	685
BATHROOM (q-pair).....	51, 118
BATHROOM PRIVILEGES (data).....	685
BATHROOM PRIVILEGES (q-pair).....	11
BEATS PER MINUTE (unit).....	1026
BED (data).....	685
BED (meth).....	963
BED (q-pair).....	607
BED EXIT ALARM (data).....	685
BED EXIT ALARM (q-pair).....	23, 452
BED POSITION.....	49
BEDREST (data).....	685
BEDREST (q-pair).....	11
BEHAVIOR.....	50
BEHAVIOR CURRENT LOCATION.....	51
BEHAVIOR OBSERVED.....	52
BEHAVIOR PRECAUTION LEVEL.....	53
BEHAVIOR PRECAUTION REASON.....	54
BILATERAL (data).....	686
BILATERAL ANKLES (loca).....	912
BILATERAL ANKLES (q-pair).....	468, 470
BILATERAL FOOT (loca).....	912
BILATERAL FOOT (q-pair).....	160
BILATERAL HAND (loca).....	912
BILATERAL HAND (q-pair).....	185
BILATERAL LEG (loca).....	913
BILATERAL LEG (q-pair).....	310, 468
BILATERAL NECK (loca).....	913
BILATERAL NECK (q-pair).....	141
BILATERAL PERIPHERALS (loca).....	913
BILATERAL PERIPHERALS (q-pair).....	440
BILATERAL WRIST (loca).....	913
BILATERAL WRIST (q-pair).....	468, 469, 573
BIPAP (meth).....	964
BIPAP (q-pair).....	551
BIPAP - RESPIRONICS (data).....	686
BIPAP - RESPIRONICS (q-pair).....	600
BIPAP - VISION (data).....	686
BIPAP - VISION (q-pair).....	600
BIPAP-CPAP (meth).....	964

BIPAP-CPAP (q-pair).....	551
BIPAP-EXPIRATORY.....	55
BIPAP-INSPIRATORY.....	56
BIPAP/CPAP (data).....	686
BIPAP/CPAP (q-pair).....	600
BIPHASIC POSITIVE AIRWAY PRESSURE (data).....	686
BIPHASIC POSITIVE AIRWAY PRESSURE (q-pair).....	586, 600
BISMUTH GAUZE (data).....	687
BISMUTH GAUZE (q-pair).....	619
BITE (data).....	687
BITE (q-pair).....	627
BITE BLOCK (data).....	687
BITE BLOCK (q-pair).....	21
BLACK (data).....	687
BLACK (q-pair).....	138, 174, 545, 575, 610, 614
BLACK/NECROTIC (data).....	687
BLACK/NECROTIC (q-pair).....	411
BLADDER (loca).....	913
BLADDER (meth).....	964
BLADDER (q-pair).....	275, 384
BLADDER ASSESSMENT.....	57
BLADDER IRRIGATION (prod).....	989
BLADDER IRRIGATION (q-pair).....	337
BLAND (data).....	688
BLAND (q-pair).....	134
BLISTERING (data).....	688
BLISTERING (q-pair).....	627
BLOOD PRESSURE.....	58
BLOOD PRESSURE CUFF SIZE.....	59
BLOOD PRODUCT (prod).....	989
BLOOD PRODUCT (q-pair).....	270
BLOOD STREAKED (data).....	688
BLOOD TRANSFUSION IV (data).....	688
BLOOD TRANSFUSION IV (q-pair).....	38
BLOODY (data).....	688
BLOODY (q-pair).....	99, 258, 413, 544
BLUE (data).....	689
BLUE (q-pair).....	614
BLUNTED (data).....	689
BLUNTED (q-pair).....	16
BODY SURFACE AREA.....	60
BONE (data).....	689
BORDERLINE REDUCED (data).....	689
BORDERLINE REDUCED (q-pair).....	306
BOTH ARMS (loca).....	914
BOTH ARMS (q-pair).....	469
BOTH HANDS (loca).....	914
BOTH HANDS (q-pair).....	141
BOTH LEGS (loca).....	914
BOTH LEGS (q-pair).....	125, 141, 470
BOWEL MOVEMENT EVENT(S).....	61
BRACHIAL (loca).....	914
BRACHIAL (q-pair).....	58, 440
BRACHIAL (qual).....	1005
BRADEN FRICTION AND SHEAR 1 PROBLEM (data).....	689

BRADEN FRICTION AND SHEAR 1 PROBLEM (q-pair).....	63
BRADEN SCALE (unit).....	1026
BRADEN SCALE ACTIVITY.....	62
BRADEN SCALE ACTIVITY 1 BEDFAST (data).....	690
BRADEN SCALE ACTIVITY 1 BEDFAST (q-pair).....	62
BRADEN SCALE ACTIVITY 2 CHAIRFAST (data).....	690
BRADEN SCALE ACTIVITY 2 CHAIRFAST (q-pair).....	62
BRADEN SCALE ACTIVITY 3 WALKS OCCASIONALLY (data).....	690
BRADEN SCALE ACTIVITY 3 WALKS OCCASIONALLY (q-pair).....	62
BRADEN SCALE ACTIVITY 4 WALKS FREQUENTLY (data).....	690
BRADEN SCALE ACTIVITY 4 WALKS FREQUENTLY (q-pair).....	62
BRADEN SCALE FRICTION AND SHEAR.....	63
BRADEN SCALE FRICTION AND SHEAR 2 POTENTIAL PROBLEM (data)	
.....	690
BRADEN SCALE FRICTION AND SHEAR 2 POTENTIAL PROBLEM (q-pair)	
.....	63
BRADEN SCALE FRICTION AND SHEAR 3 NO APPARENT PROBLEM (data)	
.....	691
BRADEN SCALE FRICTION AND SHEAR 3 NO APPARENT PROBLEM (q-pair)	
.....	63
BRADEN SCALE MOBILITY.....	64
BRADEN SCALE MOBILITY 1 COMPLETELY IMMOBILE (data).....	691
BRADEN SCALE MOBILITY 1 COMPLETELY IMMOBILE (q-pair).....	64
BRADEN SCALE MOBILITY 2 VERY LIMITED (data).....	691
BRADEN SCALE MOBILITY 2 VERY LIMITED (q-pair).....	64
BRADEN SCALE MOBILITY 3 SLIGHTLY LIMITED (data).....	691
BRADEN SCALE MOBILITY 3 SLIGHTLY LIMITED (q-pair).....	64
BRADEN SCALE MOBILITY 4 NO LIMITATION (data).....	691
BRADEN SCALE MOBILITY 4 NO LIMITATION (q-pair).....	64
BRADEN SCALE MOISTURE.....	65
BRADEN SCALE MOISTURE 1 CONSTANTLY MOIST (data).....	692
BRADEN SCALE MOISTURE 1 CONSTANTLY MOIST (q-pair).....	65
BRADEN SCALE MOISTURE 2 VERY MOIST (data).....	692
BRADEN SCALE MOISTURE 2 VERY MOIST (q-pair).....	65
BRADEN SCALE MOISTURE 3 OCCASIONALLY MOIST (data).....	692
BRADEN SCALE MOISTURE 3 OCCASIONALLY MOIST (q-pair).....	65
BRADEN SCALE MOISTURE 4 RARELY MOIST (data).....	692
BRADEN SCALE MOISTURE 4 RARELY MOIST (q-pair).....	65
BRADEN SCALE NUTRITION.....	66
BRADEN SCALE NUTRITION 1 VERY POOR (data).....	692
BRADEN SCALE NUTRITION 1 VERY POOR (q-pair).....	66
BRADEN SCALE NUTRITION 2 PROBABLY INADEQUATE (data).....	693
BRADEN SCALE NUTRITION 2 PROBABLY INADEQUATE (q-pair).....	66
BRADEN SCALE NUTRITION 3 ADEQUATE (data).....	693
BRADEN SCALE NUTRITION 3 ADEQUATE (q-pair).....	66
BRADEN SCALE NUTRITION 4 EXCELLENT (data).....	693
BRADEN SCALE NUTRITION 4 EXCELLENT (q-pair).....	66
BRADEN SCALE SENSORY PERCEPTION.....	67
BRADEN SCALE SENSORY PERCEPTION 1 COMPLETELY LIMITED (data)	
.....	693
BRADEN SCALE SENSORY PERCEPTION 1 COMPLETELY LIMITED (q-pair)	
.....	67
BRADEN SCALE SENSORY PERCEPTION 2 VERY LIMITED (data).....	693
BRADEN SCALE SENSORY PERCEPTION 2 VERY LIMITED (q-pair).....	67
BRADEN SCALE SENSORY PERCEPTION 3 SLIGHTLY LIMITED (data)	
.....	694

BRADEN SCALE SENSORY PERCEPTION 3 SLIGHTLY LIMITED (q-pair)	67
BRADEN SCALE SENSORY PERCEPTION 4 NO IMPAIRMENT (data)	694
BRADEN SCALE SENSORY PERCEPTION 4 NO IMPAIRMENT (q-pair)	67
BRADEN SCALE TOTAL	68
BREAKDOWN (data)	694
BREAKDOWN (q-pair)	511
BREATH SOUNDS	69
BREATHS PER MINUTE (unit)	1026
BRIGHT (data)	694
BRIGHT (q-pair)	16
BRISK (data)	694
BRISK (q-pair)	442
BROWN (data)	695
BROWN (q-pair)	138, 174, 433, 545, 575, 614
BRUISED (data)	695
BRUISED (q-pair)	510, 627
BRUIT (data)	695
BRUIT (q-pair)	1
BUCKS TRACTION (data)	695
BUCKS TRACTION (q-pair)	563
BULB SUCTION (data)	695
BULB SUCTION (q-pair)	616
BULLET (data)	696
BULLET (q-pair)	627
BURN (data)	696
BURN (q-pair)	627
BURNING (data)	696
BURNING (q-pair)	381
BY PROVIDER ORDER (data)	696
BY PROVIDER ORDER (q-pair)	358

- C -

CADEXOMER IODINE (data)	696
CADEXOMER IODINE (q-pair)	422
CALCIUM ALGINATE (data)	697
CALCIUM ALGINATE (q-pair)	422, 619
CALF (loca)	914
CALF (q-pair)	105
CALL HOUSE OFFICER - BP DIASTOLIC (task)	1016
CALL HOUSE OFFICER - BP SYSTOLIC (task)	1016
CALL HOUSE OFFICER - PULSE (task)	1016
CALL HOUSE OFFICER - RESPIRATORY RATE (task)	1017
CALL HOUSE OFFICER - TEMP (task)	1017
CALL LIGHT WITHIN REACH	70
CALM (data)	697
CALM (q-pair)	52
CALORIE (unit)	1026
CANE (data)	697
CANE (q-pair)	5, 568
CANES (data)	697
CANES (q-pair)	568
CAPILLARY REFILL	71

CARDIAC ABNORMAL BPM.....	72
CARDIAC ALARM PARAMETERS VERIFIED.....	73
CARDIAC BRUIT.....	74
CARDIAC CHAIR (data).....	697
CARDIAC CHAIR (q-pair).....	49
CARDIAC DIET (data).....	698
CARDIAC DIET (q-pair).....	134
CARDIAC ECTOPY.....	75
CARDIAC HEART SOUNDS.....	76
CARDIAC INDEX.....	77
CARDIAC MONITOR (meth).....	964
CARDIAC OUTPUT.....	78
CARDIAC OUTPUT INJECTATE (data).....	698
CARDIAC OUTPUT INJECTATE (prod).....	989
CARDIAC OUTPUT INJECTATE (q-pair).....	337
CARDIAC OUTPUT MONITOR (meth).....	964
CARDIAC OUTPUT MONITOR (q-pair).....	95, 307, 311, 312, 428, 430, 439, 472, 473, 548, 549, 553
CARDIAC REPOSITION (data).....	698
CARDIAC REPOSITION (q-pair).....	12
CARDIAC RHYTHM.....	79
CARDIAC RHYTHM LEAD.....	80
CARDIAC THRILL.....	81
CARDIAC THRILL TIMING.....	82
CAREGIVER (data).....	698
CAREGIVER (q-pair).....	10
CAROTID (loca).....	915
CAROTID (q-pair).....	440
CASS ENDOTRACHEAL TUBE (data).....	698
CASS ENDOTRACHEAL TUBE (q-pair).....	21
CATHETER (data).....	699
CATHETER (meth).....	965
CATHETER (q-pair).....	183
CATHETER (sche).....	1008
CELLULITIS (data).....	699
CELLULITIS (q-pair).....	627
CENTIMETER (unit).....	1026
CENTIMETERS H2O (unit).....	1027
CENTRAL ECG (data).....	699
CENTRAL ECG (q-pair).....	23
CENTRAL LINE ALARMS ON.....	83
CENTRAL LINE CHANGED.....	84
CENTRAL LINE CONDITION.....	85
CENTRAL LINE DRESSING CHANGE.....	86
CENTRAL LINE DRESSING TYPE.....	87
CENTRAL LINE INSERTION MARK.....	88
CENTRAL LINE SITE CARE.....	89
CENTRAL LINE SIZE FRENCH.....	90
CENTRAL LINE SIZE GAUGE.....	91
CENTRAL LINE TYPE.....	92
CENTRAL LINE WAVEFORM.....	93
CENTRAL LINE ZEROED.....	94
CENTRAL VENOUS CATHETER (data).....	699
CENTRAL VENOUS CATHETER (q-pair).....	302, 400
CENTRAL VENOUS PRESSURE.....	95

CEREBRAL PERFUSION PRESSURE.....	96
CERVICAL SPINAL PRECAUTIONS (data).....	699
CERVICAL SPINAL PRECAUTIONS (q-pair).....	11
CHAIR (data).....	700
CHAIR (meth).....	965
CHAIR (q-pair).....	11, 568, 607
CHANGED (data).....	700
CHANGED (q-pair).....	42, 86, 302
CHEST (loca).....	915
CHEST (q-pair).....	146, 384, 468, 471, 523
CHEST EXPANSION.....	97
CHEST PAIN.....	98
CHEST TUBE (meth).....	965
CHEST TUBE (prod).....	990
CHEST TUBE (q-pair).....	343
CHEST TUBE DRAINAGE (prod).....	990
CHEST TUBE DRAINAGE CHARACTERISTICS.....	99
CHEST TUBE DRAINAGE SUCTION TYPE.....	100
CHEST TUBE DRAINAGE SYSTEM REPLACED.....	101
CHEST TUBE DRESSING CHANGE.....	102
CHEST TUBE SITE ASSESSMENT.....	103
CHEST TUBE SUCTION (CM H2O).....	104
CHEYNE-STOKES (data).....	700
CHEYNE-STOKES (q-pair).....	447
CIRCULATION CHECK (data).....	700
CIRCULATION CHECK (q-pair).....	467
CIRCUMFERENCE/GIRTH.....	105
CIWA SCALE (unit).....	1027
CIWA: AGITATION.....	106
CIWA: ANXIETY.....	107
CIWA: AUDITORY DISTURBANCES.....	108
CIWA: HEADACHE FULLNESS IN HEAD.....	109
CIWA: NAUSEA/VOMITING.....	110
CIWA: ORIENTATION AND CLOUDING OF SENSORIUM.....	111
CIWA: PAROXYSMAL SWEATS.....	112
CIWA: TACTILE DISTURBANCES.....	113
CIWA: TOTAL.....	114
CIWA: TREMOR.....	115
CIWA: VISUAL DISTURBANCES.....	116
CLAMPED (data).....	700
CLAMPED (q-pair).....	177
CLAY (data).....	701
CLAY (q-pair).....	138, 174, 575, 614
CLEAR (data).....	701
CLEAR (q-pair).....	69, 138, 174, 258, 413, 433, 501, 545, 574, 575
CLEAR LIQUID DIET (data).....	701
CLEAR LIQUID DIET (q-pair).....	134
CLICK (data).....	701
CLICK (q-pair).....	76
CLINICAL ASSESSMENT (data).....	701
CLINICAL ASSESSMENT (q-pair).....	145
CLOSE OBSERVATION (data).....	702
CLOSE OBSERVATION (q-pair).....	53

CLOSE OBSERVATION BEHAVIOR.....	117
CLOSE OBSERVATION CURRENT LOCATION.....	118
CLOSE OBSERVATION REASON FOR PRECAUTION.....	119
CLOSED (data).....	702
CLOSED (q-pair).....	122
CLOTHED SELF (data).....	702
CLOTHED SELF (q-pair).....	4
CLOTS (data).....	702
CLOTS (q-pair).....	544
CLOUDY (data).....	702
CLOUDY (q-pair).....	258, 413, 501, 574
CLUBBING (data).....	703
CLUBBING (q-pair).....	329
CMS INTACT DISTAL TO SURGICAL SITE (data).....	703
CMS INTACT DISTAL TO SURGICAL SITE (q-pair).....	377
COCCYX (loca).....	915
COCCYX (q-pair).....	384, 421, 523, 622
CODE STATUS.....	120
COFFEE GROUND (data).....	703
COLD (data).....	703
COLD (q-pair).....	386, 523
COLLAGEN (data).....	703
COLLAGEN (q-pair).....	422
COLLOID (meth).....	965
COLLOID (prod).....	990
COLLOID (q-pair).....	270, 283
COLLOID TYPE.....	121
COLORLESS (data).....	704
COLORLESS (q-pair).....	614
COLOSTOMY (data).....	704
COLOSTOMY (q-pair).....	336
COMA SCORE EYE.....	122
COMA SCORE MOTOR.....	123
COMA SCORE VERBAL.....	124
COMBATIVE (data).....	704
COMBATIVE (q-pair).....	52
COMMODO (data).....	704
COMMODO (q-pair).....	11
COMMUNICATION BARRIER (data).....	704
COMMUNICATION BARRIER (q-pair).....	142
COMPLETE BATH (data).....	705
COMPLETE BATH (q-pair).....	4
COMPLIANCE (data).....	705
COMPLIANCE (q-pair).....	145
COMPR. HOSE SEQ. COMP. DEVICE.....	125
COMPUTED TOMOGRAPHY (meth).....	965
COMPUTED TOMOGRAPHY (q-pair).....	305, 306
CONFUSED (data).....	705
CONFUSED (q-pair).....	124, 142
CONSISTENT WITH BODY TEMPERATURE (data).....	705
CONSISTENT WITH BODY TEMPERATURE (q-pair).....	523
CONSTANT (data).....	705
CONSTANT (q-pair).....	381, 382, 628
CONSTIPATION (data).....	706
CONSTIPATION (q-pair).....	172, 544
CONSTRICTED (data).....	706

CONSTRICTED (q-pair).....	16
CONTACT (data).....	706
CONTACT (q-pair).....	292
CONTEMPLATIVE (data).....	706
CONTEMPLATIVE (q-pair).....	144
CONTINENT (data).....	706
CONTINENT (q-pair).....	183, 577
CONTINUOUS (data).....	707
CONTINUOUS (q-pair).....	82, 177, 179, 487
CONTINUOUS POSITIVE AIR PRESSURE (data).....	707
CONTINUOUS POSITIVE AIR PRESSURE (q-pair).....	586, 600
CONTINUOUS SUCTION (data).....	707
CONTINUOUS SUCTION (q-pair).....	100
CONTRAST LEFT VENTRICULOGRAPHY (meth).....	966
CONTRAST LEFT VENTRICULOGRAPHY (q-pair).....	305, 306
CONTROL MODE VENTILATION (data).....	707
CONTROL MODE VENTILATION (q-pair).....	586
CONTROLLED VENTILATOR (meth).....	966
CONTROLLED VENTILATOR (q-pair).....	448
CONTUSION (data).....	707
CONTUSION (q-pair).....	627
COOL (data).....	708
COOL (q-pair).....	523
COOL CHEM PACK (data).....	708
COOL CHEM PACK (q-pair).....	605
COOLER THAN GENERAL BODY TEMPERATURE (data).....	708
COOLER THAN GENERAL BODY TEMPERATURE (q-pair).....	523
COOLER THAN OPPOSITE BODY PART (data).....	708
COOLER THAN OPPOSITE BODY PART (q-pair).....	523
COOLING BLANKET (data).....	708
COOLING BLANKET (q-pair).....	605
COOLING OFF (data).....	709
COOLING OFF (q-pair).....	606
COOLING ON (data).....	709
COOLING ON (q-pair).....	606
COOPERATIVE (data).....	709
COOPERATIVE (q-pair).....	52
COPIUS (data).....	709
COPIUS (q-pair).....	347, 432
CORE (loca).....	915
CORE (q-pair).....	558
COUGH.....	126
CPAP.....	127
CPAP (meth).....	966
CPAP (q-pair).....	551
CRACKLES (data).....	709
CRACKLES (q-pair).....	69
CRAMPING (data).....	710
CRAMPING (q-pair).....	172, 381
CREPITUS (data).....	710
CREPITUS (q-pair).....	103
CRRT (meth).....	966
CRRT (q-pair).....	271, 342
CRUSHING (data).....	710
CRUSHING (q-pair).....	381, 627
CRUTCH (data).....	710

CRUTCH (q-pair).....	568
CRUTCHES (data).....	710
CRUTCHES (q-pair).....	5, 568
CRYING (data).....	711
CRYING (q-pair).....	385
CRYOPRECIPITATE (prod).....	990
CRYOPRECIPITATE (q-pair).....	270, 276
CRYSTALLOID (data).....	711
CRYSTALLOID (prod).....	990
CRYSTALLOID (q-pair).....	270, 283
CUBIC CENTIMETER (unit).....	1027
CUFF (meth).....	966
CUFF (q-pair).....	133, 312, 313, 554
CUFF INFLATION CHECKED.....	128
CUFF-AUTOMATED (meth).....	967
CUFF-AUTOMATED (q-pair).....	58
CUFF-MANUAL (meth).....	967
CUFF-MANUAL (q-pair).....	58
CULTURAL (data).....	711
CULTURAL (q-pair).....	143
CVVHD (meth).....	967
CVVHD (q-pair).....	271, 342
CYANOTIC (data).....	711
CYANOTIC (q-pair).....	327, 328, 510

- D -

D - DOPPLER (data).....	711
D - DOPPLER (q-pair).....	441
D10W ENTERAL NUTRITION (data).....	712
D10W ENTERAL NUTRITION (q-pair).....	149
D50W ENTERAL NUTRITION (data).....	712
D50W ENTERAL NUTRITION (q-pair).....	149
D5W ENTERAL NUTRITION (data).....	712
D5W ENTERAL NUTRITION (q-pair).....	149
D=DOPPLER (data).....	712
DAILY PLAN (data).....	712
DAILY PLAN (q-pair).....	145
DAKINS SOLUTION (prod).....	991
DAKINS SOLUTION (q-pair).....	275
DAMPENED (data).....	713
DAMPENED (q-pair).....	46, 93, 355, 358, 362
DANGLED (data).....	713
DANGLED (q-pair).....	11
DAY (unit).....	1027
DAY ROOM (data).....	713
DAY ROOM (q-pair).....	51, 118
DDD (data).....	713
DDD (q-pair).....	367
DDDR (data).....	713
DDDR (q-pair).....	367
DDI (data).....	714
DDI (q-pair).....	367
DECELERATING (data).....	714
DECELERATING (q-pair).....	601

DECREASED (data).....	714
DECREASED (q-pair).....	526
DEEP (data).....	714
DEEP (q-pair).....	381, 447
DEGREES (ANGLE) (unit).....	1027
DEGREES C (unit).....	1028
DEGREES F (unit).....	1028
DENIES HAVING SLEPT (data).....	714
DENIES HAVING SLEPT (q-pair).....	529
DENIES NEED (data).....	715
DENIES NEED (q-pair).....	144
DESCRIP (data).....	715
DEVICE RELATED PRESSURE (data).....	715
DEVICE RELATED PRESSURE (q-pair).....	419
DEVICE RELEASE (data).....	715
DEVICE RELEASE (q-pair).....	467
DIABETIC/VASCULAR ULCER (data).....	715
DIABETIC/VASCULAR ULCER (q-pair).....	627
DIAGNOSIS TYPE.....	129
DIALYSIS (data).....	716
DIALYSIS (meth).....	967
DIALYSIS (q-pair).....	92
DIALYSIS CATHETER (meth).....	967
DIALYSIS CATHETER CONDITION.....	130
DIALYSIS DRESSING TYPE.....	131
DIALYSIS RATIO.....	132
DIAPHORETIC (data).....	716
DIAPHORETIC (q-pair).....	518
DIARRHEA (data).....	716
DIARRHEA (q-pair).....	172, 544
DIASTOLIC (data).....	716
DIASTOLIC (q-pair).....	82
DIASTOLIC PRESSURE.....	133
DIET.....	134
DILATED (data).....	716
DILATED (q-pair).....	442
DILUTE DAKINS SOLUTION (data).....	717
DILUTE DAKINS SOLUTION (q-pair).....	422
DILUTE SODIUM HYPOCHLORITE SOLUTION (data).....	717
DIMINISHED (data).....	717
DIMINISHED (q-pair).....	69
DISCHARGE (data).....	717
DISCHARGE (q-pair).....	129, 145, 622
DISCHARGE ACTIVITIES.....	135
DISCHARGE MEDS PLANNED (data).....	717
DISCHARGE MEDS PLANNED (q-pair).....	135
DISCHARGE ORDER ON CHART (data).....	718
DISCHARGE ORDER ON CHART (q-pair).....	135
DISCOLORATION (data).....	718
DISCOLORATION (q-pair).....	328, 505
DISCOMFORT (data).....	718
DISCOMFORT (q-pair).....	41, 85, 130, 197, 304, 353
DISCONTINUED (data).....	718
DISCONTINUED (q-pair).....	302, 400, 487
DISEASE PROCESS (data).....	718

DISEASE PROCESS (q-pair).....	145
DISHEVELED (data).....	719
DISHEVELED (q-pair).....	30
DISORIENTED (data).....	719
DISORIENTED (q-pair).....	464
DISORIENTED TO PERSON (data).....	719
DISORIENTED TO PERSON (q-pair).....	334
DISORIENTED TO PLACE (data).....	719
DISORIENTED TO PLACE (q-pair).....	334
DISORIENTED TO SITUATION (data).....	719
DISORIENTED TO SITUATION (q-pair).....	334
DISORIENTED TO TIME (data).....	720
DISORIENTED TO TIME (q-pair).....	334
DISTAL PORT (data).....	720
DISTAL PORT (q-pair).....	303
DISTENDED (data).....	720
DISTENDED (q-pair).....	2, 331
DIVERSION (data).....	720
DIVERSION (q-pair).....	577
DO NOT PROCESS.....	136
DO NOT RESUSCITATE (data).....	720
DO NOT RESUSCITATE (q-pair).....	38, 120
DONE (data).....	721
DONE (q-pair).....	14
DONOR SITE (loca).....	915
DONOR SITE (q-pair).....	265
DOO (data).....	721
DOO (q-pair).....	367
DOPPLER (meth).....	968
DOPPLER (q-pair).....	58, 133, 440, 554
DORSALIS PEDIS (loca).....	916
DORSALIS PEDIS (q-pair).....	440
DORSI FLEXION.....	137
DOUBLE (data).....	721
DOUBLE (q-pair).....	92
DOUBLE LUMEN CATHETER (meth).....	968
DOWN (data).....	721
DRAINAGE CHARACTERISTICS.....	138
DRESSING (data).....	721
DRESSING (q-pair).....	442, 611
DRESSING CHANGED (data).....	722
DRESSING CHANGED (q-pair).....	302, 400
DRESSING WITH PRESCRIPTION RX (data).....	722
DRESSINGS CHECKED (data).....	722
DRESSINGS CHECKED (q-pair).....	135
DROPLET (data).....	722
DROPLET (q-pair).....	292
DROWSY (data).....	722
DROWSY (q-pair).....	52, 309
DRY (data).....	723
DRY (q-pair).....	42, 86, 503, 518,
	607
DRY (qual).....	1005
DRY AND INTACT (data).....	723
DRY SHAMPOO (data).....	723
DRY SHAMPOO (q-pair).....	4

DRY STERILE DRESSING (data).....	723
DULL (data).....	723
DULL (q-pair).....	381
DUMPING SYNDROME (data).....	724
DUMPING SYNDROME (q-pair).....	134
DUODENOSTOMY (data).....	724
DUODENOSTOMY (q-pair).....	147
DUSKY (data).....	724
DUSKY (q-pair).....	510
DVI (data).....	724
DVI (q-pair).....	367
DVT PROPHYLAXIS.....	139
DYE PRESENT (data).....	724
DYE PRESENT (q-pair).....	575
DYNE * SECONDS / CM(-5) * M2 (unit).....	1028
DYNE * SECONDS / CM(-5) / M2 (unit).....	1028
DYNE-SEC/CM5 (unit).....	1028
DYSPHAGIA PHASE I (data).....	725
DYSPHAGIA PHASE I (q-pair).....	134
DYSPHAGIA STAGE II (data).....	725
DYSPHAGIA STAGE II (q-pair).....	134
DYSPHAGIA STAGE III (data).....	725
DYSPHAGIA STAGE III (q-pair).....	134
DYSPHAGIA STAGE IV (data).....	725
DYSPHAGIA STAGE IV (q-pair).....	134
DYSPHAGIA STAGE V (data).....	725
DYSPHAGIA STAGE V (q-pair).....	134
DYSPNEIC (data).....	726
DYSPNEIC (q-pair).....	447
DYSURIA (data).....	726
DYSURIA (q-pair).....	577

- E -

EAGER TO LEARN (data).....	726
EAGER TO LEARN (q-pair).....	144
ECCHYMOSIS (data).....	726
ECCHYMOSIS (q-pair).....	510
ECHOCARDIOGRAPHY (meth).....	968
ECHOCARDIOGRAPHY (q-pair).....	305, 306
EDEMA AMOUNT.....	140
EDEMA SEVERITY.....	141
EDUCATION - BARRIERS.....	142
EDUCATION - BELIEFS.....	143
EDUCATION - READINESS TO LEARN.....	144
EDUCATION - TOPIC.....	145
EF RESIDUAL (prod).....	991
EF RESIDUAL (q-pair).....	337
EKG LEAD.....	146
ELASTIC WRAP (data).....	726
ELASTIC WRAP (q-pair).....	139
ELBOW PAD (data).....	727
ELBOW PAD (q-pair).....	418
ELECTRIC (data).....	727
ELECTRIC (q-pair).....	381
ELECTROCARDIOGRAM (meth).....	968

ELECTROCARDIOGRAM (q-pair).....	440, 532, 533, 534, 535, 536, 537, 538, 539, 540, 541, 542, 543, 560, 561
ELECTROCARDIOGRAM TRACE 1 (meth).....	968
ELECTROCARDIOGRAM TRACE 1 (q-pair).....	440
ELECTROCARDIOGRAM TRACE 2 (meth).....	969
ELECTROCARDIOGRAM TRACE 2 (q-pair).....	440
ELEMENTAL FORMULA (data).....	727
ELEMENTAL FORMULA (prod).....	991
ELEMENTAL FORMULA (q-pair).....	270
ELOPEMENT (data).....	727
ELOPEMENT (q-pair).....	38, 54, 119
EMERGENCY DEPARTMENT (data).....	727
EMERGENCY DEPARTMENT (q-pair).....	152, 379
EMESIS (prod).....	991
EMESIS (q-pair).....	344
EMOTIONAL (data).....	728
EMOTIONAL (q-pair).....	142
END OF LIFE CARE (data).....	728
END OF LIFE CARE (q-pair).....	419
ENTERAL FEEDING / SUPPLEMENT (meth).....	969
ENTERAL FEEDING / SUPPLEMENT (q-pair).....	270
ENTERAL FEEDING TUBE (meth).....	969
ENTERAL FEEDING TUBE (q-pair).....	270
ENTERAL FEEDING TUBE TYPE.....	147
ENTERAL NUTRITION RATE.....	148
ENTERAL NUTRITION TYPE.....	149
ENVIRONMENT (data).....	728
ENVIRONMENT (q-pair).....	145
ENZYMATIC OINTMENT (data).....	728
ENZYMATIC OINTMENT (q-pair).....	422
EPICARDIAL LEADS (data).....	728
EPICARDIAL LEADS (q-pair).....	373
EPITHELIALIZATION (data).....	729
EPITHELIALIZATION (q-pair).....	621
EQUAL (data).....	729
EQUAL (q-pair).....	310, 573
EQUIPMENT (data).....	729
EQUIPMENT (q-pair).....	145
ER EMERGENCY SEVERITY INDEX.....	150
ER TRIAGE CONDITION.....	151
ER TRIAGE DISPOSITION.....	152
ERYTHEMA (data).....	729
ERYTHEMA (q-pair).....	505
ESCAPE (data).....	729
ESCAPE (q-pair).....	54, 119
ESCHAR (data).....	730
ESCHAR (q-pair).....	621
ESTIMATED (q-pair).....	189, 259, 347, 440, 558, 607, 609
ESTIMATED (qual).....	1005
ESTIMATED BY ARM SPAN (q-pair).....	189
ESTIMATED BY ARM SPAN (qual).....	1006
ET LOCATION.....	153
ETT REPOSITIONED (data).....	730

ETT REPOSITIONED (q-pair).....	14
ETT RETAPED/SECURED (data).....	730
ETT RETAPED/SECURED (q-pair).....	14
EVENT(S) (unit).....	1029
EVENTS PER HOUR (unit).....	1029
EWALD TUBE (data).....	730
EWALD TUBE (q-pair).....	178
EXPIRATORY WHEEZES (data).....	730
EXPIRATORY WHEEZES (q-pair).....	69
EXTENDS AROUND INCISION EXTENDS INTO PERIWOUND LESS THAN 1 CM (data)	731
EXTENDS AROUND INCISION EXTENDS INTO PERIWOUND LESS THAN 1 CM (q-pair)	513
EXTENDS AROUND INCISION MORE THAN 1 CM (data).....	731
EXTENDS AROUND INCISION MORE THAN 1 CM (q-pair).....	513
EXTENDS AROUND INCISION NONE (data).....	731
EXTENDS AROUND INCISION NONE (q-pair).....	513
EXTENDS INTO PERIWOUND <1CM (data).....	731
EXTENDS INTO PERIWOUND <1CM (q-pair).....	263
EXTENDS INTO PERIWOUND>1CM (data).....	731
EXTENDS INTO PERIWOUND>1CM (q-pair).....	263
EXTENDS TO PAIN (data).....	732
EXTENDS TO PAIN (q-pair).....	123
EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND LESS THAN 1 CM (data)	732
EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND LESS THAN 1 CM (q-pair)	514
EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND MORE THAN 1 CM (data)	732
EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND MORE THAN 1 CM (q-pair)	514
EXTENSION AROUND WOUND LIMITED TO WOUND (data).....	732
EXTENSION AROUND WOUND LIMITED TO WOUND (q-pair).....	514
EXTENSION AROUND WOUND NONE (data).....	732
EXTENSION AROUND WOUND NONE (q-pair).....	514
EXTENSIVE ASSISTANCE (data).....	733
EXTENSIVE ASSISTANCE (q-pair).....	7, 8, 9
EXTERNAL CATHETER (data).....	733
EXTERNAL CATHETER (q-pair).....	182, 577
EXTREMELY DRY (data).....	733
EYE CLOSED (data).....	733
EYE CLOSED (q-pair).....	442

- F -

FACE (loca).....	916
FACE (q-pair).....	141, 384, 523
FACE TENT (meth).....	969
FACE TENT (q-pair).....	550
FACE WASH (data).....	733
FACE WASH (q-pair).....	4
FALL PRECAUTIONS.....	154
FAMILY (data).....	734
FAMILY (q-pair).....	10
FAMILY INFORMED (data).....	734
FAMILY INFORMED (q-pair).....	317, 450

FAMILY ORIENTED TO UNIT.....	155
FAMILY/FRIEND (data).....	734
FAMILY/FRIEND (q-pair).....	567
FASCIA (data).....	734
FAT RESTRICTED 50 GM DIET (data).....	734
FAT RESTRICTED 50 GM DIET (q-pair).....	134
FECAL FAT TEST DIET (data).....	735
FECAL FAT TEST DIET (q-pair).....	134
FECES (prod).....	991
FEET (loca).....	916
FEET (q-pair).....	141
FEMORAL (loca).....	916
FEMORAL (q-pair).....	58, 440
FETAL HEART TONES.....	156
FFP (prod).....	992
FFP (q-pair).....	270, 276
FFP - THAWED (prod).....	992
FFP - THAWED (q-pair).....	270, 276
FHWN1 (data).....	735
FIBER FORMULA (data).....	735
FIBRINOUS (data).....	735
FIBRINOUS (q-pair).....	621
FINGERSTICK GLUCOSE.....	157
FIRM (data).....	735
FIRM (q-pair).....	2
FIRST DEGREE AV BLOCK (data).....	736
FIRST DEGREE AV BLOCK (q-pair).....	79
FIXED (data).....	736
FIXED (q-pair).....	442
FLAT (data).....	736
FLAT (q-pair).....	2, 16, 331
FLATULENCE (data).....	736
FLATULENCE (q-pair).....	1
FLATUS.....	158
FLEXION TO PAIN (data).....	736
FLEXION TO PAIN (q-pair).....	123
FLICKER OF CONTRACTION (data).....	737
FLICKER OF CONTRACTION (q-pair).....	325
FLOTATION BED (data).....	737
FLOTATION BED (q-pair).....	530
FLOWSHEET ADDENDUM.....	159
FLUCTUANCE (data).....	737
FLUCTUATES (data).....	737
FLUCTUATES (q-pair).....	103
FLUSHED (data).....	737
FLUSHED (q-pair).....	302, 510
FOAM AIR HEEL ELEVATOR (data).....	738
FOAM AIR HEEL ELEVATOR (q-pair).....	418
FOAM DRESSING (data).....	738
FOAM DRESSING (q-pair).....	619
FOLEY (data).....	738
FOLEY (q-pair).....	182
FOOT (unit).....	1029
FOOT CARE (data).....	738
FOOT CARE (q-pair).....	4
FOOT DOWN (data).....	738

FOOT DOWN (q-pair).....	499
FOOT STRENGTH.....	160
FOOT UP (data).....	739
FOOT UP (q-pair).....	499
FOREHEAD (loca).....	916
FOREHEAD (q-pair).....	421
FORMULA - HIGH FIBER (prod).....	992
FORMULA - PROBALANCE (prod).....	992
FORMULA - PROBALANCE (q-pair).....	149
FORMULA 1.5 CAL (prod).....	992
FORMULA 1.5 CAL (q-pair).....	270
FORMULA 2.0 CAL (prod).....	993
FORMULA 2.0 CAL (q-pair).....	149, 270
FORMULA D10W ENTERAL NUTRITION (prod).....	993
FORMULA D10W ENTERAL NUTRITION (q-pair).....	149
FORMULA D50W ENTERAL NUTRITION (prod).....	993
FORMULA D50W ENTERAL NUTRITION (q-pair).....	149
FORMULA D5W ENTERAL NUTRITION (prod).....	993
FORMULA D5W ENTERAL NUTRITION (q-pair).....	149
FORMULA FIBER (prod).....	993
FORMULA FIBER (q-pair).....	270
FOUL SMELLING (data).....	739
FOUL SMELLING (q-pair).....	416, 615
FOWLER'S (data).....	739
FOWLER'S (q-pair).....	12
FRANK BLOOD (data).....	739
FRANK BLOOD (q-pair).....	544
FRENCH (unit).....	1029
FREQUENT VERBAL REMINDERS (data).....	739
FREQUENT VERBAL REMINDERS (q-pair).....	456
FRIEND (data).....	740
FRIEND (q-pair).....	10
FRONTAL (loca).....	917
FULL CODE (data).....	740
FULL CODE (q-pair).....	120
FULL LIQUID (prod).....	994
FULL LIQUID (q-pair).....	288
FULL LIQUID DIABETIC (data).....	740
FULL LIQUID DIABETIC (q-pair).....	134
FULL LIQUID DIET (data).....	740
FULL LIQUID DIET (q-pair).....	134
FULL SPINAL PRECAUTIONS (data).....	740
FULL SPINAL PRECAUTIONS (q-pair).....	11
FULL THICK SKIN LOSS INVOLVING MUSCLE, BONE, SUPPO (data)	
.....	741
FULL THICK SKIN LOSS INVOLVING MUSCLE, BONE, SUPPO (q-pair)	
.....	421
FULL THICK SKIN LOSS INVOLVING SUBCUTANEOUS TISSUE (data)	
.....	741
FULL THICK SKIN LOSS INVOLVING SUBCUTANEOUS TISSUE (q-pair)	
.....	421
FUNDAL HEIGHT.....	161

- G -

GAG REFLEX.....	162
-----------------	-----

GAIT BELT (data).....	741
GAIT BELT (q-pair).....	5
GASTRIC BYPASS CLEAR LIQUID (data).....	741
GASTRIC BYPASS CLEAR LIQUID (q-pair).....	134
GASTRIC BYPASS PUREED (data).....	741
GASTRIC BYPASS PUREED (q-pair).....	134
GASTRIC BYPASS REGULAR (data).....	742
GASTRIC BYPASS REGULAR (q-pair).....	134
GASTRIC DRAINAGE (prod).....	994
GASTRIC DRAINAGE (q-pair).....	337
GASTRIC RESIDUAL AMOUNT DISCARDED.....	163
GASTRIC RESIDUAL CHECKED.....	164
GASTRIC RESIDUAL NOT DISCARDED.....	165
GASTRIC TUBE (data).....	742
GASTRIC TUBE (q-pair).....	147, 178
GASTRIC pH.....	166
GASTROINTESTINAL (loca).....	917
GASTROINTESTINAL (q-pair).....	384
GASTROSTOMY FEEDING BUTTON (data).....	742
GASTROSTOMY FEEDING BUTTON (q-pair).....	178
GASTROSTOMY TUBE (data).....	742
GASTROSTOMY TUBE (q-pair).....	178
GAUGE (unit).....	1029
GAUZE (data).....	742
GAUZE (q-pair).....	422, 619
GAUZE AND TAPE (data).....	743
GAUZE AND TAPE (q-pair).....	43, 198
GAUZE PACKING (data).....	743
GAUZE PACKING (q-pair).....	422, 619
GAUZE WRAP (data).....	743
GAUZE WRAP (q-pair).....	422, 619
GCS - EYE OPENING.....	167
GCS - MOTOR RESPONSE.....	168
GCS - TOTAL.....	169
GCS - VERBAL RESPONSE.....	170
GENERAL UMBILICUS (loca).....	917
GENERAL UMBILICUS (q-pair).....	1
GENERALIZED (loca).....	917
GENERALIZED (q-pair).....	141, 622
GENITAL ASSESSMENT.....	171
GENITALS (loca).....	917
GENITALS (q-pair).....	622
GERI-CHAIR (data).....	743
GERI-CHAIR (q-pair).....	11, 117, 471
GI SYMPTOMS.....	172
GI TUBE (meth).....	969
GI TUBE DRAINAGE CHARACTERISTICS.....	173
GI TUBE DRAINAGE COLOR.....	174
GI TUBE PLACEMENT CONFIRMATION.....	175
GI TUBE SIZE.....	176
GI TUBE SUCTION TYPE.....	177
GI TUBE TYPE.....	178
GIVEN DISCHARGE INSTRUCTION AND VERBALIZES UNDERSTANDING (data)	743
GIVEN DISCHARGE INSTRUCTION AND VERBALIZES UNDERSTANDING (q-pair)	135

GLASGOW COMA SCALE (unit).....	1030
GLUCERNA DIABETIC (prod).....	994
GLUCERNA DIABETIC (q-pair).....	149
GLUTEN RESTRICTED (data).....	744
GLUTEN RESTRICTED (q-pair).....	134
GM-M/M2/BEAT (unit).....	1030
GNAWING (data).....	744
GNAWING (q-pair).....	381
GOOD (data).....	744
GOOD (q-pair).....	526
GRAMS PER HOUR (unit).....	1030
GRAMS PER MINUTE (unit).....	1030
GRANULATING (data).....	744
GRANULATING (q-pair).....	621
GRANULOCYTES PHERESIS (prod).....	994
GRANULOCYTES PHERESIS (q-pair).....	270, 276
GRAVITY (data).....	744
GRAVITY (q-pair).....	100, 177
GRAY (data).....	745
GRAY (q-pair).....	610
GREEN (data).....	745
GREEN (q-pair).....	138, 174, 433, 545, 575, 610, 614
GRIMACING (data).....	745
GRIMACING (q-pair).....	385
GROUP ROOM (data).....	745
GROUP ROOM (q-pair).....	51, 118
GU BLADDER IRRIGATION FREQUENCY.....	179
GU CATHETER CARE PERFORMED.....	180
GU CATHETER SIZE.....	181
GU CATHETER TYPE.....	182
GU VOIDING.....	183
GUARDING (data).....	745
GUARDING (q-pair).....	2, 385
GURNEY CHAIR (data).....	746
GURNEY CHAIR (q-pair).....	11

- H -

HAIR CONDITION.....	184
HALL (data).....	746
HALL (q-pair).....	51, 118
HAND GRIP.....	185
HARD (data).....	746
HARD (q-pair).....	2, 503
HAS NOT SLEPT (data).....	746
HAS NOT SLEPT (q-pair).....	529
HEAD (loca).....	918
HEAD (q-pair).....	105, 384, 523
HEAD CIRCUMFERENCE.....	186
HEAD DOWN (data).....	746
HEAD DOWN (q-pair).....	499
HEAD OCCIPITAL (loca).....	918
HEAD OCCIPITAL (q-pair).....	421, 622
HEAD OF BED ELEVATION.....	187
HEAD UP (data).....	747

HEAD UP (q-pair).....	499
HEARING (data).....	747
HEARING (q-pair).....	142
HEART (loca).....	918
HEART (q-pair).....	81
HEARTBURN (data).....	747
HEARTBURN (q-pair).....	172
HEAT (data).....	747
HEAT (q-pair).....	386
HEAT CHEM PACK (data).....	747
HEAT CHEM PACK (q-pair).....	605
HEAT LAMP (data).....	748
HEAT LAMP (q-pair).....	605
HEAT/MOISTURE EXCHANGER PRESENT.....	188
HEAVY (data).....	748
HEAVY (q-pair).....	381, 512
HEEL LIFT (data).....	748
HEEL LIFT (q-pair).....	418
HEEL PILLOW (data).....	748
HEEL PILLOW (q-pair).....	418
HEIGHT.....	189
HEMATOMA (data).....	748
HEMATOMA (q-pair).....	510, 627
HEMATURIA (data).....	749
HEMATURIA (q-pair).....	574
HEME (prod).....	994
HEME (q-pair).....	337
HEMODIALYSIS (meth).....	970
HEMODIALYSIS (prod).....	995
HEMODIALYSIS (q-pair).....	271, 337, 342
HEMORRHAGE (data).....	749
HEMORRHAGE (q-pair).....	328
HEMOVAC (data).....	749
HEMOVAC (meth).....	970
HEMOVAC (q-pair).....	343, 616
HEPATIC FORMULA (data).....	749
HEPATIC FORMULA (prod).....	995
HEPATIC FORMULA (q-pair).....	270
HEPLOCK (data).....	749
HEPLOCK (q-pair).....	302, 400
HIGH (data).....	750
HIGH (q-pair).....	49, 177
HIGH CALORIE (data).....	750
HIGH CALORIE (q-pair).....	134
HIGH FALL RISK (data).....	750
HIGH FALL RISK (q-pair).....	38
HIGH FIBER (data).....	750
HIGH FIBER (q-pair).....	134
HIGH FOWLER'S (data).....	750
HIGH FOWLER'S (q-pair).....	12
HIGH PROTEIN (data).....	751
HIGH PROTEIN (q-pair).....	134
HIGH PROTEIN FORMULA (prod).....	995
HIGH PROTEIN FORMULA (q-pair).....	149
HIP (loca).....	918
HIP (q-pair).....	105

HIRSUTISM (data).....	751
HIRSUTISM (q-pair).....	184
HOMANS SIGN.....	190
HOME (data).....	751
HOME (q-pair).....	379
HOME HEALTH (data).....	751
HOME HEALTH (q-pair).....	145
HOSTILE (data).....	751
HOSTILE (q-pair).....	52
HOT (data).....	752
HOT (q-pair).....	523
HOUR (unit).....	1030
HUMIDIFIER (data).....	752
HUMIDIFIER (q-pair).....	605
HYDRATION - DIET (task).....	1017
HYDRATION - HEPARIN DRIP (task).....	1017
HYDRATION - HYPERAL (task).....	1017
HYDRATION - IV FLUIDS (task).....	1018
HYDRATION - NPO (task).....	1018
HYDRATION - PERIPHERAL IV (task).....	1018
HYDRATION - PICC/CENTRAL LINE (task).....	1018
HYDRATION - SUPPLEMENT SNACK (task).....	1018
HYDRATION - TUBE FEEDING (task).....	1019
HYDROCOLLOIDAL DRESSING (data).....	752
HYDROCOLLOIDAL DRESSING (q-pair).....	422, 619, 626
HYDROFIBER (data).....	752
HYDROFIBER (q-pair).....	422
HYDROFIBER WITH SILVER (data).....	752
HYDROFIBER WITH SILVER (q-pair).....	422
HYDROGEL DRESSING (data).....	753
HYDROGEL DRESSING (q-pair).....	422, 619
HYPERACTIVE (data).....	753
HYPERACTIVE (q-pair).....	1
HYPERDYNAMIC (data).....	753
HYPERDYNAMIC (q-pair).....	46, 93, 306, 355, 362
HYPERGRANULIZATION (data).....	753
HYPERGRANULIZATION (q-pair).....	621
HYPERTONIC HYDROGEL (data).....	753
HYPERTONIC HYDROGEL (q-pair).....	619
HYPOACTIVE (data).....	754
HYPOACTIVE (q-pair).....	1
HYPOALBUMINEMIA (data).....	754
HYPOALBUMINEMIA (q-pair).....	419

- I -

I (data).....	754
I (q-pair).....	80
I:E RATIO.....	191
IABP (data).....	754
IABP (meth).....	970
IABP (q-pair).....	23, 313
IABP ASSISTED SYSTOLIC PRESSURE.....	192
IABP BAEDP.....	193
IABP BALLOON SIZE.....	194

IABP CATHETER FLUSH.....	195
IABP CIRCULATION ASSESSED.....	196
IABP CONDITION.....	197
IABP DRESSING TYPE.....	198
IABP PEAK DIASTOLIC AUGMENTED PRESSURE.....	199
IABP RATIO.....	200
IABP SIZE.....	201
IABP TIMING.....	202
IABP TIMING CHECK.....	203
IABP UAEDP.....	204
IABP UNASSISTED SYSTOLIC PRESSURE.....	205
ICE PACK (data).....	754
ICE PACK (q-pair).....	605
ICTERIC (data).....	755
ICTERIC (q-pair).....	510
ICU (data).....	755
ICU (q-pair).....	379
IDENTIFICATION (data).....	755
IDENTIFICATION (q-pair).....	38
IDIOVENTRICULAR (data).....	755
IDIOVENTRICULAR (q-pair).....	79
II (data).....	755
II (q-pair).....	80
III (data).....	756
III (q-pair).....	80
ILEAL CONDUIT (data).....	756
ILEAL CONDUIT (q-pair).....	182, 577
ILEOSTOMY (data).....	756
ILEOSTOMY (q-pair).....	336
IMAGING (data).....	756
IMAGING (q-pair).....	379
IMPAIRED LISTENING SKILLS (data).....	756
IMPAIRED LISTENING SKILLS (q-pair).....	142
INACTIVE - ARTERIAL LINE SITE.....	206
INACTIVE - ARTERIAL LINE TYPE.....	207
INACTIVE - AUDIOMETRY.....	208
INACTIVE - BISPECTRAL INDEX.....	209
INACTIVE - BODY MASS INDEX.....	210
INACTIVE - CARDIAC AV FISTULA BRUIT.....	211
INACTIVE - CATHETER TYPE.....	212
INACTIVE - CENTRAL LINE SITE.....	213
INACTIVE - CHEST TUBE DRESSING INTACT.....	214
INACTIVE - CHEST TUBE SITE.....	215
INACTIVE - COUGH AND DEEP BREATH EVENTS.....	216
INACTIVE - DATE.....	217
INACTIVE - DATE DAY.....	218
INACTIVE - DATE MONTH.....	219
INACTIVE - DATE YEAR.....	220
INACTIVE - DIALYSIS CATHETER SITE.....	221
INACTIVE - EDEMA LOCATION.....	222
INACTIVE - ENTERAL NUTRITION FORMULA.....	223
INACTIVE - GI EMESIS APPEARANCE.....	224
INACTIVE - HEARING.....	225
INACTIVE - IABP INSERTION SITE.....	226
INACTIVE - IABP SITE.....	227
INACTIVE - INCISION LOCATION.....	228

INACTIVE - INCISION SITE.....	229
INACTIVE - IV DOSE RATE.....	230
INACTIVE - IV FLUID.....	231
INACTIVE - IV INSERTION DATE.....	232
INACTIVE - IV SITE.....	233
INACTIVE - OR OUTPUT - CHEST TUBE.....	234
INACTIVE - OR OUTPUT - DRAIN.....	235
INACTIVE - OR OUTPUT - PUMP LOSS.....	236
INACTIVE - OR OUTPUT - URINE.....	237
INACTIVE - OSTOMY TYPE.....	238
INACTIVE - OXYGEN DELIVERED.....	239
INACTIVE - OXYGEN DELIVERY METHOD.....	240
INACTIVE - OXYGEN EXTRACTION RATIO.....	241
INACTIVE - PA CATHETER INSERTION DATE.....	242
INACTIVE - PA CATHETER SITE.....	243
INACTIVE - PACEMAKER DATE INSERTED.....	244
INACTIVE - PACEMAKER SITES.....	245
INACTIVE - PAIN LOCATION.....	246
INACTIVE - POST OPERATION DAY NUMBER.....	247
INACTIVE - PRESSURE ULCER SITE.....	248
INACTIVE - PRESSURE ULCER VISIBLE STRUCTURE.....	249
INACTIVE - PULSE OXIMETRY.....	250
INACTIVE - SKIN DISCHARGE COLOR.....	251
INACTIVE - SKIN DISCHARGE ODOR.....	252
INACTIVE - SKIN DRAINAGE CHARACTERISTICS.....	253
INACTIVE - VENTILATOR IN-LINE NEBULIZER.....	254
INACTIVE - WOUND SITE.....	255
INACTIVE - WOUND TOPICAL TREATMENT TYPE.....	256
INAPPROPRIATE (data).....	757
INAPPROPRIATE (q-pair).....	30, 124
INCENTIVE SPIROMETRY.....	257
INCH (unit).....	1031
INCISION DISCHARGE CHARACTERISTICS.....	258
INCISION DISCHARGE CHARACTERISTICS BLOODY (data).....	757
INCISION DISCHARGE CHARACTERISTICS CLOUDY (data).....	757
INCISION DISCHARGE CHARACTERISTICS NONE (data).....	757
INCISION DISCHARGE CHARACTERISTICS PURULENT (data).....	757
INCISION DISCHARGE CHARACTERISTICS SEROSANGUINOUS (data)	758
INCISION DISCHARGE CHARACTERISTICS SEROUS (data).....	758
INCISION DISCHARGE CHARACTERISTICS VISCOUS (data).....	758
INCISION DRAINAGE AMOUNT.....	259
INCISION DRAINAGE AMOUNT HEAVY (data).....	758
INCISION DRAINAGE AMOUNT HEAVY (q-pair).....	259
INCISION DRAINAGE AMOUNT MODERATE (data).....	758
INCISION DRAINAGE AMOUNT MODERATE (q-pair).....	259
INCISION DRAINAGE AMOUNT NONE (data).....	759
INCISION DRAINAGE AMOUNT NONE (q-pair).....	259
INCISION DRAINAGE AMOUNT SCANT (data).....	759
INCISION DRAINAGE AMOUNT SCANT (q-pair).....	259
INCISION DRAINAGE AMOUNT SMALL (data).....	759
INCISION DRAINAGE AMOUNT SMALL (q-pair).....	259
INCISION DRESSING PREVENTS ASSESSMENT.....	260
INCISION DRESSING TYPE.....	261
INCISION DRESSING TYPE ANTIMICROBIAL (data).....	759
INCISION DRESSING TYPE ANTIMICROBIAL (q-pair).....	261

INCISION DRESSING TYPE ANTIMICROBIAL DISK (data).....	759
INCISION DRESSING TYPE ANTIMICROBIAL DISK (q-pair).....	261
INCISION DRESSING TYPE BALSAM OF PERU (data).....	760
INCISION DRESSING TYPE BALSAM OF PERU (q-pair).....	261
INCISION DRESSING TYPE BISMUTH GAUZE (data).....	760
INCISION DRESSING TYPE BISMUTH GAUZE (q-pair).....	261
INCISION DRESSING TYPE CALCIUM ALGINATE (data).....	760
INCISION DRESSING TYPE CALCIUM ALGINATE (q-pair).....	261
INCISION DRESSING TYPE FOAM DRESSING (data).....	760
INCISION DRESSING TYPE FOAM DRESSING (q-pair).....	261
INCISION DRESSING TYPE GAUZE (data).....	760
INCISION DRESSING TYPE GAUZE (q-pair).....	261
INCISION DRESSING TYPE GAUZE PACKING (data).....	761
INCISION DRESSING TYPE GAUZE PACKING (q-pair).....	261
INCISION DRESSING TYPE GAUZE WRAP (data).....	761
INCISION DRESSING TYPE GAUZE WRAP (q-pair).....	261
INCISION DRESSING TYPE HYDROCOLLOIDAL (data).....	761
INCISION DRESSING TYPE HYDROCOLLOIDAL (q-pair).....	261
INCISION DRESSING TYPE HYDROGEL (data).....	761
INCISION DRESSING TYPE HYDROGEL (q-pair).....	261
INCISION DRESSING TYPE HYPERTONIC HYDROGEL (data).....	761
INCISION DRESSING TYPE HYPERTONIC HYDROGEL (q-pair).....	261
INCISION DRESSING TYPE NONE (data).....	762
INCISION DRESSING TYPE NONE (q-pair).....	261
INCISION DRESSING TYPE PRESSURE (data).....	762
INCISION DRESSING TYPE PRESSURE (q-pair).....	261
INCISION DRESSING TYPE TRANSPARENT (data).....	762
INCISION DRESSING TYPE TRANSPARENT (q-pair).....	261
INCISION EDGE.....	262
INCISION EDGE INDISTINCT/ATTACHED (data).....	762
INCISION EDGE INDISTINCT/ATTACHED (q-pair).....	262
INCISION EDGE NON-ATTACHED (data).....	762
INCISION EDGE NON-ATTACHED (q-pair).....	262
INCISION EDGE ROLLED (data).....	763
INCISION EDGE ROLLED (q-pair).....	262
INCISION EDGE TUNNELING (data).....	763
INCISION EDGE TUNNELING (q-pair).....	262
INCISION EDGE UNDERMINING (data).....	763
INCISION EDGE UNDERMINING (q-pair).....	262
INCISION ERYTHEMA.....	263
INCISION ODOR.....	264
INCISION ODOR FOUL (data).....	763
INCISION ODOR FOUL (q-pair).....	264
INCISION ODOR NONE (data).....	763
INCISION ODOR NONE (q-pair).....	264
INCISION PRESENT.....	265
INCISION STAPLES INTACT.....	266
INCISION SUTURES INTACT.....	267
INCISION TENDER TO PALPATION.....	268
INCISION WOUND CLOSURE APPROXIMATED (data).....	764
INCISION WOUND CLOSURE APPROXIMATED (q-pair).....	515
INCISION WOUND CLOSURE DEHISCENCE (data).....	764
INCISION WOUND CLOSURE DEHISCENCE (q-pair).....	515
INCISION WOUND CLOSURE GLUE (data).....	764
INCISION WOUND CLOSURE GLUE (q-pair).....	515
INCISION WOUND CLOSURE NOT CLOSED (data).....	764

INCISION WOUND CLOSURE NOT CLOSED (q-pair).....	515
INCISIONAL SITE (data).....	764
INCISIONAL SITE (loca).....	918
INCISIONAL SITE (q-pair).....	141, 384
INCOMPREHENSIBLE (data).....	765
INCOMPREHENSIBLE (q-pair).....	124
INCONTENENCE DERMATITIS (data).....	765
INCONTENENCE DERMATITIS (q-pair).....	627
INCONTINENT (data).....	765
INCONTINENT (q-pair).....	183, 577
INDEPENDENT (data).....	765
INDEPENDENT (q-pair).....	5, 7, 8, 9
INDICATES UNDERSTANDING (data).....	765
INDICATES UNDERSTANDING (q-pair).....	461, 466
INDIFFERENT (data).....	766
INDIFFERENT (q-pair).....	144
INDURATION (data).....	766
INDURATION (q-pair).....	505
INDWELLING CATHETER (data).....	766
INDWELLING CATHETER (q-pair).....	182, 577
INEFFECTIVE COPING (data).....	766
INEFFECTIVE COPING (q-pair).....	52
INFECTON (data).....	766
INFECTON (q-pair).....	145
INFECTON PREVENTION AND SCREENING (data).....	767
INFECTON PREVENTION AND SCREENING (q-pair).....	145
INFLATABLE HEEL LIFT (data).....	767
INFLATABLE HEEL LIFT (q-pair).....	418
INFUSION SET CHANGED.....	269
INITIATED (data).....	767
INITIATED (q-pair).....	487
INSERTED THROUGH EXISTING INTRODUCER (data).....	767
INSERTED THROUGH EXISTING INTRODUCER (q-pair).....	84, 357
INSERTED THROUGH NEW INTRODUCER (data).....	767
INSERTED THROUGH NEW INTRODUCER (q-pair).....	84, 357
INSPIRATORY WHEEZES (data).....	768
INSPIRATORY WHEEZES (q-pair).....	69
INTACT (data).....	768
INTACT (q-pair).....	42, 86, 483, 484, 485, 486, 511
INTAKE.....	270
INTAKE - DIALYSIS.....	271
INTAKE - EF/SUPPLEMENT CALORIES.....	272
INTAKE - ENTERAL FEEDING TUBE.....	273
INTAKE - GI TUBE.....	274
INTAKE - IRRIGATION.....	275
INTAKE - IV BLOOD PRODUCTS.....	276
INTAKE - IV CALORIES.....	277
INTAKE - IV CARDIAC OUTPUT INJECTATE.....	278
INTAKE - IV COLLOID.....	279
INTAKE - IV CRYSTALLOIDS.....	280
INTAKE - IV FLUID.....	281
INTAKE - IV MEDS.....	282
INTAKE - IV NUTRITION.....	283
INTAKE - LINE FLUSH.....	284
INTAKE - PCA VOLUME.....	285

INTAKE - PERCENT MEAL CONSUMED.....	286
INTAKE - PO CALORIES.....	287
INTAKE - PO FLUIDS.....	288
INTAKE - TPN CALORIES.....	289
INTAKE - TRANSFUSION	290
INTENTIONAL (data).....	768
INTENTIONAL (q-pair).....	326
INTERMITTENT (data).....	768
INTERMITTENT (q-pair).....	177, 179, 381, 382, 628
INTERMITTENT SLEEP (data).....	768
INTERMITTENT SLEEP (q-pair).....	529
INTERRUPTED (data).....	769
INTERRUPTED (q-pair).....	487
INTRACRANIAL PRESSURE.....	291
INTRODUCER (data).....	769
INTRODUCER (q-pair).....	92
INTRODUCER CATHETER (meth).....	970
INVASIVE (meth).....	970
INVASIVE (q-pair).....	58, 133, 312, 313, 440, 554
IRREGULAR (data).....	769
IRREGULAR (q-pair).....	447
IRRIGATION (meth).....	971
IRRITABLE (data).....	769
IRRITABLE (q-pair).....	385
ISOLATION PRECAUTIONS.....	292
IV (meth).....	971
IV (q-pair).....	270
IV DISCONTINUED (data).....	769
IV DISCONTINUED (q-pair).....	135
IV DOSE RATE - UNITS/HOUR.....	293
IV DOSE RATE - UNITS/kg/HOUR.....	294
IV DOSE RATE - g/TIME.....	295
IV DOSE RATE - mcg/TIME.....	296
IV DOSE RATE - mcg/kg/HOUR.....	297
IV DOSE RATE - mg/TIME.....	298
IV DOSE RATE - mg/kg/HOUR.....	299
IV FLOW RATE - ml/TIME.....	300
IV FLOW RATE - ml/kg/TIME.....	301
IV LINE CARE.....	302
IV PORT / NUMBER.....	303
IV PUMP (data).....	770
IV PUMP (q-pair).....	23
IV SITE CONDITION.....	304
IVPB (data).....	770

- J -

JAUNDICE (data).....	770
JAUNDICE (q-pair).....	510
JEJUNOSTOMY (data).....	770
JEJUNOSTOMY (q-pair).....	147, 178, 336
JOINT CAPSULE (data).....	770
JP DRAIN (data).....	771
JP DRAIN (meth).....	971

JP DRAIN (prod).....	995
JP DRAIN (q-pair).....	343, 616
JUICE (prod).....	995
JUICE (q-pair).....	270, 288
JUNCTIONAL (data).....	771
JUNCTIONAL (q-pair).....	79

- K -

K-PAD (data).....	771
K-PAD (q-pair).....	605
KG/M2 (unit).....	1031
KILOGRAM (unit).....	1031
KNEE (loca).....	919
KNEE (q-pair).....	105
KNEES (loca).....	919
KNEES (q-pair).....	141
KNIFE-LIKE (data).....	771
KNIFE-LIKE (q-pair).....	381
KUSSMAUL (data).....	771
KUSSMAUL (q-pair).....	447

- L -

L ABDOMEN (loca).....	919
L ABDOMEN (q-pair).....	622
L ANKLE (loca).....	919
L ANKLE (q-pair).....	265, 384, 421, 470, 523, 622
L ANTECUBITAL (loca).....	919
L ANTECUBITAL (q-pair).....	92, 130, 182, 303, 304, 421, 622
L ARM (loca).....	920
L ARM (q-pair).....	58, 92, 130, 133, 141, 146, 182, 303, 304, 384, 468, 469, 523, 554
L AXILLARY (loca).....	920
L AXILLARY (q-pair).....	41, 45, 198, 201
L BRACHIAL (loca).....	920
L BRACHIAL (q-pair).....	41, 45, 95, 368, 373
L BUNDLE BRANCH BLOCK (data).....	772
L BUNDLE BRANCH BLOCK (q-pair).....	79
L CALF (loca).....	920
L CALF (q-pair).....	265
L CAROTID (loca).....	920
L CAROTID (q-pair).....	74, 81
L CHEST (loca).....	921
L CHEST (q-pair).....	98, 384, 523, 622
L DORSALIS PEDIS (loca).....	921
L EAR (loca).....	921
L EAR (q-pair).....	265, 384, 421, 523, 622
L ELBOW (loca).....	921
L ELBOW (q-pair).....	384, 421, 523, 622

L EXTERNAL JUGULAR (loca).....921
L EXTERNAL JUGULAR (q-pair).....92, 130, 182, 303,
304, 361
L EYE (loca).....922
L EYE (q-pair).....122, 141, 442, 443,
602, 622
L FACE (loca).....922
L FACE (q-pair).....141, 265, 622
L FEMORAL (loca).....922
L FEMORAL (q-pair).....41, 45, 74, 81, 92,
95, 130, 182, 198,
201, 265, 303, 304,
361, 368, 373, 421,
441, 622
L FINGER (loca).....922
L FINGER (q-pair).....265
L FLANK (loca).....922
L FLANK (q-pair).....265
L FOOT (loca).....923
L FOOT (q-pair).....141, 160, 265, 303,
304, 384, 421, 483,
523, 622
L FOOT DOWN (data).....772
L FOOT DOWN (q-pair).....499
L FOOT UP (data).....772
L FOOT UP (q-pair).....499
L GLUTEUS (loca).....923
L GLUTEUS (q-pair).....265, 421, 622
L GROIN (loca).....923
L GROIN (q-pair).....384, 523
L HAND (loca).....923
L HAND (q-pair).....141, 185, 384, 484,
523, 622
L HEAD (loca).....923
L HEAD (q-pair).....421, 622
L HEAD DOWN (data).....772
L HEAD DOWN (q-pair).....499
L HEAD UP (data).....772
L HEAD UP (q-pair).....499
L HEEL (loca).....924
L HEEL (q-pair).....421, 622
L HIP (loca).....924
L HIP (q-pair).....265, 384, 523, 622
L INGUINAL (loca).....924
L INGUINAL (q-pair).....265
L INNER KNEE (loca).....924
L INNER KNEE (q-pair).....421, 622
L INTERNAL JUGULAR (loca).....924
L INTERNAL JUGULAR (q-pair).....92, 95, 130, 182,
303, 304, 361, 368,
373
L ISCHIAL (loca).....925
L ISCHIAL (q-pair).....421, 622
L JAW (loca).....925
L JAW (q-pair).....384, 523
L KNEE (loca).....925

L KNEE (q-pair).....	141, 265, 384, 523
L LATERAL HEEL (loca).....	925
L LATERAL HEEL (q-pair).....	421
L LATERAL MALLEOLUS (loca).....	925
L LATERAL MALLEOLUS (q-pair).....	421
L LEG (loca).....	926
L LEG (q-pair).....	58, 125, 133, 141, 146, 303, 304, 384, 468, 470, 523, 554
L LOWER ANTERIOR CHEST (loca).....	926
L LOWER ANTERIOR CHEST (q-pair).....	100, 265
L LOWER ARM (loca).....	926
L LOWER ARM (q-pair).....	265, 384, 421, 523, 622
L LOWER EXTREMITY (loca).....	926
L LOWER EXTREMITY (q-pair).....	325, 485
L LOWER LEG (loca).....	926
L LOWER LEG (q-pair).....	265, 622
L LOWER LIMB (loca).....	927
L LOWER LIMB (q-pair).....	310
L LOWER POSTERIOR CHEST (loca).....	927
L LOWER POSTERIOR CHEST (q-pair).....	100, 265
L MANDIBULAR (loca).....	927
L MANDIBULAR (q-pair).....	265
L MEDIAL HEEL (loca).....	927
L MEDIAL HEEL (q-pair).....	421
L MEDIAL MALLEOLUS (loca).....	927
L MEDIAL MALLEOLUS (q-pair).....	421
L MIDDLE ANTERIOR CHEST (loca).....	928
L MIDDLE ANTERIOR CHEST (q-pair).....	100, 265
L MIDDLE POSTERIOR CHEST (loca).....	928
L MIDDLE POSTERIOR CHEST (posi).....	985
L MIDDLE POSTERIOR CHEST (q-pair).....	100, 265
L NECK (data).....	773
L NECK (loca).....	928
L NECK (q-pair).....	265
L NEPHROSTOMY (data).....	773
L NEPHROSTOMY (q-pair).....	182, 577
L NOSE (loca).....	928
L NOSE (q-pair).....	384, 523
L OCCIPITAL (loca).....	928
L OCCIPITAL (q-pair).....	265
L ORAL (loca).....	929
L ORAL (q-pair).....	265
L ORBIT (loca).....	929
L ORBIT (q-pair).....	265
L OUTER KNEE (loca).....	929
L OUTER KNEE (q-pair).....	421, 622
L PARIETAL (loca).....	929
L PARIETAL (q-pair).....	265
L PEDAL (loca).....	929
L PEDAL (q-pair).....	441
L POST TIBIAL (loca).....	930
L POST TIBIAL (q-pair).....	441
L RADIAL (loca).....	930
L RADIAL (q-pair).....	41, 45, 441

L RIBS (loca).....	930
L RIBS (q-pair).....	421, 622
L SCAPULA (loca).....	930
L SCAPULA (q-pair).....	384, 421, 523, 622
L SHOULDER (loca).....	930
L SHOULDER (q-pair).....	384, 421, 523, 622
L STENT (data).....	773
L STENT (q-pair).....	182
L SUBCLAVIAN (loca).....	931
L SUBCLAVIAN (q-pair).....	92, 95, 130, 182, 265, 303, 304, 361, 368, 373
L TEMPORAL (loca).....	931
L TEMPORAL (q-pair).....	74, 81
L TESTICLE (loca).....	931
L TESTICLE (q-pair).....	622
L THIGH (loca).....	931
L THIGH (q-pair).....	265, 384, 523
L THORACIC ANTERIOR (loca).....	931
L THORACIC ANTERIOR (q-pair).....	265
L THORACIC LATERAL (loca).....	932
L THORACIC LATERAL (q-pair).....	265
L THORACIC POSTERIOR (loca).....	932
L THORACIC POSTERIOR (q-pair).....	265
L THORACOTOMY (loca).....	932
L THORACOTOMY (q-pair).....	265
L THUMB (loca).....	932
L THUMB (q-pair).....	265
L TOE (loca).....	932
L TOE (q-pair).....	265
L TRAPEZIUS (loca).....	933
L TRAPEZIUS (q-pair).....	265
L TROCHANTER (loca).....	933
L TROCHANTER (q-pair).....	421, 622
L UPPER ANTERIOR CHEST (loca).....	933
L UPPER ANTERIOR CHEST (q-pair).....	100, 265
L UPPER ARM (loca).....	933
L UPPER ARM (q-pair).....	92, 130, 182, 265, 303, 304, 384, 421, 523, 622
L UPPER EXTREMITY (loca).....	933
L UPPER EXTREMITY (q-pair).....	325, 486
L UPPER LEG (loca).....	934
L UPPER LEG (q-pair).....	622
L UPPER LIMB (loca).....	934
L UPPER LIMB (q-pair).....	573
L UPPER POSTERIOR CHEST (loca).....	934
L UPPER POSTERIOR CHEST (q-pair).....	100, 265
L WRIST (loca).....	934
L WRIST (q-pair).....	92, 130, 182, 303, 304, 468, 469, 622
LACERATION (data).....	773
LACERATION (q-pair).....	627
LACKS INTEREST (data).....	773
LACKS INTEREST (q-pair).....	144
LACTOSE RESTRICTED (data).....	774

LACTOSE RESTRICTED (q-pair)	134
LANGUAGE (data)	774
LANGUAGE (q-pair)	142
LARGE (data)	774
LARGE (q-pair)	341, 347, 349, 432
LARYNGEAL MASK (meth)	971
LARYNGEAL MASK (q-pair)	550
LARYNGECTOMY (data)	774
LARYNGECTOMY (q-pair)	21
LEAD (meth)	971
LEAD (q-pair)	440
LEFT (data)	774
LEFT (loca)	934
LEFT (q-pair)	12
LEFT LOWER LOBE (loca)	935
LEFT LOWER LOBE (q-pair)	69
LEFT NECK (loca)	935
LEFT NECK (q-pair)	141
LEFT UPPER LOBE (loca)	935
LEFT UPPER LOBE (q-pair)	69
LEFT VENTRICULAR EJECTION FRACTION - CATEGORICAL	306
LEFT VENTRICULAR EJECTION FRACTION - CONTINUOUS	305
LEFT VENTRICULAR STROKE WORK	307
LEG MOVEMENTS	308
LESS THAN 30 (data)	775
LESS THAN 30 (q-pair)	306
LETHARGIC (data)	775
LETHARGIC (q-pair)	52, 309
LEVEL OF CONSCIOUSNESS	309
LEVELED (data)	775
LEVELED (q-pair)	302
LG ADULT (meth)	972
LG ADULT CUFF (meth)	972
LIBERAL DIABETIC (data)	775
LIBERAL DIABETIC (q-pair)	134
LIFT SCALE (meth)	972
LIFT SCALE (q-pair)	607
LIMITED ASSISTANCE (data)	775
LIMITED ASSISTANCE (q-pair)	7, 8, 9
LIMITED ATTENTION SPAN (data)	776
LIMITED ATTENTION SPAN (q-pair)	142
LIMITED TO INCISION (data)	776
LIMITED TO INCISION (q-pair)	263, 513
LINE FLUSH (data)	776
LINE FLUSH (prod)	996
LINE FLUSH HEPARIN (prod)	996
LINE FLUSH HEPARIN (q-pair)	284
LIP SALVE (data)	776
LIP SALVE (q-pair)	14
LIPIDS (prod)	996
LIPIDS (q-pair)	283
LISTED DIET (data)	776
LISTED DIET (q-pair)	134
LITER (unit)	1031
LITERS PER MINUTE (unit)	1031
LITERS PER MINUTE PER METER2 (BSA) (unit)	1032

LOCALIZES TO PAIN (data).....	777
LOCALIZES TO PAIN (q-pair).....	123
LOOSE (data).....	777
LOOSE (q-pair).....	544
LOW (data).....	777
LOW (q-pair).....	49, 177
LOW AIRLOSS OVERLAY (data).....	777
LOW AIRLOSS OVERLAY (q-pair).....	530
LOW AIRLOSS/ALTERNATING PRESSURE COMBO (data).....	777
LOW AIRLOSS/ALTERNATING PRESSURE COMBO (q-pair).....	530
LOW CHOLESTEROL (data).....	778
LOW CHOLESTEROL (q-pair).....	134
LOW FIBER (data).....	778
LOW FIBER (q-pair).....	134
LOW PROTEIN 60 GM (data).....	778
LOW PROTEIN 60 GM (q-pair).....	134
LOW TYRAMINE (data).....	778
LOW TYRAMINE (q-pair).....	134
LOWER ABDOMEN (loca).....	935
LOWER ABDOMEN (q-pair).....	622
LOWER ARM (loca).....	935
LOWER ARM (q-pair).....	105
LOWER BACK (loca).....	936
LOWER BACK (q-pair).....	622
LOWER LIMB EXTREMITIES STRENGTH.....	310
LUMBAR DRAIN (data).....	778
LUMBAR DRAIN (q-pair).....	616
LUMBAR SPINAL PRECAUTIONS (data).....	779
LUMBAR SPINAL PRECAUTIONS (q-pair).....	11
LVSWI.....	311
LYING (posi).....	985
LYING (q-pair).....	58, 133, 440, 448, 554, 607
LYING OR SITTING (data).....	779
LYING OR SITTING (q-pair).....	117

- M -

MACERATION (data).....	779
MACERATION (q-pair).....	505
MAGNACAL (data).....	779
MAGNACAL (prod).....	996
MAGNACAL (q-pair).....	149
MAGNETIC RESONANCE IMAGING (meth).....	972
MAGNETIC RESONANCE IMAGING (q-pair).....	305, 306
MAINTAINING EFFECTIVE VENTILATION (data).....	779
MAINTAINING EFFECTIVE VENTILATION (q-pair).....	377
MAROON (data).....	780
MAROON (q-pair).....	545, 614
MASK (meth).....	972
MASK (q-pair).....	550, 551
MEAL (data).....	780
MEAL (q-pair).....	4
MEAN ARTERIAL PRESSURE.....	312
MEAN BLOOD PRESSURE.....	313
MEAN BLOOD PRESSURE ARTERIAL LINE.....	314

MEAN BLOOD PRESSURE NON-INVASIVE.....	315
MEAN BP CUFF.....	316
MECHANICAL (data).....	780
MECHANICAL (q-pair).....	134
MED RECONCILIATION DOCUMENTED (data).....	780
MED RECONCILIATION DOCUMENTED (q-pair).....	135
MEDIAL PORT (data).....	780
MEDIAL PORT (q-pair).....	303
MEDIASTINAL (loca).....	936
MEDIASTINAL (prod).....	996
MEDIASTINAL (q-pair).....	198, 201, 337
MEDICAL DEVICE SAFETY (data).....	781
MEDICAL DEVICE SAFETY (q-pair).....	54, 119
MEDICATION (data).....	781
MEDICATION (prod).....	997
MEDICATION (q-pair).....	145, 386
MEDICATION ROOM (data).....	781
MEDICATION ROOM (q-pair).....	51, 118
MEDIUM (data).....	781
MEMORY (data).....	781
MEMORY (q-pair).....	142
MENTAL HEALTH (data).....	782
MENTAL HEALTH (q-pair).....	145
MENTAL HEALTH HOLD (data).....	782
MENTAL HEALTH HOLD (q-pair).....	145
MICROGRAMS (unit).....	1032
MICROGRAMS PER HOUR (unit).....	1032
MICROGRAMS PER KILOGRAM PER HOUR (unit).....	1032
MICROGRAMS PER KILOGRAM PER MINUTE (unit).....	1032
MICROGRAMS PER MINUTE (unit).....	1033
MICU (data).....	782
MICU (q-pair).....	379
MID BACK (loca).....	936
MID BACK (q-pair).....	622
MID-STERNAL (loca).....	936
MID-STERNAL (q-pair).....	98
MIDDLE (loca).....	936
MIDLINE ABDOMINAL (loca).....	937
MIDLINE ABDOMINAL (q-pair).....	265
MIDLINE NECK (loca).....	937
MIDLINE NECK (q-pair).....	265
MILD-MODERATELY REDUCED (data).....	782
MILD-MODERATELY REDUCED (q-pair).....	306
MILLIAMPS (unit).....	1033
MILLIBAR (unit).....	1033
MILLIGRAMS (unit).....	1033
MILLIGRAMS PER HOUR (unit).....	1033
MILLIGRAMS PER KILOGRAM PER HOUR (data).....	782
MILLIGRAMS PER KILOGRAM PER HOUR (unit).....	1034
MILLIGRAMS PER KILOGRAM PER MINUTE (unit).....	1034
MILLIGRAMS PER MINUTE (unit).....	1034
MILLILITER (unit).....	1034
MILLILITERS PER BEAT (unit).....	1034
MILLILITERS PER BEAT PER METER SQUARE (BSA) (unit).....	1035
MILLILITERS PER HOUR (unit).....	1035
MILLILITERS PER METER SQUARE (BSA) (unit).....	1035

MILLILITERS PER MINUTE (unit).....	1035
MILLILITERS PER MINUTE PER METER SQUARE (unit).....	1035
MILLIMETERS (unit).....	1036
MILLIMETERS OF MERCURY (unit).....	1036
MILLIVOLTS (unit).....	1036
MILLILITERS PER KILOGRAM PER HOUR (unit).....	1036
MILLILITERS PER KILOGRAM PER MINUTE (unit).....	1036
MINUTE (unit).....	1037
MISSING (data).....	783
MOANING (data).....	783
MOANING (q-pair).....	385
MODERATE (data).....	783
MODERATE (q-pair).....	341, 347, 349, 432, 512
MODERATE-SEVERELY REDUCED (data).....	783
MODERATE-SEVERELY REDUCED (q-pair).....	306
MODERATELY REDUCED (data).....	783
MODERATELY REDUCED (q-pair).....	306
MODIFIED ALDRETE SCORE (unit).....	1037
MODULAR CUSTOMIZED FORMULA (prod).....	997
MODULAR CUSTOMIZED FORMULA (q-pair).....	149
MOIST (data).....	784
MOIST (q-pair).....	518
MONITOR (meth).....	973
MONITOR (q-pair).....	58, 77, 78, 95, 440, 448, 556
MONITORING CAMERA.....	317
MONTH (unit).....	1037
MORSE SCORE (unit).....	1037
MORSE: AMBULATORY AID DURING WALKING.....	318
MORSE: GAIT.....	319
MORSE: HISTORY OF FALLING WITHIN LAST THREE MONTHS.....	320
MORSE: IV THERAPY SALINE LOCK.....	321
MORSE: MENTAL STATUS OF PT OWN ABILITY.....	322
MORSE: SCORE TOTAL.....	323
MORSE: SECONDARY DIAGNOSIS LISTED ON CHART.....	324
MOTOR ACTIVITY.....	325
MOTOR CHECK.....	326
MOTORIZED WHEELCHAIR (data).....	784
MOTORIZED WHEELCHAIR (q-pair).....	568
MOTTLED (data).....	784
MOTTLED (q-pair).....	510
MOUTHWASH USED (data).....	784
MOUTHWASH USED (q-pair).....	14
MOVE AGAINST GRAVITY AND RESISTANCE, EASILY OVERCOME BY THE EXAMINER (data)	784
MOVE AGAINST GRAVITY AND RESISTANCE, EASILY OVERCOME BY THE EXAMINER (q-pai r).....	325
MOVE AGAINST GRAVITY BUT NO RESISTANCE (data).....	785
MOVE AGAINST GRAVITY BUT NO RESISTANCE (q-pair).....	325
MOVEMENT WITH GRAVITY ELIMINATED (data).....	785
MOVEMENT WITH GRAVITY ELIMINATED (q-pair).....	325
MUCOID (data).....	785
MUCOID (q-pair).....	544
MUCOUS MEMBRANE CONDITION.....	327
MULTIGATED NUCLEAR VENTRICULOGRAPHY (meth).....	973

MULTIGATED NUCLEAR VENTRICULOGRAPHY (q-pair).....305, 306
MULTIPODUS (data).....785
MULTIPODUS (q-pair).....418
MULTIVOLTS (unit).....1037
MURMUR (data).....785
MURMUR (q-pair).....76
MUSHROOM CATHETER (data).....786
MUSHROOM CATHETER (q-pair).....182
MUSTACHE CANNULA (meth).....973
mg/dl (unit).....1041

- N -

N/A (data).....786
N/A (q-pair).....38, 160, 185, 257,
546, 576
NAIL CARE (data).....786
NAIL CARE (q-pair).....4
NAIL COLOR.....328
NAIL SHAPE.....329
NAIL TEXTURE.....330
NAPS ONLY (data).....786
NAPS ONLY (q-pair).....529
NASAL (data).....786
NASAL (q-pair).....437
NASAL CANNULA (meth).....973
NASAL CANNULA (q-pair).....551
NASAL ENDOTRACHEAL (data).....787
NASAL ENDOTRACHEAL (q-pair).....21
NASAL INTESTINAL SMALL BORE (data).....787
NASAL INTESTINAL SMALL BORE (q-pair).....147
NASAL SMALL BOWEL FEEDING TUBE (data).....787
NASAL SMALL BOWEL FEEDING TUBE (q-pair).....178
NASAL TRUMPET (data).....787
NASAL TRUMPET (q-pair).....21
NASOGASTRIC LARGE BORE TUBE (data).....787
NASOGASTRIC LARGE BORE TUBE (q-pair).....147
NASOGASTRIC TUBE (data).....788
NASOGASTRIC TUBE (q-pair).....178
NASOGASTRIC TUBE SMALL BORE (data).....788
NAUSEA (data).....788
NAUSEA (q-pair).....172
NEAT (data).....788
NEAT (q-pair).....30
NEBULIZER (meth).....973
NEBULIZER (q-pair).....551
NECK (loca).....937
NECK (q-pair).....105, 384, 523, 622
NECK VEIN CONDITION.....331
NECROTIC (data).....788
NECROTIC (q-pair).....621
NEED FOR FLUIDS/TOILETING ASSESSED.....332
NEEDS REINFORCEMENT (data).....789
NEEDS REINFORCEMENT (q-pair).....461, 466
NEGATIVE PRESSURE (data).....789
NEGATIVE PRESSURE (q-pair).....292

NEGATIVE PRESSURE THERAPY DEVICE (data).....	789
NEGATIVE PRESSURE THERAPY DEVICE (q-pair).....	422
NEPRO RENAL FORMULA (data).....	789
NEPRO RENAL FORMULA (prod).....	997
NEPRO RENAL FORMULA (q-pair).....	149, 270
NEUTROPENIC (data).....	789
NEUTROPENIC (q-pair).....	134, 292
NEW ADMISSION EDUCATION (data).....	790
NEW ADMISSION EDUCATION (q-pair).....	145
NEW SITE (data).....	790
NEW SITE (q-pair).....	84, 357
NG RESIDUAL (prod).....	997
NG RESIDUAL (q-pair).....	337
NO (data).....	790
NO (q-pair).....	13, 15, 17, 22, 31, 32, 33, 34, 35, 36, 37, 39, 40, 44, 47, 48, 50, 57, 70, 71, 73, 83, 89, 94, 97, 101, 102, 128, 154, 155, 164, 171, 180, 188, 190, 195, 196, 203, 257, 260, 265, 266, 267, 268, 269, 332, 352, 359, 363, 366, 368, 370, 378, 380, 393, 403, 408, 409, 414, 417, 420, 423, 425, 431, 446, 453, 454, 455, 457, 458, 459, 460, 474, 481, 482, 502, 504, 506, 507, 508, 519, 520, 562, 570, 571, 579, 618, 625
NO BREAK IN EPIDERMIS (data).....	790
NO BREAK IN EPIDERMIS (q-pair).....	421
NO MOVEMENT (data).....	790
NO MOVEMENT (q-pair).....	325
NO SIGNS AND SYMPTOMS (data).....	791
NO SIGNS AND SYMPTOMS (q-pair).....	385
NON RE-BREATHER (meth).....	974
NON RE-BREATHER (q-pair).....	551
NON-ADHERENT GAUZE (data).....	791
NON-INVASIVE (meth).....	974
NON-INVASIVE (q-pair).....	58, 133, 312, 313, 440, 554
NON-PITTING (data).....	791
NON-PRODUCTIVE (data).....	791
NON-PRODUCTIVE (q-pair).....	126
NON-TENDER (data).....	791
NON-TENDER (q-pair).....	2
NONBLANCHABLE ERYTHEMA (data).....	792
NONBLANCHABLE ERYTHEMA (q-pair).....	421
NONE (data).....	792
NONE (q-pair).....	75, 122, 123, 124,

	142, 258, 263, 341, 399, 413, 416, 419, 509, 512, 615, 619
NONINVASIVE (meth).....	974
NONREACTIVE (data).....	792
NONREACTIVE (q-pair).....	442
NORMAL (data).....	792
NORMAL (q-pair).....	1, 46, 93, 103, 184, 306, 329, 330, 331, 355, 362, 518
NORMAL FOR ETHNIC GROUP (data).....	792
NORMAL FOR ETHNIC GROUP (q-pair).....	328, 510
NORMAL POWER (data).....	793
NORMAL POWER (q-pair).....	325
NORMAL SINUS (data).....	793
NORMAL STOOL (data).....	793
NORMAL STOOL (q-pair).....	544
NOSE (loca).....	937
NOSE (q-pair).....	384, 523
NOSE SIZE.....	333
NOT DONE (data).....	793
NOT DONE (q-pair).....	7, 8, 9, 14
NOT MONITORED (data).....	793
NOT MONITORED (q-pair).....	46, 93, 362
NOVASOURCE 2.0 (data).....	794
NOVASOURCE 2.0 (prod).....	997
NOVASOURCE 2.0 (q-pair).....	270
NPO DIET (data).....	794
NPO DIET (q-pair).....	134
NUCLEAR PERFUSION IMAGING (meth).....	974
NUCLEAR PERFUSION IMAGING (q-pair).....	305, 306
NUMB (data).....	794
NUMB (q-pair).....	381
NUMBNESS (data).....	794
NUMBNESS (q-pair).....	483, 484, 485, 486
NURSING INTERVENTIONS - STRAIGHT CATH (task).....	1019
NURSING INTERVENTIONS - WOUND CARE (task).....	1019
NURSING STATION (data).....	794
NURSING STATION (q-pair).....	51, 118
NUTREN 1.5 (prod).....	998
NUTREN 1.5 (q-pair).....	149
NUTRITION (data).....	795
NUTRITION (q-pair).....	145
NUTRITIONAL SUPPLEMENT (data).....	795
NUTRITIONAL SUPPLEMENT (prod).....	998
NUTRITIONAL SUPPLEMENT (q-pair).....	288

- O -

O2 SATS (meth).....	974
OBEYS COMMANDS (data).....	795
OBEYS COMMANDS (q-pair).....	123
OCCIPITAL (loca).....	937
OCCIPITAL (q-pair).....	265
OFF DURING BATHING (data).....	795
OFF DURING BATHING (q-pair).....	487

OFF WHILE AMBULATORY (data)	795
OFF WHILE AMBULATORY (q-pair)	487
OFFER PRN MEDS (data)	796
OFFER PRN MEDS (q-pair)	456
OFFERED (data)	796
OFFERED (q-pair)	14
ON WHILE SITTING IN CHAIR (data)	796
ON WHILE SITTING IN CHAIR (q-pair)	487
ONE PERSON ASSIST (data)	796
ONE PERSON ASSIST (q-pair)	5, 7, 8, 9
ONE TO ONE SUPERVISION (data)	796
ONE TO ONE SUPERVISION (q-pair)	456
OPEN SPONTANEOUSLY (data)	797
OPEN SPONTANEOUSLY (q-pair)	122
OPEN TO AIR (data)	797
OPEN TO AIR (q-pair)	422, 626
OPERATING ROOM (data)	797
OPERATING ROOM (q-pair)	379
ORAL (data)	797
ORAL (loca)	938
ORAL (q-pair)	437, 558, 622
ORAL AIRWAY (data)	797
ORAL AIRWAY (q-pair)	21
ORAL ENDOTRACHEAL (data)	798
ORAL ENDOTRACHEAL (q-pair)	21
ORAL GASTRIC TUBE (data)	798
ORAL GASTRIC TUBE (q-pair)	147, 178
ORANGE (data)	798
ORANGE (q-pair)	138, 174, 575, 614
ORIENTATION	334
ORIENTATION (data)	798
ORIENTATION (q-pair)	145
ORIENTED (data)	798
ORIENTED (q-pair)	124, 464
ORIENTED TO PERSON (data)	799
ORIENTED TO PERSON (q-pair)	334
ORIENTED TO PLACE (data)	799
ORIENTED TO PLACE (q-pair)	334
ORIENTED TO SITUATION (data)	799
ORIENTED TO SITUATION (q-pair)	334
ORIENTED TO TIME (data)	799
ORIENTED TO TIME (q-pair)	334
ORTHOPNEIC (data)	799
ORTHOPNEIC (q-pair)	447
ORTHOTIC BOOTS (data)	800
ORTHOTIC BOOTS (q-pair)	335
ORTHOTIC COLLAR (data)	800
ORTHOTIC COLLAR (q-pair)	335
ORTHOTIC KNEE IMMOBILIZER (data)	800
ORTHOTIC KNEE IMMOBILIZER (q-pair)	335
ORTHOTIC SPINE BRACE (data)	800
ORTHOTIC SPINE BRACE (q-pair)	335
ORTHOTIC TED HOSE (data)	800
ORTHOTIC TED HOSE (q-pair)	335
ORTHOTIC WEDGE PILLOW (data)	801
ORTHOTIC WEDGE PILLOW (q-pair)	335

ORTHOTICS.....	335
OSMOLITE (data).....	801
OSMOLITE (prod).....	998
OSMOLITE (q-pair).....	149
OSMOLITE HN (data).....	801
OSMOLITE HN (prod).....	998
OSMOLITE HN (q-pair).....	149
OSTOMY TYPE.....	336
OTHER (data).....	801
OTHER (meth).....	975
OTHER (q-pair).....	151, 152, 275, 305, 306, 343
OTHER (qual).....	1006
OTHER PO FLUID (prod).....	998
OTHER PO FLUID (q-pair).....	288
OUNCE (unit).....	1038
OUTPUT.....	337
OUTPUT - BLADDER SCAN.....	338
OUTPUT - BLOOD.....	339
OUTPUT - CHEST TUBE.....	340
OUTPUT - CHEST TUBE ESTIMATED.....	341
OUTPUT - DIALYSIS.....	342
OUTPUT - DRAIN.....	343
OUTPUT - EMESIS.....	344
OUTPUT - GI TUBE.....	345
OUTPUT - STOOL.....	346
OUTPUT - STOOL ESTIMATED.....	347
OUTPUT - URINE.....	348
OUTPUT - URINE ESTIMATED.....	349
OVER WIRE (data).....	801
OVER WIRE (q-pair).....	84, 357
OXIMETRIC (data).....	802
OXIMETRIC (q-pair).....	361
OXYGEN CONSUMPTION.....	350
OXYGEN SATURATION.....	351

- P -

P.O. (meth).....	975
P.O. (q-pair).....	270
PA CATHETER (meth).....	975
PA CATHETER ALARMS ON.....	352
PA CATHETER CONDITION.....	353
PA CATHETER DRESSING TYPE.....	354
PA CATHETER DYNAMIC RESPONSE TEST.....	355
PA CATHETER INSERTION MARK.....	356
PA CATHETER INSERTION TECHNIQUE.....	357
PA CATHETER REPOSITIONED/DISCONTINUED REASON.....	358
PA CATHETER SITE CARE.....	359
PA CATHETER SIZE.....	360
PA CATHETER TYPE.....	361
PA CATHETER WAVEFORM.....	362
PA CATHETER ZEROED.....	363
PA LINE (meth).....	975
PAC (data).....	802
PAC (q-pair).....	75

PACED: ATRIAL (data).....	802
PACED: ATRIAL (q-pair).....	79
PACED: AV SEQUENTIAL (data).....	802
PACED: AV SEQUENTIAL (q-pair).....	79
PACED: VENTRICULAR (data).....	802
PACED: VENTRICULAR (q-pair).....	79
PACEMAKER ATRIAL mA.....	364
PACEMAKER AV INTERVAL.....	365
PACEMAKER CAPTURE.....	366
PACEMAKER MODE.....	367
PACEMAKER ON.....	368
PACEMAKER RATE.....	369
PACEMAKER SENSE.....	370
PACEMAKER SENSITIVITY ATRIAL.....	371
PACEMAKER SENSITIVITY VENTRICULAR.....	372
PACEMAKER TYPE.....	373
PACEMAKER UNIT ID.....	374
PACEMAKER VENTRICULAR mA.....	375
PACEMAKER WIRES.....	376
PACER (data).....	803
PACER (q-pair).....	361
PACU - DISCHARGE CRITERIA.....	377
PACU - DISCHARGE CRITERIA MET.....	378
PACU - DISPOSITION.....	379
PACU - PT RECEIVED REVERSAL AGENT.....	380
PAIN (data).....	803
PAIN (q-pair).....	142, 145
PAIN CHARACTER.....	381
PAIN CONTROLLED PER PT (data).....	803
PAIN CONTROLLED PER PT (q-pair).....	377
PAIN DURATION.....	382
PAIN GOAL.....	383
PAIN LEVEL.....	384
PAIN NON-VERBAL CUES.....	385
PAIN RELIEVING FACTOR INTERVENTION.....	386
PAINAD BL 0 (data).....	803
PAINAD BL 0 (q-pair).....	387
PAINAD BL 1 (data).....	803
PAINAD BL 1 (q-pair).....	387
PAINAD BL 2 (data).....	804
PAINAD BL 2 (q-pair).....	387
PAINAD BODY LANGUAGE.....	387
PAINAD BR 0 (data).....	804
PAINAD BR 0 (q-pair).....	388
PAINAD BR 1 (data).....	804
PAINAD BR 1 (q-pair).....	388
PAINAD BR 2 (data).....	804
PAINAD BR 2 (q-pair).....	388
PAINAD BREATHING.....	388
PAINAD CO 0 (data).....	804
PAINAD CO 0 (q-pair).....	389
PAINAD CO 1 (data).....	805
PAINAD CO 1 (q-pair).....	389
PAINAD CO 2 (data).....	805
PAINAD CO 2 (q-pair).....	389
PAINAD CONSOLABILITY.....	389

PAINAD FACIAL EXPRESSION.....	390
PAINAD FE 0 (data).....	805
PAINAD FE 0 (q-pair).....	390
PAINAD FE 1 (data).....	805
PAINAD FE 1 (q-pair).....	390
PAINAD FE 2 (data).....	805
PAINAD FE 2 (q-pair).....	390
PAINAD NEGATIVE VOCALIZATION.....	391
PAINAD NV 0 (data).....	806
PAINAD NV 0 (q-pair).....	391
PAINAD NV 1 (data).....	806
PAINAD NV 1 (q-pair).....	391
PAINAD NV 2 (data).....	806
PAINAD NV 2 (q-pair).....	391
PAINAD SCALE (unit).....	1038
PAINAD TOTAL.....	392
PAINFUL EXPRESSION (data).....	806
PAINFUL EXPRESSION (q-pair).....	385
PALE (data).....	806
PALE (q-pair).....	327, 328, 510
PALPATED (meth).....	975
PALPATED (q-pair).....	58, 440
PAROXYSMAL ATRIAL TACHYCARDIA (data).....	807
PAROXYSMAL ATRIAL TACHYCARDIA (q-pair).....	79
PART THICKNESS SKIN LOSS INVOLVE EPIDERMIS/DERMIS (data)	807
PART THICKNESS SKIN LOSS INVOLVE EPIDERMIS/DERMIS (q-pair)	421
PARTIAL BATH (data).....	807
PARTIAL BATH (q-pair).....	4
PARTIAL IMPAIRED (data).....	807
PARTIAL IMPAIRED (q-pair).....	478
PARTIAL RE-BREATHING (meth).....	976
PARTIAL RE-BREATHING (q-pair).....	551
PARTIAL SIMPLE (data).....	807
PARTIAL SIMPLE (q-pair).....	478
PASSIVE RANGE OF MOTION (data).....	808
PASSIVE RANGE OF MOTION (q-pair).....	11
PATENT CLEAN (data).....	808
PATENT CLEAN (q-pair).....	41, 85, 130, 197, 302, 304, 353, 400
PATIENT (data).....	808
PATIENT (q-pair).....	10
PATIENT ABLE TO RESPOND.....	393
PATIENT INFORMED (data).....	808
PATIENT INFORMED (q-pair).....	317, 450
PATIENT ON PASS (null).....	984
PATIENT REFUSED (null).....	984
PATIENT REPORTS NOT PRESENT (data).....	808
PATIENT REPORTS NOT PRESENT (q-pair).....	158
PATIENT REPORTS PRESENT (data).....	809
PATIENT REPORTS PRESENT (q-pair).....	158
PATIENT SIGNS SAFETY CONTRACT (data).....	809
PATIENT SIGNS SAFETY CONTRACT (q-pair).....	53
PATIENT STAYS WITHIN LINE OF SITE OF NURSING STATION (data)	809

PATIENT STAYS WITHIN LINE OF SITE OF NURSING STATION (q-pair)

.....	53
PATIENT'S ROOM (data).....	809
PATIENT'S ROOM (q-pair).....	51, 118
PCA (data).....	809
PCA (prod).....	999
PCA (q-pair).....	270
PCA - CLINICAL BOLUS DOSE.....	394
PCA - DEMAND DOSE.....	395
PCA - DOSAGE DELIVERED.....	396
PCA - DOSES DELIVERED.....	397
PCA - LIMIT DOSE.....	398
PCA - LIMIT TIME.....	399
PCA - LINE CARE.....	400
PCA - LOCK OUT.....	401
PCA - MODE.....	402
PCA - PUMP CLEARED.....	403
PCA - REQUESTS.....	404
PCA - VOLUME REMAINING.....	405
PCA - VOLUME START.....	406
PCA AND BASAL RATE (data).....	810
PCA AND BASAL RATE (q-pair).....	402
PCA ONLY (data).....	810
PCA ONLY (q-pair).....	402
PEAK INSPIRATORY PRESSURE.....	407
PEDIATRIC (data).....	810
PEDIATRIC (q-pair).....	59
PEDIATRIC (qual).....	1006
PEDIATRIC CUFF (meth).....	976
PEG TUBE (data).....	810
PEG TUBE (q-pair).....	147, 178
PENDANT CANNULA (meth).....	976
PENDANT CANNULA (q-pair).....	551
PENETRATING (data).....	810
PENETRATING (q-pair).....	627
PENIS (loca).....	938
PENIS (q-pair).....	384, 523
PENROSE (meth).....	976
PENROSE (q-pair).....	343
PENROSE DRAIN (data).....	811
PENROSE DRAIN (q-pair).....	616
PERCENTAGE (unit).....	1038
PERCUSSION BED (data).....	811
PERCUSSION BED (q-pair).....	530
PERINEAL (loca).....	938
PERINEAL (q-pair).....	265, 384, 421, 622
PERINEAL CARE (data).....	811
PERINEAL CARE (q-pair).....	4
PERIORBITAL (loca).....	938
PERIORBITAL (q-pair).....	141
PERIPHERAL (loca).....	938
PERIPHERAL (q-pair).....	440
PERIRECTAL (loca).....	939
PERITONEAL DIALYSIS (meth).....	976
PERITONEAL DIALYSIS (q-pair).....	271, 342
PERMANENT (data).....	811

PERMANENT (q-pair).....	373
PERSONAL MEDICAL ALERT BRACELET (data).....	811
PERSONAL MEDICAL ALERT BRACELET (q-pair).....	38
PETROLATUM (data).....	812
PETROLATUM (q-pair).....	422
PHYSICIAN (data).....	812
PHYSICIAN (q-pair).....	567
PHYSICIAN NOTIFIED (data).....	812
PHYSICIAN NOTIFIED (q-pair).....	386
PHYSICIAN ORDER (data).....	812
PHYSICIAN ORDER (q-pair).....	462
PICC (data).....	812
PICC (q-pair).....	92, 302
PICC LINE (meth).....	977
PIN (data).....	813
PIN TRACTION (data).....	813
PIN TRACTION (q-pair).....	563
PINK (data).....	813
PINK (q-pair).....	138, 174, 327, 433, 510, 575, 610, 614
PINPOINT (data).....	813
PINPOINT (q-pair).....	442
PITTED (data).....	813
PITTED (q-pair).....	330
PITTING (data).....	814
PJC (data).....	814
PJC (q-pair).....	75
PLAIN (data).....	814
PLAIN (q-pair).....	361
PLAN OF CARE - REVISION INDICATED.....	408
PLANTAR FLEXION.....	409
PLASMAPHERESIS (prod).....	999
PLASMAPHERESIS (q-pair).....	276
PLATELETS (prod).....	999
PLATELETS (q-pair).....	270, 276
PLATELETS PHERESIS (prod).....	999
PLATELETS PHERESIS (q-pair).....	270, 276
PLATELETS PHERESIS IRRADIATED (prod).....	999
PLATELETS PHERESIS IRRADIATED (q-pair).....	270, 276
POINTS (unit).....	1038
POOR (data).....	814
POOR (q-pair).....	526
POPLITEAL (loca).....	939
POPLITEAL (q-pair).....	440
PORCINE SUBMUCOSA (data).....	814
PORCINE SUBMUCOSA (q-pair).....	422
POSITIONAL (data).....	815
POSITIONAL (q-pair).....	46, 93, 362, 381
POST CATH DIET (data).....	815
POST CATH DIET (q-pair).....	134
POST OPERATION DIAGNOSIS (data).....	815
POST PROCEDURE DIAGNOSIS (data).....	815
POST PROCEDURE DIAGNOSIS (q-pair).....	129
POST-FALL (data).....	815
POST-FALL (q-pair).....	622
POST-OPERATIVE (data).....	816

POST-OPERATIVE (q-pair).....	129, 145, 622
POST-PROCEDURE (data).....	816
POST-PROCEDURE (q-pair).....	145
POST-SEDATION (data).....	816
POST-SEDATION (q-pair).....	145
POSTERIOR TIBIAL (loca).....	939
POSTERIOR TIBIAL (q-pair).....	440
POTASSIUM MODIFIED 60 MEQ (data).....	816
POTASSIUM MODIFIED 60 MEQ (q-pair).....	134
POUND (unit).....	1038
PPN (prod).....	1000
PPN (q-pair).....	270
PR INTERVAL.....	410
PRBC (data).....	816
PRBC (prod).....	1000
PRBC (q-pair).....	270, 276
PRBC - AUTOLOGOUS (prod).....	1000
PRBC - AUTOLOGOUS (q-pair).....	270, 276
PRBC - LEUKOREduced (prod).....	1000
PRBC - LEUKOREduced (q-pair).....	270, 276
PRBC - LEUKOREduced - DIRECTED (prod).....	1000
PRBC - LEUKOREduced - DIRECTED (q-pair).....	270, 276
PRBC - LEUKOREduced - IRRADIATED (prod).....	1001
PRBC - LEUKOREduced - IRRADIATED (q-pair).....	270, 276
PRBC - LEUKOREduced - WASHED (prod).....	1001
PRBC - LEUKOREduced - WASHED (q-pair).....	270, 276
PRE-MEDICATED FOR ACTIVITY (data).....	817
PRE-MEDICATED FOR ACTIVITY (q-pair).....	386
PRE-OPERATIVE (data).....	817
PRE-OPERATIVE (q-pair).....	145
PRE-PROCEDURE (data).....	817
PRE-PROCEDURE (q-pair).....	145
PRE-SEDATION (data).....	817
PRE-SEDATION (q-pair).....	145
PRE-SURGERY (data).....	817
PRECAUTIONS - 24 HOUR CHART CHECK (task).....	1019
PRECAUTIONS - ASPIRATION (task).....	1019
PRECAUTIONS - ETOH W/D (task).....	1020
PRECAUTIONS - FALL (task).....	1020
PRECAUTIONS - PATIENT TEACHING (task).....	1020
PRECAUTIONS - VALIDATE ALARMS (task).....	1020
PREMATURE ATRIAL (data).....	818
PRESENT (data).....	818
PRESENT (q-pair).....	74, 81, 125, 137, 158, 162, 308, 569, 572
PRESSURE (data).....	818
PRESSURE (q-pair).....	198, 381, 619
PRESSURE CALIBRATION (data).....	818
PRESSURE CALIBRATION (q-pair).....	302
PRESSURE CONTROL (data).....	818
PRESSURE CONTROL (q-pair).....	586
PRESSURE DRESSING (data).....	819
PRESSURE DRESSING (q-pair).....	43, 87, 131, 354
PRESSURE ULCER (data).....	819
PRESSURE ULCER (q-pair).....	627

PRESSURE ULCER COLOR.....	411
PRESSURE ULCER DEPTH.....	412
PRESSURE ULCER DRAINAGE CHARACTER.....	413
PRESSURE ULCER DRESSING PREVENTS ASSESSMENT.....	414
PRESSURE ULCER LENGTH.....	415
PRESSURE ULCER ODOR.....	416
PRESSURE ULCER PHOTOGRAPHED.....	417
PRESSURE ULCER PREVENTION DEVICE TYPE.....	418
PRESSURE ULCER RISK FACTORS.....	419
PRESSURE ULCER SINUS/UNDERMINING.....	420
PRESSURE ULCER STAGE/SITE.....	421
PRESSURE ULCER TREATMENT.....	422
PRESSURE ULCER TUNNEL.....	423
PRESSURE ULCER TUNNEL LENGTH.....	424
PRESSURE ULCER UNDERMINING.....	425
PRESSURE ULCER WIDTH.....	426
PRESYSTOLIC (data).....	819
PRESYSTOLIC (q-pair).....	82
PREVENT ACCESS TO DEVICE (data).....	819
PREVENT ACCESS TO DEVICE (q-pair).....	456
PREVIOUS PERCUTANEOUS SITE (data).....	819
PREVIOUS PERCUTANEOUS SITE (q-pair).....	627
PRN (data).....	820
PRN (q-pair).....	179
PRN (sche).....	1008
PROCEDURE ROOM (data).....	820
PROCEDURE ROOM (q-pair).....	379
PRODUCTIVE (data).....	820
PRODUCTIVE (q-pair).....	126
PROJECTILE VOMITING (data).....	820
PROJECTILE VOMITING (q-pair).....	172
PRONE (data).....	820
PRONE (q-pair).....	12
PROSTHESIS (data).....	821
PROXIMAL PORT (data).....	821
PROXIMAL PORT (q-pair).....	303
PT DENIES NAUSEA / NO VOMITING (data).....	821
PT DENIES NAUSEA / NO VOMITING (q-pair).....	377
PT PAIN AT ACCEPTABLE LEVEL (data).....	821
PT PAIN AT ACCEPTABLE LEVEL (q-pair).....	377
PULMONARY - CHEST PHYSIO-THERAPY (task).....	1020
PULMONARY - COUGH AND DEEP BREATH (task).....	1021
PULMONARY - ET TUBE CARE (task).....	1021
PULMONARY - INCENTIVE SPIROMETER (task).....	1021
PULMONARY - OXYGEN (task).....	1021
PULMONARY - SUCTION PATIENT (task).....	1021
PULMONARY - TRACH CARE (task).....	1022
PULMONARY ARTERY DIASTOLIC.....	427
PULMONARY ARTERY MEAN PRESSURE.....	428
PULMONARY ARTERY SYSTOLIC.....	429
PULMONARY CAPILLARY WEDGE PRESSURE.....	430
PULMONARY CARE.....	431
PULMONARY SECRETION AMOUNT.....	432
PULMONARY SECRETION COLOR.....	433
PULMONARY SECRETION CONSISTENCY.....	434
PULMONARY SUCTIONING EVENT.....	435

PULMONARY SUCTIONING FREQUENCY.....	436
PULMONARY SUCTIONING SITE.....	437
PULMONARY VASCULAR RESISTANCE.....	438
PULMONARY VASCULAR RESISTANCE INDEX.....	439
PULSE.....	440
PULSE OX MONITOR (meth).....	977
PULSE OX MONITOR (q-pair).....	556
PULSE STRENGTH.....	441
PUNCTURE (data).....	821
PUNCTURE (q-pair).....	627
PUPIL REACTION.....	442
PUPIL SIZE.....	443
PUREED (data).....	822
PUREED (q-pair).....	134
PURULENT (data).....	822
PURULENT (q-pair).....	99, 138, 174, 258, 413, 574
PVC (data).....	822
PVC (q-pair).....	75
pH (data).....	909
pH (q-pair).....	175
pH (unit).....	1042

- Q -

Q2H (sche).....	1008
Q4H (sche).....	1008
Q7D (sche).....	1008
QD (sche).....	1009
QRS.....	444
QT.....	445
QUAD LUMEN CATHETER (unit).....	1039
QUIET OR RESTING (data).....	822
QUIET OR RESTING (q-pair).....	117
QUIET ROOM (data).....	822
QUIET ROOM (q-pair).....	51, 118, 386

- R -

R ABDOMEN (loca).....	939
R ABDOMEN (q-pair).....	622
R ANKLE (loca).....	939
R ANKLE (q-pair).....	265, 384, 421, 470, 523, 622
R ANTECUBITAL (loca).....	940
R ANTECUBITAL (q-pair).....	92, 130, 182, 303, 304, 421, 622
R ARM (loca).....	940
R ARM (q-pair).....	58, 92, 130, 133, 141, 146, 182, 303, 304, 384, 468, 469, 523, 554
R ARTERY (loca).....	940
R AXILLARY (loca).....	940
R AXILLARY (q-pair).....	41, 45, 198, 201
R BRACHIAL (loca).....	940

R BRACHIAL (q-pair).....	41, 45, 95, 368, 373
R BUNDLE BRANCH BLOCK (data).....	823
R BUNDLE BRANCH BLOCK (q-pair).....	79
R CALF (loca).....	941
R CALF (q-pair).....	265
R CAROTID (loca).....	941
R CAROTID (q-pair).....	74, 81
R CHEST (loca).....	941
R CHEST (q-pair).....	98, 384, 523, 622
R DORSALIS PEDIS (loca).....	941
R EAR (loca).....	941
R EAR (q-pair).....	265, 384, 421, 622
R ELBOW (loca).....	942
R ELBOW (q-pair).....	384, 421, 523, 622
R EXTERNAL JUGULAR (loca).....	942
R EXTERNAL JUGULAR (q-pair).....	92, 130, 182, 303, 304, 361
R EYE (loca).....	942
R EYE (q-pair).....	122, 141, 442, 443, 602, 622
R FACE (loca).....	942
R FACE (q-pair).....	141, 265, 622
R FEMORAL (loca).....	942
R FEMORAL (q-pair).....	41, 45, 74, 81, 92, 95, 130, 182, 198, 201, 265, 303, 304, 361, 368, 373, 421, 441, 622
R FINGER (loca).....	943
R FINGER (q-pair).....	265
R FLANK (loca).....	943
R FLANK (q-pair).....	265
R FOOT (loca).....	943
R FOOT (q-pair).....	141, 160, 265, 303, 304, 384, 421, 483, 523, 622
R FOOT DOWN (data).....	823
R FOOT DOWN (q-pair).....	499
R FOOT UP (data).....	823
R FOOT UP (q-pair).....	499
R GLUTEUS (loca).....	943
R GLUTEUS (q-pair).....	265, 421, 622
R GROIN (loca).....	943
R GROIN (q-pair).....	384, 523
R HAND (loca).....	944
R HAND (q-pair).....	141, 185, 384, 484, 523, 622
R HEAD (loca).....	944
R HEAD (q-pair).....	421, 622
R HEAD DOWN (data).....	823
R HEAD DOWN (q-pair).....	499
R HEAD UP (data).....	823
R HEAD UP (q-pair).....	499
R HEEL (loca).....	944
R HEEL (q-pair).....	421, 622

R HIP (loca).....	944
R HIP (q-pair).....	265, 384, 523, 622
R INGUINAL (loca).....	944
R INGUINAL (q-pair).....	265
R INNER KNEE (loca).....	945
R INNER KNEE (q-pair).....	421, 622
R INTERNAL JUGULAR (loca).....	945
R INTERNAL JUGULAR (q-pair).....	92, 95, 130, 182, 303, 304, 361, 368, 373
R ISCHIAL (loca).....	945
R ISCHIAL (q-pair).....	421, 622
R JAW (loca).....	945
R JAW (q-pair).....	384, 523
R KNEE (loca).....	945
R KNEE (q-pair).....	141, 265, 384, 523
R LATERAL HEEL (loca).....	946
R LATERAL HEEL (q-pair).....	421
R LATERAL MALLEOLUS (loca).....	946
R LATERAL MALLEOLUS (q-pair).....	421
R LEG (loca).....	946
R LEG (q-pair).....	58, 125, 133, 141, 146, 303, 304, 384, 468, 470, 523, 554
R LOWER ANTERIOR CHEST (loca).....	946
R LOWER ANTERIOR CHEST (q-pair).....	100, 265
R LOWER ARM (loca).....	946
R LOWER ARM (q-pair).....	265, 384, 421, 523, 622
R LOWER EXTREMITY (loca).....	947
R LOWER EXTREMITY (q-pair).....	325, 485
R LOWER LEG (loca).....	947
R LOWER LEG (q-pair).....	265, 622
R LOWER LIMB (loca).....	947
R LOWER LIMB (q-pair).....	310
R LOWER LOBE (loca).....	947
R LOWER LOBE (q-pair).....	69
R LOWER POSTERIOR CHEST (loca).....	947
R LOWER POSTERIOR CHEST (q-pair).....	100, 265
R MANDIBULAR (loca).....	948
R MANDIBULAR (q-pair).....	265
R MEDIAL HEEL (loca).....	948
R MEDIAL HEEL (q-pair).....	421
R MEDIAL MALLEOLUS (loca).....	948
R MEDIAL MALLEOLUS (q-pair).....	421
R MIDDLE ANTERIOR CHEST (loca).....	948
R MIDDLE ANTERIOR CHEST (q-pair).....	100, 265
R MIDDLE LOBE (loca).....	948
R MIDDLE LOBE (q-pair).....	69
R MIDDLE POSTERIOR CHEST (loca).....	949
R MIDDLE POSTERIOR CHEST (q-pair).....	100, 265
R NECK (data).....	824
R NECK (loca).....	949
R NECK (q-pair).....	141, 265
R NEPHROSTOMY (data).....	824
R NEPHROSTOMY (q-pair).....	182, 577

R NOSE (loca).....	949
R NOSE (q-pair).....	384, 523
R OCCIPITAL (loca).....	949
R OCCIPITAL (q-pair).....	265
R ORAL (loca).....	949
R ORAL (q-pair).....	265
R ORBIT (loca).....	950
R ORBIT (q-pair).....	265
R OUTER KNEE (loca).....	950
R OUTER KNEE (q-pair).....	421, 622
R PARIETAL (loca).....	950
R PARIETAL (q-pair).....	265
R PEDAL (loca).....	950
R PEDAL (q-pair).....	441
R POST TIBIAL (loca).....	950
R POST TIBIAL (q-pair).....	441
R RADIAL (loca).....	951
R RADIAL (q-pair).....	41, 45, 441
R RIBS (loca).....	951
R RIBS (q-pair).....	421, 622
R SCAPULA (loca).....	951
R SCAPULA (q-pair).....	384, 421, 523, 622
R SHOULDER (loca).....	951
R SHOULDER (q-pair).....	384, 421, 523, 622
R STENT (data).....	824
R STENT (q-pair).....	182
R SUBCLAVIAN (loca).....	951
R SUBCLAVIAN (q-pair).....	92, 95, 130, 182, 265, 303, 304, 361, 368, 373
R TEMPORAL (loca).....	952
R TEMPORAL (q-pair).....	74, 81
R TESTICLE (loca).....	952
R TESTICLE (q-pair).....	622
R THIGH (loca).....	952
R THIGH (q-pair).....	265, 384, 523
R THORACIC ANTERIOR (loca).....	952
R THORACIC ANTERIOR (q-pair).....	265
R THORACIC LATERAL (loca).....	952
R THORACIC LATERAL (q-pair).....	265
R THORACIC POSTERIOR (loca).....	953
R THORACIC POSTERIOR (q-pair).....	265
R THORACOTOMY (loca).....	953
R THORACOTOMY (q-pair).....	265
R THUMB (loca).....	953
R THUMB (q-pair).....	265
R TOE (loca).....	953
R TOE (q-pair).....	265
R TRAPEZIUS (loca).....	953
R TRAPEZIUS (q-pair).....	265
R TROCHANTER (loca).....	954
R TROCHANTER (q-pair).....	421, 622
R UPPER ANTERIOR CHEST (loca).....	954
R UPPER ANTERIOR CHEST (q-pair).....	100, 265
R UPPER ARM (loca).....	954
R UPPER ARM (q-pair).....	92, 130, 182, 265,

	303, 304, 384, 421, 523, 622
R UPPER EXTREMITY (loca).....	954
R UPPER EXTREMITY (q-pair).....	325, 486
R UPPER LEG (loca).....	954
R UPPER LEG (q-pair).....	622
R UPPER LIMB (loca).....	955
R UPPER LIMB (q-pair).....	573
R UPPER LOBE (loca).....	955
R UPPER LOBE (q-pair).....	69
R UPPER POSTERIOR CHEST (loca).....	955
R UPPER POSTERIOR CHEST (q-pair).....	100, 265
R WRIST (loca).....	955
R WRIST (q-pair).....	92, 130, 182, 303, 304, 468, 469, 622
RADIAL (loca).....	955
RADIAL (q-pair).....	58, 440
RADIAL (qual).....	1006
RADIATING (data).....	824
RADIATING (q-pair).....	381
RAMSEY ANXIOUS (data).....	824
RAMSEY ANXIOUS (q-pair).....	475
RAMSEY ASLEEP BRISK RESPONSE (data).....	825
RAMSEY ASLEEP BRISK RESPONSE (q-pair).....	475
RAMSEY ASLEEP NO RESPONSE (data).....	825
RAMSEY ASLEEP NO RESPONSE (q-pair).....	475
RAMSEY ASLEEP SLUGGISH (data).....	825
RAMSEY ASLEEP SLUGGISH (q-pair).....	475
RAMSEY COOPERATIVE (data).....	825
RAMSEY COOPERATIVE (q-pair).....	475
RAMSEY SEDATED (data).....	825
RAMSEY SEDATED (q-pair).....	475
RANGE (unit).....	1039
RASH/EXFOLIATION (data).....	826
RASH/EXFOLIATION (q-pair).....	545, 627
RASS AGITATED (data).....	826
RASS AGITATED (q-pair).....	476
RASS ALERT AND CALM (data).....	826
RASS ALERT AND CALM (q-pair).....	476
RASS COMBATIVE (data).....	826
RASS COMBATIVE (q-pair).....	476
RASS DEEP SEDATION (data).....	826
RASS DEEP SEDATION (q-pair).....	476
RASS DROWSY (data).....	827
RASS DROWSY (q-pair).....	476
RASS LIGHT SEDATION (data).....	827
RASS LIGHT SEDATION (q-pair).....	476
RASS MODERATE SEDATION (data).....	827
RASS MODERATE SEDATION (q-pair).....	476
RASS RESTLESS (data).....	827
RASS RESTLESS (q-pair).....	476
RASS UNAROUSABLE (data).....	827
RASS UNAROUSABLE (q-pair).....	476
RASS VERY AGITATED (data).....	828
RASS VERY AGITATED (q-pair).....	476
REACTIVE (data).....	828

REACTIVE (q-pair).....	442
RECTAL (loca).....	956
RECTAL (q-pair).....	558, 622
RECTAL TEMP PROBE (meth).....	977
RECTAL TUBE (data).....	828
RECTAL TUBE (q-pair).....	336
RECTAL TUBE PRESENT.....	446
RED (data).....	828
RED (q-pair).....	41, 85, 130, 138, 174, 197, 304, 353, 433, 545, 575, 610, 614
RED/GRANULATING (data).....	828
RED/GRANULATING (q-pair).....	411
REDDENING OF THE SKIN (data).....	829
REDDENING OF THE SKIN (q-pair).....	421
REDIRECT PATIENT ACTIVITY (data).....	829
REDIRECT PATIENT ACTIVITY (q-pair).....	456
REFLEX (data).....	829
REFLEX (q-pair).....	160, 185, 310, 326, 573
REFUSE MEDICATION (data).....	829
REFUSE MEDICATION (q-pair).....	386
REFUSE TO MOVE (data).....	829
REFUSE TO MOVE (q-pair).....	419
REFUSED (data).....	830
REFUSED (q-pair).....	14
REGISTERED NURSE (data).....	830
REGISTERED NURSE (q-pair).....	462
REGULAR (data).....	830
REGULAR (q-pair).....	134, 447
REINFORCED (data).....	830
REINFORCED (q-pair).....	42, 86
RELIGIOUS (data).....	830
RELIGIOUS (q-pair).....	143
REMOVE FROM STIM. SETTING (data).....	831
REMOVE FROM STIM. SETTING (q-pair).....	456
RENAL DIET (data).....	831
RENAL DIET (q-pair).....	134
REORIENTED TO ENVIRONMENT (data).....	831
REORIENTED TO ENVIRONMENT (q-pair).....	456
REPORTED (data).....	831
REPORTED (q-pair).....	577
REPOSITION (data).....	831
REPOSITION (q-pair).....	386, 456, 463, 487
REQUEST FOR MORE MEDICATION (data).....	832
REQUEST FOR MORE MEDICATION (q-pair).....	386
RES CRIT: AGGRESSIVE (data).....	832
RES CRIT: AGGRESSIVE (q-pair).....	449
RES CRIT: AMBULATE (data).....	832
RES CRIT: AMBULATE (q-pair).....	449
RES CRIT: COMBATIVE (data).....	832
RES CRIT: COMBATIVE (q-pair).....	449
RES CRIT: EQUIPMENT (data).....	832
RES CRIT: EQUIPMENT (q-pair).....	449
RES CRIT: HARM TO SELF (data).....	833

RES CRIT: HARM TO SELF (q-pair).....	449
RES CRIT: LINES (data).....	833
RES CRIT: LINES (q-pair).....	449
RES CRIT: UNABLE TO FOLLOW DIRECTIONS (data).....	833
RES CRIT: UNABLE TO FOLLOW DIRECTIONS (q-pair).....	449
RES EVAL: AMBULATION/EXERCISE (data).....	833
RES EVAL: AMBULATION/EXERCISE (q-pair).....	452
RES EVAL: BED/WHEELCHAIR ALARM (posi).....	985
RES EVAL: DECREASE STIMULATION (data).....	833
RES EVAL: DECREASE STIMULATION (q-pair).....	452
RES EVAL: DIVERSIONAL ACTIVITIES (data).....	834
RES EVAL: DIVERSIONAL ACTIVITIES (q-pair).....	452
RES EVAL: EDUCATION (data).....	834
RES EVAL: EDUCATION (q-pair).....	452
RES EVAL: INVOLVEMENT WITH FAMILY (data).....	834
RES EVAL: INVOLVEMENT WITH FAMILY (q-pair).....	452
RES EVAL: MEDICATION REVIEW (data).....	834
RES EVAL: MEDICATION REVIEW (q-pair).....	452
RES EVAL: REALITY ORIENTATION (data).....	834
RES EVAL: REALITY ORIENTATION (posi).....	985
RES EVAL: REALITY ORIENTATION (q-pair).....	452
RES EVAL: RELOCATE (data).....	835
RES EVAL: RELOCATE (q-pair).....	452
RES EVAL: UNABLE TO ATTEMPT (data).....	835
RES EVAL: UNABLE TO ATTEMPT (q-pair).....	452
RES LRM: DEVICE ACCESS (data).....	835
RES LRM: DEVICE ACCESS (q-pair).....	463
RES LRM: LIMIT SETTING (data).....	835
RES LRM: LIMIT SETTING (q-pair).....	463
RES LRM: PRN MEDS (data).....	835
RES LRM: PRN MEDS (q-pair).....	463
RES LRM: REDIRECT (data).....	836
RES LRM: REDIRECT (q-pair).....	463
RES LRM: REMOVE FROM SETTING (data).....	836
RES LRM: REMOVE FROM SETTING (q-pair).....	463
RES LRM: REORIENT (data).....	836
RES LRM: REORIENT (q-pair).....	463
RES LRM: SUPERVISION (data).....	836
RES LRM: SUPERVISION (q-pair).....	463
RES LRM: VERBAL DE-ESCALATION (data).....	836
RES LRM: VERBAL DE-ESCALATION (q-pair).....	463
RES LRM: VERBAL REMINDERS (data).....	837
RES LRM: VERBAL REMINDERS (q-pair).....	463
RESERVOIR CANNULA (meth).....	977
RESERVOIR CANNULA (q-pair).....	551
RESPIRATORY PATTERN.....	447
RESPIRATORY RATE.....	448
RESTLESS (data).....	837
RESTLESS (q-pair).....	52, 309, 385, 464
RESTRAINT (data).....	837
RESTRAINT (q-pair).....	465
RESTRAINT - ASSESS NEED FOR FLUIDS/TOILET (task).....	1022
RESTRAINT - CIRCULATION CHECK (task).....	1022
RESTRAINT - DEVICE RELEASE (task).....	1022
RESTRAINT - FALL PRECAUTIONS (task).....	1022
RESTRAINT - SEIZURE PRECAUTIONS (task).....	1023

RESTRAINT BILATERAL ANKLE (data).....	837
RESTRAINT BILATERAL ANKLE (q-pair).....	468
RESTRAINT BILATERAL HAND MITT (data).....	837
RESTRAINT BILATERAL HAND MITT (q-pair).....	468
RESTRAINT CRITERIA.....	449
RESTRAINT EDUCATION.....	450
RESTRAINT EDUCATION TYPE.....	451
RESTRAINT EVAL - ALT. MEASURES ATTEMPTED.....	452
RESTRAINT EVAL - ELIMINATION MONITORED.....	453
RESTRAINT EVAL - FLUID/NUTRITION INTAKE.....	454
RESTRAINT EVAL - HYGIENE/COMFORT.....	455
RESTRAINT EVAL - LEAST RESTRICTIVE MEASURES.....	456
RESTRAINT EVAL - MENTAL.....	457
RESTRAINT EVAL - PRIVACY/DIGNITY.....	458
RESTRAINT EVAL - REPOSITION/RANGE OF MOTION.....	459
RESTRAINT EVAL - SKIN INTEGRITY.....	460
RESTRAINT FAMILY ED. UNDERSTANDING.....	461
RESTRAINT GERI CHAIR (data).....	838
RESTRAINT GERI CHAIR (q-pair).....	468
RESTRAINT HARD PLASTIC/LEATHER (data).....	838
RESTRAINT HARD PLASTIC/LEATHER (q-pair).....	468, 469, 470
RESTRAINT HARD PLASTIC/LEATHER 2 POINT (data).....	838
RESTRAINT HARD PLASTIC/LEATHER 2 POINT (q-pair).....	468
RESTRAINT HARD PLASTIC/LEATHER 4 POINT (data).....	838
RESTRAINT HARD PLASTIC/LEATHER 4 POINT (q-pair).....	468
RESTRAINT INITIATED BY.....	462
RESTRAINT LEAST RESTRICTIVE MEASURES.....	463
RESTRAINT MENTAL ASSESSMENT.....	464
RESTRAINT ONE ARM (data).....	838
RESTRAINT ONE ARM (q-pair).....	468
RESTRAINT ONE LEG (data).....	839
RESTRAINT ONE LEG (q-pair).....	468
RESTRAINT ORDER TYPE.....	465
RESTRAINT PATIENT ED. UNDERSTANDING.....	466
RESTRAINT SAFETY CARE.....	467
RESTRAINT SINGLE HAND MITT (data).....	839
RESTRAINT SINGLE HAND MITT (q-pair).....	468, 469
RESTRAINT SOFT (data).....	839
RESTRAINT SOFT (q-pair).....	468, 469, 470
RESTRAINT SOFT 2 POINT (data).....	839
RESTRAINT SOFT 2 POINT (q-pair).....	468
RESTRAINT SOFT 4 POINT (data).....	839
RESTRAINT SOFT 4 POINT (q-pair).....	468
RESTRAINT SOFT ARMS + LOCKED WAIST (data).....	840
RESTRAINT SOFT ARMS + LOCKED WAIST (q-pair).....	468
RESTRAINT SOFT ARMS + WAIST (data).....	840
RESTRAINT SOFT ARMS + WAIST (q-pair).....	468
RESTRAINT STRAP (data).....	840
RESTRAINT STRAP (q-pair).....	468, 470
RESTRAINT STRAP (sche).....	1009
RESTRAINT TYPE ALL.....	468
RESTRAINT TYPE ARMS.....	469
RESTRAINT TYPE LEGS.....	470
RESTRAINT TYPE OTHER.....	471
RESTRAINT VEST (data).....	840
RESTRAINT VEST (q-pair).....	468, 471

RESTRAINT WAIST (data).....	840
RESTRAINT WAIST (q-pair).....	468, 471
RESTRAINT WAIST LOCKED (data).....	841
RESTRAINT WAIST LOCKED (q-pair).....	468
RESTRAINT WRIST (data).....	841
RESTRAINT WRIST (q-pair).....	468
RESTRAINTS BROCHURE (data).....	841
RESTRAINTS BROCHURE (q-pair).....	451
RESTRAINTS ED REMOVAL (data).....	841
RESTRAINTS ED REMOVAL (q-pair).....	451
RESTRAINTS ED TYPE (data).....	841
RESTRAINTS ED TYPE (q-pair).....	451
RESTRAINTS EXPLAINED (data).....	842
RESTRAINTS EXPLAINED (q-pair).....	451
RESTRAINTS LIMIT SETTING (data).....	842
RESTRAINTS LIMIT SETTING (q-pair).....	456
RESTRAINTS TIME LIMIT (data).....	842
RESTRAINTS TIME LIMIT (q-pair).....	451
RHONCHI (data).....	842
RHONCHI (q-pair).....	69
RHYTHM (meth).....	977
RIGHT (data).....	842
RIGHT (loca).....	956
RIGHT (q-pair).....	12
RIGHT VENTRICULAR STROKE WORK.....	472
RIGHT VENTRICULAR STROKE WORK INDEX.....	473
RIGID HEEL SPLINT (data).....	843
RIGID HEEL SPLINT (q-pair).....	418
ROOM AIR (meth).....	978
ROOM AIR15 (meth).....	978
ROOM AIR20 (meth).....	978
ROTATIONAL THERAPY (data).....	843
ROTATIONAL THERAPY (q-pair).....	530
ROUND (data).....	843
ROUND (q-pair).....	2
RUB (data).....	843
RUB (q-pair).....	76
Rh IMMUNE GLOBULIN (prod).....	1001
Rh IMMUNE GLOBULIN (q-pair).....	276

- S -

S1 (data).....	843
S1 (q-pair).....	76
S2 (data).....	844
S2 (q-pair).....	76
S3 (data).....	844
S3 (q-pair).....	76
S4 (data).....	844
S4 (q-pair).....	76
SACRAL (loca).....	956
SACRAL (q-pair).....	141, 421, 622
SAFETY (data).....	844
SAFETY (q-pair).....	145
SAFETY DEVICE WAIST BELT APPLIED.....	474
SALINE (prod).....	1001

SALINE (q-pair).....	275, 284
SALINE LOCK (data).....	844
SALINE LOCK (q-pair).....	302, 400
SALINE MOIST DRESSING (data).....	845
SALINE MOIST DRESSING (q-pair).....	422, 626
SANDOSOURCE PEPTIDE (data).....	845
SANDOSOURCE PEPTIDE (prod).....	1001
SANDOSOURCE PEPTIDE (q-pair).....	270
SAS AGITATED (data).....	845
SAS AGITATED (q-pair).....	477
SAS CALM AND COOPERATIVE (data).....	845
SAS CALM AND COOPERATIVE (q-pair).....	477
SAS DANGEROUS AGITATION (data).....	845
SAS DANGEROUS AGITATION (q-pair).....	477
SAS SEDATED (data).....	846
SAS SEDATED (q-pair).....	477
SAS UNAROUSABLE (data).....	846
SAS UNAROUSABLE (q-pair).....	477
SAS VERY AGITATED (data).....	846
SAS VERY AGITATED (q-pair).....	477
SAS VERY SEDATED (data).....	846
SAS VERY SEDATED (q-pair).....	477
SCALE (unit).....	1039
SCANT (data).....	846
SCANT (q-pair).....	341, 432, 512
SCAPHOID (data).....	847
SCAPHOID (q-pair).....	2
SCROTUM (loca).....	956
SCROTUM (q-pair).....	384, 421, 523
SECLUSION (data).....	847
SECLUSION (q-pair).....	465
SECOND (unit).....	1039
SECOND DEGREE AV BLOCK: TYPE 1 (data).....	847
SECOND DEGREE AV BLOCK: TYPE 1 (q-pair).....	79
SECOND DEGREE AV BLOCK: TYPE 2 (data).....	847
SECOND DEGREE AV BLOCK: TYPE 2 (q-pair).....	79
SEDATED (data).....	847
SEDATED (q-pair).....	52, 142, 464
SEDATION SCALE - RAMSEY.....	475
SEDATION SCALE - RASS.....	476
SEDATION SCALE - SAS.....	477
SEDIMENT (data).....	848
SEDIMENT (q-pair).....	574
SEIZURE CHARACTER.....	478
SEIZURE DURATION.....	479
SEIZURE FOCI.....	480
SEIZURE PRECAUTIONS.....	481
SEIZURES.....	482
SELF AM CARE BY PATIENT (data).....	848
SELF AM CARE BY PATIENT (q-pair).....	4
SELF CARE (data).....	848
SELF CARE (q-pair).....	12, 14
SELF CARE OF BEDSIDE ENVIRONMENT (data).....	848
SELF CARE OF BEDSIDE ENVIRONMENT (q-pair).....	4
SELF FACE WASH (data).....	848
SELF FACE WASH (q-pair).....	4

SELF FOOT CARE (data).....	849
SELF FOOT CARE (q-pair).....	4
SELF SHAVE (data).....	849
SELF SHAVE (q-pair).....	4
SELF SHOWER (data).....	849
SELF SHOWER (q-pair).....	4
SELF TRACH CARE (data).....	849
SELF TRACH CARE (q-pair).....	4
SEMIFOWLERS (posi).....	985
SEMIFOWLERS (q-pair).....	58, 440, 448
SEMIFOWLERS HIGH (data).....	849
SEMIFOWLERS HIGH (q-pair).....	49
SEMIFOWLERS LOW (data).....	850
SEMIFOWLERS LOW (q-pair).....	49
SENSATION FOOT.....	483
SENSATION HAND.....	484
SENSATION LOWER EXTREMITY.....	485
SENSATION UPPER EXTREMITY.....	486
SEQUENTIAL COMPRESSION DEVICE (data).....	850
SEQUENTIAL COMPRESSION DEVICE (q-pair).....	139
SEQUENTIAL COMPRESSION DEVICE APPLICATION.....	487
SEROSANGUINOUS (data).....	850
SEROSANGUINOUS (q-pair).....	99, 258, 413, 501
SEROUS (data).....	850
SEROUS (q-pair).....	99, 258, 413, 501
SET UP HELP ONLY (data).....	850
SET UP HELP ONLY (q-pair).....	9
SEVERELY REDUCED (data).....	851
SEVERELY REDUCED (q-pair).....	306
SEWS - AGITATION.....	488
SEWS - ANXIETY.....	489
SEWS - DISORIENTATION NAME PLACE.....	490
SEWS - HALLUCINATION ANY.....	491
SEWS - HALLUCINATION AUDITORY ONLY.....	492
SEWS - NAME ONLY.....	493
SEWS - NAUSEA.....	494
SEWS - SWEATING.....	495
SEWS - TOTAL.....	496
SEWS - TREMOR.....	497
SEWS - VITAL SIGNS.....	498
SEWS SCALE (unit).....	1039
SHALLOW (data).....	851
SHALLOW (q-pair).....	447
SHAMPOO (data).....	851
SHAMPOO (q-pair).....	4
SHARP (data).....	851
SHARP (q-pair).....	381
SHAVE (data).....	851
SHAVE (q-pair).....	4
SHOOTING (data).....	852
SHOOTING (q-pair).....	381
SHOWER (data).....	852
SHOWER (q-pair).....	4
SI UNIT (unit).....	1040
SICU (data).....	852
SICU (q-pair).....	379

SIDE RAILS.....	499
SIGMOIDOSTOMY (data).....	852
SIGMOIDOSTOMY (q-pair).....	336
SIGNIFICANT EVENT.....	500
SIGNIFICANT OTHER (data).....	852
SIGNIFICANT OTHER (q-pair).....	10
SILVADENE (data).....	853
SILVADENE (q-pair).....	422
SILVER (data).....	853
SILVER SULFADIAZINE (data).....	853
SINGLE (data).....	853
SINGLE (q-pair).....	92
SINGLE LUMEN CATHETER (meth).....	978
SINUS ARREST (data).....	853
SINUS ARREST (q-pair).....	79
SINUS ARRHYTHMIA (data).....	854
SINUS ARRHYTHMIA (q-pair).....	79
SINUS BLOCK (data).....	854
SINUS BLOCK (q-pair).....	79
SINUS BRADYCARDIA (data).....	854
SINUS BRADYCARDIA (q-pair).....	79
SINUS DRAINAGE CHARACTER.....	501
SINUS DRAINAGE PRESENT.....	502
SINUS NECROTIC TISSUE APPEARANCE.....	503
SINUS NECROTIC TISSUE PRESENT.....	504
SINUS PAUSE (data).....	854
SINUS PAUSE (q-pair).....	79
SINUS RHYTHM (data).....	854
SINUS RHYTHM (q-pair).....	79
SINUS TACHYCARDIA (data).....	855
SINUS TACHYCARDIA (q-pair).....	79
SINUS TRACT SURROUNDING TISSUE CHARACTER.....	505
SINUS TRACT ULCER.....	506
SINUS TRACT ULCER EPITHELIALIZATION.....	507
SINUS TRACT ULCER GRANULIZATION.....	508
SITTING (posi).....	986
SITTING (q-pair).....	58, 133, 440, 448, 554, 607
SKIN (ADHESIVE INDICATOR) (meth).....	978
SKIN (ADHESIVE INDICATOR) (q-pair).....	146
SKIN (loca).....	956
SKIN (q-pair).....	558
SKIN ALTERATION.....	509
SKIN COLOR.....	510
SKIN CONDITION.....	511
SKIN DRAINAGE AMOUNT.....	512
SKIN ERYTHEMA EXTENDS AROUND INCISION.....	513
SKIN ERYTHEMA EXTENSION AROUND WOUND.....	514
SKIN INCISION CLOSURE.....	515
SKIN INCISION WOUND CLOSURE OPEN WOUND (data).....	855
SKIN INCISION WOUND CLOSURE OPEN WOUND (q-pair).....	515
SKIN INCISION WOUND CLOSURE STAPLES (data).....	855
SKIN INCISION WOUND CLOSURE STAPLES (q-pair).....	515
SKIN INCISION WOUND CLOSURE STERI-STRIPS (data).....	855
SKIN INCISION WOUND CLOSURE STERI-STRIPS (q-pair).....	515
SKIN INCISION WOUND CLOSURE SUTURES (data).....	855

SKIN INCISION WOUND CLOSURE SUTURES (q-pair).....	515
SKIN INDURATION AROUND INCISION EXTENDS INTO PERIWOUND LESS THAN 1 cm (data).....	856
SKIN INDURATION AROUND INCISION EXTENDS INTO PERIWOUND LESS THAN 1 cm (q-pair).....	516
SKIN INDURATION AROUND INCISION EXTENDS INTO PERIWOUND MORE THAN 1 cm (data).....	856
SKIN INDURATION AROUND INCISION EXTENDS INTO PERIWOUND MORE THAN 1 cm (q-pair).....	516
SKIN INDURATION AROUND INCISION LIMITED TO INCISION (data)	856
SKIN INDURATION AROUND INCISION LIMITED TO INCISION (q-pair)	516
SKIN INDURATION AROUND INCISION NONE (data).....	856
SKIN INDURATION AROUND INCISION NONE (q-pair).....	516
SKIN INDURATION EXTENSION AROUND INCISION.....	516
SKIN INDURATION EXTENSION AROUND WOUND.....	517
SKIN INDURATION EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND < 1 cm (data)	856
SKIN INDURATION EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND < 1 cm (q-pair).....	517
SKIN INDURATION EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND > 1 cm (data)	857
SKIN INDURATION EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND > 1 cm (q-pair).....	517
SKIN INDURATION EXTENSION AROUND WOUND LIMITED TO WOUND (data)	857
SKIN INDURATION EXTENSION AROUND WOUND LIMITED TO WOUND (q-pair)	517
SKIN INDURATION EXTENSION AROUND WOUND NONE (data).....	857
SKIN INDURATION EXTENSION AROUND WOUND NONE (q-pair).....	517
SKIN MOISTURE.....	518
SKIN PATCHES PRESENT.....	519
SKIN RISK FACTORS PRESENT.....	520
SKIN SWELLING EXTENSION AROUND INCISION.....	521
SKIN SWELLING EXTENSION AROUND INCISION EXTENDS INTO PERIWOUND < 1 cm (data).....	857
SKIN SWELLING EXTENSION AROUND INCISION EXTENDS INTO PERIWOUND < 1 cm (q-pair).....	521
SKIN SWELLING EXTENSION AROUND INCISION EXTENDS INTO PERIWOUND > 1 cm (data).....	857
SKIN SWELLING EXTENSION AROUND INCISION EXTENDS INTO PERIWOUND > 1 cm (q-pair).....	521
SKIN SWELLING EXTENSION AROUND INCISION FLUCTUANT (data)	858
SKIN SWELLING EXTENSION AROUND INCISION FLUCTUANT (q-pair)	521
SKIN SWELLING EXTENSION AROUND INCISION LIMITED TO WOUND (data)	858
SKIN SWELLING EXTENSION AROUND INCISION LIMITED TO WOUND (q-pair)	521
SKIN SWELLING EXTENSION AROUND INCISION NONE (data).....	858
SKIN SWELLING EXTENSION AROUND INCISION NONE (q-pair).....	521
SKIN SWELLING EXTENSION AROUND WOUND.....	522
SKIN SWELLING EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND < 1 cm (data)	858

SKIN SWELLING EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND < 1 cm (q-pair)	522
SKIN SWELLING EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND > 1 cm (data)	858
SKIN SWELLING EXTENSION AROUND WOUND EXTENDS INTO PERIWOUND > 1 cm (q-pair)	522
SKIN SWELLING EXTENSION AROUND WOUND FLUCTANT (data)...	859
SKIN SWELLING EXTENSION AROUND WOUND FLUCTANT (q-pair)...	522
SKIN SWELLING EXTENSION AROUND WOUND LIMITED TO WOUND (data)	859
SKIN SWELLING EXTENSION AROUND WOUND LIMITED TO WOUND (q-pair)	522
SKIN SWELLING EXTENSION AROUND WOUND NONE (data).....	859
SKIN SWELLING EXTENSION AROUND WOUND NONE (q-pair).....	522
SKIN TEAR (data).....	859
SKIN TEAR (q-pair).....	627
SKIN TEMP.....	523
SKIN TEMPERATURE AT INCISION SITE.....	524
SKIN TEMPERATURE AT INCISION SITE CONSISTENT WITH BODY TEMPERATURE (data)	859
SKIN TEMPERATURE AT INCISION SITE CONSISTENT WITH BODY TEMPERATURE (q-pair)	524
SKIN TEMPERATURE AT INCISION SITE COOLER THAN GENERAL BODY TEMPERATURE (data)	860
a) SKIN TEMPERATURE AT INCISION SITE COOLER THAN GENERAL BODY TEMPERATURE (q-pair)	524
SKIN TEMPERATURE AT INCISION SITE COOLER THAN OPPOSITE BODY PART (data)	860
SKIN TEMPERATURE AT INCISION SITE COOLER THAN OPPOSITE BODY PART (q-pair)	524
SKIN TEMPERATURE AT INCISION SITE WARMER THAN GENERAL BODY TEMPERATURE (data)	860
a) SKIN TEMPERATURE AT INCISION SITE WARMER THAN GENERAL BODY TEMPERATURE (q-pair)	524
SKIN TEMPERATURE AT INCISION SITE WARMER THAN OPPOSITE BODY PART (data)	860
SKIN TEMPERATURE AT INCISION SITE WARMER THAN OPPOSITE BODY PART (q-pair)	524
SKIN TEMPERATURE AT WOUND SITE.....	525
SKIN TEMPERATURE AT WOUND SITE CONSISTENT WITH BODY TEMPERATURE (data)	860
SKIN TEMPERATURE AT WOUND SITE CONSISTENT WITH BODY TEMPERATURE (q-pair)	525
SKIN TEMPERATURE AT WOUND SITE COOLER THAN GENERAL BODY TEMPERATURE (data)	861
SKIN TEMPERATURE AT WOUND SITE COOLER THAN GENERAL BODY TEMPERATURE (q-pair)	525
) SKIN TEMPERATURE AT WOUND SITE COOLER THAN OPPOSITE BODY PART (data)	861
SKIN TEMPERATURE AT WOUND SITE COOLER THAN OPPOSITE BODY PART (q-pair)	525
SKIN TEMPERATURE AT WOUND SITE WARMER THAN GENERAL BODY TEMPERATURE (data)	861
SKIN TEMPERATURE AT WOUND SITE WARMER THAN GENERAL BODY TEMPERATURE (q-pair)	525
) SKIN TEMPERATURE AT WOUND SITE WARMER THAN OPPOSITE BODY PART (data)	

.....	861
SKIN TEMPERATURE AT WOUND SITE WARMER THAN OPPOSITE BODY PART (q-pair)	
.....	525
SKIN TRAUMA (data).....	861
SKIN TRAUMA (q-pair).....	627
SKIN TURGOR.....	526
SKIN WARM AND DRY COLOR APPROPRIATE (data).....	862
SKIN WARM AND DRY COLOR APPROPRIATE (q-pair).....	377
SKIN WOUND CLOSURE.....	527
SKIN WOUND CLOSURE APPROXIMATED (data).....	862
SKIN WOUND CLOSURE APPROXIMATED (q-pair).....	527
SKIN WOUND CLOSURE DEHISCENCE (data).....	862
SKIN WOUND CLOSURE DEHISCENCE (q-pair).....	527
SKIN WOUND CLOSURE GLUE (data).....	862
SKIN WOUND CLOSURE GLUE (q-pair).....	527
SKIN WOUND CLOSURE NOT CLOSED (data).....	862
SKIN WOUND CLOSURE NOT CLOSED (q-pair).....	527
SKIN WOUND CLOSURE OPEN WOUND (data).....	863
SKIN WOUND CLOSURE OPEN WOUND (q-pair).....	527
SKIN WOUND CLOSURE STAPLES (data).....	863
SKIN WOUND CLOSURE STAPLES (q-pair).....	527
SKIN WOUND CLOSURE STERI-STRIPS (data).....	863
SKIN WOUND CLOSURE STERI-STRIPS (q-pair).....	527
SKIN WOUND CLOSURE SUTURES (data).....	863
SKIN WOUND CLOSURE SUTURES (q-pair).....	527
SLEEP AMOUNT.....	528
SLEEP QUALIFIERS.....	529
SLEEPING (data).....	863
SLEEPING (q-pair).....	117, 385, 464
SLEEPS POORLY (data).....	864
SLEEPS POORLY (q-pair).....	52
SLEPT WELL (data).....	864
SLEPT WELL (q-pair).....	529
SLIGHT (data).....	864
SLIGHT (q-pair).....	349
SLING TRACTION (data).....	864
SLING TRACTION (q-pair).....	563
SLOUGH (data).....	864
SLOUGH (q-pair).....	621
SLUGGISH (data).....	865
SLUGGISH (q-pair).....	442
SM ADULT (meth).....	979
SM ADULT CUFF (meth).....	979
SMALL (data).....	865
SMALL (q-pair).....	341, 347, 349, 432,
	512
SMALL BORE FEEDING TUBE (data).....	865
SMALL BORE FEEDING TUBE (q-pair).....	178
SMEAR (data).....	865
SMEAR (q-pair).....	347
SMOKING AREA (data).....	865
SMOKING AREA (q-pair).....	51, 118
SMOKING PORCH (data).....	866
SMOKING PORCH (q-pair).....	51, 118
SMOOTH (data).....	866
SMOOTH (q-pair).....	621

SNACK (data).....	866
SNACK (q-pair).....	4
SODIUM CONTROLLED 2 GM (data).....	866
SODIUM CONTROLLED 2 GM (q-pair).....	134
SODIUM CONTROLLED 4 GM (data).....	866
SODIUM CONTROLLED 4 GM (q-pair).....	134
SOFT (data).....	867
SOFT (q-pair).....	2, 503, 544
SOFT DRINK (prod).....	1002
SOFT DRINK (q-pair).....	288
SPECIALTY BED/MATTRESS.....	530
SPECIFIC GRAVITY.....	531
SPINE CERVICAL (loca).....	957
SPINE CERVICAL (q-pair).....	265
SPINE LUMBAR (loca).....	957
SPINE LUMBAR (q-pair).....	265
SPINE THORACIC (loca).....	957
SPINE THORACIC (q-pair).....	265
SPIRITUAL (data).....	867
SPIRITUAL (q-pair).....	143, 145
SPLINTING (data).....	867
SPLINTING (q-pair).....	385
SPONTANEOUS (data).....	867
SPONTANEOUS (meth).....	979
SPONTANEOUS (q-pair).....	448, 586
SPOONING (data).....	867
SPOONING (q-pair).....	329
SQ EMPHYSEMA (data).....	868
SQ EMPHYSEMA (q-pair).....	103
SQUARE (data).....	868
SQUARE (q-pair).....	601
SQUARE METERS (unit).....	1040
SQUEEZING (data).....	868
SQUEEZING (q-pair).....	381
ST SEGMENT AVF.....	532
ST SEGMENT AVL.....	533
ST SEGMENT AVR.....	534
ST SEGMENT I.....	535
ST SEGMENT II.....	536
ST SEGMENT III.....	537
ST SEGMENT V1.....	538
ST SEGMENT V2.....	539
ST SEGMENT V3.....	540
ST SEGMENT V4.....	541
ST SEGMENT V5.....	542
ST SEGMENT V6.....	543
STABBING (data).....	868
STABBING (q-pair).....	381
STABLE (data).....	868
STABLE (q-pair).....	151
STAFF (data).....	869
STAFF (q-pair).....	567
STAGE I (data).....	869
STAGE I (q-pair).....	421
STAGE II (data).....	869
STAGE II (q-pair).....	421

STAGE III (data).....	869
STAGE III (q-pair).....	421
STAGE IV (data).....	869
STAGE IV (q-pair).....	421
STANDARD (data).....	870
STANDARD (q-pair).....	292
STANDARD PRESSURE REDUCTION (data).....	870
STANDARD PRESSURE REDUCTION (q-pair).....	530
STANDING (posi).....	986
STANDING (q-pair).....	58, 133, 440, 448, 554, 607
STANDING WEIGHT (meth).....	979
STANDING WEIGHT (q-pair).....	607
STAPLES (data).....	870
STARTED (data).....	870
STARTED (q-pair).....	302, 400
STARTED ON LEFT (data).....	870
STARTED ON LEFT (q-pair).....	480
STARTED ON RIGHT (data).....	871
STARTED ON RIGHT (q-pair).....	480
STATED (q-pair).....	189, 577, 607, 609
STATED (qual).....	1006
STEADY GATE OR BASELINE (data).....	871
STEADY GATE OR BASELINE (q-pair).....	135
STEP DOWN (data).....	871
STEP DOWN (q-pair).....	379
STERNOTOMY (loca).....	957
STINGING (data).....	871
STINGING (q-pair).....	381
STOMACH (loca).....	957
STOMACH (q-pair).....	384
STOOL CHARACTERISTICS.....	544
STOOL COLOR.....	545
STOOL NUMBER STOOLS.....	546
STRAIGHT CATHETER (data).....	871
STRAIGHT CATHETER (q-pair).....	182, 577
STRETCHER (data).....	872
STRETCHER (q-pair).....	568
STRICT (data).....	872
STRICT (q-pair).....	292
STROKE INDEX.....	547
STROKE VOLUME.....	548
STROKE VOLUME INDEX.....	549
STROKE WORK PER BEAT (unit).....	1040
STROKE WORK PER BEAT PER M2 (unit).....	1040
STRONG (data).....	872
STRONG (q-pair).....	160, 185, 310, 573
STRYKER FRAME (data).....	872
STRYKER FRAME (q-pair).....	530
STUDENT (data).....	872
STUDENT (q-pair).....	567
SUBCUTANEOUS TISSUE (data).....	873
SUBSTERNAL (loca).....	958
SUBSTERNAL (q-pair).....	384
SUCTIONED ONLY (data).....	873
SUCTIONED ONLY (q-pair).....	14

SUICIDAL (data).....	873
SUICIDAL (q-pair).....	54, 119
SUMP (data).....	873
SUMP (q-pair).....	616
SUPERVISION/SET-UP HELP (data).....	873
SUPERVISION/SET-UP HELP (q-pair).....	7, 8, 9
SUPINE (data).....	874
SUPINE (posi).....	986
SUPINE (q-pair).....	12, 58, 448
SUPPLEMENTAL OXYGEN CONCENTRATION.....	550
SUPPLEMENTAL OXYGEN FLOW RATE.....	551
SUPRAPUBIC (data).....	874
SUPRAPUBIC (q-pair).....	182, 577
SUPRAVENTRICULAR TACHYCARDIA (data).....	874
SUPRAVENTRICULAR TACHYCARDIA (q-pair).....	79
SURGICAL (data).....	874
SURGICAL (q-pair).....	627
SURGICAL DONOR (data).....	874
SURGICAL DONOR (q-pair).....	627
SURGICAL GRAFT (data).....	875
SURGICAL GRAFT (q-pair).....	627
SUSPECTED DEEP TISSUE INJURY (data).....	875
SUSPECTED DEEP TISSUE INJURY (q-pair).....	421
SUTURES (data).....	875
SWOLLEN (data).....	875
SWOLLEN (q-pair).....	41, 85, 130, 197, 304, 353
SYNCHRONOUS INTERMITTENT MANDATORY VENTILATION (data)..	875
SYNCHRONOUS INTERMITTENT MANDATORY VENTILATION (q-pair)	586
SYRINGE CHANGED (data).....	876
SYRINGE CHANGED (q-pair).....	400
SYSTEMIC VASCULAR RESISTANCE.....	552
SYSTEMIC VASCULAR RESISTANCE INDEX.....	553
SYSTOLIC (data).....	876
SYSTOLIC (q-pair).....	82
SYSTOLIC PRESSURE.....	554
SaO2%.....	555
SpO2%.....	556
SvO2%.....	557

- T -

T-PIECE (meth).....	979
T-PIECE (q-pair).....	550
TAN (data).....	876
TAN (q-pair).....	138, 174, 433, 575, 610, 614
TAPE BURN (data).....	876
TAPE BURN (q-pair).....	627
TARRY (data).....	876
TARRY (q-pair).....	544
TEA (data).....	877
TEA (q-pair).....	575
TEETH (loca).....	958
TEETH (q-pair).....	384
TEETH BRUSHED (data).....	877

TEETH BRUSHED (q-pair).....	14
TEMPERATURE.....	558
TEMPORAL (loca).....	958
TEMPORAL (q-pair).....	558
TEMPORARY (TRANSVENOUS) (data).....	877
TEMPORARY (TRANSVENOUS) (q-pair).....	373
TENACIOUS (data).....	877
TENACIOUS (q-pair).....	434
TENDER (data).....	877
TENDER (q-pair).....	2, 381
TENDON (data).....	878
TESTS/PROC. - ACCUCHECKS (task).....	1023
TESTS/PROC. - BLOOD CULTURES (task).....	1023
TESTS/PROC. - DAILY LABS (task).....	1023
TESTS/PROC. - MORNING LABS (task).....	1023
TESTS/PROC. - ZERO BALANCE HEMO. LINES (task).....	1024
THICK (data).....	878
THICK (q-pair).....	138, 173, 330, 434
THIGH (data).....	878
THIGH (loca).....	958
THIGH (q-pair).....	59, 105
THIGH CUFF (meth).....	980
THIN (data).....	878
THIN (q-pair).....	138, 173, 434
THIRD DEGREE AV BLOCK (data).....	878
THIRD DEGREE AV BLOCK (q-pair).....	79
THORACIC SPINAL PRECAUTIONS (data).....	879
THORACIC SPINAL PRECAUTIONS (q-pair).....	11
THROBBING (data).....	879
THROBBING (q-pair).....	381
THYROID (loca).....	958
THYROID (q-pair).....	265
TINGLING (data).....	879
TINGLING (q-pair).....	381, 483, 484, 485, 486
TO PAIN (data).....	879
TO PAIN (q-pair).....	122
TO SPEECH (data).....	879
TO SPEECH (q-pair).....	122
TOENAIL CARE (data).....	880
TOENAIL CARE (q-pair).....	4
TOLERATING FLUIDS/NOURISHMENT (data).....	880
TOLERATING FLUIDS/NOURISHMENT (q-pair).....	135, 377
TONIC/CLONIC (data).....	880
TONIC/CLONIC (q-pair).....	478
TONOMETRY.....	559
TOOTH (loca).....	959
TOOTH (q-pair).....	384
TOPICAL TREATMENT (data).....	880
TOPICAL TREATMENT (q-pair).....	611
TORSO (loca).....	959
TORSO (q-pair).....	141, 471
TOTAL DEPENDENCE (data).....	880
TOTAL DEPENDENCE (q-pair).....	7, 8, 9
TPN (meth).....	980
TPN (prod).....	1002

TPN (q-pair).....	270
TRACE ST SEGMENT.....	560
TRACE ST SEGMENT 2.....	561
TRACHEA MIDLINE DEVIATED.....	562
TRACHEOSTOMY (data).....	881
TRACHEOSTOMY (loca).....	959
TRACHEOSTOMY (q-pair).....	21, 421
TRACHEOSTOMY COLLAR (meth).....	980
TRACHEOSTOMY COLLAR (q-pair).....	550
TRACTION TYPE.....	563
TRACTION WEIGHT.....	564
TRAIN OF FOUR RESPONSE.....	565
TRAIN OF FOUR SMS.....	566
TRANSCUTANEOUS (data).....	881
TRANSCUTANEOUS (q-pair).....	373
TRANSFER (data).....	881
TRANSFER (q-pair).....	129, 622
TRANSFER WITH ASSIST (data).....	881
TRANSFER WITH ASSIST (q-pair).....	11
TRANSFUSION (meth).....	980
TRANSFUSION (q-pair).....	270
TRANSPARENT (data).....	881
TRANSPARENT (q-pair).....	43, 87, 131, 198, 354, 619
TRANSPARENT DRESSING (data).....	882
TRANSPARENT DRESSING (q-pair).....	422, 626
TRANSPORT FACILITATOR.....	567
TRANSPORT FACILITATOR PRESENT (data).....	882
TRANSPORT FACILITATOR PRESENT (q-pair).....	135
TRANSPORTATION MODE/DEVICE.....	568
TRANSPYLORIC (data).....	882
TRANSPYLORIC (q-pair).....	147
TRANSTRACHEAL (meth).....	980
TRANSVERSE ABDOMINAL (loca).....	959
TRANSVERSE ABDOMINAL (q-pair).....	265
TREATMENTS (data).....	882
TREATMENTS (q-pair).....	145
TREMULOUS.....	569
TRENDELENBURG (posi).....	986
TRENDELENBURG (q-pair).....	58, 448
TRIPLE (data).....	882
TRIPLE (q-pair).....	92
TRIPLE LUMEN CATHETER (meth).....	981
TUBE POSITION VERIFIED.....	570
TUBED/DYSPHAGIC (data).....	883
TUBED/DYSPHAGIC (q-pair).....	124
TUBES & DRAINS - CONDOM CATH (task).....	1024
TUBES & DRAINS - FOLEY (task).....	1024
TUBES & DRAINS - I/O CATH (task).....	1024
TUBES & DRAINS - NG TUBE (task).....	1024
TUBES & DRAINS - SUPRAPUBIC DRAIN (task).....	1025
TUBING CHANGED/LABELED (data).....	883
TUBING CHANGED/LABELED (q-pair).....	302, 400
TURN AND POSITION (data).....	883
TURN AND POSITION (q-pair).....	4
TURN AND SKIN CARE.....	571

TURNING SELF INDEPENDENTLY (data).....	883
TURNING SELF INDEPENDENTLY (q-pair).....	4
TWILL TAPE (data).....	883
TWILL TAPE (q-pair).....	7, 616
TWO PERSON ASSIST (data).....	884
TWO PERSON ASSIST (q-pair).....	5, 7, 8, 9
TYMPANIC (loca).....	959
TYMPANIC (q-pair).....	558

- U -

ULNAR (loca).....	960
ULNAR (q-pair).....	440
ULNAR DRIFT.....	572
ULTRACAL (data).....	884
ULTRACAL (prod).....	1002
ULTRACAL (q-pair).....	149
UMBILICAL (loca).....	960
UMBILICAL (q-pair).....	265
UNABLE TO ASSESS (data).....	884
UNABLE TO ASSESS (q-pair).....	142, 144
UNABLE TO INFORM (data).....	884
UNABLE TO INFORM (q-pair).....	317
UNABLE TO READ (data).....	884
UNABLE TO READ (q-pair).....	142
UNABLE TO RESPOND (data).....	885
UNABLE TO RESPOND (q-pair).....	144
UNABLE TO WEDGE (data).....	885
UNABLE TO WEDGE (q-pair).....	358
UNBEARABLE (data).....	885
UNBEARABLE (q-pair).....	381
UNCOOPERATIVE (data).....	885
UNCOOPERATIVE (q-pair).....	52
UNEQUAL (data).....	885
UNEQUAL (q-pair).....	310, 573
UNITS (unit).....	1040
UNITS PER HOUR (unit).....	1041
UNITS PER MINUTE (unit).....	1041
UNITS/kg/HOUR (unit).....	1041
UNKNOWN (data).....	886
UNKNOWN (q-pair).....	158, 480
UNSTABLE (data).....	886
UNSTABLE (q-pair).....	151, 419
UNSTAGEABLE (data).....	886
UNSTAGEABLE (q-pair).....	421
UNSUCCESSFUL (data).....	886
UNSUCCESSFUL (q-pair).....	461, 466
UP (data).....	886
UP TO CHAIR (data).....	887
UP TO CHAIR (q-pair).....	4
UPPER ABDOMEN (loca).....	960
UPPER ABDOMEN (q-pair).....	622
UPPER ARM (loca).....	960
UPPER ARM (q-pair).....	105
UPPER BACK (loca).....	960
UPPER BACK (q-pair).....	622

UPPER LIMB EXTREMITIES STRENGTH.....	573
URINE (prod).....	1002
URINE (q-pair).....	531
URINE CHARACTERISTICS.....	574
URINE COLOR.....	575
URINE NUMBER VOIDS.....	576
URINE VOIDING METHOD.....	577
UROMETER (data).....	887
UROMETER (q-pair).....	23
UROSTOMY (data).....	887
UROSTOMY (q-pair).....	182

- V -

V1 (data).....	887
V1 (q-pair).....	80
V2 (data).....	887
V2 (q-pair).....	80
V3 (data).....	888
V3 (q-pair).....	80
V4 (data).....	888
V4 (q-pair).....	80
V5 (data).....	888
V5 (q-pair).....	80
V6 (data).....	888
V6 (q-pair).....	80
VACUUM ASSISTED CLOSURE (data).....	888
VACUUM ASSISTED CLOSURE (q-pair).....	616
VASCULAR (data).....	889
VASCULAR (q-pair).....	627
VASCULAR LESION ARTERIAL (data).....	889
VASCULAR LESION ARTERIAL (q-pair).....	627
VASCULAR LESION MIXED (data).....	889
VASCULAR LESION MIXED (q-pair).....	627
VASCULAR LESION VENOUS (data).....	889
VASCULAR LESION VENOUS (q-pair).....	627
VASCULAR/NEUROPATHIC (data).....	889
VASCULAR/NEUROPATHIC (q-pair).....	627
VASOPRESSOR (data).....	890
VASOPRESSOR (q-pair).....	419
VDD (data).....	890
VDD (q-pair).....	367
VEGETARIAN DIET (data).....	890
VEGETARIAN DIET (q-pair).....	134
VENTRICULAR FIBRILLATION (data).....	890
VENTRICULAR FIBRILLATION (q-pair).....	79
VENTILATOR (data).....	890
VENTILATOR (meth).....	981
VENTILATOR (q-pair).....	23, 407, 447, 448, 550, 578, 579, 580, 581, 582, 583, 584, 585, 586, 587, 588, 589, 590, 591, 592, 593, 594, 595, 596, 597, 598, 599, 600, 601

VENTILATOR ACHIEVA P5O2 (data).....	891
VENTILATOR ACHIEVA P5O2 (q-pair).....	600
VENTILATOR ADULT STAR (data).....	891
VENTILATOR ADULT STAR (q-pair).....	600
VENTILATOR AUTO PEEP.....	578
VENTILATOR BEAR 1000 (data).....	891
VENTILATOR BEAR 1000 (q-pair).....	600
VENTILATOR CIRCUIT CHANGE.....	579
VENTILATOR DRAEGER (data).....	891
VENTILATOR DRAEGER (q-pair).....	600
VENTILATOR EXHALED VOLUME.....	580
VENTILATOR EXPIRATORY TIME.....	581
VENTILATOR FLOW RATE.....	582
VENTILATOR INSPIRATORY TIME.....	583
VENTILATOR MEAN AIR PRESSURE.....	584
VENTILATOR MINUTE VOLUME.....	585
VENTILATOR MODE.....	586
VENTILATOR OTHER (data).....	891
VENTILATOR OTHER (q-pair).....	600
VENTILATOR PATIENT EXHALE VOLUME.....	587
VENTILATOR PATIENT MINUTE VOLUME.....	588
VENTILATOR PATIENT RATE.....	589
VENTILATOR PATIENT RATE TOTAL.....	590
VENTILATOR PB 7200ae (data).....	892
VENTILATOR PB 7200ae (q-pair).....	600
VENTILATOR PB 840 (data).....	892
VENTILATOR PB 840 (q-pair).....	600
VENTILATOR PCV SET PRESSURE.....	591
VENTILATOR PEAK PRESSURE.....	592
VENTILATOR PEEP.....	593
VENTILATOR PRESSURE SUPPORT.....	594
VENTILATOR PSV SET PRESSURE.....	595
VENTILATOR RATE.....	596
VENTILATOR SENSITIVITY.....	597
VENTILATOR SERVO (data).....	892
VENTILATOR SERVO (q-pair).....	600
VENTILATOR TIDAL VOLUME ACTUAL.....	598
VENTILATOR TIDAL VOLUME SET.....	599
VENTILATOR TYPE.....	600
VENTILATOR WAVE FORM.....	601
VENTRICULAR (data).....	892
VENTRICULAR (q-pair).....	373
VENTRICULAR DRAIN (meth).....	981
VENTRICULAR DRAIN (q-pair).....	291
VENTRICULAR TACHYCARDIA (data).....	892
VENTRICULAR TACHYCARDIA (q-pair).....	79
VENTRICULAR WHIP PRESENT (data).....	893
VENTRICULAR WHIP PRESENT (q-pair).....	358
VENTRICULOSTOMY (data).....	893
VENTRICULOSTOMY (q-pair).....	616
VENTURI MASK (meth).....	981
VENTURI MASK (q-pair).....	551
VERBAL DE-ESCALATION (data).....	893
VERBAL DE-ESCALATION (q-pair).....	456
VERY LOW BMI (data).....	893
VERY LOW BMI (q-pair).....	419

VIOLET (data).....	893
VIOLET (q-pair).....	614
VIP PORT (data).....	894
VIP PORT (q-pair).....	361
VISCOUS (data).....	894
VISCOUS (q-pair).....	258, 413
VISION CORRECTED.....	602
VISION UNCORRECTED.....	603
VISUAL IMPAIRMENT (data).....	894
VISUAL IMPAIRMENT (q-pair).....	142
VITAL ENTERAL NUTRITION TYPE (data).....	894
VITAL ENTERAL NUTRITION TYPE (prod).....	1002
VITAL ENTERAL NUTRITION TYPE (q-pair).....	149
VITAL HN ENTERAL NUTRITION (data).....	894
VITAL HN ENTERAL NUTRITION (prod).....	1003
VITAL HN ENTERAL NUTRITION (q-pair).....	149
VITAL SIGNS STABLE (data).....	895
VITAL SIGNS STABLE (q-pair).....	135
VOIDED (data).....	895
VOIDED (q-pair).....	135, 377
VOLUNTEER (data).....	895
VOLUNTEER (q-pair).....	567
VOMITING (data).....	895
VOMITING (q-pair).....	172
VT/VC.....	604
VVI (data).....	895
VVI (q-pair).....	367

- W -

WAFFLE (data).....	896
WAFFLE (q-pair).....	418
WAIST (loca).....	961
WAIST (q-pair).....	468, 471
WAITING ROOM (data).....	896
WAITING ROOM (q-pair).....	152, 379
WALKER (data).....	896
WALKER (q-pair).....	5, 568
WALKING OR PACING (data).....	896
WALKING OR PACING (q-pair).....	117
WANDERER (data).....	896
WANDERER (q-pair).....	38
WARD (data).....	897
WARD (q-pair).....	379
WARM (data).....	897
WARM (q-pair).....	523
WARM HUMIDITY WITH O2 (data).....	897
WARM HUMIDITY WITH O2 (q-pair).....	605
WARMER THAN GENERAL BODY TEMP (data).....	897
WARMER THAN GENERAL BODY TEMP (q-pair).....	523
WARMER THAN OPPOSITE BODY PART (data).....	897
WARMER THAN OPPOSITE BODY PART (q-pair).....	523
WARMING BLANKET (data).....	898
WARMING BLANKET (q-pair).....	605
WARMING BOTTLE (data).....	898
WARMING BOTTLE (q-pair).....	605

WARMING COOLING DEVICE.....	605
WARMING COOLING DEVICE APPLICATION.....	606
WARMING OFF (data).....	898
WARMING OFF (q-pair).....	606
WARMING ON (data).....	898
WARMING ON (q-pair).....	606
WATER (prod).....	1003
WATER (q-pair).....	270, 288
WATERSEAL (data).....	898
WATERSEAL (q-pair).....	100
WAVEFORM REFLECTS ATRIAL PLACEMENT (data).....	899
WAVEFORM REFLECTS ATRIAL PLACEMENT (q-pair).....	358
WBC - WASHED (prod).....	1003
WBC - WASHED (q-pair).....	270, 276
WEAK (data).....	899
WEAK (q-pair).....	160, 185, 310, 573
WEDGED (data).....	899
WEDGED (q-pair).....	358
WEIGHT.....	607
WEIGHT AT ADMISSION.....	608
WEIGHT DRY.....	609
WHEELCHAIR (data).....	899
WHEELCHAIR (q-pair).....	568
WHEELCHAIR EXIT (posi).....	986
WHEELCHAIR EXIT ALARM (data).....	899
WHEELCHAIR EXIT ALARM (q-pair).....	23, 452
WHEELCHAIR SCALE (meth).....	981
WHEELCHAIR SCALE (q-pair).....	607
WHEEZES (data).....	900
WHITE (data).....	900
WHITE (q-pair).....	138, 174, 433, 575, 610, 614
WHOLE BLOOD (prod).....	1003
WHOLE BLOOD (q-pair).....	270
WITH ACTIVITY (meth).....	982
WITH ACTIVITY (q-pair).....	58, 440, 448
WITH AMBULATION (meth).....	982
WITH CAST OR BRACE (meth).....	982
WITH CAST OR BRACE (q-pair).....	607
WITH PHYSICIAN (data).....	900
WITH PHYSICIAN (q-pair).....	117
WITH PROSTHESIS (meth).....	982
WITH PROSTHESIS (q-pair).....	607
WITH THERAPIST (data).....	900
WITH THERAPIST (q-pair).....	117
WITHDRAWN (data).....	900
WITHDRAWN (q-pair).....	52
WITHDRAWS TO PAIN (data).....	901
WITHDRAWS TO PAIN (q-pair).....	123
WITHOUT ASSISTANCE (data).....	901
WITHOUT ASSISTANCE (q-pair).....	5
WITHOUT PROSTHESIS (meth).....	982
WITHOUT PROSTHESIS (q-pair).....	607
WOUND (data).....	901
WOUND (q-pair).....	509
WOUND BED COLOR.....	610

WOUND CARE TYPE.....	611
WOUND DEPTH.....	612
WOUND DISCHARGE CHARACTERISTICS.....	613
WOUND DISCHARGE CHARACTERISTICS BLOODY (data).....	901
WOUND DISCHARGE CHARACTERISTICS BLOODY (q-pair).....	613
WOUND DISCHARGE CHARACTERISTICS CLOUDY (data).....	901
WOUND DISCHARGE CHARACTERISTICS CLOUDY (q-pair).....	613
WOUND DISCHARGE CHARACTERISTICS NONE (data).....	902
WOUND DISCHARGE CHARACTERISTICS NONE (q-pair).....	613
WOUND DISCHARGE CHARACTERISTICS PURULENT (data).....	902
WOUND DISCHARGE CHARACTERISTICS PURULENT (q-pair).....	613
WOUND DISCHARGE CHARACTERISTICS SEROSANGUINOUS (data)..	902
WOUND DISCHARGE CHARACTERISTICS SEROSANGUINOUS (q-pair)	613
WOUND DISCHARGE CHARACTERISTICS SEROUS (data).....	902
WOUND DISCHARGE CHARACTERISTICS SEROUS (q-pair).....	613
WOUND DISCHARGE CHARACTERISTICS VISCOUS (data).....	902
WOUND DISCHARGE CHARACTERISTICS VISCOUS (q-pair).....	613
WOUND DISCHARGE COLOR.....	614
WOUND DISCHARGE ODOR.....	615
WOUND DRAIN TYPE.....	616
WOUND DRAINAGE (prod).....	1003
WOUND DRAINAGE (q-pair).....	337
WOUND DRAINAGE AMOUNT.....	617
WOUND DRAINAGE AMOUNT HEAVY (data).....	903
WOUND DRAINAGE AMOUNT HEAVY (q-pair).....	617
WOUND DRAINAGE AMOUNT MODERATE (data).....	903
WOUND DRAINAGE AMOUNT MODERATE (q-pair).....	617
WOUND DRAINAGE AMOUNT NONE (data).....	903
WOUND DRAINAGE AMOUNT NONE (q-pair).....	617
WOUND DRAINAGE AMOUNT SCANT (data).....	903
WOUND DRAINAGE AMOUNT SCANT (q-pair).....	617
WOUND DRAINAGE AMOUNT SMALL (data).....	903
WOUND DRAINAGE AMOUNT SMALL (q-pair).....	617
WOUND DRESSING CHANGE.....	618
WOUND DRESSING TYPE.....	619
WOUND DRESSING TYPE ANTIMICROBIAL (data).....	904
WOUND DRESSING TYPE CALCIUM ALGINATE (data).....	904
WOUND DRESSING TYPE GAUZE (data).....	904
WOUND DRESSING TYPE GAUZE PACKING (data).....	904
WOUND DRESSING TYPE GAUZE WRAP (data).....	904
WOUND DRESSING TYPE HYDROCOLLOIDAL (data).....	905
WOUND DRESSING TYPE HYDROGEL (data).....	905
WOUND DRESSING TYPE NONE (data).....	905
WOUND DRESSING TYPE PRESSURE (data).....	905
WOUND DRESSING TYPE TRANSPARENT (data).....	905
WOUND EDGE.....	620
WOUND EDGE INDISTINCT/ATTACHED (data).....	906
WOUND EDGE INDISTINCT/ATTACHED (q-pair).....	620
WOUND EDGE NON-ATTACHED (data).....	906
WOUND EDGE NON-ATTACHED (q-pair).....	620
WOUND EDGE ROLLED (data).....	906
WOUND EDGE ROLLED (q-pair).....	620
WOUND EDGE TUNNELING (data).....	906
WOUND EDGE TUNNELING (q-pair).....	620
WOUND EDGE UNDERMINING (data).....	906
WOUND EDGE UNDERMINING (q-pair).....	620

WOUND HEALING STATUS.....	621
WOUND IDENTIFIED/SITE.....	622
WOUND IRRIGATION (prod).....	1004
WOUND LENGTH.....	623
WOUND ODOR.....	624
WOUND ODOR FOUL (data).....	907
WOUND ODOR FOUL (q-pair).....	624
WOUND ODOR NONE (data).....	907
WOUND ODOR NONE (q-pair).....	624
WOUND PHOTOGRAPHED.....	625
WOUND TREATMENT.....	626
WOUND TYPE.....	627
WOUND VAC DEVICE (data).....	907
WOUND VAC DEVICE (q-pair).....	616
WOUND VACUUM (prod).....	1004
WOUND VACUUM (q-pair).....	337
WOUND VACUUM ASSISTED CLOSURE.....	628
WOUND WIDTH.....	629
WRIST (loca).....	961
WRIST (q-pair).....	105

- X -

X-RAY (data).....	907
X-RAY (q-pair).....	175
XENADERM (data).....	907
XENADERM (q-pair).....	422
XEROFORM (data).....	908
XEROFORM (q-pair).....	422, 626

- Y -

YEAR (unit).....	1041
YELLOW (data).....	908
YELLOW (q-pair).....	138, 174, 433, 503, 545, 575, 610, 614
YELLOW/FIBRINOUS (data).....	908
YELLOW/FIBRINOUS (q-pair).....	411
YELLOW/GREEN (data).....	908
YELLOW/GREEN (q-pair).....	545
YES (data).....	908
YES (q-pair).....	13, 15, 17, 22, 31, 32, 33, 34, 35, 36, 37, 39, 40, 44, 47, 48, 50, 57, 70, 71, 73, 83, 89, 94, 97, 101, 102, 128, 154, 155, 164, 171, 180, 188, 190, 195, 196, 203, 257, 260, 265, 266, 267, 268, 269, 332, 352, 359, 363, 366, 368, 370, 378, 380, 393, 403, 408, 409, 414, 417, 420, 423, 425, 431, 446,

453, 454, 455, 457,
458, 459, 460, 474,
481, 482, 502, 504,
506, 507, 508, 519,
520, 562, 570, 571,
579, 618, 625