

# **EVENT CAPTURE GRAPHICAL USER INTERFACE**

## **RELEASE NOTES**



**Software Version 2.0**

**Patch EC\*2\*116**

**December 2012**

**Department of Veterans Affairs  
Office of Information and Technology (OIT)**



## Table of Contents

1	Defect Fixes with Remedy Tickets .....	1
1.1	Splash Screen Does Not Close When Launched in CPRS [732662].....	1
1.2	Associated Clinic Cannot Be Changed on Data Entry by Patient Screen When New Procedure and Category Are Set [726618 and duplicates] .....	1



# 1 Defect Fixes with Remedy Tickets

The EC\*2\*116 patch addresses two issues identified in the Event Capture GUI module since the release of patch EC\*2\*112 in June 2012.

## 1.1 Splash Screen Does Not Close When Launched in CPRS [732662]

### Issue

The splash screen does not close when launched from CPRS. If the Event Capture application is minimized and restored, the splash screen stays on top of all EC forms.

### Solution

The splash screen is closed in the CPRS specific Event Capture launch code.

## 1.2 Associated Clinic Cannot Be Changed on Data Entry by Patient Screen When New Procedure and Category Are Set [726618 and duplicates]

### Issue

Data Entry by Patient: When a recent visit is selected as a basis for adding a new procedure, there are multiple categories available. When the user selects both a new category and a new procedure, the associated clinic is set to the default value. The user is unable to change the associated clinic to the desired value because the interface is disabled.

### Solution

The code was removed that caused the associated clinic to be incorrectly set.