



# **VISUAL IMPAIRMENT SERVICE TEAM (VIST)**

## **Installation Guide**

(Revised)

Version 4.0

June 1998

Department of Veterans Affairs  
Veterans Health Administration  
Office of the Chief Information Officer



# Table of Contents


- Introduction..... 1**
- Pre Installation..... 1
- Options to be Deleted during Install ..... 2
- Templates to be Deleted during Install..... 2
- Routines to Delete ..... 4
- Routines..... 4**
- Installation..... 4**
- Sample KIDS Install ..... 4
- Post Installation ..... 7**
- System Configuration..... 7
- Resources Requirements..... 7
- Routine Mapping..... 7
- Journaling Globals ..... 7
- Files ..... 8**
- Security Keys..... 8
- Technical Assistance ..... 8



## Introduction

It is recommended that this software be loaded into a test account before going into production.

It is *very important* to read these instructions completely before installing this software! Read these instructions through at least once before attempting to install.

 During the installation of VIST V. 4.0, all data entered using the class III Blind Rehabilitation software will be moved into the new files created for VIST V. 4.0. This conversion time will vary depending on the amount of patient information that was entered into the class III files.

This version (4.0) of the Visual Impairment Service Team software package can only be run with a standard M operating system. It also requires the following VA software applications:

<b><u>Package</u></b>	<b><u>Minimum Version Needed</u></b>
VA FileMan	21.0
Kernel	8.0
Patient Information Management System (PIMS)	5.3
Fee Basis	3.5

The above software is not included in this package and must be installed before this package is completely functional.

## Pre Installation

- (1) Performing a system backup before running the KIDS install is highly recommended. We also recommend that you install the Visual Impairment Service Team (VIST) package with no VIST users on the system.
- (2) System Configuration: ^ANRV is a new global and will need to be journaled after installation. There are no specific guidelines for system configuration or global placement for the VIST software.

- (3) Installation time will take 10 to 15 minutes. The times for the installation will depend on the amount of data in your system and whether you are using a VAX or PC. Installation will take longer on a VAX system.
- (4) To begin installing the VIST package use the *Kernel Installation & Distribution System* [XPD MAIN] option. Load the **ANRV\_4\_0.KID** distribution file using the option *Load a Distribution* [XPD LOAD DISTRIBUTION]. Next, you can verify routine checksums using the option *Verify Checksums in Transport Global* [XPD PRINT CHECKSUM]. Finally, install the package using the *Install Package(s)* [XPD INSTALL BUILD] option. An example of the installation process is included in this section.

### Options to be Deleted during Install

Option Name	Menu Text
ANRV ASSESSMENT INQ	Individual Assessment - Old

### Templates to be Deleted during Install

Input	File
ANRV EDIT CHECKLIST	623061
ANRV EDIT LETTER	623033
ANRV VIST PARM	623050
ANRVACE	623036
ANRVEID	623158
ANRVRE	623158
ANRVRRRE	623160
Print	File
ANRV ADDRESS LIST	623158
ANRV ADDRESS LIST HDR	623158
ANRV AGE LIST	623158
ANRV ASSESSMENT INQ	623158
ANRV BIRTH LIST	623158
ANRV CHECKLIST	623061
ANRV CITY LIST	623158
ANRV CLAIM REC INQ	623036
ANRV COUNTY LIST	623158
ANRV DECEASED LIST	623158
ANRV ELIG AMIS LIST	623158
ANRV EYE DIAG INQ	623158
ANRV EYE DIAG LIST	623158
ANRV EYE DIAG NAR	623158
ANRV EYE DIAG PRINT	623158

<b>Print (cont'd)</b>	<b>File</b>
ANRV FEE PT	623158
ANRV INC AMIS	623158
ANRV MAIL LABELS	623158
ANRV ML BY CITY	623158
ANRV ML BY COUNTY	623158
ANRV POS LIST	623158
ANRV PRINT LETTER	623158
ANRV REFERRAL LIST	623158
ANRV REVIEW INQ	623158
ANRV ROSTER A/R	623158
ANRV STATE LIST	623158
ANRV ZIP CODE LIST	623158
ANRVACP	623036
ANRVAP	623158
ANRVED	623158
ANRVFVD	623158
ANRVID	623158
ANRVPR	623160
ANRVREV	623158
ANRVRP1	623158
ANVRRL	623160
<b>Sort</b>	<b>File</b>
ANRV ADDRESS LIST	623158
ANRV AGE LIST	623158
ANRV BIRTH LIST	623158
ANRV CHECKLIST	623061
ANRV CITY LIST	623158
ANRV COUNTY LIST	623158
ANRV DECEASED LIST	623158
ANRV ELIG AMIS LIST	623158
ANRV EYE DIAG LIST	623158
ANRV EYE DIAG NAR	623158
ANRV EYE DIAG PRINT	623158
ANRV FEE PT	623158
ANRV INC AMIS	623158
ANRV MAIL LABELS	623158
ANRV ML BY CITY	623158
ANRV ML BY COUNTY	623158
ANRV ML BY STATE	623158
ANRV POS LIST	623158

Sort (cont'd)	File
ANRV REFERRAL LIST	623158
ANRV STATE LIST	623158
ANRV ZIP CODE LIST	623158
ANRVACP	623036
ANRVAD	623158
ANRVAP	623158
ANRVED	623158
ANRVFVD	623158
ANRVID	623158
ANRVPR	623160
ANRVREV1	623158
ANRVRP1	623158
ANRVRRL	623160

## Routines to Delete

The ANRVPRE, ANRV PST, ANRVZ\*, and ANRV1\* routines will be deleted.

## Routines

The following is a list of routines you will see for VIST when you load the new routine set. The first line of each routine contains a brief description of the general function of the routine. Use the Kernel option XU FIRST LINE PRINT (*First Line Routine Print*) to print a list of just the first line of each ANRV\* routine.

ANRVAM1	ANRVAM2	ANRVAP	ANRVDEL	ANRVLET
ANRVML	ANRVML1	ANRVML2	ANRVML3	ANRVML4
ANRVPOST	ANRVPR	ANRVPR1	ANRVPR2	ANRVRP
ANRVRP2	ANRVRP3	ANRVRP4	ANRVRP5	ANRVRP6
ANRVRP7	ANRVRP8	ANRVRP9	ANRVRRL	ANRVSITE

## Installation

### Sample KIDS Install

```
>D ^XUP
```

```
Setting up programmer environment
Terminal Type set to: C-VT100
```

```
Select OPTION NAME: XPD MAIN           Kernel Installation & Distribution
System
```



Select Kernel Installation & Distribution System Option: INSTALLation

- 1 Load a Distribution
- 2 Verify Checksums in Transport Global
- 3 Print Transport Global
- 4 Compare Transport Global to Current System
- 5 Backup a Transport Global
- 6 Install Package(s)  
Restart Install of Package(s)  
Unload a Distribution

Select Installation Option: LOAD a Distribution  
Enter a Host File: ANRV\_4\_0.KID

KIDS Distribution saved on Jun 10, 1998@18:31:57  
Comment: Visual Impairment Service Team (VIST) 4.0

This Distribution contains Transport Globals for the following Package(s):  
Visual Impairment Service Team (VIST) 4.0

Want to Continue with Load? YES// <Enter>  
Loading Distribution...

Visual Impairment Service Team (VIST) 4.0  
Use INSTALL NAME: Visual Impairment Service Team (VIST) 4.0 to install this  
Distribution.

**Note:** At this point, you can run the *Verify Checksum in Transport Global* [XPD PRINT CHECKSUM] option. If any of the routine checksums failed, you should not proceed with the installation process. Contact your distributing CIO Field Office.

Select Installation Option: INSTAll Package(s)  
Select INSTALL NAME: Visual Impairment Service Team (VIST) 4.0  
Loaded from Distribution 6/10/98@12:58:45  
=> VIST 4.0 ;Created on Mar 23, 1998@12:57:50

This Distribution was loaded on Jun 10, 1998@12:58:45 with header of  
VIST TEST INSTALL ;Created on Jun 10, 1998@12:57:50  
It consisted of the following Install(s):  
Visual Impairment Service Team (VIST) 4.0

Visual Impairment Service Team (VIST) 4.0

Install Questions for Visual Impairment Service Team (VIST) 4.0

Incoming Files:

2040 VIST ROSTER

- 2041 VIST PARAMETERS
- 2041.5 VIST EYE DIAGNOSIS (including data)
- 2041.6 VIST CHECKLIST OPTIONS (including data)
- 2041.7 VIST BENEFITS AND SERVICES CHECKLIST (including data)
- 2042 VIST REFERRAL FACILITY (including data)
- 2042.5 VIST REFERRAL ROSTER
- 2043 VIST LETTER (including data)
- 2043.5 VARO CLAIMS
- 2044 VIST LOCAL BENEFITS AND SERVICES (including data)

Want to DISABLE Scheduled Options, Menu Options, and Protocols? YES// **NO**  
 Want to MOVE routines to other CPUs? NO// **<ENTER>**

Enter the Device you want to print the Install messages.  
 You can queue the install by enter a 'Q' at the device prompt.  
 Enter a '^' to abort the install.

DEVICE: ***[Select Print Device]***

Install Started for Visual Impairment Service Team (VIST) 4.0 :  
 Jun 10, 1998@14:28:59

Installing Routines:.....  
 Jun 10, 1998@14:29:04

Installing Data Dictionaries: .....  
 Jun 10, 1998@14:29:20

Installing Data:  
 Jun 10, 1998@14:29:23

Installing PACKAGE COMPONENTS:

Installing FUNCTION..

Installing PRINT TEMPLATE.....  
 .....

```
Installing SORT TEMPLATE.....
Installing INPUT TEMPLATE.....
Installing OPTION.....
      Jun 10, 1998@14:30:25

Running Post-Install Routine: ^ANRVPOST.

Updating Routine file.....

Updating KIDS files.....

Visual Impairment Service Team (VIST) 4.0 Installed.
      Jun 10, 1998@14:30:33..

Install Message sent #26535

  An Option(s) has been added by this build.
  Electing to rebuild menus now is highly recommended.

  Do you want to rebuild Menu's now? Y// NO [It is recommended that you respond "No"
  to this prompt.]
```

## Post Installation



During the installation of VIST V. 4.0, all data entered using the class III Blind Rehabilitation software will be moved into the new files created for VIST V. 4.0. This conversion time will vary depending on the amount of patient information that was entered into the class III files.

### System Configuration

^ANRV is a new global and will need to be journaled after installation. There are no specific guidelines for system configuration or global placement for the Visual Impairment Service Team software.

### Resources Requirements

This version (4.0) of VIST contains approximately 25 routines including all ANRV\* routines and compiled templates that take up approximately 69 K disk space.

### Routine Mapping

It is not necessary to map routines.

### Journaling Globals

It is recommended that the ^ANRV global be journaled after installation.

## Files

This package requires the 10 files listed below. Information about the files can be obtained by using the VA FileMan to generate a list of file attributes.

The Data Dictionaries (DDs) are considered part of the online documentation for this software application. Use VA FileMan option *List File Attributes* [DILIST], under *Data Dictionary Utilities* [DI DDU], to print the DDs. The following are the files for which you should print DDs:

FILE #	NAME	DATE DD	SEC. CODE	COMES W/FILE	SITE DATA	RSLV PTS	OVER RIDE
2040	VIST ROSTER	YES	NO	NO			
2041	VIST PARAMETERS	YES	NO	NO			
2041.5	VIST EYE DIAGNOSIS	YES	NO	YES	MERG	NO	NO
2041.6	VIST CHECKLIST OPTIONS	YES	NO	YES	MERG	NO	NO
2041.7	VIST BENEFITS AND SERVICES CHECKLIST	YES	NO	NO			
2042	VIST REFERRAL FACILITY	YES	NO	YES	MERG	NO	NO
2042.5	VIST REFERRAL ROSTER	YES	NO	NO			
2043	VIST LETTER	YES	NO	YES	MERG	NO	NO
2043.5	VARO CLAIMS	YES	NO	NO			
2044	VIST LOCAL BENEFITS AND SERVICES	YES	NO	YES	MERG	NO	NO

The namespace for the VIST package is ANRV.

### Security Keys

There are **no** security keys for the VIST software package.

### Technical Assistance

If you experience problems of a technical nature with this software package, you are encouraged to contact the CUSTOMER SERVICE NATIONAL HELP DESK.